Elabscience®

DAB Detection Kit

Cat. No: E-IR-R101

Size: 1 mL/3 mL/6mL

Cat	Products	1 mL	3 mL	6 mL	Storage
E-IR-R101A	DAB Concentrate (20×)	1 mL	3 mL	6 mL	2~8°C
E-IR-R101B	DAB Substrate	20 mL	60 mL	120 mL	2~8°C
Manual		One Copy			

Introduction

Elabscience® DAB Detection kit can be used for immunohistochemical development. DAB (3,3N-Diaminobenzidine Tertrahydrochloride) is a common substrate of horseradish peroxidase (HRP), under the catalysis of horseradish peroxidase, DAB will produce a brown precipitate.

Instructions

DAB Concentrate (20×) is a concentrate solution, dilute the DAB Concentrate (20×) with DAB Substrate to 1× DAB working solution before use. For example: add 10 μ L of DAB Concentrate (20×) to 190 μ L of DAB Substrate and mix fully.

Note: The prepared 1× DAB working solution should be used within 4 hours, and be protected from light.

Storage

The kit can be store at 2~8°C with shading light for 12 months. Avoid freezing and thawing cycles, put the reagent back in the refrigerator immediately after use, the reagent is valid for 6 months after opening.

Cautions

- 1. This kit is for research use only.
- 2. For your safety and health, please wear the lab coat and disposable gloves before the experiments.
- It is necessary to control the color developing time of DAB, usually incubate at room temperature for 5~10 minutes. The positive signal is brownish yellow, be careful not to develop too deep. Then rinse the slides with distilled water to stop color development.
- 4. Please dry the liquid around the slice during the experiment to prevent weakening the dyeing intensity.
- 5. The positive and negative controls should be set up during the experiment to prevent possible false positive and false negative results.