

## Dulbecco's Phosphate Buffer (DPBS) Ingredient Formula

Ingredient (mg/L) / Cat. No.	PB180329	PB180330	PB180331	PB180332	PB180333
<b>Inorganic Salt</b>					
CaCl <sub>2</sub> (Anhydrous)	–	50	50	50	–
MgCl <sub>2</sub> •6H <sub>2</sub> O	–	100	100	100	–
KCl	200	200	200	200	200
KH <sub>2</sub> PO <sub>4</sub>	200	200	200	200	200
NaCl	8000	8600	8600	8600	8000
Na <sub>2</sub> HPO <sub>4</sub> •12H <sub>2</sub> O	2886	1150	1150	1150	2886
<b>Others</b>					
D-Glucose	–	–	–	1000	–
Sodium Pyruvate	–	–	–	36	–
Phenol Red	–	10	–	–	10

Pricella<sup>®</sup>  
by Elabscience

Pricella<sup>®</sup>  
by Elabscience