

## Normal Rabbit Blocking Buffer (Ready-to-Use)

Cat. No: E-IR-R128

Size: 5 mL/ 10 mL/ 50 mL

Cat	Products	5 mL	10 mL	50 mL	Storage
E-IR-R128	Normal Rabbit Blocking Buffer (Ready-to-Use)	5 mL	10 mL	50 mL	2~8°C

### Introduction

This Normal Rabbit Blocking Buffer is suitable for sample sealing in immunohistochemical or immunofluorescence experiments. It can block the non-specific binding sites of proteins in the samples to be tested, and also block the endogenous Fc fragment binding sites in the samples to block the binding of antibodies to the Fc receptors in the samples, thus reducing the background and false positive.

When selecting the blocking buffer, it should be noted that the blocking buffer should not contain the target protein, and the source of the blocking buffer is different from that of the primary antibody, which can be blocked with the same source of the secondary antibody.

### Instructions

This product is ready to use Normal Rabbit Blocking Buffer and can be used directly.

### Storage

Store at 2~8°C for 12 months.

### Cautions

1. This product couldn't be used for clinical diagnosis or treatment, food or medicine, and can't be stored in residence.
2. For your safety and health, please wear the lab coat and disposable gloves before the experiments.