




Elabscience®

www.elabscience.com

**ELISA Kits for
Cardiovascular Research**

Catalogue

Blood

Angiogenesis

Atherosclerosis

Vasculature

Lipids Lipoproteins

Heart

Featured Publications



● Cardiovascular

The essential components of human cardiovascular system are heart, blood and blood vessels. It includes the pulmonary circulation, a loop through the lungs where blood is oxygenated; and the systemic circulation, a loop through the rest of the body to provide oxygenated blood. The systemic circulation can also be seen to function in two parts—a macrocirculation and a microcirculation. The cardiovascular system of human is closed, meaning that the blood never leaves the network of blood vessels. In contrast, oxygen and nutrients diffuse across the blood vessel layers and enter interstitial fluid, which carries oxygen and nutrients to the target cells, and carbon dioxide and wastes in the opposite direction.

Cardiovascular disease is the leading cause of death globally. Many biologists are trying to understand the cardiovascular system problems during disease development and solve them. Elabscience offers various and quality ELISA kits for the study associated with blood, angiogenesis, atherosclerosis, vasculature, heart, lipids/lipoproteins and hypoxia.



Blood

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0001	Human AGP1(Alpha-1-Acid Glycoprotein 1) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H0002	Human ACE(Angiotensin I Converting Enzyme) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0040	Human TF(Tissue Factor) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0055	Human Cys-C(Cystatin C) ELISA Kit	6.25--400ng/mL	3.75ng/mL	Sandwich
E-EL-H0070	Human α FP(Alpha-Fetoprotein) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0152	Human CP(Ceruloplasmin) ELISA Kit	0.63--40 μ g/mL	0.38 μ g/mL	Sandwich
E-EL-H0165	Human TFPI(Tissue Factor Pathway Inhibitor) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0166	Human TM(Thrombomodulin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0234	Human α 1ACT(Alpha1 Antichymotrypsin) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-H0760	Human F9(Coagulation Factor IX) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-H0762	Human F13(Coagulation Factor X III) ELISA Kit	25--1600pg/mL	15pg/mL	Sandwich
E-EL-H0763	Human F12(Coagulation Factor XII) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0766	Human F13A1(Coagulation Factor X III A1 Polypeptide) ELISA Kit	18.75--1200pg/mL	11.25pg/mL	Sandwich
E-EL-H0767	Human F13B(Coagulation Factor X III B Polypeptide) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H0967	Human α2-M(Alpha-2 Macroglobulin) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H1507	Human HP(Haptoglobin) ELISA Kit	12.5--800ng/mL	7.5ng/mL	Sandwich
E-EL-H1577	Human PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1588	Human TPO(Thrombopoietin) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1788	Human TAT(Thrombin/Antithrombin Complex) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1792	Human PT(Prothrombin) ELISA Kit	9.38--600ng/mL	5.63ng/mL	Sandwich
E-EL-H1855	Human KLK1(Kallikrein 1) ELISA Kit	18.75--1200pg/mL	11.25pg/mL	Sandwich
E-EL-H2102	Human Plg(Plasminogen) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H2106	Human tPA(Plasminogen Activator, Tissue) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-H2168	Human VWF(Von Willebrand Factor) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H2172	Human HO1(Heme Oxygenase 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H2183	Human SAA(Serum Amyloid A) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-H2193	Human FG(Fibrinogen) ELISA Kit	9.38--600ng/mL	5.63ng/mL	Sandwich
E-EL-H2216	Human PF4(Platelet Factor 4) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H2255	Human KNG1(Kininogen 1) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H2342	Human TTR(Transthyretin) ELISA Kit	15.63--1000ng/mL	9.38ng/mL	Sandwich
E-EL-H5200	Human LTF/LF(Lactoferrin) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0091	Mouse ACE(Angiotensin I Converting Enzyme) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0304	Mouse F II (coagulation factor II) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-M0305	Mouse GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0389	Mouse Cys-C(Cystatin C) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0498	Mouse FG(Fibrinogen) ELISA Kit	15.63--1000ng/mL	9.38ng/mL	Sandwich
E-EL-M0613	Mouse HP(Haptoglobin) ELISA Kit	78.13--5000ng/mL	46.88ng/mL	Sandwich
E-EL-M0640	Mouse TPO(Thrombopoietin) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0651	Mouse HO1(Heme Oxygenase 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0826	Mouse NAP-2(Neutrophil Activating Protein-2) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-M0926	Mouse PF4(Platelet Factor 4) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M1046	Mouse SAA(Serum Amyloid A) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-M1138	Mouse TAT(Thrombin/Antithrombin Complex) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M2428	Mouse α 1-AT(Alpha 1-Antitrypsin) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-R0182	Rat CP(Ceruloplasmin) ELISA Kit	6.25--400ng/mL	3.75ng/mL	Sandwich
E-EL-R0473	Rat HP(Haptoglobin) ELISA Kit	15.63--1000ng/mL	9.38ng/mL	Sandwich
E-EL-R0661	Rat NAP-2(Neutrophil Activating Protein-2) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R0748	Rat Plg(Plasminogen) ELISA Kit	7.81--500ng/mL	4.69ng/mL	Sandwich
E-EL-R1081	Rat α 1-AGP(Alpha-1-Acid Glycoprotein) ELISA Kit	6.25--400ng/mL	3.75ng/mL	Sandwich
E-EL-R1125	Rat FG(Fibrinogen) ELISA Kit	62.5--4000ng/mL	37.5ng/mL	Sandwich
E-EL-R1139	Rat α 2-M(Alpha-2 Macroglobulin) ELISA Kit	9.38--600ng/mL	5.63ng/mL	Sandwich



Angiogenesis

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0007	Human ANG1(Angiopoietin 1) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0008	Human ANG2(Angiopoietin 2) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-H0012	Human BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0048	Human IL-8(Interleukin 8) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0062	Human ENG(Endoglin) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0066	Human EPO(Erythropoietin) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0071	Human AFGF/FGF1(Acidic Fibroblast Growth Factor 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0074	Human FGF21(Fibroblast Growth Factor 21) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0075	Human Flt3L(FMS-like Tyrosine Kinase 3 Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0111	Human VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0132	Human MMP-10(Matrix Metalloproteinase 10) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0134	Human MMP-13(Matrix Metalloproteinase 13) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0150	Human IL-12(Interleukin 12) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0151	Human IL-12 p40(Interleukin 12 p40) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0184	Human TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0334	Human ANGPTL1(Angiopoietin Like Protein 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0335	Human ANGPTL2(Angiopoietin Like Protein 2) ELISA Kit	0.47--30ng/mL	0.28ng/mL	Sandwich
E-EL-H0336	Human ANGPTL3(Angiopoietin Like Protein 3) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0337	Human ANGPTL4(Angiopoietin Like Protein 4) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0338	Human ANGPTL6(Angiopoietin Like Protein 6) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0483	Human bFGF/FGF2(Basic Fibroblast Growth Factor) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H0552	Human FGF4(Fibroblast Growth Factor 4) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0701	Human CCL16(Chemokine C-C-Motif Ligand 16) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0985	Human PDGFR β (Platelet Derived Growth Factor Receptor Beta) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H1087	Human VEGFR1/FLT1(Vascular Endothelial Growth Factor Receptor 1) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H1117	Human FGF9(Fibroblast Growth Factor 9) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H1347	Human OPN(Osteopontin) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1446	Human MMP-3(Matrix Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H1449	Human MMP-7(Matrix Metalloproteinase 7) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H1450	Human MMP-8(Matrix Metalloproteinase 8) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H1451	Human MMP-9(Matrix Metalloproteinase 9) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1530	Human MIF(Macrophage Migration Inhibitory Factor) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H1555	Human PGF(Placental Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1575	Human PDGF-AA(Platelet Derived Growth Factor AA) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1576	Human PDGF-AB(Platelet Derived Growth Factor AB) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1586	Human TGF- α (Transforming Growth Factor Alpha) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1600	Human VEGF-C(Vascular Endothelial Growth Factor C) ELISA Kit	109.38--7000pg/mL	65.63pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H1603	Human VEGFR-2/KDR(Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H2247	Human OSM(Oncostatin M) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H2312	Human TWEAK(Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0027	Mouse EPO(Erythropoietin) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0098	Mouse ANG2(Angiopoietin 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0626	Mouse MMP-3(Matrix Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-M0629	Mouse OPN(Osteopontin) ELISA Kit	39.06--2500pg/mL	23.44pg/mL	Sandwich
E-EL-M0641	Mouse TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-M0784	Mouse MMP-8(Matrix Metalloproteinase 8) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M1232	Mouse VEGFR3/Flt4(Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0540	Rat TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R0623	Rat MMP-8(Matrix Metalloproteinase 8) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0702	Rat OPN(Osteopontin) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich



Atherosclerosis

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0004	Human ADP/Acrp30(Adiponectin) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0020	Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0023	Human RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0036	Human sCD163(Soluble Cluster of Differentiation 163) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0043	Human CRP(C-Reactive Protein) ELISA Kit	0.39--25ng/mL	0.23ng/mL	Sandwich
E-EL-H0044	Human CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0046	Human ENA-78(Epithelial Neutrophil Activating Peptide 78) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0054	Human CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0055	Human Cys-C(Cystatin C) ELISA Kit	6.25--400ng/mL	3.75ng/mL	Sandwich
E-EL-H0064	Human ET-1(Endothelin 1) ELISA Kit	1.25--80pg/mL	0.75pg/mL	Sandwich
E-EL-H0087	Human IGFBP-3(Insulin-like Growth Factor Binding Protein 3) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0097	Human MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0113	Human LEP(Leptin) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0184	Human TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0295	Human AGER(Advanced Glycosylation End Product Specific Receptor) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0422	Human ANXA5(Annexin A5) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0442	Human IGFBP-1(Insulin-like Growth Factor Binding Protein 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0463	Human ApoA4(Apolipoprotein A4) ELISA Kit	4.69--300ng/mL	2.81ng/mL	Sandwich
E-EL-H0464	Human ApoB(Apolipoprotein B) ELISA Kit	78.13--5000ng/mL	46.88ng/mL	Sandwich
E-EL-H0466	Human ApoC2(Apolipoprotein C2) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H0651	Human CT1(Cardiotrophin 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0876	Human sSELE(Soluble E-Selectin) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0951	Human TLR-2(Toll-like Receptor 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0999	Human CYSLTR1(Cysteinyl Leukotriene Receptor 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H1050	Human CYSLTR2(Cysteinyl Leukotriene Receptor 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1213	Human RETN(Resistin) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1387	Human RANKL(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H1446	Human MMP-3(Matrix Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H1450	Human MMP-8(Matrix Metalloproteinase 8) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H1451	Human MMP-9(Matrix Metalloproteinase 9) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1574	Human PTX 3/TSG-14(Pentraxin 3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1642	Human ApoB48(Apolipoprotein B48) ELISA Kit	15.63--1000ng/mL	9.38ng/mL	Sandwich
E-EL-H1967	Human MyD88(Myeloid Differentiation Factor 88) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H2166	Human sVCAM-1/CD106(soluble Vascular Cell Adhesion Molecule 1) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H2172	Human HO1(Heme Oxygenase 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H2183	Human SAA(Serum Amyloid A) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-H2193	Human FG(Fibrinogen) ELISA Kit	9.38--600ng/mL	5.63ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H2214	Human PDGFR α (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H2216	Human PF4(Platelet Factor 4) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H2247	Human OSM(Oncostatin M) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H2484	Human CD200(Cluster of Differentiation 200) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H5587	Human VCAM-1/CD106(Vascular Cell Adhesion Molecule 1) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-M0002	Mouse ADP/Acrp30(Adiponectin) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0006	Mouse MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0009	Mouse RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0039	Mouse LEP(Leptin) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-M0149	Mouse AGER(Advanced Glycosylation End Product Specific Receptor) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-M0267	Mouse CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0270	Mouse CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0282	Mouse Gas6(Growth Arrest Specific Protein 6) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0389	Mouse Cys-C(Cystatin C) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0457	Mouse ET-1(Endothelin 1) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0471	Mouse ENA-78(Epithelial Neutrophil Activating Peptide 78) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0474	Mouse SELE(E-selectin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0515	Mouse FK(Fractalkine) ELISA Kit	0.39--25ng/mL	0.23ng/mL	Sandwich
E-EL-M0626	Mouse MMP-3(Matrix Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-M0627	Mouse MMP-9(Matrix Metalloproteinase-9) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0635	Mouse RETN(Resistin) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0641	Mouse TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-M0644	Mouse RANkL(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0647	Mouse sVCAM-1/CD106(soluble Vascular Cell Adhesion Molecule 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-M0651	Mouse HO1(Heme Oxygenase 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0784	Mouse MMP-8(Matrix Metalloproteinase 8) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0859	Mouse OSM(Oncostatin M) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0926	Mouse PF4(Platelet Factor 4) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0929	Mouse PDGFsR- α (Platelet-Derived Growth Factor Soluble Receptor α) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M1046	Mouse SAA(Serum Amyloid A) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-M2445	Mouse MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0022	Rat CRP(C-Reactive Protein) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0167	Rat ET-1(Endothelin 1) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-R0329	Rat ADP/Acrp30(Adiponectin) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-R0370	Rat ENA-78(Epithelial Neutrophil Activating Peptide 78) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0385	Rat CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R0540	Rat TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R0582	Rat LEP(Leptin) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R0601	Rat MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0619	Rat MMP-3(Matrix Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R0623	Rat MMP-8(Matrix Metalloproteinase 8) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0624	Rat MMP-9(Matrix Metalloproteinase 9) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-R0633	Rat MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R0643	Rat AGER(Advanced Glycosylation End Product Specific Receptor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0692	Rat OSM(Oncostatin M) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0757	Rat PDGFsR- α (Platelet-Derived Growth Factor Soluble Receptor α) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R0845	Rat RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0881	Rat SAA(Serum Amyloid A) ELISA Kit	0.16--10 μ g/mL	0.1 μ g/mL	Sandwich
E-EL-R1002	Rat PTX-3/TSG-14(Pentraxin 3) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich



Vasculature

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0002	Human ACE(Angiotensin I Converting Enzyme) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0062	Human ENG(Endoglin) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H0272	Human ADH(Antidiuretic Hormone) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Competitive
E-EL-H0326	Human Ang- II (Angiotensin II) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Competitive
E-EL-H0327	Human Ang- I (Angiotensin I) ELISA Kit	125--8000pg/mL	75pg/mL	Competitive
E-EL-H0596	Human BK(Bradykinin) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Competitive
E-EL-H0598	Human BNP(Brain Natriuretic Peptide) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Competitive
E-EL-H1003	Human GUS β (β -glucuronidase) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H1574	Human PTX 3/TSG-14(Pentraxin 3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1964	Human Anti-MPO(Anti-Myeloperoxidase Antibody) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H2060	Human GT(Gastrin) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Competitive



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H2172	Human HO1(Heme Oxygenase 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H2358	Human 11,12-DHETs(11,12-dihydroxyeicosatrienoic acids) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Competitive
E-EL-H5309	Human sENG(Soluble Endoglin) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-M0091	Mouse ACE(Angiotensin I Converting Enzyme) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0202	Mouse BK(Bradykinin) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Competitive
E-EL-M0204	Mouse BNP(Brain Natriuretic Peptide) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Competitive
E-EL-M0567	Mouse GUSb(Glucuronidase, Beta) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0651	Mouse HO1(Heme Oxygenase 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M1072	Mouse ENG(Endoglin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M2611	Mouse Ang- I (Angiotensin I) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Competitive
E-EL-M2612	Mouse Ang- II (Angiotensin II) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Competitive
E-EL-M2669	Mouse GT(Gastrin) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Competitive
E-EL-R0124	Rat BK(Bradykinin) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Competitive



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0126	Rat BNP(Brain Natriuretic Peptide) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Competitive
E-EL-R0432	Rat GUSb(Glucuronidase, Beta) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0472	Rat GT(Gastrin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Competitive
E-EL-R0488	Rat HO1(Heme Oxygenase 1) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R1002	Rat PTX-3/TSG-14(Pentraxin 3) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-R1429	Rat ANG I (Angiotensin I) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Competitive
E-EL-R1430	Rat ANG II (Angiotensin II) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Competitive
E-EL-R2401	Rat ACE(Angiotensin I Converting Enzyme) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R2434	Rat ENG(Endoglin) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich



Lipids/Lipoproteins

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0001	Human AGP1(Alpha-1-Acid Glycoprotein 1) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H0056	Human Cyt-C(Cytochrome C) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0125	Human ApoA1(Apolipoprotein A1) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H0337	Human ANGPTL4(Angiopietin Like Protein 4) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0462	Human ApoA2(Apolipoprotein A2) ELISA Kit	62.5--4000ng/mL	37.5ng/mL	Sandwich
E-EL-H0464	Human ApoB(Apolipoprotein B) ELISA Kit	78.13--5000ng/mL	46.88ng/mL	Sandwich
E-EL-H0465	Human ApoC1(Apolipoprotein C1) ELISA Kit	9.38--600ng/mL	5.63ng/mL	Sandwich
E-EL-H0466	Human ApoC2(Apolipoprotein C2) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H0467	Human ApoC3(Apolipoprotein C3) ELISA Kit	7.81--500ng/mL	4.69ng/mL	Sandwich
E-EL-H0470	Human ApoE(Apolipoprotein E) ELISA Kit	23.44--1500ng/mL	14.06ng/mL	Sandwich
E-EL-H0472	Human ApoH(Apolipoprotein H) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0698	Human CHEM(Chemerin) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-M0002	Mouse ADP/Acrp30(Adiponectin) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0015	Mouse CHEM(Chemerin) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-M1046	Mouse SAA(Serum Amyloid A) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-R0329	Rat ADP/Acrp30(Adiponectin) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-R2435	Rat FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich



Heart

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0142	Human MYO(Myoglobin) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H0327	Human Ang- I (Angiotensin I) ELISA Kit	125--8000pg/mL	75pg/mL	Competitive
E-EL-H0937	Human PIAS3(Protein Inhibitor of Activated STAT3) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H1577	Human PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1655	Human mTOR(Mammalian Target of Rapamycin) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich
E-EL-H1867	Human HSF1(Heat Shock Factor 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H5316	Human HDAC1(Histone Deacetylase 1) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-M0084	Mouse MYO(Myoglobin) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-M0236	Mouse CT1(Cardiotrophin 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0614	Mouse HSF1(Heat Shock Factor 1) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-M1057	Mouse STAT3(Signal Transducer And Activator Of Transcription 3) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0022	Rat CRP(C-Reactive Protein) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0053	Rat MYO(Myoglobin) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-R1226	Rat STAT3(Signal Transducer and Activator of Transcription 3) ELISA Kit	0.16--10ng/mL	0.1ng/mL	Sandwich



Featured Publications

Cat. No.	Targets	Species	Publications
E-EL-H2340 E-EL-H0965 E-EL-H2447 E-EL-H5200 E-EL-H1163	4-HNE(4-Hydroxynonenal) α 1-AT(Alpha 1-Antitrypsin) PI3/Elafin(Peptidase Inhibitor 3, Skin Derived) LTF/LF(Lactoferrin)ELN(Elastin)	Human	Liu R, Chen L, Wu W, et al. Neutrophil Serine Proteases And Their Endogenous Inhibitors In Coronary Artery Ectasia Patients[J]. Anatolian Journal of Cardiology, 2015, 16(1): 23-28.
E-EL-H2164 E-EL-H0483	VEGF-B(Vascular Endothelial Cell Growth Factor B) bFGF/FGF2(Basic Fibroblast Growth Factor)	Human	Ahmed F, Mohammed H. Effect Of Coronary Artery Bypass Grafting On The Plasma Levels Of VEGF, bFGF, Copper And NO As Angiogenic Biomarkers In Chronic Ischemic Heart Disease Patients[J]. Annals of Clinical and Laboratory Research, 2016, 4(3): 107-113.
E-EL-H2102 E-EL-H0965 E-EL-H0968 E-EL-H2340	Plg(Plasminogen) α 1-AT(Alpha 1-Antitrypsin) α 2-AP(α 2-Antiplasmin) 4-HNE(4-Hydroxynonenal)	Human	Wu W, Liu R, Chen L, et al. Disequilibrium Of Blood Coagulation And Fibrinolytic System In Patients With Coronary Artery Ectasia[J]. Medicine, 2016.
E-EL-H1779 E-EL-H1445 E-EL-H1451	PLAU/uPA(Urokinase-Type Plasminogen Activator) MMP-2(Matrix Metalloproteinase 2) MMP-9(Matrix Metalloproteinase 9)	Human	Wang Z, Yin H, Wang T, et al. Downregulation Of Foxm1 Inhibits Proliferation, Invasion And Angiogenesis Of Hela Cells In Vitro And In Vivo[J]. Bone, 2014, 45: 2355-2364.
E-EL-H1586	TGF- α (Transforming Growth Factor Alpha)	Human	Chen Y, Bai X, Qi Z, et al. The Hepatitis B Virus X Protein Promotes Pancreatic Cancer Through Modulation Of The PI3K/AKT Signaling Pathway[J]. Cancer Letters, 2016, 380(1): 98-105.



Cat. No.	Targets	Species	Publications
E-EL-H1539 E-EL-H0951 E-EL-H1655	TLR4(Toll-Like Receptor 4) TLR-2(Toll-like Receptor 2) mTOR(Mammalian Target of Rapamycin)	Human	Kang M, Min K, Jang J, et al. Involvement Of Immune Responses In The Efficacy Of Cord Blood Cell Therapy For Cerebral Palsy[J]. Journal of Back & Musculoskeletal Rehabilitation, 2015, 24(19): 2259-2268.
E-EL-H1507	HP(Haptoglobin)	Human	Esnafoğlu E, Cırrık S, Ayyıldız S N, et al. Increased Serum Zonulin Levels As An Intestinal Permeability Marker In Autistic Subjects[J]. The Journal of Pediatrics, 2017.
E-EL-H1507	HP(Haptoglobin)	Human	Cangemi R, Della V P, Calvieri C, et al. Low-Grade Endotoxemia And Clotting Activation In The Early Phase Of Pneumonia[J]. Respirology, 2016, 21(8): 1465-1471.
E-EL-H1507	HP(Haptoglobin)	Human	Cangemi R, Pignatelli P, Carnevale R, et al. Low-Grade Endotoxemia, Gut Permeability And Platelet Activation In Community-Acquired Pneumonia[J]. Journal of Infection, 2016, 73(2): 107-114.
E-EL-H1507	HP(Haptoglobin)	Human	Raparelli V, Basili S, Carnevale R, et al. Low-Grade Endotoxemia And Platelet Activation In Cirrhosis[J]. Hepatology, 2017, 65(2): 571-581.
E-EL-H1451 E-EL-H2340	MMP-9(Matrix Metalloproteinase 9) 4-HNE(4-Hydroxynonenal)	Human	Wang D, Ma Y, Yang X, et al. Hypermethylation Of The Keap1 Gene Inactivates Its Function, Promotes Nrf2 Nuclear Accumulation, And Is Involved In Arsenite-Induced Human Keratinocyte Transformation[J]. Free Radical Biology and Medicine, 2015, 89: 209-219.



Cat. No.	Targets	Species	Publications
E-EL-H1451	MMP-9(Matrix Metalloproteinase 9)	Human	Mohamad S A, Gafar H S, Elhady M, et al. Cord Blood Matrix Metalloproteinase-9 As Predictor Of Respiratory Outcome In Preterm Infants[J]. Clinical Medicine and Diagnostics, 2015, 5(4): 53-59.
E-EL-H1451	MMP-9(Matrix Metalloproteinase 9)	Human	Yang X, Wang D, Ma Y, et al. Continuous Activation Of Nrf2 And Its Target Antioxidant Enzymes Leads To Arsenite-Induced Malignant Transformation Of Human Bronchial Epithelial Cells[J]. Toxicology & Applied Pharmacology, 2015, 289(2): 231-239.
E-EL-H0698	CHEM(Chemerin)	Human	Calvet J, Orellana C, Gratacós J, et al. Synovial Fluid Adipokines Are Associated With Clinical Severity In Knee Osteoarthritis: A Cross-Sectional Study In Female Patients With Joint Effusion[J]. Arthritis Research & Therapy, 2016.
E-EL-H0598	BNP(Brain Natriuretic Peptide)	Human	Gupta R K, Himashree G, Singh K, et al. Elevated Pulmonary Artery Pressure And Brain Natriuretic Peptide In High Altitude Pulmonary Edema Susceptible Non-Mountaineers[J]. Scientific Reports, 2016.
E-EL-H0463	ApoA4(Apolipoprotein A4)	Human	Luczak M, Formanowicz D, Łukasz Marczak, et al. Itraq-Based Proteomic Analysis Of Plasma Reveals Abnormalities In Lipid Metabolism Proteins In Chronic Kidney Disease-Related Atherosclerosis[J]. Scientific Reports, 2016.
E-EL-H0427 E-EL-H0148 E-EL-H1642 E-EL-H0058	IAA(Insulin Autoantibody) GLP-1(Glucagon Like Peptide 1) ApoB48(Apolipoprotein B48) DPP4(Dipeptidyl Peptidase IV)	Human	Carnevale R, Loffredo L, Ben M D, et al. Extra Virgin Olive Oil Improves Post-Prandial Glycemic And Lipid Profile In Patients With Impaired Fasting Glucose[J]. Clinical Nutrition, 2016.



Cat. No.	Targets	Species	Publications
E-EL-H0315 E-EL-H0055 E-EL-H2188 E-EL-H0463 E-EL-H2110 E-EL-H2111 E-EL-H0387	AMBP(Alpha-1-Microglobulin/Bikunin Precursor) Cys-C(Cystatin C)BMG/β2-MG(Beta-2-Microglobulin) ApoA4(Apolipoprotein A4) FGα(Fibrinogen Alpha) FGβ(Fibrinogen Beta) FGγ(Fibrinogen Gamma chain)	Human	Luczak M, Suszynska-Zajczyk J, Marczak L, et al. Label-Free Quantitative Proteomics Reveals Differences In Molecular Mechanism Of Atherosclerosis Related And Non-Related To Chronic Kidney Disease[J]. International Journal of Molecular Sciences, 2016.
E-EL-H0150	IL-12(Interleukin 12)	Human	Jiang L, Hu S, Tian D, et al. A-Mannosidase I Protein Expression In Peripheral Blood Mononuclear Cells Is Upregulated During Hepatitis B Virus Infection[J]. Viral Immunology, 2016, 29(1): 33-39.
E-EL-H0119 E-EL-H0300 E-EL-H0002 E-EL-H0326 E-EL-H2164 E-EL-H1277	REN(Renin) AGT(Angiotensinogen) ACE(Angiotensin I Converting Enzyme) Ang- II (Angiotensin II) VEGF-B(Vascular Endothelial Cell Growth Factor B) HIF-1α(Hypoxia Inducible Factor 1 Alpha)	Human	Nath M, Chandra P, Halder N, et al. Involvement Of Renin-Angiotensin System In Retinopathy Of Prematurity - A Possible Target For Therapeutic Intervention[J]. Plos One, 2016.
E-EL-H0111 E-EL-H1451	VEGF-A(Vascular Endothelial Cell Growth Factor A) MMP-9(Matrix Metalloproteinase 9)	Human	Zhang F, Lu Y, Chen Q, et al. Identification Of NCK1 As A Novel Downstream Effector Of STAT3 In Colorectal Cancer Metastasis And Angiogenesis[J]. Cellular Signalling, 2017.
E-EL-H0111	VEGF-A(Vascular Endothelial Cell Growth Factor A)	Human	He Y H, Li M F, Zhang X Y, et al. NLRC5 Promotes Cell Proliferation Via Regulating The AKT/VEGF-A Signaling Pathway In Hepatocellular Carcinoma[J]. Toxicology, 2016, 359: 47-57.



Cat. No.	Targets	Species	Publications
E-EL-H0110 E-EL-H2443 E-EL-H2172	TGF- β 1(Transforming Growth Factor Beta 1) IL-35(Interleukin 35) HO1(Heme Oxygenase 1)	Human	Li Y, Dong L, Wei Y, et al. Overexpression Of The Transcription Factor Foxp3 In Lung Adenocarcinoma Sustains Malignant Character By Promoting G1/S Transition Gene CCND1[J]. Tumor Biology, 2016.
E-EL-H0109 E-EL-H0150 E-EL-H0149	TNF- α (Tumor Necrosis Factor Alpha) IL-12(Interleukin 12) IL-1 β (Interleukin 1 Beta)	Human	Jiang Y, Kou J, Han X, et al. Ros-Dependent Activation Of Autophagy Through The PI3K/Akt/mTOR Pathway Is Induced By Hydroxysafflor Yellow A-Sonodynamic Therapy In THP-1 Macrophages[J]. Oxidative Medicine and Cellular Longevity, 2017.
E-EL-H0108 E-EL-H0109 E-EL-H0150 E-EL-H0102	IFN- γ (Interferon Gamma) TNF- α (Tumor Necrosis Factor Alpha) IL-12(Interleukin 12) IL-6(Interleukin 6)	Human	He W, Ren Y, Wang X, et al. C Reactive Protein And Enzymatically Modified LDL Cooperatively Promote Dendritic Cell-Mediated T Cell Activation[J]. Cardiovascular Pathology, 2017, 29: 1-6.
E-EL-H0102 E-EL-H0048 E-EL-H0149	IL-6(Interleukin 6) IL-8(Interleukin 8) IL-1 β (Interleukin 1 Beta)	Human	Tian P, Li B, He C, et al. Antidiabetic (Type 2) Effects Of Lactobacillus G15 And Q14 In Rats Through Regulation Of Intestinal Permeability And Microbiota[J]. Food & Function, 2016.
E-EL-H0088 E-EL-H0048	IL-1 α (Interleukin 1 Alpha) IL-8(Interleukin 8)	Human	Tang Y, Jin X, Xiang Y, et al. The lncRNA MALAT1 Protects The Endothelium Against ox-LDL-Induced Dysfunction Via Upregulating The Expression Of The miR-22-3p Target Genes CXCR2 And AKT[J]. Febs Letters, 2015, 589(20): 3189-3196.
E-EL-H0071	AFGF/FGF1(Acidic Fibroblast Growth Factor 1)	Human	Bai Y P, Shang K, Chen H, et al. FGF-1/-3/FGFR4 Signaling In Cancer-Associated Fibroblasts Promotes Tumor Progression In Colon Cancer Through Erk And MMP-7[J]. Cancer Science, 2015, 106(10): 1278-1287.



Cat. No.	Targets	Species	Publications
E-EL-H0064 E-EL-H1976	ET-1(Endothelin 1) PLAP/ALPP(Placental Alkaline Phosphatase)	Human	Uyanikoglu H, Turp A B, Hilali N G, et al. Serum Endothelin-1 And Placental Alkaline Phosphatase Levels In Placenta Percreta And Normal Pregnancies[J]. The Journal of Maternal-Fetal & Neonatal Medicine, 2017: 1-6.
E-EL-H0056	Cyt-C(Cytochrome C)	Human	Emokpae M A, Chima H N, Ahmed M. Seminal Plasma Caspase 3, Cytochrome C And Total Antioxidant Capacity In Oligospermic Males And Association With Sperm Indices[J]. Journal of Experimental and Integrative Medicine, 2016.
E-EL-H0008	ANG2(Angiopoietin 2)	Human	El-Asrar M A, Elbarbary N S, Ismail E A R, et al. Circulating Angiopoietin-2 Levels In Children And Adolescents With Type 1 Diabetes Mellitus: Relation To Carotid And Aortic Intima-Media Thickness[J]. Angiogenesis, 2016, 19(3): 421-431.
E-EL-H0004	ADP/Acrp30(Adiponectin)	Human	Shanaki M, Fadaei R, Moradi N, et al. The Circulating CTRP13 In Type 2 Diabetes And Non-Alcoholic Fatty Liver Patients[J]. Plos One, 2016.
E-EL-H0004	ADP/Acrp30(Adiponectin)	Human	Gao S, Cheng Y, Zhao L, et al. The Relationships Of Irisin With Bone Mineral Density And Body Composition In PCOS Patients[J]. Diabetes/Metabolism Research & Reviews, 2015, 32: 421-428.
E-EL-M1272 E-EL-M0160	ALT(Alanine Transaminase) AST(Aspartate Aminotransferase)	Mouse	Wu T, Qi X, Liu Y, et al. Dietary Supplementation With Purified Mulberry (Morus Australis, Poir) Anthocyanins Suppresses Body Weight Gain In High-Fat Diet Fed C57Bl/6 Mice[J]. Food Chemistry, 2013, 141(1): 482-487.



Cat. No.	Targets	Species	Publications
E-EL-R2466 E-EL-R0329 E-EL-R0022 E-EL-R0015 E-EL-R0019	INS(Insulin) ADP/Acrp30(Adiponectin) CRP(C-Reactive Protein) IL-6(Interleukin 6) TNF- α (Tumor Necrosis Factor Alpha)	Rat	Zhang Y, Imam M U, Ismail M, et al. High Fat Diet-Induced Inflammation And Oxidative Stress Are Attenuated By N-Acetylneuraminic Acid In Rats[J]. Journal of Biomedical Science, 2015.
E-EL-R2466 E-EL-R0329	INS(Insulin) ADP/Acrp30(Adiponectin)	Rat	Morakinyo A O, Samuel T A, Adekunbi D A, et al. Niacin Improves Adiponectin Secretion, Glucose Tolerance And Insulin Sensitivity In Diet-Induced Obese Rats[J]. Egyptian Journal of Basic & Applied Sciences, 2015, 2(4): 261-267.
E-EL-R1430 E-EL-R0030	ANG II (Angiotensin II) REN(Renin)	Rat	Hohl M. Modulation Of The Sympathetic Nervous System By Renal Denervation Prevents Reduction Of Aortic Distensibility In Atherosclerosis Prone ApoE-Deficient Rats[J]. Journal of Translational Medicine, 2016.
E-EL-R1327 E-EL-R0338 E-EL-R0151 E-EL-R0520 E-EL-R0633	CKMB(Creatine Kinase MB Isoenzyme) D-LDH(D-Lactate Dehydrogenase) cTnT/TNNT2(Troponin T Type 2, Cardiac) NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) MCP-1(Monocyte Chemotactic Protein 1)	Rat	Lin X, Wu M, Liu B, et al. Candesartan Ameliorates Acute Myocardial Infarction In Rats Through Inducible Nitric Oxide Synthase, Nuclear Factor- κ B, Monocyte Chemoattractant Protein-1, Activator Protein-1 And Restoration Of Heat Shock Protein 72[J]. Molecular Medicine Reports, 2015, 12(6): 8193-8200.
E-EL-R0710 E-EL-R0965 E-EL-R0966 E-EL-R0582 E-EL-R1079	OxLDL(Oxidized Low Density Lipoprotein) TXB2(Thromboxane B2) TXB3(Thromboxanes B3) LEP(Leptin) VWF(Von Willebrand Factor)	Rat	Zhang Y, Imam M U, Ismail M, et al. Edible Bird'S Nest Attenuates Procoagulation Effects Of High-Fat Diet In Rats[J]. Drug Design Development and Therapy, 2015, 9 : 3951-3959.



Cat. No.	Targets	Species	Publications
E-EL-R0633 E-EL-R0019 E-EL-R0015 E-EL-R0016	MCP-1(Monocyte Chemotactic Protein 1) TNF- α (Tumor Necrosis Factor Alpha) IL-6(Interleukin 6) IL-10(Interleukin 10)	Rat	Wu D, Luo N, Wang L, et al. Hydrogen Sulfide Ameliorates Chronic Renal Failure In Rats By Inhibiting Apoptosis And Inflammation Through ROS/MAPK And NF- κ B Signaling Pathways[J]. Scientific Reports (Nature Publisher Group), 2017.
E-EL-R0624 E-EL-R0618	MMP-9(Matrix Metalloproteinase 9) MMP-2(Matrix Metalloproteinase 2)	Rat	Liu Y, Zhong G, Zhao Z. Effect Of Chronic Left Ventricular Unloading On Myocardial Remodeling: Multimodal Assessment Of Two Heterotopic Heart Transplantation Techniques[J]. Journal of Heart and Lung Transplantation, 2014, 34(4): 594-603.
E-EL-R0624	MMP-9(Matrix Metalloproteinase 9)	Rat	Lu H, Liu X, Nan Z, et al. Neuroprotective Effects Of Brain-Derived Neurotrophic Factor And Noggin-Modified Bone Mesenchymal Stem Cells In Focal Cerebral Ischemia In Rats[J]. Journal of Stroke & Cerebrovascular Diseases, 2016, 25(2): 410-418.
E-EL-R0582 E-EL-R0710 E-EL-R1079	LEP(Leptin) OxLDL(Oxidized Low Density Lipoprotein) VWF(Von Willebrand Factor)	Rat	Zhang Y, Imam M U, Ismail M, et al. N-Acetylneuraminic Acid Attenuates Hypercoagulation On High Fat Diet-Induced Hyperlipidemic Rats[J]. Food & Nutrition Research, 2015, 59: 29046-29055.
E-EL-R0582 E-EL-0041	LEP(Leptin) 8-epi-PGF 2α (8-Epi-Prostaglandin F2 Alpha)	Rat&Univer- sal	Yida Z, Imam Mu, Ismail M, et al. Edible Bird'S Nest Prevents High Fat Diet-Induced Insulin Resistance In Rats[J]. Journal of Diabetes Research, 2015.



Cat. No.	Targets	Species	Publications
E-EL-R0582	LEP(Leptin)	Rat	Adamu H A, Imam M U, Ooi D J, et al. In Utero Exposure To Germinated Brown Rice And Its Oryzanol-Rich Extract Attenuated High Fat Diet-Induced Insulin Resistance In F1 Generation Of Rats[J]. BMC Complementary and Alternative Medicine, 2017, 17: 67-75.
E-EL-R0582	LEP(Leptin)	Rat	López Y R, Péreztorres I, Zúñigamuñoz A, et al. Effect Of Glycine On Adipocyte Hypertrophy In A Metabolic Syndrome Rat Model[J]. Current Drug Delivery, 2016, 13(1): 158-169.
E-EL-R0912 E-EL-R2492 E-EL-R0872 E-EL-R0329 E-EL-R0689	LDL-C(low density lipoprotein cholesterol) CHEM(Chemerin) HDL-C(high density lipoprotein cholesterol) ADP/Acrp30(Adiponectin) ITLN1(Intelectin 1/Omentin)	Rat	Wu Q M, Ni H X, Lu X. Changes Of Adipocytokine Expression After Diabetic Rats Received Sitagliptin And The Molecular Mechanism [J]. Asian Pacific Journal of Tropical Medicine, 2016.
E-EL-R0391 E-EL-R0026 E-EL-R0582	FSH(Follicle Stimulating Hormone) LH(Luteinizing Hormone) LEP(Leptin)	Rat	Yang M, Yan W J, Yin T L, et al. Curcumin Ameliorates High-Fat Diet-Induced Spermatogenesis Dysfunction[J]. Molecular Medicine Reports, 2016, 14(4): 3588-3594.
E-EL-R0182 E-EL-R0473 E-EL-R0791	CP(Ceruloplasmin) HP(Haptoglobin) PTGS1(Prostaglandin Endoperoxide Synthase 1)	Rat	Hou D, Zhu R, Sun Z, et al. Serum Proteomics Analysis In Rats Of Immunosuppression Induced By Chronic Stress[J]. Scandinavian Journal of Immunology, 2016, 84(3): 165-173.
E-EL-R0167 E-EL-R0046 E-EL-R1061	ET-1(Endothelin 1) sICAM-1/CD54(Soluble Intercellular Adhesion Molecule 1) VCAM-1/CD106(Vascular Cell Adhesion Molecule 1)	Rat	Zhang K, Wang D, Jiang Z, et al. Endothelial Dysfunction Correlates With Decompression Bubbles In Rats[J]. Scientific Reports, 2016.



Cat. No.	Targets	Species	Publications
E-EL-R0167	ET-1(Endothelin 1)	Rat	Kurt A, Tumkaya L, Kalkan Y, et al. Is Adalimumab Protective In Ischemia-Reperfusion Injury In Lung[J]. Iranian Journal of Basic Medical Sciences, 2015, 18(11): 1093-1099.
E-EL-R0167	ET-1(Endothelin 1)	Rat	Damiani E, Ince C, Orlando F, et al. Effects Of The Infusion Of 4% Or 20% Human Serum Albumin On The Skeletal Muscle Microcirculation In Endotoxemic Rats[J]. Plos One, 2016.
E-EL-R0167	ET-1(Endothelin 1)	Rat	Kurt A, Tumkaya L, Turut H, et al. Protective Effects Of Infliximab On Lung Injury Induced By Methotrexate[J]. Archivos De Bronconeumología, 2015, 51(11): 551-557.
E-EL-R0151 E-EL-R0019 E-EL-R0012 E-EL-R0015 E-EL-R0608 E-EL-R1430	cTnT/TNNT2(Troponin T Type 2, Cardiac) TNF- α (Tumor Necrosis Factor Alpha) IL-1 β (Interleukin 1 Beta) IL-6(Interleukin 6) MIF(Macrophage Migration Inhibitory Factor) ANG II (Angiotensin II)	Rat	Zhou Q, Pan X, Wang L, et al. The Protective Role Of Neuregulin-1: A Potential Therapy For Sepsis-Induced Cardiomyopathy[J]. European Journal of Pharmacology, 2016, 788: 234-240.
E-EL-R0022 E-EL-R0015 E-EL-R0019	CRP(C-Reactive Protein) IL-6(Interleukin 6) TNF- α (Tumor Necrosis Factor Alpha)	Rat	Zhang Y, Imam M U, Ismail M, et al. Edible Bird'S Nest Attenuates High Fat Diet-Induced Oxidative Stress And Inflammation Via Regulation Of Hepatic Antioxidant And Inflammatory Genes[J]. BMC Complementary and Alternative Medicine, 2015, 15: 310-316.
E-EL-R0014 E-EL-R0016 E-EL-R0022	IL-4(Interleukin 4) IL-10(Interleukin 10) CRP(C-Reactive Protein)	Rat	Barquilla P C, Pagano E S, Jiménez-Ortega V, et al. Melatonin Normalizes Clinical And Biochemical Parameters Of Mild Inflammation In Diet-Induced Metabolic Syndrome In Rats[J]. Journal of Pineal Research, 2014, 57(3): 280-90.



Cat. No.	Targets	Species	Publications
E-EL-0065 E-EL-0072 E-EL-R0391 E-EL-R0026 E-EL-R0582	E2(Estradiol) T(Testosterone) FSH(Follicle Stimulating Hormone) LH(Luteinizing Hormone) LEP(Leptin)	Rat	Mu Y, Yan W J, Yin T L, et al. Curcumin Ameliorates High-Fat Diet-Induced Spermatogenesis Dysfunction[J]. Molecular Medicine Reports, 2016, 14(4): 3588-3594.
E-EL-M2614 E-EL-M0002 E-EL-M0039	INS(Insulin) ADP/Acrp30(Adiponectin) LEP(Leptin)	Mouse	Chirag A. Patel, Jalaram H. Thakkar, Priyanka Patel, et al. Lipase Inhibitory And Anti-Obesity Effect Of The Hydro Alcoholic Extract Of Pinus Sylvestris Bark[J]. World Journal of Pharmacy and Pharmaceutical Sciences, 2015, 4(4): 1213-1229.
E-EL-M0780 E-EL-M0627 E-EL-M2442	MMP-2(Matrix Metalloproteinase 2) MMP-9(Matrix Metalloproteinase-9) EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor)	Mouse	Dai X, Ji Y, Jiang P, et al. Marsdenia Tenacissima Extract Suppresses Tumor Growth And Angiogenesis In A20 Mouse Lymphoma[J]. Oncology Letters, 2017, 13(5): 2897-2902.
E-EL-M0644	RANKL(Receptor Activator of Nuclear Factor Kappa B Ligand)	Mouse	Sartori M, Vincenzi F, Ravani A, et al. RAW 264.7 Co-Cultured With Ultra-High Molecular Weight Polyethylene Particles Spontaneously Differentiate Into Osteoclasts: An In Vitro Model Of Periprosthetic Osteolysis[J]. Journal of Biomedical Materials Research Part A, 2017, 105(2): 510-520.
E-EL-M0389 E-EL-M1272	Cys-C(Cystatin C) ALT(Alanine Transaminase)	Mouse	Wu Q, Huang L, Li Z, et al. The Potential Application Of Raw Cadmium Sulfide Nanoparticles As Ct Photographic Developer[J]. Nanoscale Research Letters, 2016, 11: 232-239.
E-EL-M0204	BNP(Brain Natriuretic Peptide)	Mouse	Diané A, Nikolic N, Rudecki A P, et al. Pacap Is Essential For The Adaptive Thermogenic Response Of Brown Adipose Tissue To Cold Exposure[J]. Journal of Endocrinology, 2014, 222(3): 327-339.



Cat. No.	Targets	Species	Publications
E-EL-M0037 E-EL-M0730 E-EL-M0077 E-EL-M0647	IL-1 β (Interleukin 1 Beta) IL-18(Interleukin 18) sICAM-1/CD54(Soluble Intercellular Adhesion Molecule 1) sVCAM-1/CD106(soluble Vascular Cell Adhesion Molecule 1)	Mouse	Zhang J L, Xia L Y, Zhang F, et al. A Novel Mechanism Of Diabetic Vascular Endothelial Dysfunction: Hypoadiponectinemia-Induced NLRP3 Inflammasome Activation[J]. Biochimica Et Biophysica Acta-Molecular Basis of Disease, 2017.
E-EL-M0018 E-EL-M0267	GRO α /CXCL1(Growth Regulated Oncogene Alpha) CX3CL1(Chemokine C-X3-C-Motif Ligand 1)	Mouse	Wang H, Hong L J, Huang J Y, et al. P2RX7 Sensitizes Mac-1/ICAM-1-Dependent Leukocyte-Endothelial Adhesion And Promotes Neurovascular Injury During Septic Encephalopathy[J]. Cell Research, 2015, 25(6): 674-690.
E-EL-H0125	ApoA1(Apolipoprotein A1)	Human	Reddy M N, Dungdung R, Pilankatta R. Proteomic Analysis Reveals The Enhancement Of Human Serum Apolipoprotein A-1 (Apo A-1) In Individuals Infected With Multiple Dengue Virus Serotypes[J]. Tropical Medicine & International Health, 2017.
E-EL-H0483	bFGF/FGF2(Basic Fibroblast Growth Factor)	Human	Esnafoğlu E, Ayyıldız S N. Decreased Levels Of Serum Fibroblast Growth Factor-2 In Children With Autism Spectrum Disorder[J]. Psychiatry Research, 2017.
E-EL-H0113 E-EL-H1763 E-EL-H1705	LEP(Leptin) VF(Visfatin) TGF- β (Transforming Growth Factor Beta)	Human	Partan R U, Hidayat R, Mukti M. The Relationships Among Leptin, Resistin, Visfatin, Transforming Growth Factor (TGF)- β Level And The Severity Grade Of Osteoarthritis In Osteoarthritis Patients With Obesity[J]. Journal of Research in Medical and Dental Science, 2017, 5(2): 10-14.
E-EL-H0698	CHEM(Chemerin)	Human	Gashlan H M. Relationship Between Levels Of Retinol Binding Protein 4, Vaspın And Chemerin And Insulin Resistance In Gestational Diabetes Mellitus[J]. International Journal of Pharmaceutical Research & Allied Sciences, 2017, 6(1): 236-250.



Cat. No.	Targets	Species	Publications
E-EL-H0113	LEP(Leptin)	Human	Yeboah F A, Ngala R A, Bawah A T, et al. Adiposity And Hyperleptinemia During The First Trimester Among Pregnant Women With Preeclampsia[J]. International Journal of Women's Health, 2017, 9: 449-454.
E-EL-H0113 E-EL-H0093	LEP(Leptin) LEPR(Leptin Receptor)	Human	Tu X, Yu C, Gao M, et al. Lepr Gene Polymorphism And Plasma Soluble Leptin Receptor Levels Are Associated With Polycystic Ovary Syndrome In Han Chinese Women[J]. Personalized Medicine, 2017, 14(4): 299-307.
E-EL-H0048 E-EL-H0109 E-EL-H2172 E-EL-H1640 E-EL-R0580 E-EL-R0504 E-EL-R0019 E-EL-R0488 E-EL-R0902	IL-8(Interleukin 8) TNF- α (Tumor Necrosis Factor Alpha) HO1(Heme Oxygenase 1) PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) LRP-1(Low-Density Lipoprotein-Receptor-Related Protein) HDL(High Density Lipoprotein) TNF- α (Tumor Necrosis Factor Alpha) HO1(Heme Oxygenase 1) sPECAM-1(Soluble Platelet Endothelial Cell Adhesion Molecule 1)	Human&Rat	Yan Z, Wang W, Wu Y, et al. Zinc Oxide Nanoparticle-Induced Atherosclerotic Alterations In Vitro And In Vivo[J]. International Journal of Nanomedicine, 2017, 12: 4433-4442.
E-EL-R0167 E-EL-R0361 E-EL-R0039 E-EL-R0236	ET-1(Endothelin 1) EG-VEGF(Endocrine Gland Derived Vascular-Endothelial Growth Factor) LN(Laminin) COL4(Collagen Type IV)	Rat	Li G, Peng Y, Zhao T, et al. Plumbagin Alleviates Capillarization Of Hepatic Sinusoids In Vitro By Downregulating ET-1, VEGF, LN, And Type IV Collagen[J]. Biomed Research International, 2017.



Cat. No.	Targets	Species	Publications
E-EL-H2193	FG(Fibrinogen)	Human	Yılmaz N, Yılmaz M, Sirin B, et al. The Relationship Between Levels Of Plasma-Soluble Urokinase Plasminogen Activator Receptor (suPAR) And Presence Of Migraine Attack And Aura[J]. Journal of Receptors and Signal Transduction, 2017.
E-EL-H1919	GHRL(Ghrelin)	Human	Wu J, Wang K, Xu J, et al. Associations Between Serum Ghrelin And Knee Symptoms, Joint Structures And Cartilage Or Bone Biomarkers In Patients With Knee Osteoarthritis[J]. Osteoarthritis and Cartilage, 2017.
E-EL-H0654	COMP(Cartilage Oligomeric Matrix Protein)		
E-EL-H0835	CTX I (Cross Linked C-telopeptide of Type I Collagen)		
E-EL-H0836	NTX I (Cross Linked N-telopeptide of Type I Collagen)		
E-EL-H0183	P III NP(N-Terminal Procollagen III Propeptide)		
E-EL-H0132	MMP-3(Matrix Metalloproteinase 3)		
E-EL-H0134	MMP-10(Matrix Metalloproteinase 10) MMP-13(Matrix Metalloproteinase 13)		







Our Worldwide Distributors





Brochure of Elisa kits & Antibodies



A microscopic view of several red blood cells (erythrocytes) against a dark red, blurred background. The cells are biconcave discs, appearing as reddish-brown, oval shapes with a lighter center. They are scattered across the frame, with some in sharp focus and others blurred in the background.

Tel: 1-240-252-7368 Fax: 1-240-252-7376
Website: www.elabscience.com
Email: orders@elabscience.com
Address: 14780 memorial drive, suite 216, houston, TX77079