

Tel:240-252-7368(USA) Fax: 240-252-7376(USA) techsupport@elabscience.com Website: www.elabscience.com

# Biotin Anti-Human/Mouse/Rat CD47 Antibody[MIAP410]

Catalog No. E-AB-F1016B Reactivity Human, Mouse, Rat

Storage Store at 2~8°C, Avoid freeze / thaw cycles Applications FCM

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## **Antigen Information**

Alternate Names Leukocyte surface antigen CD47,Cd47,Integrin-associated protein,IAP

**Uniprot ID** Q08722,Q61735,P97829

Background CD47, also known as Integrin-Associated Protein (IAP), is a membrane protein of about 50 kD

with an IgV-like extracelluluar domain, a five membrane-spanning segment and a short terminal cytoplasmic region. It is widely expressed on many cell types and often associated with beta 3 integrins. It has been reported that CD47 functions as a self marker. Red cells lacking CD47 were rapidly cleared from the bloodstream by splenic macrophages. By binding to SIRPα, CD47 controls hemostatic innate immune functions, such as phagocytosis and cell trafficking.

#### **Product Details**

 Form
 Liquid

 Concentration
 0.5 mg/mL

 Size
 25µg/100µg

 Clone No.
 MIAP410

 Host
 Mouse

 $\begin{tabular}{lll} \textbf{Isotype} & Mouse IgG1, \kappa \\ \textbf{Reactivity} & Human, Mouse, Rat \\ \end{tabular}$ 

**Application** FCM

**Isotype Control**Biotin Mouse IgG1, κ Isotype Control[MOPC-21] [Product E-AB-F09793B]

**Storage Buffer** Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.

Shipping
Biological ice pack at 4 °C
Stability & Storage
Keep as concentrated solution.
Store at 2~8°C .Do not freeze.

This product is guaranteed up to one year from purchase.



Tel:240-252-7368(USA) Fax: 240-252-7376(USA) techsupport@elabscience.com Website: www.elabscience.com

### Recommended usage

Each lot of this antibody is quality control tested by flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is  $\leq 1.0 \,\mu g$  per  $10^6$  cells in 100  $\mu L$  volume or 100  $\mu L$  of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

#### **Related Information**

- 1. Sample Preparation for Flow Cytometry <a href="https://www.elabscience.com/List-detail-5594.html">https://www.elabscience.com/List-detail-5594.html</a>
- 2. Staining Cell Surface Targets for Flow Cytometry <a href="https://www.elabscience.com/List-detail-5568.html">https://www.elabscience.com/List-detail-5568.html</a>
- 3. Flow Cytometry Troubleshooting Tips <a href="https://www.elabscience.com/List-detail-5593.html">https://www.elabscience.com/List-detail-5593.html</a>
- 4. How to select the appropriate detection channel through the spectrogram? https://www.elabscience.com/Listdetail-459742.html