# **Elabscience**®

## APC/Cyanine7 Rat IgM, κ Isotype Control[RTK2118]

Catalog No.E-AB-F09772NStorageStore at 2~8°C, Avoid freeze / thaw cycles

**Important Note:** Centrifuge before opening to ensure complete recovery of vial contents.

#### **Product Details**

| Form                | Liquid   |
|---------------------|--|
| Size                | 20 Tests/50 Tests/100 Tests/200 Tests  |
| Clone No.           | RTK2118  |
| Host                | Rat  |
| Isotype             | IgM, κ   |
| Application         | FCM  |
| Storage Buffer      | PBS with 0.05% Proclin300, 1% BSA  |
| Shipping            | Biological ice pack at 4 °C  |
| Stability & Storage | Keep as concentrated solution.   |
|                     | Store at 2~8°C and protected from prolonged exposure to light.Do not freeze. |
|                     | This product is guaranteed up to one year from purchase.                     |

### Fluorophore

Conjugation: APC/Cyanine7

### **Recommended usage**

Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. The amount of the reagent is suggested to be used 2  $\mu$ L or 5  $\mu$ L of antibody per test (million cells in 100  $\mu$ l staining volume or per 100  $\mu$ l of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

### **Related Information**

- 1. Sample Preparation for Flow Cytometry <u>https://www.elabscience.com/List-detail-5594.html</u>
- 2. Staining Intracellular Antigens for Flow Cytometry https://www.elabscience.com/List-detail-5570.html
- 3. Flow Cytometry Troubleshooting Tips https://www.elabscience.com/List-detail-5593.html
- 4. How to select the appropriate detection channel through the spectrogram? <u>https://www.elabscience.com/List-detail-459742.html</u>

For Research Use Only Thank you for your recent purchase

If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications:Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA-Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFCM: Intracellular Staining for Flow Cytometry; WB-Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation