Elabscience®

Anti-Ebola virus EBOV(subtype Bundibugyo, strain Uganda 2007) GP1/Glycoprotein Polyclonal Antibody

E-AB-V1104

| Application | WB | Host | Rabbit |
|-------------|---|------|--------|
| Storage | Store at -20°C. Avoid freeze / thaw cycles. | | |

Important Note Centrifuge before opening to ensure complete recovery of vial contents.

Product Details

| Immunogen | Recombinant EBOV (subtype Bundibugyo, strain Uganda 2007) GP1 / Glycoprotein Protein | |
|---------------------|--|--|
| | (His Tag) | |
| Isotype | IgG | |
| Host | Rabbit | |
| Reactivity | Ebola virus | |
| Dilution | WB 1:1000-1:5000 | |
| Storage Buffer | 0.2 μm filtered solution in PBS | |
| Stability & Storage | Ships on ice packs. Store at -20°C | |
| Description | This antibody was produced in rabbits immunized with purified Recombinant EBOV (subtype | |
| | Bundibugyo, strain Uganda 2007) GP1 / Glycoprotein Protein (His Tag). And the antibody was | |
| | purified by Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) GP1 / | |

Antigen Infomation

Background

Alternate Names Glycoprotein,GP

The fourth gene of the EBOV genome encodes a 16-kDa envelope-attached glycoprotein (GP) and a 11 kDa secreted glycoprotein (sGP). Both GP and sGP have an identical 295-residue Nterminus, however, they have different C-terminal sequences. Recently, great attention has been paid to GP for vaccines design and entry inhibitors isolation. GP is a class I fusion protein which assembles as trimers on viral surface and plays an important role in virus entry and attachment. Mature GP is a disulfide-linked heterodimer formed by two subunits, GP1 and GP2, which are generated from the proteolytical process of GP precursor (pre-GP) by cellular furin during virus assembly. The GP1 subunit contains a mucin domain and a receptor-binding domain (RBD); the GP2 subunit has a fusion peptide, a helical heptad-repeat (HR) region, a transmembrane (TM) domain, and a 4-residue cytoplasmic tail. The RBD of GP1 mediates the interaction of EBOV with cellular receptor (e.g. DC-SIGN/LSIGN, TIM-1, hMGL, NPC1, β-integrins, folate receptor- α , and Tyro3 family receptors), of which TIM1 and NPC1 are essential for EBOV entry; the mucin domain having N- and O-linked glycans enhances the viral attachment to cellular hMGL, and participates in shielding key neutralization epitopes, which helps the virus evades immune elimination. There are large conformation changes of GP2 during membrane fusion, which enhance the insertion of fusion loop into cellular membrane and facilitate the release of viral nucleocapsid core to cytoplasm.


For Research Use Only Thank you for your recent purchase. If you would like to learn more about antibodies,please visit www.elabscience.com. Focus on your research Service for life science

Applications:Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA-Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFCM: Intracellular Staining for Flow Cytometry; WB-Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation

Elabscience®

Tel:240-252-7368(USA) Fax: 240-252-7376(USA) techsupport@elabscience.com Website: www.elabscience.com

Images



Western Blot analysis of Recombinant EBOV (subtype Bundibugyo, strain Uganda 2007) GP1 / Glycoprotein Protein (His Tag)(PKSV030136 with 50ng and 10ng) using Anti-Ebola virus EBOV(subtype Bundibugyo, strain Uganda 2007) GP1/Glycoprotein Polyclonal Antibody at dilution of 1:1000.

For Research Use Only

Thank you for your recent purchase. If you would like to learn more about antibodies, please visit www.elabscience.com. Focus on your research Service for life science

Applications:Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA-Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFCM: Intracellular Staining for Flow Cytometry; WB-Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation