

Cat. No.	Product Name	Reactivity	Size	Detection Range	Price
E-EL-H0111	VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$495
E-EL-H1763	VF (Visfatin) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H2168	VWF (Von Willebrand Factor) ELISA Kit	Human	96 Tests	1.56-100 ng/mL	\$609
E-EL-H5560	Zonulin ELISA Kit	Human	96 Tests	0.78-50 ng/mL	\$609
E-EL-M0079	ACTH (Adrenocorticotropic Hormone) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$495
E-EL-M0002	ADP/Acrp30 (Adiponectin) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$495
E-EL-M0337	C3a (Complement Component 3a) ELISA Kit	Mouse	96 Tests	0.16-10 ng/mL	\$609
E-EL-M0355	CKMB (Creatine Kinase MB Isoenzyme) ELISA Kit	Mouse	96 Tests	78.13-5000 pg/mL	\$609
E-EL-M0325	COL1 (Collagen Type I) ELISA Kit	Mouse	96 Tests	0.31-20 ng/mL	\$609
E-EL-0161	CORT (Corticosterone) ELISA Kit	Mouse	96 Tests	3.12-200 ng/mL	\$495
E-EL-M0053	CRP (C-Reactive Protein) ELISA Kit	Mouse	96 Tests	0.03-25 ng/mL	\$609
E-EL-M0340	CTGF (Connective Tissue Growth Factor) ELISA Kit	Mouse	96 Tests	0.16-10 ng/mL	\$609
E-EL-M3023	CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	Mouse	96 Tests	31.25-2000 pg/mL	\$609
E-EL-M0511	FSH (Follicle Stimulating Hormone) ELISA Kit	Mouse	96 Tests	1.56-100 ng/mL	\$495
E-EL-M0676	HMGB-1 (High Mobility Group Protein B1) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$609
E-EL-M0033	IFN-β (Interferon Beta) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$495
E-EL-M0048	IFN-γ (Interferon Gamma) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$495
E-EL-M0690	IgA (Immunoglobulin A) ELISA Kit	Mouse	96 Tests	6.25-400 ng/mL	\$609
E-EL-M3034	IgE (Immunoglobulin E) ELISA Kit	Mouse	96 Tests	125-8000 ng/mL	\$495
E-EL-M0692	IgG (Immunoglobulin G) ELISA Kit	Mouse	96 Tests	1.56-100 ng/mL	\$495
E-EL-M3036	IgM (Immunoglobulin M) ELISA Kit	Mouse	96 Tests	3.13-200 ng/mL	\$495
E-EL-M0037	IL-1β (Interleukin 1 Beta) ELISA Kit	Mouse	96 Tests	7.81-500 pg/mL	\$495
E-EL-M0042	IL-2 (Interleukin 2) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$495
E-EL-M0043	IL-4 (Interleukin 4) ELISA Kit	Mouse	96 Tests	31.25-2000 pg/mL	\$495
E-EL-M0722	IL-5 (Interleukin 5) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$609
E-EL-M0044	IL-6 (Interleukin 6) ELISA Kit	Mouse	96 Tests	31.25-2000 pg/mL	\$495
E-EL-M0046	IL-10 (Interleukin 10) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$495
E-EL-M3062	IL-12 (Interleukin 12) ELISA Kit	Mouse	96 Tests	0.31-20 pg/mL	\$495
E-EL-M0047	IL-17A (Interleukin 17A) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$495
E-EL-M0730	IL-18 (Interleukin 18) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$609
E-EL-M0731	IL-23 (Interleukin 23) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$609
E-EL-M1382	INS (Insulin) ELISA Kit	Mouse	96 Tests	3.13-200 pg/mL	\$609
E-EL-M2614	INS (Insulin) ELISA Kit	Mouse	96 Tests	0.31-20 ng/mL	\$609
E-EL-M3008	LEP (Leptin) ELISA Kit	Mouse	96 Tests	0.31-20 ng/mL	\$495
E-EL-M3053	LH (Luteinizing Hormone) ELISA Kit	Mouse	96 Tests	0.31-20 ng/mL	\$495
E-EL-M0746	LTF/LF (Lactoferrin) ELISA Kit	Mouse	96 Tests	62.5-4000 pg/mL	\$495
E-EL-M0792	MAU (Microalbuminuria) ELISA Kit	Mouse	96 Tests	0.16-10 ng/mL	\$495
E-EL-M3001	MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	Mouse	96 Tests	31.25-2000 pg/mL	\$609
E-EL-M0788	MT (Melatonin) ELISA Kit	Mouse	96 Tests	7.81-500 pg/mL	\$495
E-EL-M0828	NGAL (Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	Mouse	96 Tests	39.06-2500 pg/mL	\$609
E-EL-M0834	NT-proBNP (N-terminal pro-Brain Natriuretic Peptide) ELISA Kit	Mouse	96 Tests	15.63-1000 pg/mL	\$609
E-EL-M0864	OC/BGP (Osteocalcin) ELISA Kit	Mouse	96 Tests	0.39-25 ng/mL	\$495
E-EL-M0233	P I NP (Procollagen I N-Terminal Propeptide) ELISA Kit	Mouse	96 Tests	93.75-6000 pg/mL	\$609
E-EL-M1040	sIgA (Secretory Immunoglobulin A) ELISA Kit	Mouse	96 Tests	0.16-10 ng/mL	\$609
E-EL-M3063	TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit	Mouse	96 Tests	7.81-500 pg/mL	\$495

The Products are for Research Use Only

Cat. No.	Product Name	Reactivity	Size	Detection Range	Price
E-EL-M1203	TNNI3/cTn-I (Troponin I Type 3, Cardiac) ELISA Kit	Mouse	96 Tests	7.81-500 pg/mL	\$495
E-EL-M1292	VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit	Mouse	96 Tests	31.25-2000 pg/mL	\$609
E-EL-M1247	VWF (Von Willebrand Factor) ELISA Kit	Mouse	96 Tests	1.25-80 ng/mL	\$609
E-EL-R0048	ACTH (Adrenocorticotropic Hormone) ELISA Kit	Rat	96 Tests	15.63-1000 pg/mL	\$495
E-EL-R3012	ADP/Acrp30 (Adiponectin) ELISA Kit	Rat	96 Tests	1.56-100 ng/mL	\$495
E-EL-R1235	BDNF (Brain Derived Neurotrophic Factor) ELISA Kit	Rat	96 Tests	31.25-2000 pg/mL	\$495
E-EL-R0255	C3a (Complement Component 3a) ELISA Kit	Rat	96 Tests	0.63-40 ng/mL	\$609
E-EL-0160	CORT (Corticosterone) ELISA Kit	Rat/Chicken	96 Tests	0.39-25 ng/mL	\$495
E-EL-R0506	CRP (C-Reactive Protein) ELISA Kit	Rat	96 Tests	0.16-10 ng/mL	\$495
E-EL-0152	E2 (Estradiol) ELISA Kit	Rat/Porcine	96 Tests	1.56-100 pg/mL	\$495
E-EL-R0391	FSH (Follicle Stimulating Hormone) ELISA Kit	Rat	96 Tests	3.13-200 ng/mL	\$495
E-EL-R0009	IFN-γ (Interferon Gamma) ELISA Kit	Rat	96 Tests	31.25-2000 pg/mL	\$495
E-EL-R3001	IGF-1 (Insulin-like Growth Factor 1) ELISA Kit	Rat	96 Tests	3.13-200 ng/mL	\$609
E-EL-R0518	IgG (Immunoglobulin G) ELISA Kit	Rat	96 Tests	3.13-200 ng/mL	\$495
E-EL-R0012	IL-1β (Interleukin 1 Beta) ELISA Kit	Rat	96 Tests	31.25-2000 pg/mL	\$495
E-EL-R0013	IL-2 (Interleukin 2) ELISA Kit	Rat	96 Tests	31.25-2000 pg/mL	\$495
E-EL-R0014	IL-4 (Interleukin 4) ELISA Kit	Rat	96 Tests	15.63-1000 pg/mL	\$495
E-EL-R0015	IL-6 (Interleukin 6) ELISA Kit	Rat	96 Tests	12.5-800 pg/mL	\$495
E-EL-R0016	IL-10 (Interleukin 10) ELISA Kit	Rat	96 Tests	31.25-2000 pg/mL	\$495
E-EL-R0567	IL-18 (Interleukin 18) ELISA Kit	Rat	96 Tests	15.63-1000 pg/mL	\$609
E-EL-R3034	INS (Insulin) ELISA Kit	Rat	96 Tests	6.25-400 pg/mL	\$609
E-EL-R0582	LEP (Leptin) ELISA Kit	Rat	96 Tests	0.16-10 ng/mL	\$495
E-EL-R0026	LH (Luteinizing Hormone) ELISA Kit	Rat	96 Tests	1.56-100 mIU/mL	\$495
E-EL-R0633	MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	Rat	96 Tests	0.16-10 ng/mL	\$609
E-EL-R0031	MT (Melatonin) ELISA Kit	Rat	96 Tests	15.63-1000 pg/mL	\$495
E-EL-R3006	PRL (Prolactin) ELISA Kit	Rat	96 Tests	3.13-200 ng/mL	\$609
E-EL-R2856	TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit	Rat	96 Tests	15.63-1000 pg/mL	\$495
E-EL-R1253	TNNI3/cTn-I (Troponin I Type 3, Cardiac) ELISA Kit	Rat	96 Tests	31.25-2000 pg/mL	\$495
E-EL-R0976	TSH (Thyroid Stimulating Hormone) ELISA Kit	Rat	96 Tests	1.25-80 ng/mL	\$495
E-EL-R1079	VWF (Von Willebrand Factor) ELISA Kit	Rat	96 Tests	0.16-10 ng/mL	\$609



Elabscience Bionovation Inc.

 Toll-free : 1-888-852-8623  Tel : 1-832-243-6086  Fax : 1-832-243-6017
 Web : www.elabscience.com  Email : orders@elabscience.com techsupport@elabscience.com

EFO98-2410v2

Elabscience®

Featured ELISA Products

A Reliable Research Partner in Life Science and Medicine



www.elabscience.com

Elabscience Bionovation Inc.



Elabscience® ELISA Kits have been strictly quality-controlled to ensure the accuracy and consistency of results.



Elabscience® ELISA Kits’ good performance and reproducibility brought us well reputation as a trustworthy supplier.



CiteAb recognized Elabscience® as the "ELISA KIT SUPPLIER TO WATCH IN 2019".



Elabscience® ELISA products have been cited in over 10,000 SCI papers, and featured in top-tier journals like *Nature Medicine, Nature, Cell, Immunity, Signal Transduction and Targeted Therapy*, etc.

Featured ELISA Kits

Cat. No.	Product Name	Reactivity	Size	Detection Range	Price
E-EL-0028	8-OHdG (8-Hydroxydeoxyguanosine) ELISA Kit	Universal	96 Tests	1.56-100 ng/mL	\$495
E-EL-0015	25-HVD3 (25-Hydroxy Vitamin D3) ELISA Kit	Universal	96 Tests	3.13-200 ng/mL	\$495
E-EL-0125	A-CoA (Acetyl Coenzyme A) ELISA Kit	Universal	96 Tests	0.31-20 ng/mL	\$495
E-EL-0070	ALD (Aldosterone) ELISA Kit	Universal	96 Tests	31.25-2000 pg/mL	\$495
E-EL-0056	Camp (Cyclic adenosine monophosphate) ELISA Kit	Universal	96 Tests	1.56-100 ng/mL	\$495
E-EL-0046	DA (Dopamine) ELISA Kit	Universal	96 Tests	31.25-2000 pg/mL	\$495
E-EL-0031	DHT (Dihydrotestosterone) ELISA Kit	Universal	96 Tests	39.06-2500 pg/mL	\$495
E-EL-0016	DHVD3 (1,25-Dihydroxyvitamin D3) ELISA Kit	Universal	96 Tests	7.81-500 pg/mL	\$495
E-EL-0009	FA/VB9 (Folic Acid/Vitamin B9) ELISA Kit	Universal	96 Tests	1.56-100 ng/mL	\$495
E-EL-0079	ft3 (Free Triiodothyronine) ELISA Kit	Universal	96 Tests	1.56-100 pg/mL	\$495
E-EL-0122	ft4 (Free Thyroxine) ELISA Kit	Universal	96 Tests	1.56-100 pg/mL	\$495
E-EL-0154	Pg (Progesterone) ELISA Kit	Universal	96 Tests	0.31-20 ng/mL	\$495
E-EL-0033	ST/5-HT (Serotonin/5-Hydroxytryptamine) ELISA Kit	Universal	96 Tests	15.63-1000 ng/mL	\$495
E-EL-0155	T (Testosterone) ELISA Kit	Universal	96 Tests	0.31-20 ng/mL	\$495
E-EL-0162	TGF-β1 (Transforming Growth Factor Beta 1) ELISA Kit	Universal	96 Tests	0.16-10 ng/mL	\$495
E-EL-H0281	ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit	Human	96 Tests	0.39-25 ng/mL	\$609
E-EL-H0137	ACTH (Adrenocorticotrophic Hormone) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$495
E-EL-H6084	ApoB48 (Apolipoprotein B48) ELISA Kit	Human	96 Tests	1.56-100 ng/mL	\$609
E-EL-H0470	ApoE (Apolipoprotein E) ELISA Kit	Human	96 Tests	23.44-1500 ng/mL	\$495
E-EL-H0542	Aβ1-40 (Amyloid Beta 1-40) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$495
E-EL-H0543	Aβ1-42 (Amyloid Beta 1-42) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$495
E-EL-H0010	BDNF (Brain Derived Neurotrophic Factor) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$495
E-EL-H6042	bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$609
E-EL-H0818	C3a (Complement Component 3a) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H0190	C5a (Complement Component 5a) ELISA Kit	Human	96 Tests	78.13-5000 pg/mL	\$609
E-EL-H2357	CALP (Calprotectin) ELISA Kit	Human	96 Tests	1.56-100 ng/mL	\$609
E-EL-H0698	CHEM (Chemerin) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H0869	COL1α1 (Collagen Type I Alpha 1) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H0781	COL7 (Collagen Type VII) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$495
E-EL-0157	Cortisol ELISA Kit	Human	96 Tests	6.25-400 ng/mL	\$495
E-EL-H0043	CRP (C-Reactive Protein) ELISA Kit	Human	96 Tests	0.39-25 ng/mL	\$495
E-EL-H0835	CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609

The Products are for Research Use Only

Cat. No.	Product Name	Reactivity	Size	Detection Range	Price
E-EL-H3643	Cys-C (Cystatin C) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$495
E-EL-H6145	D2D (D-Dimer) ELISA Kit	Human	96 Tests	6.25-400 ng/mL	\$609
E-EL-H0996	DEFβ2/DEFB2 (Defensin Beta 2) ELISA Kit	Human	96 Tests	62.5-4000 pg/mL	\$609
E-EL-Q150	E2 (Estradiol) ELISA Kit	Human/Mouse/Monkey	96 Tests	9.37-600 pg/mL	\$495
E-EL-H0059	EGF (Epidermal Growth Factor) ELISA Kit	Human	96 Tests	3.91-250 pg/mL	\$495
E-EL-H3640	EPO (Erythropoietin) ELISA Kit	Human	96 Tests	1.56-100 mIU/mL	\$495
E-EL-H1557	ESM1 (Endothelial Cell Specific Molecule 1) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H0064	ET-1 (Endothelin 1) ELISA Kit	Human	96 Tests	1.25-80 pg/mL	\$495
E-EL-H0760	F9 (Coagulation Factor IX) ELISA Kit	Human	96 Tests	1.25-80 ng/mL	\$495
E-EL-H1745	Fcy (Fc Fragment of IgG) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H0168	FE (Ferritin) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$495
E-EL-H0386	FETUA (Fetuin A) ELISA Kit	Human	96 Tests	9.38-600 ng/mL	\$609
E-EL-H0074	FGF21 (Fibroblast Growth Factor 21) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$609
E-EL-H1116	FGF23 (Fibroblast Growth Factor 23) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H0075	Flt3L (FMS-like Tyrosine Kinase 3 Ligand) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H0179	FN (Fibronectin) ELISA Kit	Human	96 Tests	1.56-100 ng/mL	\$609
E-EL-H1470	GAL3 (Galectin 3) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H0080	GDF15 (Growth Differentiation Factor 15) ELISA Kit	Human	96 Tests	23.44-1500 pg/mL	\$609
E-EL-H6093	GFAP (Glial Fibrillary Acidic Protein) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H0177	GH (Growth Hormone) ELISA Kit	Human	96 Tests	78.13-5000 pg/mL	\$495
E-EL-H6025	GLP-1 (Glucagon Like Peptide 1) ELISA Kit	Human	96 Tests	1.56-100 pg/mL	\$495
E-EL-H1617	GzmB (Granzyme B) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H6080	HBsAg (Hepatitis B surface antigen) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H0084	HGF (Hepatocyte Growth Factor) ELISA Kit	Human	96 Tests	125-8000 pg/mL	\$495
E-EL-H6066	HIF-1α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	Human	96 Tests	62.5-4000 pg/mL	\$495
E-EL-H1554	HMGB-1 (High Mobility Group Protein B1) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$609
E-EL-H2364	HS (Heparan Sulfate) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H5134	hs-CRP (high-sensitivity C-Reactive Protein) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$495
E-EL-H1863	HSP-70/HSPA9 (Heat Shock 70 kDa Protein 9) ELISA Kit	Human	96 Tests	0.78-50 ng/mL	\$495
E-EL-H6114	ICAM-1/CD54 (intercellular adhesion molecule 1) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H2162	IDO (Indoleamine-2,3-Dioxygenase) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H0159	IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H0085	IFN-β (Interferon Beta) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$495
E-EL-H0108	IFN-γ (Interferon Gamma) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$495
E-EL-H0086	IGF-1 (Insulin-like Growth Factor 1) ELISA Kit	Human	96 Tests	1.56-100 ng/mL	\$609
E-EL-H0169	IgG (Immunoglobulin G) ELISA Kit	Human	96 Tests	1.56-100 ng/mL	\$495
E-EL-H0088	IL-1α (Interleukin 1 Alpha) ELISA Kit	Human	96 Tests	3.91-250 pg/mL	\$495
E-EL-H0149	IL-1β (Interleukin 1 Beta) ELISA Kit	Human	96 Tests	7.81-500 pg/mL	\$495
E-EL-H0099	IL-2 (Interleukin 2) ELISA Kit	Human	96 Tests	7.81-500 pg/mL	\$495
E-EL-H0101	IL-4 (Interleukin 4) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$495
E-EL-H6156	IL-6 (Interleukin 6) ELISA Kit	Human	96 Tests	1.56-100 pg/mL	\$495
E-EL-H6008	IL-8 (Interleukin 8) ELISA Kit	Human	96 Tests	7.81-500 pg/mL	\$495
E-EL-H6154	IL-10 (Interleukin 10)ELISA Kit	Human	96 Tests	1.56-100 pg/mL	\$495
E-EL-H0105	IL-17A (Interleukin 17A) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$495
E-EL-H0253	IL-18 (Interleukin 18) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H0107	IL-23 (Interleukin 23) ELISA Kit	Human	96 Tests	39.06-2500 pg/mL	\$609
E-EL-H2402	IL-33 (Interleukin 33) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609

The Products are for Research Use Only

Cat. No.	Product Name	Reactivity	Size	Detection Range	Price
E-EL-H2443	IL-35 (Interleukin 35) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H0313	INHb (Inhibin B) ELISA Kit	Human	96 Tests	7.81-500 pg/mL	\$609
E-EL-H2665	INS (Insulin) ELISA Kit	Human	96 Tests	0.78-50 μIU/mL	\$609
E-EL-H0050	IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit	Human	96 Tests	7.81-500 pg/mL	\$495
E-EL-H6029	KIM-1 (Kidney Injury Molecule 1) ELISA Kit	Human	96 Tests	7.81-500 pg/mL	\$495
E-EL-H6017	LEP (Leptin) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$495
E-EL-H2438	LL-37 (Antibacterial Protein LL-37) ELISA Kit	Human	96 Tests	1.56-100 ng/mL	\$609
E-EL-H5200	LTF/LF (Lactoferrin) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$495
E-EL-H6005	MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	Human	96 Tests	62.5-4000 pg/mL	\$609
E-EL-H1446	MMP-3 (Matrix Metalloproteinase 3) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$495
E-EL-H1450	MMP-8 (Matrix Metalloproteinase 8) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$495
E-EL-H6075	MMP-9 (Matrix Metalloproteinase 9) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$495
E-EL-H2016	MT (Melatonin) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$495
E-EL-H0616	MUC1 (Mucin 1) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$495
E-EL-H2279	MUC5AC (Mucin 5 Subtype AC) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H1205	NGF (Nerve Growth Factor) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H1047	NSE (Neuron Specific Enolase) ELISA Kit	Human	96 Tests	2.34-150 ng/mL	\$495
E-EL-H6126	NT-proBNP (N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H1343	OC/BGP (Osteocalcin) ELISA Kit	Human	96 Tests	1.25-80 ng/mL	\$495
E-EL-H1341	OPG (Osteoprotegerin) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$495
E-EL-H1347	OPN (Osteopontin) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H6021	OxLDL (Oxidized Low Density Lipoprotein) ELISA Kit	Human	96 Tests	62.5-4000 pg/mL	\$495
E-EL-H0185	P I NP (Procollagen I N-Terminal Propeptide) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H0183	PIIINP (N-Terminal Procollagen III Propeptide) ELISA Kit	Human	96 Tests	23.44-1500 pg/mL	\$495
E-EL-H2199	PAF (Platelet Activating Factor) ELISA Kit	Human	96 Tests	78.13-5000 pg/mL	\$609
E-EL-H2104	PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H1492	PCT (Procalcitonin) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$495
E-EL-H1577	PDGF-BB (Platelet Derived Growth Factor BB) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$609
E-EL-H1547	PD-L1 (Programmed Cell Death Protein 1 Ligand 1) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H1555	PGF (Placental Growth Factor) ELISA Kit	Human	96 Tests	1.96-125 pg/mL	\$609
E-EL-H5442	PIVKA- II (Protein Induced by Vitamin K Absence or Antagonist- II) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H1558	Rac1 (Ras Related C3 Botulinum Toxin Substrate 1) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H1213	RETN (Resistin) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$609
E-EL-H1297	S100B (S100 Calcium Binding Protein B) ELISA Kit	Human	96 Tests	31.25-2000 pg/mL	\$495
E-EL-H2183	SAA (Serum Amyloid A) ELISA Kit	Human	96 Tests	1.25-80 ng/mL	\$609
E-EL-H2376	sC5b-9 (Soluble Terminal Complement Complex C5b-9) ELISA Kit	Human	96 Tests	31.25-2000 ng/mL	\$609
E-EL-H0052	SDF-1/CXCL12 (Stromal Cell Derived Factor 1) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H1275	slgA (Secretory Immunoglobulin A) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H0983	SNCα(Synuclein Alpha) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H1266	SP-A1 (Pulmonary surfactant-associated protein A1) ELISA Kit	Human	96 Tests	15.63-1000 pg/mL	\$609
E-EL-H6082	sST2 (Soluble ST2) ELISA Kit	Human	96 Tests	0.31-20 ng/mL	\$609
E-EL-H6085	sTfR1 (Soluble Transferrin Receptor 1) ELISA Kit	Human	96 Tests	6.25-400 ng/mL	\$609
E-EL-H2584	suPAR (Soluble Urokinase-type Plasminogen Activator Receptor) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$609
E-EL-H0184	TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	Human	96 Tests	0.16-10 ng/mL	\$495
E-EL-H0109	TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit	Human	96 Tests	7.81-500 pg/mL	\$495
E-EL-H5587	VCAM-1/CD106 (Vascular Cell Adhesion Molecule 1) ELISA Kit	Human	96 Tests	1.56-100 ng/mL	\$609

The Products are for Research Use Only