

Top 200 Best-selling ELISA Kits from Elabscience®

—Focus on hot research targets, contribute to your successful experiments!

Why Choose Elabscience® ELISA Kits?

👍 Strict Quality Control

- Elabscience® ensures the stability and high performance of its ELISA products through a rigorous three-step quality inspection of raw materials, semi-finished products, and finished products.

🏆 Authoritative Recognition

- CiteAb recognized Elabscience® as the "ELISA KIT SUPPLIER TO WATCH IN 2019".
- Elabscience® ELISA products have been cited in over 10,000 SCI articles, featured in top-tier journals like *Nature Medicine*, *Nature*, *Cell*, *Immunity*, *Molecular Cancer*, etc.

💡 Efficient Service

- Ready-to-use products with guaranteed stock, shipped within 2 to 3 business days.
- Professional technical support responds within 8 h and provides solutions within 24 h.



ELISA Kits Features



- ✓ **High Precision:** Both inter and intra CV < 10%
- ✓ **High Sensitivity:** pg level
- ✓ **Multiple Species:** Human, Mouse, Rat, Rabbit, Monkey, Porcine, Chicken, etc.
- ✓ **High Specificity:** Cross-reactivity <10%
- ✓ **Multiple Sample Validation:** Serum, Plasma, Cell culture supernatant, Cell lysate, Tissue homogenate, etc.
- ✓ **Multiple Specifications:** 24 T/48 T/96 T*5 plates /96 T*10 plates

ELISA Product Series Comparison



Items	Traditional ELISA Kits	Quickey Pro™ ELISA Kits	MS ELISA Kits	High Sensitivity ELISA Kits
Benifit	Various Protein Targets	Minimum Experiment Time	Less Sample Volume	Higher Sensitivity
Sample Volume	100/50 µL	50 µL	25 µL	100 µL
Assay Time	3.5/2.5 h	90 min	3.5 h	3.5 h
Operation Steps	3/2 Steps	1 Step	3 Steps	3 Steps
Sensitivity	pg/mL	pg/mL	pg/mL	fg/mL
Storage	- 20°C/4°C	4°C	- 20°C/4°C	- 20°C/4°C

Top 200 Best-selling ELISA Kits (partial items)



Cat.No.	Product Name	Detection Range	Sensitivity
E-EL-0056	cAMP(Cyclic adenosine monophosphate) ELISA Kit	1.56-100 ng/mL	0.94 ng/mL
E-EL-0033	ST/5-HT(Serotonin/5-Hydroxytryptamine) ELISA Kit	15.63-1000 ng/mL	9.38 ng/mL
E-EL-0155	T(Testosterone) ELISA Kit	0.31-20 ng/mL	0.17 ng/mL
E-EL-0162	TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit	0.16-10 ng/mL	0.10 ng/mL
E-EL-H0108	Human IFN-γ(Interferon Gamma) ELISA Kit	15.63-1000 pg/mL	9.38 pg/mL
E-EL-H0149	Human IL-1β(Interleukin 1 Beta) ELISA Kit	7.81-500 pg/mL	4.69 pg/mL
E-EL-H6156	Human IL-6(Interleukin 6) ELISA Kit	1.56-100 pg/mL	0.94 pg/mL
E-EL-H0109	Human TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	7.81-500 pg/mL	4.69 pg/mL
E-EL-M0048	Mouse IFN-γ(Interferon Gamma) ELISA Kit	15.63-1000 pg/mL	9.38 pg/mL
E-EL-M0037	Mouse IL-1β(Interleukin 1 Beta) ELISA Kit	7.81-500 pg/mL	4.69 pg/mL
E-EL-M0044	Mouse IL-6(Interleukin 6) ELISA Kit	31.25-2000 pg/mL	18.75 pg/mL
E-EL-M0046	Mouse IL-10(Interleukin 10) ELISA Kit	15.63-1000 pg/mL	9.38 pg/mL
E-EL-M0730	Mouse IL-18(Interleukin 18) ELISA Kit	15.63-1000 pg/mL	9.38 pg/mL
E-EL-M1382	Mouse INS(Insulin) ELISA Kit	3.13-200 pg/mL	1.88 pg/mL
E-EL-M0792	Mouse MAU(Microalbuminuria) ELISA Kit	0.16-10 ng/mL	0.10 ng/mL
E-EL-M3063	Mouse TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	7.81-500 pg/mL	4.69 pg/mL
E-EL-M0676	Mouse HMGB-1(High Mobility Group Protein B1) ELISA Kit	15.63-1000 pg/mL	9.38 pg/mL
E-EL-R0012	Rat IL-1β(Interleukin 1 Beta) ELISA Kit	31.25-2000 pg/mL	18.75 pg/mL
E-EL-R0015	Rat IL-6(Interleukin 6) ELISA Kit	12.5-800 pg/mL	7.50 pg/mL
E-EL-R2856	Rat TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	15.63-1000 pg/mL	9.38 pg/mL

For more ELISA hot targets, please contact your local distributors or visit our website: www.elabscience.com

Part of High IF Publications

Targets	Cat.No.	Literature Information	Research Areas	IF
EPI NA/NE Mouse CORT Mouse GC Mouse GHRL Mouse GLP-1 Mouse OXA	E-EL-0045 E-EL-0047 E-OSEL-M0001 E-EL-M0555 E-EL-M0551 E-EL-M3012 E-EL-M0860	Meng J J, Shen J W, Li G, et al. Light modulates glucose metabolism by a retina-hypothalamus-brown adipose tissue axis[J]. <i>Cell</i> , 2023,186(2): 398-412.	Neuroscience Metabolism	66.85
cGMP Mouse ADH Mouse ALB Mouse BMG/β2-MG	E-EL-0083 E-EL-M0106 E-EL-M3032 E-EL-M2411	Xu W, Li G, Chen Y, et al. A novel antidiuretic hormone governs tumour-induced renal dysfunction[J]. <i>Nature</i> , 2023, 624, 425-432.	Cancer Metabolism	64.8
Human F8(Coagulation Factor VIII)	E-EL-H3641	Kaffe E, Roulis M, Zhao J , et al. Humanized mouse liver reveals endothelial control of essential hepatic metabolic functions[J]. <i>Cell</i> , 2023,186(18): 2793-3809.	Metabolism	64.5
Mouse ADH	E-EL-M0106	Bi Q Q , Wang C, Cheng G, et al. Microglia-derived PDGFB promotes neuronal potassium currents to suppress basal sympathetic tonicity and limit hypertension[J]. <i>Immunity</i> , 2022, 55(8): 1466-1482.	Neuroscience Cardiovascular	43.474

For more published information, please visit our website: www.elabscience.com

Elabscience Bionovation Inc.

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Email: techsupport@elabscience.com

Web: www.elabscience.com