

## Recombinant Monkeypox Virus Antigen (His tag)

Cat. No.: E-ELC-080      Size: 1mg

### Product Information

<b>Concentration:</b>	6.93mg/mL
<b>Source:</b>	Mammalian
<b>Protein Construction:</b>	A DNA sequence encoding the Monkey Virus Antigen was expressed with a polyhistidine tag.
<b>Application:</b>	Lateral Flow (Rapid Test) Other applications have not been validated yet.
<b>Formulation:</b>	Liquid
<b>Purity:</b>	> 90 % as determined by SDS-PAGE
<b>Molecular Mass:</b>	15.33 kDa
<b>Storage Buffer:</b>	20 mM PBS
<b>Storage:</b>	Can be stable for up to 12 months when stored at -20°C. Avoid repeated freeze-thaw cycles.