

Mouse anti-human IgG monoclonal

Cat. No.: E-ELC-003 Size: 1mg/5mg/10mg

Product Information

- Concentration:** 3.98mg/mL (Might be different among batches, refer to the manual for delivered product.)
- Host:** Mouse
- Specificity:** Specific to human IgG heavy chain
- Application:** Lateral Flow (Rapid Test)
Other applications should be validated by the end users
- Formulation:** Liquid
- Purified Method:** Protein A purified
- Purity:** > 95 %
- Cross Reactivity:** This antibody has been validated to not bind with human IgA or IgM
- Storage Buffer:** 0.02M PBS, pH 7.4, with 0.02% sodium azide
- Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.