Elabscience®

Tel:240-252-7368(USA)
Fax:240-252-7376(USA)
www.elabscience.com
E-mail:techsupport@elabscience.com
Elabscience Biotechnology Inc.

Mouse anti-human IgM monoclonal

Cat. No.: E-ELC-004 Size: 1mg/5mg/10mg

Product Information

Concentration: 2mg/mL (Might be different among batches, refer to the manual for

delivered product.)

Host: Mouse

Specificity: Specific to human IgM heavy chain

Application: Lateral Flow (Rapid Test)

Other applications should be validated by the end users

Formulation: Liquid

Purified Method: Protein A purified

Purity: > 99 %

Cross Reactivity: This antibody has been validated to not bind with human IgG or IgA

Storage Buffer: 0.02M PBS, pH 7.4, with 0.02% sodium azide

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.