

## Goat anti-chicken IgY polyclonal

Cat. No.: E-ELC-005      Size: 1mg/5mg/10mg

### Product Information

- Concentration:** Might be different among batches, refer to the manual for delivered product.)
- Host:** Goat
- Specificity:** Specific to chicken IgY
- Application:** Lateral Flow (Rapid Test)  
Other applications should be validated by the end users
- Formulation:** Liquid
- Purified Method:** Affinity purified
- Purity:** > 90 %
- Cross Reactivity:** This antibody has been validated to not bind with human IgG, rat IgG and mouse IgG
- Storage Buffer:** 10mM PBS, pH 7.4, with 0.02% sodium azide
- Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.