Elabscience®

Tel:240-252-7368(USA)
Fax:240-252-7376(USA)
www.elabscience.com
E-mail:techsupport@elabscience.com
Elabscience Biotechnology Inc.

Chicken IgY(with colloidal gold conjugated)

Cat. No.: E-ELC-006 Size: 1mL/5mL/10mL

Product Information

Conjugation: Colloidal gold

Source: Egg yolks

Application: Lateral Flow (Rapid Test)

Formulation: Liquid

Purity: > 90 % as determined by SDS-PAGE

Storage Buffer: 10mM PBS, pH 7.4, with 0.02% sodium azide

Storage: Store at 2-8°C for 2 months

The Chicken IgY has been conjugated with colloidal gold and is able to be applied as the control in COVID-19 Rapid Test Card development. It is recommended to use Elabscience's related products for this assay. Below is the corresponding product list.

Cat.NO.	Product name
E-ELC-002	Recombinant SARS-CoV-2 N Protein (with colloidal gold conjugated)
E-ELC-003	Mouse anti-human IgG monoclonal
E-ELC-004	Mouse anti-human IgM monoclonal
E-ELC-005	Goat anti-chicken IgY polyclonal
E-ELC-008	Buffer for Rapid Test development(02)