

Buffers for Rapid Test development(02)

Cat. No: E-ELC-008 Size: 1000/5000/10000 Tests

Product Information

| | Description | Storage |
|---------------|---|---------------|
| Buffer | Diluent for antibodies coated to NC filters | 4°C, 3 months |

It is recommended to use Elabscience's related products for this assay. Below is the corresponding product list.

| Cat.NO. | Product name |
|-----------|---|
| E-ELC-002 | Recombinant SARS-CoV-2 N Protein (with colloidal gold conjugated) |
| E-ELC-003 | Mouse anti-human IgG monoclonal |
| E-ELC-004 | Mouse anti-human IgM monoclonal |
| E-ELC-005 | Goat anti-chicken IgY polyclonal |
| E-ELC-006 | Chicken IgY(with colloidal gold conjugated) |