Elabscience®

Tel:240-252-7368(USA)
Fax:240-252-7376(USA)
www.elabscience.com
E-mail:techsupport@elabscience.com
Elabscience Biotechnology Inc.

Chicken IgY

Cat. No.: E-ELC-009 Size: 1mg/5mg/10mg

Product Information

Concentration Might be different among batches, refer to the manual for delivered product.

Source: Egg yolks

Application: Lateral Flow (Rapid Test)

Other applications have not been validated yet.

Formulation: Liquid

Purity: > 90 % as determined by SDS-PAGE

Storage Buffer: 10mM PBS, pH 7.4, with 0.02% sodium azide

Storage: Can be stable for up to 12 months when stored at -20 to -80℃. Avoid

repeated freeze-thaw cycles.