(FOR RESEARCH USE ONLY, DO NOT USE IT IN CLINICAL DIAGNOSTICS!)

Catalog No: E-TSEL-H0014

Product size: 96T/48T/24T/96T*5

QuicKey Human CLU(Clusterin) ELISA Kit

This manual must be read attentively and completely before using this product.

If you have any problems, please contact our Technical Service Center for help (info in the header of each page).

Tel: Fax: 1-832-243-6086 1-832-243-6017

Email:

techsupport@elabscience.com

Website:

www.elabscience.com

Please refer to specific expiry date from label on the side of box.

Please kindly provide us with the lot number (on the outside of the box) of the kit for more

efficient service.

-1-

QuicKey Series

Get more sensitive and precise results with saving at least 1h comparing to traditional ELISA Kits. The new developed technology in house will help to accelerate your science research in a more efficient way.

Intended use

This ELISA kit applies to the in vitro quantitative determination of Human CLU concentrations in serum, plasma, urine, saliva. Please consult technical support for the applicability if other biological fluids need to be tested.

Specification

- •Sensitivity: 0.94ng/mL.
- •Detection Range: 1.56-100ng/mL
- Specificity: This kit recognizes Human CLU in samples. No significant cross-reactivity or interference between Human CLU and analogues was observed.
- Repeatability: Coefficient of variation is < 10%.

Background

Clusterin is a ubiquitous protein that is arousing increasing interest owing to its widespread diffusion and multifunctional role. It is composed of two 40 kD subments (NA1, NA2) encoded by a single gene, which hold together by disulphide bonds. The mature dimer has a molecular mass of 80 kD and is micro-heterogeneous in electrical charge due to glycosylation. This protein has been differently named as apolipoprotein J (ApoJ), sulfated glycoprotein 2 (SGP-2), glycoprotein III (GP- III), testosterone repressed prostate message (TRPM-2), SP-40, 40cytolysis inhibitor (CL1) and others [1].It is a component of the high-density lipoprotein (HDL) complex with putative functions in the transport of lipids, apoptosis, protection of cells from stress and, most important, in the regulation of complement activity. It has been postulated that clusterin is involved in a variety of human diseases including Alzheimer disease, brain injury following ischemia, nephrosis and myocarditis [2].

- Wilson M R, Easterbrooksmith S B. Clusterin is a secreted mammalian chaperone. [J]. Trends in Biochemical Sciences, 2000, 25(3):95-98.
- Mclaughlin L, Zhu G, Mistry M, et al. Apolipoprotein J/clusterin limits the severity of murine autoimmune myocarditis. [J]. Journal of Clinical Investigation, 2000, 106(9):1105-1113.

Test principle

This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human CLU. Samples (or Standards) and biotinylated detection antibody specific for Human CLU are added to the micro ELISA plate wells. Human CLU would combine with the specific antibody. Then Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human CLU, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 ±2 nm. The OD value is proportional to the concentration of Human CLU. You can calculate the concentration of Human CLU in the samples by comparing the OD of the samples to the standard curve.

Kit components & Storage

An unopened kit can be stored at $2-8^{\circ}$ C for six months. After test, the unused wells and reagents should be stored according to the table below.

Item	Specifications	Storage	
Micro ELISA Plate (Dismountable)	96T: 8 wells ×12 strips 48T: 8 wells ×6 strips 24T: 8 wells ×3 strips 96T*5: 5 plates, 96T	2-8°C, 1 month	
Reference Standard	96T: 2 vials 48T/24T: 1 vial 96T*5: 10 vials	2-8°C, use the reconstituted standard within 24h	
Reference Standard & Sample Diluent	96T/48T/24T: 2 vials, 20 mL 96T*5: 10 vials, 20 mL		
Biotinylated Detection Ab Working Solution	96T/48T/24T: 1 vial, 6 mL 96T*5: 5 vials, 6 mL	2-8°C (Protect from light)	
HRP Conjugate Diluent	96T/48T/24T: 1 vial, 14 mL 96T*5: 5 vials, 14 mL		
Concentrated Wash Buffer(25×)	96T/48T/24T: 1 vial, 30 mL 96T*5: 5 vials, 30 mL		
Concentrated HRP Conjugate (100×)	96T: 1 vial, 120 μL 48T/24T: 1 vial, 60 μL 96T*5: 5 vials, 120 μL		
Substrate Reagent	96T/48T/24T: 1 vial, 10 mL 96T*5: 5 vials, 10 mL		
Stop Solution	96T/48T/24T: 1 vial, 10 mL 96T*5: 5 vials, 10 mL	2-8℃	
Plate Sealer	96T/48T/24T: 5 pieces 96T*5: 25 pieces		
Product Description	1 copy		
Certificate of Analysis	1 copy		

Note: All reagent bottle caps must be tightened to prevent evaporation and microbial pollution. The volume of reagents in partial shipments is a little more than the volume marked on the label, please use accurate measuring equipment instead of directly pouring into the vial(s).

Other supplies required

Note

- Please wear lab coats, eye protection and latex gloves for protection. Please perform
 the experiment following the national security protocols of biological laboratories,
 especially when detecting blood samples or other bodily fluids.
- 2. A freshly opened ELISA Plate may appear to have a water-like substance, which is normal and will not have any impact on the experimental results. Return the unused wells to the foil pouch provided in the kit, store it according to the conditions suggested in the above table
- 3. Do not reuse the reconstituted standard and HRP conjugate working solution. The unspent biotinylated detection Ab working solution and other stock solutions should be stored according to the storage conditions suggested in the above table.
- 4. The microplate reader should be able to be installed with a filter that can detect the wave length at 450±2 nm. The optical density should be within 0~3.5.
- 5. Do not mix or use components with those from other lots.
- Change pipette tips in between adding of each standard level, between sample adding and between reagent adding. Also, use separate reservoirs for each reagent.

Sample collection

(More detailed information please view our website: http://www.elabscience.com/List-detail-253.html)

Serum: Allow samples to clot for 1 hour at room temperature or overnight at 2-8°C before centrifugation for 15 min at $1000 \times g$ at $2 \sim 8$ °C. Collect the supernatant to carry out the assay. Blood collection tubes should be disposable and be non-endotoxin.

Plasma: Collect plasma using EDTA or heparin as an anticoagulant. Centrifuge samples for 15 min at 1000×g at 2~8°C within 30 min of collection. Collect the supernatant to carry out the assay. Hemolysed samples are not suitable for ELISA assay!

Saliva: Remove particulates by centrifugation for 10 minutes at $4000 \times g$ at $2 \sim 8$ °C. Collect the supernatant to carry out the assay. Recommend to use fresh saliva samples.


Urine: Use a sterile container to collect urine samples. Remove particulates by centrifugation for 15 minutes at $1000 \times g$ at $2 \sim 8$ °C.Collect the supernatant to carry out the assay.

Note for sample

- Samples should be assayed within 7 days when stored at 2-8℃, otherwise samples
 must be divided up and stored at -20℃ (≤1 month) or -80℃ (≤3 months). Avoid
 repeated freeze-thaw cycles.
- Please predict the concentration before assaying. If the sample concentration is not within the range of the standard curve, users must determine the optimal sample dilutions for their particular experiments.
- Some recombinant protein may not be detected due to a mismatching with the coated antibody or detection antibody.

Reagent preparation

- 1. Bring all reagents to room temperature ($18\sim25^{\circ}$ C) before use. Follow the Microplate reader manual for set-up and preheat it for 15 min before OD measurement.
- Wash Buffer: Dilute 30 mL of Concentrated Wash Buffer with 720 mL of deionized
 or distilled water to prepare 750 mL of Wash Buffer. Note: if crystals have formed in
 the concentrate, warm it in a 40°C water bath and mix it gently until the crystals have
 completely dissolved.
- 3. Standard working solution: Centrifuge the standard at 10,000×g for 1 min. Add 1.0 mL of Reference Standard &Sample Diluent, let it stand for 10 min and invert it gently several times. After it dissolves fully, mix it thoroughly with a pipette. This reconstitution produces a working solution of 100ng/mL. Then make serial dilutions as needed. The recommended dilution gradient is as follows: 100、50、25、12.5、6.25、3.13、1.56、0ng/mL.
 - Dilution method: Take 7 EP tubes, add $500\mu L$ of Reference Standard & Sample Diluent to each tube. Pipette $500\mu L$ of the 100ng/mL working solution to the first tube and mix up to produce a 50ng/mL working solution. Pipette $500\mu L$ of the solution from the former tube into the latter one according to this step. The illustration below is for reference. Note: the last tube is regarded as a blank. Don't pipette solution into it from the former tube.



4. **HRP Conjugate working solution:** HRP Conjugate is HRP conjugated avidin.Calculate the required amount before the experiment (100 μL/well). In preparation, slightly more than calculated should be prepared. Centrifuge the Concentrated HRP Conjugate at 800×g for 1 min, then dilute the 100×Concentrated HRP Conjugate to 1×working solution with HRP Conjugate Diluent.

Assay procedure (A brief assay procedure is on the 12th page)

- 1. Add the Standard working solution to the first two columns: Each concentration of the solution is added in duplicate, to one well each, side by side (50 μL for each well). Add the samples to the other wells (50 μL for each well). Immediately add 50 μL of Biotinylated Detection Ab working solution to each well. Cover the plate with the sealer provided in the kit. Incubate for 90 min at 37 °C. Note: solutions should be added to the bottom of the micro ELISA plate well, avoid touching the inside wall and causing foaming as much as possible.
- 2. Aspirate or decant the solution from each well, add 350μL of wash buffer to each well. Soak for 1~2 min and aspirate or decant the solution from each well and pat it dry against clean absorbent paper. Repeat this wash step 3 times. Note: a microplate washer can be used in this step and other wash steps.
- Add 100µL of HRP Conjugate working solution to each well. Cover with the Plate sealer. Incubate for 30 min at 37 °C.
- Aspirate or decant the solution from each well, repeat the wash process for 5 times as conducted in step 2.
- 5. Add 90μL of Substrate Reagent to each well. Cover with a new plate sealer. Incubate for about 15 min at 37 °C. Protect the plate from light. Note: the reaction time can be shortened or extended according to the actual color change, but not more than 30min.
- Add 50μLof Stop Solution to each well. Note: adding the stop solution should be done
 in the same order as the substrate solution.
- 7. Determine the optical density (OD value) of each well at once with a micro-plate reader

set to 450 nm.

Calculation of results


Average the duplicate readings for each standard and samples, then subtract the average zero standard optical density. Plot a four parameter logistic curve on log-log graph paper, with standard concentration on the x-axis and OD values on the y-axis.

If the samples have been diluted, the concentration calculated from the standard curve must be multiplied by the dilution factor. If the OD of the sample surpasses the upper limit of the standard curve, you should re-test it with an appropriate dilution. The actual concentration is the calculated concentration multiplied by the dilution factor.

Typical data

As the OD values of the standard curve may vary according to the conditions of the actual assay performance (e.g. operator, pipetting technique, washing technique or temperature effects), the operator should establish a standard curve for each test. Typical standard curve and data is provided below for reference only.

Concentration(ng/mL)	100	50	25	12.5	6.25	3.13	1.56	0
OD	2.460	1.804	1.113	0.589	0.376	0.227	0.162	0.071
Corrected OD	2.389	1.733	1.042	0.518	0.305	0.156	0.091	



Sample values

Serum/Plasma/Urine/Saliva—Samples from apparently healthy volunteers were evaluated for the presence of Human CLU in this assay.

Sample Type	Source	Range	Dilution Factor
Serum (n=14)	Healthy human	60.58-346.42μg/mL	4000-32000
EDTA plasma (n=11)	Healthy human	158.60-381.21μg/mL	8000-32000
Saliva (n=10)	Healthy human	30.12-614.25ng/mL	2-50
Urine (n=12)	Healthy human	25.99-138.36ng/mL	2-10

Precision

Intra-assay Precision (Precision within an assay): 3 samples with low, mid range and high level Human CLU were tested 20 times on one plate, respectively.

Inter-assay Precision (Precision between assays): 3 samples with low, mid range and high level Human CLU were tested on 3 different plates, 20 replicates in each plate.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
n	20	20	20	20	20	20
Mean((ng/mL)	4.57	10.92	48.94	4.6	12.45	41.92
Standard deviation	0.24	0.51	1.99	0.24	0.62	1.95
CV (%)	5.25	4.67	4.07	5.22	4.98	4.65

Recovery

The recovery of Human CLU spiked at three different levels in samples throughout the range of the assay was evaluated in various matrices.

Sample Type	Range (%)	Average Recovery (%)
Serum (n=14)	90-101	97
EDTA plasma (n=11)	88-98	95
Saliva (n=10)	83-104	91
Urine (n=12)	83-96	91

Linearity


Samples were spiked with high concentrations of Human CLU and diluted with Reference Standard & Sample Diluent to produce samples with values within the range of the assay.

		Serum (n=4)	EDTA plasma (n=4)	Saliva (n=4)	Urine (n=4)
1:2	Range (%)	88-96	90-96	87-120	83-97
1.2	Average (%)	91	93	106	91
1:4	Range (%)	89-98	85-94	87-114	98-112
1.4	Average (%)	95	90	100	105
1.0	Range (%)	94-107	92-100	86-115	93-106
1:8	Average (%)	100	96	96	101
1:16	Range (%)	101-110	99-112	87-110	94-107
	Average (%)	105	105	99	102

Troubleshooting

Troubleshooting Problem causes solutions					
Problem		0.0-1-00			
	Inaccurate pipetting	Check pipettes.			
Poor standard curve	Improper standard dilution	Ensure briefly spin the vial of standard and dissolve the powder thoroughly by gentle mixing.			
	Wells are not completely aspirated	Completely aspirate wells in between steps.			
	Insufficient incubation time	Ensure sufficient incubation time.			
	Incorrect assay temperature	Use recommended incubation temperature. Bring substrate to room temperature before use.			
Low signal	Inadequate reagent volumes	Check pipettes and ensure correct preparation.			
	Improper dilution	propulation.			
	HRP conjugate inactive or TMB failure	Mix HRP conjugate and TMB, rapid coloring.			
Deep color but low value	Plate reader setting is not optimal	Verify the wavelength and filter setting on the Microplate reader.			
Large CV	Inaccurate pipetting	Check pipettes			
	Concentration of target protein is too high	Use recommended dilution factor.			
High background	Plate is insufficiently washed	Review the manual for proper wash. If using a plate washer, check that all ports are unobstructed.			
	Contaminated wash buffer	Prepare fresh wash buffer.			
Low sensitivity	Improper storage of the ELISA kit	All the reagents should be stored according to the instructions.			
Low sensitivity	Stop solution is not added	Stop solution should be added to each well before measurement.			

SUMMARY



Declaration

- Limited by current conditions and scientific technology, we can't conduct comprehensive identification and analysis on all the raw material provided. So there might be some qualitative and technical risks for users using the kit.
- The final experimental results will be closely related to the validity of products, operational skills of the operators and the experimental environments. Please make sure that sufficient samples are available.