
(FOR RESEARCH USE ONLY. DO NOT USE IT IN CLINICAL DIAGNOSIS !)

TBST Buffer (10 ×)

Catalog No: E-BC-R335

Sizes: 100 mL/ 200 mL/ 500 mL

Cat	Products	100 mL	200 mL	500 mL	Storage
E-BC-R335	TBST Buffer (10 ×)	100 mL	200 mL	500 mL	RT
Manual			1 copy		

This manual must be read attentively and completely before using this product.

If you have any problems, please contact our Technical Service Center for help.

Phone: 240-252-7368(USA) 240-252-7376(USA)

Email: techsupport@elabscience.com

Website: www.elabscience.com

Please kindly provide us the lot number (on the outside of the box) of the kit for more efficient service.

Introduction

TBST buffer is used in Western Blot to wash away non-specific binding on the membrane, this buffer contains 0.1% Tween20. Tween20 is a non-ionic detergent that increases the buffer's ability to elute and reduce non-specific binding of antibodies. This product is supplied with 10 × concentrate and easy for transportation.

Instructions

This product is a 10 × concentrate solution, dilute with pure water to 1 × working solution before use. For example: take 100 mL TBST buffer (10 ×), dilute with pure water to 1 L.

1 × TBST Buffer Components

145.4 mM NaCl, 10 mM Tris-base, 0.1% (v/v) Tween20, pH7.5.

Storage

Store at room temperature for 12 months.

Cautions

1. If the product crystallizes, please use water bath to dissolve at 37 °C, and then dilute it after it is completely dissolved.
2. Keep the product sealed to prevent from pollution.
3. For your safety and health, please wear the lab coat and disposable gloves before the experiments.