Blocking Powder

Cat. No: E-IR-R108 Size: 5 g / 10 g / 50 g



Product Content

| Cat. | Products | 5 g | 10 g | 50 g | Storage |
|-----------|-----------------|-----|------|------|---------|
| E-IR-R108 | Blocking Powder | 5 g | 10 g | 50 g | 2~8 ℃ |

Introduction

The product has high purity and is light yellow white powder.

This product is widely used in immune related experiments, such as Western Blotting, ELISA and immunostaining, dilution of primary antibody or secondary antibody, etc. In addition, the product can also be used as a standard for quantitative detection of protein concentration.

The solubility of the product in PBS is more than 50 mg/mL, it can be prepared into a transparent 5% BSA blocking buffer.

Storage

Store at $2 \sim 8 \text{ } \text{C}$ for 12 months.

Cautions

- 1. This product couldn't be used for clinical diagnosis or treatment, food or medicine, and can't be stored in residence.
- 2. For your safety and health, please wear the lab coat and disposable gloves before the experiments.

Tel: 1-832-243-6086 Email: techsupport@elabscience.com