

# Normal Goat Blocking Buffer



Cat. No: E-IR-R111

Size: 5 mL / 10 mL

## Product Content

Cat.	Products	5 mL	10 mL	Storage
E-IR-R111	Normal Goat Blocking Buffer	5 mL	10 mL	2~8 °C

## Introduction

This product is normal goat Blocking Buffer, added stabilizer, As a blocking buffer, It can effectively reduce the nonspecific staining of IHC/IF experiment, so as to reduce the background.

## Instructions

The product is concentrated normal goat Blocking Buffer. Please dilute with PBS to 1× working solution before use.

For example: take 10 μL Normal Goat Blocking Buffer, add it to 90 μL PBS and mix.

## Storage

Store at 2~8 °C for 12 months.

## Cautions

1. This product couldn't be used for clinical diagnosis or treatment, food or medicine, and can't be stored in residence.
2. For your safety and health, please wear the lab coat and disposable gloves before the experiments.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017

Rev. V1.3