


Elabscience®

Focus on lab research Service for life science

2015 Catalog

For CLIA Kits

www.elabscience.com



About us

Elabscience is a high-tech biotech company, specializing in reagents of immunodiagnostic technology. It offers more than 20,000 products, including ELISA kits, CLIA kits, Antibodies, and Proteins.

Besides, Elabscience has established seven technical service platforms , including ELISA Customized Service, Immunoassay service, Molecular Biology Service, Protein Expression and Purification Service, Antibody Preparation Service, Cell Culture Service and Whole Research Project Service, to provide one-stop experimental service with professional knowledge, facilities and attitude.




Products

Elabscience has a portfolio of 20,000 products at present. Its website keeps updated by adding hundreds of products every month, so please keep close attention to our website to get your needed products.

Elabscience's products cover following categories:

- ▲ ELISA kits
- ▲ CLIA kits
- ▲ Monoclonal antibodies
- ▲ Polyclonal antibodies
- ▲ Proteins
- ▲ cDNA



Technical support

In Elabscience, we have experienced experts in the field of research and development, endeavoring to provide customers all-round technical support.

Any quality problems of products, we would replace the products or refund.

Any technical support in the process of operation, you can ask for help through Tel, Email, or livechat at the same time, you

can also loge in our technical forum to check the related articles of your experiment.



Tel: 86-27-87385095

Fax: 86-27- 87645690

Email: techsupport@elabscience.com

Technical Forum:

<http://www.elabscience.com/index.php/resources/show/tid/49>

Products:

Cat.NO	Product Name
E-CL-H0001	Human α 1-AGP (Alpha-1-Acid Glycoprotein) CLIA Kit
E-CL-H0002	Human ACE (Angiotensin I Converting Enzyme) CLIA Kit
E-CL-H0003	Human ACVA (Activin A) CLIA Kit
E-CL-H0004	Human ADP/Acrp30 (Adiponectin) CLIA Kit
E-CL-H0005	Human AgRP (Agouti Related Protein) CLIA Kit
E-CL-H0006	Human ANG (Angiogenin) CLIA Kit
E-CL-H0007	Human ANG1 (Angiopoietin 1) CLIA Kit
E-CL-H0008	Human ANG2 (Angiopoietin 2) CLIA Kit
E-CL-H0009	Human BAFF/CD257 (B-Cell Activating Factor) CLIA Kit
E-CL-H0010	Human BDNF (Brain Derived Neurotrophic Factor) CLIA Kit
E-CL-H0011	Human BMP-2 (Bone Morphogenetic Protein 2) CLIA Kit
E-CL-H0012	Human BMP-4 (Bone Morphogenetic Protein 4) CLIA Kit
E-CL-H0013	Human BMP-7 (Bone Morphogenetic Protein 7) CLIA Kit
E-CL-H0014	Human E-Cad (E-Cadherin) CLIA Kit
E-CL-H0015	Human CA9 (Carbonic Anhydrase IX) CLIA Kit
E-CL-H0016	Human CASP1 (Caspase 1) CLIA Kit
E-CL-H0017	Human CASP3 (Caspase 3) CLIA Kit
E-CL-H0018	Human CTSB (Cathepsin B) CLIA Kit

Cat.NO	Product Name
E-CL-H0019	Human CTSV (Cathepsin V) CLIA Kit
E-CL-H0020	Human MCP-1 (Monocyte Chemotactic Protein 1) CLIA Kit
E-CL-H0021	Human MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) CLIA Kit
E-CL-H0022	Human MIP-1 β (Macrophage Inflammatory Protein 1 Beta) CLIA Kit
E-CL-H0023	Human RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) CLIA Kit
E-CL-H0024	Human MCP-3 (Monocyte Chemotactic Protein 3) CLIA Kit
E-CL-H0025	Human ECF/CCL11 (Eosinophil Chemotactic Factor) CLIA Kit
E-CL-H0026	Human TARC (Thymus Activation Regulated Chemokine) CLIA Kit
E-CL-H0027	Human MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) CLIA Kit
E-CL-H0028	Human SLC (Secondary Lymphoid Tissue Chemokine) CLIA Kit
E-CL-H0029	Human MDC (Macrophage Derived Chemokine) CLIA Kit
E-CL-H0030	Human MPIF2 (Myeloid Progenitor Inhibitory Factor 2) CLIA Kit
E-CL-H0031	Human MIP4 α (Macrophage Inflammatory Protein 4 Alpha) CLIA Kit
E-CL-H0032	Human CTACK (Cutaneous T-cell Attracting Chemokine) CLIA Kit
E-CL-H0033	Human sCD14 (Soluble Cluster of Differentiation14) CLIA Kit
E-CL-H0034	Human Fc ϵ RII/CD23 (Receptor II for the Fc Region of Immunoglobulin E) CLIA Kit
E-CL-H0035	Human CD40L/TNFSF5 (Cluster of Differentiation 40 Ligand) CLIA Kit
E-CL-H0036	Human CD163 (Cluster of Differentiation 163) CLIA Kit
E-CL-H0037	Human GP39 (Glycoprotein 39, Cartilage) CLIA Kit

Cat.NO	Product Name
E-CL-H0038	Human CLU (Clusterin) CLIA Kit
E-CL-H0039	Human CNTF (Ciliary Neurotrophic Factor) CLIA Kit
E-CL-H0040	Human TF (Tissue Factor) CLIA Kit
E-CL-H0041	Human CFD (Complement Factor D) CLIA Kit
E-CL-H0042	Human CRN (Corin) CLIA Kit
E-CL-H0043	Human CRP (C-Reactive Protein) CLIA Kit
E-CL-H0044	Human CX3CL1 (Chemokine C-X3-C-Motif Ligand 1) CLIA Kit
E-CL-H0045	Human GRO α /CXCL1 (Growth Regulated Oncogene Alpha) CLIA Kit
E-CL-H0046	Human ENA-78 (Epithelial Neutrophil Activating Peptide 78) CLIA Kit
E-CL-H0047	Human GCP-2 (Granulocyte Chemotactic Protein 2) CLIA Kit
E-CL-H0048	Human IL-8 (Interleukin 8) CLIA Kit
E-CL-H0049	Human AAHA (Anti-Apolipoprotein Antibody) CLIA Kit
E-CL-H0050	Human IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) CLIA Kit
E-CL-H0051	Human I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) CLIA Kit
E-CL-H0052	Human SDF-1 (Stromal Cell Derived Factor 1) CLIA Kit
E-CL-H0053	Human BLC (B-Lymphocyte Chemoattractant) CLIA Kit
E-CL-H0054	Human CXCL16 (Chemokine C-X-C-Motif Ligand 16) CLIA Kit
E-CL-H0055	Human Cys-C (Cystatin C) CLIA Kit
E-CL-H0056	Human Cyt-C (Cytochrome C) CLIA Kit

Cat.NO	Product Name
E-CL-H0057	Human DKK1 (Dickkopf Related Protein 1) CLIA Kit
E-CL-H0058	Human DPP4 (Dipeptidyl Peptidase IV) CLIA Kit
E-CL-H0059	Human EGF (Epidermal Growth Factor) CLIA Kit
E-CL-H0060	Human EGFR (Epidermal Growth Factor Receptor) CLIA Kit
E-CL-H0061	Human EG-VEGF (Endocrine Gland Derived Vascular Endothelial Growth Factor)CLIA Kit
E-CL-H0062	Human ENG (Endoglin) CLIA Kit
E-CL-H0063	Human ES (Endostatin) CLIA Kit
E-CL-H0064	Human ET-1 (Endothelin 1) CLIA Kit
E-CL-H0065	Human EPCR (Endothelial Protein C Receptor) CLIA Kit
E-CL-H0066	Human EPO (Erythropoietin) CLIA Kit
E-CL-H0067	Human FAS/CD95 (Factor Related Apoptosis) CLIA Kit
E-CL-H0068	Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) CLIA Kit
E-CL-H0069	Human FABP1 (Fatty Acid Binding Protein 1, Liver) CLIA Kit
E-CL-H0070	Human α FP (Alpha-Fetoprotein) CLIA Kit
E-CL-H0071	Human AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) CLIA Kit
E-CL-H0072	Human AAP (Alanine Aminopeptidase) CLIA Kit
E-CL-H0073	Human FGF19 (Fibroblast Growth Factor 19) CLIA Kit
E-CL-H0074	Human FGF21 (Fibroblast Growth Factor 21) CLIA Kit
E-CL-H0075	Human Flt3L (FMS Like Tyrosine Kinase 3 Ligand) CLIA Kit

Cat.NO	Product Name
E-CL-H0076	Human FST (Follistatin) CLIA Kit
E-CL-H0077	Human Hepc (Hepcidin) CLIA Kit
E-CL-H0078	Human Gas6 (Growth Arrest Specific Protein 6) CLIA Kit
E-CL-H0079	Human G-CSF (Granulocyte Colony Stimulating Factor 3) CLIA Kit
E-CL-H0080	Human GDF15 (Growth Differentiation Factor 15) CLIA Kit
E-CL-H0081	Human GM-CSF (Granulocyte Macrophage Colony Stimulating Factor) CLIA Kit
E-CL-H0082	Human sgp130 (soluble Glucoprotein 130) CLIA Kit
E-CL-H0083	Human HGF (Hepatocyte Growth Factor) CLIA Kit
E-CL-H0084	Human IFN- β (Interferon Beta) CLIA Kit
E-CL-H0085	Human IGF-1 (Insulin-Like Growth Factor 1) CLIA Kit
E-CL-H0086	Human IGFBP-3 (Insulin Like Growth Factor Binding Protein 3) CLIA Kit
E-CL-H0087	Human IL-1 α (Interleukin 1 Alpha) CLIA Kit
E-CL-H0088	Human IL-1ra (Interleukin 1 Receptor Antagonist) CLIA Kit
E-CL-H0089	Human IL-1R2 (Interleukin 1 Receptor Type II) CLIA Kit
E-CL-H0090	Human PSA/KLK3 (Prostate Specific Antigen) CLIA Kit
E-CL-H0091	Human FGF7 (Fibroblast Growth Factor 7) CLIA Kit
E-CL-H0092	Human LEPR (Leptin Receptor) CLIA Kit
E-CL-H0093	Human LIF (Leukemia Inhibitory Factor) CLIA Kit
E-CL-H0094	Human TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14)CLIA Kit

Cat.NO	Product Name
E-CL-H0095	Human NGAL (Neutrophil Gelatinase Associated Lipocalin) CLIA Kit
E-CL-H0096	Human MCSF (Macrophage Colony Stimulating Factor 1) CLIA Kit
E-CL-H0097	Human Pro-MMP-1 (pro-Matrix Metalloproteinase 1) CLIA Kit
E-CL-H0098	Human IL-2 (Interleukin 2) CLIA Kit
E-CL-H0099	Human IL-3 (Interleukin 3) CLIA Kit
E-CL-H0100	Human IL-4 (Interleukin 4) CLIA Kit
E-CL-H0101	Human IL-6 (Interleukin 6) CLIA Kit
E-CL-H0102	Human IL-10 (Interleukin 10) CLIA Kit
E-CL-H0103	Human IL-13 (Interleukin 13) CLIA Kit
E-CL-H0104	Human IL-17 (Interleukin 17) CLIA Kit
E-CL-H0105	Human IL-22 (Interleukin 22) CLIA Kit
E-CL-H0106	Human IL-23 (Interleukin 23) CLIA Kit
E-CL-H0107	Human IFN- γ (Interferon Gamma) CLIA Kit
E-CL-H0108	Human TNF- α (Tumor Necrosis Factor Alpha) CLIA Kit
E-CL-H0109	Human TGF- β 1 (Transforming Growth Factor Beta 1) CLIA Kit
E-CL-H0110	Human VEGF (Vascular Endothelial Cell Growth Factor) CLIA Kit
E-CL-H0111	Human EREG (Epiregulin) CLIA Kit
E-CL-H0112	Human LEP (Leptin) CLIA Kit
E-CL-H0113	Human Bcl-2 (B-cell Leukemia/Lymphoma 2) CLIA Kit

Cat.NO	Product Name
E-CL-H0114	Human MAU (Microalbuminuria) CLIA Kit
E-CL-H0115	Human LH (Luteinizing Hormone) CLIA Kit
E-CL-H0116	Human REN (Renin) CLIA Kit
E-CL-H0117	Human CTSC (Cathepsin C) CLIA Kit
E-CL-H0118	Human C-P (C-Peptide) CLIA Kit
E-CL-H0119	Human T (Testosterone) CLIA Kit
E-CL-H0120	Human OxLDL (Oxidized Low Density Lipoprotein) CLIA Kit
E-CL-H0121	Human ApoA1 (Apolipoprotein A1) CLIA Kit
E-CL-H0122	Human LN (Laminin) CLIA Kit
E-CL-H0123	Human COL1 (Collagen Type I) CLIA Kit
E-CL-H0124	Human MMP-25 (Matrix Metalloproteinase 25) CLIA Kit
E-CL-H0125	Human MMP-10 (Matrix Metalloproteinase 10) CLIA Kit
E-CL-H0126	Human MMP-12 (Matrix Metalloproteinase 12) CLIA Kit
E-CL-H0127	Human MMP-13 (Matrix Metalloproteinase 13) CLIA Kit
E-CL-H0128	Human sICAM-1/CD54 (Soluble Intercellular Adhesion Molecule 1) CLIA Kit
E-CL-H0129	Human CT (Calcitonin) CLIA Kit
E-CL-H0130	Human ACTH (Adrenocorticotrophic Hormone) CLIA Kit
E-CL-H0131	Human DA (Dopamine) CLIA Kit
E-CL-H0132	Human ACAT2 (Acetyl Coenzyme A Acetyltransferase 2, cytosolic) CLIA Kit

Cat.NO	Product Name
E-CL-H0133	Human PRL (Prolactin) CLIA Kit
E-CL-H0134	Human MYO (Myoglobin) CLIA Kit
E-CL-H0135	Human Tn-T (Troponin T) CLIA Kit
E-CL-H0136	Human Tn-I (Troponin I) CLIA Kit
E-CL-H0137	Human PAPPA (Pregnancy Associated Plasma Protein A) CLIA Kit
E-CL-H0138	Human ALP (Alkaline Phosphatase) CLIA Kit
E-CL-H0139	Human CSTA (Cystatin A) CLIA Kit
E-CL-H0140	Human GLP-1 (Glucagon Like Peptide 1) CLIA Kit
E-CL-H0141	Human IL-1 β (Interleukin 1 Beta) CLIA Kit
E-CL-H0142	Human IL-12 (Interleukin 12) CLIA Kit
E-CL-H0143	Human IL-12 p40 (Interleukin 12 p40) CLIA Kit
E-CL-H0144	Human CP (Ceruloplasmin) CLIA Kit
E-CL-H0145	Human α -CTx (Alpha Crosslaps) CLIA Kit
E-CL-H0146	Human CD42 (Cluster of Differentiation 42) CLIA Kit
E-CL-H0147	Human HLA-B27 (Human Leukocyte Antigen B27) CLIA Kit
E-CL-H0148	Human FGFR2 (Fibroblast Growth Factor Receptor 2) CLIA Kit
E-CL-H0149	Human IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) CLIA Kit
E-CL-H0150	Human LP-a (Lipoprotein a) CLIA Kit
E-CL-H0151	Human MBP (Myelin Basic Protein) CLIA Kit

Cat.NO	Product Name
E-CL-H0152	Human TFPI (Tissue Factor Pathway Inhibitor) CLIA Kit
E-CL-H0153	Human TM (Thrombomodulin) CLIA Kit
E-CL-H0154	Human C4 (Complement Component 4) CLIA Kit
E-CL-H0155	Human FE (Ferritin) CLIA Kit
E-CL-H0156	Human IgG (Immunoglobulin G) CLIA Kit
E-CL-H0157	Human TRF (Transferrin) CLIA Kit
E-CL-H0158	Human FSTL1 (Follistatin Like Protein 1) CLIA Kit
E-CL-H0159	Human LOX (Lysyl Oxidase) CLIA Kit
E-CL-H0160	Human CG (Chorionic Gonadotropin) CLIA Kit
E-CL-H0161	Human GH (Growth Hormone) CLIA Kit
E-CL-H0162	Human COL4 (Collagen Type IV) CLIA Kit
E-CL-H0163	Human FN (Fibronectin) CLIA Kit
E-CL-H0164	Human IL-9 (Interleukin 9) CLIA Kit
E-CL-H0165	Human PCI (Procollagen I) CLIA Kit
E-CL-H0166	Human PCIII (Procollagen III) CLIA Kit
E-CL-H0167	Human PIIINP (N-Terminal Procollagen III Propeptide) CLIA Kit
E-CL-H0168	Human TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) CLIA Kit
E-CL-H0169	Human PINP (Procollagen I N-terminal Propeptide) CLIA Kit
E-CL-H0170	Human KIM-1 (Kidney Injury Molecule 1) CLIA Kit

Cat.NO	Product Name
E-CL-H0171	Human ADCY2 (Adenylate Cyclase 2) CLIA Kit
E-CL-H0172	Human AGGF1 (Angiogenic Factor with G Patch and FHA Domains 1) CLIA Kit
E-CL-H0173	Human C4a (Complement Component 4a) CLIA Kit
E-CL-H0174	Human C5a (Complement Component 5a) CLIA Kit
E-CL-H0175	Human IL-5 (Interleukin 5) CLIA Kit
E-CL-H0176	Human IL-6R (Interleukin 6 Receptor) CLIA Kit
E-CL-H0177	Human FcγR I (Receptor I for the Fc Region of Immunoglobulin G) CLIA Kit
E-CL-H0178	Human NCAD (Neural Cadherin) CLIA Kit
E-CL-H0179	Human PICP (Procollagen I C-terminal Propeptide) CLIA Kit
E-CL-H0180	Human PDGFRL (Platelet Derived Growth Factor Receptor Like Protein) CLIA Kit
E-CL-H0181	Human LUM (Lumican) CLIA Kit
E-CL-H0182	Human PRR4 (Proline Rich Protein 4, Lacrimal) CLIA Kit
E-CL-H0183	Human PSAP (Prosaposin) CLIA Kit
E-CL-H0184	Human PTGES (Prostaglandin E Synthase, Microsomal) CLIA Kit
E-CL-H0185	Human ISG15 (Interferon Stimulated Gene 15) CLIA Kit
E-CL-H0186	Human AIMP1 (Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1) CLIA Kit
E-CL-H0187	Human 11-DH-TXB2 (11-Dehydro-thromboxane B2) CLIA Kit
E-CL-H0188	Human 12-HETE (12-Hydroxyeicosatetraenoic Acid) CLIA Kit
E-CL-H0189	Human 15-LO (Arachidonate 15-Lipoxygenase) CLIA Kit

Cat.NO	Product Name
E-CL-H0190	Human BACE1 (Beta-site APP Cleaving Enzyme 1) CLIA Kit
E-CL-H0191	Human IL-32 (Interleukin 32) CLIA Kit
E-CL-H0192	Human TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) CLIA Kit
E-CL-H0193	Human 26S-PSM (26S Proteasome) CLIA Kit
E-CL-H0194	Human 3-NT (3-Nitrotyrosine) CLIA Kit
E-CL-H0195	Human 5-HIAA (5-Hydroxyindolacetic Acid) CLIA Kit
E-CL-H0196	Human IL-15 (Interleukin 15) CLIA Kit
E-CL-H0197	Human 5-LO (Arachidonate 5-Lipoxygenase) CLIA Kit
E-CL-H0198	Human 5-NT (5-Nucleotidase) CLIA Kit
E-CL-H0199	Human 8-epi-PGF2 α (8-Epi-Prostaglandin F2 Alpha) CLIA Kit
E-CL-H0200	Human 5-HETE (5-Hydroxyeicosatetraenoic Acid) CLIA Kit
E-CL-H0201	Human AAA (Anti-Albumin Antibody) CLIA Kit
E-CL-H0202	Human α 1ACT (Alpha1 Antichymotrypsin) CLIA Kit
E-CL-H0203	Human IL-16 (Interleukin 16) CLIA Kit
E-CL-H0204	Human AP17 (Apelin 17) CLIA Kit
E-CL-H0205	Human AR (Amphiregulin) CLIA Kit
E-CL-H0206	Human INHA (Inhibin A) CLIA Kit
E-CL-H0207	Human ARO (Aromatase) CLIA Kit
E-CL-H0208	Human BLVRB (Biliverdin Reductase B) CLIA Kit

Cat.NO	Product Name
E-CL-H0209	Human ACA-IgA (Anti-Cardiolipin Antibody IgA) CLIA Kit
E-CL-H0210	Human AChE (Acetylcholinesterase) CLIA Kit
E-CL-H0211	Human IL-18 (Interleukin 18) CLIA Kit
E-CL-H0212	Human IL-19 (Interleukin 19) CLIA Kit
E-CL-H0213	Human ACTG1 (Actin Gamma 1) CLIA Kit
E-CL-H0214	Human ACTN2 (Actinin Alpha 2) CLIA Kit
E-CL-H0215	Human ACTN3 (Actinin Alpha 3) CLIA Kit
E-CL-H0216	Human ACT α 2 (Actin Alpha 2, Smooth Muscle) CLIA Kit
E-CL-H0217	Human ACT α 1 (Actin Alpha 1, Skeletal Muscle) CLIA Kit
E-CL-H0218	Human ADA (Adenosine Deaminase) CLIA Kit
E-CL-H0219	Human ADAM10 (A Disintegrin And Metalloprotease 10) CLIA Kit
E-CL-H0220	Human ADAM8 (A Disintegrin And Metalloprotease 8) CLIA Kit
E-CL-H0221	Human ADAM9 (A Disintegrin And Metalloprotease 9) CLIA Kit
E-CL-H0222	Human ADAMTS4 (ADAM with Thrombospondin Type 1 Motif 4) CLIA Kit
E-CL-H0223	Human ADCY1 (Adenylate Cyclase 1, Brain) CLIA Kit
E-CL-H0224	Human ADH1 (Alcohol Dehydrogenase 1) CLIA Kit
E-CL-H0225	Human ADH (Antidiuretic Hormone) CLIA Kit
E-CL-H0226	Human ADIPOR1 (Adiponectin Receptor 1) CLIA Kit
E-CL-H0227	Human ADIPOR2 (Adiponectin Receptor 2) CLIA Kit

Cat.NO	Product Name
E-CL-H0228	Human ADM (Adrenomedullin) CLIA Kit
E-CL-H0229	Human ADRA1A (Adrenergic Receptor alpha-1A) CLIA Kit
E-CL-H0230	Human ADRP (Adipose Differentiation Related Protein) CLIA Kit
E-CL-H0231	Human IL-20 (Interleukin 20) CLIA Kit
E-CL-H0232	Human DLK-1 (Delta Like 1 Homolog) CLIA Kit
E-CL-H0233	Human ACE2 (Angiotensin I Converting Enzyme 2) CLIA Kit
E-CL-H0234	Human BMP-3 (Bone Morphogenetic Protein 3) CLIA Kit
E-CL-H0235	Human CAPS (Calcyphosine) CLIA Kit
E-CL-H0236	Human FABP4 (Fatty Acid Binding Protein 4, Adipocyte) CLIA Kit
E-CL-H0237	Human AFP-L3 (Alpha-Fetoprotein Lens Culinaris Agglutinin 3) CLIA Kit
E-CL-H0238	Human FUCA (Alpha-L-Fucosidase, Tissue) CLIA Kit
E-CL-H0239	Human CCNB (Cyclin-B) CLIA Kit
E-CL-H0240	Human AGC (Aggrecan) CLIA Kit
E-CL-H0241	Human AGER (Advanced Glycosylation End Product Specific Receptor) CLIA Kit
E-CL-H0242	Human CD3d (Cluster of Differentiation 3d) CLIA Kit
E-CL-H0243	Human AGR2 (Anterior Gradient Protein 2) CLIA Kit
E-CL-H0244	Human AGRN (Agrin) CLIA Kit
E-CL-H0245	Human AGT (Angiotensinogen) CLIA Kit
E-CL-H0246	Human AIF1 (Allograft Inflammatory Factor 1) CLIA Kit

Cat.NO	Product Name
E-CL-H0247	Human ALCAM (Activated Leukocyte Cell Adhesion Molecule) CLIA Kit
E-CL-H0248	Human ALD (Aldosterone) CLIA Kit
E-CL-H0249	Human ALDH (Aldehyde Dehydrogenase, Mitochondrial) CLIA Kit
E-CL-H0250	Human ALDOA (Aldolase A, Fructose Bisphosphate) CLIA Kit
E-CL-H0251	Human Alk-Smase (Alkaline Sphingomyelinase) CLIA Kit
E-CL-H0252	Human ALT (Alanine Transaminase) CLIA Kit
E-CL-H0253	Human INHB (Inhibin B) CLIA Kit
E-CL-H0254	Human CDK6 (Cyclin Dependent Kinase 6) CLIA Kit
E-CL-H0255	Human AMBP (Alpha-1-Microglobulin/Bikunin Precursor) CLIA Kit
E-CL-H0256	Human CDK8 (Cyclin Dependent Kinase 8) CLIA Kit
E-CL-H0257	Human AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) CLIA Kit
E-CL-H0258	Human AMY1 (Amylase Alpha 1, Salivary) CLIA Kit
E-CL-H0259	Human IAPP (Islet Amyloid Polypeptide) CLIA Kit
E-CL-H0260	Human CDKN3 (Cyclin Dependent Kinase Inhibitor 3) CLIA Kit
E-CL-H0261	Human ANG (Angiostatin) CLIA Kit
E-CL-H0262	Human ANG II (Angiotensin II) CLIA Kit
E-CL-H0263	Human ANG I (Angiotensin I) CLIA Kit
E-CL-H0264	Human ANG II R1/AGTR1 (Angiotensin II Receptor 1) CLIA Kit
E-CL-H0265	Human ANG II R2/AGTR2 (Angiotensin II Receptor 2) CLIA Kit

Cat.NO	Product Name
E-CL-H0266	Human ANGPT4 (Angiopoietin 4) CLIA Kit
E-CL-H0267	Human ANGPTL1 (Angiopoietin Like Protein 1) CLIA Kit
E-CL-H0268	Human ANGPTL2 (Angiopoietin Like Protein 2) CLIA Kit
E-CL-H0269	Human ANGPTL3 (Angiopoietin Like Protein 3) CLIA Kit
E-CL-H0270	Human ANGPTL4 (Angiopoietin Like Protein 4) CLIA Kit
E-CL-H0271	Human ANGPTL6 (Angiopoietin Like Protein 6) CLIA Kit
E-CL-H0272	Human ANG-R-Tie1 (Angiopoietin Receptor Tie1) CLIA Kit
E-CL-H0273	Human ANG-R-Tie2 (Angiopoietin Receptor Tie2) CLIA Kit
E-CL-H0274	Human ANNA1(Anti-Neuronal Nuclear Antibody 1) CLIA Kit
E-CL-H0275	Human MAG (Myelin Associated Glycoprotein) CLIA Kit
E-CL-H0276	Human CDNF (Cerebral Dopamine Neurotrophic Factor) CLIA Kit
E-CL-H0277	Human FGFR1 (Fibroblast Growth Factor Receptor 1) CLIA Kit
E-CL-H0278	Human CNN1 (Calponin 1, Basic, Smooth Muscle) CLIA Kit
E-CL-H0279	Human Anti-AChR (Anti-Acetylcholine Receptor Antibody) CLIA Kit
E-CL-H0280	Human CNTNAP2 (Contactin Associated Protein Like Protein 2) CLIA Kit
E-CL-H0281	Human ACHAP (Acetylcholinesterase Associated Protein) CLIA Kit
E-CL-H0282	Human ACO1 (Aconitase 1) CLIA Kit
E-CL-H0283	Human CPE (Carboxypeptidase E) CLIA Kit
E-CL-H0284	Human Anti-GM1 (Anti-Ganglioside M1 Antibody) CLIA Kit

Cat.NO	Product Name
E-CL-H0285	Human CG α (Chorionic Gonadotropin Alpha Polypeptide) CLIA Kit
E-CL-H0286	Human ANG III (Angiotensin III) CLIA Kit
E-CL-H0287	Human RARS (Arginyl tRNA Synthetase, cytoplasmic) CLIA Kit
E-CL-H0288	Human FGFR3 (Fibroblast Growth Factor Receptor 3) CLIA Kit
E-CL-H0289	Human FGFR4 (Fibroblast Growth Factor Receptor 4) CLIA Kit
E-CL-H0290	Human DPP10 (Dipeptidyl Peptidase X) CLIA Kit
E-CL-H0291	Human CPR (Cytochrome P450 Reductase) CLIA Kit
E-CL-H0292	Human CRLF1 (Cytokine Receptor Like Factor 1) CLIA Kit
E-CL-H0293	Human CSF2R β (Colony Stimulating Factor 2 Receptor Beta) CLIA Kit
E-CL-H0294	Human CA3 (Carbonic Anhydrase III, Muscle Specific) CLIA Kit
E-CL-H0295	Human AIRA (Anti-Insulin Receptor Antibody) CLIA Kit
E-CL-H0296	Human CYP1A1 (Cytochrome P450, family 1, subfamily A, polypeptide 1) CLIA Kit
E-CL-H0297	Human CYP1A2 (Cytochrome P450, family 1, subfamily A, polypeptide 2) CLIA Kit
E-CL-H0298	Human DCLK1 (Doublecortin Like Kinase 1) CLIA Kit
E-CL-H0299	Human DPYD (Dihydropyrimidine Dehydrogenase) CLIA Kit
E-CL-H0300	Human DSE (Dermatan Sulfate Epimerase) CLIA Kit
E-CL-H0301	Human EPYC (Epiphykan) CLIA Kit
E-CL-H0302	Human FABP8 (Fatty Acid Binding Protein 8, Myelin) CLIA Kit
E-CL-H0303	Human FADS2 (Fatty Acid Desaturase 2) CLIA Kit

Cat.NO	Product Name
E-CL-H0304	Human FAIM3 (Fas Apoptotic Inhibitory Molecule 3) CLIA Kit
E-CL-H0305	Human FASN (Fatty Acid Synthase) CLIA Kit
E-CL-H0306	Human FBLN4 (Fibulin 4) CLIA Kit
E-CL-H0307	Human CYP24A1 (Cytochrome P450, family 24, subfamily A, polypeptide 1) CLIA Kit
E-CL-H0308	Human FCN1 (Ficolin 1) CLIA Kit
E-CL-H0309	Human AMA IgA (Anti-Myelin antibody IgA) CLIA Kit
E-CL-H0310	Human FcγBP (Fc Fragment of IgG Binding Protein) CLIA Kit
E-CL-H0311	Human FcγR3A (Fc Fragment of IgG Low Affinity IIIa Receptor) CLIA Kit
E-CL-H0312	Human CYP1B1 (Cytochrome P450, family 1, subfamily B, polypeptide 1) CLIA Kit
E-CL-H0313	Human FcεR II (Receptor II for the Fc Fragment of IgE) CLIA Kit
E-CL-H0314	Human FETUA (Fetuin A) CLIA Kit
E-CL-H0315	Human FGγ (Fibrinogen Gamma chain) CLIA Kit
E-CL-H0316	Human FLOT2 (Flotillin 2) CLIA Kit
E-CL-H0317	Human GC (Glycocalicin) CLIA Kit
E-CL-H0318	Human GLRX3 (Glutaredoxin 3) CLIA Kit
E-CL-H0319	Human GIF (Gastric Intrinsic Factor) CLIA Kit
E-CL-H0320	Human C4BPβ (C4 Binding Protein Beta) CLIA Kit
E-CL-H0321	Human CA4 (Carbonic Anhydrase IV) CLIA Kit
E-CL-H0322	Human C4BPα (C4 Binding Protein Alpha) CLIA Kit

Cat.NO	Product Name
E-CL-H0323	Human GPLD1 (Glycosylphosphatidylinositol Specific Phospholipase D1) CLIA Kit
E-CL-H0324	Human PRELP (Proline Arginine Rich End Leucine Rich Repeat Protein) CLIA Kit
E-CL-H0325	Human TARS (Threonyl tRNA Synthetase, cytoplasmic) CLIA Kit
E-CL-H0326	Human GTF2A1 (General Transcription Factor IIA, Polypeptide 1) CLIA Kit
E-CL-H0327	Human GUCA2A (Guanylate Cyclase Activator 2A) CLIA Kit
E-CL-H0328	Human INSL3 (Insulin Like Protein 3) CLIA Kit
E-CL-H0329	Human INSL5 (Insulin Like Protein 5) CLIA Kit
E-CL-H0330	Human iPLA2 (Phospholipidase A2, Calcium Independent) CLIA Kit
E-CL-H0331	Human GUCA2B (Guanylate Cyclase Activator 2B) CLIA Kit
E-CL-H0332	Human GZMK (Granzyme K) CLIA Kit
E-CL-H0333	Human ALOXE3 (Arachidonate Lipoxygenase 3) CLIA Kit
E-CL-H0334	Human CYP26A1 (Cytochrome P450, family 26, subfamily A, polypeptide 1) CLIA Kit
E-CL-H0335	Human HB (Hemoglobin) CLIA Kit
E-CL-H0336	Human HCCS (Holocytochrome C Synthase) CLIA Kit
E-CL-H0337	Human CAMK2 α (Calcium/Calmodulin-Dependent Protein Kinase II Alpha) CLIA Kit
E-CL-H0338	Human BAG3 (Bcl2 Associated Athanogene 3) CLIA Kit
E-CL-H0339	Human ARSB (Arylsulfatase B) CLIA Kit
E-CL-H0340	Human HRG (Histidine Rich Glycoprotein) CLIA Kit
E-CL-H0341	Human HSPG (Heparan Sulfate Proteoglycan) CLIA Kit

Cat.NO	Product Name
E-CL-H0342	Human ANXA5 (Annexin A5) CLIA Kit
E-CL-H0343	Human HTRA1 (HtrA Serine Peptidase 1) CLIA Kit
E-CL-H0344	Human IL10R β (Interleukin 10 Receptor Beta) CLIA Kit
E-CL-H0345	Human IGF1R (Insulin Like Growth Factor 1 Receptor) CLIA Kit
E-CL-H0346	Human HYOU1 (Hypoxia Up Regulated 1) CLIA Kit
E-CL-H0347	Human IAA (Insulin Autoantibody) CLIA Kit
E-CL-H0348	Human ICAM4 (Intercellular Adhesion Molecule 4) CLIA Kit
E-CL-H0349	Human IDE (Insulin Degrading Enzyme) CLIA Kit
E-CL-H0350	Human Anti-Surv (Anti-Survivin Antibody) CLIA Kit
E-CL-H0351	Human TA (Anti-Teichoic Acid Antibody) CLIA Kit
E-CL-H0352	Human AT (AntiThrombin) CLIA Kit
E-CL-H0353	Human LIPH (Lipase H) CLIA Kit
E-CL-H0354	Human LT β R (Lymphotoxin Beta Receptor) CLIA Kit
E-CL-H0355	Human IGF-2 (Insulin Like Growth Factor 2) CLIA Kit
E-CL-H0356	Human TPO-Ab (Anti-Thyroid-Peroxidase Antibody) CLIA Kit
E-CL-H0357	Human IGF2BP3 (Insulin Like Growth Factor 2 mRNA Binding Protein 3) CLIA Kit
E-CL-H0358	Human TTN (Titin) CLIA Kit
E-CL-H0359	Human IGF2R (Insulin Like Growth Factor 2 Receptor) CLIA Kit
E-CL-H0360	Human IGFALS (Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) CLIA Kit

Cat.NO	Product Name
E-CL-H0361	Human IGFBP-1 (Insulin Like Growth Factor Binding Protein 1) CLIA Kit
E-CL-H0362	Human IGFBP-5 (Insulin Like Growth Factor Binding Protein 5) CLIA Kit
E-CL-H0363	Human IGFBP-4 (Insulin Like Growth Factor Binding Protein 4) CLIA Kit
E-CL-H0364	Human IGFBP-2 (Insulin Like Growth Factor Binding Protein 2) CLIA Kit
E-CL-H0365	Human IGFBP-7 (Insulin Like Growth Factor Binding Protein 7) CLIA Kit
E-CL-H0366	Human ANXA2 (Annexin A2) CLIA Kit
E-CL-H0367	Human IGFBP-6 (Insulin Like Growth Factor Binding Protein 6) CLIA Kit
E-CL-H0368	Human INSR (Insulin Receptor) CLIA Kit
E-CL-H0369	Human APC (Activated Protein C) CLIA Kit
E-CL-H0370	Human APLN (Apelin) CLIA Kit
E-CL-H0371	Human AP13 (Apelin 13) CLIA Kit
E-CL-H0372	Human AP36 (Apelin 36) CLIA Kit
E-CL-H0373	Human Apo A1-Ab (Anti-Apolipoprotein A1 Antibody) CLIA Kit
E-CL-H0374	Human LCN1 (Lipocalin 1) CLIA Kit
E-CL-H0375	Human ApoA2 (Apolipoprotein A2) CLIA Kit
E-CL-H0376	Human ApoA4 (Apolipoprotein A4) CLIA Kit
E-CL-H0377	Human ApoB (Apolipoprotein B) CLIA Kit
E-CL-H0378	Human ApoC1 (Apolipoprotein C1) CLIA Kit
E-CL-H0379	Human ApoC2 (Apolipoprotein C2) CLIA Kit

Cat.NO	Product Name
E-CL-H0380	Human ApoC3 (Apolipoprotein C3) CLIA Kit
E-CL-H0381	Human ApoC4 (Apolipoprotein C4) CLIA Kit
E-CL-H0382	Human ApoD (Apolipoprotein D) CLIA Kit
E-CL-H0383	Human ApoE (Apolipoprotein E) CLIA Kit
E-CL-H0384	Human ApoF (Apolipoprotein F) CLIA Kit
E-CL-H0385	Human ApoH (Apolipoprotein H) CLIA Kit
E-CL-H0386	Human ApoM (Apolipoprotein M) CLIA Kit
E-CL-H0387	Human CA6 (Carbonic Anhydrase VI) CLIA Kit
E-CL-H0388	Human ApoA5 (Apolipoprotein A5) CLIA Kit
E-CL-H0389	Human ApoB100 (Apolipoprotein B100) CLIA Kit
E-CL-H0390	Human LIPC (Lipase, Hepatic) CLIA Kit
E-CL-H0391	Human AIF (Apoptosis Inducing Factor) CLIA Kit
E-CL-H0392	Human LIPF (Lipase, Gastric) CLIA Kit
E-CL-H0393	Human APAF1 (Apoptosis Protease Activating Factor 1) CLIA Kit
E-CL-H0394	Human bFGF/FGF2 (Basic Fibroblast Growth Factor) CLIA Kit
E-CL-H0395	Human AP (Aprotinin) CLIA Kit
E-CL-H0396	Human ITIH4 (Inter Alpha-Globulin Inhibitor H4) CLIA Kit
E-CL-H0397	Human AQP-1 (Aquaporin 1) CLIA Kit
E-CL-H0398	Human AQP-2 (Aquaporin 2) CLIA Kit

Cat.NO	Product Name
E-CL-H0399	Human AQP-3 (Aquaporin 3) CLIA Kit
E-CL-H0400	Human AQP-4 (Aquaporin 4) CLIA Kit
E-CL-H0401	Human AR (Androgen Receptor) CLIA Kit
E-CL-H0402	Human AA (Arachidonic Acid) CLIA Kit
E-CL-H0403	Human Arg (Arginase) CLIA Kit
E-CL-H0404	Human KLK6 (Kallikrein 6) CLIA Kit
E-CL-H0405	Human KLK7 (Kallikrein 7) CLIA Kit
E-CL-H0406	Human ARHGDI α (Rho GDP Dissociation Inhibitor Alpha) CLIA Kit
E-CL-H0407	Human ARR β 2 (Arrestin Beta 2) CLIA Kit
E-CL-H0408	Human MANF (Mesencephalic Astrocyte Derived Neurotrophic Factor) CLIA Kit
E-CL-H0409	Human MAP2K1/MKK1 (Mitogen Activated Protein Kinase Kinase 1) CLIA Kit
E-CL-H0410	Human AhR (Aryl Hydrocarbon Receptor) CLIA Kit
E-CL-H0411	Human ARSA (Arylsulfatase A) CLIA Kit
E-CL-H0412	Human ASD (Androstenedione) CLIA Kit
E-CL-H0413	Human ASGR2 (Asialoglycoprotein Receptor 2) CLIA Kit
E-CL-H0414	Human ASGR1 (Asialoglycoprotein Receptor 1) CLIA Kit
E-CL-H0415	Human ASM (Acid Sphingomyelinase) CLIA Kit
E-CL-H0416	Human ASP (Acylation Stimulating Protein) CLIA Kit
E-CL-H0417	Human MAPK14 (Mitogen Activated Protein Kinase 14) CLIA Kit

Cat.NO	Product Name
E-CL-H0418	Human AST (Aspartate Aminotransferase) CLIA Kit
E-CL-H0419	Human ASPN (Asporin) CLIA Kit
E-CL-H0420	Human ADMA (Asymmetrical Dimethylarginine) CLIA Kit
E-CL-H0421	Human ATM (Ataxia Telangiectasia Mutated) CLIA Kit
E-CL-H0422	Human ATXN1 (Ataxin 1) CLIA Kit
E-CL-H0423	Human ABCG1 (ATP Binding Cassette Transporter G1) CLIA Kit
E-CL-H0424	Human ABCA1(ATP-binding cassette transporter A1) CLIA Kit
E-CL-H0425	Human ABCG2 (ATP Binding Cassette Transporter G2) CLIA Kit
E-CL-H0426	Human ANP (Atrial Natriuretic Peptide) CLIA Kit
E-CL-H0427	Human MGL (Monoacylglycerol Lipase) CLIA Kit
E-CL-H0428	Human AIRE (Autoimmune Regulator) CLIA Kit
E-CL-H0429	Human AXIN2 (Axis Inhibition Protein 2) CLIA Kit
E-CL-H0430	Human AZU1 (Azurocidin 1) CLIA Kit
E-CL-H0431	Human A β 40 (Amyloid Beta 40) CLIA Kit
E-CL-H0432	Human A β 42 (Amyloid Beta 42) CLIA Kit
E-CL-H0433	Human BAFFR (B Cell Activating Factor Receptor) CLIA Kit
E-CL-H0434	Human BCL2L1/BCL-X (Bcl-2 Like Protein 1/Bcl2 Associated X Protein) CLIA Kit
E-CL-H0435	Human BPI (Bactericidal/Permeability Increasing Protein) CLIA Kit
E-CL-H0436	Human FGF4 (Fibroblast Growth Factor 4) CLIA Kit

Cat.NO	Product Name
E-CL-H0437	Human FGF6 (Fibroblast Growth Factor 6) CLIA Kit
E-CL-H0438	Human LDHA (Lactate Dehydrogenase A) CLIA Kit
E-CL-H0439	Human BAK1 (Bcl2 Antagonist/Killer 1) CLIA Kit
E-CL-H0440	Human Bax (Bcl-2 Associated X Protein) CLIA Kit
E-CL-H0441	Human BECN1 (Beclin 1)CLIA Kit
E-CL-H0442	Human AQP-8 (Aquaporin 8) CLIA Kit
E-CL-H0443	Human β 2-GP (Beta2 Glycoprotein) CLIA Kit
E-CL-H0444	Human bTC (Betacellulin) CLIA Kit
E-CL-H0445	Human β -EPR (Beta-Endorphin Receptor) CLIA Kit
E-CL-H0446	Human β -EP (Beta-Endorphin) CLIA Kit
E-CL-H0447	Human Bid (BH3 Interacting Domain Death Agonist) CLIA Kit
E-CL-H0448	Human MMRN1 (Multimerin 1) CLIA Kit
E-CL-H0449	Human MFAP4 (Microfibrillar Associated Protein 4) CLIA Kit
E-CL-H0450	Human Bcl2L11 (Bcl-2 Like Protein 11) CLIA Kit
E-CL-H0451	Human BMG/ β 2-MG (Beta-2-Microglobulin) CLIA Kit
E-CL-H0452	Human BALP (Bone Alkaline Phosphatase) CLIA Kit
E-CL-H0453	Human BMP-1 (Bone Morphogenetic Protein 1) CLIA Kit
E-CL-H0454	Human BMP-10 (Bone Morphogenetic Protein 10) CLIA Kit
E-CL-H0455	Human BMP-15 (Bone Morphogenetic Protein 15) CLIA Kit

Cat.NO	Product Name
E-CL-H0456	Human BMP-6 (Bone Morphogenetic Protein 6) CLIA Kit
E-CL-H0457	Human BSP (Bone Sialoprotein) CLIA Kit
E-CL-H0458	Human MOG (Myelin Oligodendrocyte Glycoprotein) CLIA Kit
E-CL-H0459	Human MT1E (Metallothionein 1E) CLIA Kit
E-CL-H0460	Human BK (Bradykinin) CLIA Kit
E-CL-H0461	Human BNP (Brain Natriuretic Peptide) CLIA Kit
E-CL-H0462	Human MSM β (Microseminoprotein Beta) CLIA Kit
E-CL-H0463	Human BRAK/CXCL14 (Breast and Kidney Expressed Chemokine) CLIA Kit
E-CL-H0464	Human BRCA1 (Breast Cancer Susceptibility Protein 1) CLIA Kit
E-CL-H0465	Human BRCA2 (Breast Cancer Susceptibility Protein 2) CLIA Kit
E-CL-H0466	Human BCAN (Brevican) CLIA Kit
E-CL-H0467	Human MT1M (Metallothionein 1M) CLIA Kit
E-CL-H0468	Human β TG/PBP/CXCL7/NAP2 (Thromboglobulin, Beta) CLIA Kit
E-CL-H0469	Human MT3 (Metallothionein 3) CLIA Kit
E-CL-H0470	Human CNP (C-type Natriuretic Peptide) CLIA Kit
E-CL-H0471	Human MUC13 (Mucin 13) CLIA Kit
E-CL-H0472	Human CA724 (Tumor Marker CA724) CLIA Kit
E-CL-H0473	Human MUC1 (Mucin 1) CLIA Kit
E-CL-H0474	Human CaN (Calcineurin) CLIA Kit

Cat.NO	Product Name
E-CL-H0475	Human CGRP (Calcitonin Gene Related Peptide) CLIA Kit
E-CL-H0476	Human CACYBP (Calcyclin Binding Protein) CLIA Kit
E-CL-H0477	Human CALD (Caldesmon) CLIA Kit
E-CL-H0478	Human MSR1 (Macrophage Scavenger Receptor 1) CLIA Kit
E-CL-H0479	Human CAM (Calmodulin) CLIA Kit
E-CL-H0480	Human CRT (Calreticulin) CLIA Kit
E-CL-H0481	Human CR (Calretinin) CLIA Kit
E-CL-H0482	Human CALU (Calumenin) CLIA Kit
E-CL-H0483	Human MUC17 (Mucin 17) CLIA Kit
E-CL-H0484	Human MUC2 (Mucin 2) CLIA Kit
E-CL-H0485	Human MUC20 (Mucin 20) CLIA Kit
E-CL-H0486	Human CNR1 (Cannabinoid Receptor 1, Brain) CLIA Kit
E-CL-H0487	Human CA125 (Carbohydrate Antigen 125) CLIA Kit
E-CL-H0488	Human CA19-9 (Carbohydrate antigen19-9) CLIA Kit
E-CL-H0489	Human CA1 (Carbonic Anhydrase I) CLIA Kit
E-CL-H0490	Human CA2 (Carbonic Anhydrase II) CLIA Kit
E-CL-H0491	Human HbCO (Carboxyhemoglobin) CLIA Kit
E-CL-H0492	Human CEA (Carcinoembryonic Antigen) CLIA Kit
E-CL-H0493	Human cTnT/TNNT2 (Troponin T Type 2, Cardiac) CLIA Kit

Cat.NO	Product Name
E-CL-H0494	Human IL-7 (Interleukin 7) CLIA Kit
E-CL-H0495	Human cTn-I/TNNI3 (cardiac Troponin I)CLIA Kit
E-CL-H0496	Human CT-1 (Cardiotrophin 1) CLIA Kit
E-CL-H0497	Human CLCF1 (Cardiotrophin Like Cytokine Factor 1) CLIA Kit
E-CL-H0498	Human PAR1 (Protease Activated Receptor 1) CLIA Kit
E-CL-H0499	Human COMP (Cartilage Oligomeric Matrix Protein) CLIA Kit
E-CL-H0500	Human CASP14 (Caspase 14) CLIA Kit
E-CL-H0501	Human CASP8 (Caspase 8) CLIA Kit
E-CL-H0502	Human CASP4 (Caspase 4) CLIA Kit
E-CL-H0503	Human CASP7 (Caspase 7) CLIA Kit
E-CL-H0504	Human CASP 9 (Caspase 9) CLIA Kit
E-CL-H0505	Human CTNNβ1 (Catenin, Beta 1) CLIA Kit
E-CL-H0506	Human CTSD (Cathepsin D) CLIA Kit
E-CL-H0507	Human CTSK (Cathepsin K) CLIA Kit
E-CL-H0508	Human CTSL (Cathepsin L) CLIA Kit
E-CL-H0509	Human Cav-1 (Caveolin-1) CLIA Kit
E-CL-H0510	Human CD59 (Protectin) CLIA Kit
E-CL-H0511	Human CA5B (Carbonic Anhydrase VB, Mitochondrial) CLIA Kit
E-CL-H0512	Human CFL1(Cofilin 1, Non-Muscle) CLIA Kit

Cat.NO	Product Name
E-CL-H0513	Human MUC21 (Mucin 21) CLIA Kit
E-CL-H0514	Human CHEM (Chemerin) CLIA Kit
E-CL-H0515	Human CCL14 (Chemokine C-C-Motif Ligand 14) CLIA Kit
E-CL-H0516	Human CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) CLIA Kit
E-CL-H0517	Human TERT (Telomerase Reverse Transcriptase) CLIA Kit
E-CL-H0518	Human MUC3 (Mucin 3) CLIA Kit
E-CL-H0519	Human MUC4 (Mucin 4) CLIA Kit
E-CL-H0520	Human MUC6 (Mucin 6) CLIA Kit
E-CL-H0521	Human MYBPC3 (Myosin Binding Protein C, Cardiac) CLIA Kit
E-CL-H0522	Human MYH10 (Myosin Heavy Chain 10, Non Muscle) CLIA Kit
E-CL-H0523	Human MYH9 (Myosin Heavy Chain 9, Non Muscle) CLIA Kit
E-CL-H0524	Human MYL9 (Myosin Light Chain 9, Regulatory) CLIA Kit
E-CL-H0525	Human CCKAR (Cholecystokinin A Receptor) CLIA Kit
E-CL-H0526	Human CCK-8 (Cholecystokinin 8) CLIA Kit
E-CL-H0527	Human CETP (Cholesteryl Ester Transfer Protein) CLIA Kit
E-CL-H0528	Human CH (Cholesterol) CLIA Kit
E-CL-H0529	Human CENPF (Centromere Protein F) CLIA Kit
E-CL-H0530	Human CHRM2 (Cholinergic Receptor, Muscarinic 2) CLIA Kit
E-CL-H0531	Human CS (Chondroitin Sulfate) CLIA Kit

Cat.NO	Product Name
E-CL-H0532	Human CG β (Chorionic Gonadotrophin Beta) CLIA Kit
E-CL-H0533	Human CA12 (Carbonic Anhydrase XII) CLIA Kit
E-CL-H0534	Human CgA (Chromogranin A/Pancreastatin) CLIA Kit
E-CL-H0535	Human NEF3 (Neurofilament 3) CLIA Kit
E-CL-H0536	Human NEFL (Neurofilament, Light Polypeptide) CLIA Kit
E-CL-H0537	Human NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) CLIA Kit
E-CL-H0538	Human CLDN1 (Claudin 1) CLIA Kit
E-CL-H0539	Human CD99 (Cluster of Differentiation 99) CLIA Kit
E-CL-H0540	Human CD4 (Cluster of Differentiation 4) CLIA Kit
E-CL-H0541	Human NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) CLIA Kit
E-CL-H0542	Human CLDN3 (Claudin 3) CLIA Kit
E-CL-H0543	Human NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) CLIA Kit
E-CL-H0544	Human F2 (Coagulation Factor II) CLIA Kit
E-CL-H0545	Human F7 (Coagulation Factor VII) CLIA Kit
E-CL-H0546	Human F9 (Coagulation Factor IX) CLIA Kit
E-CL-H0547	Human F10 (Coagulation Factor X) CLIA Kit
E-CL-H0548	Human F13 (Coagulation Factor XIII) CLIA Kit
E-CL-H0549	Human F12 (Coagulation Factor XII) CLIA Kit
E-CL-H0550	Human F11 (Coagulation Factor XI) CLIA Kit

Cat.NO	Product Name
E-CL-H0551	Human F13A1 (Coagulation Factor XIII A1 Polypeptide) CLIA Kit
E-CL-H0552	Human COL15 (Collagen Type XV) CLIA Kit
E-CL-H0553	Human COL2 (Collagen Type II) CLIA Kit
E-CL-H0554	Human COL3 (Collagen Type III) CLIA Kit
E-CL-H0555	Human COL6 (Collagen Type VI) CLIA Kit
E-CL-H0556	Human COL1 α 2 (Collagen Type I Alpha 2) CLIA Kit
E-CL-H0557	Human COL2 α 1 (Collagen Type II Alpha 1) CLIA Kit
E-CL-H0558	Human COL3 α 1 (Collagen Type III Alpha 1) CLIA Kit
E-CL-H0559	Human COL4 α 1 (Collagen Type IV Alpha 1) CLIA Kit
E-CL-H0560	Human COL6 α 1 (Collagen Type VI Alpha 1) CLIA Kit
E-CL-H0561	Human COL7 (Collagen Type VII) CLIA Kit
E-CL-H0562	Human COL8 α 1 (Collagen Type VIII Alpha 1) CLIA Kit
E-CL-H0563	Human COL10 (Collagen Type X) CLIA Kit
E-CL-H0564	Human COL11 α 1 (Collagen Type XI Alpha 1) CLIA Kit
E-CL-H0565	Human COL17 (Collagen Type XVII) CLIA Kit
E-CL-H0566	Human Collagenase I CLIA Kit
E-CL-H0567	Human NPYR1 (Neuropeptide Y Receptor Y1) CLIA Kit
E-CL-H0568	Human CLDN9 (Claudin 9) CLIA Kit
E-CL-H0569	Human PIINP (Procollagen II N-Terminal ProPeptide) CLIA Kit

Cat.NO	Product Name
E-CL-H0570	Human C1INH (Complement 1 Inhibitor) CLIA Kit
E-CL-H0571	Human C1q (Complement 1q) CLIA Kit
E-CL-H0572	Human C3 (Complement Component 3) CLIA Kit
E-CL-H0573	Human CLDN4 (Claudin 4) CLIA Kit
E-CL-H0574	Human NPYR2 (Neuropeptide Y Receptor Y2) CLIA Kit
E-CL-H0575	Human C4b (Complement Component 4b) CLIA Kit
E-CL-H0576	Human C5c (Complement C5 Convertase) CLIA Kit
E-CL-H0577	Human NPYR5 (Neuropeptide Y Receptor Y5) CLIA Kit
E-CL-H0578	Human ConA (Concanavalin A) CLIA Kit
E-CL-H0579	Human CTGF (Connective Tissue Growth Factor) CLIA Kit
E-CL-H0580	Human CX26 (Connexin 26) CLIA Kit
E-CL-H0581	Human CHGB (Chromogranin B) CLIA Kit
E-CL-H0582	Human DEF α 1 (Defensin Alpha 1, Neutrophil) CLIA Kit
E-CL-H0583	Human CTXI (Cross Linked C-telopeptide of Type I Collagen) CLIA Kit
E-CL-H0584	Human NTXI (Cross Linked N-telopeptide of Type I Collagen) CLIA Kit
E-CL-H0585	Human CTX-II (Cross Linked C-telopeptide of Type II Collagen) CLIA Kit
E-CL-H0586	Human LDHB (Lactate Dehydrogenase B) CLIA Kit
E-CL-H0587	Human P III CP (Procollagen III C-Terminal ProPeptide) CLIA Kit
E-CL-H0588	Human DKK4 (Dickkopf Related Protein 4) CLIA Kit

Cat.NO	Product Name
E-CL-H0589	Human LDH (Lactate Dehydrogenase) CLIA Kit
E-CL-H0590	Human CPP (Copeptin) CLIA Kit
E-CL-H0591	Human CLDN11 (Claudin 11) CLIA Kit
E-CL-H0592	Human CXCR3 (Chemokine CXC Receptor 3) CLIA Kit
E-CL-H0593	Human SCG2 (Secretogranin II) CLIA Kit
E-CL-H0594	Human MYCBP (C-Myc Binding Protein) CLIA Kit
E-CL-H0595	Human MYC (V-Myc Myelocytomatosis Viral Oncogene Homolog) CLIA Kit
E-CL-H0596	Human LDHC (Lactate Dehydrogenase C) CLIA Kit
E-CL-H0597	Human EphA2 (Tyrosine protein kinase receptor A2) CLIA Kit
E-CL-H0598	Human DKK2 (Dickkopf Related Protein 2) CLIA Kit
E-CL-H0599	Human D2D (D-Dimer) CLIA Kit
E-CL-H0600	Human F8 (Coagulation Factor VIII) CLIA Kit
E-CL-H0601	Human PA2G4 (Proliferation Associated Protein 2G4) CLIA Kit
E-CL-H0602	Human COL1 α 1 (Collagen Type I Alpha 1) CLIA Kit
E-CL-H0603	Human COL6 α 3 (Collagen Type VI Alpha 3) CLIA Kit
E-CL-H0604	Human COL13 (Collagen Type XIII) CLIA Kit
E-CL-H0605	Human COL14 (Collagen Type XIV) CLIA Kit
E-CL-H0606	Human CTX-III (Cross Linked C-telopeptide of Type I Collagen) CLIA Kit
E-CL-H0607	Human EFNA1 (Ephrin A1) CLIA Kit

Cat.NO	Product Name
E-CL-H0608	Human SELE (E-Selectin) CLIA Kit
E-CL-H0609	Human GATA3 (GATA Binding Protein 3) CLIA Kit
E-CL-H0610	Human GN β 1 (G Protein Beta 1) CLIA Kit
E-CL-H0611	Human HHIP (Hedgehog Interacting Protein) CLIA Kit
E-CL-H0612	Human H-Cadherin (Cadherin, Heart) CLIA Kit
E-CL-H0613	Human NRG-4 (Neuregulin 4) CLIA Kit
E-CL-H0614	Human Fc α R (Fc Fragment of IgA Receptor) CLIA Kit
E-CL-H0615	Human SELL (L-Selectin) CLIA Kit
E-CL-H0616	Human HB μ (Hemoglobin, Mu) CLIA Kit
E-CL-H0617	Human NAGase (N-Acetyl Beta-D-Glucosaminidase) CLIA Kit
E-CL-H0618	Human NRG-2 (Neuregulin 2) CLIA Kit
E-CL-H0619	Human NRG-3 (Neuregulin 3) CLIA Kit
E-CL-H0620	Human NT-ProBNP (N-Terminal Pro-Brain Natriuretic Peptide) CLIA Kit
E-CL-H0621	Human TP53 (Tumor Protein p53) CLIA Kit
E-CL-H0622	Human PDCN (Podocin) CLIA Kit
E-CL-H0623	Human Preptin CLIA Kit
E-CL-H0624	Human P-Cadherin (Cadherin, Placental) CLIA Kit
E-CL-H0625	Human P-GP (Permeability Glycoprotein) CLIA Kit
E-CL-H0626	Human SP (Substance P) CLIA Kit

Cat.NO	Product Name
E-CL-H0627	Human SELP (P-Selectin) CLIA Kit
E-CL-H0628	Human Rac-GAP1 (Rac-GTPase Activating Protein 1) CLIA Kit
E-CL-H0629	Human Rock-2 (Rho Associated Coiled Coil Containing Protein Kinase 2) CLIA Kit
E-CL-H0630	Human RUNX2 (Runt-Related Transcription Factor 2) CLIA Kit
E-CL-H0631	Human COL9α3 (Collagen Type IX Alpha 3) CLIA Kit
E-CL-H0632	Human S-100 (Soluble Protein-100) CLIA Kit
E-CL-H0633	Human Salusin α (Salusin Alpha) CLIA Kit
E-CL-H0634	Human SHH (Hedgehog Homolog, Sonic) CLIA Kit
E-CL-H0635	Human Slit2 (Slit Homolog 2) CLIA Kit
E-CL-H0636	Human Smad1/MADH1 (Mothers Against Decapentaplegic Homolog 1) CLIA Kit
E-CL-H0637	Human Smad4/MADH4 (Mothers Against Decapentaplegic Homolog 4)CLIA Kit
E-CL-H0638	Human Smad7/MADH7 (Mothers Against Decapentaplegic Homolog 7) CLIA Kit
E-CL-H0639	Human PCDHβ16 (Protocadherin Beta 16) CLIA Kit
E-CL-H0640	Human COL5α2 (Collagen Type V Alpha 2) CLIA Kit
E-CL-H0641	Human THP (Tamm–Horsfall Glycoprotein) CLIA Kit
E-CL-H0642	Human TDP43 (Tar DNA Binding Protein 43KDa)CLIA Kit
E-CL-H0643	Human TAF (TATA Box Binding Protein/TBP-Associated Factors) CLIA Kit
E-CL-H0644	Human TAF13 (TATA Box Binding Protein Associated Factor 13) CLIA Kit
E-CL-H0645	Human MAPτ (Microtubule Associated Protein Tau/Tau Protein) CLIA Kit

Cat.NO	Product Name
E-CL-H0646	Human TIEG1 (TGF Beta Inducible Early Response Gene 1) CLIA Kit
E-CL-H0647	Human TLR2 (Toll-Like Receptor 2) CLIA Kit
E-CL-H0648	Human TLR9 (Toll-Like Receptor 9) CLIA Kit
E-CL-H0649	Human LAT (Linker For Activation of T-cell) CLIA Kit
E-CL-H0650	Human β -CTx (Beta Crosslaps) CLIA Kit
E-CL-H0651	Human C3c (Complement C3 Convertase) CLIA Kit
E-CL-H0652	Human PCDH15 (Protocadherin 15) CLIA Kit
E-CL-H0653	Human α 1BG (Alpha-1-B-Glycoprotein) CLIA Kit
E-CL-H0654	Human α 1-AT (Alpha 1-Antitrypsin) CLIA Kit
E-CL-H0655	Human PLC γ 2 (Phospholipase C Gamma 2, Phosphatidylinositol Specific) CLIA Kit
E-CL-H0656	Human α 2-M (Alpha-2 Macroglobulin) CLIA Kit
E-CL-H0657	Human α 2-PI (α 2-Plasmin Inhibitor) CLIA Kit
E-CL-H0658	Human IDU α (α -L-Iduronidase) CLIA Kit
E-CL-H0659	Human NAGLU (N-Acetyl Alpha-D-Glucosaminidase) CLIA Kit
E-CL-H0660	Human GL α (Galactosidase Alpha) CLIA Kit
E-CL-H0661	Human α -SMA (Alpha-Smooth Muscle Actin) CLIA Kit
E-CL-H0662	Human α -La (Alpha-Lactalbumin) CLIA Kit
E-CL-H0663	Human SNC α (Synuclein Alpha) CLIA Kit
E-CL-H0664	Human PDGFR β (Platelet Derived Growth Factor Receptor Beta) CLIA Kit

Cat.NO	Product Name
E-CL-H0665	Human β 2-GP1 Ab IgA (β 2-Glycoprotein 1 antibody IgA) CLIA Kit
E-CL-H0666	Human PIICP (Procollagen II C-Terminal ProPeptide) CLIA Kit
E-CL-H0667	Human PCDHyA2 (Protocadherin Gamma A2) CLIA Kit
E-CL-H0668	Human GL β (Galactosidase Beta) CLIA Kit
E-CL-H0669	Human MAN1A1 (Mannosidase Alpha Class 1A Member 1) CLIA Kit
E-CL-H0670	Human DEFB1 (Defensin Beta 1) CLIA Kit
E-CL-H0671	Human DEF β 2/DEFB2 (Defensin Beta 2) CLIA Kit
E-CL-H0672	Human β -actin (Actin Beta) CLIA Kit
E-CL-H0673	Human CYSLTR1 (Cysteinyl Leukotriene Receptor 1) CLIA Kit
E-CL-H0674	Human LACT β (β -Lactamase) CLIA Kit
E-CL-H0675	Human IRAK3 (Interleukin 1 Receptor Associated Kinase 3) CLIA Kit
E-CL-H0676	Human GUS β (β -glucuronidase) CLIA Kit
E-CL-H0677	Human SYCN (Syncollin) CLIA Kit
E-CL-H0678	Human TUB β (β -Tubulin) CLIA Kit
E-CL-H0679	Human β -TG (β -Thromboglobulin) CLIA Kit
E-CL-H0680	Human Mly/CXCL9 (Monocyte Interferon Gamma Inducing Factor) CLIA Kit
E-CL-H0681	Human PLA2G10 (Phospholipase A2, Group X) CLIA Kit
E-CL-H0682	Human γ GT1 (Gamma Glutamyltransferase 1) CLIA Kit
E-CL-H0683	Human POMC (Pro-Opiomelanocortin) CLIA Kit

Cat.NO	Product Name
E-CL-H0684	Human OXA (Orexin A) CLIA Kit
E-CL-H0685	Human PIK3C β (Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) CLIA Kit
E-CL-H0686	Human OIT3 (Oncoprotein Induced Transcript 3) CLIA Kit
E-CL-H0687	Human PLD (Phospholipase D) CLIA Kit
E-CL-H0688	Human PIK3C σ (Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) CLIA Kit
E-CL-H0689	Human PRDX5 (Peroxiredoxin 5) CLIA Kit
E-CL-H0690	Human IL1R1 (Interleukin 1 Receptor Type I) CLIA Kit
E-CL-H0691	Human PLA2G12B (Phospholipase A2, Group XII B) CLIA Kit
E-CL-H0692	Human IL-2R α (Interleukin 2 Receptor Alpha) CLIA Kit
E-CL-H0693	Human IL-2R β (Interleukin 2 Receptor Beta) CLIA Kit
E-CL-H0694	Human LTB4 (Leukotriene B4) CLIA Kit
E-CL-H0695	Human LTC4 (Leukotriene C4) CLIA Kit
E-CL-H0696	Human LTD4 (Leukotriene D4) CLIA Kit
E-CL-H0697	Human LILRB4 (Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) CLIA Kit
E-CL-H0698	Human LIFR (Leukemia Inhibitory Factor Receptor) CLIA Kit
E-CL-H0699	Human GARS (Glycyl tRNA Synthetase) CLIA Kit
E-CL-H0700	Human EPHA1 (Ephrin Type A Receptor 1) CLIA Kit
E-CL-H0701	Human GP4/CD36 (Platelet Membrane Glycoprotein IV) CLIA Kit
E-CL-H0702	Human NSE (Neuron-Specific Enolase) CLIA Kit

Cat.NO	Product Name
E-CL-H0703	Human PLA2G4D (Phospholipase A2, Group IV D) CLIA Kit
E-CL-H0704	Human GAL1 (Galectin 1) CLIA Kit
E-CL-H0705	Human GAL12 (Galectin 12) CLIA Kit
E-CL-H0706	Human GAL14 (Galectin 14) CLIA Kit
E-CL-H0707	Human GAL2 (Galectin 2) CLIA Kit
E-CL-H0708	Human EPHA3 (Ephrin Type A Receptor 3) CLIA Kit
E-CL-H0709	Human GAL4 (Galectin 4) CLIA Kit
E-CL-H0710	Human GAL8 (Galectin 8) CLIA Kit
E-CL-H0711	Human GAL7 (Galectin 7) CLIA Kit
E-CL-H0712	Human GAL9 (Galectin 9) CLIA Kit
E-CL-H0713	Human PLA2R1 (Phospholipase A2 Receptor 1) CLIA Kit
E-CL-H0714	Human PLB (Phospholipase B) CLIA Kit
E-CL-H0715	Human cPLA2 (Cytosolic Phospholipase A2) CLIA Kit
E-CL-H0716	Human EFNA3 (Ephrin A3) CLIA Kit
E-CL-H0717	Human PLD2 (Phospholipase D2) CLIA Kit
E-CL-H0718	Human PYD (Pyridinoline) CLIA Kit
E-CL-H0719	Human PDH β (Pyruvate Dehydrogenase Beta) CLIA Kit
E-CL-H0720	Human OCLN (Occludin) CLIA Kit
E-CL-H0721	Human SEMA7A (Semaphorin 7A) CLIA Kit

Cat.NO	Product Name
E-CL-H0722	Human PLVAP (Plasmalemma Vesicle Associated Protein) CLIA Kit
E-CL-H0723	Human IGF1R3 (Insulin Like Growth Factor 1 Receptor 3) CLIA Kit
E-CL-H0724	Human SPJ (Surfactant Associated Protein J) CLIA Kit
E-CL-H0725	Human mIgM (Membrane Immumoglobulin M) CLIA Kit
E-CL-H0726	Human GPRC5A (G Protein Coupled Receptor, Family C, Group 5, Member A) CLIA Kit
E-CL-H0727	Human FABP5 (Fatty Acid Binding Protein 5, Epidermal) CLIA Kit
E-CL-H0728	Human VEGFR1/FLT1 (Vascular Endothelial Growth Factor Receptor 1) CLIA Kit
E-CL-H0729	Human ET (Epitestosterone) CLIA Kit
E-CL-H0730	Human PKM2 (Pyruvate Kinase, Muscle) CLIA Kit
E-CL-H0731	Human PDH α (Pyruvate Dehydrogenase Alpha) CLIA Kit
E-CL-H0732	Human PDP (Pyruvate Dehydrogenase Phosphatase) CLIA Kit
E-CL-H0733	Human VIM (Vimentin) CLIA Kit
E-CL-H0734	Human VTN (Vitronectin) CLIA Kit
E-CL-H0735	Human PPOX (Protoporphyrinogen Oxidase) CLIA Kit
E-CL-H0736	Human PRB1 (Basic Salivary Proline Rich Protein 1) CLIA Kit
E-CL-H0737	Human FOXC1 (Forkhead Box Protein C1) CLIA Kit
E-CL-H0738	Human FOXO1 (Forkhead Box Protein O1) CLIA Kit
E-CL-H0739	Human FOXO3 (Forkhead Box Protein O3) CLIA Kit
E-CL-H0740	Human FOXP3 (Forkhead Box Protein P3) CLIA Kit

Cat.NO	Product Name
E-CL-H0741	Human EFNA4 (Ephrin A4) CLIA Kit
E-CL-H0742	Human TFF3 (Trefoil Factor 3, Intestinal) CLIA Kit
E-CL-H0743	Human TAGLN (Transgelin) CLIA Kit
E-CL-H0744	Human TEP1 (Telomerase Associated Protein 1) CLIA Kit
E-CL-H0745	Human PRB2 (Basic Salivary Proline Rich Protein 2) CLIA Kit
E-CL-H0746	Human SOD1 (Superoxide Dismutase 1, Soluble) CLIA Kit
E-CL-H0747	Human FGF23 (Fibroblast Growth Factor 23) CLIA Kit
E-CL-H0748	Human FGF9 (Fibroblast Growth Factor 9) CLIA Kit
E-CL-H0749	Human PRF1 (Perforin 1/Pore-Forming Protein) CLIA Kit
E-CL-H0750	Human FSHR (Follicle Stimulating Hormone Receptor) CLIA Kit
E-CL-H0751	Human TGM1 (Transglutaminase 1, Keratinocyte) CLIA Kit
E-CL-H0752	Human E2 (Estradiol) CLIA Kit
E-CL-H0753	Human EPHA10 (Ephrin Type A Receptor 10) CLIA Kit
E-CL-H0754	Human ER β (Estrogen Receptor Beta) CLIA Kit
E-CL-H0755	Human ER α (Estrogen Receptor Alpha) CLIA Kit
E-CL-H0756	Human EFNA5 (Ephrin A5) CLIA Kit
E-CL-H0757	Human PRDX2 (Peroxisome oxidoreductin 2) CLIA Kit
E-CL-H0758	Human TSH (Thyroid Stimulating Hormone) CLIA Kit
E-CL-H0759	Human TRH (Thyrotropin Releasing Hormone) CLIA Kit

Cat.NO	Product Name
E-CL-H0760	Human FSH (Follicle-Stimulating Hormone) CLIA Kit
E-CL-H0761	Human CRH (Corticotropin Releasing Hormone) CLIA Kit
E-CL-H0762	Human GHRH (Growth Hormone Releasing Hormone) CLIA Kit
E-CL-H0763	Human δ SIP- α (Delta Sleep Inducing Peptide α) CLIA Kit
E-CL-H0764	Human GnRH (Gonadotropin Releasing Hormone) CLIA Kit
E-CL-H0765	Human SCT (Secretin) CLIA Kit
E-CL-H0766	Human FABP7 (Fatty Acid Binding Protein 7, Brain) CLIA Kit
E-CL-H0767	Human MCP-2/CCL8 (Monocyte Chemotactic Protein 2) CLIA Kit
E-CL-H0768	Human MCP-4/CCL13 (Monocyte Chemotactic Protein 4) CLIA Kit
E-CL-H0769	Human GPER1 (G Protein Coupled Estrogen Receptor 1) CLIA Kit
E-CL-H0770	Human PRDX1 (Peroxisomal oxidoreductase 1) CLIA Kit
E-CL-H0771	Human ELN (Elastin) CLIA Kit
E-CL-H0772	Human ELA4 (Elastase 4) CLIA Kit
E-CL-H0773	Human PROC (Protein C) CLIA Kit
E-CL-H0774	Human PDIA2 (Protein Disulfide Isomerase A2) CLIA Kit
E-CL-H0775	Human PROS (Protein S) CLIA Kit
E-CL-H0776	Human PROZ (Protein Z) CLIA Kit
E-CL-H0777	Human PDIA3 (Protein Disulfide Isomerase A3) CLIA Kit
E-CL-H0778	Human PRSS2 (Protease, Serine, 2) CLIA Kit

Cat.NO	Product Name
E-CL-H0779	Human PKC β 1 (Protein Kinase C Beta 1) CLIA Kit
E-CL-H0780	Human PKC ϵ (Protein Kinase C Epsilon) CLIA Kit
E-CL-H0781	Human PRSS23 (Protease, Serine 23) CLIA Kit
E-CL-H0782	Human PRSS1 (Protease, Serine, 1) CLIA Kit
E-CL-H0783	Human PAR2 (Protease Activated Receptor 2) CLIA Kit
E-CL-H0784	Human PCI (Protein C Inhibitor) CLIA Kit
E-CL-H0785	Human PDI (Protein Disulfide Isomerase) CLIA Kit
E-CL-H0786	Human PDIA5 (Protein Disulfide Isomerase A5) CLIA Kit
E-CL-H0787	Human NGF (Nerve Growth Factor) CLIA Kit
E-CL-H0788	Human PRSS33 (Protease, Serine 33) CLIA Kit
E-CL-H0789	Human LDL (Low Density Lipoprotein) CLIA Kit
E-CL-H0790	Human LDLR (Low Density Lipoprotein Receptor) CLIA Kit
E-CL-H0791	Human LRP6 (Low Density Lipoprotein Receptor Related Protein) CLIA Kit
E-CL-H0792	Human RETN (Resistin) CLIA Kit
E-CL-H0793	Human RPS6K β 1 (Ribosomal Protein S6 Kinase Beta 1) CLIA Kit
E-CL-H0794	Human APP (Amyloid Precursor Protein) CLIA Kit
E-CL-H0795	Human D1R/DRD1 (Dopamine Receptor D1) CLIA Kit
E-CL-H0796	Human D2R/DRD2 (Dopamine Receptor D2) CLIA Kit
E-CL-H0797	Human D β H (Dopamine- β -Hydroxylase) CLIA Kit

Cat.NO	Product Name
E-CL-H0798	Human DDC (Dopamine Decarboxylase) CLIA Kit
E-CL-H0799	Human DAT (Dopamine Transporter) CLIA Kit
E-CL-H0800	Human DT (Dopachrome Tautomerase) CLIA Kit
E-CL-H0801	Human VS (Versican/PG-M) CLIA Kit
E-CL-H0802	Human PIGR/SC (Polymeric Immunoglobulin Receptor/Membrane Secretory Component) CLIA Kit
E-CL-H0803	Human PYY (Peptide YY) CLIA Kit
E-CL-H0804	Human PTN (Pleiotrophin) CLIA Kit
E-CL-H0805	Human DAO (Diamine Oxidase) CLIA Kit
E-CL-H0806	Human DLD (Dihydrolipoyl Dehydrogenase) CLIA Kit
E-CL-H0807	Human NRP1 (Neuropilin 1) CLIA Kit
E-CL-H0808	Human DNM1 (Dynamin 1) CLIA Kit
E-CL-H0809	Human DNM2 (Dynamin 2) CLIA Kit
E-CL-H0810	Human fFDFT1 (Farnesyl Diphosphate Farnesyltransferase 1) CLIA Kit
E-CL-H0811	Human rT3 (Reverse Triiodothyronine) CLIA Kit
E-CL-H0812	Human Ub (Ubiquitin) CLIA Kit
E-CL-H0813	Human UBD (Ubiquitin D) CLIA Kit
E-CL-H0814	Human UBE1L (Ubiquitin Activating Enzyme E1 Like Protein) CLIA Kit
E-CL-H0815	Human UBE2C (Ubiquitin Conjugating Enzyme E2C) CLIA Kit
E-CL-H0816	Human UBE3A (Ubiquitin Protein Ligase E3A) CLIA Kit

Cat.NO	Product Name
E-CL-H0817	Human NNE (Non-Neuronal Enolase) CLIA Kit
E-CL-H0818	Human SDC3 (Syndecan 3) CLIA Kit
E-CL-H0819	Human TPS (Tryptase) CLIA Kit
E-CL-H0820	Human FBLN5 (Fibulin 5) CLIA Kit
E-CL-H0821	Human REG1 β (Regenerating Islet Derived Protein 1 Beta) CLIA Kit
E-CL-H0822	Human SPA (Pulmonary Surfactant Associated Protein A) CLIA Kit
E-CL-H0823	Human SPB (Pulmonary Surfactant Associated Protein B) CLIA Kit
E-CL-H0824	Human SPC (Pulmonary Surfactant Associated Protein C) CLIA Kit
E-CL-H0825	Human SPD (Pulmonary Surfactant Associated Protein D) CLIA Kit
E-CL-H0826	Human PARC/CCL18 (Pulmonary Activation Regulated Chemokine) CLIA Kit
E-CL-H0827	Human PRSS8 (Protease, Serine 8) CLIA Kit
E-CL-H0828	Human SPARCL1 (SPARC Like Protein 1) CLIA Kit
E-CL-H0829	Human SEMA3A (Semaphorin 3A) CLIA Kit
E-CL-H0830	Human sIgA (Secretory Immunoglobulin A) CLIA Kit
E-CL-H0831	Human SLPI (Secretory Leukocyte Protease Inhibitor) CLIA Kit
E-CL-H0832	Human HIF-1 α (Hypoxia Inducible Factor 1 Alpha) CLIA Kit
E-CL-H0833	Human SAP (Serum Amyloid P Component) CLIA Kit
E-CL-H0834	Human PAR4 (Protease Activated Receptor 4) CLIA Kit
E-CL-H0835	Human GPR120 (G Protein Coupled Receptor 120) CLIA Kit

Cat.NO	Product Name
E-CL-H0836	Human REG3α (Regenerating Islet Derived Protein 3 Alpha) CLIA Kit
E-CL-H0837	Human Cyr61 (Cysteine Rich Protein, Angiogenic Inducer 61) CLIA Kit
E-CL-H0838	Human LRG1 (Leucine Rich Alpha-2-Glycoprotein 1) CLIA Kit
E-CL-H0839	Human REG3γ (Regenerating Islet Derived Protein 3 Gamma) CLIA Kit
E-CL-H0840	Human S100A8 (S100 Calcium Binding Protein A8) CLIA Kit
E-CL-H0841	Human S100A9 (S100 Calcium Binding Protein A9) CLIA Kit
E-CL-H0842	Human GPR37 (G Protein Coupled Receptor 37) CLIA Kit
E-CL-H0843	Human S100A11 (S100 Calcium Binding Protein A11) CLIA Kit
E-CL-H0844	Human S100A12 (S100 Calcium Binding Protein A12) CLIA Kit
E-CL-H0845	Human S100A13 (S100 Calcium Binding Protein A13) CLIA Kit
E-CL-H0846	Human S100A6 (S100 Calcium Binding Protein A6) CLIA Kit
E-CL-H0847	Human S100A7 (S100 Calcium Binding Protein A7) CLIA Kit
E-CL-H0848	Human S100B (S100 Calcium Binding Protein B) CLIA Kit
E-CL-H0849	Human SDC1 (Syndecan 1) CLIA Kit
E-CL-H0850	Human SCF (Stem Cell Factor) CLIA Kit
E-CL-H0851	Human SCFR/c-Kit (Stem Cell Factor Receptor) CLIA Kit
E-CL-H0852	Human GAL (Galanin) CLIA Kit
E-CL-H0853	Human MBP/MBL (Mannose Binding Protein/Mannose Binding Lectin) CLIA Kit
E-CL-H0854	Human MR (Mannose Receptor) CLIA Kit

Cat.NO	Product Name
E-CL-H0855	Human PGAM1 (Phosphoglycerate Mutase 1) CLIA Kit
E-CL-H0856	Human PK (Pyruvate Kinase, Liver And RBC) CLIA Kit
E-CL-H0857	Human PYGL (Glycogen Phosphorylase, Liver) CLIA Kit
E-CL-H0858	Human SFRP2 (Secreted Frizzled Related Protein 2) CLIA Kit
E-CL-H0859	Human GP-73 (Golgi Protein 73) CLIA Kit
E-CL-H0860	Human SIAE (Sialic Acid Acetylesterase) CLIA Kit
E-CL-H0861	Human ST2 (Syntenin 2) CLIA Kit
E-CL-H0862	Human SIGLEC1 (Sialic Acid Binding Ig Like Lectin 1) CLIA Kit
E-CL-H0863	Human MHB (Methemoglobin) CLIA Kit
E-CL-H0864	Human SIGLEC10 (Sialic Acid Binding Ig Like Lectin 10) CLIA Kit
E-CL-H0865	Human SIGLEC2 (Sialic Acid Binding Ig Like Lectin 2) CLIA Kit
E-CL-H0866	Human PADI4 (Peptidyl Arginine Deiminase Type IV) CLIA Kit
E-CL-H0867	Human SIGLEC8 (Sialic Acid Binding Ig Like Lectin 8) CLIA Kit
E-CL-H0868	Human Smad3/MADH3 (Mothers Against Decapentaplegic Homolog 3) CLIA Kit
E-CL-H0869	Human GDH/GLDH (Glutamate dehydrogenase) CLIA Kit
E-CL-H0870	Human GAD2 (Glutamate Decarboxylase 2) CLIA Kit
E-CL-H0871	Human GAD-Ab (Anti-Glutamic Acid Decarboxylase antibody) CLIA Kit
E-CL-H0872	Human GSTα1 (Glutathione S Transferase Alpha 1) CLIA Kit
E-CL-H0873	Human GSTθ1/GSTt1 (Glutathione S Transferase Theta 1) CLIA Kit

Cat.NO	Product Name
E-CL-H0874	Human GSTθ2/GSTt2 (Glutathione S Transferase Theta 2) CLIA Kit
E-CL-H0875	Human GSTω1 (Glutathione S Transferase Omega 1) CLIA Kit
E-CL-H0876	Human HDAC2 (Histone Deacetylase 2) CLIA Kit
E-CL-H0877	Human C1QTNF1 (C1q and Tumor Necrosis Factor Related Protein 1) CLIA Kit
E-CL-H0878	Human GR (Glutathione Reductase) CLIA Kit
E-CL-H0879	Human GSTpi (Glutathione S Transferases Pi) CLIA Kit
E-CL-H0880	Human OPG (Osteoprotegerin) CLIA Kit
E-CL-H0881	Human OC/BGP (Osteocalcin) CLIA Kit
E-CL-H0882	Human ELA3B (Elastase 3B) CLIA Kit
E-CL-H0883	Human MUSK (Muscle Skeletal Receptor Tyrosine Kinase) CLIA Kit
E-CL-H0884	Human GSTα5 (Glutathione S Transferase Alpha 5) CLIA Kit
E-CL-H0885	Human OPN (Osteopontin) CLIA Kit
E-CL-H0886	Human MPIF-1/CCL23 (Myeloid Progenitor Inhibitory Factor 1) CLIA Kit
E-CL-H0887	Human SOX9 (Sex Determining Region Y Box Protein 9) CLIA Kit
E-CL-H0888	Human GREM1 (Gremlin 1) CLIA Kit
E-CL-H0889	Human ON (Osteonectin) CLIA Kit
E-CL-H0890	Human OP (Oligomeric Protein) CLIA Kit
E-CL-H0891	Human PAPPA2 (Pappalysin 2) CLIA Kit
E-CL-H0892	Human IgA (Immunoglobulin A) CLIA Kit

Cat.NO	Product Name
E-CL-H0893	Human ALDOC (Aldolase C, Fructose Bisphosphate) CLIA Kit
E-CL-H0894	Human PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) CLIA Kit
E-CL-H0895	Human PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) CLIA Kit
E-CL-H0896	Human PPAR- γ (Peroxisome Proliferator Activated Receptor Gamma) CLIA Kit
E-CL-H0897	Human LPO (Lactoperoxidase) ELSIA Kit
E-CL-H0898	Human RNASE3/ECP (Ribonuclease A3/Eosinophil Cationic Protein) CLIA Kit
E-CL-H0899	Human RI (Ribonuclease Inhibitor) CLIA Kit
E-CL-H0900	Human NF κ B (Nuclear Factor Kappa B) CLIA Kit
E-CL-H0901	Human RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) CLIA Kit
E-CL-H0902	Human NF- κ B p65 (Nuclear Factor Kappa B p65) CLIA Kit
E-CL-H0903	Human SERPINA1/A1AT (serpin peptidase inhibitor clade A (alpha-1 antiproteinase antitrypsin) member 1) CLIA Kit
E-CL-H0904	Human THRSP (Thyroid Hormone Responsive) CLIA Kit
E-CL-H0905	Human EPOR (Erythropoietin Receptor) CLIA Kit
E-CL-H0906	Human EXT1 (Exostoses 1) CLIA Kit
E-CL-H0907	Human GST κ 1 (Glutathione S Transferase Kappa 1) CLIA Kit
E-CL-H0908	Human LTA4H (Leukotriene A4 Hydrolase) CLIA Kit
E-CL-H0909	Human GST α 4 (Glutathione S Transferase Alpha 4) CLIA Kit
E-CL-H0910	Human SNC γ (Synuclein Gamma) CLIA Kit
E-CL-H0911	Human TLN (Talin) CLIA Kit

Cat.NO	Product Name
E-CL-H0912	Human Sec23b (protein transport protein Sec23B) CLIA Kit
E-CL-H0913	Human LHRH/GnRH (Luteinizing Hormone Releasing Hormone) CLIA Kit
E-CL-H0914	Human FABP6 (Fatty Acid Binding Protein 6, Ileal) CLIA Kit
E-CL-H0915	Human TMPRSS2 (Transmembrane Protease, Serine 2) CLIA Kit
E-CL-H0916	Human TNNI1 (Troponin I Type 1, Slow Skeletal) CLIA Kit
E-CL-H0917	Human FSCN (Fascin) CLIA Kit
E-CL-H0918	Human TNNC2 (Troponin C Type 2) CLIA Kit
E-CL-H0919	Human FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) CLIA Kit
E-CL-H0920	Human CKM (Creatine Kinase, Muscle) CLIA Kit
E-CL-H0921	Human CKMB (Creatine Kinase MB Isoenzyme) CLIA Kit
E-CL-H0922	Human DMD (Dystrophin) CLIA Kit
E-CL-H0923	Human MSTN (Myostatin) CLIA Kit
E-CL-H0924	Human PUMA (P53 Upregulated Modulator of Apoptosis) CLIA Kit
E-CL-H0925	Human MGP (Matrix Gla Protein) CLIA Kit
E-CL-H0926	Human MMP-1 (Matrix Metalloproteinase 1) CLIA Kit
E-CL-H0927	Human IgA2 (Immunoglobulin A2) CLIA Kit
E-CL-H0928	Human MMP-11 (Matrix Metalloproteinase 11) CLIA Kit
E-CL-H0929	Human MGST1 (Microsomal Glutathione S Transferase 1) CLIA Kit
E-CL-H0930	Human MMP-2 (Matrix Metalloproteinase 2) CLIA Kit

Cat.NO	Product Name
E-CL-H0931	Human MMP-3 (Matrix Metalloproteinase 3) CLIA Kit
E-CL-H0932	Human MMP-17 (Matrix Metalloproteinase 17) CLIA Kit
E-CL-H0933	Human MMP-24 (Matrix Metalloproteinase 24) CLIA Kit
E-CL-H0934	Human MMP-7 (Matrix Metalloproteinase 7) CLIA Kit
E-CL-H0935	Human MMP-8 (Matrix Metalloproteinase 8) CLIA Kit
E-CL-H0936	Human MMP-9 (Matrix Metalloproteinase 9) CLIA Kit
E-CL-H0937	Human ELA2B (Elastase 2B) CLIA Kit
E-CL-H0938	Human TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) CLIA Kit
E-CL-H0939	Human TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) CLIA Kit
E-CL-H0940	Human TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) CLIA Kit
E-CL-H0941	Human LGALS3BP (Lectin Galactoside Binding, Soluble 3 Binding Protein) CLIA Kit
E-CL-H0942	Human ST2/IL-1 R4 (Stromelysin-2) CLIA Kit
E-CL-H0943	Human TNNI2 (Troponin I Type 2, Fast Skeletal) CLIA Kit
E-CL-H0944	Human TP/TXA2R (Thromboxane Receptor) CLIA Kit
E-CL-H0945	Human TDGF1 (Teratocarcinoma Derived Growth Factor 1) CLIA Kit
E-CL-H0946	Human AICAR (AICA Riboside) CLIA Kit
E-CL-H0947	Human KLK8 (Kallikrein 8) CLIA Kit
E-CL-H0948	Human GAL3 (Galectin 3) CLIA Kit
E-CL-H0949	Human MTX (Methotrexate) CLIA Kit

Cat.NO	Product Name
E-CL-H0950	Human fMet (Formylmethionine) CLIA Kit
E-CL-H0951	Human FPR2 (Formyl Peptide Receptor 2) CLIA Kit
E-CL-H0952	Human PTH (Parathyroid Hormone) CLIA Kit
E-CL-H0953	Human PTHrP (Parathyroid Hormone Related Protein) CLIA Kit
E-CL-H0954	Human TPP1 (Tripeptidyl Peptidase I) CLIA Kit
E-CL-H0955	Human CD276 (Cluster of Differentiation 276) CLIA Kit
E-CL-H0956	Human TBG (Thyroxine Binding Globulin) CLIA Kit
E-CL-H0957	Human PAH (Phenylalanine Hydroxylase) CLIA Kit
E-CL-H0958	Human TG (Thyroglobulin) CLIA Kit
E-CL-H0959	Human T4 (Thyroxine) CLIA Kit
E-CL-H0960	Human GJβ1 (Gap Junction Protein Beta 1) CLIA Kit
E-CL-H0961	Human PCT (Procalcitonin) CLIA Kit
E-CL-H0962	Human PCPE1 (Procollagen C Proteinase Enhancer 1) CLIA Kit
E-CL-H0963	Human PYD (Pyridinoline Crosslinks) CLIA Kit
E-CL-H0964	Human GDNF (Glial Cell Line Derived Neurotrophic Factor) CLIA Kit
E-CL-H0965	Human CDSN (Corneodesmosin) CLIA Kit
E-CL-H0966	Human CNTN1 (Contactin 1) CLIA Kit
E-CL-H0967	Human Des (Desmin) CLIA Kit
E-CL-H0968	Human HP (Haptoglobin) CLIA Kit

Cat.NO	Product Name
E-CL-H0969	Human HPR (Haptoglobin Related Protein) CLIA Kit
E-CL-H0970	Human TNNT1 (Troponin T Type 1, Slow Skeletal) CLIA Kit
E-CL-H0971	Human TNPO2 (Transportin 2) CLIA Kit
E-CL-H0972	Human UCP-2 (Uncoupling Protein 2, Mitochondrial) CLIA Kit
E-CL-H0973	Human MT1 (Metallothionein 1) CLIA Kit
E-CL-H0974	Human MT2 (Metallothionein 2) CLIA Kit
E-CL-H0975	Human TJP1 (Tight Junction Protein 1) CLIA Kit
E-CL-H0976	Human CLC (Charcot Leyden Crystal Protein) CLIA Kit
E-CL-H0977	Human ACVB (Activin B) CLIA Kit
E-CL-H0978	Human CRY α B (Crystallin Alpha B) CLIA Kit
E-CL-H0979	Human MIP-3 β /ELC/CCL19 (Macrophage Inflammatory Protein 3 β) CLIA Kit
E-CL-H0980	Human MIF (Macrophage Migration Inhibitory Factor) CLIA Kit
E-CL-H0981	Human PDCD1 (Programmed Cell Death Protein 1) CLIA Kit
E-CL-H0982	Human AP (Allopregnanolone) CLIA Kit
E-CL-H0983	Human CRHBP (Corticotropin Releasing Hormone Binding Protein) CLIA Kit
E-CL-H0984	Human α HSP (Alpha-Hemoglobin Stabilizing Protein) CLIA Kit
E-CL-H0985	Human MMP-28 (Matrix Metalloproteinase 28) CLIA Kit
E-CL-H0986	Human TLR4 (Toll-Like Receptor 4) CLIA Kit
E-CL-H0987	Human SPTAN1 (Alpha-Fodrin) CLIA Kit

Cat.NO	Product Name
E-CL-H0988	Human α ZGP1 (Alpha-2-Glycoprotein 1, Zine Binding) CLIA Kit
E-CL-H0989	Human MSLN (Mesothelin) CLIA Kit
E-CL-H0990	Human SOST (Sclerostin) CLIA Kit
E-CL-H0991	Human IRF5 (Interferon Regulatory Factor 5) CLIA Kit
E-CL-H0992	Human SIRT1 (Sirtuin 1) CLIA Kit
E-CL-H0993	Human PDCD1LG1 (Programmed Cell Death Protein 1 Ligand 1) CLIA Kit
E-CL-H0994	Human τ PK1 (Tau-Protein Kinase) CLIA Kit
E-CL-H0995	Human ACP5 (Tartrate Resistant Acid Phosphatase 5) CLIA Kit
E-CL-H0996	Human VWFCP/ADAMTS13 (Von Willebrand Factor Cleaving Protease) CLIA Kit
E-CL-H0997	Human HMGB-1 (High mobility group protein B1) CLIA Kit
E-CL-H0998	Human PGF (Placental Growth Factor) CLIA Kit
E-CL-H0999	Human FOXC2 (Forkhead Box Protein C2) CLIA Kit
E-CL-H1000	Human ESM1 (Endothelial Cell Specific Molecule 1) CLIA Kit
E-CL-H1001	Human Rac1 (Ras Related C3 Botulinum Toxin Substrate 1) CLIA Kit
E-CL-H1002	Human LILRB2 (Leukocyte Immunoglobulin Like Receptor Subfamily B Member 2) CLIA Kit
E-CL-H1003	Human OTR (Oxytocin Receptor) CLIA Kit
E-CL-H1004	Human PDIA6 (Protein Disulfide Isomerase A6) CLIA Kit
E-CL-H1005	Human NFE2L2 (Nuclear Factor, Erythroid Derived 2 Like 2) CLIA Kit
E-CL-H1006	Human DNP-Ab (Deoxyribonucleic Protein Antibody) CLIA Kit

Cat.NO	Product Name
E-CL-H1007	Human Pro-MMP-10 (Pro-Matrix Metalloproteinase 10) CLIA Kit
E-CL-H1008	Human Pro-MMP-13 (Pro-Matrix Metalloproteinase 13) CLIA Kit
E-CL-H1009	Human PTX 3/TSG-14 (Pentraxin 3, Long) CLIA Kit
E-CL-H1010	Human PDGF-AA (Platelet Derived Growth Factor-AA) CLIA Kit
E-CL-H1011	Human PDGF-AB (Platelet Derived Growth Factor AB) CLIA Kit
E-CL-H1012	Human PDGF-BB (Platelet Derived Growth Factor-BB) CLIA Kit
E-CL-H1013	Human PGRN (Progranulin) CLIA Kit
E-CL-H1014	Human PCSK9 (Proprotein Convertase Subtilisin/Kexin Type 9) CLIA Kit
E-CL-H1015	Human RBP4 (Retinol Binding Protein 4, Plasma) CLIA Kit
E-CL-H1016	Human RLN2 (Relaxin-2) CLIA Kit
E-CL-H1017	Human RBP1 (Retinol Binding Protein 1, Cellular) CLIA Kit
E-CL-H1018	Human Surv (Survivin) CLIA Kit
E-CL-H1019	Human RBP3 (Retinol Binding Protein 3, Interstitial) CLIA Kit
E-CL-H1020	Human TGF- α (Transforming Growth Factor α) CLIA Kit
E-CL-H1021	Human TGF- β 2 (Transforming Growth Factor β 2) CLIA Kit
E-CL-H1022	Human TPO (Thrombopoietin) CLIA Kit
E-CL-H1023	Human TSP-1 (Thrombospondin-1) CLIA Kit
E-CL-H1024	Human TSP-2 (Thrombospondin-2) CLIA Kit
E-CL-H1025	Human TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) CLIA Kit

Cat.NO	Product Name
E-CL-H1026	Human TREM-1 (Triggering Receptor Expressed on Myeloid Cells-1) CLIA Kit
E-CL-H1027	Human TSLP (Thymic Stromal Lymphopoietin) CLIA Kit
E-CL-H1028	Human PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) CLIA Kit
E-CL-H1029	Human VEGF-C (Vascular Endothelial Growth Factor C) CLIA Kit
E-CL-H1030	Human VEGF-D (Vascular Endothelial Growth Factor D) CLIA Kit
E-CL-H1031	Human VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) CLIA Kit
E-CL-H1032	Human DBP (Vitamin D Binding Protein) CLIA Kit
E-CL-H1033	Human INS (Insulin) CLIA Kit
E-CL-H1034	Human EMMPRIN/CD147 (Extracellular Matrix Metalloproteinase Inducer) CLIA Kit
E-CL-H1035	Human CYB5A (Cytochrome b5) CLIA Kit
E-CL-H1036	Human CG β (Chorionic Gonadotropin Beta Polypeptide) CLIA Kit
E-CL-H1037	Human GzmB (Granzyme B) CLIA Kit
E-CL-H1038	Human GNLY (Granulysin) CLIA Kit
E-CL-H1039	Human GRN (Granulin) CLIA Kit
E-CL-H1040	Human CR2/CD21 (Complement Receptor 2) CLIA Kit
E-CL-H1041	Human CD30L (Cluster of Differentiation 30 Lligand) CLIA Kit
E-CL-H1042	Human cADPRH/CD38 (Cyclic ADP Ribose Hydrolase) CLIA Kit
E-CL-H1043	Human LAB7-2/CD86 (B-Lymphocyte Activation Antigen B7-2) CLIA Kit
E-CL-H1044	Human EGFR2 (Epidermal Growth Factor Receptor 2) CLIA Kit

Cat.NO	Product Name
E-CL-H1045	Human LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) CLIA Kit
E-CL-H1046	Human PECAM1/CD31 (Platelet/Endothelial Cell Adhesion Molecule 1) CLIA Kit
E-CL-H1047	Human sTfR1 (Soluble Transferrin Receptor1) CLIA Kit
E-CL-H1048	Human Tyk-2 (Tyrosine Kinase 2) CLIA Kit
E-CL-H1049	Human TH (Tyrosine Hydroxylase) CLIA Kit
E-CL-H1050	Human mTOR (Mammalian Target of Rapamycin) CLIA Kit
E-CL-H1051	Human PSPN (Persephin) CLIA Kit
E-CL-H1052	Human TGFβI (Transforming Growth Factor Beta Induced Protein) CLIA Kit
E-CL-H1053	Human HLA-A (Leukocyte Antigen A) CLIA Kit
E-CL-H1054	Human PKCη (Protein Kinase C Eta) CLIA Kit
E-CL-H1055	Human MHCG (Major Histocompatibility Complex Class I G) CLIA Kit
E-CL-H1056	Human TYR (Tyrosinase) CLIA Kit
E-CL-H1057	Human FGL1 (Fibrinogen Like Protein 1) CLIA Kit
E-CL-H1058	Human ORM2 (Orosomucoid 2) CLIA Kit
E-CL-H1059	Human PRG4 (Proteoglycan 4) CLIA Kit
E-CL-H1060	Human CX40 (Connexin 40) CLIA Kit
E-CL-H1061	Human CX43 (Connexin 43) CLIA Kit
E-CL-H1062	Human SDC4 (Syndecan 4) CLIA Kit
E-CL-H1063	Human LAP (Leucine Aminopeptidase) CLIA Kit

Cat.NO	Product Name
E-CL-H1064	Human LT β (Lymphotoxin Beta) CLIA Kit
E-CL-H1065	Human LFA-3 (Lymphocyte Function Associated Antigen 3) CLIA Kit
E-CL-H1066	Human Lptn/XCL1 (Lymphotactin) CLIA Kit
E-CL-H1067	Human PFKM (Phosphofructokinase, Muscle) CLIA Kit
E-CL-H1068	Human A-CoA (Acetyl Coenzyme A) CLIA Kit
E-CL-H1069	Human PCK1 (Phosphoenolpyruvate Carboxykinase 1) CLIA Kit
E-CL-H1070	Human PLAA (Phospholipase A2 Activating Protein) CLIA Kit
E-CL-H1071	Human PLSCR1 (Phospholipid Scramblase 1) CLIA Kit
E-CL-H1072	Human PLSCR2 (Phospholipid Scramblase 2) CLIA Kit
E-CL-H1073	Human PI3K2 α (Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) CLIA Kit
E-CL-H1074	Human GPC1 (Glypican 1) CLIA Kit
E-CL-H1075	Human GPC2 (Glypican 2) CLIA Kit
E-CL-H1076	Human GPC3 (Glypican 3) CLIA Kit
E-CL-H1077	Human GPC4 (Glypican 4) CLIA Kit
E-CL-H1078	Human PS (Phosphatidylserine) CLIA Kit
E-CL-H1079	Human PLTP (Phospholipid Transfer Protein) CLIA Kit
E-CL-H1080	Human PLC γ 1 (Phospholipase C Gamma 1) CLIA Kit
E-CL-H1081	Human SCCA1 (Squamous Cell Carcinoma Antigen 1) CLIA Kit
E-CL-H1082	Human SCCA2 (Squamous Cell Carcinoma Antigen 2) CLIA Kit

Cat.NO	Product Name
E-CL-H1083	Human DS (Dermatan Sulfate) CLIA Kit
E-CL-H1084	Human Trx (Thioredoxin) CLIA Kit
E-CL-H1085	Human TrxR1 (Thioredoxin Reductase 1) CLIA Kit
E-CL-H1086	Human TBP2 (Thioredoxin Binding Protein 2) CLIA Kit
E-CL-H1087	Human OVA sIgE (Ovalbumin Specific IgE) CLIA Kit
E-CL-H1088	Human FcαR/CD89 (Receptor I for the Fc Fragment of IgA) CLIA Kit
E-CL-H1089	Human FcγR II/CD32 (Receptor II for the Fc Fragment of IgG) CLIA Kit
E-CL-H1090	Human FcγR3B (Fc Fragment of IgG Low Affinity IIIb Receptor) CLIA Kit
E-CL-H1091	Human Fcγ (Fc Fragment of IgG) CLIA Kit
E-CL-H1092	Human IgG1 (Immunoglobulin G1) CLIA Kit
E-CL-H1093	Human IgG2 (Immunoglobulin G2) CLIA Kit
E-CL-H1094	Human IgG3 (Immunoglobulin G3) CLIA Kit
E-CL-H1095	Human IgG4 (Immunoglobulin G4) CLIA Kit
E-CL-H1096	Human MCP/CD46 (Membrane Cofactor Protein) CLIA Kit
E-CL-H1097	Human TSP-4 (Thrombospondin-4) CLIA Kit
E-CL-H1098	Human TCC C5b-9 (Terminal Complement Complex C5b-9) CLIA Kit
E-CL-H1099	Human SMN1 (Survival of Motor Neuron 1, Telomeric) CLIA Kit
E-CL-H1100	Human ETRA (Endothelin Receptor A) CLIA Kit
E-CL-H1101	Human ECE1 (Endothelin Converting Enzyme 1) CLIA Kit

Cat.NO	Product Name
E-CL-H1102	Human REPIN1 (Replication Initiator 1) CLIA Kit
E-CL-H1103	Human EL/LIPG (Endothelial Lipase) CLIA Kit
E-CL-H1104	Human VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) CLIA Kit
E-CL-H1105	Human VF (Visfatin) CLIA Kit
E-CL-H1106	Human CFTR (Cystic Fibrosis Transmembrane Conductance Regulator) CLIA Kit
E-CL-H1107	Human NGB (Neuroglobin) CLIA Kit
E-CL-H1108	Human PGD2S (Prostaglandin D2 Synthase, Brain) CLIA Kit
E-CL-H1109	Human OCT (Ornithine Carbamoyl Transferase) CLIA Kit
E-CL-H1110	Human ODC (Ornithine Decarboxylase) CLIA Kit
E-CL-H1111	Human PLAU/uPA(Urokinase-Type Plasminogen Activator) CLIA Kit
E-CL-H1112	Human ETRB (Endothelin Receptor B) CLIA Kit
E-CL-H1113	Human GS (Gelsolin) CLIA Kit
E-CL-H1114	Human TAT (Thrombin/Antithrombin Complex) CLIA Kit
E-CL-H1115	Human F1+2 (Prothrombin Fragment 1+2) CLIA Kit
E-CL-H1116	Human KERA (Keratocan) CLIA Kit
E-CL-H1117	Human VCL (Vinculin) CLIA Kit
E-CL-H1118	Human DEFα5 (Defensin Alpha 5, Paneth Cell Specific) CLIA Kit
E-CL-H1119	Human DPT (Dermatopontin) CLIA Kit
E-CL-H1120	Human SYK (Spleen Tyrosine Kinase) CLIA Kit

Cat.NO	Product Name
E-CL-H1121	Human SMM (Smooth Muscle Myosin) CLIA Kit
E-CL-H1122	Human PDPN (Podoplanin) CLIA Kit
E-CL-H1123	Human OSCAR (Osteoclast Associated Receptor) CLIA Kit
E-CL-H1124	Human IgM (Immunoglobulin M) CLIA Kit
E-CL-H1125	Human G6PD (Glucose 6 Phosphate Dehydrogenase) CLIA Kit
E-CL-H1126	Human GPI (Glucose 6 Phosphate Isomerase) CLIA Kit
E-CL-H1127	Human GCK (Glucokinase) CLIA Kit
E-CL-H1128	Human GOD (Glucose Oxidase) CLIA Kit
E-CL-H1129	Human GLUT1 (Glucose Transporter 1) CLIA Kit
E-CL-H1130	Human GLUT3 (Glucose Transporter 3) CLIA Kit
E-CL-H1131	Human GLUT4 (Glucose Transporter 4) CLIA Kit
E-CL-H1132	Human PA (Prealbumin) CLIA Kit
E-CL-H1133	Human PK2 (Prokineticin 2) CLIA Kit
E-CL-H1134	Human PKR1 (Prokineticin Receptor 1) CLIA Kit
E-CL-H1135	Human PKR2 (Prokineticin Receptor 2) CLIA Kit
E-CL-H1136	Human PGDS (Prostaglandin D Synthase) CLIA Kit
E-CL-H1137	Human PGE1 (Prostaglandin E1) CLIA Kit
E-CL-H1138	Human EP2 (Prostaglandin E Receptor 2) CLIA Kit
E-CL-H1139	Human PGF2 α (Prostaglandin F2 α) CLIA Kit

Cat.NO	Product Name
E-CL-H1140	Human PTGIS (Prostaglandin I Synthase) CLIA Kit
E-CL-H1141	Human PTGS1 (Prostaglandin Endoperoxide Synthase 1) CLIA Kit
E-CL-H1142	Human PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) CLIA Kit
E-CL-H1143	Human ACP3/PAP (Prostatic Acid Phosphatase 3) CLIA Kit
E-CL-H1144	Human NT-ProANP (N-Terminal Pro Atrial Natriuretic Peptide) CLIA Kit
E-CL-H1145	Human DSP (Desmoplakin) CLIA Kit
E-CL-H1146	Human SM/SPH (Sphingomyelin) CLIA Kit
E-CL-H1147	Human HJV/HFE2 (Hemojuvelin/Hemochromatosis Type 2) CLIA Kit
E-CL-H1148	Human NA/NE (Noradrenaline/Norepinephrine) CLIA Kit
E-CL-H1149	Human HSP β 6/HSP-20 (Heat Shock Protein Beta 6) CLIA Kit
E-CL-H1150	Human HSP-27/HSPB1 (Heat Shock Protein 27) CLIA Kit
E-CL-H1151	Human HSP-40 (Heat Shock Protein 40) CLIA Kit
E-CL-H1152	Human HSP-60/HSPD1 (Heat Shock Protein 60) CLIA Kit
E-CL-H1153	Human HSP-70 (Heat Shock Protein 70) CLIA Kit
E-CL-H1154	Human HSP-90 (Heat Shock Protein 90) CLIA Kit
E-CL-H1155	Human HSF2 (Heat Shock Transcription Factor 2) CLIA Kit
E-CL-H1156	Human HSF1 (Heat Shock Transcription Factor 1) CLIA Kit
E-CL-H1157	Human LZM (Lysozyme) CLIA Kit
E-CL-H1158	Human T3 (Triiodothyronine) CLIA Kit

Cat.NO	Product Name
E-CL-H1159	Human TFF1 (Trefoil Factor 1) CLIA Kit
E-CL-H1160	Human TFF2 (Trefoil Factor 2) CLIA Kit
E-CL-H1161	Human PEDF (Pigment Epithelium Derived Factor) CLIA Kit
E-CL-H1162	Human SDH (Sorbitol Dehydrogenase) CLIA Kit
E-CL-H1163	Human Ep-CAM/CD326 (Epithelial Cell Adhesion Molecule) CLIA Kit
E-CL-H1164	Human BGN (Biglycan) CLIA Kit
E-CL-H1165	Human NRG-1 (Neuregulin 1) CLIA Kit
E-CL-H1166	Human NKB (Neurokinin B) CLIA Kit
E-CL-H1167	Human GFAP (Glial Fibrillary Acidic Protein) CLIA Kit
E-CL-H1168	Human NPY (Neuropeptide Y) CLIA Kit
E-CL-H1169	Human NCAM/CD56 (Neural Cell Adhesion Molecule) CLIA Kit
E-CL-H1170	Human L1CAM/CD171 (L1-Cell Adhesion Molecule) CLIA Kit
E-CL-H1171	Human NT-3 (Neurotrophin-3) CLIA Kit
E-CL-H1172	Human NT-4 (Neurotrophin 4) CLIA Kit
E-CL-H1173	Human NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) CLIA Kit
E-CL-H1174	Human NAIP (Neuronal Apoptosis Inhibitory Protein) CLIA Kit
E-CL-H1175	Human NPHN (Nephrin) CLIA Kit
E-CL-H1176	Human EPI (Epinephrine/Adrenaline) CLIA Kit
E-CL-H1177	Human GROβ (Growth Regulated Oncogene Beta) CLIA Kit

Cat.NO	Product Name
E-CL-H1178	Human GROγ (Growth Regulated Oncogene Gamma) CLIA Kit
E-CL-H1179	Human GDF1 (Growth Differentiation Factor 1) CLIA Kit
E-CL-H1180	Human GDF11 (Growth Differentiation Factor 11) CLIA Kit
E-CL-H1181	Human GDF2 (Growth Differentiation Factor 2) CLIA Kit
E-CL-H1182	Human GDF3 (Growth Differentiation Factor 3) CLIA Kit
E-CL-H1183	Human GDF5 (Growth Differentiation Factor 5) CLIA Kit
E-CL-H1184	Human GDF9 (Growth Differentiation Factor 9) CLIA Kit
E-CL-H1185	Human GH2 (Growth Hormone 2) CLIA Kit
E-CL-H1186	Human GHR/GHBP (Growth Hormone Receptor) CLIA Kit
E-CL-H1187	Human GHITM (Growth Hormone Inducible Transmembrane Protein) CLIA Kit
E-CL-H1188	Human GAP43 (Growth Associated Protein 43) CLIA Kit
E-CL-H1189	Human SST (Somatostatin) CLIA Kit
E-CL-H1190	Human Grb10 (Growth Factor Receptor Bound Protein 10) CLIA Kit
E-CL-H1191	Human Grb2 (Growth Factor Receptor Bound Protein 2) CLIA Kit
E-CL-H1192	Human OXB (Orexin B) CLIA Kit
E-CL-H1193	Human CYPA (Cyclophilin A) CLIA Kit
E-CL-H1194	Human CYPB (Cyclophilin B) CLIA Kit
E-CL-H1195	Human EAR2 (Eosinophil Associated Ribonuclease A family Member 2) CLIA Kit
E-CL-H1196	Human NE/ELA2 (Elastase 2, Neutrophil) CLIA Kit

Cat.NO	Product Name
E-CL-H1197	Human TNK2 (Tyrosine Kinase, Non Receptor 2) CLIA Kit
E-CL-H1198	Human DAF (Decay Accelerating Factor) CLIA Kit
E-CL-H1199	Human STK39 (Serine/Threonine Kinase 39) CLIA Kit
E-CL-H1200	Human TACR1 (Tachykinin Receptor 1) CLIA Kit
E-CL-H1201	Human TACR2 (Tachykinin Receptor 2) CLIA Kit
E-CL-H1202	Human Anti-MPO (Anti-Myeloperoxidase Antibody) CLIA Kit
E-CL-H1203	Human MYD88 (Myeloid Differentiation Factor 88) CLIA Kit
E-CL-H1204	Human PR3 (Proteinase 3) CLIA Kit
E-CL-H1205	Human PP13 (Placental Protein13) CLIA Kit
E-CL-H1206	Human PLAP/ALPP (Placental Alkaline Phosphatase) CLIA Kit
E-CL-H1207	Human FETUB (Fetuin B) CLIA Kit
E-CL-H1208	Human PGRP-S (Peptidoglycan Recognition Protein-S) CLIA Kit
E-CL-H1209	Human GP5 (Glycoprotein V, Platelet) CLIA Kit
E-CL-H1210	Human GP6 (Glycoprotein VI, Platelet) CLIA Kit
E-CL-H1211	Human GSP (Glycosylated Serum Protein) CLIA Kit
E-CL-H1212	Human HbA1c (Glycosylated Hemoglobin/Hemoglobin A1c) CLIA Kit
E-CL-H1213	Human GPI (Glycophosphatidylinositol) CLIA Kit
E-CL-H1214	Human GR α (Glucocorticoid Receptor Alpha) CLIA Kit
E-CL-H1215	Human GR β (Glucocorticoid Receptor Beta) CLIA Kit

Cat.NO	Product Name
E-CL-H1216	Human GSK3 β (Glycogen Synthase Kinase 3 Beta) CLIA Kit
E-CL-H1217	Human GSK3 α (Glycogen Synthase Kinase 3 Alpha) CLIA Kit
E-CL-H1218	Human FTL (Ferritin, Light Polypeptide) CLIA Kit
E-CL-H1219	Human FTH (Ferritin, Heavy Polypeptide) CLIA Kit
E-CL-H1220	Human MRP1 (Multidrug Resistance Associated Protein 1) CLIA Kit
E-CL-H1221	Human SYNPO (Synaptopodin) CLIA Kit
E-CL-H1222	Human MT (Melatonin) CLIA Kit
E-CL-H1223	Human DHEA (Dehydroepiandrosterone) CLIA Kit
E-CL-H1224	Human DPD (Deoxypyridinoline) CLIA Kit
E-CL-H1225	Human DNase-I (Deoxyribonuclease I) CLIA Kit
E-CL-H1226	Human SA (Sialic Acid) CLIA Kit
E-CL-H1227	Human SIGLEC3 (Sialic Acid Binding Ig Like Lectin 3) CLIA Kit
E-CL-H1228	Human ITLN1 (Intelectin 1/Omentin) CLIA Kit
E-CL-H1229	Human STMN1 (Stathmin 1) CLIA Kit
E-CL-H1230	Human DCX (Doublecortin) CLIA Kit
E-CL-H1231	Human MAP-2 (Microtubule Associated Protein 2) CLIA Kit
E-CL-H1232	Human VDR (Vitamin D Receptor) CLIA Kit
E-CL-H1233	Human GKN1 (Gastrokine 1) CLIA Kit
E-CL-H1234	Human GKN2 (Gastrokine 2) CLIA Kit

Cat.NO	Product Name
E-CL-H1235	Human MTL (Motilin) CLIA Kit
E-CL-H1236	Human VIPR1 (Vasoactive Intestinal Peptide Receptor) CLIA Kit
E-CL-H1237	Human GT-34 (Gastrin 34) CLIA Kit
E-CL-H1238	Human ProGRP (Pro-Gastrin Releasing Peptide) CLIA Kit
E-CL-H1239	Human GIP (Gastric Inhibitory Polypeptide) CLIA Kit
E-CL-H1240	Human PTD ((Pentosidine) CLIA Kit
E-CL-H1241	Human CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) CLIA Kit
E-CL-H1242	Human CK-13/KRT13 (Cytokeratin 13) CLIA Kit
E-CL-H1243	Human CK-17/KRT17 (Cytokeratin 17) CLIA Kit
E-CL-H1244	Human CK-18/KRT18 (Cytokeratin 18) CLIA Kit
E-CL-H1245	Human CK-19/KRT19 (Cytokeratin 19) CLIA Kit
E-CL-H1246	Human CK-20/KRT20 (Cytokeratin 20) CLIA Kit
E-CL-H1247	Human CYFRA21-1 (Cytokeratin Fragment Antigen 21-1) CLIA Kit
E-CL-H1248	Human CYP7A1 (Cytochrome P450 7A1) CLIA Kit
E-CL-H1249	Human ECM1 (Extracellular Matrix Protein 1) CLIA Kit
E-CL-H1250	Human MEPE (Matrix Extracellular Phosphoglycoprotein) CLIA Kit
E-CL-H1251	Human ERK1 (Extracellular Signal Regulated Kinase 1) CLIA Kit
E-CL-H1252	Human ERK2 (Extracellular Signal Regulated Kinase 2) CLIA Kit
E-CL-H1253	Human SOCS3 (Suppressors of Cytokine Signaling 3) CLIA Kit

Cat.NO	Product Name
E-CL-H1254	Human CCND1 (Cyclin-D1) CLIA Kit
E-CL-H1255	Human CCND2 (Cyclin-D2) CLIA Kit
E-CL-H1256	Human FLN β (Filamin B, Beta) CLIA Kit
E-CL-H1257	Human FLNC (Filamin C, Gamma) CLIA Kit
E-CL-H1258	Human FMOD (Fibromodulin) CLIA Kit
E-CL-H1259	Human PAP (Plasmin-Antiplasmin Complex) CLIA Kit
E-CL-H1260	Human Plg (Plasminogen) CLIA Kit
E-CL-H1261	Human PAI1 (Plasminogen Activator Inhibitor 1) CLIA Kit
E-CL-H1262	Human PAI2 (Plasminogen Activator Inhibitor 2)CLIA Kit
E-CL-H1263	Human tPA (Plasminogen Activator, Tissue) CLIA Kit
E-CL-H1264	Human TAFI (Thrombin Activatable Fibrinolysis Inhibitor) CLIA Kit
E-CL-H1265	Human FG α (Fibrinogen Alpha) CLIA Kit
E-CL-H1266	Human FG β (Fibrinogen Beta) CLIA Kit
E-CL-H1267	Human FGL2 (Fibrinogen Like Protein 2) CLIA Kit
E-CL-H1268	Human KRT10 (Keratin, Type I Cytoskeletal 10) CLIA Kit
E-CL-H1269	Human SBEM (Small Breast Epithelial Mucin) CLIA Kit
E-CL-H1270	Human STAT3 (Signal Transducer and Activator of Transcription 3) CLIA Kit
E-CL-H1271	Human STAT4 (Signal Transducer and Activator of Transcription 4) CLIA Kit
E-CL-H1272	Human SOX18 (Sex Determining Region Y Box Protein 18) CLIA Kit

Cat.NO	Product Name
E-CL-H1273	Human SHBG (Sex Hormone Binding Globulin) CLIA Kit
E-CL-H1274	Human TS/TYMS (Thymidylate Synthetase) CLIA Kit
E-CL-H1275	Human TP (Thymidine Phosphorylase) CLIA Kit
E-CL-H1276	Human TMSβ10 (Thymosin Beta 10) CLIA Kit
E-CL-H1277	Human TP5 (Thymopentin) CLIA Kit
E-CL-H1278	Human OLFM4 (Olfactomedin 4) CLIA Kit
E-CL-H1279	Human VIP (Vasoactive Intestinal Peptide) CLIA Kit
E-CL-H1280	Human VE-Cadherin (Vascular Endothelial Cadherin) CLIA Kit
E-CL-H1281	Human VEGF165 (Vascular Endothelial Growth Factor165) CLIA Kit
E-CL-H1282	Human GSTμ2/GSTm2 (Glutathione S Transferase Mu 2, Muscle) CLIA Kit
E-CL-H1283	Human SUMO2 (Small Ubiquitin Related Modifier 2) CLIA Kit
E-CL-H1284	Human IgE (Immunoglobulin E) CLIA Kit
E-CL-H1285	Human VEGF-B (Vascular Endothelial Cell Growth Factor B) CLIA Kit
E-CL-H1286	Human VEGFR3/Flt4 (Vascular Endothelial Cell Growth Factor Receptor 3) CLIA Kit
E-CL-H1287	Human sVCAM-1/CD106 (soluble Vasculuar Cell Adhesion Molecule 1) CLIA Kit
E-CL-H1288	Human VWF (Von Willebrand Factor) CLIA Kit
E-CL-H1289	Human HbC (Hemoglobin C) CLIA Kit
E-CL-H1290	Human HBβ (Hemoglobin Beta) CLIA Kit
E-CL-H1291	Human HO1 (Heme Oxygenase 1) CLIA Kit

Cat.NO	Product Name
E-CL-H1292	Human HO2 (Heme Oxygenase 2) CLIA Kit
E-CL-H1293	Human SEPP1 (Selenoprotein P, Plasma) CLIA Kit
E-CL-H1294	Human ALB (Albumin) CLIA Kit
E-CL-H1295	Human SAA (Serum Amyloid A) CLIA Kit
E-CL-H1296	Human CR1 (Complement Receptor type 1) CLIA Kit
E-CL-H1297	Human MFG-E8 (Milk Fat Globule EGF Factor 8) CLIA Kit
E-CL-H1298	Human SERT (Serotonin Transporter) CLIA Kit
E-CL-H1299	Human TXB2 (Thromboxane B2) CLIA Kit
E-CL-H1300	Human FG (Fibrinogen) CLIA Kit
E-CL-H1301	Human FDP (Fibrinogen Degradation Product) CLIA Kit
E-CL-H1302	Human FPA (Fibrinopeptide A) CLIA Kit
E-CL-H1303	Human FPB (Fibrinopeptide B) CLIA Kit
E-CL-H1304	Human PAF (Platelet Activating Factor) CLIA Kit
E-CL-H1305	Human GPIIb/IIIa (Platelet Membrane Glycoprotein IIb/IIIa) CLIA Kit
E-CL-H1306	Human PKP1 (Plakophilin 1) CLIA Kit
E-CL-H1307	Human PFKP (Phosphofructokinase, Platelet) CLIA Kit
E-CL-H1308	Human PDGFC (Platelet Derived Growth Factor C) CLIA Kit
E-CL-H1309	Human PDGFD (Platelet Derived Growth Factor D) CLIA Kit
E-CL-H1310	Human PDGFR α (Platelet Derived Growth Factor Receptor Alpha) CLIA Kit

Cat.NO	Product Name
E-CL-H1311	Human PF4 (Platelet Factor 4) CLIA Kit
E-CL-H1312	Human PF4V1 (Platelet Factor 4 Variant 1) CLIA Kit
E-CL-H1313	Human NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) CLIA Kit
E-CL-H1314	Human OLAB (Anti-Oxidized Low Density Lipoprotein Antibody) CLIA Kit
E-CL-H1315	Human TRY (Trypsin) CLIA Kit
E-CL-H1316	Human TAP (Trypsinogen Activation Peptide) CLIA Kit
E-CL-H1317	Human PI (Proinsulin) CLIA Kit
E-CL-H1318	Human AMY2 (Amylase Alpha 2, Pancreatic) CLIA Kit
E-CL-H1319	Human PP (Pancreatic Polypeptide) CLIA Kit
E-CL-H1320	Human GC (Glucagon) CLIA Kit
E-CL-H1321	Human GLP-2 (Glucagon Like Peptide 2) CLIA Kit
E-CL-H1322	Human PTF1 α (Pancreas Specific Transcription Factor 1 α) CLIA Kit
E-CL-H1323	Human REG1 α (Regenerating Islet Derived Protein 1 Alpha) CLIA Kit
E-CL-H1324	Human REG4 (Regenerating Islet Derived Protein 4) CLIA Kit
E-CL-H1325	Human SULT1A1 (Sulfotransferase 1A1) CLIA Kit
E-CL-H1326	Human PL (Pancreatic Lipase) CLIA Kit
E-CL-H1327	Human OSM (Oncostatin M) CLIA Kit
E-CL-H1328	Human DCN (Decorin) CLIA Kit
E-CL-H1329	Human PHB (Prohibitin) CLIA Kit

Cat.NO	Product Name
E-CL-H1330	Human CRY1 (Cryptochrome 1) CLIA Kit
E-CL-H1331	Human KNG1 (Kininogen 1) CLIA Kit
E-CL-H1332	Human F-TESTO (Free Testosterone) CLIA Kit
E-CL-H1333	Human fT4 (Free Thyroxine) CLIA Kit
E-CL-H1334	Human fT3 (Free Triiodothyronine) CLIA Kit
E-CL-H1335	Human FBN1 (Fibrillin 1) CLIA Kit
E-CL-H1336	Human Pg (Progesterone) CLIA Kit
E-CL-H1337	Human PGR (Progesterone Receptor) CLIA Kit
E-CL-H1338	Human EGR1 (Early Growth Response Protein 1) CLIA Kit
E-CL-H1339	Human NFYC (Nuclear Transcription Factor Y Gamma) CLIA Kit
E-CL-H1340	Human MUC5AC (Mucin 5 Subtype AC) CLIA Kit
E-CL-H1341	Human MUC5B (Mucin 5 Subtype B) CLIA Kit
E-CL-H1342	Human MUC7 (Mucin 7) CLIA Kit
E-CL-H1343	Human LpPLA2 (Phospholipase A2, Lipoprotein Associated) CLIA Kit
E-CL-H1344	Human LPL (Lipoprotein Lipase) CLIA Kit
E-CL-H1345	Human LBP (Lipopolysaccharide Binding Protein) CLIA Kit
E-CL-H1346	Human LXA4 (Lipoxin A4) CLIA Kit
E-CL-H1347	Human OSMR (Oncostatin M Receptor) CLIA Kit
E-CL-H1348	Human MK (Midkine) CLIA Kit

Cat.NO	Product Name
E-CL-H1349	Human PON1 (Paraoxonase 1) CLIA Kit
E-CL-H1350	Human HNP1-3 (Neutrophil Peptide 1-3) CLIA Kit
E-CL-H1351	Human NAP3 (Neutrophil Activating Protein 3) CLIA Kit
E-CL-H1352	Human TP63 (Tumor Protein p63) CLIA Kit
E-CL-H1353	Human TNF- β (Tumor Necrosis Factor Beta) CLIA Kit
E-CL-H1354	Human CDK1 (Cyclin Dependent Kinase 1) CLIA Kit
E-CL-H1355	Human CDK2 (Cyclin Dependent Kinase 2) CLIA Kit
E-CL-H1356	Human CDK4 (Cyclin Dependent Kinase 4) CLIA Kit
E-CL-H1357	Human CDK5 (Cyclin Dependent Kinase 5) CLIA Kit
E-CL-H1358	Human CDK7 (Cyclin Dependent Kinase 7) CLIA Kit
E-CL-H1359	Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) CLIA Kit
E-CL-H1360	Human Ntn1 (Netrin 1) CLIA Kit
E-CL-H1361	Human Ntn4 (Netrin 4) CLIA Kit
E-CL-H1362	Human MHCC/HLA-C (Major Histocompatibility Complex Class I C) CLIA Kit
E-CL-H1363	Human MHCE/HLA-E (Major Histocompatibility Complex Class I E) CLIA Kit
E-CL-H1364	Human TGM2 (Transglutaminase 2, Tissue) CLIA Kit
E-CL-H1365	Human TGM3 (Transglutaminase 3, Epidermal) CLIA Kit
E-CL-H1366	Human TGF- β 3 (Transforming Growth Factor β 3) CLIA Kit
E-CL-H1367	Human TTR (Transthyretin) CLIA Kit

Cat.NO	Product Name
E-CL-H1368	Human TFR/CD71 (Transferrin Receptor) CLIA Kit
E-CL-H1369	Human TFR2 (Transferrin Receptor 2) CLIA Kit
E-CL-H1370	Human TNPO1 (Transportin 1) CLIA Kit
E-CL-H1371	Human Pax (Paxillin) CLIA Kit
E-CL-H1372	Human PCX (Podocalyxin) CLIA Kit
E-CL-H1373	Human TPA (Tissue Polypeptide Antigen) CLIA Kit
E-CL-H1374	Human TPS (Tissue Polypeptide Specific Antigen) CLIA Kit
E-CL-H1375	Human TFPI2 (Tissue Factor Pathway Inhibitor 2) CLIA Kit
E-CL-H1376	Human SEPCR/PROCR (Soluble Endothelial Protein C Receptor) CLIA Kit
E-CL-H1377	Human proPSA (precursor of Prostate Specific Antigen) CLIA Kit
E-CL-H1378	Human UBA52 (Ubiquitin A 52 Residue Ribosomal Protein Fusion Product 1) CLIA Kit
E-CL-H1379	Human NES1 (Nesfatin 1) CLIA Kit
E-CL-H1380	Human PON3 (Paraoxonase 3) CLIA Kit
E-CL-H1381	Human DDT (D-Dopachrome Tautomerase) CLIA Kit
E-CL-H1382	Human UCHL1 (Ubiquitin Carboxyl Terminal Hydrolase L1) CLIA Kit
E-CL-H1383	Human UBR4 (Ubiquitin Protein Ligase E3 Component N-Recognin 4) CLIA Kit
E-CL-H1384	Human MX1 (Myxovirus Resistance 1) CLIA Kit
E-CL-H1385	Human SAA2 (Serum Amyloid A2) CLIA Kit
E-CL-H1386	Human UBE2A (Ubiquitin Conjugating Enzyme E2A) CLIA Kit

Cat.NO	Product Name
E-CL-H1387	Human SOD3 (Superoxide Dismutase 3, Extracellular) CLIA Kit
E-CL-H1388	Human BCAR1 (Breast Cancer Anti Estrogen Resistance 1) CLIA Kit
E-CL-H1389	Human PLA2G2D (Phospholipase A2, Group II D) CLIA Kit
E-CL-H1390	Human PLA2G2A (Phospholipase A2, Group II A) CLIA Kit
E-CL-H1391	Human UBAP2 (Ubiquitin Associated Protein 2) CLIA Kit
E-CL-H1392	Human MASP1 (Mannan Associated Serine Protease 1) CLIA Kit
E-CL-H1393	Human MASP2 (Mannan Associated Serine Protease 2) CLIA Kit
E-CL-H1394	Human USP8 (Ubiquitin Specific Peptidase 8) CLIA Kit
E-CL-H1395	Human DPP6 (Dipeptidyl Peptidase VI) CLIA Kit
E-CL-H1396	Human MRC2 (Mannose Receptor C Type 2) CLIA Kit
E-CL-H1397	Human MRC1 (Mannose Receptor C Type 1) CLIA Kit
E-CL-H1398	Human MRC1 (Mannose Receptor C Type 1 Like Protein 1) CLIA Kit
E-CL-H1399	Human GADL1 (Glutamate Decarboxylase Like Protein 1) CLIA Kit
E-CL-H1400	Human GAD1 (Glutamate Decarboxylase 1) CLIA Kit
E-CL-H1401	Human PEX2 (Peroxisomal Biogenesis Factor 2) CLIA Kit
E-CL-H1402	Human PPAR- δ (Peroxisome Proliferator Activated Receptor Delta) CLIA Kit
E-CL-H1403	Human NUP62 (Nucleoporin 62kDa) CLIA Kit
E-CL-H1404	Human PCNA (Proliferating Cell Nuclear Antigen) CLIA Kit
E-CL-H1405	Human RANK (Receptor Activator of Nuclear Factor Kappa B) CLIA Kit

Cat.NO	Product Name
E-CL-H1406	Human NUMA1 (Nuclear Mitotic Apparatus Protein 1) CLIA Kit
E-CL-H1407	Human IL-33 (Interleukin 33) CLIA Kit
E-CL-H1408	Human MCAM (Melanoma Cell Adhesion Molecule) CLIA Kit
E-CL-H1409	Human MIA1 (Melanoma Inhibitory Activity Protein 1) CLIA Kit
E-CL-H1410	Human MFI2 (Melanotransferrin) CLIA Kit
E-CL-H1411	Human SPTβN4 (Spectrin Beta, Non Erythrocytic 4) CLIA Kit
E-CL-H1412	Human IκBα (Inhibitory Subunit of NF Kappa B Alpha) CLIA Kit
E-CL-H1413	Human RARα (Retinoic Acid Receptor Alpha) CLIA Kit
E-CL-H1414	Human RNASEN (Ribonuclease III, Nuclear) CLIA Kit
E-CL-H1415	Human MCSP (Melanoma Associated Chondroitin Sulfate Proteoglycan) CLIA Kit
E-CL-H1416	Human αMSH (Alpha-Melanocyte Stimulating Hormone) CLIA Kit
E-CL-H1417	Human EPP (Erythrocyte Protoporphyrin) CLIA Kit
E-CL-H1418	Human THRα (Thyroid Hormone Receptor Alpha) CLIA Kit
E-CL-H1419	Human NFRκB (Nuclear Factor Related To Kappa B Binding Protein) CLIA Kit
E-CL-H1420	Human TNKS1 (Tankyrase 1) CLIA Kit
E-CL-H1421	Human TNKS2 (Tankyrase 2) CLIA Kit
E-CL-H1422	Human PARP4 (Poly ADP Ribose Polymerase 4) CLIA Kit
E-CL-H1423	Human Neopterin CLIA Kit
E-CL-H1424	Human ME20M (Melanoma Associated ME20) CLIA Kit

Cat.NO	Product Name
E-CL-H1425	Human PXR (Pregnane X Receptor) CLIA Kit
E-CL-H1426	Human RELB (V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) CLIA Kit
E-CL-H1427	Human IL-24 (Interleukin 24) CLIA Kit
E-CL-H1428	Human AIM2 (Absent In Melanoma 2) CLIA Kit
E-CL-H1429	Human TSPAN30/CD63 (Tetraspanin 30) CLIA Kit
E-CL-H1430	Human BANP (BTG3 Associated Nuclear Protein) CLIA Kit
E-CL-H1431	Human PRKG1 (Protein Kinase, cGMP Dependent Type I) CLIA Kit
E-CL-H1432	Human PDE4D (Phosphodiesterase 4D, cAMP Specific) CLIA Kit
E-CL-H1433	Human CKB (Creatine Kinase, Brain) CLIA Kit
E-CL-H1434	Human LOXL4 (Lysyl Oxidase Like Protein 4) CLIA Kit
E-CL-H1435	Human TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) CLIA Kit
E-CL-H1436	Human IL-35 (Interleukin 35) CLIA Kit
E-CL-H1437	Human TK1 (Thymidine Kinase 1, Soluble) CLIA Kit
E-CL-H1438	Human 12-LO (Arachidonate 12-Lipoxygenase) CLIA Kit
E-CL-H1439	Human VIL1 (Villin 1) CLIA Kit
E-CL-H1440	Human PI3/Elafin (Peptidase Inhibitor 3, Skin Derived) CLIA Kit
E-CL-H1441	Human TTPα (Alpha Tocopherol Transfer Protein) CLIA Kit
E-CL-H1442	Human FNDC5 (Fibronectin Type III Domain Containing Protein 5) CLIA Kit
E-CL-H1443	Human IL-21 (Interleukin 21) CLIA Kit

Cat.NO	Product Name
E-CL-H1444	Human UCRP (Ubiquitin Cross Reactive Protein) CLIA Kit
E-CL-H1445	Human IFN- α (Interferon Alpha) CLIA Kit
E-CL-H1446	Human DAG1 (Dystrophin Associated Glycoprotein 1) CLIA Kit
E-CL-H1447	Human HABP1/C1QBP (Hyaluronan Binding Protein 1) CLIA Kit
E-CL-H1448	Human Anti-CRP (Anti-C Reactive Protein Antibody) CLIA Kit
E-CL-H1449	Human PAEP (Progestagen-Associated Endometrial Protein) CLIA Kit
E-CL-H1450	Human IL-11 (Interleukin 11) CLIA Kit
E-CL-H1451	Human hs-CRP (high-sensitivity C-reactive protein) CLIA Kit
E-CL-R0001	Rat ACV-A (Activin A) CLIA Kit
E-CL-R0002	Rat BMP-2 (Bone Morphogenetic Protein 2) CLIA Kit
E-CL-R0003	Rat GRO α /CXCL1 (Growth Regulated Oncogene Alpha) CLIA Kit
E-CL-R0004	Rat ELN (Elastin) CLIA Kit
E-CL-R0005	Rat MIP-3 (Cytokine-induced neutrophil chemoattractant 3) CLIA Kit
E-CL-R0006	Rat Cyt-C (Cytochrome C) CLIA Kit
E-CL-R0007	Rat EPO (Erythropoietin) CLIA Kit
E-CL-R0008	Rat GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) CLIA Kit
E-CL-R0009	Rat IFN- γ (Interferon Gamma) CLIA Kit
E-CL-R0010	Rat IGF-1 (Insulin-Like Growth Factor 1) CLIA Kit
E-CL-R0011	Rat IL-1 α (Interleukin 1 Alpha) CLIA Kit

Cat.NO	Product Name
E-CL-R0012	Rat IL-1 β (Interleukin 1 Beta) CLIA Kit
E-CL-R0013	Rat IL-2 (Interleukin 2) CLIA Kit
E-CL-R0014	Rat IL-4 (Interleukin 4) CLIA Kit
E-CL-R0015	Rat IL-6 (Interleukin 6) CLIA Kit
E-CL-R0016	Rat IL-10 (Interleukin 10) CLIA Kit
E-CL-R0017	Rat ANP (Atrial Natriuretic Peptide) CLIA Kit
E-CL-R0018	Rat EPI (Epinephrine) CLIA Kit
E-CL-R0019	Rat TNF- α (Tumor Necrosis Factor Alpha) CLIA Kit
E-CL-R0020	Rat VEGF (Vascular Endothelial cell Growth Factor) CLIA Kit
E-CL-R0021	Rat CRP (C-Reactive Protein) CLIA Kit
E-CL-R0022	Rat INS (Insulin) CLIA Kit
E-CL-R0023	Rat MAU (Microalbuminuria) CLIA Kit
E-CL-R0024	Rat LH (Luteinizing Hormone) CLIA Kit
E-CL-R0025	Rat GH (Growth Hormone) CLIA Kit
E-CL-R0026	Rat REN (Renin) CLIA Kit
E-CL-R0027	Rat MT (Melatonin) CLIA Kit
E-CL-R0028	Rat C-P (C-peptide) CLIA Kit
E-CL-R0029	Rat T (Testosterone) CLIA Kit
E-CL-R0030	Rat ApoB100 (Apolipoprotein B100) CLIA Kit

Cat.NO	Product Name
E-CL-R0031	Rat LN (Laminin) CLIA Kit
E-CL-R0032	Rat COL-1 α 1(Collagen I α 1) CLIA Kit
E-CL-R0033	Rat MMP-10 (Matrix Metalloproteinase 10) CLIA Kit
E-CL-R0034	Rat MMP-12 (Matrix Metalloproteinase 12) CLIA Kit
E-CL-R0035	Rat MMP-13 (Matrix Metalloproteinase 13) CLIA Kit
E-CL-R0036	Rat sICAM-1/CD54 (Intercellular Adhesion Molecule-1) CLIA Kit
E-CL-R0037	Rat CT (Calcitonin) CLIA Kit
E-CL-R0038	Rat ACTH (Adrenocorticotrophic Hormone) CLIA Kit
E-CL-R0039	Rat DA (Dopamine) CLIA Kit
E-CL-R0040	Rat OPG (Osteoprotegerin) CLIA Kit
E-CL-R0041	Rat PRL (Prolactin) CLIA Kit
E-CL-R0042	Rat MYO (Myoglobin) CLIA Kit
E-CL-R0043	Rat Tn-T (Troponin T) CLIA Kit
E-CL-R0044	Rat Tn-I (Troponin I) CLIA Kit
E-CL-R0045	Rat PAPPA (Pregnancy Associated Plasma Protein A) CLIA Kit
E-CL-R0046	Rat ALP (Alkaline Phosphatase) CLIA Kit
E-CL-R0047	Rat NSE (Neuron Specific Enolase) CLIA Kit
E-CL-R0048	Rat GLP-1 (Glucagon Like Peptide 1) CLIA Kit
E-CL-R0049	Rat SP-A (Pulmonary Surfactant-Associated Protein A) CLIA Kit

Cat.NO	Product Name
E-CL-R0050	Rat NE (Neutrophil Elastase) CLIA Kit
E-CL-R0051	Rat IL-12/p40 (Interleukin 12 p40) CLIA Kit
E-CL-R0052	Rat IL-12 (Interleukin 12) CLIA Kit
E-CL-R0053	Rat MBP/MBL (Mannma Binding Protein/Mannan Binding Lectin) CLIA Kit
E-CL-R0054	Rat AQP-3 (Aquaporin 3) CLIA Kit
E-CL-R0055	Rat AQP-4 (Aquaporin 4) CLIA Kit
E-CL-R0056	Rat Arg (Arginase) CLIA Kit
E-CL-R0057	Rat AVPR1A (Arginine Vasopressin Receptor 1A) CLIA Kit
E-CL-R0058	Rat AhR (Aryl Hydrocarbon Receptor) CLIA Kit
E-CL-R0059	Rat AST (Aspartate Aminotransferase) CLIA Kit
E-CL-R0060	Rat SP-B (Pulmonary Surfatcant-Associated Protein B) CLIA Kit
E-CL-R0061	Rat ATXN1 (Ataxin 1) CLIA Kit
E-CL-R0062	Rat ABCA1 (ATP Binding Cassette Transporter A1) CLIA Kit
E-CL-R0063	Rat TGF- β 1 (Transforming Growth Factor β 1) CLIA Kit
E-CL-R0064	Rat FGF4 (Fibroblast Growth Factor 4) CLIA Kit
E-CL-R0065	Rat FGF6 (Fibroblast Growth Factor 6) CLIA Kit
E-CL-R0066	Rat FGF9 (Fibroblast Growth Factor 9) CLIA Kit
E-CL-R0067	Rat bFGF/FGF2 (Basic Fibroblast Growth Factor) CLIA Kit
E-CL-R0068	Rat BCL3 (B-Cell Leukemia/Lymphoma 3) CLIA Kit

Cat.NO	Product Name
E-CL-R0069	Rat BCL-2 (B-Cell Leukemia/Lymphoma 2) CLIA Kit
E-CL-R0070	Rat BAX (BCL-2 Associated X Protein) CLIA Kit
E-CL-R0071	Rat BAD (BCL-4 Associated Death Promoter) CLIA Kit
E-CL-R0072	Rat BTC (Beta Cellulin) CLIA Kit
E-CL-R0073	Rat BACE1 (Beta Site APP Cleaving Enzyme 1) CLIA Kit
E-CL-R0074	Rat β -EPR (Beta-Endorphin Receptor) CLIA Kit
E-CL-R0075	Rat β -EP (Beta-Endorphin) CLIA Kit
E-CL-R0076	Rat BID (BH3 Interacting Domain Death Agonist) CLIA Kit
E-CL-R0077	Rat BTA (Bladder Tumor Antigen) CLIA Kit
E-CL-R0078	Rat BLC (B-Lymphocyte Chemoattractant) CLIA Kit
E-CL-R0079	Rat BALP (Bone Alkaline Phosphatase) CLIA Kit
E-CL-R0080	Rat BMP-1 (Bone Morphogenetic Protein 1) CLIA Kit
E-CL-R0081	Rat BMP-3 (Bone Morphogenetic Protein 3) CLIA Kit
E-CL-R0082	Rat BMP-4 (Bone Morphogenetic Protein 4) CLIA Kit
E-CL-R0083	Rat BMP-5 (Bone Morphogenetic Protein 5) CLIA Kit
E-CL-R0084	Rat BMP-6 (Bone Morphogenetic Protein 6) CLIA Kit
E-CL-R0085	Rat BMP-7 (Bone Morphogenetic Protein 7) CLIA Kit
E-CL-R0086	Rat BMPR-1A (Bone Morphogenetic Protein Receptor 1A) CLIA Kit
E-CL-R0087	Rat BMPR- II (Bone Morphogenetic Protein Receptor II) CLIA Kit

Cat.NO	Product Name
E-CL-R0088	Rat ALP-B (Bone-Specific Alkphase B) CLIA Kit
E-CL-R0089	Rat BK (Bradykinin) CLIA Kit
E-CL-R0090	Rat ANG2 (Angiopoietin 2) CLIA Kit
E-CL-R0091	Rat BNP (Brain Natriuretic Peptide) CLIA Kit
E-CL-R0092	Rat BRAK (Breast And Kidney Expressed Chemokine) CLIA Kit
E-CL-R0093	Rat VE-Cadherin (Vascular Endothelial Cadherin) CLIA Kit
E-CL-R0094	Rat NCAD (Neural Cadherin) CLIA Kit
E-CL-R0095	Rat CAN (Calcineurin) CLIA Kit
E-CL-R0096	Rat CGRP (Calcitonin Gene Related Peptide) CLIA Kit
E-CL-R0097	Rat β -Lg (Beta-Lactoglobulin) CLIA Kit
E-CL-R0098	Rat CASR (Calcium Sensing Receptor) CLIA Kit
E-CL-R0099	Rat CAMK2 (Calcium/Calmodulin-Dependent Protein Kinase- II) CLIA Kit
E-CL-R0100	Rat CAM (Calmodulin) CLIA Kit
E-CL-R0101	Rat CRT (Calreticulin) CLIA Kit
E-CL-R0102	Rat CR (Calretinin) CLIA Kit
E-CL-R0103	Rat CALU (Calumenin) CLIA Kit
E-CL-R0104	Rat CDT (Carbohydrate-Deficient Transferrin) CLIA Kit
E-CL-R0105	Rat CA1 (Carbonic Anhydrase I) CLIA Kit
E-CL-R0106	Rat CA2 (Carbonic Anhydrase II) CLIA Kit

Cat.NO	Product Name
E-CL-R0107	Rat P I CP (Carboxyterminal Propeptide of Type I Procollagen) CLIA Kit
E-CL-R0108	Rat CEA (Carcinoembryonic Antigen) CLIA Kit
E-CL-R0109	Rat cTnT/TNNT2 (Troponin T Type 2, Cardiac) CLIA Kit
E-CL-R0110	Rat RANK (Receptor Activator of Nuclear Factor Kappa B) CLIA Kit
E-CL-R0111	Rat α FP (Alpha-Fetoprotein) CLIA Kit
E-CL-R0112	Rat CT1 (Cardiotrophin 1) CLIA Kit
E-CL-R0113	Rat CLCF1 (Cardiotrophin Like Cytokine Factor 1) CLIA Kit
E-CL-R0114	Rat 6-keto-PGF1 α (6-keto-prostaglandin F1 α) CLIA Kit
E-CL-R0115	Rat CACT (Carnitine-Acylcarnitine Translocase) CLIA Kit
E-CL-R0116	Rat COMP (Cartilage Oligomeric Matrix Protein) CLIA Kit
E-CL-R0117	Rat CASP3 (Caspase 3) CLIA Kit
E-CL-R0118	Rat CASP7 (Caspase 7) CLIA Kit
E-CL-R0119	Rat CASP9 (Caspase 9) CLIA Kit
E-CL-R0120	Rat CASP14 (Caspase 14) CLIA Kit
E-CL-R0121	Rat ET-1 (Endothelin 1) CLIA Kit
E-CL-R0122	Rat CTSD (Cathepsin D) CLIA Kit
E-CL-R0123	Rat CTSK (Cathepsin K) CLIA Kit
E-CL-R0124	Rat CTSL (Cathepsin L) CLIA Kit
E-CL-R0125	Rat CAV-1 (Caveolin-1) CLIA Kit

Cat.NO	Product Name
E-CL-R0126	Rat CP (Ceruloplasmin) CLIA Kit
E-CL-R0127	Rat c-fos CLIA Kit
E-CL-R0128	Rat CCL1 (Chemokine C-C-Motif Ligand 1) CLIA Kit
E-CL-R0129	Rat CCRL1(Chemokine C-C-Motif Receptor Like Protein 1) CLIA Kit
E-CL-R0130	Rat SCD (Stearoyl Coenzyme A Desaturase) CLIA Kit
E-CL-R0131	Rat CXCL16 (Chemokine C-X-C-Motif Ligand 16) CLIA Kit
E-CL-R0132	Rat CKLF (Chemokine Like Factor) CLIA Kit
E-CL-R0133	Rat CKLFSF4 (Chemokine Like Factor Superfamily 4) CLIA Kit
E-CL-R0134	Rat CKLFSF5 (Chemokine Like Factor Superfamily 5) CLIA Kit
E-CL-R0135	Rat CKLFSF6 (Chemokine Like Factor Superfamily 6) CLIA Kit
E-CL-R0136	Rat CCKAR (Cholecystokinin A Receptor) CLIA Kit
E-CL-R0137	Rat CETP (Cholest Erolester Transfer Protein) CLIA Kit
E-CL-R0138	Rat AGT (Angiotensinogen) CLIA Kit
E-CL-R0139	Rat CH (Cholesterol) CLIA Kit
E-CL-R0140	Rat CHRM2 (Cholinergic Receptor, Muscarinic 2) CLIA Kit
E-CL-R0141	Rat CG (Chorionic Gonadotrophin) CLIA Kit
E-CL-R0142	Rat CHGA (Chromogranin A) CLIA Kit
E-CL-R0143	Rat CNTF (Ciliary Neurotrophic Factor) CLIA Kit
E-CL-R0144	Rat CS (Citrate Synthase) CLIA Kit

Cat.NO	Product Name
E-CL-R0145	Rat CD3 (Cluster of Differentiation 3) CLIA Kit
E-CL-R0146	Rat CD30L (Cluster Of Differentiation 30 Ligand) CLIA Kit
E-CL-R0147	Rat CD40L (Cluster Of Differentiation 40 Ligand) CLIA Kit
E-CL-R0148	Rat CD8 (Cluster of Differentiation 8) CLIA Kit
E-CL-R0149	Rat CD97 (Cluster of Differentiation 97) CLIA Kit
E-CL-R0150	Rat CLU (Clusterin) CLIA Kit
E-CL-R0151	Rat F II (Coagulation Factor II)CLIA Kit
E-CL-R0152	Rat F III (Coagulation Factor III) CLIA Kit
E-CL-R0153	Rat F X (Coagulation Factor X) CLIA Kit
E-CL-R0154	Rat F X III (Coagulation Factor X III) CLIA Kit
E-CL-R0155	Rat F XII (Coagulation Factor XII) CLIA Kit
E-CL-R0156	Rat FIX (Coagulation Factor IX) CLIA Kit
E-CL-R0157	Rat FV (Coagulation Factor V) CLIA Kit
E-CL-R0158	Rat CFL1 (Cofilin 1, Non-Muscle) CLIA Kit
E-CL-R0159	Rat COL1 (Collagen Type I) CLIA Kit
E-CL-R0160	Rat COL2 (Collagen Type II) CLIA Kit
E-CL-R0161	Rat COL3 (Collagen Type III) CLIA Kit
E-CL-R0162	Rat COL4 (Collagen Type IV) CLIA Kit
E-CL-R0163	Rat COL1α2 (Collagen Type I Alpha 2) CLIA Kit

Cat.NO	Product Name
E-CL-R0164	Rat COL4 α 1 (Collagen Type IV Alpha 1) CLIA Kit
E-CL-R0165	Rat COL9 α 1 (Collagen Type IX Alpha 1) CLIA Kit
E-CL-R0166	Rat COL10 (Collagen Type X) CLIA Kit
E-CL-R0167	Rat Collagenase I CLIA Kit
E-CL-R0168	Rat Collagenase II CLIA Kit
E-CL-R0169	Rat OC/BGP (Osteocalcin) CLIA Kit
E-CL-R0170	Rat BSP (Bone Sialoprotein) CLIA Kit
E-CL-R0171	Rat GCSFR (Colony Stimulating Factor Receptor, Granulocyte) CLIA Kit
E-CL-R0172	Rat TMS β 4 (Thymosin Beta 4) CLIA Kit
E-CL-R0173	Rat CTGF (Connective Tissue Growth Factor) CLIA Kit
E-CL-R0174	Rat CX26 (Connexin 26) CLIA Kit
E-CL-R0175	Rat CX40 (Connexin 40) CLIA Kit
E-CL-R0176	Rat CX43 (Connexin 43) CLIA Kit
E-CL-R0177	Rat CNTN1 (Contactin 1) CLIA Kit
E-CL-R0178	Rat CORT (Corticosterone) CLIA Kit
E-CL-R0179	Rat CRH (Corticotropin Releasing Hormone) CLIA Kit
E-CL-R0180	Rat CKMT1B (Creatine Kinase 1B, Mitochondrial) CLIA Kit
E-CL-R0181	Rat CKM (Creatine Kinase, Muscle) CLIA Kit
E-CL-R0182	Rat Cr (Creatinine) CLIA Kit

Cat.NO	Product Name
E-CL-R0183	Rat NTXI (Cross Linked N-Telopeptide of Type I Collagen) CLIA Kit
E-CL-R0184	Rat Anti-CIIAb (Anti-Type II Collagen Antibody) CLIA Kit
E-CL-R0185	Rat ICTP (Cross-linked C-Telopeptide of Type I Collagen) CLIA Kit
E-CL-R0186	Rat CASP8 (Caspase 8) CLIA Kit
E-CL-R0187	Rat CTX III (Cross-linked C-Telopeptides Of Type III Collagen) CLIA Kit
E-CL-R0188	Rat CNP (C-Type Natriuretic Peptide)CLIA Kit
E-CL-R0189	Rat CTACK (Cutaneous T-Cell Attracting Chemokine) CLIA Kit
E-CL-R0190	Rat CXCR3 (CXC-Chemokine Receptor 3) CLIA Kit
E-CL-R0191	Rat TPS (Tryptase) CLIA Kit
E-CL-R0192	Rat CREB (Cyclic AMP Response Element Binding Protein) CLIA Kit
E-CL-R0193	Rat CCND1 (Cyclin D1) CLIA Kit
E-CL-R0194	Rat CCND2 (Cyclin D2) CLIA Kit
E-CL-R0195	Rat CCND3 (Cyclin D3) CLIA Kit
E-CL-R0196	Rat CYPA (Cyclophilin A) CLIA Kit
E-CL-R0197	Rat CYPB (Cyclophilin B) CLIA Kit
E-CL-R0198	Rat β -actin (Beta actin) CLIA Kit
E-CL-R0199	Rat Cys-C (Cystatin C) CLIA Kit
E-CL-R0200	Rat CYSLTR2 (Cysteinyl Leukotriene Receptor 2) CLIA Kit
E-CL-R0201	Rat FGF7 (Fibroblast Growth Factor 7) CLIA Kit

Cat.NO	Product Name
E-CL-R0202	Rat CYP2E1 (Cytochrome P450, family 2, subfamily E, polypeptide 1) CLIA Kit
E-CL-R0203	Rat CYP450 (Cytochrome P450) CLIA Kit
E-CL-R0204	Rat Cpla2 (Cytosolic Phospholipase A2) CLIA Kit
E-CL-R0205	Rat CVL (Cytovillin) CLIA Kit
E-CL-R0206	Rat DCN (Decorin) CLIA Kit
E-CL-R0207	Rat CA9 (Carbonic Anhydrase IX) CLIA Kit
E-CL-R0208	Rat DHEA-S7 (Dehydroepiandrosterone S7) CLIA Kit
E-CL-R0209	Rat DHEA-S (Dehydroepiandrosterone Sulfate) CLIA Kit
E-CL-R0210	Rat DHEA (Dehydroepiandrosterone) CLIA Kit
E-CL-R0211	Rat DPD (Deoxypyridinoline) CLIA Kit
E-CL-R0212	Rat DNASE1 (Deoxyribonuclease I) CLIA Kit
E-CL-R0213	Rat ADP/Acrp30 (Adiponectin) CLIA Kit
E-CL-R0214	Rat Des (Desmin) CLIA Kit
E-CL-R0215	Rat DAO (Diamine Oxidase) CLIA Kit
E-CL-R0216	Rat DLD (Dihydrolipoyl Dehydrogenase) CLIA Kit
E-CL-R0217	Rat DPP IV (Dipeptidyl Peptidase IV) CLIA Kit
E-CL-R0218	Rat D-LDH (D-Lactate Dehydrogenase) CLIA Kit
E-CL-R0219	Rat DRD1 (Dopamine Receptor D1) CLIA Kit
E-CL-R0220	Rat DRD2 (Dopamine Receptor D2) CLIA Kit

Cat.NO	Product Name
E-CL-R0221	Rat DRD3 (Dopamine Receptor D3) CLIA Kit
E-CL-R0222	Rat DRD4 (Dopamine Receptor D4) CLIA Kit
E-CL-R0223	Rat DAT (Dopamine Transporter) CLIA Kit
E-CL-R0224	Rat TM (Thrombin) CLIA Kit
E-CL-R0225	Rat DMD (Dystrophin) CLIA Kit
E-CL-R0226	Rat E-Cad (E-Cadherin) CLIA Kit
E-CL-R0227	Rat E2 (Estradiol) CLIA Kit
E-CL-R0228	Rat EGR1 (Early Growth Response Protein 1) CLIA Kit
E-CL-R0229	Rat EGR2 (Early Growth Response Protein 2) CLIA Kit
E-CL-R0230	Rat AChE (Acetylcholinesterase) CLIA Kit
E-CL-R0231	Rat MPIF2 (Myeloid Progenitor Inhibitory Factor 2) CLIA Kit
E-CL-R0232	Rat EG-VEGF (Endocrine Gland Vascular Endothelial Growth Factor) CLIA Kit
E-CL-R0233	Rat ALB (Albumin) CLIA Kit
E-CL-R0234	Rat ERN1 (Endoplasmic Reticulum To Nucleus Signalling 1) CLIA Kit
E-CL-R0235	Rat ES (Endostatin) CLIA Kit
E-CL-R0236	Rat EL (Endothelial lipase) CLIA Kit
E-CL-R0237	Rat NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) CLIA Kit
E-CL-R0238	Rat ECF/CCL11 (Eosinophil Chemotactic Factor) CLIA Kit
E-CL-R0239	Rat EGF (Epidermal Growth Factor) CLIA Kit

Cat.NO	Product Name
E-CL-R0240	Rat ENA-78 (Epithelial Neutrophil Activating Peptide 78) CLIA Kit
E-CL-R0241	Rat CASP1 (Caspase 1) CLIA Kit
E-CL-R0242	Rat ERK1(Extracellular Signal Regulated Kinase 1) CLIA Kit
E-CL-R0243	Rat FAS/CD95 (Factor Related Apoptosis) CLIA Kit
E-CL-R0244	Rat ACS (Fatty Acyl-CoA synthetase) CLIA Kit
E-CL-R0245	Rat Fga (Fibrinogen Alpha) CLIA Kit
E-CL-R0246	Rat FDP (Fibrinogen Degradation Product) CLIA Kit
E-CL-R0247	Rat FRS3 (Fibroblast Growth Factor Receptor Substrate 3) CLIA Kit
E-CL-R0248	Rat FGF-10 (Fibroblast Growth Factor-10) CLIA Kit
E-CL-R0249	Rat FST (Follistatin) CLIA Kit
E-CL-R0250	Rat CX3CL1 (Chemokine C-X3-C-Motif Ligand 1) CLIA Kit
E-CL-R0251	Rat CG β (Chorionic Gonadotropin Beta Polypeptide) CLIA Kit
E-CL-R0252	Rat ATGL (Adipose Triglyceride Lipase) CLIA Kit
E-CL-R0253	Rat F-TESTO (Free Testosterone) CLIA Kit
E-CL-R0254	Rat fT4 (Free Thyroxine) CLIA Kit
E-CL-R0255	Rat FSH (Follicle Stimulating Hormone) CLIA Kit
E-CL-R0256	Rat FUCA (Alpha-L-Fucosidase, Tissue) CLIA Kit
E-CL-R0257	Rat GNB1 (G Protein Beta 1) CLIA Kit
E-CL-R0258	Rat GL β (Galactosidase, Beta) CLIA Kit

Cat.NO	Product Name
E-CL-R0259	Rat FASN (Fatty Acid Synthase) CLIA Kit
E-CL-R0260	Rat GAL (Galanin) CLIA Kit
E-CL-R0261	Rat GAL1 (Galectin 1) CLIA Kit
E-CL-R0262	Rat GAL2 (Galectin 2) CLIA Kit
E-CL-R0263	Rat GAL3 (Galectin 3) CLIA Kit
E-CL-R0264	Rat GAL4 (Galectin 4) CLIA Kit
E-CL-R0265	Rat GAL7 (Galectin 7) CLIA Kit
E-CL-R0266	Rat GAL8 (Galectin 8) CLIA Kit
E-CL-R0267	Rat GAL9 (Galectin 9) CLIA Kit
E-CL-R0268	Rat γ -GT1 (Gamma Glutamyltransferase 1) CLIA Kit
E-CL-R0269	Rat GIP (Gastric Inhibitory Polypeptide) CLIA Kit
E-CL-R0270	Rat GT-34 (Gastrin 34) CLIA Kit
E-CL-R0271	Rat CA199 (Gastrointestinalcancer Marker-CA199) CLIA Kit
E-CL-R0272	Rat GKN2 (Gastrokine 2) CLIA Kit
E-CL-R0273	Rat GATA3 (GATA Binding Protein 3) CLIA Kit
E-CL-R0274	Rat GATA4 (GATA Binding Protein 4) CLIA Kit
E-CL-R0275	Rat GS (Gelsolin) CLIA Kit
E-CL-R0276	Rat G-6-Pase/G6PT/GSD (Glucose 6 Phosphate) CLIA Kit
E-CL-R0277	Rat GMFB (Glia Maturation Factor, Beta) CLIA Kit

Cat.NO	Product Name
E-CL-R0278	Rat GDNF (Glial Cell Line Derived Neurotrophic Factor) CLIA Kit
E-CL-R0279	Rat GR (Glucocorticoid Receptor) CLIA Kit
E-CL-R0280	Rat GR α (Glucocorticoid Receptor α) CLIA Kit
E-CL-R0281	Rat GR β (Glucocorticoid Receptor β) CLIA Kit
E-CL-R0282	Rat GC (Glucagon) CLIA Kit
E-CL-R0283	Rat GCK (Glucokinase) CLIA Kit
E-CL-R0284	Rat gp130 (Glucoprotein 130) CLIA Kit
E-CL-R0285	Rat G6PD (Glucose 6 Phosphate Dehydrogenase) CLIA Kit
E-CL-R0286	Rat GLUT3 (Glucose Transporter 3) CLIA Kit
E-CL-R0287	Rat GLUT4 (Glucose Transporter 4) CLIA Kit
E-CL-R0288	Rat GIP (Glucose-Dependent Insulin-Releasing Polypeptide) CLIA Kit
E-CL-R0289	Rat GUSb (Glucuronidase, Beta) CLIA Kit
E-CL-R0290	Rat GCLM (Glutamate Cysteine Ligase, Modifier Subunit) CLIA Kit
E-CL-R0291	Rat GAD2 (Glutamate Decarboxylase 2, Acid) CLIA Kit
E-CL-R0292	Rat GDH (Glutamate Dehydrogenase) CLIA Kit
E-CL-R0293	Rat CCK-8 (Cholecystokinin 8) CLIA Kit
E-CL-R0294	Rat GSTA1 (Glutathione S Transferase Alpha 1) CLIA Kit
E-CL-R0295	Rat GSTO1 (Glutathione S Transferase Omega 1) CLIA Kit
E-CL-R0296	Rat GSTT1 (Glutathione S Transferase Theta 1) CLIA Kit

Cat.NO	Product Name
E-CL-R0297	Rat GSTT2 (Glutathione S Transferase Theta 2) CLIA Kit
E-CL-R0298	Rat HbA1C (Glycosylated Hemoglobin/Hemoglobin A1c) CLIA Kit
E-CL-R0299	Rat PYGB (Glycogen Phosphorylase, Brain) CLIA Kit
E-CL-R0300	Rat PYGM (Glycogen Phosphorylase, Muscle) CLIA Kit
E-CL-R0301	Rat GSK3A (Glycogen Synthase Kinase 3 Alpha) CLIA Kit
E-CL-R0302	Rat GSK (Glycogen Synthase Kinase) CLIA Kit
E-CL-R0303	Rat GnRHR (Gonadotropin-Releasing Hormone Receptor) CLIA Kit
E-CL-R0304	Rat GnRH (Gonadotropin-Releasing Hormone) CLIA Kit
E-CL-R0305	Rat GRN (Granulin) CLIA Kit
E-CL-R0306	Rat GCP-2 (Granulocyte Chemotactic Protein 2) CLIA Kit
E-CL-R0307	Rat G-CSF (Granulocyte Colony Stimulating Factor) CLIA Kit
E-CL-R0308	Rat Gzms-B(Granzymes B) CLIA Kit
E-CL-R0309	Rat GREM1 (GremLin 1) CLIA Kit
E-CL-R0310	Rat GDF1 (Growth Differentiation Factor 1) CLIA Kit
E-CL-R0311	Rat GDF2 (Growth Differentiation Factor 2) CLIA Kit
E-CL-R0312	Rat GDF3 (Growth Differentiation Factor 3) CLIA Kit
E-CL-R0313	Rat GDF5 (Growth Differentiation Factor 5) CLIA Kit
E-CL-R0314	Rat GDF11 (Growth Differentiation Factor 11) CLIA Kit
E-CL-R0315	Rat GDF15 (Growth Differentiation Factor 15) CLIA Kit

Cat.NO	Product Name
E-CL-R0316	Rat GRB10 (Growth Factor Receptor Bound Protein 10) CLIA Kit
E-CL-R0317	Rat GRB14 (Growth Factor Receptor Bound Protein 14) CLIA Kit
E-CL-R0318	Rat GHRH (Growth Hormone Releasing Hormone) CLIA Kit
E-CL-R0319	Rat GT (Gastrin) CLIA Kit
E-CL-R0320	Rat HP (Haptoglobin) CLIA Kit
E-CL-R0321	Rat HSPA4 (Heat Shock 70kDa Protein 4) CLIA Kit
E-CL-R0322	Rat HSF1 (Heat Shock Factor 1) CLIA Kit
E-CL-R0323	Rat HSP-27 (Heat Shock Protein 27) CLIA Kit
E-CL-R0324	Rat HSP-40 (Heat Shock Protein 40) CLIA Kit
E-CL-R0325	Rat HSP-60 (Heat Shock Protein 60) CLIA Kit
E-CL-R0326	Rat HSP-70 (Heat Shock Protein 70) CLIA Kit
E-CL-R0327	Rat HSP-90 (Heat Shock Protein 90) CLIA Kit
E-CL-R0328	Rat HSP GP96 (Heat Shock Protein Glycoprotein 96) CLIA Kit
E-CL-R0329	Rat HSP47 (Heat Shock Protein) CLIA Kit
E-CL-R0330	Rat HSF2 (Heat Shock Transcription Factor 2) CLIA Kit
E-CL-R0331	Rat SHH (Hedgehog Homolog, Sonic) CLIA Kit
E-CL-R0332	Rat Hp-IgG (Helicobacter Pylori IgG) CLIA Kit
E-CL-R0333	Rat HO-1 (Heme Oxygenase 1) CLIA Kit
E-CL-R0334	Rat HO-2 (heme oxygenase 2) CLIA Kit

Cat.NO	Product Name
E-CL-R0335	Rat HPX (Hemopexin) CLIA Kit
E-CL-R0336	Rat HSPG (Heparan Sulphate Proteoglycan) CLIA Kit
E-CL-R0337	Rat HPA (Heparanase) CLIA Kit
E-CL-R0338	Rat HBEGF (Heparin Binding Epidermal Growth Factor Like Growth Factor) CLIA Kit
E-CL-R0339	Rat HC II (Heparin Cofactor II) CLIA Kit
E-CL-R0340	Rat LIPC (Lipase, Hepatic) CLIA Kit
E-CL-R0341	Rat HGF (Hepatocyte Growth Factor) CLIA Kit
E-CL-R0342	Rat HNF1a (Hepatocyte Nuclear Factor 1 Alpha) CLIA Kit
E-CL-R0343	Rat HDGF (Hepatoma Derived Growth Factor) CLIA Kit
E-CL-R0344	Rat Hpc (Hepcidin) CLIA Kit
E-CL-R0345	Rat HK1 (Hexokinase 1) CLIA Kit
E-CL-R0346	Rat HK (Hexokinase) CLIA Kit
E-CL-R0347	Rat HDL (High Density Lipoprotein) CLIA Kit
E-CL-R0348	Rat HMGB-1 (High Mobility Group Protein B1) CLIA Kit
E-CL-R0349	Rat hs-CRP (High Sensitivity C-Reactive Protein) CLIA Kit
E-CL-R0350	Rat COL1 (Collagen Type I) CLIA Kit
E-CL-R0351	Rat HABP1 (Hyaluronan Binding Protein 1) CLIA Kit
E-CL-R0352	Rat HIF-1 α (Hypoxia Inducible Factor 1 Alpha) CLIA Kit
E-CL-R0353	Rat HMG-CoA (3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase) CLIA Kit

Cat.NO	Product Name
E-CL-R0354	Rat IgA (Immunoglobulin A) CLIA Kit
E-CL-R0355	Rat IgE (Immunoglobulin E) CLIA Kit
E-CL-R0356	Rat IgG (Immunoglobulin G) CLIA Kit
E-CL-R0357	Rat IgM (Immunoglobulin M) CLIA Kit
E-CL-R0358	Rat NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) CLIA Kit
E-CL-R0359	Rat INHA (Inhibin A) CLIA Kit
E-CL-R0360	Rat ADH (Antidiuretic Hormone) CLIA Kit
E-CL-R0361	Rat INHbC (Inhibin Beta C) CLIA Kit
E-CL-R0362	Rat INHBP (Inhibin Binding Protein) CLIA Kit
E-CL-R0363	Rat IKB α (Inhibitory Subunit of NF- κ B α)CLIA Kit
E-CL-R0364	Rat IGFBP5 (Insulin Like Growth Factor Binding Protein 5) CLIA Kit
E-CL-R0365	Rat IGFBP6 (Insulin Like Growth Factor Binding Protein 6) CLIA Kit
E-CL-R0366	Rat INSR- β (Insulin Receptor Beta Subunit) CLIA Kit
E-CL-R0367	Rat IGF-2 (Insulin Like Growth Factor 2) CLIA Kit
E-CL-R0368	Rat IGFBP-1 (Insulin-Like Growth Factor Binding Protein 1) CLIA Kit
E-CL-R0369	Rat IGFBP-2 (Insulin-Like Growth Factor Binding Protein 2) CLIA Kit
E-CL-R0370	Rat IGFBP-3 (Insulin-Like Growth Factor Binding Protein 3) CLIA Kit
E-CL-R0371	Rat IGFBP-4 (Insulin-Like Growth Factor Binding Protein 4) CLIA Kit
E-CL-R0372	Rat I-PTH (Intact Parathormone) CLIA Kit

Cat.NO	Product Name
E-CL-R0373	Rat PDGF-AB (Platelet-Derived Growth Factor AB) CLIA Kit
E-CL-R0374	Rat PDGF-BB (Platelet-Derived Growth Factor-BB) CLIA Kit
E-CL-R0375	Rat SP (Substance P) CLIA Kit
E-CL-R0376	Rat TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) CLIA Kit
E-CL-R0377	Rat α -GST(α -Glutathione S-Transferases) CLIA Kit
E-CL-R0378	Rat ICAM-2 (Intercellular Adhesion Molecule 2) CLIA Kit
E-CL-R0379	Rat ICAM-3 (Intercellular Adhesion Molecule 3) CLIA Kit
E-CL-R0380	Rat IFN- α (Interferon α) CLIA Kit
E-CL-R0381	Rat IFN- β (Interferon Beta) CLIA Kit
E-CL-R0382	Rat IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) CLIA Kit
E-CL-R0383	Rat IFI16/p16 (Interferon-Inducible Protein 16) CLIA Kit
E-CL-R0384	Rat INF- α -Ab (Interferon- α Antibody) CLIA Kit
E-CL-R0385	Rat PEPCK (Phosphoenolpyruvate Carboxykinase) CLIA Kit
E-CL-R0386	Rat IL-1ra (Interleukin 1 Receptor Antagonist) CLIA Kit
E-CL-R0387	Rat LCAT (Lecithin-Cholesterol Acyltransferase) CLIA Kit
E-CL-R0388	Rat ALD (Aldosterone) CLIA Kit
E-CL-R0389	Rat IL-2R α (Interleukin-2 Receptor- α) CLIA Kit
E-CL-R0390	Rat IL-3 (Interleukin 3) CLIA Kit
E-CL-R0391	Rat FE (Ferritin) CLIA Kit

Cat.NO	Product Name
E-CL-R0392	Rat IL-5 (Interleukin 5) CLIA Kit
E-CL-R0393	Rat IL-7 (Interleukin 7) CLIA Kit
E-CL-R0394	Rat IL-8 (Interleukin 8) CLIA Kit
E-CL-R0395	Rat IL-9 (Interleukin 9) CLIA Kit
E-CL-R0396	Rat IL-11 (Interleukin 11) CLIA Kit
E-CL-R0397	Rat IL-13 (Interleukin 13) CLIA Kit
E-CL-R0398	Rat IL-15 (Interleukin 15) CLIA Kit
E-CL-R0399	Rat IL-16 (Interleukin 16) CLIA Kit
E-CL-R0400	Rat IL-17 (Interleukin 17) CLIA Kit
E-CL-R0401	Rat IL-18 (Interleukin 18) CLIA Kit
E-CL-R0402	Rat IL-21 (Interleukin 21) CLIA Kit
E-CL-R0403	Rat IL-23 (Interleukin 23) CLIA Kit
E-CL-R0404	Rat IL-27 (Interleukin 27) CLIA Kit
E-CL-R0405	Rat IL-32 (Interleukin 32) CLIA Kit
E-CL-R0406	Rat IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) CLIA Kit
E-CL-R0407	Rat MUC2 (Mucin 2) CLIA Kit
E-CL-R0408	Rat AR (Amphiregulin) CLIA Kit
E-CL-R0409	Rat KIM-1 (Kidney Injury Molecule 1) CLIA Kit
E-CL-R0410	Rat LDH (Lactate Dehydrogenase) CLIA Kit

Cat.NO	Product Name
E-CL-R0411	Rat LTF/LF (Lactoferrin) CLIA Kit
E-CL-R0412	Rat FN (Fibronectin) CLIA Kit
E-CL-R0413	Rat LDL (Low Density Lipoprotein) CLIA Kit
E-CL-R0414	Rat LRP-1 (Low-Density Lipoprotein-Receptor-Related Protein) CLIA Kit
E-CL-R0415	Rat LDL-IC (Low Density Lipoprotein Immunity Complex) CLIA Kit
E-CL-R0416	Rat LEP (Leptin) CLIA Kit
E-CL-R0417	Rat LIFR (Leukemia Inhibitory Factor Receptor) CLIA Kit
E-CL-R0418	Rat LIF (Leukemia Inhibitory Factor) CLIA Kit
E-CL-R0419	Rat LT-B4 (Leukotriene B4) CLIA Kit
E-CL-R0420	Rat LT-C4 (Leukotriene C4) CLIA Kit
E-CL-R0421	Rat LT-E4 (Leukotriene E4) CLIA Kit
E-CL-R0422	Rat LBP (Lipolysaccharide Binding Protein) CLIA Kit
E-CL-R0423	Rat LPL (Lipoprotein Lipase) CLIA Kit
E-CL-R0424	Rat Lp-a (Lipoprotein a) CLIA Kit
E-CL-R0425	Rat Lp-PL-A2 (Lipoprotein-Associated Phospholipase A2) CLIA Kit
E-CL-R0426	Rat LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) CLIA Kit
E-CL-R0427	Rat SELL (L-Selectin) CLIA Kit
E-CL-R0428	Rat LHRH (Luteinizing Hormone-Releasing Hormone) CLIA Kit
E-CL-R0429	Rat Lptn/XCL1 (Lymphotactin) CLIA Kit

Cat.NO	Product Name
E-CL-R0430	Rat M-CSF (Macrophage Colony-Stimulating Factor) CLIA Kit
E-CL-R0431	Rat MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) CLIA Kit
E-CL-R0432	Rat MIP-1 β (Macrophage Inflammatory Protein 1 Beta) CLIA Kit
E-CL-R0433	Rat PYGL (Glycogen Phosphorylase, Liver) CLIA Kit
E-CL-R0434	Rat MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) CLIA Kit
E-CL-R0435	Rat MIP-3 β (Macrophage Inflammatory Protein 3 β) CLIA Kit
E-CL-R0436	Rat MIF (Macrophage Migration Inhibitory Factor) CLIA Kit
E-CL-R0437	Rat MDC/CCL22 (Macrophage-Derived Chemokine) CLIA Kit
E-CL-R0438	Rat MHC I /AgB I /H-1 I (Major Histocompatibility Complex I) CLIA Kit
E-CL-R0439	Rat MHC II /AgB II /H-1 II (Major Histocompatibility Complex II) CLIA Kit
E-CL-R0440	Rat MHC III /AgB III /H-1 III (Major Histocompatibility Complex III) CLIA Kit
E-CL-R0441	Rat RETN (Resistin) CLIA Kit
E-CL-R0442	Rat MMP-1 (Matrix Metalloproteinase 1) CLIA Kit
E-CL-R0443	Rat MMP-2 (Matrix Metalloproteinase 2) CLIA Kit
E-CL-R0444	Rat MMP-3 (Matrix Metalloproteinase 3) CLIA Kit
E-CL-R0445	Rat MMP-17 (Matrix Metalloproteinase 17) CLIA Kit
E-CL-R0446	Rat MMP-24 (Matrix Metalloproteinase 24) CLIA Kit
E-CL-R0447	Rat MMP-7 (Matrix Metalloproteinase 7) CLIA Kit
E-CL-R0448	Rat MMP-8 (Matrix Metalloproteinase 8) CLIA Kit

Cat.NO	Product Name
E-CL-R0449	Rat MMP-9 (Matrix Metalloproteinase 9) CLIA Kit
E-CL-R0450	Rat ANG (Angiogenin) CLIA Kit
E-CL-R0451	Rat ANG1 (Angiopoietin 1) CLIA Kit
E-CL-R0452	Rat SLC (Secondary Lymphoid Tissue Chemokine) CLIA Kit
E-CL-R0453	Rat COL2 (Collagen Type II) CLIA Kit
E-CL-R0454	Rat MT 1 (Metallothionein1) CLIA Kit
E-CL-R0455	Rat MHB (Methemoglobin) CLIA Kit
E-CL-R0456	Rat Methylase CLIA Kit
E-CL-R0457	Rat MALB (Microalbumin) CLIA Kit
E-CL-R0458	Rat MCP-1 (Monocyte Chemotactic Protein 1) CLIA Kit
E-CL-R0459	Rat MCP-2 (Monocyte Chemotactic Protein 2) CLIA Kit
E-CL-R0460	Rat MCP-3 (Monocyte Chemotactic Protein 3) CLIA Kit
E-CL-R0461	Rat MCP-4 (Monocyte Chemotactic Protein 4) CLIA Kit
E-CL-R0462	Rat MTL (Motilin) CLIA Kit
E-CL-R0463	Rat AMH (Anti-Mullerian Hormone) CLIA Kit
E-CL-R0464	Rat M-ACHR (Muscarinic AcetylcholineReceptor) CLIA Kit
E-CL-R0465	Rat MBP (Myelin Basic Protein) CLIA Kit
E-CL-R0466	Rat AGER (Advanced Glycosylation End Product Specific Receptor) CLIA Kit
E-CL-R0467	Rat OXB (Orexin B) CLIA Kit

Cat.NO	Product Name
E-CL-R0468	Rat MYS (Myosin) CLIA Kit
E-CL-R0469	Rat NAGase (N-Acetyl Beta-D-Glucosaminidase) CLIA Kit
E-CL-R0470	Rat NN-T4 (Neonatal Thyroxine) CLIA Kit
E-CL-R0471	Rat NGF (Nerve Growth Factor) CLIA Kit
E-CL-R0472	Rat NGB (Neuroglobin) CLIA Kit
E-CL-R0473	Rat NPY (Neuropeptide Y) CLIA Kit
E-CL-R0474	Rat NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) CLIA Kit
E-CL-R0475	Rat NT-3 (Neurotrophin 3) CLIA Kit
E-CL-R0476	Rat NT-4 (Neurotrophin 4) CLIA Kit
E-CL-R0477	Rat NAP3 (Neutrophil Activating Protein 3) CLIA Kit
E-CL-R0478	Rat NAP-2 (Neutrophil Activating Protein-2) CLIA Kit
E-CL-R0479	Rat NGAL (Neutrophil Gelatinase Associated Lipocalin) CLIA Kit
E-CL-R0480	Rat NOX1 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) CLIA Kit
E-CL-R0481	Rat NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) CLIA Kit
E-CL-R0482	Rat NADPH (Nicotinamide Adenine Dinucleotide Phosphate) CLIA Kit
E-CL-R0483	Rat N-ACHR (Nicotinic Acetylcholine Receptor) CLIA Kit
E-CL-R0484	Rat NA/NE (Noradrenaline/Norepinephrine) CLIA Kit
E-CL-R0485	Rat NT-proBNP (N-Terminal Pro-Brain Natriuretic Peptide) CLIA Kit
E-CL-R0486	Rat PIIINP (N-Terminal Procollagen III Propeptide) CLIA Kit

Cat.NO	Product Name
E-CL-R0487	Rat NFIL3 (Nuclear Factor, Interleukin 3 Regulated) CLIA Kit
E-CL-R0488	Rat NF-κB (Nuclear Factor Kappa B) CLIA Kit
E-CL-R0489	Rat NF-κB p65(Nuclear Factor-κB P65) CLIA Kit
E-CL-R0490	Rat OBTF4 (Octamer Binding Transcription Factor 4) CLIA Kit
E-CL-R0491	Rat ITLN1 (Intelectin 1/Omentin) CLIA Kit
E-CL-R0492	Rat OSM (Oncostatin-M) CLIA Kit
E-CL-R0493	Rat OXA (Orexin A) CLIA Kit
E-CL-R0494	Rat ODC (Ornithine Decarboxylase) CLIA Kit
E-CL-R0495	Rat GROβ (Growth Regulated Oncogene Beta) CLIA Kit
E-CL-R0496	Rat GROγ (Growth Regulated Oncogene Gamma) CLIA Kit
E-CL-R0497	Rat ODF (Osteoclast Differentiation Factor) CLIA Kit
E-CL-R0498	Rat ON (Osteonectin) CLIA Kit
E-CL-R0499	Rat OPN (Osteopontin) CLIA Kit
E-CL-R0500	Rat OPGL (Osteoprotegerin Ligand) CLIA Kit
E-CL-R0501	Rat CA125 (Carbohydrate Antigen 125) CLIA Kit
E-CL-R0502	Rat OxLDL (Oxidized Lowdensity Lipoprotein) CLIA Kit
E-CL-R0503	Rat PUMA (P53 Upregulated Modulator of Apoptosis) CLIA Kit
E-CL-R0504	Rat PTH (Parathyroid Hormone) CLIA Kit
E-CL-R0505	Rat PTHR2 (Parathyroid Hormone Receptor 2) CLIA Kit

Cat.NO	Product Name
E-CL-R0506	Rat PTHRP (Parathyroid Hormone Related Protein) CLIA Kit
E-CL-R0507	Rat PTD (Pentosidine) CLIA Kit
E-CL-R0508	Rat PP (Pepsin) CLIA Kit
E-CL-R0509	Rat PYY (Peptide YY) CLIA Kit
E-CL-R0510	Rat PRF1 (Perforin 1) CLIA Kit
E-CL-R0511	Rat PMP3 (Peroxisomal Membrane Protein 3) CLIA Kit
E-CL-R0512	Rat PPAR- γ (Peroxisome Proliferator-Activated Receptor γ) CLIA Kit
E-CL-R0513	Rat P-PKC (Phospho Protein Kinase C) CLIA Kit
E-CL-R0514	Rat PERK (Phospho-Extracellular Signal-Regulated Kinase) CLIA Kit
E-CL-R0515	Rat PI3K2a (Phosphoinositide 3 Kinase 2a) CLIA Kit
E-CL-R0516	Rat pACCCase (Phosphorylated Acetyl-CoA Carboxylase) CLIA Kit
E-CL-R0517	Rat AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) CLIA Kit
E-CL-R0518	Rat PI3K (Phosphotylinosital 3 Kinase) CLIA Kit
E-CL-R0519	Rat PICP (Procollagen I C-terminal Propeptide) CLIA Kit
E-CL-R0520	Rat PGF (Placental Growth Factor) CLIA Kit
E-CL-R0521	Rat PAP (Plasmin-Antiplasmin Complex) CLIA Kit
E-CL-R0522	Rat PAI-1 (Plasminogen Activator Inhibitor 1) CLIA Kit
E-CL-R0523	Rat PAI-2 (Plasminogen Activator Inhibitor 2) CLIA Kit
E-CL-R0524	Rat PAI (Plasminogen Activator Inhibitor) CLIA Kit

Cat.NO	Product Name
E-CL-R0525	Rat tPA (Plasminogen Activator, Tissue) CLIA Kit
E-CL-R0526	Rat Plg (Plasminogen) CLIA Kit
E-CL-R0527	Rat PAF (Platelet Activating Factor) CLIA Kit
E-CL-R0528	Rat PA-IgG/M/A (Platelet Antibodies IgG/M/A) CLIA Kit
E-CL-R0529	Rat PDGF (platelet-derived growth factor) CLIA Kit
E-CL-R0530	Rat PDGFA (Platelet Derived Growth Factor Subunit A) CLIA Kit
E-CL-R0531	Rat PDGFB (Platelet Derived Growth Factor Subunit B) CLIA Kit
E-CL-R0532	Rat PDGF-C (Platelet Derived Growth Factor C) CLIA Kit
E-CL-R0533	Rat PDGF-D (Platelet Derived Growth Factor D) CLIA Kit
E-CL-R0534	Rat PDGFRL (Platelet Derived Growth Factor Receptor Like Protein) CLIA Kit
E-CL-R0535	Rat PDGFsR- α (Platelet-Derived Growth Factor Soluble Receptor α) CLIA Kit
E-CL-R0536	Rat PF3 (Platelet Factor 3) CLIA Kit
E-CL-R0537	Rat PF-4 (Platelet Factor-4) CLIA Kit
E-CL-R0538	Rat GP- II b III a (Platelet Membrane Glycoprotein II b III a) CLIA Kit
E-CL-R0539	Rat PECAM1 (Platelet/Endothelial Cell Adhesion Molecule) CLIA Kit
E-CL-R0540	Rat PTN (Pleiotrophin) CLIA Kit
E-CL-R0541	Rat PCX (Podocalyxin) CLIA Kit
E-CL-R0542	Rat PC4 (Positive Cofactor 4) CLIA Kit
E-CL-R0543	Rat PA/PLB (Pre-Albumin) CLIA Kit

Cat.NO	Product Name
E-CL-R0544	Rat PRO-ANP (Pro Atrial Natriuretic Peptide) CLIA Kit
E-CL-R0545	Rat PCT (Procalcitonin) CLIA Kit
E-CL-R0546	Rat P I NP (Procollagen I N-Terminal Peptide) CLIA Kit
E-CL-R0547	Rat P III NT (Procollagen III N-Terminal Peptide)CLIA Kit
E-CL-R0548	Rat Pg (Progesterone) CLIA Kit
E-CL-R0549	Rat PGR (Progesterone Receptor) CLIA Kit
E-CL-R0550	Rat FSTL1 (Follistatin Like Protein 1) CLIA Kit
E-CL-R0551	Rat PHB (Prohibitin) CLIA Kit
E-CL-R0552	Rat PI (Proinsulin) CLIA Kit
E-CL-R0553	Rat PK2 (Prokineticin 2) CLIA Kit
E-CL-R0554	Rat PKR1 (Prokineticin Receptor 1) CLIA Kit
E-CL-R0555	Rat PCNA (Proliferating Cell Nuclear Antigen) CLIA Kit
E-CL-R0556	Rat GZMK (Granzyme K) CLIA Kit
E-CL-R0557	Rat PRCP (Prolylcarboxypeptidase) CLIA Kit
E-CL-R0558	Rat EP2 (Prostaglandin E Receptor 2) CLIA Kit
E-CL-R0559	Rat PTGES (Prostaglandin E Synthase, Microsomal) CLIA Kit
E-CL-R0560	Rat PGE1 (Prostaglandin E1)CLIA Kit
E-CL-R0561	Rat NRG-1 (Neuregulin 1) CLIA Kit
E-CL-R0562	Rat PTGS1 (Prostaglandin Endoperoxide Synthase 1) CLIA Kit

Cat.NO	Product Name
E-CL-R0563	Rat PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) CLIA Kit
E-CL-R0564	Rat PGF2 α (Prostaglandin F2 Alpha) CLIA Kit
E-CL-R0565	Rat PSA (Prostate Specific Antigen) CLIA Kit
E-CL-R0566	Rat PAR2 (Protease Activated Receptor 2) CLIA Kit
E-CL-R0567	Rat PRSS1 (Protease, Serine, 1) CLIA Kit
E-CL-R0568	Rat PROC (Protein C) CLIA Kit
E-CL-R0569	Rat PROCR (Protein C Receptor, Endothelial) CLIA Kit
E-CL-R0570	Rat PKB (Protein Kinase B) CLIA Kit
E-CL-R0571	Rat PKCB1 (Protein Kinase C Beta 1) CLIA Kit
E-CL-R0572	Rat PKCD (Protein Kinase C Delta) CLIA Kit
E-CL-R0573	Rat PKC (Protein Kinase C) CLIA Kit
E-CL-R0574	Rat PP (Protein Phosphatase) CLIA Kit
E-CL-R0575	Rat PR3 (Proteinase 3) CLIA Kit
E-CL-R0576	Rat SELP (P-Selectin) CLIA Kit
E-CL-R0577	Rat SPC (Pulmonary Surfactant Associated Protein C) CLIA Kit
E-CL-R0578	Rat SPD (Pulmonary Surfactant Associated Protein D) CLIA Kit
E-CL-R0579	Rat PDHA (Pyruvate Dehydrogenase Alpha) CLIA Kit
E-CL-R0580	Rat PDK4 (Pyruvate Dehydrogenase Kinase Isozyme 4) CLIA Kit
E-CL-R0581	Rat PK (Pyruvate Kinase, Liver And RBC) CLIA Kit

Cat.NO	Product Name
E-CL-R0582	Rat M2-PK (Pyruvate Kinase,M2) CLIA Kit
E-CL-R0583	Rat NRG-2 (Neuregulin 2) CLIA Kit
E-CL-R0584	Rat RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) CLIA Kit
E-CL-R0585	Rat REG1 α (Regenerating Islet Derived Protein 1 Alpha) CLIA Kit
E-CL-R0586	Rat RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) CLIA Kit
E-CL-R0587	Rat RBP4 (Retinol Binding Protein 4, Plasma) CLIA Kit
E-CL-R0588	Rat RNASET2 (Ribonuclease T2) CLIA Kit
E-CL-R0589	Rat RUNX2 (Runt Related Transcription Factor 2) CLIA Kit
E-CL-R0590	Rat S100A11 (S100 Calcium Binding Protein A11) CLIA Kit
E-CL-R0591	Rat S100B (S100 Calcium Binding Protein B) CLIA Kit
E-CL-R0592	Rat S100A6 (S100 Calcium-Binding Protein A6) CLIA Kit
E-CL-R0593	Rat FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) CLIA Kit
E-CL-R0594	Rat SPLA2 (Secreted Phospholipase A2) CLIA Kit
E-CL-R0595	Rat SCT (Secretin) CLIA Kit
E-CL-R0596	Rat sIgA (Secretory Immunoglobulin A)CLIA Kit
E-CL-R0597	Rat SLPI (Secretory Leukocyte Peptidase Inhibitor) CLIA Kit
E-CL-R0598	Rat SELE (Selectin, Endothelium) CLIA Kit
E-CL-R0599	Rat SEPP1 (Selenoprotein P) CLIA Kit
E-CL-R0600	Rat SEMA7A (Semaphorin 7A) CLIA Kit

Cat.NO	Product Name
E-CL-R0601	Rat STK (Serine/Threonine Protein Kinase) CLIA Kit
E-CL-R0602	Rat SAA (Serum Amyloid A) CLIA Kit
E-CL-R0603	Rat SHBG (Sex Hormone-Binding Globulin) CLIA Kit
E-CL-R0604	Rat sCD14 (Soluble Cluster of Differentiation 14) CLIA Kit
E-CL-R0605	Rat EPCR (Endothelial Protein C Receptor) CLIA Kit
E-CL-R0606	Rat SELE (E-Selectin) CLIA Kit
E-CL-R0607	Rat sFAS/sAPO-1 (Soluble Factor-Related Apoptosis) CLIA Kit
E-CL-R0608	Rat IL-2SR (Soluble Interleukin 2 Receptor) CLIA Kit
E-CL-R0609	Rat IL-6R (Interleukin 6 Receptor) CLIA Kit
E-CL-R0610	Rat IL-1SR I (Soluble Interleukin-1 Receptor I) CLIA Kit
E-CL-R0611	Rat IL-1SR II (Soluble Interleukin-1 Receptor II) CLIA Kit
E-CL-R0612	Rat IL-2SR α (Soluble Interleukin-2 Receptor) CLIA Kit
E-CL-R0613	Rat S-100 (Soluble Protein-100) CLIA Kit
E-CL-R0614	Rat TLR-2 (Toll-Like Receptor 2) CLIA Kit
E-CL-R0615	Rat SVCAM-1 (Soluble Vascular Cell Adhesion Molecule 1) CLIA Kit
E-CL-R0616	Rat FLT1/VEGFR1 (Vascular Endothelial Growth Factor Receptor 1) CLIA Kit
E-CL-R0617	Rat SS (Somatostatin) CLIA Kit
E-CL-R0618	Rat SCF (Stem Cell Factor) CLIA Kit
E-CL-R0619	Rat SDF-1 (Stromal Cell Derived Factor 1) CLIA Kit

Cat.NO	Product Name
E-CL-R0620	Rat SDF-1 β (Stromal Cell Derived Factor 1 β) CLIA Kit
E-CL-R0621	Rat SOCS3 (Suppressors Of Cytokine Signaling 3) CLIA Kit
E-CL-R0622	Rat SYP (Synaptophysin) CLIA Kit
E-CL-R0623	Rat SNAP-25 (Synaptosome Associated Protein 25) CLIA Kit
E-CL-R0624	Rat SDC4 (Syndecan 4) CLIA Kit
E-CL-R0625	Rat TACR1 (Tachykinin Receptor 1) CLIA Kit
E-CL-R0626	Rat TACR2 (Tachykinin Receptor 2) CLIA Kit
E-CL-R0627	Rat THP (Tamm–Horsfall Glycoprotein) CLIA Kit
E-CL-R0628	Rat TRACP-5b (Tartrate-Resistant Acid Phosphatase 5b) CLIA Kit
E-CL-R0629	Rat MAP τ (Microtubule Associated Protein Tau) CLIA Kit
E-CL-R0630	Rat TEK (Tyrosine Kinase, Endothelial) CLIA Kit
E-CL-R0631	Rat TE (Telomerase) CLIA Kit
E-CL-R0632	Rat TNC (Tenascin C) CLIA Kit
E-CL-R0633	Rat TDGF1 (Teratocarcinoma Derived Growth Factor 1) CLIA Kit
E-CL-R0634	Rat TrxR1 (Thioredoxin Reductase 1) CLIA Kit
E-CL-R0635	Rat TRX (Thioredoxin) CLIA Kit
E-CL-R0636	Rat TAFI (Thrombin Activatable Fibrinolysis Inhibitor) CLIA Kit
E-CL-R0637	Rat TAT (Thrombin-Antithrombin Complex) CLIA Kit
E-CL-R0638	Rat TM (Thrombomodulin) CLIA Kit

Cat.NO	Product Name
E-CL-R0639	Rat TPO (Thrombopoietin) CLIA Kit
E-CL-R0640	Rat TSP-1 (Thrombospondin 1) CLIA Kit
E-CL-R0641	Rat TXB2 (Thromboxane B2) CLIA Kit
E-CL-R0642	Rat TSLP (Thymic Stromal Lymphopoietin) CLIA Kit
E-CL-R0643	Rat TARC (Thymus Activation Regulated Chemokine) CLIA Kit
E-CL-R0644	Rat TECK (Thymus Expressed Chemokine) CLIA Kit
E-CL-R0645	Rat TG (Thyroglobulin) CLIA Kit
E-CL-R0646	Rat GPC3 (Glypican 3) CLIA Kit
E-CL-R0647	Rat TSH (Thyroid Stimulating Hormone) CLIA Kit
E-CL-R0648	Rat TRH (Thyrotropin-Releasing Hormone) CLIA Kit
E-CL-R0649	Rat TAB (Thyroxine Antibody) CLIA Kit
E-CL-R0650	Rat T4 (Thyroxine) CLIA Kit
E-CL-R0651	Rat TFPI (Tissue Factor Pathway Inhibitor) CLIA Kit
E-CL-R0652	Rat TF (Tissue Factor) CLIA Kit
E-CL-R0653	Rat TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) CLIA Kit
E-CL-R0654	Rat TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) CLIA Kit
E-CL-R0655	Rat TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) CLIA Kit
E-CL-R0656	Rat TPA (Tissue Polypeptide Antigen) CLIA Kit
E-CL-R0657	Rat TLR4 (Toll-Like Receptor 4) CLIA Kit

Cat.NO	Product Name
E-CL-R0658	Rat TLR9 (Toll-Like Receptor 9) CLIA Kit
E-CL-R0659	Rat 5-LO (Arachidonate 5-Lipoxygenase) CLIA Kit
E-CL-R0660	Rat TFR2 (Transferrin Receptor Protein 2) CLIA Kit
E-CL-R0661	Rat TFR (Transferrin Receptor) CLIA Kit
E-CL-R0662	Rat TGF- α (Transforming Growth Factor Alpha) CLIA Kit
E-CL-R0663	Rat GP4/CD36 (Platelet Membrane Glycoprotein IV) CLIA Kit
E-CL-R0664	Rat TGF- β 2 (Transforming Growth Factor Beta 2) CLIA Kit
E-CL-R0665	Rat TGF- β 3 (Transforming Growth Factor Beta 3) CLIA Kit
E-CL-R0666	Rat TGF β I (Transforming Growth Factor Beta Induced) CLIA Kit
E-CL-R0667	Rat TTR (Transthyretin) CLIA Kit
E-CL-R0668	Rat TFF2 (Trefoil Factor 2) CLIA Kit
E-CL-R0669	Rat TFF3 (Trefoil Factor 3, Intestinal) CLIA Kit
E-CL-R0670	Rat TREM-1 (Triggering Receptor Expresses on Myeloid Cells-1) CLIA Kit
E-CL-R0671	Rat INHB (Inhibin B) CLIA Kit
E-CL-R0672	Rat TRY (Trypsin) CLIA Kit
E-CL-R0673	Rat TAP (Trypsinogen Activation Peptide) CLIA Kit
E-CL-R0674	Rat TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) CLIA Kit
E-CL-R0675	Rat TNF- β (Tumor Necrosis Factor β) CLIA Kit
E-CL-R0676	Rat TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) CLIA Kit

Cat.NO	Product Name
E-CL-R0677	Rat TP53 (Tumor Protein 53) CLIA Kit
E-CL-R0678	Rat TSGF (Tumor Specific Growth Factor/Tumor Supplied Group of Factor) CLIA Kit
E-CL-R0679	Rat TYR (Tyrosinase) CLIA Kit
E-CL-R0680	Rat TIE-2 (Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2) CLIA Kit
E-CL-R0681	Rat U-T3 (Ultrasensitivity Triiodothyronine) CLIA Kit
E-CL-R0682	Rat PLAU/uPA(Urokinase-Type Plasminogen Activator) CLIA Kit
E-CL-R0683	Rat PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) CLIA Kit
E-CL-R0684	Rat VEGF-B (Vascular Endothelial Cell Growth Factor B) CLIA Kit
E-CL-R0685	Rat VEGF-C (Vascular Endothelial Cell Growth Factor C) CLIA Kit
E-CL-R0686	Rat VEGF-D (Vascular Endothelial Cell Growth Factor D) CLIA Kit
E-CL-R0687	Rat VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) CLIA Kit
E-CL-R0688	Rat VCAM-1 (Vascular Cell Adhesion Molecule 1) CLIA Kit
E-CL-R0689	Rat VIP (Vasoactive Intestinal Peptide) CLIA Kit
E-CL-R0690	Rat VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) CLIA Kit
E-CL-R0691	Rat VF (Visfatin) CLIA Kit
E-CL-R0692	Rat DBP (Vitamin D Binding Protein) CLIA Kit
E-CL-R0693	Rat ADAMTS13/VWFPC (Von Willebrand Factor Cleaving Protease) CLIA Kit
E-CL-R0694	Rat VWF (Von Willebrand Factor) CLIA Kit
E-CL-R0695	Rat α 1-AGP (α 1-Acid Glycoprotein) CLIA Kit

Cat.NO	Product Name
E-CL-R0696	Rat α 1-MG (α 1-Microglobulin) CLIA Kit
E-CL-R0697	Rat α -Glu (Alpha Glucosidase) CLIA Kit
E-CL-R0698	Rat β 2-GP1 IgA/G/M(β 2-Glycoprotein 1 Antibody IgA/G/M) CLIA Kit
E-CL-R0699	Rat BMG/ β 2-MG (Beta-2-Microglobulin) CLIA Kit
E-CL-R0700	Rat β -TG (β -Thromboglobulin) CLIA Kit
E-CL-R0701	Rat pr (phospho Tau Protein) CLIA Kit
E-CL-R0702	Rat TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) CLIA Kit
E-CL-R0703	Rat T3 (Triiodothyronine) CLIA Kit
E-CL-R0704	Rat MDH1 (Malate Dehydrogenase 1) CLIA Kit
E-CL-R0705	Rat LEPR (Leptin Receptor) CLIA Kit
E-CL-R0706	Rat PCIII (Procollagen III) CLIA Kit
E-CL-R0707	Rat TRF (Transferrin) CLIA Kit
E-CL-R0708	Rat VLDL (Very Low Density Lipoprotein) CLIA Kit
E-CL-R0709	Rat SDH (Succinate Dehydrogenase) CLIA Kit
E-CL-R0710	Rat ApoA1 (Apolipoprotein A1) CLIA Kit
E-CL-R0711	Rat ICD (Isocitrate Dehydrogenase) CLIA Kit
E-CL-R0712	Rat α -KGDHC (α Ketoglutarate Dehydrogenase) CLIA Kit
E-CL-R0713	Rat ND1 (NADH Dehydrogenase 1) CLIA Kit
E-CL-R0714	Rat 6-PFK (Phosphofructokinase) CLIA Kit

Cat.NO	Product Name
E-CL-R0715	Rat ApoA2 (Apolipoprotein A2) CLIA Kit
E-CL-R0716	Rat 8-epi-PGF2 α (8-Epi-Prostaglandin F2 Alpha) CLIA Kit
E-CL-R0717	Rat ApoB (Apolipoprotein B) CLIA Kit
E-CL-R0718	Rat ApoE (Apolipoprotein E) CLIA Kit
E-CL-R0719	Rat ALT (Alanine Transaminase) CLIA Kit
E-CL-R0720	Rat BDNF (Brain Derived Neurotrophic Factor) CLIA Kit
E-CL-R0721	Rat cTn-I/TNNI3 (Cardiac Troponin I) CLIA Kit
E-CL-R0722	Rat CK-MB (Creatine Kinase MB Isoenzyme) CLIA Kit
E-CL-R0723	Rat A β 40 (Amyloid Beta 40) CLIA Kit
E-CL-R0724	Rat A β 42 (Amyloid Beta 42) CLIA Kit
E-CL-R0725	Rat CYP7A1 (Cholesterol 7-Alpha-Hydroxylase) CLIA Kit
E-CL-R0726	Rat α -CTx (Alpha Crosslaps) CLIA Kit
E-CL-R0727	Rat β -CTx (Beta Crosslaps) CLIA Kit
E-CL-R0728	Rat IL-2R β (Interleukin-2 Receptor- β) CLIA Kit
E-CL-R0729	Rat PGA (Pepsinogen A) CLIA Kit
E-CL-R0730	Rat PGC (PepsinogenC) CLIA Kit
E-CL-R0731	Rat NT-Pro-ANP (N-Terminal Pro-Atrial Natriuretic Peptide) CLIA Kit
E-CL-R0732	Rat PINP (Procollagen Type I N-Terminal Propeptide) CLIA Kit
E-CL-R0733	Rat Rt3 (Reverse Triiodothyronine) CLIA Kit

Cat.NO	Product Name
E-CL-R0734	Rat TGM2 (Transglutaminase 2, Tissue) CLIA Kit
E-CL-R0735	Rat 11-DH-TXB2 (11-Dehydrothromboxane B2) CLIA Kit
E-CL-R0736	Rat SOD1(Superoxide Dismutase 1, Soluble) CLIA Kit
E-CL-R0737	Rat SOD2 (Superoxide Dismutase 2, Mitochondrial) CLIA Kit
E-CL-R0738	Rat SOD3 (Superoxide Dismutase 3, Extracellular) CLIA Kit
E-CL-R0739	Rat SOD4 (Superoxide Dismutase 4, Copper Chaperone) CLIA Kit
E-CL-R0740	Rat GFAP (Glial Fibrillary Acidic Protein) CLIA Kit
E-CL-R0741	Rat ANG I (Angiotensin I) CLIA Kit
E-CL-R0742	Rat ANG II (Angiotensin II) CLIA Kit
E-CL-R0743	Rat PS1 (Presenilin 1) CLIA Kit
E-CL-R0744	Rat AgRP (Agouti Related Protein) CLIA Kit
E-CL-R0745	Rat TH (Tyrosine Hydroxylase) CLIA Kit
E-CL-R0746	Rat NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) CLIA Kit
E-CL-R0747	Rat CTXI (Cross Linked C-telopeptide of Type I Collagen) CLIA Kit
E-CL-R0748	Rat RTN4 (Reticulon 4) CLIA Kit
E-CL-R0749	Rat FASL (Factor Related Apoptosis Ligand) CLIA Kit
E-CL-M0001	Mouse ACV-A (Activin A) CLIA Kit
E-CL-M0002	Mouse ADP/Acrp30 (Adiponectin) CLIA Kit
E-CL-M0003	Mouse ANG3 (Angiopoietin 3) CLIA Kit

Cat.NO	Product Name
E-CL-M0004	Mouse BAFF/CD257 (B-cell Activating Factor) CLIA Kit
E-CL-M0005	Mouse LOX (Lysyl Oxidase) CLIA Kit
E-CL-M0006	Mouse MCP-1 (Monocyte Chemotactic Protein 1) CLIA Kit
E-CL-M0007	Mouse MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) CLIA Kit
E-CL-M0008	Mouse MIP-1 β (Macrophage Inflammatory Protein 1 Beta)CLIA Kit
E-CL-M0009	Mouse RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) CLIA Kit
E-CL-M0010	Mouse TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) CLIA Kit
E-CL-M0011	Mouse TARC (Thymus Activation Regulated Chemokine) CLIA Kit
E-CL-M0012	Mouse MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) CLIA Kit
E-CL-M0013	Mouse MDC/CCL22 (Macrophage-Derived Chemokine) CLIA Kit
E-CL-M0014	Mouse CHE (Chemerin) CLIA Kit
E-CL-M0015	Mouse CHI3L1 (Chitinase 3-like 1) CLIA Kit
E-CL-M0016	Mouse IL-1R1 (Interleukin 1 Receptor Type I) CLIA Kit
E-CL-M0017	Mouse GRO α /CXCL1 (Growth Regulated Oncogene Alpha) CLIA Kit
E-CL-M0018	Mouse GRO β (Growth Regulated Oncogene Beta) CLIA Kit
E-CL-M0019	Mouse CXCL9/MIG(Monokine induced by interferon-gamma)CLIA Kit
E-CL-M0020	Mouse IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) CLIA Kit
E-CL-M0021	Mouse Cyt C (Cytochrome C) CLIA Kit
E-CL-M0022	Mouse DKK1 (Dickkopf Related Protein 1)CLIA Kit

Cat.NO	Product Name
E-CL-M0023	Mouse EGF (Epidermal growth factor) CLIA Kit
E-CL-M0024	Mouse ESM1 (Endothelial Cell Specific Molecule 1) CLIA Kit
E-CL-M0025	Mouse EPO (Erythropoietin) CLIA Kit
E-CL-M0026	Mouse FASL (Factor Related Apoptosis Ligand) CLIA Kit
E-CL-M0027	Mouse FGF21(Fibroblast Growth Factor 21) CLIA Kit
E-CL-M0028	Mouse Flt3L (FMS Like Tyrosine Kinase 3 Ligand) CLIA Kit
E-CL-M0029	Mouse G-CSF (Granulocyte Colony Stimulating Factor) CLIA Kit
E-CL-M0030	Mouse GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) CLIA Kit
E-CL-M0031	Mouse IFN- β (Interferon Beta) CLIA Kit
E-CL-M0032	Mouse IGF-1 (Insulin-like growth factor 1) CLIA Kit
E-CL-M0033	Mouse IGFBP-3 (Insulin-like Growth Factor Binding Protein 3) CLIA Kit
E-CL-M0034	Mouse IL-1 α (Interleukin 1 Alpha) CLIA Kit
E-CL-M0035	Mouse IL-1 β (Interleukin 1 Beta) CLIA Kit
E-CL-M0036	Mouse IL-1ra (Interleukin 1 Receptor Antagonist) CLIA Kit
E-CL-M0037	Mouse LEP (Leptin) CLIA Kit
E-CL-M0038	Mouse LIF (Leukemia Inhibitory Factor) CLIA Kit
E-CL-M0039	Mouse M-CSF (Macrophage Colony-Stimulating Factor) CLIA Kit
E-CL-M0040	Mouse IL-2 (Interleukin 2) CLIA Kit
E-CL-M0041	Mouse IL-4 (Interleukin 4) CLIA Kit

Cat.NO	Product Name
E-CL-M0042	Mouse IL-6 (Interleukin 6) CLIA Kit
E-CL-M0043	Mouse IL-8 (Interleukin 8) CLIA Kit
E-CL-M0044	Mouse IL-10 (Interleukin 10) CLIA Kit
E-CL-M0045	Mouse IL-17 (Interleukin 17) CLIA Kit
E-CL-M0046	Mouse IFN- γ (Interferon Gamma) CLIA Kit
E-CL-M0047	Mouse TNF- α (Tumor Necrosis Factor Alpha) CLIA Kit
E-CL-M0048	Mouse VEGF (Vascular Endothelial cell Growth Factor) CLIA Kit
E-CL-M0049	Mouse TGF- β 1 (Transforming Growth Factor β 1) CLIA Kit
E-CL-M0050	Mouse CRP (C-Reactive Protein) CLIA Kit
E-CL-M0051	Mouse INS (Insulin) CLIA Kit
E-CL-M0052	Mouse I-TAC (Interferon Inducible T-Cell Alpha Chemoattractant) CLIA Kit
E-CL-M0053	Mouse LH (Luteinizing Hormone) CLIA Kit
E-CL-M0054	Mouse GH (Growth Hormone) CLIA Kit
E-CL-M0055	Mouse REN (Renin) CLIA Kit
E-CL-M0056	Mouse CASP8 (Caspase 8) CLIA Kit
E-CL-M0057	Mouse T (Testosterone) CLIA Kit
E-CL-M0058	Mouse OxLDL (Oxidized Low-Density Lipoprotein) CLIA Kit
E-CL-M0059	Mouse A β 40 (Amyloid Beta 40) CLIA Kit
E-CL-M0060	Mouse A β 42 (Amyloid Beta 42) CLIA Kit

Cat.NO	Product Name
E-CL-M0061	Mouse TPS (Tryptase) CLIA Kit
E-CL-M0062	Mouse MMP-10 (Matrix Metalloproteinase 10) CLIA Kit
E-CL-M0063	Mouse MMP-12 (Matrix Metalloproteinase 12) CLIA Kit
E-CL-M0064	Mouse MMP-13 (Matrix Metalloproteinase 13) CLIA Kit
E-CL-M0065	Mouse sICAM-1/CD54 (Soluble Intercellular Adhesion Molecule 1) CLIA Kit
E-CL-M0066	Mouse ANG (Angiogenin) CLIA Kit
E-CL-M0067	Mouse ACTH (Adrenocorticotrophic Hormone) CLIA Kit
E-CL-M0068	Mouse DA (Dopamine) CLIA Kit
E-CL-M0069	Mouse OPG (Osteoprotegerin) CLIA Kit
E-CL-M0070	Mouse PRL (Prolactin) CLIA Kit
E-CL-M0071	Mouse MYO (Myoglobin) CLIA Kit
E-CL-M0072	Mouse Tn-T (Troponin T) CLIA Kit
E-CL-M0073	Mouse Tn-I (Troponin I) CLIA Kit
E-CL-M0074	Mouse PAPPA (Pregnancy Associated Plasma Protein A) CLIA Kit
E-CL-M0075	Mouse ALP (Alkaline Phosphatase) CLIA Kit
E-CL-M0076	Mouse NSE (Neuron-Specific Enolase) CLIA Kit
E-CL-M0077	Mouse GLP-1 (Glucagon Like Peptide 1) CLIA Kit
E-CL-M0078	Mouse ACE (Angiotensin I Converting Enzyme) CLIA Kit
E-CL-M0079	Mouse ANGPTL2 (Angiopoietin Like Protein 2) CLIA Kit

Cat.NO	Product Name
E-CL-M0080	Mouse ANGPTL4 (Angiopoietin Like Protein 4) CLIA Kit
E-CL-M0081	Mouse AGT (Angiotensinogen) CLIA Kit
E-CL-M0082	Mouse ANG1 (Angiopoietin 1) CLIA Kit
E-CL-M0083	Mouse ANG2 (Angiopoietin 2) CLIA Kit
E-CL-M0084	Mouse ANXA5 (Annexin A5) CLIA Kit
E-CL-M0085	Mouse ANXA2(Annexin A2) CLIA Kit
E-CL-M0086	Mouse ADH (Antidiuretic Hormone) CLIA Kit
E-CL-M0087	Mouse NRG-1 (Neuregulin 1) CLIA Kit
E-CL-M0088	Mouse AECA (anti-Endothelial Cell Antibody) CLIA Kit
E-CL-M0089	Mouse anti-IgE receptor antibody CLIA Kit
E-CL-M0090	Mouse AMH (Anti-Mullerian Hormone) CLIA Kit
E-CL-M0091	Mouse ANA (anti-Nuclear Antibody) CLIA Kit
E-CL-M0092	Mouse CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) CLIA Kit
E-CL-M0093	Mouse HDAC2 (Histone Deacetylase 2) CLIA Kit
E-CL-M0094	Mouse NRG-4 (Neuregulin 4) CLIA Kit
E-CL-M0095	Mouse AT- III (anti-Thrombin III) CLIA Kit
E-CL-M0096	Mouse APO-A1 (Apolipoprotein A1) CLIA Kit
E-CL-M0097	Mouse APO-A2 (Apolipoprotein A2) CLIA Kit
E-CL-M0098	Mouse ApoB (Apolipoprotein B) CLIA Kit

Cat.NO	Product Name
E-CL-M0099	Mouse APO-C1 (Apolipoprotein C1) CLIA Kit
E-CL-M0100	Mouse APO-C4 (Apolipoprotein C4) CLIA Kit
E-CL-M0101	Mouse APO-E (Apolipoprotein E) CLIA Kit
E-CL-M0102	Mouse APO-F (Apolipoprotein F) CLIA Kit
E-CL-M0103	Mouse ApoB100 (Apolipoprotein B100) CLIA Kit
E-CL-M0104	Mouse AIF (Apoptosis Inducing Factor) CLIA Kit
E-CL-M0105	Mouse ASK1 (Apoptosis Signal Regulating Kinase 1) CLIA Kit
E-CL-M0106	Mouse APAF1 (Apoptotic Peptidase Activating Factor 1) CLIA Kit
E-CL-M0107	Mouse SLC (Secondary Lymphoid Tissue Chemokine) CLIA Kit
E-CL-M0108	Mouse AQP-1 (Aquaporin 1) CLIA Kit
E-CL-M0109	Mouse GROγ (Growth Regulated Oncogene Gamma) CLIA Kit
E-CL-M0110	Mouse AQP-2 (Aquaporin 2) CLIA Kit
E-CL-M0111	Mouse AGER (Advanced Glycosylation End Product Specific Receptor) CLIA Kit
E-CL-M0112	Mouse AQP-3 (Aquaporin 3) CLIA Kit
E-CL-M0113	Mouse AQP-4 (Aquaporin 4) CLIA Kit
E-CL-M0114	Mouse AQP-5 (Aquaporin 5) CLIA Kit
E-CL-M0115	Mouse AQP-9 (Aquaporin 9) CLIA Kit
E-CL-M0116	Mouse Arg (Arginase) CLIA Kit
E-CL-M0117	Mouse AhR (Aryl Hydrocarbon Receptor) CLIA Kit

Cat.NO	Product Name
E-CL-M0118	Mouse ARNT2 (Aryl Hydrocarbon Receptor Nuclear Translocator 2) CLIA Kit
E-CL-M0119	Mouse ASGR1 (Asialoglycoprotein Receptor 1) CLIA Kit
E-CL-M0120	Mouse AST (Aspartate amino-transferase) CLIA Kit
E-CL-M0121	Mouse ATM (Ataxia Telangiectasia Mutated) CLIA Kit
E-CL-M0122	Mouse ATXN1 (Ataxin 1) CLIA Kit
E-CL-M0123	Mouse ABCA1 (ATP Binding Cassette Transporter A1) CLIA Kit
E-CL-M0124	Mouse ANP (Atrial Natriuretic Peptide) CLIA Kit
E-CL-M0125	Mouse AXIN2 (Axis Inhibition Protein 2) CLIA Kit
E-CL-M0126	Mouse bFGF/FGF2 (Basic Fibroblast Growth Factor) CLIA Kit
E-CL-M0127	Mouse FGF4 (Fibroblast Growth Factor 4) CLIA Kit
E-CL-M0128	Mouse FGF6 (Fibroblast Growth Factor 6) CLIA Kit
E-CL-M0129	Mouse BAFFR (B-Cell Activation Factor Receptor) CLIA Kit
E-CL-M0130	Mouse Bcl-2 (B-Cell Leukemia/Lymphoma 2) CLIA Kit
E-CL-M0131	Mouse Bax (Bcl-2 Associated X Protein) CLIA Kit
E-CL-M0132	Mouse BMF (Bcl-2 Modifying Factor) CLIA Kit
E-CL-M0133	Mouse BTC (Beta Cellulin) CLIA Kit
E-CL-M0134	Mouse BACE1 (Beta Site APP Cleaving Enzyme 1) CLIA Kit
E-CL-M0135	Mouse β -EPR (Beta-Endorphin Receptor) CLIA Kit
E-CL-M0136	Mouse β -EP (Beta-Endorphin) CLIA Kit

Cat.NO	Product Name
E-CL-M0137	Mouse Bid (BH3 Interacting Domain Death Agonist) CLIA Kit
E-CL-M0138	Mouse IL-25 (Interleukin 25) CLIA Kit
E-CL-M0139	Mouse BLC (B-Lymphocyte Chemoattractant) CLIA Kit
E-CL-M0140	Mouse BMP-1 (Bone Morphogenetic Protein 1) CLIA Kit
E-CL-M0141	Mouse BMP-2 (Bone Morphogenetic Protein 2) CLIA Kit
E-CL-M0142	Mouse BMP-4 (Bone Morphogenetic Protein 4) CLIA Kit
E-CL-M0143	Mouse BMP-6 (Bone Morphogenetic Protein 6) CLIA Kit
E-CL-M0144	Mouse BMP-7 (Bone Morphogenetic Protein 7) CLIA Kit
E-CL-M0145	Mouse BMP8B (Bone Morphogenetic Protein 8B) CLIA Kit
E-CL-M0146	Mouse BALP (Bone Alkaline Phosphatase) CLIA Kit
E-CL-M0147	Mouse CASP1 (Caspase 1) CLIA Kit
E-CL-M0148	Mouse BK (Bradykinin) CLIA Kit
E-CL-M0149	Mouse BDNF (Brain Derived Neurotrophic Factor) CLIA Kit
E-CL-M0150	Mouse BNP (Brain Natriuretic Peptide) CLIA Kit
E-CL-M0151	Mouse BRAK (Breast And Kidney Expressed Chemokine) CLIA Kit
E-CL-M0152	Mouse Ca-ATPase (Calcium ATPase) CLIA Kit
E-CL-M0153	Mouse VE-Cadherin (Vascular Endothelial Cadherin) CLIA Kit
E-CL-M0154	Mouse E-Cad (E-Cadherin) CLIA Kit
E-CL-M0155	Mouse NCAD (Neural Cadherin) CLIA Kit

Cat.NO	Product Name
E-CL-M0156	Mouse CaN (Calcineurin) CLIA Kit
E-CL-M0157	Mouse CGRP (Calcitonin Gene Related Peptide) CLIA Kit
E-CL-M0158	Mouse CT (Calcitonin) CLIA Kit
E-CL-M0159	Mouse CACYBP (Calcyclin Binding Protein) CLIA Kit
E-CL-M0160	Mouse CAST (Calpastatin) CLIA Kit
E-CL-M0161	Mouse CRT (Calreticulin) CLIA Kit
E-CL-M0162	Mouse CALU (Calumenin) CLIA Kit
E-CL-M0163	Mouse CA9 (Carbonic Anhydrase IX) CLIA Kit
E-CL-M0164	Mouse CA1 (Carbonic Anhydrase I) CLIA Kit
E-CL-M0165	Mouse CA2 (Carbonic Anhydrase 2) CLIA Kit
E-CL-M0166	Mouse PICP (Procollagen I C-terminal Propeptide) CLIA Kit
E-CL-M0167	Mouse CEA (Carcinoembryonic Antign) CLIA Kit
E-CL-M0168	Mouse PINP (Procollagen Type I N-Terminal Propeptide) CLIA Kit
E-CL-M0169	Mouse CLS (Cardiolipin Synthase) CLIA Kit
E-CL-M0170	Mouse CLCF1 (Cardiotrophin Like Cytokine Factor 1) CLIA Kit
E-CL-M0171	Mouse CT-1 (Cardiotrophln 1) CLIA Kit
E-CL-M0172	Mouse COMP (Cartilage Oligomeric Matrix Protein) CLIA Kit
E-CL-M0173	Mouse CASP3 (Caspase 3) CLIA Kit
E-CL-M0174	Mouse CASP7 (Caspase 7) CLIA Kit

Cat.NO	Product Name
E-CL-M0175	Mouse CASP9 (Caspase 9) CLIA Kit
E-CL-M0176	Mouse CASP10 (Caspase 10) CLIA Kit
E-CL-M0177	Mouse CASP11 (Caspase 11) CLIA Kit
E-CL-M0178	Mouse CASP14 (Caspase 14) CLIA Kit
E-CL-M0179	Mouse CTNNb1 (Catenin, Beta 1) CLIA Kit
E-CL-M0180	Mouse CTSD (Cathepsin D) CLIA Kit
E-CL-M0181	Mouse CTSK (Cathepsin K) CLIA Kit
E-CL-M0182	Mouse CTSL (Cathepsin L) CLIA Kit
E-CL-M0183	Mouse CAV1 (Caveolin 1) CLIA kit
E-CL-M0184	Mouse CP (Ceruloplasmin) CLIA Kit
E-CL-M0185	Mouse CCL1 (Chemokine C-C-Motif Ligand 1) CLIA Kit
E-CL-M0186	Mouse CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) CLIA Kit
E-CL-M0187	Mouse CCRL1 (Chemokine C-C-Motif Receptor Like Protein 1) CLIA Kit
E-CL-M0188	Mouse CX3CL1 (Chemokine C-X3-C-Motif Ligand 1) CLIA Kit
E-CL-M0189	Mouse CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) CLIA Kit
E-CL-M0190	Mouse CXCL16 (Chemokine C-X-C-Motif Ligand 16) CLIA Kit
E-CL-M0191	Mouse CCKAR (Cholecystokinin A Receptor) CLIA Kit
E-CL-M0192	Mouse CETP (cholest erolester transfer protein) CLIA Kit
E-CL-M0193	Mouse Gas6 (Growth Arrest Specific Protein 6) CLIA Kit

Cat.NO	Product Name
E-CL-M0194	Mouse CS (chondroitin sulfate) CLIA Kit
E-CL-M0195	Mouse CNTF (Ciliary Neurotrophic Factor) CLIA Kit
E-CL-M0196	Mouse CC17 (Clara Cell Protein) CLIA Kit
E-CL-M0197	Mouse CD3 (Cluster of differentiation 3) CLIA Kit
E-CL-M0198	Mouse CD4 (Cluster Of Differentiation 4) CLIA Kit
E-CL-M0199	Mouse CD19 (Cluster of differentiation 19) CLIA Kit
E-CL-M0200	Mouse CD30L (Cluster of Differentiation 30 Ligand) CLIA Kit
E-CL-M0201	Mouse CD40L (Cluster of Differentiation 40 Ligand) CLIA Kit
E-CL-M0202	Mouse CLU (Clusterin) CLIA Kit
E-CL-M0203	Mouse F II (coagulation factor II) CLIA Kit
E-CL-M0204	Mouse GP4/CD36 (Platelet Membrane Glycoprotein IV) CLIA Kit
E-CL-M0205	Mouse FV (Coagulation Factor V) CLIA Kit
E-CL-M0206	Mouse FVIII (Coagulation Factor VIII) CLIA Kit
E-CL-M0207	Mouse FIX (Coagulation Factor IX) CLIA Kit
E-CL-M0208	Mouse FX (coagulation factor X) CLIA Kit
E-CL-M0209	Mouse FXII (Coagulation Factor XII) CLIA Kit
E-CL-M0210	Mouse F XIII B (Coagulation Factor XIII B Polypeptide) CLIA Kit
E-CL-M0211	Mouse CFL1 (Cofilin 1, Non-Muscle) CLIA Kit
E-CL-M0212	Mouse COL1 (Collagen Type I) CLIA Kit

Cat.NO	Product Name
E-CL-M0213	Mouse COL3 (Collagen Type III) CLIA Kit
E-CL-M0214	Mouse COL4 (Collagen Type IV) CLIA Kit
E-CL-M0215	Mouse COL1 α 2 (Collagen Type I Alpha 2) CLIA Kit
E-CL-M0216	Mouse COL4 α 1 (Collagen Type IV Alpha 1) CLIA Kit
E-CL-M0217	Mouse COL10 (Collagen Type X) CLIA Kit
E-CL-M0218	Mouse COL1 (Collagen Type I) CLIA Kit
E-CL-M0219	Mouse GPC3 (Glypican 3) CLIA Kit
E-CL-M0220	Mouse CTGF (Connective Tissue Growth Factor) CLIA Kit
E-CL-M0221	Mouse CX26 (Connexin 26) CLIA Kit
E-CL-M0222	Mouse CX40 (Connexin 40) CLIA Kit
E-CL-M0223	Mouse CX43 (Connexin 43) CLIA Kit
E-CL-M0224	Mouse CNTN1 (Contactin 1) CLIA Kit
E-CL-M0225	Mouse CORT (Corticosterone) CLIA Kit
E-CL-M0226	Mouse CRH (Corticotropin Releasing Hormone) CLIA Kit
E-CL-M0227	Mouse C-P (C-Peptide) CLIA Kit
E-CL-M0228	Mouse CK-MB (Creatine Kinase MB Isoenzyme) CLIA Kit
E-CL-M0229	Mouse CKMT1A (Creatine Kinase, Mitochondrial 1A) CLIA Kit
E-CL-M0230	Mouse CKMT1B (Creatine Kinase, Mitochondrial 1B) CLIA Kit
E-CL-M0231	Mouse CKM (Creatine Kinase, Muscle) CLIA Kit

Cat.NO	Product Name
E-CL-M0232	Mouse NTXI (Cross Linked N-Telopeptide of Type I Collagen) CLIA Kit
E-CL-M0233	Mouse α -CTx (Alpha Crosslaps) CLIA Kit
E-CL-M0234	Mouse CICP (cross-linked C-telopeptides of Type I collagen) CLIA Kit
E-CL-M0235	Mouse CRYaB (Crystallin Alpha B) CLIA Kit
E-CL-M0236	Mouse CTXI (Cross Linked C-telopeptide of Type I Collagen) CLIA Kit
E-CL-M0237	Mouse PIIICP (c-terminal Procollagen III propeptide) CLIA Kit
E-CL-M0238	Mouse CNP (C-type natriuretic peptide) CLIA Kit
E-CL-M0239	Mouse CTACK (Cutaneous T-Cell Attracting Chemokine) CLIA Kit
E-CL-M0240	Mouse β -CTx (Beta Crosslaps) CLIA Kit
E-CL-M0241	Mouse CXCR3 (CXC-chemokine receptor 3) CLIA Kit
E-CL-M0242	Mouse COL2 (Collagen Type II) CLIA Kit
E-CL-M0243	Mouse CDKN1B (Cyclin Dependent Kinase Inhibitor 1B) CLIA Kit
E-CL-M0244	Mouse CCND1 (Cyclin-D1) CLIA Kit
E-CL-M0245	Mouse CCND2 (Cyclin-D2) CLIA Kit
E-CL-M0246	Mouse CYPA (cyclophilin A) CLIA Kit
E-CL-M0247	Mouse CYPB (Cyclophilin B) CLIA Kit
E-CL-M0248	Mouse Cys-C (Cystatin C) CLIA Kit
E-CL-M0249	Mouse CYB561 (Cytochrome b-561) CLIA Kit
E-CL-M0250	Mouse Cyt-C (Cytochrome C) CLIA Kit

Cat.NO	Product Name
E-CL-M0251	Mouse CYP7A1 (Cytochrome P450 7A1) CLIA Kit
E-CL-M0252	Mouse CTLA4 (Cytotoxic T-Lymphocyte Associated Antigen 4) CLIA Kit
E-CL-M0253	Mouse DCN (Decorin) CLIA Kit
E-CL-M0254	Mouse DEFA5 (Defensin Alpha 5, Paneth Cell Specific) CLIA Kit
E-CL-M0255	Mouse DHEA (Dehydroepiandrosterone) CLIA Kit
E-CL-M0256	Mouse DPD (Deoxypyridinoline) CLIA Kit
E-CL-M0257	Mouse DNASE1 (Deoxyribonuclease I) CLIA Kit
E-CL-M0258	Mouse Des (Desmin) CLIA Kit
E-CL-M0259	Mouse DAO (Diamine Oxidase) CLIA Kit
E-CL-M0260	Mouse DKK2 (Dickkopf Related Protein 2) CLIA Kit
E-CL-M0261	Mouse DLD (Dihydrolipoyl Dehydrogenase) CLIA Kit
E-CL-M0262	Mouse DLAT (Dihydrolipoyl Transacetylase) CLIA Kit
E-CL-M0263	Mouse D-LDH (D-Lactate Dehydrogenase) CLIA Kit
E-CL-M0264	Mouse D1R (Dopamine D1 receptor) CLIA Kit
E-CL-M0265	Mouse DAT (Dopamine Transporter) CLIA Kit
E-CL-M0266	Mouse DNM1 (Dynamin 1) CLIA Kit
E-CL-M0267	Mouse DMD (Dystrophin) CLIA Kit
E-CL-M0268	Mouse E2 (Estradiol) CLIA Kit
E-CL-M0269	Mouse EGR1 (Early Growth Response Protein 1) CLIA Kit

Cat.NO	Product Name
E-CL-M0270	Mouse EGR2 (Early Growth Response Protein 2) CLIA Kit
E-CL-M0271	Mouse NE (Neutrophil Elastase) CLIA Kit
E-CL-M0272	Mouse ELA4 (Elastase 4) CLIA Kit
E-CL-M0273	Mouse ELN (Elastin) CLIA Kit
E-CL-M0274	Mouse EMILIN1 (Elastin Microfibril Interface Located Protein 1) CLIA Kit
E-CL-M0275	Mouse ERN1 (Endoplasmic Reticulum To Nucleus Signalling 1) CLIA Kit
E-CL-M0276	Mouse ES (Endostatin) CLIA Kit
E-CL-M0277	Mouse EL (Endothelial lipase) CLIA Kit
E-CL-M0278	Mouse NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) CLIA Kit
E-CL-M0279	Mouse ET-1 (Endothelin 1) CLIA Kit
E-CL-M0280	Mouse ECE1 (Endothelin Converting Enzyme 1) CLIA Kit
E-CL-M0281	Mouse Ear 14 (Eosinophil-associated, ribonuclease A family, member 14) CLIA Kit
E-CL-M0282	Mouse ECF/CCL11 (Eosinophil Chemotactic Factor) CLIA Kit
E-CL-M0283	Mouse MPIF2 (Myeloid Progenitor Inhibitory Factor 2) CLIA Kit
E-CL-M0284	Mouse EPI (Epinephrine/Adrenaline) CLIA Kit
E-CL-M0285	Mouse EPR (Epiregulin) CLIA Kit
E-CL-M0286	Mouse EPCAM (Epithelial Cell Adhesion Molecule) CLIA Kit
E-CL-M0287	Mouse ENA-78 (Epithelial Neutrophil Activating Peptide 78) CLIA Kit
E-CL-M0288	Mouse SELE (E-Selectin) CLIA Kit

Cat.NO	Product Name
E-CL-M0289	Mouse ER α (Estrogen Receptor Alpha) CLIA Kit
E-CL-M0290	Mouse ECM1 (Extracellular Matrix Protein 1) CLIA Kit
E-CL-M0291	Mouse FVIII Δ g (Factor VIII-related Antigen) CLIA Kit
E-CL-M0292	Mouse FAS/CD95 (Factor Related Apoptosis) CLIA Kit
E-CL-M0293	Mouse FABP5 (Fatty Acid Binding Protein 5, Epidermal) CLIA Kit
E-CL-M0294	Mouse FABP6 (Fatty Acid Binding Protein 6, Ileal) CLIA Kit
E-CL-M0295	Mouse ER β (Estrogen Receptor Beta) CLIA Kit
E-CL-M0296	Mouse FE (Ferritin) CLIA Kit
E-CL-M0297	Mouse FBN1 (Fibrillin 1) CLIA Kit
E-CL-M0298	Mouse FGA (Fibrinogen Alpha chain) CLIA Kit
E-CL-M0299	Mouse FGB (Fibrinogen Beta chain) CLIA Kit
E-CL-M0300	Mouse FDP (Fibrinogen Degradation Product) CLIA Kit
E-CL-M0301	Mouse FGL2 (Fibrinogen Like Protein 2) CLIA Kit
E-CL-M0302	Mouse Fbg (Fibrinogen) CLIA Kit
E-CL-M0303	Mouse FGF7 (Fibroblast Growth Factor 7) CLIA Kit
E-CL-M0304	Mouse FGF9 (Fibroblast Growth Factor 9) CLIA Kit
E-CL-M0305	Mouse FRS2 (Fibroblast Growth Factor Receptor Substrate 2) CLIA Kit
E-CL-M0306	Mouse FMOD (Fibromodulin) CLIA Kit
E-CL-M0307	Mouse FN (Fibronectin) CLIA Kit

Cat.NO	Product Name
E-CL-M0308	Mouse FSH (follicle-stimulating hormone) CLIA Kit
E-CL-M0309	Mouse FST (Follistatin) CLIA Kit
E-CL-M0310	Mouse FPR2 (Formyl Peptide Receptor 2) CLIA Kit
E-CL-M0311	Mouse FK (Fractalkine) CLIA Kit
E-CL-M0312	Mouse F-TESTO (Free Testosterone) CLIA Kit
E-CL-M0313	Mouse fT4 (Free Thyroxine) CLIA Kit
E-CL-M0314	Mouse fT3 (Free Triiodothyronine Index) CLIA Kit
E-CL-M0315	Mouse GNB1 (G Protein Beta 1) CLIA Kit
E-CL-M0316	Mouse GLB (Galactosidase, Beta) CLIA Kit
E-CL-M0317	Mouse GAL1 (Galectin 1) CLIA Kit
E-CL-M0318	Mouse GAL2 (Galectin 2) CLIA Kit
E-CL-M0319	Mouse GAL3 (Galectin 3) CLIA Kit
E-CL-M0320	Mouse GAL4 (Galectin 4) CLIA Kit
E-CL-M0321	Mouse GAL6 (Galectin 6) CLIA Kit
E-CL-M0322	Mouse GAL7 (Galectin 7) CLIA Kit
E-CL-M0323	Mouse GAL8 (Galectin 8) CLIA Kit
E-CL-M0324	Mouse GAL9 (Galectin 9) CLIA Kit
E-CL-M0325	Mouse GAL12 (Galectin 12) CLIA Kit
E-CL-M0326	Mouse gGT1 (Gamma Glutamyltransferase 1) CLIA Kit

Cat.NO	Product Name
E-CL-M0327	Mouse GIP (Gastric Inhibitory Polypeptide) CLIA Kit
E-CL-M0328	Mouse GT-34 (Gastrin 34) CLIA Kit
E-CL-M0329	Mouse GKN2 (Gastrokine 2) CLIA Kit
E-CL-M0330	Mouse GS (Gelsolin) CLIA Kit
E-CL-M0331	Mouse GTF2H1 (General Transcription Factor IIH, Polypeptide 1) CLIA Kit
E-CL-M0332	Mouse GDNF (Glial Cell Line Derived Neurotrophic Factor) CLIA Kit
E-CL-M0333	Mouse GFAP (glial fibrillary acidic protein) CLIA Kit
E-CL-M0334	Mouse GC (Glucagon) CLIA Kit
E-CL-M0335	Mouse GR α (glucocorticoid receptor α) CLIA Kit
E-CL-M0336	Mouse GR β (glucocorticoid receptor β) CLIA Kit
E-CL-M0337	Mouse GCK (Glucokinase) CLIA Kit
E-CL-M0338	Mouse gp130 (Glucoprotein 130) CLIA Kit
E-CL-M0339	Mouse G6PD (Glucose 6 Phosphate Dehydrogenase) CLIA Kit
E-CL-M0340	Mouse GLUT4 (Glucose Transporter 4) CLIA Kit
E-CL-M0341	Mouse GPI (Glucose 6 Phosphate Isomerase) CLIA Kit
E-CL-M0342	Mouse GUSb (Glucuronidase, Beta) CLIA Kit
E-CL-M0343	Mouse GDH/GLDH (Glutamate dehydrogenase) CLIA Kit
E-CL-M0344	Mouse OXB (Orexin B) CLIA Kit
E-CL-M0345	Mouse HB (Hemoglobin) CLIA Kit

Cat.NO	Product Name
E-CL-M0346	Mouse GR (Glutathione Reductase) CLIA Kit
E-CL-M0347	Mouse GSTa1 (Glutathione S Transferase Alpha 1) CLIA Kit
E-CL-M0348	Mouse GSTo1 (Glutathione S Transferase Omega 1) CLIA Kit
E-CL-M0349	Mouse GSTt1 (Glutathione S Transferase Theta 1) CLIA Kit
E-CL-M0350	Mouse GSTt2 (Glutathione S Transferase Theta 2) CLIA Kit
E-CL-M0351	Mouse HbA1C (Glycosylated Hemoglobin/Hemoglobin A1c) CLIA Kit
E-CL-M0352	Mouse PYGL (Glycogen Phosphorylase, Liver) CLIA Kit
E-CL-M0353	Mouse GSK3a (Glycogen Synthase Kinase 3 Alpha) CLIA Kit
E-CL-M0354	Mouse GPC4 (Glypican 4) CLIA Kit
E-CL-M0355	Mouse GnRH (Gonadotropin-Releasing Hormone) CLIA Kit
E-CL-M0356	Mouse GRN (Granulin) CLIA Kit
E-CL-M0357	Mouse GCP-2 (Granulocyte Chemotactic Protein 2) CLIA Kit
E-CL-M0358	Mouse Gzms-A (granzyme A) CLIA Kit
E-CL-M0359	Mouse Gzms-B (granzyme B) CLIA Kit
E-CL-M0360	Mouse GREM1 (GremLin 1) CLIA Kit
E-CL-M0361	Mouse GDF1 (Growth Differentiation Factor 1) CLIA Kit
E-CL-M0362	Mouse GDF2 (Growth Differentiation Factor 2) CLIA Kit
E-CL-M0363	Mouse GDF3 (Growth Differentiation Factor 3) CLIA Kit
E-CL-M0364	Mouse GDF5 (Growth Differentiation Factor 5) CLIA Kit

Cat.NO	Product Name
E-CL-M0365	Mouse GDF7 (Growth Differentiation Factor 7) CLIA Kit
E-CL-M0366	Mouse GDF9 (Growth Differentiation Factor 9) CLIA Kit
E-CL-M0367	Mouse GDF11 (Growth Differentiation Factor 11) CLIA Kit
E-CL-M0368	Mouse GDF15 (Growth Differentiation Factor 15) CLIA Kit
E-CL-M0369	Mouse GH2 (Growth Hormone 2) CLIA Kit
E-CL-M0370	Mouse GHRH (Growth Hormone Releasing Hormone) CLIA Kit
E-CL-M0371	Mouse CAM (Calmodulin) CLIA Kit
E-CL-M0372	Mouse HP (Haptoglobin) CLIA Kit
E-CL-M0373	Mouse HSF1 (Heat Shock Factor 1) CLIA Kit
E-CL-M0374	Mouse HSP-60 (Heat Shock Protein 60) CLIA Kit
E-CL-M0375	Mouse HSP-70 (Heat Shock Protein 70) CLIA Kit
E-CL-M0376	Mouse HSP-90 (Heat Shock Protein 90) CLIA Kit
E-CL-M0377	Mouse HSF2 (Heat Shock Transcription Factor 2) CLIA Kit
E-CL-M0378	Mouse SHH (Hedgehog Homolog, Sonic) CLIA Kit
E-CL-M0379	Mouse HHIP (Hedgehog Interacting Protein) CLIA Kit
E-CL-M0380	Mouse MMP-3 (Matrix Metalloproteinase 3) CLIA Kit
E-CL-M0381	Mouse MMP-9 (Matrix Metalloproteinase-9) CLIA Kit
E-CL-M0382	Mouse OPN (Osteopontin) CLIA Kit
E-CL-M0383	Mouse PDGF-BB (Platelet-Derived Growth Factor-BB) CLIA Kit

Cat.NO	Product Name
E-CL-M0384	Mouse POSTN/OSF-2 (Periostin) CLIA Kit
E-CL-M0385	Mouse PCSK9 (Proprotein Convertase Subtilisin/Kexin Type 9) CLIA Kit
E-CL-M0386	Mouse RETN (Resistin) CLIA Kit
E-CL-M0387	Mouse SCF (Stem Cell Factor) CLIA Kit
E-CL-M0388	Mouse SELL (L-Selectin) CLIA Kit
E-CL-M0389	Mouse SELP (P-Selectin) CLIA Kit
E-CL-M0390	Mouse SP (Substance P) CLIA Kit
E-CL-M0391	Mouse TPO (Thrombopoietin) CLIA Kit
E-CL-M0392	Mouse TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) CLIA Kit
E-CL-M0393	Mouse TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) CLIA Kit
E-CL-M0394	Mouse RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) CLIA Kit
E-CL-M0395	Mouse TREM-1 (Triggering Receptor Expresses on Myeloid Cells-1) CLIA Kit
E-CL-M0396	Mouse TSLP (Thymic Stromal Lymphopoietin) CLIA Kit
E-CL-M0397	Mouse sVCAM-1 (soluble vascular cell adhesion molecule 1) CLIA Kit
E-CL-M0398	Mouse FLT1/VEGFR1 (Vascular Endothelial Growth Factor Receptor 1) CLIA Kit
E-CL-M0399	Mouse VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) CLIA Kit
E-CL-M0400	Mouse HO-1 (Heme Oxygenase 1) CLIA Kit
E-CL-M0401	Mouse HO-2 (Heme Oxygenase 2, Decycling) CLIA Kit
E-CL-M0402	Mouse RANK (Receptor Activator of Nuclear Factor Kappa B) CLIA Kit

Cat.NO	Product Name
E-CL-M0403	Mouse HJV (Hemojuvelin) CLIA Kit
E-CL-M0404	Mouse HPA (Heparanase) CLIA Kit
E-CL-M0405	Mouse HC II (Heparin Cofactor II) CLIA Kit
E-CL-M0406	Mouse HB-EGF (Heparin-binding Epidermal Growth Factor-like Growth Factor) CLIA Kit
E-CL-M0407	Mouse LIPC (Lipase, Hepatic) CLIA Kit
E-CL-M0408	Mouse HGF (Hepatocyte Growth Factor) CLIA Kit
E-CL-M0409	Mouse HNF1a (Hepatocyte Nuclear Factor 1 Alpha) CLIA Kit
E-CL-M0410	Mouse HNF1b (Hepatocyte Nuclear Factor 1 Beta) CLIA Kit
E-CL-M0411	Mouse HNF4a (Hepatocyte Nuclear Factor 4 Alpha) CLIA Kit
E-CL-M0412	Mouse HNF6a (Hepatocyte Nuclear Factor 6 Alpha) CLIA Kit
E-CL-M0413	Mouse HNF6b (Hepatocyte Nuclear Factor 6 Beta) CLIA Kit
E-CL-M0414	Mouse HDGF (Hepatoma Derived Growth Factor)CLIA Kit
E-CL-M0415	Mouse Hpc (Hepcidin) CLIA Kit
E-CL-M0416	Mouse HNRPA1 (Heterogeneous Nuclear Ribonucleoprotein A1) CLIA Kit
E-CL-M0417	Mouse HK1 (Hexokinase 1) CLIA Kit
E-CL-M0418	Mouse HK (Hexokinase) CLIA Kit
E-CL-M0419	Mouse HMGB-1 (High mobility group protein B1) CLIA Kit
E-CL-M0420	Mouse hs-CRP (high sensitivity C-Reactive Protein) CLIA Kit
E-CL-M0421	Mouse histon-H2b CLIA Kit

Cat.NO	Product Name
E-CL-M0422	Mouse HABP1 (Hyaluronan Binding Protein 1) CLIA Kit
E-CL-M0423	Mouse HAase (Hyaluronidase) CLIA Kit
E-CL-M0424	Mouse HIF-1 α (Hypoxia Inducible Factor 1 Alpha) CLIA Kit
E-CL-M0425	Mouse IgA (Immunoglobulin A) CLIA Kit
E-CL-M0426	Mouse IgE (Immunoglobulin E) CLIA Kit
E-CL-M0427	Mouse IgG (Immunoglobulin G) CLIA Kit
E-CL-M0428	Mouse IgG1 (Immunoglobulin G1) CLIA Kit
E-CL-M0429	Mouse IgG2a (Immunoglobulin G2a) CLIA Kit
E-CL-M0430	Mouse IgM (Immunoglobulin M) CLIA Kit
E-CL-M0431	Mouse NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) CLIA Kit
E-CL-M0432	Mouse INHA (Inhibin A) CLIA Kit
E-CL-M0433	Mouse INHB (Inhibin B) CLIA Kit
E-CL-M0434	Mouse INHbA (Inhibin Beta A) CLIA Kit
E-CL-M0435	Mouse INHBP (Inhibin Binding Protein) CLIA Kit
E-CL-M0436	Mouse IAA (Insulin Autoantibodies) CLIA Kit
E-CL-M0437	Mouse INSL5 (Insulin Like Protein 5) CLIA Kit
E-CL-M0438	Mouse INSR- β (Insulin Receptor Beta Subunit) CLIA Kit
E-CL-M0439	Mouse IGF-2 (Insulin Like Growth Factor 2) CLIA Kit
E-CL-M0440	Mouse IGFBP-1 (Insulin-Like Growth Factor Binding Protein 1) CLIA Kit

Cat.NO	Product Name
E-CL-M0441	Mouse IGFBP-2 (Insulin-Like Growth Factor Binding Protein 2) CLIA Kit
E-CL-M0442	Mouse IGFBP-4 (Insulin-Like Growth Factor Binding Protein 4) CLIA Kit
E-CL-M0443	Mouse IGFBP-6 (Insulin-Like Growth Factor Binding Protein 6) CLIA Kit
E-CL-M0444	Mouse I-PTH (intact Parathormone) CLIA Kit
E-CL-M0445	Mouse ICAM-1 (intercellular adhesion molecule 1) CLIA Kit
E-CL-M0446	Mouse ICAM-2 (intercellular adhesion molecule 2) CLIA Kit
E-CL-M0447	Mouse ICAM-3 (intercellular adhesion molecule 3) CLIA Kit
E-CL-M0448	Mouse IFN- α (Interferon α) CLIA Kit
E-CL-M0449	Mouse FOXO1 (Forkhead Box Protein O1) CLIA Kit
E-CL-M0450	Mouse IFI16/p16 (Interferon-inducible Protein 16) CLIA Kit
E-CL-M0451	Mouse IFN- α -Ab (Interferon- α Antibody) CLIA Kit
E-CL-M0452	Mouse UCP-1 (Uncoupling Protein 1, Mitochondrial) CLIA Kit
E-CL-M0453	Mouse IL-34 (Interleukin 34) CLIA Kit
E-CL-M0454	Mouse IL-3 (Interleukin 3) CLIA Kit
E-CL-M0455	Mouse IL-5 (Interleukin 5) CLIA Kit
E-CL-M0456	Mouse IL-7 (Interleukin 7) CLIA Kit
E-CL-M0457	Mouse IL-9 (Interleukin 9) CLIA Kit
E-CL-M0458	Mouse IL-11 (Interleukin 11) CLIA Kit
E-CL-M0459	Mouse IL-12 (Interleukin 12) CLIA Kit

Cat.NO	Product Name
E-CL-M0460	Mouse IL-13 (Interleukin 13) CLIA Kit
E-CL-M0461	Mouse IL-15 (Interleukin 15) CLIA Kit
E-CL-M0462	Mouse IL-16 (Interleukin 16) CLIA Kit
E-CL-M0463	Mouse IL-18 (Interleukin 18) CLIA Kit
E-CL-M0464	Mouse IL-23 (Interleukin 23) CLIA Kit
E-CL-M0465	Mouse IL-27 (Interleukin 27) CLIA Kit
E-CL-M0466	Mouse IL-35 (Interleukin 35) CLIA Kit
E-CL-M0467	Mouse IL-2R (Interleukin-2 Receptor) CLIA Kit
E-CL-M0468	Mouse IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) CLIA Kit
E-CL-M0469	Mouse IMA (Ischemia Modified Albumin) CLIA Kit
E-CL-M0470	Mouse IL-32 (Interleukin 32) CLIA Kit
E-CL-M0471	Mouse IL18BP (Interleukin 18 Binding Protein) CLIA Kit
E-CL-M0472	Mouse AR (Amphiregulin) CLIA Kit
E-CL-M0473	Mouse KIM-1 (Kidney Injury Molecule 1) CLIA Kit
E-CL-M0474	Mouse LDH (Lactate Dehydrogenase) CLIA Kit
E-CL-M0475	Mouse LTF/LF (Lactoferrin) CLIA Kit
E-CL-M0476	Mouse LN (Laminin) CLIA Kit
E-CL-M0477	Mouse LEPR (Leptin Receptor) CLIA Kit
E-CL-M0478	Mouse LIFR (Leukemia Inhibitory Factor Receptor) CLIA Kit

Cat.NO	Product Name
E-CL-M0479	Mouse LT-B4 (Leukotriene B4) CLIA Kit
E-CL-M0480	Mouse LT-C4 (Leukotriene C4) CLIA Kit
E-CL-M0481	Mouse LAT (Linker for activation of T cell) CLIA Kit
E-CL-M0482	Mouse LPL (Lipoprotein Lipase) CLIA Kit
E-CL-M0483	Mouse Lp-a (Lipoprotein a) CLIA Kit
E-CL-M0484	Mouse Lp-PL-A2 (Lipoprotein-associated Phospholipase A2) CLIA Kit
E-CL-M0485	Mouse Lptn/XCL1 (Lymphotactin) CLIA Kit
E-CL-M0486	Mouse LZM (Lysozyme) CLIA Kit
E-CL-M0487	Mouse MIP-3 β (Macrophage Inflammatory Protein 3 β) CLIA Kit
E-CL-M0488	Mouse MIF (Macrophage Migration Inhibitory Factor) CLIA Kit
E-CL-M0489	Mouse MHCI/H-2 I (Major Histocompatibility Complex I) CLIA Kit
E-CL-M0490	Mouse MMP-1 (Matrix Metalloproteinase 1) CLIA Kit
E-CL-M0491	Mouse MMP-2 (Matrix Metalloproteinase 2) CLIA Kit
E-CL-M0492	Mouse MMP-17 (Matrix Metalloproteinase 17) CLIA Kit
E-CL-M0493	Mouse MMP-24 (Matrix Metalloproteinase 24) CLIA Kit
E-CL-M0494	Mouse MMP-7 (Matrix Metalloproteinase 7) CLIA Kit
E-CL-M0495	Mouse MMP-8 (Matrix Metalloproteinase 8) CLIA Kit
E-CL-M0496	Mouse MT (Melatonin) CLIA Kit
E-CL-M0497	Mouse MHB (Methemoglobin) CLIA Kit

Cat.NO	Product Name
E-CL-M0498	Mouse MCP-2 (Monocyte Chemotactic Protein 2) CLIA Kit
E-CL-M0499	Mouse MCP-3 (Monocyte Chemotactic Protein 3) CLIA Kit
E-CL-M0500	Mouse MTL (Motilin) CLIA Kit
E-CL-M0501	Mouse MUC5AC (Mucin-5 subtype AC) CLIA Kit
E-CL-M0502	Mouse MUC5B (Mucin-5 Subtype B) CLIA Kit
E-CL-M0503	Mouse CHRM2 (Cholinergic Receptor, Muscarinic 2) CLIA Kit
E-CL-M0504	Mouse MBP (Myelin Basic Protein) CLIA Kit
E-CL-M0505	Mouse NAGase (N-Acetyl Beta-D-Glucosaminidase) CLIA Kit
E-CL-M0506	Mouse NGF (Nerve growth factor) CLIA Kit
E-CL-M0507	Mouse NGB (Neuroglobin) CLIA Kit
E-CL-M0508	Mouse Slit2 (Slit Homolog 2) CLIA Kit
E-CL-M0509	Mouse NPY (Neuropeptide Y) CLIA Kit
E-CL-M0510	Mouse NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) CLIA Kit
E-CL-M0511	Mouse NT-3 (Neurotrophin 3) CLIA Kit
E-CL-M0512	Mouse NT-4 (Neurotrophin 4) CLIA Kit
E-CL-M0513	Mouse NAP3 (Neutrophil Activating Protein 3) CLIA Kit
E-CL-M0514	Mouse NGAL (Neutrophil Gelatinase Associated Lipocalin) CLIA Kit
E-CL-M0515	Mouse NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) CLIA Kit
E-CL-M0516	Mouse NADPH (Nicotinamide Adenine Dinucleotide Phosphate) CLIA Kit

Cat.NO	Product Name
E-CL-M0517	Mouse NA/NE (Noradrenaline/Norepinephrine) CLIA Kit
E-CL-M0518	Mouse NT-proBNP (N-terminal pro-Brain Natriuretic Peptide) CLIA Kit
E-CL-M0519	Mouse P III NP (N-terminal procollagen III propeptide) CLIA Kit
E-CL-M0520	Mouse NFkB (Nuclear Factor Kappa B) CLIA Kit
E-CL-M0521	Mouse NFIL3 (Nuclear Factor, Interleukin 3 Regulated) CLIA Kit
E-CL-M0522	Mouse NRIP1 (Nuclear Receptor Interacting Protein 1) CLIA Kit
E-CL-M0523	Mouse NUP160 (Nucleoporin 160kDa) CLIA Kit
E-CL-M0524	Mouse OLFM4 (Olfactomedin 4) CLIA Kit
E-CL-M0525	Mouse ITLN1 (Intelectin 1/Omentin) CLIA Kit
E-CL-M0526	Mouse OSM (Oncostatin-M) CLIA Kit
E-CL-M0527	Mouse OXA (Orexin A) CLIA Kit
E-CL-M0528	Mouse ODC (Ornithine Decarboxylase) CLIA Kit
E-CL-M0529	Mouse OC/BGP (Osteocalcin) CLIA Kit
E-CL-M0530	Mouse ON (Osteonectin) CLIA Kit
E-CL-M0531	Mouse OVA sIgE (Ovalbumin Specific IgE) CLIA Kit
E-CL-M0532	Mouse CA125 (Carbohydrate Antigen 125) CLIA Kit
E-CL-M0533	Mouse OxLDL (Oxidized Lowdensity Lipoprotein) CLIA Kit
E-CL-M0534	Mouse MAPK P38 (P38 Mitogen-activated Protein Kinase) CLIA Kit
E-CL-M0535	Mouse PP (Pancreatic Polypeptide) CLIA Kit

Cat.NO	Product Name
E-CL-M0536	Mouse PTH (Parathormone) CLIA Kit
E-CL-M0537	Mouse PTHrP (Parathyroid Hormone Related Protein) CLIA Kit
E-CL-M0538	Mouse PTD (Pentosidine) CLIA Kit
E-CL-M0539	Mouse PP (Pepsin) CLIA Kit
E-CL-M0540	Mouse PGC (Pepsinogen C) CLIA Kit
E-CL-M0541	Mouse PRF1 (Perforin 1) CLIA Kit
E-CL-M0542	Mouse PMP3 (Peroxisomal Membrane Protein 3) CLIA Kit
E-CL-M0543	Mouse PPAR- γ (Peroxisome Proliferator-activated receptor γ) CLIA Kit
E-CL-M0544	Mouse PPAR- α (Peroxisome Proliferators-activator Receptors alpha) CLIA Kit
E-CL-M0545	Mouse PI Ab-IgG/IgM (phosphatidylinositol antibody IgG/IgM) CLIA Kit
E-CL-M0546	Mouse PS (Phosphatidylserine) CLIA Kit
E-CL-M0547	Mouse PEDF (Pigment Epithelium Derived Factor) CLIA Kit
E-CL-M0548	Mouse PACAP (Pituitary Adenylate Cyclase Activating Polypeptide) CLIA Kit
E-CL-M0549	Mouse PGF (Placental Growth Factor) CLIA Kit
E-CL-M0550	Mouse PAP (Plasmin-Antiplasmin Complex) CLIA Kit
E-CL-M0551	Mouse PAI-1 (Plasminogen Activator Inhibitor 1) CLIA Kit
E-CL-M0552	Mouse PAI-2 (Plasminogen Activator Inhibitor 2) CLIA Kit
E-CL-M0553	Mouse tPA (Plasminogen Activator, Tissue) CLIA Kit
E-CL-M0554	Mouse PAF (Platelet Activating Factor) CLIA Kit

Cat.NO	Product Name
E-CL-M0555	Mouse PA-IgG/M/A (Platelet Antibodies IgG/M/A) CLIA Kit
E-CL-M0556	Mouse PDGFC (Platelet Derived Growth Factor C) CLIA Kit
E-CL-M0557	Mouse PDGFD (Platelet Derived Growth Factor D) CLIA Kit
E-CL-M0558	Mouse PDGFRL (Platelet Derived Growth Factor Receptor Like Protein) CLIA Kit
E-CL-M0559	Mouse PDGFB (Platelet Derived Growth Factor Subunit B) CLIA Kit
E-CL-M0560	Mouse PF4 (Platelet Factor 4) CLIA Kit
E-CL-M0561	Mouse GP II b III a (Platelet Membrane Glycoprotein II b III a) CLIA Kit
E-CL-M0562	Mouse PECAM1 (Platelet/Endothelial Cell Adhesion Molecule 1) CLIA Kit
E-CL-M0563	Mouse PDGFsR- α (Platelet-Derived Growth Factor Soluble Receptor α) CLIA Kit
E-CL-M0564	Mouse PTN (Pleiotrophin) CLIA Kit
E-CL-M0565	Mouse PCX (Podocalyxin) CLIA Kit
E-CL-M0566	Mouse PAPP-A (pregnancy associated plasma protein A) CLIA Kit
E-CL-M0567	Mouse Pro-ANP (Pro Atrial Natriuretic Peptide) CLIA Kit
E-CL-M0568	Mouse PCT (Procalcitonin) CLIA Kit
E-CL-M0569	Mouse P III NP (Procollagen III N-terminal Peptide) CLIA Kit
E-CL-M0570	Mouse Pg (Progesterone) CLIA Kit
E-CL-M0571	Mouse PHB (Prohibitin) CLIA Kit
E-CL-M0572	Mouse PI (Proinsulin) CLIA Kit
E-CL-M0573	Mouse PK2 (Prokineticin 2) CLIA Kit

Cat.NO	Product Name
E-CL-M0574	Mouse PKR1 (Prokineticin Receptor 1) CLIA Kit
E-CL-M0575	Mouse POMC (Proopiomelanocortin) CLIA Kit
E-CL-M0576	Mouse PGD2S (Prostaglandin D2 Synthase) CLIA Kit
E-CL-M0577	Mouse PTGES (Prostaglandin E Synthase, Microsomal) CLIA Kit
E-CL-M0578	Mouse PGE1 (Prostaglandin E1) CLIA Kit
E-CL-M0579	Mouse PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) CLIA Kit
E-CL-M0580	Mouse PSA (Prostate Specific Antigen) CLIA Kit
E-CL-M0581	Mouse PAR2 (Protease Activated Receptor 2) CLIA Kit
E-CL-M0582	Mouse PRSS1 (Protease, Serine, 1) CLIA Kit
E-CL-M0583	Mouse PROCR (Protein C Receptor, Endothelial) CLIA Kit
E-CL-M0584	Mouse PDIA3 (Protein Disulfide Isomerase A3) CLIA Kit
E-CL-M0585	Mouse PDI (Protein Disulfide Isomerase) CLIA Kit
E-CL-M0586	Mouse PIAS1 (Protein Inhibitor Of Activated STAT 1) CLIA Kit
E-CL-M0587	Mouse PARC (Pulmonary Activation Regulated Chemokine) CLIA Kit
E-CL-M0588	Mouse SPA (Pulmonary Surfactant Associated Protein A) CLIA Kit
E-CL-M0589	Mouse SPD (Pulmonary Surfactant Associated Protein D) CLIA Kit
E-CL-M0590	Mouse PDHa (Pyruvate Dehydrogenase Alpha) CLIA Kit
E-CL-M0591	Mouse PDP (Pyruvate Dehydrogenase Phosphatase) CLIA Kit
E-CL-M0592	Mouse M2-PK (Pyruvate Kinase) CLIA Kit

Cat.NO	Product Name
E-CL-M0593	Mouse PK (Pyruvate Kinase) CLIA Kit
E-CL-M0594	Mouse REG4 (Regenerating Islet Derived Protein 4) CLIA Kit
E-CL-M0595	Mouse RBP4 (Retinol Binding Protein 4, Plasma) CLIA Kit
E-CL-M0596	Mouse RNASET2 (Ribonuclease T2) CLIA Kit
E-CL-M0597	Mouse RUNX2 (Runt Related Transcription Factor 2) CLIA Kit
E-CL-M0598	Mouse S100A11 (S100 Calcium Binding Protein A11) CLIA Kit
E-CL-M0599	Mouse S100A7 (S100 Calcium Binding Protein A7) CLIA Kit
E-CL-M0600	Mouse S100B (S100 Calcium Binding Protein B) CLIA Kit
E-CL-M0601	Mouse S100A6 (S100 Calcium Binding Protein A6) CLIA Kit
E-CL-M0602	Mouse SCT (Secretin) CLIA Kit
E-CL-M0603	Mouse slgA (Secretory Immunoglobulin A) CLIA Kit
E-CL-M0604	Mouse SELE (Selectin, Endothelium) CLIA Kit
E-CL-M0605	Mouse SEPP1 (Selenoprotein P) CLIA Kit
E-CL-M0606	Mouse SEMA7A (Semaphorin 7A) CLIA Kit
E-CL-M0607	Mouse SAA (Serum amyloid A) CLIA Kit
E-CL-M0608	Mouse SESN1 (Sestrin 1) CLIA Kit
E-CL-M0609	Mouse SESN2 (Sestrin 2) CLIA Kit
E-CL-M0610	Mouse SESN3 (Sestrin 3) CLIA Kit
E-CL-M0611	Mouse SHBG (Sex Hormone-Binding Globulin) CLIA Kit

Cat.NO	Product Name
E-CL-M0612	Mouse ENG (Endoglin) CLIA Kit
E-CL-M0613	Mouse EPCR (Endothelial Protein C Receptor) CLIA Kit
E-CL-M0614	Mouse IL-2sR α /CD25 (Soluble Interleukin-2 Receptor) CLIA Kit
E-CL-M0615	Mouse IL-2sR β /CD122 (Soluble Interleukin-2 Receptor) CLIA Kit
E-CL-M0616	Mouse TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) CLIA Kit
E-CL-M0617	Mouse SS (Somatostatin) CLIA Kit
E-CL-M0618	Mouse SCCA2 (Squamous Cell Carcinoma Antigen 2) CLIA Kit
E-CL-M0619	Mouse SCFR (Stem Cell Factor Receptor) CLIA Kit
E-CL-M0620	Mouse SDF1 (Stromal Cell Derived Factor 1) CLIA Kit
E-CL-M0621	Mouse SDF4 (Stromal Cell Derived Factor 4) CLIA Kit
E-CL-M0622	Mouse SOCS3 (Suppressors Of Cytokine Signaling 3) CLIA Kit
E-CL-M0623	Mouse SYP (Synaptophysin) CLIA Kit
E-CL-M0624	Mouse SYT1 (Synaptotagmin 1) CLIA Kit
E-CL-M0625	Mouse SDC4 (Syndecan 4) CLIA Kit
E-CL-M0626	Mouse SNCa (Synuclein, Alpha) CLIA Kit
E-CL-M0627	Mouse TACR1 (Tachykinin Receptor 1) CLIA Kit
E-CL-M0628	Mouse TACR2 (Tachykinin Receptor 2) CLIA Kit
E-CL-M0629	Mouse TLN (Talin) CLIA Kit
E-CL-M0630	Mouse THP (Tamm–Horsfall Glycoprotein) CLIA Kit

Cat.NO	Product Name
E-CL-M0631	Mouse TRACP-5b (Tartrate-Resistant Acid Phosphatase 5b) CLIA Kit
E-CL-M0632	Mouse MAP τ (Microtubule Associated Protein Tau/Tau Protein) CLIA Kit
E-CL-M0633	Mouse TDG1 (Teratocarcinoma Derived Growth Factor 1) CLIA Kit
E-CL-M0634	Mouse TBP2 (Thioredoxin Binding Protein 2) CLIA Kit
E-CL-M0635	Mouse TrxR1 (Thioredoxin Reductase 1) CLIA Kit
E-CL-M0636	Mouse Trx (Thioredoxin) CLIA Kit
E-CL-M0637	Mouse TAFI (Thrombin activatable fibrinolysis inhibitor) CLIA Kit
E-CL-M0638	Mouse TSP-1 (Thrombin Sensitive Protein 1) CLIA Kit
E-CL-M0639	Mouse TAT (Thrombin-Antithrombin Complex) CLIA Kit
E-CL-M0640	Mouse TM (Thrombomodulin) CLIA Kit
E-CL-M0641	Mouse Thbs-1 (Thrombospondin 1) CLIA Kit
E-CL-M0642	Mouse TXB2 (Thromboxane B2) CLIA Kit
E-CL-M0643	Mouse TpP (Thrombus Precursor Protein) CLIA Kit
E-CL-M0644	Mouse TG (Thyroglobulin) CLIA Kit
E-CL-M0645	Mouse TSH (Thyroid Stimulating Hormone) CLIA Kit
E-CL-M0646	Mouse TRH (Thyrotropin-releasing Hormone) CLIA Kit
E-CL-M0647	Mouse TBG (Thyroxine Binding Globulin) CLIA Kit
E-CL-M0648	Mouse T4 (Thyroxine) CLIA Kit
E-CL-M0649	Mouse TJP1 (Tight Junction Protein 1) CLIA Kit

Cat.NO	Product Name
E-CL-M0650	Mouse TFPI (Tissue Factor Pathway Inhibitor) CLIA Kit
E-CL-M0651	Mouse TF (Tissue factor) CLIA Kit
E-CL-M0652	Mouse TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) CLIA Kit
E-CL-M0653	Mouse TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) CLIA Kit
E-CL-M0654	Mouse TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) CLIA Kit
E-CL-M0655	Mouse TPA (Tissue Polypeptide Antigen) CLIA Kit
E-CL-M0656	Mouse TLR-9 (Toll-like Receptor 9) CLIA Kit
E-CL-M0657	Mouse TRF (Transferrin) CLIA Kit
E-CL-M0658	Mouse TFR2 (Transferrin Receptor Protein 2) CLIA Kit
E-CL-M0659	Mouse TGF β 1 (Transforming Growth Factor Beta Induced) CLIA Kit
E-CL-M0660	Mouse TGF- α (Transforming Growth Factor α) CLIA Kit
E-CL-M0661	Mouse TGF β 2 (Transforming Growth Factor Beta 2) CLIA Kit
E-CL-M0662	Mouse TGF- β 3 (Transforming Growth Factor β 3) CLIA Kit
E-CL-M0663	Mouse TNPO1 (Transportin 1) CLIA Kit
E-CL-M0664	Mouse TTR (Transthyretin) CLIA Kit
E-CL-M0665	Mouse TFF2 (Trefoil Factor 2) CLIA Kit
E-CL-M0666	Mouse TFF3 (Trefoil Factor 3, Intestinal) CLIA Kit
E-CL-M0667	Mouse T3 (Triiodothyronine) CLIA Kit
E-CL-M0668	Mouse TNNI3/cTn-I (Troponin I Type 3, Cardiac) CLIA Kit

Cat.NO	Product Name
E-CL-M0669	Mouse TRY (Trypsin) CLIA Kit
E-CL-M0670	Mouse TAP (Trypsinogen Activation Peptide) CLIA Kit
E-CL-M0671	Mouse TNF- β (Tumor necrosis factor β) CLIA Kit
E-CL-M0672	Mouse Tumor necrosis factor-related apoptosis-inducing ligand CLIA Kit
E-CL-M0673	Mouse TP53 (Tumor Protein 53) CLIA Kit
E-CL-M0674	Mouse TRY (Tyrosinase) CLIA Kit
E-CL-M0675	Mouse Ub (Ubiquitin) CLIA Kit
E-CL-M0676	Mouse PLAU/uPA(Urokinase-Type Plasminogen Activator) CLIA Kit
E-CL-M0677	Mouse Cr (Creatinine) CLIA Kit
E-CL-M0678	Mouse VEGF-B (Vascular Endothelial cell Growth Factor B) CLIA Kit
E-CL-M0679	Mouse VEGF-C (Vascular Endothelial cell Growth Factor C) CLIA Kit
E-CL-M0680	Mouse VEGF-D (Vascular Endothelial cell Growth Factor D) CLIA Kit
E-CL-M0681	Mouse VEGFR3/Flt4 (Vascular Endothelial Cell Growth Factor Receptor 3) CLIA Kit
E-CL-M0682	Mouse VIP (Vasoactive Intestinal Peptide) CLIA Kit
E-CL-M0683	Mouse FOS (V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) CLIA Kit
E-CL-M0684	Mouse VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) CLIA Kit
E-CL-M0685	Mouse VF (Visfatin) CLIA Kit
E-CL-M0686	Mouse vWF (von Willebrand Factor) CLIA Kit
E-CL-M0687	Mouse α 1-AGP (α 1-Acid glycoprotein) CLIA Kit

Cat.NO	Product Name
E-CL-M0688	Mouse α 1-MG (α 1-Microglobulin) CLIA Kit
E-CL-M0689	Mouse α 2-PI (α 2-plasmin inhibitor) CLIA Kit
E-CL-M0690	Mouse β GAL (β -Galactosidase) CLIA Kit
E-CL-M0691	Mouse β -TG (β -Thromboglobulin) CLIA Kit
E-CL-M0692	Mouse ALT (Alanine Transaminase) CLIA Kit
E-CL-M0693	Mouse IL-21 (Interleukin 21) CLIA Kit
E-CL-M0694	Mouse α HBDH (α -Hydroxybutyrate Dehydrogenase) CLIA Kit
E-CL-M0695	Mouse GLUT1 (Glucose Transporter 1) CLIA Kit
E-CL-M0696	Mouse 5-LO (Arachidonate 5-Lipoxygenase) CLIA Kit
E-CL-M0697	Mouse NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) CLIA Kit
E-CL-M0698	Mouse BFP (Brain Finger Protein) CLIA Kit
E-CL-M0699	Mouse S100A8 (S100 Calcium Binding Protein A8) CLIA Kit
E-CL-M0700	Mouse S100A9 (S100 Calcium Binding Protein A9) CLIA Kit
E-CL-M0701	Mouse S100A10 (S100 Calcium-binding Protein A10) CLIA Kit
E-CL-M0702	Mouse PLA2G2D (Phospholipase A2, Group IID) CLIA Kit
E-CL-M0703	Mouse cPLA2 (Phospholipase A2, Cytosolic) CLIA Kit
E-CL-M0704	Mouse iPLA2 (Phospholipase A2, Calcium Independent) CLIA Kit
E-CL-M0705	Mouse SAA2 (Serum Amyloid A2) CLIA Kit
E-CL-M0706	Mouse AgRP (Agouti Related Protein) CLIA Kit

Cat.NO	Product Name
E-CL-M0707	Mouse PGF2 α (Prostaglandin F2 Alpha) CLIA Kit
E-CL-M0708	Mouse AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) CLIA Kit
E-CL-M0709	Mouse LDL (Low Density Lipoprotein) CLIA Kit
E-CL-M0710	Mouse FASN (Fatty Acid Synthase) CLIA Kit
E-CL-M0711	Mouse HDL (High Density Lipoprotein) CLIA Kit
E-CL-M0712	Mouse cTnT/TNNT2 (Troponin T Type 2, Cardiac) CLIA Kit
E-CL-M0713	Mouse PYY (Peptide YY) CLIA Kit
E-CL-M0714	Mouse SOD1 (Superoxide Dismutase 1, Soluble) CLIA Kit
E-CL-M0715	Mouse CYP2E1 (Cytochrome P450, family 2, subfamily E, polypeptide 1) CLIA Kit
E-CL-0001	CsA (Cyclosporine A) CLIA Kit
E-CL-0002	VA (Vitamin A) CLIA Kit
E-CL-0003	FA/VB9 (Folic Acid/Vitamin B9) CLIA Kit
E-CL-0004	VB12 (Vitamin B12) CLIA Kit
E-CL-0005	VD2 (Vitamin D2) CLIA Kit
E-CL-0006	VD3 (Vitamin D3) CLIA Kit
E-CL-0007	25-HVD3 (25-Hydroxy Vitamin D3) CLIA Kit
E-CL-0008	PGI2 (Prostacyclin) CLIA Kit
E-CL-0009	GSH (Glutathione) CLIA Kit
E-CL-0010	8-OHdG (8-Hydroxydeoxyguanosine) CLIA Kit

Cat.NO	Product Name
E-CL-0011	OT (Oxytocin) CLIA Kit
E-CL-0012	Cortisol CLIA Kit
E-CL-0013	DHT (Dihydrotestosterone) CLIA Kit
E-CL-0014	HIS (Histamine) CLIA Kit
E-CL-0015	ST/5-HT (Serotonin/5-Hydroxytryptamine) CLIA Kit
E-CL-0016	PGE2 (Prostaglandin E2) CLIA Kit
E-CL-0017	β-Nph (Beta-Naphthol) CLIA Kit
E-CL-0018	HA (Hyaluronic Acid) CLIA Kit
E-CL-0019	E3 (Estriol) CLIA Kit
E-CL-0020	F-E3 (Free Estriol)CLIA Kit
E-CL-0021	17-OHP (17-Hydroxyprogesterone) CLIA Kit