

www.elabscience.com

Elabscience®

Focus on lab research Service for life science



Catalog₂₀₁₅₋₂₀₁₆

for ELISA/CLIA KITS

About us

Elabscience is a high-tech biotech company, specializing in reagents of immunodiagnostic technology. It offers more than 20,000 products, including ELISA kits, CLIA kits, Antibodies, and Proteins.

Besides, Elabscience has established seven technical service platforms , including ELISA Customized Service, Immunoassay service, Molecular Biology Service, Protein Expression and Purification Service, Antibody Preparation Service, Cell Culture Service and Whole Research Project Service, to provide one-stop experimental service with professional knowledge, facilities and attitude.




Products

Elabscience has a portfolio of 20,000 products at present. Its website keeps updated by adding hundreds of products every month, so please keep close attention to our website to get your needed products.

Elabscience's products cover following categories:

- ▲ ELISA kits
- ▲ CLIA kits
- ▲ Monoclonal antibodies
- ▲ Polyclonal antibodies
- ▲ Proteins
- ▲ cDNA



Technical support

In Elabscience, we have experienced experts in the field of research and development, endeavoring to provide customers all-round technical support.

Any quality problems of products, we would replace the products or refund.

Any technical support in the process of operation, you can ask for help through Tel, Email, or livechat at the same time, you

can also loge in our technical forum to check the related articles of your experiment.



Tel: 86-27-87385095

Fax: 86-27- 87645690

Email: techsupport@elabscience.com


Technical Forum:

<http://www.elabscience.com/index.php/resources/show/tid/49>

ELISA principle


Enzyme-linked immunosorbent assay (ELISA), is a popular format of a "wet-lab" type analytic biochemistry assay that uses one sub-type of heterogeneous, solid-phase **enzyme immunoassay (EIA)** to detect the presence of a substance, usually an antigen, in a liquid sample or wet sample.

As a heterogenous assay, ELISA separates some component of the analytical reaction mixture by adsorbing certain components onto a solid phase which is physically immobilized. In ELISA, a liquid sample is added onto a stationary solid phase with special binding properties and is followed by multiple liquid reagents that are sequentially added, incubated and washed followed by some optical change (e.g. color development by the product of an enzymatic reaction) in the final liquid in the well from which the quantity of the analyte is measured. The sensitivity of detection depends on amplification of the signal during the analytic reactions. Since enzyme reactions are very well known amplification processes, the signal is generated by enzymes which are linked to the detection reagents in fixed proportions to allow accurate quantification - thus the name "enzyme linked".




ELISA kits

Two major types of ELISA formats we used



(1) Sandwich-ELISA



(2) Competitive-ELISA

Which type of ELISA method will be used is determined by the characteristics of the analyte itself and corresponding antibodies we could obtain.

When you order ELISA kits from Elabscience, we will provide manuals for you to get detailed information about the kits you wanted. If you didn't get the manuals, please contact with the local distributor or our market department directly.

Species covering:

- ▲ Human
- ▲ Rat
- ▲ Mouse
- ▲ Guinea pig
- ▲ Rabbit
- ▲ Porcine
- ▲ Canine
- ▲ Chicken
- ▲ Monkey
- ▲ Sheep



- Compatible with serum, plasma, cell supernatant and urine analysis and tissue homogenates
- Specifically react with both native and recombinant antigens

ELISA sample preparation

Serum: Allow samples to clot for 2 hours at room temperature or overnight at 4°C before centrifugation for 15 minutes at 1000×g. Collect the supernatant and carry out the assay immediately. Blood collection tubes should be disposable, non-pyrogenic, and non-endotoxin.

Plasma: Collect plasma using EDTA or heparin as an anticoagulant. Centrifuge samples for 15 minutes at 1000×g at 2 - 8°C within 30 minutes of collection. Collect the supernatant and carry out the assay immediately. Hemolysis samples are not suitable for ELISA assay!

Tissue homogenates: You'd better get detailed references from other literatures before assay aiming at different tissue types. For general information, hemolysis blood may affect the result, so you should mince the tissues to small pieces and rinse them in ice-cold PBS (0.02M, pH=7.4) to remove excess blood thoroughly. Tissue pieces should be weighed and then homogenized in PBS (the volume depends on the weight of the tissue) with a glass homogenizer on ice. To further break the cells, you can sonicate the suspension with an ultrasonic cell disrupter or subject it to freeze-thaw cycles. The homogenates are then centrifugated for 5 minutes at 5000×g to get the supernate.

Cell culture supernatant: Pipette cell culture media into a centrifuge tube and centrifuge at 1,200 rpm for 10 min at 4°C. Immediately aliquot supernatant and store samples at -80°C.

Cell extract: Place cell culture plates on ice. Aspirate medium and gently wash cells once with cold PBS. Aspirate PBS and add 1.0 ml complete extraction buffer per 100mm plate. Scrape cells to collect in tilted plate and remove to pre-chilled tube. Vortex briefly and incubate on ice for 5-30 min. Centrifuge at 12,000 x rpm for 10 min at 4°C to get rid of insoluble contents. Aliquot supernatant (this is the soluble cell extract) to clean, chilled tubes on ice and store samples at -80°C.

Milk: Collect samples and centrifuge at 1,000 x g for 10 min at 4°C.

Urine: Collect samples and centrifuge at 1,000 x g for 10 min at 4°C.

Note:

1. Samples should be used within 7 days when stored at 2-8°C , otherwise samples must be divided and stored at -20°C (≤1 month) or -80°C (≤6 months) to avoid the loss of bioactivity and contamination. Avoid repeated freeze-thaw cycles.

2. Please take the samples to room temperature (18-25°C) without extra heating before performing the assay.

3. Please predict the concentration before assaying. If the sample concentration is not within the range of the standard curve, users must determine the optimal sample dilutions for their particular experiments.

Troubleshooting for ELISA

Problem	Causes	Solutions
Poor standard curve	Inaccurate pipetting	Check pipettes
	Improper standard dilution	Ensure briefly spin the vial of standard and dissolve the powder thoroughly by a gentle mix.
	Wells not completely aspirated	Completely aspirate wells between steps.
Low signal	Too brief incubation times	Ensure sufficient incubation time;
	Incorrect assay temperature	Use recommended incubation temperature. Bring substrate to room temperature before use.
	Inadequate reagent volumes	Check pipettes and ensure correct preparation
	Improper dilution	
Deep color but low value	Plate reader settings not optimal	Verify the wavelength and filter setting in the plate reader.
		Open the Plate Reader ahead to pre-heat
Large CV	Inaccurate pipetting	Check pipettes
High background	Concentration of detector too high	Use recommended dilution factor.
	Plate is insufficiently washed	Review the manual for proper wash. If using a plate washer, check that all ports are unobstructed.
	Contaminated wash buffer	Make fresh wash buffer
Low sensitivity	Improper storage of the ELISA kit	All the reagents should be stored according to the instructions
	Stop solution not added	Stop solution should be added to each well before measurement

Products:

Cat.NO	Product Name
E-EL-H0001	Human α 1-AGP (Alpha-1-Acid Glycoprotein) ELISA Kit
E-EL-H0002	Human ACE (Angiotensin I Converting Enzyme) ELISA Kit
E-EL-H0003	Human ACVA (Activin A) ELISA Kit
E-EL-H0004	Human ADP/Acrp30 (Adiponectin) ELISA Kit
E-EL-H0005	Human AgRP (Agouti Related Protein) ELISA Kit
E-EL-H0006	Human ANG (Angiogenin) ELISA Kit
E-EL-H0007	Human ANG1 (Angiopoietin 1) ELISA Kit
E-EL-H0008	Human ANG2 (Angiopoietin 2) ELISA Kit
E-EL-H0009	Human BAFF/CD257 (B-Cell Activating Factor) ELISA Kit
E-EL-H0010	Human BDNF (Brain Derived Neurotrophic Factor) ELISA Kit
E-EL-H0011	Human BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit
E-EL-H0012	Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit
E-EL-H0013	Human BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit
E-EL-H0014	Human E-Cad (E-Cadherin) ELISA Kit
E-EL-H0015	Human CA9 (Carbonic Anhydrase IX) ELISA Kit
E-EL-H0016	Human CASP1 (Caspase 1) ELISA Kit
E-EL-H0017	Human CASP3 (Caspase 3) ELISA Kit
E-EL-H0018	Human CTSB (Cathepsin B) ELISA Kit

Cat.NO	Product Name
E-EL-H0019	Human CTSV (Cathepsin V) ELISA Kit
E-EL-H0020	Human MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit
E-EL-H0021	Human MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
E-EL-H0022	Human MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit
E-EL-H0023	Human RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
E-EL-H0024	Human MCP-3 (Monocyte Chemotactic Protein 3) ELISA Kit
E-EL-H0025	Human ECF/CCL11 (Eosinophil Chemotactic Factor) ELISA Kit
E-EL-H0026	Human TARC (Thymus Activation Regulated Chemokine) ELISA Kit
E-EL-H0027	Human MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
E-EL-H0028	Human SLC (Secondary Lymphoid Tissue Chemokine) ELISA Kit
E-EL-H0029	Human MDC (Macrophage Derived Chemokine) ELISA Kit
E-EL-H0030	Human MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
E-EL-H0031	Human MIP4 α (Macrophage Inflammatory Protein 4 Alpha) ELISA Kit
E-EL-H0032	Human CTACK (Cutaneous T-cell Attracting Chemokine) ELISA Kit
E-EL-H0033	Human sCD14 (Soluble Cluster of Differentiation14) ELISA Kit
E-EL-H0034	Human Fc ϵ RII/CD23 (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
E-EL-H0035	Human CD40L/TNFSF5 (Cluster of Differentiation 40 Ligand) ELISA Kit
E-EL-H0036	Human CD163 (Cluster of Differentiation 163) ELISA Kit
E-EL-H0037	Human GP39 (Glycoprotein 39, Cartilage) ELISA Kit
E-EL-H0038	Human CLU (Clusterin) ELISA Kit

Cat.NO	Product Name
E-EL-H0039	Human CNTF (Ciliary Neurotrophic Factor) ELISA Kit
E-EL-H0040	Human TF (Tissue Factor) ELISA Kit
E-EL-H0041	Human CFD (Complement Factor D) ELISA Kit
E-EL-H0042	Human CRN (Corin) ELISA Kit
E-EL-H0043	Human CRP (C-Reactive Protein) ELISA Kit
E-EL-H0044	Human CX3CL1 (Chemokine C-X3-C-Motif Ligand 1) ELISA Kit
E-EL-H0045	Human GRO α /CXCL1 (Growth Regulated Oncogene Alpha) ELISA Kit
E-EL-H0046	Human ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit
E-EL-H0047	Human GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit
E-EL-H0048	Human IL-8 (Interleukin 8) ELISA Kit
E-EL-H0049	Human AAHA (Anti-Apolipoprotein Antibody) ELISA Kit
E-EL-H0050	Human IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
E-EL-H0051	Human I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit
E-EL-H0052	Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit
E-EL-H0053	Human BLC (B-Lymphocyte Chemoattractant) ELISA Kit
E-EL-H0054	Human CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit
E-EL-H0055	Human Cys-C (Cystatin C) ELISA Kit
E-EL-H0056	Human Cyt-C (Cytochrome C) ELISA Kit
E-EL-H0057	Human DKK1 (Dickkopf Related Protein 1) ELISA Kit
E-EL-H0058	Human DPP4 (Dipeptidyl Peptidase IV) ELISA Kit

Cat.NO	Product Name
E-EL-H0059	Human EGF (Epidermal Growth Factor) ELISA Kit
E-EL-H0060	Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit
E-EL-H0061	Human EG-VEGF (Endocrine Gland Derived Vascular Endothelial Growth Factor)ELISA Kit
E-EL-H0062	Human ENG (Endoglin) ELISA Kit
E-EL-H0063	Human ES (Endostatin) ELISA Kit
E-EL-H0064	Human ET-1 (Endothelin 1) ELISA Kit
E-EL-H0065	Human EPCR (Endothelial Protein C Receptor) ELISA Kit
E-EL-H0066	Human EPO (Erythropoietin) ELISA Kit
E-EL-H0067	Human FAS/CD95 (Factor Related Apoptosis) ELISA Kit
E-EL-H0068	Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit
E-EL-H0069	Human FABP1 (Fatty Acid Binding Protein 1, Liver) ELISA Kit
E-EL-H0070	Human α FP (Alpha-Fetoprotein) ELISA Kit
E-EL-H0071	Human AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit
E-EL-H0072	Human AAP (Alanine Aminopeptidase) ELISA Kit
E-EL-H0073	Human FGF19 (Fibroblast Growth Factor 19) ELISA Kit
E-EL-H0074	Human FGF21 (Fibroblast Growth Factor 21) ELISA Kit
E-EL-H0075	Human Flt3L (FMS Like Tyrosine Kinase 3 Ligand) ELISA Kit
E-EL-H0076	Human FST (Follistatin) ELISA Kit
E-EL-H0077	Human Hcp (Hepcidin) ELISA Kit
E-EL-H0078	Human Gas6 (Growth Arrest Specific Protein 6) ELISA Kit

Cat.NO	Product Name
E-EL-H0079	Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit
E-EL-H0080	Human GDF15 (Growth Differentiation Factor 15) ELISA Kit
E-EL-H0081	Human GM-CSF (Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit
E-EL-H0082	Human sgp130 (soluble Glucoprotein 130) ELISA Kit
E-EL-H0084	Human HGF (Hepatocyte Growth Factor) ELISA Kit
E-EL-H0085	Human IFN- β (Interferon Beta)ELISA Kit
E-EL-H0086	Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit
E-EL-H0087	Human IGFBP-3 (Insulin Like Growth Factor Binding Protein 3) ELISA Kit
E-EL-H0088	Human IL-1 α (Interleukin 1 Alpha) ELISA Kit
E-EL-H0089	Human IL-1ra (Interleukin 1 Receptor Antagonist) ELISA Kit
E-EL-H0090	Human IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit
E-EL-H0091	Human PSA/KLK3 (Prostate Specific Antigen) ELISA Kit
E-EL-H0092	Human FGF7 (Fibroblast Growth Factor 7) ELISA Kit
E-EL-H0093	Human LEPR (Leptin Receptor) ELISA Kit
E-EL-H0094	Human LIF (Leukemia Inhibitory Factor) ELISA Kit
E-EL-H0095	Human TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14)ELISA Kit
E-EL-H0096	Human NGAL (Neutrophil Gelatinase Associated Lipocalin) ELISA Kit
E-EL-H0097	Human MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit
E-EL-H0098	Human Pro-MMP-1 (pro-Matrix Metalloproteinase 1) ELISA Kit
E-EL-H0099	Human IL-2 (Interleukin 2) ELISA Kit

Cat.NO	Product Name
E-EL-H0100	Human IL-3 (Interleukin 3) ELISA Kit
E-EL-H0101	Human IL-4 (Interleukin 4) ELISA Kit
E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit
E-EL-H0103	Human IL-10 (Interleukin 10) ELISA Kit
E-EL-H0104	Human IL-13(Interleukin 13)ELISA Kit
E-EL-H0105	Human IL-17 (Interleukin 17) ELISA Kit
E-EL-H0106	Human IL-22 (Interleukin 22) ELISA Kit
E-EL-H0107	Human IL-23 (Interleukin 23) ELISA Kit
E-EL-H0108	Human IFN- γ (Interferon Gamma) ELISA Kit
E-EL-H0109	Human TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-H0110	Human TGF- β 1 (Transforming Growth Factor Beta 1) ELISA Kit
E-EL-H0111	Human VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit
E-EL-H0112	Human EREG (Epiregulin) ELISA Kit
E-EL-H0113	Human LEP (Leptin) ELISA Kit
E-EL-H0114	Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit
E-EL-H0115	Human MAU (Microalbuminuria) ELISA Kit
E-EL-H0116	Human LH (Luteinizing Hormone) ELISA Kit
E-EL-H0119	Human REN (Renin) ELISA Kit
E-EL-H0120	Human CTSC (Cathepsin C) ELISA Kit
E-EL-H0121	Human C-P (C-Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-H0124	Human OxLDL (Oxidized Low Density Lipoprotein) ELISA Kit
E-EL-H0125	Human ApoA1 (Apolipoprotein A1) ELISA Kit
E-EL-H0128	Human LN (Laminin) ELISA Kit
E-EL-H0130	Human COL1 (Collagen Type I) ELISA Kit
E-EL-H0131	Human MMP-25 (Matrix Metalloproteinase 25) ELISA Kit
E-EL-H0132	Human MMP-10 (Matrix Metalloproteinase 10) ELISA Kit
E-EL-H0133	Human MMP-12 (Matrix Metalloproteinase 12) ELISA Kit
E-EL-H0134	Human MMP-13 (Matrix Metalloproteinase 13) ELISA Kit
E-EL-H0135	Human sICAM-1/CD54 (Soluble Intercellular Adhesion Molecule 1) ELISA Kit
E-EL-H0136	Human CT (Calcitonin) ELISA Kit
E-EL-H0137	Human ACTH (Adrenocorticotrophic Hormone) ELISA Kit
E-EL-H0139	Human ACAT2 (Acetyl Coenzyme A Acetyltransferase 2, cytosolic) ELISA Kit
E-EL-H0141	Human PRL (Prolactin) ELISA Kit
E-EL-H0142	Human MYO (Myoglobin) ELISA Kit
E-EL-H0143	Human Tn-T (Troponin T) ELISA Kit
E-EL-H0144	Human Tn-I (Troponin I) ELISA Kit
E-EL-H0145	Human PAPPA (Pregnancy Associated Plasma Protein A) ELISA Kit
E-EL-H0146	Human ALP (Alkaline Phosphatase) ELISA Kit
E-EL-H0147	Human CSTA (Cystatin A) ELISA Kit
E-EL-H0148	Human GLP-1 (Glucagon Like Peptide 1) ELISA Kit

Cat.NO	Product Name
E-EL-H0149	Human IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-H0150	Human IL-12 (Interleukin 12) ELISA Kit
E-EL-H0151	Human IL-12 p40 (Interleukin 12 p40) ELISA Kit
E-EL-H0152	Human CP (Ceruloplasmin) ELISA Kit
E-EL-H0154	Human α -CTx (Alpha Crosslaps) ELISA Kit
E-EL-H0155	Human CD42 (Cluster of Differentiation 42) ELISA Kit
E-EL-H0157	Human HLA-B27 (Human Leukocyte Antigen B27) ELISA Kit
E-EL-H0158	Human FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit
E-EL-H0159	Human IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) ELISA Kit
E-EL-H0160	Human LP-a (Lipoprotein a) ELISA Kit
E-EL-H0161	Human MBP (Myelin Basic Protein) ELISA Kit
E-EL-H0164	Human TE (Telomerase) ELISA Kit
E-EL-H0165	Human TFPI (Tissue Factor Pathway Inhibitor) ELISA Kit
E-EL-H0166	Human TM (Thrombomodulin) ELISA Kit
E-EL-H0167	Human C4 (Complement Component 4) ELISA Kit
E-EL-H0168	Human FE (Ferritin) ELISA Kit
E-EL-H0169	Human IgG (Immunoglobulin G) ELISA Kit
E-EL-H0171	Human TF (Transferrin) ELISA Kit
E-EL-H0173	Human FSTL1 (Follistatin Like Protein 1) ELISA Kit
E-EL-H0174	Human LOX (Lysyl Oxidase) ELISA Kit

Cat.NO	Product Name
E-EL-H0175	Human CG (Chorionic Gonadotropin) ELISA Kit
E-EL-H0177	Human GH (Growth Hormone) ELISA Kit
E-EL-H0178	Human COL4 (Collagen Type IV) ELISA Kit
E-EL-H0179	Human FN (Fibronectin) ELISA Kit
E-EL-H0180	Human IL-9 (Interleukin 9) ELISA Kit
E-EL-H0181	Human PCI (Procollagen I) ELISA Kit
E-EL-H0182	Human PCIII (Procollagen III) ELISA Kit
E-EL-H0183	Human PIIINP (N-Terminal Procollagen III Propeptide) ELISA Kit
E-EL-H0184	Human TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
E-EL-H0185	Human PINP(Procollagen I N-terminal Propeptide)ELISA Kit
E-EL-H0186	Human KIM-1 (Kidney Injury Molecule 1) ELISA Kit
E-EL-H0187	Human ADCY2 (Adenylate Cyclase 2) ELISA Kit
E-EL-H0188	Human AGGF1 (Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit
E-EL-H0189	Human C4a (Complement Component 4a) ELISA Kit
E-EL-H0190	Human C5a (Complement Component 5a) ELISA Kit
E-EL-H0191	Human IL-5 (Interleukin 5) ELISA Kit
E-EL-H0192	Human IL-6R (Interleukin 6 Receptor) ELISA Kit
E-EL-H0194	Human FcγR I (Receptor I for the Fc Region of Immunoglobulin G) ELISA Kit
E-EL-H0195	Human NCAD (Neural Cadherin) ELISA Kit
E-EL-H0196	Human PICP (Procollagen I C-terminal Propeptide) ELISA Kit

Cat.NO	Product Name
E-EL-H0197	Human PDGFR α (Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
E-EL-H0198	Human LUM (Lumican) ELISA Kit
E-EL-H0199	Human PRR4 (Proline Rich Protein 4, Lacrimal) ELISA Kit
E-EL-H0200	Human PSAP (Prosaposin) ELISA Kit
E-EL-H0201	Human PTGES (Prostaglandin E Synthase, Microsomal) ELISA Kit
E-EL-H0203	Human ISG15 (Interferon Stimulated Gene 15) ELISA Kit
E-EL-H0204	Human AIMP1 (Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1) ELISA Kit
E-EL-H0207	Human 11-DH-TXB2 (11-Dehydro-thromboxane B2) ELISA Kit
E-EL-H0208	Human 12-HETE (12-Hydroxyeicosatetraenoic Acid) ELISA Kit
E-EL-H0209	Human 15-LO (Arachidonate 15-Lipoxygenase) ELISA Kit
E-EL-H0210	Human 16- α OHE-1 (16- α Hydroxyestrogen 1) ELISA Kit
E-EL-H0214	Human BACE1 (Beta-site APP Cleaving Enzyme 1) ELISA Kit
E-EL-H0216	Human IL-32(Interleukin 32)ELISA Kit
E-EL-H0217	Human TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
E-EL-H0219	Human 26S-PSM (26S Proteasome) ELISA Kit
E-EL-H0221	Human 5-HIAA (5-Hydroxyindolacetic Acid) ELISA Kit
E-EL-H0222	Human IL-15 (Interleukin 15) ELISA Kit
E-EL-H0223	Human 5-LO (Arachidonate 5-Lipoxygenase) ELISA Kit
E-EL-H0224	Human 5-NT (5-Nucleotidase) ELISA Kit
E-EL-H0231	Human AAA (Anti-Actin Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-H0232	Human 5-HETE (5-Hydroxyeicosatetraenoic Acid) ELISA Kit
E-EL-H0233	Human AAA (Anti-Albumin Antibody) ELISA Kit
E-EL-H0234	Human α 1ACT (Alpha1 Antichymotrypsin) ELISA Kit
E-EL-H0235	Human IL-16 (Interleukin 16) ELISA Kit
E-EL-H0236	Human AP17 (Apelin 17) ELISA Kit
E-EL-H0237	Human AR (Amphiregulin) ELISA Kit
E-EL-H0239	Human INHA (Inhibin A) ELISA Kit
E-EL-H0240	Human ARO (Aromatase) ELISA Kit
E-EL-H0241	Human BLVRB (Biliverdin Reductase B) ELISA Kit
E-EL-H0242	Human ACA-IgA (Anti-Cardiolipin Antibody IgA) ELISA Kit
E-EL-H0243	Human ACA-IgG (Anti-Cardiolipin Antibody IgG) ELISA Kit
E-EL-H0244	Human ACA-IgM (Anti-Cardiolipin Antibody IgM) ELISA Kit
E-EL-H0245	Human ACD (Adrenocortical Dysplasia Homolog) ELISA Kit
E-EL-H0248	Human AChE (Acetylcholinesterase) ELISA Kit
E-EL-H0249	Human AChRab (Acetylcholine Receptor Antibody) ELISA Kit
E-EL-H0250	Human ACO-2 (Aconitase 2) ELISA Kit
E-EL-H0251	Human ACP (Acid Phosphatase) ELISA Kit
E-EL-H0252	Human ACP1 (Acid Phosphatase 1) ELISA Kit
E-EL-H0253	Human IL-18 (Interleukin 18) ELISA Kit
E-EL-H0254	Human IL-19 (Interleukin 19) ELISA Kit

Cat.NO	Product Name
E-EL-H0256	Human ACTG1 (Actin Gamma 1) ELISA Kit
E-EL-H0257	Human ACTN2 (Actinin Alpha 2) ELISA Kit
E-EL-H0258	Human ACTN3 (Actinin Alpha 3) ELISA Kit
E-EL-H0259	Human ACT α 2 (Actin Alpha 2, Smooth Muscle) ELISA Kit
E-EL-H0260	Human ACT α 1 (Actin Alpha 1, Skeletal Muscle) ELISA Kit
E-EL-H0262	Human ADA (Adenosine Deaminase) ELISA Kit
E-EL-H0263	Human ADAM10 (A Disintegrin And Metalloprotease 10) ELISA Kit
E-EL-H0264	Human ADAM8 (A Disintegrin And Metalloprotease 8) ELISA Kit
E-EL-H0265	Human ADAM9 (A Disintegrin And Metalloprotease 9) ELISA Kit
E-EL-H0266	Human ADAMTS4 (ADAM with Thrombospondin Type 1 Motif 4) ELISA Kit
E-EL-H0267	Human ADCY1 (Adenylate Cyclase 1, Brain) ELISA Kit
E-EL-H0271	Human ADH1 (Alcohol Dehydrogenase 1) ELISA Kit
E-EL-H0272	Human ADH (Antidiuretic Hormone) ELISA Kit
E-EL-H0273	Human ADIPOR1 (Adiponectin Receptor 1) ELISA Kit
E-EL-H0274	Human ADIPOR2 (Adiponectin Receptor 2) ELISA Kit
E-EL-H0275	Human ADM (Adrenomedullin) ELISA Kit
E-EL-H0277	Human ADRA1A (Adrenergic Receptor alpha-1A) ELISA Kit
E-EL-H0278	Human ADRP (Adipose Differentiation Related Protein) ELISA Kit
E-EL-H0279	Human IL-20 (Interleukin 20) ELISA Kit
E-EL-H0280	Human DLK-1 (Delta Like 1 Homolog) ELISA Kit

Cat.NO	Product Name
E-EL-H0281	Human ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit
E-EL-H0282	Human BMP-3 (Bone Morphogenetic Protein 3) ELISA Kit
E-EL-H0283	Human AFA (Anti-Filaggrin Antibody) ELISA Kit
E-EL-H0284	Human CAPS (Calcyphosine) ELISA Kit
E-EL-H0285	Human FABP4 (Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit
E-EL-H0286	Human AFAP1 (Actin Filament Associated Protein 1) ELISA Kit
E-EL-H0289	Human AFP-L3 (Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit
E-EL-H0290	Human FUCA (Alpha-L-Fucosidase, Tissue) ELISA Kit
E-EL-H0291	Human AGA (Anti-Gliadin Antibody) ELISA Kit
E-EL-H0293	Human CCNB (Cyclin-B) ELISA Kit
E-EL-H0294	Human AGC (Aggrecan) ELISA Kit
E-EL-H0295	Human AGER (Advanced Glycosylation End Product Specific Receptor) ELISA Kit
E-EL-H0297	Human CD3d (Cluster of Differentiation 3d) ELISA Kit
E-EL-H0298	Human AGR2 (Anterior Gradient Protein 2) ELISA Kit
E-EL-H0299	Human AGRN (Agrin) ELISA Kit
E-EL-H0300	Human AGT (Angiotensinogen) ELISA Kit
E-EL-H0302	Human AIF1 (Allograft Inflammatory Factor 1) ELISA Kit
E-EL-H0303	Human AIG1 (Androgen Induced Protein 1) ELISA Kit
E-EL-H0304	Human AKAP1 (A Kinase Anchor Protein 1) ELISA Kit
E-EL-H0305	Human ALCAM (Activated Leukocyte Cell Adhesion Molecule) ELISA Kit

Cat.NO	Product Name
E-EL-H0308	Human ALDH (Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit
E-EL-H0309	Human ALDOA (Aldolase A, Fructose Bisphosphate) ELISA Kit
E-EL-H0310	Human ALK (Anaplastic Lymphoma Kinase) ELISA Kit
E-EL-H0311	Human Alk-Smase (Alkaline Sphingomyelinase) ELISA Kit
E-EL-H0312	Human ALT (Alanine Transaminase) ELISA Kit
E-EL-H0313	Human INHB (Inhibin B) ELISA Kit
E-EL-H0314	Human CDK6 (Cyclin Dependent Kinase 6) ELISA Kit
E-EL-H0315	Human AMBP (Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit
E-EL-H0316	Human CDK8 (Cyclin Dependent Kinase 8) ELISA Kit
E-EL-H0317	Human AMH (Anti-Müllerian Hormone) ELISA Kit
E-EL-H0318	Human AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) ELISA Kit
E-EL-H0320	Human AMY1 (Amylase Alpha 1, Salivary) ELISA Kit
E-EL-H0322	Human IAPP (Islet Amyloid Polypeptide) ELISA Kit
E-EL-H0323	Human CDKN3 (Cyclin Dependent Kinase Inhibitor 3) ELISA Kit
E-EL-H0324	Human ANG (Angiostatin) ELISA Kit
E-EL-H0326	Human Ang-II (Angiotensin II) ELISA Kit
E-EL-H0327	Human Ang-I (Angiotensin I) ELISA Kit
E-EL-H0329	Human ANG II R1/AGTR1 (Angiotensin II Receptor 1) ELISA Kit
E-EL-H0330	Human ANG II R1-Ab (Angiotensin II Receptor 1 Antibody) ELISA Kit
E-EL-H0331	Human ANG II R2/AGTR2 (Angiotensin II Receptor 2) ELISA Kit

Cat.NO	Product Name
E-EL-H0333	Human ANGPT4 (Angiopoietin 4) ELISA Kit
E-EL-H0334	Human ANGPTL1 (Angiopoietin Like Protein 1) ELISA Kit
E-EL-H0335	Human ANGPTL2 (Angiopoietin Like Protein 2) ELISA Kit
E-EL-H0336	Human ANGPTL3 (Angiopoietin Like Protein 3) ELISA Kit
E-EL-H0337	Human ANGPTL4 (Angiopoietin Like Protein 4) ELISA Kit
E-EL-H0338	Human ANGPTL6 (Angiopoietin Like Protein 6) ELISA Kit
E-EL-H0339	Human ANG-R-Tie1 (Angiopoietin Receptor Tie1) ELISA Kit
E-EL-H0340	Human ANG-R-Tie2 (Angiopoietin Receptor Tie2) ELISA Kit
E-EL-H0341	Human Ank-B (Ankyrin B) ELISA Kit
E-EL-H0342	Human ANNA1(Anti-Neuronal Nuclear Antibody 1) ELISA Kit
E-EL-H0343	Human MAG (Myelin Associated Glycoprotein) ELISA Kit
E-EL-H0344	Human CDNF (Cerebral Dopamine Neurotrophic Factor) ELISA Kit
E-EL-H0345	Human FGFR1 (Fibroblast Growth Factor Receptor 1) ELISA Kit
E-EL-H0346	Human CNN1 (Calponin 1, Basic, Smooth Muscle) ELISA Kit
E-EL-H0347	Human Anti-AChR (Anti-Acetylcholine Receptor Antibody) ELISA Kit
E-EL-H0348	Human CNTNAP2 (Contactin Associated Protein Like Protein 2) ELISA Kit
E-EL-H0349	Human ACHAP (Acetylcholinesterase Associated Protein) ELISA Kit
E-EL-H0350	Human ACO1 (Aconitase 1) ELISA Kit
E-EL-H0351	Human CPE (Carboxypeptidase E) ELISA Kit
E-EL-H0352	Human Anti-GM1 (Anti-Ganglioside M1 Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-H0353	Human CG α (Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit
E-EL-H0354	Human Ang-III (Angiotensin III) ELISA Kit
E-EL-H0355	Human RARS (Arginyl tRNA Synthetase, cytoplasmic) ELISA Kit
E-EL-H0356	Human FGFR3 (Fibroblast Growth Factor Receptor 3) ELISA Kit
E-EL-H0357	Human FGFR4 (Fibroblast Growth Factor Receptor 4) ELISA Kit
E-EL-H0358	Human DPP10 (Dipeptidyl Peptidase X) ELISA Kit
E-EL-H0359	Human CPR (Cytochrome P450 Reductase) ELISA Kit
E-EL-H0360	Human CRLF1 (Cytokine Receptor Like Factor 1) ELISA Kit
E-EL-H0361	Human CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit
E-EL-H0362	Human CA3 (Carbonic Anhydrase III, Muscle Specific) ELISA Kit
E-EL-H0363	Human AIRA(Anti-Insulin Receptor Antibody)ELISA Kit
E-EL-H0364	Human CYP1A1 (Cytochrome P450, family 1, subfamily A, polypeptide 1) ELISA Kit
E-EL-H0365	Human CYP1A2 (Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit
E-EL-H0366	Human AKA (Anti-Keratin Antibody) ELISA Kit
E-EL-H0367	Human DCLK1 (Doublecortin Like Kinase 1) ELISA Kit
E-EL-H0368	Human DPYD (Dihydropyrimidine Dehydrogenase) ELISA Kit
E-EL-H0369	Human DSE (Dermatan Sulfate Epimerase) ELISA Kit
E-EL-H0370	Human EPYC (Epiphycan) ELISA Kit
E-EL-H0371	Human FABP8 (Fatty Acid Binding Protein 8, Myelin) ELISA Kit
E-EL-H0372	Human FADS2 (Fatty Acid Desaturase 2) ELISA Kit

Cat.NO	Product Name
E-EL-H0374	Human FAIM3 (Fas Apoptotic Inhibitory Molecule 3) ELISA Kit
E-EL-H0375	Human FASN(Fatty Acid Synthase)ELISA Kit
E-EL-H0376	Human FBLN4 (Fibulin 4) ELISA Kit
E-EL-H0377	Human CYP24A1 (Cytochrome P450, family 24, subfamily A, polypeptide 1) ELISA Kit
E-EL-H0378	Human FCN1 (Ficolin 1) ELISA Kit
E-EL-H0379	Human AMA IgA (Anti-Myelin antibody IgA) ELISA Kit
E-EL-H0380	Human MAG Ab (Anti-Myelin Associated Glycoprotein Antibody) ELISA Kit
E-EL-H0381	Human FcγBP (Fc Fragment of IgG Binding Protein) ELISA Kit
E-EL-H0382	Human FcγR3A (Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit
E-EL-H0383	Human ANEA (Anti-Neutrophil Elastase Antibody) ELISA Kit
E-EL-H0384	Human CYP1B1 (Cytochrome P450, family 1, subfamily B, polypeptide 1) ELISA Kit
E-EL-H0385	Human FcεR II (Receptor II for the Fc Fragment of IgE) ELISA Kit
E-EL-H0386	Human FETUA (Fetuin A) ELISA Kit
E-EL-H0387	Human FGγ (Fibrinogen Gamma chain) ELISA Kit
E-EL-H0388	Human FLOT2 (Flotillin 2) ELISA Kit
E-EL-H0389	Human GC (Glycocalicin) ELISA Kit
E-EL-H0390	Human OJ/IleRS (Anti-OJ-Antibody) ELISA Kit
E-EL-H0391	Human GLRX3 (Glutaredoxin 3) ELISA Kit
E-EL-H0392	Human GIF (Gastric Intrinsic Factor) ELISA Kit
E-EL-H0393	Human AP62A (Anti-P62 Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-H0394	Human C4BP β (C4 Binding Protein Beta) ELISA Kit
E-EL-H0395	Human CA4 (Carbonic Anhydrase IV) ELISA Kit
E-EL-H0396	Human C4BP α (C4 Binding Protein Alpha) ELISA Kit
E-EL-H0397	Human GPLD1 (Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Kit
E-EL-H0398	Human PRELP (Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit
E-EL-H0399	Human APSA (Anti-Phosphatidyl Serine Antibody) ELISA Kit
E-EL-H0400	Human APA (Anti-Phospholipid Antibody) ELISA Kit
E-EL-H0401	Human Anti-AlaRS (Anti-Alanyl-tRNA Synthetase/Anti-PL12-Antibody) ELISA Kit
E-EL-H0403	Human TARS (Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit
E-EL-H0404	Human GTF2A1 (General Transcription Factor IIA, Polypeptide 1) ELISA Kit
E-EL-H0405	Human GUCA2A (Guanylate Cyclase Activator 2A) ELISA Kit
E-EL-H0406	Human INSL3 (Insulin Like Protein 3) ELISA Kit
E-EL-H0407	Human INSL5 (Insulin Like Protein 5) ELISA Kit
E-EL-H0409	Human iPLA2 (Phospholipidase A2, Calcium Independent) ELISA Kit
E-EL-H0410	Human GUCA2B (Guanylate Cyclase Activator 2B) ELISA Kit
E-EL-H0411	Human GZMK (Granzyme K) ELISA Kit
E-EL-H0412	Human GUK1 (Guanylate Kinase 1) ELISA Kit
E-EL-H0413	Human ALOXE3 (Arachidonate Lipoxygenase 3) ELISA Kit
E-EL-H0414	Human CYP26A1 (Cytochrome P450, family 26, subfamily A, polypeptide 1) ELISA Kit
E-EL-H0415	Human HB (Hemoglobin) ELISA Kit

Cat.NO	Product Name
E-EL-H0416	Human HCCS (Holocytochrome C Synthase) ELISA Kit
E-EL-H0417	Human CAMK2 α (Calcium/Calmodulin-Dependent Protein Kinase II Alpha) ELISA Kit
E-EL-H0418	Human BAG3 (Bcl2 Associated Athanogene 3) ELISA Kit
E-EL-H0419	Human ARSB (Arylsulfatase B) ELISA Kit
E-EL-H0420	Human HRG (Histidine Rich Glycoprotein) ELISA Kit
E-EL-H0421	Human HSPG (Heparan Sulfate Proteoglycan) ELISA Kit
E-EL-H0422	Human ANXA5 (Annexin A5) ELISA Kit
E-EL-H0423	Human HTRA1 (HtrA Serine Peptidase 1) ELISA Kit
E-EL-H0424	Human IL10R β (Interleukin 10 Receptor Beta) ELISA Kit
E-EL-H0425	Human IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit
E-EL-H0426	Human HYOU1 (Hypoxia Up Regulated 1) ELISA Kit
E-EL-H0427	Human IAA (Insulin Autoantibody) ELISA Kit
E-EL-H0428	Human ICAM4 (Intercellular Adhesion Molecule 4) ELISA Kit
E-EL-H0429	Human IDE (Insulin Degrading Enzyme) ELISA Kit
E-EL-H0430	Human Anti-Surv (Anti-Survivin Antibody) ELISA Kit
E-EL-H0431	Human TA (Anti-Teichoic Acid Antibody) ELISA Kit
E-EL-H0432	Human AT (AntiThrombin) ELISA Kit
E-EL-H0433	Human LIPH (Lipase H) ELISA Kit
E-EL-H0434	Human LT β R (Lymphotoxin Beta Receptor) ELISA Kit
E-EL-H0435	Human IGF-2 (Insulin Like Growth Factor 2) ELISA Kit

Cat.NO	Product Name
E-EL-H0436	Human ATGA/TGAB (Anti-Thyroid-Globulin Antibody) ELISA Kit
E-EL-H0437	Human TPO-Ab (Anti-Thyroid-Peroxidase Antibody) ELISA Kit
E-EL-H0438	Human IGF2BP3 (Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit
E-EL-H0439	Human TTN (Titin) ELISA Kit
E-EL-H0440	Human IGF2R (Insulin Like Growth Factor 2 Receptor) ELISA Kit
E-EL-H0441	Human IGFALS (Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) ELISA Kit
E-EL-H0442	Human IGFBP-1 (Insulin Like Growth Factor Binding Protein 1) ELISA Kit
E-EL-H0443	Human IGFBP-5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit
E-EL-H0444	Human IGFBP-4 (Insulin Like Growth Factor Binding Protein 4) ELISA Kit
E-EL-H0445	Human Anti-CIIAb (Anti-Type II Collagen Antibody) ELISA Kit
E-EL-H0446	Human IGFBP-2 (Insulin Like Growth Factor Binding Protein 2) ELISA Kit
E-EL-H0447	Human IGFBP-7 (Insulin Like Growth Factor Binding Protein 7) ELISA Kit
E-EL-H0448	Human ANXA2 (Annexin A2) ELISA Kit
E-EL-H0450	Human IGFBP-6 (Insulin Like Growth Factor Binding Protein 6) ELISA Kit
E-EL-H0452	Human INSR (Insulin Receptor) ELISA Kit
E-EL-H0453	Human APC (Activated Protein C) ELISA Kit
E-EL-H0456	Human APLN (Apelin) ELISA Kit
E-EL-H0458	Human AP13 (Apelin 13) ELISA Kit
E-EL-H0459	Human AP36 (Apelin 36) ELISA Kit
E-EL-H0460	Human Apo A1-Ab (Anti-Apolipoprotein A1 Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-H0461	Human LCN1(Lipocalin 1)ELISA Kit
E-EL-H0462	Human ApoA2 (Apolipoprotein A2) ELISA Kit
E-EL-H0463	Human ApoA4 (Apolipoprotein A4) ELISA Kit
E-EL-H0464	Human ApoB (Apolipoprotein B) ELISA Kit
E-EL-H0465	Human ApoC1 (Apolipoprotein C1) ELISA Kit
E-EL-H0466	Human ApoC2 (Apolipoprotein C2) ELISA Kit
E-EL-H0467	Human ApoC3 (Apolipoprotein C3) ELISA Kit
E-EL-H0468	Human ApoC4 (Apolipoprotein C4) ELISA Kit
E-EL-H0469	Human ApoD (Apolipoprotein D) ELISA Kit
E-EL-H0470	Human ApoE (Apolipoprotein E) ELISA Kit
E-EL-H0471	Human ApoF (Apolipoprotein F) ELISA Kit
E-EL-H0472	Human ApoH (Apolipoprotein H) ELISA Kit
E-EL-H0473	Human ApoM (Apolipoprotein M) ELISA Kit
E-EL-H0474	Human CA6 (Carbonic Anhydrase VI) ELISA Kit
E-EL-H0475	Human ApoA5 (Apolipoprotein A5) ELISA Kit
E-EL-H0476	Human ApoB100 (Apolipoprotein B100) ELISA Kit
E-EL-H0477	Human AATF (Apoptosis Antagonizing Transcription Factor) ELISA Kit
E-EL-H0478	Human LIPC (Lipase, Hepatic) ELISA Kit
E-EL-H0479	Human AIF (Apoptosis Inducing Factor) ELISA Kit
E-EL-H0480	Human LIPF (Lipase, Gastric) ELISA Kit

Cat.NO	Product Name
E-EL-H0481	Human APAF1 (Apoptosis Protease Activating Factor 1) ELISA Kit
E-EL-H0482	Human ASK-1 (Apoptosis Signal Regulating Kinase 1) ELISA Kit
E-EL-H0483	Human bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit
E-EL-H0484	Human AP (Aprotinin) ELISA Kit
E-EL-H0485	Human APT (Abnormal Prothrombin) ELISA Kit
E-EL-H0486	Human ITIH4 (Inter Alpha-Globulin Inhibitor H4) ELISA Kit
E-EL-H0487	Human AQP-1 (Aquaporin 1) ELISA Kit
E-EL-H0488	Human AQP-2 (Aquaporin 2) ELISA Kit
E-EL-H0489	Human AQP-3 (Aquaporin 3) ELISA Kit
E-EL-H0490	Human AQP-4 (Aquaporin 4) ELISA Kit
E-EL-H0493	Human AR (Androgen Receptor) ELISA Kit
E-EL-H0497	Human ARG1 (Arginase I) ELISA Kit
E-EL-H0498	Human AVPR1A (Arginine Vasopressin Receptor 1A) ELISA Kit
E-EL-H0499	Human KLK6 (Kallikrein 6) ELISA Kit
E-EL-H0500	Human KLK7 (Kallikrein 7) ELISA Kit
E-EL-H0501	Human ARHGDI α (Rho GDP Dissociation Inhibitor Alpha) ELISA Kit
E-EL-H0502	Human ARR β 1 (Arrestin Beta 1) ELISA Kit
E-EL-H0503	Human ARR β 2 (Arrestin Beta 2) ELISA Kit
E-EL-H0504	Human MANF (Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit
E-EL-H0505	Human MAP2K1 (Mitogen Activated Protein Kinase Kinase 1) ELISA Kit

Cat.NO	Product Name
E-EL-H0506	Human AhR (Aryl Hydrocarbon Receptor) ELISA Kit
E-EL-H0507	Human ARSA (Arylsulfatase A) ELISA Kit
E-EL-H0509	Human ASGR2 (Asialoglycoprotein Receptor 2) ELISA Kit
E-EL-H0510	Human ASGR1 (Asialoglycoprotein Receptor 1) ELISA Kit
E-EL-H0511	Human ASM (Acid Sphingomyelinase) ELISA Kit
E-EL-H0513	Human MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit
E-EL-H0514	Human AST (Aspartate Aminotransferase) ELISA Kit
E-EL-H0515	Human ASPN (Asporin) ELISA Kit
E-EL-H0517	Human ATA (Angiotensinase A) ELISA Kit
E-EL-H0518	Human ATM (Ataxia Telangiectasia Mutated) ELISA Kit
E-EL-H0519	Human ATXN1 (Ataxin 1) ELISA Kit
E-EL-H0520	Human ATXN2 (Ataxin 2) ELISA Kit
E-EL-H0521	Human ATF1 (Activating Transcription Factor 1) ELISA Kit
E-EL-H0522	Human ATF2 (Activating Transcription Factor 2) ELISA Kit
E-EL-H0523	Human ATF3 (Activating Transcription Factor 3) ELISA Kit
E-EL-H0524	Human ATF4 (Activating Transcription Factor 4) ELISA Kit
E-EL-H0525	Human ATF5 (Activating Transcription Factor 5) ELISA Kit
E-EL-H0526	Human ATF6 (Activating Transcription Factor 6) ELISA Kit
E-EL-H0527	Human ATF7 (Activating Transcription Factor 7) ELISA Kit
E-EL-H0528	Human ABCG1 (ATP Binding Cassette Transporter G1) ELISA Kit

Cat.NO	Product Name
E-EL-H0529	Human ABCA1(ATP Binding Cassette Transporter A1) ELISA Kit
E-EL-H0530	Human ABCG2 (ATP Binding Cassette Transporter G2) ELISA Kit
E-EL-H0532	Human ANP (Atrial Natriuretic Peptide) ELISA Kit
E-EL-H0533	Human bFGF-2 (Basic Fibroblast Growth Factor 2) ELISA Kit
E-EL-H0535	Human MGL (Monoacylglycerol Lipase) ELISA Kit
E-EL-H0536	Human AIRE (Autoimmune Regulator) ELISA Kit
E-EL-H0538	Human AXIN2 (Axis Inhibition Protein 2) ELISA Kit
E-EL-H0540	Human AZU1 (Azurocidin 1) ELISA Kit
E-EL-H0541	Human MCARS1 (Microspherule Protein 1) ELISA Kit
E-EL-H0542	Human A β 40 (Amyloid Beta 40) ELISA Kit
E-EL-H0543	Human A β 42 (Amyloid Beta 42) ELISA Kit
E-EL-H0544	Human BAFFR (B Cell Activating Factor Receptor) ELISA Kit
E-EL-H0545	Human BCDF (B-Cell Differentiation Factor) ELISA Kit
E-EL-H0547	Human BCR (B-Cell Receptor) ELISA Kit
E-EL-H0548	Human BlysR (B Lymphocyte Stimulator Receptor) ELISA Kit
E-EL-H0549	Human BCL2L1/BCL-X (Bcl-2 Like Protein 1) ELISA Kit
E-EL-H0550	Human BPI (Bactericidal/Permeability Increasing Protein) ELISA Kit
E-EL-H0552	Human FGF4 (Fibroblast Growth Factor 4) ELISA Kit
E-EL-H0553	Human FGF6 (Fibroblast Growth Factor 6) ELISA Kit
E-EL-H0556	Human LDHA (Lactate Dehydrogenase A) ELISA Kit

Cat.NO	Product Name
E-EL-H0557	Human Bcl-10 (B-cell Leukemia/Lymphoma 10) ELISA Kit
E-EL-H0558	Human Bcl-3 (B-cell Leukemia/Lymphoma 3) ELISA Kit
E-EL-H0559	Human BLNK (B-Cell Linker Protein) ELISA Kit
E-EL-H0560	Human BAK1 (Bcl2 Antagonist/Killer 1) ELISA Kit
E-EL-H0561	Human BAD (Bcl2 Associated Death Promoter) ELISA Kit
E-EL-H0562	Human Bax (Bcl-2 Associated X Protein) ELISA Kit
E-EL-H0563	Human BMF (Bcl-2 Modifying Factor) ELISA Kit
E-EL-H0564	Human BECN1 (Beclin 1)ELISA Kit
E-EL-H0565	Human SEMA3B (Semaphorin 3B) ELISA Kit
E-EL-H0568	Human β 2-GP (Beta2 Glycoprotein) ELISA Kit
E-EL-H0570	Human bTC (Betacellulin) ELISA Kit
E-EL-H0571	Human β -EPR (Beta-Endorphin Receptor) ELISA Kit
E-EL-H0572	Human β -EP (Beta-Endorphin) ELISA Kit
E-EL-H0575	Human BACE2 (β -site APP Cleaving Enzyme 2) ELISA Kit
E-EL-H0577	Human BID (BH3 Interacting Domain Death Agonist) ELISA Kit
E-EL-H0578	Human MMRN1 (Multimerin 1) ELISA Kit
E-EL-H0579	Human BTA (Bladder Tumor Antigen) ELISA Kit
E-EL-H0581	Human MFAP4 (Microfibrillar Associated Protein 4) ELISA Kit
E-EL-H0582	Human Bcl2L11 (Bcl-2 Like Protein 11) ELISA Kit
E-EL-H0584	Human BALP (Bone Alkaline Phosphatase) ELISA Kit

Cat.NO	Product Name
E-EL-H0585	Human BMP-1 (Bone Morphogenetic Protein 1) ELISA Kit
E-EL-H0586	Human BMP-10 (Bone Morphogenetic Protein 10) ELISA Kit
E-EL-H0587	Human BMP-15 (Bone Morphogenetic Protein 15) ELISA Kit
E-EL-H0588	Human BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit
E-EL-H0589	Human BMPR1A (Bone Morphogenetic Protein Receptor, type IA) ELISA Kit
E-EL-H0590	Human BMPR2 (Bone Morphogenetic Protein Receptor II) ELISA Kit
E-EL-H0592	Human BSP (Bone Sialoprotein) ELISA Kit
E-EL-H0594	Human MOG (Myelin Oligodendrocyte Glycoprotein) ELISA Kit
E-EL-H0595	Human MT1E (Metallothionein 1E) ELISA Kit
E-EL-H0596	Human BK (Bradykinin) ELISA Kit
E-EL-H0597	Human BFP (Brain Finger Protein) ELISA Kit
E-EL-H0598	Human BNP (Brain Natriuretic Peptide) ELISA Kit
E-EL-H0599	Human MSM β (Microseminoprotein Beta) ELISA Kit
E-EL-H0600	Human BRAK/CXCL14 (Breast and Kidney Expressed Chemokine) ELISA Kit
E-EL-H0601	Human BRCA1 (Breast Cancer Susceptibility Protein 1) ELISA Kit
E-EL-H0602	Human BRCA2 (Breast Cancer Susceptibility Protein 2) ELISA Kit
E-EL-H0603	Human BCAN (Brevican) ELISA Kit
E-EL-H0605	Human MT1M (Metallothionein 1M) ELISA Kit
E-EL-H0607	Human MT3 (Metallothionein 3) ELISA Kit
E-EL-H0608	Human CNP (C-type Natriuretic Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-H0609	Human C1q-Ab (Anti-Complement 1q Antibody) ELISA Kit
E-EL-H0610	Human C4BP (C4 Binding Protein) ELISA Kit
E-EL-H0612	Human MUC13 (Mucin 13) ELISA Kit
E-EL-H0613	Human CA724 (Tumor Marker CA724) ELISA Kit
E-EL-H0615	Human CNR-1 (Cadherin Related Neuronal Receptor 1) ELISA Kit
E-EL-H0616	Human MUC1 (Mucin 1) ELISA Kit
E-EL-H0617	Human CABIN1 (Calcineurin Binding Protein 1) ELISA Kit
E-EL-H0618	Human CaN (Calcineurin) ELISA Kit
E-EL-H0619	Human CGRP (Calcitonin Gene Related Peptide) ELISA Kit
E-EL-H0620	Human CTR (Calcitonin Receptor) ELISA Kit
E-EL-H0621	Human CaSR (Calcium Sensing Receptor) ELISA Kit
E-EL-H0622	Human CACYBP (Calcyclin Binding Protein) ELISA Kit
E-EL-H0623	Human CALD (Caldesmon) ELISA Kit
E-EL-H0624	Human MSR1 (Macrophage Scavenger Receptor 1) ELISA Kit
E-EL-H0625	Human CAM (Calmodulin) ELISA Kit
E-EL-H0626	Human CNX (Calnexin) ELISA Kit
E-EL-H0627	Human CRT (Calreticulin) ELISA Kit
E-EL-H0628	Human CR (Calretinin) ELISA Kit
E-EL-H0629	Human CALU (Calumenin) ELISA Kit
E-EL-H0630	Human CAM-Ab (Anti-Calmodulin Specific Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-H0631	Human MUC17 (Mucin 17) ELISA Kit
E-EL-H0632	Human MUC2 (Mucin 2) ELISA Kit
E-EL-H0633	Human MUC20 (Mucin 20) ELISA Kit
E-EL-H0634	Human CNR1 (Cannabinoid Receptor 1, Brain) ELISA Kit
E-EL-H0635	Human CAP1 (Adenylate Cyclase Associated Protein 1) ELISA Kit
E-EL-H0636	Human CA125 (Carbohydrate Antigen 125) ELISA Kit
E-EL-H0637	Human CA19-9 (Carbohydrate antigen19-9) ELISA Kit
E-EL-H0638	Human CDT (Carbohydrate Deficient Transferrin) ELISA Kit
E-EL-H0639	Human CA1 (Carbonic Anhydrase I) ELISA Kit
E-EL-H0640	Human CA2 (Carbonic Anhydrase II) ELISA Kit
E-EL-H0642	Human HbCO (Carboxyhemoglobin) ELISA Kit
E-EL-H0643	Human CEA (Carcinoembryonic Antigen) ELISA Kit
E-EL-H0645	Human ANKRD1 (Ankyrin Repeat Domain 1) ELISA Kit
E-EL-H0646	Human cTnT/TNNT2 (Troponin T Type 2, Cardiac) ELISA Kit
E-EL-H0647	Human CMLC-1 (Cardiac Myosin Light Chains 1) ELISA Kit
E-EL-H0648	Human IL-7 (Interleukin 7) ELISA Kit
E-EL-H0649	Human cTn-I/TNNI3 (cardiac Troponin I)ELISA Kit
E-EL-H0650	Human CLS (Cardiolipin Synthase) ELISA Kit
E-EL-H0651	Human CT-1 (Cardiotrophin 1) ELISA Kit
E-EL-H0652	Human CLCF1 (Cardiotrophin Like Cytokine Factor 1) ELISA Kit

Cat.NO	Product Name
E-EL-H0653	Human PAR1 (Protease Activated Receptor 1) ELISA Kit
E-EL-H0654	Human COMP (Cartilage Oligomeric Matrix Protein) ELISA Kit
E-EL-H0655	Human CASP10 (Caspase 10) ELISA Kit
E-EL-H0656	Human CASP11 (Caspase 11) ELISA Kit
E-EL-H0657	Human CASP12 (Caspase 12) ELISA Kit
E-EL-H0658	Human CASP14 (Caspase 14) ELISA Kit
E-EL-H0659	Human CASP8 (Caspase 8) ELISA Kit
E-EL-H0660	Human CASP4 (Caspase 4) ELISA Kit
E-EL-H0661	Human CASP7 (Caspase 7) ELISA Kit
E-EL-H0662	Human CAD (Caspase Activated Deoxyribonuclease) ELISA Kit
E-EL-H0663	Human CASP 9 (Caspase 9) ELISA Kit
E-EL-H0665	Human CTNN α 1 (Catenin, Alpha 1) ELISA Kit
E-EL-H0666	Human CTNN β 1 (Catenin, Beta 1) ELISA Kit
E-EL-H0669	Human CTSD (Cathepsin D) ELISA Kit
E-EL-H0670	Human CTSK (Cathepsin K) ELISA Kit
E-EL-H0671	Human CTSL (Cathepsin L) ELISA Kit
E-EL-H0672	Human CDX2 (Caudal Type Homeobox Transcription Factor 2) ELISA Kit
E-EL-H0673	Human Cav-1 (Caveolin-1) ELISA Kit
E-EL-H0674	Human C/EBP α (CCAAT/Enhancer Binding Protein Alpha) ELISA Kit
E-EL-H0675	Human C/EBP β (CCAAT/Enhancer Binding Protein Beta) ELISA Kit

Cat.NO	Product Name
E-EL-H0676	Human C/EBPδ (CCAAT/Enhancer Binding Protein, Delta) ELISA Kit
E-EL-H0677	Human C/EBPε (CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit
E-EL-H0678	Human CEBPγ (CCAAT/Enhancer Binding Protein Gamma) ELISA Kit
E-EL-H0679	Human CCR1 (CC-Chemokine Receptor 1) ELISA Kit
E-EL-H0682	Human CD59 (Protectin) ELISA Kit
E-EL-H0685	Human CDC23 (Cell Division Cycle Protein 23) ELISA Kit
E-EL-H0687	Human CENPA (Centromere Protein A) ELISA Kit
E-EL-H0688	Human CENPE (Centromere Protein E) ELISA Kit
E-EL-H0689	Human CENPH (Centromere Protein H) ELISA Kit
E-EL-H0690	Human CENPI (Centromere Protein I) ELISA Kit
E-EL-H0691	Human CENPB (Centromere Protein B) ELISA Kit
E-EL-H0692	Human CEP110 (Centrosomal Protein 110kDa) ELISA Kit
E-EL-H0693	Human CA5B (Carbonic Anhydrase VB, Mitochondrial) ELISA Kit
E-EL-H0694	Human CFL1(Cofilin 1, Non-Muscle) ELISA Kit
E-EL-H0695	Human V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog (FOS) ELISA Kit
E-EL-H0696	Human MUC21 (Mucin 21) ELISA Kit
E-EL-H0697	Human CPN10 (Chaperonin 10) ELISA Kit
E-EL-H0698	Human CHEM(Chemerin) ELISA Kit
E-EL-H0699	Human CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit
E-EL-H0700	Human CCL14 (Chemokine C-C-Motif Ligand 14) ELISA Kit

Cat.NO	Product Name
E-EL-H0701	Human CCL16 (Chemokine C-C-Motif Ligand 16) ELISA Kit
E-EL-H0702	Human CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit
E-EL-H0703	Human CCL4L1 (Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit
E-EL-H0704	Human CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit
E-EL-H0705	Human CCRL1 (Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
E-EL-H0706	Human TERT (Telomerase Reverse Transcriptase) ELISA Kit
E-EL-H0707	Human CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
E-EL-H0708	Human MUC3 (Mucin 3) ELISA Kit
E-EL-H0709	Human CKLFSF4 (Chemokine Like Factor Superfamily 4) ELISA Kit
E-EL-H0710	Human CKLFSF5 (Chemokine Like Factor Superfamily 5) ELISA Kit
E-EL-H0711	Human CKLFSF6 (Chemokine Like Factor Superfamily 6) ELISA Kit
E-EL-H0712	Human CKLFSF7 (Chemokine Like Factor Superfamily 7) ELISA Kit
E-EL-H0713	Human CKLFSF8 (Chemokine Like Factor Superfamily 8) ELISA Kit
E-EL-H0714	Human CKLF (Chemokine Like Factor) ELISA Kit
E-EL-H0715	Human MUC4 (Mucin 4) ELISA Kit
E-EL-H0716	Human MUC6 (Mucin 6) ELISA Kit
E-EL-H0717	Human MYBPC3 (Myosin Binding Protein C, Cardiac) ELISA Kit
E-EL-H0718	Human MYH10 (Myosin Heavy Chain 10, Non Muscle) ELISA Kit
E-EL-H0719	Human MYH9 (Myosin Heavy Chain 9, Non Muscle) ELISA Kit
E-EL-H0720	Human MYL9 (Myosin Light Chain 9, Regulatory) ELISA Kit

Cat.NO	Product Name
E-EL-H0721	Human CCKAR (Cholecystokinin A Receptor) ELISA Kit
E-EL-H0722	Human CCKBR (Cholecystokinin B Receptor) ELISA Kit
E-EL-H0723	Human CCK-8 (Cholecystokinin 8) ELISA Kit
E-EL-H0725	Human CETP (Cholesteryl Ester Transfer Protein) ELISA Kit
E-EL-H0727	Human CENPF (Centromere Protein F) ELISA Kit
E-EL-H0729	Human CHAc (Choline acetylase) ELISA Kit
E-EL-H0730	Human PC/CPG (Choline Phosphoglyceride) ELISA Kit
E-EL-H0731	Human CHRM2 (Cholinergic Receptor, Muscarinic 2) ELISA Kit
E-EL-H0733	Human CS (Chondroitin Sulfate) ELISA Kit
E-EL-H0734	Human CG β (Chorionic Gonadotrophin Beta) ELISA Kit
E-EL-H0735	Human CA12 (Carbonic Anhydrase XII) ELISA Kit
E-EL-H0736	Human CBX3 (Chromobox Homolog 3) ELISA Kit
E-EL-H0737	Human CHD3 (Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
E-EL-H0738	Human CHD5 (Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
E-EL-H0739	Human CgA (Chromogranin A/Pancreastatin) ELISA Kit
E-EL-H0740	Human NEF3 (Neurofilament 3) ELISA Kit
E-EL-H0741	Human NEFL (Neurofilament, Light Polypeptide) ELISA Kit
E-EL-H0742	Human NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) ELISA Kit
E-EL-H0743	Human C-jun (V-Jun Sarcoma Virus 17 Oncogene Homolog) ELISA Kit
E-EL-H0744	Human CC16 (Clara Cell 16kD protein) ELISA Kit

Cat.NO	Product Name
E-EL-H0745	Human CLDN1 (Claudin 1) ELISA Kit
E-EL-H0746	Human CLDN14 (Claudin 14) ELISA Kit
E-EL-H0747	Human CD109 (Cluster of Differentiation 109) ELISA Kit
E-EL-H0748	Human CD30 (Cluster of Differentiation 30) ELISA Kit
E-EL-H0750	Human CD97 (Cluster of Differentiation 97) ELISA Kit
E-EL-H0751	Human CD99 (Cluster of Differentiation 99) ELISA Kit
E-EL-H0752	Human CD4 (Cluster of Differentiation 4) ELISA Kit
E-EL-H0753	Human NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) ELISA Kit
E-EL-H0754	Human CLDN3 (Claudin 3) ELISA Kit
E-EL-H0755	Human NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit
E-EL-H0756	Human c-myc (c-myc Oncogene Product) ELISA Kit
E-EL-H0758	Human F7 (Coagulation Factor VII) ELISA Kit
E-EL-H0760	Human F9 (Coagulation Factor IX) ELISA Kit
E-EL-H0761	Human F10 (Coagulation Factor X) ELISA Kit
E-EL-H0762	Human F13 (Coagulation Factor XIII) ELISA Kit
E-EL-H0763	Human F12 (Coagulation Factor XII) ELISA Kit
E-EL-H0764	Human F5 (Coagulation Factor V/Plasma Factor V) ELISA Kit
E-EL-H0765	Human F11 (Coagulation Factor XI) ELISA Kit
E-EL-H0766	Human F13A1 (Coagulation Factor XIII A1 Polypeptide) ELISA Kit
E-EL-H0767	Human F13B (Coagulation Factor XIII B Polypeptide) ELISA Kit

Cat.NO	Product Name
E-EL-H0769	Human CFL2 (Cofilin 2, Muscle) ELISA Kit
E-EL-H0771	Human CLA (Collagen Autoantibody) ELISA Kit
E-EL-H0772	Human COL15 (Collagen Type XV) ELISA Kit
E-EL-H0773	Human COL2 (Collagen Type II) ELISA Kit
E-EL-H0774	Human COL3 (Collagen Type III) ELISA Kit
E-EL-H0775	Human COL6 (Collagen Type VI) ELISA Kit
E-EL-H0776	Human COL1 α 2 (Collagen Type I Alpha 2) ELISA Kit
E-EL-H0777	Human COL2 α 1 (Collagen Type II Alpha 1) ELISA Kit
E-EL-H0778	Human COL3 α 1 (Collagen Type III Alpha 1) ELISA Kit
E-EL-H0779	Human COL4 α 1 (Collagen Type IV Alpha 1) ELISA Kit
E-EL-H0780	Human COL6 α 1 (Collagen Type VI Alpha 1) ELISA Kit
E-EL-H0781	Human COL7 (Collagen Type VII) ELISA Kit
E-EL-H0782	Human COL8 α 1 (Collagen Type VIII Alpha 1) ELISA Kit
E-EL-H0783	Human COL10 (Collagen Type X) ELISA Kit
E-EL-H0784	Human COL11 α 1 (Collagen Type XI Alpha 1) ELISA Kit
E-EL-H0785	Human COL17 (Collagen Type XVII) ELISA Kit
E-EL-H0786	Human Collagenase I ELISA Kit
E-EL-H0787	Human Collagenase II ELISA Kit
E-EL-H0788	Human NPYR1 (Neuropeptide Y Receptor Y1) ELISA Kit
E-EL-H0790	Human CLL1 (Collectin Liver 1) ELISA Kit

Cat.NO	Product Name
E-EL-H0791	Human COLEC11 (Collectin Sub-Family Member 11) ELISA Kit
E-EL-H0792	Human CCA (Colon Cancer Antigen) ELISA Kit
E-EL-H0793	Human CCSA-2 (Colon Cancer Specific Antigen-2) ELISA Kit
E-EL-H0794	Human CCSA-3 (Colon Cancer Specific Antigen-3) ELISA Kit
E-EL-H0795	Human CCSA-4 (Colon Cancer Specific Antigen-4) ELISA Kit
E-EL-H0796	Human CSF2R α (Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit
E-EL-H0797	Human CLDN9 (Claudin 9) ELISA Kit
E-EL-H0798	Human PIINP (Procollagen II N-Terminal ProPeptide) ELISA Kit
E-EL-H0799	Human GCSFR (Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
E-EL-H0802	Human C1INH (Complement 1 Inhibitor) ELISA Kit
E-EL-H0803	Human C1q (Complement 1q) ELISA Kit
E-EL-H0804	Human C3 (Complement Component 3) ELISA Kit
E-EL-H0805	Human CLDN4 (Claudin 4) ELISA Kit
E-EL-H0806	Human C7 (Complement Component 7) ELISA Kit
E-EL-H0807	Human C1qR1 (Complement Component 1q Receptor) ELISA Kit
E-EL-H0808	Human C1r (Complement Component 1R) ELISA Kit
E-EL-H0809	Human C2 (Complement Component 2) ELISA Kit
E-EL-H0810	Human C5 (Complement Component 5) ELISA Kit
E-EL-H0811	Human C8b (Complement Component 8b) ELISA Kit
E-EL-H0812	Human C8g (Complement Component 8g) ELISA Kit

Cat.NO	Product Name
E-EL-H0813	Human C9 (Complement Component 9) ELISA Kit
E-EL-H0814	Human CFB (Complement Factor B) ELISA Kit
E-EL-H0815	Human CFH (Complement Factor H) ELISA Kit
E-EL-H0816	Human CFI (Complement Factor I) ELISA Kit
E-EL-H0817	Human CFP (Complement Factor P) ELISA Kit
E-EL-H0818	Human C3a (Complement Component 3a) ELISA Kit
E-EL-H0819	Human C3bR (Complement Fragment 3b Receptor) ELISA Kit
E-EL-H0820	Human C3b (Complement Fragment 3b) ELISA Kit
E-EL-H0821	Human C3c (Complement Fragment 3c) ELISA Kit
E-EL-H0822	Human NPYR2 (Neuropeptide Y Receptor Y2) ELISA Kit
E-EL-H0823	Human C4b (Complement Component 4b) ELISA Kit
E-EL-H0824	Human C5c (Complement C5 Convertase) ELISA Kit
E-EL-H0825	Human C5b (Complement Fragment 5b) ELISA Kit
E-EL-H0826	Human NPYR5 (Neuropeptide Y Receptor Y5) ELISA Kit
E-EL-H0827	Human ConA (Concanavalin A) ELISA Kit
E-EL-H0828	Human CTGF (Connective Tissue Growth Factor) ELISA Kit
E-EL-H0829	Human CTAPIII (Connective Tissue Activating Peptide III) ELISA Kit
E-EL-H0830	Human CX26 (Connexin 26) ELISA Kit
E-EL-H0831	Human CX31 (Connexin 31) ELISA Kit
E-EL-H0832	Human CHGB (Chromogranin B) ELISA Kit

Cat.NO	Product Name
E-EL-H0833	Human DEF α 1 (Defensin Alpha 1, Neutrophil) ELISA Kit
E-EL-H0835	Human CTXI (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
E-EL-H0836	Human NTXI (Cross Linked N-telopeptide of Type I Collagen) ELISA Kit
E-EL-H0837	Human CTX-II(Cross Linked C-telopeptide of Type II Collagen)ELISA Kit
E-EL-H0838	Human LDHB (Lactate Dehydrogenase B) ELISA Kit
E-EL-H0839	Human HELIX- II (Type II Collagen Helical Peptide) ELISA Kit
E-EL-H0840	Human HMGCS1 (3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit
E-EL-H0841	Human P III CP (Procollagen III C-Terminal ProPeptide) ELISA Kit
E-EL-H0842	Human DKK4 (Dickkopf Related Protein 4) ELISA Kit
E-EL-H0844	Human MTHFR (Methylenetetrahydrofolate Reductase) ELISA Kit
E-EL-H0845	Human PRKA α 1 (Protein Kinase, AMP Activated Alpha 1) ELISA Kit
E-EL-H0846	Human NFKBIE (Nuclear Factor Of Kappa Light Polypeptide Gene Enhancer In B-cells Inhibitor Epsilon) ELISA Kit
E-EL-H0847	Human PTPRB (Protein Tyrosine Phosphatase Receptor Type B) ELISA Kit
E-EL-H0848	Human BF (Factor B) ELISA Kit
E-EL-H0849	Human CREB (Cyclic AMP Response Element Binding Protein) ELISA Kit
E-EL-H0851	Human CPP (Copeptin) ELISA Kit
E-EL-H0852	Human CLDN11 (Claudin 11) ELISA Kit
E-EL-H0853	Human CXCR1 (CXC-Chemokine Receptor 1) ELISA Kit
E-EL-H0854	Human CXCR3 (Chemokine CXC Receptor 3) ELISA Kit
E-EL-H0855	Human SCG2 (Secretogranin II) ELISA Kit

Cat.NO	Product Name
E-EL-H0856	Human MYCBP (C-Myc Binding Protein) ELISA Kit
E-EL-H0857	Human MYC (V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit
E-EL-H0858	Human CLEC2C (C-Type Lectin Domain Family 2 Member C) ELISA Kit
E-EL-H0859	Human CLEC4E (C-Type Lectin Domain Family 4 Member E) ELISA Kit
E-EL-H0860	Human LDHC (Lactate Dehydrogenase C) ELISA Kit
E-EL-H0861	Human EphA2 (Tyrosine protein kinase receptor A2) ELISA Kit
E-EL-H0862	Human DKK2 (Dickkopf Related Protein 2) ELISA Kit
E-EL-H0863	Human DBP (DNA Binding Protein) ELISA Kit
E-EL-H0864	Human D2D (D-Dimer) ELISA Kit
E-EL-H0866	Human D-LDH (D-Lactate Dehydrogenase) ELISA Kit
E-EL-H0867	Human F8 (Coagulation Factor VIII) ELISA Kit
E-EL-H0868	Human PA2G4 (Proliferation Associated Protein 2G4) ELISA Kit
E-EL-H0869	Human COL1 α 1 (Collagen Type I Alpha 1) ELISA Kit
E-EL-H0870	Human COL6 α 3 (Collagen Type VI Alpha 3) ELISA Kit
E-EL-H0871	Human COL13 (Collagen Type XIII) ELISA Kit
E-EL-H0872	Human COL14 (Collagen Type XIV) ELISA Kit
E-EL-H0873	Human CTX-III (Cross Linked C-telopeptide of Type III Collagen) ELISA Kit
E-EL-H0874	Human EJ/GlyRS (Anti-EJ-Antibody) ELISA Kit
E-EL-H0875	Human EFNA1 (Ephrin A1) ELISA Kit
E-EL-H0876	Human SELE (E-Selectin) ELISA Kit

Cat.NO	Product Name
E-EL-H0877	Human F-Box Protein 2(FBXO2) ELISA Kit
E-EL-H0878	Human Flt3 (FMS Like Tyrosine Kinase 3) ELISA Kit
E-EL-H0879	Human GATA1 (GATA Binding Protein 1) ELISA Kit
E-EL-H0880	Human GATA2 (GATA Binding Protein 2) ELISA Kit
E-EL-H0881	Human GATA3 (GATA Binding Protein 3) ELISA Kit
E-EL-H0882	Human GATA4 (GATA Binding Protein 4) ELISA Kit
E-EL-H0883	Human GATA5 (GATA Binding Protein 5) ELISA Kit
E-EL-H0884	Human GDI1 (GDP Dissociation Inhibitor 1) ELISA Kit
E-EL-H0885	Human GAP (GTPase Activating Protein) ELISA Kit
E-EL-H0886	Human GNβ1 (G Protein Beta 1) ELISA Kit
E-EL-H0887	Human RGS4 (Regulator of G Protein Signaling 4) ELISA Kit
E-EL-H0888	Human HHIP (Hedgehog Interacting Protein) ELISA Kit
E-EL-H0889	Human H-Cadherin (Cadherin, Heart) ELISA Kit
E-EL-H0890	Human NRG-4 (Neuregulin 4) ELISA Kit
E-EL-H0891	Human FcαR (Fc Fragment of IgA Receptor) ELISA Kit
E-EL-H0892	Human PTPRK (Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit
E-EL-H0895	Human SELL (L-Selectin) ELISA Kit
E-EL-H0896	Human MXD1 (MAX Dimerization Protein 1) ELISA Kit
E-EL-H0897	Human HBμ (Hemoglobin, Mu) ELISA Kit
E-EL-H0898	Human NAGase (N-Acetyl Beta-D-Glucosaminidase) ELISA Kit

Cat.NO	Product Name
E-EL-H0899	Human NRG-2 (Neuregulin 2) ELISA Kit
E-EL-H0900	Human NRG-3 (Neuregulin 3) ELISA Kit
E-EL-H0901	Human NLR (NOD-Like Receptor9) ELISA Kit
E-EL-H0902	Human NT-ProBNP (N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
E-EL-H0903	Human Ext-N (N-Terminal Extein) ELISA Kit
E-EL-H0904	Human N-MID-OT (N-MID Osteocalcin) ELISA Kit
E-EL-H0905	Human PTPRN (Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
E-EL-H0906	Human Ras-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit
E-EL-H0907	Human GAS (N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
E-EL-H0908	Human PTPRO (Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
E-EL-H0909	Human P27 (P27 Protein) ELISA Kit
E-EL-H0910	Human TP53 (Tumor Protein p53) ELISA Kit
E-EL-H0911	Human PDAP1 (Pdgfa Associated Protein 1) ELISA Kit
E-EL-H0912	Human PDCN (Podocin) ELISA Kit
E-EL-H0913	Human Preptin ELISA Kit
E-EL-H0914	Human P-Cadherin (Cadherin, Placental) ELISA Kit
E-EL-H0915	Human P-GP (Permeability Glycoprotein) ELISA Kit
E-EL-H0917	Human SELP (P-Selectin) ELISA Kit
E-EL-H0918	Human RAPGEF1 (Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit
E-EL-H0919	Human Rac-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-H0920	Human RHOB (Ras Homolog Gene Family, Member B) ELISA Kit
E-EL-H0921	Human RCOR3 (Rest Corepressor 3) ELISA Kit
E-EL-H0922	Human RHOA (Ras Homolog Gene Family, Member A) ELISA Kit
E-EL-H0923	Human Rock-2 (Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit
E-EL-H0924	Human RND1 (Rho Family GTPase 1) ELISA Kit
E-EL-H0925	Human RUNX2 (Runt-Related Transcription Factor 2) ELISA Kit
E-EL-H0926	Human COL9α3 (Collagen Type IX Alpha 3) ELISA Kit
E-EL-H0927	Human S-100 (Soluble Protein-100) ELISA Kit
E-EL-H0928	Human Salusin α (Salusin Alpha) ELISA Kit
E-EL-H0929	Human SHC/SLP-76 (SH2-Containing Protein) ELISA Kit
E-EL-H0930	Human SHH (Hedgehog Homolog, Sonic) ELISA Kit
E-EL-H0931	Human Slit2 (Slit Homolog 2) ELISA Kit
E-EL-H0932	Human Smad1/MADH1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit
E-EL-H0933	Human Smad4/MADH4 (Mothers Against Decapentaplegic Homolog 4) ELISA Kit
E-EL-H0934	Human Smad7/MADH7 (Mothers Against Decapentaplegic Homolog 7) ELISA Kit
E-EL-H0935	Human Sp100 (Sp100 Nuclear Antigen) ELISA Kit
E-EL-H0936	Human PIAS1 (Protein Inhibitor of Activated STAT 1) ELISA Kit
E-EL-H0937	Human PIAS3 (Protein Inhibitor of Activated STAT3) ELISA Kit
E-EL-H0938	Human PCDHβ16 (Protocadherin Beta 16) ELISA Kit
E-EL-H0939	Human COL5α2 (Collagen Type V Alpha 2) ELISA Kit

Cat.NO	Product Name
E-EL-H0940	Human SAG (S Antigen) ELISA Kit
E-EL-H0941	Human THP (Tamm–Horsfall Glycoprotein) ELISA Kit
E-EL-H0942	Human TDP43 (Tar DNA Binding Protein 43KDa)ELISA Kit
E-EL-H0943	Human TBP (TATA Binding Protein) ELISA Kit
E-EL-H0944	Human TAF (TATA Box Binding Protein/TBP-Associated Factors) ELISA Kit
E-EL-H0945	Human TAF1 (TATA Box Binding Protein Associated Factor 1) ELISA Kit
E-EL-H0946	Human TAF2 (TATA Box Binding Protein Associated Factor 2) ELISA Kit
E-EL-H0947	Human TAF13 (TATA Box Binding Protein Associated Factor 13) ELISA Kit
E-EL-H0948	Human MAP τ (Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
E-EL-H0949	Human TIEG1 (TGF Beta Inducible Early Response Gene 1) ELISA Kit
E-EL-H0951	Human TLR2 (Toll-Like Receptor 2) ELISA Kit
E-EL-H0952	Human TLR9 (Toll-Like Receptor 9) ELISA Kit
E-EL-H0953	Human TU3A (TU3A Protein) ELISA Kit
E-EL-H0954	Human PAR3 (Protease Activated Receptor 3) ELISA Kit
E-EL-H0955	Human LAT (Linker For Activation of T-cell) ELISA Kit
E-EL-H0956	Human TALLA-1 (T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit
E-EL-H0957	Human TCGF- III (Tcell Growth Factor III) ELISA Kit
E-EL-H0958	Human TCR (Tcell Receptor) ELISA Kit
E-EL-H0959	Human Tgtp (T-cell Specific Gtpase) ELISA Kit
E-EL-H0960	Human β -CTx(Beta Crosslaps)ELISA Kit

Cat.NO	Product Name
E-EL-H0961	Human C3c (Complement C3 Convertase) ELISA Kit
E-EL-H0962	Human XCR1 (XC-Chemokine Receptor 1) ELISA Kit
E-EL-H0963	Human PCDH15 (Protocadherin 15) ELISA Kit
E-EL-H0964	Human α 1BG (Alpha-1-B-Glycoprotein) ELISA Kit
E-EL-H0965	Human α 1-AT (Alpha 1-Antitrypsin) ELISA Kit
E-EL-H0966	Human PLC γ 2 (Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit
E-EL-H0967	Human α 2-M (Alpha-2 Macroglobulin) ELISA Kit
E-EL-H0968	Human α 2-AP (α 2-Antiplasmin) ELISA Kit
E-EL-H0970	Human IDU α (Alpha-L-Iduronidase) ELISA Kit
E-EL-H0971	Human NAGLU (N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit
E-EL-H0972	Human GL α (Galactosidase Alpha) ELISA Kit
E-EL-H0975	Human α -GST (α -glutathione S-transferases) ELISA Kit
E-EL-H0976	Human α -SCA (α -sarcomeric Actin) ELISA Kit
E-EL-H0977	Human α Actin (Alpha-Actin) ELISA Kit
E-EL-H0979	Human α -SMA (Alpha-Smooth Muscle Actin) ELISA Kit
E-EL-H0980	Human α -Glu (Alpha Glucosidase) ELISA Kit
E-EL-H0981	Human α HBDH (α -Hydroxybutyrate Dehydrogenase) ELISA Kit
E-EL-H0982	Human α -La (Alpha-Lactalbumin) ELISA Kit
E-EL-H0983	Human SNC α (Synuclein Alpha) ELISA Kit
E-EL-H0985	Human PDGFR β (Platelet Derived Growth Factor Receptor Beta) ELISA Kit

Cat.NO	Product Name
E-EL-H0986	Human β 2-GP1 Ab IgA (β 2-Glycoprotein 1 antibody IgA) ELISA Kit
E-EL-H0987	Human PIICP (Procollagen II C-Terminal ProPeptide) ELISA Kit
E-EL-H0988	Human β 2-GP1 Ab IgG (β 2-Glycoprotein 1 Antibody IgG) ELISA Kit
E-EL-H0989	Human β 2-GP1 Ab IgM (β 2-Glycoprotein 1 Antibody IgM) ELISA Kit
E-EL-H0990	Human PCDHyA2 (Protocadherin Gamma A2) ELISA Kit
E-EL-H0991	Human GL β (Galactosidase Beta) ELISA Kit
E-EL-H0992	Human MAN1A1 (Mannosidase Alpha Class 1A Member 1) ELISA Kit
E-EL-H0993	Human β -LPH (β -Lipotropic Hormone) ELISA Kit
E-EL-H0994	Human β APP (β Amyloid Precursor Protein) ELISA Kit
E-EL-H0995	Human DEFB1 (Defensin Beta 1) ELISA Kit
E-EL-H0996	Human DEF β 2/DEFB2 (Defensin Beta 2) ELISA Kit
E-EL-H0997	Human β Manase (β Mannosidase) ELISA Kit
E-EL-H0998	Human β -actin (Actin Beta) ELISA Kit
E-EL-H0999	Human CYSLTR1 (Cysteinyl Leukotriene Receptor 1) ELISA Kit
E-EL-H1000	Human LACT β (β -Lactamase) ELISA Kit
E-EL-H1001	Human β -Glu (Beta Glucosidase) ELISA Kit
E-EL-H1002	Human IRAK3 (Interleukin 1 Receptor Associated Kinase 3) ELISA Kit
E-EL-H1003	Human GUS β (β -glucuronidase) ELISA Kit
E-EL-H1005	Human SYCN (Syncollin) ELISA Kit
E-EL-H1006	Human TUB β (β -Tubulin) ELISA Kit

Cat.NO	Product Name
E-EL-H1007	Human β -TG (β -Thromboglobulin) ELISA Kit
E-EL-H1009	Human Mly/CXCL9 (Monocyte Interferon Gamma Inducing Factor) ELISA Kit
E-EL-H1010	Human γ -ECS (γ -Glutamyl Systeine Synthetase) ELISA Kit
E-EL-H1011	Human PLA2G10 (Phospholipase A2, Group X) ELISA Kit
E-EL-H1012	Human γ GT1 (Gamma Glutamyltransferase 1) ELISA Kit
E-EL-H1014	Human POMC(Pro-Opiomelanocortin)ELISA Kit
E-EL-H1015	Human OXA (Orexin A) ELISA Kit
E-EL-H1016	Human PIK3C β (Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit
E-EL-H1017	Human OIT3 (Oncoprotein Induced Transcript 3) ELISA Kit
E-EL-H1019	Human PIK3C δ (Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit
E-EL-H1020	Human OBTF1 (Octamer Binding Transcription Factor 1) ELISA Kit
E-EL-H1021	Human OBTF2 (Octamer Binding Transcription Factor 2) ELISA Kit
E-EL-H1022	Human OCT4 (Octamer Binding Transcription Factor 4) ELISA Kit
E-EL-H1023	Human OTF2A (Octamer Transcription Factor 2A) ELISA Kit
E-EL-H1024	Human OTF2B (Octamer Transcription Factor 2B) ELISA Kit
E-EL-H1025	Human TR (Targeted Ribonuclease) ELISA Kit
E-EL-H1026	Human PRDX5 (Peroxisome Proliferator-Induced Protein 5) ELISA Kit
E-EL-H1028	Human IL1R1 (Interleukin 1 Receptor Type I) ELISA Kit
E-EL-H1029	Human PLA2G12B (Phospholipase A2, Group XII B) ELISA Kit
E-EL-H1030	Human IL-2R α (Interleukin 2 Receptor Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-H1031	Human IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit
E-EL-H1032	Human NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
E-EL-H1035	Human LTC4 (Leukotriene C4) ELISA Kit
E-EL-H1036	Human LTD4 (Leukotriene D4) ELISA Kit
E-EL-H1037	Human LT-E4 (Leukotriene E4) ELISA Kit
E-EL-H1038	Human HLE (Leukocyte Exastase)ELISA Kit
E-EL-H1039	Human LR (Leukoregulin) ELISA Kit
E-EL-H1040	Human LCA/CD45(Leukocyte Common Antigen) ELISA Kit
E-EL-H1041	Human LILRB4 (Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit
E-EL-H1042	Human LIFR (Leukemia Inhibitory Factor Receptor) ELISA Kit
E-EL-H1043	Human GARS (Glycyl tRNA Synthetase) ELISA Kit
E-EL-H1044	Human EPHA1 (Ephrin Type A Receptor 1) ELISA Kit
E-EL-H1045	Human SPOP (Speckle Type POZ Protein) ELISA Kit
E-EL-H1046	Human GP4/CD36 (Platelet Membrane Glycoprotein IV) ELISA Kit
E-EL-H1047	Human NSE (Neuron-Specific Enolase) ELISA Kit
E-EL-H1048	Human PLA2G4D (Phospholipase A2, Group IV D) ELISA Kit
E-EL-H1050	Human CYSLTR2 (Cysteinyl Leukotriene Receptor 2) ELISA Kit
E-EL-H1051	Human GAL1 (Galectin 1) ELISA Kit
E-EL-H1052	Human GAL12 (Galectin 12) ELISA Kit
E-EL-H1053	Human GAL14 (Galectin 14) ELISA Kit

Cat.NO	Product Name
E-EL-H1054	Human GAL2 (Galectin 2) ELISA Kit
E-EL-H1055	Human EPHA3 (Ephrin Type A Receptor 3) ELISA Kit
E-EL-H1056	Human GAL4 (Galectin 4) ELISA Kit
E-EL-H1057	Human GAL8 (Galectin 8) ELISA Kit
E-EL-H1058	Human GAL7 (Galectin 7) ELISA Kit
E-EL-H1059	Human GAL9 (Galectin 9) ELISA Kit
E-EL-H1060	Human Gal-6S (Galactose-6-Sulphatase) ELISA Kit
E-EL-H1061	Human UBC (Urinary Bladder Cancer Antigen) ELISA Kit
E-EL-H1062	Human PLA2R1 (Phospholipase A2 Receptor 1) ELISA Kit
E-EL-H1063	Human PLB (Phospholipase B) ELISA Kit
E-EL-H1064	Human cPLA2 (Cytosolic Phospholipase A2) ELISA Kit
E-EL-H1066	Human NP (Nandrolone Phenylpropionate) ELISA Kit
E-EL-H1067	Human EFNA3 (Ephrin A3) ELISA Kit
E-EL-H1068	Human PLD2 (Phospholipase D2) ELISA Kit
E-EL-H1069	Human PYD (Pyridinoline) ELISA Kit
E-EL-H1071	Human PDXK (Pyridoxal Kinase) ELISA Kit
E-EL-H1072	Human PDH β (Pyruvate Dehydrogenase Beta) ELISA Kit
E-EL-H1073	Human OCLN (Occludin) ELISA Kit
E-EL-H1074	Human SEMA7A (Semaphorin 7A) ELISA Kit
E-EL-H1075	Human PLVAP (Plasmalemma Vesicle Associated Protein) ELISA Kit

Cat.NO	Product Name
E-EL-H1076	Human MN (Metanephrine) ELISA Kit
E-EL-H1077	Human IGF1R3 (Insulin Like Growth Factor 1 Receptor 3) ELISA Kit
E-EL-H1078	Human HAGH (Hydroxyacylglutathione Hydrolase) ELISA Kit
E-EL-H1079	Human SPJ (Surfactant Associated Protein J) ELISA Kit
E-EL-H1080	Human MIgA (Membrane Immumoglobulin A) ELISA Kit
E-EL-H1081	Human MIgD (Membrane Immumoglobulin D) ELISA Kit
E-EL-H1082	Human mIgM (Membrane Immumoglobulin M) ELISA Kit
E-EL-H1083	Human EK (Eepidermal Keratin) ELISA Kit
E-EL-H1084	Human GPRC5A (G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit
E-EL-H1085	Human CAPF (Cuticular Active Peptide Factor) ELISA Kit
E-EL-H1086	Human FABP5 (Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
E-EL-H1087	Human VEGFR1/FLT1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit
E-EL-H1088	Human ET (Epitestosterone) ELISA Kit
E-EL-H1089	Human PKM2 (Pyruvate Kinase, Muscle) ELISA Kit
E-EL-H1090	Human PDH E1 (Pyruvate Dehydrogenase E1) ELISA Kit
E-EL-H1091	Human PDH α (Pyruvate Dehydrogenase Alpha) ELISA Kit
E-EL-H1092	Human PDK4(Pyruvate Dehydrogenase Kinase Isozyme 4)ELISA Kit
E-EL-H1093	Human PDP (Pyruvate Dehydrogenase Phosphatase) ELISA Kit
E-EL-H1094	Human VIM (Vimentin) ELISA Kit
E-EL-H1095	Human VTN (Vitronectin) ELISA Kit

Cat.NO	Product Name
E-EL-H1096	Human PPOX (Protoporphyrinogen Oxidase) ELISA Kit
E-EL-H1097	Human PRB1 (Basic Salivary Proline Rich Protein 1) ELISA Kit
E-EL-H1098	Human FOXC1 (Forkhead Box Protein C1) ELISA Kit
E-EL-H1099	Human FOXO1 (Forkhead Box Protein O1) ELISA Kit
E-EL-H1100	Human FOXO2 (Forkhead Box Protein O2) ELISA Kit
E-EL-H1101	Human FOXO3 (Forkhead Box Protein O3) ELISA Kit
E-EL-H1102	Human FOXO4 (Forkhead Box Protein O4) ELISA Kit
E-EL-H1103	Human FOXP1 (Forkhead Box Protein P1) ELISA Kit
E-EL-H1104	Human FOXP3 (Forkhead Box Protein P3) ELISA Kit
E-EL-H1105	Human LATS (Long Acting Thyroid Stimulator) ELISA Kit
E-EL-H1106	Human EFNA4 (Ephrin A4) ELISA Kit
E-EL-H1108	Human TFF3 (Trefoil Factor 3, Intestinal) ELISA Kit
E-EL-H1110	Human TAGLN (Transgelin) ELISA Kit
E-EL-H1111	Human TEP1 (Telomerase Associated Protein 1) ELISA Kit
E-EL-H1112	Human PRB2 (Basic Salivary Proline Rich Protein 2) ELISA Kit
E-EL-H1113	Human SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit
E-EL-H1114	Human OGP (Osteogenic Growth Peptide) ELISA Kit
E-EL-H1115	Human MPF (Maturation Promoting Factor) ELISA Kit
E-EL-H1116	Human FGF23(Fibroblast Growth Factor 23)ELISA Kit
E-EL-H1117	Human FGF9 (Fibroblast Growth Factor 9) ELISA Kit

Cat.NO	Product Name
E-EL-H1118	Human FGFBP1 (Fibroblast Growth Factor Binding Protein 1) ELISA Kit
E-EL-H1119	Human FGFBP2 (Fibroblast Growth Factor Binding Protein 2) ELISA Kit
E-EL-H1120	Human FRS2 (Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit
E-EL-H1121	Human FRS3 (Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
E-EL-H1122	Human VLA (Very Late Appearing Antigen) ELISA Kit
E-EL-H1123	Human PRF1 (Perforin 1/Pore-Forming Protein) ELISA Kit
E-EL-H1124	Human FSHR (Follicle Stimulating Hormone Receptor) ELISA Kit
E-EL-H1125	Human VMA (Vanillylmandelic Acid) ELISA Kit
E-EL-H1126	Human TGM1 (Transglutaminase 1, Keratinocyte) ELISA Kit
E-EL-H1127	Human PACAPR (Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit
E-EL-H1129	Human EPHA10 (Ephrin Type A Receptor 10) ELISA Kit
E-EL-H1132	Human ER β (Estrogen Receptor Beta) ELISA Kit
E-EL-H1133	Human ER α (Estrogen Receptor Alpha) ELISA Kit
E-EL-H1134	Human EIPS2 (Estrogen Induced Protein pS2) ELISA Kit
E-EL-H1136	Human EFNA5 (Ephrin A5) ELISA Kit
E-EL-H1137	Human PRDX2 (Peroxiredoxin 2) ELISA Kit
E-EL-H1138	Human MAPKAPK3 (Mitogen Activated Protein Kinase Activated Protein Kinase 3) ELISA Kit
E-EL-H1139	Human TSA (Thyroid Stimulating Antibody) ELISA Kit
E-EL-H1140	Human TEF (Thyrotroph Embryonic Factor) ELISA Kit
E-EL-H1141	Human TSH (Thyroid Stimulating Hormone) ELISA Kit

Cat.NO	Product Name
E-EL-H1142	Human TRH (Thyrotropin Releasing Hormone) ELISA Kit
E-EL-H1143	Human FSH (Follicle-Stimulating Hormone) ELISA Kit
E-EL-H1144	Human CRH (Corticotropin Releasing Hormone) ELISA Kit
E-EL-H1146	Human GHRH (Growth Hormone Releasing Hormone) ELISA Kit
E-EL-H1147	Human δ SIP- α (Delta Sleep Inducing Peptide α) ELISA Kit
E-EL-H1148	Human GasR (Gastrin Receptor) ELISA Kit
E-EL-H1149	Human GnRH(Gonadotropin Releasing Hormone)ELISA Kit
E-EL-H1150	Human SCTR (Secretin Receptor) ELISA Kit
E-EL-H1151	Human SCT (Secretin) ELISA Kit
E-EL-H1154	Human PYGB (Glycogen Phosphorylase, Brain) ELISA Kit
E-EL-H1155	Human FABP7 (Fatty Acid Binding Protein 7, Brain) ELISA Kit
E-EL-H1157	Human MAOA (Type A Monoamine Oxidase) ELISA Kit
E-EL-H1158	Human MCP-2/CCL8 (Monocyte Chemotactic Protein 2) ELISA Kit
E-EL-H1159	Human MCP-4/CCL13 (Monocyte Chemotactic Protein 4) ELISA Kit
E-EL-H1160	Human GPER1 (G Protein Coupled Estrogen Receptor 1) ELISA Kit
E-EL-H1161	Human PRDX1 (Peroxisome oxidoreductase 1) ELISA Kit
E-EL-H1162	Human SREBP (Sterol Regulatory Element Binding Proteins) ELISA Kit
E-EL-H1163	Human ELN (Elastin) ELISA Kit
E-EL-H1165	Human ELA4 (Elastase 4) ELISA Kit
E-EL-H1166	Human EMILIN 1 (Elastin Microfibril Interface Located Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-H1167	Human PROC (Protein C) ELISA Kit
E-EL-H1168	Human PDIA2 (Protein Disulfide Isomerase A2) ELISA Kit
E-EL-H1169	Human PROS (Protein S) ELISA Kit
E-EL-H1170	Human PROZ (Protein Z) ELISA Kit
E-EL-H1171	Human PDIA3 (Protein Disulfide Isomerase A3) ELISA Kit
E-EL-H1172	Human PRSS2 (Protease, Serine, 2) ELISA Kit
E-EL-H1174	Human PKB (Protein Kinase B) ELISA Kit
E-EL-H1175	Human PKBG (Protein Kinase B Gamma) ELISA Kit
E-EL-H1176	Human PKC (Protein Kinase C) ELISA Kit
E-EL-H1177	Human PICK1 (Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
E-EL-H1178	Human PKC β 1 (Protein Kinase C Beta 1) ELISA Kit
E-EL-H1179	Human PKC γ (Protein Kinase C Gamma) ELISA Kit
E-EL-H1180	Human PKC δ 1 (Protein Kinase C Delta) ELISA Kit
E-EL-H1181	Human PKC ϵ (Protein Kinase C Epsilon) ELISA Kit
E-EL-H1182	Human PKC ζ (Protein Kinase C Zeta) ELISA Kit
E-EL-H1183	Human PKC θ (Protein Kinase C Theta)ELISA Kit
E-EL-H1184	Human PRSS23 (Protease, Serine 23) ELISA Kit
E-EL-H1185	Human PKD1 (Protein Kinase D1) ELISA Kit
E-EL-H1186	Human PKD2 (Protein Kinase D2)ELISA Kit
E-EL-H1187	Human PKN1 (Protein Kinase N1) ELISA Kit

Cat.NO	Product Name
E-EL-H1188	Human PKI α (Protein Kinase Inhibitor Alpha) ELISA Kit
E-EL-H1189	Human PKI β (Protein Kinase Inhibitor Beta) ELISA Kit
E-EL-H1190	Human PKI γ (Protein Kinase Inhibitor Gamma) ELISA Kit
E-EL-H1191	Human PG (Proteoglycan) ELISA Kit
E-EL-H1192	Human PTK (Protein Tyrosine Kinase) ELISA Kit
E-EL-H1193	Human PTPRJ (Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit
E-EL-H1194	Human PP (Protein Phosphatase) ELISA Kit
E-EL-H1195	Human PPP1R1A (Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit
E-EL-H1196	Human PRSS1 (Protease, Serine, 1) ELISA Kit
E-EL-H1197	Human PRSS12 (Protease, Serine 12) ELISA Kit
E-EL-H1198	Human PR3Ab (Proteinase 3 Antibody) ELISA Kit
E-EL-H1200	Human PAR2 (Protease Activated Receptor 2) ELISA Kit
E-EL-H1201	Human PCI (Protein C Inhibitor) ELISA Kit
E-EL-H1202	Human PLP (Proteolipid Protein) ELISA Kit
E-EL-H1203	Human PDI (Protein Disulfide Isomerase) ELISA Kit
E-EL-H1204	Human PDIA5 (Protein Disulfide Isomerase A5) ELISA Kit
E-EL-H1205	Human NGF (Nerve Growth Factor) ELISA Kit
E-EL-H1206	Human PRSS33 (Protease, Serine 33) ELISA Kit
E-EL-H1208	Human CK-LMW (Cytokeratin Low Molecular Weight) ELISA Kit
E-EL-H1209	Human LMP7/PSMB9 (Low Molecular-weight Protein/proteasome Beta-type Subunit) ELISA Kit

Cat.NO	Product Name
E-EL-H1210	Human LDL (Low Density Lipoprotein) ELISA Kit
E-EL-H1211	Human LDLR (Low Density Lipoprotein Receptor) ELISA Kit
E-EL-H1212	Human LRP6 (Low Density Lipoprotein Receptor Related Protein) ELISA Kit
E-EL-H1213	Human RETN (Resistin) ELISA Kit
E-EL-H1214	Human RPS6K β 1 (Ribosomal Protein S6 Kinase Beta 1) ELISA Kit
E-EL-H1215	Human ETF β (Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit
E-EL-H1216	Human APP (Amyloid Precursor Protein) ELISA Kit
E-EL-H1218	Human M-AChR (Muscarinic Acetylcholine Receptor) ELISA Kit
E-EL-H1219	Human M-AChR M1 (Muscarinic Acetylcholine Receptor M1) ELISA Kit
E-EL-H1220	Human M-AChR M2 (Muscarinic Acetylcholine Receptor M2) ELISA Kit
E-EL-H1221	Human M-AChR M3 (Muscarinic Acetylcholine Receptor M3) ELISA Kit
E-EL-H1224	Human D1R/DRD1 (Dopamine Receptor D1) ELISA Kit
E-EL-H1225	Human D2R/DRD2 (Dopamine Receptor D2) ELISA Kit
E-EL-H1226	Human D β H (Dopamine- β -Hydroxylase) ELISA Kit
E-EL-H1227	Human D3R/DRD3 (Dopamine Receptor D3) ELISA Kit
E-EL-H1228	Human D4R/DRD4 (Dopamine Receptor D4) ELISA Kit
E-EL-H1229	Human D5R/DRD5 (Dopamine Receptor D5) ELISA Kit
E-EL-H1230	Human DDC (Dopamine Decarboxylase) ELISA Kit
E-EL-H1231	Human DAT (Dopamine Transporter) ELISA Kit
E-EL-H1232	Human DT (Dopachrome Tautomerase) ELISA Kit

Cat.NO	Product Name
E-EL-H1233	Human VS (Versican/PG-M) ELISA Kit
E-EL-H1235	Human PHSA (Poly Serum Antigen) ELISA Kit
E-EL-H1236	Human PIGR/SC (Polymeric Immunoglobulin Receptor/Membrane Secretory Component) ELISA Kit
E-EL-H1237	Human PYY (Peptide YY) ELISA Kit
E-EL-H1238	Human PTN (Pleiotrophin) ELISA Kit
E-EL-H1240	Human MRP (Multidrug Resistance-Associated Protein) ELISA Kit
E-EL-H1241	Human DAO (Diamine Oxidase) ELISA Kit
E-EL-H1242	Human UGT2B4 (Uridine Diphosphate Glucuronosyltransferase 2 Family Polypeptide B4) ELISA Kit
E-EL-H1243	Human DLD (Dihydrolipoyl Dehydrogenase) ELISA Kit
E-EL-H1244	Human DLAT (Dihydrolipoyl Transacetylase) ELISA Kit
E-EL-H1245	Human DPYSL2 (Dihydropyrimidinase Like Protein 2) ELISA Kit
E-EL-H1246	Human DPYSL3 (Dihydropyrimidinase Like Protein 3) ELISA Kit
E-EL-H1247	Human NRP1 (Neuropilin 1) ELISA Kit
E-EL-H1248	Human DNM1 (Dynamin 1) ELISA Kit
E-EL-H1249	Human DNM2 (Dynamin 2) ELISA Kit
E-EL-H1250	Human fFDFT1 (Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit
E-EL-H1251	Human rT3 (Reverse Triiodothyronine) ELISA Kit
E-EL-H1252	Human Ub (Ubiquitin) ELISA Kit
E-EL-H1253	Human UBD (Ubiquitin D) ELISA Kit
E-EL-H1254	Human DUB (Deubiquitinating Enzyme) ELISA Kit

Cat.NO	Product Name
E-EL-H1255	Human UBE1L (Ubiquitin Activating Enzyme E1 Like Protein) ELISA Kit
E-EL-H1256	Human UBE2C (Ubiquitin Conjugating Enzyme E2C) ELISA Kit
E-EL-H1257	Human UBE3A (Ubiquitin Protein Ligase E3A) ELISA Kit
E-EL-H1260	Human NNE (Non-Neuronal Enolase) ELISA Kit
E-EL-H1261	Human SDC3 (Syndecan 3) ELISA Kit
E-EL-H1262	Human TPS (Tryptase) ELISA Kit
E-EL-H1263	Human FBLN5 (Fibulin 5) ELISA Kit
E-EL-H1264	Human DR-70TM (Tumor Marker DR-70 For Lung Cancer) ELISA Kit
E-EL-H1265	Human REG1 β (Regenerating Islet Derived Protein 1 Beta) ELISA Kit
E-EL-H1266	Human SPA (Pulmonary Surfactant Associated Protein A) ELISA Kit
E-EL-H1267	Human SPB (Pulmonary Surfactant Associated Protein B) ELISA Kit
E-EL-H1268	Human SPC (Pulmonary Surfactant Associated Protein C) ELISA Kit
E-EL-H1269	Human SPD (Pulmonary Surfactant Associated Protein D) ELISA Kit
E-EL-H1270	Human PARC/CCL18 (Pulmonary Activation Regulated Chemokine) ELISA Kit
E-EL-H1271	Human LRP (Lung Resistance-Related Protein) ELISA Kit
E-EL-H1272	Human PRSS8 (Protease, Serine 8) ELISA Kit
E-EL-H1273	Human SPARCL1 (SPARC Like Protein 1) ELISA Kit
E-EL-H1274	Human SEMA3A (Semaphorin 3A) ELISA Kit
E-EL-H1275	Human sIgA (Secretory Immunoglobulin A) ELISA Kit
E-EL-H1276	Human SLPI (Secretory Leukocyte Protease Inhibitor) ELISA Kit

Cat.NO	Product Name
E-EL-H1277	Human HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit
E-EL-H1278	Human SAA3 (Serum Amyloid A3) ELISA Kit
E-EL-H1279	Human SAP (Serum Amyloid P Component) ELISA Kit
E-EL-H1280	Human PAR4 (Protease Activated Receptor 4) ELISA Kit
E-EL-H1281	Human PRCP (Prolylcarboxypeptidase) ELISA Kit
E-EL-H1282	Human RPA-70 (replication protein A) ELISA Kit
E-EL-H1283	Human GPR120 (G Protein Coupled Receptor 120) ELISA Kit
E-EL-H1284	Human REG3 α (Regenerating Islet Derived Protein 3 Alpha) ELISA Kit
E-EL-H1285	Human Cyr61 (Cysteine Rich Protein, Angiogenic Inducer 61) ELISA Kit
E-EL-H1286	Human SPARC (Secreted Protein Acidic And Rich In Cysteine) ELISA Kit
E-EL-H1287	Human LRG1 (Leucine Rich Alpha-2-Glycoprotein 1) ELISA Kit
E-EL-H1288	Human REG3 γ (Regenerating Islet Derived Protein 3 Gamma) ELISA Kit
E-EL-H1289	Human S100A8 (S100 Calcium Binding Protein A8) ELISA Kit
E-EL-H1290	Human S100A9 (S100 Calcium Binding Protein A9) ELISA Kit
E-EL-H1291	Human GPR37 (G Protein Coupled Receptor 37) ELISA Kit
E-EL-H1292	Human S100A11 (S100 Calcium Binding Protein A11) ELISA Kit
E-EL-H1293	Human S100A12 (S100 Calcium Binding Protein A12) ELISA Kit
E-EL-H1294	Human S100A13 (S100 Calcium Binding Protein A13) ELISA Kit
E-EL-H1295	Human S100A6 (S100 Calcium Binding Protein A6) ELISA Kit
E-EL-H1296	Human S100A7 (S100 Calcium Binding Protein A7) ELISA Kit

Cat.NO	Product Name
E-EL-H1297	Human S100B (S100 Calcium Binding Protein B) ELISA Kit
E-EL-H1298	Human SDC1 (Syndecan 1) ELISA Kit
E-EL-H1299	Human SCF (Stem Cell Factor) ELISA Kit
E-EL-H1300	Human SCFR/c-Kit (Stem Cell Factor Receptor) ELISA Kit
E-EL-H1301	Human GAL (Galanin) ELISA Kit
E-EL-H1302	Human GALR1 (Galanin Receptor 1) ELISA Kit
E-EL-H1303	Human MASP (Mannan Binding Lectin Associated Serine Protease) ELISA Kit
E-EL-H1305	Human MBP/MBL (Mannose Binding Protein/Mannose Binding Lectin) ELISA Kit
E-EL-H1306	Human MBLR (Mannose Binding Lectin Receptor) ELISA Kit
E-EL-H1308	Human PGAM1 (Phosphoglycerate Mutase 1) ELISA Kit
E-EL-H1309	Human PHC (Primary Hepatic Carcinoma) ELISA Kit
E-EL-H1310	Human PK (Pyruvate Kinase, Liver And RBC) ELISA Kit
E-EL-H1311	Human PYGL (Glycogen Phosphorylase, Liver) ELISA Kit
E-EL-H1312	Human SFRP2 (Secreted Frizzled Related Protein 2) ELISA Kit
E-EL-H1313	Human GP-73 (Golgi Protein 73) ELISA Kit
E-EL-H1314	Human CK-HMW (Cytokeratin High Molecular Weight) ELISA Kit
E-EL-H1315	Human SIAE (Sialic Acid Acetylesterase) ELISA Kit
E-EL-H1316	Human ST2 (Syntenin 2) ELISA Kit
E-EL-H1317	Human SIGLEC1 (Sialic Acid Binding Ig Like Lectin 1) ELISA Kit
E-EL-H1318	Human MHB(Methemoglobin)ELISA Kit

Cat.NO	Product Name
E-EL-H1319	Human TR (Testosterone Receptor) ELISA Kit
E-EL-H1320	Human SIGLEC10 (Sialic Acid Binding Immunoglobulin Like Lectin 10) ELISA Kit
E-EL-H1321	Human SIGLEC2 (Sialic Acid Binding Immunoglobulin Like Lectin 2) ELISA Kit
E-EL-H1322	Human PADI4 (Peptidyl Arginine Deiminase Type IV) ELISA Kit
E-EL-H1323	Human SIGLEC8 (Sialic Acid Binding Ig Like Lectin 8) ELISA Kit
E-EL-H1324	Human Smad3/MADH3 (Mothers Against Decapentaplegic Homolog 3) ELISA Kit
E-EL-H1325	Human OFQ/N (Orphanin FQ/Nociceptin) ELISA Kit
E-EL-H1326	Human GCLM (Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit
E-EL-H1327	Human GDH/GLDH (Glutamate dehydrogenase) ELISA Kit
E-EL-H1328	Human GAD (Glutamic Acid Decarboxylase) ELISA Kit
E-EL-H1329	Human GAD2 (Glutamate Decarboxylase 2) ELISA Kit
E-EL-H1330	Human GAD-Ab (Anti-Glutamic Acid Decarboxylase antibody) ELISA Kit
E-EL-H1332	Human GSTs (Glutathione S-Transferases) ELISA Kit
E-EL-H1333	Human GSTα1 (Glutathione S Transferase Alpha 1) ELISA Kit
E-EL-H1334	Human GSTθ1/GSTt1 (Glutathione S Transferase Theta 1) ELISA Kit
E-EL-H1335	Human GSTθ2/GSTt2 (Glutathione S Transferase Theta 2) ELISA Kit
E-EL-H1336	Human GSTω1 (Glutathione S Transferase Omega 1) ELISA Kit
E-EL-H1337	Human HDAC2 (Histone Deacetylase 2) ELISA Kit
E-EL-H1338	Human C1QTNF1 (C1q and Tumor Necrosis Factor Related Protein 1) ELISA Kit
E-EL-H1339	Human GR (Glutathione Reductase) ELISA Kit

Cat.NO	Product Name
E-EL-H1340	Human GSTpi (Glutathione S Transferases Pi) ELISA Kit
E-EL-H1341	Human OPG (Osteoprotegerin) ELISA Kit
E-EL-H1343	Human OC/BGP (Osteocalcin) ELISA Kit
E-EL-H1344	Human ELA3B (Elastase 3B) ELISA Kit
E-EL-H1345	Human MUSK (Muscle Skeletal Receptor Tyrosine Kinase) ELISA Kit
E-EL-H1346	Human GSTα5 (Glutathione S Transferase Alpha 5) ELISA Kit
E-EL-H1347	Human OPN (Osteopontin) ELISA Kit
E-EL-H1348	Human MPIF-1/CCL23 (Myeloid Progenitor Inhibitory Factor 1) ELISA Kit
E-EL-H1349	Human SOX9 (Sex Determining Region Y Box Protein 9) ELISA Kit
E-EL-H1350	Human GREM1 (Gremlin 1) ELISA Kit
E-EL-H1351	Human ON (Osteonectin) ELISA Kit
E-EL-H1352	Human SF1 (Steroidogenic Factor 1) ELISA Kit
E-EL-H1353	Human OP (Oligomeric Protein) ELISA Kit
E-EL-H1354	Human PAPPA2 (Pappalysin 2) ELISA Kit
E-EL-H1355	Human IgA (Immunoglobulin A) ELISA Kit
E-EL-H1356	Human P-CK (Pan-Cytokeratin) ELISA Kit
E-EL-H1357	Human ALDOC (Aldolase C, Fructose Bisphosphate) ELISA Kit
E-EL-H1358	Human PMP3 (Peroxisomal Membrane Protein 3) ELISA Kit
E-EL-H1359	Human PPARγC1α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit
E-EL-H1360	Human PPAR-α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-H1361	Human PPAR-γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
E-EL-H1362	Human LPO (Lactoperoxidase) ELISA Kit
E-EL-H1363	Human NXF1 (Nuclear RNA Export Factor 1) ELISA Kit
E-EL-H1364	Human NMP-22 (Nuclear Matrix Protein 22) ELISA Kit
E-EL-H1365	Human NUP107 (Nucleoporin 107kDa) ELISA Kit
E-EL-H1366	Human NUP133 (Nucleoporin 133kDa) ELISA Kit
E-EL-H1367	Human NUP153 (Nucleoporin 153kDa) ELISA Kit
E-EL-H1368	Human NUP155 (Nucleoporin 155kDa) ELISA Kit
E-EL-H1369	Human NUP160 (Nucleoporin 160kDa) ELISA Kit
E-EL-H1370	Human NUP188 (Nucleoporin 188kDa) ELISA Kit
E-EL-H1371	Human NUP205 (Nucleoporin 205kDa) ELISA Kit
E-EL-H1372	Human NUP214 (Nucleoporin 214kDa) ELISA Kit
E-EL-H1373	Human NUP35 (Nucleoporin 35kDa) ELISA Kit
E-EL-H1374	Human NUP37 (Nucleoporin 37kDa) ELISA Kit
E-EL-H1375	Human NUP50 (Nucleoporin 50kDa) ELISA Kit
E-EL-H1376	Human NUP85 (Nucleoporin 85kDa) ELISA Kit
E-EL-H1377	Human NUP98 (Nucleoporin 98kDa) ELISA Kit
E-EL-H1378	Human NRIP1 (Nuclear Receptor Interacting Protein 1) ELISA Kit
E-EL-H1379	Human RNASE3/ECP (Ribonuclease A3/Eosinophil Cationic Protein) ELISA Kit
E-EL-H1380	Human RNASEH (Ribonuclease H)ELISA Kit

Cat.NO	Product Name
E-EL-H1381	Human RNASET2 (Ribonuclease T2) ELISA Kit
E-EL-H1382	Human RI (Ribonuclease Inhibitor) ELISA Kit
E-EL-H1383	Human RNASE (Ribonuclease) ELISA Kit
E-EL-H1384	Human RPL10 (Ribosomal Protein L10) ELISA Kit
E-EL-H1385	Human RPL26 (Ribosomal Protein L26) ELISA Kit
E-EL-H1386	Human NFκB (Nuclear Factor Kappa B) ELISA Kit
E-EL-H1387	Human RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
E-EL-H1388	Human NF-κB p65 (Nuclear Factor Kappa B p65) ELISA Kit
E-EL-H1389	Human MART/Melan-A (Melanoma Marker A) ELISA Kit
E-EL-H1390	Human MAGE (Melanoma-associated Antigen) ELISA Kit
E-EL-H1391	Human MMSAM (Melanoma Metastasis Surface Adhesion Molecule) ELISA Kit
E-EL-H1392	Human βMSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-H1394	Human MLZE (Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit
E-EL-H1396	Human ESF (Erythropoiesis Stimulating Factor) ELISA Kit
E-EL-H1397	Human EMP (Erythrocyte Membrane Protein) ELISA Kit
E-EL-H1398	Human THRSP (Thyroid Hormone Responsive) ELISA Kit
E-EL-H1399	Human EPOR (Erythropoietin Receptor) ELISA Kit
E-EL-H1400	Human GRIK2 (Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit
E-EL-H1401	Human EXT1 (Exostoses 1) ELISA Kit
E-EL-H1402	Human GSTκ1 (Glutathione S Transferase Kappa 1) ELISA Kit

Cat.NO	Product Name
E-EL-H1403	Human LTA4H (Leukotriene A4 Hydrolase) ELISA Kit
E-EL-H1404	Human GSTα4 (Glutathione S Transferase Alpha 4) ELISA Kit
E-EL-H1405	Human SNCy (Synuclein Gamma) ELISA Kit
E-EL-H1406	Human TXLNa (Taxilin Alpha) ELISA Kit
E-EL-H1407	Human TXLNb (Taxilin Beta) ELISA Kit
E-EL-H1408	Human TLN (Talin) ELISA Kit
E-EL-H1411	Human Sec23b (protein transport protein Sec23B) ELISA Kit
E-EL-H1412	Human MUSKab (Muscle Skeletal Receptor Tyrosine Kinase antibody) ELISA Kit
E-EL-H1414	Human COX-2R (Cyclooxygenase-2 Receptor) ELISA Kit
E-EL-H1415	Human RNF4 (Ring Finger Protein 4) ELISA Kit
E-EL-H1416	Human LHRH/GnRH (Luteinizing Hormone Releasing Hormone) ELISA Kit
E-EL-H1417	Human FABP6 (Fatty Acid Binding Protein 6, Ileal) ELISA Kit
E-EL-H1418	Human TMPRSS2 (Transmembrane Protease, Serine 2) ELISA Kit
E-EL-H1419	Human F XII a (Activated Coagulation Factor XII) ELISA Kit
E-EL-H1420	Human TNNI1 (Troponin I Type 1, Slow Skeletal) ELISA Kit
E-EL-H1421	Human FSCN (Fascin) ELISA Kit
E-EL-H1422	Human TNNC2 (Troponin C Type 2) ELISA Kit
E-EL-H1423	Human TNC (Tenascin C) ELISA Kit
E-EL-H1424	Human TNR (Tenascin R) ELISA Kit
E-EL-H1425	Human TNX (Tenascin X) ELISA Kit

Cat.NO	Product Name
E-EL-H1426	Human IgA1 (Immunoglobulin A1) ELISA Kit
E-EL-H1427	Human MYO (Myosin) ELISA Kit
E-EL-H1428	Human MYL1 (Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit
E-EL-H1429	Human MLCK/MYLK (Myosin Light Chain Kinase) ELISA Kit
E-EL-H1430	Human MYH1 (Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit
E-EL-H1431	Human FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
E-EL-H1432	Human PYGM (Glycogen Phosphorylase, Muscle) ELISA Kit
E-EL-H1433	Human CKM (Creatine Kinase, Muscle) ELISA Kit
E-EL-H1434	Human CKMB(Creatine Kinase MB Isoenzyme)ELISA Kit
E-EL-H1435	Human MYOG (Myogenin) ELISA Kit
E-EL-H1436	Human DMD (Dystrophin) ELISA Kit
E-EL-H1437	Human MSTN (Myostatin) ELISA Kit
E-EL-H1439	Human PUMA (P53 Upregulated Modulator of Apoptosis) ELISA Kit
E-EL-H1440	Human MGP (Matrix Gla Protein) ELISA Kit
E-EL-H1441	Human MMP-1 (Matrix Metalloproteinase 1) ELISA Kit
E-EL-H1442	Human IgA2 (Immunoglobulin A2) ELISA Kit
E-EL-H1443	Human MMP-11 (Matrix Metalloproteinase 11) ELISA Kit
E-EL-H1444	Human MGST1 (Microsomal Glutathione S Transferase 1) ELISA Kit
E-EL-H1445	Human MMP-2 (Matrix Metalloproteinase 2) ELISA Kit
E-EL-H1446	Human MMP-3 (Matrix Metalloproteinase 3) ELISA Kit

Cat.NO	Product Name
E-EL-H1447	Human MMP-17 (Matrix Metalloproteinase 17) ELISA Kit
E-EL-H1448	Human MMP-24 (Matrix Metalloproteinase 24) ELISA Kit
E-EL-H1449	Human MMP-7 (Matrix Metalloproteinase 7) ELISA Kit
E-EL-H1450	Human MMP-8 (Matrix Metalloproteinase 8) ELISA Kit
E-EL-H1451	Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit
E-EL-H1452	Human ELA2B (Elastase 2B) ELISA Kit
E-EL-H1453	Human TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
E-EL-H1454	Human TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
E-EL-H1455	Human TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) ELISA Kit
E-EL-H1456	Human LGALS3BP (Lectin Galactoside Binding, Soluble 3 Binding Protein) ELISA Kit
E-EL-H1457	Human MAT (Matrilysin) ELISA Kit
E-EL-H1461	Human TNNI2 (Troponin I Type 2, Fast Skeletal) ELISA Kit
E-EL-H1462	Human TP/TXA2R (Thromboxane Receptor) ELISA Kit
E-EL-H1463	Human SDF-2 (Stromal Cell Derived Factor 2) ELISA Kit
E-EL-H1464	Human SDF2L1 (Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
E-EL-H1465	Human SDF-4 (Stromal Cell Derived Factor 4) ELISA Kit
E-EL-H1466	Human TDGF1 (Teratocarcinoma Derived Growth Factor 1) ELISA Kit
E-EL-H1467	Human AICAR (AICA Riboside) ELISA Kit
E-EL-H1468	Human KLK8 (Kallikrein 8) ELISA Kit
E-EL-H1469	Human VLDLR (Very Low Density Lipoprotein Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-H1470	Human GAL3 (Galectin 3) ELISA Kit
E-EL-H1471	Human MTX (Methotrexate) ELISA Kit
E-EL-H1472	Human Methylase ELISA Kit
E-EL-H1473	Human PGD2-MOX (Prostaglandin D2-Methoxime) ELISA Kit
E-EL-H1475	Human FPR2 (Formyl Peptide Receptor 2) ELISA Kit
E-EL-H1476	Human PTH (Parathyroid Hormone) ELISA Kit
E-EL-H1477	Human PTHR2 (Parathyroid Hormone Receptor 2) ELISA Kit
E-EL-H1478	Human PTHrP (Parathyroid Hormone Related Protein) ELISA Kit
E-EL-H1479	Human TPP1 (Tripeptidyl Peptidase I) ELISA Kit
E-EL-H1480	Human PLP (Parathyroid Hormone-Like Protein) ELISA Kit
E-EL-H1482	Human TSI (Thyroid Stimulating Immunoglobulin) ELISA Kit
E-EL-H1483	Human THAA (Thyroid Hormone Autoantibody) ELISA Kit
E-EL-H1484	Human CD276 (Cluster of Differentiation 276) ELISA Kit
E-EL-H1485	Human TBG (Thyroxine Binding Globulin) ELISA Kit
E-EL-H1486	Human PAH (Phenylalanine Hydroxylase) ELISA Kit
E-EL-H1487	Human TG (Thyroglobulin) ELISA Kit
E-EL-H1489	Human TAB (Thyroxine Antibody) ELISA Kit
E-EL-H1490	Human Thyronine ELISA Kit
E-EL-H1491	Human GJβ1 (Gap Junction Protein Beta 1) ELISA Kit
E-EL-H1492	Human PCT (Procalcitonin) ELISA Kit

Cat.NO	Product Name
E-EL-H1493	Human PCPE1 (Procollagen C Proteinase Enhancer 1) ELISA Kit
E-EL-H1495	Human GDNF (Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
E-EL-H1496	Human KAF (Keratinocyte Autocrine Factor) ELISA Kit
E-EL-H1497	Human CDSN (Corneodesmosin) ELISA Kit
E-EL-H1498	Human CNTN1 (Contactin 1) ELISA Kit
E-EL-H1499	Human CNTN2 (Contactin 2) ELISA Kit
E-EL-H1500	Human CNTN3 (Contactin 3) ELISA Kit
E-EL-H1501	Human CNTN4 (Contactin 4) ELISA Kit
E-EL-H1502	Human CNTNAP1 (Contactin Associated Protein 1) ELISA Kit
E-EL-H1503	Human CNTNAP4 (Contactin Associated Protein 4) ELISA Kit
E-EL-H1504	Human ELMO2 (Engulfment And Cell Motility 2) ELISA Kit
E-EL-H1505	Human Des (Desmin) ELISA Kit
E-EL-H1506	Human SSRP1 (Structure Specific Recognition Protein 1) ELISA Kit
E-EL-H1507	Human HP (Haptoglobin) ELISA Kit
E-EL-H1508	Human HPR (Haptoglobin Related Protein) ELISA Kit
E-EL-H1509	Human TNNT1 (Troponin T Type 1, Slow Skeletal) ELISA Kit
E-EL-H1510	Human TNPO2 (Transportin 2) ELISA Kit
E-EL-H1511	Human UCP-2 (Uncoupling Protein 2, Mitochondrial) ELISA Kit
E-EL-H1512	Human HLTF (Helicase Like Transcription Factor) ELISA Kit
E-EL-H1513	Human UPP2 (Uridine Phosphorylase 2) ELISA Kit

Cat.NO	Product Name
E-EL-H1514	Human MT1 (Metallothionein 1) ELISA Kit
E-EL-H1515	Human MT2 (Metallothionein 2) ELISA Kit
E-EL-H1516	Human TJP1 (Tight junction protein ZO-1) ELISA Kit
E-EL-H1517	Human CLC (Charcot Leyden Crystal Protein) ELISA Kit
E-EL-H1518	Human ACVB (Activin B) ELISA Kit
E-EL-H1519	Human CRYαB (Crystallin Alpha B) ELISA Kit
E-EL-H1521	Human Sp17 (Sperm Protein 17) ELISA Kit
E-EL-H1522	Human SSFA2 (Sperm Specific Antigen 2) ELISA Kit
E-EL-H1523	Human MSP (Macrophage Stimulating Protein) ELISA Kit
E-EL-H1524	Human MAF (Macrophage Activating Factor) ELISA Kit
E-EL-H1525	Human MCSFR (Macrophage Colony Stimulating Factor Receptor) ELISA Kit
E-EL-H1526	Human MCF (Macrophage Chemotatic Factor) ELISA Kit
E-EL-H1528	Human MIP-3β/ELC/CCL19 (Macrophage Inflammatory Protein 3β) ELISA Kit
E-EL-H1529	Human MIP-5 (Macrophage Inflammatory Protein 5) ELISA Kit
E-EL-H1530	Human MIF (Macrophage Migration Inhibitory Factor) ELISA Kit
E-EL-H1531	Human FZD8 (Frizzled Homolog 8) ELISA Kit
E-EL-H1534	Human PDCD1 (Programmed Cell Death Protein 1) ELISA Kit
E-EL-H1535	Human AP (Allopregnanolone) ELISA Kit
E-EL-H1536	Human CRHBP (Corticotropin Releasing Hormone Binding Protein) ELISA Kit
E-EL-H1537	Human αHSP (Alpha-Hemoglobin Stabilizing Protein) ELISA Kit

Cat.NO	Product Name
E-EL-H1538	Human MMP-28 (Matrix Metalloproteinase 28) ELISA Kit
E-EL-H1539	Human TLR4 (Toll-Like Receptor 4) ELISA Kit
E-EL-H1540	Human SPTAN1 (Alpha-Fodrin) ELISA Kit
E-EL-H1541	Human AZGP1 (Zinc-alpha-2-glycoprotein) ELISA Kit
E-EL-H1542	Human TSHRab (Thyroid Stimulating Hormone Receptor Antidoby) ELISA Kit
E-EL-H1543	Human MSLN (Mesothelin) ELISA Kit
E-EL-H1544	Human SOST (Sclerostin) ELISA Kit
E-EL-H1545	Human IRF5 (Interferon Regulatory Factor 5) ELISA Kit
E-EL-H1546	Human SIRT1 (Sirtuin 1) ELISA Kit
E-EL-H1547	Human PDCD1LG1 (Programmed Cell Death Protein 1 Ligand 1) ELISA Kit
E-EL-H1548	Human LSP (Liver Specific Lipoprotein) ELISA Kit
E-EL-H1549	Human AGPA (Anti-Glucoprotein 210 Antibody) ELISA Kit
E-EL-H1550	Human τ PK1 (Tau-Protein Kinase) ELISA Kit
E-EL-H1551	Human ACP5 (Tartrate Resistant Acid Phosphatase 5) ELISA Kit
E-EL-H1552	Human VWFCP/ADAMTS13 (Von Willebrand Factor Cleaving Protease) ELISA Kit
E-EL-H1553	Human MTR (5-Methyltetrahydrofolate Homocysteine Methyltransferase) ELISA Kit
E-EL-H1554	Human HMGB-1(High mobility group protein B1)ELISA Kit
E-EL-H1555	Human PGF (Placental Growth Factor) ELISA Kit
E-EL-H1556	Human FOXC2 (Forkhead Box Protein C2) ELISA Kit
E-EL-H1557	Human ESM1 (Endothelial Cell Specific Molecule 1) ELISA Kit

Cat.NO	Product Name
E-EL-H1558	Human Rac1 (Ras Related C3 Botulinum Toxin Substrate 1) ELISA Kit
E-EL-H1559	Human GHAb (Anti-Growth Hormone Antibody) ELISA Kit
E-EL-H1560	Human LILRB2 (Leukocyte Immunoglobulin Like Receptor Subfamily B Member 2) ELISA Kit
E-EL-H1561	Human OTR (Oxytocin Receptor) ELISA Kit
E-EL-H1562	Human GP (Anti-Glycoprotein Antibody) ELISA Kit
E-EL-H1563	Human PDIA6 (Protein Disulfide Isomerase A6) ELISA Kit
E-EL-H1564	Human NFE2L2 (Nuclear Factor, Erythroid Derived 2 Like 2) ELISA Kit
E-EL-H1565	Human HDL (High Density Lipoprotein) ELISA Kit
E-EL-H1566	Human DNP-Ab (Deoxyribonucleic Protein Antibody) ELISA Kit
E-EL-H1567	Human HBsAg (Hepatitis B surface antigen) ELISA Kit
E-EL-H1571	Human Pro-MMP-10 (Pro-Matrix Metalloproteinase 10) ELISA Kit
E-EL-H1572	Human Pro-MMP-13 (Pro-Matrix Metalloproteinase 13) ELISA Kit
E-EL-H1573	Human SIRPA (Signal Regulatory Protein Alpha) ELISA Kit
E-EL-H1574	Human PTX 3/TSG-14 (Pentraxin 3, Long) ELISA Kit
E-EL-H1575	Human PDGF-AA (Platelet Derived Growth Factor-AA) ELISA Kit
E-EL-H1576	Human PDGF-AB (Platelet Derived Growth Factor AB) ELISA Kit
E-EL-H1577	Human PDGF-BB (Platelet Derived Growth Factor-BB) ELISA Kit
E-EL-H1579	Human PCSK9 (Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit
E-EL-H1581	Human RBP4 (Retinol Binding Protein 4, Plasma) ELISA Kit
E-EL-H1582	Human RLN2 (Relaxin-2) ELISA Kit

Cat.NO	Product Name
E-EL-H1583	Human RBP1 (Retinol Binding Protein 1, Cellular) ELISA Kit
E-EL-H1584	Human Surv (Survivin) ELISA Kit
E-EL-H1585	Human RBP3 (Retinol Binding Protein 3, Interstitial) ELISA Kit
E-EL-H1586	Human TGF- α (Transforming Growth Factor α) ELISA Kit
E-EL-H1587	Human TGF- β 2 (Transforming Growth Factor β 2) ELISA Kit
E-EL-H1588	Human TPO (Thrombopoietin) ELISA Kit
E-EL-H1589	Human TSP-1 (Thrombospondin-1) ELISA Kit
E-EL-H1590	Human TSP-2 (Thrombospondin-2) ELISA Kit
E-EL-H1592	Human IL7R (Interleukin 7 Receptor) ELISA Kit
E-EL-H1593	Human TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit
E-EL-H1594	Human ODF (Osteoclast Differentiation Factor) ELISA Kit
E-EL-H1595	Human PLD1 (Phospholipase D1) ELISA Kit
E-EL-H1596	Human sTREM-1 (soluble Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
E-EL-H1597	Human TREM-1 (Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
E-EL-H1598	Human TSLP (Thymic Stromal Lymphopoietin) ELISA Kit
E-EL-H1599	Human PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit
E-EL-H1600	Human VEGF-C (Vascular Endothelial Growth Factor C) ELISA Kit
E-EL-H1601	Human VEGF-D (Vascular Endothelial Growth Factor D) ELISA Kit
E-EL-H1603	Human VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
E-EL-H1604	Human DBP (Vitamin D Binding Protein) ELISA Kit

Cat.NO	Product Name
E-EL-H1606	Human EMMPRIN/CD147 (Extracellular Matrix Metalloproteinase Inducer) ELISA Kit
E-EL-H1607	Human CD34 (Cluster Of Differentiation 34) ELISA Kit
E-EL-H1609	Human snRNP/Sm (Small Nuclear Ribonucleoprotein) ELISA Kit
E-EL-H1610	Human CYB5A (Cytochrome b5) ELISA Kit
E-EL-H1611	Human CG β (Chorionic Gonadotropin Beta Polypeptide) ELISA Kit
E-EL-H1612	Human sOX40L (soluble OX40L) ELISA Kit
E-EL-H1613	Human TAP (Transporter Associated with Antigen Processing) ELISA Kit
E-EL-H1614	Human BRAF (B-Raf Proto Oncogene Serine/Threonine Protein Kinase) ELISA Kit
E-EL-H1615	Human sST2 (soluble ST2) ELISA Kit
E-EL-H1616	Human GzmA (Granzyme A) ELISA Kit
E-EL-H1617	Human GzmB (Granzyme B) ELISA Kit
E-EL-H1618	Human GNLY (Granulysin) ELISA Kit
E-EL-H1619	Human GRN (Granulin) ELISA Kit
E-EL-H1620	Human S100P (S100 Calcium Binding Protein P) ELISA Kit
E-EL-H1621	Human CR2/CD21 (Complement Receptor 2) ELISA Kit
E-EL-H1622	Human CD28 (Cluster of Differentiation 28) ELISA Kit
E-EL-H1623	Human CD30L (Cluster of Differentiation 30 Ligand) ELISA Kit
E-EL-H1624	Human cADPRH/CD38 (Cyclic ADP Ribose Hydrolase) ELISA Kit
E-EL-H1625	Human F-PSA (Free Prostate-Specific Antigen) ELISA Kit
E-EL-H1626	Human LAB7-2/CD86 (B-Lymphocyte Activation Antigen B7-2) ELISA Kit

Cat.NO	Product Name
E-EL-H1627	Human CGb (Chorionic Gonadotropin Beta Polypeptide) ELISA Kit
E-EL-H1628	Human LGALS3 (Lectin-Galactoside Binding Soluble 3) ELISA Kit
E-EL-H1629	Human EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit
E-EL-H1630	Human CLDN5 (Claudin 5) ELISA Kit
E-EL-H1631	Human GRIN2B (Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 2B) ELISA Kit
E-EL-H1632	Human GRIN1 (Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 1) ELISA Kit
E-EL-H1633	Human GRIN2A (Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 2A) ELISA Kit
E-EL-H1634	Human FXN (Frataxin) ELISA Kit
E-EL-H1635	Human LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit
E-EL-H1636	Human ACTN4 (Actinin Alpha 4) ELISA Kit
E-EL-H1638	Human sAC (Soluble Adenylate Cyclase) ELISA Kit
E-EL-H1639	Human sFMC (Soluble Fibrin Monomer Complex) ELISA Kit
E-EL-H1640	Human PECAM1/CD31 (Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit
E-EL-H1641	Human COL18α1 (Collagen Type XVIII Alpha 1) ELISA Kit
E-EL-H1642	Human ApoB48 (Apolipoprotein B48) ELISA Kit
E-EL-H1643	Human sTNF-1 (Soluble Tumor Necrosis Factor 1) ELISA Kit
E-EL-H1644	Human MITF (Microphthalmia Associated Transcription Factor) ELISA Kit
E-EL-H1645	Human BMP-5 (Bone Morphogenetic Protein 5) ELISA Kit
E-EL-H1646	Human sTfR1 (Soluble Transferrin Receptor1) ELISA Kit
E-EL-H1647	Human IL-12 p35 (Interleukin 12 p35) ELISA Kit

Cat.NO	Product Name
E-EL-H1648	Human IL-25 (Interleukin 25) ELISA Kit
E-EL-H1649	Human VacA (Vacuolating Cytotoxin A) ELISA Kit
E-EL-H1650	Human IL-34 (Interleukin 34) ELISA Kit
E-EL-H1651	Human RLN1 (Relaxin-1) ELISA Kit
E-EL-H1652	Human RLN3 (Relaxin-3) ELISA Kit
E-EL-H1653	Human Tyk-2 (Tyrosine Kinase 2) ELISA Kit
E-EL-H1654	Human TH (Tyrosine Hydroxylase) ELISA Kit
E-EL-H1655	Human mTOR (Mammalian Target of Rapamycin) ELISA Kit
E-EL-H1656	Human PSPN (Persephin) ELISA Kit
E-EL-H1657	Human TGFβ1 (Transforming Growth Factor Beta Induced Protein) ELISA Kit
E-EL-H1658	Human TSC22 (Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
E-EL-H1659	Human HLA-A (Leukocyte Antigen A) ELISA Kit
E-EL-H1660	Human PKCη (Protein Kinase C Eta) ELISA Kit
E-EL-H1661	Human UCP-1 (Uncoupling Protein 1, Mitochondrial) ELISA Kit
E-EL-H1662	Human IL18BP (Interleukin 18 Binding Protein) ELISA Kit
E-EL-H1663	Human MHCG (Major Histocompatibility Complex Class I G) ELISA Kit
E-EL-H1665	Human TYR (Tyrosinase) ELISA Kit
E-EL-H1666	Human NFκB2 (Nuclear Factor Kappa B2) ELISA Kit
E-EL-H1667	Human FGL1 (Fibrinogen Like Protein 1) ELISA Kit
E-EL-H1668	Human FBLN2 (Fibulin 2) ELISA Kit

Cat.NO	Product Name
E-EL-H1669	Human PDGFA (Platelet Derived Growth Factor Subunit A) ELISA Kit
E-EL-H1670	Human PDGFB (Platelet Derived Growth Factor Subunit B) ELISA Kit
E-EL-H1671	Human ORM2 (Orosomucoid 2) ELISA Kit
E-EL-H1672	Human CG (Cryoglobulin) ELISA Kit
E-EL-H1673	Human FBLN3 (Fibulin 3) ELISA Kit
E-EL-H1674	Human GM-CSF Ab (Anti-Granulocyte Macrophage Colony Stimulating Factor Antibody) ELISA Kit
E-EL-H1675	Human PRG4 (Proteoglycan 4) ELISA Kit
E-EL-H1676	Human GSANA (Granulocyte Specific Anti-Nuclear Antibody) ELISA Kit
E-EL-H1677	Human FBLN1 (Fibulin 1) ELISA Kit
E-EL-H1678	Human Cx (Connexin) ELISA Kit
E-EL-H1679	Human CX37 (Connexin 37) ELISA Kit
E-EL-H1680	Human CX40 (Connexin 40) ELISA Kit
E-EL-H1681	Human CX43 (Connexin 43) ELISA Kit
E-EL-H1682	Human SDC4 (Syndecan 4) ELISA Kit
E-EL-H1684	Human LAP (Leucine Aminopeptidase) ELISA Kit
E-EL-H1686	Human LT β (Lymphotoxin Beta) ELISA Kit
E-EL-H1687	Human LFA-1 (Lymphocyte Function Associated Antigen 1) ELISA Kit
E-EL-H1688	Human LFA-2/CD2 (Lymphocyte Function Associated Antigen 2) ELISA Kit
E-EL-H1689	Human LFA-3 (Lymphocyte Function Associated Antigen 3) ELISA Kit
E-EL-H1690	Human Ly6a (Lymphocyte Antigen 6 Complex Locus A) ELISA Kit

Cat.NO	Product Name
E-EL-H1691	Human Lptn/XCL1 (Lymphotactin) ELISA Kit
E-EL-H1694	Human PGAM2 (Phosphoglycerate Mutase 2) ELISA Kit
E-EL-H1695	Human PFKM (Phosphofructokinase, Muscle) ELISA Kit
E-EL-H1696	Human P-PKC (Phospho Protein Kinase C) ELISA Kit
E-EL-H1697	Human plkB α (Phospho-Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
E-EL-H1698	Human pERK(Phospho Extracellular Signal Regulated Kinase)ELISA Kit
E-EL-H1699	Human A-CoA (Acetyl Coenzyme A) ELISA Kit
E-EL-H1700	Human PI3K (Phosphoinositide-3-Kinase) ELISA Kit
E-EL-H1701	Human PGM (Phosphoglucomutase) ELISA Kit
E-EL-H1702	Human PCK1 (Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
E-EL-H1703	Human PLAA (Phospholipase A2 Activating Protein) ELISA Kit
E-EL-H1704	Human PLA1 (Phospholipase A1) ELISA Kit
E-EL-H1706	Human GLO1 (Glyoxalase I)ELISA Kit
E-EL-H1707	Human PLSCR1 (Phospholipid Scramblase 1) ELISA Kit
E-EL-H1708	Human PLSCR2 (Phospholipid Scramblase 2) ELISA Kit
E-EL-H1709	Human PI3K2 α (Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit
E-EL-H1710	Human GPC1 (Glypican 1) ELISA Kit
E-EL-H1711	Human GPC2 (Glypican 2) ELISA Kit
E-EL-H1712	Human GPC3 (Glypican 3) ELISA Kit
E-EL-H1713	Human GPC4 (Glypican 4) ELISA Kit

Cat.NO	Product Name
E-EL-H1715	Human PS (Phosphatidylserine) ELISA Kit
E-EL-H1716	Human PLTP(Phospholipid Transfer Protein)ELISA Kit
E-EL-H1717	Human PLC (Phospholipase C) ELISA Kit
E-EL-H1718	Human PLC γ 1 (Phospholipase C Gamma 1) ELISA Kit
E-EL-H1719	Human PLC β 1/PI-PLC (Phospholipase C Beta 1, Phosphoinositide Specific) ELISA Kit
E-EL-H1720	Human FUR (Furin) ELISA Kit
E-EL-H1721	Human SCCA1 (Squamous Cell Carcinoma Antigen 1) ELISA Kit
E-EL-H1722	Human SCCA2 (Squamous Cell Carcinoma Antigen 2) ELISA Kit
E-EL-H1723	Human pPLA2 (Phospholipase A2, Pancreas) ELISA Kit
E-EL-H1725	Human DS (Dermatan Sulfate) ELISA Kit
E-EL-H1726	Human MS(Melatonin Sulfate)ELISA Kit
E-EL-H1727	Human Trx (Thioredoxin) ELISA Kit
E-EL-H1728	Human TrxR1 (Thioredoxin Reductase 1) ELISA Kit
E-EL-H1729	Human TBP2 (Thioredoxin Binding Protein 2) ELISA Kit
E-EL-H1730	Human IAP (immunosuppressive acidic protein) ELISA Kit
E-EL-H1732	Human OVA slgE (Ovalbumin Specific IgE) ELISA Kit
E-EL-H1733	Human OVA slgG (Ovalbumin Specific IgG) ELISA Kit
E-EL-H1739	Human MPA (Mycophenolic Acid) ELISA Kit
E-EL-H1741	Human Fc α R/CD89 (Receptor I for the Fc Fragment of IgA) ELISA Kit
E-EL-H1742	Human Fc ϵ R I (Receptor I for the Fc Fragment of IgE) ELISA Kit

Cat.NO	Product Name
E-EL-H1743	Human FcγR II/CD32 (Receptor II for the Fc Fragment of IgG) ELISA Kit
E-EL-H1744	Human FcγR3B (Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit
E-EL-H1745	Human Fcγ (Fc Fragment of IgG) ELISA Kit
E-EL-H1746	Human IgG1 (Immunoglobulin G1) ELISA Kit
E-EL-H1747	Human IgG2 (Immunoglobulin G2) ELISA Kit
E-EL-H1748	Human IgG3 (Immunoglobulin G3) ELISA Kit
E-EL-H1749	Human IgG4 (Immunoglobulin G4) ELISA Kit
E-EL-H1750	Human MCP/CD46(Membrane Cofactor Protein)ELISA Kit
E-EL-H1751	Human TSP-4 (Thrombospondin-4) ELISA Kit
E-EL-H1752	Human Ponticulin ELISA Kit
E-EL-H1753	Human TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit
E-EL-H1754	Human TDT (Terminal Deoxynucleotidyl Transferase) ELISA Kit
E-EL-H1755	Human SMN1 (Survival of Motor Neuron 1, Telomeric) ELISA Kit
E-EL-H1757	Human EM2 (Endomorphin 2) ELISA Kit
E-EL-H1758	Human ETRA (Endothelin Receptor A) ELISA Kit
E-EL-H1759	Human ECE1 (Endothelin Converting Enzyme 1) ELISA Kit
E-EL-H1760	Human REPIN1 (Replication Initiator 1) ELISA Kit
E-EL-H1761	Human EL/LIPG (Endothelial Lipase) ELISA Kit
E-EL-H1762	Human VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
E-EL-H1763	Human VF (Visfatin) ELISA Kit

Cat.NO	Product Name
E-EL-H1764	Human ERN1 (Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit
E-EL-H1766	Human CFTR (Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit
E-EL-H1768	Human NGB (Neuroglobin) ELISA Kit
E-EL-H1769	Human PGD2S (Prostaglandin D2 Synthase, Brain) ELISA Kit
E-EL-H1770	Human MAdCAM1 (Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit
E-EL-H1771	Human FAK (Focal Adhesion Kinase) ELISA Kit
E-EL-H1772	Human OCT (Ornithine Carbamoyl Transferase) ELISA Kit
E-EL-H1773	Human ODC (Ornithine Decarboxylase) ELISA Kit
E-EL-H1774	Human GEF (Guanylate Exchange Factor) ELISA Kit
E-EL-H1775	Human GBP4 (Guanylate Binding Protein 4) ELISA Kit
E-EL-H1776	Human PIBF (Progesterone-induced-blocking factor) ELISA Kit
E-EL-H1778	Human UGCG (UDP-glucose ceramide glucosyltransferase) ELISA Kit
E-EL-H1779	Human PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit
E-EL-H1781	Human UPP1 (Uridine Phosphorylase 1) ELISA Kit
E-EL-H1782	Human UreC (Urease Related Protein C) ELISA Kit
E-EL-H1783	Human ETRB (Endothelin Receptor B) ELISA Kit
E-EL-H1786	Human GS (Gelsolin) ELISA Kit
E-EL-H1787	Human TM (Thrombin)ELISA Kit
E-EL-H1788	Human TAT (Thrombin/Antithrombin Complex) ELISA Kit
E-EL-H1789	Human TSP1 (Thrombin Sensitive Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-H1790	Human TR (Thrombin Receptor) ELISA Kit
E-EL-H1792	Human PT (Prothrombin) ELISA Kit
E-EL-H1793	Human F1+2 (Prothrombin Fragment 1+2) ELISA Kit
E-EL-H1794	Human KERA (Keratocan) ELISA Kit
E-EL-H1795	Human VCL (Vinculin) ELISA Kit
E-EL-H1798	Human DEF α 5 (Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
E-EL-H1799	Human FSA (Fetal Sulfoslycoprotein Antigen) ELISA Kit
E-EL-H1800	Human EED (Embryonic Ectoderm Development) ELISA Kit
E-EL-H1801	Human CTAGE (Cutaneous T-Cell Lymphocyte Associated Antigen) ELISA Kit
E-EL-H1802	Human DPT (Dermatopontin) ELISA Kit
E-EL-H1803	Human SRF/IF (Skin Reactive Factor/Inflammatory Factor) ELISA Kit
E-EL-H1804	Human CLA/CTAGE (Cutaneous Lymphocyte Associated Antigen) ELISA Kit
E-EL-H1808	Human SYK (Spleen Tyrosine Kinase) ELISA Kit
E-EL-H1809	Human SMM (Smooth Muscle Myosin) ELISA Kit
E-EL-H1810	Human PDPN (Podoplanin) ELISA Kit
E-EL-H1812	Human OSF (Osteoclast Stimulating Factor 1) ELISA Kit
E-EL-H1813	Human OSCAR (Osteoclast Associated Receptor) ELISA Kit
E-EL-H1814	Human IgM (Immunoglobulin M) ELISA Kit
E-EL-H1816	Human G6PD (Glucose 6 Phosphate Dehydrogenase) ELISA Kit
E-EL-H1817	Human GPI (Glucose 6 Phosphate Isomerase) ELISA Kit

Cat.NO	Product Name
E-EL-H1818	Human GCK (Glucokinase) ELISA Kit
E-EL-H1819	Human GKRP (Glucokinase Regulatory Protein) ELISA Kit
E-EL-H1820	Human GOD (Glucose Oxidase) ELISA Kit
E-EL-H1822	Human GLUT1 (Glucose Transporter 1) ELISA Kit
E-EL-H1823	Human GLUT2 (Glucose Transporter 2) ELISA Kit
E-EL-H1824	Human GLUT3 (Glucose Transporter 3) ELISA Kit
E-EL-H1825	Human GLUT4 (Glucose Transporter 4) ELISA Kit
E-EL-H1826	Human GTC (General Transcription Complex) ELISA Kit
E-EL-H1827	Human GTF2B (General Transcription Factor IIB) ELISA Kit
E-EL-H1828	Human GTF2F1 (General Transcription Factor IIF, Polypeptide 1) ELISA Kit
E-EL-H1830	Human PK2 (Prokineticin 2) ELISA Kit
E-EL-H1831	Human PKR1 (Prokineticin Receptor 1) ELISA Kit
E-EL-H1832	Human PKR2 (Prokineticin Receptor 2) ELISA Kit
E-EL-H1833	Human PCP (Procollagen C Peptidase) ELISA Kit
E-EL-H1835	Human Cys-C (high-sensitivity Cystatin C) ELISA Kit
E-EL-H1836	Human PGDS (Prostaglandin D Synthase) ELISA Kit
E-EL-H1838	Human EP2 (Prostaglandin E Receptor 2) ELISA Kit
E-EL-H1839	Human EP3 (Prostaglandin E Receptor 3) ELISA Kit
E-EL-H1840	Human PGF (Prostaglandin F) ELISA Kit
E-EL-H1841	Human PGF2 α (Prostaglandin F2 α) ELISA Kit

Cat.NO	Product Name
E-EL-H1842	Human FP (Prostaglandin F Receptor) ELISA Kit
E-EL-H1843	Human PTGIS (Prostaglandin I Synthase) ELISA Kit
E-EL-H1844	Human PTGES2 (Prostaglandin E Synthase 2) ELISA Kit
E-EL-H1845	Human PTGS1 (Prostaglandin Endoperoxide Synthase 1) ELISA Kit
E-EL-H1846	Human PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit
E-EL-H1847	Human ACP3/PAP (Prostatic Acid Phosphatase 3) ELISA Kit
E-EL-H1848	Human NT-ProANP (N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit
E-EL-H1849	Human Dyn (Big Dynorphin) ELISA Kit
E-EL-H1850	Human DSP (Desmoplakin) ELISA Kit
E-EL-H1851	Human DSG-E1 (Desmogleins 1) ELISA Kit
E-EL-H1852	Human SM/SPH (Sphingomyelin) ELISA Kit
E-EL-H1853	Human ERCC1 (Excision Repair Cross Complementing Rodent Repair Deficiency Complementation 1) ELISA Kit
E-EL-H1854	Human HJV/HFE2 (Hemojuvelin/Hemochromatosis Type 2) ELISA Kit
E-EL-H1855	Human KLK1 (Kallikrein 1) ELISA Kit
E-EL-H1856	Human LXR α (Liver X Receptor Alpha) ELISA Kit
E-EL-H1857	Human dDAVP (1-Desamino-8D-Arginine Vasopressin) ELISA Kit
E-EL-H1859	Human HSP β 6/HSP-20 (Heat Shock Protein Beta 6) ELISA Kit
E-EL-H1860	Human HSP-27/HSPB1 (Heat Shock Protein 27) ELISA Kit
E-EL-H1861	Human HSP-40 (Heat Shock Protein 40) ELISA Kit
E-EL-H1862	Human HSP-60/HSPD1 (Heat Shock Protein 60) ELISA Kit

Cat.NO	Product Name
E-EL-H1863	Human HSP-70 (Heat Shock Protein 70) ELISA Kit
E-EL-H1864	Human HSP-90 (Heat Shock Protein 90) ELISA Kit
E-EL-H1865	Human HSP gp96 (Heat Shock Protein glycoprotein 96) ELISA Kit
E-EL-H1866	Human HSF2 (Heat Shock Transcription Factor 2) ELISA Kit
E-EL-H1867	Human HSF1 (Heat Shock Transcription Factor 1) ELISA Kit
E-EL-H1868	Human PSG1/SP1 (Pregnancy Specific Beta-1-Glycoprotein 1) ELISA Kit
E-EL-H1869	Human LZM (Lysozyme) ELISA Kit
E-EL-H1870	Human TERF2 (Telomeric Repeat Binding Factor 2) ELISA Kit
E-EL-H1872	Human TFF1 (Trefoil Factor 1) ELISA Kit
E-EL-H1873	Human TFF2 (Trefoil Factor 2) ELISA Kit
E-EL-H1874	Human WARS (Tryptophanyl tRNA Synthetase) ELISA Kit
E-EL-H1875	Human PEDF(Pigment Epithelium Derived Factor)ELISA Kit
E-EL-H1876	Human SDH (Sorbitol Dehydrogenase) ELISA Kit
E-EL-H1877	Human CCR2 (Chemokine C-C-Motif Receptor 2) ELISA Kit
E-EL-H1878	Human EMA (Epithelial Membrane Antigen) ELISA Kit
E-EL-H1880	Human Ep-CAM/CD326 (Epithelial Cell Adhesion Molecule) ELISA Kit
E-EL-H1881	Human BGN (Biglycan) ELISA Kit
E-EL-H1883	Human NRG-1 (Neuregulin 1) ELISA Kit
E-EL-H1884	Human NKA (Neurokinin A) ELISA Kit
E-EL-H1885	Human NKB (Neurokinin B) ELISA Kit

Cat.NO	Product Name
E-EL-H1886	Human NT (Neurotensin) ELISA Kit
E-EL-H1887	Human GMFB (Glia Maturation Factor, Beta) ELISA Kit
E-EL-H1888	Human GFAP (Glial Fibrillary Acidic Protein) ELISA Kit
E-EL-H1889	Human CL (Cardiolipin) ELISA Kit
E-EL-H1890	Human HK1 (Hexokinase 1) ELISA Kit
E-EL-H1891	Human NF (Neurofilament) ELISA Kit
E-EL-H1892	Human P2 (Myelin Protein 2) ELISA Kit
E-EL-H1893	Human NPY (Neuropeptide Y) ELISA Kit
E-EL-H1894	Human NCAM/CD56 (Neural Cell Adhesion Molecule) ELISA Kit
E-EL-H1895	Human L1CAM/CD171 (L1-Cell Adhesion Molecule) ELISA Kit
E-EL-H1896	Human NT-3 (Neurotrophin-3) ELISA Kit
E-EL-H1897	Human NT-4 (Neurotrophin 4) ELISA Kit
E-EL-H1898	Human NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
E-EL-H1899	Human NTRK3 (Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
E-EL-H1900	Human NAIP (Neuronal Apoptosis Inhibitory Protein) ELISA Kit
E-EL-H1901	Human NPHN (Nephrin) ELISA Kit
E-EL-H1903	Human GTEAGE (Glomerular Tissue Extract Advanced Glycosylation End Products) ELISA Kit
E-EL-H1904	Human GROβ (Growth Regulated Oncogene Beta) ELISA Kit
E-EL-H1905	Human GROγ (Growth Regulated Oncogene Gamma) ELISA Kit
E-EL-H1906	Human GDF1 (Growth Differentiation Factor 1) ELISA Kit

Cat.NO	Product Name
E-EL-H1907	Human GDF10 (Growth Differentiation Factor 10) ELISA Kit
E-EL-H1908	Human GDF11 (Growth Differentiation Factor 11) ELISA Kit
E-EL-H1909	Human GDF2 (Growth Differentiation Factor 2) ELISA Kit
E-EL-H1910	Human GDF3 (Growth Differentiation Factor 3) ELISA Kit
E-EL-H1911	Human GDF5 (Growth Differentiation Factor 5) ELISA Kit
E-EL-H1912	Human GDF6 (Growth Differentiation Factor 6) ELISA Kit
E-EL-H1913	Human GDF7 (Growth Differentiation Factor 7) ELISA Kit
E-EL-H1914	Human GDF9 (Growth Differentiation Factor 9) ELISA Kit
E-EL-H1915	Human HK2 (Hexokinase 2) ELISA Kit
E-EL-H1916	Human GH2 (Growth Hormone 2) ELISA Kit
E-EL-H1917	Human GHR/GHBP (Growth Hormone Binding Protein) ELISA Kit
E-EL-H1918	Human OSBPL1A (Oxysterol Binding Protein Like Protein 1A) ELISA Kit
E-EL-H1919	Human GHRL (Ghrelin) ELISA Kit
E-EL-H1920	Human ARG2 (Arginase II) ELISA Kit
E-EL-H1921	Human GHITM (Growth Hormone Inducible Transmembrane Protein) ELISA Kit
E-EL-H1922	Human GAP43 (Growth Associated Protein 43) ELISA Kit
E-EL-H1923	Human SST (Somatostatin) ELISA Kit
E-EL-H1924	Human SSTR1 (Somatostatin Receptor 1) ELISA Kit
E-EL-H1925	Human Grb10 (Growth Factor Receptor Bound Protein 10) ELISA Kit
E-EL-H1926	Human Grb14 (Growth Factor Receptor Bound Protein 14) ELISA Kit

Cat.NO	Product Name
E-EL-H1927	Human Grb2 (Growth Factor Receptor Bound Protein 2) ELISA Kit
E-EL-H1928	Human Grb7 (Growth Factor Receptor Bound Protein 7) ELISA Kit
E-EL-H1929	Human OXB (Orexin B) ELISA Kit
E-EL-H1931	Human RLR (Retinoic Acid Inducible Gene/RIG-Like Receptor) ELISA Kit
E-EL-H1932	Human RB1 (Retinoblastoma Protein 1) ELISA Kit
E-EL-H1933	Human PRB (Retinoblastoma Tumor Suppressor Protein) ELISA Kit
E-EL-H1934	Human CYPA (Cyclophilin A) ELISA Kit
E-EL-H1935	Human CYPB (Cyclophilin B) ELISA Kit
E-EL-H1936	Human CYPD (Cyclophilin D) ELISA Kit
E-EL-H1937	Human NRP2 (Neuropilin 2) ELISA Kit
E-EL-H1938	Human CS (Citrate Synthase) ELISA Kit
E-EL-H1939	Human OGG1 (8-Oxoguanine Glycosylase 1) ELISA Kit
E-EL-H1940	Human ECFR (Eosinophil Chemotactic Factor Receptor) ELISA Kit
E-EL-H1941	Human EAR1 (Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit
E-EL-H1942	Human EAR12 (Eosinophil Associated Ribonuclease A family Member 12) ELISA Kit
E-EL-H1943	Human EAR2 (Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit
E-EL-H1944	Human EAR3 (Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit
E-EL-H1945	Human NOG (Noggin) ELISA Kit
E-EL-H1946	Human NE/ELA2 (Elastase 2, Neutrophil) ELISA Kit
E-EL-H1947	Human PLN (Phospholamban) ELISA Kit

Cat.NO	Product Name
E-EL-H1948	Human RIPK2 (Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit
E-EL-H1949	Human TNK2 (Tyrosine Kinase, Non Receptor 2) ELISA Kit
E-EL-H1950	Human DAF (Decay Accelerating Factor) ELISA Kit
E-EL-H1952	Human STK39 (Serine/Threonine Kinase 39) ELISA Kit
E-EL-H1953	Human PRSS (Serine Protease) ELISA Kit
E-EL-H1954	Human TNNI3K (Serine/threonine-protein kinase TNNI3K) ELISA Kit
E-EL-H1955	Human MAP2K6/MKK6 (Mitogen Activated Protein Kinase Kinase 6) ELISA Kit
E-EL-H1956	Human MAPK (Mitogen Activated Protein Kinase) ELISA Kit
E-EL-H1957	Human FHA (Filamentous Hemagglutinin) ELISA Kit
E-EL-H1958	Human DAP6 (Death Associated Protein 6) ELISA Kit
E-EL-H1959	Human RXFP1 (Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
E-EL-H1961	Human TACR1 (Tachykinin Receptor 1) ELISA Kit
E-EL-H1962	Human TACR2 (Tachykinin Receptor 2) ELISA Kit
E-EL-H1963	Human TACR3 (Tachykinin Receptor 3) ELISA Kit
E-EL-H1964	Human Anti-MPO (Anti-Myeloperoxidase Antibody) ELISA Kit
E-EL-H1965	Human PLP1 (Proteolipid Protein 1, Myelin) ELISA Kit
E-EL-H1967	Human MyD88 (Myeloid Differentiation Factor 88) ELISA Kit
E-EL-H1968	Human POT1 (Protection of Telomeres Protein 1) ELISA Kit
E-EL-H1969	Human HBF (Fetal Hemoglobin) ELISA Kit
E-EL-H1970	Human PR3 (Proteinase 3) ELISA Kit

Cat.NO	Product Name
E-EL-H1971	Human PL (Placental Lactogen) ELISA Kit
E-EL-H1972	Human PP (Placental Protein) ELISA Kit
E-EL-H1973	Human PP13(Placental Protein13)ELISA Kit
E-EL-H1974	Human PP14 (Placental Protein14) ELISA Kit
E-EL-H1975	Human PRI (Placental Ribonuclease Inhibitor) ELISA Kit
E-EL-H1976	Human PLAP/ALPP(Placental Alkaline Phosphatase)ELISA Kit
E-EL-H1977	Human CPT1A (Carnitine Palmitoyltransferase 1A, Liver) ELISA Kit
E-EL-H1978	Human CPT1B (Carnitine Palmitoyltransferase 1B, Muscle) ELISA Kit
E-EL-H1979	Human FETUB (Fetuin B) ELISA Kit
E-EL-H1980	Human OSBPL8 (Oxysterol Binding Protein Like Protein 8) ELISA Kit
E-EL-H1981	Human CFHR3 (Complement Factor H Related Protein 3) ELISA Kit
E-EL-H1983	Human PGRP-S (Peptidoglycan Recognition Protein-S) ELISA Kit
E-EL-H1984	Human PPIL1 (Peptidylprolyl Isomerase/Cyclophilin Like 1) ELISA Kit
E-EL-H1989	Human OB (Obestatin) ELISA Kit
E-EL-H1990	Human Gp49a (Glycoprotein 49A) ELISA Kit
E-EL-H1991	Human GP5 (Glycoprotein V, Platelet) ELISA Kit
E-EL-H1992	Human GP6 (Glycoprotein VI, Platelet) ELISA Kit
E-EL-H1993	Human GA (Glycated Albumin) ELISA Kit
E-EL-H1994	Human GSP (Glycosylated Serum Protein) ELISA Kit
E-EL-H1995	Human HbA1c (Hemoglobin A1c) ELISA Kit

Cat.NO	Product Name
E-EL-H1996	Human GlyCAM1 (Glycosylation Dependent Cell Adhesion Molecule 1) ELISA Kit
E-EL-H1997	Human GPI (Glycophosphatidylinositol) ELISA Kit
E-EL-H1998	Human GR α (Glucocorticoid Receptor Alpha) ELISA Kit
E-EL-H1999	Human GR β (Glucocorticoid Receptor Beta) ELISA Kit
E-EL-H2000	Human GSK3 β (Glycogen Synthase Kinase 3 Beta) ELISA Kit
E-EL-H2001	Human GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit
E-EL-H2002	Human A-GHRL (Acylated Ghrelin) ELISA Kit
E-EL-H2004	Human COX (Cytochrome C Oxidase) ELISA Kit
E-EL-H2005	Human GLTP (Glycolipid Transfer Protein) ELISA Kit
E-EL-H2006	Human SMAF (Specific Macrophage Arming Factor) ELISA Kit
E-EL-H2008	Human CFHR5 (Complement Factor H Related Protein 5) ELISA Kit
E-EL-H2009	Human FTL (Ferritin, Light Polypeptide) ELISA Kit
E-EL-H2010	Human FTH (Ferritin, Heavy Polypeptide) ELISA Kit
E-EL-H2011	Human PIR (Pirin) ELISA Kit
E-EL-H2012	Human MRP1 (Multidrug Resistance Associated Protein 1) ELISA Kit
E-EL-H2013	Human SYNPO (Synaptopodin) ELISA Kit
E-EL-H2014	Human SYP (Synaptophysin) ELISA Kit
E-EL-H2015	Human SAP25 (Synaptosome Associated Protein 25) ELISA Kit
E-EL-H2016	Human MT (Melatonin) ELISA Kit
E-EL-H2017	Human DHEA (Dehydroepiandrosterone) ELISA Kit

Cat.NO	Product Name
E-EL-H2018	Human DHEA-S7 (Dehydroepiandrosterone S7) ELISA Kit
E-EL-H2019	Human DHEA-S (Dehydroepiandrosterone Sulfate) ELISA Kit
E-EL-H2020	Human DPD (Deoxypyridinoline) ELISA Kit
E-EL-H2021	Human DNase-I (Deoxyribonuclease I) ELISA Kit
E-EL-H2022	Human TERF1 (Telomeric Repeat Binding Factor 1) ELISA Kit
E-EL-H2024	Human SA (Sialic Acid) ELISA Kit
E-EL-H2025	Human SIGLEC12 (Sialic Acid Binding Ig Like Lectin 12) ELISA Kit
E-EL-H2026	Human SIGLEC3 (Sialic Acid Binding Ig Like Lectin 3) ELISA Kit
E-EL-H2027	Human Extein ELISA Kit
E-EL-H2028	Human ITLN1 (Intelectin 1/Omentin) ELISA Kit
E-EL-H2029	Human TUB α 3C (Tubulin Alpha 3c)ELISA Kit
E-EL-H2030	Human CFHR4 (Complement Factor H Related Protein 4) ELISA Kit
E-EL-H2031	Human TPPP (Tubulin Polymerization Promoting Protein) ELISA Kit
E-EL-H2032	Human STMN1 (Stathmin 1) ELISA Kit
E-EL-H2033	Human DCX (Doublecortin) ELISA Kit
E-EL-H2034	Human MAP-2 (Microtubule Associated Protein 2) ELISA Kit
E-EL-H2035	Human MTF (Microtransferrin) ELISA Kit
E-EL-H2043	Human VDR (Vitamin D Receptor) ELISA Kit
E-EL-H2046	Human CFHR2 (Complement Factor H Related Protein 2) ELISA Kit
E-EL-H2047	Human UST2 (Urotensin 2) ELISA Kit

Cat.NO	Product Name
E-EL-H2049	Human CFHR1 (Complement Factor H Related Protein 1) ELISA Kit
E-EL-H2050	Human PP (Pepsin) ELISA Kit
E-EL-H2051	Human PGA (Pepsinogen A) ELISA Kit
E-EL-H2052	Human PGC (Pepsinogen C) ELISA Kit
E-EL-H2053	Human GKN1 (Gastrokine 1) ELISA Kit
E-EL-H2054	Human GKN2 (Gastrokine 2) ELISA Kit
E-EL-H2055	Human MTL (Motilin) ELISA Kit
E-EL-H2056	Human VIPR1 (Vasoactive Intestinal Peptide Receptor) ELISA Kit
E-EL-H2057	Human GT-34 (Gastrin 34) ELISA Kit
E-EL-H2059	Human ProGRP (Pro-Gastrin Releasing Peptide) ELISA Kit
E-EL-H2060	Human GT (Gastrin) ELISA Kit
E-EL-H2061	Human GIP (Gastric Inhibitory Polypeptide) ELISA Kit
E-EL-H2062	Human GM (Gastric Mucin) ELISA Kit
E-EL-H2063	Human MFGE8 (Milk Fat Globule EGF Factor 8) ELISA Kit
E-EL-H2064	Human PTD (Pentosidine) ELISA Kit
E-EL-H2065	Human SEP1 (Selenoprotein 1) ELISA Kit
E-EL-H2066	Human MMP23B (Matrix Metalloproteinase 23B) ELISA Kit
E-EL-H2068	Human CAPA (Cytotoxin Associated Protein A) ELISA Kit
E-EL-H2069	Human CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) ELISA Kit
E-EL-H2070	Human CK-13/KRT13 (Cytokeratin 13) ELISA Kit

Cat.NO	Product Name
E-EL-H2071	Human CK-17/KRT17 (Cytokeratin 17) ELISA Kit
E-EL-H2072	Human CK-18/KRT18 (Cytokeratin 18) ELISA Kit
E-EL-H2074	Human CK18-M65 (Cytokeratin 18-M65) ELISA Kit
E-EL-H2075	Human CK-19/KRT19 (Cytokeratin 19) ELISA Kit
E-EL-H2076	Human CK-20/KRT20 (Cytokeratin 20) ELISA Kit
E-EL-H2077	Human CYFRA21-1 (Cytokeratin Fragment Antigen 21-1) ELISA Kit
E-EL-H2079	Human CVL/EZR (Cytovillin/Ezrin) ELISA Kit
E-EL-H2080	Human CYB561 (Cytochrome b-561) ELISA Kit
E-EL-H2081	Human CYP7A1 (Cytochrome P450 7A1) ELISA Kit
E-EL-H2082	Human CYP21B (Cytochrome P450 21B) ELISA Kit
E-EL-H2083	Human ECM1 (Extracellular Matrix Protein 1) ELISA Kit
E-EL-H2084	Human MEPE (Matrix Extracellular Phosphoglycoprotein) ELISA Kit
E-EL-H2085	Human ERK (Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-H2086	Human ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit
E-EL-H2087	Human ERK2 (Extracellular Signal Regulated Kinase 2) ELISA Kit
E-EL-H2088	Human ERK3 (Extracellular Signal Regulated Kinase 3) ELISA Kit
E-EL-H2089	Human SOCS3 (Suppressors of Cytokine Signaling 3) ELISA Kit
E-EL-H2090	Human CYTH1 (Cytohesin 1) ELISA Kit
E-EL-H2091	Human CYTIP (Cytohesin Interacting Protein) ELISA Kit
E-EL-H2092	Human CCND1 (Cyclin-D1) ELISA Kit

Cat.NO	Product Name
E-EL-H2093	Human CCND2 (Cyclin-D2) ELISA Kit
E-EL-H2094	Human CCND3 (Cyclin-D3) ELISA Kit
E-EL-H2095	Human FLN α (Filamin A, Alpha) ELISA Kit
E-EL-H2096	Human FLN β (Filamin B, Beta) ELISA Kit
E-EL-H2097	Human FLNC (Filamin C, Gamma) ELISA Kit
E-EL-H2098	Human FMOD (Fibromodulin) ELISA Kit
E-EL-H2099	Human FLRT2 (Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit
E-EL-H2100	Human Fbn (Fibrinolysin) ELISA Kit
E-EL-H2101	Human PAP (Plasmin-Antiplasmin Complex) ELISA Kit
E-EL-H2102	Human Plg (Plasminogen) ELISA Kit
E-EL-H2103	Human PLGA (Plasminogen Activator) ELISA Kit
E-EL-H2104	Human PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit
E-EL-H2105	Human PAI2 (Plasminogen Activator Inhibitor 2)ELISA Kit
E-EL-H2106	Human tPA (Plasminogen Activator, Tissue) ELISA Kit
E-EL-H2107	Human TAFI (Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit
E-EL-H2108	Human FBL (Fibrillarlin) ELISA Kit
E-EL-H2109	Human GHR (Growth Hormone Receptor) ELISA Kit
E-EL-H2110	Human FG α (Fibrinogen Alpha) ELISA Kit
E-EL-H2111	Human FG β (Fibrinogen Beta) ELISA Kit
E-EL-H2112	Human FGL2 (Fibrinogen Like Protein 2) ELISA Kit

Cat.NO	Product Name
E-EL-H2113	Human FRA (Fibronectin Related Antigen) ELISA Kit
E-EL-H2114	Human FSP (Fibroblast Surface Protein) ELISA Kit
E-EL-H2115	Human SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-H2116	Human CKMT1A (Creatine Kinase, Mitochondrial 1A) ELISA Kit
E-EL-H2117	Human CKMT1B (Creatine Kinase, Mitochondrial 1B) ELISA Kit
E-EL-H2118	Human KRT10 (Keratin, Type I Cytoskeletal 10) ELISA Kit
E-EL-H2119	Human SNRPD1 (Small Nuclear Ribonucleoprotein Polypeptide D1) ELISA Kit
E-EL-H2120	Human SNRPC (Small Nuclear Ribonucleoprotein Polypeptide C) ELISA Kit
E-EL-H2121	Human SBEM (Small Breast Epithelial Mucin) ELISA Kit
E-EL-H2122	Human CMR (Costimulatory Molecules Receptor) ELISA Kit
E-EL-H2124	Human TINF2 (TERF1 Interacting Nuclear Factor 2) ELISA Kit
E-EL-H2125	Human TRAb (Thyroid Stimulating Hormone Receptor Antidoby) ELISA Kit
E-EL-H2127	Human SRP9 (Signal Recognition Particle 9kDa) ELISA Kit
E-EL-H2128	Human NPTX1 (Neuronal Pentraxin 1) ELISA Kit
E-EL-H2129	Human KISS1 (Kisspeptin 1) ELISA Kit
E-EL-H2130	Human STAT1 (Signal Transducer and Activator of Transcription 1) ELISA Kit
E-EL-H2131	Human STAT2 (Signal Transducer and Activator of Transcription 2) ELISA Kit
E-EL-H2132	Human STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit
E-EL-H2133	Human STAT4 (Signal Transducer and Activator of Transcription 4) ELISA Kit
E-EL-H2134	Human STAT5A (Signal Transducer and Activator of Transcription 5A) ELISA Kit

Cat.NO	Product Name
E-EL-H2135	Human STAT6 (Signal Transducer and Activator of Transcription 6) ELISA Kit
E-EL-H2136	Human STAT6B (Signal Transducer and Activator of Transcription 6B) ELISA Kit
E-EL-H2137	Human SOX1 (Sex Determining Region Y Box Protein 1) ELISA Kit
E-EL-H2138	Human SOX18 (Sex Determining Region Y Box Protein 18) ELISA Kit
E-EL-H2139	Human SOX2 (Sex Determining Region Y Box Protein 2) ELISA Kit
E-EL-H2140	Human SOX3 (Sex Determining Region Y Box Protein 3) ELISA Kit
E-EL-H2141	Human SHBG (Sex Hormone Binding Globulin) ELISA Kit
E-EL-H2142	Human TS/TYMS (Thymidylate Synthetase) ELISA Kit
E-EL-H2143	Human TLA (Thymus Leukemia Antigen) ELISA Kit
E-EL-H2144	Human TECK/CCL25 (Thymus Expressed Chemokine) ELISA Kit
E-EL-H2145	Human TI-Ag (Thymus Independent Antigen) ELISA Kit
E-EL-H2146	Human TP (Thymidine Phosphorylase) ELISA Kit
E-EL-H2148	Human TMSβ10 (Thymosin Beta 10) ELISA Kit
E-EL-H2149	Human TMSB4 (Thymosin Beta 4) ELISA Kit
E-EL-H2151	Human TP5 (Thymopentin) ELISA Kit
E-EL-H2152	Human TD-Ag (Thymus Dependent Antigen) ELISA Kit
E-EL-H2153	Human OLFM4 (Olfactomedin 4) ELISA Kit
E-EL-H2154	Human NPTX2 (Neuronal Pentraxin 2) ELISA Kit
E-EL-H2155	Human VIP (Vasoactive Intestinal Peptide) ELISA Kit
E-EL-H2156	Human VPI (Vasoactive Peptidase Inhibitors) ELISA Kit

Cat.NO	Product Name
E-EL-H2157	Human VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit
E-EL-H2158	Human VEGF165 (Vascular Endothelial Growth Factor165) ELISA Kit
E-EL-H2159	Human GST μ 2/GST μ 2 (Glutathione S Transferase Mu 2, Muscle) ELISA Kit
E-EL-H2160	Human SUMO2 (Small Ubiquitin Related Modifier 2) ELISA Kit
E-EL-H2161	Human IgE (Immunoglobulin E) ELISA Kit
E-EL-H2162	Human IDO (Indoleamine-2,3-Dioxygenase) ELISA Kit
E-EL-H2163	Human OGN (Osteoglycin) ELISA Kit
E-EL-H2164	Human VEGF-B (Vascular Endothelial Cell Growth Factor B) ELISA Kit
E-EL-H2165	Human VEGFR3/Flt4 (Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit
E-EL-H2166	Human sVCAM-1/CD106 (soluble Vascular Cell Adhesion Molecule 1) ELISA Kit
E-EL-H2167	Human ES-Ab (Anti-Endostatin Antibody) ELISA Kit
E-EL-H2168	Human VWF (Von Willebrand Factor) ELISA Kit
E-EL-H2169	Human HbC (Hemoglobin C) ELISA Kit
E-EL-H2170	Human HbH (Hemoglobin H) ELISA Kit
E-EL-H2171	Human HB β (Hemoglobin Beta) ELISA Kit
E-EL-H2172	Human HO1 (Heme Oxygenase 1) ELISA Kit
E-EL-H2173	Human HO2 (Heme Oxygenase 2) ELISA Kit
E-EL-H2175	Human PC (Plasma Anticoagulant Protein C) ELISA Kit
E-EL-H2176	Human PS (Plasma Anticoagulant Protein S) ELISA Kit
E-EL-H2177	Human SEPP1(Selenoprotein P, Plasma)ELISA Kit

Cat.NO	Product Name
E-EL-H2179	Human HA (Hemagglutinin) ELISA Kit
E-EL-H2180	Human SGK2 (Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit
E-EL-H2181	Human ALB (Albumin) ELISA Kit
E-EL-H2182	Human KRAS (V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog) ELISA Kit
E-EL-H2183	Human SAA (Serum Amyloid A) ELISA Kit
E-EL-H2184	Human SEMA4C (Semaphorin 4C) ELISA Kit
E-EL-H2185	Human CR1 (Complement Receptor type 1) ELISA Kit
E-EL-H2186	Human MFG-E8 (Milk Fat Globule EGF Factor 8) ELISA Kit
E-EL-H2187	Human SERT (Serotonin Transporter) ELISA Kit
E-EL-H2188	Human BMG/β2-MG (Beta-2-Microglobulin) ELISA Kit
E-EL-H2189	Human TPP (Thrombus Precursor Protein) ELISA Kit
E-EL-H2191	Human TXB2 (Thromboxane B2) ELISA Kit
E-EL-H2192	Human STC2 (Stanniocalcin 2) ELISA Kit
E-EL-H2193	Human FG (Fibrinogen) ELISA Kit
E-EL-H2194	Human FDP (Fibrinogen Degradation Product) ELISA Kit
E-EL-H2195	Human FPA (Fibrinopeptide A) ELISA Kit
E-EL-H2196	Human FPB(Fibrinopeptide B)ELISA Kit
E-EL-H2197	Human SSTR2 (Somatostatin Receptor 2) ELISA Kit
E-EL-H2198	Human PHLDA2 (Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit
E-EL-H2199	Human PAF (Platelet Activating Factor) ELISA Kit

Cat.NO	Product Name
E-EL-H2200	Human PAFAHIB3 (Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit
E-EL-H2201	Human HAase (Hyaluronidase) ELISA Kit
E-EL-H2202	Human GP II b/ III a (Platelet Membrane Glycoprotein II b/ III a complex) ELISA Kit
E-EL-H2203	Human ITGB3 (Integrin Beta 3) ELISA Kit
E-EL-H2204	Human PKP1 (Plakophilin 1) ELISA Kit
E-EL-H2205	Human PKP2 (Plakophilin 2) ELISA Kit
E-EL-H2206	Human ANGPTL8 (Angiopoietin Like Protein 8) ELISA Kit
E-EL-H2207	Human PAC3 (Platelet Associated Complement 3) ELISA Kit
E-EL-H2208	Human PA-IgM (Platelet Associated Antibody) ELISA Kit
E-EL-H2209	Human PAIg (Platelet Associated Immunoglobulins) ELISA Kit
E-EL-H2210	Human PFKP (Phosphofructokinase, Platelet) ELISA Kit
E-EL-H2211	Human PDGF (Platelet-Derived Growth Factor) ELISA Kit
E-EL-H2212	Human PDGFC (Platelet Derived Growth Factor C) ELISA Kit
E-EL-H2213	Human PDGFD (Platelet Derived Growth Factor D) ELISA Kit
E-EL-H2214	Human PDGFR α (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit
E-EL-H2215	Human PF3 (Platelet Factor 3) ELISA Kit
E-EL-H2216	Human PF4 (Platelet Factor 4) ELISA Kit
E-EL-H2217	Human PF4V1 (Platelet Factor 4 Variant 1) ELISA Kit
E-EL-H2218	Human GYPE (Glycophorin E) ELISA Kit
E-EL-H2219	Human SSTR4 (Somatostatin Receptor 4) ELISA Kit

Cat.NO	Product Name
E-EL-H2221	Human NADPH (Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
E-EL-H2222	Human NOX1 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
E-EL-H2223	Human NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
E-EL-H2224	Human NOX5 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit
E-EL-H2225	Human EloA (Elongin A) ELISA Kit
E-EL-H2227	Human OLAB (Anti-Oxidized Low Density Lipoprotein Antibody) ELISA Kit
E-EL-H2228	Human SSTR3 (Somatostatin Receptor 3) ELISA Kit
E-EL-H2231	Human TRY (Trypsin) ELISA Kit
E-EL-H2232	Human Try2 (Trypsinogen II) ELISA Kit
E-EL-H2233	Human TAP (Trypsinogen Activation Peptide) ELISA Kit
E-EL-H2234	Human PI (Proinsulin) ELISA Kit
E-EL-H2235	Human AMY2 (Amylase Alpha 2, Pancreatic) ELISA Kit
E-EL-H2236	Human PP (Pancreatic Polypeptide) ELISA Kit
E-EL-H2237	Human GC (Glucagon) ELISA Kit
E-EL-H2238	Human GLP-2 (Glucagon Like Peptide 2) ELISA Kit
E-EL-H2239	Human JAK2 (Janus Kinase 2) ELISA Kit
E-EL-H2240	Human CA242 (Pancreatic Carcinoma Markers-CA242) ELISA Kit
E-EL-H2241	Human PTF1 α (Pancreas Specific Transcription Factor 1 α) ELISA Kit
E-EL-H2242	Human REG1 α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit
E-EL-H2243	Human REG4 (Regenerating Islet Derived Protein 4) ELISA Kit

Cat.NO	Product Name
E-EL-H2244	Human SULT1A1 (Sulfotransferase 1A1) ELISA Kit
E-EL-H2245	Human PL (Pancreatic Lipase) ELISA Kit
E-EL-H2246	Human PTPN11 (Protein Tyrosine Phosphatase, Non Receptor Type 11) ELISA Kit
E-EL-H2247	Human OSM (Oncostatin M) ELISA Kit
E-EL-H2248	Human DCN (Decorin) ELISA Kit
E-EL-H2249	Human PHB (Prohibitin) ELISA Kit
E-EL-H2250	Human CRY1 (Cryptochrome 1) ELISA Kit
E-EL-H2251	Human EL (EuglobulinLysis) ELISA Kit
E-EL-H2252	Human SEMA4D (Semaphorin 4D) ELISA Kit
E-EL-H2253	Human SSTR5 (Somatostatin Receptor 5) ELISA Kit
E-EL-H2254	Human FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit
E-EL-H2255	Human KNG1 (Kininogen 1) ELISA Kit
E-EL-H2256	Human S100A2 (S100 Calcium Binding Protein A2) ELISA Kit
E-EL-H2257	Human F-TESTO (Free Testosterone) ELISA Kit
E-EL-H2258	Human fT4 (Free Thyroxine) ELISA Kit
E-EL-H2259	Human VAP-1 (Vascular Adhesion Protein 1) ELISA Kit
E-EL-H2260	Human fT3 (Free Triiodothyronine) ELISA Kit
E-EL-H2261	Human f-Hb (Free Haemoglobin) ELISA Kit
E-EL-H2262	Human PTPN6 (Protein Tyrosine Phosphatase, Non Receptor Type 6) ELISA Kit
E-EL-H2263	Human MGMT (O-6-Methylguanine DNA Methyltransferase) ELISA Kit

Cat.NO	Product Name
E-EL-H2264	Human PCDH1 (Protocadherin 1) ELISA Kit
E-EL-H2265	Human PCDHβ2 (Protocadherin Beta 2) ELISA Kit
E-EL-H2266	Human FBN1 (Fibrillin 1) ELISA Kit
E-EL-H2267	Human FBN2 (Fibrillin 2) ELISA Kit
E-EL-H2268	Human FBN3 (Fibrillin 3) ELISA Kit
E-EL-H2269	Human Pg (Progesterone) ELISA Kit
E-EL-H2270	Human PAQR7 (Progestin and AdipoQ Receptor Family Member VII) ELISA Kit
E-EL-H2271	Human PGR (Progesterone Receptor) ELISA Kit
E-EL-H2272	Human sCD44 (Soluble Cluster of Differentiation 44) ELISA Kit
E-EL-H2273	Human PS2 (Presenilin 2) ELISA Kit
E-EL-H2274	Human EGR1 (Early Growth Response Protein 1) ELISA Kit
E-EL-H2275	Human EGR2 (Early Growth Response Protein 2) ELISA Kit
E-EL-H2276	Human EGR3 (Early Growth Response Protein 3) ELISA Kit
E-EL-H2277	Human EGR4 (Early Growth Response Protein 4) ELISA Kit
E-EL-H2278	Human NFYC (Nuclear Transcription Factor Y Gamma) ELISA Kit
E-EL-H2279	Human MUC5AC (Mucin 5 Subtype AC) ELISA Kit
E-EL-H2280	Human MUC5B (Mucin 5 Subtype B) ELISA Kit
E-EL-H2281	Human MUC7 (Mucin 7) ELISA Kit
E-EL-H2282	Human MEC/CCL28 (Mucosae Associated Epithelia Chemokine) ELISA Kit
E-EL-H2283	Human eIF3a (Eukaryotic Translation Initiation Factor 3a) ELISA Kit

Cat.NO	Product Name
E-EL-H2284	Human PC4 (Positive Cofactor 4) ELISA Kit
E-EL-H2286	Human LpPLA2 (Phospholipase A2, Lipoprotein Associated) ELISA Kit
E-EL-H2287	Human LPL (Lipoprotein Lipase) ELISA Kit
E-EL-H2289	Human LBP (Lipopolysaccharide Binding Protein) ELISA Kit
E-EL-H2290	Human IL4I1 (Interleukin 4 Induced Protein 1) ELISA Kit
E-EL-H2291	Human GYS1 (Glycogen Synthase 1, Muscle) ELISA Kit
E-EL-H2292	Human GYS2 (Glycogen Synthase 2, Liver) ELISA Kit
E-EL-H2296	Human OSMR (Oncostatin M Receptor) ELISA Kit
E-EL-H2297	Human MK (Midkine) ELISA Kit
E-EL-H2298	Human PON1 (Paraoxonase 1) ELISA Kit
E-EL-H2299	Human HNP1-3 (Neutrophil Peptide 1-3) ELISA Kit
E-EL-H2300	Human NAP3 (Neutrophil Activating Protein 3) ELISA Kit
E-EL-H2301	Human NAP (Neutrophil Alkaline Phosphatase) ELISA Kit
E-EL-H2302	Human PFN1 (Profilin 1) ELISA Kit
E-EL-H2303	Human ADAM22 (A Disintegrin And Metalloprotease 22) ELISA Kit
E-EL-H2304	Human TP63 (Tumor Protein p63) ELISA Kit
E-EL-H2305	Human TACE (TNF α Converting Enzyme) ELISA Kit
E-EL-H2306	Human TNF- β (Tumor Necrosis Factor Beta) ELISA Kit
E-EL-H2307	Human TL1A (Tumor Necrosis Factor Related Ligand 1A) ELISA Kit
E-EL-H2308	Human CCT8 (T-complex protein 1 subunit theta) ELISA Kit

Cat.NO	Product Name
E-EL-H2309	Human CCT8L1 (T-complex 1 protein subunit theta-like 1) ELISA Kit
E-EL-H2310	Human CCT5 (T-complex 1 protein subunit epsilon) ELISA Kit
E-EL-H2311	Human CCT6A (T-complex 1 protein subunit zeta) ELISA Kit
E-EL-H2312	Human TWEAK (Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit
E-EL-H2313	Human TSGF (Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
E-EL-H2315	Human TSTA (Tumor Specific Transplantation Antigen) ELISA Kit
E-EL-H2316	Human TAA (Tumor Associated Antigen) ELISA Kit
E-EL-H2318	Human NEB (Noradrenaline Bitartas) ELISA Kit
E-EL-H2319	Human RAG-1 (Recombination Activating Gene 1) ELISA Kit
E-EL-H2320	Human RAG-2 (Recombination Activating Gene 2) ELISA Kit
E-EL-H2321	Human CDK1 (Cyclin Dependent Kinase 1) ELISA Kit
E-EL-H2322	Human CDK2 (Cyclin Dependent Kinase 2) ELISA Kit
E-EL-H2323	Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit
E-EL-H2324	Human CDK5 (Cyclin Dependent Kinase 5) ELISA Kit
E-EL-H2325	Human CDK7 (Cyclin Dependent Kinase 7) ELISA Kit
E-EL-H2326	Human CDK9 (Cyclin Dependent Kinase 9) ELISA Kit
E-EL-H2327	Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
E-EL-H2328	Human Ntn1(Netrin 1)ELISA Kit
E-EL-H2329	Human Ntn4 (Netrin 4) ELISA Kit
E-EL-H2330	Human MVP (Major Vault Protein) ELISA Kit

Cat.NO	Product Name
E-EL-H2331	Human MHCC/HLA-C (Major Histocompatibility Complex Class I C) ELISA Kit
E-EL-H2332	Human MHCE/HLA-E (Major Histocompatibility Complex Class I E) ELISA Kit
E-EL-H2333	Human AMACR(Alpha-methylacyl-CoA racemase) ELISA Kit
E-EL-H2334	Human CCNA1(Cyclin-A1) ELISA Kit
E-EL-H2335	Human RELN (Reelin) ELISA Kit
E-EL-H2336	Human TGM2 (Transglutaminase 2, Tissue) ELISA Kit
E-EL-H2337	Human TGM3 (Transglutaminase 3, Epidermal) ELISA Kit
E-EL-H2338	Human IL-27 (Interleukin 27) ELISA Kit
E-EL-H2339	Human TGF-β3 (Transforming Growth Factor β3) ELISA Kit
E-EL-H2340	Human 4-HNE (4-Hydroxynonenal)ELISA Kit
E-EL-H2341	Human EDN (EDNEosinophil-Derived Neurotoxin) ELISA Kit
E-EL-H2342	Human TTR (Transthyretin) ELISA Kit
E-EL-H2343	Human TCF20 (Transcription Factor 20) ELISA Kit
E-EL-H2344	Human GAPDH (Glyceraldehyde-3-Phosphate Dehydrogenase) ELISA Kit
E-EL-H2345	Human TFR/CD71 (Transferrin Receptor) ELISA Kit
E-EL-H2346	Human TFR2 (Transferrin Receptor 2) ELISA Kit
E-EL-H2347	Human F-TESTO (Free Testosterone) ELISA Kit
E-EL-H2348	Human TF (Transfer Factor) ELISA Kit
E-EL-H2349	Human TNPO1 (Transportin 1) ELISA Kit
E-EL-H2350	Human Pax (Paxillin) ELISA Kit

Cat.NO	Product Name
E-EL-H2355	Human FPN (Ferroportin) ELISA Kit
E-EL-H2356	Human HEPH (Hephaestin) ELISA Kit
E-EL-H2357	Human CALP (Calprotectin) ELISA Kit
E-EL-H2358	Human 11,12-DHETs(11,12-dihydroxyeicosatrienoic acids) ELISA
E-EL-H2359	Human CD8 (Cluster of Differentiation 8) ELISA Kit
E-EL-H2360	Human PCX (Podocalyxin) ELISA Kit
E-EL-H2361	Human TPA (Tissue Polypeptide Antigen) ELISA Kit
E-EL-H2362	Human TPS (Tissue Polypeptide Specific Antigen) ELISA Kit
E-EL-H2363	Human SCARB1 (Scavenger receptor class B member 1) ELISA Kit
E-EL-H2364	Human HS (heparan sulfate) ELISA Kit
E-EL-H2365	Human HPA (Heparanase) ELISA Kit
E-EL-H2366	Human TFPI2(Tissue Factor Pathway Inhibitor 2)ELISA Kit
E-EL-H2367	Human DKK3 (Dickkopf Related Protein 3) ELISA Kit
E-EL-H2368	Human iC3b (Complement iC3b) ELISA Kit
E-EL-H2369	Human SEPCR/PROCR (Soluble Endothelial Protein C Receptor) ELISA Kit
E-EL-H2370	Human proPSA (precursor of Prostate Specific Antigen) ELISA Kit
E-EL-H2371	Human UBA52 (Ubiquitin A 52 Residue Ribosomal Protein Fusion Product 1) ELISA Kit
E-EL-H2372	Human PON2 (Paraoxonase 2) ELISA Kit
E-EL-H2373	Human NES1 (Nesfatin 1) ELISA Kit
E-EL-H2374	Human PON3 (Paraoxonase 3) ELISA Kit

Cat.NO	Product Name
E-EL-H2375	Human DDT (D-Dopachrome Tautomerase) ELISA Kit
E-EL-H2376	Human SC5b9 (soluble terminal complement complex) ELISA Kit
E-EL-H2377	Human UCHL1 (Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
E-EL-H2378	Human UBR4 (Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit
E-EL-H2379	Human MX1 (Myxovirus Resistance 1) ELISA Kit
E-EL-H2380	Human SAA2 (Serum Amyloid A2) ELISA Kit
E-EL-H2381	Human UBE2A (Ubiquitin Conjugating Enzyme E2A) ELISA Kit
E-EL-H2382	Human SOD3 (Superoxide Dismutase 3, Extracellular) ELISA Kit
E-EL-H2383	Human BCAR1 (Breast Cancer Anti Estrogen Resistance 1) ELISA Kit
E-EL-H2384	Human PLA2G2D (Phospholipase A2, Group II D) ELISA Kit
E-EL-H2385	Human PLA2G2A (Phospholipase A2, Group II A) ELISA Kit
E-EL-H2386	Human UBAP2 (Ubiquitin Associated Protein 2) ELISA Kit
E-EL-H2387	Human MASP1 (Mannan Associated Serine Protease 1) ELISA Kit
E-EL-H2388	Human MASP2 (Mannan Associated Serine Protease 2) ELISA Kit
E-EL-H2389	Human USP8 (Ubiquitin Specific Peptidase 8) ELISA Kit
E-EL-H2390	Human DPP6 (Dipeptidyl Peptidase VI) ELISA Kit
E-EL-H2391	Human MRC2 (Mannose Receptor C Type 2) ELISA Kit
E-EL-H2392	Human MRC1 (Mannose Receptor C Type 1) ELISA Kit
E-EL-H2393	Human MRC1 (Mannose Receptor C Type 1 Like Protein 1) ELISA Kit
E-EL-H2394	Human GADL1 (Glutamate Decarboxylase Like Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-H2395	Human GAD1 (Glutamate Decarboxylase 1) ELISA Kit
E-EL-H2396	Human PEX2 (Peroxisomal Biogenesis Factor 2) ELISA Kit
E-EL-H2397	Human PPAR- δ (Peroxisome Proliferator Activated Receptor Delta) ELISA Kit
E-EL-H2398	Human NUP62 (Nucleoporin 62kDa) ELISA Kit
E-EL-H2399	Human PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit
E-EL-H2400	Human RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
E-EL-H2401	Human NUMA1 (Nuclear Mitotic Apparatus Protein 1) ELISA Kit
E-EL-H2402	Human IL-33 (Interleukin 33) ELISA Kit
E-EL-H2403	Human MCAM (Melanoma Cell Adhesion Molecule) ELISA Kit
E-EL-H2404	Human MIA1 (Melanoma Inhibitory Activity Protein 1) ELISA Kit
E-EL-H2405	Human MFI2 (Melanotransferrin) ELISA Kit
E-EL-H2406	Human SPT β N4 (Spectrin Beta, Non Erythrocytic 4) ELISA Kit
E-EL-H2407	Human I κ B α (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
E-EL-H2408	Human I κ B β (Inhibitory Subunit of NF Kappa B Beta) ELISA Kit
E-EL-H2409	Human HNRPA2B1 (Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit
E-EL-H2410	Human HNRPU (Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit
E-EL-H2411	Human RAR α (Retinoic Acid Receptor Alpha) ELISA Kit
E-EL-H2412	Human RNASEN (Ribonuclease III, Nuclear) ELISA Kit
E-EL-H2413	Human MCSP (Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit
E-EL-H2414	Human α MSH (Alpha-Melanocyte Stimulating Hormone) ELISA Kit

Cat.NO	Product Name
E-EL-H2415	Human EP (Erythrocyte Protoporphyrin) ELISA Kit
E-EL-H2416	Human THR α (Thyroid Hormone Receptor Alpha) ELISA Kit
E-EL-H2417	Human NFR κ B (Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit
E-EL-H2418	Human TNKS1 (Tankyrase 1) ELISA Kit
E-EL-H2419	Human TNKS2 (Tankyrase 2) ELISA Kit
E-EL-H2420	Human PARP4 (Poly ADP Ribose Polymerase 4) ELISA Kit
E-EL-H2421	Human Neopterin ELISA Kit
E-EL-H2422	Human ME20M (Melanoma Associated ME20) ELISA Kit
E-EL-H2423	Human PXR (Pregnane X Receptor) ELISA Kit
E-EL-H2424	Human ALOX15B (Arachidonate-15-Lipoxygenase, Type B) ELISA Kit
E-EL-H2425	Human RELB (V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
E-EL-H2426	Human IL-24 (Interleukin 24) ELISA Kit
E-EL-H2427	Human AIM2 (Absent In Melanoma 2) ELISA Kit
E-EL-H2428	Human TSPAN30/CD63 (Tetraspanin 30) ELISA Kit
E-EL-H2430	Human BANP (BTG3 Associated Nuclear Protein) ELISA Kit
E-EL-H2431	Human PRKG1 (Protein Kinase, cGMP Dependent Type I) ELISA Kit
E-EL-H2432	Human PDE4D (Phosphodiesterase 4D, cAMP Specific) ELISA Kit
E-EL-H2433	Human CKB (Creatine Kinase, Brain) ELISA Kit
E-EL-H2434	Human Anti-LKM-1 IgG (Anti-Liver-Kidney Microsomes-1 IgG) ELISA Kit
E-EL-H2435	Human LOXL4 (Lysyl Oxidase Like Protein 4) ELISA Kit

Cat.NO	Product Name
E-EL-H2436	Human TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
E-EL-H2437	Human ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-H2438	Human LL-37(Antibacterial Protein LL-37)ELISA Kit
E-EL-H2439	Human Tα1 (Thymosin-α1) ELISA Kit
E-EL-H2440	Human PS1 (Presenilin 1) ELISA Kit
E-EL-H2441	Human PACAP-38 (Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
E-EL-H2442	Human PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-H2443	Human IL-35 (Interleukin 35) ELISA Kit
E-EL-H2444	Human TK1 (Thymidine Kinase 1, Soluble) ELISA Kit
E-EL-H2445	Human 12-LO (Arachidonate 12-Lipoxygenase) ELISA Kit
E-EL-H2446	Human VIL1 (Villin 1) ELISA Kit
E-EL-H2447	Human PI3/Elafin (Peptidase Inhibitor 3, Skin Derived) ELISA Kit
E-EL-H2448	Human TTPα (Alpha Tocopherol Transfer Protein) ELISA Kit
E-EL-H2449	Human FNDC5 (Fibronectin Type III Domain Containing Protein 5) ELISA Kit
E-EL-H2450	Human IL-21 (Interleukin 21) ELISA Kit
E-EL-H2451	Human SSADH (Succinate-semialdehyde dehydrogenase, mitochondrial) ELISA Kit
E-EL-H2452	Human POSTN/OSF-2 (Periostin) ELISA Kit
E-EL-H2453	Human MMP-15 (Matrix Metalloproteinase 15) ELISA Kit
E-EL-H2454	Human MMP-14 (Matrix Metalloproteinase 14) ELISA Kit
E-EL-H2455	Human VLDL (Very Low Density Lipoprotein) ELISA Kit

Cat.NO	Product Name
E-EL-H2456	Human LCAT (Lecithin-Cholesterol Acyltransferase) ELISA Kit
E-EL-H2457	Human CYP2E1 (Cytochrome P450, family 2, subfamily E, polypeptide 1) ELISA Kit
E-EL-H2531	Human UCRP (Ubiquitin Cross Reactive Protein) ELISA Kit
E-EL-H2532	Human IFN- α (Interferon Alpha) ELISA Kit
E-EL-H2533	Human DAG1 (Dystrophin Associated Glycoprotein 1) ELISA Kit
E-EL-H2534	Human HABP1/C1QBP (Hyaluronan Binding Protein 1) ELISA Kit
E-EL-H2535	Human RTN4 (Reticulon 4) ELISA Kit
E-EL-H2568	Human Anti-CRP (Anti-C Reactive Protein Antibody) ELISA Kit
E-EL-H2569	Human VEGF (Vascular Endothelial Cell Growth Factor) ELISA Kit
E-EL-H2570	Human CLDN (Claudin) ELISA Kit
E-EL-H2571	Human IL-37 (Interleukin 37) ELISA Kit
E-EL-H2572	Human sAGER (Soluble Advanced Glycosylation End Product Specific Receptor) ELISA Kit
E-EL-H2578	Human CXCR2(CXC-Chemokine Receptor 2) ELISA Kit
E-EL-H2579	Human MEP1 α (Meprin A Alpha) ELISA Kit
E-EL-H2580	Human IL26 (Interleukin 26) ELISA Kit
E-EL-H2581	Human SGPP1 (Sphingosine-1-Phosphate Phosphatase 1) ELISA Kit
E-EL-H2582	Human SPHK1 (Sphingosine Kinase 1) ELISA Kit
E-EL-H2583	Human S1P (Sphingosine 1 phosphate) ELISA Kit
E-EL-H2584	Human suPAR(soluble urokinase-type plasminogen activator receptor) ELISA Kit
E-EL-H2585	Human ICAM-1 (intercellular adhesion molecule 1) ELISA Kit

Cat.NO	Product Name
E-EL-H4184	Human PAEP (Progestagen-Associated Endometrial Protein) ELISA Kit
E-EL-H5022	Human IL-11 (Interleukin 11) ELISA Kit
E-EL-H5105	Human sdLDL (Small Dense Low Density Lipoprotein) ELISA Kit
E-EL-H5134	Human hs-CRP (high-sensitivity C-reactive protein) ELISA Kit
E-EL-H5135	Human TSHR (Thyroid Stimulating Hormone Receptor) ELISA Kit
E-EL-H5200	Human LTF(Lactoferrin)ELISA Kit
E-EL-H5300	Human TERF2IP (Telomeric Repeat Binding Factor 2 Interacting Protein) ELISA Kit
E-EL-H5301	Human FNDC5(Fibronectin Type III Domain Containing Protein 5) ELISA Kit
E-EL-H5302	Human Anti-dsDNA (Anti-Double Stranded DNA Antibody) ELISA Kit
E-EL-H5303	Human DBT (Dihydrolipoamide branched chain transacylase E2) ELISA Kit
E-EL-H5307	Human AD (Adropin) ELISA Kit
E-EL-H5308	Human Salusin β (Salusin Beta) ELISA Kit
E-EL-H5309	Human ENG (Soluble Endoglin) ELISA Kit
E-EL-H5310	Human c-MET/HGFR (hepatocyte growth factor receptor) ELISA Kit
E-EL-H5311	Human RECK(Reversion-inducing-cysteine-rich protein with kazal motifs)ELISA Kit
E-EL-H5312	Human CD44v5 (cluster of differentiation 44 variant 5)ELISA Kit
E-EL-H5313	Human TPH1(Tryptophan Hydroxylase 1)ELISA Kit
E-EL-H5314	Human pMAPT /pTAU(phosphorylated microtubule-associated protein tau)ELISA Kit
E-EL-H5315	Human HSD17b10 (17-Beta-Hydroxysteroid Dehydrogenase Type 10) ELISA Kit
E-EL-H5316	Human HDAC1 (Histone Deacetylase 1) ELISA Kit

Cat.NO	Product Name
E-EL-H5317	Human MICA (major histocompatibility complex class I chain-related gene A) ELISA Kit
E-EL-H5401	Human KISS(Kisspeptin)ELISA Kit
E-EL-H5402	Human ENG(Endoglin)ELISA Kit
E-EL-H5405	Human SCUBE1 (Signal Peptide, CUB and EGF-like Domain-containing Protein 1)ELISA Kit
E-EL-H5407	Human GSSG(Oxidized Glutathione) ELISA Kit
E-EL-H5408	Human CAT(Catalase) ELISA Kit
E-EL-H5409	Human TBA(Total bile acid)ELISA Kit
E-EL-H5410	Human GPX1 (Glutathione Peroxidase 1) ELISA Kit
E-EL-H5411	Human Placenta growth factor(PGF)ELISA Kit
E-EL-H5412	Human sFLT-1(Receptor for Vascular Endothelial cell Growth Factor VEGF)ELISA Kit
E-EL-H5413	Human PSMA(Prostate specific membrane antigen)ELISA Kit
E-EL-H5414	Human SORT1(Sortilin 1) ELISA Kit
E-EL-H5415	Human NID2 (Nidogen 2) ELISA Kit
E-EL-H5416	Human ICA(islet cell antibody)ELISA Kit
E-EL-H5417	Human Peth(Phosphatidylethanol)ELISA Kit
E-EL-H5418	Human MAC (Membrane Attack Complex) ELISA Kit
E-EL-H5419	Human 25(OH)D(25-Hydroxy Vitamin D) ELISA Kit
E-EL-H5420	Human HPF4(Heparin/Platelet Factor 4 Antibodies)ELISA Kit
E-EL-H5421	Human GRIA1(Glutamate receptor 1)ELISA Kit
E-EL-H5422	Human IMA(Ischemia Modified Albumin)ELISA Kit

Cat.NO	Product Name
E-EL-H5423	Human MeCP2(Methyl CpG Binding Protein 2) Elisa Kit
E-EL-H5424	Human MBD1(Methyl CpG Binding Domain Protein 1) Elisa Kit
E-EL-H5425	Human MBD2(Methyl CpG Binding Domain Protein 2)Elisa Kit
E-EL-H5426	Human DNMT1(DNA Methyltransferase 1) ELISA Kit
E-EL-H5427	Human DNMT2(DNA Methyltransferase 2)ELISA Kit
E-EL-H5428	Human DNMT3A(DNA Methyltransferase 3A)ELISA Kit
E-EL-H5429	Human DNMT3B(DNA Methyltransferase 3B)ELISA Kit
E-EL-H5430	Human ROMO1(Reactive oxygen species modulator 1)ELISA Kit
E-EL-H5431	Human CTSS(Cathepsin S)ELISA Kit
E-EL-H5432	Human Ki67P(Ki-67 Protein)ELISA Kit
E-EL-H5433	Human HE4(epididymal protein 4)ELISA Kit
E-EL-H5434	Human FGF18 (Fibroblast Growth Factor 18)ELISA Kit
E-EL-H5435	Human TACSTD2(Tumor-associated calcium signal transducer 2) ELISA Kit
E-EL-H5436	Human ROS(Reactive OXyen Species) ELISA Kit
E-EL-H5437	Human TG (Triglyceride) ELISA Kit
E-EL-H5438	Human PDCD1LG2 (Programmed Cell Death Protein 1 Ligand 2) ELISA Kit
E-EL-H5439	Human INS (Insulin) ELISA Kit
E-EL-H5440	Human PEBP4(Phosphatidylethanolamine Binding Protein 4) ELISA Kit
E-EL-H5441	Human MR-ProADM(Mid-regional pro-adrenomedullin)ELISA Kit
E-EL-H5442	Human PIVKA-II(Protein Induced by Vitamin K Absence or Antagonist- II)ELISA Kit

Cat.NO	Product Name
E-EL-H5443	Human MGO(Methylglyoxal) ELISA Kit
E-EL-H5444	Human IAPPI(slet Amyloid Polypeptide)ELISA Kit
E-EL-H5445	Human PRGN(progranulin)ELISA Kit
E-EL-H5446	Human C4(Complement 4)ELISA Kit
E-EL-H5447	Human SFRP4(secreted frizzled-related protein 4)ELISA Kit
E-EL-H5448	Human FCN3(Ficolin 3)ELISA Kit
E-EL-H5449	Human PDH (Pyruvate Dehydrogenase) ELISA Kit
E-EL-H5450	Human HLA-G(leukocyte antigen G)ELISA Kit
E-EL-H5451	Human KL(Klotho)ELISA Kit
E-EL-H5453	Human ARL1(Aldose Reductase Like Protein 1)ELISA Kit
E-EL-H5454	Human PTP1B(protein tyrosine phosphatase 1B)ELISA Kit
E-EL-H5455	Human BIRC7(Baculoviral IAP repeat-containing protein 7)ELISA Kit
E-EL-H5456	Human CIC(circulating immune complex)ELISA Kit
E-EL-H5457	Human C3d(Complement Fragment 3d)ELISA Kit
E-EL-H5458	Human Hsp60 IgA(Heat Shock Protein 60 IgA) ELISA Kit
E-EL-H5459	Human Hsp60 IgG(Heat Shock Protein 60 IgG) ELISA Kit
E-EL-H5460	Human Hsp60 IgM(Heat Shock Protein 60 IgM) ELISA Kit
E-EL-H5461	Human FGF8 (Fibroblast Growth Factor 8)ELISA Kit
E-EL-H5462	Human FGF10(Fibroblast Growth Factor 10)ELISA Kit
E-EL-H5463	Human ABCC5 (ATP Binding Cassette Transporter C5)ELISA Kit

Cat.NO	Product Name
E-EL-H5464	Human N-ERC/Mesothelin ELISA Kit
E-EL-H5465	Human SAM(S-Adenosyl Methionine)ELISA Kit
E-EL-H5466	Human LAG3(Lymphocyte Activation Gene 3)ELISA Kit
E-EL-H5467	Human HGS(Hepatocyte growth factor-regulated tyrosine kinase substrate)ELISA Kit
E-EL-H5468	Human LEFTY1(Left/Right Determination Factor 1)ELISA Kit
E-EL-H5469	Human IL31(Interleukin 31)ELISA Kit
E-EL-H5470	Human SCGB2A2(Mammaglobin-A)ELISA Kit
E-EL-H5471	Human DNP-Ab (Deoxyribonucleic Protein Antibody) ELISA Kit
E-EL-H5472	Human HBsAb (Hepatitis B surface antibody) ELISA Kit
E-EL-H5473	Human CD44v9 (cluster of differentiation 44 variant 9)ELISA Kit
E-EL-H5474	Human BLVRA(Biliverdin Reductase A)ELISA Kit
E-EL-H5475	Human ICTP(C-telopeptide of type I collagen) ELISA Kit
E-EL-H5476	Human NE (Neutrophil Elastase) ELISA Kit
E-EL-H5478	Human BPn(biopyrrin)ELISA Kit
E-EL-H5479	Human SLC16A8(Solute Carrier Family 16, Member 8)ELISA Kit
E-EL-H5480	Human IER3(Immediate Early Response 3)ELISA Kit
E-EL-H5481	Human ADAMTS9(A Disintegrin And Metalloproteinase With Thrombospondin 9)ELISA Kit
E-EL-H5482	Human B3GALT1 (Beta-1,3-glucosyltransferase)ELISA Kit
E-EL-H5483	Human HES1(Transcription factor HES-1)ELISA Kit
E-EL-H5484	Human CSNK2A1 (Casein Kinase 2, alpha 1 Polypeptide) ELISA Kit

Cat.NO	Product Name
E-EL-H5485	Human CSNK2B (Casein Kinase 2, beta) ELISA Kit
E-EL-H5486	Human TE(telomerase)ELISA Kit
E-EL-H5487	Human 2, 3-DPG(2,3-Disphosphoglycerate)ELISA Kit
E-EL-H5488	Human SMRP(Soluble Mesothelin Related Peptide)ELISA Kit
E-EL-H5489	Human CCR4(Chemokine C-C-Motif Receptor)ELISA Kit
E-EL-H5490	Human CCR4(Chemokine C-X-C-Motif Receptor)ELISA Kit
E-EL-H5491	Human C-ERC/Mesothelin ELISA Kit
E-EL-H5493	Human ucOC(Undercarboxylated Osteocalcin)ELISA Kit
E-EL-H5494	Human Pro-ADM(Proadrenomedullin)ELISA Kit
E-EL-H5495	Human CECR1 (Cat Eye Syndrome Chromosome Region, Candidate 1) ELISA Kit
E-EL-H5496	Human PROM1(Prominin-1)ELISA Kit
E-EL-H5497	Human Hepc25 (Hepcidin 25) ELISA Kit
E-EL-H5498	Human HOXA10(Homeobox protein Hox-A10)ELISA Kit
E-EL-H5499	Human HOXB4(Homeobox protein Hox-B4)ELISA Kit
E-EL-H5500	Human PL(Placental Lactogen)ELISA Kit
E-EL-H5501	Human TF(Thymic Factor)ELISA Kit
E-EL-H5502	Human MC5R(Melanocortin 5 Receptor)ELISA Kit
E-EL-H5503	Human QSOX1(Sulfhydryl oxidase 1)ELISA Kit
E-EL-H5504	Human ATRN (Attractin) ELISA Kit
E-EL-H5505	Human SCG3(Secretogranin III)ELISA Kit

Cat.NO	Product Name
E-EL-H5507	Human IL-28A (Interleukin 28A) ELISA Kit
E-EL-H5508	Human IL-28B (Interleukin 28B) ELISA Kit
E-EL-H5509	Human IL-29 (Interleukin 29) ELISA Kit
E-EL-H5510	Human C4d (Complement Fragment 4d) ELISA Kit
E-EL-H5511	Human VNN1(Vanin 1)ELISA Kit
E-EL-H5512	Human ANXA1(Annexin A1) ELISA Kit
E-EL-H5513	Human CF6(coupling factor 6)ELISA Kit
E-EL-H5514	Human BCL2A1(Bcl-2-related protein A1)ELISA Kit
E-EL-H5515	Human TIMD4(T-cell immunoglobulin and mucin domain-containing protein 4)ELISA Kit
E-EL-H5516	Human VSNL1(Visinin Like Protein 1)ELISA Kit
E-EL-H5517	Human IAP(Integrin Associated Protein)ELISA Kit
E-EL-H5518	Human Ang1-7(Angiotensin 1-7)ELISA Kit
E-EL-H5519	Human MC Ab(Melanocyte Antibody)ELISA Kit
E-EL-H5520	Human ANXA3 (Annexin A3) ELISA Kit
E-EL-H5521	Human AKAP12(A Kinase Anchor Protein 12) ELISA Kit
E-EL-H5522	Human C1QTNF2(C1q and Tumor Necrosis Factor Related Protein 2) ELISA Kit
E-EL-H5523	Human C1QTNF3(C1q and Tumor Necrosis Factor Related Protein 3) ELISA Kit
E-EL-H5524	Human C1QTNF9(C1q and Tumor Necrosis Factor Related Protein 9) ELISA Kit
E-EL-H5525	Human FAP α (Fibroblast Activation Protein Alpha)ELISA Kit
E-EL-H5526	Human GDN(Glia Derived Nexin)ELISA Kit

Cat.NO	Product Name
E-EL-H5527	Human PLOD3(Procollagen Lysine-2-Oxoglutarate-5-Dioxygenase 3) ELISA Kit
E-EL-H5528	Human CLIP(Corticotropin Like Intermediate Lobe Peptide)ELISA Kit
E-EL-H5529	Human γ MSH (gama-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-H5530	Human SESN1(Sestrin 1)ELISA Kit
E-EL-H5531	Human LAMA4(Laminin Alpha 4)ELISA Kit
E-EL-H5532	Human CSNK2A2 (Casein Kinase 2, alpha 2 Polypeptide) ELISA Kit
E-EL-H5533	Human STUB1(E3 ubiquitin-protein ligase CHIP)ELISA Kit
E-EL-H5534	Human pNF-H (Phosphorylated Neurofilament H) ELISA Kit
E-EL-H5535	Human DHFR(Dihydrofolate Reductase)ELISA Kit
E-EL-H5536	Human sSELP (Soluble P-Selectin) ELISA Kit
E-EL-M0001	Mouse ACV-A (Activin A) ELISA Kit
E-EL-M0002	Mouse ADP/Acrp30 (Adiponectin) ELISA Kit
E-EL-M0003	Mouse ANG3 (Angiopoietin 3) ELISA Kit
E-EL-M0004	Mouse BAFF/CD257 (B-cell Activating Factor) ELISA Kit
E-EL-M0005	Mouse LOX (Lysyl Oxidase) ELISA Kit
E-EL-M0006	Mouse MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit
E-EL-M0007	Mouse MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
E-EL-M0008	Mouse MIP-1 β (Macrophage Inflammatory Protein 1 Beta)ELISA Kit
E-EL-M0009	Mouse RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
E-EL-M0010	Mouse TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit

Cat.NO	Product Name
E-EL-M0011	Mouse CCL12/MCP-5 (Monocyte Chemotactic Protein 5) ELISA Kit
E-EL-M0012	Mouse TARC (Thymus Activation Regulated Chemokine) ELISA Kit
E-EL-M0013	Mouse MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
E-EL-M0014	Mouse MDC/CCL22 (Macrophage-Derived Chemokine) ELISA Kit
E-EL-M0015	Mouse CHEM (Chemerin) ELISA Kit
E-EL-M0016	Mouse CHI3L1 (Chitinase 3-like 1) ELISA Kit
E-EL-M0017	Mouse IL-1R1 (Interleukin 1 Receptor Type I) ELISA Kit
E-EL-M0018	Mouse GRO α /CXCL1 (Growth Regulated Oncogene Alpha) ELISA Kit
E-EL-M0019	Mouse GRO β (Growth Regulated Oncogene Beta) ELISA Kit
E-EL-M0020	Mouse CXCL9/MIG(Monokine induced by interferon-gamma)ELISA Kit
E-EL-M0021	Mouse IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
E-EL-M0023	Mouse Cyt C (Cytochrome C) ELISA Kit
E-EL-M0024	Mouse DKK1 (Dickkopf Related Protein 1)ELISA Kit
E-EL-M0025	Mouse EGF (Epidermal growth factor) ELISA Kit
E-EL-M0026	Mouse ESM1 (Endothelial Cell Specific Molecule 1) ELISA Kit
E-EL-M0027	Mouse EPO (Erythropoietin) ELISA Kit
E-EL-M0028	Mouse FASL (Factor Related Apoptosis Ligand) ELISA Kit
E-EL-M0029	Mouse FGF21(Fibroblast Growth Factor 21) ELISA Kit
E-EL-M0030	Mouse Flt3L (FMS Like Tyrosine Kinase 3 Ligand) ELISA Kit
E-EL-M0031	Mouse G-CSF (Granulocyte Colony Stimulating Factor) ELISA Kit

Cat.NO	Product Name
E-EL-M0032	Mouse GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
E-EL-M0033	Mouse IFN- β (Interferon Beta)ELISA Kit
E-EL-M0034	Mouse IGF-1 (Insulin-like growth factor 1) ELISA Kit
E-EL-M0035	Mouse IGFBP-3 (Insulin-like Growth Factor Binding Protein 3) ELISA Kit
E-EL-M0036	Mouse IL-1 α (Interleukin 1 Alpha)ELISA Kit
E-EL-M0037	Mouse IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-M0038	Mouse IL-1ra (Interleukin 1 Receptor Antagonist) ELISA Kit
E-EL-M0039	Mouse LEP (Leptin) ELISA Kit
E-EL-M0040	Mouse LIF (Leukemia Inhibitory Factor) ELISA Kit
E-EL-M0042	Mouse IL-2 (Interleukin 2) ELISA Kit
E-EL-M0043	Mouse IL-4 (Interleukin 4) ELISA Kit
E-EL-M0044	Mouse IL-6 (Interleukin 6) ELISA Kit
E-EL-M0045	Mouse IL-24 (Interleukin 24) ELISA Kit
E-EL-M0046	Mouse IL-10 (Interleukin 10) ELISA Kit
E-EL-M0047	Mouse IL-17 (Interleukin 17) ELISA Kit
E-EL-M0048	Mouse IFN- γ (Interferon Gamma) ELISA Kit
E-EL-M0049	Mouse TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-M0050	Mouse VEGF (Vascular Endothelial cell Growth Factor) ELISA Kit
E-EL-M0051	Mouse TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit
E-EL-M0053	Mouse CRP (C-Reactive Protein) ELISA Kit

Cat.NO	Product Name
E-EL-M0056	Mouse I-TAC (Interferon Inducible T-Cell Alpha Chemoattractant) ELISA Kit
E-EL-M0057	Mouse LH (Luteinizing Hormone) ELISA Kit
E-EL-M0060	Mouse GH (Growth Hormone) ELISA Kit
E-EL-M0061	Mouse REN (Renin) ELISA Kit
E-EL-M0062	Mouse CIC (Circulating Immune Complexes) ELISA Kit
E-EL-M0063	Mouse CASP8 (Caspase 8) ELISA Kit
E-EL-M0066	Mouse OxLDL (Oxidized Low-Density Lipoprotein) ELISA Kit
E-EL-M0067	Mouse A β 40 (Amyloid Beta 40) ELISA Kit
E-EL-M0068	Mouse A β 42 (Amyloid Beta 42) ELISA Kit
E-EL-M0072	Mouse TPS (Tryptase) ELISA Kit
E-EL-M0074	Mouse MMP-10 (Matrix Metalloproteinase 10) ELISA Kit
E-EL-M0075	Mouse MMP-12 (Matrix Metalloproteinase 12) ELISA Kit
E-EL-M0076	Mouse MMP-13 (Matrix Metalloproteinase 13) ELISA Kit
E-EL-M0077	Mouse sICAM-1/CD54 (Soluble Intercellular Adhesion Molecule 1) ELISA Kit
E-EL-M0078	Mouse ANG (Angiogenin) ELISA Kit
E-EL-M0079	Mouse ACTH (Adrenocorticotrophic Hormone) ELISA Kit
E-EL-M0081	Mouse OPG (Osteoprotegerin) ELISA Kit
E-EL-M0083	Mouse PRL (Prolactin) ELISA Kit
E-EL-M0084	Mouse MYO (Myoglobin) ELISA Kit
E-EL-M0085	Mouse Tn-T (Troponin T) ELISA Kit

Cat.NO	Product Name
E-EL-M0086	Mouse Tn-I (Troponin I) ELISA Kit
E-EL-M0087	Mouse PAPPA(Pregnancy Associated Plasma Protein A)ELISA Kit
E-EL-M0088	Mouse ALP (Alkaline Phosphatase) ELISA Kit
E-EL-M0089	Mouse NSE (Neuron-Specific Enolase) ELISA Kit
E-EL-M0090	Mouse GLP-1 (Glucagon Like Peptide 1)ELISA Kit
E-EL-M0091	Mouse ACE (Angiotensin I Converting Enzyme) ELISA Kit
E-EL-M0092	Mouse ANGPTL2 (Angiopoietin Like Protein 2) ELISA Kit
E-EL-M0093	Mouse ANGPTL4 (Angiopoietin Like Protein 4) ELISA Kit
E-EL-M0094	Mouse ANG-R-Tie1 (Angiopoietin Receptor Tie1) ELISA Kit
E-EL-M0095	Mouse ANG-R-Tie2 (Angiopoietin Receptor Tie2) ELISA Kit
E-EL-M0096	Mouse AGT (Angiotensinogen) ELISA Kit
E-EL-M0097	Mouse ANG1 (Angiopoietin 1) ELISA Kit
E-EL-M0098	Mouse ANG2 (Angiopoietin 2) ELISA Kit
E-EL-M0099	Mouse ANXA5 (Annexin A5) ELISA Kit
E-EL-M0100	Mouse ANXA2(Annexin A2) ELISA Kit
E-EL-M0101	Mouse α -Fodrin IgG/IgA (anti-alpha-Fodrin IgG/IgA) ELISA Kit
E-EL-M0102	Mouse ACA-IgA (anti-Cardiolipin Antibody IgA) ELISA Kit
E-EL-M0103	Mouse ACA-IgG (anti-Cardiolipin Antibody IgG) ELISA Kit
E-EL-M0104	Mouse ACA-IgM (anti-Cardiolipin Antibody IgM) ELISA Kit
E-EL-M0105	Mouse ACA (anti-Centrol and Centrosome Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-M0106	Mouse ADH (Antidiuretic Hormone) ELISA Kit
E-EL-M0107	Mouse NRG-1 (Neuregulin 1) ELISA Kit
E-EL-M0108	Mouse EMA IgA (anti-Endomysial Antibody IgA) ELISA Kit
E-EL-M0109	Mouse AECA (anti-Endothelial Cell Antibody) ELISA Kit
E-EL-M0110	Mouse GP210 (anti-Glucoprotein 210) ELISA Kit
E-EL-M0111	Mouse IgERA(anti-IgE receptor antibody)ELISA Kit
E-EL-M0113	Mouse AMH (Anti-Mullerian Hormone) ELISA Kit
E-EL-M0116	Mouse OMgp-Ab (anti-Oligodendrocyte-Myelin Glycoprotein Antibody) ELISA Kit
E-EL-M0117	Mouse CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
E-EL-M0118	Mouse CDKN1A (Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit
E-EL-M0120	Mouse HDAC2 (Histone Deacetylase 2) ELISA Kit
E-EL-M0121	Mouse NRG-4 (Neuregulin 4) ELISA Kit
E-EL-M0122	Mouse AT- III (anti-Thrombin III) ELISA Kit
E-EL-M0123	Mouse ATR (anti-Thrombin Receptor) ELISA Kit
E-EL-M0124	Mouse TGAB (anti-Thyroid-Globulin Antibody) ELISA Kit
E-EL-M0125	Mouse TPO-Ab (anti-Thyroid-Peroxidase Antibody) ELISA Kit
E-EL-M0126	Mouse AT (anti-Trypsin) ELISA Kit
E-EL-M0129	Mouse ACF (Apobec 1 Complementation Factor) ELISA Kit
E-EL-M0130	Mouse APO-A1(Apolipoprotein A1)ELISA Kit
E-EL-M0132	Mouse ApoB (Apolipoprotein B) ELISA Kit

Cat.NO	Product Name
E-EL-M0133	Mouse APO-C1 (Apolipoprotein C1) ELISA Kit
E-EL-M0134	Mouse APO-C4 (Apolipoprotein C4) ELISA Kit
E-EL-M0135	Mouse APO-E (Apolipoprotein E) ELISA Kit
E-EL-M0136	Mouse APO-F (Apolipoprotein F) ELISA Kit
E-EL-M0137	Mouse APO-H (Apolipoprotein H) ELISA Kit
E-EL-M0138	Mouse ApoB100 (Apolipoprotein B100) ELISA Kit
E-EL-M0139	Mouse AATF (Apoptosis Antagonizing Transcription Factor) ELISA Kit
E-EL-M0140	Mouse AATK (Apoptosis Associated Tyrosine Kinase) ELISA Kit
E-EL-M0141	Mouse AIF (Apoptosis Inducing Factor) ELISA Kit
E-EL-M0142	Mouse ASK1 (Apoptosis Signal Regulating Kinase 1) ELISA Kit
E-EL-M0143	Mouse APAF1 (Apoptotic Peptidase Activating Factor 1) ELISA Kit
E-EL-M0144	Mouse AQP-0 (Aquaporin 0) ELISA Kit
E-EL-M0145	Mouse SLC (Secondary Lymphoid Tissue Chemokine) ELISA Kit
E-EL-M0146	Mouse AQP-1 (Aquaporin 1) ELISA Kit
E-EL-M0147	Mouse GRO γ (Growth Regulated Oncogene Gamma) ELISA Kit
E-EL-M0148	Mouse AQP-2 (Aquaporin 2) ELISA Kit
E-EL-M0149	Mouse AGER (Advanced Glycosylation End Product Specific Receptor) ELISA Kit
E-EL-M0150	Mouse AQP-3 (Aquaporin 3) ELISA Kit
E-EL-M0151	Mouse AQP-4 (Aquaporin 4) ELISA Kit
E-EL-M0154	Mouse Arg (Arginase) ELISA Kit

Cat.NO	Product Name
E-EL-M0155	Mouse AVPR1A (Arginine Vasopressin Receptor 1A) ELISA Kit
E-EL-M0156	Mouse AhR (Aryl Hydrocarbon Receptor) ELISA Kit
E-EL-M0157	Mouse AIP (Aryl Hydrocarbon Receptor Interacting Protein) ELISA Kit
E-EL-M0158	Mouse ARNT2 (Aryl Hydrocarbon Receptor Nuclear Translocator 2) ELISA Kit
E-EL-M0159	Mouse ASGR1 (Asialoglycoprotein Receptor 1) ELISA Kit
E-EL-M0160	Mouse AST (Aspartate Aminotransferase) ELISA Kit
E-EL-M0161	Mouse ATM (Ataxia Telangiectasia Mutated) ELISA Kit
E-EL-M0162	Mouse ATXN1 (Ataxin 1) ELISA Kit
E-EL-M0163	Mouse ATXN2 (Ataxin 2) ELISA Kit
E-EL-M0164	Mouse ABCA1 (ATP Binding Cassette Transporter A1) ELISA Kit
E-EL-M0165	Mouse ABCG1 (ATP Binding Cassette Transporter G1) ELISA Kit
E-EL-M0166	Mouse ANP (Atrial Natriuretic Peptide) ELISA Kit
E-EL-M0167	Mouse AXIN2 (Axis Inhibition Protein 2) ELISA Kit
E-EL-M0170	Mouse bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit
E-EL-M0171	Mouse FGF4 (Fibroblast Growth Factor 4) ELISA Kit
E-EL-M0172	Mouse FGF6 (Fibroblast Growth Factor 6) ELISA Kit
E-EL-M0174	Mouse BAFFR (B-Cell Activation Factor Receptor) ELISA Kit
E-EL-M0175	Mouse Bcl-2 (B-Cell Leukemia/Lymphoma 2) ELISA Kit
E-EL-M0176	Mouse Bcl-3 (B-Cell CLL/Lymphoma 3) ELISA Kit
E-EL-M0177	Mouse BAK1 (BCL2 Antagonist/Killer 1) ELISA Kit

Cat.NO	Product Name
E-EL-M0178	Mouse Bax (Bcl-2 Associated X Protein) ELISA Kit
E-EL-M0179	Mouse BMF (Bcl-2 Modifying Factor) ELISA Kit
E-EL-M0180	Mouse BAD (Bcl-2 Associated Death Promoter) ELISA Kit
E-EL-M0181	Mouse BTC (Betacellulin) ELISA Kit
E-EL-M0182	Mouse BACE1 (Beta Site APP Cleaving Enzyme 1) ELISA Kit
E-EL-M0183	Mouse β -EPR (Beta-Endorphin Receptor) ELISA Kit
E-EL-M0184	Mouse β -EP (Beta-Endorphin) ELISA Kit
E-EL-M0186	Mouse BID (BH3 Interacting Domain Death Agonist) ELISA Kit
E-EL-M0187	Mouse IL-25 (Interleukin 25) ELISA Kit
E-EL-M0189	Mouse BTA (Bladder Tumor Antigen) ELISA Kit
E-EL-M0190	Mouse BLC (B-Lymphocyte Chemoattractant) ELISA Kit
E-EL-M0192	Mouse BMP-1 (Bone Morphogenetic Protein 1) ELISA Kit
E-EL-M0193	Mouse BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit
E-EL-M0194	Mouse BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit
E-EL-M0195	Mouse BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit
E-EL-M0196	Mouse BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit
E-EL-M0197	Mouse BMP8B (Bone Morphogenetic Protein 8B) ELISA Kit
E-EL-M0198	Mouse BMPR1A (Bone Morphogenetic Protein Receptor 1A) ELISA Kit
E-EL-M0199	Mouse BMPR2 (Bone Morphogenetic Protein Receptor II) ELISA Kit
E-EL-M0200	Mouse BALP (Bone Alkaline Phosphatase) ELISA Kit

Cat.NO	Product Name
E-EL-M0201	Mouse CASP1 (Caspase 1) ELISA Kit
E-EL-M0202	Mouse BK (Bradykinin) ELISA Kit
E-EL-M0203	Mouse BDNF (Brain Derived Neurotrophic Factor) ELISA Kit
E-EL-M0204	Mouse BNP (Brain Natriuretic Peptide) ELISA Kit
E-EL-M0206	Mouse BRAK (Breast And Kidney Expressed Chemokine) ELISA Kit
E-EL-M0208	Mouse CA724 (Tumor Marker CA724) ELISA Kit
E-EL-M0209	Mouse Ca-ATPase (Calcium ATPase) ELISA Kit
E-EL-M0210	Mouse VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit
E-EL-M0211	Mouse E-Cad (E-Cadherin) ELISA Kit
E-EL-M0212	Mouse NCAD (Neural Cadherin) ELISA Kit
E-EL-M0213	Mouse CABIN1 (Calcineurin Binding Protein 1) ELISA Kit
E-EL-M0214	Mouse CaN (Calcineurin) ELISA Kit
E-EL-M0215	Mouse CGRP (Calcitonin Gene Related Peptide) ELISA Kit
E-EL-M0216	Mouse CTR (Calcitonin Receptor) ELISA Kit
E-EL-M0217	Mouse CT (Calcitonin) ELISA Kit
E-EL-M0218	Mouse CCB (Calcium Channel Blockers) ELISA Kit
E-EL-M0219	Mouse CaSR (Calcium Sensing Receptor) ELISA Kit
E-EL-M0220	Mouse CACYBP (Calcyclin Binding Protein) ELISA Kit
E-EL-M0221	Mouse CIAP (Calf Intestinal Alkaline Phosphatase) ELISA Kit
E-EL-M0222	Mouse CNX (Calnexin) ELISA Kit

Cat.NO	Product Name
E-EL-M0223	Mouse CAST (Calpastatin) ELISA Kit
E-EL-M0224	Mouse CRT (Calreticulin) ELISA Kit
E-EL-M0225	Mouse CR (Calretinin) ELISA Kit
E-EL-M0226	Mouse CALU (Calumenin) ELISA Kit
E-EL-M0227	Mouse CA9 (Carbonic Anhydrase IX) ELISA Kit
E-EL-M0228	Mouse CA1 (Carbonic Anhydrase I) ELISA Kit
E-EL-M0229	Mouse CA2 (Carbonic Anhydrase 2) ELISA Kit
E-EL-M0231	Mouse PICP (Procollagen I C-terminal Propeptide) ELISA Kit
E-EL-M0232	Mouse CEA (Carcinoembryonic Antign) ELISA Kit
E-EL-M0233	Mouse PINP (Procollagen Type I N-Terminal Propeptide) ELISA Kit
E-EL-M0234	Mouse CLS (Cardiolipin Synthase) ELISA Kit
E-EL-M0235	Mouse CLCF1 (Cardiotrophin Like Cytokine Factor 1) ELISA Kit
E-EL-M0236	Mouse CT-1 (Cardiotrophin 1) ELISA Kit
E-EL-M0237	Mouse COMP (Cartilage Oligomeric Matrix Protein) ELISA Kit
E-EL-M0238	Mouse CASP3 (Caspase 3) ELISA Kit
E-EL-M0239	Mouse CASP4 (Caspase 4) ELISA Kit
E-EL-M0240	Mouse CASP7 (Caspase 7) ELISA Kit
E-EL-M0241	Mouse CASP9 (Caspase 9) ELISA Kit
E-EL-M0242	Mouse CASP10 (Caspase 10) ELISA Kit
E-EL-M0243	Mouse CASP11 (Caspase 11) ELISA Kit

Cat.NO	Product Name
E-EL-M0244	Mouse CASP14 (Caspase 14) ELISA Kit
E-EL-M0245	Mouse CAD (Caspase Activated Dnase) ELISA Kit
E-EL-M0247	Mouse CTNNA1 (Catenin, Alpha 1) ELISA Kit
E-EL-M0248	Mouse CTNNB1 (Catenin, Beta 1) ELISA Kit
E-EL-M0249	Mouse CTSD (Cathepsin D) ELISA Kit
E-EL-M0250	Mouse CTSK (Cathepsin K) ELISA Kit
E-EL-M0251	Mouse CTSL (Cathepsin L) ELISA Kit
E-EL-M0252	Mouse CDX2 (Caudal Type Homeobox Transcription Factor 2) ELISA Kit
E-EL-M0253	Mouse CAV1 (Caveolin 1) ELISA Kit
E-EL-M0254	Mouse CEBPD (CCAAT/Enhancer Binding Protein Delta) ELISA Kit
E-EL-M0255	Mouse CDC24 (Cell Division Cycle 24) ELISA Kit
E-EL-M0256	Mouse CENPA (Centromere Protein A) ELISA Kit
E-EL-M0257	Mouse CENPB (Centromere Protein B) ELISA Kit
E-EL-M0258	Mouse CENPE (Centromere Protein E) ELISA Kit
E-EL-M0259	Mouse CENPH (Centromere Protein H) ELISA Kit
E-EL-M0260	Mouse CENPI (Centromere Protein I) ELISA Kit
E-EL-M0261	Mouse CEP110 (Centrosomal Protein 110kDa) ELISA Kit
E-EL-M0262	Mouse CP (Ceruloplasmin) ELISA Kit
E-EL-M0263	Mouse CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit
E-EL-M0264	Mouse CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-M0265	Mouse CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit
E-EL-M0266	Mouse CCRL1 (Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
E-EL-M0267	Mouse CX3CL1 (Chemokine C-X3-C-Motif Ligand 1) ELISA Kit
E-EL-M0268	Mouse CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
E-EL-M0269	Mouse CXCL15 (Chemokine C-X-C-Motif Ligand 15) ELISA Kit
E-EL-M0270	Mouse CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit
E-EL-M0271	Mouse CKLF (Chemokine Like Factor) ELISA Kit
E-EL-M0272	Mouse CKLFSF4 (Chemokine Like Factor Superfamily 4) ELISA Kit
E-EL-M0273	Mouse CKLFSF5 (Chemokine Like Factor Superfamily 5) ELISA Kit
E-EL-M0274	Mouse CKLFSF6 (Chemokine Like Factor Superfamily 6) ELISA Kit
E-EL-M0275	Mouse CKLFSF7 (Chemokine Like Factor Superfamily 7) ELISA Kit
E-EL-M0276	Mouse CKLFSF8 (Chemokine Like Factor Superfamily 8) ELISA Kit
E-EL-M0277	Mouse CCKAR (Cholecystokinin A Receptor) ELISA Kit
E-EL-M0278	Mouse CCKBR (Cholecystokinin B Receptor) ELISA Kit
E-EL-M0280	Mouse CETP (cholest erolester transfer protein) ELISA Kit
E-EL-M0281	Mouse PC/CPG (choline phosphoglyceride) ELISA Kit
E-EL-M0282	Mouse Gas6 (Growth Arrest Specific Protein 6) ELISA Kit
E-EL-M0283	Mouse CS (chondroitin sulfate) ELISA Kit
E-EL-M0285	Mouse CG (Chorionic Gonadotrophin) ELISA Kit
E-EL-M0286	Mouse CBX3 (Chromobox Homolog 3) ELISA Kit

Cat.NO	Product Name
E-EL-M0287	Mouse CHD3 (Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
E-EL-M0288	Mouse CHD5 (Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
E-EL-M0289	Mouse CNTF (Ciliary Neurotrophic Factor) ELISA Kit
E-EL-M0290	Mouse CROP (Cisplatin Resistance Associated Overexpressed Protein) ELISA Kit
E-EL-M0291	Mouse C-jun ELISA Kit
E-EL-M0292	Mouse CC17 (Clara Cell Protein) ELISA Kit
E-EL-M0293	Mouse CD3 (Cluster of differentiation 3) ELISA Kit
E-EL-M0294	Mouse CD4 (Cluster Of Differentiation 4) ELISA Kit
E-EL-M0295	Mouse CD19 (Cluster of differentiation 19) ELISA Kit
E-EL-M0296	Mouse CD30 (Cluster of differentiation 30) ELISA Kit
E-EL-M0297	Mouse CD30L (Cluster of Differentiation 30 Ligand) ELISA Kit
E-EL-M0298	Mouse CD40 (Cluster of Differentiation 40) ELISA Kit
E-EL-M0299	Mouse CD40L (Cluster of Differentiation 40 Ligand) ELISA Kit
E-EL-M0300	Mouse CD97 (Cluster of Differentiation 97) ELISA Kit
E-EL-M0301	Mouse CD100 (Cluster of Differentiation 100) ELISA Kit
E-EL-M0302	Mouse CD109 (Cluster of Differentiation 109) ELISA Kit
E-EL-M0303	Mouse CLU (Clusterin) ELISA Kit
E-EL-M0304	Mouse F II (coagulation factor II) ELISA Kit
E-EL-M0305	Mouse GP4/CD36 (Platelet Membrane Glycoprotein IV) ELISA Kit
E-EL-M0306	Mouse FV (Coagulation Factor V) ELISA Kit

Cat.NO	Product Name
E-EL-M0307	Mouse FVIII (Coagulation Factor VIII) ELISA Kit
E-EL-M0308	Mouse FIX(Coagulation Factor IX)ELISA Kit
E-EL-M0309	Mouse FX (coagulation factor X) ELISA Kit
E-EL-M0310	Mouse FXII (Coagulation Factor XII) ELISA Kit
E-EL-M0311	Mouse F13B (Coagulation Factor XIII B Polypeptide) ELISA Kit
E-EL-M0312	Mouse CFL1 (Cofilin 1, Non-Muscle) ELISA Kit
E-EL-M0313	Mouse CFL2 (Cofilin 2, Muscle) ELISA Kit
E-EL-M0315	Mouse COL2 (Collagen Type II) ELISA Kit
E-EL-M0316	Mouse COL3 (Collagen Type III) ELISA Kit
E-EL-M0317	Mouse COL4 (Collagen Type IV) ELISA Kit
E-EL-M0318	Mouse COL1α2 (Collagen Type I Alpha 2) ELISA Kit
E-EL-M0319	Mouse COL4α1 (Collagen Type IV Alpha 1) ELISA Kit
E-EL-M0320	Mouse COL10 (Collagen Type X) ELISA Kit
E-EL-M0322	Mouse Collagenase II ELISA Kit
E-EL-M0323	Mouse CCA (Colon Cancer Antigen) ELISA Kit
E-EL-M0324	Mouse CSF2Ra (Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit
E-EL-M0325	Mouse COL1 (Collagen Type I) ELISA Kit
E-EL-M0326	Mouse DCLK1 (Doublecortin Like Kinase 1) ELISA Kit
E-EL-M0327	Mouse GCSFR (Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
E-EL-M0328	Mouse GPC3 (Glypican 3) ELISA Kit

Cat.NO	Product Name
E-EL-M0329	Mouse C1INH (Complement 1 Inhibitor) ELISA Kit
E-EL-M0330	Mouse C3 (Complement 3) ELISA Kit
E-EL-M0331	Mouse C4 (Complement Component 4) ELISA Kit
E-EL-M0332	Mouse C1qR1 (Complement Component 1Q Receptor1) ELISA Kit
E-EL-M0333	Mouse C5 (Complement Component 5) ELISA Kit
E-EL-M0334	Mouse CFB (Complement Factor B) ELISA Kit
E-EL-M0335	Mouse CFD (Complement Factor D) ELISA Kit
E-EL-M0336	Mouse CFH (Complement Factor H) ELISA Kit
E-EL-M0337	Mouse C3a (Complement Component 3a) ELISA Kit
E-EL-M0338	Mouse C4a (Complement Component 4a) ELISA Kit
E-EL-M0339	Mouse C5a (Complement Component 5a) ELISA Kit
E-EL-M0340	Mouse CTGF (Connective Tissue Growth Factor) ELISA Kit
E-EL-M0341	Mouse CX26 (Connexin 26) ELISA Kit
E-EL-M0342	Mouse CX31 (Connexin 31) ELISA Kit
E-EL-M0343	Mouse CX37 (Connexin 37) ELISA Kit
E-EL-M0344	Mouse CX40 (Connexin 40) ELISA Kit
E-EL-M0345	Mouse CX43 (Connexin 43) ELISA Kit
E-EL-M0346	Mouse CNTN1 (Contactin 1) ELISA Kit
E-EL-M0347	Mouse CNTN2 (Contactin 2) ELISA Kit
E-EL-M0348	Mouse CNTN3 (Contactin 3) ELISA Kit

Cat.NO	Product Name
E-EL-M0349	Mouse CORT (Corticosterone) ELISA Kit
E-EL-M0350	Mouse SIRT1 (NAD-dependent protein deacetylase sirtuin-1) ELISA Kit
E-EL-M0351	Mouse CRH (Corticotropin Releasing Hormone) ELISA Kit
E-EL-M0354	Mouse C-P (C-Peptide) ELISA Kit
E-EL-M0355	Mouse CK-MB (Creatine Kinase MB Isoenzyme) ELISA Kit
E-EL-M0356	Mouse CKMT1A (Creatine Kinase, Mitochondrial 1A) ELISA Kit
E-EL-M0357	Mouse CKMT1B (Creatine Kinase, Mitochondrial 1B) ELISA Kit
E-EL-M0358	Mouse CKM (Creatine Kinase, Muscle) ELISA Kit
E-EL-M0359	Mouse AGEs (Advanced Glycation End Products) ELISA Kit
E-EL-M0360	Mouse NTXI (Cross Linked N-Telopeptide of Type I Collagen) ELISA Kit
E-EL-M0361	Mouse α -CTx (Alpha Crosslaps) ELISA Kit
E-EL-M0363	Mouse CICP (cross-linked C-telopeptides of Type I collagen) ELISA Kit
E-EL-M0364	Mouse CRY1 (Cryptochrome 1) ELISA Kit
E-EL-M0365	Mouse CRYaB (Crystallin Alpha B) ELISA Kit
E-EL-M0366	Mouse CTXI (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
E-EL-M0367	Mouse PIIICP (c-terminal Procollagen III propeptide) ELISA Kit
E-EL-M0368	Mouse CTX-2 (Cytotoxin 2) ELISA Kit
E-EL-M0369	Mouse CLEC2C (C-Type Lectin Domain Family 2, Member C) ELISA Kit
E-EL-M0370	Mouse CNP (C-type natriuretic peptide) ELISA Kit
E-EL-M0371	Mouse CTACK (Cutaneous T-Cell Attracting Chemokine) ELISA Kit

Cat.NO	Product Name
E-EL-M0372	Mouse β -CTx (Beta Crosslaps) ELISA Kit
E-EL-M0373	Mouse CXCR3 (CXC-chemokine receptor 3) ELISA Kit
E-EL-M0374	Mouse COL2 (Collagen Type II) ELISA Kit
E-EL-M0375	Mouse CREB (Cyclic AMP Response Element Binding Protein) ELISA Kit
E-EL-M0376	Mouse c-myc (c-myc Oncogene Product) ELISA Kit
E-EL-M0377	Mouse CDK5 (Cyclin Dependent Kinase 5) ELISA Kit
E-EL-M0378	Mouse CDK9 (Cyclin Dependent Kinase 9) ELISA Kit
E-EL-M0380	Mouse CDKN1B (Cyclin Dependent Kinase Inhibitor 1B) ELISA Kit
E-EL-M0381	Mouse CCND1 (Cyclin-D1) ELISA Kit
E-EL-M0382	Mouse CCND2 (Cyclin-D2) ELISA Kit
E-EL-M0383	Mouse CCND3 (Cyclin-D3) ELISA Kit
E-EL-M0384	Mouse CYPA(Cyclophilin A)ELISA Kit
E-EL-M0385	Mouse CYPB (Cyclophilin B) ELISA Kit
E-EL-M0386	Mouse CYPD (Cyclophilin D) ELISA Kit
E-EL-M0389	Mouse Cys-C (Cystatin C) ELISA Kit
E-EL-M0390	Mouse CYSLTR2 (Cysteinyl Leukotriene Receptor 2) ELISA Kit
E-EL-M0391	Mouse CYB561 (Cytochrome b-561) ELISA Kit
E-EL-M0392	Mouse FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit
E-EL-M0393	Mouse Cyt-C (Cytochrome C) ELISA Kit
E-EL-M0394	Mouse CYP7A1 (Cytochrome P450 7A1) ELISA Kit

Cat.NO	Product Name
E-EL-M0395	Mouse CCO (cytochrom-oxidase) ELISA Kit
E-EL-M0396	Mouse CYTH1 (Cytohesin 1) ELISA Kit
E-EL-M0397	Mouse CYTIP (Cytohesin Interacting Protein) ELISA Kit
E-EL-M0398	Mouse CTLA4 (Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit
E-EL-M0400	Mouse D2D (D-Dimer) ELISA Kit
E-EL-M0401	Mouse DAP6 (Death Associated Protein 6) ELISA Kit
E-EL-M0402	Mouse DCN (Decorin) ELISA Kit
E-EL-M0403	Mouse DEFA5 (Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
E-EL-M0405	Mouse DHEA-S7 (Dehydroepiandrosterone S7) ELISA Kit
E-EL-M0406	Mouse DHEA-S (Dehydroepiandrosterone sulfate) ELISA Kit
E-EL-M0407	Mouse DHEA (Dehydroepiandrosterone) ELISA Kit
E-EL-M0408	Mouse DPD (Deoxypyridinoline) ELISA Kit
E-EL-M0409	Mouse DNASE1 (Deoxyribonuclease I) ELISA Kit
E-EL-M0410	Mouse VAP-1 (Vascular Adhesion Protein 1) ELISA Kit
E-EL-M0411	Mouse Des (Desmin) ELISA Kit
E-EL-M0412	Mouse DAO (Diamine Oxidase) ELISA Kit
E-EL-M0413	Mouse H3 (Histone H3) ELISA Kit
E-EL-M0414	Mouse DKK2 (Dickkopf Related Protein 2) ELISA Kit
E-EL-M0415	Mouse DLD (Dihydrolipoyl Dehydrogenase) ELISA Kit
E-EL-M0416	Mouse DLAT (Dihydrolipoyl Transacetylase) ELISA Kit

Cat.NO	Product Name
E-EL-M0418	Mouse DPPIV (dipeptidyl peptidase IV) ELISA Kit
E-EL-M0419	Mouse D-LDH (D-Lactate Dehydrogenase) ELISA Kit
E-EL-M0420	Mouse D1R (Dopamine D1 receptor) ELISA Kit
E-EL-M0421	Mouse DRD3 (Dopamine Receptor D3) ELISA Kit
E-EL-M0422	Mouse DRD4 (Dopamine Receptor D4) ELISA Kit
E-EL-M0423	Mouse DAT (Dopamine Transporter) ELISA Kit
E-EL-M0424	Mouse GAPDH (Glyceraldehyde-3-Phosphate Dehydrogenase) ELISA Kit
E-EL-M0425	Mouse DNM1 (Dynamin 1) ELISA Kit
E-EL-M0426	Mouse DMD (Dystrophin) ELISA Kit
E-EL-M0427	Mouse SCARB1 (Scavenger receptor class B member 1) ELISA Kit
E-EL-M0430	Mouse EGR1 (Early Growth Response Protein 1) ELISA Kit
E-EL-M0431	Mouse EGR2 (Early Growth Response Protein 2) ELISA Kit
E-EL-M0432	Mouse EGR3 (Early Growth Response Protein 3) ELISA Kit
E-EL-M0433	Mouse EGR4 (Early Growth Response Protein 4) ELISA Kit
E-EL-M0443	Mouse DKK3 (Dickkopf Related Protein 3) ELISA Kit
E-EL-M0444	Mouse NE (Neutrophil Elastase) ELISA Kit
E-EL-M0445	Mouse ELA4 (Elastase 4) ELISA Kit
E-EL-M0446	Mouse ELN (Elastin) ELISA Kit
E-EL-M0447	Mouse EMILIN1 (Elastin Microfibril Interface Located Protein 1) ELISA Kit
E-EL-M0448	Mouse ELOA (Elongin A) ELISA Kit

Cat.NO	Product Name
E-EL-M0449	Mouse EED (Embryonic Ectoderm Development) ELISA Kit
E-EL-M0452	Mouse ERN1 (Endoplasmic Reticulum To Nucleus Signalling 1) ELISA Kit
E-EL-M0453	Mouse ES (Endostatin) ELISA Kit
E-EL-M0454	Mouse EL (Endothelial lipase) ELISA Kit
E-EL-M0456	Mouse NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit
E-EL-M0457	Mouse ET-1 (Endothelin 1) ELISA Kit
E-EL-M0458	Mouse ECE1 (Endothelin Converting Enzyme 1) ELISA Kit
E-EL-M0459	Mouse ETRA (Endothelin Receptor A) ELISA Kit
E-EL-M0461	Mouse ELMO2 (Engulfment And Cell Motility 2) ELISA Kit
E-EL-M0462	Mouse Ear 14 (Eosinophil-associated, ribonuclease A family, member 14) ELISA Kit
E-EL-M0463	Mouse ECFR (Eosinophil chemotactic factor receptor) ELISA Kit
E-EL-M0464	Mouse ECF/CCL11 (Eosinophil Chemotactic Factor) ELISA Kit
E-EL-M0466	Mouse MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
E-EL-M0470	Mouse EPCAM (Epithelial Cell Adhesion Molecule) ELISA Kit
E-EL-M0471	Mouse ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit
E-EL-M0472	Mouse EPOR (Erythropoietin receptor) ELISA Kit
E-EL-M0473	Mouse UCP-2 (Uncoupling Protein 2, Mitochondrial) ELISA Kit
E-EL-M0474	Mouse SELE (E-Selectin) ELISA Kit
E-EL-M0475	Mouse ER (Estradiol Receptor) ELISA Kit
E-EL-M0476	Mouse ER α (Estrogen Receptor Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-M0478	Mouse EIF3A (Eukaryotic Translation Initiation Factor 3A) ELISA Kit
E-EL-M0479	Mouse ECM1 (Extracellular Matrix Protein 1) ELISA Kit
E-EL-M0480	Mouse ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit
E-EL-M0481	Mouse ERK2 (Extracellular Signal Regulated Kinase 2) ELISA Kit
E-EL-M0482	Mouse FVIIIa (Factor VIII-related Antigen) ELISA Kit
E-EL-M0483	Mouse BF (factor B) ELISA Kit
E-EL-M0484	Mouse FAS/CD95 (Factor Related Apoptosis) ELISA Kit
E-EL-M0485	Mouse FSCN (Fascin) ELISA Kit
E-EL-M0487	Mouse FABP5 (Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
E-EL-M0488	Mouse FABP6 (Fatty Acid Binding Protein 6, Ileal) ELISA Kit
E-EL-M0489	Mouse FBXO2 (F-Box Protein 2) ELISA Kit
E-EL-M0490	Mouse ER β (Estrogen Receptor Beta) ELISA Kit
E-EL-M0491	Mouse FE (Ferritin) ELISA Kit
E-EL-M0492	Mouse FBL (Fibrillarlin) ELISA Kit
E-EL-M0493	Mouse FBN1 (Fibrillin 1) ELISA Kit
E-EL-M0494	Mouse FGA (Fibrinogen Alpha chain) ELISA Kit
E-EL-M0495	Mouse FGB (Fibrinogen Beta chain) ELISA Kit
E-EL-M0496	Mouse FDP (Fibrinogen Degradation Product) ELISA Kit
E-EL-M0497	Mouse FGL2 (Fibrinogen Like Protein 2) ELISA Kit
E-EL-M0498	Mouse Fbg (Fibrinogen) ELISA Kit

Cat.NO	Product Name
E-EL-M0499	Mouse FGF7 (Fibroblast Growth Factor 7) ELISA Kit
E-EL-M0500	Mouse FGF9 (Fibroblast Growth Factor 9) ELISA Kit
E-EL-M0501	Mouse FGFBP1 (Fibroblast Growth Factor Binding Protein 1) ELISA Kit
E-EL-M0502	Mouse FRS2 (Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit
E-EL-M0503	Mouse FRS3 (Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
E-EL-M0504	Mouse FMOD (Fibromodulin) ELISA Kit
E-EL-M0505	Mouse FLRT2 (Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit
E-EL-M0506	Mouse FN (Fibronectin) ELISA Kit
E-EL-M0507	Mouse FLNa (Filamin A) ELISA Kit
E-EL-M0508	Mouse Flt3 (FMS-like Tyrosine Kinase 3) ELISA Kit
E-EL-M0509	Mouse FAK (Focal Adhesion Kinase) ELISA Kit
E-EL-M0511	Mouse FSH (follicle-stimulating hormone) ELISA Kit
E-EL-M0512	Mouse FST (Follistatin) ELISA Kit
E-EL-M0513	Mouse FOXP3 (Forkhead Box Protein P3) ELISA Kit
E-EL-M0514	Mouse FPR2 (Formyl Peptide Receptor 2) ELISA Kit
E-EL-M0515	Mouse FK (Fractalkine) ELISA Kit
E-EL-M0516	Mouse GAD-Ab-IgM (Glutamic Acid Decarboxylase Autoantibody IgM) ELISA Kit
E-EL-M0518	Mouse F-TESTO (Free Testosterone) ELISA Kit
E-EL-M0519	Mouse fT4 (Free Thyroxine) ELISA Kit
E-EL-M0520	Mouse fT3 (Free Triiodothyronine) ELISA Kit

Cat.NO	Product Name
E-EL-M0521	Mouse FLI1 (Friend Leukemia Virus Integration 1) ELISA Kit
E-EL-M0522	Mouse FUCA (Alpha-L-Fucosidase, Tissue) ELISA Kit
E-EL-M0523	Mouse GNB1 (G Protein Beta 1) ELISA Kit
E-EL-M0524	Mouse Gal-6S (galactose-6-sulphate sulphotase) ELISA Kit
E-EL-M0525	Mouse GLB (Galactosidase, Beta) ELISA Kit
E-EL-M0526	Mouse GALR1 (Galanin Receptor 1) ELISA Kit
E-EL-M0527	Mouse GAL1 (Galectin 1) ELISA Kit
E-EL-M0528	Mouse GAL2 (Galectin 2) ELISA Kit
E-EL-M0529	Mouse GAL3 (Galectin 3) ELISA Kit
E-EL-M0530	Mouse GAL4 (Galectin 4) ELISA Kit
E-EL-M0531	Mouse GAL6 (Galectin 6) ELISA Kit
E-EL-M0532	Mouse GAL7 (Galectin 7) ELISA Kit
E-EL-M0533	Mouse GAL8 (Galectin 8) ELISA Kit
E-EL-M0534	Mouse GAL9 (Galectin 9) ELISA Kit
E-EL-M0535	Mouse GAL12 (Galectin 12) ELISA Kit
E-EL-M0536	Mouse gGT1 (Gamma Glutamyltransferase 1) ELISA Kit
E-EL-M0538	Mouse GJB1 (Gap Junction Protein Beta 1) ELISA Kit
E-EL-M0539	Mouse GIP (Gastric Inhibitory Polypeptide) ELISA Kit
E-EL-M0540	Mouse GT-34 (Gastrin 34) ELISA Kit
E-EL-M0541	Mouse CA199 (Gastrointestinalcancer Marker-CA199) ELISA Kit

Cat.NO	Product Name
E-EL-M0542	Mouse GKN2 (Gastrokine 2) ELISA Kit
E-EL-M0543	Mouse GATA1 (GATA Binding Protein 1) ELISA Kit
E-EL-M0544	Mouse GATA2 (GATA Binding Protein 2) ELISA Kit
E-EL-M0545	Mouse GATA3 (GATA Binding Protein 3) ELISA Kit
E-EL-M0546	Mouse GATA4 (GATA Binding Protein 4) ELISA Kit
E-EL-M0547	Mouse GATA5 (GATA Binding Protein 5) ELISA Kit
E-EL-M0548	Mouse GS (Gelsolin) ELISA Kit
E-EL-M0549	Mouse GTFIIB (General Transcription Factor IIB) ELISA Kit
E-EL-M0550	Mouse GTF2H1 (General Transcription Factor IIH, Polypeptide 1) ELISA Kit
E-EL-M0551	Mouse GHRL (Ghrelin) ELISA Kit
E-EL-M0552	Mouse GMFB (Glia Maturation Factor, Beta) ELISA Kit
E-EL-M0553	Mouse GDNF (Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
E-EL-M0554	Mouse GFAP (glial fibrillary acidic protein) ELISA Kit
E-EL-M0555	Mouse GC (Glucagon) ELISA Kit
E-EL-M0556	Mouse GR (glucocorticoid receptor) ELISA Kit
E-EL-M0557	Mouse GR α (glucocorticoid receptor α) ELISA Kit
E-EL-M0558	Mouse GR β (glucocorticoid receptor β) ELISA Kit
E-EL-M0559	Mouse GCK (Glucokinase) ELISA Kit
E-EL-M0560	Mouse gp130 (Glucoprotein 130) ELISA Kit
E-EL-M0561	Mouse G6PD (Glucose 6 Phosphate Dehydrogenase) ELISA Kit

Cat.NO	Product Name
E-EL-M0562	Mouse GLUT2 (Glucose Transporter 2) ELISA Kit
E-EL-M0563	Mouse GLUT3 (Glucose Transporter 3) ELISA Kit
E-EL-M0564	Mouse GLUT4 (Glucose Transporter 4) ELISA Kit
E-EL-M0565	Mouse GPI (Glucose 6 Phosphate Isomerase) ELISA Kit
E-EL-M0567	Mouse GUSb (Glucuronidase, Beta) ELISA Kit
E-EL-M0568	Mouse GCLM (Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit
E-EL-M0569	Mouse GAD2 (Glutamate Decarboxylase 2, Acid) ELISA Kit
E-EL-M0570	Mouse GDH/GLDH (Glutamate dehydrogenase) ELISA Kit
E-EL-M0571	Mouse GAD-Ab (Glutamic Acid Decarboxylase Autoantibody) ELISA Kit
E-EL-M0572	Mouse OXB (Orexin B) ELISA Kit
E-EL-M0573	Mouse HB(Hemoglobin)ELISA Kit
E-EL-M0574	Mouse GR (Glutathione Reductase) ELISA Kit
E-EL-M0575	Mouse GSTa1 (Glutathione S Transferase Alpha 1) ELISA Kit
E-EL-M0576	Mouse GSTo1 (Glutathione S Transferase Omega 1) ELISA Kit
E-EL-M0577	Mouse GSTt1 (Glutathione S Transferase Theta 1) ELISA Kit
E-EL-M0578	Mouse GSTt2 (Glutathione S Transferase Theta 2) ELISA Kit
E-EL-M0579	Mouse HbA1C (Glycosylated Hemoglobin/Hemoglobin A1c) ELISA Kit
E-EL-M0583	Mouse PYGB (Glycogen Phosphorylase, Brain) ELISA Kit
E-EL-M0584	Mouse PYGL (Glycogen Phosphorylase, Liver) ELISA Kit
E-EL-M0585	Mouse PYGM (Glycogen Phosphorylase, Muscle) ELISA Kit

Cat.NO	Product Name
E-EL-M0586	Mouse GSK3a (Glycogen Synthase Kinase 3 Alpha) ELISA Kit
E-EL-M0587	Mouse GSK (glycogen synthase kinase) ELISA Kit
E-EL-M0588	Mouse GP5 (Glycoprotein V, Platelet) ELISA Kit
E-EL-M0589	Mouse GPC4 (Glypican 4) ELISA Kit
E-EL-M0591	Mouse GRN (Granulin) ELISA Kit
E-EL-M0592	Mouse GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit
E-EL-M0593	Mouse Gzms-A (granzyme A) ELISA Kit
E-EL-M0594	Mouse Gzms-B (granzyme B) ELISA Kit
E-EL-M0595	Mouse GREM1 (GremLin 1) ELISA Kit
E-EL-M0596	Mouse GDF1 (Growth Differentiation Factor 1) ELISA Kit
E-EL-M0597	Mouse GDF2 (Growth Differentiation Factor 2) ELISA Kit
E-EL-M0598	Mouse GDF3 (Growth Differentiation Factor 3) ELISA Kit
E-EL-M0599	Mouse GDF5 (Growth Differentiation Factor 5) ELISA Kit
E-EL-M0600	Mouse GDF7 (Growth Differentiation Factor 7) ELISA Kit
E-EL-M0601	Mouse GDF9 (Growth Differentiation Factor 9) ELISA Kit
E-EL-M0602	Mouse GDF10 (Growth Differentiation Factor 10) ELISA Kit
E-EL-M0603	Mouse GDF11(Growth Differentiation Factor 11)ELISA Kit
E-EL-M0604	Mouse GDF15 (Growth Differentiation Factor 15) ELISA Kit
E-EL-M0605	Mouse Grb2 (Growth Factor Receptor Bound Protein 2) ELISA Kit
E-EL-M0606	Mouse Grb7 (Growth Factor Receptor Bound Protein 7) ELISA Kit

Cat.NO	Product Name
E-EL-M0607	Mouse Grb10 (Growth Factor Receptor Bound Protein 10) ELISA Kit
E-EL-M0608	Mouse Grb14 (Growth Factor Receptor Bound Protein 14) ELISA Kit
E-EL-M0609	Mouse GH2 (Growth Hormone 2) ELISA Kit
E-EL-M0610	Mouse GHRH (Growth Hormone Releasing Hormone) ELISA Kit
E-EL-M0611	Mouse CAM (Calmodulin) ELISA Kit
E-EL-M0613	Mouse HP (Haptoglobin) ELISA Kit
E-EL-M0614	Mouse HSF1 (Heat Shock Factor 1) ELISA Kit
E-EL-M0615	Mouse HSP-20 (Heat Shock Protein 20) ELISA Kit
E-EL-M0616	Mouse HSP-27 (heat shock protein 27) ELISA Kit
E-EL-M0617	Mouse HSP-40 (Heat Shock Protein 40) ELISA Kit
E-EL-M0618	Mouse HSP-60 (Heat Shock Protein 60) ELISA Kit
E-EL-M0619	Mouse HSP-70 (Heat Shock Protein 70) ELISA Kit
E-EL-M0620	Mouse HSP-90 (Heat Shock Protein 90) ELISA Kit
E-EL-M0621	Mouse HSP gp96 (Heat Shock Protein glycoprotein 96) ELISA Kit
E-EL-M0622	Mouse HSF2 (Heat Shock Transcription Factor 2) ELISA Kit
E-EL-M0623	Mouse SHH (Hedgehog Homolog, Sonic) ELISA Kit
E-EL-M0624	Mouse HHIP (Hedgehog Interacting Protein) ELISA Kit
E-EL-M0625	Mouse HLTF (Helicase Like Transcription Factor) ELISA Kit
E-EL-M0626	Mouse MMP-3 (Matrix Metalloproteinase 3) ELISA Kit
E-EL-M0627	Mouse MMP-9 (Matrix Metalloproteinase-9) ELISA Kit

Cat.NO	Product Name
E-EL-M0628	Mouse Pro-MMP-9 (Pro-Matrix Metalloproteinase-9) ELISA Kit
E-EL-M0629	Mouse OPN (Osteopontin) ELISA Kit
E-EL-M0630	Mouse PDGF-AA (Platelet-Derived Growth Factor-AA) ELISA Kit
E-EL-M0631	Mouse PDGF-AB (Platelet-Derived Growth Factor AB) ELISA Kit
E-EL-M0632	Mouse PDGF-BB (Platelet-Derived Growth Factor-BB) ELISA Kit
E-EL-M0633	Mouse POSTN/OSF-2 (Periostin) ELISA Kit
E-EL-M0634	Mouse PCSK9 (Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit
E-EL-M0635	Mouse RETN (Resistin) ELISA Kit
E-EL-M0636	Mouse SCF (Stem Cell Factor) ELISA Kit
E-EL-M0637	Mouse SELL (L-Selectin) ELISA Kit
E-EL-M0638	Mouse SELP (P-Selectin) ELISA Kit
E-EL-M0640	Mouse TPO (Thrombopoietin) ELISA Kit
E-EL-M0641	Mouse TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
E-EL-M0642	Mouse ADIPOR1 (Adiponectin Receptor 1) ELISA Kit
E-EL-M0643	Mouse TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
E-EL-M0644	Mouse RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
E-EL-M0645	Mouse TREM-1 (Triggering Receptor Expresses on Myeloid Cells-1) ELISA Kit
E-EL-M0646	Mouse TSLP (Thymic Stromal Lymphopoietin) ELISA Kit
E-EL-M0647	Mouse sVCAM-1 (soluble vascular cell adhesion molecule 1) ELISA Kit
E-EL-M0648	Mouse FLT1/VEGFR1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-M0649	Mouse VEGFR-2/KDR (Vasculoar Endothelial Growth Factor Receptor 2) ELISA Kit
E-EL-M0651	Mouse HO-1 (Heme Oxygenase 1) ELISA Kit
E-EL-M0652	Mouse HO-2 (Heme Oxygenase 2, Decycling) ELISA Kit
E-EL-M0653	Mouse RANk (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
E-EL-M0654	Mouse HJV (Hemojuvelin) ELISA Kit
E-EL-M0655	Mouse HPA (Heparanase) ELISA Kit
E-EL-M0656	Mouse ALB (Albumin) ELISA Kit
E-EL-M0657	Mouse HC II (Heparin Cofactor II) ELISA Kit
E-EL-M0658	Mouse HB-EGF (Heparin-binding Epidermal Growth Factor-like Growth Factor) ELISA Kit
E-EL-M0659	Mouse LIPC (Lipase, Hepatic) ELISA Kit
E-EL-M0664	Mouse HGF (Hepatocyte Growth Factor) ELISA Kit
E-EL-M0665	Mouse HNF1a (Hepatocyte Nuclear Factor 1 Alpha) ELISA Kit
E-EL-M0666	Mouse HNF1b (Hepatocyte Nuclear Factor 1 Beta) ELISA Kit
E-EL-M0667	Mouse HNF4a (Hepatocyte Nuclear Factor 4 Alpha) ELISA Kit
E-EL-M0668	Mouse HNF6a (Hepatocyte Nuclear Factor 6 Alpha) ELISA Kit
E-EL-M0669	Mouse HNF6b (Hepatocyte Nuclear Factor 6 Beta) ELISA Kit
E-EL-M0670	Mouse HDGF (Hepatoma Derived Growth Factor) ELISA Kit
E-EL-M0671	Mouse Hpc (Hepcidin) ELISA Kit
E-EL-M0672	Mouse HNRPA1 (Heterogeneous Nuclear Ribonucleoprotein A1) ELISA Kit
E-EL-M0673	Mouse HK1 (Hexokinase 1) ELISA Kit

Cat.NO	Product Name
E-EL-M0674	Mouse HK (Hexokinase) ELISA Kit
E-EL-M0676	Mouse HMGB-1 (High mobility group protein B1) ELISA Kit
E-EL-M0677	Mouse hs-CRP (high sensitivity C-Reactive Protein) ELISA Kit
E-EL-M0678	Mouse HDAC (Histone Deacetylase) ELISA Kit
E-EL-M0679	Mouse histon-H2b ELISA Kit
E-EL-M0681	Mouse Bcl10 (B-Cell CLL/Lymphoma 10) ELISA Kit
E-EL-M0682	Mouse HABP (Hya]uronate binding protein) ELISA Kit
E-EL-M0683	Mouse HABP1 (Hyaluronan Binding Protein 1) ELISA Kit
E-EL-M0684	Mouse HAase (Hyaluronidase) ELISA Kit
E-EL-M0687	Mouse HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit
E-EL-M0688	Mouse ID(deiodinase) ELISA Kit
E-EL-M0689	Mouse IgG-Fc ELISA Kit
E-EL-M0690	Mouse IgA (Immunoglobulin A) ELISA Kit
E-EL-M0691	Mouse IgE (Immunoglobulin E) ELISA Kit
E-EL-M0692	Mouse IgG (Immunoglobulin G) ELISA Kit
E-EL-M0693	Mouse IgG1 (Immunoglobulin G1) ELISA Kit
E-EL-M0694	Mouse IgG2a (Immunoglobulin G2a) ELISA Kit
E-EL-M0695	Mouse IgM (Immunoglobulin M) ELISA Kit
E-EL-M0696	Mouse NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) ELISA Kit
E-EL-M0697	Mouse INHA (Inhibin A) ELISA Kit

Cat.NO	Product Name
E-EL-M0698	Mouse INHB (Inhibin B) ELISA Kit
E-EL-M0699	Mouse INHbA (Inhibin Beta A) ELISA Kit
E-EL-M0700	Mouse INHBP (Inhibin Binding Protein) ELISA Kit
E-EL-M0701	Mouse IAA (Insulin Autoantibodies) ELISA Kit
E-EL-M0702	Mouse INSL5 (Insulin Like Protein 5) ELISA Kit
E-EL-M0703	Mouse INSR- β (Insulin Receptor Beta Subunit) ELISA Kit
E-EL-M0704	Mouse IGF-2 (Insulin Like Growth Factor 2) ELISA Kit
E-EL-M0705	Mouse IGFBP-1 (Insulin-Like Growth Factor Binding Protein 1) ELISA Kit
E-EL-M0706	Mouse IGFBP-2 (Insulin-Like Growth Factor Binding Protein 2) ELISA Kit
E-EL-M0707	Mouse IGFBP-4 (Insulin-Like Growth Factor Binding Protein 4) ELISA Kit
E-EL-M0708	Mouse IGFBP-6 (Insulin-Like Growth Factor Binding Protein 6) ELISA Kit
E-EL-M0709	Mouse I-PTH (intact Parathormone) ELISA Kit
E-EL-M0710	Mouse ICAM-1 (intercellular adhesion molecule 1) ELISA Kit
E-EL-M0711	Mouse ICAM-2 (intercellular adhesion molecule 2) ELISA Kit
E-EL-M0712	Mouse ICAM-3(intercellular adhesion molecule 3)ELISA Kit
E-EL-M0713	Mouse IFN- α (Interferon Alpha) ELISA Kit
E-EL-M0714	Mouse FOXO1 (Forkhead Box Protein O1) ELISA Kit
E-EL-M0715	Mouse IFI16/p16 (Interferon-inducible Protein 16) ELISA Kit
E-EL-M0716	Mouse IFN- α -Ab (Interferon- α Antibody) ELISA Kit
E-EL-M0717	Mouse UCP-1 (Uncoupling Protein 1, Mitochondrial) ELISA Kit

Cat.NO	Product Name
E-EL-M0718	Mouse FIL1d (Interleukin 1 Delta) ELISA Kit
E-EL-M0719	Mouse FIL1h (Interleukin 1 Eta) ELISA Kit
E-EL-M0720	Mouse IL-34 (Interleukin 34) ELISA Kit
E-EL-M0721	Mouse IL-3 (Interleukin 3) ELISA Kit
E-EL-M0722	Mouse IL-5 (Interleukin 5) ELISA Kit
E-EL-M0723	Mouse IL-7 (Interleukin 7) ELISA Kit
E-EL-M0724	Mouse IL-9 (Interleukin 9) ELISA Kit
E-EL-M0725	Mouse IL-11 (Interleukin 11) ELISA Kit
E-EL-M0726	Mouse IL-12 (Interleukin 12) ELISA Kit
E-EL-M0727	Mouse IL-13 (Interleukin 13) ELISA Kit
E-EL-M0728	Mouse IL-15 (Interleukin 15) ELISA Kit
E-EL-M0729	Mouse IL-16 (Interleukin 16) ELISA Kit
E-EL-M0730	Mouse IL-18 (Interleukin 18) ELISA Kit
E-EL-M0731	Mouse IL-23 (Interleukin 23) ELISA Kit
E-EL-M0732	Mouse IL-27 (Interleukin 27) ELISA Kit
E-EL-M0733	Mouse IL-35 (Interleukin 35) ELISA Kit
E-EL-M0734	Mouse IL-2R (Interleukin-2 Receptor) ELISA Kit
E-EL-M0735	Mouse IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) ELISA Kit
E-EL-M0736	Mouse IMA (Ischemia Modified Albumin) ELISA Kit
E-EL-M0737	Mouse IL-32 (Interleukin 32) ELISA Kit

Cat.NO	Product Name
E-EL-M0739	Mouse IL18BP (Interleukin 18 Binding Protein) ELISA Kit
E-EL-M0740	Mouse AR (Amphiregulin) ELISA Kit
E-EL-M0744	Mouse KIM-1 (Kidney Injury Molecule 1) ELISA Kit
E-EL-M0746	Mouse LTF/LF (Lactoferrin) ELISA Kit
E-EL-M0747	Mouse LN (Laminin) ELISA Kit
E-EL-M0748	Mouse LDL-IC (Low Density Lipoprotein Immune Complex) ELISA Kit
E-EL-M0749	Mouse LEPR (Leptin Receptor) ELISA Kit
E-EL-M0750	Mouse Lep IgG (Leptospira IgG) ELISA Kit
E-EL-M0751	Mouse LIFR (Leukemia Inhibitory Factor Receptor) ELISA Kit
E-EL-M0753	Mouse LT-C4 (Leukotriene C4) ELISA Kit
E-EL-M0754	Mouse LAT (Linker for activation of T cell) ELISA Kit
E-EL-M0755	Mouse LPL (Lipoprotein Lipase) ELISA Kit
E-EL-M0756	Mouse Lp-a (Lipoprotein a) ELISA Kit
E-EL-M0757	Mouse Lp-PL-A2 (Lipoprotein-associated Phospholipase A2) ELISA Kit
E-EL-M0758	Mouse PAL (L-Phenylalanine ammonia-lyase) ELISA Kit
E-EL-M0760	Mouse LHRH (Luteinizing Hormone-Releasing Hormone) ELISA Kit
E-EL-M0762	Mouse LFA-3 (lymphocyte function associated antigen 3) ELISA Kit
E-EL-M0763	Mouse Lptn/XCL1 (Lymphotactin) ELISA Kit
E-EL-M0764	Mouse LZM (Lysozyme) ELISA Kit
E-EL-M0765	Mouse BSP (Bone Sialoprotein) ELISA Kit

Cat.NO	Product Name
E-EL-M0766	Mouse MIP-1 δ (Macrophage Inflammatory Protein 1 δ) ELISA Kit
E-EL-M0767	Mouse FPN (Ferroportin) ELISA Kit
E-EL-M0768	Mouse ITGB3 (Integrin Beta 3) ELISA Kit
E-EL-M0769	Mouse MIP-3 β (Macrophage Inflammatory Protein 3 β) ELISA Kit
E-EL-M0770	Mouse MIP-5 (Macrophage Inflammatory Protein 5) ELISA Kit
E-EL-M0771	Mouse MIF (Macrophage Migration Inhibitory Factor) ELISA Kit
E-EL-M0772	Mouse H4 (Histone H4) ELISA Kit
E-EL-M0773	Mouse MHCI/H-2 I (Major Histocompatibility Complex I) ELISA Kit
E-EL-M0774	Mouse MHC- II /H-2 II (Major Histocompatibility Complex- II) ELISA Kit
E-EL-M0775	Mouse MHC- III /H-2 III (Major Histocompatibility Complex- III) ELISA Kit
E-EL-M0776	Mouse HAT (Histone Acetyltransferase) ELISA Kit
E-EL-M0778	Mouse CA15-3 (Carbohydrate Antigen 15-3) ELISA Kit
E-EL-M0779	Mouse MMP-1 (Matrix Metalloproteinase 1) ELISA Kit
E-EL-M0780	Mouse MMP-2 (Matrix Metalloproteinase 2) ELISA Kit
E-EL-M0781	Mouse MMP-17 (Matrix Metalloproteinase 17) ELISA Kit
E-EL-M0782	Mouse MMP-24 (Matrix Metalloproteinase 24) ELISA Kit
E-EL-M0783	Mouse MMP-7 (Matrix Metalloproteinase 7) ELISA Kit
E-EL-M0784	Mouse MMP-8 (Matrix Metalloproteinase 8) ELISA Kit
E-EL-M0788	Mouse MT (Melatonin) ELISA Kit
E-EL-M0789	Mouse mIgA (membrane IgA) ELISA Kit

Cat.NO	Product Name
E-EL-M0790	Mouse MHB (Methemoglobin) ELISA Kit
E-EL-M0791	Mouse Methylase ELISA Kit
E-EL-M0792	Mouse MAU (Microalbuminuria) ELISA Kit
E-EL-M0793	Mouse MAO (Monoamine Oxidase) ELISA Kit
E-EL-M0794	Mouse MCP-2 (Monocyte Chemotactic Protein 2) ELISA Kit
E-EL-M0795	Mouse MCP-3 (Monocyte Chemotactic Protein 3) ELISA Kit
E-EL-M0796	Mouse MCP-4 (Monocyte Chemotactic Protein 4) ELISA Kit
E-EL-M0797	Mouse Smad4 (Mothers Against Decapentaplegic Homolog 4) ELISA Kit
E-EL-M0798	Mouse MTL (Motilin) ELISA Kit
E-EL-M0799	Mouse MUC5AC (Mucin-5 subtype AC) ELISA Kit
E-EL-M0800	Mouse MUC5B (Mucin-5 Subtype B) ELISA Kit
E-EL-M0801	Mouse MUC7 (Mucin-7) ELISA Kit
E-EL-M0802	Mouse CHRM1 (Cholinergic Receptor, Muscarinic 1) ELISA Kit
E-EL-M0803	Mouse CHRM2 (Cholinergic Receptor, Muscarinic 2) ELISA Kit
E-EL-M0805	Mouse MBP (Myelin Basic Protein) ELISA Kit
E-EL-M0806	Mouse ApoM (Apolipoprotein M) ELISA Kit
E-EL-M0808	Mouse MYS (Myosin) ELISA Kit
E-EL-M0810	Mouse GAS (N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
E-EL-M0812	Mouse NAGase (N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
E-EL-M0813	Mouse ATPase (ATPase, Na ⁺ /K ⁺ Transporting) ELISA Kit

Cat.NO	Product Name
E-EL-M0814	Mouse NN-T4 (Neonatal Thyroxine) ELISA Kit
E-EL-M0815	Mouse NGF (Nerve growth factor) ELISA Kit
E-EL-M0816	Mouse NCAM/CD56 (Neural Cell Adhesion Molecule) ELISA Kit
E-EL-M0817	Mouse NF (Neurofilament Protein) ELISA Kit
E-EL-M0818	Mouse NGB (Neuroglobin) ELISA Kit
E-EL-M0819	Mouse Slit2 (Slit Homolog 2) ELISA Kit
E-EL-M0820	Mouse NPY (Neuropeptide Y) ELISA Kit
E-EL-M0821	Mouse NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
E-EL-M0822	Mouse NTRK3 (Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
E-EL-M0824	Mouse NT-4 (Neurotrophin 4) ELISA Kit
E-EL-M0825	Mouse NAP3 (Neutrophil Activating Protein 3) ELISA Kit
E-EL-M0826	Mouse NAP-2 (Neutrophil Activating Protein-2) ELISA Kit
E-EL-M0828	Mouse NGAL (Neutrophil Gelatinase Associated Lipocalin) ELISA Kit
E-EL-M0829	Mouse NOX1 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
E-EL-M0830	Mouse NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
E-EL-M0831	Mouse NADPH (Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
E-EL-M0834	Mouse NT-proBNP(N-terminal pro-Brain Natriuretic Peptide)ELISA Kit
E-EL-M0836	Mouse NFkB (Nuclear Factor Kappa B) ELISA Kit
E-EL-M0837	Mouse NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
E-EL-M0838	Mouse NF-κB p65 (Nuclear Factor-κB p65) ELISA Kit

Cat.NO	Product Name
E-EL-M0839	Mouse NRIP1 (Nuclear Receptor Interacting Protein 1) ELISA Kit
E-EL-M0840	Mouse NXF1 (Nuclear RNA Export Factor 1) ELISA Kit
E-EL-M0841	Mouse NUP35 (Nucleoporin 35kDa) ELISA Kit
E-EL-M0842	Mouse NUP37 (Nucleoporin 37kDa) ELISA Kit
E-EL-M0843	Mouse NUP50 (Nucleoporin 50kDa) ELISA Kit
E-EL-M0844	Mouse NUP85 (Nucleoporin 85kDa) ELISA Kit
E-EL-M0845	Mouse NUP98 (Nucleoporin 98kDa) ELISA Kit
E-EL-M0846	Mouse NUP107 (Nucleoporin 107kDa) ELISA Kit
E-EL-M0847	Mouse NUP133 (Nucleoporin 133kDa) ELISA Kit
E-EL-M0848	Mouse NUP153 (Nucleoporin 153kDa) ELISA Kit
E-EL-M0849	Mouse NUP155 (Nucleoporin 155kDa) ELISA Kit
E-EL-M0850	Mouse NUP160 (Nucleoporin 160kDa) ELISA Kit
E-EL-M0851	Mouse NUP188 (Nucleoporin 188kDa) ELISA Kit
E-EL-M0852	Mouse NUP205 (Nucleoporin 205kDa) ELISA Kit
E-EL-M0853	Mouse NUP214 (Nucleoporin 214kDa) ELISA Kit
E-EL-M0854	Mouse OBTF1 (Octamer Binding Transcription Factor 1) ELISA Kit
E-EL-M0855	Mouse OBTF2 (Octamer Binding Transcription Factor 2) ELISA Kit
E-EL-M0856	Mouse OLFM4 (Olfactomedin 4) ELISA Kit
E-EL-M0857	Mouse ITLN1 (Intelectin 1/Omentin) ELISA Kit
E-EL-M0858	Mouse OIT3 (Oncoprotein Induced Transcript 3) ELISA Kit

Cat.NO	Product Name
E-EL-M0859	Mouse OSM (Oncostatin-M) ELISA Kit
E-EL-M0860	Mouse OXA (Orexin A) ELISA Kit
E-EL-M0861	Mouse ODC (Ornithine Decarboxylase) ELISA Kit
E-EL-M0862	Mouse OFQ/N (Orphanin FQ/Nociceptin) ELISA Kit
E-EL-M0863	Mouse IκBα (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
E-EL-M0864	Mouse OC/BGP (Osteocalcin) ELISA Kit
E-EL-M0865	Mouse OSCAR (Osteoclast Associated Receptor) ELISA Kit
E-EL-M0866	Mouse ODF (Osteoclast Differentiation Factor) ELISA Kit
E-EL-M0867	Mouse ON (Osteonectin) ELISA Kit
E-EL-M0868	Mouse OPGL (Osteoprotegerin Ligand) ELISA Kit
E-EL-M0869	Mouse OVA sIgE (Ovalbumin Specific IgE) ELISA Kit
E-EL-M0870	Mouse CA125 (Carbohydrate Antigen 125) ELISA Kit
E-EL-M0871	Mouse pIκBα (Phospho-Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
E-EL-M0872	Mouse OLab (Oxidized Lowdensity Lipoprotein Antibody) ELISA Kit
E-EL-M0876	Mouse MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit
E-EL-M0878	Mouse PP (Pancreatic Polypeptide) ELISA Kit
E-EL-M0880	Mouse PTH (Parathyroid Hormone) ELISA Kit
E-EL-M0881	Mouse PTHR2 (Parathyroid Hormone Receptor 2) ELISA Kit
E-EL-M0882	Mouse PTHrP (Parathyroid Hormone Related Protein) ELISA Kit
E-EL-M0883	Mouse Pax (Paxillin) ELISA Kit

Cat.NO	Product Name
E-EL-M0884	Mouse PDAP1 (PDGFA Associated Protein 1) ELISA Kit
E-EL-M0885	Mouse PTD (Pentosidine) ELISA Kit
E-EL-M0886	Mouse PP (Pepsin) ELISA Kit
E-EL-M0887	Mouse PGA (Pepsinogen A) ELISA Kit
E-EL-M0888	Mouse PGC (Pepsinogen C) ELISA Kit
E-EL-M0889	Mouse PIN4 (Peptidyl Prolyl Cis/Trans Isomerase NIMA Interacting 4 Protein) ELISA Kit
E-EL-M0890	Mouse PRF1 (Perforin 1) ELISA Kit
E-EL-M0892	Mouse PMP3 (Peroxisomal Membrane Protein 3) ELISA Kit
E-EL-M0893	Mouse PPAR- γ (Peroxisome Proliferator-activated receptor γ) ELISA Kit
E-EL-M0894	Mouse PPAR- α (Peroxisome Proliferators-activator Receptors alpha) ELISA Kit
E-EL-M0895	Mouse PT-IgA (Pertussis virus IgA) ELISA Kit
E-EL-M0896	Mouse PI Ab-IgG/IgM (phosphatidylinositol antibody IgG/IgM) ELISA Kit
E-EL-M0897	Mouse PS (Phosphatidylserine) ELISA Kit
E-EL-M0898	Mouse P-PKC (Phospho Protein Kinase C) ELISA Kit
E-EL-M0899	Mouse pERK (Phospho-Extracellular Signal-Regulated Kinase) ELISA Kit
E-EL-M0900	Mouse PI3K2a (Phosphoinositide 3 Kinase 2a) ELISA Kit
E-EL-M0901	Mouse PLA1 (Phospholipase A1) ELISA Kit
E-EL-M0902	Mouse PLAA (Phospholipase A1 Activating Protein) ELISA Kit
E-EL-M0903	Mouse LDLR (Low Density Lipoprotein Receptor) ELISA Kit
E-EL-M0904	Mouse PLC γ 1 (Phospholipase C Gamma 1) ELISA Kit

Cat.NO	Product Name
E-EL-M0906	Mouse PI3K (Phosphotylinosital 3 kinase) ELISA Kit
E-EL-M0907	Mouse PEDF (Pigment Epithelium Derived Factor) ELISA Kit
E-EL-M0908	Mouse PIR (Pirin) ELISA Kit
E-EL-M0909	Mouse PACAPR1 (Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit
E-EL-M0910	Mouse PACAP (Pituitary Adenylate Cyclase Activating Polypeptide) ELISA Kit
E-EL-M0911	Mouse PGF (Placental Growth Factor) ELISA Kit
E-EL-M0912	Mouse PKP2 (Plakophilin 2) ELISA Kit
E-EL-M0913	Mouse PAP (Plasmin-Antiplasmin Complex) ELISA Kit
E-EL-M0914	Mouse PAI-1 (Plasminogen Activator Inhibitor 1) ELISA Kit
E-EL-M0915	Mouse PAI-2 (Plasminogen Activator Inhibitor 2) ELISA Kit
E-EL-M0917	Mouse tPA (Plasminogen Activator, Tissue) ELISA Kit
E-EL-M0918	Mouse PAF (Platelet Activating Factor) ELISA Kit
E-EL-M0920	Mouse PDGFC (Platelet Derived Growth Factor C) ELISA Kit
E-EL-M0921	Mouse PDGFD (Platelet Derived Growth Factor D) ELISA Kit
E-EL-M0922	Mouse PDGFRL (Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
E-EL-M0923	Mouse PDGFA (Platelet Derived Growth Factor Subunit A) ELISA Kit
E-EL-M0924	Mouse PDGFB (Platelet Derived Growth Factor Subunit B) ELISA Kit
E-EL-M0925	Mouse PF3 (Platelet Factor 3) ELISA Kit
E-EL-M0926	Mouse PF4 (Platelet Factor 4) ELISA Kit
E-EL-M0927	Mouse GP II b III a (Platelet Membrane Glycoprotein II b III a) ELISA Kit

Cat.NO	Product Name
E-EL-M0928	Mouse PECAM1 (Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit
E-EL-M0929	Mouse PDGFR- α (Platelet-Derived Growth Factor Soluble Receptor α) ELISA Kit
E-EL-M0930	Mouse PDGF (Platelet-Derived Growth Factor) ELISA Kit
E-EL-M0931	Mouse PHLDA2 (Pleckstrin Homology Like Domain Family A Member 2) ELISA Kit
E-EL-M0932	Mouse PTN (Pleiotrophin) ELISA Kit
E-EL-M0933	Mouse PCX (Podocalyxin) ELISA Kit
E-EL-M0934	Mouse PVRL1 (Poliovirus Receptor Related Protein 1) ELISA Kit
E-EL-M0935	Mouse PVRL2 (Poliovirus Receptor Related Protein 2) ELISA Kit
E-EL-M0936	Mouse PVRL3 (Poliovirus Receptor Related Protein 3) ELISA Kit
E-EL-M0937	Mouse PC4 (activated RNA polymerase II transcription cofactor 4) ELISA Kit
E-EL-M0940	Mouse PAPP-A (pregnancy associated plasma protein A) ELISA Kit
E-EL-M0941	Mouse Pro-ANP (Pro Atrial Natriuretic Peptide) ELISA Kit
E-EL-M0944	Mouse PGR (Progesterone receptor) ELISA Kit
E-EL-M0945	Mouse Pg (Progesterone) ELISA Kit
E-EL-M0946	Mouse PHB (Prohibitin) ELISA Kit
E-EL-M0947	Mouse PI (Proinsulin) ELISA Kit
E-EL-M0948	Mouse PK2 (Prokineticin 2) ELISA Kit
E-EL-M0949	Mouse PKR1 (Prokineticin Receptor 1) ELISA Kit
E-EL-M0950	Mouse GPX1 (Glutathione Peroxidase 1) ELISA Kit
E-EL-M0951	Mouse PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit

Cat.NO	Product Name
E-EL-M0952	Mouse POMC (Proopiomelanocortin) ELISA Kit
E-EL-M0953	Mouse PGD2S (Prostaglandin D2 Synthase) ELISA Kit
E-EL-M0954	Mouse EP2 (Prostaglandin E Receptor 2) ELISA Kit
E-EL-M0955	Mouse PTGES2 (Prostaglandin E Synthase 2) ELISA Kit
E-EL-M0956	Mouse PTGES (Prostaglandin E Synthase, Microsomal) ELISA Kit
E-EL-M0958	Mouse PTGS1 (Prostaglandin Endoperoxide Synthase 1) ELISA Kit
E-EL-M0959	Mouse PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit
E-EL-M0960	Mouse PG-F (Prostaglandin F) ELISA Kit
E-EL-M0961	Mouse PSA (Prostate Specific Antigen) ELISA Kit
E-EL-M0962	Mouse PAP (Prostatic Acid Phosphatase) ELISA Kit
E-EL-M0963	Mouse PAR2 (Protease Activated Receptor 2) ELISA Kit
E-EL-M0964	Mouse PRSS1 (Protease, Serine, 1) ELISA Kit
E-EL-M0965	Mouse PROC (Protein C) ELISA Kit
E-EL-M0967	Mouse PDIA3 (Protein Disulfide Isomerase A3) ELISA Kit
E-EL-M0968	Mouse PDI (Protein Disulfide Isomerase) ELISA Kit
E-EL-M0969	Mouse PIAS1 (Protein Inhibitor Of Activated STAT 1) ELISA Kit
E-EL-M0970	Mouse PICK1 (Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
E-EL-M0971	Mouse PKBg (Protein Kinase B Gamma) ELISA Kit
E-EL-M0972	Mouse PKB (Protein Kinase B) ELISA Kit
E-EL-M0973	Mouse PKCb1 (Protein Kinase C Beta 1) ELISA Kit

Cat.NO	Product Name
E-EL-M0974	Mouse PKCd (Protein Kinase C Delta) ELISA Kit
E-EL-M0975	Mouse PKCe (Protein Kinase C Epsilon) ELISA Kit
E-EL-M0976	Mouse PKCg (Protein Kinase C Gamma) ELISA Kit
E-EL-M0977	Mouse PKCi (Protein Kinase C Iota) ELISA Kit
E-EL-M0978	Mouse PKCq (Protein Kinase C Theta) ELISA Kit
E-EL-M0979	Mouse PKCz (Protein Kinase C Zeta) ELISA Kit
E-EL-M0980	Mouse PKD1 (Protein Kinase D1) ELISA Kit
E-EL-M0981	Mouse PKD2 (Protein Kinase D2) ELISA Kit
E-EL-M0982	Mouse PKIa (Protein Kinase Inhibitor Alpha) ELISA Kit
E-EL-M0983	Mouse PKIb (Protein Kinase Inhibitor Beta) ELISA Kit
E-EL-M0984	Mouse PKIg (Protein Kinase Inhibitor Gamma) ELISA Kit
E-EL-M0985	Mouse PKN1 (Protein Kinase N1) ELISA Kit
E-EL-M0986	Mouse PRKAa1 (Protein Kinase, AMP Activated Alpha 1) ELISA Kit
E-EL-M0987	Mouse PP (Protein Phosphatase) ELISA Kit
E-EL-M0988	Mouse PTPRN (Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
E-EL-M0989	Mouse PTPRO (Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
E-EL-M0990	Mouse PCDH1 (Protocadherin 1) ELISA Kit
E-EL-M0991	Mouse PARC (Pulmonary Activation Regulated Chemokine) ELISA Kit
E-EL-M0992	Mouse SPA (Pulmonary Surfactant Associated Protein A) ELISA Kit
E-EL-M0993	Mouse SPD (Pulmonary Surfactant Associated Protein D) ELISA Kit

Cat.NO	Product Name
E-EL-M0996	Mouse PDHa (Pyruvate Dehydrogenase Alpha) ELISA Kit
E-EL-M0997	Mouse PDK4 (Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit
E-EL-M0998	Mouse PDP (Pyruvate Dehydrogenase Phosphatase) ELISA Kit
E-EL-M0999	Mouse PDH E1 (Pyruvate Dehydrogenase-E1) ELISA Kit
E-EL-M1000	Mouse PK (Pyruvate Kinase, Liver And RBC) ELISA Kit
E-EL-M1001	Mouse M2-PK (Pyruvate Kinase) ELISA Kit
E-EL-M1002	Mouse PK (Pyruvate Kinase) ELISA Kit
E-EL-M1003	Mouse Qa-2 (Qa Lymphocyte Antigen 2 Region) ELISA Kit
E-EL-M1004	Mouse RAPGEF2 (Rap Guanine Nucleotide Exchange Factor 2) ELISA Kit
E-EL-M1005	Mouse RasGAP (Ras GTPase Activating Protein) ELISA Kit
E-EL-M1008	Mouse REG4 (Regenerating Islet Derived Protein 4) ELISA Kit
E-EL-M1009	Mouse RGS4 (Regulator Of G Protein Signaling 4) ELISA Kit
E-EL-M1010	Mouse RXFP1 (Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
E-EL-M1012	Mouse RSV-IgA (Respiratory syncytial virus IgA) ELISA Kit
E-EL-M1013	Mouse RSV-IgD (Respiratory syncytial virus IgD) ELISA Kit
E-EL-M1014	Mouse RSV-IgG (Respiratory syncytial virus IgG) ELISA Kit
E-EL-M1015	Mouse RSV-IgM (Respiratory syncytial virus IgM) ELISA Kit
E-EL-M1016	Mouse RBP4 (Retinol Binding Protein 4, Plasma) ELISA Kit
E-EL-M1019	Mouse RNASEA (Ribonuclease A) ELISA Kit
E-EL-M1020	Mouse RNASEH (Ribonuclease H) ELISA Kit

Cat.NO	Product Name
E-EL-M1021	Mouse RNASET2 (Ribonuclease T2) ELISA Kit
E-EL-M1022	Mouse RPL10 (Ribosomal Protein L10) ELISA Kit
E-EL-M1023	Mouse RPL26 (Ribosomal Protein L26) ELISA Kit
E-EL-M1024	Mouse RNF4 (Ring Finger Protein 4) ELISA Kit
E-EL-M1030	Mouse RUNX2 (Runt Related Transcription Factor 2) ELISA Kit
E-EL-M1031	Mouse S100A11 (S100 Calcium Binding Protein A11) ELISA Kit
E-EL-M1032	Mouse S100A7 (S100 Calcium Binding Protein A7) ELISA Kit
E-EL-M1033	Mouse S100B (S100 Calcium Binding Protein B) ELISA Kit
E-EL-M1035	Mouse S100A6 (S100 Calcium Binding Protein A6) ELISA Kit
E-EL-M1036	Mouse S100A12 (S100 Calcium Binding Protein A12) ELISA Kit
E-EL-M1037	Mouse SDH (Sorbitol Dehydrogenase) ELISA Kit
E-EL-M1038	Mouse sPLA2 (Secreted Phospholipase A2) ELISA Kit
E-EL-M1039	Mouse SCT (Secretin) ELISA Kit
E-EL-M1040	Mouse sIgA (Secretory Immunoglobulin A) ELISA Kit
E-EL-M1042	Mouse SEPP1 (Selenoprotein P) ELISA Kit
E-EL-M1043	Mouse SEMA7A (Semaphorin 7A) ELISA Kit
E-EL-M1045	Mouse STK (serine/threonine protein kinase) ELISA Kit
E-EL-M1046	Mouse SAA (Serum amyloid A) ELISA Kit
E-EL-M1047	Mouse SESN1 (Sestrin 1) ELISA Kit
E-EL-M1048	Mouse SESN2 (Sestrin 2) ELISA Kit

Cat.NO	Product Name
E-EL-M1049	Mouse SESN3 (Sestrin 3) ELISA Kit
E-EL-M1050	Mouse SOX3 (Sex Determining Region Y Box Protein 3) ELISA Kit
E-EL-M1051	Mouse SHBG (Sex Hormone-Binding Globulin) ELISA Kit
E-EL-M1052	Mouse SIGLEC3 (Sialic Acid Binding Ig Like Lectin 3) ELISA Kit
E-EL-M1053	Mouse SIGLEC12 (Sialic Acid Binding Ig Like Lectin 12) ELISA Kit
E-EL-M1054	Mouse SRP10 (Signal Recognition Particle 9kDa) ELISA Kit
E-EL-M1055	Mouse STAT1 (Signal Transducer And Activator Of Transcription 1) ELISA Kit
E-EL-M1056	Mouse STAT2 (Signal Transducer And Activator Of Transcription 2) ELISA Kit
E-EL-M1057	Mouse STAT3 (Signal Transducer And Activator Of Transcription 3) ELISA Kit
E-EL-M1058	Mouse STAT4 (Signal Transducer And Activator Of Transcription 4) ELISA Kit
E-EL-M1059	Mouse STAT5A (Signal Transducer And Activator Of Transcription 5A) ELISA Kit
E-EL-M1060	Mouse STAT6B (Signal Transducer And Activator Of Transcription 5B) ELISA Kit
E-EL-M1061	Mouse STAT6 (Signal Transducer And Activator Of Transcription 6) ELISA Kit
E-EL-M1062	Mouse Smad1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit
E-EL-M1063	Mouse OCT (Ornithine Carbamoyl Transferase) ELISA Kit
E-EL-M1064	Mouse ACTN2 (Actinin Alpha 2) ELISA Kit
E-EL-M1066	Mouse Smad7 (Mothers Against Decapentaplegic Homolog 7) ELISA Kit
E-EL-M1067	Mouse SNRPD1 (Small Nuclear Ribonucleoprotein D1) ELISA Kit
E-EL-M1071	Mouse sCD86 (Soluble Cluster of differentiation 86) ELISA Kit
E-EL-M1072	Mouse ENG (Endoglin) ELISA Kit

Cat.NO	Product Name
E-EL-M1073	Mouse EPCR (Endothelial Protein C Receptor) ELISA Kit
E-EL-M1074	Mouse IL-1sR I (Soluble Interleukin-1 Receptor I) ELISA Kit
E-EL-M1075	Mouse IL-1sR II (Soluble Interleukin-1 Receptor II) ELISA Kit
E-EL-M1076	Mouse IL-2sR α /CD25 (Soluble Interleukin-2 Receptor) ELISA Kit
E-EL-M1077	Mouse IL-2sR β /CD122 (Soluble Interleukin-2 Receptor) ELISA Kit
E-EL-M1078	Mouse sMHC-1 (soluble myosin heavy chain 1) ELISA Kit
E-EL-M1079	Mouse sMHC-2 (soluble myosin heavy chain 2) ELISA Kit
E-EL-M1081	Mouse sPECAM-1 (Soluble Platelet Endothelial Cell Adhesion Molecule 1) ELISA Kit
E-EL-M1082	Mouse S-100 (Soluble protein-100) ELISA Kit
E-EL-M1084	Mouse TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit
E-EL-M1085	Mouse SSTR1 (Somatostatin Receptor 1) ELISA Kit
E-EL-M1086	Mouse SST (Somatostatin) ELISA Kit
E-EL-M1087	Mouse SPOP (Speckle Type POZ Protein) ELISA Kit
E-EL-M1088	Mouse SSFA2 (Sperm Specific Antigen 2) ELISA Kit
E-EL-M1090	Mouse SCCA2 (Squamous Cell Carcinoma Antigen 2) ELISA Kit
E-EL-M1091	Mouse SCFR (Stem Cell Factor Receptor) ELISA Kit
E-EL-M1092	Mouse SF1 (Steroidogenic Factor 1) ELISA Kit
E-EL-M1093	Mouse SREBP (Sterol Regulatory Element Binding Protein) ELISA Kit
E-EL-M1094	Mouse SDF1 (Stromal Cell Derived Factor 1) ELISA Kit
E-EL-M1095	Mouse SDF-1 β (Stromal Cell Derived Factor 1 Beta) ELISA Kit

Cat.NO	Product Name
E-EL-M1096	Mouse SDF2 (Stromal Cell Derived Factor 2) ELISA Kit
E-EL-M1097	Mouse SDF2L1 (Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
E-EL-M1098	Mouse SDF4 (Stromal Cell Derived Factor 4) ELISA Kit
E-EL-M1099	Mouse SSRP1 (Structure Specific Recognition Protein 1) ELISA Kit
E-EL-M1102	Mouse SOCS3 (Suppressors Of Cytokine Signaling 3) ELISA Kit
E-EL-M1104	Mouse SMN1 (Survival Of Motor Neuron 1, Telomeric) ELISA Kit
E-EL-M1105	Mouse SYP (Synaptophysin) ELISA Kit
E-EL-M1106	Mouse SNAP-25 (Synaptosome Associated Protein 25) ELISA Kit
E-EL-M1107	Mouse SYT1 (Synaptotagmin 1) ELISA Kit
E-EL-M1108	Mouse SDC4 (Syndecan 4) ELISA Kit
E-EL-M1109	Mouse SNCa (Synuclein, Alpha) ELISA Kit
E-EL-M1111	Mouse TACR1 (Tachykinin Receptor 1) ELISA Kit
E-EL-M1112	Mouse TACR2 (Tachykinin Receptor 2) ELISA Kit
E-EL-M1113	Mouse TLN (Talin) ELISA Kit
E-EL-M1114	Mouse THP (Tamm–Horsfall Glycoprotein) ELISA Kit
E-EL-M1115	Mouse TDP43 (TAR DNA Binding Protein 43) ELISA Kit
E-EL-M1116	Mouse TRACP-5b (Tartrate-Resistant Acid Phosphatase 5b) ELISA Kit
E-EL-M1117	Mouse TBP (TATA Binding Protein) ELISA Kit
E-EL-M1118	Mouse TAF1 (TATA Box Binding Protein Associated Factor 1) ELISA Kit
E-EL-M1119	Mouse TAF2 (TATA Box Binding Protein Associated Factor 2) ELISA Kit

Cat.NO	Product Name
E-EL-M1120	Mouse TAF4 (TATA Box Binding Protein Associated Factor 4) ELISA Kit
E-EL-M1121	Mouse MAP τ (Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
E-EL-M1122	Mouse TXLNa (Taxilin Alpha) ELISA Kit
E-EL-M1123	Mouse TXLNb (Taxilin Beta) ELISA Kit
E-EL-M1124	Mouse TEK (Tyrosine Kinase, Endothelial) ELISA Kit
E-EL-M1125	Mouse TE (telomerase) ELISA Kit
E-EL-M1126	Mouse TNC (Tenascin C) ELISA Kit
E-EL-M1127	Mouse TNX (Tenascin X) ELISA Kit
E-EL-M1128	Mouse TDG1 (Teratocarcinoma Derived Growth Factor 1) ELISA Kit
E-EL-M1129	Mouse TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit
E-EL-M1130	Mouse TIEG1 (TGF-beta-inducible early response gene-1) ELISA Kit
E-EL-M1131	Mouse TBP2 (Thioredoxin Binding Protein 2) ELISA Kit
E-EL-M1133	Mouse TrxR1(Thioredoxin Reductase 1)ELISA Kit
E-EL-M1134	Mouse Trx (Thioredoxin) ELISA Kit
E-EL-M1135	Mouse TAFI (Thrombin activatable fibrinolysis inhibitor) ELISA Kit
E-EL-M1136	Mouse TR (Thrombin Receptor) ELISA Kit
E-EL-M1137	Mouse TSP-1 (Thrombospondin 1) ELISA Kit
E-EL-M1138	Mouse TAT (Thrombin-Antithrombin Complex) ELISA Kit
E-EL-M1140	Mouse TM (Thrombomodulin) ELISA Kit
E-EL-M1141	Mouse Thrombopoietin antibody ELISA Kit

Cat.NO	Product Name
E-EL-M1142	Mouse Thrombopoietin Receptor antibody ELISA Kit
E-EL-M1144	Mouse TXB2 (Thromboxane B2) ELISA Kit
E-EL-M1145	Mouse TXB4 (Thromboxanes B4) ELISA Kit
E-EL-M1146	Mouse TpP (Thrombus Precursor Protein) ELISA Kit
E-EL-M1147	Mouse RELN (Reelin) ELISA Kit
E-EL-M1148	Mouse TECK (Thymus Expressed Chemokine) ELISA Kit
E-EL-M1149	Mouse TG (Thyroglobulin) ELISA Kit
E-EL-M1150	Mouse TSAb (Thyroid Stimulating Hormone Antibody) ELISA Kit
E-EL-M1152	Mouse TRAb (Thyroid Stimulating Hormone Receptor Antidoby) ELISA Kit
E-EL-M1153	Mouse TSH (Thyroid Stimulating Hormone) ELISA Kit
E-EL-M1156	Mouse TEF (Thyrotroph Embryonic Factor) ELISA Kit
E-EL-M1157	Mouse TRH (Thyrotropin-releasing Hormone) ELISA Kit
E-EL-M1158	Mouse TAb (Thyroxine Antibody) ELISA Kit
E-EL-M1159	Mouse TBG (Thyroxine Binding Globulin) ELISA Kit
E-EL-M1161	Mouse TJP1 (Tight Junction Protein 1) ELISA Kit
E-EL-M1162	Mouse TFPI (Tissue Factor Pathway Inhibitor) ELISA Kit
E-EL-M1163	Mouse TF (Tissue factor) ELISA Kit
E-EL-M1164	Mouse TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
E-EL-M1165	Mouse TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
E-EL-M1166	Mouse TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) ELISA Kit

Cat.NO	Product Name
E-EL-M1167	Mouse TPA (Tissue Polypeptide Antigen) ELISA Kit
E-EL-M1168	Mouse t-PA (Tissue-type Plasiminogen Actilyse) ELISA Kit
E-EL-M1169	Mouse TLR9 (Toll-like Receptor 9) ELISA Kit
E-EL-M1171	Mouse FBXO32 (F-Box Protein 32) ELISA Kit
E-EL-M1172	Mouse CD8 (Cluster of Differentiation 8) ELISA Kit
E-EL-M1174	Mouse MuRF1 (muscle-specific RING-finger protein 1) ELISA Kit
E-EL-M1181	Mouse TCF20 (Transcription Factor 20) ELISA Kit
E-EL-M1183	Mouse RELB (V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
E-EL-M1184	Mouse TRF (Transferrin) ELISA Kit
E-EL-M1185	Mouse TFR2 (Transferrin Receptor Protein 2) ELISA Kit
E-EL-M1186	Mouse TFR (Transferrin Receptor) ELISA Kit
E-EL-M1187	Mouse TGFb1 (Transforming Growth Factor Beta Induced) ELISA Kit
E-EL-M1190	Mouse TGF- α (Transforming Growth Factor α) ELISA Kit
E-EL-M1191	Mouse TGFb2 (Transforming Growth Factor Beta 2) ELISA Kit
E-EL-M1192	Mouse TGF- β 3 (Transforming Growth Factor β 3) ELISA Kit
E-EL-M1193	Mouse TGF- β 4 (Transforming Growth Factor β 4) ELISA Kit
E-EL-M1195	Mouse TG (Transglutaminase) ELISA Kit
E-EL-M1196	Mouse TMED1 (Transmembrane Emp24 Domain Containing Protein 1) ELISA Kit
E-EL-M1197	Mouse TNPO1 (Transportin 1) ELISA Kit
E-EL-M1198	Mouse TTR (Transthyretin) ELISA Kit

Cat.NO	Product Name
E-EL-M1199	Mouse TFF2 (Trefoil Factor 2) ELISA Kit
E-EL-M1200	Mouse TFF3 (Trefoil Factor 3, Intestinal) ELISA Kit
E-EL-M1203	Mouse TNNI3/cTn-I (Troponin I Type 3, Cardiac) ELISA Kit
E-EL-M1204	Mouse TRY (Trypsin) ELISA Kit
E-EL-M1205	Mouse TAP (Trypsinogen Activation Peptide) ELISA Kit
E-EL-M1207	Mouse TUBa3C (Tubulin Alpha 3c) ELISA Kit
E-EL-M1208	Mouse TUBb (Tubulin Beta) ELISA Kit
E-EL-M1209	Mouse TPPP (Tubulin Polymerization Promoting Protein) ELISA Kit
E-EL-M1210	Mouse TNF- β (Tumor necrosis factor β) ELISA Kit
E-EL-M1211	Mouse TRAIL R1 (Tumor Necrosis Factor-related Apoptosis-inducing Ligand 1) ELISA Kit
E-EL-M1212	Mouse TRAIL R3 (Tumor Necrosis Factor-related Apoptosis-inducing Ligand 3) ELISA Kit
E-EL-M1213	Mouse TRAIL R4 (Tumor Necrosis Factor-related Apoptosis-inducing Ligand 4) ELISA Kit
E-EL-M1215	Mouse TP53 (Tumor Protein 53) ELISA Kit
E-EL-M1216	Mouse TSGF (Tumor Specific Growth Factor/tumor supplied group of factor) ELISA Kit
E-EL-M1217	Mouse TRY (Tyrosinase) ELISA Kit
E-EL-M1219	Mouse Ub (Ubiquitin) ELISA Kit
E-EL-M1221	Mouse U-T4 (Ultrasensitivity Thyroxine) ELISA Kit
E-EL-M1222	Mouse U-T3 (Ultrasensitivity Triiodothyronine) ELISA Kit
E-EL-M1226	Mouse PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit
E-EL-M1228	Mouse Cr (Creatinine) ELISA Kit

Cat.NO	Product Name
E-EL-M1229	Mouse VEGF-B (Vascular Endothelial cell Growth Factor B) ELISA Kit
E-EL-M1230	Mouse VEGF-C (Vascular Endothelial cell Growth Factor C) ELISA Kit
E-EL-M1231	Mouse VEGF-D (Vascular Endothelial cell Growth Factor D) ELISA Kit
E-EL-M1232	Mouse VEGFR3/Flt4 (Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit
E-EL-M1233	Mouse VCAM-1 (Vasculuar cell adhesion molecule 1) ELISA Kit
E-EL-M1234	Mouse VIP (Vasoactive Intestinal Peptide) ELISA Kit
E-EL-M1235	Mouse FOS (V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) ELISA Kit
E-EL-M1236	Mouse VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
E-EL-M1237	Mouse VF (Visfatin) ELISA Kit
E-EL-M1247	Mouse vWF (von Willebrand Factor) ELISA Kit
E-EL-M1249	Mouse α Manase (α mannosidase) ELISA Kit
E-EL-M1250	Mouse α 1-AGP (α 1-Acid glycoprotein) ELISA Kit
E-EL-M1251	Mouse α 1-MG(α 1-Microglobulin)ELISA Kit
E-EL-M1252	Mouse α 2-AP (α 2-Antiplasmin) ELISA Kit
E-EL-M1253	Mouse α 2-PI (α 2-plasmin inhibitor) ELISA Kit
E-EL-M1254	Mouse Gal (α -galactoyl) ELISA Kit
E-EL-M1255	Mouse α -Glu (Alpha Glucosidase) ELISA Kit
E-EL-M1256	Mouse α -GST (α -glutathione S-transferases) ELISA Kit
E-EL-M1258	Mouse IDUA (Alpha-L-iduronidase) ELISA Kit
E-EL-M1259	Mouse α NAG (α N-acetylglucosaminidase) ELISA Kit

Cat.NO	Product Name
E-EL-M1260	Mouse β Manase (β mannosidase) ELISA Kit
E-EL-M1261	Mouse β 2-GP1 IgA/G/M (β 2-glycoprotein 1 antibody IgA/G/M) ELISA Kit
E-EL-M1264	Mouse β GAL (β -Galactosidase) ELISA Kit
E-EL-M1265	Mouse β -Glu (Beta Glucosidase) ELISA Kit
E-EL-M1266	Mouse β GD (β -glucuronidase) ELISA Kit
E-EL-M1267	Mouse β -Hex A (β -hexosaminidase A) ELISA Kit
E-EL-M1269	Mouse β -TG (β -Thromboglobulin) ELISA Kit
E-EL-M1270	Mouse GLTP (Glycolipid Transfer Protein) ELISA Kit
E-EL-M1271	Mouse GDF6 (Growth Differentiation Factor 6) ELISA Kit
E-EL-M1272	Mouse ALT (Alanine Transaminase) ELISA Kit
E-EL-M1273	Mouse IL-21 (Interleukin 21) ELISA Kit
E-EL-M1274	Mouse PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit
E-EL-M1275	Mouse α HBDH (α -Hydroxybutyrate Dehydrogenase) ELISA Kit
E-EL-M1277	Mouse GLUT1 (Glucose Transporter 1) ELISA Kit
E-EL-M1278	Mouse 5-LO (Arachidonate 5-Lipoxygenase) ELISA Kit
E-EL-M1279	Mouse 12-LO (Arachidonate 12-Lipoxygenase) ELISA Kit
E-EL-M1280	Mouse 15-LO (Arachidonate 15-Lipoxygenase) ELISA Kit
E-EL-M1281	Mouse NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) ELISA Kit
E-EL-M1342	Mouse BFP (Brain Finger Protein) ELISA Kit
E-EL-M1343	Mouse S100A8 (S100 Calcium Binding Protein A8) ELISA Kit

Cat.NO	Product Name
E-EL-M1344	Mouse S100A9 (S100 Calcium Binding Protein A9) ELISA Kit
E-EL-M1345	Mouse S100A10 (S100 Calcium-binding Protein A10) ELISA Kit
E-EL-M1346	Mouse PLA2G2D (Phospholipase A2, Group IID) ELISA Kit
E-EL-M1347	Mouse cPLA2 (Phospholipase A2, Cytosolic) ELISA Kit
E-EL-M1348	Mouse iPLA2 (Phospholipase A2, Calcium Independent) ELISA Kit
E-EL-M1349	Mouse SAA2 (Serum Amyloid A2) ELISA Kit
E-EL-M1350	Mouse RTN4 (Reticulon 4) ELISA Kit
E-EL-M1352	Mouse ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-M1353	Mouse β MSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-M1354	Mouse Ta1 (Thymosin- α 1) ELISA Kit
E-EL-M1355	Mouse PS1 (Presenilin 1) ELISA Kit
E-EL-M1356	Mouse AgRP (Agouti Related Protein) ELISA Kit
E-EL-M1357	Mouse PACAP-38 (Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
E-EL-M1358	Mouse PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-M1360	Mouse PGF2 α (Prostaglandin F2 Alpha) ELISA Kit
E-EL-M1361	Mouse AMPK(Phosphorylated Adenosine Monophosphate Activated Protein Kinase)ELISA Kit
E-EL-M1362	Mouse G6PC (Glucose-6-Phosphatase, Catalytic) ELISA Kit
E-EL-M1363	Mouse LDL (Low Density Lipoprotein) ELISA Kit
E-EL-M1372	Mouse FASN (Fatty Acid Synthase) ELISA Kit
E-EL-M1402	Mouse HDL (High Density Lipoprotein) ELISA Kit

Cat.NO	Product Name
E-EL-M1801	Mouse cTnT/TNNT2 (Troponin T Type 2, Cardiac) ELISA Kit
E-EL-M2375	Mouse PYY (Peptide YY) ELISA Kit
E-EL-M2398	Mouse SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit
E-EL-M2399	Mouse IgG2b (Immunoglobulin G2b) ELISA Kit
E-EL-M2400	Mouse IgG2c (Immunoglobulin G2c) ELISA Kit
E-EL-M2401	Mouse IgG3 (Immunoglobulin G3) ELISA Kit
E-EL-M2402	Mouse FABP1 (Fatty Acid Binding Protein 1, Liver) ELISA Kit
E-EL-M2403	Mouse FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
E-EL-M2404	Mouse FABP4 (Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit
E-EL-M2405	Mouse α FP (Alpha-Fetoprotein) ELISA Kit
E-EL-M2406	Mouse CSTA (Cystatin A) ELISA Kit
E-EL-M2407	Mouse P-Cadherin (Cadherin, Placental) ELISA Kit
E-EL-M2408	Mouse PCDH15 (Protocadherin 15) ELISA Kit
E-EL-M2409	Mouse PDPN (Podoplanin) ELISA Kit
E-EL-M2410	Mouse NF κ B2 (Nuclear Factor Kappa B2) ELISA Kit
E-EL-M2411	Mouse BMG/ β 2-MG (Beta-2-Microglobulin) ELISA Kit
E-EL-M2412	Mouse NPHN (Nephrin) ELISA Kit
E-EL-M2413	Mouse AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit
E-EL-M2415	Mouse FGF23 (Fibroblast Growth Factor 23) ELISA Kit
E-EL-M2416	Mouse FGFR1 (Fibroblast Growth Factor Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-M2417	Mouse TLR4 (Toll-Like Receptor 4) ELISA Kit
E-EL-M2418	Mouse TLR-2 (Toll-like Receptor 2) ELISA Kit
E-EL-M2419	Mouse PCT (Procalcitonin) ELISA Kit
E-EL-M2420	Mouse TH (Tyrosine Hydroxylase) ELISA Kit
E-EL-M2421	Mouse ADM (Adrenomedullin) ELISA Kit
E-EL-M2422	Mouse APLN (Apelin) ELISA Kit
E-EL-M2423	Mouse CTSB (Cathepsin B) ELISA Kit
E-EL-M2424	Mouse SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-M2425	Mouse SOD3 (Superoxide Dismutase 3, Extracellular) ELISA Kit
E-EL-M2426	Mouse MMP-11 (Matrix Metalloproteinase 11) ELISA Kit
E-EL-M2427	Mouse EREG (Epiregulin) ELISA Kit
E-EL-M2428	Mouse α 1-AT (Alpha 1-Antitrypsin) ELISA Kit
E-EL-M2429	Mouse MMP-14 (Matrix Metalloproteinase 14) ELISA Kit
E-EL-M2430	Mouse CCK-8 (Cholecystokinin 8) ELISA Kit
E-EL-M2431	Mouse sCD14 (Soluble Cluster of Differentiation 14) ELISA Kit
E-EL-M2432	Mouse Fc ϵ RII/CD23 (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
E-EL-M2433	Mouse PCK1 (Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
E-EL-M2434	Mouse ACT α 2 (Actin Alpha 2, Smooth Muscle) ELISA Kit
E-EL-M2435	Mouse SOST (Sclerostin) ELISA Kit
E-EL-M2436	Mouse LCAT (Lecithin-Cholesterol Acyltransferase) ELISA Kit

Cat.NO	Product Name
E-EL-M2437	Mouse PIIINP (N-Terminal Procollagen III Propeptide) ELISA Kit
E-EL-M2438	Mouse NT-3 (Neurotrophin-3) ELISA Kit
E-EL-M2439	Mouse CD163 (Cluster of Differentiation 163) ELISA Kit
E-EL-M2440	Mouse DPP4 (Dipeptidyl Peptidase IV) ELISA Kit
E-EL-M2441	Mouse EGFR (Epidermal Growth Factor Receptor) ELISA Kit
E-EL-M2442	Mouse EG-VEGF (Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit
E-EL-M2443	Mouse IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit
E-EL-M2444	Mouse TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit
E-EL-M2445	Mouse MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit
E-EL-M2446	Mouse IL-22 (Interleukin 22) ELISA Kit
E-EL-M2447	Mouse ApoA1 (Apolipoprotein A1) ELISA Kit
E-EL-M2448	Mouse PL (Pancreatic Lipase) ELISA Kit
E-EL-M2449	Mouse IRF5 (Interferon Regulatory Factor 5)
E-EL-M2450	Mouse IL-12 p35 (Interleukin 12 p35) ELISA Kit
E-EL-M2451	Mouse IL-12/p40 (Interleukin 12 p40) ELISA Kit
E-EL-M2452	Mouse FTH (Ferritin, Heavy Polypeptide) ELISA Kit
E-EL-M2453	Mouse IL-6R (Interleukin 6 Receptor) ELISA Kit
E-EL-M2454	Mouse IL-19 (Interleukin 19) ELISA Kit
E-EL-M2455	Mouse IL-20 (Interleukin 20) ELISA Kit
E-EL-M2456	Mouse ALCAM (Activated Leukocyte Cell Adhesion Molecule) ELISA Kit

Cat.NO	Product Name
E-EL-M2457	Mouse IAPP (Islet Amyloid Polypeptide) ELISA Kit
E-EL-M2458	Mouse ANGPTL3 (Angiopoietin Like Protein 3) ELISA Kit
E-EL-M2459	Mouse FETUA (Fetuin A) ELISA Kit
E-EL-M2460	Mouse SDC1 (Syndecan 1) ELISA Kit
E-EL-M2461	Mouse IGFBP5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit
E-EL-M2462	Mouse ApoA2 (Apolipoprotein A2) ELISA Kit
E-EL-M2463	Mouse ApoC3 (Apolipoprotein C3) ELISA Kit
E-EL-M2567	Mouse CYP2E1 (Cytochrome P450, family 2, subfamily E, polypeptide 1) ELISA Kit
E-EL-M2600	Mouse MOG (Myelin Oligodendrocyte Glycoprotein) ELISA Kit
E-EL-M2601	Mouse AZGP1 (Zinc-alpha-2-glycoprotein) ELISA Kit
E-EL-M2602	Mouse LRP-1 (Low-Density Lipoprotein-Receptor-Related Protein) ELISA Kit
E-EL-M2603	Mouse TG (Triglyceride) ELISA Kit
E-EL-M2604	Mouse MUC1 (Mucin 1) ELISA Kit
E-EL-M2605	Mouse CAT (Catalase) ELISA Kit
E-EL-M2606	Mouse Klotho ELISA Kit
E-EL-M2607	Mouse Nrf2(Nuclear factor erythroid 2-related factor 2)ELISA Kit
E-EL-M2608	Mouse FCHOL(Free cholesterol)ELISA Kit
E-EL-M2609	Mouse IL8Ra (Interleukin 8 Receptor Alpha)ELISA Kit
E-EL-M2610	Mouse IL8Rb(Interleukin 8 Receptor Beta)ELISA Kit
E-EL-M2611	Mouse Ang- I (Angiotensin I) ELISA Kit

Cat.NO	Product Name
E-EL-M2612	Mouse Ang-II (Angiotensin II) ELISA Kit
E-EL-M2614	Mouse INS(Insulin) ELISA Kit
E-EL-M2615	Mouse UCHL1 (Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
E-EL-M2616	Mouse PRGN(progranulin)ELISA Kit
E-EL-M2617	Mouse HCY (Homocysteine) ELISA Kit
E-EL-M2618	Mouse MAOA(Amine oxidaseA)ELISA Kit
E-EL-M2619	Mouse MAOB(Amine oxidaseB)ELISA Kit
E-EL-M2620	Mouse Slc18a2(Solute Carrier Family 18,Member 2)ELISA Kit
E-EL-M2621	Mouse Slc18a1(Solute Carrier Family 18,Member 1)ELISA Kit
E-EL-M2622	Mouse S1PR1(Sphingosine 1 Phosphate Receptor 1)ELISA Kit
E-EL-M2623	Mouse S1PR2(Sphingosine 1 Phosphate Receptor 2)ELISA Kit
E-EL-M2624	Mouse S1PR3(Sphingosine 1 Phosphate Receptor 3)ELISA Kit
E-EL-M2625	Mouse S1PR4(Sphingosine 1 Phosphate Receptor 4)ELISA Kit
E-EL-M2626	Mouse S1PR5(Sphingosine 1 Phosphate Receptor 5)ELISA Kit
E-EL-M2627	Mouse CHGA (Chromogranin A) ELISA Kit
E-EL-M2628	Mouse PIBF(progesterone induced blocking factor)ELISA Kit
E-EL-M2629	Mouse DEFb14(beta-defensin 14)ELISA Kit
E-EL-M2630	Mouse A-CoA(acetyl-CoA)ELISA Kit
E-EL-M2631	Mouse AK1(Adenylate Kinase 1)ELISA Kit
E-EL-M2632	Mouse MPO(Myeloperoxidase)ELISA Kit

Cat.NO	Product Name
E-EL-M2633	Mouse GLP-2(Glucagon Like Peptide 2)ELISA Kit
E-EL-M2634	Mouse MLCK/MYLK(Myosin light chain kinase)ELISA Kit
E-EL-M2635	Mouse α MSH (Alpha-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-M2636	Mouse γ MSH (gamma-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-M2637	Mouse AChE (Acetylcholinesterase) ELISA Kit
E-EL-M2638	Mouse MT1 (Metallothionein 1) ELISA Kit
E-EL-M2639	Mouse AMT(Aminomethyltransferase)ELISA Kit
E-EL-M2640	Mouse RHOA (Ras Homolog Gene Family, Member A) ELISA Kit
E-EL-M2641	Mouse TFPI2(Tissue Factor Pathway Inhibitor 2) ELISA Kit
E-EL-M2642	Mouse IL-33 (Interleukin 33) ELISA Kit
E-EL-M2643	Mouse CYP1A2(Cytochrome P450 1A2)ELISA Kit
E-EL-M2644	Mouse FVII(Coagulation Factor VII)ELISA Kit
E-EL-M2646	Muose GDN(Glia Derived Nexin)ELISA Kit
E-EL-M2647	Mouse sIL-6R(Soluble Interleukin 6 Receptor)ELISA Kit
E-EL-M2648	Mouse sgp130(Soluble Glucoprotein 130)ELISA Kit
E-EL-M2649	Mouse CVL/EZR (Cytovillin/Ezrin) ELISA Kit
E-EL-M2650	Mouse sCD146(Soluble Cluster of Differentiation 146)ELISA Kit
E-EL-M2651	Mouse STUB1(E3 ubiquitin-protein ligase CHIP)ELISA Kit
E-EL-R0001	Rat ACV-A (Activin A) ELISA Kit
E-EL-R0002	Rat BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit

Cat.NO	Product Name
E-EL-R0003	Rat GRO α /CXCL1 (Growth Regulated Oncogene Alpha) ELISA Kit
E-EL-R0004	Rat ELN (Elastin) ELISA Kit
E-EL-R0005	Rat MIP-3 (Cytokine-induced neutrophil chemoattractant 3) ELISA Kit
E-EL-R0006	Rat Cyt-C (Cytochrome C) ELISA Kit
E-EL-R0007	Rat EPO (Erythropoietin) ELISA Kit
E-EL-R0008	Rat GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
E-EL-R0009	Rat IFN- γ (Interferon Gamma) ELISA Kit
E-EL-R0010	Rat IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit
E-EL-R0011	Rat IL-1 α (Interleukin 1 Alpha) ELISA Kit
E-EL-R0012	Rat IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-R0013	Rat IL-2 (Interleukin 2) ELISA Kit
E-EL-R0014	Rat IL-4 (Interleukin 4) ELISA Kit
E-EL-R0015	Rat IL-6 (Interleukin 6) ELISA Kit
E-EL-R0016	Rat IL-10 (Interleukin 10) ELISA Kit
E-EL-R0017	Rat ANP (Atrial Natriuretic Peptide) ELISA Kit
E-EL-R0019	Rat TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-R0020	Rat VEGF (Vascular Endothelial cell Growth Factor) ELISA Kit
E-EL-R0022	Rat CRP (C-Reactive Protein) ELISA Kit
E-EL-R0025	Rat MAU (Microalbuminuria) ELISA Kit
E-EL-R0026	Rat LH (Luteinizing Hormone) ELISA Kit

Cat.NO	Product Name
E-EL-R0029	Rat GH (Growth Hormone) ELISA Kit
E-EL-R0030	Rat REN (Renin) ELISA Kit
E-EL-R0031	Rat MT (Melatonin) ELISA Kit
E-EL-R0032	Rat C-P (C-peptide) ELISA Kit
E-EL-R0035	Rat LOX (Lysyl Oxidase) ELISA Kit
E-EL-R0037	Rat ApoB100 (Apolipoprotein B100) ELISA Kit
E-EL-R0039	Rat LN (Laminin) ELISA Kit
E-EL-R0041	Rat COL1 α 1 (Collagen Type I Alpha 1) ELISA Kit
E-EL-R0043	Rat MMP-10 (Matrix Metalloproteinase 10) ELISA Kit
E-EL-R0044	Rat MMP-12 (Matrix Metalloproteinase 12) ELISA Kit
E-EL-R0045	Rat MMP-13 (Matrix Metalloproteinase 13) ELISA Kit
E-EL-R0046	Rat sICAM-1/CD54 (Intercellular Adhesion Molecule-1) ELISA Kit
E-EL-R0047	Rat CT (Calcitonin) ELISA Kit
E-EL-R0048	Rat ACTH (Adrenocorticotrophic Hormone) ELISA Kit
E-EL-R0050	Rat OPG (Osteoprotegerin) ELISA Kit
E-EL-R0052	Rat PRL (Prolactin) ELISA Kit
E-EL-R0053	Rat MYO (Myoglobin) ELISA Kit
E-EL-R0054	Rat Tn-T (Troponin T) ELISA Kit
E-EL-R0055	Rat Tn-I (Troponin I) ELISA Kit
E-EL-R0056	Rat PAPPA (Pregnancy Associated Plasma Protein A) ELISA Kit

Cat.NO	Product Name
E-EL-R0057	Rat ALP (Alkaline Phosphatase) ELISA Kit
E-EL-R0058	Rat NSE (Neuron Specific Enolase) ELISA Kit
E-EL-R0059	Rat GLP-1 (Glucagon Like Peptide 1) ELISA Kit
E-EL-R0060	Rat SP-A (Pulmonary Surfactant-Associated Protein A) ELISA Kit
E-EL-R0061	Rat Smad7 (Mothers Against Decapentaplegic Homolog 7) ELISA Kit
E-EL-R0062	Rat NE (Neutrophil Elastase) ELISA Kit
E-EL-R0063	Rat IL-12/p40 (Interleukin 12 p40) ELISA Kit
E-EL-R0064	Rat IL-12 (Interleukin 12) ELISA Kit
E-EL-R0065	Rat MBP/MBL (Mannan Binding Protein/Mannan Binding Lectin) ELISA Kit
E-EL-R0066	Rat AQP-3 (Aquaporin 3) ELISA Kit
E-EL-R0067	Rat AQP-4 (Aquaporin 4) ELISA Kit
E-EL-R0070	Rat Arg (Arginase) ELISA Kit
E-EL-R0071	Rat AVPR1A (Arginine Vasopressin Receptor 1A) ELISA Kit
E-EL-R0072	Rat AhR (Aryl Hydrocarbon Receptor) ELISA Kit
E-EL-R0073	Rat AIP (Aryl Hydrocarbon Receptor Interacting Protein) ELISA Kit
E-EL-R0074	Rat ARNT2 (Aryl Hydrocarbon Receptor Nuclear Translocator 2) ELISA Kit
E-EL-R0075	Rat ASGR1 (Asialoglycoprotein Receptor 1) ELISA Kit
E-EL-R0076	Rat AST (Aspartate Aminotransferase) ELISA Kit
E-EL-R0077	Rat SP-B (Pulmonary Surfactant-Associated Protein B) ELISA Kit
E-EL-R0078	Rat TBX2 (T-Box Protein 2) ELISA Kit

Cat.NO	Product Name
E-EL-R0079	Rat ATM (Ataxia Telangiectasia Mutated) ELISA Kit
E-EL-R0080	Rat ATXN1 (Ataxin 1) ELISA Kit
E-EL-R0081	Rat ATXN2 (Ataxin 2) ELISA Kit
E-EL-R0082	Rat ABCA1 (ATP Binding Cassette Transporter A1) ELISA Kit
E-EL-R0083	Rat ABCG1 (ATP Binding Cassette Transporter G1) ELISA Kit
E-EL-R0084	Rat TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit
E-EL-R0085	Rat AXIN2 (Axis Inhibition Protein 2) ELISA Kit
E-EL-R0088	Rat FGF4 (Fibroblast Growth Factor 4) ELISA Kit
E-EL-R0089	Rat FGF6 (Fibroblast Growth Factor 6) ELISA Kit
E-EL-R0090	Rat FGF9 (Fibroblast Growth Factor 9) ELISA Kit
E-EL-R0091	Rat bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit
E-EL-R0092	Rat BAFFR (B-Cell Activation Factor Receptor) ELISA Kit
E-EL-R0093	Rat BCL10 (B-Cell Leukemia/Lymphoma 10) ELISA Kit
E-EL-R0094	Rat BCL3 (B-Cell Leukemia/Lymphoma 3) ELISA Kit
E-EL-R0096	Rat BCL-2 (B-Cell Leukemia/Lymphoma 2) ELISA Kit
E-EL-R0097	Rat BAK1 (BCL2 Antagonist/Killer 1) ELISA Kit
E-EL-R0098	Rat BAX (BCL-2 Associated X Protein) ELISA Kit
E-EL-R0099	Rat BMF (BCL-2 Modifying Factor) ELISA Kit
E-EL-R0100	Rat BAD (BCL-2 Associated Death Promoter) ELISA Kit
E-EL-R0101	Rat BTC (Betacellulin) ELISA Kit

Cat.NO	Product Name
E-EL-R0102	Rat BACE1 (Beta Site APP Cleaving Enzyme 1) ELISA Kit
E-EL-R0103	Rat BACE2 (Beta-Site APP-Cleaving Enzyme 2) ELISA Kit
E-EL-R0104	Rat β -EPR (Beta-Endorphin Receptor) ELISA Kit
E-EL-R0105	Rat β -EP (Beta-Endorphin) ELISA Kit
E-EL-R0107	Rat BID (BH3 Interacting Domain Death Agonist) ELISA Kit
E-EL-R0110	Rat BTA (Bladder Tumor Antigen) ELISA Kit
E-EL-R0112	Rat BLC (B-Lymphocyte Chemoattractant) ELISA Kit
E-EL-R0113	Rat BALP (Bone Alkaline Phosphatase) ELISA Kit
E-EL-R0114	Rat BMP-1 (Bone Morphogenetic Protein 1) ELISA Kit
E-EL-R0115	Rat BMP-3 (Bone Morphogenetic Protein 3) ELISA Kit
E-EL-R0116	Rat BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit
E-EL-R0117	Rat BMP-5 (Bone Morphogenetic Protein 5) ELISA Kit
E-EL-R0118	Rat BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit
E-EL-R0119	Rat BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit
E-EL-R0120	Rat BMPR1A (Bone Morphogenetic Protein Receptor, type IA) ELISA Kit
E-EL-R0121	Rat BMPR2 (Bone Morphogenetic Protein Receptor II) ELISA Kit
E-EL-R0124	Rat BK (Bradykinin) ELISA Kit
E-EL-R0125	Rat ANG2 (Angiopoietin 2) ELISA Kit
E-EL-R0126	Rat BNP (Brain Natriuretic Peptide) ELISA Kit
E-EL-R0128	Rat BRAK (Breast And Kidney Expressed Chemokine) ELISA Kit

Cat.NO	Product Name
E-EL-R0130	Rat VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit
E-EL-R0132	Rat NCAD (Neural Cadherin) ELISA Kit
E-EL-R0133	Rat CABIN1 (Calcineurin Binding Protein 1) ELISA Kit
E-EL-R0134	Rat CAN (Calcineurin) ELISA Kit
E-EL-R0135	Rat CGRP (Calcitonin Gene Related Peptide) ELISA Kit
E-EL-R0136	Rat β -Lg (Beta-Lactoglobulin) ELISA Kit
E-EL-R0137	Rat CASR (Calcium Sensing Receptor) ELISA Kit
E-EL-R0138	Rat CAMK2 (Calcium/Calmodulin-Dependent Protein Kinase- II) ELISA Kit
E-EL-R0139	Rat Alkaline Phosphatase, Intestinal (ALPI) ELISA Kit
E-EL-R0140	Rat CAM (Calmodulin) ELISA Kit
E-EL-R0141	Rat CNX (Calnexin) ELISA Kit
E-EL-R0142	Rat CRT (Calreticulin) ELISA Kit
E-EL-R0143	Rat CR (Calretinin) ELISA Kit
E-EL-R0144	Rat CALU (Calumenin) ELISA Kit
E-EL-R0145	Rat CDT (Carbohydrate-Deficient Transferrin) ELISA Kit
E-EL-R0146	Rat CA1 (Carbonic Anhydrase I) ELISA Kit
E-EL-R0147	Rat CA2 (Carbonic Anhydrase II) ELISA Kit
E-EL-R0149	Rat PICP (Carboxyterminal Propeptide of Type I Procollagen) ELISA Kit
E-EL-R0150	Rat CEA (Carcinoembryonic Antigen) ELISA Kit
E-EL-R0151	Rat cTnT/TNNT2 (Troponin T Type 2, Cardiac) ELISA Kit

Cat.NO	Product Name
E-EL-R0152	Rat RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
E-EL-R0153	Rat α FP(Alpha-Fetoprotein)ELISA Kit
E-EL-R0154	Rat CLS (Cardiolipin Synthase) ELISA Kit
E-EL-R0155	Rat CT1 (Cardiotrophin 1) ELISA Kit
E-EL-R0156	Rat CLCF1 (Cardiotrophin Like Cytokine Factor 1) ELISA Kit
E-EL-R0158	Rat CACT (Carnitine-Acylcarnitine Translocase) ELISA Kit
E-EL-R0159	Rat COMP (Cartilage Oligomeric Matrix Protein) ELISA Kit
E-EL-R0160	Rat CASP3 (Caspase 3) ELISA Kit
E-EL-R0161	Rat CASP4 (Caspase 4) ELISA Kit
E-EL-R0162	Rat CASP7 (Caspase 7) ELISA Kit
E-EL-R0163	Rat CASP9 (Caspase 9) ELISA Kit
E-EL-R0164	Rat CASP10 (Caspase 10) ELISA Kit
E-EL-R0165	Rat CASP12 (Caspase 12) ELISA Kit
E-EL-R0166	Rat CASP14 (Caspase 14) ELISA Kit
E-EL-R0167	Rat ET-1 (Endothelin 1) ELISA Kit
E-EL-R0169	Rat CTNNA1 (Catenin, Alpha 1) ELISA Kit
E-EL-R0170	Rat CTNNB (Catenin, Beta) ELISA Kit
E-EL-R0171	Rat CTSD (Cathepsin D) ELISA Kit
E-EL-R0172	Rat CTSK (Cathepsin K) ELISA Kit
E-EL-R0173	Rat CTSL (Cathepsin L) ELISA Kit

Cat.NO	Product Name
E-EL-R0174	Rat CDX2 (Caudal Type Homeobox Transcription Factor 2) ELISA Kit
E-EL-R0175	Rat CAV-1 (Caveolin-1) ELISA Kit
E-EL-R0176	Rat C/EBPδ (CCAAT/Enhancer Binding Protein data, C/EBPδ) ELISA Kit
E-EL-R0177	Rat CENPA (Centromere Protein A) ELISA Kit
E-EL-R0178	Rat CENPB (Centromere Protein B) ELISA Kit
E-EL-R0179	Rat CENPE (Centromere Protein E) ELISA Kit
E-EL-R0180	Rat CENP1 (Centromere Protein 1) ELISA Kit
E-EL-R0181	Rat CEP110 (Centrosomal Protein 110kDa) ELISA Kit
E-EL-R0182	Rat CP (Ceruloplasmin) ELISA Kit
E-EL-R0183	Rat V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog (FOS) ELISA Kit
E-EL-R0184	Rat CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit
E-EL-R0185	Rat CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit
E-EL-R0186	Rat CCRL1 (Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
E-EL-R0187	Rat SCD (Stearoyl Coenzyme A Desaturase) ELISA Kit
E-EL-R0188	Rat CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
E-EL-R0189	Rat CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit
E-EL-R0190	Rat CKLF (Chemokine Like Factor) ELISA Kit
E-EL-R0191	Rat CKLFSF4 (Chemokine Like Factor Superfamily 4) ELISA Kit
E-EL-R0192	Rat CKLFSF5 (Chemokine Like Factor Superfamily 5) ELISA Kit
E-EL-R0193	Rat CKLFSF6 (Chemokine Like Factor Superfamily 6) ELISA Kit

Cat.NO	Product Name
E-EL-R0194	Rat CKLFSF7 (Chemokine Like Factor Superfamily 7) ELISA Kit
E-EL-R0195	Rat CKLFSF10 (Chemokine Like Factor Superfamily 10) ELISA Kit
E-EL-R0196	Rat CCKAR (Cholecystokinin A Receptor) ELISA Kit
E-EL-R0198	Rat CETP (Cholest Erolester Transfer Protein) ELISA Kit
E-EL-R0199	Rat AGT (Angiotensinogen) ELISA Kit
E-EL-R0201	Rat CHRM2 (Cholinergic Receptor, Muscarinic 2) ELISA Kit
E-EL-R0203	Rat CG (Chorionic Gonadotrophin) ELISA Kit
E-EL-R0204	Rat CHD3 (Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
E-EL-R0205	Rat CHD5 (Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
E-EL-R0206	Rat CHGA (Chromogranin A) ELISA Kit
E-EL-R0207	Rat CNTF (Ciliary Neurotrophic Factor) ELISA Kit
E-EL-R0208	Rat CROP (Cisplatin Resistance Associated Overexpressed Protein) ELISA Kit
E-EL-R0209	Rat CS (Citrate Synthase) ELISA Kit
E-EL-R0210	Rat c-jun ELISA Kit
E-EL-R0211	Rat JNK (C-Jun N-Terminal Kinases) ELISA Kit
E-EL-R0212	Rat CC17 (Clara Cell Protein) ELISA Kit
E-EL-R0213	Rat CD109 (Cluster Of Differentiation 109) ELISA Kit
E-EL-R0214	Rat CD3 (Cluster of Differentiation 3) ELISA Kit
E-EL-R0215	Rat CD30L (Cluster Of Differentiation 30 Ligand) ELISA Kit
E-EL-R0216	Rat CD30 (Cluster of differentiation 30) ELISA Kit

Cat.NO	Product Name
E-EL-R0217	Rat CD40L (Cluster Of Differentiation 40 Ligand) ELISA Kit
E-EL-R0218	Rat CD6 (Cluster of Differentiation 6) ELISA Kit
E-EL-R0219	Rat CD8 (Cluster of Differentiation 8) ELISA Kit
E-EL-R0220	Rat CD97 (Cluster of Differentiation 97) ELISA Kit
E-EL-R0221	Rat CD44 (Cluster of Differentiation) ELISA Kit
E-EL-R0222	Rat CLU (Clusterin) ELISA Kit
E-EL-R0223	Rat F II (Coagulation Factor II)ELISA Kit
E-EL-R0224	Rat F III (Coagulation Factor III) ELISA Kit
E-EL-R0226	Rat F X (Coagulation Factor X) ELISA Kit
E-EL-R0227	Rat F X III (Coagulation Factor X III) ELISA Kit
E-EL-R0228	Rat F XII (Coagulation Factor XII) ELISA Kit
E-EL-R0229	Rat FIX (Coagulation Factor IX) ELISA Kit
E-EL-R0230	Rat FV (Coagulation Factor V) ELISA Kit
E-EL-R0231	Rat CFL1 (Cofilin 1, Non-Muscle) ELISA Kit
E-EL-R0232	Rat CFL2 (Cofilin 2, Muscle) ELISA Kit
E-EL-R0233	Rat COL1 (Collagen Type I) ELISA Kit
E-EL-R0234	Rat COL2 (Collagen Type II) ELISA Kit
E-EL-R0235	Rat COL3 (Collagen Type III) ELISA Kit
E-EL-R0236	Rat COL4 (Collagen Type IV) ELISA Kit
E-EL-R0237	Rat COL1α2 (Collagen Type I Alpha 2) ELISA Kit

Cat.NO	Product Name
E-EL-R0238	Rat COL4α1 (Collagen Type IV Alpha 1) ELISA Kit
E-EL-R0239	Rat COL9α1 (Collagen Type IX Alpha 1) ELISA Kit
E-EL-R0240	Rat COL10 (Collagen Type X) ELISA Kit
E-EL-R0241	Rat Collagenase I ELISA Kit
E-EL-R0242	Rat Collagenase II ELISA Kit
E-EL-R0243	Rat OC/BGP (Osteocalcin) ELISA Kit
E-EL-R0244	Rat BSP (Bone Sialoprotein) ELISA Kit
E-EL-R0245	Rat GCSFR (Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
E-EL-R0246	Rat C1INH (Complement 1 Inhibitor) ELISA Kit
E-EL-R0247	Rat TMSβ4 (Thymosin Beta 4) ELISA Kit
E-EL-R0248	Rat C4 (Complement Component 4) ELISA Kit
E-EL-R0249	Rat C1R (Complement Component 1Q Receptor) ELISA Kit
E-EL-R0250	Rat C3 (Complement Component 3) ELISA Kit
E-EL-R0251	Rat C5 (Complement Component 5) ELISA Kit
E-EL-R0252	Rat C9 (Complement Component 9) ELISA Kit
E-EL-R0253	Rat CFB (Complement Factor B) ELISA Kit
E-EL-R0254	Rat CFH (Complement Factor H) ELISA Kit
E-EL-R0255	Rat C3a (Complement Component 3a) ELISA Kit
E-EL-R0256	Rat C4a (Complement Component 4a) ELISA Kit
E-EL-R0257	Rat C5a (Complement Component 5a) ELISA Kit

Cat.NO	Product Name
E-EL-R0259	Rat CTGF (Connective Tissue Growth Factor) ELISA Kit
E-EL-R0260	Rat CX26 (Connexin 26) ELISA Kit
E-EL-R0261	Rat CX31 (Connexin 31) ELISA Kit
E-EL-R0262	Rat CX37 (Connexin 37) ELISA Kit
E-EL-R0263	Rat CX40 (Connexin 40) ELISA Kit
E-EL-R0264	Rat CX43 (Connexin 43) ELISA Kit
E-EL-R0265	Rat CNTN1 (Contactin 1) ELISA Kit
E-EL-R0266	Rat CNTN2 (Contactin 2) ELISA Kit
E-EL-R0267	Rat CNTN3 (Contactin 3) ELISA Kit
E-EL-R0268	Rat CNTN4(Contactin 4) ELISA Kit
E-EL-R0269	Rat CORT (Corticosterone) ELISA Kit
E-EL-R0270	Rat CRH (Corticotropin Releasing Hormone) ELISA Kit
E-EL-R0273	Rat CKMT1B (Creatine Kinase 1B, Mitochondrial) ELISA Kit
E-EL-R0274	Rat CKM (Creatine Kinase, Muscle) ELISA Kit
E-EL-R0275	Rat Cr (Creatinine) ELISA Kit
E-EL-R0276	Rat NTXI (Cross Linked N-Telopeptide of Type I Collagen) ELISA Kit
E-EL-R0277	Rat Anti-CIIAb (Anti-Type II Collagen Antibody) ELISA Kit
E-EL-R0279	Rat ICTP (Cross-linked C-Telopeptide of Type I Collagen) ELISA Kit
E-EL-R0280	Rat CASP8 (Caspase 8) ELISA Kit
E-EL-R0281	Rat CTX III (Cross-linked C-Telopeptides Of Type III Collagen) ELISA Kit

Cat.NO	Product Name
E-EL-R0282	Rat CRY1 (Cryptochrome 1) ELISA Kit
E-EL-R0283	Rat CLEC2C (C-Type Lectin Domain Family 2 Member C) ELISA Kit
E-EL-R0284	Rat CNP (C-Type Natriuretic Peptide) ELISA Kit
E-EL-R0285	Rat CTACK (Cutaneous T-Cell Attracting Chemokine) ELISA Kit
E-EL-R0287	Rat CXCR3 (CXC-Chemokine Receptor 3) ELISA Kit
E-EL-R0288	Rat TPS (Tryptase) ELISA Kit
E-EL-R0289	Rat CREB (Cyclic AMP Response Element Binding Protein) ELISA Kit
E-EL-R0291	Rat CDK5(Cyclin Dependent Kinase 5) ELISA Kit
E-EL-R0292	Rat CDK9 (Cyclin Dependent Kinase 9) ELISA Kit
E-EL-R0293	Rat CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
E-EL-R0294	Rat CCND1 (Cyclin D1) ELISA Kit
E-EL-R0295	Rat CCND2 (Cyclin D2) ELISA Kit
E-EL-R0296	Rat CCND3 (Cyclin D3) ELISA Kit
E-EL-R0298	Rat CYPA (Cyclophilin A) ELISA Kit
E-EL-R0299	Rat CYPB (Cyclophilin B) ELISA Kit
E-EL-R0300	Rat CYPD (Cyclophilin D) ELISA Kit
E-EL-R0303	Rat β -actin (Beta actin) ELISA Kit
E-EL-R0304	Rat Cys-C (Cystatin C) ELISA Kit
E-EL-R0305	Rat CYSLTR2 (Cysteinyl Leukotriene Receptor 2) ELISA Kit
E-EL-R0306	Rat CYB561 (Cytochrome b-561) ELISA Kit

Cat.NO	Product Name
E-EL-R0307	Rat FGF7 (Fibroblast Growth Factor 7) ELISA Kit
E-EL-R0308	Rat CYP2E1 (Cytochrome P450, family 2, subfamily E, polypeptide 1) ELISA Kit
E-EL-R0309	Rat CYP450 (Cytochrome P450) ELISA Kit
E-EL-R0310	Rat CDKN1A (Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit
E-EL-R0311	Rat CYTH1 (Cytohesin 1) ELISA Kit
E-EL-R0312	Rat CYTIP (Cytohesin Interacting Protein) ELISA Kit
E-EL-R0313	Rat Cpla2 (Cytosolic Phospholipase A2) ELISA Kit
E-EL-R0314	Rat CTLA4 (Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit
E-EL-R0315	Rat CVL (Cytovillin) ELISA Kit
E-EL-R0317	Rat D2D (D-Dimer) ELISA Kit
E-EL-R0318	Rat DC-SIGN/CD209 (DC Specific Intercellular Adhesion Molecule 3-Grabbing Nonintegrin) ELISA Kit
E-EL-R0319	Rat DAP6 (Death Associated Protein 6) ELISA Kit
E-EL-R0320	Rat DAF (Decay Accelerating Factor) ELISA Kit
E-EL-R0321	Rat DCN (Decorin) ELISA Kit
E-EL-R0322	Rat D-5 (Defensin-5) ELISA Kit
E-EL-R0323	Rat CA9 (Carbonic Anhydrase IX) ELISA Kit
E-EL-R0324	Rat DHEA-S7 (Dehydroepiandrosterone S7) ELISA Kit
E-EL-R0325	Rat DHEA-S(Dehydroepiandrosterone Sulfate)ELISA Kit
E-EL-R0326	Rat DHEA (Dehydroepiandrosterone) ELISA Kit
E-EL-R0327	Rat DPD (Deoxypyridinoline) ELISA Kit

Cat.NO	Product Name
E-EL-R0328	Rat DNASE1 (Deoxyribonuclease I) ELISA Kit
E-EL-R0329	Rat ADP/Acrp30 (Adiponectin) ELISA Kit
E-EL-R0330	Rat Des (Desmin) ELISA Kit
E-EL-R0331	Rat DAO (Diamine Oxidase) ELISA Kit
E-EL-R0332	Rat DKK1 (Dickkopf Related Protein 1) ELISA Kit
E-EL-R0333	Rat DKK2 (Dickkopf Related Protein 2) ELISA Kit
E-EL-R0334	Rat DLD (Dihydrolipoyl Dehydrogenase) ELISA Kit
E-EL-R0335	Rat DLAT (Dihydrolipoyl Transacetylase) ELISA Kit
E-EL-R0337	Rat DPP IV (Dipeptidyl Peptidase IV) ELISA Kit
E-EL-R0338	Rat D-LDH (D-Lactate Dehydrogenase) ELISA Kit
E-EL-R0339	Rat DRD1 (Dopamine Receptor D1) ELISA Kit
E-EL-R0340	Rat DRD2 (Dopamine Receptor D2) ELISA Kit
E-EL-R0341	Rat DRD3 (Dopamine Receptor D3) ELISA Kit
E-EL-R0342	Rat DRD4 (Dopamine Receptor D4) ELISA Kit
E-EL-R0343	Rat DAT (Dopamine Transporter) ELISA Kit
E-EL-R0344	Rat TM (Thrombin) ELISA Kit
E-EL-R0345	Rat DNM1 (Dynamin 1) ELISA Kit
E-EL-R0346	Rat DMD (Dystrophin) ELISA Kit
E-EL-R0347	Rat E-Cad (E-Cadherin) ELISA Kit
E-EL-R0350	Rat EGR1 (Early Growth Response Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-R0351	Rat EGR2 (Early Growth Response Protein 2) ELISA Kit
E-EL-R0352	Rat EGR3 (Early Growth Response Protein 3) ELISA Kit
E-EL-R0353	Rat GLUT1 (Glucose Transporter 1) ELISA Kit
E-EL-R0354	Rat GLUT2 (Glucose Transporter 2) ELISA Kit
E-EL-R0355	Rat AChE (Acetylcholinesterase) ELISA Kit
E-EL-R0356	Rat MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
E-EL-R0357	Rat ELA4 (Elastase 4) ELISA Kit
E-EL-R0358	Rat EMILIN 1 (Elastin Microfibril Interface Located Protein 1) ELISA Kit
E-EL-R0359	Rat ELOA (Elongin A) ELISA Kit
E-EL-R0360	Rat EED (Embryonic Ectoderm Development) ELISA Kit
E-EL-R0361	Rat EG-VEGF (Endocrine Gland Vascular Endothelial Growth Factor) ELISA Kit
E-EL-R0362	Rat ALB (Albumin) ELISA Kit
E-EL-R0363	Rat ERN1 (Endoplasmic Reticulum To Nucleus Signalling 1) ELISA Kit
E-EL-R0364	Rat ES (Endostatin) ELISA Kit
E-EL-R0365	Rat EL (Endothelial lipase) ELISA Kit
E-EL-R0367	Rat NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit
E-EL-R0368	Rat ECF/CCL11 (Eosinophil Chemotactic Factor) ELISA Kit
E-EL-R0369	Rat EGF (Epidermal Growth Factor) ELISA Kit
E-EL-R0370	Rat ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit
E-EL-R0371	Rat CASP1 (Caspase 1) ELISA Kit

Cat.NO	Product Name
E-EL-R0372	Rat ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit
E-EL-R0373	Rat FAS/CD95 (Factor Related Apoptosis) ELISA Kit
E-EL-R0375	Rat ACS (Fatty Acyl-CoA synthetase) ELISA Kit
E-EL-R0376	Rat Fga (Fibrinogen Alpha) ELISA Kit
E-EL-R0377	Rat FDP (Fibrinogen Degradation Product) ELISA Kit
E-EL-R0379	Rat FRS3 (Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
E-EL-R0380	Rat FGF-10 (Fibroblast Growth Factor-10) ELISA Kit
E-EL-R0381	Rat FLRT2 (Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit
E-EL-R0382	Rat FAK (Focal Adhesion Kinase) ELISA Kit
E-EL-R0384	Rat FST (Follistatin) ELISA Kit
E-EL-R0385	Rat CX3CL1 (Chemokine C-X3-C-Motif Ligand 1) ELISA Kit
E-EL-R0386	Rat CG β (Chorionic Gonadotropin Beta Polypeptide) ELISA Kit
E-EL-R0387	Rat ATGL (Adipose Triglyceride Lipase) ELISA Kit
E-EL-R0389	Rat F-TESTO (Free Testosterone) ELISA Kit
E-EL-R0390	Rat fT4 (Free Thyroxine) ELISA Kit
E-EL-R0391	Rat FSH (Follicle Stimulating Hormone) ELISA Kit
E-EL-R0392	Rat FUCA (Alpha-L-Fucosidase, Tissue) ELISA Kit
E-EL-R0393	Rat GNB1 (G Protein Beta 1) ELISA Kit
E-EL-R0394	Rat GL β (Galactosidase, Beta) ELISA Kit
E-EL-R0395	Rat FASN (Fatty Acid Synthase) ELISA Kit

Cat.NO	Product Name
E-EL-R0396	Rat GAL (Galanin) ELISA Kit
E-EL-R0397	Rat GAL1 (Galectin 1) ELISA Kit
E-EL-R0398	Rat GAL2 (Galectin 2) ELISA Kit
E-EL-R0399	Rat GAL3(Galectin 3)ELISA Kit
E-EL-R0400	Rat GAL4 (Galectin 4) ELISA Kit
E-EL-R0401	Rat GAL7 (Galectin 7) ELISA Kit
E-EL-R0402	Rat GAL8 (Galectin 8) ELISA Kit
E-EL-R0403	Rat GAL9 (Galectin 9) ELISA Kit
E-EL-R0404	Rat γGT1 (Gamma Glutamyltransferase 1) ELISA Kit
E-EL-R0406	Rat GJb1(Gap Junction Protein Beta 1) ELISA Kit
E-EL-R0407	Rat GIP (Gastric Inhibitory Polypeptide) ELISA Kit
E-EL-R0408	Rat GT-34 (Gastrin 34) ELISA Kit
E-EL-R0409	Rat CA199 (Gastrointestinalcancer Marker-CA199) ELISA Kit
E-EL-R0410	Rat GKN2 (Gastrokine 2) ELISA Kit
E-EL-R0411	Rat GATA1 (GATA Binding Protein 1) ELISA Kit
E-EL-R0412	Rat GATA2 (GATA Binding Protein 2) ELISA Kit
E-EL-R0413	Rat GATA3 (GATA Binding Protein 3) ELISA Kit
E-EL-R0414	Rat GATA4 (GATA Binding Protein 4) ELISA Kit
E-EL-R0415	Rat GATA5 (GATA Binding Protein 5) ELISA Kit
E-EL-R0416	Rat GS (Gelsolin) ELISA Kit

Cat.NO	Product Name
E-EL-R0417	Rat GTF II B (General Transcription Factor IIB) ELISA Kit
E-EL-R0418	Rat G-6-Pase/G6PT/GSD (Glucose 6 Phosphate) ELISA Kit
E-EL-R0419	Rat GMFB (Glia Maturation Factor, Beta) ELISA Kit
E-EL-R0420	Rat GDNF (Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
E-EL-R0421	Rat GTR (Glucocorticoid Induced Tumor Necrosis Factor Receptor) ELISA Kit
E-EL-R0422	Rat GR (Glucocorticoid Receptor) ELISA Kit
E-EL-R0423	Rat GR α (Glucocorticoid Receptor α) ELISA Kit
E-EL-R0424	Rat GR β (Glucocorticoid Receptor β) ELISA Kit
E-EL-R0425	Rat GC (Glucagon) ELISA Kit
E-EL-R0426	Rat GCK (Glucokinase) ELISA Kit
E-EL-R0427	Rat gp130 (Glucoprotein 130) ELISA Kit
E-EL-R0428	Rat G6PD (Glucose 6 Phosphate Dehydrogenase) ELISA Kit
E-EL-R0429	Rat GLUT3 (Glucose Transporter 3) ELISA Kit
E-EL-R0430	Rat GLUT4 (Glucose Transporter 4) ELISA Kit
E-EL-R0432	Rat GUSb (Glucuronidase, Beta) ELISA Kit
E-EL-R0433	Rat GCLM (Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit
E-EL-R0434	Rat GAD2 (Glutamate Decarboxylase 2, Acid) ELISA Kit
E-EL-R0435	Rat GDH (Glutamate Dehydrogenase) ELISA Kit
E-EL-R0436	Rat GRIK2 (Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit
E-EL-R0437	Rat CCK-8 (Cholecystokinin 8) ELISA Kit

Cat.NO	Product Name
E-EL-R0438	Rat GSTA1 (Glutathione S Transferase Alpha 1) ELISA Kit
E-EL-R0439	Rat GSTO1 (Glutathione S Transferase Omega 1) ELISA Kit
E-EL-R0440	Rat GSTT1 (Glutathione S Transferase Theta 1) ELISA Kit
E-EL-R0441	Rat GSTT2 (Glutathione S Transferase Theta 2) ELISA Kit
E-EL-R0442	Rat HbA1C (Glycosylated Hemoglobin/Hemoglobin A1c) ELISA Kit
E-EL-R0445	Rat PYGB (Glycogen Phosphorylase, Brain) ELISA Kit
E-EL-R0446	Rat PYGM (Glycogen Phosphorylase, Muscle) ELISA Kit
E-EL-R0447	Rat GSK3A (Glycogen Synthase Kinase 3 Alpha) ELISA Kit
E-EL-R0448	Rat GSK (Glycogen Synthase Kinase) ELISA Kit
E-EL-R0449	Rat GP5 (Glycoprotein V, Platelet) ELISA Kit
E-EL-R0450	Rat GnRHR (Gonadotropin-Releasing Hormone Receptor) ELISA Kit
E-EL-R0451	Rat GnRH (Gonadotropin-Releasing Hormone) ELISA Kit
E-EL-R0452	Rat GRN (Granulin) ELISA Kit
E-EL-R0453	Rat GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit
E-EL-R0454	Rat G-CSF (Granulocyte Colony Stimulating Factor) ELISA Kit
E-EL-R0455	Rat Gzms-A (Granzymes A) ELISA Kit
E-EL-R0456	Rat Gzms-B (Granzymes B) ELISA Kit
E-EL-R0457	Rat GREM1 (GremLin 1) ELISA Kit
E-EL-R0458	Rat GDF1 (Growth Differentiation Factor 1) ELISA Kit
E-EL-R0459	Rat GDF2 (Growth Differentiation Factor 2) ELISA Kit

Cat.NO	Product Name
E-EL-R0460	Rat GDF3 (Growth Differentiation Factor 3) ELISA Kit
E-EL-R0461	Rat GDF5 (Growth Differentiation Factor 5) ELISA Kit
E-EL-R0462	Rat GDF9 (Growth Differentiation Factor 9) ELISA Kit
E-EL-R0463	Rat GDF11(Growth Differentiation Factor 11)ELISA Kit
E-EL-R0464	Rat GDF15 (Growth Differentiation Factor 15) ELISA Kit
E-EL-R0465	Rat GRB10 (Growth Factor Receptor Bound Protein 10) ELISA Kit
E-EL-R0466	Rat GRB14 (Growth Factor Receptor Bound Protein 14) ELISA Kit
E-EL-R0467	Rat GRB2 (Growth Factor Receptor Bound Protein 2) ELISA Kit
E-EL-R0468	Rat GH2 (Growth Hormone 2) ELISA Kit
E-EL-R0469	Rat GHRH (Growth Hormone Releasing Hormone) ELISA Kit
E-EL-R0470	Rat TG (Triglyceride) ELISA Kit
E-EL-R0472	Rat GT (Gastrin) ELISA Kit
E-EL-R0473	Rat HP (Haptoglobin) ELISA Kit
E-EL-R0474	Rat HSPA4 (Heat Shock 70kDa Protein 4) ELISA Kit
E-EL-R0475	Rat HSF1 (Heat Shock Factor 1) ELISA Kit
E-EL-R0476	Rat HSP-27 (Heat Shock Protein 27) ELISA Kit
E-EL-R0477	Rat HSP-40 (Heat Shock Protein 40) ELISA Kit
E-EL-R0478	Rat HSP-60 (Heat Shock Protein 60) ELISA Kit
E-EL-R0479	Rat HSP-70 (Heat Shock Protein 70) ELISA Kit
E-EL-R0480	Rat HSP-90 (Heat Shock Protein 90) ELISA Kit

Cat.NO	Product Name
E-EL-R0481	Rat HSP GP96 (Heat Shock Protein Glycoprotein 96) ELISA Kit
E-EL-R0482	Rat HSP47 (Heat Shock Protein) ELISA Kit
E-EL-R0483	Rat HSF2 (Heat Shock Transcription Factor 2) ELISA Kit
E-EL-R0484	Rat SHH (Hedgehog Homolog, Sonic) ELISA Kit
E-EL-R0485	Rat HHIP (Hedgehog Interacting Protein) ELISA Kit
E-EL-R0486	Rat HLTF (Helicase Like Transcription Factor) ELISA Kit
E-EL-R0487	Rat Hp-IgG (Helicobacter Pylori IgG) ELISA Kit
E-EL-R0488	Rat HO-1 (Heme Oxygenase 1) ELISA Kit
E-EL-R0489	Rat HO-2 (Heme Oxygenase 2) ELISA Kit
E-EL-R0490	Rat HPX (Hemopexin) ELISA Kit
E-EL-R0491	Rat HSPG (Heparan Sulphate Proteoglycan) ELISA Kit
E-EL-R0492	Rat HPA (Heparanase) ELISA Kit
E-EL-R0493	Rat HBEGF (Heparin Binding Epidermal Growth Factor Like Growth Factor) ELISA Kit
E-EL-R0494	Rat HC II (Heparin Cofactor II) ELISA Kit
E-EL-R0495	Rat LIPC (Lipase, Hepatic) ELISA Kit
E-EL-R0496	Rat HGF (Hepatocyte Growth Factor) ELISA Kit
E-EL-R0497	Rat HNF1a (Hepatocyte Nuclear Factor 1 Alpha) ELISA Kit
E-EL-R0498	Rat HNF4a (Hepatocyte Nuclear Factor 4 Alpha) ELISA Kit
E-EL-R0499	Rat HDGF (Hepatoma Derived Growth Factor) ELISA Kit
E-EL-R0500	Rat Hpc (Hepcidin) ELISA Kit

Cat.NO	Product Name
E-EL-R0501	Rat HNRPA1 (Heterogeneous Nuclear Ribonucleoprotein A1) ELISA Kit
E-EL-R0502	Rat HK1 (Hexokinase 1) ELISA Kit
E-EL-R0503	Rat HK (Hexokinase) ELISA Kit
E-EL-R0504	Rat HDL (High Density Lipoprotein) ELISA Kit
E-EL-R0505	Rat HMGB-1 (High Mobility Group Protein B1) ELISA Kit
E-EL-R0506	Rat hs-CRP (High Sensitivity C-Reactive Protein) ELISA Kit
E-EL-R0507	Rat HAT1 (Histone Acetyltransferase 1) ELISA Kit
E-EL-R0508	Rat Histon-H2b ELISA Kit
E-EL-R0509	Rat COL1 (Collagen Type I) ELISA Kit
E-EL-R0510	Rat HABP1 (Hyaluronan Binding Protein 1) ELISA Kit
E-EL-R0513	Rat HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit
E-EL-R0514	Rat IDS (Iduronate Sulfatase) ELISA Kit
E-EL-R0515	Rat HMG-CoAR(3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase)ELISA Kit
E-EL-R0516	Rat IgA (Immunoglobulin A) ELISA Kit
E-EL-R0517	Rat IgE (Immunoglobulin E) ELISA Kit
E-EL-R0518	Rat IgG (Immunoglobulin G) ELISA Kit
E-EL-R0519	Rat IgM (Immunoglobulin M) ELISA Kit
E-EL-R0520	Rat NOS2/iNOS(Nitric Oxide Synthase 2, Inducible)ELISA Kit
E-EL-R0521	Rat INHA (Inhibin A) ELISA Kit
E-EL-R0522	Rat ADH (Antidiuretic Hormone) ELISA Kit

Cat.NO	Product Name
E-EL-R0523	Rat INHbC (Inhibin Beta C) ELISA Kit
E-EL-R0524	Rat INHBP (Inhibin Binding Protein) ELISA Kit
E-EL-R0525	Rat IKB α (Inhibitory Subunit of NF- κ B α)ELISA Kit
E-EL-R0526	Rat IGFBP5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit
E-EL-R0527	Rat IGFBP6 (Insulin Like Growth Factor Binding Protein 6) ELISA Kit
E-EL-R0528	Rat ISR- β (Insulin Receptor Beta Subunit) ELISA Kit
E-EL-R0529	Rat IL-35 (Interleukin 35) ELISA Kit
E-EL-R0530	Rat IGF-2 (Insulin Like Growth Factor 2) ELISA Kit
E-EL-R0531	Rat IGFBP-1 (Insulin-Like Growth Factor Binding Protein 1) ELISA Kit
E-EL-R0532	Rat IGFBP-2 (Insulin-Like Growth Factor Binding Protein 2) ELISA Kit
E-EL-R0533	Rat IGFBP-3 (Insulin-Like Growth Factor Binding Protein 3) ELISA Kit
E-EL-R0534	Rat IGFBP-4 (Insulin-Like Growth Factor Binding Protein 4) ELISA Kit
E-EL-R0535	Rat I-PTH (Intact Parathormone) ELISA Kit
E-EL-R0536	Rat PDGF-AB (Platelet-Derived Growth Factor AB) ELISA Kit
E-EL-R0537	Rat PDGF-BB (Platelet-Derived Growth Factor-BB) ELISA Kit
E-EL-R0540	Rat TIMP-1(Tissue Inhibitors of Metalloproteinase 1)ELISA Kit
E-EL-R0541	Rat α -GST(α -Glutathione S-Transferases) ELISA Kit
E-EL-R0542	Rat ICAM-2 (Intercellular Adhesion Molecule 2) ELISA Kit
E-EL-R0543	Rat ICAM-3 (Intercellular Adhesion Molecule 3) ELISA Kit
E-EL-R0544	Rat IFN- α (Interferon Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-R0545	Rat IFN- β (Interferon Beta) ELISA Kit
E-EL-R0546	Rat IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
E-EL-R0547	Rat IFI16/p16 (Interferon-Inducible Protein 16) ELISA Kit
E-EL-R0548	Rat INF- α -Ab (Interferon- α Antibody) ELISA Kit
E-EL-R0549	Rat PEPCK (Phosphoenolpyruvate Carboxykinase) ELISA Kit
E-EL-R0550	Rat FIL1d (Interleukin 1Delta) ELISA Kit
E-EL-R0551	Rat FIL1h (Interleukin 1 Eta) ELISA Kit
E-EL-R0552	Rat IL-1ra (Interleukin 1 Receptor Antagonist) ELISA Kit
E-EL-R0553	Rat LCAT (Lecithin-Cholesterol Acyltransferase) ELISA Kit
E-EL-R0554	Rat ALD (Aldosterone) ELISA Kit
E-EL-R0555	Rat IL-2R α (Interleukin-2 Receptor- α) ELISA Kit
E-EL-R0556	Rat IL-3 (Interleukin 3) ELISA Kit
E-EL-R0557	Rat FE (Ferritin) ELISA Kit
E-EL-R0558	Rat IL-5 (Interleukin 5) ELISA Kit
E-EL-R0559	Rat IL-7 (Interleukin 7) ELISA Kit
E-EL-R0561	Rat IL-9 (Interleukin 9) ELISA Kit
E-EL-R0562	Rat IL-11 (Interleukin 11) ELISA Kit
E-EL-R0563	Rat IL-13 (Interleukin 13) ELISA Kit
E-EL-R0564	Rat IL-15 (Interleukin 15) ELISA Kit
E-EL-R0565	Rat IL-16 (Interleukin 16) ELISA Kit

Cat.NO	Product Name
E-EL-R0566	Rat IL-17 (Interleukin 17) ELISA Kit
E-EL-R0567	Rat IL-18 (Interleukin 18) ELISA Kit
E-EL-R0568	Rat IL-21 (Interleukin 21) ELISA Kit
E-EL-R0569	Rat IL-23 (Interleukin 23) ELISA Kit
E-EL-R0570	Rat IL-27 (Interleukin 27) ELISA Kit
E-EL-R0571	Rat IL-32 (Interleukin 32) ELISA Kit
E-EL-R0572	Rat IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) ELISA Kit
E-EL-R0573	Rat MUC2 (Mucin 2) ELISA Kit
E-EL-R0574	Rat AR (Amphiregulin) ELISA Kit
E-EL-R0575	Rat KIM-1 (Kidney Injury Molecule 1) ELISA Kit
E-EL-R0577	Rat LTF/LF (Lactoferrin) ELISA Kit
E-EL-R0578	Rat FN (Fibronectin) ELISA Kit
E-EL-R0579	Rat LDL (Low Density Lipoprotein) ELISA Kit
E-EL-R0580	Rat LRP-1 (Low-Density Lipoprotein-Receptor-Related Protein) ELISA Kit
E-EL-R0581	Rat LDL-IC (Low Density Lipoprotein Immunity Complex) ELISA Kit
E-EL-R0582	Rat LEP (Leptin) ELISA Kit
E-EL-R0583	Rat Lep IgG (Leptospira IgG) ELISA Kit
E-EL-R0584	Rat LIFR (Leukemia Inhibitory Factor Receptor) ELISA Kit
E-EL-R0585	Rat LIF (Leukemia Inhibitory Factor) ELISA Kit
E-EL-R0587	Rat LT-C4 (Leukotriene C4) ELISA Kit

Cat.NO	Product Name
E-EL-R0588	Rat LT-E4 (Leukotriene E4) ELISA Kit
E-EL-R0589	Rat LBP (Lipolysaccharide Binding Protein) ELISA Kit
E-EL-R0590	Rat LPL (Lipoprotein Lipase) ELISA Kit
E-EL-R0591	Rat Lp-a (Lipoprotein a) ELISA Kit
E-EL-R0592	Rat Lp-PL-A2 (Lipoprotein-Associated Phospholipase A2) ELISA Kit
E-EL-R0594	Rat LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit
E-EL-R0596	Rat SELL (L-Selectin) ELISA Kit
E-EL-R0597	Rat LHRH (Luteinizing Hormone-Releasing Hormone) ELISA Kit
E-EL-R0599	Rat LFA-3 (Lymphocyte Function Associated Antigen 3) ELISA Kit
E-EL-R0600	Rat Lptn/XCL1 (Lymphotactin) ELISA Kit
E-EL-R0601	Rat M-CSF (Macrophage Colony-Stimulating Factor) ELISA Kit
E-EL-R0602	Rat MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
E-EL-R0603	Rat MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit
E-EL-R0604	Rat PYGL (Glycogen Phosphorylase, Liver) ELISA Kit
E-EL-R0605	Rat MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
E-EL-R0606	Rat MIP-3 β (Macrophage Inflammatory Protein 3 β) ELISA Kit
E-EL-R0607	Rat MIP-5 (Macrophage Inflammatory Protein 5) ELISA Kit
E-EL-R0608	Rat MIF (Macrophage Migration Inhibitory Factor) ELISA Kit
E-EL-R0609	Rat MDC/CCL22 (Macrophage-Derived Chemokine) ELISA Kit
E-EL-R0611	Rat MHC I /AgB I /H-1 I (Major Histocompatibility Complex I) ELISA Kit

Cat.NO	Product Name
E-EL-R0612	Rat MHC II /AgB II /H-1 II (Major Histocompatibility Complex II) ELISA Kit
E-EL-R0613	Rat MHC III /AgB III /H-1 III (Major Histocompatibility Complex III) ELISA Kit
E-EL-R0614	Rat RETN (Resistin) ELISA Kit
E-EL-R0615	Rat CA15-3 (Carbohydrate Antigen 15-3) ELISA Kit
E-EL-R0616	Rat MAT (Matrilysin) ELISA Kit
E-EL-R0617	Rat MMP-1 (Matrix Metalloproteinase 1) ELISA Kit
E-EL-R0618	Rat MMP-2(Matrix Metalloproteinase 2)ELISA Kit
E-EL-R0619	Rat MMP-3 (Matrix Metalloproteinase 3) ELISA Kit
E-EL-R0620	Rat MMP-17 (Matrix Metalloproteinase 17) ELISA Kit
E-EL-R0621	Rat MMP-24 (Matrix Metalloproteinase 24) ELISA Kit
E-EL-R0622	Rat MMP-7 (Matrix Metalloproteinase 7) ELISA Kit
E-EL-R0623	Rat MMP-8 (Matrix Metalloproteinase 8) ELISA Kit
E-EL-R0624	Rat MMP-9(Matrix Metalloproteinase 9)ELISA Kit
E-EL-R0625	Rat ANG (Angiogenin) ELISA Kit
E-EL-R0626	Rat ANG1 (Angiopoietin 1) ELISA Kit
E-EL-R0627	Rat SLC (Secondary Lymphoid Tissue Chemokine) ELISA Kit
E-EL-R0629	Rat MT1(Metallothionein 1)ELISA Kit
E-EL-R0630	Rat MHB (Methemoglobin) ELISA Kit
E-EL-R0631	Rat Methylase ELISA Kit
E-EL-R0633	Rat MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-R0634	Rat MCP-2 (Monocyte Chemotactic Protein 2) ELISA Kit
E-EL-R0635	Rat MCP-3 (Monocyte Chemotactic Protein 3) ELISA Kit
E-EL-R0636	Rat MCP-4 (Monocyte Chemotactic Protein 4) ELISA Kit
E-EL-R0637	Rat Smad1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit
E-EL-R0638	Rat Smad4 (Mothers Against Decapentaplegic Homolog 4) ELISA Kit
E-EL-R0639	Rat MTL (Motilin) ELISA Kit
E-EL-R0640	Rat AMH (Anti-Mullerian Hormone) ELISA Kit
E-EL-R0641	Rat M-ACHR (Muscarinic Acetylcholine Receptor) ELISA Kit
E-EL-R0642	Rat MBP (Myelin Basic Protein) ELISA Kit
E-EL-R0643	Rat AGER (Advanced Glycosylation End Product Specific Receptor) ELISA Kit
E-EL-R0644	Rat OXB (Orexin B) ELISA Kit
E-EL-R0645	Rat PTEN (Phosphatase and Tensin Homolog) ELISA Kit
E-EL-R0646	Rat MYS (Myosin) ELISA Kit
E-EL-R0647	Rat NAGase (N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
E-EL-R0648	Rat GAS (N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
E-EL-R0651	Rat NN-T4 (Neonatal Thyroxine) ELISA Kit
E-EL-R0652	Rat NGF (Nerve Growth Factor) ELISA Kit
E-EL-R0653	Rat NCAM/CD56 (Neural Cell Adhesion Molecule) ELISA Kit
E-EL-R0654	Rat NGB (Neuroglobin) ELISA Kit
E-EL-R0655	Rat NPY (Neuropeptide Y) ELISA Kit

Cat.NO	Product Name
E-EL-R0656	Rat NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
E-EL-R0657	Rat NTRK3 (Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
E-EL-R0658	Rat NT-3 (Neurotrophin 3) ELISA Kit
E-EL-R0659	Rat NT-4 (Neurotrophin 4) ELISA Kit
E-EL-R0660	Rat NAP3 (Neutrophil Activating Protein 3) ELISA Kit
E-EL-R0661	Rat NAP-2 (Neutrophil Activating Protein-2) ELISA Kit
E-EL-R0662	Rat NGAL (Neutrophil Gelatinase Associated Lipocalin) ELISA Kit
E-EL-R0663	Rat NOX1 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
E-EL-R0664	Rat NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
E-EL-R0665	Rat NADPH (Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
E-EL-R0666	Rat N-ACHR (Nicotinic Acetylcholine Receptor) ELISA Kit
E-EL-R0670	Rat NT-proBNP (N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
E-EL-R0671	Rat PIIINP (N-Terminal Procollagen III Propeptide) ELISA Kit
E-EL-R0672	Rat NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
E-EL-R0673	Rat NF-κB (Nuclear Factor Kappa B) ELISA Kit
E-EL-R0674	Rat NF-κB p65 (Nuclear Factor-κB P65) ELISA Kit
E-EL-R0675	Rat NRIP1 (Nuclear Receptor Interacting Protein 1) ELISA Kit
E-EL-R0676	Rat NXF1 (Nuclear RNA Export Factor 1) ELISA Kit
E-EL-R0677	Rat NUP107 (Nucleoporin 107kDa) ELISA Kit
E-EL-R0678	Rat NUP153 (Nucleoporin 153kDa) ELISA Kit

Cat.NO	Product Name
E-EL-R0679	Rat NUP155 (Nucleoporin 155kDa) ELISA Kit
E-EL-R0680	Rat NUP160 (Nucleoporin 160kDa) ELISA Kit
E-EL-R0681	Rat NUP188 (Nucleoporin 188kDa) ELISA Kit
E-EL-R0682	Rat NUP205 (Nucleoporin 205kDa) ELISA Kit
E-EL-R0683	Rat NUP214 (Nucleoporin 214kDa) ELISA Kit
E-EL-R0684	Rat NUP35 (Nucleoporin 35kDa) ELISA Kit
E-EL-R0685	Rat NUP37 (Nucleoporin 37kDa) ELISA Kit
E-EL-R0686	Rat NUP85 (Nucleoporin 85kDa) ELISA Kit
E-EL-R0687	Rat NUP198 (Nucleoporin 198kDa) ELISA Kit
E-EL-R0688	Rat OCT4 (Octamer Binding Transcription Factor 4) ELISA Kit
E-EL-R0689	Rat ITLN1 (Intelectin 1/Omentin) ELISA Kit
E-EL-R0690	Rat OIT3 (Oncoprotein Induced Transcript 3) ELISA Kit
E-EL-R0691	Rat OSMR (Oncostatin M Receptor) ELISA Kit
E-EL-R0692	Rat OSM (Oncostatin-M) ELISA Kit
E-EL-R0693	Rat OXA (Orexin A) ELISA Kit
E-EL-R0694	Rat ODC (Ornithine Decarboxylase) ELISA Kit
E-EL-R0695	Rat OFQ/N (Orphanin FQ/Nociceptin) ELISA Kit
E-EL-R0696	Rat GROβ (Growth Regulated Oncogene Beta) ELISA Kit
E-EL-R0697	Rat GROγ (Growth Regulated Oncogene Gamma) ELISA Kit
E-EL-R0698	Rat OSCAR (Osteoclast Associated Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-R0699	Rat ODF (Osteoclast Differentiation Factor) ELISA Kit
E-EL-R0700	Rat OSF (Osteoclast Stimulating Factor 1) ELISA Kit
E-EL-R0701	Rat ON (Osteonectin) ELISA Kit
E-EL-R0702	Rat OPN (Osteopontin) ELISA Kit
E-EL-R0703	Rat OPGL (Osteoprotegerin Ligand) ELISA Kit
E-EL-R0705	Rat OVA SIgM (Ovalbumin Specific IgM) ELISA Kit
E-EL-R0706	Rat OVA sIgE (Ovalbumin specific Immunoglobulin E) ELISA Kit
E-EL-R0707	Rat CA125 (Carbohydrate Antigen 125) ELISA Kit
E-EL-R0708	Rat Tryptophan Hydroxylase 1(TPH1) ELISA Kit
E-EL-R0709	Rat OLab (Oxidized Lowdensity Lipoprotein Antibody) ELISA Kit
E-EL-R0710	Rat OxLDL (Oxidized Lowdensity Lipoprotein) ELISA Kit
E-EL-R0711	Rat MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit
E-EL-R0712	Rat PUMA (P53 Upregulated Modulator of Apoptosis) ELISA Kit
E-EL-R0713	Rat Smad3 (Mothers Against Decapentaplegic Homolog 3) ELISA Kit
E-EL-R0714	Rat PTH (Parathyroid Hormone) ELISA Kit
E-EL-R0715	Rat PTHR2 (Parathyroid Hormone Receptor 2) ELISA Kit
E-EL-R0716	Rat PTHRP (Parathyroid Hormone Related Protein) ELISA Kit
E-EL-R0717	Rat PDAP1 (Pdgfa Associated Protein 1) ELISA Kit
E-EL-R0718	Rat PTD (Pentosidine) ELISA Kit
E-EL-R0719	Rat PP (Pepsin) ELISA Kit

Cat.NO	Product Name
E-EL-R0720	Rat PYY (Peptide YY) ELISA Kit
E-EL-R0721	Rat PRF1 (Perforin 1) ELISA Kit
E-EL-R0722	Rat Hp-IgM (Helicobacter Pylori IgM) ELISA Kit
E-EL-R0723	Rat PMP3 (Peroxisomal Membrane Protein 3) ELISA Kit
E-EL-R0724	Rat PPAR- γ (Peroxisome Proliferator-Activated Receptor γ) ELISA Kit
E-EL-R0725	Rat PPAR- α (Peroxisome Proliferators Activator Receptors α) ELISA Kit
E-EL-R0727	Rat PI Ab-IgG/IgM (Phosphatidylinositol Antibody IgG/IgM) ELISA Kit
E-EL-R0728	Rat P-PKC (Phospho Protein Kinase C) ELISA Kit
E-EL-R0729	Rat PERK (Phospho-Extracellular Signal-Regulated Kinase) ELISA Kit
E-EL-R0730	Rat PIKBA (Phospho-Inhibitory Subunit Of NF-KBA) ELISA Kit
E-EL-R0731	Rat PI3K2 α (Phosphoinositide 3 Kinase 2 Alpha) ELISA Kit
E-EL-R0732	Rat PLAA (Phospholipase A1 Activating Protein) ELISA Kit
E-EL-R0733	Rat PLA1 (Phospholipase A1) ELISA Kit
E-EL-R0735	Rat PLC γ 1 (Phospholipase C Gamma 1) ELISA Kit
E-EL-R0736	Rat PLC (Phospholipase C) ELISA Kit
E-EL-R0737	Rat pACCase (Phosphorylated Acetyl-CoA Carboxylase) ELISA Kit
E-EL-R0738	Rat AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) ELISA Kit
E-EL-R0739	Rat PI3K (Phosphotylinosital 3 Kinase) ELISA Kit
E-EL-R0740	Rat PICP (Procollagen I C-terminal Propeptide) ELISA Kit
E-EL-R0742	Rat PGF (Placental Growth Factor) ELISA Kit

Cat.NO	Product Name
E-EL-R0743	Rat PAP (Plasmin-Antiplasmin Complex) ELISA Kit
E-EL-R0744	Rat PAI-1 (Plasminogen Activator Inhibitor 1) ELISA Kit
E-EL-R0745	Rat PAI-2 (Plasminogen Activator Inhibitor 2) ELISA Kit
E-EL-R0747	Rat tPA (Plasminogen Activator, Tissue) ELISA Kit
E-EL-R0748	Rat Plg (Plasminogen) ELISA Kit
E-EL-R0749	Rat PAF (Platelet Activating Factor) ELISA Kit
E-EL-R0751	Rat PDGF (platelet-derived growth factor) ELISA Kit
E-EL-R0752	Rat PDGFA (Platelet Derived Growth Factor Subunit A) ELISA Kit
E-EL-R0753	Rat PDGFB (Platelet Derived Growth Factor Subunit B) ELISA Kit
E-EL-R0754	Rat PDGF-C (Platelet Derived Growth Factor C) ELISA Kit
E-EL-R0755	Rat PDGF-D (Platelet Derived Growth Factor D) ELISA Kit
E-EL-R0756	Rat PDGFRL (Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
E-EL-R0757	Rat PDGFsR- α (Platelet-Derived Growth Factor Soluble Receptor α) ELISA Kit
E-EL-R0758	Rat PF3 (Platelet Factor 3) ELISA Kit
E-EL-R0759	Rat PF-4 (Platelet Factor-4) ELISA Kit
E-EL-R0760	Rat GP- II b III a (Platelet Membrane Glycoprotein II b III a) ELISA Kit
E-EL-R0761	Rat PECAM1 (Platelet/Endothelial Cell Adhesion Molecule) ELISA Kit
E-EL-R0762	Rat PHLDA2 (Pleckstrin Homology Like Domain Family A Member 2) ELISA Kit
E-EL-R0763	Rat PTN (Pleiotrophin) ELISA Kit
E-EL-R0764	Rat PCX (Podocalyxin) ELISA Kit

Cat.NO	Product Name
E-EL-R0765	Rat PVRL1 (Poliovirus Receptor Related Protein 1) ELISA Kit
E-EL-R0766	Rat PVRL2 (Poliovirus Receptor Related Protein 2) ELISA Kit
E-EL-R0767	Rat PC4 (Positive Cofactor 4) ELISA Kit
E-EL-R0768	Rat PA/PLB (Pre-Albumin) ELISA Kit
E-EL-R0769	Rat FXI (Coagulation Factor XI) ELISA Kit
E-EL-R0770	Rat PRO-ANP (Pro Atrial Natriuretic Peptide) ELISA Kit
E-EL-R0771	Rat PCT (Procalcitonin) ELISA Kit
E-EL-R0772	Rat PINP (Procollagen I N-Terminal Peptide) ELISA Kit
E-EL-R0773	Rat P III NT (Procollagen III N-Terminal Peptide) ELISA Kit
E-EL-R0775	Rat Pg (Progesterone) ELISA Kit
E-EL-R0776	Rat PGR (Progesterone Receptor) ELISA Kit
E-EL-R0777	Rat FSTL1 (Follistatin Like Protein 1) ELISA Kit
E-EL-R0778	Rat PHB (Prohibitin) ELISA Kit
E-EL-R0779	Rat PI (Proinsulin) ELISA Kit
E-EL-R0780	Rat PK2 (Prokineticin 2) ELISA Kit
E-EL-R0781	Rat PKR1 (Prokineticin Receptor 1) ELISA Kit
E-EL-R0782	Rat PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit
E-EL-R0783	Rat GZMK (Granzyme K) ELISA Kit
E-EL-R0784	Rat PRCP (Prolylcarboxypeptidase) ELISA Kit
E-EL-R0785	Rat PGD2S (Prostaglandin D2 Synthase, Brain) ELISA Kit

Cat.NO	Product Name
E-EL-R0786	Rat EP2 (Prostaglandin E Receptor 2) ELISA Kit
E-EL-R0787	Rat PTGES2 (Prostaglandin E Synthase 2) ELISA Kit
E-EL-R0788	Rat PTGES (Prostaglandin E Synthase, Microsomal) ELISA Kit
E-EL-R0790	Rat NRG-1 (Neuregulin 1) ELISA Kit
E-EL-R0791	Rat PTGS1 (Prostaglandin Endoperoxide Synthase 1) ELISA Kit
E-EL-R0792	Rat PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit
E-EL-R0793	Rat FP (Prostaglandin F Receptor) ELISA Kit
E-EL-R0794	Rat PG-F (Prostaglandin F) ELISA Kit
E-EL-R0795	Rat PGF2 α (Prostaglandin F2 Alpha) ELISA Kit
E-EL-R0796	Rat PSA (Prostate Specific Antigen) ELISA Kit
E-EL-R0797	Rat PAP (Prostatic Acid Phosphatase) ELISA Kit
E-EL-R0798	Rat PAR2 (Protease Activated Receptor 2) ELISA Kit
E-EL-R0799	Rat PRSS1 (Protease, Serine, 1) ELISA Kit
E-EL-R0800	Rat PROC (Protein C) ELISA Kit
E-EL-R0801	Rat PROCR (Protein C Receptor, Endothelial) ELISA Kit
E-EL-R0802	Rat PDIA3 (Protein Disulfide Isomerase A3) ELISA Kit
E-EL-R0803	Rat PDI (Protein Disulfide Isomerase) ELISA Kit
E-EL-R0804	Rat PIAS1 (Protein Inhibitor Of Activated STAT 1) ELISA Kit
E-EL-R0805	Rat PICK1 (Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
E-EL-R0806	Rat PKBG (Protein Kinase B Gamma) ELISA Kit

Cat.NO	Product Name
E-EL-R0807	Rat PKB (Protein Kinase B) ELISA Kit
E-EL-R0808	Rat PKCB1 (Protein Kinase C Beta 1) ELISA Kit
E-EL-R0809	Rat PKCD (Protein Kinase C Delta) ELISA Kit
E-EL-R0810	Rat PKCE (Protein Kinase C Epsilon) ELISA Kit
E-EL-R0811	Rat PKCG (Protein Kinase C Gamma) ELISA Kit
E-EL-R0812	Rat PKCI (Protein Kinase C Iota) ELISA Kit
E-EL-R0813	Rat PKCq (Protein Kinase C Theta) ELISA Kit
E-EL-R0814	Rat PKCZ (Protein Kinase C Zeta) ELISA Kit
E-EL-R0815	Rat PKC (Protein Kinase C) ELISA Kit
E-EL-R0816	Rat PKD1 (Protein Kinase D1) ELISA Kit
E-EL-R0817	Rat PKD2 (Protein Kinase D2) ELISA Kit
E-EL-R0818	Rat PKIA (Protein Kinase Inhibitor Alpha) ELISA Kit
E-EL-R0819	Rat PKIB (Protein Kinase Inhibitor Beta) ELISA Kit
E-EL-R0820	Rat PKIG (Protein Kinase Inhibitor Gamma) ELISA Kit
E-EL-R0821	Rat PKN1 (Protein Kinase N1) ELISA Kit
E-EL-R0822	Rat PRKAA1 (Protein Kinase, AMP Activated Alpha 1) ELISA Kit
E-EL-R0823	Rat PP (Protein Phosphatase) ELISA Kit
E-EL-R0824	Rat PTPRN (Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
E-EL-R0825	Rat PTPRO (Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
E-EL-R0826	Rat PR3 (Proteinase 3) ELISA Kit

Cat.NO	Product Name
E-EL-R0827	Rat PCDH1 (Protocadherin 1) ELISA Kit
E-EL-R0828	Rat SELP (P-Selectin) ELISA Kit
E-EL-R0829	Rat GRAP2 (GRB2-Related Adapter protein 2) ELISA Kit
E-EL-R0830	Rat SPC (Pulmonary Surfactant Associated Protein C) ELISA Kit
E-EL-R0831	Rat SPD (Pulmonary Surfactant Associated Protein D) ELISA Kit
E-EL-R0834	Rat PDHA (Pyruvate Dehydrogenase Alpha) ELISA Kit
E-EL-R0835	Rat PDK4 (Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit
E-EL-R0836	Rat PDP (Pyruvate Dehydrogenase Phosphatase) ELISA Kit
E-EL-R0837	Rat PK (Pyruvate Kinase, Liver And RBC) ELISA Kit
E-EL-R0838	Rat M2-PK (Pyruvate Kinase,M2) ELISA Kit
E-EL-R0839	Rat GAP (Ras Gtpase Activating Protein) ELISA Kit
E-EL-R0840	Rat NRG-2 (Neuregulin 2) ELISA Kit
E-EL-R0841	Rat RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
E-EL-R0842	Rat GHRL (Ghrelin) ELISA Kit
E-EL-R0843	Rat REG1 α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit
E-EL-R0844	Rat REG4 (Regenerating Islet Derived Protein 4) ELISA Kit
E-EL-R0845	Rat RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
E-EL-R0846	Rat RGS4 (Regulator Of G Protein Signaling 4) ELISA Kit
E-EL-R0847	Rat RXFP1 (Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
E-EL-R0848	Rat RSV-IgA (Respiratory Syncytial Virus IgA) ELISA Kit

Cat.NO	Product Name
E-EL-R0849	Rat RSV- IgD (Respiratory Syncytial Virus IgD) ELISA Kit
E-EL-R0850	Rat RSV- IgG (Respiratory Syncytial Virus IgG) ELISA Kit
E-EL-R0851	Rat RSV- IgM (Respiratory Syncytial Virus IgM) ELISA Kit
E-EL-R0852	Rat RBP4 (Retinol Binding Protein 4, Plasma) ELISA Kit
E-EL-R0854	Rat Mly/CXCL9 (Monocyte Interferon Gamma Inducing Factor) ELISA Kit
E-EL-R0855	Rat RNASEA (Ribonuclease A) ELISA Kit
E-EL-R0856	Rat RNASEH (Ribonuclease H) ELISA Kit
E-EL-R0857	Rat RNASET2 (Ribonuclease T2) ELISA Kit
E-EL-R0858	Rat RPL10 (Ribosomal Protein L10) ELISA Kit
E-EL-R0859	Rat RPL26 (Ribosomal Protein L26) ELISA Kit
E-EL-R0860	Rat RNF4 (Ring Finger Protein 4) ELISA Kit
E-EL-R0865	Rat TACE (TNF α Converting Enzyme) ELISA Kit
E-EL-R0866	Rat RUNX2 (Runt Related Transcription Factor 2) ELISA Kit
E-EL-R0867	Rat S100A11 (S100 Calcium Binding Protein A11) ELISA Kit
E-EL-R0868	Rat S100B (S100 Calcium Binding Protein B) ELISA Kit
E-EL-R0870	Rat S100A6 (S100 Calcium-Binding Protein A6) ELISA Kit
E-EL-R0871	Rat FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
E-EL-R0872	Rat HDL-C (high density lipoprotein cholesterol) ELISA Kit
E-EL-R0873	Rat SPLA2 (Secreted Phospholipase A2) ELISA Kit
E-EL-R0874	Rat SCT (Secretin) ELISA Kit

Cat.NO	Product Name
E-EL-R0875	Rat sIgA (Secretory Immunoglobulin A) ELISA Kit
E-EL-R0876	Rat SLPI (Secretory Leukocyte Peptidase Inhibitor) ELISA Kit
E-EL-R0877	Rat SELE (Selectin, Endothelium) ELISA Kit
E-EL-R0878	Rat SEPP1 (Selenoprotein P) ELISA Kit
E-EL-R0879	Rat SEMA7A (Semaphorin 7A) ELISA Kit
E-EL-R0880	Rat STK (Serine/Threonine Protein Kinase) ELISA Kit
E-EL-R0881	Rat SAA (Serum Amyloid A) ELISA Kit
E-EL-R0882	Rat SOX2 (Sex Determining Region Y Box Protein 2) ELISA Kit
E-EL-R0883	Rat SHBG (Sex Hormone-Binding Globulin) ELISA Kit
E-EL-R0884	Rat STAT5B (Signal Transducer And Activator Of Transcription 5B) ELISA Kit
E-EL-R0885	Rat ACTN2 (Actinin Alpha 2) ELISA Kit
E-EL-R0886	Rat SLIT2 (Slit Homolog 2) ELISA Kit
E-EL-R0887	Rat SNRPD1 (Small Nuclear Ribonucleoprotein D1) ELISA Kit
E-EL-R0888	Rat sCD14 (Soluble Cluster of Differentiation 14) ELISA Kit
E-EL-R0889	Rat SCD30L (Soluble Cluster of Differentiation 30 Ligand) ELISA Kit
E-EL-R0891	Rat SCD86 (Soluble Cluster of Differentiation 41 Ligand) ELISA Kit
E-EL-R0892	Rat EPCR(Endothelial Protein C Receptor)ELISA Kit
E-EL-R0893	Rat SELE (E-Selectin) ELISA Kit
E-EL-R0894	Rat sFAS/sAPO-1 (Soluble Factor-Related Apoptosis) ELISA Kit
E-EL-R0895	Rat IL-2SR (Soluble Interleukin 2 Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-R0896	Rat IL-6R (Interleukin 6 Receptor) ELISA Kit
E-EL-R0897	Rat IL-1SR I (Soluble Interleukin-1 Receptor I) ELISA Kit
E-EL-R0898	Rat IL-1SR II (Soluble Interleukin-1 Receptor II) ELISA Kit
E-EL-R0899	Rat IL-2SR α (Soluble Interleukin-2 Receptor) ELISA Kit
E-EL-R0900	Rat IL-2SR β (Soluble Interleukin-2 Receptor) ELISA Kit
E-EL-R0902	Rat SPECAM-1 (Soluble Platelet Endothelial Cell Adhesion Molecule 1) ELISA Kit
E-EL-R0903	Rat S-100 (Soluble Protein-100) ELISA Kit
E-EL-R0907	Rat TLR-2 (Toll-Like Receptor 2) ELISA Kit
E-EL-R0908	Rat TLR-6 (Toll-Like Receptor 6) ELISA Kit
E-EL-R0910	Rat SVCAM-1 (Soluble Vascular Cell Adhesion Molecule 1) ELISA Kit
E-EL-R0911	Rat FLT1/VEGFR1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit
E-EL-R0912	Rat LDL-C (low density lipoprotein cholesterol) ELISA Kit
E-EL-R0913	Rat SSTR1 (Somatostatin Receptor 1) ELISA Kit
E-EL-R0914	Rat SST (Somatostatin) ELISA Kit
E-EL-R0916	Rat SSFA2 (Sperm Specific Antigen 2) ELISA Kit
E-EL-R0918	Rat SCFR (Stem Cell Factor Receptor) ELISA Kit
E-EL-R0919	Rat SCF (Stem Cell Factor) ELISA Kit
E-EL-R0920	Rat SF1 (Steroidogenic Factor 1) ELISA Kit
E-EL-R0921	Rat SREBP (Sterol Regulatory Element Binding Proteins) ELISA Kit
E-EL-R0922	Rat SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit

Cat.NO	Product Name
E-EL-R0923	Rat SDF-1 β (Stromal Cell Derived Factor 1 β) ELISA Kit
E-EL-R0924	Rat SDF2 (Stromal Cell Derived Factor 2) ELISA Kit
E-EL-R0925	Rat SDF2L1 (Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
E-EL-R0926	Rat SDF4 (Stromal Cell Derived Factor 4) ELISA Kit
E-EL-R0927	Rat SSRP1 (Structure Specific Recognition Protein 1) ELISA Kit
E-EL-R0930	Rat SOCS3 (Suppressors Of Cytokine Signaling 3) ELISA Kit
E-EL-R0931	Rat SMN1 (Survival of Motor Neuron 1, Telomeric) ELISA Kit
E-EL-R0932	Rat SYP (Synaptophysin) ELISA Kit
E-EL-R0933	Rat SNAP-25 (Synaptosome Associated Protein 25) ELISA Kit
E-EL-R0934	Rat SDC4 (Syndecan 4) ELISA Kit
E-EL-R0935	Rat ChAT (Choline Acetyltransferase) ELISA Kit
E-EL-R0936	Rat TACR1 (Tachykinin Receptor 1) ELISA Kit
E-EL-R0937	Rat TACR2 (Tachykinin Receptor 2) ELISA Kit
E-EL-R0938	Rat THP (Tamm–Horsfall Glycoprotein) ELISA Kit
E-EL-R0939	Rat TRACP-5b (Tartrate-Resistant Acid Phosphatase 5b) ELISA Kit
E-EL-R0940	Rat TBP (Tata Binding Protein) ELISA Kit
E-EL-R0941	Rat TAF1 (Tata Box Binding Protein Associated Factor 1) ELISA Kit
E-EL-R0942	Rat TAF4 (Tata Box Binding Protein Associated Factor 2) ELISA Kit
E-EL-R0943	Rat MAP τ (Microtubule Associated Protein Tau) ELISA Kit
E-EL-R0944	Rat TXLNA (Taxilin Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-R0945	Rat TXLNB (Taxilin Beta) ELISA Kit
E-EL-R0946	Rat TEK (Tyrosine Kinase, Endothelial) ELISA Kit
E-EL-R0947	Rat TE (Telomerase) ELISA Kit
E-EL-R0948	Rat TNC (Tenascin C) ELISA Kit
E-EL-R0949	Rat TNX (Tenascin X) ELISA Kit
E-EL-R0950	Rat TDGF1 (Teratocarcinoma Derived Growth Factor 1) ELISA Kit
E-EL-R0951	Rat TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit
E-EL-R0952	Rat TIEG1 (TGF-Beta-Inducible Early Response Gene-1) ELISA Kit
E-EL-R0953	Rat TrxR1 (Thioredoxin Reductase 1) ELISA Kit
E-EL-R0954	Rat TRX (Thioredoxin) ELISA Kit
E-EL-R0955	Rat TAFI (Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit
E-EL-R0956	Rat TR (Thrombin Receptor) ELISA Kit
E-EL-R0957	Rat MSRA (Methionine Sulfoxide Reductase A) ELISA Kit
E-EL-R0958	Rat TAT (Thrombin-Antithrombin Complex) ELISA Kit
E-EL-R0959	Rat MSRB1 (Methionine Sulfoxide Reductase B1) ELISA Kit
E-EL-R0960	Rat TM (Thrombomodulin) ELISA Kit
E-EL-R0961	Rat TPO-Ab(Thrombopoietin Antibody)ELISA Kit
E-EL-R0962	Rat TPOR(Thrombopoietin Receptor)ELISA Kit
E-EL-R0963	Rat TPO (Thrombopoietin) ELISA Kit
E-EL-R0964	Rat TSP-1(Thrombospondin 1)ELISA Kit

Cat.NO	Product Name
E-EL-R0965	Rat TXB2 (Thromboxane B2) ELISA Kit
E-EL-R0966	Rat TXB3 (Thromboxanes B3) ELISA Kit
E-EL-R0967	Rat TPP (Thrombus Precursor Protein) ELISA Kit
E-EL-R0968	Rat TSLP (Thymic Stromal Lymphopoietin) ELISA Kit
E-EL-R0969	Rat TARC (Thymus Activation Regulated Chemokine) ELISA Kit
E-EL-R0970	Rat TECK (Thymus Expressed Chemokine) ELISA Kit
E-EL-R0971	Rat TG (Thyroglobulin) ELISA Kit
E-EL-R0972	Rat HDAC2 (Histone Deacetylase 2) ELISA Kit
E-EL-R0974	Rat GPC3 (Glypican 3) ELISA Kit
E-EL-R0975	Rat HDAC1 (Histone Deacetylase 1) ELISA Kit
E-EL-R0976	Rat TSH (Thyroid Stimulating Hormone) ELISA Kit
E-EL-R0977	Rat COL3 α 1 (Collagen Type III Alpha 1) ELISA Kit
E-EL-R0978	Rat TEF (Thyrotroph Embryonic Factor) ELISA Kit
E-EL-R0979	Rat TRH (Thyrotropin-Releasing Hormone) ELISA Kit
E-EL-R0980	Rat TAB (Thyroxine Antibody) ELISA Kit
E-EL-R0982	Rat TJP1 (Tight Junction Protein 1) ELISA Kit
E-EL-R0983	Rat TFPI (Tissue Factor Pathway Inhibitor) ELISA Kit
E-EL-R0984	Rat TF (Tissue Factor) ELISA Kit
E-EL-R0985	Rat TIMP-2(Tissue Inhibitors of Metalloproteinase 2)ELISA Kit
E-EL-R0986	Rat TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit

Cat.NO	Product Name
E-EL-R0987	Rat TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) ELISA Kit
E-EL-R0988	Rat TPA (Tissue Polypeptide Antigen) ELISA Kit
E-EL-R0989	Rat SERP1 (Stress-associated Endoplasmic Reticulum Protein 1) ELISA Kit
E-EL-R0990	Rat TLR4 (Toll-Like Receptor 4) ELISA Kit
E-EL-R0991	Rat TLR9 (Toll-Like Receptor 9) ELISA Kit
E-EL-R0995	Rat c-myc (c-myc Oncogene Product) ELISA Kit
E-EL-R0996	Rat SDC1 (Syndecan 1) ELISA Kit
E-EL-R0999	Rat 5-LO (Arachidonate 5-Lipoxygenase) ELISA Kit
E-EL-R1000	Rat 12-LO (Arachidonate 12-Lipoxygenase) ELISA Kit
E-EL-R1001	Rat 15-LO (Arachidonate 15-Lipoxygenase) ELISA Kit
E-EL-R1002	Rat PTX3 (Pentraxin 3) ELISA Kit
E-EL-R1003	Rat 20S-PSM (20S-Proteasome) ELISA Kit
E-EL-R1004	Rat CK-18/KRT18 (Cytokeratin 18) ELISA Kit
E-EL-R1005	Rat TCF20 (Transcription factor 20) ELISA Kit
E-EL-R1006	Rat ERR α (Estrogen Related Receptor Alpha) ELISA Kit
E-EL-R1007	Rat α 2-AP (α 2-Antiplasmin) ELISA Kit
E-EL-R1008	Rat TFR2 (Transferrin Receptor Protein 2) ELISA Kit
E-EL-R1009	Rat TFR (Transferrin Receptor) ELISA Kit
E-EL-R1010	Rat IL-24 (Interleukin 24) ELISA Kit
E-EL-R1011	Rat TGF- β R (Transforming Growth Factor Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-R1012	Rat TGF- α (Transforming Growth Factor Alpha) ELISA Kit
E-EL-R1014	Rat GP4/CD36 (Platelet Membrane Glycoprotein IV) ELISA Kit
E-EL-R1015	Rat TGF- β 2 (Transforming Growth Factor Beta 2) ELISA Kit
E-EL-R1016	Rat TGF- β 3 (Transforming Growth Factor Beta 3) ELISA Kit
E-EL-R1017	Rat TGFBI (Transforming Growth Factor Beta Induced Protein) ELISA Kit
E-EL-R1018	Rat TSC22 (Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
E-EL-R1019	Rat TG (Transglutaminase) ELISA Kit
E-EL-R1020	Rat TMED1(Transmembrane Emp24 Domain Containing Protein 1) ELISA Kit
E-EL-R1021	Rat TNPO1(Transportin 1) ELISA Kit
E-EL-R1022	Rat TTR (Transthyretin) ELISA Kit
E-EL-R1023	Rat TFF2 (Trefoil Factor 2) ELISA Kit
E-EL-R1024	Rat TFF3 (Trefoil Factor 3, Intestinal) ELISA Kit
E-EL-R1025	Rat TREM-1 (Triggering Receptor Expresses on Myeloid Cells-1) ELISA Kit
E-EL-R1027	Rat INHB (Inhibin B) ELISA Kit
E-EL-R1029	Rat TRY (Trypsin) ELISA Kit
E-EL-R1030	Rat TAP (Trypsinogen Activation Peptide) ELISA Kit
E-EL-R1033	Rat TUBB (Tubulin Beta) ELISA Kit
E-EL-R1034	Rat TPPP (Tubulin Polymerization Promoting Protein) ELISA Kit
E-EL-R1036	Rat TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
E-EL-R1037	Rat TNF- β (Tumor Necrosis Factor β) ELISA Kit

Cat.NO	Product Name
E-EL-R1038	Rat TRAIL R1 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand 1) ELISA Kit
E-EL-R1039	Rat TRAIL R3 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand 3) ELISA Kit
E-EL-R1040	Rat TRAIL R4 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand 4) ELISA Kit
E-EL-R1041	Rat TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit
E-EL-R1042	Rat TP53 (Tumor Protein 53) ELISA Kit
E-EL-R1043	Rat TSGF (Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
E-EL-R1045	Rat TYR (Tyrosinase) ELISA Kit
E-EL-R1046	Rat TRKB (Tyrosine Kinase B) ELISA Kit
E-EL-R1047	Rat TIE-2 (Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2) ELISA Kit
E-EL-R1049	Rat U-T4 (Ultrasensitivity Thyroxine) ELISA Kit
E-EL-R1050	Rat U-T3 (Ultrasensitivity Triiodothyronine) ELISA Kit
E-EL-R1052	Rat NFE2L2 (Nuclear Factor, Erythroid Derived 2 Like Protein 2) ELISA Kit
E-EL-R1053	Rat PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit
E-EL-R1055	Rat PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit
E-EL-R1056	Rat VEGF-B (Vascular Endothelial Cell Growth Factor B) ELISA Kit
E-EL-R1057	Rat VEGF-C (Vascular Endothelial Cell Growth Factor C) ELISA Kit
E-EL-R1058	Rat VEGF-D (Vascular Endothelial Cell Growth Factor D) ELISA Kit
E-EL-R1060	Rat VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
E-EL-R1061	Rat VCAM-1 (Vascular Cell Adhesion Molecule 1) ELISA Kit
E-EL-R1062	Rat VIP (Vasoactive Intestinal Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-R1063	Rat VLA (Very Late Appearing Antigen) ELISA Kit
E-EL-R1064	Rat FOS (V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) ELISA Kit
E-EL-R1065	Rat HRAS (V-Ha-Ras Harvey Rat Sarcoma Viral Oncogene Homolog) ELISA Kit
E-EL-R1066	Rat VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
E-EL-R1067	Rat VF (Visfatin) ELISA Kit
E-EL-R1072	Rat DBP (Vitamin D Binding Protein) ELISA Kit
E-EL-R1078	Rat ADAMTS13/VWFPCP(Von Willebrand Factor Cleaving Protease)ELISA Kit
E-EL-R1079	Rat VWF (Von Willebrand Factor) ELISA Kit
E-EL-R1080	Rat XIAP (X-Linked Inhibitor Of Apoptosis Protein) ELISA Kit
E-EL-R1081	Rat α 1-AGP (α 1-Acid Glycoprotein) ELISA Kit
E-EL-R1082	Rat α 1-MG (α 1-Microglobulin) ELISA Kit
E-EL-R1083	Rat α -Glu (Alpha Glucosidase) ELISA Kit
E-EL-R1084	Rat β 2-GP1 IgA/G/M(β 2-Glycoprotein 1 Antibody IgA/G/M) ELISA Kit
E-EL-R1085	Rat BMG/ β 2-MG (Beta-2-Microglobulin) ELISA Kit
E-EL-R1087	Rat β -TG (β -Thromboglobulin) ELISA Kit
E-EL-R1088	Rat GLTP (Glycolipid Transfer Protein) ELISA Kit
E-EL-R1089	Rat GDF6 (Growth Differentiation Factor 6) ELISA Kit
E-EL-R1090	Rat p τ (phospho Tau Protein) ELISA Kit
E-EL-R1091	Rat TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
E-EL-R1093	Rat MDH1 (Malate Dehydrogenase 1) ELISA Kit

Cat.NO	Product Name
E-EL-R1094	Rat Ta1 (Thymosin- α 1) ELISA Kit
E-EL-R1095	Rat LEPR (Leptin Receptor) ELISA Kit
E-EL-R1096	Rat PCIII (Procollagen III) ELISA Kit
E-EL-R1097	Rat fT3 (Free Triiodothyronine) ELISA Kit
E-EL-R1098	Rat 5-HIAA (5-Hydroxyindolacetic Acid) ELISA Kit
E-EL-R1099	Rat FLAP (5-Lipoxygenase Activating Protein) ELISA Kit
E-EL-R1100	Rat PON1 (Paraoxonase 1) ELISA Kit
E-EL-R1101	Rat DNMT1 (DNA Methyltransferase 1) ELISA Kit
E-EL-R1102	Rat SIRT1 (NAD-dependent protein deacetylase sirtuin-1) ELISA Kit
E-EL-R1103	Rat SCUBE1 (Signal Peptide, CUB and EGF-like Domain-containing Protein 1) ELISA Kit
E-EL-R1104	Rat FNDC5 (Fibronectin type III domain-containing protein 5) ELISA Kit
E-EL-R1105	Rat GNMT (Glycine N-methyltransferase) ELISA Kit
E-EL-R1124	Rat HSD11 β 1 (11-Beta-Hydroxysteroid Dehydrogenase Type 1) ELISA Kit
E-EL-R1125	Rat FG (Fibrinogen) ELISA Kit
E-EL-R1156	Rat MRC1 (Mannose Receptor C Type 1) ELISA Kit
E-EL-R1157	Rat PEMT (Phosphatidylethanolamine N-methyltransferase) ELISA Kit
E-EL-R1179	Rat TRF (Transferrin) ELISA Kit
E-EL-R1203	Rat VLDL (Very Low Density Lipoprotein) ELISA Kit
E-EL-R1205	Rat SDH (Succinate Dehydrogenase) ELISA Kit
E-EL-R1206	Rat F1+2 (Prothrombin Fragment 1+2) ELISA Kit

Cat.NO	Product Name
E-EL-R1207	Rat VASP (Vasodilator Stimulated Phosphoprotein) ELISA Kit
E-EL-R1208	Rat DNMT3A (DNA Methyltransferase 3A) ELISA Kit
E-EL-R1209	Rat VAP-1 (Vascular Adhesion Protein 1) ELISA Kit
E-EL-R1210	Rat ApoA1 (Apolipoprotein A1) ELISA Kit
E-EL-R1211	Rat ICD (Isocitrate Dehydrogenase) ELISA Kit
E-EL-R1212	Rat α -KGDHC (α Ketoglutarate Dehydrogenase) ELISA Kit
E-EL-R1213	Rat ND1 (NADH Dehydrogenase 1) ELISA Kit
E-EL-R1214	Rat 6-PFK (Phosphofructokinase) ELISA Kit
E-EL-R1215	Rat ApoA2 (Apolipoprotein A2) ELISA Kit
E-EL-R1217	Rat SNCA (Synuclein Alpha) ELISA Kit
E-EL-R1218	Rat ApoB (Apolipoprotein B) ELISA Kit
E-EL-R1219	Rat AMACR(Alpha-methylacyl-CoA racemase) ELISA Kit
E-EL-R1221	Rat CCNA1(Cyclin-A1) ELISA Kit
E-EL-R1222	Rat RELN (Reelin) ELISA Kit
E-EL-R1223	Rat IL-27 (Interleukin 27) ELISA Kit
E-EL-R1224	Rat GPIb (Platelet Membrane Glycoprotein Ib) ELISA Kit
E-EL-R1225	Rat GP6 (Glycoprotein VI, Platelet) ELISA Kit
E-EL-R1226	Rat STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit
E-EL-R1227	Rat JAK2 (Janus Kinase 2) ELISA Kit
E-EL-R1228	Rat SOCS3 (Suppressors of Cytokine Signaling 3) ELISA Kit

Cat.NO	Product Name
E-EL-R1229	Rat FPN (Ferroportin) ELISA Kit
E-EL-R1230	Rat ApoE (Apolipoprotein E) ELISA Kit
E-EL-R1231	Rat DNMT3B (DNA Methyltransferase 3B) ELISA Kit
E-EL-R1232	Rat ALT (Alanine Transaminase) ELISA Kit
E-EL-R1235	Rat BDNF (Brain Derived Neurotrophic Factor) ELISA Kit
E-EL-R1253	Rat cTn-I/TNNI3 (Cardiac Troponin I) ELISA Kit
E-EL-R1327	Rat CK-MB (Creatine Kinase MB Isoenzyme) ELISA Kit
E-EL-R1401	Rat A β 40 (Amyloid Beta 40) ELISA Kit
E-EL-R1402	Rat A β 42 (Amyloid Beta 42) ELISA Kit
E-EL-R1403	Rat CYP7A1 (Cholesterol 7-Alpha-Hydroxylase) ELISA Kit
E-EL-R1404	Rat α -CTx (Alpha Crosslaps) ELISA Kit
E-EL-R1405	Rat β -CTx(Beta Crosslaps)ELISA Kit
E-EL-R1406	Rat GALR1 (Galanin Receptor1) ELISA Kit
E-EL-R1407	Rat GALR2 (Galanin Receptor2) ELISA Kit
E-EL-R1408	Rat GALR3 (Galanin Receptor3) ELISA Kit
E-EL-R1409	Rat GALR4 (Galanin Receptor4) ELISA Kit
E-EL-R1410	Rat IL-2R β (Interleukin-2 Receptor- β) ELISA Kit
E-EL-R1411	Rat PGA (Pepsinogen A) ELISA Kit
E-EL-R1412	Rat PGC (PepsinogenC) ELISA Kit
E-EL-R1413	Rat NT-Pro-ANP (N-Terminal Pro-Atrial Natriuretic Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-R1414	Rat PINP (Procollagen Type I N-Terminal Propeptide) ELISA Kit
E-EL-R1420	Rat Rt3 (Reverse Triiodothyronine) ELISA Kit
E-EL-R1421	Rat TGM2 (Transglutaminase 2, Tissue) ELISA Kit
E-EL-R1423	Rat 11-DH-TXB2 (11-Dehydrothromboxane B2) ELISA Kit
E-EL-R1424	Rat SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit
E-EL-R1425	Rat SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-R1426	Rat SOD3 (Superoxide Dismutase 3, Extracellular) ELISA Kit
E-EL-R1427	Rat SOD4 (Superoxide Dismutase 4, Copper Chaperone) ELISA Kit
E-EL-R1428	Rat GFAP (Glial Fibrillary Acidic Protein) ELISA Kit
E-EL-R1429	Rat ANG I (Angiotensin I) ELISA Kit
E-EL-R1430	Rat ANG II (Angiotensin II) ELISA Kit
E-EL-R1431	Rat ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-R1432	Rat β MSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-R1433	Rat PS1 (Presenilin 1) ELISA Kit
E-EL-R1434	Rat AgRP (Agouti Related Protein) ELISA Kit
E-EL-R1435	Rat PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38)ELISA Kit
E-EL-R1436	Rat PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-R1437	Rat TH (Tyrosine Hydroxylase) ELISA Kit
E-EL-R1438	Rat NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) ELISA Kit
E-EL-R1456	Rat CTXI (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit

Cat.NO	Product Name
E-EL-R1503	Rat RTN4 (Reticulon 4) ELISA Kit
E-EL-R2373	Rat FASL (Factor Related Apoptosis Ligand) ELISA Kit
E-EL-R2376	Rat F VII (Coagulation Factor VII) ELISA Kit
E-EL-R2377	Rat F VIII (Coagulation Factor VIII) ELISA Kit
E-EL-R2378	Rat AQP-2 (Aquaporin 2) ELISA Kit
E-EL-R2379	Rat ApoC3 (Apolipoprotein C3) ELISA Kit
E-EL-R2380	Rat AGRN (Agrin) ELISA Kit
E-EL-R2386	Rat EGFR (Epidermal Growth Factor Receptor) ELISA Kit
E-EL-R2387	Rat PACAP (Pituitary Adenylate Cyclase Activating Polypeptide) ELISA Kit
E-EL-R2388	Rat 5-HETE (5-Hydroxyeicosatetraenoic Acid) ELISA Kit
E-EL-R2389	Rat CALP (Calprotectin) ELISA Kit
E-EL-R2391	Rat ANXA5 (Annexin A5) ELISA Kit
E-EL-R2392	Rat GPIb-IX-V (Platelet Membrane Glycoprotein GPIb-IX-V Complex) ELISA Kit
E-EL-R2393	Rat SCARB1 (Scavenger receptor class B member 1) ELISA Kit
E-EL-R2394	Rat IL-25 (Interleukin 25) ELISA Kit
E-EL-R2395	Rat POSTN/OSF-2 (Periostin) ELISA Kit
E-EL-R2396	Rat BAFF/CD257 (B-Cell Activating Factor) ELISA Kit
E-EL-R2397	Rat NFκB2 (Nuclear Factor Kappa B2) ELISA Kit
E-EL-R2398	Rat RELB (V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
E-EL-R2399	Rat PDPN (Podoplanin) ELISA Kit

Cat.NO	Product Name
E-EL-R2400	Rat PCT (Procalcitonin) ELISA Kit
E-EL-R2401	Rat ACE (Angiotensin I Converting Enzyme) ELISA Kit
E-EL-R2402	Rat FMOD (Fibromodulin) ELISA Kit
E-EL-R2403	Rat NNE (Non-Neuronal Enolase) ELISA Kit
E-EL-R2404	Rat cADPRH/CD38 (Cyclic ADP Ribose Hydrolase) ELISA Kit
E-EL-R2405	Rat ADM (Adrenomedullin) ELISA Kit
E-EL-R2406	Rat NPHN (Nephrin) ELISA Kit
E-EL-R2407	Rat AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit
E-EL-R2408	Rat FGF21(Fibroblast Growth Factor 21) ELISA Kit
E-EL-R2409	Rat FGF19 (Fibroblast Growth Factor 19) ELISA Kit
E-EL-R2410	Rat FGF23 (Fibroblast Growth Factor 23) ELISA Kit
E-EL-R2411	Rat FGFR1 (Fibroblast Growth Factor Receptor 1) ELISA Kit
E-EL-R2412	Rat PEDF (Pigment Epithelium Derived Factor) ELISA Kit
E-EL-R2413	Rat APLN (Apelin) ELISA Kit
E-EL-R2414	Rat CTSB (Cathepsin B) ELISA Kit
E-EL-R2415	Rat IgG2a (Immunoglobulin G2a) ELISA Kit
E-EL-R2416	Rat IgG2b (Immunoglobulin G2b) ELISA Kit
E-EL-R2417	Rat IgG2c (Immunoglobulin G2c) ELISA Kit
E-EL-R2418	Rat MMP-11 (Matrix Metalloproteinase 11) ELISA Kit
E-EL-R2419	Rat EREG (Epiregulin) ELISA Kit

Cat.NO	Product Name
E-EL-R2420	Rat α 1-AT (Alpha 1-Antitrypsin) ELISA Kit
E-EL-R2421	Rat MMP-14 (Matrix Metalloproteinase 14) ELISA Kit
E-EL-R2422	Rat Fc ϵ RII/CD23 (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
E-EL-R2423	Rat IL-34 (Interleukin 34) ELISA Kit
E-EL-R2424	Rat PCK1 (Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
E-EL-R2425	Rat ACT α 2 (Actin Alpha 2, Smooth Muscle) ELISA Kit
E-EL-R2426	Rat IL18BP (Interleukin 18 Binding Protein) ELISA Kit
E-EL-R2427	Rat SOST (Sclerostin) ELISA Kit
E-EL-R2428	Rat FOXP3 (Forkhead Box Protein P3) ELISA Kit
E-EL-R2429	Rat NT-3 (Neurotrophin-3) ELISA Kit
E-EL-R2430	Rat COL2 α 1 (Collagen Type II Alpha 1) ELISA Kit
E-EL-R2431	Rat CD163 (Cluster of Differentiation 163) ELISA Kit
E-EL-R2432	Rat CFD (Complement Factor D) ELISA Kit
E-EL-R2433	Rat I-TAC (Interferon Inducible T-Cell Alpha Chemoattractant) ELISA Kit
E-EL-R2434	Rat ENG (Endoglin) ELISA Kit
E-EL-R2435	Rat FABP1 (Fatty Acid Binding Protein 1, Liver) ELISA Kit
E-EL-R2436	Rat FABP4 (Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit
E-EL-R2437	Rat IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit
E-EL-R2438	Rat IL-1R1 (Interleukin 1 Receptor Type I) ELISA Kit
E-EL-R2439	Rat MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit

Cat.NO	Product Name
E-EL-R2440	Rat IL-22 (Interleukin 22) ELISA Kit
E-EL-R2441	Rat PL (Pancreatic Lipase) ELISA Kit
E-EL-R2442	Rat IL-12 p35 (Interleukin 12 p35) ELISA Kit
E-EL-R2443	Rat CALD (Caldesmon) ELISA Kit
E-EL-R2444	Rat FTH (Ferritin, Heavy Polypeptide) ELISA Kit
E-EL-R2445	Rat IL-19 (Interleukin 19) ELISA Kit
E-EL-R2446	Rat IL-20 (Interleukin 20) ELISA Kit
E-EL-R2447	Rat ALCAM (Activated Leukocyte Cell Adhesion Molecule) ELISA Kit
E-EL-R2448	Rat IAPP (Islet Amyloid Polypeptide) ELISA Kit
E-EL-R2449	Rat ANGPTL3 (Angiopoietin Like Protein 3) ELISA Kit
E-EL-R2450	Rat MAG (Myelin Associated Glycoprotein) ELISA Kit
E-EL-R2451	Rat FETUA (Fetuin A) ELISA Kit
E-EL-R2452	Rat DKK3 (Dickkopf Related Protein 3) ELISA Kit
E-EL-R2453	Rat ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit
E-EL-R2455	Rat IDE (Insulin Degrading Enzyme) ELISA Kit
E-EL-R2456	Rat CAT(Catalase) ELISA Kit
E-EL-R2457	Rat TBA(Total bile acid) Kit
E-EL-R2458	Rat Neopterin ELISA Kit
E-EL-R2459	Rat CD4(Cluster of Differentiation 4)ELISA Kit
E-EL-R2460	Rat AGEs (advanced glycation end products)ELISA Kit

Cat.NO	Product Name
E-EL-R2461	Rat GRIA1(Glutamate receptor 1)ELISA Kit
E-EL-R2462	Rat AKR1B1(Aldose reductase)ELISA Kit
E-EL-R2463	Rat SDH (Sorbitol Dehydrogenase) ELISA Kit
E-EL-R2464	Rat GYS2 (Glycogen Synthase 2, Liver) ELISA Kit
E-EL-R2465	Rat GMP-140 (Alpha-granular Membrane Protein) ELISA Kit
E-EL-R2466	Rat INS (Insulin) ELISA Kit
E-EL-R2467	Rat CYP1A1 (Cytochrome P450, family 1, subfamily A, polypeptide 1) ELISA Kit
E-EL-R2468	Rat CYP1A2 (Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit
E-EL-R2469	Rat CYP1B1 (Cytochrome P450, family 1, subfamily B, polypeptide 1) ELISA Kit
E-EL-R2470	Rat CYP2C9 (Cytochrome P450, family 2, subfamily C, polypeptide 9) ELISA Kit
E-EL-R2471	Rat CYP2C19(Cytochrome P450, family 2, subfamily C, polypeptide 19) ELISA Kit
E-EL-R2472	Rat CYP3A4 (Cytochrome P450, family 3, subfamily A, polypeptide 4) ELISA Kit
E-EL-R2473	Rat CYP3A5 (Cytochrome P450, family 3, subfamily A, polypeptide 5) ELISA Kit
E-EL-R2474	Rat ALAD(Aminolevulinate Delta Dehydratase)ELISA Kit
E-EL-R2475	Rat LIPE(Lipase, Hormone Sensitive)ELISA Kit
E-EL-R2476	Rat cGMP(Cyclic GMP)ELISA Kit
E-EL-R2477	Rat ucOC(Undercarboxylated Osteocalcin)ELISA Kit
E-EL-R2478	Rat UCHL1 (Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
E-EL-R2479	Rat TDP43 (TAR DNA Binding Protein 43) ELISA Kit
E-EL-R2480	Rat PRGN(progranulin)ELISA Kit

Cat.NO	Product Name
E-EL-R2481	Rat LPO(Lipid Peroxide)ELISA Kit
E-EL-R2482	Rat 8-OHdG(8-Hydroxydeoxyguanosine)ELISA Kit
E-EL-R2483	Rat sCD80(Soluble Cluster of Differentiation 80)ELISA Kit
E-EL-R2484	Rat ESM1(Endothelial Cell Specific Molecule 1)ELISA Kit
E-EL-R2485	Rat Pro-ADM(Proadrenomedullin)ELISA Kit
E-EL-R2486	Rat IMA(Ischemia Modified Albumin)ELISA Kit
E-EL-R2487	Rat PCSK9(Proprotein convertase subtilisin/kexin type 9)ELISA Kit
E-EL-R2488	Rat 8-iso-PG(8-iso prostaglandin)ELISA Kit
E-EL-R2489	Rat γ -Secretase ELISA Kit
E-EL-R2490	Rat APP(Amyloid Precursor Protein)ELISA Kit
E-EL-R2491	Rat GPX1 (Glutathione Peroxidase 1) ELISA Kit
E-EL-R2492	Rat CHEM(Chemerin)ELISA Kit
E-EL-R2493	Rat PKC α (Protein Kinase C α)ELISA Kit
E-EL-R2494	Rat GS(Glutamine synthetase)ELISA Kit
E-EL-R2495	Rat SDF-1 α (Stromal Cell Derived Factor 1 α) ELISA Kit
E-EL-R2496	Rat dGHR(Desacyl Ghrelin) ELISA Kit
E-EL-R2497	Rat esRAGE(Endogenous secretory Advanced Glycosylation End Product Specific Receptor)ELISA Kit
E-EL-R2498	Rat AT1R(angiotensin II receptor type 1)ELISA Kit
E-EL-R2499	Rat ADRb1(Adrenergic Receptor Beta 1)ELISA Kit
E-EL-R2500	Rat ADR β K1(Adrenergic Receptor Beta Kinase 1)ELISA Kit

Cat.NO	Product Name
E-EL-R2501	Rat MIP-2 (Macrophage Inflammatory Protein 2) ELISA Kit
E-EL-R2502	Rat CLDN5(Claudin-5)ELISA Kit
E-EL-R2503	Rat OCLN(Occludin)ELISA Kit
E-EL-R2504	Rat F-actin ELISA Kit
E-EL-R2505	Rat EP (Erythrocyte Protoporphyrin) ELISA Kit
E-EL-R2506	Rat PDCD1LG1 (Programmed Cell Death Protein 1 Ligand 1) ELISA Kit
E-EL-R2507	Rat ADORA1(Adenosine A1 Receptor)ELISA Kit
E-EL-R2508	Rat ADORA2a(Adenosine A2a Receptor)ELISA Kit
E-EL-R2509	Rat ADORA2b(Adenosine A2b Receptor)ELISA Kit
E-EL-R2510	Rat ADORA3(Adenosine A3 Receptor)ELISA Kit
E-EL-R2511	Rat VMA (Vanillylmandelic Acid) ELISA Kit
E-EL-R2512	Rat ACC(Acetyl-CoA Carboxylase)ELISA Kit
E-EL-R2513	Rat PC(Pyruvate Carboxylase)ELISA Kit
E-EL-R2514	Rat NES(Nesfatin)ELISA Kit
E-EL-R2515	Rat PLA2G2A(Phospholipase A2,Group IIA)ELISA Kit
E-EL-R2516	Rat HSD11B1(Corticosteroid 11-beta-dehydrogenase isozyme 1)ELISA Kit
E-EL-R2517	Rat SEMA3A(Semaphorin-3A)ELISA Kit
E-EL-R2518	Rat CAPN2(Calpain-2 catalytic subunit)ELISA Kit
E-EL-R2519	Rat β -Hex A (β -hexosaminidase A) ELISA Kit
E-EL-R2520	Rat β -Hex B (β -hexosaminidase B) ELISA Kit

Cat.NO	Product Name
E-EL-R2521	Rat CD28 (Cluster of Differentiation 28) ELISA Kit
E-EL-R2522	Rat PDCD1LG2(Programmed Cell Death Protein 1 Ligand 2) ELISA Kit
E-EL-R2523	Rat ELA2(Elastase 2, Neutrophil) ELISA Kit
E-EL-R2524	Rat ICOS(Inducible T-cell costimulator)ELISA Kit
E-EL-R2525	Rat ICOSL(Glnducible T-Cell Co Stimulator Ligand)ELISA Kit
E-EL-R2526	Rat TPO(Thyroid Peroxidase) ELISA Kit
E-EL-R2527	Rat CF6(Coupling Factor 6)ELISA Kit
E-EL-R2528	Rat ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit
E-EL-R2529	Rat ERK1/2(Extracellular Signal Regulated Kinase 1/2) ELISA Kit
E-EL-R2530	Rat KISS1 (Kisspeptin 1) ELISA Kit
E-EL-R2531	Rat PIINP(Procollagen II N-Terminal Propeptide)ELISA Kit
E-EL-R2532	Rat AChRab (Acetylcholine Receptor Antibody) ELISA Kit
E-EL-R2533	Rat HYOU1 (Hypoxia Up Regulated 1) ELISA Kit
E-EL-R2534	Rat GDN(Glia Derived Nexin)ELISA Kit
E-EL-R2535	Rat ACP (Acid Phosphatase) ELISA Kit
E-EL-R2536	Rat NEFL (Neurofilament, Light Polypeptide) ELISA Kit
E-EL-R2537	Rat CAPN1(calpain 1) ELISA Kit
E-EL-R2538	Rat AZU/HBP(Azurocidin/Heparin Binding Protein)ELISA Kit
E-EL-R2539	Rat CA19-9 (Carbohydrate antigen19-9) ELISA Kit
E-EL-R2540	Rat AsAb(Anti-sperm Antibody)ELISA Kit

Cat.NO	Product Name
E-EL-P0001	Porcine IL-8 (Interleukin 8) ELISA Kit
E-EL-P0002	Porcine IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-P0003	Porcine IFN- γ (Interferon Gamma) ELISA Kit
E-EL-P0004	Porcine IgG (Immunoglobulin G) ELISA Kit
E-EL-P0005	Porcine IL-10 (Interleukin-10) ELISA Kit
E-EL-P0006	Porcine IL-12 p40 (Interleukin 12 p40) ELISA Kit
E-EL-P0007	Porcine IL-18 (Interleukin 18) ELISA Kit
E-EL-P0008	Porcine IL-23 p40 (Interleukin 23 p40) ELISA Kit
E-EL-P0009	Porcine TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit
E-EL-P0010	Porcine TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-P0013	Porcine ALB (Albumin) ELISA Kit
E-EL-P0014	Porcine HP (Haptoglobin) ELISA Kit
E-EL-P0015	Porcine CRP (C-Reactive Protein) ELISA Kit
E-EL-P0016	Porcine PLG (Plasminogen) ELISA Kit
E-EL-P0017	Porcine FBG (Fibrinogen) ELISA Kit
E-EL-P0018	Porcine MYO (Myoglobin) ELISA Kit
E-EL-P0019	Porcine FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
E-EL-P0020	Porcine IL-6 (Interleukin 6) ELISA Kit
E-EL-P0021	Porcine PRL (Prolactin) ELISA Kit
E-EL-P0022	Porcine ANP (Atrial Natriuretic Peptide) ELISA Kit

E-EL-P0023	Porcine BNP (Brain Natriuretic Peptide) ELISA Kit
E-EL-P0024	Porcine FN (Fibronectin) ELISA Kit
E-EL-P0025	Porcine DPT (Dermatopontin) ELISA Kit
E-EL-P0027	Porcine CFP (Complement Factor P) ELISA Kit
E-EL-P0028	Porcine CFTR (Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit
E-EL-P0030	Porcine TDGF1 (Teratocarcinoma Derived Growth Factor 1) ELISA Kit
E-EL-P0031	Porcine PON3 (Paraoxonase 3) ELISA Kit
E-EL-P0032	Porcine MMP-17 (Matrix Metalloproteinase 17) ELISA Kit
E-EL-P0033	Porcine CPR (Cytochrome P450 Reductase) ELISA Kit
E-EL-P0034	Porcine GDF9 (Growth Differentiation Factor 9) ELISA Kit
E-EL-P0035	Porcine ACP3/PAP (Prostatic Acid Phosphatase 3) ELISA Kit
E-EL-P0036	Porcine CCKBR (Cholecystokinin B Receptor) ELISA Kit
E-EL-P0037	Porcine TGM1 (Transglutaminase 1, Keratinocyte) ELISA Kit
E-EL-P0038	Porcine PDE4D (Phosphodiesterase 4D, cAMP Specific) ELISA Kit
E-EL-P0039	Porcine NUP107 (Nucleoporin 107kDa) ELISA Kit
E-EL-P0040	Porcine SPC (Pulmonary Surfactant Associated Protein C) ELISA Kit
E-EL-P0041	Porcine PRB2 (Basic Salivary Proline Rich Protein 2) ELISA Kit
E-EL-P0042	Porcine BFP (Brain Finger Protein) ELISA Kit
E-EL-P0043	Porcine GLRX3 (Glutaredoxin 3) ELISA Kit
E-EL-P0044	Porcine FPR2 (Formyl Peptide Receptor 2) ELISA Kit
E-EL-P0048	Porcine LpPLA2 (Phospholipase A2, Lipoprotein Associated) ELISA Kit

E-EL-P0049	Porcine CNTF (Ciliary Neurotrophic Factor) ELISA Kit
E-EL-P0050	Porcine PDGF-AA (Platelet Derived Growth Factor-AA) ELISA Kit
E-EL-P0051	Porcine C/EBP β (CCAAT/Enhancer Binding Protein Beta) ELISA Kit
E-EL-P0052	Porcine ATF-2 (Activating Transcription Factor 2) ELISA Kit
E-EL-P0053	Porcine MUC17 (Mucin 17) ELISA Kit
E-EL-P0054	Porcine ERK3 (Extracellular Signal Regulated Kinase 3) ELISA Kit
E-EL-P0055	Porcine PLAU/uPA (Urokinase-Type Plasminogen Activator) ELISA Kit
E-EL-P0056	Porcine VWF (Von Willebrand Factor) ELISA Kit
E-EL-P0057	Porcine CG α (Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit
E-EL-P0058	Porcine GFAP (Glial Fibrillary Acidic Protein) ELISA Kit
E-EL-P0059	Porcine FGF-7/KGF (Fibroblast Growth Factor 7) ELISA Kit
E-EL-P0060	Porcine MUC13 (Mucin 13) ELISA Kit
E-EL-P0061	Porcine NT-ProBNP (N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
E-EL-P0063	Porcine LEP (Leptin) ELISA Kit
E-EL-P0065	Porcine PAC3 (Platelet Associated Complement 3) ELISA Kit
E-EL-P0066	Porcine PSA/KLK3 (Prostate Specific Antigen) ELISA Kit
E-EL-P0068	Porcine GAL1 (Galectin 1) ELISA Kit
E-EL-P0069	Porcine COL6 (Collagen Type VI) ELISA Kit
E-EL-P0070	Porcine CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) ELISA Kit
E-EL-P0071	Porcine NP (Nandrolone Phenylpropionate) ELISA Kit
E-EL-P0072	Porcine Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit

E-EL-P0073	Porcine F13 (Coagulation Factor X III) ELISA Kit
E-EL-P0074	Porcine CTAPIII (Connective Tissue Activating Peptide III) ELISA Kit
E-EL-P0075	Porcine GH (Growth Hormone) ELISA Kit
E-EL-P0076	Porcine PP13 (Placental Protein13) ELISA Kit
E-EL-P0077	Porcine ApoB100 (Apolipoprotein B100) ELISA Kit
E-EL-P0078	Porcine CRT (Calreticulin) ELISA Kit
E-EL-P0080	Porcine PRI (Placental Ribonuclease Inhibitor) ELISA Kit
E-EL-P0082	Porcine Methylase ELISA Kit
E-EL-P0084	Porcine NAGLU (N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit
E-EL-P0086	Porcine IL-19 (Interleukin 19) ELISA Kit
E-EL-P0087	Porcine PXR (Pregnane X Receptor) ELISA Kit
E-EL-P0088	Porcine C4 (Complement Component 4) ELISA Kit
E-EL-P0089	Porcine AGT (Angiotensinogen) ELISA Kit
E-EL-P0090	Porcine PS (Plasma Anticoagulant Protein S) ELISA Kit
E-EL-P0091	Porcine PIINP (Procollagen II N-Terminal ProPeptide) ELISA Kit
E-EL-P0092	Porcine FcεR II (Receptor II for the Fc Fragment of IgE) ELISA Kit
E-EL-P0093	Porcine HSF2 (Heat Shock Transcription Factor 2) ELISA Kit
E-EL-P0094	Porcine HSP-70/HSPA9 (Heat Shock Protein 70) ELISA Kit
E-EL-P0095	Porcine MYOG (Myogenin) ELISA Kit
E-EL-P0096	Porcine CR (Calretinin) ELISA Kit
E-EL-P0097	Porcine BK (Bradykinin) ELISA Kit

E-EL-P0099	Porcine NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) ELISA Kit
E-EL-P0100	Porcine TF (Tissue Factor) ELISA Kit
E-EL-P0101	Porcine SCT (Secretin) ELISA Kit
E-EL-P0102	Porcine ACTg1 (Actin Gamma 1) ELISA Kit
E-EL-P0103	Porcine CYP7A1 (Cytochrome P450 7A1) ELISA Kit
E-EL-P0104	Porcine C/EBP α (CCAAT/Enhancer Binding Protein Alpha) ELISA Kit
E-EL-P0105	Porcine PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit
E-EL-P0106	Porcine ADCY3 (Adenylate Cyclase 3) ELISA Kit
E-EL-P0107	Porcine PGD2S (Prostaglandin D2 Synthase, Brain) ELISA Kit
E-EL-P0108	Porcine MMP-9 (Matrix Metalloproteinase 9) ELISA Kit
E-EL-P0109	Porcine CAM-Ab (Anti-Calmodulin Specific Antibody) ELISA Kit
E-EL-P0110	Porcine MN (Metanephrine) ELISA Kit
E-EL-P0111	Porcine FTL (Ferritin, Light Polypeptide) ELISA Kit
E-EL-P0113	Porcine H-Cadherin (Cadherin, Heart) ELISA Kit
E-EL-P0114	Porcine OLAB (Anti-Oxidized Low Density Lipoprotein Antibody) ELISA Kit
E-EL-P0115	Porcine LATS (Long Acting Thyroid Stimulator) ELISA Kit
E-EL-P0116	Porcine PK (Pyruvate Kinase, Liver And RBC) ELISA Kit
E-EL-P0117	Porcine OTF2A (Octamer Transcription Factor 2A) ELISA Kit
E-EL-P0118	Porcine PMP3 (Peroxisomal Membrane Protein 3) ELISA Kit
E-EL-P0119	Porcine 15-LO (Arachidonate 15-Lipoxygenase) ELISA Kit
E-EL-P0121	Porcine BMP-15 (Bone Morphogenetic Protein 15) ELISA Kit

E-EL-P0122	Porcine DPP6 (Dipeptidyl Peptidase VI) ELISA Kit
E-EL-P0123	Porcine ATXN1 (Ataxin 1) ELISA Kit
E-EL-P0124	Porcine EAR3 (Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit
E-EL-P0125	Porcine CLDN (Claudin) ELISA Kit
E-EL-P0126	Porcine PDIA2 (Protein Disulfide Isomerase A2) ELISA Kit
E-EL-P0127	Porcine C1q (Complement 1q) ELISA Kit
E-EL-P0128	Porcine PGRP-S (Peptidoglycan Recognition Protein-S) ELISA Kit
E-EL-P0129	Porcine D2R/DRD2 (Dopamine Receptor D2) ELISA Kit
E-EL-P0130	Porcine SDF-2 (Stromal Cell Derived Factor 2) ELISA Kit
E-EL-P0131	Porcine CLA (Collagen Autoantibody) ELISA Kit
E-EL-P0132	Porcine Gas6 (Growth Arrest Specific Protein 6) ELISA Kit
E-EL-P0133	Porcine SIAE (Sialic Acid Acetyltransferase) ELISA Kit
E-EL-P0134	Porcine SLPI (Secretory Leukocyte Protease Inhibitor) ELISA Kit
E-EL-P0135	Porcine dDAVP (1-Desamino-8D-Arginine Vasopressin) ELISA Kit
E-EL-P0136	Porcine DAF (Decay Accelerating Factor) ELISA Kit
E-EL-P0137	Porcine ARSA (Arylsulfatase A) ELISA Kit
E-EL-P0138	Porcine CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
E-EL-P0139	Porcine GUK1 (Guanylate Kinase 1) ELISA Kit
E-EL-P0140	Porcine ADRP (Adipose Differentiation Related Protein) ELISA Kit
E-EL-P0141	Porcine HSP-27/HSPB1 (Heat Shock Protein 27) ELISA Kit
E-EL-P0142	Porcine PI3K2 α (Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit

E-EL-P0143	Porcine PTHR2 (Parathyroid Hormone Receptor 2) ELISA Kit
E-EL-P0144	Porcine GRN (Granulin) ELISA Kit
E-EL-P0145	Porcine GATA5 (GATA Binding Protein 5) ELISA Kit
E-EL-P0146	Porcine TH (Tyrosine Hydroxylase) ELISA Kit
E-EL-P0147	Porcine GMP-140 (Alpha-granular Membrane Protein) ELISA Kit
E-EL-P0148	Porcine EED (Embryonic Ectoderm Development) ELISA Kit
E-EL-P0149	Porcine CENPF (Centromere Protein F) ELISA Kit
E-EL-P0150	Porcine D4R/DRD4 (Dopamine Receptor D4) ELISA Kit
E-EL-P0151	Porcine CAMK2 α (Calcium/Calmodulin Dependent Protein Kinase II Alpha) ELISA Kit
E-EL-P0153	Porcine DSG-E1 (Desmogleins 1) ELISA Kit
E-EL-P0154	Porcine PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit
E-EL-P0155	Porcine MMP-12 (Matrix Metalloproteinase 12) ELISA Kit
E-EL-P0156	Porcine NEF3 (Neurofilament 3) ELISA Kit
E-EL-P0157	Porcine NGAL (Neutrophil Gelatinase Associated Lipocalin) ELISA Kit
E-EL-P0158	Porcine IL-17 (Interleukin 17) ELISA Kit
E-EL-P0161	Porcine MYL9 (Myosin Light Chain 9, Regulatory) ELISA Kit
E-EL-P0162	Porcine PK2 (Prokineticin 2) ELISA Kit
E-EL-P0164	Porcine C/EBP ϵ (CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit
E-EL-P0165	Porcine LDL (Low Density Lipoprotein) ELISA Kit
E-EL-P0168	Porcine CASP4 (Caspase 4) ELISA Kit
E-EL-P0169	Porcine RIPK2 (Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit

E-EL-P0171	Porcine PDGFA (Platelet Derived Growth Factor Subunit A) ELISA Kit
E-EL-P0172	Porcine TDT (Terminal Deoxynucleotidyl Transferase) ELISA Kit
E-EL-P0173	Porcine PLN (Phospholamban) ELISA Kit
E-EL-P0174	Porcine I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit
E-EL-P0175	Porcine CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit
E-EL-P0176	Porcine GIP(Gastric Inhibitory Polypeptide)ELISA Kit
E-EL-P0177	Porcine F2 (Coagulation Factor II) ELISA Kit
E-EL-P0179	Porcine ALDOA (Aldolase A, Fructose Bisphosphate) ELISA Kit
E-EL-P0180	Porcine LIPH (Lipase H) ELISA Kit
E-EL-P0181	Porcine FLN β (Filamin B, Beta) ELISA Kit
E-EL-P0182	Porcine Fc γ R II/CD32 (Receptor II for the Fc Fragment of IgG) ELISA Kit
E-EL-P0183	Porcine ANG (Angiostatin) ELISA Kit
E-EL-P0184	Porcine BDNF (Brain Derived Neurotrophic Factor) ELISA Kit
E-EL-P0185	Porcine ACTH (Adrenocorticotrophic Hormone) ELISA Kit
E-EL-P0187	Porcine Extein ELISA Kit
E-EL-P0188	Porcine UCRP (Ubiquitin Cross Reactive Protein) ELISA Kit
E-EL-P0189	Porcine UCP-2 (Uncoupling Protein 2, Mitochondrial) ELISA Kit
E-EL-P0190	Porcine C1INH (Complement 1 Inhibitor) ELISA Kit
E-EL-P0191	Porcine sFMC (Soluble Fibrin Monomer Complex) ELISA Kit
E-EL-P0192	Porcine ER β (Estrogen Receptor Beta) ELISA Kit
E-EL-P0193	Porcine SIGLEC1 (Sialic Acid Binding Ig Like Lectin 1) ELISA Kit

E-EL-P0194	Porcine C1r (Complement Component 1R) ELISA Kit
E-EL-P0195	Porcine UGCG (UDP-glucose ceramide glucosyltransferase) ELISA Kit
E-EL-P0196	Porcine BACE1 (Beta-site APP Cleaving Enzyme 1) ELISA Kit
E-EL-P0197	Porcine DAT (Dopamine Transporter) ELISA Kit
E-EL-P0198	Porcine EAR1 (Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit
E-EL-P0199	Porcine BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit
E-EL-P0200	Porcine DCLK1 (Doublecortin Like Kinase 1) ELISA Kit
E-EL-P0201	Porcine SOCS3 (Suppressors of Cytokine Signaling 3) ELISA Kit
E-EL-P0202	Porcine Smad7/MADH7 (Mothers Against Decapentaplegic Homolog 7) ELISA Kit
E-EL-P0203	Porcine SIGLEC10 (Sialic Acid Binding Ig Like Lectin 10) ELISA Kit
E-EL-P0204	Porcine TAFI (Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit
E-EL-P0205	Porcine FASN (Fatty Acid Synthase) ELISA Kit
E-EL-P0206	Porcine α -CTx (Alpha Crosslaps) ELISA Kit
E-EL-P0207	Porcine CK (Creatine Kinase) ELISA Kit
E-EL-P0208	Porcine SNCy (Synuclein Gamma) ELISA Kit
E-EL-P0209	Porcine EPYC (Epiphycan) ELISA Kit
E-EL-P0210	Porcine F12 (Coagulation Factor XII) ELISA Kit
E-EL-P0211	Porcine NT-ProANP (N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit
E-EL-P0212	Porcine NT-3 (Neurotrophin-3) ELISA Kit
E-EL-P0213	Porcine LFA-3 (Lymphocyte Function Associated Antigen 3) ELISA Kit
E-EL-P0214	Porcine RNF4 (Ring Finger Protein 4) ELISA Kit

E-EL-P0216	Porcine GOD (Glucose Oxidase) ELISA Kit
E-EL-P0217	Porcine CD59 (Protectin) ELISA Kit
E-EL-P0219	Porcine Gal-6S (Galactose-6-Sulphatase) ELISA Kit
E-EL-P0222	Porcine LTβR (Lymphotoxin Beta Receptor) ELISA Kit
E-EL-P0223	Porcine MT (Melatonin) ELISA Kit
E-EL-P0225	Porcine PEX2 (Peroxisomal Biogenesis Factor 2) ELISA Kit
E-EL-P0226	Porcine ANGPTL2 (Angiopoietin Like Protein 2) ELISA Kit
E-EL-P0227	Porcine TNNI1 (Troponin I Type 1, Slow Skeletal) ELISA Kit
E-EL-P0228	Porcine MMP-25 (Matrix Metalloproteinase 25) ELISA Kit
E-EL-P0230	Porcine GKRP (Glucokinase Regulatory Protein) ELISA Kit
E-EL-P0231	Porcine EL (EuglobulinLysis) ELISA Kit
E-EL-P0232	Porcine ACA-IgA (Anti-Cardiolipin Antibody IgA) ELISA Kit
E-EL-P0233	Porcine SELP (P-Selectin) ELISA Kit
E-EL-P0234	Porcine FAS/CD95/TNFRSF6 (Factor Related Apoptosis) ELISA Kit
E-EL-P0235	Porcine PGF2α (Prostaglandin F2α) ELISA Kit
E-EL-P0236	Porcine FcεR I (Receptor I for the Fc Fragment of IgE) ELISA Kit
E-EL-P0238	Porcine FBLN5 (Fibulin 5) ELISA Kit
E-EL-P0239	Porcine APSA (Anti-Phosphatidyl Serine Antibody) ELISA Kit
E-EL-P0240	Porcine D3R/DRD3 (Dopamine Receptor D3) ELISA Kit
E-EL-P0241	Porcine BTA (Bladder Tumor Antigen) ELISA Kit
E-EL-P0242	Porcine PYGB (Glycogen Phosphorylase, Brain) ELISA Kit

E-EL-P0244	Porcine CHD5 (Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
E-EL-P0245	Porcine TLR2 (Toll-Like Receptor 2) ELISA Kit
E-EL-P0246	Porcine MTF (Microtransferrin) ELISA Kit
E-EL-P0247	Porcine RI (Ribonuclease Inhibitor) ELISA Kit
E-EL-P0248	Porcine CaSR (Calcium Sensing Receptor) ELISA Kit
E-EL-P0249	Porcine MTX (Methotrexate) ELISA Kit
E-EL-P0250	Porcine α -SMA (α -Smooth muscle actin) ELISA Kit
E-EL-P0251	Porcine CD3d (Cluster of Differentiation 3d) ELISA Kit
E-EL-P0252	Porcine TNPO2 (Transportin 2) ELISA Kit
E-EL-P0253	Porcine Tn-I (Troponin I) ELISA Kit
E-EL-P0254	Porcine ATF-3 (Activating Transcription Factor 3) ELISA Kit
E-EL-P0255	Porcine PKD1 (Protein Kinase D1) ELISA Kit
E-EL-P0256	Porcine POMC (Pro-Opiomelanocortin) ELISA Kit
E-EL-P0258	Porcine GST μ 2/GST μ 2 (Glutathione S Transferase Mu 2, Muscle) ELISA Kit
E-EL-P0259	Porcine GTC (General Transcription Complex) ELISA Kit
E-EL-P0260	Porcine OPGL (Osteoprotegerin Ligand) ELISA Kit
E-EL-P0261	Porcine GPRC5A (G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit
E-EL-P0262	Porcine iPLA2 (Phospholipidase A2, Calcium Independent) ELISA Kit
E-EL-P0263	Porcine GJ β 1 (Gap Junction Protein Beta 1) ELISA Kit
E-EL-P0264	Porcine AR (Amphiregulin) ELISA Kit
E-EL-P0265	Porcine CgB (Chromogranin B) ELISA Kit

E-EL-P0266	Porcine Gp49a (Glycoprotein 49A) ELISA Kit
E-EL-P0267	Porcine ACT α 1 (Actin Alpha 1, Skeletal Muscle) ELISA Kit
E-EL-P0268	Porcine GPLD1 (Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Kit
E-EL-P0269	Porcine PMN Elastase (Polymorphonuclear Elastase) ELISA Kit
E-EL-P0270	Porcine PPAR- δ (Peroxisome Proliferator Activated Receptor Delta) ELISA Kit
E-EL-P0271	Porcine ADAM9 (A Disintegrin And Metalloprotease 9) ELISA Kit
E-EL-P0273	Porcine OVA sIgG (Ovalbumin Specific IgG) ELISA Kit
E-EL-P0274	Porcine BMF (Bcl-2 Modifying Factor) ELISA Kit
E-EL-P0276	Porcine VMA (Vanillylmandelic Acid) ELISA Kit
E-EL-P0278	Porcine SHBG (Sex Hormone Binding Globulin) ELISA Kit
E-EL-P0279	Porcine TGM2 (Transglutaminase 2, Tissue) ELISA Kit
E-EL-P0280	Porcine GPR120 (G Protein Coupled Receptor 120) ELISA Kit
E-EL-P0281	Porcine FBN2 (Fibrillin 2) ELISA Kit
E-EL-P0282	Porcine CLDN14 (Claudin 14) ELISA Kit
E-EL-P0284	Porcine ANGPTL4 (Angiopoietin Like Protein 4) ELISA Kit
E-EL-P0285	Porcine BPI (Bactericidal/Permeability Increasing Protein) ELISA Kit
E-EL-P0286	Porcine IgE (Immunoglobulin E) ELISA Kit
E-EL-P0288	Porcine ECM1 (Extracellular Matrix Protein 1) ELISA Kit
E-EL-P0289	Porcine MASP2 (Mannan Associated Serine Protease 2) ELISA Kit
E-EL-P0290	Porcine IL-15 (Interleukin 15) ELISA Kit
E-EL-P0291	Porcine LP-a (Lipoprotein a) ELISA Kit

E-EL-P0293	Porcine SCFR/c-Kit (Stem Cell Factor Receptor) ELISA Kit
E-EL-P0294	Porcine DCN (Decorin) ELISA Kit
E-EL-P0295	Porcine CLDN11 (Claudin 11) ELISA Kit
E-EL-P0296	Porcine DKK4 (Dickkopf Related Protein 4) ELISA Kit
E-EL-P0297	Porcine S100A8 (S100 Calcium Binding Protein A8/Calgranulin A) ELISA Kit
E-EL-P0298	Porcine Slit2 (Slit Homolog 2) ELISA Kit
E-EL-P0299	Porcine PTHrP (Parathyroid Hormone Related Protein) ELISA Kit
E-EL-P0300	Porcine Pg (Progesterone) ELISA Kit
E-EL-P0303	Porcine Collagenase I ELISA Kit
E-EL-P0304	Porcine ITIH4 (Inter Alpha-Globulin Inhibitor H4) ELISA Kit
E-EL-P0305	Porcine GDF3 (Growth Differentiation Factor 3) ELISA Kit
E-EL-P0306	Porcine GS (Gelsolin) ELISA Kit
E-EL-P0307	Porcine OBTF4 (Octamer Binding Transcription Factor 4) ELISA Kit
E-EL-P0308	Porcine CNTNAP2 (Contactin Associated Protein Like Protein 2) ELISA Kit
E-EL-P0309	Porcine INSL3 (Insulin Like Protein 3) ELISA Kit
E-EL-P0311	Porcine CCRL1 (Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
E-EL-P0312	Porcine α -GST (α -glutathione S-transferases) ELISA Kit
E-EL-P0313	Porcine CAPA (Cytotoxin Associated Protein A) ELISA Kit
E-EL-P0314	Porcine ANG II R1-Ab (Angiotensin II Receptor 1 Antibody) ELISA Kit
E-EL-P0315	Porcine GAS (N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
E-EL-P0316	Porcine PYY (Peptide YY) ELISA Kit

E-EL-P0317	Porcine TAF13 (TATA Box Binding Protein Associated Factor 13) ELISA Kit
E-EL-P0318	Porcine ASK-1 (Apoptosis Signal Regulating Kinase 1) ELISA Kit
E-EL-P0319	Porcine MIA1 (Melanoma Inhibitory Activity Protein 1) ELISA Kit
E-EL-P0320	Porcine Neopterin ELISA Kit
E-EL-P0321	Porcine SEMA3A (Semaphorin 3A) ELISA Kit
E-EL-P0323	Porcine VDR (Vitamin D Receptor) ELISA Kit
E-EL-P0324	Porcine SPD (Pulmonary Surfactant Associated Protein D) ELISA Kit
E-EL-P0325	Porcine LEPR (Leptin Receptor) ELISA Kit
E-EL-P0326	Porcine DEF α 5 (Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
E-EL-P0327	Porcine TE (Telomerase) ELISA Kit
E-EL-P0328	Porcine ORM2 (Orosomucoid 2) ELISA Kit
E-EL-P0329	Porcine VEGF-C (Vascular Endothelial Growth Factor C) ELISA Kit
E-EL-P0331	Porcine c-fos ELISA Kit
E-EL-P0333	Porcine VF (Visfatin) ELISA Kit
E-EL-P0334	Porcine ADIPOR1 (Adiponectin Receptor 1) ELISA Kit
E-EL-P0335	Porcine GTF2B (General Transcription Factor IIB) ELISA Kit
E-EL-P0336	Porcine SYK (Spleen Tyrosine Kinase) ELISA Kit
E-EL-P0337	Porcine OGP (Osteogenic Growth Peptide) ELISA Kit
E-EL-P0338	Porcine MMP-2 (Matrix Metalloproteinase 2) ELISA Kit
E-EL-P0339	Porcine GDF6 (Growth Differentiation Factor 6) ELISA Kit
E-EL-P0340	Porcine LRP6 (Low Density Lipoprotein Receptor Related Protein) ELISA Kit

E-EL-P0341	Porcine ApoM (Apolipoprotein M) ELISA Kit
E-EL-P0342	Porcine TR (Thrombin Receptor) ELISA Kit
E-EL-P0343	Porcine ASPN (Asporin) ELISA Kit
E-EL-P0344	Porcine UGT2B4 (Uridine Diphosphate Glucuronosyltransferase 2 Family Polypeptide B4) ELISA Kit
E-EL-P0345	Porcine SERT (Serotonin Transporter) ELISA Kit
E-EL-P0346	Porcine GLUT2 (Glucose Transporter 2) ELISA Kit
E-EL-P0347	Porcine RAPGEF1 (Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit
E-EL-P0348	Porcine I κ B β (Inhibitory Subunit of NF Kappa B Beta) ELISA Kit
E-EL-P0349	Porcine MAP-2 (Microtubule Associated Protein 2) ELISA Kit
E-EL-P0350	Porcine 6-K-PG (6-Keto-prostaglandin) ELISA Kit
E-EL-P0351	Porcine EPP (Erythrocyte Protoporphyrin) ELISA Kit
E-EL-P0352	Porcine NT-4 (Neurotrophin 4) ELISA Kit
E-EL-P0353	Porcine NUP205 (Nucleoporin 205kDa) ELISA Kit
E-EL-P0354	Porcine TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
E-EL-P0356	Porcine CNTN2 (Contactin 2) ELISA Kit
E-EL-P0357	Porcine ALD (Aldolase) ELISA Kit
E-EL-P0358	Porcine IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit
E-EL-P0360	Porcine AQP-2 (Aquaporin 2) ELISA Kit
E-EL-P0362	Porcine ACP1 (Acid Phosphatase 1) ELISA Kit
E-EL-P0363	Porcine EloA (Elongin A) ELISA Kit
E-EL-P0364	Porcine PRSS23 (Protease, Serine 23) ELISA Kit

E-EL-P0365	Porcine OBTF1 (Octamer Binding Transcription Factor 1) ELISA Kit
E-EL-P0366	Porcine CTSC (Cathepsin C) ELISA Kit
E-EL-P0367	Porcine NF (Neurofilament) ELISA Kit
E-EL-P0369	Porcine PFKM (Phosphofructokinase, Muscle) ELISA Kit
E-EL-P0370	Porcine Trx (Thioredoxin) ELISA Kit
E-EL-P0372	Porcine FTH (Ferritin, Heavy Polypeptide) ELISA Kit
E-EL-P0373	Porcine GPI (Glucose 6 Phosphate Isomerase) ELISA Kit
E-EL-P0374	Porcine MAP2K6/MKK6 (Mitogen Activated Protein Kinase Kinase 6) ELISA Kit
E-EL-P0375	Porcine CLCF1 (Cardiotrophin Like Cytokine Factor 1) ELISA Kit
E-EL-P0376	Porcine Ank-B (Ankyrin B) ELISA Kit
E-EL-P0377	Porcine CD97 (Cluster of Differentiation 97) ELISA Kit
E-EL-P0379	Porcine DPYSL2 (Dihydropyrimidinase Like Protein 2) ELISA Kit
E-EL-P0380	Porcine CACYBP (Calcyclin Binding Protein) ELISA Kit
E-EL-P0382	Porcine cADPRH/CD38 (Cyclic ADP Hydrolase) ELISA Kit
E-EL-P0383	Porcine ATF-1 (Activating Transcription Factor 1) ELISA Kit
E-EL-P0384	Porcine TACR2 (Tachykinin Receptor 2) ELISA Kit
E-EL-P0385	Porcine CR2/CD21 (Complement Receptor 2) ELISA Kit
E-EL-P0386	Porcine sCKR (Soluble Cytokine Receptor) ELISA Kit
E-EL-P0387	Porcine GPR37 (G Protein Coupled Receptor 37) ELISA Kit
E-EL-P0388	Porcine FGF19 (Fibroblast Growth Factor 19) ELISA Kit
E-EL-P0389	Porcine MCP-2/CCL8 (Monocyte Chemotactic Protein 2) ELISA Kit

E-EL-P0390	Porcine PAP (Plasmin-Antiplasmin Complex) ELISA Kit
E-EL-P0391	Porcine PLA2G4D (Phospholipase A2, Group IV D) ELISA Kit
E-EL-P0392	Porcine ACHAP (Acetylcholinesterase Associated Protein) ELISA Kit
E-EL-P0393	Porcine MBLR (Mannose Binding Lectin Receptor) ELISA Kit
E-EL-P0394	Porcine MAT (Matrilysin) ELISA Kit
E-EL-P0395	Porcine MYH1 (Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit
E-EL-P0397	Porcine KNG1 (Kininogen 1) ELISA Kit
E-EL-P0398	Porcine P-PKC (Phospho Protein Kinase C) ELISA Kit
E-EL-P0399	Porcine MCT (Mast Cell Tryptase) ELISA Kit
E-EL-P0400	Porcine CLS (Cardiolipin Synthase) ELISA Kit
E-EL-P0402	Porcine GzmA (Granzyme A) ELISA Kit
E-EL-P0403	Porcine MUC3 (Mucin 3) ELISA Kit
E-EL-P0404	Porcine PLA1 (Phospholipase A1) ELISA Kit
E-EL-P0405	Porcine SSFA2 (Sperm Specific Antigen 2) ELISA Kit
E-EL-P0406	Porcine LFA-2/CD2 (Lymphocyte Function Associated Antigen 2) ELISA Kit
E-EL-P0407	Porcine MAG Ab (Anti-Myelin Associated Glycoprotein Antibody) ELISA Kit
E-EL-P0408	Porcine MSLN (Mesothelin) ELISA Kit
E-EL-P0409	Porcine CYP1A2 (Cytochrome P450 1A2) ELISA Kit
E-EL-P0410	Porcine ATF6 (Activating Transcription Factor 6) ELISA Kit
E-EL-P0411	Porcine MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit
E-EL-P0412	Porcine REG1 β (Regenerating Islet Derived Protein 1 Beta) ELISA Kit

E-EL-P0413	Porcine TAP (Trypsinogen Activation Peptide) ELISA Kit
E-EL-P0414	Porcine ACA-IgG (Anti-Cardiolipin Antibody IgG) ELISA Kit
E-EL-P0415	Porcine AGA (Anti-Glycolipid Antibody) ELISA Kit
E-EL-P0416	Porcine PAH (Phenylalanine Hydroxylase) ELISA Kit
E-EL-P0417	Porcine DNM1 (Dynamin 1) ELISA Kit
E-EL-P0418	Porcine TAGLN (Transgelin) ELISA Kit
E-EL-P0419	Porcine ATGA/TGAB (Anti-Thyroid-Globulin Antibody) ELISA Kit
E-EL-P0421	Porcine FCN1 (Ficolin 1) ELISA Kit
E-EL-P0422	Porcine KAF (Keratinocyte Autocrine Factor) ELISA Kit
E-EL-P0423	Porcine β Manase ELISA Kit
E-EL-P0424	Porcine PKBG (Protein Kinase B Gamma) ELISA Kit
E-EL-P0425	Porcine MVP (Major Vault Protein) ELISA Kit
E-EL-P0426	Porcine VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
E-EL-P0427	Porcine ApoB (Apolipoprotein B) ELISA Kit
E-EL-P0428	Porcine ER α (Estrogen Receptor Alpha) ELISA Kit
E-EL-P0429	Porcine RNASEH (Ribonuclease H) ELISA Kit
E-EL-P0430	Porcine sgp130 (soluble Glucoprotein 130) ELISA Kit
E-EL-P0431	Porcine OSM (Oncostatin M) ELISA Kit
E-EL-P0432	Porcine LUM (Lumican) ELISA Kit
E-EL-P0433	Porcine MCRS1 (Microspherule Protein 1) ELISA Kit
E-EL-P0434	Porcine UPP2 (Uridine Phosphorylase 2) ELISA Kit

E-EL-P0435	Porcine SEMA3B (Semaphorin 3B) ELISA Kit
E-EL-P0436	Porcine FE (Ferritin) ELISA Kit
E-EL-P0437	Porcine snRNP/Sm (Small Nuclear Ribonucleoprotein) ELISA Kit
E-EL-P0438	Porcine CNX (Calnexin) ELISA Kit
E-EL-P0439	Porcine ARRβ1 (Arrestin Beta 1) ELISA Kit
E-EL-P0440	Porcine GAL4 (Galectin 4) ELISA Kit
E-EL-P0441	Porcine AVPR1A (Arginine Vasopressin Receptor 1A) ELISA Kit
E-EL-P0442	Porcine Bid (BH3 Interacting Domain Death Agonist) ELISA Kit
E-EL-P0443	Porcine C3c (Complement Fragment 3c) ELISA Kit
E-EL-P0444	Porcine TI-Ag (Thymus Independent Antigen) ELISA Kit
E-EL-P0446	Porcine plkBα (Phospho-Inhibitory Subunit Of NF Kappa B Alpha) ELISA Kit
E-EL-P0447	Porcine ERK (Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-P0448	Porcine PAPPA2 (Pappalysin 2) ELISA Kit
E-EL-P0450	Porcine PP (Protein Phosphatase) ELISA Kit
E-EL-P0451	Porcine CETP (Cholesteryl Ester Transfer Protein) ELISA Kit
E-EL-P0452	Porcine TACR3 (Tachykinin Receptor 3) ELISA Kit
E-EL-P0453	Porcine α Actin (Alpha-Actin) ELISA Kit
E-EL-P0454	Porcine PLP (Parathyroid Hormone-Like Protein) ELISA Kit
E-EL-P0455	Porcine IL-12 (Interleukin 12) ELISA Kit
E-EL-P0456	Porcine RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
E-EL-P0457	Porcine STAT2 (Signal Transducer and Activator of Transcription 2) ELISA Kit

E-EL-P0458	Porcine OCT (Ornithine Carbamoyl Transferase) ELISA Kit
E-EL-P0459	Porcine MAF (Macrophage Activating Factor) ELISA Kit
E-EL-P0460	Porcine APC (Activated Protein C) ELISA Kit
E-EL-P0461	Porcine S100B (S100 Calcium Binding Protein B) ELISA Kit
E-EL-P0462	Porcine TECK/CCL25 (Thymus Expressed Chemokine) ELISA Kit
E-EL-P0463	Porcine CCND2 (Cyclin-D2) ELISA Kit
E-EL-P0464	Porcine TSP-1 (Thrombospondin-1) ELISA Kit
E-EL-P0465	Porcine FSHR (Follicle Stimulating Hormone Receptor) ELISA Kit
E-EL-P0466	Porcine MOG (Myelin Oligodendrocyte Glycoprotein) ELISA Kit
E-EL-P0468	Porcine Dyn (Big Dynorphin) ELISA Kit
E-EL-P0469	Porcine DPP10 (Dipeptidyl Peptidase X) ELISA Kit
E-EL-P0470	Porcine CFI (Complement Factor I) ELISA Kit
E-EL-P0472	Porcine FOXP3 (Forkhead Box Protein P3) ELISA Kit
E-EL-P0473	Porcine RHOA (Ras Homolog Gene Family, Member A) ELISA Kit
E-EL-P0474	Porcine NUP160 (Nucleoporin 160kDa) ELISA Kit
E-EL-P0475	Porcine TAF1 (TATA Box Binding Protein Associated Factor 1) ELISA Kit
E-EL-P0476	Porcine PTGS1 (Prostaglandin Endoperoxide Synthase 1) ELISA Kit
E-EL-P0478	Porcine CX26 (Connexin 26) ELISA Kit
E-EL-P0479	Porcine CK-17/KRT17 (Cytokeratin 17) ELISA Kit
E-EL-P0481	Porcine TMPRSS2 (Transmembrane Protease, Serine 2) ELISA Kit
E-EL-P0482	Porcine PIR (Pirin) ELISA Kit

E-EL-P0483	Porcine GP (Anti-Glycoprotein Antibody) ELISA Kit
E-EL-P0484	Porcine F7 (Coagulation Factor VII) ELISA Kit
E-EL-P0485	Porcine TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit
E-EL-P0486	Porcine FPB (Fibrinopeptide B) ELISA Kit
E-EL-P0488	Porcine NAP3 (Neutrophil Activating Protein 3) ELISA Kit
E-EL-P0489	Porcine AQP-1 (Aquaporin 1) ELISA Kit
E-EL-P0491	Porcine EAR2 (Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit
E-EL-P0492	Porcine CASP11 (Caspase 11) ELISA Kit
E-EL-P0493	Porcine NUP98 (Nucleoporin 98kDa) ELISA Kit
E-EL-P0494	Porcine FRS2 (Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit
E-EL-P0495	Porcine ATF7 (Activating Transcription Factor 7) ELISA Kit
E-EL-P0496	Porcine SRP9 (Signal Recognition Particle 9kDa) ELISA Kit
E-EL-P0498	Porcine NUP35 (Nucleoporin 35kDa) ELISA Kit
E-EL-P0499	Porcine CSF2R α (Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit
E-EL-P0501	Porcine NOX5 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit
E-EL-P0502	Porcine AGA (Anti-Gliadin Antibody) ELISA Kit
E-EL-P0503	Porcine CABIN1 (Calcineurin Binding Protein 1) ELISA Kit
E-EL-P0504	Porcine WARS (Tryptophanyl tRNA Synthetase) ELISA Kit
E-EL-P0505	Porcine MCAM (Melanoma Cell Adhesion Molecule) ELISA Kit
E-EL-P0506	Porcine SPT β N4 (Spectrin Beta, Non Erythrocytic 4) ELISA Kit
E-EL-P0507	Porcine GAP (GTPase Activating Protein) ELISA Kit

E-EL-P0508	Porcine ELN (Elastin) ELISA Kit
E-EL-P0509	Porcine FOXO1 (Forkhead Box Protein O1) ELISA Kit
E-EL-P0510	Porcine DPYD (Dihydropyrimidine Dehydrogenase) ELISA Kit
E-EL-P0511	Porcine ANG (Angiogenin) ELISA Kit
E-EL-P0512	Porcine ATM (Ataxia Telangiectasia Mutated) ELISA Kit
E-EL-P0514	Porcine HSPG (Heparan Sulfate Proteoglycan) ELISA Kit
E-EL-P0515	Porcine GKN2 (Gastrokeine 2) ELISA Kit
E-EL-P0516	Porcine IGFBP-7 (Insulin Like Growth Factor Binding Protein 7) ELISA Kit
E-EL-P0517	Porcine RHOB (Ras Homolog Gene Family, Member B) ELISA Kit
E-EL-P0518	Porcine TSH (Thyroid Stimulating Hormone) ELISA Kit
E-EL-P0521	Porcine LR (Leukoregulin) ELISA Kit
E-EL-P0522	Porcine NFYC (Nuclear Transcription Factor Y Gamma) ELISA Kit
E-EL-P0523	Porcine PAFAHIB3 (Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit
E-EL-P0524	Porcine RPS6K β 1 (Ribosomal Protein S6 Kinase Beta 1) ELISA Kit
E-EL-P0525	Porcine FBN3 (Fibrillin 3) ELISA Kit
E-EL-P0527	Porcine GH2 (Growth Hormone 2) ELISA Kit
E-EL-P0528	Porcine MYD88 (Myeloid Differentiation Factor 88) ELISA Kit
E-EL-P0529	Porcine PTGES (Prostaglandin E Synthase, Microsomal) ELISA Kit
E-EL-P0530	Porcine CDK7 (Cyclin Dependent Kinase 7) ELISA Kit
E-EL-P0531	Porcine FGL2 (Fibrinogen Like Protein 2) ELISA Kit
E-EL-P0533	Porcine EPCR (Endothelial Protein C Receptor) ELISA Kit

E-EL-P0534	Porcine C3c (Complement C3 Convertase) ELISA Kit
E-EL-P0535	Porcine GC (Glucagon) ELISA Kit
E-EL-P0536	Porcine LZM (Lysozyme) ELISA Kit
E-EL-P0537	Porcine CNR-1 (Cadherin Related Neuronal Receptor1) ELISA Kit
E-EL-P0538	Porcine EPHA3 (Ephrin Type A Receptor 3) ELISA Kit
E-EL-P0539	Porcine PEDF (Pigment Epithelium Derived Factor) ELISA Kit
E-EL-P0540	Porcine PCX (Podocalyxin) ELISA Kit
E-EL-P0541	Porcine GAD (Glutamic Acid Decarboxylase) ELISA Kit
E-EL-P0542	Porcine DHEA-S7 (Dehydroepiandrosterone S7) ELISA Kit
E-EL-P0543	Porcine PTD (Pentosidine) ELISA Kit
E-EL-P0544	Porcine TNNI2 (Troponin I Type 2, Fast Skeletal) ELISA Kit
E-EL-P0545	Porcine CASP1 (Caspase 1) ELISA Kit
E-EL-P0546	Porcine DCX (Doublecortin) ELISA Kit
E-EL-P0547	Porcine CKB (Creatine Kinase, Brain) ELISA Kit
E-EL-P0548	Porcine CNTN4 (Contactin 4) ELISA Kit
E-EL-P0549	Porcine CA2 (Carbonic Anhydrase II) ELISA Kit
E-EL-P0550	Porcine KERA (Keratocan) ELISA Kit
E-EL-P0551	Porcine TU3A (TU3A Protein) ELISA Kit
E-EL-P0552	Porcine ANG-R-Tie2 (Angiopoietin Receptor Tie2) ELISA Kit
E-EL-P0553	Porcine DS (Dermatan Sulfate) ELISA Kit
E-EL-P0554	Porcine PAR1 (Protease Activated Receptor 1) ELISA Kit

E-EL-P0555	Porcine rT3 (Reverse Triiodothyronine) ELISA Kit
E-EL-P0556	Porcine PTPRN (Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
E-EL-P0557	Porcine C5 (Complement Component 5) ELISA Kit
E-EL-P0559	Porcine PON1 (Paraoxonase 1) ELISA Kit
E-EL-P0560	Porcine β -LPH (β -Lipotropic Hormone) ELISA Kit
E-EL-P0562	Porcine FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit
E-EL-P0563	Porcine MYC (V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit
E-EL-P0564	Porcine FMOD (Fibromodulin) ELISA Kit
E-EL-P0565	Porcine TS/TYMS (Thymidylate Synthetase) ELISA Kit
E-EL-P0566	Porcine PSAP (Prosaposin) ELISA Kit
E-EL-P0567	Porcine CD23/Fc ϵ RII (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
E-EL-P0568	Porcine PAIg (Platelet Associated Immunoglobulins) ELISA Kit
E-EL-P0569	Porcine RAG-2 (Recombination Activating Gene 2) ELISA Kit
E-EL-P0570	Porcine Ext-N (N-Terminal Extein) ELISA Kit
E-EL-P0571	Porcine cTn-I/TNNI3 (cardiac Troponin I) ELISA Kit
E-EL-P0572	Porcine FLN α (Filamin A, Alpha) ELISA Kit
E-EL-P0574	Porcine NT (Neurotensin) ELISA Kit
E-EL-P0575	Porcine GATA3 (GATA Binding Protein 3) ELISA Kit
E-EL-P0576	Porcine PLB (Phospholipase B) ELISA Kit
E-EL-P0577	Porcine GHAb (Anti-Growth Hormone Antibody) ELISA Kit
E-EL-P0578	Porcine GATA1 (GATA Binding Protein 1) ELISA Kit

E-EL-P0579	Porcine HPR (Haptoglobin Related Protein) ELISA Kit
E-EL-P0580	Porcine OXB (Orexin B) ELISA Kit
E-EL-P0581	Porcine CTACK (Cutaneous T-cell Attracting Chemokine) ELISA Kit
E-EL-P0583	Porcine AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit
E-EL-P0584	Porcine Ponticulin ELISA Kit
E-EL-P0585	Porcine ALK (Anaplastic Lymphoma Kinase) ELISA Kit
E-EL-P0586	Porcine TD-Ag (Thymus Dependent Antigen) ELISA Kit
E-EL-P0587	Porcine PKM2 (Pyruvate Kinase, Muscle) ELISA Kit
E-EL-P0588	Porcine TACE (TNF α Converting Enzyme) ELISA Kit
E-EL-P0589	Porcine MAO (Monoamine Oxidase) ELISA Kit
E-EL-P0590	Porcine MYL1 (Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit
E-EL-P0591	Porcine TFF1 (Trefoil Factor 1) ELISA Kit
E-EL-P0593	Porcine PAF (Platelet Activating Factor) ELISA Kit
E-EL-P0594	Porcine f-Hb (Free Haemoglobin) ELISA Kit
E-EL-P0596	Porcine TR (Targeted Ribonuclease) ELISA Kit
E-EL-P0597	Porcine PDGFR α (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit
E-EL-P0598	Porcine NCAD (Neural Cadherin) ELISA Kit
E-EL-P0601	Porcine CGRP (Calcitonin Gene Related Peptide) ELISA Kit
E-EL-P0602	Porcine GAL8 (Galectin 8) ELISA Kit
E-EL-P0603	Porcine CASP10 (Caspase 10) ELISA Kit
E-EL-P0604	Porcine PTN (Pleiotrophin) ELISA Kit

E-EL-P0605	Porcine IL-16 (Interleukin 16) ELISA Kit
E-EL-P0606	Porcine GSTpi (Glutathione S Transferases Pi) ELISA Kit
E-EL-P0607	Porcine pERK (Phospho Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-P0608	Porcine TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
E-EL-P0609	Porcine Cx (Connexin) ELISA Kit
E-EL-P0610	Porcine β 2-GP (Beta2 Glycoprotein) ELISA Kit
E-EL-P0611	Porcine α MSH (Alpha-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-P0612	Porcine CASP 9 (Caspase 9) ELISA Kit
E-EL-P0614	Porcine TTR (Transthyretin) ELISA Kit
E-EL-P0615	Porcine CALU (Calumenin) ELISA Kit
E-EL-P0616	Porcine FGF4 (Fibroblast Growth Factor 4) ELISA Kit
E-EL-P0618	Porcine AST (Aspartate Aminotransferase) ELISA Kit
E-EL-P0619	Porcine IGF-2 (Insulin Like Growth Factor 2) ELISA Kit
E-EL-P0620	Porcine MYH9 (Myosin Heavy Chain 9, Non Muscle) ELISA Kit
E-EL-P0621	Porcine MK (Midkine) ELISA Kit
E-EL-P0622	Porcine CYSLTR1 (Cysteinyl Leukotriene Receptor 1) ELISA Kit
E-EL-P0623	Porcine TSPAN30 (Tetraspanin 30) ELISA Kit
E-EL-P0624	Porcine RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
E-EL-P0625	Porcine GAL9 (Galectin 9) ELISA Kit
E-EL-P0626	Porcine SA (Sialic Acid) ELISA Kit
E-EL-P0627	Porcine OPN (Osteopontin) ELISA Kit

E-EL-P0628	Porcine CEP110 (Centrosomal Protein 110kDa) ELISA Kit
E-EL-P0629	Porcine PYGM (Glycogen Phosphorylase, Muscle) ELISA Kit
E-EL-P0630	Porcine PKD2 (Protein Kinase D2) ELISA Kit
E-EL-P0631	Porcine COL1α2 (Collagen Type I Alpha 2) ELISA Kit
E-EL-P0633	Porcine ESF (Erythropoiesis Stimulating Factor) ELISA Kit
E-EL-P0634	Porcine GKN1 (Gastrokine 1) ELISA Kit
E-EL-P0635	Porcine CKLFSF7 (Chemokine Like Factor Superfamily 7) ELISA Kit
E-EL-P0636	Porcine PDGFRβ (Platelet Derived Growth Factor Receptor Beta) ELISA Kit
E-EL-P0637	Porcine GCK (Glucokinase) ELISA Kit
E-EL-P0638	Porcine IGFBP-3 (Insulin Like Growth Factor Binding Protein 3) ELISA Kit
E-EL-P0639	Porcine PKIγ (Protein Kinase Inhibitor Gamma) ELISA Kit
E-EL-P0640	Porcine STAT1 (Signal Transducer and Activator of Transcription 1) ELISA Kit
E-EL-P0641	Porcine REG3α (Regenerating Islet Derived Protein 3 Alpha) ELISA Kit
E-EL-P0642	Porcine GHRH (Growth Hormone Releasing Hormone) ELISA Kit
E-EL-P0643	Porcine AChE (Acetylcholinesterase) ELISA Kit
E-EL-P0644	Porcine SMAF (Specific Macrophage Arming Factor) ELISA Kit
E-EL-P0645	Porcine Fcγ (Fc Fragment of IgG) ELISA Kit
E-EL-P0646	Porcine NCAM/CD56 (Neural Cell Adhesion Molecule) ELISA Kit
E-EL-P0647	Porcine MCSP (Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit
E-EL-P0648	Porcine DLK-1 (Delta Like 1 Homolog) ELISA Kit
E-EL-P0649	Porcine PLGA (Plasminogen Activator) ELISA Kit

E-EL-P0650	Porcine Ly6a (Lymphocyte Antigen 6 Complex Locus A) ELISA Kit
E-EL-P0651	Porcine BAG3 (Bcl2 Associated Athanogene 3) ELISA Kit
E-EL-P0652	Porcine PRELP (Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit
E-EL-P0653	Porcine LTC4 (Leukotriene C4) ELISA Kit
E-EL-P0654	Porcine EFNA5 (Ephrin A5) ELISA Kit
E-EL-P0656	Porcine MSR1 (Macrophage Scavenger Receptor 1) ELISA Kit
E-EL-P0657	Porcine UBE3A (Ubiquitin Protein Ligase E3A) ELISA Kit
E-EL-P0658	Porcine CEBPγ (CCAAT/Enhancer Binding Protein Gamma) ELISA Kit
E-EL-P0659	Porcine DUB (Deubiquitinating Enzyme) ELISA Kit
E-EL-P0660	Porcine PRDX5 (Peroxisiredoxin 5) ELISA Kit
E-EL-P0662	Porcine NUP133 (Nucleoporin 133kDa) ELISA Kit
E-EL-P0664	Porcine NAP (Neutrophil Alkaline Phosphatase) ELISA Kit
E-EL-P0665	Porcine ECF/CCL11 (Eosinophil Chemotactic Factor) ELISA Kit
E-EL-P0666	Porcine LPL (Lipoprotein Lipase) ELISA Kit
E-EL-P0667	Porcine TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
E-EL-P0668	Porcine NUP50 (Nucleoporin 50kDa) ELISA Kit
E-EL-P0669	Porcine Grb7 (Growth Factor Receptor Bound Protein 7) ELISA Kit
E-EL-P0670	Porcine PF3 (Platelet Factor 3) ELISA Kit
E-EL-P0671	Porcine COLEC11 (Collectin Sub-Family Member 11) ELISA Kit
E-EL-P0672	Porcine ABCG1 (ATP Binding Cassette Transporter G1) ELISA Kit
E-EL-P0673	Porcine GUSβ (β-glucuronidase) ELISA Kit

E-EL-P0674	Porcine PKP1 (Plakophilin 1) ELISA Kit
E-EL-P0675	Porcine GP-73 (Golgi Protein 73) ELISA Kit
E-EL-P0677	Porcine CHD3 (Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
E-EL-P0679	Porcine ITLN1 (Intelectin 1/Omentin) ELISA Kit
E-EL-P0680	Porcine SDF-1/CXCL12 (Stromal Cell Derived Factor 1) ELISA Kit
E-EL-P0681	Porcine ERN1 (Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit
E-EL-P0682	Porcine TBP2 (Thioredoxin Binding Protein 2) ELISA Kit
E-EL-P0683	Porcine PKP2 (Plakophilin 2) ELISA Kit
E-EL-P0684	Porcine MCF (Macrophage Chemotatic Factor) ELISA Kit
E-EL-P0685	Porcine AGGF1 (Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit
E-EL-P0686	Porcine GSTω1 (Glutathione S Transferase Omega 1) ELISA Kit
E-EL-P0687	Porcine FcγBP (Fc Fragment of IgG Binding Protein) ELISA Kit
E-EL-P0689	Porcine COL7 (Collagen Type VII) ELISA Kit
E-EL-P0690	Porcine PDGF-BB (Platelet Derived Growth Factor-BB) ELISA Kit
E-EL-P0691	Porcine TARC (Thymus Activation Regulated Chemokine) ELISA Kit
E-EL-P0693	Porcine TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
E-EL-P0694	Porcine Tgtp (T-cell Specific Gtpase) ELISA Kit
E-EL-P0695	Porcine EMA (Epithelial Membrane Antigen) ELISA Kit
E-EL-P0696	Porcine Des (Desmin) ELISA Kit
E-EL-P0697	Porcine LDHC (Lactate Dehydrogenase C) ELISA Kit
E-EL-P0698	Porcine TNK2 (Tyrosine Kinase, Non Receptor 2) ELISA Kit

E-EL-P0699	Porcine MMSAM (Melanoma Metastasis Surface Adhesion Molecule) ELISA Kit
E-EL-P0700	Porcine FABP5 (Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
E-EL-P0701	Porcine GIF (Gastric Intrinsic Factor) ELISA Kit
E-EL-P0703	Porcine DHEA-S (Dehydroepiandrosterone Sulfate) ELISA Kit
E-EL-P0704	Porcine LDHB (Lactate Dehydrogenase B) ELISA Kit
E-EL-P0705	Porcine D-LDH (D-Lactate Dehydrogenase) ELISA Kit
E-EL-P0706	Porcine ASB (Arylsulfatase B) ELISA Kit
E-EL-P0707	Porcine TPA (Tissue Polypeptide Antigen) ELISA Kit
E-EL-P0708	Porcine FcαR/CD89 (Receptor I for the Fc Fragment of IgA) ELISA Kit
E-EL-P0709	Porcine TMSβ4 (Thymosin Beta 4) ELISA Kit
E-EL-P0710	Porcine PTGES2 (Prostaglandin E Synthase 2) ELISA Kit
E-EL-P0712	Porcine ON (Osteonectin) ELISA Kit
E-EL-P0713	Porcine TUBα3C (Tubulin Alpha 3c) ELISA Kit
E-EL-P0714	Porcine CYSLTR2 (Cysteinyl Leukotriene Receptor 2) ELISA Kit
E-EL-P0715	Porcine FGFR4 (Fibroblast Growth Factor Receptor 4) ELISA Kit
E-EL-P0716	Porcine PROC (Protein C) ELISA Kit
E-EL-P0717	Porcine PICP (Procollagen I C-terminal Propeptide) ELISA Kit
E-EL-P0718	Porcine CRH (Corticotropin Releasing Hormone) ELISA Kit
E-EL-P0719	Porcine GHITM (Growth Hormone Inducible Transmembrane Protein) ELISA Kit
E-EL-P0720	Porcine SCTR (Secretin Receptor) ELISA Kit
E-EL-P0721	Porcine Grb10 (Growth Factor Receptor Bound Protein 10) ELISA Kit

E-EL-P0722	Porcine VPI (Vasoactive Peptidase Inhibitors) ELISA Kit
E-EL-P0723	Porcine TPO-Ab (Anti-Thyroid-Peroxidase Antibody) ELISA Kit
E-EL-P0724	Porcine HNRPA2B1 (Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit
E-EL-P0725	Porcine DLAT (Dihydrolipoyl Transacetylase) ELISA Kit
E-EL-P0726	Porcine PTPRB (Protein Tyrosine Phosphatase Receptor Type B) ELISA Kit
E-EL-P0727	Porcine IL-2R α (Interleukin 2 Receptor Alpha) ELISA Kit
E-EL-P0728	Porcine MlgD (Membrane IgD) ELISA Kit
E-EL-P0729	Porcine CCND1 (Cyclin-D1) ELISA Kit
E-EL-P0730	Porcine ANEA (Anti-Neutrophil Elastase Antibody) ELISA Kit
E-EL-P0731	Porcine IGFBP-1 (Insulin Like Growth Factor Binding Protein 1) ELISA Kit
E-EL-P0732	Porcine EFNA3 (Ephrin A3) ELISA Kit
E-EL-P0734	Porcine CASP3 (Caspase 3) ELISA Kit
E-EL-P0735	Porcine AQP-3 (Aquaporin 3) ELISA Kit
E-EL-P0736	Porcine PIICP (Procollagen II C-Terminal ProPeptide) ELISA Kit
E-EL-P0737	Porcine USP8 (Ubiquitin Specific Peptidase 8) ELISA Kit
E-EL-P0738	Porcine COL13 (Collagen Type XIII) ELISA Kit
E-EL-P0739	Porcine SPB (Pulmonary Surfactant Associated Protein B) ELISA Kit
E-EL-P0742	Porcine CDX2 (Caudal Type Homeobox Transcription Factor 2) ELISA Kit
E-EL-P0744	Porcine bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit
E-EL-P0745	Porcine c-jun ELISA Kit
E-EL-P0746	Porcine 6-OHDA (6-Hydroxydopamine) ELISA Kit

E-EL-P0747	Porcine CXCR3 (Chemokine CXC Receptor 3) ELISA Kit
E-EL-P0748	Porcine PDGFD (Platelet Derived Growth Factor D) ELISA Kit
E-EL-P0749	Porcine EMMPRIN/CD147 (Extracellular Matrix Metalloproteinase Inducer) ELISA Kit
E-EL-P0750	Porcine GRO γ (Growth Regulated Oncogene Gamma) ELISA Kit
E-EL-P0751	Porcine UFC (Urinary Free Cortisol) ELISA Kit
E-EL-P0753	Porcine PLC β 1/PI-PLC (Phospholipase C Beta 1, Phosphoinositide Specific) ELISA Kit
E-EL-P0754	Porcine TERT (Telomerase Reverse Transcriptase) ELISA Kit
E-EL-P0755	Porcine PDI (Protein Disulfide Isomerase) ELISA Kit
E-EL-P0758	Porcine MGST1 (Microsomal Glutathione S Transferase 1) ELISA Kit
E-EL-P0759	Porcine CALD (Caldesmon) ELISA Kit
E-EL-P0760	Porcine NRG-1 (Neuregulin 1) ELISA Kit
E-EL-P0761	Porcine GDF5 (Growth Differentiation Factor 5) ELISA Kit
E-EL-P0762	Porcine NOX1 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
E-EL-P0764	Porcine NPYR5 (Neuropeptide Y Receptor Y5) ELISA Kit
E-EL-P0765	Porcine LOXL4 (Lysyl Oxidase Like Protein 4) ELISA Kit
E-EL-P0766	Porcine SDF-4 (Stromal Cell Derived Factor 4) ELISA Kit
E-EL-P0769	Porcine CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit
E-EL-P0770	Porcine TTN (Titin) ELISA Kit
E-EL-P0772	Porcine GATA4 (GATA Binding Protein 4) ELISA Kit
E-EL-P0773	Porcine SEPCR/PROCR (Soluble Endothelial Protein C Receptor) ELISA Kit
E-EL-P0775	Porcine CCA (Colon Cancer Antigen) ELISA Kit

E-EL-P0776	Porcine CYP21B (Cytochrome P450 21B) ELISA Kit
E-EL-P0777	Porcine α HSP (Alpha-Hemoglobin Stabilizing Protein) ELISA Kit
E-EL-P0778	Porcine ApoC4 (Apolipoprotein C4) ELISA Kit
E-EL-P0779	Porcine BLC (B-Lymphocyte Chemoattractant) ELISA Kit
E-EL-P0780	Porcine GSANA (Granulocyte Specific Anti-Nuclear Antibody) ELISA Kit
E-EL-P0781	Porcine IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit
E-EL-P0782	Porcine NGB (Neuroglobin) ELISA Kit
E-EL-P0783	Porcine BMP-1 (Bone Morphogenetic Protein 1) ELISA Kit
E-EL-P0784	Porcine MLCK/MYLK (Myosin Light Chain Kinase) ELISA Kit
E-EL-P0785	Porcine CFH (Complement Factor H) ELISA Kit
E-EL-P0786	Porcine PLD2 (Phospholipase D2) ELISA Kit
E-EL-P0787	Porcine PTF1 α (Pancreas Specific Transcription Factor 1 α) ELISA Kit
E-EL-P0788	Porcine AP (Aprotinin) ELISA Kit
E-EL-P0789	Porcine GPC1 (Glypican 1) ELISA Kit
E-EL-P0790	Porcine SEP1 (Selenoprotein 1) ELISA Kit
E-EL-P0791	Porcine AXIN2 (Axis Inhibition Protein 2) ELISA Kit
E-EL-P0792	Porcine DMD (Dystrophin) ELISA Kit
E-EL-P0793	Porcine FGFBP1 (Fibroblast Growth Factor Binding Protein 1) ELISA Kit
E-EL-P0794	Porcine Anti-AlaRS (Anti-Alanyl-tRNA Synthetase/Anti-PL12-Antibody) ELISA Kit
E-EL-P0795	Porcine CLDN4 (Claudin 4) ELISA Kit
E-EL-P0796	Porcine C4BP α (C4 Binding Protein Alpha) ELISA Kit

E-EL-P0797	Porcine CTR (Calcitonin Receptor) ELISA Kit
E-EL-P0799	Porcine CTNN β 1 (Catenin, Beta 1) ELISA Kit
E-EL-P0800	Porcine CLL1 (Collectin Liver 1) ELISA Kit
E-EL-P0801	Porcine GL α (Galactosidase Alpha) ELISA Kit
E-EL-P0802	Porcine PROZ (Protein Z) ELISA Kit
E-EL-P0804	Porcine β -EP (Beta-Endorphin) ELISA Kit
E-EL-P0808	Porcine CKMT1B (Creatine Kinase, Mitochondrial 1B) ELISA Kit
E-EL-P0810	Porcine PIK3C α (Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit
E-EL-P0811	Porcine APP (Amyloid Precursor Protein) ELISA Kit
E-EL-P0812	Porcine SPARCL1 (SPARC Like Protein 1) ELISA Kit
E-EL-P0813	Porcine PAEP (Progestagen-Associated Endometrial Protein) ELISA Kit
E-EL-P0814	Porcine NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
E-EL-P0815	Porcine BAFF/CD257 (B-Cell Activating Factor) ELISA Kit
E-EL-P0816	Porcine FSCN (Fascin) ELISA Kit
E-EL-P0817	Porcine TXB2 (Thromboxane B2) ELISA Kit
E-EL-P0818	Porcine PCDH15 (Protocadherin 15) ELISA Kit
E-EL-P0819	Porcine SELE (E-Selectin) ELISA Kit
E-EL-P0820	Porcine MUC5AC (Mucin 5 Subtype AC) ELISA Kit
E-EL-P0821	Porcine BLNK (B-Cell Linker Protein) ELISA Kit
E-EL-P0822	Porcine TNKS2 (Tankyrase 2) ELISA Kit
E-EL-P0823	Porcine NPYR2 (Neuropeptide Y Receptor Y2) ELISA Kit

E-EL-P0824	Porcine CDK4 (Cyclin Dependent Kinase 4) ELISA Kit
E-EL-P0825	Porcine NNE (Non-Neuronal Enolase) ELISA Kit
E-EL-P0826	Porcine COMP (Cartilage Oligomeric Matrix Protein) ELISA Kit
E-EL-P0827	Porcine BAK1 (Bcl2 Antagonist/Killer 1) ELISA Kit
E-EL-P0828	Porcine PGF (Placental Growth Factor) ELISA Kit
E-EL-P0830	Porcine SEPP1 (Selenoprotein P) ELISA Kit
E-EL-P0831	Porcine FLOT2 (Flotillin 2) ELISA Kit
E-EL-P0833	Porcine MAN1A1 (Mannosidase Alpha Class 1A Member 1) ELISA Kit
E-EL-P0834	Porcine ATA (Angiotensinase A) ELISA Kit
E-EL-P0835	Porcine CYFRA21-1 (Cytokeratin Fragment Antigen 21-1) ELISA Kit
E-EL-P0837	Porcine GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit
E-EL-P0838	Porcine PDGF (Platelet-Derived Growth Factor) ELISA Kit
E-EL-P0839	Porcine BANP (BTG3 Associated Nuclear Protein) ELISA Kit
E-EL-P0841	Porcine SRF/IF (Skin Reactive Factor/Inflammatory Factor) ELISA Kit
E-EL-P0843	Porcine LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit
E-EL-P0844	Porcine CT-1 (Cardiotrophin 1) ELISA Kit
E-EL-P0845	Porcine MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit
E-EL-P0846	Porcine TFPI (Tissue Factor Pathway Inhibitor) ELISA Kit
E-EL-P0847	Porcine HA (Hemagglutinin) ELISA Kit
E-EL-P0848	Porcine EphA2 (Tyrosine protein kinase receptor A2) ELISA Kit
E-EL-P0849	Porcine XCR1 (XC-Chemokine Receptor 1) ELISA Kit

E-EL-P0851	Porcine AGPA (Anti-Glucoprotein 210 Antibody) ELISA Kit
E-EL-P0852	Porcine PCDHyA2 (Protocadherin Gamma A2) ELISA Kit
E-EL-P0853	Porcine IGF-1(Insulin Like Growth Factor 1)ELISA Kit
E-EL-P0855	Porcine AP13 (Apelin 13) ELISA Kit
E-EL-P0856	Porcine SPJ (Surfactant Associated Protein J) ELISA Kit
E-EL-P0857	Porcine ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit
E-EL-P0858	Porcine ACD (Adrenocortical Dysplasia Homolog) ELISA Kit
E-EL-P0859	Porcine FGF23 (Fibroblast Growth Factor 23) ELISA Kit
E-EL-P0860	Porcine THR α (Thyroid Hormone Receptor Alpha) ELISA Kit
E-EL-P0861	Porcine S100A9 (S100 Calcium Binding Protein A9) ELISA Kit
E-EL-P0862	Porcine Try2 (Trypsinogen II) ELISA Kit
E-EL-P0863	Porcine C5c (Complement C5 Convertase) ELISA Kit
E-EL-P0864	Porcine SNRPD1 (Small Nuclear Ribonucleoprotein Polypeptide D1) ELISA Kit
E-EL-P0865	Porcine OSF (Osteoclast Stimulating Factor 1) ELISA Kit
E-EL-P0866	Porcine GZMK (Granzyme K) ELISA Kit
E-EL-P0867	Porcine MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
E-EL-P0868	Porcine CRN (Corin) ELISA Kit
E-EL-P0869	Porcine D5R/DRD5 (Dopamine Receptor D5) ELISA Kit
E-EL-P0871	Porcine PI3K (Phosphoinositide-3-Kinase) ELISA Kit
E-EL-P0872	Porcine C4b (Complement Component 4b) ELISA Kit
E-EL-P0873	Porcine TF (Transfer Factor) ELISA Kit

E-EL-P0874	Porcine SCF (Stem Cell Factor) ELISA Kit
E-EL-P0876	Porcine GLUT3 (Glucose Transporter 3) ELISA Kit
E-EL-P0877	Porcine ANGPT4 (Angiopoietin 4) ELISA Kit
E-EL-P0878	Porcine EXT1 (Exostoses 1) ELISA Kit
E-EL-P0879	Porcine FRS3 (Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
E-EL-P0882	Porcine EG-VEGF (Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit
E-EL-P0884	Porcine CKMT1A (Creatine Kinase, Mitochondrial 1A) ELISA Kit
E-EL-P0885	Porcine ET-1 (Endothelin 1) ELISA Kit
E-EL-P0887	Porcine CKLFSF5 (Chemokine Like Factor Superfamily 5) ELISA Kit
E-EL-P0888	Porcine HYOU1 (Hypoxia Up Regulated 1) ELISA Kit
E-EL-P0889	Porcine P III CP (Procollagen III C-Terminal ProPeptide) ELISA Kit
E-EL-P0890	Porcine PCPE1 (Procollagen C Proteinase Enhancer 1) ELISA Kit
E-EL-P0891	Porcine CA153 (Carbohydrate Antigen 153) ELISA Kit
E-EL-P0892	Porcine PRSS8 (Protease, Serine 8) ELISA Kit
E-EL-P0893	Porcine PDK4 (Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit
E-EL-P0894	Porcine BCAN (Brevican) ELISA Kit
E-EL-P0895	Porcine TM (Thrombin) ELISA Kit
E-EL-P0896	Porcine PACAPR (Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit
E-EL-P0897	Porcine KLK6 (Kallikrein 6) ELISA Kit
E-EL-P0898	Porcine TM (Thrombomodulin) ELISA Kit
E-EL-P0899	Porcine FAIM3 (Fas Apoptotic Inhibitory Molecule 3) ELISA Kit

E-EL-P0900	Porcine CLA/CTAGE (Cutaneous Lymphocyte Associated Antigen) ELISA Kit
E-EL-P0901	Porcine GAL7 (Galectin 7) ELISA Kit
E-EL-P0902	Porcine NTXI (Cross Linked N-telopeptide of Type I Collagen) ELISA Kit
E-EL-P0903	Porcine GPC3 (Glypican 3) ELISA Kit
E-EL-P0904	Porcine PPOX (Protoporphyrinogen Oxidase) ELISA Kit
E-EL-P0905	Porcine ANG2 (Angiopoietin 2) ELISA Kit
E-EL-P0906	Porcine ABCG2 (ATP Binding Cassette Transporter G2) ELISA Kit
E-EL-P0907	Porcine ADAMTS4 (ADAM with Thrombospondin Type 1 Motif4) ELISA Kit
E-EL-P0908	Porcine LILRB4 (Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit
E-EL-P0909	Porcine PROS (Protein S) ELISA Kit
E-EL-P0910	Porcine TGF- β 2 (Transforming Growth Factor β 2) ELISA Kit
E-EL-P0911	Porcine EK (Eepidermal Keratin) ELISA Kit
E-EL-P0912	Porcine TP/TXA2R (Thromboxane Receptor) ELISA Kit
E-EL-P0913	Porcine STAT6B (Signal Transducer and Activator of Transcription 6B) ELISA Kit
E-EL-P0914	Porcine PTX 3/TSG-14 (Pentraxin 3, Long) ELISA Kit
E-EL-P0915	Porcine PFKP (Phosphofructokinase, Platelet) ELISA Kit
E-EL-P0916	Porcine GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit
E-EL-P0917	Porcine TAB (Thyroxine Antibody) ELISA Kit
E-EL-P0918	Porcine ARO (Aromatase) ELISA Kit
E-EL-P0919	Porcine CT (Calcitonin) ELISA Kit
E-EL-P0920	Porcine UPP1 (Uridine Phosphorylase 1) ELISA Kit

E-EL-P0921	Porcine Grb14 (Growth Factor Receptor Bound Protein 14) ELISA Kit
E-EL-P0922	Porcine PDH β (Pyruvate Dehydrogenase Beta) ELISA Kit
E-EL-P0923	Porcine FBLN4 (Fibulin 4) ELISA Kit
E-EL-P0924	Porcine LCN1 (Lipocalin 1) ELISA Kit
E-EL-P0925	Porcine NMP-22 (Nuclear Matrix Protein 22) ELISA Kit
E-EL-P0926	Porcine CRY α B (Crystallin Alpha B) ELISA Kit
E-EL-P0927	Porcine CASP14 (Caspase 14) ELISA Kit
E-EL-P0928	Porcine RXFP1 (Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
E-EL-P0931	Porcine GPIIb/IIIa (Platelet Membrane Glycoprotein IIb/IIIa) ELISA Kit
E-EL-P0932	Porcine PIGR/SC (Polymeric Immunoglobulin Receptor/Membrane Secretory Component) ELISA Kit
E-EL-P0934	Porcine CTNN α 1 (Catenin, Alpha 1) ELISA Kit
E-EL-P0935	Porcine ECFR (Eosinophil Chemotactic Factor Receptor) ELISA Kit
E-EL-P0936	Porcine MEC/CCL28 (Mucosae Associated Epithelia Chemokine) ELISA Kit
E-EL-P0938	Porcine LIFR (Leukemia Inhibitory Factor Receptor) ELISA Kit
E-EL-P0940	Porcine CLC (Charcot Leyden Crystal Protein) ELISA Kit
E-EL-P0941	Porcine MAG (Myelin Associated Glycoprotein) ELISA Kit
E-EL-P0942	Porcine 11-DH-TXB2 (11-Dehydro-thromboxane B2) ELISA Kit
E-EL-P0943	Porcine AFA (Anti-Filaggrin Antibody) ELISA Kit
E-EL-P0947	Porcine Flt3 (FMS Like Tyrosine Kinase 3) ELISA Kit
E-EL-P0948	Porcine PCP (Procollagen C Peptidase) ELISA Kit
E-EL-P0949	Porcine α HBDH (α -Hydroxybutyrate Dehydrogenase) ELISA Kit

E-EL-P0950	Porcine M-AChR M2 (Muscarinic Acetylcholine Receptor M2) ELISA Kit
E-EL-P0951	Porcine CA5B (Carbonic Anhydrase VB, Mitochondrial) ELISA Kit
E-EL-P0952	Porcine AP62A (Anti-P62 Antibody) ELISA Kit
E-EL-P0953	Porcine CENPI (Centromere Protein I) ELISA Kit
E-EL-P0954	Porcine RPL10 (Ribosomal Protein L10) ELISA Kit
E-EL-P0955	Porcine G6PD (Glucose 6 Phosphate Dehydrogenase) ELISA Kit
E-EL-P0956	Porcine MTHFR (Methylenetetrahydrofolate Reductase) ELISA Kit
E-EL-P0957	Porcine TFF2 (Trefoil Factor 2) ELISA Kit
E-EL-P0959	Porcine SOX3 (Sex Determining Region Y Box Protein 3) ELISA Kit
E-EL-P0960	Porcine C2 (Complement Component 2) ELISA Kit
E-EL-P0961	Porcine FcγR3A (Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit
E-EL-P0962	Porcine PPAR-α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
E-EL-P0964	Porcine STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit
E-EL-P0965	Porcine NFKBIE (Nuclear Factor Of Kappa Light Polypeptide Gene Enhancer In B-cells Inhibitor Epsilon) ELISA Kit
E-EL-P0967	Porcine UreC (Urease Related Protein C) ELISA Kit
E-EL-P0968	Porcine N-MID-OT (N-MID Osteocalcin) ELISA Kit
E-EL-P0969	Porcine GDF15 (Growth Differentiation Factor 15) ELISA Kit
E-EL-P0970	Porcine SGK2 (Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit
E-EL-P0971	Porcine CA19-9 (Carbohydrate antigen19-9) ELISA Kit
E-EL-P0972	Porcine F-Box Protein 2(FBXO2) ELISA Kit
E-EL-P0973	Porcine TIEG1 (TGF Beta Inducible Early Response Gene 1) ELISA Kit

E-EL-P0974	Porcine DBP (Vitamin D Binding Protein) ELISA Kit
E-EL-P0975	Porcine FDP (Fibrinogen Degradation Product) ELISA Kit
E-EL-P0977	Porcine TAF2 (TATA Box Binding Protein Associated Factor 2) ELISA Kit
E-EL-P0978	Porcine HO2 (Heme Oxygenase 2) ELISA Kit
E-EL-P0979	Porcine GDF11 (Growth Differentiation Factor 11) ELISA Kit
E-EL-P0980	Porcine CaN (Calcineurin) ELISA Kit
E-EL-P0981	Porcine TMSβ10 (Thymosin Beta 10) ELISA Kit
E-EL-P0982	Porcine BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit
E-EL-P0983	Porcine LAT (Linker For Activation of T-cell) ELISA Kit
E-EL-P0984	Porcine SAP (Serum Amyloid P Component) ELISA Kit
E-EL-P0985	Porcine TSC22 (Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
E-EL-P0986	Porcine FOXP1 (Forkhead Box Protein P1) ELISA Kit
E-EL-P0987	Porcine ApoA1 (Apolipoprotein A1) ELISA Kit
E-EL-P0988	Porcine TNR (Tenascin R) ELISA Kit
E-EL-P0989	Porcine EPHA1 (Ephrin Type A Receptor 1) ELISA Kit
E-EL-P0990	Porcine RAG-1 (Recombination Activating Gene 1) ELISA Kit
E-EL-P0991	Porcine MCP/CD46 (Membrane Cofactor Protein) ELISA Kit
E-EL-P0992	Porcine GLP-1(Glucagon Like Peptide 1)ELISA Kit
E-EL-P0993	Porcine NRG-3 (Neuregulin 3) ELISA Kit
E-EL-P0994	Porcine TP53 (Tumor Protein p53) ELISA Kit
E-EL-P0995	Porcine NPY (Neuropeptide Y) ELISA Kit

E-EL-P0996	Porcine TPP1 (Tripeptidyl Peptidase I) ELISA Kit
E-EL-P0998	Porcine ATF-5 (Activating Transcription Factor 5) ELISA Kit
E-EL-P0999	Porcine IFN- α (Interferon Alpha) ELISA Kit
E-EL-P1000	Porcine UCHL1 (Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
E-EL-P1001	Porcine PPIL1 (Peptidylprolyl Isomerase/Cyclophilin Like 1) ELISA Kit
E-EL-P1002	Porcine PG (Proteoglycan) ELISA Kit
E-EL-P1003	Porcine FBN1 (Fibrillin 1) ELISA Kit
E-EL-P1004	Porcine AhR (Aryl Hydrocarbon Receptor) ELISA Kit
E-EL-P1005	Porcine IDE (Insulin Degrading Enzyme) ELISA Kit
E-EL-P1006	Porcine LIPC (Lipase, Hepatic) ELISA Kit
E-EL-P1007	Porcine C3b (Complement Fragment 3b) ELISA Kit
E-EL-P1008	Porcine HB β (Hemoglobin Beta) ELISA Kit
E-EL-P1009	Porcine α 2-PI (α 2-Plasmin Inhibitor) ELISA Kit
E-EL-P1010	Porcine ANG1 (Angiopoietin 1) ELISA Kit
E-EL-P1011	Porcine ERCC1 (Excision Repair Cross Complementing Rodent Repair Deficiency Complementation 1) ELISA Kit
E-EL-P1012	Porcine CLDN3 (Claudin 3) ELISA Kit
E-EL-P1015	Porcine GCSFR (Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
E-EL-P1017	Porcine PTPRJ (Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit
E-EL-P1018	Porcine TEF (Thyrotroph Embryonic Factor) ELISA Kit
E-EL-P1019	Porcine F XII a (Activated Coagulation Factor XII) ELISA Kit
E-EL-P1022	Porcine GR α (Glucocorticoid Receptor Alpha) ELISA Kit

E-EL-P1023	Porcine PLSCR1 (Phospholipid Scramblase 1) ELISA Kit
E-EL-P1024	Porcine MRC1 (Mannose Receptor C Type 1 Like Protein 1) ELISA Kit
E-EL-P1025	Porcine PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit
E-EL-P1026	Porcine CA724 (Tumor Marker CA724) ELISA Kit
E-EL-P1028	Porcine FOXC1 (Forkhead Box Protein C1) ELISA Kit
E-EL-P1029	Porcine TFPI2 (Tissue Factor Pathway Inhibitor 2) ELISA Kit
E-EL-P1030	Porcine MART/Melan-A (Melanoma Marker A) ELISA Kit
E-EL-P1031	Porcine PP (Pancreatic Polypeptide) ELISA Kit
E-EL-P1032	Porcine PTPRO (Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
E-EL-P1034	Porcine α 2-AP (α 2-Antiplasmin) ELISA Kit
E-EL-P1035	Porcine COL9 α 3 (Collagen Type IX Alpha 3) ELISA Kit
E-EL-P1036	Porcine IL-6R (Interleukin 6 Receptor) ELISA Kit
E-EL-P1037	Porcine MIP-3 β /ELC/CCL19 (Macrophage Inflammatory Protein 3 β) ELISA Kit
E-EL-P1038	Porcine ACE (Angiotensin I Converting Enzyme) ELISA Kit
E-EL-P1039	Porcine HLA-A (Leukocyte Antigen A) ELISA Kit
E-EL-P1040	Porcine CYTH1 (Cytohesin 1) ELISA Kit
E-EL-P1041	Porcine INSR (Insulin Receptor) ELISA Kit
E-EL-P1042	Porcine GP39 (Glycoprotein 39, Cartilage) ELISA Kit
E-EL-P1043	Porcine OSMR (Oncostatin M Receptor) ELISA Kit
E-EL-P1045	Porcine RND1 (Rho Family GTPase 1) ELISA Kit
E-EL-P1046	Porcine SSRP1 (Structure Specific Recognition Protein 1) ELISA Kit

E-EL-P1047	Porcine PAR3 (Protease Activated Receptor 3) ELISA Kit
E-EL-P1049	Porcine Bax (Bcl-2 Associated X Protein) ELISA Kit
E-EL-P1051	Porcine NUP155 (Nucleoporin 155kDa) ELISA Kit
E-EL-P1053	Porcine TLA (Thymus Leukemia Antigen) ELISA Kit
E-EL-P1054	Porcine MCP-4/CCL13 (Monocyte Chemotactic Protein 4) ELISA Kit
E-EL-P1055	Porcine COL6α3 (Collagen Type VI Alpha 3) ELISA Kit
E-EL-P1056	Porcine 5-NT (5-Nucleotidase) ELISA Kit
E-EL-P1057	Porcine FG (Fibrinogen) ELISA Kit
E-EL-P1058	Porcine MIP-5 (Macrophage Inflammatory Protein 5) ELISA Kit
E-EL-P1060	Porcine βTG/PBP/CXCL7/NAP2 (Thromboglobulin, Beta) ELISA Kit
E-EL-P1061	Porcine FABP4 (Fatty Acid-binding Protein 4, Adipocyte) ELISA Kit
E-EL-P1062	Porcine PKCθ (Protein Kinase C Theta) ELISA Kit
E-EL-P1063	Porcine OSCAR (Osteoclast Associated Receptor) ELISA Kit
E-EL-P1064	Porcine CA4 (Carbonic Anhydrase IV) ELISA Kit
E-EL-P1065	Porcine TWEAK (Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit
E-EL-P1066	Porcine NPYR1 (Neuropeptide Y Receptor Y1) ELISA Kit
E-EL-P1067	Porcine GLTP (Glycolipid Transfer Protein) ELISA Kit
E-EL-P1068	Porcine CYP1B1 (Cytochrome P450 1B1) ELISA Kit
E-EL-P1070	Porcine ANXA5 (Annexin A5) ELISA Kit
E-EL-P1071	Porcine FST (Follistatin) ELISA Kit
E-EL-P1072	Porcine NRIP1 (Nuclear Receptor Interacting Protein 1) ELISA Kit

E-EL-P1073	Porcine c-myc (c-myc Oncogene Product) ELISA Kit
E-EL-P1074	Porcine A β 1-42 (Amyloid Beta Peptide 1-42) ELISA Kit
E-EL-P1075	Porcine CD109 (Cluster of Differentiation 109) ELISA Kit
E-EL-P1076	Porcine CDK2 (Cyclin Dependent Kinase 2) ELISA Kit
E-EL-P1077	Porcine MRP1 (Multidrug Resistance Associated Protein 1) ELISA Kit
E-EL-P1078	Porcine CPP (Copeptin) ELISA Kit
E-EL-P1079	Porcine AT (AntiThrombin) ELISA Kit
E-EL-P1080	Porcine α Manase (α Mannosidase) ELISA Kit
E-EL-P1081	Porcine DAG1 (Dystrophin Associated Glycoprotein 1) ELISA Kit
E-EL-P1083	Porcine SYNPO (Synaptopodin) ELISA Kit
E-EL-P1084	Porcine ALDOC (Aldolase C, Fructose Bisphosphate) ELISA Kit
E-EL-P1085	Porcine EREG (Epiregulin) ELISA Kit
E-EL-P1086	Porcine IGF2R (Insulin Like Growth Factor 2 Receptor) ELISA Kit
E-EL-P1088	Porcine NRP1 (Neuropilin 1) ELISA Kit
E-EL-P1089	Porcine Fc α R (Fc Fragment of IgA Receptor) ELISA Kit
E-EL-P1090	Porcine IGF2BP3 (Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit
E-EL-P1091	Porcine TRY (Trypsin) ELISA Kit
E-EL-P1093	Porcine ET (Epitestosterone) ELISA Kit
E-EL-P1094	Porcine GAL (Galanin) ELISA Kit
E-EL-P1095	Porcine EIPS2 (Estrogen Induced Protein pS2) ELISA Kit
E-EL-P1096	Porcine TACR1 (Tachykinin Receptor 1) ELISA Kit

E-EL-P1097	Porcine PRF1 (Perforin 1/Pore-Forming Protein) ELISA Kit
E-EL-P1098	Porcine NADPH (Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
E-EL-P1099	Porcine ACVA (Activin A) ELISA Kit
E-EL-P1100	Porcine CDC23 (Cell Division Cycle Protein 23) ELISA Kit
E-EL-P1101	Porcine ST2 (Syntenin 2) ELISA Kit
E-EL-P1102	Porcine TSP-4 (Thrombospondin-4) ELISA Kit
E-EL-P1103	Porcine α -Glu (Alpha Glucosidase) ELISA Kit
E-EL-P1104	Porcine GST θ 1/GSTt1 (Glutathione S Transferase Theta 1) ELISA Kit
E-EL-P1105	Porcine NFRkB (Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit
E-EL-P1106	Porcine VEGF (Vascular Endothelial Cell Growth Factor) ELISA Kit
E-EL-P1107	Porcine PGA (Pepsinogen A) ELISA Kit
E-EL-P1108	Porcine GM-CSF Ab (Anti-Granulocyte Macrophage Colony Stimulating Factor Antibody) ELISA Kit
E-EL-P1109	Porcine IL-23 (Interleukin 23) ELISA Kit
E-EL-P1110	Porcine NGF (Nerve Growth Factor) ELISA Kit
E-EL-P1111	Porcine CENPB (Centromere Protein B) ELISA Kit
E-EL-P1113	Porcine MUC20 (Mucin 20) ELISA Kit
E-EL-P1114	Porcine GTF2A1 (General Transcription Factor IIA, Polypeptide 1) ELISA Kit
E-EL-P1116	Porcine EL/LIPG (Endothelial Lipase) ELISA Kit
E-EL-P1118	Porcine sAC (Soluble Adenylate Cyclase) ELISA Kit
E-EL-P1119	Porcine IL-4 (Interleukin 4) ELISA Kit
E-EL-P1120	Porcine MANF (Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit

E-EL-P1121	Porcine VCAM-1/CD106 (soluble Vascular Cell Adhesion Molecule 1) ELISA Kit
E-EL-P1122	Porcine S100A12 (S100 Calcium Binding Protein A12) ELISA Kit
E-EL-P1123	Porcine TBG (Thyroxine Binding Globulin) ELISA Kit
E-EL-P1124	Porcine CDK5 (Cyclin Dependent Kinase 5) ELISA Kit
E-EL-P1125	Porcine VIP (Vasoactive Intestinal Peptide) ELISA Kit
E-EL-P1126	Porcine AIG1 (Androgen Induced Protein 1) ELISA Kit
E-EL-P1127	Porcine Rac-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit
E-EL-P1128	Porcine AGR2 (Anterior Gradient Protein 2) ELISA Kit
E-EL-P1129	Porcine EMP (Erythrocyte Membrane Protein) ELISA Kit
E-EL-P1130	Porcine sCD14 (Soluble Cluster of Differentiation 14) ELISA Kit
E-EL-P1132	Porcine NAIP (Neuronal Apoptosis Inhibitory Protein) ELISA Kit
E-EL-P1133	Porcine PP (Placental Protein) ELISA Kit
E-EL-P1134	Porcine MUC2 (Mucin 2) ELISA Kit
E-EL-P1136	Porcine AIRE (Autoimmune Regulator) ELISA Kit
E-EL-P1137	Porcine GRIK2 (Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit
E-EL-P1138	Porcine PC/CPG (Choline Phosphoglyceride) ELISA Kit
E-EL-P1139	Porcine TXLNb (Taxilin Beta) ELISA Kit
E-EL-P1140	Porcine AMY1 (Amylase Alpha 1, Salivary) ELISA Kit
E-EL-P1142	Porcine C4a (Complement Component 4a) ELISA Kit
E-EL-P1143	Porcine β 2-GP1 Ab IgG (β 2-Glycoprotein 1 Antibody IgG) ELISA Kit
E-EL-P1144	Porcine GREM1 (Gremlin 1) ELISA Kit

E-EL-P1145	Porcine ELMO2 (Engulfment And Cell Motility 2) ELISA Kit
E-EL-P1146	Porcine MAGE (Melanoma-associated Antigen) ELISA Kit
E-EL-P1147	Porcine MR (Mannose Receptor) ELISA Kit
E-EL-P1148	Porcine PRDX2 (Peroxisiredoxin 2) ELISA Kit
E-EL-P1149	Porcine STK39 (Serine/Threonine Kinase 39) ELISA Kit
E-EL-P1150	Porcine PCK1 (Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
E-EL-P1151	Porcine MRP (Multidrug Resistance-Associated Protein) ELISA Kit
E-EL-P1152	Porcine HSP-60/HSPD1 (Heat Shock Protein 60) ELISA Kit
E-EL-P1157	Porcine PDXK (Pyridoxal Kinase) ELISA Kit
E-EL-P1158	Porcine β -Glu (Beta Glucosidase) ELISA Kit
E-EL-P1159	Porcine P2 (Myelin Protein 2) ELISA Kit
E-EL-P1160	Porcine AMA IgA (Anti-Myelin antibody IgA) ELISA Kit
E-EL-P1162	Porcine Bcl-10 (B-cell Leukemia/Lymphoma 10) ELISA Kit
E-EL-P1163	Porcine NE/ELA2 (Elastase 2, Neutrophil) ELISA Kit
E-EL-P1164	Porcine PCDH β 2 (Protocadherin Beta 2) ELISA Kit
E-EL-P1165	Porcine CKLFSF8 (Chemokine Like Factor Superfamily 8) ELISA Kit
E-EL-P1166	Porcine MEPE (Matrix Extracellular Phosphoglycoprotein) ELISA Kit
E-EL-P1167	Porcine GEF (Guanylate Exchange Factor) ELISA Kit
E-EL-P1168	Porcine CREB (Cyclic AMP Response Element Binding Protein) ELISA Kit
E-EL-P1169	Porcine FZD8 (Frizzled Homolog 8) ELISA Kit
E-EL-P1170	Porcine SIGLEC3 (Sialic Acid Binding Ig Like Lectin 3) ELISA Kit

E-EL-P1171	Porcine GATA2 (GATA Binding Protein 2) ELISA Kit
E-EL-P1172	Porcine GMFβ (Glia Maturation Factor, Beta) ELISA Kit
E-EL-P1173	Porcine RUNX2(Runt-Related Transcription Factor 2)ELISA Kit
E-EL-P1174	Porcine γGT1 (Gamma Glutamyltransferase 1) ELISA Kit
E-EL-P1175	Porcine TGF-α (Transforming Growth Factor α) ELISA Kit
E-EL-P1176	Porcine EM2 (Endomorphin 2) ELISA Kit
E-EL-P1177	Porcine TL1A (Tumor Necrosis Factor Related Ligand 1A) ELISA Kit
E-EL-P1178	Porcine ETFβ (Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit
E-EL-P1179	Porcine NUP37 (Nucleoporin 37kDa) ELISA Kit
E-EL-P1182	Porcine β-CTx (Beta Crosslaps) ELISA Kit
E-EL-P1183	Porcine SEMA4C (Semaphorin 4C) ELISA Kit
E-EL-P1184	Porcine AKAP1 (A Kinase Anchor Protein 1) ELISA Kit
E-EL-P1185	Porcine ADCY4 (Adenylate Cyclase 4) ELISA Kit
E-EL-P1187	Porcine CMLC-1 (Cardiac Myosin Light Chains 1) ELISA Kit
E-EL-P1188	Porcine HB (Hemoglobin) ELISA Kit
E-EL-P1189	Porcine EGR1 (Early Growth Response Protein 1) ELISA Kit
E-EL-P1190	Porcine APA (Anti-Phospholipid Antibody) ELISA Kit
E-EL-P1191	Porcine PRKG1 (Protein Kinase, cGMP Dependent Type I) ELISA Kit
E-EL-P1192	Porcine HbA1c (Glycated Hemoglobin) ELISA Kit
E-EL-P1194	Porcine αZGP1 (Alpha-2-Glycoprotein 1, Zine Binding) ELISA Kit
E-EL-P1195	Porcine SPOP (Speckle Type POZ Protein) ELISA Kit

E-EL-P1196	Porcine CDNF (Cerebral Dopamine Neurotrophic Factor) ELISA Kit
E-EL-P1198	Porcine VIM (Vimentin) ELISA Kit
E-EL-P1199	Porcine EGF (Epidermal Growth Factor) ELISA Kit
E-EL-P1201	Porcine Pro-MMP-10 (Pro-Matrix Metalloproteinase 10) ELISA Kit
E-EL-P1202	Porcine NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit
E-EL-P1203	Porcine sICAM-1 (Soluble Intercellular Adhesion Molecule-1) ELISA Kit
E-EL-P1204	Porcine BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit
E-EL-P1205	Porcine Salusin α (Salusin Alpha) ELISA Kit
E-EL-P1206	Porcine DNM2 (Dynamin 2) ELISA Kit
E-EL-P1207	Porcine PLSCR2 (Phospholipid Scramblase 2) ELISA Kit
E-EL-P1208	Porcine APT (Abnormal Prothrombin) ELISA Kit
E-EL-P1209	Porcine GM (Gastric Mucin) ELISA Kit
E-EL-P1210	Porcine MT1 (Metallothionein 1) ELISA Kit
E-EL-P1211	Porcine FLNC (Filamin C, Gamma) ELISA Kit
E-EL-P1212	Porcine ADCY1 (Adenylate Cyclase 1, Brain) ELISA Kit
E-EL-P1214	Porcine LAP (Leucine Aminopeptidase) ELISA Kit
E-EL-P1215	Porcine PARC/CCL18 (Pulmonary Activation Regulated Chemokine) ELISA Kit
E-EL-P1216	Porcine Anti-CIIAb (Anti-Type II Collagen Antibody) ELISA Kit
E-EL-P1217	Porcine P-GP (Permeability Glycoprotein) ELISA Kit
E-EL-P1218	Porcine P-CK (Pan-Cytokeratin) ELISA Kit
E-EL-P1219	Porcine IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit

E-EL-P1220	Porcine MT2 (Metallothionein 2) ELISA Kit
E-EL-P1223	Porcine CCL4L1 (Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit
E-EL-P1224	Porcine CD163 (Cluster of Differentiation 163) ELISA Kit
E-EL-P1225	Porcine CDKN3 (Cyclin Dependent Kinase Inhibitor 3) ELISA Kit
E-EL-P1226	Porcine MSM β (Microseminoprotein Beta) ELISA Kit
E-EL-P1227	Porcine PDGFRL (Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
E-EL-P1228	Porcine CENPA (Centromere Protein A) ELISA Kit
E-EL-P1229	Porcine HELIX- II (Type II Collagen Helical Peptide) ELISA Kit
E-EL-P1230	Porcine PDCN (Podocin) ELISA Kit
E-EL-P1232	Porcine PRG4 (Proteoglycan 4) ELISA Kit
E-EL-P1233	Porcine TSP-2 (Thrombospondin-2) ELISA Kit
E-EL-P1234	Porcine CVNPF (Cobra Venom Neuronal Protective Factor) ELISA Kit
E-EL-P1235	Porcine DLD (Dihydrolipoyl Dehydrogenase) ELISA Kit
E-EL-P1236	Porcine VCL (Vinculin) ELISA Kit
E-EL-P1237	Porcine CKLF (Chemokine Like Factor) ELISA Kit
E-EL-P1238	Porcine CK-LMW (Cytokeratin Low Molecular Weight) ELISA Kit
E-EL-P1239	Porcine ANXA2 (Annexin A2) ELISA Kit
E-EL-P1240	Porcine ALP (Alkaline Phosphatase) ELISA Kit
E-EL-P1241	Porcine HNRPU (Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit
E-EL-P1242	Porcine CNR1 (Cannabinoid Receptor 1, Brain) ELISA Kit
E-EL-P1243	Porcine AMY2 (Amylase Alpha 2, Pancreatic) ELISA Kit

E-EL-P1244	Porcine AIRA (Anti-Insulin Receptor Antibody) ELISA Kit
E-EL-P1247	Porcine LTA4H (Leukotriene A4 Hydrolase) ELISA Kit
E-EL-P1248	Porcine GNβ1 (G Protein Beta 1) ELISA Kit
E-EL-P1249	Porcine PCDHβ16 (Protocadherin Beta 16) ELISA Kit
E-EL-P1250	Porcine CDSN (Corneodesmosin) ELISA Kit
E-EL-P1252	Porcine CCNB (Cyclin-B) ELISA Kit
E-EL-P1253	Porcine PKIβ (Protein Kinase Inhibitor Beta) ELISA Kit
E-EL-P1255	Porcine MFG-E8 (Milk Fat Globule EGF Factor 8) ELISA Kit
E-EL-P1256	Porcine ELA2B (Elastase 2B) ELISA Kit
E-EL-P1257	Porcine GAD-Ab (Anti-Glutamic Acid Decarboxylase antibody) ELISA Kit
E-EL-P1258	Porcine SF1 (Steroidogenic Factor 1) ELISA Kit
E-EL-P1259	Porcine CA242 (Pancreatic Carcinoma Markers-CA242) ELISA Kit
E-EL-P1261	Porcine HJV/HFE2 (Hemojuvelin/Hemochromatosis Type 2) ELISA Kit
E-EL-P1263	Porcine TRH (Thyrotropin Releasing Hormone) ELISA Kit
E-EL-P1264	Porcine CTGF (Connective Tissue Growth Factor) ELISA Kit
E-EL-P1265	Porcine SIGLEC12 (Sialic Acid Binding Ig Like Lectin 12) ELISA Kit
E-EL-P1266	Porcine AZU1 (Azurocidin 1) ELISA Kit
E-EL-P1267	Porcine M-AChR (Muscarinic Acetylcholine Receptor) ELISA Kit
E-EL-P1268	Porcine A-CoA (Acetyl Coenzyme A) ELISA Kit
E-EL-P1269	Porcine GAL12 (Galectin 12) ELISA Kit
E-EL-P1270	Porcine PRSS12 (Protease, Serine 12) ELISA Kit

E-EL-P1271	Porcine ADH1 (Alcohol Dehydrogenase 1) ELISA Kit
E-EL-P1272	Porcine TCR (Tcell Receptor) ELISA Kit
E-EL-P1273	Porcine IgA (Immunoglobulin A) ELISA Kit
E-EL-P1274	Porcine PRKAA1 (Protein Kinase, AMP Activated Alpha 1) ELISA Kit
E-EL-P1275	Porcine EP2 (Prostaglandin E Receptor 2) ELISA Kit
E-EL-P1277	Porcine RNASEN (Ribonuclease III, Nuclear) ELISA Kit
E-EL-P1278	Porcine CLEC4E (C-Type Lectin Domain Family 4 Member E) ELISA Kit
E-EL-P1279	Porcine NKA (Neurokinin A) ELISA Kit
E-EL-P1280	Porcine FABP6 (Fatty Acid Binding Protein 6, Ileal) ELISA Kit
E-EL-P1281	Porcine MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit
E-EL-P1282	Porcine CK-HMW (Cytokeratin High Molecular Weight) ELISA Kit
E-EL-P1283	Porcine HLTF (Helicase Like Transcription Factor) ELISA Kit
E-EL-P1286	Porcine fFDT1 (Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit
E-EL-P1288	Porcine Apo A1-Ab (Anti-Apolipoprotein A1 antibody) ELISA Kit
E-EL-P1289	Porcine COL1α1 (Collagen Type I Alpha 1) ELISA Kit
E-EL-P1290	Porcine TrxR1 (Thioredoxin Reductase 1) ELISA Kit
E-EL-P1291	Porcine GnRH (Gonadotropin Releasing Hormone) ELISA Kit
E-EL-P1292	Porcine CA1 (Carbonic Anhydrase I) ELISA Kit
E-EL-P1293	Porcine GROα/CXCL1/NAP3 (Growth Regulated Oncogene Alpha) ELISA Kit
E-EL-P1294	Porcine ELA4 (Elastase 4) ELISA Kit
E-EL-P1295	Porcine FLRT2 (Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit

E-EL-P1296	Porcine CDK1 (Cyclin Dependent Kinase 1) ELISA Kit
E-EL-P1297	Porcine ES (Endostatin) ELISA Kit
E-EL-P1298	Porcine CK-18/KRT18 (Cytokeratin 18) ELISA Kit
E-EL-P1300	Porcine MUSK (Muscle Skeletal Receptor Tyrosine Kinase) ELISA Kit
E-EL-P1302	Porcine EPHA10 (Ephrin Type A Receptor 10) ELISA Kit
E-EL-P1303	Porcine F8 (Coagulation Factor VIII) ELISA Kit
E-EL-P1304	Porcine MCP-3 (Monocyte Chemotactic Protein 3) ELISA Kit
E-EL-P1305	Porcine PKCε (Protein Kinase C Epsilon) ELISA Kit
E-EL-P1306	Porcine CD28 (Cluster of Differentiation 28) ELISA Kit
E-EL-P1307	Porcine SDH (Sorbitol Dehydrogenase) ELISA Kit
E-EL-P1308	Porcine PRCP (Prolylcarboxypeptidase) ELISA Kit
E-EL-P1309	Porcine C8b (Complement Component 8b) ELISA Kit
E-EL-P1310	Porcine TGM3 (Transglutaminase 3, Epidermal) ELISA Kit
E-EL-P1311	Porcine DAO (Diamine Oxidase) ELISA Kit
E-EL-P1312	Porcine M-AChR M3 (Muscarinic Acetylcholine Receptor M3) ELISA Kit
E-EL-P1314	Porcine 26S-PSM (26S Proteasome) ELISA Kit
E-EL-P1315	Porcine SOX9(Sex Determining Region Y Box Protein 9)ELISA Kit
E-EL-P1316	Porcine BF (Factor B) ELISA Kit
E-EL-P1317	Porcine PA2G4 (Proliferation Associated Protein 2G4) ELISA Kit
E-EL-P1318	Porcine TSTA (Tumor Specific Transplantation Antigen) ELISA Kit
E-EL-P1320	Porcine CSF3/GCSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit

E-EL-P1321	Porcine GDNF (Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
E-EL-P1322	Porcine EGFR (Epidermal Growth Factor Receptor) ELISA Kit
E-EL-P1323	Porcine BMP-3 (Bone Morphogenetic Protein 3) ELISA Kit
E-EL-P1324	Porcine MASP (Mannan Binding Lectin Associated Serine Protease) ELISA Kit
E-EL-P1327	Porcine MIF (Macrophage Migration Inhibitory Factor) ELISA Kit
E-EL-P1328	Porcine SDC1 (Syndecan 1) ELISA Kit
E-EL-P1329	Porcine RNASE (Ribonuclease) ELISA Kit
E-EL-P1330	Porcine ANG II R1/AGTR1 (Angiotensin II Receptor 1) ELISA Kit
E-EL-P1333	Porcine UBD (Ubiquitin D) ELISA Kit
E-EL-P1334	Porcine ApoD (Apolipoprotein D) ELISA Kit
E-EL-P1336	Porcine MXD1 (MAX Dimerization Protein 1) ELISA Kit
E-EL-P1337	Porcine MMP-8 (Matrix Metalloproteinase 8) ELISA Kit
E-EL-P1339	Porcine FABP7 (Fatty Acid Binding Protein 7, Brain) ELISA Kit
E-EL-P1340	Porcine PLA2G2A (Phospholipase A2, Group II A) ELISA Kit
E-EL-P1341	Porcine CPE (Carboxypeptidase E) ELISA Kit
E-EL-P1342	Porcine UBC (Urinary Bladder Cancer Antigen) ELISA Kit
E-EL-P1343	Porcine VEGFR3/Flt4 (Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit
E-EL-P1345	Porcine Ep-CAM/CD326 (Epithelial Cell Adhesion Molecule) ELISA Kit
E-EL-P1346	Porcine FAK (Focal Adhesion Kinase) ELISA Kit
E-EL-P1347	Porcine COL17 (Collagen Type XVII) ELISA Kit
E-EL-P1349	Porcine MMP-11 (Matrix Metalloproteinase 11) ELISA Kit

E-EL-P1350	Porcine IGFBP-5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit
E-EL-P1351	Porcine PTGIS (Prostaglandin I Synthase) ELISA Kit
E-EL-P1352	Porcine TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) ELISA Kit
E-EL-P1353	Porcine CX40 (Connexin 40) ELISA Kit
E-EL-P1354	Porcine COL10 (Collagen Type X) ELISA Kit
E-EL-P1355	Porcine FGF6 (Fibroblast Growth Factor 6) ELISA Kit
E-EL-P1356	Porcine C1qR1 (Complement Component 1q Receptor) ELISA Kit
E-EL-P1358	Porcine NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) ELISA Kit
E-EL-P1359	Porcine C5b (Complement Fragment 5b) ELISA Kit
E-EL-P1360	Porcine GSTs (Glutathione S-Transferases) ELISA Kit
E-EL-P1361	Porcine Ras-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit
E-EL-P1362	Porcine TNPO1 (Transportin 1) ELISA Kit
E-EL-P1363	Porcine UBAP2 (Ubiquitin Associated Protein 2) ELISA Kit
E-EL-P1364	Porcine IL-13 (Interleukin 13) ELISA Kit
E-EL-P1366	Porcine EJ/GlyRS (Anti-EJ-Antibody) ELISA Kit
E-EL-P1367	Porcine AATF (Apoptosis Antagonizing Transcription Factor) ELISA Kit
E-EL-P1368	Porcine cPLA2 (Cytosolic Phospholipase A2) ELISA Kit
E-EL-P1369	Porcine BECN1 (Beclin 1) ELISA Kit
E-EL-P1370	Porcine C4BP (C4 Binding Protein) ELISA Kit
E-EL-P1371	Porcine PLP (Proteolipid Protein) ELISA Kit
E-EL-P1372	Porcine SPTAN1 (Alpha-Fodrin) ELISA Kit

E-EL-P1373	Porcine GUCA2A (Guanylate Cyclase Activator 2A) ELISA Kit
E-EL-P1374	Porcine CK-13/KRT13 (Cytokeratin 13) ELISA Kit
E-EL-P1375	Porcine FOXO4 (Forkhead Box Protein O4) ELISA Kit
E-EL-P1376	Porcine CX43 (Connexin 43) ELISA Kit
E-EL-P1377	Porcine C3 (Complement Component 3) ELISA Kit
E-EL-P1379	Porcine EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit
E-EL-P1380	Porcine SAA2 (Serum Amyloid A2) ELISA Kit
E-EL-P1381	Porcine SIGLEC8 (Sialic Acid Binding Ig Like Lectin 8) ELISA Kit
E-EL-P1382	Porcine FGF21(Fibroblast Growth Factor 21) ELISA Kit
E-EL-P1383	Porcine PA-IgM (Platelet Associated Antibody) ELISA Kit
E-EL-P1384	Porcine TRF (Transferrin) ELISA Kit
E-EL-P1385	Porcine CX3CL1 (Chemokine C-X3-C-Motif Ligand) ELISA Kit
E-EL-P1386	Porcine CASP12 (Caspase 12) ELISA Kit
E-EL-P1387	Porcine TGF-β3 (Transforming Growth Factor β3) ELISA Kit
E-EL-P1389	Porcine NF-κB p65 (Nuclear Factor-κB p65) ELISA Kit
E-EL-P1390	Porcine OCLN (Occludin) ELISA Kit
E-EL-P1391	Porcine MlgA (Membrane IgA) ELISA Kit
E-EL-P1393	Porcine Cyt-C (Cytochrome C) ELISA Kit
E-EL-P1394	Porcine ApoE (Apolipoprotein E) ELISA Kit
E-EL-P1395	Porcine BMG/β2-MG (Beta-2-Microglobulin) ELISA Kit
E-EL-P1397	Porcine CD99 (Cluster of Differentiation 99) ELISA Kit

E-EL-P1398	Porcine MIP-1 β (Macrophage Inflammatory Protein 1 Beta)ELISA Kit
E-EL-P1399	Porcine ACTN3 (Actinin Alpha 3) ELISA Kit
E-EL-P1400	Porcine IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
E-EL-P1402	Porcine CNTN3 (Contactin 3) ELISA Kit
E-EL-P1403	Porcine CAPF (Cuticular Active Peptide Factor) ELISA Kit
E-EL-P1404	Porcine NRG-2 (Neuregulin 2) ELISA Kit
E-EL-P1405	Porcine FG β (Fibrinogen Beta) ELISA Kit
E-EL-P1406	Porcine NUP153 (Nucleoporin 153kDa) ELISA Kit
E-EL-P1407	Porcine SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-P1409	Porcine CYP26A1 (Cytochrome P450 26A1) ELISA Kit
E-EL-P1410	Porcine Grb2 (Growth Factor Receptor Bound Protein 2) ELISA Kit
E-EL-P1412	Porcine ADIPOR2 (Adiponectin Receptor 2) ELISA Kit
E-EL-P1413	Porcine HO1 (Heme Oxygenase 1) ELISA Kit
E-EL-P1414	Porcine PS (Phosphatidylserine) ELISA Kit
E-EL-P1415	Porcine KLK7 (Kallikrein 7) ELISA Kit
E-EL-P1416	Porcine SCG2 (Secretogranin II) ELISA Kit
E-EL-P1417	Porcine RAR α (Retinoic Acid Receptor Alpha) ELISA Kit
E-EL-P1419	Porcine ACA-IgM (Anti-Cardiolipin Antibody IgM) ELISA Kit
E-EL-P1420	Porcine γ -ECS (γ -Glutamyl Systeine Synthetase) ELISA Kit
E-EL-P1421	Porcine PDPN (Podoplanin) ELISA Kit
E-EL-P1422	Porcine MT1E (Metallothionein 1E) ELISA Kit

E-EL-P1423	Porcine GR β (Glucocorticoid Receptor Beta) ELISA Kit
E-EL-P1424	Porcine CLEC2C (C-Type Lectin Domain Family 2 Member C) ELISA Kit
E-EL-P1425	Porcine PDGF-AB (Platelet Derived Growth Factor AB) ELISA Kit
E-EL-P1426	Porcine Anti-GM1 (Anti-Ganglioside M1 Antibody) ELISA Kit
E-EL-P1427	Porcine PHC (Primary Hepatic Carcinoma) ELISA Kit
E-EL-P1428	Porcine TNC (Tenascin C) ELISA Kit
E-EL-P1429	Porcine GP4 (Platelet Membrane Glycoprotein IV) ELISA Kit
E-EL-P1430	Porcine DKK2 (Dickkopf Related Protein 2) ELISA Kit
E-EL-P1431	Porcine PLTP (Phospholipid Transfer Protein) ELISA Kit
E-EL-P1432	Porcine NUP85 (Nucleoporin 85kDa) ELISA Kit
E-EL-P1434	Porcine TLN (Talin) ELISA Kit
E-EL-P1435	Porcine LH (Luteinizing Hormone) ELISA Kit
E-EL-P1436	Porcine EFNA1 (Ephrin A1) ELISA Kit
E-EL-P1437	Porcine IL-1 α (Interleukin 1 Alpha)ELISA Kit
E-EL-P1438	Porcine CAP1 (Adenylate Cyclase Associated Protein 1) ELISA Kit
E-EL-P1439	Porcine CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit
E-EL-P1440	Porcine DSP (Desmoplakin) ELISA Kit
E-EL-P1441	Porcine HCCS (Holocytochrome C Synthase) ELISA Kit
E-EL-P1442	Porcine GNLY (Granulysin) ELISA Kit
E-EL-P1443	Porcine ApoA2 (Apolipoprotein A2) ELISA Kit
E-EL-P1444	Porcine TA (Anti-Teichoic Acid Antibody) ELISA Kit

E-EL-P1445	Porcine GasR (Gastrin Receptor) ELISA Kit
E-EL-P1446	Porcine PRSS1 (Protease, Serine, 1) ELISA Kit
E-EL-P1447	Porcine REG1α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit
E-EL-P1448	Porcine GPC4 (Glypican 4) ELISA Kit
E-EL-P1449	Porcine CRHBP (Corticotropin Releasing Hormone Binding Protein) ELISA Kit
E-EL-P1450	Porcine REG3γ (Regenerating Islet Derived Protein 3 Gamma) ELISA Kit
E-EL-P1451	Porcine TNX (Tenascin X) ELISA Kit
E-EL-P1452	Porcine PGR (Progesterone Receptor) ELISA Kit
E-EL-P1453	Porcine CCL16 (Chemokine C-C-Motif Ligand 16) ELISA Kit
E-EL-P1455	Porcine CYPD (Cyclophilin D) ELISA Kit
E-EL-P1456	Porcine DAP6 (Death Associated Protein 6) ELISA Kit
E-EL-P1457	Porcine GPI (Glycophosphatidylinositol) ELISA Kit
E-EL-P1458	Porcine Bcl-3 (B-cell Leukemia/Lymphoma 3) ELISA Kit
E-EL-P1459	Porcine TAP (Transporter Associated with Antigen Processing) ELISA Kit
E-EL-P1461	Porcine PKN1 (Protein Kinase N1) ELISA Kit
E-EL-P1462	Porcine RGS4 (Regulator of G Protein Signaling 4) ELISA Kit
E-EL-P1463	Porcine Pro-MMP-13 (Pro-Matrix Metalloproteinase 13) ELISA Kit
E-EL-P1464	Porcine CNTNAP1 (Contactin Associated Protein 1) ELISA Kit
E-EL-P1465	Porcine TSGF (Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
E-EL-P1466	Porcine SSTR1 (Somatostatin Receptor 1) ELISA Kit
E-EL-P1467	Porcine MYO (Myosin) ELISA Kit

E-EL-P1468	Porcine DR-70TM (Tumor Marker DR-70 For Lung Cancer) ELISA Kit
E-EL-P1469	Porcine VTN (Vitronectin) ELISA Kit
E-EL-P1470	Porcine MRC1 (Mannose Receptor C Type 1) ELISA Kit
E-EL-P1473	Porcine OLFM4 (Olfactomedin 4) ELISA Kit
E-EL-P1475	Porcine COL14 (Collagen Type XIV) ELISA Kit
E-EL-P1476	Porcine CRLF1 (Cytokine Receptor Like Factor 1) ELISA Kit
E-EL-P1477	Porcine GYPE (Glycophorin E) ELISA Kit
E-EL-P1478	Porcine CNTNAP4 (Contactin Associated Protein 4) ELISA Kit
E-EL-P1479	Porcine CENPE (Centromere Protein E) ELISA Kit
E-EL-P1480	Porcine AP (Allopregnanolone) ELISA Kit
E-EL-P1481	Porcine PDH α (Pyruvate Dehydrogenase Alpha) ELISA Kit
E-EL-P1482	Porcine FG α (Fibrinogen Alpha) ELISA Kit
E-EL-P1483	Porcine GDF1 (Growth Differentiation Factor 1) ELISA Kit
E-EL-P1486	Porcine CFB (Complement Factor B) ELISA Kit
E-EL-P1487	Porcine CYB561 (Cytochrome b-561) ELISA Kit
E-EL-P1488	Porcine MAOA (Type A Monoamine Oxidase) ELISA Kit
E-EL-P1489	Porcine EPOR (Erythropoietin Receptor) ELISA Kit
E-EL-P1490	Porcine PKR2 (Prokineticin Receptor 2) ELISA Kit
E-EL-P1491	Porcine TAA (Tumor Associated Antigen) ELISA Kit
E-EL-P1492	Porcine VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit
E-EL-P1493	Porcine PIK3C β (Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit

E-EL-P1494	Porcine PDAP1 (Pdgfa Associated Protein 1) ELISA Kit
E-EL-P1495	Porcine ODF (Osteoclast Differentiation Factor) ELISA Kit
E-EL-P1496	Porcine MMP-1 (Matrix Metalloproteinase 1) ELISA Kit
E-EL-P1500	Porcine SDC4 (Syndecan 4) ELISA Kit
E-EL-P1501	Porcine S100A13 (S100 Calcium Binding Protein A13) ELISA Kit
E-EL-P1503	Porcine MMP-3 (Matrix Metalloproteinase 3) ELISA Kit
E-EL-P1504	Porcine PYGL (Glycogen Phosphorylase, Liver) ELISA Kit
E-EL-P1505	Porcine IGFBP-4 (Insulin Like Growth Factor Binding Protein 4) ELISA Kit
E-EL-P1507	Porcine LHRH/GnRH (Luteinizing Hormone Releasing Hormone) ELISA Kit
E-EL-P1508	Porcine AIF (Apoptosis Inducing Factor) ELISA Kit
E-EL-P1509	Porcine EPO (Erythropoietin) ELISA Kit
E-EL-P1510	Porcine PHSA (Poly Serum Antigen) ELISA Kit
E-EL-P1511	Porcine FOXO2 (Forkhead Box Protein O2) ELISA Kit
E-EL-P1512	Porcine NES1 (Nesfatin 1) ELISA Kit
E-EL-P1513	Porcine PAPPA (Pregnancy Associated Plasma Protein A/ Pappalysin-1) ELISA Kit
E-EL-P1514	Porcine CALLA (Common Acute Lymphocytic Leukaemia Antigen) ELISA Kit
E-EL-P1515	Porcine GLβ (Galactosidase Beta) ELISA Kit
E-EL-P1516	Porcine CKLFSF4 (Chemokine Like Factor Superfamily 4) ELISA Kit
E-EL-P1517	Porcine PKC (Protein Kinase C) ELISA Kit
E-EL-P1519	Porcine NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
E-EL-P1520	Porcine EFNA4 (Ephrin A4) ELISA Kit

E-EL-P1521	Porcine FGF9 (Fibroblast Growth Factor 9) ELISA Kit
E-EL-P1522	Porcine MPIF-1/CCL23 (Myeloid Progenitor Inhibitory Factor 1) ELISA Kit
E-EL-P1523	Porcine LBP (Lipopolysaccharide Binding Protein) ELISA Kit
E-EL-P1524	Porcine IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) ELISA Kit
E-EL-P1525	Porcine NPHN (Nephrin) ELISA Kit
E-EL-P1526	Porcine FETUA (Fetuin A) ELISA Kit
E-EL-P1527	Porcine TARS (Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit
E-EL-P1528	Porcine sIgA (Secretory Immunoglobulin A) ELISA Kit
E-EL-P1529	Porcine CCSA-4 (Colon Cancer Specific Antigen-4) ELISA Kit
E-EL-P1530	Porcine C-P (C-Peptide) ELISA Kit
E-EL-P1531	Porcine REG4 (Regenerating Islet Derived Protein 4) ELISA Kit
E-EL-P1532	Porcine PS2 (Presenilin 2) ELISA Kit
E-EL-P1533	Porcine GC (Glycocalicin) ELISA Kit
E-EL-P1534	Porcine IL-22 (Interleukin 22) ELISA Kit
E-EL-P1535	Porcine β APP (β Amyloid Precursor Protein) ELISA Kit
E-EL-P1536	Porcine SOX18 (Sex Determining Region Y Box Protein 18) ELISA Kit
E-EL-P1537	Porcine Arg (Arginase) ELISA Kit
E-EL-P1538	Porcine BGN (Biglycan) ELISA Kit
E-EL-P1540	Porcine UBE1L (Ubiquitin Activating Enzyme E1 Like Protein) ELISA Kit
E-EL-P1541	Porcine PADI4 (Peptidyl Arginine Deiminase Type IV) ELISA Kit
E-EL-P1542	Porcine TUB β (β -Tubulin) ELISA Kit

E-EL-P1545	Porcine RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
E-EL-P1546	Porcine PP14 (Placental Protein14) ELISA Kit
E-EL-P1547	Porcine IRAK3 (Interleukin 1 Receptor Associated Kinase 3) ELISA Kit
E-EL-P1549	Porcine PAI2 (Plasminogen Activator Inhibitor 2) ELISA Kit
E-EL-P1550	Porcine PSPN (Persephin) ELISA Kit
E-EL-P1551	Porcine PGC (Pepsinogen C) ELISA Kit
E-EL-P1553	Porcine CCL14 (Chemokine C-C-Motif Ligand 14) ELISA Kit
E-EL-P1554	Porcine RARS (Arginyl tRNA Synthetase) ELISA Kit
E-EL-P1555	Porcine APAF-1 (Apoptosis Protease Activating Factor-1) ELISA Kit
E-EL-P1557	Porcine MPF (Maturation Promoting Factor) ELISA Kit
E-EL-P1559	Porcine ES-Ab (Anti-Endostatin Antibody) ELISA Kit
E-EL-P1560	Porcine AAA (Anti-Albumin Antibody) ELISA Kit
E-EL-P1561	Porcine FGL1 (Fibrinogen Like Protein 1) ELISA Kit
E-EL-P1562	Porcine DDC (Dopamine Decarboxylase) ELISA Kit
E-EL-P1563	Porcine MPA (Mycophenolic Acid) ELISA Kit
E-EL-P1564	Porcine PCI (Procollagen I) ELISA Kit
E-EL-P1565	Porcine CA3 (Carbonic Anhydrase III, Muscle Specific) ELISA Kit
E-EL-P1566	Porcine PUMA (P53 Upregulated Modulator of Apoptosis) ELISA Kit
E-EL-P1567	Porcine PDIA6 (Protein Disulfide Isomerase A6) ELISA Kit
E-EL-P1568	Porcine MTL (Motilin) ELISA Kit
E-EL-P1569	Porcine TFR2 (Transferrin Receptor 2) ELISA Kit

E-EL-P1570	Porcine SELL (L-Selectin) ELISA Kit
E-EL-P1571	Porcine CTXI (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
E-EL-P1572	Porcine τ PK1 (Tau-Protein Kinase) ELISA Kit
E-EL-P1574	Porcine PL (Pancreatic Lipase) ELISA Kit
E-EL-P1578	Porcine GMCSF/CSF2(Granulocyte Macrophage Colony Stimulating Factor 2) ELISA Kit
E-EL-P1579	Porcine VS (Versican/PG-M) ELISA Kit
E-EL-P1582	Porcine GAL2 (Galectin 2) ELISA Kit
E-EL-P1583	Porcine HMGCS1 (3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit
E-EL-P1584	Porcine NEB (Noradrenaline Bitartas) ELISA Kit
E-EL-P1585	Porcine PF4V1 (Platelet Factor 4 Variant 1) ELISA Kit
E-EL-P1586	Porcine FGy (Fibrinogen Gamma) ELISA Kit
E-EL-P1587	Porcine BMPR1A (Bone Morphogenetic Protein Receptor, type IA) ELISA Kit
E-EL-P1588	Porcine AKA (Anti-Keratin Antibody) ELISA Kit
E-EL-P1589	Porcine LRG1 (Leucine Rich Alpha-2-Glycoprotein 1) ELISA Kit
E-EL-P1590	Porcine IAPP (Islet Amyloid Polypeptide) ELISA Kit
E-EL-P1591	Porcine TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
E-EL-P1592	Porcine CTAGE (Cutaneous T-Cell Lymphocyte Associated Antigen) ELISA Kit
E-EL-P1593	Porcine PGM (Phosphoglucomutase) ELISA Kit
E-EL-P1594	Porcine ANG-R-Tie1 (Angiopoietin Receptor Tie1) ELISA Kit
E-EL-P1595	Porcine RETN (Resistin) ELISA Kit
E-EL-P1597	Porcine PKC δ 1 (Protein Kinase C Delta) ELISA Kit

E-EL-P1598	Porcine PSG1/SP1 (Pregnancy Specific Beta-1-Glycoprotein 1) ELISA Kit
E-EL-P1600	Porcine ApoA4 (Apolipoprotein A4) ELISA Kit
E-EL-P1602	Porcine MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
E-EL-P1603	Porcine KLK8 (Kallikrein 8) ELISA Kit
E-EL-P1604	Porcine C8g (Complement Component 8g) ELISA Kit
E-EL-P1605	Porcine FSTL1 (Follistatin Like Protein 1) ELISA Kit
E-EL-P1606	Porcine CNN1 (Calponin 1, Basic, Smooth Muscle) ELISA Kit
E-EL-P1607	Porcine FRA (Fibronectin Related Antigen) ELISA Kit
E-EL-P1608	Porcine IL-2 (Interleukin 2) ELISA Kit
E-EL-P1609	Porcine ASGR1 (Asialoglycoprotein Receptor 1) ELISA Kit
E-EL-P1611	Porcine CDK8 (Cyclin Dependent Kinase 8) ELISA Kit
E-EL-P1612	Porcine TP (Thymidine Phosphorylase) ELISA Kit
E-EL-P1615	Porcine PCT (Procalcitonin) ELISA Kit
E-EL-P1616	Porcine MS (Melatonin Sulfate) ELISA Kit
E-EL-P1617	Porcine CNTN1 (Contactin 1) ELISA Kit
E-EL-P1618	Porcine TNKS1 (Tankyrase 1) ELISA Kit
E-EL-P1619	Porcine PDGFB (Platelet Derived Growth Factor Subunit B) ELISA Kit
E-EL-P1620	Porcine HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit
E-EL-P1621	Porcine sTNF-1 (Soluble Tumor Necrosis Factor 1) ELISA Kit
E-EL-P1623	Porcine MUC7 (Mucin 7) ELISA Kit
E-EL-P1624	Porcine CCR1 (CC-Chemokine Receptor 1) ELISA Kit

E-EL-P1625	Porcine IGFBP-2 (Insulin Like Growth Factor Binding Protein 2) ELISA Kit
E-EL-P1626	Porcine SAP25 (Synaptosome Associated Protein 25) ELISA Kit
E-EL-P1627	Porcine ABCA1(ATP-binding cassette transporter A1) ELISA Kit
E-EL-P1629	Porcine Ntn4 (Netrin 4) ELISA Kit
E-EL-P1630	Porcine EMILIN 1 (Elastin Microfibril Interface Located Protein 1) ELISA Kit
E-EL-P1631	Porcine BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit
E-EL-P1632	Porcine C3bR (Complement Fragment 3b Receptor) ELISA Kit
E-EL-P1633	Porcine TR (Testosterone Receptor) ELISA Kit
E-EL-P1634	Porcine CD4 (Cluster of Differentiation 4) ELISA Kit
E-EL-P1635	Porcine BCL2L1/BCL-X (Bcl-2 Like Protein 1) ELISA Kit
E-EL-P1636	Porcine PL (Placental Lactogen) ELISA Kit
E-EL-P1637	Porcine HRG (Histidine Rich Glycoprotein) ELISA Kit
E-EL-P1638	Porcine ANGPTL3 (Angiopoietin Like Protein 3) ELISA Kit
E-EL-P1640	Porcine TALLA-1 (T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit
E-EL-P1642	Porcine PGAM2 (Phosphoglycerate Mutase 2) ELISA Kit
E-EL-P1643	Porcine SM/SPH (Sphingomyelin) ELISA Kit
E-EL-P1644	Porcine AR (Androgen Receptor) ELISA Kit
E-EL-P1645	Porcine SMM (Smooth Muscle Myosin) ELISA Kit
E-EL-P1647	Porcine ELA3B (Elastase 3B) ELISA Kit
E-EL-P1648	Porcine CA6 (Carbonic Anhydrase VI) ELISA Kit
E-EL-P1649	Porcine LSP (Liver Specific Lipoprotein) ELISA Kit

E-EL-P1650	Porcine CCKAR (Cholecystokinin A Receptor) ELISA Kit
E-EL-P1651	Porcine DSE (Dermatan Sulfate Epimerase) ELISA Kit
E-EL-P1652	Porcine CLDN1 (Claudin 1) ELISA Kit
E-EL-P1654	Porcine NUP214 (Nucleoporin 214kDa) ELISA Kit
E-EL-P1655	Porcine β -TG (β -Thromboglobulin) ELISA Kit
E-EL-P1656	Porcine UBR4 (Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit
E-EL-P1658	Porcine CTSL (Cathepsin L) ELISA Kit
E-EL-P1659	Porcine NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
E-EL-P1661	Porcine CASP7 (Caspase 7) ELISA Kit
E-EL-P1663	Porcine GT-34 (Gastrin 34) ELISA Kit
E-EL-P1665	Porcine Ub (Ubiquitin) ELISA Kit
E-EL-P1666	Porcine PRSS (Serine Protease) ELISA Kit
E-EL-P1667	Porcine EGR3 (Early Growth Response Protein 3) ELISA Kit
E-EL-P1668	Porcine PICK1 (Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
E-EL-P1669	Porcine CK-20/KRT20 (Cytokeratin 20) ELISA Kit
E-EL-P1670	Porcine F-TESTO (Free Testosterone) ELISA Kit
E-EL-P1671	Porcine LFA-1 (Lymphocyte Function Associated Antigen 1) ELISA Kit
E-EL-P1674	Porcine PC4 (Positive Cofactor 4) ELISA Kit
E-EL-P1675	Porcine UBE2A (Ubiquitin Conjugating Enzyme E2A) ELISA Kit
E-EL-P1676	Porcine PR3 (Proteinase 3) ELISA Kit
E-EL-P1677	Porcine FSP (Fibroblast Surface Protein) ELISA Kit

E-EL-P1678	Porcine ADM (Adrenomedullin) ELISA Kit
E-EL-P1679	Porcine CCSA-3 (Colon Cancer Specific Antigen-3) ELISA Kit
E-EL-P1680	Porcine DPYSL3 (Dihydropyrimidinase Like Protein 3) ELISA Kit
E-EL-P1681	Porcine PKC γ (Protein Kinase C Gamma) ELISA Kit
E-EL-P1682	Porcine IGFALS (Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) ELISA Kit
E-EL-P1683	Porcine MT1M (Metallothionein 1M) ELISA Kit
E-EL-P1684	Porcine BLVRB (Biliverdin Reductase B) ELISA Kit
E-EL-P1685	Porcine Anti-AChR (Anti-Acetylcholine Receptor Antibody) ELISA Kit
E-EL-P1686	Porcine GP5 (Glycoprotein V, Platelet) ELISA Kit
E-EL-P1688	Porcine ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit
E-EL-P1689	Porcine LGALS3BP (Lectin Galactoside Binding, Soluble 3 Binding Protein) ELISA Kit
E-EL-P1690	Porcine STMN1 (Stathmin 1) ELISA Kit
E-EL-P1691	Porcine CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit
E-EL-P1692	Porcine CYP24A1 (Cytochrome P450 24A1) ELISA Kit
E-EL-P1693	Porcine LDHA (Lactate Dehydrogenase A) ELISA Kit
E-EL-P1694	Porcine MSP (Macrophage Stimulating Protein) ELISA Kit
E-EL-P1695	Porcine COL4 (Collagen Type IV) ELISA Kit
E-EL-P1696	Porcine IL-20 (Interleukin 20) ELISA Kit
E-EL-P1697	Porcine PDP (Pyruvate Dehydrogenase Phosphatase) ELISA Kit
E-EL-P1698	Porcine F10 (Coagulation Factor X) ELISA Kit
E-EL-P1699	Porcine fT4 (Free Thyroxine) ELISA Kit

E-EL-P1701	Porcine LTD4 (Leukotriene D4) ELISA Kit
E-EL-P1702	Porcine COX-2R (Cyclooxygenase-2 Receptor) ELISA Kit
E-EL-P1703	Porcine α 1ACT (Alpha1 Antichymotrypsin) ELISA Kit
E-EL-P1704	Porcine CCND3 (Cyclin-D3) ELISA Kit
E-EL-P1705	Porcine GAD1 (Glutamate Decarboxylase 1) ELISA Kit
E-EL-P1706	Porcine MAPKAPK3 (Mitogen Activated Protein Kinase Activated Protein Kinase 3) ELISA Kit
E-EL-P1708	Porcine P-Cadherin (Cadherin, Placental) ELISA Kit
E-EL-P1709	Porcine BAFFR (B Cell Activating Factor Receptor) ELISA Kit
E-EL-P1711	Porcine LDLR (Low Density Lipoprotein Receptor) ELISA Kit
E-EL-P1712	Porcine SYP (Synaptophysin) ELISA Kit
E-EL-P1713	Porcine TCGF- III (Tcell Growth Factor III) ELISA Kit
E-EL-P1714	Porcine LT-E4 (Leukotriene E4) ELISA Kit
E-EL-P1715	Porcine CA9 (Carbonic Anhydrase IX) ELISA Kit
E-EL-P1717	Porcine BCDF (B-Cell Differentiation Factor) ELISA Kit
E-EL-P1718	Porcine VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
E-EL-P1719	Porcine HBF (Fetal Hemoglobin) ELISA Kit
E-EL-P1720	Porcine GARS (Glycyl tRNA Synthetase) ELISA Kit
E-EL-P1721	Porcine ADH (Antidiuretic Hormone) ELISA Kit
E-EL-P1722	Porcine SIGLEC2 (Sialic Acid Binding Ig Like Lectin 2) ELISA Kit
E-EL-P1724	Porcine CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
E-EL-P1725	Porcine TJP1 (Tight Junction Protein ZO-1) ELISA Kit

E-EL-P1726	Porcine TSLP (Thymic Stromal Lymphopoietin) ELISA Kit
E-EL-P1727	Porcine RPL26 (Ribosomal Protein L26) ELISA Kit
E-EL-P1728	Porcine FHA (Filamentous Hemagglutinin) ELISA Kit
E-EL-P1729	Porcine GAD2 (Glutamate Decarboxylase 2) ELISA Kit
E-EL-P1730	Porcine FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit
E-EL-P1731	Porcine REN (Renin) ELISA Kit
E-EL-P1732	Porcine AAA (Anti-Actin Antibody) ELISA Kit
E-EL-P1733	Porcine EGR2 (Early Growth Response Protein 2) ELISA Kit
E-EL-P1734	Porcine VEGF-B (Vascular Endothelial Cell Growth Factor B) ELISA Kit
E-EL-P1735	Porcine GSTθ2/GSTt2 (Glutathione S Transferase Theta 2) ELISA Kit
E-EL-P1736	Porcine CR1 (Complement Receptor type 1) ELISA Kit
E-EL-P1737	Porcine MYH10 (Myosin Heavy Chain 10, Non Muscle) ELISA Kit
E-EL-P1738	Porcine GSTκ1 (Glutathione S Transferase Kappa 1) ELISA Kit
E-EL-P1739	Porcine CYP1A1 (Cytochrome P450 1A1) ELISA Kit
E-EL-P1740	Porcine SNRPC (Small Nuclear Ribonucleoprotein Polypeptide C) ELISA Kit
E-EL-P1741	Porcine E-Cad (E-Cadherin) ELISA Kit
E-EL-P1742	Porcine DBP (DNA Binding Protein) ELISA Kit
E-EL-P1743	Porcine CDK9 (Cyclin Dependent Kinase 9) ELISA Kit
E-EL-P1744	Porcine STAT6 (Signal Transducer and Activator of Transcription 6) ELISA Kit
E-EL-P1745	Porcine AChRab (Acetylcholine Receptor Antibody) ELISA Kit
E-EL-P1747	Porcine MT3 (Metallothionein 3) ELISA Kit

E-EL-P1748	Porcine MBP/MBL (Mannose Binding Protein/Mannose Binding Lectin) ELISA Kit
E-EL-P1749	Porcine MYCBP (C-Myc Binding Protein) ELISA Kit
E-EL-P1750	Porcine HTRA1 (HtrA Serine, Peptidase 1) ELISA Kit
E-EL-P1751	Porcine AP 17 (Apelin 17) ELISA Kit
E-EL-P1752	Porcine BCAR1 (Breast Cancer Anti Estrogen Resistance 1) ELISA Kit
E-EL-P1753	Porcine Alk-Smase (Alkaline Sphingomyelinase) ELISA Kit
E-EL-P1754	Porcine COL5α2 (Collagen Type V Alpha 2) ELISA Kit
E-EL-P1755	Porcine IgG3 (Immunoglobulin G3) ELISA Kit
E-EL-P1756	Porcine C9 (Complement Component 9) ELISA Kit
E-EL-P1757	Porcine IκBα (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
E-EL-P1758	Porcine Collagenase II ELISA Kit
E-EL-P1759	Porcine TPO (Thrombopoietin) ELISA Kit
E-EL-P1760	Porcine Fbn (Fibrinolysin) ELISA Kit
E-EL-P1761	Porcine OJ/IleRS (Anti-OJ-Antibody) ELISA Kit
E-EL-P1762	Porcine FSA (Fetal Sulfoslycoprotein Antigen) ELISA Kit
E-EL-P1763	Porcine TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit
E-EL-P1764	Porcine βMSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-P1765	Porcine CA125 (Carbohydrate Antigen 125) ELISA Kit
E-EL-P1766	Porcine ANGPTL6 (Angiopoietin Like Protein 6) ELISA Kit
E-EL-P1768	Porcine F5 (Coagulation Factor V) ELISA Kit
E-EL-P1769	Porcine FBL (Fibrillarlin) ELISA Kit

E-EL-P1770	Porcine SYCN (Syncollin) ELISA Kit
E-EL-P1771	Porcine CHRM2 (Cholinergic Receptor, Muscarinic 2) ELISA Kit
E-EL-P1773	Porcine C1q-Ab (Anti-Complement 1q Antibody) ELISA Kit
E-EL-P1774	Porcine DKK1 (Dickkopf Related Protein 1) ELISA Kit
E-EL-P1775	Porcine S100A7 (S100 Calcium Binding Protein A7) ELISA Kit
E-EL-P1776	Porcine CAD (Caspase Activated Deoxyribonuclease) ELISA Kit
E-EL-P1777	Porcine SDF2L1 (Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
E-EL-P1778	Porcine CVL/EZR (Cytovillin/Ezrin) ELISA Kit
E-EL-P1779	Porcine LN (Laminin) ELISA Kit
E-EL-P1780	Porcine STAT5A (Signal Transducer and Activator of Transcription 5A) ELISA Kit
E-EL-P1781	Porcine GPC2 (Glypican 2) ELISA Kit
E-EL-P1782	Porcine CD30L (Cluster of Differentiation 30 Lligand) ELISA Kit
E-EL-P1783	Porcine LRP (Lung Resistance-Related Protein) ELISA Kit
E-EL-P1784	Porcine SAA (Serum Amyloid A) ELISA Kit
E-EL-P1785	Porcine NXF1 (Nuclear RNA Export Factor 1) ELISA Kit
E-EL-P1786	Porcine GLUT4 (Glucose Transporter 4) ELISA Kit
E-EL-P1787	Porcine CA12 (Carbonic Anhydrase XII) ELISA Kit
E-EL-P1788	Porcine PLA2G2D (Phospholipase A2, Group II D) ELISA Kit
E-EL-P1789	Porcine ALCAM (Activated Leukocyte Cell Adhesion Molecule) ELISA Kit
E-EL-P1790	Porcine GLP-2 (Glucagon Like Peptide 2) ELISA Kit
E-EL-P1791	Porcine RLR (Retinoic Acid Inducible Gene/RIG-Like Receptor) ELISA Kit

E-EL-P1792	Porcine BAD (Bcl2 Associated Death Promoter) ELISA Kit
E-EL-P1793	Porcine PKCζ (Protein Kinase C Zeta) ELISA Kit
E-EL-P1794	Porcine C5a (Complement Component 5a) ELISA Kit
E-EL-P1795	Porcine APLN (Apelin) ELISA Kit
E-EL-P1796	Porcine Cys-C (Cystatin C) ELISA Kit
E-EL-P1797	Porcine TP63 (Tumor Protein p63) ELISA Kit
E-EL-P1799	Porcine β-actin (Actin Beta) ELISA Kit
E-EL-P1800	Porcine cTnT/TNNT2 (Troponin T Type 2, Cardiac) ELISA Kit
E-EL-P1801	Porcine SOX2 (Sex Determining Region Y Box Protein 2) ELISA Kit
E-EL-P1802	Porcine CHAc (Choline acetylase) ELISA Kit
E-EL-P1803	Porcine MLZE (Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit
E-EL-P1805	Porcine ANKRD1 (Ankyrin Repeat Domain 1) ELISA Kit
E-EL-P1806	Porcine COL11α1 (Collagen Type XI Alpha 1) ELISA Kit
E-EL-P1807	Porcine PRB (Retinoblastoma Tumor Suppressor Protein) ELISA Kit
E-EL-P1808	Porcine TNF-β (Tumor Necrosis Factor Beta) ELISA Kit
E-EL-P1809	Porcine MAP2K1 (Mitogen Activated Protein Kinase Kinase 1) ELISA Kit
E-EL-P1810	Porcine FADS2 (Fatty Acid Desaturase 2) ELISA Kit
E-EL-P1811	Porcine Smad3/MADH3 (Mothers Against Decapentaplegic Homolog 3) ELISA Kit
E-EL-P1812	Porcine PDGFC (Platelet Derived Growth Factor C) ELISA Kit
E-EL-P1813	Porcine 12-HETE (12-Hydroxyeicosatetraenoic Acid) ELISA Kit
E-EL-P1814	Porcine DPD (Deoxypyridinoline) ELISA Kit

E-EL-P1815	Porcine Mly/CXCL9 (Monocyte Interferon Gamma Inducing Factor) ELISA Kit
E-EL-P1817	Porcine Pro-MMP-1 (pro-Matrix Metalloproteinase 1) ELISA Kit
E-EL-P1818	Porcine Tyk-2 (Tyrosine Kinase 2) ELISA Kit
E-EL-P1819	Porcine IFN- β (Interferon Beta) ELISA Kit
E-EL-P1820	Porcine eIF3a (Eukaryotic Translation Initiation Factor 3a) ELISA Kit
E-EL-P1821	Porcine RB1 (Retinoblastoma Protein 1) ELISA Kit
E-EL-P1822	Porcine COL6 α 1 (Collagen Type VI Alpha 1) ELISA Kit
E-EL-P1823	Porcine CDT (Carbohydrate Deficient Transferrin) ELISA Kit
E-EL-P1825	Porcine MMRN1 (Multimerin 1) ELISA Kit
E-EL-P1829	Porcine DPP4 (Dipeptidyl Peptidase IV) ELISA Kit
E-EL-P1830	Porcine BlySR (B Lymphocyte Stimulator Receptor) ELISA Kit
E-EL-P1832	Porcine FABP1 (Fatty Acid Binding Protein 1, Liver) ELISA Kit
E-EL-P1833	Porcine TPS (Tissue Polypeptide Specific Antigen) ELISA Kit
E-EL-P1834	Porcine INSL5 (Insulin Like Protein 5) ELISA Kit
E-EL-P1835	Porcine AIM2 (Absent In Melanoma 2) ELISA Kit
E-EL-P1836	Porcine PCI (Protein C Inhibitor) ELISA Kit
E-EL-P1837	Porcine CAM (Calmodulin) ELISA Kit
E-EL-P1838	Porcine BACE2 (β -site APP Cleaving Enzyme 2) ELISA Kit
E-EL-P1839	Porcine IAA (Insulin Autoantibody) ELISA Kit
E-EL-P1840	Porcine PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit
E-EL-P1841	Porcine PKC β 1 (Protein Kinase C Beta 1) ELISA Kit

E-EL-P1842	Porcine MAP τ (Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
E-EL-P1844	Porcine ACTN2 (Actinin Alpha 2) ELISA Kit
E-EL-P1845	Porcine IL-5 (Interleukin 5) ELISA Kit
E-EL-P1847	Porcine GST α 1 (Glutathione S Transferase Alpha 1) ELISA Kit
E-EL-P1848	Porcine LMP7/PSMB9 (Low Molecular-weight Protein/proteasome Beta-type Subunit) ELISA Kit
E-EL-P1849	Porcine tPA (Plasminogen Activator, Tissue) ELISA Kit
E-EL-P1850	Porcine TFR/CD71 (Transferrin Receptor) ELISA Kit
E-EL-P1851	Porcine NKB (Neurokinin B) ELISA Kit
E-EL-P1852	Porcine PPP1R1A (Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit
E-EL-P1853	Porcine NUP188 (Nucleoporin 188kDa) ELISA Kit
E-EL-P1854	Porcine GAL14 (Galectin 14) ELISA Kit
E-EL-P1855	Porcine ADAM8 (A Disintegrin And Metalloprotease 8) ELISA Kit
E-EL-P1856	Porcine VIPR1 (Vasoactive Intestinal Peptide Receptor) ELISA Kit
E-EL-P1857	Porcine TNNC2 (Troponin C Type 2) ELISA Kit
E-EL-P1858	Porcine PDIA5 (Protein Disulfide Isomerase A5) ELISA Kit
E-EL-P1859	Porcine MGP (Matrix Gla Protein) ELISA Kit
E-EL-P1860	Porcine IGFBP-6 (Insulin Like Growth Factor Binding Protein 6) ELISA Kit
E-EL-P1862	Porcine Pax (Paxillin) ELISA Kit
E-EL-P1865	Porcine CC16 (Clara Cell 16KD Protein) ELISA Kit
E-EL-P1867	Porcine GCLM (Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit
E-EL-P1868	Porcine HbC (Hemoglobin C) ELISA Kit

E-EL-P1869	Porcine GDI1 (GDP Dissociation Inhibitor 1) ELISA Kit
E-EL-P1870	Porcine AMBP (Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit
E-EL-P1871	Porcine CTSD (Cathepsin D) ELISA Kit
E-EL-P1872	Porcine GROβ (Growth Regulated Oncogene Beta) ELISA Kit
E-EL-P1874	Porcine THRSP (Thyroid Hormone Responsive) ELISA Kit
E-EL-P1875	Porcine EGR4 (Early Growth Response Protein 4) ELISA Kit
E-EL-P1876	Porcine Flt3L (FMS Like Tyrosine Kinase 3 Ligand) ELISA Kit
E-EL-P1877	Porcine GDH/GLDH (Glutamate dehydrogenase) ELISA Kit
E-EL-P1878	Porcine MAdCAM1 (Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit
E-EL-P1879	Porcine SEMA7A (Semaphorin 7A) ELISA Kit
E-EL-P1880	Porcine ETRB (Endothelin Receptor B) ELISA Kit
E-EL-P1881	Porcine IL-11 (Interleukin 11) ELISA Kit
E-EL-P1882	Porcine LAB7-2/CD86 (B-Lymphocyte Activation Antigen B7-2) ELISA Kit
E-EL-P1883	Porcine MYBPC3 (Myosin Binding Protein C, Cardiac) ELISA Kit
E-EL-P1884	Porcine Surv (Survivin) ELISA Kit
E-EL-P1885	Porcine UBE2C (Ubiquitin Conjugating Enzyme E2C) ELISA Kit
E-EL-P1887	Porcine MRC2 (Mannose Receptor C Type 2) ELISA Kit
E-EL-P1888	Porcine F9 (Coagulation Factor IX) ELISA Kit
E-EL-P1889	Porcine STAT4 (Signal Transducer and Activator of Transcription 4) ELISA Kit
E-EL-P1890	Porcine MAPK (Mitogen Activated Protein Kinase) ELISA Kit
E-EL-P1891	Porcine TEP1 (Telomerase Associated Protein 1) ELISA Kit

E-EL-P1892	Porcine OTF2B (Octamer Transcription Factor 2B) ELISA Kit
E-EL-P1893	Porcine ANGPTL1 (Angiopoietin Like Protein 1) ELISA Kit
E-EL-P1894	Porcine 16- α OHE-1 (16- α Hydroxyestrogen 1) ELISA Kit
E-EL-P1895	Porcine THAA (Thyroid Hormone Autoantibody) ELISA Kit
E-EL-P1896	Porcine ODC (Ornithine Decarboxylase) ELISA Kit
E-EL-P1897	Porcine Ntn1 (Netrin 1) ELISA Kit
E-EL-P1899	Porcine GR (Glutathione Reductase) ELISA Kit
E-EL-P1900	Porcine PLC γ 1 (Phospholipase C Gamma 1) ELISA Kit
E-EL-P1901	Porcine CYPB (Cyclophilin B) ELISA Kit
E-EL-P1902	Porcine ME20M (Melanoma Associated ME20) ELISA Kit
E-EL-P1903	Porcine ANG II R2/AGTR2 (Angiotensin II Receptor 2) ELISA Kit
E-EL-P1904	Porcine DT (Dopachrome Tautomerase) ELISA Kit
E-EL-P1905	Porcine AGC (Aggrecan) ELISA Kit
E-EL-P1906	Porcine CKLFSF6 (Chemokine Like Factor Superfamily 6) ELISA Kit
E-EL-P1907	Porcine FGF2BP2 (Fibroblast Growth Factor Binding Protein 2) ELISA Kit
E-EL-P1908	Porcine PCSK9 (Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit
E-EL-P1909	Porcine LIF (Leukemia Inhibitory Factor) ELISA Kit
E-EL-P1910	Porcine Anti-Surv (Anti-Survivin Antibody) ELISA Kit
E-EL-P1911	Porcine Rock-2 (Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit
E-EL-P1912	Porcine CPN10 (Chaperonin 10) ELISA Kit
E-EL-P1913	Porcine GTF2F1 (General Transcription Factor IIF, Polypeptide 1) ELISA Kit

E-EL-P1914	Porcine SHC/SLP-76 (SH2-Containing Protein) ELISA Kit
E-EL-P1916	Porcine TFF3 (Trefoil Factor 3, Intestinal) ELISA Kit
E-EL-P1917	Porcine CCSA-2 (Colon Cancer Specific Antigen-2) ELISA Kit
E-EL-P1919	Porcine IL-3 (Interleukin 3) ELISA Kit
E-EL-P1920	Porcine C4BP β (C4 Binding Protein Beta) ELISA Kit
E-EL-P1921	Porcine α 1-AT (Alpha 1-Antitrypsin) ELISA Kit
E-EL-P1922	Porcine GLUT1 (Glucose Transporter 1) ELISA Kit
E-EL-P1923	Porcine CRY1 (Cryptochrome 1) ELISA Kit
E-EL-P1924	Porcine ETRA (Endothelin Receptor A) ELISA Kit
E-EL-P1926	Porcine GDF10 (Growth Differentiation Factor 10) ELISA Kit
E-EL-P1928	Porcine D1R/DRD1 (Dopamine Receptor D1) ELISA Kit
E-EL-P1929	Porcine 5-LO (Arachidonate 5-Lipoxygenase) ELISA Kit
E-EL-P1930	Porcine ANNA1(Anti-Neuronal Nuclear Antibody 1) ELISA Kit
E-EL-P1932	Porcine CLU (Clusterin) ELISA Kit
E-EL-P1933	Porcine MX1 (Myxovirus Resistance 1) ELISA Kit
E-EL-P1934	Porcine TBP (TATA Binding Protein) ELISA Kit
E-EL-P1935	Porcine TPP (Thrombus Precursor Protein) ELISA Kit
E-EL-P1936	Porcine TSA (Thyroid Stimulating Antibody) ELISA Kit
E-EL-P1937	Porcine FETUB (Fetuin B) ELISA Kit
E-EL-P1938	Porcine RNASET2 (Ribonuclease T2) ELISA Kit
E-EL-P1939	Porcine DHEA (Dehydroepiandrosterone) ELISA Kit

E-EL-P1940	Porcine fT3 (Free Triiodothyronine) ELISA Kit
E-EL-P1942	Porcine BALP (Bone Alkaline Phosphatase) ELISA Kit
E-EL-P1943	Porcine HSP-90 (Heat Shock Protein 90) ELISA Kit
E-EL-P1944	Porcine ICAM-4 (Intercellular Adhesion Molecule-4) ELISA Kit
E-EL-P1945	Porcine ApoC3 (Apolipoprotein C3) ELISA Kit
E-EL-P1946	Porcine ANG II (Angiotensin II) ELISA Kit
E-EL-P1947	Porcine GTEAGE (Glomerular Tissue Extract Advanced Glycosylation End Products) ELISA Kit
E-EL-P1948	Porcine NAGase (N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
E-EL-P1949	Porcine ADAM10 (A Disintegrin And Metalloprotease 10) ELISA Kit
E-EL-P1951	Porcine DNase-I (Deoxyribonuclease I) ELISA Kit
E-EL-P1952	Porcine COL3 (Collagen Type III) ELISA Kit
E-EL-P1954	Porcine δ SIP- α (Delta Sleep Inducing Peptide α) ELISA Kit
E-EL-P1955	Porcine MASP1 (Mannan Associated Serine Protease 1) ELISA Kit
E-EL-P1956	Porcine C/EBP δ (CCAAT/Enhancer Binding Protein, Delta) ELISA Kit
E-EL-P1957	Porcine CNP (C-type Natriuretic Peptide) ELISA Kit
E-EL-P1958	Porcine PHLDA2 (Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit
E-EL-P1959	Porcine MUC5B (Mucin 5 Subtype B) ELISA Kit
E-EL-P1960	Porcine COL15 (Collagen Type XV) ELISA Kit
E-EL-P1962	Porcine LT β (Lymphotoxin Beta) ELISA Kit
E-EL-P1964	Porcine FPA (Fibrinopeptide A) ELISA Kit
E-EL-P1965	Porcine HLE (Leukocyte Exastase) ELISA Kit

E-EL-P1966	Porcine MBP (Myelin Basic Protein) ELISA Kit
E-EL-P1967	Porcine PIAS1 (Protein Inhibitor of Activated STAT 1) ELISA Kit
E-EL-P1968	Porcine ASGR2 (Asialoglycoprotein Receptor 2) ELISA Kit
E-EL-P1969	Porcine EAR12 (Eosinophil Associated Ribonuclease A family Member 12) ELISA Kit
E-EL-P1970	Porcine CEA (Carcinoembryonic Antigen) ELISA Kit
E-EL-P1971	Porcine NFκB (Nuclear Factor Kappa B) ELISA Kit
E-EL-P1972	Porcine Bcl2L11 (Bcl-2 Like Protein 11) ELISA Kit
E-EL-P1973	Porcine ADRA1A (Adrenergic Receptor alpha-1A) ELISA Kit
E-EL-P1974	Porcine Aβ1-40 (Amyloid Beta Peptide 1-40) ELISA Kit
E-EL-P1975	Porcine AMH (Anti-Mullerian Hormone) ELISA Kit
E-EL-P1976	Porcine CCP (Complement Control Protein) ELISA Kit
E-EL-P1977	Porcine CYTIP (Cytohesin Interacting Protein) ELISA Kit
E-EL-P1978	Porcine PCDH1 (Protocadherin 1) ELISA Kit
E-EL-P1979	Porcine PARP4 (Poly ADP Ribose Polymerase 4) ELISA Kit
E-EL-P1980	Porcine DEFα1 (Defensin Alpha 1, Neutrophil) ELISA Kit
E-EL-P1981	Porcine SAG (S Antigen) ELISA Kit
E-EL-P1983	Porcine ACO1 (Aconitase 1) ELISA Kit
E-EL-P1984	Porcine CENPH (Centromere Protein H) ELISA Kit
E-EL-P1985	Porcine BSP (Bone Sialoprotein) ELISA Kit
E-EL-P1987	Porcine ACAT2 (Acetyl Coenzyme A Acetyltransferase 2) ELISA Kit
E-EL-P1988	Porcine CMR (Costimulatory Molecules Receptor) ELISA Kit

E-EL-P1989	Porcine α 2-M (Alpha-2 Macroglobulin) ELISA Kit
E-EL-P1990	Porcine PLC γ 2 (Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit
E-EL-P1991	Porcine Fc γ R3B (Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit
E-EL-P1993	Porcine HB μ (Hemoglobin, Mu) ELISA Kit
E-EL-P1994	Porcine CBX3 (Chromobox Homolog 3) ELISA Kit
E-EL-P1995	Porcine GSP (Glycosylated Serum Protein) ELISA Kit
E-EL-P1996	Porcine OXA (Orexin A) ELISA Kit
E-EL-P1997	Porcine ANG I (Angiotensin I) ELISA Kit
E-EL-P1998	Porcine SPARC (Secreted Protein Acidic And Rich In Cysteine) ELISA Kit
E-EL-P1999	Porcine α 1BG (Alpha-1-B-Glycoprotein) ELISA Kit
E-EL-P2000	Porcine Anti-MPO (Anti-Myeloperoxidase Antibody) ELISA Kit
E-EL-P2001	Porcine OxLDL (Oxidized Low Density Lipoprotein) ELISA Kit
E-EL-P2003	Porcine KIM-1(Kidney Injury Molecule 1)ELISA Kit
E-EL-P2004	Porcine P27 (P27 Protein) ELISA Kit
E-EL-P2005	Porcine PF4 (Platelet Factor 4) ELISA Kit
E-EL-P2007	Porcine CS (Chondroitin Sulfate) ELISA Kit
E-EL-P2008	Porcine PC (Plasma Anticoagulant Protein C) ELISA Kit
E-EL-P2009	Porcine α 1-AGP (Alpha-1-Acid Glycoprotein) ELISA Kit
E-EL-P2010	Porcine PLC (Phospholipase C) ELISA Kit
E-EL-P2011	Porcine TYR (Tyrosinase) ELISA Kit
E-EL-P2013	Porcine IL-9 (Interleukin 9) ELISA Kit

E-EL-P2014	Porcine ARR β 2 (Arrestin Beta 2) ELISA Kit
E-EL-P2015	Porcine ARHGD1 α (Rho GDP Dissociation Inhibitor Alpha) ELISA Kit
E-EL-P2016	Porcine ANG II R2-Ab (Angiotensin II Receptor 2 Antibody) ELISA Kit
E-EL-P2017	Porcine NLR (NOD-Like Receptor9) ELISA Kit
E-EL-P2018	Porcine M-AChR M1 (Muscarinic Acetylcholine Receptor M1) ELISA Kit
E-EL-P2019	Porcine FGFR3 (Fibroblast Growth Factor Receptor 3) ELISA Kit
E-EL-P2020	Porcine NEFL (Neurofilament, Light Polypeptide) ELISA Kit
E-EL-P2021	Porcine TCF20 (Transcription Factor 20) ELISA Kit
E-EL-P2022	Porcine GzmB (Granzyme B) ELISA Kit
E-EL-P2023	Porcine CG (Cryoglobulin) ELISA Kit
E-EL-P2024	Porcine LCA/CD45(Leukocyte Common Antigen) ELISA Kit
E-EL-P2025	Porcine Sp100 (Sp100 Nuclear Antigen) ELISA Kit
E-EL-P2026	Porcine FSH (Follicle-Stimulating Hormone) ELISA Kit
E-EL-P2027	Porcine TDP43 (Tar DNA Binding Protein 43KDa) ELISA Kit
E-EL-P2028	Porcine CK-19/KRT19 (Cytokeratin 19) ELISA Kit
E-EL-P2029	Porcine AGER (Advanced Glycosylation End Product Specific Receptor) ELISA Kit
E-EL-P2030	Porcine LGALS3 (Lectin-Galactoside Binding Soluble 3) ELISA Kit
E-EL-P2031	Porcine BRAK/CXCL14 (Breast and Kidney Expressed Chemokine) ELISA Kit
E-EL-P2032	Porcine OTR (Oxytocin Receptor) ELISA Kit
E-EL-P2033	Porcine COL4 α 1 (Collagen Type IV Alpha 1) ELISA Kit
E-EL-P2034	Porcine PKR1 (Prokineticin Receptor 1) ELISA Kit

E-EL-P2035	Porcine OBTF2 (Octamer Binding Transcription Factor 2) ELISA Kit
E-EL-P2036	Porcine BMP-10 (Bone Morphogenetic Protein 10) ELISA Kit
E-EL-P2039	Porcine TGFβ1 (Transforming Growth Factor Beta Induced Protein) ELISA Kit
E-EL-P2040	Porcine PTPRK (Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit
E-EL-P2041	Porcine HbH (Hemoglobin H) ELISA Kit
E-EL-P2042	Porcine BCR (B-Cell Receptor) ELISA Kit
E-EL-P2044	Porcine Cav-1 (Caveolin-1) ELISA Kit
E-EL-P2046	Porcine HSF1 (Heat Shock Transcription Factor 1) ELISA Kit
E-EL-P2047	Porcine CSTA (Cystatin A) ELISA Kit
E-EL-P2048	Porcine GHR (Growth Hormone Receptor) ELISA Kit
E-EL-P2049	Porcine S100A11 (S100 Calcium Binding Protein A11) ELISA Kit
E-EL-P2050	Porcine C7 (Complement Component 7) ELISA Kit
E-EL-P2051	Porcine FABP8 (Fatty Acid Binding Protein 8, Myelin) ELISA Kit
E-EL-P2052	Porcine MlgM (Membrane IgM) ELISA Kit
E-EL-P2054	Porcine ConA (Concanavalin A) ELISA Kit
E-EL-P2055	Porcine SAA3 (Serum Amyloid A3) ELISA Kit
E-EL-P2056	Porcine SDC3 (Syndecan 3) ELISA Kit
E-EL-P2057	Porcine PLP1 (Proteolipid Protein 1, Myelin) ELISA Kit
E-EL-P2058	Porcine SREBP (Sterol Regulatory Element Binding Proteins) ELISA Kit
E-EL-P2059	Porcine SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit
E-EL-P2060	Porcine GPER1 (G Protein Coupled Estrogen Receptor 1) ELISA Kit

E-EL-P2061	Porcine GALR1 (Galanin Receptor 1) ELISA Kit
E-EL-P2063	Porcine MFAP4 (Microfibrillar Associated Protein 4) ELISA Kit
E-EL-P2064	Porcine PDIA3 (Protein Disulfide Isomerase A3) ELISA Kit
E-EL-P2065	Porcine PRDX1 (Peroxisome oxidoreductin 1) ELISA Kit
E-EL-P2066	Porcine PLAP/ALPP (Placental Alkaline Phosphatase) ELISA Kit
E-EL-P2067	Porcine CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit
E-EL-P2068	Porcine NSE (Neuron-Specific Enolase) ELISA Kit
E-EL-P2069	Porcine GSK3 β (Glycogen Synthase Kinase 3 Beta) ELISA Kit
E-EL-P2070	Porcine PLA2G12B (Phospholipase A2, Group XII B) ELISA Kit
E-EL-P2072	Porcine TREM-1 (Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
E-EL-P2073	Porcine MFI2 (Melanotransferrin) ELISA Kit
E-EL-P2074	Porcine ACP5 (Tartrate Resistant Acid Phosphatase 5) ELISA Kit
E-EL-P2075	Porcine Thyronine ELISA Kit
E-EL-P2076	Porcine ADP/Acrp30(Adiponectin)ELISA Kit
E-EL-P2077	Porcine AGRN (Agrin) ELISA Kit
E-EL-P2078	Porcine AP36 (Apelin 36) ELISA Kit
E-EL-P2079	Porcine IL-1ra (Interleukin 1 Receptor Antagonist) ELISA Kit
E-EL-P2080	Porcine GBP4 (Guanylate Binding Protein 4) ELISA Kit
E-EL-P2082	Porcine MMP-10 (Matrix Metalloproteinase 10) ELISA Kit
E-EL-P2083	Porcine CTSB (Cathepsin B) ELISA Kit
E-EL-P2084	Porcine NTRK3 (Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit

E-EL-P2085	Porcine GUCA2B (Guanylate Cyclase Activator 2B) ELISA Kit
E-EL-P2086	Porcine GlyCAM1 (Glycosylation Dependent Cell Adhesion Molecule 1) ELISA Kit
E-EL-P2087	Porcine SLC/CCL21 (Secondary Lymphoid Tissue Chemokine) ELISA Kit
E-EL-P2088	Porcine PRSS33 (Protease, Serine 33) ELISA Kit
E-EL-P2089	Porcine PKI α (Protein Kinase Inhibitor Alpha) ELISA Kit
E-EL-P2090	Porcine PAQR7 (Progesterin and AdipoQ Receptor Family Member VII) ELISA Kit
E-EL-P2091	Porcine GST α 5 (Glutathione S Transferase Alpha 5) ELISA Kit
E-EL-P2093	Porcine AAP (Alanine Aminopeptidase) ELISA Kit
E-EL-P2094	Porcine COL3 α 1 (Collagen Type III Alpha 1) ELISA Kit
E-EL-P2095	Porcine ATXN2 (Ataxin 2) ELISA Kit
E-EL-P2096	Porcine ASM (Acid Sphingomyelinase) ELISA Kit
E-EL-P2097	Porcine ACT α 2 (Actin Alpha 2, Smooth Muscle) ELISA Kit
E-EL-P2098	Porcine F13A1 (Coagulation Factor XIII A1 Polypeptide) ELISA Kit
E-EL-P2099	Porcine MSTN (Myostatin) ELISA Kit
E-EL-P2100	Porcine SST (Somatostatin) ELISA Kit
E-EL-P2101	Porcine GDF2 (Growth Differentiation Factor 2) ELISA Kit
E-EL-P2102	Porcine CKMB (Creatine Kinase MB Isoenzyme) ELISA Kit
E-EL-P2103	Porcine ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit
E-EL-P2104	Porcine PKC η (Protein Kinase C Eta) ELISA Kit
E-EL-P2105	Porcine PGF (Prostaglandin F) ELISA Kit
E-EL-P2106	Porcine C3a (Complement Component 3a) ELISA Kit

E-EL-P2108	Porcine SOD3 (Superoxide Dismutase 3, Extracellular) ELISA Kit
E-EL-P2109	Porcine CYPA (Cyclophilin A) ELISA Kit
E-EL-P2110	Porcine RPA-70 (replication protein A) ELISA Kit
E-EL-P2112	Porcine FGFR1 (Fibroblast Growth Factor Receptor 1) ELISA Kit
E-EL-P2113	Porcine LIPF (Lipase, Gastric) ELISA Kit
E-EL-P2114	Porcine MUC6 (Mucin 6) ELISA Kit
E-EL-P2115	Porcine CX31 (Connexin 31) ELISA Kit
E-EL-P2117	Porcine DDT (D-Dopachrome Tautomerase) ELISA Kit
E-EL-P2118	Porcine ProGRP (Pro-Gastrin Releasing Peptide) ELISA Kit
E-EL-P2119	Porcine PR3Ab (Proteinase 3 Antibody) ELISA Kit
E-EL-P2120	Porcine Lptn/XCL1 (Lymphotactin) ELISA Kit
E-EL-P2121	Porcine FUCA (Alpha-L-Fucosidase, Tissue) ELISA Kit
E-EL-P2122	Porcine α FP (Alpha-Fetoprotein) ELISA Kit
E-EL-P2123	Porcine CFD (Complement Factor D) ELISA Kit
E-EL-P2124	Porcine NRG-4 (Neuregulin 4) ELISA Kit
E-EL-P2125	Porcine ADCY2 (Adenylate Cyclase 2) ELISA Kit
E-EL-P2126	Porcine PGD2-MOX (Prostaglandin D2-Methoxime) ELISA Kit
E-EL-P2127	Porcine OFQ/N (Orphanin FQ/Nociceptin) ELISA Kit
E-EL-P2128	Porcine RCOR3 (Rest Corepressor 3) ELISA Kit
E-EL-P2129	Porcine CD30 (Cluster of Differentiation 30) ELISA Kit
E-EL-P2131	Porcine PTH (Parathyroid Hormone) ELISA Kit

E-EL-P2132	Porcine β -EPR (Beta-Endorphin Receptor) ELISA Kit
E-EL-P2133	Porcine HHIP (Hedgehog Interacting Protein) ELISA Kit
E-EL-P2134	Porcine α -La (Alpha-Lactalbumin) ELISA Kit
E-EL-P2136	Porcine MF/MPF (Maturation Promoting Factor) ELISA Kit
E-EL-P2137	Porcine FOXO3 (Forkhead Box Protein O3) ELISA Kit
E-EL-P2138	Porcine PLVAP (Plasmalemma Vesicle Associated Protein) ELISA Kit
E-EL-P2140	Porcine Gal (α -Galactoyl) ELISA Kit
E-EL-P2141	Porcine PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit
E-EL-P2142	Porcine Tn-T (Troponin T) ELISA Kit
E-EL-P2143	Porcine HSP-40 (Heat Shock Protein 40) ELISA Kit
E-EL-P2144	Porcine PRB1 (Basic Salivary Proline Rich Protein 1) ELISA Kit
E-EL-P2145	Porcine AFP-L3 (Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit
E-EL-P2146	Porcine PGDS (Prostaglandin D Synthase) ELISA Kit
E-EL-P2147	Porcine PAR2 (Protease Activated Receptor 2) ELISA Kit
E-EL-P2148	Porcine EP3 (Prostaglandin E Receptor 3) ELISA Kit
E-EL-P2149	Porcine SBEM (Small Breast Epithelial Mucin) ELISA Kit
E-EL-P2150	Porcine BMPR2 (Bone Morphogenetic Protein Receptor II) ELISA Kit
E-EL-P2151	Porcine NUP62 (Nucleoporin 62kDa) ELISA Kit
E-EL-P2152	Porcine ECE1 (Endothelin Converting Enzyme 1) ELISA Kit
E-EL-P2153	Porcine AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) ELISA Kit
E-EL-P2154	Porcine CD42 (Cluster of Differentiation 42) ELISA Kit

E-EL-P2155	Porcine MUC4 (Mucin 4) ELISA Kit
E-EL-P2156	Porcine NUMA1 (Nuclear Mitotic Apparatus Protein 1) ELISA Kit
E-EL-P2157	Porcine SOX1 (Sex Determining Region Y Box Protein 1) ELISA Kit
E-EL-P2158	Porcine MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
E-EL-P2159	Porcine ApoA5 (Apolipoprotein A5) ELISA Kit
E-EL-P2160	Porcine 5-HIAA (5-Hydroxyindolacetic Acid) ELISA Kit
E-EL-P2161	Porcine RELB (V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
E-EL-P2162	Porcine VacA (Vacuolating Cytotoxin A) ELISA Kit
E-EL-P2163	Porcine F1+2 (Prothrombin Fragment 1+2) ELISA Kit
E-EL-P2164	Porcine IDU α (Alpha-L-Iduronidase) ELISA Kit
E-EL-P2165	Porcine GDF7 (Growth Differentiation Factor 7) ELISA Kit
E-EL-P2166	Porcine Cyr61 (Cysteine Rich Protein, Angiogenic Inducer 61) ELISA Kit
E-EL-P2167	Porcine GA (Glycated Albumin) ELISA Kit
E-EL-P2168	Porcine TAT (Thrombin/Antithrombin Complex) ELISA Kit
E-EL-P2169	Porcine ENG/CD105 (Endoglin) ELISA Kit
E-EL-P2170	Porcine CP (Ceruloplasmin) ELISA Kit
E-EL-P2172	Porcine AIF1 (Allograft Inflammatory Factor 1) ELISA Kit
E-EL-P2173	Porcine L1CAM/CD171 (L1-Cell Adhesion Molecule) ELISA Kit
E-EL-P2174	Porcine CTSK (Cathepsin K) ELISA Kit
E-EL-P2175	Porcine GADL1 (Glutamate Decarboxylase Like Protein 1) ELISA Kit
E-EL-P2176	Porcine GAP43 (Growth Associated Protein 43) ELISA Kit

E-EL-P2178	Porcine OC/BGP (Osteocalcin/Bone Gla Protein) ELISA Kit
E-EL-P2179	Porcine DβH (Dopamine-β-Hydroxylase) ELISA Kit
E-EL-P2180	Porcine MDC/CCL22 (Macrophage Derived Chemokine) ELISA Kit
E-EL-P2182	Porcine Preptin ELISA Kit
E-EL-P2184	Porcine TP5 (Thymopentin) ELISA Kit
E-EL-P2185	Porcine SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit
E-EL-P2186	Porcine TG (Thyroglobulin) ELISA Kit
E-EL-P2187	Porcine PT (Prothrombin) ELISA Kit
E-EL-P2188	Porcine SNCα (Synuclein Alpha) ELISA Kit
E-EL-P2190	Porcine GSTα4 (Glutathione S Transferase Alpha 4) ELISA Kit
E-EL-P2191	Porcine S100A6 (S100 Calcium Binding Protein A6) ELISA Kit
E-EL-P2192	Porcine CDK6 (Cyclin Dependent Kinase 6) ELISA Kit
E-EL-P2193	Porcine PP(Pepsin)ELISA Kit
E-EL-P2194	Porcine SUMO2 (Small Ubiquitin Related Modifier 2) ELISA Kit
E-EL-P2195	Porcine Smad4/MADH4 (Mothers Against Decapentaplegic Homolog 4) ELISA Kit
E-EL-P2196	Porcine PRSS2 (Protease, Serine, 2) ELISA Kit
E-EL-P2197	Porcine TNNT1 (Troponin T Type 1, Slow Skeletal) ELISA Kit
E-EL-P2198	Porcine OPG (Osteoprotegerin) ELISA Kit
E-EL-P2199	Porcine PKB (Protein Kinase B) ELISA Kit
E-EL-P2200	Porcine ERK2 (Extracellular Signal Regulated Kinase 2) ELISA Kit
E-EL-P2202	Porcine Sp17 (Sperm Protein 17) ELISA Kit

E-EL-P2203	Porcine PLA2G10 (Phospholipase A2, Group X) ELISA Kit
E-EL-P2205	Porcine TLR9 (Toll-Like Receptor 9) ELISA Kit
E-EL-P2206	Porcine CLDN9 (Claudin 9) ELISA Kit
E-EL-P2207	Porcine PHB (Prohibitin) ELISA Kit
E-EL-P2208	Porcine β 2-GP1 Ab IgM (β 2-Glycoprotein 1 Antibody IgM) ELISA Kit
E-EL-P2209	Porcine F13B (Coagulation Factor XIII B Polypeptide) ELISA Kit
E-EL-P2210	Porcine ATF-4 (Activating Transcription Factor 4) ELISA Kit
E-EL-P2211	Porcine CgA (Chromogranin A/Pancreastatin) ELISA Kit
E-EL-P2212	Porcine ALDH (Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit
E-EL-P2214	Porcine PPAR- γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
E-EL-P2215	Porcine VEGFR1/Flt1 (Vasculoar Endothelial Growth Factor Receptor 1) ELISA Kit
E-EL-P2216	Porcine PLA2R1 (Phospholipase A2 Receptor 1) ELISA Kit
E-EL-P2217	Porcine Smad1/MADH1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit
E-EL-P2218	Porcine TAF (TATA Box Binding Protein/TBP-Associated Factors) ELISA Kit
E-EL-P2219	Porcine OIT3 (Oncoprotein Induced Transcript 3) ELISA Kit
E-EL-P2222	Porcine TSI (Thyroid Stimulating Immunoglobulin) ELISA Kit
E-EL-P2223	Porcine VLA (Very Late Appearing Antigen) ELISA Kit
E-EL-P2224	Porcine sTfR1 (Soluble Transferrin Receptor1) ELISA Kit
E-EL-P2225	Porcine FP (Prostaglandin F Receptor) ELISA Kit
E-EL-P2226	Porcine SHH (Hedgehog Homolog, Sonic) ELISA Kit
E-EL-P2227	Porcine GP6 (Glycoprotein VI, Platelet) ELISA Kit

E-EL-P2229	Porcine TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit
E-EL-P2230	Porcine GAL3 (Galectin 3) ELISA Kit
E-EL-P2232	Porcine HGF (Hepatocyte Growth Factor) ELISA Kit
E-EL-P2234	Porcine PIAS3 (Protein Inhibitor of Activated STAT3) ELISA Kit
E-EL-P2236	Porcine MCSFR (Macrophage Colony Stimulating Factor Receptor) ELISA Kit
E-EL-P2237	Porcine D2D (D-Dimer) ELISA Kit
E-EL-P2238	Porcine β 2-GP1 Ab IgA (β 2-Glycoprotein 1 antibody IgA) ELISA Kit
E-EL-P2239	Porcine CX37 (Connexin 37) ELISA Kit
E-EL-P2240	Porcine VLDLR (Very Low Density Lipoprotein Receptor) ELISA Kit
E-EL-P2241	Porcine bTC (Betacellulin) ELISA Kit
E-EL-P2242	Porcine PLAA (Phospholipase A2 Activating Protein) ELISA Kit
E-EL-P2243	Porcine PDH E1 (Pyruvate Dehydrogenase E1) ELISA Kit
E-EL-P2244	Porcine AQP-4 (Aquaporin 4) ELISA Kit
E-EL-P2245	Porcine SFRP2 (Secreted Frizzled Related Protein 2) ELISA Kit
E-EL-P2246	Porcine α -SCA (α -sarcomeric Actin) ELISA Kit
E-EL-P2247	Porcine Hepc (Hepcidin) ELISA Kit
E-EL-P2248	Porcine PAR4 (Protease Activated Receptor 4) ELISA Kit
E-EL-P2249	Porcine CD40L/TNFSF5 (Cluster of Differentiation 40 Ligand) ELISA Kit
E-EL-P2251	Porcine SPA (Pulmonary Surfactant Associated Protein A) ELISA Kit
E-EL-P2252	Porcine β -Lg (β -Lactoglobulin) ELISA Kit
E-EL-P2253	Porcine CXCR1 (CXC-Chemokine Receptor 1) ELISA Kit

E-EL-P2254	Porcine RBP4 (Retinol Binding Protein 4, Plasma) ELISA Kit
E-EL-P2255	Porcine TPPP (Tubulin Polymerization Promoting Protein) ELISA Kit
E-EL-P2256	Porcine TXLNa (Taxilin Alpha) ELISA Kit
E-EL-P2257	Porcine MMP-13 (Matrix Metalloproteinase 13) ELISA Kit
E-EL-P2258	Porcine CCK-33 (Cholecystokinin 33) ELISA Kit
E-EL-P2259	Porcine CCK-39 (Cholecystokinin 39) ELISA Kit
E-EL-P2260	Porcine CCK-58 (Cholecystokinin 58) ELISA Kit
E-EL-P2261	Porcine ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-P2262	Porcine Ta1 (Thymosin- α 1) ELISA Kit
E-EL-P2263	Porcine PS1 (Presenilin 1) ELISA Kit
E-EL-P2264	Porcine PACAP-38 (Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
E-EL-P2265	Porcine PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-P2267	Porcine GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
E-EL-P2268	Porcine HMGB-1(High mobility group protein B1)ELISA Kit
E-EL-P2269	Porcine IgM (Immunoglobulin M) ELISA Kit
E-EL-P2270	Porcine ApoA2 (Apolipoprotein A2) ELISA Kit
E-EL-P2272	Porcine INS (Insulin) ELISA Kit
E-EL-P2273	Porcine VMA(vanillylmandelic acid)ELISA kit
E-EL-MK0001	Monkey ApoC3(Apolipoprotein C3)ELISA Kit
E-EL-MK0002	Monkey IFN- γ (Interferon Gamma) ELISA Kit
E-EL-MK0003	Monkey sTNF-1 (Soluble Tumor Necrosis Factor 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0004	Monkey TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-MK0005	Monkey IL-12 (Interleukin 12) ELISA Kit
E-EL-MK0006	Monkey IL-2 (Interleukin 2) ELISA Kit
E-EL-MK0007	Monkey IL-10 (Interleukin-10) ELISA Kit
E-EL-MK0008	Monkey IL-6 (Interleukin 6) ELISA Kit
E-EL-MK0009	Monkey IL-8(Interleukin 8) ELISA Kit
E-EL-MK0010	Monkey sICAM-1 (Soluble Intercellular Adhesion Molecule-1) ELISA Kit
E-EL-MK0011	Monkey SELE (E-Selectin) ELISA Kit
E-EL-MK0013	Monkey SELP (P-Selectin) ELISA Kit
E-EL-MK0014	Monkey TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
E-EL-MK0016	Monkey IgG (Immunoglobulin G) ELISA Kit
E-EL-MK0017	Monkey CRP (C-reactive protein) ELISA Kit
E-EL-MK0018	Monkey MYO (Myoglobin) ELISA Kit
E-EL-MK0019	Monkey HP (Haptoglobin) ELISA Kit
E-EL-MK0020	Monkey SAA (Serum Amyloid A) ELISA Kit
E-EL-MK0021	Monkey IgE (Immunoglobulin E) ELISA Kit
E-EL-MK0022	Monkey Cys-C (Cystatin C) ELISA Kit
E-EL-MK0023	Monkey cTn-I/TNNI3 (cardiac Troponin I) ELISA Kit
E-EL-MK0024	Monkey ADM (Adrenomedullin) ELISA Kit
E-EL-MK0025	Monkey E-Cad (E-Cadherin) ELISA Kit

Cat.NO	Product Name
E-EL-MK0026	Monkey β 2-GP1 Ab IgM (β 2-Glycoprotein 1 Antibody IgM) ELISA Kit
E-EL-MK0027	Monkey ARO (Aromatase) ELISA Kit
E-EL-MK0028	Monkey Rock-2 (Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit
E-EL-MK0029	Monkey Bcl2L11 (Bcl-2 Like Protein 11) ELISA Kit
E-EL-MK0030	Monkey CRY α B (Crystallin Alpha B) ELISA Kit
E-EL-MK0031	Monkey COL6 α 3 (Collagen Type VI Alpha 3) ELISA Kit
E-EL-MK0032	Monkey SHC/SLP-76 (SH2-Containing Protein) ELISA Kit
E-EL-MK0033	Monkey CNN1 (Calponin 1, Basic, Smooth Muscle) ELISA Kit
E-EL-MK0035	Monkey CNP (C-type Natriuretic Peptide) ELISA Kit
E-EL-MK0036	Monkey CFP (Complement Factor P) ELISA Kit
E-EL-MK0037	Monkey ApoF (Apolipoprotein F) ELISA Kit
E-EL-MK0038	Monkey NOX5 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit
E-EL-MK0039	Monkey DR-70TM (Tumor Marker DR-70 For Lung Cancer) ELISA Kit
E-EL-MK0040	Monkey MCF (Macrophage Chemotatic Factor) ELISA Kit
E-EL-MK0041	Monkey CKMT1A (Creatine Kinase, Mitochondrial 1A) ELISA Kit
E-EL-MK0042	Monkey CRT (Calreticulin) ELISA Kit
E-EL-MK0043	Monkey CPP (Copeptin) ELISA Kit
E-EL-MK0044	Monkey ApoB100 (Apolipoprotein B100) ELISA Kit
E-EL-MK0045	Monkey CX3CL1 (Chemokine C-X3-C-Motif Ligand) ELISA Kit
E-EL-MK0046	Monkey TACR1 (Tachykinin Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0047	Monkey HO1 (Heme Oxygenase 1) ELISA Kit
E-EL-MK0048	Monkey EMILIN 1 (Elastin Microfibril Interface Located Protein 1) ELISA Kit
E-EL-MK0050	Monkey RAPGEF1 (Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit
E-EL-MK0051	Monkey NGB (Neuroglobin) ELISA Kit
E-EL-MK0052	Monkey FLN β (Filamin B, Beta) ELISA Kit
E-EL-MK0053	Monkey BAD (Bcl2 Associated Death Promoter) ELISA Kit
E-EL-MK0054	Monkey HO2 (Heme Oxygenase 2) ELISA Kit
E-EL-MK0055	Monkey PIAS3 (Protein Inhibitor of Activated STAT3) ELISA Kit
E-EL-MK0056	Monkey GH (Growth Hormone) ELISA Kit
E-EL-MK0057	Monkey GRIK2 (Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit
E-EL-MK0058	Monkey FOXC1 (Forkhead Box Protein C1) ELISA Kit
E-EL-MK0059	Monkey SUMO2 (Small Ubiquitin Related Modifier 2) ELISA Kit
E-EL-MK0061	Monkey PDCN (Podocin) ELISA Kit
E-EL-MK0063	Monkey CCSA-3 (Colon Cancer Specific Antigen-3) ELISA Kit
E-EL-MK0064	Monkey MT1M (Metallothionein 1M) ELISA Kit
E-EL-MK0065	Monkey TSP-2 (Thrombospondin-2) ELISA Kit
E-EL-MK0066	Monkey dDAVP (1-Desamino-8D-Arginine Vasopressin) ELISA Kit
E-EL-MK0067	Monkey TARC (Thymus Activation Regulated Chemokine) ELISA Kit
E-EL-MK0068	Monkey ACO1 (Aconitase 1) ELISA Kit
E-EL-MK0069	Monkey α -EP (α -Endomorphin) ELISA Kit

Cat.NO	Product Name
E-EL-MK0070	Monkey ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit
E-EL-MK0071	Monkey SPTAN1 (Alpha-Fodrin) ELISA Kit
E-EL-MK0072	Monkey Methylase ELISA Kit
E-EL-MK0073	Monkey DPD (Deoxypyridinoline) ELISA Kit
E-EL-MK0074	Monkey NAP3 (Neutrophil Activating Protein 3) ELISA Kit
E-EL-MK0076	Monkey KERA (Keratocan) ELISA Kit
E-EL-MK0077	Monkey VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit
E-EL-MK0078	Monkey CABIN1 (Calcineurin Binding Protein 1) ELISA Kit
E-EL-MK0081	Monkey GHAb (Anti-Growth Hormone Antibody) ELISA Kit
E-EL-MK0082	Monkey PGF (Placental Growth Factor) ELISA Kit
E-EL-MK0083	Monkey ANGPTL1 (Angiopoietin Like Protein 1) ELISA Kit
E-EL-MK0084	Monkey MLZE (Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit
E-EL-MK0086	Monkey RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
E-EL-MK0088	Monkey Lptn/XCL1 (Lymphotactin) ELISA Kit
E-EL-MK0089	Monkey TARS (Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit
E-EL-MK0090	Monkey AT (AntiThrombin) ELISA Kit
E-EL-MK0091	Monkey NUP155 (Nucleoporin 155kDa) ELISA Kit
E-EL-MK0094	Monkey CENPB (Centromere Protein B) ELISA Kit
E-EL-MK0095	Monkey FOXP3 (Forkhead Box Protein P3) ELISA Kit
E-EL-MK0096	Monkey F5 (Coagulation Factor V) ELISA Kit

Cat.NO	Product Name
E-EL-MK0097	Monkey PKCδ1 (Protein Kinase C Delta) ELISA Kit
E-EL-MK0099	Monkey IL-19 (Interleukin 19) ELISA Kit
E-EL-MK0100	Monkey GLRX3 (Glutaredoxin 3) ELISA Kit
E-EL-MK0101	Monkey PIINP (Procollagen II N-Terminal ProPeptide) ELISA Kit
E-EL-MK0102	Monkey NCAM/CD56 (Neural Cell Adhesion Molecule) ELISA Kit
E-EL-MK0103	Monkey CK (Creatine Kinase) ELISA Kit
E-EL-MK0104	Monkey ALK (Anaplastic Lymphoma Kinase) ELISA Kit
E-EL-MK0106	Monkey TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
E-EL-MK0107	Monkey MVP (Major Vault Protein) ELISA Kit
E-EL-MK0108	Monkey SYNPO (Synaptopodin) ELISA Kit
E-EL-MK0109	Monkey GSTκ1 (Glutathione S Transferase Kappa 1) ELISA Kit
E-EL-MK0110	Monkey CAPF (Cuticular Active Peptide Factor) ELISA Kit
E-EL-MK0111	Monkey Collagenase I ELISA Kit
E-EL-MK0112	Monkey VMA (Vanillylmandelic Acid) ELISA Kit
E-EL-MK0113	Monkey β Manase ELISA Kit
E-EL-MK0114	Monkey PLAP/ALPP (Placental Alkaline Phosphatase) ELISA Kit
E-EL-MK0115	Monkey IAA (Insulin Autoantibody) ELISA Kit
E-EL-MK0116	Monkey MPF (Maturation Promoting Factor) ELISA Kit
E-EL-MK0117	Monkey PP14 (Placental Protein14) ELISA Kit
E-EL-MK0118	Monkey MAPK (Mitogen Activated Protein Kinase) ELISA Kit

Cat.NO	Product Name
E-EL-MK0119	Monkey GMP-140 (Alpha-granular Membrane Protein) ELISA Kit
E-EL-MK0120	Monkey DAG1 (Dystrophin Associated Glycoprotein 1) ELISA Kit
E-EL-MK0121	Monkey PAEP (Progestagen-Associated Endometrial Protein) ELISA Kit
E-EL-MK0123	Monkey CX37 (Connexin 37) ELISA Kit
E-EL-MK0124	Monkey C/EBP ϵ (CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit
E-EL-MK0125	Monkey MAP2K6/MKK6 (Mitogen Activated Protein Kinase Kinase 6) ELISA Kit
E-EL-MK0127	Monkey PDK4 (Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit
E-EL-MK0128	Monkey γ GT1 (Gamma Glutamyltransferase 1) ELISA Kit
E-EL-MK0129	Monkey AP (Allopregnanolone) ELISA Kit
E-EL-MK0130	Monkey CYB561 (Cytochrome b-561) ELISA Kit
E-EL-MK0131	Monkey SSRP1 (Structure Specific Recognition Protein 1) ELISA Kit
E-EL-MK0132	Monkey Cyr61 (Cysteine Rich Protein, Angiogenic Inducer 61) ELISA Kit
E-EL-MK0133	Monkey ATF6 (Activating Transcription Factor 6) ELISA Kit
E-EL-MK0134	Monkey THRSP (Thyroid Hormone Responsive) ELISA Kit
E-EL-MK0135	Monkey LZM (Lysozyme) ELISA Kit
E-EL-MK0136	Monkey VEGFR1/Flt1 (Vasculoar Endothelial Growth Factor Receptor 1) ELISA Kit
E-EL-MK0137	Monkey AMA IgA (Anti-Myelin antibody IgA) ELISA Kit
E-EL-MK0138	Monkey FBLN5 (Fibulin 5) ELISA Kit
E-EL-MK0139	Monkey OXB (Orexin B) ELISA Kit
E-EL-MK0140	Monkey C1q-Ab (Anti-Complement 1q Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-MK0141	Monkey CNTF (Ciliary Neurotrophic Factor) ELISA Kit
E-EL-MK0142	Monkey IL-22 (Interleukin 22) ELISA Kit
E-EL-MK0143	Monkey PRB2 (Basic Salivary Proline Rich Protein 2) ELISA Kit
E-EL-MK0144	Monkey CDSN (Corneodesmosin) ELISA Kit
E-EL-MK0146	Monkey TNR (Tenascin R) ELISA Kit
E-EL-MK0147	Monkey LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit
E-EL-MK0148	Monkey S-100 (Soluble Protein-100) ELISA Kit
E-EL-MK0149	Monkey PARC/CCL18 (Pulmonary Activation Regulated Chemokine) ELISA Kit
E-EL-MK0150	Monkey SIGLEC10 (Sialic Acid Binding Ig Like Lectin 10) ELISA Kit
E-EL-MK0151	Monkey slgA (Secretory Immunoglobulin A) ELISA Kit
E-EL-MK0152	Monkey AGGF1 (Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit
E-EL-MK0153	Monkey SNRPD1 (Small Nuclear Ribonucleoprotein Polypeptide D1) ELISA Kit
E-EL-MK0154	Monkey ECE1 (Endothelin Converting Enzyme 1) ELISA Kit
E-EL-MK0156	Monkey FBG (Fibrinogen) ELISA Kit
E-EL-MK0157	Monkey SEMA4C (Semaphorin 4C) ELISA Kit
E-EL-MK0158	Monkey ANG II R1/AGTR1 (Angiotensin II Receptor 1) ELISA Kit
E-EL-MK0159	Monkey TNX (Tenascin X) ELISA Kit
E-EL-MK0160	Monkey PKM2 (Pyruvate Kinase, Muscle) ELISA Kit
E-EL-MK0161	Monkey LIFR (Leukemia Inhibitory Factor Receptor) ELISA Kit
E-EL-MK0162	Monkey TRH (Thyrotropin Releasing Hormone) ELISA Kit

Cat.NO	Product Name
E-EL-MK0163	Monkey TSGF (Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
E-EL-MK0164	Monkey CNTN4 (Contactin 4) ELISA Kit
E-EL-MK0165	Monkey TRF (Transferrin) ELISA Kit
E-EL-MK0167	Monkey P-Cadherin (Cadherin, Placental) ELISA Kit
E-EL-MK0168	Monkey VIPR1 (Vasoactive Intestinal Peptide Receptor) ELISA Kit
E-EL-MK0169	Monkey NCAD (Neural Cadherin) ELISA Kit
E-EL-MK0170	Monkey C9 (Complement Component 9) ELISA Kit
E-EL-MK0171	Monkey MMP-10 (Matrix Metalloproteinase 10) ELISA Kit
E-EL-MK0172	Monkey FSTL1 (Follistatin Like Protein 1) ELISA Kit
E-EL-MK0173	Monkey PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit
E-EL-MK0175	Monkey NT-4 (Neurotrophin 4) ELISA Kit
E-EL-MK0176	Monkey MUC20 (Mucin 20) ELISA Kit
E-EL-MK0177	Monkey CLS (Cardiolipin Synthase) ELISA Kit
E-EL-MK0179	Monkey RARS (Arginyl tRNA Synthetase) ELISA Kit
E-EL-MK0180	Monkey CTGF (Connective Tissue Growth Factor) ELISA Kit
E-EL-MK0181	Monkey LIF (Leukemia Inhibitory Factor) ELISA Kit
E-EL-MK0182	Monkey COL4α1 (Collagen Type IV Alpha 1) ELISA Kit
E-EL-MK0183	Monkey PMP3 (Peroxisomal Membrane Protein 3) ELISA Kit
E-EL-MK0184	Monkey EK (Eepidermal Keratin) ELISA Kit
E-EL-MK0185	Monkey cADPRH/CD38 (Cyclic ADP Hydrolase) ELISA Kit

Cat.NO	Product Name
E-EL-MK0187	Monkey NF (Neurofilament) ELISA Kit
E-EL-MK0188	Monkey VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
E-EL-MK0189	Monkey AGA (Anti-Gliadin Antibody) ELISA Kit
E-EL-MK0190	Monkey IGF-2 (Insulin Like Growth Factor 2) ELISA Kit
E-EL-MK0191	Monkey STAT6B (Signal Transducer and Activator of Transcription 6B) ELISA Kit
E-EL-MK0192	Monkey C4a (Complement Component 4a) ELISA Kit
E-EL-MK0193	Monkey EFNA3 (Ephrin A3) ELISA Kit
E-EL-MK0194	Monkey NAGLU (N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit
E-EL-MK0195	Monkey SPOP (Speckle Type POZ Protein) ELISA Kit
E-EL-MK0196	Monkey MMP-12 (Matrix Metalloproteinase 12) ELISA Kit
E-EL-MK0197	Monkey CDK8 (Cyclin Dependent Kinase 8) ELISA Kit
E-EL-MK0198	Monkey EFNA5 (Ephrin A5) ELISA Kit
E-EL-MK0200	Monkey TNNI2 (Troponin I Type 2, Fast Skeletal) ELISA Kit
E-EL-MK0201	Monkey TP (Thymidine Phosphorylase) ELISA Kit
E-EL-MK0202	Monkey CCA (Colon Cancer Antigen) ELISA Kit
E-EL-MK0203	Monkey NUP37 (Nucleoporin 37kDa) ELISA Kit
E-EL-MK0205	Monkey GnRH (Gonadotropin Releasing Hormone) ELISA Kit
E-EL-MK0206	Monkey PKC (Protein Kinase C) ELISA Kit
E-EL-MK0207	Monkey STAT1 (Signal Transducer and Activator of Transcription 1) ELISA Kit
E-EL-MK0209	Monkey ATF-3 (Activating Transcription Factor 3) ELISA Kit

Cat.NO	Product Name
E-EL-MK0210	Monkey β -LPH (β -Lipotropic Hormone) ELISA Kit
E-EL-MK0211	Monkey GSK3 β (Glycogen Synthase Kinase 3 Beta) ELISA Kit
E-EL-MK0212	Monkey SPARCL1 (SPARC Like Protein 1) ELISA Kit
E-EL-MK0214	Monkey COL6 (Collagen Type VI) ELISA Kit
E-EL-MK0215	Monkey C1qR1 (Complement Component 1q Receptor) ELISA Kit
E-EL-MK0216	Monkey TTN (Titin) ELISA Kit
E-EL-MK0217	Monkey CYP1A2 (Cytochrome P450 1A2) ELISA Kit
E-EL-MK0219	Monkey PRSS33 (Protease, Serine 33) ELISA Kit
E-EL-MK0221	Monkey γ -ECS (γ -Glutamyl Systeine Synthetase) ELISA Kit
E-EL-MK0222	Monkey SOX18 (Sex Determining Region Y Box Protein 18) ELISA Kit
E-EL-MK0223	Monkey 5-HIAA (5-Hydroxyindolacetic Acid) ELISA Kit
E-EL-MK0224	Monkey MMP-1 (Matrix Metalloproteinase 1) ELISA Kit
E-EL-MK0225	Monkey COL5 α 2 (Collagen Type V Alpha 2) ELISA Kit
E-EL-MK0226	Monkey BlysR (B Lymphocyte Stimulator Receptor) ELISA Kit
E-EL-MK0227	Monkey GDNF (Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
E-EL-MK0228	Monkey LUM (Lumican) ELISA Kit
E-EL-MK0229	Monkey ERK3 (Extracellular Signal Regulated Kinase 3) ELISA Kit
E-EL-MK0230	Monkey HbA1c (Glycosylated Hemoglobin/Hemoglobin A1c) ELISA Kit
E-EL-MK0231	Monkey CLDN4 (Claudin 4) ELISA Kit
E-EL-MK0232	Monkey GPX1 (Glutathione Peroxidase 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0233	Monkey TPO-Ab (Anti-Thyroid-Peroxidase Antibody) ELISA Kit
E-EL-MK0234	Monkey ACD (Adrenocortical Dysplasia Homolog) ELISA Kit
E-EL-MK0235	Monkey LP-a (Lipoprotein a) ELISA Kit
E-EL-MK0236	Monkey RPL10 (Ribosomal Protein L10) ELISA Kit
E-EL-MK0237	Monkey MBP (Myelin Basic Protein) ELISA Kit
E-EL-MK0238	Monkey PRSS23 (Protease, Serine 23) ELISA Kit
E-EL-MK0239	Monkey BMP-10 (Bone Morphogenetic Protein 10) ELISA Kit
E-EL-MK0240	Monkey PP (Protein Phosphatase) ELISA Kit
E-EL-MK0241	Monkey NFRkB (Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit
E-EL-MK0242	Monkey MAdCAM1 (Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit
E-EL-MK0243	Monkey SRP9 (Signal Recognition Particle 9kDa) ELISA Kit
E-EL-MK0244	Monkey UBE2A (Ubiquitin Conjugating Enzyme E2A) ELISA Kit
E-EL-MK0245	Monkey ECF/CCL11 (Eosinophil Chemotactic Factor) ELISA Kit
E-EL-MK0246	Monkey CDK1 (Cyclin Dependent Kinase 1) ELISA Kit
E-EL-MK0247	Monkey CC16 (Clara Cell 16KD Protein) ELISA Kit
E-EL-MK0248	Monkey ANG I (Angiotensin I) ELISA Kit
E-EL-MK0249	Monkey ECM1 (Extracellular Matrix Protein 1) ELISA Kit
E-EL-MK0250	Monkey CTACK (Cutaneous T-cell Attracting Chemokine) ELISA Kit
E-EL-MK0251	Monkey CCKAR (Cholecystokinin A Receptor) ELISA Kit
E-EL-MK0252	Monkey CLDN9 (Claudin 9) ELISA Kit

Cat.NO	Product Name
E-EL-MK0253	Monkey NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) ELISA Kit
E-EL-MK0254	Monkey PRF1 (Perforin 1/Pore-Forming Protein) ELISA Kit
E-EL-MK0255	Monkey LATS (Long Acting Thyroid Stimulator) ELISA Kit
E-EL-MK0256	Monkey FSHR (Follicle Stimulating Hormone Receptor) ELISA Kit
E-EL-MK0257	Monkey ACP3/PAP (Prostatic Acid Phosphatase 3) ELISA Kit
E-EL-MK0258	Monkey PDH β (Pyruvate Dehydrogenase Beta) ELISA Kit
E-EL-MK0259	Monkey Neopterin ELISA Kit
E-EL-MK0260	Monkey sFMC (Soluble Fibrin Monomer Complex) ELISA Kit
E-EL-MK0261	Monkey FAS/CD95/TNFRSF6 (Factor Related Apoptosis) ELISA Kit
E-EL-MK0262	Monkey LDHC (Lactate Dehydrogenase C) ELISA Kit
E-EL-MK0264	Monkey MYC (V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit
E-EL-MK0265	Monkey NAIP (Neuronal Apoptosis Inhibitory Protein) ELISA Kit
E-EL-MK0266	Monkey AP62A (Anti-P62 Antibody) ELISA Kit
E-EL-MK0267	Monkey PTN (Pleiotrophin) ELISA Kit
E-EL-MK0268	Monkey FBLN4 (Fibulin 4) ELISA Kit
E-EL-MK0269	Monkey PAR2 (Protease Activated Receptor 2) ELISA Kit
E-EL-MK0270	Monkey FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit
E-EL-MK0271	Monkey BRAK/CXCL14 (Breast and Kidney Expressed Chemokine) ELISA Kit
E-EL-MK0273	Monkey snRNP/Sm (Small Nuclear Ribonucleoprotein) ELISA Kit
E-EL-MK0274	Monkey AXIN2 (Axis Inhibition Protein 2) ELISA Kit

Cat.NO	Product Name
E-EL-MK0277	Monkey GTF2F1 (General Transcription Factor IIF, Polypeptide 1) ELISA Kit
E-EL-MK0278	Monkey CDK2 (Cyclin Dependent Kinase 2) ELISA Kit
E-EL-MK0280	Monkey PHLDA2 (Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit
E-EL-MK0281	Monkey DUB (Deubiquitinating Enzyme) ELISA Kit
E-EL-MK0282	Monkey HSPG (Heparan Sulfate Proteoglycan) ELISA Kit
E-EL-MK0283	Monkey IGFBP-1 (Insulin Like Growth Factor Binding Protein 1) ELISA Kit
E-EL-MK0286	Monkey PON3 (Paraoxonase 3) ELISA Kit
E-EL-MK0287	Monkey PAPPA2 (Pappalysin 2) ELISA Kit
E-EL-MK0288	Monkey MYL1 (Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit
E-EL-MK0289	Monkey CA153 (Carbohydrate Antigen 153) ELISA Kit
E-EL-MK0290	Monkey CYSLTR2 (Cysteinyl Leukotriene Receptor 2) ELISA Kit
E-EL-MK0292	Monkey GAD-Ab (Anti-Glutamic Acid Decarboxylase antibody) ELISA Kit
E-EL-MK0293	Monkey CENPA (Centromere Protein A) ELISA Kit
E-EL-MK0294	Monkey PTPRJ (Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit
E-EL-MK0295	Monkey PIIINP (N-Terminal Procollagen III Propeptide) ELISA Kit
E-EL-MK0296	Monkey PTPRN (Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
E-EL-MK0297	Monkey TGF-β2 (Transforming Growth Factor β2) ELISA Kit
E-EL-MK0298	Monkey CDC23 (Cell Division Cycle Protein 23) ELISA Kit
E-EL-MK0299	Monkey LILRB4 (Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit
E-EL-MK0300	Monkey sgp130 (soluble Glucoprotein 130) ELISA Kit

Cat.NO	Product Name
E-EL-MK0301	Monkey MS (Melatonin Sulfate) ELISA Kit
E-EL-MK0302	Monkey CAM (Calmodulin) ELISA Kit
E-EL-MK0303	Monkey AP (Aprotinin) ELISA Kit
E-EL-MK0304	Monkey TR (Targeted Ribonuclease) ELISA Kit
E-EL-MK0306	Monkey IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit
E-EL-MK0307	Monkey CA3 (Carbonic Anhydrase III, Muscle Specific) ELISA Kit
E-EL-MK0308	Monkey MEC/CCL28 (Mucosae Associated Epithelia Chemokine) ELISA Kit
E-EL-MK0310	Monkey PCDH γ A2 (Protocadherin Gamma A2) ELISA Kit
E-EL-MK0311	Monkey C5 (Complement Component 5) ELISA Kit
E-EL-MK0312	Monkey RGS4 (Regulator of G Protein Signaling 4) ELISA Kit
E-EL-MK0313	Monkey CDT (Carbohydrate Deficient Transferrin) ELISA Kit
E-EL-MK0314	Monkey CYSLTR1 (Cysteinyl Leukotriene Receptor 1) ELISA Kit
E-EL-MK0315	Monkey NUP107 (Nucleoporin 107kDa) ELISA Kit
E-EL-MK0316	Monkey REG1 β (Regenerating Islet Derived Protein 1 Beta) ELISA Kit
E-EL-MK0317	Monkey NLR (NOD-Like Receptor9) ELISA Kit
E-EL-MK0319	Monkey CaN (Calcineurin) ELISA Kit
E-EL-MK0320	Monkey VCAM-1/CD106 (soluble Vascular Cell Adhesion Molecule 1) ELISA Kit
E-EL-MK0321	Monkey DKK2 (Dickkopf Related Protein 2) ELISA Kit
E-EL-MK0323	Monkey GNLY (Granulysin) ELISA Kit
E-EL-MK0324	Monkey HJV/HFE2 (Hemojuvelin/Hemochromatosis Type 2) ELISA Kit

Cat.NO	Product Name
E-EL-MK0325	Monkey Rac-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit
E-EL-MK0326	Monkey ASPN (Asporin) ELISA Kit
E-EL-MK0327	Monkey EJ/GlyRS (Anti-EJ-Antibody) ELISA Kit
E-EL-MK0328	Monkey DPYSL2 (Dihydropyrimidinase Like Protein 2) ELISA Kit
E-EL-MK0330	Monkey DHEA-S (Dehydroepiandrosterone Sulfate) ELISA Kit
E-EL-MK0331	Monkey FGF4 (Fibroblast Growth Factor 4) ELISA Kit
E-EL-MK0332	Monkey PP13 (Placental Protein13) ELISA Kit
E-EL-MK0334	Monkey CLEC4E (C-Type Lectin Domain Family 4 Member E) ELISA Kit
E-EL-MK0336	Monkey TGF- β 3 (Transforming Growth Factor β 3) ELISA Kit
E-EL-MK0337	Monkey CNX (Calnexin) ELISA Kit
E-EL-MK0338	Monkey F1+2 (Prothrombin Fragment 1+2) ELISA Kit
E-EL-MK0340	Monkey GSANA (Granulocyte Specific Anti-Nuclear Antibody) ELISA Kit
E-EL-MK0341	Monkey sAC (Soluble Adenylate Cyclase) ELISA Kit
E-EL-MK0342	Monkey TWEAK (Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit
E-EL-MK0343	Monkey EloA (Elongin A) ELISA Kit
E-EL-MK0344	Monkey IL-2R α (Interleukin 2 Receptor Alpha) ELISA Kit
E-EL-MK0345	Monkey NUP214 (Nucleoporin 214kDa) ELISA Kit
E-EL-MK0346	Monkey MUC7 (Mucin 7) ELISA Kit
E-EL-MK0347	Monkey NAGase (N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
E-EL-MK0348	Monkey MAU (Microalbuminuria) ELISA Kit

Cat.NO	Product Name
E-EL-MK0349	Monkey I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit
E-EL-MK0350	Monkey CCP (Complement Control Protein) ELISA Kit
E-EL-MK0351	Monkey MUC3 (Mucin 3) ELISA Kit
E-EL-MK0352	Monkey TPA (Tissue Polypeptide Antigen) ELISA Kit
E-EL-MK0353	Monkey TG (Thyroglobulin) ELISA Kit
E-EL-MK0354	Monkey α MSH (Alpha-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-MK0355	Monkey ABCA1(ATP-binding cassette transporter A1) ELISA Kit
E-EL-MK0356	Monkey CX43 (Connexin 43) ELISA Kit
E-EL-MK0357	Monkey MANF (Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit
E-EL-MK0358	Monkey LT β (Lymphotoxin Beta) ELISA Kit
E-EL-MK0359	Monkey FDP (Fibrinogen Degradation Product) ELISA Kit
E-EL-MK0361	Monkey AQP-4 (Aquaporin 4) ELISA Kit
E-EL-MK0362	Monkey EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit
E-EL-MK0363	Monkey LGALS3BP (Lectin Galactoside Binding, Soluble 3 Binding Protein) ELISA Kit
E-EL-MK0364	Monkey NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
E-EL-MK0365	Monkey CREB (Cyclic AMP Response Element Binding Protein) ELISA Kit
E-EL-MK0366	Monkey β -TG (β -Thromboglobulin) ELISA Kit
E-EL-MK0367	Monkey ANP (Atrial Natriuretic Peptide) ELISA Kit
E-EL-MK0368	Monkey TACE (TNF α Converting Enzyme) ELISA Kit
E-EL-MK0369	Monkey PKD2 (Protein Kinase D2) ELISA Kit

Cat.NO	Product Name
E-EL-MK0370	Monkey BECN1 (Beclin 1) ELISA Kit
E-EL-MK0371	Monkey RHOB (Ras Homolog Gene Family, Member B) ELISA Kit
E-EL-MK0373	Monkey GARS (Glycyl tRNA Synthetase) ELISA Kit
E-EL-MK0374	Monkey GDF15 (Growth Differentiation Factor 15) ELISA Kit
E-EL-MK0376	Monkey CYPA (Cyclophilin A) ELISA Kit
E-EL-MK0377	Monkey MART/Melan-A (Melanoma Marker A) ELISA Kit
E-EL-MK0378	Monkey EMA (Epithelial Membrane Antigen) ELISA Kit
E-EL-MK0379	Monkey COL11α1 (Collagen Type XI Alpha 1) ELISA Kit
E-EL-MK0380	Monkey PLB (Phospholipase B) ELISA Kit
E-EL-MK0381	Monkey CRN (Corin) ELISA Kit
E-EL-MK0383	Monkey BLC (B-Lymphocyte Chemoattractant) ELISA Kit
E-EL-MK0384	Monkey PRCP (Prolylcarboxypeptidase) ELISA Kit
E-EL-MK0385	Monkey D2D (D-Dimer) ELISA Kit
E-EL-MK0386	Monkey GLUT2 (Glucose Transporter 2) ELISA Kit
E-EL-MK0387	Monkey ESF (Erythropoiesis Stimulating Factor) ELISA Kit
E-EL-MK0388	Monkey RB1 (Retinoblastoma Protein 1) ELISA Kit
E-EL-MK0389	Monkey CLDN3 (Claudin 3) ELISA Kit
E-EL-MK0390	Monkey IGFBP-3 (Insulin Like Growth Factor Binding Protein 3) ELISA Kit
E-EL-MK0391	Monkey PROZ (Protein Z) ELISA Kit
E-EL-MK0392	Monkey PRSS8 (Protease, Serine 8) ELISA Kit

Cat.NO	Product Name
E-EL-MK0394	Monkey EPCR (Endothelial Protein C Receptor) ELISA Kit
E-EL-MK0395	Monkey SFRP2 (Secreted Frizzled Related Protein 2) ELISA Kit
E-EL-MK0396	Monkey TXB2 (Thromboxane B2) ELISA Kit
E-EL-MK0397	Monkey IgA (Immunoglobulin A) ELISA Kit
E-EL-MK0399	Monkey SSFA2 (Sperm Specific Antigen 2) ELISA Kit
E-EL-MK0400	Monkey SIGLEC8 (Sialic Acid Binding Ig Like Lectin 8) ELISA Kit
E-EL-MK0401	Monkey FP (Prostaglandin F Receptor) ELISA Kit
E-EL-MK0402	Monkey FGFBP1 (Fibroblast Growth Factor Binding Protein 1) ELISA Kit
E-EL-MK0403	Monkey Cx (Connexin) ELISA Kit
E-EL-MK0405	Monkey PGA (Pepsinogen A) ELISA Kit
E-EL-MK0406	Monkey mIgM (Membrane IgM) ELISA Kit
E-EL-MK0407	Monkey LPL (Lipoprotein Lipase) ELISA Kit
E-EL-MK0408	Monkey ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit
E-EL-MK0409	Monkey c-jun ELISA Kit
E-EL-MK0410	Monkey Fcγ (Fc Fragment of IgG) ELISA Kit
E-EL-MK0411	Monkey ATF-2 (Activating Transcription Factor 2) ELISA Kit
E-EL-MK0412	Monkey NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit
E-EL-MK0413	Monkey MAG (Myelin Associated Glycoprotein) ELISA Kit
E-EL-MK0414	Monkey NADPH (Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
E-EL-MK0415	Monkey Ras-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0416	Monkey EPYC (Epiphycan) ELISA Kit
E-EL-MK0418	Monkey TS/TYMS (Thymidylate Synthetase) ELISA Kit
E-EL-MK0419	Monkey ASGR2 (Asialoglycoprotein Receptor 2) ELISA Kit
E-EL-MK0420	Monkey ARRβ1 (Arrestin Beta 1) ELISA Kit
E-EL-MK0421	Monkey IL-17 (Interleukin 17) ELISA Kit
E-EL-MK0422	Monkey IGFBP-5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit
E-EL-MK0423	Monkey NFYC (Nuclear Transcription Factor Y Gamma) ELISA Kit
E-EL-MK0424	Monkey MUC6 (Mucin 6) ELISA Kit
E-EL-MK0425	Monkey TLR2 (Toll-Like Receptor 2) ELISA Kit
E-EL-MK0426	Monkey UCHL1 (Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
E-EL-MK0427	Monkey TSC22 (Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
E-EL-MK0428	Monkey Collagenase II ELISA Kit
E-EL-MK0429	Monkey AQP-3 (Aquaporin 3) ELISA Kit
E-EL-MK0430	Monkey SDF-4 (Stromal Cell Derived Factor 4) ELISA Kit
E-EL-MK0431	Monkey SPB (Pulmonary Surfactant Associated Protein B) ELISA Kit
E-EL-MK0433	Monkey ANXA5 (Annexin A5) ELISA Kit
E-EL-MK0434	Monkey CRLF1 (Cytokine Receptor Like Factor 1) ELISA Kit
E-EL-MK0435	Monkey PAR1 (Protease Activated Receptor 1) ELISA Kit
E-EL-MK0436	Monkey PLTP (Phospholipid Transfer Protein) ELISA Kit
E-EL-MK0437	Monkey GATA1 (GATA Binding Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0438	Monkey TU3A (TU3A Protein) ELISA Kit
E-EL-MK0441	Monkey PDGFA (Platelet Derived Growth Factor Subunit A) ELISA Kit
E-EL-MK0442	Monkey OJ/IleRS (Anti-OJ-Antibody) ELISA Kit
E-EL-MK0443	Monkey COL3α1 (Collagen Type III Alpha 1) ELISA Kit
E-EL-MK0444	Monkey PR3 (Proteinase 3) ELISA Kit
E-EL-MK0445	Monkey MEPE (Matrix Extracellular Phosphoglycoprotein) ELISA Kit
E-EL-MK0446	Monkey CASP1 (Caspase 1) ELISA Kit
E-EL-MK0447	Monkey FLRT2 (Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit
E-EL-MK0448	Monkey Anti-GM1 (Anti-Ganglioside M1 Antibody) ELISA Kit
E-EL-MK0449	Monkey NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
E-EL-MK0451	Monkey POMC (Pro-Opiomelanocortin) ELISA Kit
E-EL-MK0452	Monkey DPYSL3 (Dihydropyrimidinase Like Protein 3) ELISA Kit
E-EL-MK0453	Monkey ATF-4 (Activating Transcription Factor 4) ELISA Kit
E-EL-MK0454	Monkey ConA (Concanavalin A) ELISA Kit
E-EL-MK0456	Monkey EPOR (Erythropoietin Receptor) ELISA Kit
E-EL-MK0457	Monkey PROS (Protein S) ELISA Kit
E-EL-MK0458	Monkey Des (Desmin) ELISA Kit
E-EL-MK0459	Monkey SDF-1/CXCL12 (Stromal Cell Derived Factor 1) ELISA Kit
E-EL-MK0460	Monkey CXCR3 (Chemokine CXC Receptor 3) ELISA Kit
E-EL-MK0461	Monkey D5R/DRD5 (Dopamine Receptor D5) ELISA Kit

Cat.NO	Product Name
E-EL-MK0462	Monkey FABP8 (Fatty Acid Binding Protein 8, Myelin) ELISA Kit
E-EL-MK0463	Monkey TH (Tyrosine Hydroxylase) ELISA Kit
E-EL-MK0464	Monkey CEBP γ (CCAAT/Enhancer Binding Protein Gamma) ELISA Kit
E-EL-MK0465	Monkey GDI1 (GDP Dissociation Inhibitor 1) ELISA Kit
E-EL-MK0466	Monkey MFAP4 (Microfibrillar Associated Protein 4) ELISA Kit
E-EL-MK0467	Monkey LpPLA2 (Phospholipase A2, Lipoprotein Associated) ELISA Kit
E-EL-MK0468	Monkey AP36 (Apelin 36) ELISA Kit
E-EL-MK0469	Monkey α 1-AT (Alpha 1-Antitrypsin) ELISA Kit
E-EL-MK0470	Monkey DPYD (Dihydropyrimidine Dehydrogenase) ELISA Kit
E-EL-MK0472	Monkey DS (Dermatan Sulfate) ELISA Kit
E-EL-MK0473	Monkey OGP (Osteogenic Growth Peptide) ELISA Kit
E-EL-MK0474	Monkey GSP (Glycosylated Serum Protein) ELISA Kit
E-EL-MK0476	Monkey APAF-1 (Apoptosis Protease Activating Factor-1) ELISA Kit
E-EL-MK0477	Monkey A β 1-40 (Amyloid Beta Peptide 1-40) ELISA Kit
E-EL-MK0478	Monkey C3bR (Complement Fragment 3b Receptor) ELISA Kit
E-EL-MK0479	Monkey Dyn (Big Dynorphin) ELISA Kit
E-EL-MK0480	Monkey PLSCR2 (Phospholipid Scramblase 2) ELISA Kit
E-EL-MK0481	Monkey HYOU1 (Hypoxia Up Regulated 1) ELISA Kit
E-EL-MK0482	Monkey AIG1 (Androgen Induced Protein 1) ELISA Kit
E-EL-MK0484	Monkey COL15 (Collagen Type XV) ELISA Kit

Cat.NO	Product Name
E-EL-MK0485	Monkey CD40L/TNFSF5 (Cluster of Differentiation 40 Ligand) ELISA Kit
E-EL-MK0487	Monkey EIPS2 (Estrogen Induced Protein pS2) ELISA Kit
E-EL-MK0488	Monkey F12 (Coagulation Factor XII) ELISA Kit
E-EL-MK0489	Monkey COLEC11 (Collectin Sub-Family Member 11) ELISA Kit
E-EL-MK0490	Monkey TAF2 (TATA Box Binding Protein Associated Factor 2) ELISA Kit
E-EL-MK0491	Monkey Grb2 (Growth Factor Receptor Bound Protein 2) ELISA Kit
E-EL-MK0492	Monkey DBP (Vitamin D Binding Protein) ELISA Kit
E-EL-MK0493	Monkey ATXN2 (Ataxin 2) ELISA Kit
E-EL-MK0494	Monkey CCL16 (Chemokine C-C-Motif Ligand 16) ELISA Kit
E-EL-MK0496	Monkey TCR (Tcell Receptor) ELISA Kit
E-EL-MK0497	Monkey TR (Testosterone Receptor) ELISA Kit
E-EL-MK0498	Monkey MCAM (Melanoma Cell Adhesion Molecule) ELISA Kit
E-EL-MK0501	Monkey C3 (Complement Component 3) ELISA Kit
E-EL-MK0502	Monkey MF/MPF (Maturation Promoting Factor) ELISA Kit
E-EL-MK0503	Monkey β -Glu (Beta Glucosidase) ELISA Kit
E-EL-MK0504	Monkey CNTNAP2 (Contactin Associated Protein Like Protein 2) ELISA Kit
E-EL-MK0505	Monkey ST2 (Syntenin 2) ELISA Kit
E-EL-MK0507	Monkey PPAR- γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
E-EL-MK0508	Monkey PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit
E-EL-MK0509	Monkey MUC21 (Mucin 21) ELISA Kit

Cat.NO	Product Name
E-EL-MK0510	Monkey TI-Ag (Thymus Independent Antigen) ELISA Kit
E-EL-MK0511	Monkey α HBDH (α -Hydroxybutyrate Dehydrogenase) ELISA Kit
E-EL-MK0512	Monkey LH (Luteinizing Hormone) ELISA Kit
E-EL-MK0513	Monkey PI (Proinsulin) ELISA Kit
E-EL-MK0514	Monkey PAFAHIB3 (Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit
E-EL-MK0515	Monkey NFKBIE (Nuclear Factor Of Kappa Light Polypeptide Gene Enhancer In B-cells Inhibitor Epsilon) ELISA Kit
E-EL-MK0516	Monkey CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit
E-EL-MK0517	Monkey GAP (GTPase Activating Protein) ELISA Kit
E-EL-MK0518	Monkey Bcl-3 (B-cell Leukemia/Lymphoma 3) ELISA Kit
E-EL-MK0519	Monkey C1r (Complement Component 1R) ELISA Kit
E-EL-MK0520	Monkey AFP-L3 (Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit
E-EL-MK0521	Monkey LT β R (Lymphotoxin Beta Receptor) ELISA Kit
E-EL-MK0522	Monkey STAT6 (Signal Transducer and Activator of Transcription 6) ELISA Kit
E-EL-MK0523	Monkey DPP10 (Dipeptidyl Peptidase X) ELISA Kit
E-EL-MK0524	Monkey TEP1 (Telomerase Associated Protein 1) ELISA Kit
E-EL-MK0525	Monkey AKAP1 (A Kinase Anchor Protein 1) ELISA Kit
E-EL-MK0526	Monkey FETUA (Fetuin A) ELISA Kit
E-EL-MK0527	Monkey D1R/DRD1 (Dopamine Receptor D1) ELISA Kit
E-EL-MK0528	Monkey TMS β 4 (Thymosin Beta 4) ELISA Kit
E-EL-MK0529	Monkey GLUT1 (Glucose Transporter 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0530	Monkey MASP2 (Mannan Associated Serine Protease 2) ELISA Kit
E-EL-MK0531	Monkey TMSβ10 (Thymosin Beta 10) ELISA Kit
E-EL-MK0532	Monkey TrxR1 (Thioredoxin Reductase 1) ELISA Kit
E-EL-MK0533	Monkey FMOD (Fibromodulin) ELISA Kit
E-EL-MK0535	Monkey ACTH (Adrenocorticotrophic Hormone) ELISA Kit
E-EL-MK0536	Monkey PT (Prothrombin) ELISA Kit
E-EL-MK0538	Monkey PDGFRα (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit
E-EL-MK0539	Monkey PP (Placental Protein) ELISA Kit
E-EL-MK0540	Monkey FTL (Ferritin, Light Polypeptide) ELISA Kit
E-EL-MK0543	Monkey TGF-β1 (Transforming Growth Factor β1) ELISA Kit
E-EL-MK0544	Monkey MYH10 (Myosin Heavy Chain 10, Non Muscle) ELISA Kit
E-EL-MK0545	Monkey GLα (Galactosidase Alpha) ELISA Kit
E-EL-MK0546	Monkey FRS3 (Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
E-EL-MK0547	Monkey M-AChR M3 (Muscarinic Acetylcholine Receptor M3) ELISA Kit
E-EL-MK0548	Monkey PAI2 (Plasminogen Activator Inhibitor 2)ELISA Kit
E-EL-MK0549	Monkey CA6 (Carbonic Anhydrase VI) ELISA Kit
E-EL-MK0550	Monkey FPB (Fibrinopeptide B) ELISA Kit
E-EL-MK0551	Monkey CK-17/KRT17 (Cytokeratin 17) ELISA Kit
E-EL-MK0552	Monkey plkBα (Phospho-Inhibitory Subunit Of NF Kappa B Alpha) ELISA Kit
E-EL-MK0553	Monkey MAN1A1 (Mannosidase Alpha Class 1A Member 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0554	Monkey MLCK/MYLK (Myosin Light Chain Kinase) ELISA Kit
E-EL-MK0555	Monkey CTSC (Cathepsin C) ELISA Kit
E-EL-MK0556	Monkey SIGLEC12 (Sialic Acid Binding Ig Like Lectin 12) ELISA Kit
E-EL-MK0557	Monkey AMY1 (Amylase Alpha 1, Salivary) ELISA Kit
E-EL-MK0558	Monkey iPLA2 (Phospholipidase A2, Calcium Independent) ELISA Kit
E-EL-MK0559	Monkey CNTN3 (Contactin 3) ELISA Kit
E-EL-MK0560	Monkey CCSA-4 (Colon Cancer Specific Antigen-4) ELISA Kit
E-EL-MK0562	Monkey IL-1ra (Interleukin 1 Receptor Antagonist) ELISA Kit
E-EL-MK0563	Monkey ApoA5 (Apolipoprotein A5) ELISA Kit
E-EL-MK0564	Monkey TF (Transfer Factor) ELISA Kit
E-EL-MK0565	Monkey ABCG1 (ATP Binding Cassette Transporter G1) ELISA Kit
E-EL-MK0566	Monkey PCT (Procalcitonin) ELISA Kit
E-EL-MK0569	Monkey FHA (Filamentous Hemagglutinin) ELISA Kit
E-EL-MK0570	Monkey SMM (Smooth Muscle Myosin) ELISA Kit
E-EL-MK0571	Monkey OTR (Oxytocin Receptor) ELISA Kit
E-EL-MK0572	Monkey EGR1 (Early Growth Response Protein 1) ELISA Kit
E-EL-MK0573	Monkey ORM2 (Orosomucoid 2) ELISA Kit
E-EL-MK0574	Monkey AGR2 (Anterior Gradient Protein 2) ELISA Kit
E-EL-MK0575	Monkey ACE (Angiotensin I Converting Enzyme) ELISA Kit
E-EL-MK0576	Monkey CD28 (Cluster of Differentiation 28) ELISA Kit

Cat.NO	Product Name
E-EL-MK0577	Monkey IgG2 (Immunoglobulin G2) ELISA Kit
E-EL-MK0578	Monkey FABP5 (Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
E-EL-MK0579	Monkey PK2 (Prokineticin 2) ELISA Kit
E-EL-MK0580	Monkey CD3d (Cluster of Differentiation 3d) ELISA Kit
E-EL-MK0581	Monkey c-myc (c-myc Oncogene Product) ELISA Kit
E-EL-MK0582	Monkey STK39 (Serine/Threonine Kinase 39) ELISA Kit
E-EL-MK0583	Monkey PDGF (Platelet-Derived Growth Factor) ELISA Kit
E-EL-MK0584	Monkey Pg (Progesterone) ELISA Kit
E-EL-MK0585	Monkey TSTA (Tumor Specific Transplantation Antigen) ELISA Kit
E-EL-MK0587	Monkey EPCR (Endothelial Protein C Receptor) ELISA Kit
E-EL-MK0588	Monkey PYGB (Glycogen Phosphorylase, Brain) ELISA Kit
E-EL-MK0589	Monkey GSTα5 (Glutathione S Transferase Alpha 5) ELISA Kit
E-EL-MK0590	Monkey ALDH (Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit
E-EL-MK0591	Monkey HPR (Haptoglobin Related Protein) ELISA Kit
E-EL-MK0592	Monkey GPC2 (Glypican 2) ELISA Kit
E-EL-MK0593	Monkey CD42 (Cluster of Differentiation 42) ELISA Kit
E-EL-MK0594	Monkey MCT (Mast Cell Tryptase) ELISA Kit
E-EL-MK0596	Monkey RNF4 (Ring Finger Protein 4) ELISA Kit
E-EL-MK0597	Monkey CK-HMW (Cytokeratin High Molecular Weight) ELISA Kit
E-EL-MK0599	Monkey CTSK (Cathepsin K) ELISA Kit

Cat.NO	Product Name
E-EL-MK0600	Monkey ADRA1A (Adrenergic Receptor alpha-1A) ELISA Kit
E-EL-MK0601	Monkey SLPI (Secretory Leukocyte Protease Inhibitor) ELISA Kit
E-EL-MK0602	Monkey IL-20 (Interleukin 20) ELISA Kit
E-EL-MK0603	Monkey IL-13 (Interleukin 13) ELISA Kit
E-EL-MK0604	Monkey P2 (Myelin Protein 2) ELISA Kit
E-EL-MK0605	Monkey MUC13 (Mucin 13) ELISA Kit
E-EL-MK0606	Monkey VF (Visfatin) ELISA Kit
E-EL-MK0607	Monkey PRELP (Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit
E-EL-MK0608	Monkey DMD (Dystrophin) ELISA Kit
E-EL-MK0609	Monkey AMH (Anti-Mullerian Hormone) ELISA Kit
E-EL-MK0612	Monkey ITLN1 (Intelectin 1/Omentin) ELISA Kit
E-EL-MK0615	Monkey CTAGE (Cutaneous T-Cell Lymphocyte Associated Antigen) ELISA Kit
E-EL-MK0616	Monkey ATF-5 (Activating Transcription Factor 5) ELISA Kit
E-EL-MK0617	Monkey MGST1 (Microsomal Glutathione S Transferase 1) ELISA Kit
E-EL-MK0619	Monkey MMP-25 (Matrix Metalloproteinase 25) ELISA Kit
E-EL-MK0620	Monkey OCLN (Occludin) ELISA Kit
E-EL-MK0622	Monkey TNC (Tenascin C) ELISA Kit
E-EL-MK0623	Monkey EPP (Erythrocyte Protoporphyrin) ELISA Kit
E-EL-MK0624	Monkey IgA1 (Immunoglobulin A1) ELISA Kit
E-EL-MK0625	Monkey MSTN (Myostatin) ELISA Kit

Cat.NO	Product Name
E-EL-MK0627	Monkey EGFR (Epidermal Growth Factor Receptor) ELISA Kit
E-EL-MK0628	Monkey UFC (Urinary Free Cortisol) ELISA Kit
E-EL-MK0629	Monkey PDGFC (Platelet Derived Growth Factor C) ELISA Kit
E-EL-MK0630	Monkey GPR120 (G Protein Coupled Receptor 120) ELISA Kit
E-EL-MK0631	Monkey CFB (Complement Factor B) ELISA Kit
E-EL-MK0632	Monkey COL17 (Collagen Type XVII) ELISA Kit
E-EL-MK0633	Monkey CTSB (Cathepsin B) ELISA Kit
E-EL-MK0634	Monkey IGF2BP3 (Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit
E-EL-MK0635	Monkey NEB (Noradrenaline Bitartas) ELISA Kit
E-EL-MK0636	Monkey APLN (Apelin) ELISA Kit
E-EL-MK0637	Monkey BCDF (B-Cell Differentiation Factor) ELISA Kit
E-EL-MK0638	Monkey BAG3 (Bcl2 Associated Athanogene 3) ELISA Kit
E-EL-MK0640	Monkey PCDH1 (Protocadherin 1) ELISA Kit
E-EL-MK0641	Monkey GAD (Glutamic Acid Decarboxylase) ELISA Kit
E-EL-MK0642	Monkey CX31 (Connexin 31) ELISA Kit
E-EL-MK0643	Monkey CYP7A1 (Cytochrome P450 7A1) ELISA Kit
E-EL-MK0644	Monkey Sp17 (Sperm Protein 17) ELISA Kit
E-EL-MK0645	Monkey Gal-6S (Galactose-6-Sulphatase) ELISA Kit
E-EL-MK0646	Monkey CASP12 (Caspase 12) ELISA Kit
E-EL-MK0647	Monkey GATA3 (GATA Binding Protein 3) ELISA Kit

Cat.NO	Product Name
E-EL-MK0648	Monkey ANGPTL3 (Angiopoietin Like Protein 3) ELISA Kit
E-EL-MK0649	Monkey C2 (Complement Component 2) ELISA Kit
E-EL-MK0650	Monkey MYH9 (Myosin Heavy Chain 9, Non Muscle) ELISA Kit
E-EL-MK0651	Monkey SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit
E-EL-MK0652	Monkey SOX2 (Sex Determining Region Y Box Protein 2) ELISA Kit
E-EL-MK0653	Monkey CKLFSF4 (Chemokine Like Factor Superfamily 4) ELISA Kit
E-EL-MK0654	Monkey TFR2 (Transferrin Receptor 2) ELISA Kit
E-EL-MK0655	Monkey PKB (Protein Kinase B) ELISA Kit
E-EL-MK0656	Monkey PEDF (Pigment Epithelium Derived Factor) ELISA Kit
E-EL-MK0657	Monkey PYGM (Glycogen Phosphorylase, Muscle) ELISA Kit
E-EL-MK0658	Monkey UPP1 (Uridine Phosphorylase 1) ELISA Kit
E-EL-MK0659	Monkey ADH (Antidiuretic Hormone) ELISA Kit
E-EL-MK0660	Monkey PLA2G2A (Phospholipase A2, Group II A) ELISA Kit
E-EL-MK0661	Monkey BCR (B-Cell Receptor) ELISA Kit
E-EL-MK0662	Monkey CA4 (Carbonic Anhydrase IV) ELISA Kit
E-EL-MK0663	Monkey C3c (Complement C3 Convertase) ELISA Kit
E-EL-MK0664	Monkey AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit
E-EL-MK0665	Monkey UGT2B4 (Uridine Diphosphate Glucuronosyltransferase 2 Family Polypeptide B4) ELISA Kit
E-EL-MK0666	Monkey FGF6 (Fibroblast Growth Factor 6) ELISA Kit
E-EL-MK0667	Monkey PGR (Progesterone Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-MK0669	Monkey NPYR5 (Neuropeptide Y Receptor Y5) ELISA Kit
E-EL-MK0670	Monkey BCL2L1/BCL-X (Bcl-2 Like Protein 1) ELISA Kit
E-EL-MK0671	Monkey MUC2 (Mucin 2) ELISA Kit
E-EL-MK0672	Monkey ES (Endostatin) ELISA Kit
E-EL-MK0673	Monkey OCT (Ornithine Carbamoyl Transferase) ELISA Kit
E-EL-MK0674	Monkey 5-LO (Arachidonate 5-Lipoxygenase) ELISA Kit
E-EL-MK0675	Monkey CENPE (Centromere Protein E) ELISA Kit
E-EL-MK0676	Monkey KLK7 (Kallikrein 7) ELISA Kit
E-EL-MK0677	Monkey TNF- β (Tumor Necrosis Factor Beta) ELISA Kit
E-EL-MK0678	Monkey DNM2 (Dynamin 2) ELISA Kit
E-EL-MK0679	Monkey ANG (Angiogenin) ELISA Kit
E-EL-MK0680	Monkey PAF (Platelet Activating Factor) ELISA Kit
E-EL-MK0681	Monkey PUMA (P53 Upregulated Modulator of Apoptosis) ELISA Kit
E-EL-MK0682	Monkey Sp100 (Sp100 Nuclear Antigen) ELISA Kit
E-EL-MK0683	Monkey HSP-90 (Heat Shock Protein 90) ELISA Kit
E-EL-MK0685	Monkey ECFR (Eosinophil Chemotactic Factor Receptor) ELISA Kit
E-EL-MK0686	Monkey fT4 (Free Thyroxine) ELISA Kit
E-EL-MK0687	Monkey Bcl-10 (B-cell Leukemia/Lymphoma 10) ELISA Kit
E-EL-MK0688	Monkey RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
E-EL-MK0689	Monkey ADCY4 (Adenylate Cyclase 4) ELISA Kit

Cat.NO	Product Name
E-EL-MK0690	Monkey FcγR II/CD32 (Receptor II for the Fc Fragment of IgG) ELISA Kit
E-EL-MK0691	Monkey PADI4 (Peptidyl Arginine Deiminase Type IV) ELISA Kit
E-EL-MK0692	Monkey AhR (Aryl Hydrocarbon Receptor) ELISA Kit
E-EL-MK0693	Monkey CSF3/GCSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit
E-EL-MK0695	Monkey CCND2 (Cyclin-D2) ELISA Kit
E-EL-MK0696	Monkey LBP (Lipopolysaccharide Binding Protein) ELISA Kit
E-EL-MK0697	Monkey PTH (Parathyroid Hormone) ELISA Kit
E-EL-MK0698	Monkey PDXK (Pyridoxal Kinase) ELISA Kit
E-EL-MK0699	Monkey PLC (Phospholipase C) ELISA Kit
E-EL-MK0701	Monkey APSA (Anti-Phosphatidyl Serine Antibody) ELISA Kit
E-EL-MK0702	Monkey C3c (Complement Fragment 3c) ELISA Kit
E-EL-MK0703	Monkey PRG4 (Proteoglycan 4) ELISA Kit
E-EL-MK0704	Monkey C4 (Complement Component 4) ELISA Kit
E-EL-MK0705	Monkey BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit
E-EL-MK0706	Monkey GSTα4 (Glutathione S Transferase Alpha 4) ELISA Kit
E-EL-MK0707	Monkey MAPKAPK3 (Mitogen Activated Protein Kinase Activated Protein Kinase 3) ELISA Kit
E-EL-MK0708	Monkey GDF10 (Growth Differentiation Factor 10) ELISA Kit
E-EL-MK0709	Monkey GAL1 (Galectin 1) ELISA Kit
E-EL-MK0710	Monkey CT-1 (Cardiotrophin 1) ELISA Kit
E-EL-MK0711	Monkey LFA-1 (Lymphocyte Function Associated Antigen 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0712	Monkey ATF-1 (Activating Transcription Factor 1) ELISA Kit
E-EL-MK0713	Monkey SIGLEC3 (Sialic Acid Binding Ig Like Lectin 3) ELISA Kit
E-EL-MK0714	Monkey EPO (Erythropoietin) ELISA Kit
E-EL-MK0715	Monkey SELL (L-Selectin) ELISA Kit
E-EL-MK0717	Monkey SAP (Serum Amyloid P Component) ELISA Kit
E-EL-MK0719	Monkey SDF2L1 (Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
E-EL-MK0720	Monkey C4BP (C4 Binding Protein) ELISA Kit
E-EL-MK0721	Monkey CASP10 (Caspase 10) ELISA Kit
E-EL-MK0722	Monkey STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit
E-EL-MK0723	Monkey HB μ (Hemoglobin, Mu) ELISA Kit
E-EL-MK0725	Monkey CNTNAP1 (Contactin Associated Protein 1) ELISA Kit
E-EL-MK0726	Monkey Ntn1 (Netrin 1) ELISA Kit
E-EL-MK0727	Monkey HGF (Hepatocyte Growth Factor) ELISA Kit
E-EL-MK0728	Monkey GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit
E-EL-MK0729	Monkey PC/CPG (Choline Phosphoglyceride) ELISA Kit
E-EL-MK0730	Monkey MYO (Myosin) ELISA Kit
E-EL-MK0731	Monkey G6PD (Glucose 6 Phosphate Dehydrogenase) ELISA Kit
E-EL-MK0732	Monkey MDC/CCL22 (Macrophage Derived Chemokine) ELISA Kit
E-EL-MK0733	Monkey REG3 α (Regenerating Islet Derived Protein 3 Alpha) ELISA Kit
E-EL-MK0734	Monkey TGF β 1 (Transforming Growth Factor Beta Induced Protein) ELISA Kit

Cat.NO	Product Name
E-EL-MK0735	Monkey BFP (Brain Finger Protein) ELISA Kit
E-EL-MK0737	Monkey IL-9 (Interleukin 9) ELISA Kit
E-EL-MK0738	Monkey DT (Dopachrome Tautomerase) ELISA Kit
E-EL-MK0739	Monkey M-AChR M1 (Muscarinic Acetylcholine Receptor M1) ELISA Kit
E-EL-MK0740	Monkey IDU α (Alpha-L-Iduronidase) ELISA Kit
E-EL-MK0741	Monkey EMMPRIN/CD147 (Extracellular Matrix Metalloproteinase Inducer) ELISA Kit
E-EL-MK0742	Monkey β -EPR (Beta-Endorphin Receptor) ELISA Kit
E-EL-MK0743	Monkey GST μ 2/GSTm2 (Glutathione S Transferase Mu 2, Muscle) ELISA Kit
E-EL-MK0744	Monkey NRG-2 (Neuregulin 2) ELISA Kit
E-EL-MK0745	Monkey MAOA (Type A Monoamine Oxidase) ELISA Kit
E-EL-MK0747	Monkey MCRS1 (Microspherule Protein 1) ELISA Kit
E-EL-MK0748	Monkey S100A12 (S100 Calcium Binding Protein A12) ELISA Kit
E-EL-MK0749	Monkey ICAM-4 (Intercellular Adhesion Molecule-4) ELISA Kit
E-EL-MK0750	Monkey EG-VEGF (Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit
E-EL-MK0751	Monkey PSA/KLK3 (Prostate Specific Antigen) ELISA Kit
E-EL-MK0752	Monkey CX26 (Connexin 26) ELISA Kit
E-EL-MK0753	Monkey P-CK (Pan-Cytokeratin) ELISA Kit
E-EL-MK0754	Monkey C/EBP δ (CCAAT/Enhancer Binding Protein, Delta) ELISA Kit
E-EL-MK0755	Monkey MCSFR (Macrophage Colony Stimulating Factor Receptor) ELISA Kit
E-EL-MK0756	Monkey CAPA (Cytotoxin Associated Protein A) ELISA Kit

Cat.NO	Product Name
E-EL-MK0757	Monkey OBTF4 (Octamer Binding Transcription Factor 4) ELISA Kit
E-EL-MK0758	Monkey NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
E-EL-MK0759	Monkey PCDH15 (Protocadherin 15) ELISA Kit
E-EL-MK0760	Monkey SIGLEC1 (Sialic Acid Binding Ig Like Lectin 1) ELISA Kit
E-EL-MK0761	Monkey DAP6 (Death Associated Protein 6) ELISA Kit
E-EL-MK0762	Monkey VIM (Vimentin) ELISA Kit
E-EL-MK0763	Monkey PPIL1 (Peptidylprolyl Isomerase/Cyclophilin Like 1) ELISA Kit
E-EL-MK0764	Monkey CTNN α 1 (Catenin, Alpha 1) ELISA Kit
E-EL-MK0765	Monkey HbC (Hemoglobin C) ELISA Kit
E-EL-MK0766	Monkey CD99 (Cluster of Differentiation 99) ELISA Kit
E-EL-MK0767	Monkey EGR4 (Early Growth Response Protein 4) ELISA Kit
E-EL-MK0768	Monkey TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit
E-EL-MK0769	Monkey Ub (Ubiquitin) ELISA Kit
E-EL-MK0770	Monkey GST θ 1/GSTt1 (Glutathione S Transferase Theta 1) ELISA Kit
E-EL-MK0771	Monkey GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit
E-EL-MK0773	Monkey PDGFB (Platelet Derived Growth Factor Subunit B) ELISA Kit
E-EL-MK0774	Monkey GIP (Gastric Inhibitory Polypeptide) ELISA Kit
E-EL-MK0776	Monkey SPARC (Secreted Protein Acidic And Rich In Cysteine) ELISA Kit
E-EL-MK0777	Monkey PTGIS (Prostaglandin I Synthase) ELISA Kit
E-EL-MK0778	Monkey D β H (Dopamine- β -Hydroxylase) ELISA Kit

Cat.NO	Product Name
E-EL-MK0779	Monkey PC4 (Positive Cofactor 4) ELISA Kit
E-EL-MK0780	Monkey CK-18/KRT18 (Cytokeratin 18) ELISA Kit
E-EL-MK0781	Monkey BF (Factor B) ELISA Kit
E-EL-MK0782	Monkey MRC1 (Mannose Receptor C Type 1 Like Protein 1) ELISA Kit
E-EL-MK0783	Monkey CPN10 (Chaperonin 10) ELISA Kit
E-EL-MK0784	Monkey PTD (Pentosidine) ELISA Kit
E-EL-MK0785	Monkey Smad3/MADH3 (Mothers Against Decapentaplegic Homolog 3) ELISA Kit
E-EL-MK0787	Monkey GR β (Glucocorticoid Receptor Beta) ELISA Kit
E-EL-MK0788	Monkey CD163 (Cluster of Differentiation 163) ELISA Kit
E-EL-MK0789	Monkey CCKBR (Cholecystokinin B Receptor) ELISA Kit
E-EL-MK0790	Monkey KAF (Keratinocyte Autocrine Factor) ELISA Kit
E-EL-MK0791	Monkey GDH/GLDH (Glutamate dehydrogenase) ELISA Kit
E-EL-MK0792	Monkey PCPE1 (Procollagen C Proteinase Enhancer 1) ELISA Kit
E-EL-MK0793	Monkey LN (Laminin) ELISA Kit
E-EL-MK0794	Monkey IGFBP-4 (Insulin Like Growth Factor Binding Protein 4) ELISA Kit
E-EL-MK0795	Monkey CENPF (Centromere Protein F) ELISA Kit
E-EL-MK0796	Monkey AP 17 (Apelin 17) ELISA Kit
E-EL-MK0797	Monkey F2 (Coagulation Factor II) ELISA Kit
E-EL-MK0798	Monkey CNR-1 (Cadherin Related Neuronal Receptor1) ELISA Kit
E-EL-MK0799	Monkey CT (Calcitonin) ELISA Kit

Cat.NO	Product Name
E-EL-MK0801	Monkey GPIIb/IIIa (Platelet Membrane Glycoprotein IIb/IIIa) ELISA Kit
E-EL-MK0802	Monkey PFKP (Phosphofructokinase, Platelet) ELISA Kit
E-EL-MK0803	Monkey BSP (Bone Sialoprotein) ELISA Kit
E-EL-MK0804	Monkey α 1BG (Alpha-1-B-Glycoprotein) ELISA Kit
E-EL-MK0805	Monkey Tgtp (T-cell Specific Gtpase) ELISA Kit
E-EL-MK0806	Monkey DAO (Diamine Oxidase) ELISA Kit
E-EL-MK0807	Monkey IGFBP-7 (Insulin Like Growth Factor Binding Protein 7) ELISA Kit
E-EL-MK0809	Monkey CDNF (Cerebral Dopamine Neurotrophic Factor) ELISA Kit
E-EL-MK0810	Monkey PTGES (Prostaglandin E Synthase, Microsomal)ELISA Kit
E-EL-MK0811	Monkey Gp49a (Glycoprotein 49A) ELISA Kit
E-EL-MK0812	Monkey BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit
E-EL-MK0813	Monkey GZMK (Granzyme K) ELISA Kit
E-EL-MK0814	Monkey TMPRSS2 (Transmembrane Protease, Serine 2) ELISA Kit
E-EL-MK0815	Monkey MYD88 (Myeloid Differentiation Factor 88) ELISA Kit
E-EL-MK0816	Monkey ALCAM (Activated Leukocyte Cell Adhesion Molecule) ELISA Kit
E-EL-MK0817	Monkey BMPR2 (Bone Morphogenetic Protein Receptor II) ELISA Kit
E-EL-MK0818	Monkey TNNC2 (Troponin C Type 2) ELISA Kit
E-EL-MK0819	Monkey TSH (Thyroid Stimulating Hormone) ELISA Kit
E-EL-MK0820	Monkey AChRab (Acetylcholine Receptor Antibody) ELISA Kit
E-EL-MK0821	Monkey FGF9 (Fibroblast Growth Factor 9) ELISA Kit

Cat.NO	Product Name
E-EL-MK0822	Monkey TPS (Tissue Polypeptide Specific Antigen) ELISA Kit
E-EL-MK0823	Monkey IgG1 (Immunoglobulin G1) ELISA Kit
E-EL-MK0824	Monkey PDH α (Pyruvate Dehydrogenase Alpha) ELISA Kit
E-EL-MK0825	Monkey Anti-AChR (Anti-Acetylcholine Receptor Antibody) ELISA Kit
E-EL-MK0826	Monkey ARHGD1 α (Rho GDP Dissociation Inhibitor Alpha) ELISA Kit
E-EL-MK0827	Monkey CDK6 (Cyclin Dependent Kinase 6) ELISA Kit
E-EL-MK0828	Monkey CLL1 (Collectin Liver 1) ELISA Kit
E-EL-MK0829	Monkey bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit
E-EL-MK0830	Monkey α -SMA (α -Smooth muscle actin) ELISA Kit
E-EL-MK0831	Monkey GR (Glutathione Reductase) ELISA Kit
E-EL-MK0832	Monkey LTC4 (Leukotriene C4) ELISA Kit
E-EL-MK0833	Monkey MRC2 (Mannose Receptor C Type 2) ELISA Kit
E-EL-MK0834	Monkey RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
E-EL-MK0835	Monkey PKI α (Protein Kinase Inhibitor Alpha) ELISA Kit
E-EL-MK0836	Monkey SAA3 (Serum Amyloid A3) ELISA Kit
E-EL-MK0837	Monkey DAF (Decay Accelerating Factor) ELISA Kit
E-EL-MK0838	Monkey PTHrP (Parathyroid Hormone Related Protein) ELISA Kit
E-EL-MK0839	Monkey MTX (Methotrexate) ELISA Kit
E-EL-MK0840	Monkey GL β (Galactosidase Beta) ELISA Kit
E-EL-MK0841	Monkey FRS2 (Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit

Cat.NO	Product Name
E-EL-MK0842	Monkey TSP-4 (Thrombospondin-4) ELISA Kit
E-EL-MK0843	Monkey CRHBP (Corticotropin Releasing Hormone Binding Protein) ELISA Kit
E-EL-MK0844	Monkey IgA2 (Immunoglobulin A2) ELISA Kit
E-EL-MK0845	Monkey MK (Midkine) ELISA Kit
E-EL-MK0846	Monkey Tyk-2 (Tyrosine Kinase 2) ELISA Kit
E-EL-MK0847	Monkey MN (Metanephrine) ELISA Kit
E-EL-MK0848	Monkey GlyCAM1 (Glycosylation Dependent Cell Adhesion Molecule 1) ELISA Kit
E-EL-MK0849	Monkey PP (Pepsin) ELISA Kit
E-EL-MK0850	Monkey HSF2 (Heat Shock Transcription Factor 2) ELISA Kit
E-EL-MK0851	Monkey SPJ (Surfactant Associated Protein J) ELISA Kit
E-EL-MK0852	Monkey PLA2G4D (Phospholipase A2, Group IV D) ELISA Kit
E-EL-MK0853	Monkey MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
E-EL-MK0856	Monkey MUSK (Muscle Skeletal Receptor Tyrosine Kinase) ELISA Kit
E-EL-MK0858	Monkey CD30 (Cluster of Differentiation 30) ELISA Kit
E-EL-MK0859	Monkey TE (Telomerase) ELISA Kit
E-EL-MK0860	Monkey ATF7 (Activating Transcription Factor 7) ELISA Kit
E-EL-MK0861	Monkey TNNI1 (Troponin I Type 1, Slow Skeletal) ELISA Kit
E-EL-MK0862	Monkey GPER1 (G Protein Coupled Estrogen Receptor 1) ELISA Kit
E-EL-MK0863	Monkey PRB1 (Basic Salivary Proline Rich Protein 1) ELISA Kit
E-EL-MK0864	Monkey Bid (BH3 Interacting Domain Death Agonist) ELISA Kit

Cat.NO	Product Name
E-EL-MK0867	Monkey PI3K (Phosphoinositide-3-Kinase) ELISA Kit
E-EL-MK0868	Monkey Smad4/MADH4 (Mothers Against Decapentaplegic Homolog 4)ELISA Kit
E-EL-MK0869	Monkey AMBP (Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit
E-EL-MK0870	Monkey LEPR (Leptin Receptor) ELISA Kit
E-EL-MK0871	Monkey Ly6a (Lymphocyte Antigen 6 Complex Locus A) ELISA Kit
E-EL-MK0874	Monkey VWF (Von Willebrand Factor) ELISA Kit
E-EL-MK0875	Monkey EP2 (Prostaglandin E Receptor 2) ELISA Kit
E-EL-MK0878	Monkey SAP25 (Synaptosome Associated Protein 25) ELISA Kit
E-EL-MK0879	Monkey MMP-8 (Matrix Metalloproteinase 8) ELISA Kit
E-EL-MK0880	Monkey BK (Bradykinin) ELISA Kit
E-EL-MK0881	Monkey ACHAP (Acetylcholinesterase Associated Protein) ELISA Kit
E-EL-MK0882	Monkey PLP (Parathyroid Hormone-Like Protein) ELISA Kit
E-EL-MK0883	Monkey PRSS1 (Protease, Serine, 1) ELISA Kit
E-EL-MK0884	Monkey SDH (Sorbitol Dehydrogenase) ELISA Kit
E-EL-MK0885	Monkey PGF2 α (Prostaglandin F2 α) ELISA Kit
E-EL-MK0886	Monkey RND1 (Rho Family GTPase 1) ELISA Kit
E-EL-MK0887	Monkey TAP (Transporter Associated with Antigen Processing) ELISA Kit
E-EL-MK0888	Monkey NA/NE (Noradrenaline/Norepinephrine) ELISA Kit
E-EL-MK0889	Monkey TPO (Thrombopoietin) ELISA Kit
E-EL-MK0890	Monkey BACE1 (Beta-site APP Cleaving Enzyme 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0891	Monkey Cyt-C (Cytochrome C) ELISA Kit
E-EL-MK0892	Monkey LRP (Lung Resistance-Related Protein) ELISA Kit
E-EL-MK0893	Monkey ASB (Arylsulfatase B) ELISA Kit
E-EL-MK0895	Monkey ANEA (Anti-Neutrophil Elastase Antibody) ELISA Kit
E-EL-MK0896	Monkey UBR4 (Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit
E-EL-MK0897	Monkey OLFM4 (Olfactomedin 4) ELISA Kit
E-EL-MK0898	Monkey PAP (Plasmin-Antiplasmin Complex) ELISA Kit
E-EL-MK0899	Monkey CTAPIII (Connective Tissue Activating Peptide III) ELISA Kit
E-EL-MK0900	Monkey TFF2 (Trefoil Factor 2) ELISA Kit
E-EL-MK0901	Monkey HCCS (Holocytochrome C Synthase) ELISA Kit
E-EL-MK0902	Monkey FGFR4 (Fibroblast Growth Factor Receptor 4) ELISA Kit
E-EL-MK0903	Monkey MGP (Matrix Gla Protein) ELISA Kit
E-EL-MK0904	Monkey C1INH (Complement 1 Inhibitor) ELISA Kit
E-EL-MK0905	Monkey GLP-1 (Glucagon Like Peptide 1) ELISA Kit
E-EL-MK0906	Monkey KIM-1 (Kidney Injury Molecule 1) ELISA Kit
E-EL-MK0907	Monkey ACA-IgM (Anti-Cardiolipin Antibody IgM) ELISA Kit
E-EL-MK0908	Monkey F13 (Coagulation Factor X III) ELISA Kit
E-EL-MK0909	Monkey NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) ELISA Kit
E-EL-MK0910	Monkey PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit
E-EL-MK0911	Monkey UPP2 (Uridine Phosphorylase 2) ELISA Kit

Cat.NO	Product Name
E-EL-MK0912	Monkey GPR37 (G Protein Coupled Receptor 37) ELISA Kit
E-EL-MK0913	Monkey UBC (Urinary Bladder Cancer Antigen) ELISA Kit
E-EL-MK0914	Monkey ATXN1 (Ataxin 1) ELISA Kit
E-EL-MK0915	Monkey COL13 (Collagen Type XIII) ELISA Kit
E-EL-MK0916	Monkey IL-4 (Interleukin 4) ELISA Kit
E-EL-MK0917	Monkey CD109 (Cluster of Differentiation 109) ELISA Kit
E-EL-MK0918	Monkey GS (Gelsolin) ELISA Kit
E-EL-MK0919	Monkey CK-20/KRT20 (Cytokeratin 20) ELISA Kit
E-EL-MK0920	Monkey cPLA2 (Cytosolic Phospholipase A2) ELISA Kit
E-EL-MK0921	Monkey CLEC2C (C-Type Lectin Domain Family 2 Member C) ELISA Kit
E-EL-MK0922	Monkey PGM (Phosphoglucomutase) ELISA Kit
E-EL-MK0923	Monkey f-Hb (Free Haemoglobin) ELISA Kit
E-EL-MK0924	Monkey SLC/CCL21 (Secondary Lymphoid Tissue Chemokine) ELISA Kit
E-EL-MK0925	Monkey C3b (Complement Fragment 3b) ELISA Kit
E-EL-MK0926	Monkey GC (Glycocalicin) ELISA Kit
E-EL-MK0927	Monkey TYR (Tyrosinase) ELISA Kit
E-EL-MK0928	Monkey DSP (Desmoplakin) ELISA Kit
E-EL-MK0929	Monkey CYP24A1 (Cytochrome P450 24A1) ELISA Kit
E-EL-MK0930	Monkey C/EBP β (CCAAT/Enhancer Binding Protein Beta) ELISA Kit
E-EL-MK0931	Monkey CASP4 (Caspase 4) ELISA Kit

Cat.NO	Product Name
E-EL-MK0932	Monkey IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-MK0933	Monkey ETF β (Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit
E-EL-MK0934	Monkey TDGF1 (Teratocarcinoma Derived Growth Factor 1) ELISA Kit
E-EL-MK0936	Monkey LMP7/PSMB9 (Low Molecular-weight Protein/proteasome Beta-type Subunit) ELISA Kit
E-EL-MK0937	Monkey CALLA (Common Acute Lymphocytic Leukaemia Antigen) ELISA Kit
E-EL-MK0938	Monkey CD59 (Protectin) ELISA Kit
E-EL-MK0939	Monkey AChE (Acetylcholinesterase) ELISA Kit
E-EL-MK0941	Monkey MPIF-1/CCL23 (Myeloid Progenitor Inhibitory Factor 1) ELISA Kit
E-EL-MK0942	Monkey MMP-24 (Matrix Metalloproteinase 24) ELISA Kit
E-EL-MK0943	Monkey FGL2 (Fibrinogen Like Protein 2) ELISA Kit
E-EL-MK0944	Monkey GAL2 (Galectin 2) ELISA Kit
E-EL-MK0945	Monkey GATA5 (GATA Binding Protein 5) ELISA Kit
E-EL-MK0946	Monkey TGM3 (Transglutaminase 3, Epidermal) ELISA Kit
E-EL-MK0947	Monkey FG β (Fibrinogen Beta) ELISA Kit
E-EL-MK0950	Monkey TBG (Thyroxine Binding Globulin) ELISA Kit
E-EL-MK0952	Monkey PARP4 (Poly ADP Ribose Polymerase 4) ELISA Kit
E-EL-MK0954	Monkey KLK8 (Kallikrein 8) ELISA Kit
E-EL-MK0955	Monkey Tn-T (Troponin T) ELISA Kit
E-EL-MK0956	Monkey CTR (Calcitonin Receptor) ELISA Kit
E-EL-MK0957	Monkey CKLFSF8 (Chemokine Like Factor Superfamily 8) ELISA Kit

Cat.NO	Product Name
E-EL-MK0958	Monkey MT2 (Metallothionein 2) ELISA Kit
E-EL-MK0959	Monkey CAD (Caspase Activated Deoxyribonuclease) ELISA Kit
E-EL-MK0960	Monkey ANG II R2-Ab (Angiotensin II Receptor 2 Antibody) ELISA Kit
E-EL-MK0961	Monkey EREG (Epiregulin) ELISA Kit
E-EL-MK0962	Monkey MAT (Matrilysin) ELISA Kit
E-EL-MK0964	Monkey TFPI (Tissue Factor Pathway Inhibitor) ELISA Kit
E-EL-MK0965	Monkey PROC (Protein C) ELISA Kit
E-EL-MK0967	Monkey BLNK (B-Cell Linker Protein) ELISA Kit
E-EL-MK0968	Monkey TPP1 (Tripeptidyl Peptidase I) ELISA Kit
E-EL-MK0969	Monkey CgA (Chromogranin A/Pancreastatin) ELISA Kit
E-EL-MK0970	Monkey IL-18 (Interleukin 18) ELISA Kit
E-EL-MK0971	Monkey GPI (Glucose 6 Phosphate Isomerase) ELISA Kit
E-EL-MK0972	Monkey GzmB (Granzyme B) ELISA Kit
E-EL-MK0973	Monkey MIF (Macrophage Migration Inhibitory Factor) ELISA Kit
E-EL-MK0974	Monkey RAG-2 (Recombination Activating Gene 2) ELISA Kit
E-EL-MK0975	Monkey Thyronine ELISA Kit
E-EL-MK0977	Monkey PLAA (Phospholipase A2 Activating Protein) ELISA Kit
E-EL-MK0979	Monkey OTF2A (Octamer Transcription Factor 2A) ELISA Kit
E-EL-MK0980	Monkey TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
E-EL-MK0981	Monkey EXT1 (Exostoses 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK0982	Monkey ANKRD1 (Ankyrin Repeat Domain 1) ELISA Kit
E-EL-MK0983	Monkey ApoM (Apolipoprotein M) ELISA Kit
E-EL-MK0984	Monkey HNRPU (Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit
E-EL-MK0985	Monkey NGF (Nerve Growth Factor) ELISA Kit
E-EL-MK0986	Monkey CFD (Complement Factor D) ELISA Kit
E-EL-MK0987	Monkey MAGE (Melanoma-associated Antigen) ELISA Kit
E-EL-MK0988	Monkey MCP-2/CCL8 (Monocyte Chemotactic Protein 2) ELISA Kit
E-EL-MK0989	Monkey PIICP (Procollagen II C-Terminal ProPeptide) ELISA Kit
E-EL-MK0990	Monkey DHEA (Dehydroepiandrosterone) ELISA Kit
E-EL-MK0991	Monkey PDAP1 (Pdgfa Associated Protein 1) ELISA Kit
E-EL-MK0992	Monkey β 2-GP1 Ab IgG (β 2-Glycoprotein 1 Antibody IgG) ELISA Kit
E-EL-MK0993	Monkey GLUT4 (Glucose Transporter 4) ELISA Kit
E-EL-MK0994	Monkey tPA (Plasminogen Activator, Tissue) ELISA Kit
E-EL-MK0995	Monkey N-MID-OT (N-MID Osteocalcin) ELISA Kit
E-EL-MK0998	Monkey ADAM8 (A Disintegrin And Metalloprotease 8) ELISA Kit
E-EL-MK0999	Monkey SMAF (Specific Macrophage Arming Factor) ELISA Kit
E-EL-MK1000	Monkey SEPP1 (Selenoprotein P) ELISA Kit
E-EL-MK1001	Monkey PKD1 (Protein Kinase D1) ELISA Kit
E-EL-MK1002	Monkey COL6 α 1 (Collagen Type VI Alpha 1) ELISA Kit
E-EL-MK1003	Monkey CG α (Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit

Cat.NO	Product Name
E-EL-MK1004	Monkey TNPO1 (Transportin 1) ELISA Kit
E-EL-MK1005	Monkey TDT (Terminal Deoxynucleotidyl Transferase) ELISA Kit
E-EL-MK1006	Monkey PC (Plasma Anticoagulant Protein C) ELISA Kit
E-EL-MK1007	Monkey PGC (Pepsinogen C) ELISA Kit
E-EL-MK1008	Monkey UBE3A (Ubiquitin Protein Ligase E3A) ELISA Kit
E-EL-MK1009	Monkey ADAM10 (A Disintegrin And Metalloprotease 10) ELISA Kit
E-EL-MK1010	Monkey TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
E-EL-MK1011	Monkey C8g (Complement Component 8g) ELISA Kit
E-EL-MK1012	Monkey CLC (Charcot Leyden Crystal Protein) ELISA Kit
E-EL-MK1013	Monkey NUP153 (Nucleoporin 153kDa) ELISA Kit
E-EL-MK1014	Monkey GTF2A1 (General Transcription Factor IIA, Polypeptide 1) ELISA Kit
E-EL-MK1017	Monkey D3R/DRD3 (Dopamine Receptor D3) ELISA Kit
E-EL-MK1018	Monkey HRG (Histidine Rich Glycoprotein) ELISA Kit
E-EL-MK1019	Monkey MAP τ (Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
E-EL-MK1020	Monkey CR (Calretinin) ELISA Kit
E-EL-MK1021	Monkey PDGFD (Platelet Derived Growth Factor D) ELISA Kit
E-EL-MK1022	Monkey HMGCS1 (3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit
E-EL-MK1023	Monkey ACTN3 (Actinin Alpha 3) ELISA Kit
E-EL-MK1024	Monkey GPI (Glycophosphatidylinositol) ELISA Kit
E-EL-MK1025	Monkey HSP-60/HSPD1 (Heat Shock Protein 60) ELISA Kit

Cat.NO	Product Name
E-EL-MK1026	Monkey TBP2 (Thioredoxin Binding Protein 2) ELISA Kit
E-EL-MK1027	Monkey CP (Ceruloplasmin) ELISA Kit
E-EL-MK1028	Monkey SREBP (Sterol Regulatory Element Binding Proteins) ELISA Kit
E-EL-MK1029	Monkey AKA (Anti-Keratin Antibody) ELISA Kit
E-EL-MK1030	Monkey SRF/IF (Skin Reactive Factor/Inflammatory Factor) ELISA Kit
E-EL-MK1031	Monkey DLD (Dihydrolipoyl Dehydrogenase) ELISA Kit
E-EL-MK1032	Monkey OxLDL (Oxidized Low Density Lipoprotein) ELISA Kit
E-EL-MK1034	Monkey APT (Abnormal Prothrombin) ELISA Kit
E-EL-MK1035	Monkey CEP110 (Centrosomal Protein 110kDa) ELISA Kit
E-EL-MK1036	Monkey Hpc (Hepcidin) ELISA Kit
E-EL-MK1037	Monkey α HSP (Alpha-Hemoglobin Stabilizing Protein) ELISA Kit
E-EL-MK1038	Monkey IDE (Insulin Degrading Enzyme) ELISA Kit
E-EL-MK1039	Monkey TAF (TATA Box Binding Protein/TBP-Associated Factors) ELISA Kit
E-EL-MK1041	Monkey GCLM (Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit
E-EL-MK1042	Monkey HB β (Hemoglobin Beta) ELISA Kit
E-EL-MK1043	Monkey AGA (Anti-Glycolipid Antibody) ELISA Kit
E-EL-MK1044	Monkey CHRM2 (Cholinergic Receptor, Muscarinic 2) ELISA Kit
E-EL-MK1045	Monkey Pro-MMP-10 (Pro-Matrix Metalloproteinase 10) ELISA Kit
E-EL-MK1046	Monkey β -actin (Actin Beta) ELISA Kit
E-EL-MK1048	Monkey GDF3 (Growth Differentiation Factor 3) ELISA Kit

Cat.NO	Product Name
E-EL-MK1050	Monkey PG (Proteoglycan) ELISA Kit
E-EL-MK1051	Monkey Ponticulin ELISA Kit
E-EL-MK1052	Monkey OPN (Osteopontin) ELISA Kit
E-EL-MK1053	Monkey CTX-III (Cross Linked C-telopeptide of Type III Collagen) ELISA Kit
E-EL-MK1054	Monkey CSF2R α (Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit
E-EL-MK1057	Monkey CACYBP (Calcyclin Binding Protein) ELISA Kit
E-EL-MK1058	Monkey BMP-1 (Bone Morphogenetic Protein 1) ELISA Kit
E-EL-MK1059	Monkey TXLNb (Taxilin Beta) ELISA Kit
E-EL-MK1060	Monkey TACR3 (Tachykinin Receptor 3) ELISA Kit
E-EL-MK1061	Monkey MTHFR (Methylenetetrahydrofolate Reductase) ELISA Kit
E-EL-MK1062	Monkey Fc α R/CD89 (Receptor I for the Fc Fragment of IgA) ELISA Kit
E-EL-MK1063	Monkey DCLK1 (Doublecortin Like Kinase 1) ELISA Kit
E-EL-MK1064	Monkey β MSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-MK1065	Monkey PKC γ (Protein Kinase C Gamma) ELISA Kit
E-EL-MK1066	Monkey EGF (Epidermal Growth Factor) ELISA Kit
E-EL-MK1067	Monkey α ZGP1 (Alpha-2-Glycoprotein 1, Zine Binding) ELISA Kit
E-EL-MK1068	Monkey RUNX2 (Runt-Related Transcription Factor 2) ELISA Kit
E-EL-MK1069	Monkey CKLFSF7 (Chemokine Like Factor Superfamily 7) ELISA Kit
E-EL-MK1070	Monkey ApoH (Apolipoprotein H) ELISA Kit
E-EL-MK1071	Monkey SCT (Secretin) ELISA Kit

Cat.NO	Product Name
E-EL-MK1072	Monkey VEGF-B (Vascular Endothelial Cell Growth Factor B) ELISA Kit
E-EL-MK1073	Monkey GLUT3 (Glucose Transporter 3) ELISA Kit
E-EL-MK1074	Monkey KLK6 (Kallikrein 6) ELISA Kit
E-EL-MK1075	Monkey HLE (Leukocyte Exostase) ELISA Kit
E-EL-MK1076	Monkey pERK (Phospho Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-MK1077	Monkey PCI (Protein C Inhibitor) ELISA Kit
E-EL-MK1078	Monkey EM2 (Endomorphin 2) ELISA Kit
E-EL-MK1080	Monkey SNRPC (Small Nuclear Ribonucleoprotein Polypeptide C) ELISA Kit
E-EL-MK1081	Monkey ME20M (Melanoma Associated ME20) ELISA Kit
E-EL-MK1082	Monkey CKMB (Creatine Kinase MB Isoenzyme) ELISA Kit
E-EL-MK1083	Monkey VLDLR (Very Low Density Lipoprotein Receptor) ELISA Kit
E-EL-MK1084	Monkey CCR1 (CC-Chemokine Receptor 1) ELISA Kit
E-EL-MK1085	Monkey VEGF (Vascular Endothelial Cell Growth Factor) ELISA Kit
E-EL-MK1086	Monkey CDKN3 (Cyclin Dependent Kinase Inhibitor 3) ELISA Kit
E-EL-MK1087	Monkey S100B (S100 Calcium Binding Protein B) ELISA Kit
E-EL-MK1088	Monkey β 2-GP1 Ab IgA (β 2-Glycoprotein 1 antibody IgA) ELISA Kit
E-EL-MK1089	Monkey EAR1 (Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit
E-EL-MK1091	Monkey SNCy (Synuclein Gamma) ELISA Kit
E-EL-MK1092	Monkey MUC17 (Mucin 17) ELISA Kit
E-EL-MK1093	Monkey RAG-1 (Recombination Activating Gene 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK1094	Monkey CCND1 (Cyclin-D1) ELISA Kit
E-EL-MK1095	Monkey ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit
E-EL-MK1096	Monkey RPS6K β 1 (Ribosomal Protein S6 Kinase Beta 1) ELISA Kit
E-EL-MK1097	Monkey GR α (Glucocorticoid Receptor Alpha) ELISA Kit
E-EL-MK1098	Monkey FBN2 (Fibrillin 2) ELISA Kit
E-EL-MK1099	Monkey HBF (Fetal Hemoglobin) ELISA Kit
E-EL-MK1100	Monkey CA724 (Tumor Marker CA724) ELISA Kit
E-EL-MK1101	Monkey SIGLEC2 (Sialic Acid Binding Ig Like Lectin 2) ELISA Kit
E-EL-MK1103	Monkey UBE1L (Ubiquitin Activating Enzyme E1 Like Protein) ELISA Kit
E-EL-MK1104	Monkey MFG-E8 (Milk Fat Globule EGF Factor 8) ELISA Kit
E-EL-MK1105	Monkey NPY (Neuropeptide Y) ELISA Kit
E-EL-MK1106	Monkey SM/SPH (Sphingomyelin) ELISA Kit
E-EL-MK1108	Monkey MT (Melatonin) ELISA Kit
E-EL-MK1109	Monkey RETN (Resistin) ELISA Kit
E-EL-MK1110	Monkey MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
E-EL-MK1111	Monkey P III CP (Procollagen III C-Terminal ProPeptide) ELISA Kit
E-EL-MK1112	Monkey PK (Pyruvate Kinase, Liver And RBC) ELISA Kit
E-EL-MK1113	Monkey GSTs (Glutathione S-Transferases) ELISA Kit
E-EL-MK1114	Monkey NUP205 (Nucleoporin 205kDa) ELISA Kit
E-EL-MK1115	Monkey fFDFT1 (Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK1116	Monkey IL-24 (Interleukin 24) ELISA Kit
E-EL-MK1117	Monkey DEFα5 (Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
E-EL-MK1118	Monkey C7 (Complement Component 7) ELISA Kit
E-EL-MK1119	Monkey MASP (Mannan Binding Lectin Associated Serine Protease) ELISA Kit
E-EL-MK1121	Monkey PDIA6 (Protein Disulfide Isomerase A6) ELISA Kit
E-EL-MK1124	Monkey PIR (Pirin) ELISA Kit
E-EL-MK1125	Monkey BLVRB (Biliverdin Reductase B) ELISA Kit
E-EL-MK1126	Monkey ADAMTS4 (ADAM with Thrombospondin Type 1 Motif4) ELISA Kit
E-EL-MK1127	Monkey TPPP (Tubulin Polymerization Promoting Protein) ELISA Kit
E-EL-MK1128	Monkey TREM-1 (Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
E-EL-MK1129	Monkey CgB (Chromogranin B) ELISA Kit
E-EL-MK1130	Monkey AZU1 (Azurocidin 1) ELISA Kit
E-EL-MK1131	Monkey ATA (Angiotensinase A) ELISA Kit
E-EL-MK1132	Monkey AR (Amphiregulin) ELISA Kit
E-EL-MK1133	Monkey HbH (Hemoglobin H) ELISA Kit
E-EL-MK1134	Monkey FG (Fibrinogen) ELISA Kit
E-EL-MK1135	Monkey GP5 (Glycoprotein V, Platelet) ELISA Kit
E-EL-MK1136	Monkey SPA (Pulmonary Surfactant Associated Protein A) ELISA Kit
E-EL-MK1137	Monkey CDK4 (Cyclin Dependent Kinase 4) ELISA Kit
E-EL-MK1138	Monkey CKMT1B (Creatine Kinase, Mitochondrial 1B) ELISA Kit

Cat.NO	Product Name
E-EL-MK1139	Monkey PDIA5 (Protein Disulfide Isomerase A5) ELISA Kit
E-EL-MK1140	Monkey TPO (Thyroid Peroxidase) ELISA Kit
E-EL-MK1141	Monkey PSG1/SP1 (Pregnancy Specific Beta-1-Glycoprotein 1) ELISA Kit
E-EL-MK1142	Monkey MBP/MBL (Mannose Binding Protein/Mannose Binding Lectin) ELISA Kit
E-EL-MK1143	Monkey Anti-CIIAb (Anti-Type II Collagen Antibody) ELISA Kit
E-EL-MK1144	Monkey SDC1 (Syndecan 1) ELISA Kit
E-EL-MK1145	Monkey OSM (Oncostatin M) ELISA Kit
E-EL-MK1146	Monkey NUMA1 (Nuclear Mitotic Apparatus Protein 1) ELISA Kit
E-EL-MK1147	Monkey SHBG (Sex Hormone Binding Globulin) ELISA Kit
E-EL-MK1148	Monkey NPYR1 (Neuropeptide Y Receptor Y1) ELISA Kit
E-EL-MK1150	Monkey MIgD (Membrane IgD) ELISA Kit
E-EL-MK1151	Monkey BALP (Bone Alkaline Phosphatase) ELISA Kit
E-EL-MK1152	Monkey GAD1 (Glutamate Decarboxylase 1) ELISA Kit
E-EL-MK1153	Monkey NUP188 (Nucleoporin 188kDa) ELISA Kit
E-EL-MK1154	Monkey FTH (Ferritin, Heavy Polypeptide) ELISA Kit
E-EL-MK1155	Monkey GALR1 (Galanin Receptor 1) ELISA Kit
E-EL-MK1156	Monkey PCDHβ2 (Protocadherin Beta 2) ELISA Kit
E-EL-MK1157	Monkey OIT3 (Oncoprotein Induced Transcript 3) ELISA Kit
E-EL-MK1158	Monkey SPC (Pulmonary Surfactant Associated Protein C) ELISA Kit
E-EL-MK1159	Monkey AGC (Aggrecan) ELISA Kit

Cat.NO	Product Name
E-EL-MK1160	Monkey HSP-27/HSPB1 (Heat Shock Protein 27) ELISA Kit
E-EL-MK1161	Monkey LAP (Leucine Aminopeptidase) ELISA Kit
E-EL-MK1162	Monkey TNNT1 (Troponin T Type 1, Slow Skeletal) ELISA Kit
E-EL-MK1163	Monkey NKA (Neurokinins A) ELISA Kit
E-EL-MK1164	Monkey HSP-40 (Heat Shock Protein 40) ELISA Kit
E-EL-MK1165	Monkey LSP (Liver Specific Lipoprotein) ELISA Kit
E-EL-MK1166	Monkey MT1E (Metallothionein 1E) ELISA Kit
E-EL-MK1167	Monkey LFA-3 (Lymphocyte Function Associated Antigen 3) ELISA Kit
E-EL-MK1168	Monkey DPT (Dermatopontin) ELISA Kit
E-EL-MK1169	Monkey C-P (C-Peptide) ELISA Kit
E-EL-MK1170	Monkey FcεR II (Receptor II for the Fc Fragment of IgE) ELISA Kit
E-EL-MK1173	Monkey SEMA3B (Semaphorin 3B) ELISA Kit
E-EL-MK1174	Monkey αFP (Alpha-Fetoprotein) ELISA Kit
E-EL-MK1176	Monkey NUP85 (Nucleoporin 85kDa) ELISA Kit
E-EL-MK1177	Monkey NT (Neurotensin) ELISA Kit
E-EL-MK1178	Monkey ADCY3 (Adenylate Cyclase 3) ELISA Kit
E-EL-MK1179	Monkey TFR/CD71 (Transferrin Receptor) ELISA Kit
E-EL-MK1180	Monkey PS2 (Presenilin 2) ELISA Kit
E-EL-MK1181	Monkey MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit
E-EL-MK1182	Monkey Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit

Cat.NO	Product Name
E-EL-MK1183	Monkey EGR2 (Early Growth Response Protein 2) ELISA Kit
E-EL-MK1184	Monkey TNKS2 (Tankyrase 2) ELISA Kit
E-EL-MK1185	Monkey FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit
E-EL-MK1186	Monkey PTPRO (Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
E-EL-MK1187	Monkey GP (Anti-Glycoprotein Antibody) ELISA Kit
E-EL-MK1188	Monkey 11-DH-TXB2 (11-Dehydro-thromboxane B2) ELISA Kit
E-EL-MK1189	Monkey ELA4 (Elastase 4) ELISA Kit
E-EL-MK1190	Monkey IL-1 α (Interleukin 1 Alpha) ELISA Kit
E-EL-MK1191	Monkey SDF-2 (Stromal Cell Derived Factor 2) ELISA Kit
E-EL-MK1192	Monkey GC (Glucagon) ELISA Kit
E-EL-MK1193	Monkey BMG/ β 2-MG (Beta-2-Microglobulin) ELISA Kit
E-EL-MK1194	Monkey GDF6 (Growth Differentiation Factor 6) ELISA Kit
E-EL-MK1195	Monkey SAG (S Antigen) ELISA Kit
E-EL-MK1196	Monkey SDC4 (Syndecan 4) ELISA Kit
E-EL-MK1197	Monkey SYK (Spleen Tyrosine Kinase) ELISA Kit
E-EL-MK1198	Monkey GAL12 (Galectin 12) ELISA Kit
E-EL-MK1199	Monkey PCI (Procollagen I) ELISA Kit
E-EL-MK1201	Monkey ACA-IgG (Anti-Cardiolipin Antibody IgG) ELISA Kit
E-EL-MK1202	Monkey α 2-PI (α 2-Plasmin Inhibitor) ELISA Kit
E-EL-MK1204	Monkey EPHA1 (Ephrin Type A Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK1205	Monkey ADP/Acrp30 (Adiponectin) ELISA Kit
E-EL-MK1206	Monkey OLAB (Anti-Oxidized Low Density Lipoprotein Antibody) ELISA Kit
E-EL-MK1208	Monkey GHITM (Growth Hormone Inducible Transmembrane Protein) ELISA Kit
E-EL-MK1209	Monkey CASP11 (Caspase 11) ELISA Kit
E-EL-MK1210	Monkey CCND3 (Cyclin-D3) ELISA Kit
E-EL-MK1211	Monkey Alk-Smase (Alkaline Sphingomyelinase) ELISA Kit
E-EL-MK1212	Monkey FGα (Fibrinogen Alpha) ELISA Kit
E-EL-MK1213	Monkey LDL (Low Density Lipoprotein) ELISA Kit
E-EL-MK1214	Monkey 16-α OHE-1 (16-α Hydroxyestrogen 1) ELISA Kit
E-EL-MK1215	Monkey TCF20 (Transcription Factor 20) ELISA Kit
E-EL-MK1216	Monkey EMP (Erythrocyte Membrane Protein) ELISA Kit
E-EL-MK1217	Monkey LEP (Leptin) ELISA Kit
E-EL-MK1219	Monkey PRDX2 (Peroxiredoxin 2) ELISA Kit
E-EL-MK1220	Monkey GAD2 (Glutamate Decarboxylase 2) ELISA Kit
E-EL-MK1222	Monkey H-Cadherin (Cadherin, Heart) ELISA Kit
E-EL-MK1224	Monkey SEP1 (Selenoprotein 1) ELISA Kit
E-EL-MK1226	Monkey TLR9 (Toll-Like Receptor 9) ELISA Kit
E-EL-MK1227	Monkey MT1 (Metallothionein 1) ELISA Kit
E-EL-MK1228	Monkey ADIPOR2 (Adiponectin Receptor 2) ELISA Kit
E-EL-MK1229	Monkey TM (Thrombomodulin) ELISA Kit

Cat.NO	Product Name
E-EL-MK1230	Monkey PP (Pancreatic Polypeptide) ELISA Kit
E-EL-MK1231	Monkey NUP98 (Nucleoporin 98kDa) ELISA Kit
E-EL-MK1232	Monkey GA (Glycated Albumin) ELISA Kit
E-EL-MK1233	Monkey PI3K2 α (Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit
E-EL-MK1234	Monkey ATGA/TGAB (Anti-Thyroid-Globulin Antibody) ELISA Kit
E-EL-MK1235	Monkey CTX-II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit
E-EL-MK1236	Monkey PTF1 α (Pancreas Specific Transcription Factor 1 α) ELISA Kit
E-EL-MK1237	Monkey P-PKC (Phospho Protein Kinase C) ELISA Kit
E-EL-MK1238	Monkey THAA (Thyroid Hormone Autoantibody) ELISA Kit
E-EL-MK1239	Monkey SDC3 (Syndecan 3) ELISA Kit
E-EL-MK1240	Monkey REG3 γ (Regenerating Islet Derived Protein 3 Gamma) ELISA Kit
E-EL-MK1241	Monkey FZD8 (Frizzled Homolog 8) ELISA Kit
E-EL-MK1243	Monkey EphA2 (Tyrosine protein kinase receptor A2) ELISA Kit
E-EL-MK1244	Monkey CFH (Complement Factor H) ELISA Kit
E-EL-MK1245	Monkey MMP-9 (Matrix Metalloproteinase 9) ELISA Kit
E-EL-MK1246	Monkey NTRK3 (Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
E-EL-MK1247	Monkey PDH E1 (Pyruvate Dehydrogenase E1) ELISA Kit
E-EL-MK1248	Monkey CR2/CD21 (Complement Receptor 2) ELISA Kit
E-EL-MK1249	Monkey PLA2R1 (Phospholipase A2 Receptor 1) ELISA Kit
E-EL-MK1250	Monkey NSE (Neuron-Specific Enolase) ELISA Kit

Cat.NO	Product Name
E-EL-MK1251	Monkey CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit
E-EL-MK1252	Monkey CVNPF (Cobra Venom Neuronal Protective Factor) ELISA Kit
E-EL-MK1254	Monkey MMP-2 (Matrix Metalloproteinase 2) ELISA Kit
E-EL-MK1255	Monkey SF1 (Steroidogenic Factor 1) ELISA Kit
E-EL-MK1256	Monkey EPHA10 (Ephrin Type A Receptor 10) ELISA Kit
E-EL-MK1257	Monkey PKR1 (Prokineticin Receptor 1) ELISA Kit
E-EL-MK1258	Monkey ODC (Ornithine Decarboxylase) ELISA Kit
E-EL-MK1259	Monkey OPGL (Osteoprotegerin Ligand) ELISA Kit
E-EL-MK1261	Monkey PAR4 (Protease Activated Receptor 4) ELISA Kit
E-EL-MK1262	Monkey APA (Anti-Phospholipid Antibody) ELISA Kit
E-EL-MK1264	Monkey ANG II R2/AGTR2 (Angiotensin II Receptor 2) ELISA Kit
E-EL-MK1265	Monkey MCP-4/CCL13 (Monocyte Chemotactic Protein 4) ELISA Kit
E-EL-MK1266	Monkey SGK2 (Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit
E-EL-MK1268	Monkey EED (Embryonic Ectoderm Development) ELISA Kit
E-EL-MK1270	Monkey AIF1 (Allograft Inflammatory Factor 1) ELISA Kit
E-EL-MK1271	Monkey CASP14 (Caspase 14) ELISA Kit
E-EL-MK1272	Monkey DKK4 (Dickkopf Related Protein 4) ELISA Kit
E-EL-MK1273	Monkey AAA (Anti-Albumin Antibody) ELISA Kit
E-EL-MK1274	Monkey S100A8 (S100 Calcium Binding Protein A8/Calgranulin A) ELISA Kit
E-EL-MK1275	Monkey ACTα2 (Actin Alpha 2, Smooth Muscle) ELISA Kit

Cat.NO	Product Name
E-EL-MK1276	Monkey OPG (Osteoprotegerin) ELISA Kit
E-EL-MK1277	Monkey TSLP (Thymic Stromal Lymphopoietin) ELISA Kit
E-EL-MK1278	Monkey BPI (Bactericidal/Permeability Increasing Protein) ELISA Kit
E-EL-MK1279	Monkey Flt3 (FMS Like Tyrosine Kinase 3) ELISA Kit
E-EL-MK1280	Monkey GNβ1 (G Protein Beta 1) ELISA Kit
E-EL-MK1281	Monkey TAF13 (TATA Box Binding Protein Associated Factor 13) ELISA Kit
E-EL-MK1282	Monkey ANG III (Angiotensin III) ELISA Kit
E-EL-MK1283	Monkey MTL (Motilin) ELISA Kit
E-EL-MK1284	Monkey CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) ELISA Kit
E-EL-MK1285	Monkey SCTR (Secretin Receptor) ELISA Kit
E-EL-MK1286	Monkey GP-73 (Golgi Protein 73) ELISA Kit
E-EL-MK1287	Monkey M-AChR (Muscarinic Acetylcholine Receptor) ELISA Kit
E-EL-MK1288	Monkey GHR (Growth Hormone Receptor) ELISA Kit
E-EL-MK1289	Monkey ANGPTL2 (Angiopoietin Like Protein 2) ELISA Kit
E-EL-MK1291	Monkey AIM2 (Absent In Melanoma 2) ELISA Kit
E-EL-MK1292	Monkey GP4 (Platelet Membrane Glycoprotein IV) ELISA Kit
E-EL-MK1294	Monkey PPAR-δ (Peroxisome Proliferator Activated Receptor Delta) ELISA Kit
E-EL-MK1295	Monkey MUC5AC (Mucin 5 Subtype AC) ELISA Kit
E-EL-MK1298	Monkey UGCG (UDP-glucose ceramide glucosyltransferase) ELISA Kit
E-EL-MK1299	Monkey FABP1 (Fatty Acid Binding Protein 1, Liver) ELISA Kit

Cat.NO	Product Name
E-EL-MK1300	Monkey GUS β (β -glucuronidase) ELISA Kit
E-EL-MK1301	Monkey PLA2G12B (Phospholipase A2, Group XII B) ELISA Kit
E-EL-MK1302	Monkey FE (Ferritin) ELISA Kit
E-EL-MK1303	Monkey UreC (Urease Related Protein C) ELISA Kit
E-EL-MK1305	Monkey GTEAGE (Glomerular Tissue Extract Advanced Glycosylation End Products) ELISA Kit
E-EL-MK1306	Monkey DLAT (Dihydrolipoyl Transacetylase) ELISA Kit
E-EL-MK1307	Monkey TF (Tissue Factor) ELISA Kit
E-EL-MK1308	Monkey Anti-MPO (Anti-Myeloperoxidase Antibody) ELISA Kit
E-EL-MK1309	Monkey α 1-AGP (Alpha-1-Acid Glycoprotein) ELISA Kit
E-EL-MK1310	Monkey TEF (Thyrotroph Embryonic Factor) ELISA Kit
E-EL-MK1312	Monkey TNK2 (Tyrosine Kinase, Non Receptor 2) ELISA Kit
E-EL-MK1313	Monkey GMCSF/CSF2(Granulocyte Macrophage Colony Stimulating Factor 2) ELISA Kit
E-EL-MK1314	Monkey PTGS1 (Prostaglandin Endoperoxide Synthase 1) ELISA Kit
E-EL-MK1315	Monkey Pro-MMP-1 (pro-Matrix Metalloproteinase 1) ELISA Kit
E-EL-MK1317	Monkey PA (Prealbumin) ELISA Kit
E-EL-MK1318	Monkey α -CTx (Alpha Crosslaps) ELISA Kit
E-EL-MK1319	Monkey ADRP (Adipose Differentiation Related Protein) ELISA Kit
E-EL-MK1320	Monkey ARSA (Arylsulfatase A) ELISA Kit
E-EL-MK1322	Monkey ES-Ab (Anti-Endostatin Antibody) ELISA Kit
E-EL-MK1324	Monkey GAL7 (Galectin 7) ELISA Kit

Cat.NO	Product Name
E-EL-MK1325	Monkey CAMK2 α (Calcium/Calmodulin Dependent Protein Kinase II Alpha) ELISA Kit
E-EL-MK1326	Monkey DLK-1 (Delta Like 1 Homolog) ELISA Kit
E-EL-MK1328	Monkey MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit
E-EL-MK1330	Monkey MR (Mannose Receptor) ELISA Kit
E-EL-MK1333	Monkey α 2-M (Alpha-2 Macroglobulin) ELISA Kit
E-EL-MK1334	Monkey TUB α 3C (Tubulin Alpha 3c)ELISA Kit
E-EL-MK1335	Monkey GATA2 (GATA Binding Protein 2) ELISA Kit
E-EL-MK1336	Monkey INSR (Insulin Receptor) ELISA Kit
E-EL-MK1337	Monkey MMP-17 (Matrix Metalloproteinase 17) ELISA Kit
E-EL-MK1338	Monkey CaSR (Calcium Sensing Receptor) ELISA Kit
E-EL-MK1340	Monkey CD97 (Cluster of Differentiation 97) ELISA Kit
E-EL-MK1341	Monkey IAPP (Islet Amyloid Polypeptide) ELISA Kit
E-EL-MK1342	Monkey ANG1 (Angiopoietin 1) ELISA Kit
E-EL-MK1343	Monkey IL-11 (Interleukin 11) ELISA Kit
E-EL-MK1344	Monkey OSMR (Oncostatin M Receptor) ELISA Kit
E-EL-MK1345	Monkey TP63 (Tumor Protein p63) ELISA Kit
E-EL-MK1346	Monkey GT-34 (Gastrin 34) ELISA Kit
E-EL-MK1349	Monkey NEF3 (Neurofilament 3) ELISA Kit
E-EL-MK1351	Monkey GAS (N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
E-EL-MK1353	Monkey BANP (BTG3 Associated Nuclear Protein) ELISA Kit

Cat.NO	Product Name
E-EL-MK1354	Monkey FcγR3A (Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit
E-EL-MK1355	Monkey MSMβ (Microseminoprotein Beta) ELISA Kit
E-EL-MK1356	Monkey PA2G4 (Proliferation Associated Protein 2G4) ELISA Kit
E-EL-MK1357	Monkey Ntn4 (Netrin 4) ELISA Kit
E-EL-MK1358	Monkey fT3 (Free Triiodothyronine) ELISA Kit
E-EL-MK1359	Monkey BNP (Brain Natriuretic Peptide) ELISA Kit
E-EL-MK1360	Monkey LCN2 (Lipocalin 2) ELISA Kit
E-EL-MK1361	Monkey Smad7/MADH7 (Mothers Against Decapentaplegic Homolog 7) ELISA Kit
E-EL-MK1362	Monkey BMF (Bcl-2 Modifying Factor) ELISA Kit
E-EL-MK1365	Monkey C1q (Complement 1q) ELISA Kit
E-EL-MK1367	Monkey CTSD (Cathepsin D) ELISA Kit
E-EL-MK1368	Monkey FBN3 (Fibrillin 3) ELISA Kit
E-EL-MK1369	Monkey BDNF(Brain Derived Neurotrophic Factor) ELISA Kit
E-EL-MK1370	Monkey sTfR1 (Soluble Transferrin Receptor1) ELISA Kit
E-EL-MK1371	Monkey FABP6 (Fatty Acid Binding Protein 6, Ileal) ELISA Kit
E-EL-MK1372	Monkey CKLFSF5 (Chemokine Like Factor Superfamily 5) ELISA Kit
E-EL-MK1373	Monkey MMP-7 (Matrix Metalloproteinase 7) ELISA Kit
E-EL-MK1374	Monkey RCOR3 (Rest Corepressor 3) ELISA Kit
E-EL-MK1375	Monkey D-LDH (D-Lactate Dehydrogenase) ELISA Kit
E-EL-MK1377	Monkey LIPC (Lipase, Hepatic) ELISA Kit

Cat.NO	Product Name
E-EL-MK1378	Monkey GDF2 (Growth Differentiation Factor 2) ELISA Kit
E-EL-MK1379	Monkey PRDX5 (Peroxisomal Peroxidase 5) ELISA Kit
E-EL-MK1380	Monkey VCL (Vinculin) ELISA Kit
E-EL-MK1381	Monkey STAT4 (Signal Transducer and Activator of Transcription 4) ELISA Kit
E-EL-MK1382	Monkey MYCBP (C-Myc Binding Protein) ELISA Kit
E-EL-MK1383	Monkey AVPR1A (Arginine Vasopressin Receptor 1A) ELISA Kit
E-EL-MK1385	Monkey BCAR1 (Breast Cancer Anti Estrogen Resistance 1) ELISA Kit
E-EL-MK1386	Monkey CFI (Complement Factor I) ELISA Kit
E-EL-MK1387	Monkey NNE (Non-Neuronal Enolase) ELISA Kit
E-EL-MK1388	Monkey AIRA (Anti-Insulin Receptor Antibody) ELISA Kit
E-EL-MK1389	Monkey PDPN (Podoplanin) ELISA Kit
E-EL-MK1390	Monkey VacA (Vacuolating Cytotoxin A) ELISA Kit
E-EL-MK1391	Monkey FUCA (Alpha-L-Fucosidase, Tissue) ELISA Kit
E-EL-MK1392	Monkey PYY (Peptide YY) ELISA Kit
E-EL-MK1393	Monkey RNASE (Ribonuclease) ELISA Kit
E-EL-MK1395	Monkey PON1 (Paraoxonase 1) ELISA Kit
E-EL-MK1396	Monkey TLN (Talin) ELISA Kit
E-EL-MK1397	Monkey PAPPA (Pregnancy Associated Plasma Protein A/ Pappalysin-1) ELISA Kit
E-EL-MK1398	Monkey SERT (Serotonin Transporter) ELISA Kit
E-EL-MK1399	Monkey PAR3 (Protease Activated Receptor 3) ELISA Kit

Cat.NO	Product Name
E-EL-MK1400	Monkey RELB (V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
E-EL-MK1401	Monkey DDC (Dopamine Decarboxylase) ELISA Kit
E-EL-MK1402	Monkey Smad1/MADH1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit
E-EL-MK1403	Monkey LCN1 (Lipocalin 1) ELISA Kit
E-EL-MK1404	Monkey CD30L (Cluster of Differentiation 30 Ligand) ELISA Kit
E-EL-MK1406	Monkey ADA (Adenosine Deaminase) ELISA Kit
E-EL-MK1407	Monkey TXLNa (Taxilin Alpha) ELISA Kit
E-EL-MK1408	Monkey PLD2 (Phospholipase D2) ELISA Kit
E-EL-MK1409	Monkey INSL5 (Insulin Like Protein 5) ELISA Kit
E-EL-MK1410	Monkey FSH (Follicle-Stimulating Hormone) ELISA Kit
E-EL-MK1411	Monkey CASP 9 (Caspase 9) ELISA Kit
E-EL-MK1412	Monkey DPP4 (Dipeptidyl Peptidase IV) ELISA Kit
E-EL-MK1414	Monkey ANXA2 (Annexin A2) ELISA Kit
E-EL-MK1416	Monkey NPHN (Nephrin) ELISA Kit
E-EL-MK1417	Monkey GSTω1 (Glutathione S Transferase Omega 1) ELISA Kit
E-EL-MK1418	Monkey IL-23 p40(Interleukin 23) ELISA Kit
E-EL-MK1419	Monkey EAR12 (Eosinophil Associated Ribonuclease A family Member 12) ELISA Kit
E-EL-MK1420	Monkey IL-3 (Interleukin 3) ELISA Kit
E-EL-MK1421	Monkey IGFBP-6 (Insulin Like Growth Factor Binding Protein 6) ELISA Kit
E-EL-MK1422	Monkey COL1α1 (Collagen Type I Alpha 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK1424	Monkey CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit
E-EL-MK1425	Monkey DCN (Decorin) ELISA Kit
E-EL-MK1426	Monkey CCNB (Cyclin-B) ELISA Kit
E-EL-MK1427	Monkey ANG-R-Tie1 (Angiopoietin Receptor Tie1) ELISA Kit
E-EL-MK1428	Monkey CS (Chondroitin Sulfate) ELISA Kit
E-EL-MK1429	Monkey GDF1 (Growth Differentiation Factor 1) ELISA Kit
E-EL-MK1430	Monkey MAO (Monoamine Oxidase) ELISA Kit
E-EL-MK1432	Monkey ELN (Elastin) ELISA Kit
E-EL-MK1433	Monkey Cav-1 (Caveolin-1) ELISA Kit
E-EL-MK1434	Monkey LCA/CD45(Leukocyte Common Antigen) ELISA Kit
E-EL-MK1436	Monkey PKP2 (Plakophilin 2) ELISA Kit
E-EL-MK1437	Monkey RPA-70 (replication protein A) ELISA Kit
E-EL-MK1438	Monkey GP39 (Glycoprotein 39, Cartilage) ELISA Kit
E-EL-MK1439	Monkey NKB (Neurokinin B) ELISA Kit
E-EL-MK1440	Monkey GMFβ (Glia Maturation Factor, Beta) ELISA Kit
E-EL-MK1441	Monkey NPYR2 (Neuropeptide Y Receptor Y2) ELISA Kit
E-EL-MK1442	Monkey MRP1 (Multidrug Resistance Associated Protein 1) ELISA Kit
E-EL-MK1443	Monkey STMN1 (Stathmin 1) ELISA Kit
E-EL-MK1444	Monkey CALU (Calumenin) ELISA Kit
E-EL-MK1445	Monkey F XII a (Activated Coagulation Factor XII) ELISA Kit

Cat.NO	Product Name
E-EL-MK1446	Monkey 5-NT (5-Nucleotidase) ELISA Kit
E-EL-MK1447	Monkey ApoC4 (Apolipoprotein C4) ELISA Kit
E-EL-MK1448	Monkey PA-IgM (Platelet Associated Antibody) ELISA Kit
E-EL-MK1449	Monkey ELA3B (Elastase 3B) ELISA Kit
E-EL-MK1450	Monkey TA (Anti-Teichoic Acid Antibody) ELISA Kit
E-EL-MK1451	Monkey GADL1 (Glutamate Decarboxylase Like Protein 1) ELISA Kit
E-EL-MK1452	Monkey Apo A1-Ab (Anti-Apolipoprotein A1 antibody) ELISA Kit
E-EL-MK1453	Monkey GSTpi (Glutathione S Transferases Pi) ELISA Kit
E-EL-MK1454	Monkey PS (Phosphatidylserine) ELISA Kit
E-EL-MK1455	Monkey CA19-9 (Carbohydrate antigen19-9) ELISA Kit
E-EL-MK1456	Monkey GYPE (Glycophorin E) ELISA Kit
E-EL-MK1457	Monkey TFF1 (Trefoil Factor 1) ELISA Kit
E-EL-MK1460	Monkey DNase-I (Deoxyribonuclease I) ELISA Kit
E-EL-MK1461	Monkey GRN (Granulin) ELISA Kit
E-EL-MK1463	Monkey NT-proBNP (N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
E-EL-MK1464	Monkey CDK5 (Cyclin Dependent Kinase 5) ELISA Kit
E-EL-MK1466	Monkey NTXI (Cross Linked N-telopeptide of Type I Collagen) ELISA Kit
E-EL-MK1467	Monkey LT-E4 (Leukotriene E4) ELISA Kit
E-EL-MK1468	Monkey TP/TXA2R (Thromboxane Receptor) ELISA Kit
E-EL-MK1470	Monkey CALD (Caldesmon) ELISA Kit

Cat.NO	Product Name
E-EL-MK1471	Monkey PRSS2 (Protease, Serine, 2) ELISA Kit
E-EL-MK1472	Monkey F-TESTO (Free Testosterone) ELISA Kit
E-EL-MK1473	Monkey δ SIP- α (Delta Sleep Inducing Peptide α) ELISA Kit
E-EL-MK1474	Monkey PLP (Proteolipid Protein) ELISA Kit
E-EL-MK1475	Monkey PTGES2 (Prostaglandin E Synthase 2) ELISA Kit
E-EL-MK1476	Monkey CCL4L1 (Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit
E-EL-MK1478	Monkey DPP6 (Dipeptidyl Peptidase VI) ELISA Kit
E-EL-MK1479	Monkey ANNA1 (Anti-Neuronal Nuclear Antibody 1) ELISA Kit
E-EL-MK1480	Monkey MIA1 (Melanoma Inhibitory Activity Protein 1) ELISA Kit
E-EL-MK1481	Monkey NE/ELA2 (Elastase 2, Neutrophil) ELISA Kit
E-EL-MK1483	Monkey SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit
E-EL-MK1484	Monkey PRSS12 (Protease, Serine 12) ELISA Kit
E-EL-MK1485	Monkey OSF (Osteoclast Stimulating Factor 1) ELISA Kit
E-EL-MK1486	Monkey PKI γ (Protein Kinase Inhibitor Gamma) ELISA Kit
E-EL-MK1487	Monkey MBLR (Mannose Binding Lectin Receptor) ELISA Kit
E-EL-MK1488	Monkey PDGFR α (Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
E-EL-MK1489	Monkey BAFF/CD257/TNFSF13B (B-Cell Activating Factor) ELISA Kit
E-EL-MK1490	Monkey TECK/CCL25 (Thymus Expressed Chemokine) ELISA Kit
E-EL-MK1491	Monkey ApoA1 (Apolipoprotein A1) ELISA Kit
E-EL-MK1492	Monkey FABP7 (Fatty Acid Binding Protein 7, Brain) ELISA Kit

Cat.NO	Product Name
E-EL-MK1493	Monkey PDGFR β (Platelet Derived Growth Factor Receptor Beta) ELISA Kit
E-EL-MK1494	Monkey ACA-IgA (Anti-Cardiolipin Antibody IgA) ELISA Kit
E-EL-MK1495	Monkey α 2-AP (α 2-Antiplasmin) ELISA Kit
E-EL-MK1496	Monkey FOXO4 (Forkhead Box Protein O4) ELISA Kit
E-EL-MK1498	Monkey GATA4 (GATA Binding Protein 4) ELISA Kit
E-EL-MK1499	Monkey GKRP (Glucokinase Regulatory Protein) ELISA Kit
E-EL-MK1500	Monkey NP (Nandrolone Phenylpropionate) ELISA Kit
E-EL-MK1501	Monkey GAP43 (Growth Associated Protein 43) ELISA Kit
E-EL-MK1502	Monkey RI (Ribonuclease Inhibitor) ELISA Kit
E-EL-MK1503	Monkey EAR2 (Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit
E-EL-MK1504	Monkey ABCG2 (ATP Binding Cassette Transporter G2) ELISA Kit
E-EL-MK1506	Monkey UBAP2 (Ubiquitin Associated Protein 2) ELISA Kit
E-EL-MK1507	Monkey IL-23 (Interleukin 23) ELISA Kit
E-EL-MK1509	Monkey RXFP1 (Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
E-EL-MK1510	Monkey REN (Renin) ELISA Kit
E-EL-MK1513	Monkey PLP1 (Proteolipid Protein 1, Myelin) ELISA Kit
E-EL-MK1514	Monkey CHD3 (Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
E-EL-MK1515	Monkey CK-13/KRT13 (Cytokeratin 13) ELISA Kit
E-EL-MK1516	Monkey PKC η (Protein Kinase C Eta) ELISA Kit
E-EL-MK1517	Monkey β APP (β Amyloid Precursor Protein) ELISA Kit

Cat.NO	Product Name
E-EL-MK1519	Monkey ET-1 (Endothelin 1) ELISA Kit
E-EL-MK1520	Monkey GEF (Guanylate Exchange Factor) ELISA Kit
E-EL-MK1522	Monkey F7 (Coagulation Factor VII) ELISA Kit
E-EL-MK1523	Monkey ADH1 (Alcohol Dehydrogenase 1) ELISA Kit
E-EL-MK1524	Monkey Anti-AlaRS (Anti-Alanyl-tRNA Synthetase/Anti-PL12-Antibody) ELISA Kit
E-EL-MK1525	Monkey FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
E-EL-MK1526	Monkey AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) ELISA Kit
E-EL-MK1527	Monkey F9 (Coagulation Factor IX) ELISA Kit
E-EL-MK1528	Monkey ENG/CD105 (Endoglin) ELISA Kit
E-EL-MK1529	Monkey Tn-I (Troponin I) ELISA Kit
E-EL-MK1530	Monkey PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit
E-EL-MK1531	Monkey APC (Activated Protein C) ELISA Kit
E-EL-MK1532	Monkey GAL (Galanin) ELISA Kit
E-EL-MK1534	Monkey FLNC (Filamin C, Gamma) ELISA Kit
E-EL-MK1535	Monkey PCK1 (Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
E-EL-MK1536	Monkey CLA (Collagen Autoantibody) ELISA Kit
E-EL-MK1537	Monkey UBD (Ubiquitin D) ELISA Kit
E-EL-MK1538	Monkey WARS (Tryptophanyl tRNA Synthetase) ELISA Kit
E-EL-MK1540	Monkey Fc α R (Fc Fragment of IgA Receptor) ELISA Kit
E-EL-MK1541	Monkey FGF21(Fibroblast Growth Factor 21) ELISA Kit

Cat.NO	Product Name
E-EL-MK1542	Monkey LR (Leukoregulin) ELISA Kit
E-EL-MK1543	Monkey VEGFR3/Flt4 (Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit
E-EL-MK1544	Monkey AGER (Advanced Glycosylation End Product Specific Receptor) ELISA Kit
E-EL-MK1545	Monkey PRSS (Serine Protease) ELISA Kit
E-EL-MK1547	Monkey MFI2 (Melanotransferrin) ELISA Kit
E-EL-MK1548	Monkey GBP4 (Guanylate Binding Protein 4) ELISA Kit
E-EL-MK1551	Monkey FGL1 (Fibrinogen Like Protein 1) ELISA Kit
E-EL-MK1552	Monkey BMPR1A (Bone Morphogenetic Protein Receptor, type IA) ELISA Kit
E-EL-MK1553	Monkey TGM1 (Transglutaminase 1, Keratinocyte) ELISA Kit
E-EL-MK1556	Monkey GHRH (Growth Hormone Releasing Hormone) ELISA Kit
E-EL-MK1557	Monkey CG (Cryoglobulin) ELISA Kit
E-EL-MK1558	Monkey Grb10 (Growth Factor Receptor Bound Protein 10) ELISA Kit
E-EL-MK1559	Monkey AP 13 (Apelin 13) ELISA Kit
E-EL-MK1560	Monkey NUP62 (Nucleoporin 62kDa) ELISA Kit
E-EL-MK1562	Monkey FG γ (Fibrinogen Gamma) ELISA Kit
E-EL-MK1563	Monkey BACE2 (β -site APP Cleaving Enzyme 2) ELISA Kit
E-EL-MK1564	Monkey AGT (Angiotensinogen) ELISA Kit
E-EL-MK1565	Monkey CA2 (Carbonic Anhydrase II) ELISA Kit
E-EL-MK1567	Monkey SCF (Stem Cell Factor) ELISA Kit
E-EL-MK1568	Monkey FPR2 (Formyl Peptide Receptor 2) ELISA Kit

Cat.NO	Product Name
E-EL-MK1569	Monkey EPHA3 (Ephrin Type A Receptor 3) ELISA Kit
E-EL-MK1570	Monkey ACTg1 (Actin Gamma 1) ELISA Kit
E-EL-MK1571	Monkey PICP (Procollagen I C-terminal Propeptide) ELISA Kit
E-EL-MK1573	Monkey VIP (Vasoactive Intestinal Peptide) ELISA Kit
E-EL-MK1574	Monkey C4BP α (C4 Binding Protein Alpha) ELISA Kit
E-EL-MK1575	Monkey CETP (Cholesteryl Ester Transfer Protein) ELISA Kit
E-EL-MK1576	Monkey CYPB (Cyclophilin B) ELISA Kit
E-EL-MK1577	Monkey IL-12 p40 (Interleukin 12 p40) ELISA Kit
E-EL-MK1578	Monkey ApoA2 (Apolipoprotein A2) ELISA Kit
E-EL-MK1579	Monkey IFN- β (Interferon Beta) ELISA Kit
E-EL-MK1580	Monkey DBP (DNA Binding Protein) ELISA Kit
E-EL-MK1581	Monkey RNASET2 (Ribonuclease T2) ELISA Kit
E-EL-MK1582	Monkey COMP (Cartilage Oligomeric Matrix Protein) ELISA Kit
E-EL-MK1583	Monkey CVL/EZR (Cytovillin/Ezrin) ELISA Kit
E-EL-MK1584	Monkey OSCAR (Osteoclast Associated Receptor) ELISA Kit
E-EL-MK1586	Monkey GAL14 (Galectin 14) ELISA Kit
E-EL-MK1587	Monkey PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
E-EL-MK1588	Monkey TALLA-1 (T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit
E-EL-MK1589	Monkey TFF3 (Trefoil Factor 3, Intestinal) ELISA Kit
E-EL-MK1591	Monkey OBTF2 (Octamer Binding Transcription Factor 2) ELISA Kit

Cat.NO	Product Name
E-EL-MK1592	Monkey FRA (Fibronectin Related Antigen) ELISA Kit
E-EL-MK1593	Monkey ETRB (Endothelin Receptor B) ELISA Kit
E-EL-MK1594	Monkey TBP (TATA Binding Protein) ELISA Kit
E-EL-MK1595	Monkey PKC ζ (Protein Kinase C Zeta) ELISA Kit
E-EL-MK1596	Monkey EL/LIPG (Endothelial Lipase) ELISA Kit
E-EL-MK1597	Monkey CASP7 (Caspase 7) ELISA Kit
E-EL-MK1598	Monkey ASK-1 (Apoptosis Signal Regulating Kinase 1) ELISA Kit
E-EL-MK1599	Monkey AQP-1 (Aquaporin 1) ELISA Kit
E-EL-MK1601	Monkey bTC (Betacellulin) ELISA Kit
E-EL-MK1602	Monkey PIK3C β (Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit
E-EL-MK1603	Monkey NRIP1 (Nuclear Receptor Interacting Protein 1) ELISA Kit
E-EL-MK1604	Monkey PCX (Podocalyxin) ELISA Kit
E-EL-MK1605	Monkey α Actin (Alpha-Actin) ELISA Kit
E-EL-MK1606	Monkey MMP-13 (Matrix Metalloproteinase 13) ELISA Kit
E-EL-MK1607	Monkey I κ B α (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
E-EL-MK1609	Monkey IL-15 (Interleukin 15) ELISA Kit
E-EL-MK1610	Monkey COL4 (Collagen Type IV) ELISA Kit
E-EL-MK1611	Monkey S100A13 (S100 Calcium Binding Protein A13) ELISA Kit
E-EL-MK1612	Monkey C5a(Complement Component 5a)ELISA Kit
E-EL-MK1613	Monkey S100A6 (S100 Calcium Binding Protein A6) ELISA Kit

Cat.NO	Product Name
E-EL-MK1614	Monkey CCSA-2 (Colon Cancer Specific Antigen-2) ELISA Kit
E-EL-MK1615	Monkey ALDOA (Aldolase A, Fructose Bisphosphate) ELISA Kit
E-EL-MK1616	Monkey INSL3 (Insulin Like Protein 3) ELISA Kit
E-EL-MK1617	Monkey COL3 (Collagen Type III) ELISA Kit
E-EL-MK1618	Monkey SEMA3A (Semaphorin 3A) ELISA Kit
E-EL-MK1620	Monkey AMY2 (Amylase Alpha 2, Pancreatic) ELISA Kit
E-EL-MK1621	Monkey TAP (Trypsinogen Activation Peptide) ELISA Kit
E-EL-MK1622	Monkey SIAE (Sialic Acid Acetyltransferase) ELISA Kit
E-EL-MK1623	Monkey SST (Somatostatin) ELISA Kit
E-EL-MK1624	Monkey ACT α 1 (Actin Alpha 1, Skeletal Muscle) ELISA Kit
E-EL-MK1625	Monkey MMRN1 (Multimerin 1) ELISA Kit
E-EL-MK1626	Monkey CA5B (Carbonic Anhydrase VB, Mitochondrial) ELISA Kit
E-EL-MK1627	Monkey IFN- α (Interferon Alpha) ELISA Kit
E-EL-MK1628	Monkey 12-HETE (12-Hydroxyeicosatetraenoic Acid) ELISA Kit
E-EL-MK1629	Monkey MT3 (Metallothionein 3) ELISA Kit
E-EL-MK1630	Monkey GH2 (Growth Hormone 2) ELISA Kit
E-EL-MK1631	Monkey BMP-3 (Bone Morphogenetic Protein 3) ELISA Kit
E-EL-MK1633	Monkey TP5 (Thymopentin) ELISA Kit
E-EL-MK1634	Monkey α -GST (α -glutathione S-transferases) ELISA Kit
E-EL-MK1635	Monkey F10 (Coagulation Factor X) ELISA Kit

Cat.NO	Product Name
E-EL-MK1636	Monkey TGF- α (Transforming Growth Factor α) ELISA Kit
E-EL-MK1637	Monkey SOD3 (Superoxide Dismutase 3, Extracellular) ELISA Kit
E-EL-MK1638	Monkey CFTR (Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit
E-EL-MK1639	Monkey GasR (Gastrin Receptor) ELISA Kit
E-EL-MK1640	Monkey GP6 (Glycoprotein VI, Platelet) ELISA Kit
E-EL-MK1641	Monkey LHRH/GnRH (Luteinizing Hormone Releasing Hormone) ELISA Kit
E-EL-MK1642	Monkey TRY (Trypsin) ELISA Kit
E-EL-MK1644	Monkey LDHB (Lactate Dehydrogenase B) ELISA Kit
E-EL-MK1645	Monkey TNKS1 (Tankyrase 1) ELISA Kit
E-EL-MK1646	Monkey MRC1 (Mannose Receptor C Type 1) ELISA Kit
E-EL-MK1647	Monkey Gas6 (Growth Arrest Specific Protein 6) ELISA Kit
E-EL-MK1648	Monkey NF κ B (Nuclear Factor Kappa B) ELISA Kit
E-EL-MK1649	Monkey OTF2B (Octamer Transcription Factor 2B) ELISA Kit
E-EL-MK1650	Monkey CAP1 (Adenylate Cyclase Associated Protein 1) ELISA Kit
E-EL-MK1651	Monkey SOCS3 (Suppressors of Cytokine Signaling 3) ELISA Kit
E-EL-MK1652	Monkey CYTIP (Cytohesin Interacting Protein) ELISA Kit
E-EL-MK1653	Monkey TCGF- III (Tcell Growth Factor III) ELISA Kit
E-EL-MK1654	Monkey FGFBP2 (Fibroblast Growth Factor Binding Protein 2) ELISA Kit
E-EL-MK1655	Monkey ACP1 (Acid Phosphatase 1) ELISA Kit
E-EL-MK1656	Monkey CA9 (Carbonic Anhydrase IX) ELISA Kit

Cat.NO	Product Name
E-EL-MK1657	Monkey CKB (Creatine Kinase, Brain) ELISA Kit
E-EL-MK1658	Monkey IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) ELISA Kit
E-EL-MK1659	Monkey UBE2C (Ubiquitin Conjugating Enzyme E2C) ELISA Kit
E-EL-MK1660	Monkey CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit
E-EL-MK1661	Monkey CG β (Chorionic Gonadotrophin Beta)ELISA Kit
E-EL-MK1662	Monkey MIP-3 β /ELC/CCL19 (Macrophage Inflammatory Protein 3 β) ELISA Kit
E-EL-MK1663	Monkey ALP (Alkaline Phosphatase) ELISA Kit
E-EL-MK1664	Monkey FOXO1 (Forkhead Box Protein O1) ELISA Kit
E-EL-MK1665	Monkey ApoC2 (Apolipoprotein C2) ELISA Kit
E-EL-MK1666	Monkey PTHR2 (Parathyroid Hormone Receptor 2) ELISA Kit
E-EL-MK1667	Monkey PGAM2 (Phosphoglycerate Mutase 2) ELISA Kit
E-EL-MK1668	Monkey MIgA (Membrane IgA) ELISA Kit
E-EL-MK1669	Monkey MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit
E-EL-MK1670	Monkey PAIg (Platelet Associated Immunoglobulins) ELISA Kit
E-EL-MK1671	Monkey VS (Versican/PG-M) ELISA Kit
E-EL-MK1672	Monkey BCAN (Brevican) ELISA Kit
E-EL-MK1673	Monkey HHIP (Hedgehog Interacting Protein) ELISA Kit
E-EL-MK1674	Monkey PRKG1 (Protein Kinase, cGMP Dependent Type I) ELISA Kit
E-EL-MK1675	Monkey GLTP (Glycolipid Transfer Protein) ELISA Kit
E-EL-MK1676	Monkey PLN (Phospholamban) ELISA Kit

Cat.NO	Product Name
E-EL-MK1677	Monkey CA12 (Carbonic Anhydrase XII) ELISA Kit
E-EL-MK1678	Monkey α -SCA (α -sarcomeric Actin) ELISA Kit
E-EL-MK1679	Monkey SAA2 (Serum Amyloid A2) ELISA Kit
E-EL-MK1680	Monkey AGRN (Agrin) ELISA Kit
E-EL-MK1681	Monkey TACR2 (Tachykinin Receptor 2) ELISA Kit
E-EL-MK1682	Monkey GREM1 (Gremlin 1) ELISA Kit
E-EL-MK1684	Monkey Try2 (Trypsinogen II) ELISA Kit
E-EL-MK1685	Monkey MX1 (Myxovirus Resistance 1) ELISA Kit
E-EL-MK1686	Monkey HSP-70/HSPA9 (Heat Shock Protein 70) ELISA Kit
E-EL-MK1687	Monkey LAT (Linker For Activation of T-cell) ELISA Kit
E-EL-MK1688	Monkey I κ B β (Inhibitory Subunit of NF Kappa B Beta) ELISA Kit
E-EL-MK1689	Monkey ANGPT4 (Angiopoietin 4) ELISA Kit
E-EL-MK1690	Monkey PAQR7 (Progestin and AdipoQ Receptor Family Member VII) ELISA Kit
E-EL-MK1691	Monkey β TG/PBP/CXCL7/NAP2 (Thromboglobulin, Beta) ELISA Kit
E-EL-MK1692	Monkey NUP133 (Nucleoporin 133kDa) ELISA Kit
E-EL-MK1693	Monkey PGRP-S (Peptidoglycan Recognition Protein-S) ELISA Kit
E-EL-MK1694	Monkey DKK1 (Dickkopf Related Protein 1) ELISA Kit
E-EL-MK1695	Monkey MSP (Macrophage Stimulating Protein) ELISA Kit
E-EL-MK1696	Monkey PSPN (Persephin) ELISA Kit
E-EL-MK1697	Monkey PL (Pancreatic Lipase) ELISA Kit

Cat.NO	Product Name
E-EL-MK1699	Monkey L-LDH (L-Lactate Dehydrogenase) ELISA Kit
E-EL-MK1700	Monkey SYCN (Syncollin) ELISA Kit
E-EL-MK1701	Monkey ANG II (Angiotensin II) ELISA Kit
E-EL-MK1702	Monkey NMP-22 (Nuclear Matrix Protein 22) ELISA Kit
E-EL-MK1703	Monkey GAL8 (Galectin 8) ELISA Kit
E-EL-MK1704	Monkey TLA (Thymus Leukemia Antigen) ELISA Kit
E-EL-MK1706	Monkey CYP1B1 (Cytochrome P450 1B1) ELISA Kit
E-EL-MK1707	Monkey NRP1 (Neuropilin 1) ELISA Kit
E-EL-MK1708	Monkey SPTβN4 (Spectrin Beta, Non Erythrocytic 4) ELISA Kit
E-EL-MK1709	Monkey VEGF-D (Vascular Endothelial Growth Factor D) ELISA Kit
E-EL-MK1710	Monkey Pro-MMP-13 (Pro-Matrix Metalloproteinase 13) ELISA Kit
E-EL-MK1711	Monkey BMP-15 (Bone Morphogenetic Protein 15) ELISA Kit
E-EL-MK1712	Monkey NES1 (Nesfatin 1) ELISA Kit
E-EL-MK1713	Monkey IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit
E-EL-MK1714	Monkey FPA (Fibrinopeptide A) ELISA Kit
E-EL-MK1715	Monkey CBX3 (Chromobox Homolog 3) ELISA Kit
E-EL-MK1716	Monkey TAB (Thyroxine Antibody) ELISA Kit
E-EL-MK1717	Monkey PDIA3 (Protein Disulfide Isomerase A3) ELISA Kit
E-EL-MK1718	Monkey F13A1 (Coagulation Factor XIII A1 Polypeptide) ELISA Kit
E-EL-MK1719	Monkey MIP-5 (Macrophage Inflammatory Protein 5) ELISA Kit

Cat.NO	Product Name
E-EL-MK1720	Monkey PGD2-MOX (Prostaglandin D2-Methoxime) ELISA Kit
E-EL-MK1721	Monkey FSCN (Fascin) ELISA Kit
E-EL-MK1722	Monkey COL10 (Collagen Type X) ELISA Kit
E-EL-MK1723	Monkey MAP2K1 (Mitogen Activated Protein Kinase Kinase 1) ELISA Kit
E-EL-MK1724	Monkey DEF α 1 (Defensin Alpha 1, Neutrophil) ELISA Kit
E-EL-MK1725	Monkey TP53 (Tumor Protein p53) ELISA Kit
E-EL-MK1726	Monkey TDP43 (Tar DNA Binding Protein 43KDa)ELISA Kit
E-EL-MK1727	Monkey GST α 1 (Glutathione S Transferase Alpha 1) ELISA Kit
E-EL-MK1728	Monkey PCIII (Procollagen III) ELISA Kit
E-EL-MK1729	Monkey LRP6 (Low Density Lipoprotein Receptor Related Protein) ELISA Kit
E-EL-MK1730	Monkey NOX1 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
E-EL-MK1731	Monkey ACVA (Activin A) ELISA Kit
E-EL-MK1732	Monkey SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-MK1733	Monkey C8b (Complement Component 8b) ELISA Kit
E-EL-MK1734	Monkey GM-CSF Ab (Anti-Granulocyte Macrophage Colony Stimulating Factor Antibody) ELISA Kit
E-EL-MK1735	Monkey LFA-2/CD2 (Lymphocyte Function Associated Antigen 2) ELISA Kit
E-EL-MK1736	Monkey ON (Osteonectin) ELISA Kit
E-EL-MK1737	Monkey α -La (Alpha-Lactalbumin) ELISA Kit
E-EL-MK1739	Monkey sCKR (Soluble Cytokine Receptor) ELISA Kit
E-EL-MK1740	Monkey Slit2 (Slit Homolog 2) ELISA Kit

Cat.NO	Product Name
E-EL-MK1741	Monkey F11 (Coagulation Factor XI) ELISA Kit
E-EL-MK1742	Monkey ApoB (Apolipoprotein B) ELISA Kit
E-EL-MK1743	Monkey EFNA4 (Ephrin A4) ELISA Kit
E-EL-MK1744	Monkey CKLFSF6 (Chemokine Like Factor Superfamily 6) ELISA Kit
E-EL-MK1746	Monkey FABP4 (Fatty Acid-binding Protein 4, Adipocyte) ELISA Kit
E-EL-MK1747	Monkey PAH (Phenylalanine Hydroxylase) ELISA Kit
E-EL-MK1748	Monkey SBEM (Small Breast Epithelial Mucin) ELISA Kit
E-EL-MK1749	Monkey GPLD1 (Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Kit
E-EL-MK1750	Monkey AIRE (Autoimmune Regulator) ELISA Kit
E-EL-MK1751	Monkey NRG-1 (Neuregulin 1) ELISA Kit
E-EL-MK1752	Monkey ADIPOR1 (Adiponectin Receptor 1) ELISA Kit
E-EL-MK1753	Monkey MYBPC3 (Myosin Binding Protein C, Cardiac) ELISA Kit
E-EL-MK1754	Monkey CG(Chorionic Gonadotropin)ELISA Kit
E-EL-MK1755	Monkey HB (Hemoglobin) ELISA Kit
E-EL-MK1756	Monkey HLTF (Helicase Like Transcription Factor) ELISA Kit
E-EL-MK1757	Monkey Flt3L(FMS Like Tyrosine Kinase 3 Ligand) ELISA Kit
E-EL-MK1758	Monkey MYL9 (Myosin Light Chain 9, Regulatory) ELISA Kit
E-EL-MK1759	Monkey HELIX- II (Type II Collagen Helical Peptide) ELISA Kit
E-EL-MK1760	Monkey Preptin ELISA Kit
E-EL-MK1761	Monkey PDI (Protein Disulfide Isomerase) ELISA Kit

Cat.NO	Product Name
E-EL-MK1762	Monkey C5c (Complement C5 Convertase) ELISA Kit
E-EL-MK1763	Monkey LOXL4 (Lysyl Oxidase Like Protein 4) ELISA Kit
E-EL-MK1764	Monkey FGFR1 (Fibroblast Growth Factor Receptor 1) ELISA Kit
E-EL-MK1765	Monkey RNASEN (Ribonuclease III, Nuclear) ELISA Kit
E-EL-MK1766	Monkey FGF23 (Fibroblast Growth Factor 23) ELISA Kit
E-EL-MK1767	Monkey ApoC1 (Apolipoprotein C1) ELISA Kit
E-EL-MK1768	Monkey ELA2B (Elastase 2B) ELISA Kit
E-EL-MK1769	Monkey GFAP (Glial Fibrillary Acidic Protein) ELISA Kit
E-EL-MK1770	Monkey GUK1 (Guanylate Kinase 1) ELISA Kit
E-EL-MK1771	Monkey TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14)ELISA Kit
E-EL-MK1773	Monkey PLG (Plasminogen) ELISA Kit
E-EL-MK1774	Monkey PS (Plasma Anticoagulant Protein S) ELISA Kit
E-EL-MK1775	Monkey Gal (α -Galactoyl) ELISA Kit
E-EL-MK1776	Monkey AIF (Apoptosis Inducing Factor) ELISA Kit
E-EL-MK1778	Monkey PKC ϵ (Protein Kinase C Epsilon) ELISA Kit
E-EL-MK1779	Monkey IL-6R (Interleukin 6 Receptor) ELISA Kit
E-EL-MK1780	Monkey ITIH4 (Inter Alpha-Globulin Inhibitor H4) ELISA Kit
E-EL-MK1781	Monkey RBP4 (Retinol Binding Protein 4, Plasma) ELISA Kit
E-EL-MK1783	Monkey FADS2 (Fatty Acid Desaturase 2) ELISA Kit
E-EL-MK1785	Monkey CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK1786	Monkey CLDN14 (Claudin 14) ELISA Kit
E-EL-MK1787	Monkey CASP3 (Caspase 3) ELISA Kit
E-EL-MK1788	Monkey CD23/FcεRII (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
E-EL-MK1790	Monkey PCP (Procollagen C Peptidase)ELISA Kit
E-EL-MK1791	Monkey HIF-1α(Hypoxia Inducible Factor 1 Alpha)ELISA Kit
E-EL-MK1792	Monkey TSP1 (Thrombin Sensitive Protein 1) ELISA Kit
E-EL-MK1793	Monkey NRG-3 (Neuregulin 3) ELISA Kit
E-EL-MK1794	Monkey MYOG (Myogenin) ELISA Kit
E-EL-MK1795	Monkey GTF2B (General Transcription Factor IIB) ELISA Kit
E-EL-MK1797	Monkey C4BPβ (C4 Binding Protein Beta) ELISA Kit
E-EL-MK1798	Monkey GCSFR (Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
E-EL-MK1799	Monkey MMP-3 (Matrix Metalloproteinase 3) ELISA Kit
E-EL-MK1800	Monkey GAL9 (Galectin 9) ELISA Kit
E-EL-MK1801	Monkey OC/BGP (Osteocalcin/Bone Gla Protein) ELISA Kit
E-EL-MK1802	Monkey GIF (Gastric Intrinsic Factor) ELISA Kit
E-EL-MK1804	Monkey USP8 (Ubiquitin Specific Peptidase 8) ELISA Kit
E-EL-MK1805	Monkey PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit
E-EL-MK1806	Monkey LDLR (Low Density Lipoprotein Receptor) ELISA Kit
E-EL-MK1807	Monkey ERN1 (Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit
E-EL-MK1809	Monkey TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit

Cat.NO	Product Name
E-EL-MK1810	Monkey ETRA (Endothelin Receptor A) ELISA Kit
E-EL-MK1812	Monkey CNTN2 (Contactin 2) ELISA Kit
E-EL-MK1813	Monkey MAG Ab (Anti-Myelin Associated Glycoprotein Antibody) ELISA Kit
E-EL-MK1814	Monkey MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit
E-EL-MK1815	Monkey ProGRP (Pro-Gastrin Releasing Peptide) ELISA Kit
E-EL-MK1816	Monkey IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
E-EL-MK1818	Monkey NUP160 (Nucleoporin 160kDa) ELISA Kit
E-EL-MK1819	Monkey IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit
E-EL-MK1820	Monkey DSG-E1 (Desmogleins 1) ELISA Kit
E-EL-MK1821	Monkey AR (Androgen Receptor) ELISA Kit
E-EL-MK1822	Monkey MYH1 (Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit
E-EL-MK1823	Monkey Extein ELISA Kit
E-EL-MK1824	Monkey A β 1-42 (Amyloid Beta Peptide 1-42) ELISA Kit
E-EL-MK1825	Monkey DSE (Dermatan Sulfate Epimerase) ELISA Kit
E-EL-MK1826	Monkey LDHA (Lactate Dehydrogenase A) ELISA Kit
E-EL-MK1827	Monkey RIPK2 (Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit
E-EL-MK1828	Monkey CAM-Ab (Anti-Calmodulin Specific Antibody) ELISA Kit
E-EL-MK1829	Monkey MPA (Mycophenolic Acid) ELISA Kit
E-EL-MK1830	Monkey TNPO2 (Transportin 2) ELISA Kit
E-EL-MK1831	Monkey VTN (Vitronection) ELISA Kit

Cat.NO	Product Name
E-EL-MK1833	Monkey Grb7 (Growth Factor Receptor Bound Protein 7) ELISA Kit
E-EL-MK1834	Monkey SHH (Hedgehog Homolog, Sonic) ELISA Kit
E-EL-MK1835	Monkey PLGA (Plasminogen Activator) ELISA Kit
E-EL-MK1836	Monkey LRG1 (Leucine Rich Alpha-2-Glycoprotein 1) ELISA Kit
E-EL-MK1837	Monkey OVA sIgG (Ovalbumin Specific IgG) ELISA Kit
E-EL-MK1838	Monkey FBL (Fibrillarin) ELISA Kit
E-EL-MK1839	Monkey OFQ/N (Orphanin FQ/Nociceptin) ELISA Kit
E-EL-MK1840	Monkey PLVAP (Plasmalemma Vesicle Associated Protein) ELISA Kit
E-EL-MK1841	Monkey NRG-4 (Neuregulin 4) ELISA Kit
E-EL-MK1842	Monkey CDK7 (Cyclin Dependent Kinase 7) ELISA Kit
E-EL-MK1843	Monkey STAT5A (Signal Transducer and Activator of Transcription 5A) ELISA Kit
E-EL-MK1844	Monkey FCN1 (Ficolin 1) ELISA Kit
E-EL-MK1845	Monkey FST (Follistatin) ELISA Kit
E-EL-MK1847	Monkey HLA-A (Leukocyte Antigen A) ELISA Kit
E-EL-MK1848	Monkey C5b (Complement Fragment 5b) ELISA Kit
E-EL-MK1849	Monkey ASM (Acid Sphingomyelinase) ELISA Kit
E-EL-MK1852	Monkey L1CAM/CD171 (L1-Cell Adhesion Molecule) ELISA Kit
E-EL-MK1853	Monkey AATF (Apoptosis Antagonizing Transcription Factor) ELISA Kit
E-EL-MK1854	Monkey DDT (D-Dopachrome Tautomerase) ELISA Kit
E-EL-MK1855	Monkey LIPH (Lipase H) ELISA Kit

Cat.NO	Product Name
E-EL-MK1856	Monkey PRKAA1 (Protein Kinase, AMP Activated Alpha 1) ELISA Kit
E-EL-MK1857	Monkey ANGPTL6 (Angiopoietin Like Protein 6) ELISA Kit
E-EL-MK1858	Monkey CNR1 (Cannabinoid Receptor 1, Brain) ELISA Kit
E-EL-MK1859	Monkey CD4 (Cluster of Differentiation 4) ELISA Kit
E-EL-MK1860	Monkey Pax (Paxillin) ELISA Kit
E-EL-MK1861	Monkey TPP (Thrombus Precursor Protein) ELISA Kit
E-EL-MK1862	Monkey CLDN1 (Claudin 1) ELISA Kit
E-EL-MK1863	Monkey PLSCR1 (Phospholipid Scramblase 1) ELISA Kit
E-EL-MK1864	Monkey KNG1 (Kininogen 1) ELISA Kit
E-EL-MK1865	Monkey TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
E-EL-MK1866	Monkey CX40 (Connexin 40) ELISA Kit
E-EL-MK1868	Monkey PTX 3/TSG-14 (Pentraxin 3, Long) ELISA Kit
E-EL-MK1869	Monkey SSTR1 (Somatostatin Receptor 1) ELISA Kit
E-EL-MK1870	Monkey Grb14 (Growth Factor Receptor Bound Protein 14) ELISA Kit
E-EL-MK1871	Monkey SCG2 (Secretogranin II) ELISA Kit
E-EL-MK1872	Monkey P27 (P27 Protein) ELISA Kit
E-EL-MK1873	Monkey MTF (Microtransferrin) ELISA Kit
E-EL-MK1874	Monkey TAGLN (Transgelin) ELISA Kit
E-EL-MK1875	Monkey D4R/DRD4 (Dopamine Receptor D4) ELISA Kit
E-EL-MK1876	Monkey IGF-1 (Insulin Like Growth Factor 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK1877	Monkey PLC γ 1 (Phospholipase C Gamma 1) ELISA Kit
E-EL-MK1878	Monkey ALDOC (Aldolase C, Fructose Bisphosphate) ELISA Kit
E-EL-MK1879	Monkey eIF3a (Eukaryotic Translation Initiation Factor 3a) ELISA Kit
E-EL-MK1880	Monkey Trx (Thioredoxin) ELISA Kit
E-EL-MK1881	Monkey α 1ACT (Alpha1 Antichymotrypsin) ELISA Kit
E-EL-MK1882	Monkey GDF5 (Growth Differentiation Factor 5) ELISA Kit
E-EL-MK1883	Monkey ADAM9 (A Disintegrin And Metalloprotease 9) ELISA Kit
E-EL-MK1884	Monkey IL-5 (Interleukin 5) ELISA Kit
E-EL-MK1885	Monkey PDGF-AB (Platelet Derived Growth Factor AB) ELISA Kit
E-EL-MK1886	Monkey GCK (Glucokinase) ELISA Kit
E-EL-MK1887	Monkey NT-3 (Neurotrophin-3) ELISA Kit
E-EL-MK1888	Monkey PPP1R1A (Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit
E-EL-MK1889	Monkey NAP (Neutrophil Alkaline Phosphatase) ELISA Kit
E-EL-MK1890	Monkey HSF1 (Heat Shock Transcription Factor 1) ELISA Kit
E-EL-MK1891	Monkey M γ /CXCL9 (Monocyte Interferon Gamma Inducing Factor) ELISA Kit
E-EL-MK1892	Monkey NEFL (Neurofilament, Light Polypeptide) ELISA Kit
E-EL-MK1893	Monkey PDGF-BB (Platelet Derived Growth Factor-BB) ELISA Kit
E-EL-MK1894	Monkey GM (Gastric Mucin) ELISA Kit
E-EL-MK1895	Monkey TD-Ag (Thymus Dependent Antigen) ELISA Kit
E-EL-MK1896	Monkey MCSP (Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit

Cat.NO	Product Name
E-EL-MK1897	Monkey CRH (Corticotropin Releasing Hormone) ELISA Kit
E-EL-MK1898	Monkey PCDH β 16 (Protocadherin Beta 16) ELISA Kit
E-EL-MK1899	Monkey AQP-2 (Aquaporin 2) ELISA Kit
E-EL-MK1900	Monkey cTnT/TNNT2 (Troponin T Type 2, Cardiac) ELISA Kit
E-EL-MK1901	Monkey Ank-B (Ankyrin B) ELISA Kit
E-EL-MK1902	Monkey MCP-3 (Monocyte Chemotactic Protein 3) ELISA Kit
E-EL-MK1903	Monkey FLOT2 (Flotillin 2) ELISA Kit
E-EL-MK1904	Monkey BAFFR (B Cell Activating Factor Receptor) ELISA Kit
E-EL-MK1905	Monkey PKI β (Protein Kinase Inhibitor Beta) ELISA Kit
E-EL-MK1906	Monkey FOXO3 (Forkhead Box Protein O3) ELISA Kit
E-EL-MK1907	Monkey COL9 α 3 (Collagen Type IX Alpha 3) ELISA Kit
E-EL-MK1909	Monkey MMSAM (Melanoma Metastasis Surface Adhesion Molecule) ELISA Kit
E-EL-MK1910	Monkey CDX2 (Caudal Type Homeobox Transcription Factor 2) ELISA Kit
E-EL-MK1911	Monkey FGFR3 (Fibroblast Growth Factor Receptor 3) ELISA Kit
E-EL-MK1912	Monkey COX-2R (Cyclooxygenase-2 Receptor) ELISA Kit
E-EL-MK1914	Monkey C3a (Complement Component 3a) ELISA Kit
E-EL-MK1915	Monkey PCSK9 (Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit
E-EL-MK1916	Monkey F8 (Coagulation Factor VIII) ELISA Kit
E-EL-MK1917	Monkey ODF (Osteoclast Differentiation Factor) ELISA Kit
E-EL-MK1918	Monkey F-Box Protein 2(FBXO2) ELISA Kit

Cat.NO	Product Name
E-EL-MK1919	Monkey CNTNAP4 (Contactin Associated Protein 4) ELISA Kit
E-EL-MK1920	Monkey PHC (Primary Hepatic Carcinoma) ELISA Kit
E-EL-MK1921	Monkey LTD4 (Leukotriene D4) ELISA Kit
E-EL-MK1922	Monkey GOD (Glucose Oxidase) ELISA Kit
E-EL-MK1923	Monkey PMN Elastase (Polymorphonuclear Elastase) ELISA Kit
E-EL-MK1925	Monkey DHEA-S7 (Dehydroepiandrosterone S7) ELISA Kit
E-EL-MK1928	Monkey ALT (Alanine Transaminase) ELISA Kit
E-EL-MK1929	Monkey SYP (Synaptophysin) ELISA Kit
E-EL-MK1930	Monkey C4b (Complement Component 4b) ELISA Kit
E-EL-MK1931	Monkey PIK3C α (Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit
E-EL-MK1934	Monkey CRY1 (Cryptochrome 1) ELISA Kit
E-EL-MK1935	Monkey CMLC-1 (Cardiac Myosin Light Chains 1) ELISA Kit
E-EL-MK1936	Monkey β 2-GP (Beta2 Glycoprotein) ELISA Kit
E-EL-MK1938	Monkey TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) ELISA Kit
E-EL-MK1940	Monkey ER α (Estrogen Receptor Alpha) ELISA Kit
E-EL-MK1941	Monkey ApoE (Apolipoprotein E) ELISA Kit
E-EL-MK1942	Monkey AFA (Anti-Filaggrin Antibody) ELISA Kit
E-EL-MK1944	Monkey PR3Ab (Proteinase 3 Antibody) ELISA Kit
E-EL-MK1946	Monkey TL1A (Tumor Necrosis Factor Related Ligand 1A) ELISA Kit
E-EL-MK1947	Monkey SEMA7A (Semaphorin 7A) ELISA Kit

Cat.NO	Product Name
E-EL-MK1948	Monkey PKBG (Protein Kinase B Gamma) ELISA Kit
E-EL-MK1950	Monkey ACP5 (Tartrate Resistant Acid Phosphatase 5) ELISA Kit
E-EL-MK1951	Monkey GAL3 (Galectin 3) ELISA Kit
E-EL-MK1952	Monkey MUC1 (Mucin 1) ELISA Kit
E-EL-MK1953	Monkey CLDN11 (Claudin 11) ELISA Kit
E-EL-MK1954	Monkey ET (Epitestosterone) ELISA Kit
E-EL-MK1956	Monkey IGFBP-2 (Insulin Like Growth Factor Binding Protein 2) ELISA Kit
E-EL-MK1957	Monkey TUB β (β -Tubulin) ELISA Kit
E-EL-MK1958	Monkey TAF1 (TATA Box Binding Protein Associated Factor 1) ELISA Kit
E-EL-MK1960	Monkey Ext-N (N-Terminal Extein) ELISA Kit
E-EL-MK1961	Monkey COL1 α 2 (Collagen Type I Alpha 2) ELISA Kit
E-EL-MK1962	Monkey TFPI2 (Tissue Factor Pathway Inhibitor 2) ELISA Kit
E-EL-MK1963	Monkey AAP (Alanine Aminopeptidase) ELISA Kit
E-EL-MK1965	Monkey PKR2 (Prokineticin Receptor 2) ELISA Kit
E-EL-MK1966	Monkey NUP50 (Nucleoporin 50kDa) ELISA Kit
E-EL-MK1967	Monkey AAA (Anti-Actin Antibody) ELISA Kit
E-EL-MK1968	Monkey ARR β 2 (Arrestin Beta 2) ELISA Kit
E-EL-MK1969	Monkey Arg (Arginase) ELISA Kit
E-EL-MK1970	Monkey τ PK1 (Tau-Protein Kinase) ELISA Kit
E-EL-MK1972	Monkey GDF9 (Growth Differentiation Factor 9) ELISA Kit

Cat.NO	Product Name
E-EL-MK1973	Monkey OBTF1 (Octamer Binding Transcription Factor 1) ELISA Kit
E-EL-MK1974	Monkey Bax (Bcl-2 Associated X Protein) ELISA Kit
E-EL-MK1975	Monkey TM (Thrombin)ELISA Kit
E-EL-MK1976	Monkey PRI (Placental Ribonuclease Inhibitor) ELISA Kit
E-EL-MK1978	Monkey ER β (Estrogen Receptor Beta) ELISA Kit
E-EL-MK1979	Monkey GUCA2B (Guanylate Cyclase Activator 2B) ELISA Kit
E-EL-MK1980	Monkey CPR (Cytochrome P450 Reductase) ELISA Kit
E-EL-MK1981	Monkey PAC3 (Platelet Associated Complement 3) ELISA Kit
E-EL-MK1982	Monkey c-fos ELISA Kit
E-EL-MK1983	Monkey ASGR1 (Asialoglycoprotein Receptor 1) ELISA Kit
E-EL-MK1985	Monkey EP3 (Prostaglandin E Receptor 3) ELISA Kit
E-EL-MK1986	Monkey ERCC1 (Excision Repair Cross Complementing Rodent Repair Deficiency Complementation 1) ELISA Kit
E-EL-MK1987	Monkey ALB (Albumin) ELISA Kit
E-EL-MK1989	Monkey PIAS1 (Protein Inhibitor of Activated STAT 1) ELISA Kit
E-EL-MK1990	Monkey REG4 (Regenerating Islet Derived Protein 4) ELISA Kit
E-EL-MK1991	Monkey PRL(Prolactin)ELISA Kit
E-EL-MK1992	Monkey Neo-T4 (Neonatal Thyroxine) ELISA Kit
E-EL-MK1993	Monkey HNRPA2B1 (Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit
E-EL-MK1994	Monkey LIPF (Lipase, Gastric) ELISA Kit
E-EL-MK1995	Monkey ApoA4 (Apolipoprotein A4) ELISA Kit

Cat.NO	Product Name
E-EL-MK1996	Monkey S100A9 (S100 Calcium Binding Protein A9) ELISA Kit
E-EL-MK1997	Monkey PTPRB (Protein Tyrosine Phosphatase Receptor Type B) ELISA Kit
E-EL-MK1998	Monkey ADCY2 (Adenylate Cyclase 2) ELISA Kit
E-EL-MK1999	Monkey BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit
E-EL-MK2001	Monkey ANGPTL4 (Angiopoietin Like Protein 4) ELISA Kit
E-EL-MK2003	Monkey S100A11 (S100 Calcium Binding Protein A11) ELISA Kit
E-EL-MK2004	Monkey CMR (Costimulatory Molecules Receptor) ELISA Kit
E-EL-MK2005	Monkey OXA (Orexin A) ELISA Kit
E-EL-MK2008	Monkey HTRA1 (HtrA Serine, Peptidase 1) ELISA Kit
E-EL-MK2009	Monkey SPD (Pulmonary Surfactant Associated Protein D) ELISA Kit
E-EL-MK2010	Monkey VPI (Vasoactive Peptidase Inhibitors) ELISA Kit
E-EL-MK2012	Monkey CXCR1 (CXC-Chemokine Receptor 1) ELISA Kit
E-EL-MK2013	Monkey FOXO2 (Forkhead Box Protein O2) ELISA Kit
E-EL-MK2014	Monkey CYP21B (Cytochrome P450 21B) ELISA Kit
E-EL-MK2015	Monkey Salusin α (Salusin Alpha) ELISA Kit
E-EL-MK2016	Monkey CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
E-EL-MK2017	Monkey CPE (Carboxypeptidase E) ELISA Kit
E-EL-MK2018	Monkey DCX (Doublecortin) ELISA Kit
E-EL-MK2019	Monkey IRAK3 (Interleukin 1 Receptor Associated Kinase 3) ELISA Kit
E-EL-MK2021	Monkey PLA2G2D (Phospholipase A2, Group II D) ELISA Kit

Cat.NO	Product Name
E-EL-MK2022	Monkey GAL4 (Galectin 4) ELISA Kit
E-EL-MK2023	Monkey PKP1 (Plakophilin 1) ELISA Kit
E-EL-MK2025	Monkey 15-LO (Arachidonate 15-Lipoxygenase) ELISA Kit
E-EL-MK2026	Monkey HA (Hemagglutinin) ELISA Kit
E-EL-MK2027	Monkey rT3 (Reverse Triiodothyronine) ELISA Kit
E-EL-MK2029	Monkey VLA (Very Late Appearing Antigen) ELISA Kit
E-EL-MK2031	Monkey SOX9 (Sex Determining Region Y Box Protein 9) ELISA Kit
E-EL-MK2033	Monkey MRP (Multidrug Resistance-Associated Protein) ELISA Kit
E-EL-MK2034	Monkey GLP-2 (Glucagon Like Peptide 2) ELISA Kit
E-EL-MK2035	Monkey PKN1 (Protein Kinase N1) ELISA Kit
E-EL-MK2036	Monkey RHOA (Ras Homolog Gene Family, Member A) ELISA Kit
E-EL-MK2037	Monkey BGN (Biglycan) ELISA Kit
E-EL-MK2038	Monkey PHB (Prohibitin) ELISA Kit
E-EL-MK2039	Monkey RAR α (Retinoic Acid Receptor Alpha) ELISA Kit
E-EL-MK2040	Monkey FBN1 (Fibrillin 1) ELISA Kit
E-EL-MK2041	Monkey RLR (Retinoic Acid Inducible Gene/RIG-Like Receptor) ELISA Kit
E-EL-MK2042	Monkey CYTH1 (Cytohesin 1) ELISA Kit
E-EL-MK2043	Monkey Anti-Surv (Anti-Survivin Antibody) ELISA Kit
E-EL-MK2044	Monkey ALD (Aldolase) ELISA Kit
E-EL-MK2045	Monkey TJP1 (Tight Junction Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK2046	Monkey PFKM (Phosphofructokinase, Muscle) ELISA Kit
E-EL-MK2047	Monkey VEGF-C (Vascular Endothelial Growth Factor C) ELISA Kit
E-EL-MK2048	Monkey PPOX (Protoporphyrinogen Oxidase) ELISA Kit
E-EL-MK2049	Monkey ANG-R-Tie2 (Angiopoietin Receptor Tie2) ELISA Kit
E-EL-MK2050	Monkey TGM2 (Transglutaminase 2, Tissue) ELISA Kit
E-EL-MK2051	Monkey M-AChR M2 (Muscarinic Acetylcholine Receptor M2) ELISA Kit
E-EL-MK2052	Monkey MCP/CD46 (Membrane Cofactor Protein) ELISA Kit
E-EL-MK2053	Monkey MASP1 (Mannan Associated Serine Protease 1) ELISA Kit
E-EL-MK2055	Monkey ACAT2 (Acetyl Coenzyme A Acetyltransferase 2) ELISA Kit
E-EL-MK2056	Monkey PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit
E-EL-MK2058	Monkey CNTN1 (Contactin 1) ELISA Kit
E-EL-MK2059	Monkey COL7 (Collagen Type VII) ELISA Kit
E-EL-MK2060	Monkey A-CoA (Acetyl Coenzyme A) ELISA Kit
E-EL-MK2061	Monkey EL (EuglobulinLysis) ELISA Kit
E-EL-MK2062	Monkey 26S-PSM (26S Proteasome) ELISA Kit
E-EL-MK2063	Monkey GPC3 (Glypican 3) ELISA Kit
E-EL-MK2064	Monkey PICK1 (Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
E-EL-MK2065	Monkey GzmA (Granzyme A) ELISA Kit
E-EL-MK2066	Monkey sCD14 (Soluble Cluster of Differentiation14) ELISA Kit
E-EL-MK2067	Monkey FGF-7/KGF (Fibroblast Growth Factor 7) ELISA Kit

Cat.NO	Product Name
E-EL-MK2068	Monkey GROβ (Growth Regulated Oncogene Beta) ELISA Kit
E-EL-MK2069	Monkey GJβ1 (Gap Junction Protein Beta 1) ELISA Kit
E-EL-MK2070	Monkey MAP-2 (Microtubule Associated Protein 2) ELISA Kit
E-EL-MK2072	Monkey VEGF165 (Vascular Endothelial Growth Factor165) ELISA Kit
E-EL-MK2073	Monkey TERT (Telomerase Reverse Transcriptase) ELISA Kit
E-EL-MK2074	Monkey UCRP (Ubiquitin Cross Reactive Protein) ELISA Kit
E-EL-MK2075	Monkey GPC1 (Glypican 1) ELISA Kit
E-EL-MK2076	Monkey S100A7 (S100 Calcium Binding Protein A7) ELISA Kit
E-EL-MK2077	Monkey SCFR/c-Kit (Stem Cell Factor Receptor) ELISA Kit
E-EL-MK2078	Monkey HSPβ6/HSP-20 (Heat Shock Protein Beta 6) ELISA Kit
E-EL-MK2079	Monkey CLU (Clusterin) ELISA Kit
E-EL-MK2080	Monkey DNM1 (Dynamin 1) ELISA Kit
E-EL-MK2081	Monkey PLCγ2 (Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit
E-EL-MK2082	Monkey FAIM3 (Fas Apoptotic Inhibitory Molecule 3) ELISA Kit
E-EL-MK2083	Monkey EGR3 (Early Growth Response Protein 3) ELISA Kit
E-EL-MK2084	Monkey IL-16 (Interleukin 16) ELISA Kit
E-EL-MK2085	Monkey IgG3 (Immunoglobulin G3) ELISA Kit
E-EL-MK2086	Monkey NUP35 (Nucleoporin 35kDa) ELISA Kit
E-EL-MK2087	Monkey PLCβ1/PI-PLC (Phospholipase C Beta 1, Phosphoinositide Specific) ELISA Kit
E-EL-MK2088	Monkey MSR1 (Macrophage Scavenger Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-MK2089	Monkey SOX1 (Sex Determining Region Y Box Protein 1) ELISA Kit
E-EL-MK2090	Monkey ELMO2 (Engulfment And Cell Motility 2) ELISA Kit
E-EL-MK2091	Monkey FcεR I (Receptor I for the Fc Fragment of IgE) ELISA Kit
E-EL-MK2092	Monkey CK-LMW (Cytokeratin Low Molecular Weight) ELISA Kit
E-EL-MK2094	Monkey PYGL (Glycogen Phosphorylase, Liver) ELISA Kit
E-EL-MK2095	Monkey PRB (Retinoblastoma Tumor Suppressor Protein) ELISA Kit
E-EL-MK2096	Monkey GPC4 (Glypican 4) ELISA Kit
E-EL-MK2097	Monkey PL (Placental Lactogen) ELISA Kit
E-EL-MK2098	Monkey AGPA (Anti-Glucoprotein 210 Antibody) ELISA Kit
E-EL-MK2099	Monkey PXR (Pregnane X Receptor) ELISA Kit
E-EL-MK2100	Monkey SNCα (Synuclein Alpha) ELISA Kit
E-EL-MK2101	Monkey FETUB (Fetuin B) ELISA Kit
E-EL-MK2102	Monkey PDGF-AA (Platelet Derived Growth Factor-AA) ELISA Kit
E-EL-MK2103	Monkey DAT (Dopamine Transporter) ELISA Kit
E-EL-MK2104	Monkey FcγR3B (Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit
E-EL-MK2105	Monkey GROα/CXCL1/NAP3 (Growth Regulated Oncogene Alpha) ELISA Kit
E-EL-MK2106	Monkey FSP (Fibroblast Surface Protein) ELISA Kit
E-EL-MK2107	Monkey PLA2G10 (Phospholipase A2, Group X) ELISA Kit
E-EL-MK2108	Monkey COL14 (Collagen Type XIV) ELISA Kit
E-EL-MK2109	Monkey ERK (Extracellular Signal Regulated Kinase) ELISA Kit

Cat.NO	Product Name
E-EL-MK2110	Monkey GSTθ2/GSTt2 (Glutathione S Transferase Theta 2) ELISA Kit
E-EL-MK2111	Monkey GDF11 (Growth Differentiation Factor 11) ELISA Kit
E-EL-MK2112	Monkey Fbn (Fibrinolysin) ELISA Kit
E-EL-MK2113	Monkey GROγ (Growth Regulated Oncogene Gamma) ELISA Kit
E-EL-MK2114	Monkey PIGR/SC (Polymeric Immunoglobulin Receptor/Membrane Secretory Component) ELISA Kit
E-EL-MK2115	Monkey α-Glu (Alpha Glucosidase) ELISA Kit
E-EL-MK2116	Monkey β-CTx (Beta Crosslaps) ELISA Kit
E-EL-MK2117	Monkey FSA (Fetal Sulfoslycoprotein Antigen) ELISA Kit
E-EL-MK2118	Monkey CA242 (Pancreatic Carcinoma Markers-CA242) ELISA Kit
E-EL-MK2119	Monkey PKCβ1 (Protein Kinase C Beta 1) ELISA Kit
E-EL-MK2121	Monkey SOX3 (Sex Determining Region Y Box Protein 3) ELISA Kit
E-EL-MK2122	Monkey MXD1 (MAX Dimerization Protein 1) ELISA Kit
E-EL-MK2123	Monkey STAT2 (Signal Transducer and Activator of Transcription 2) ELISA Kit
E-EL-MK2124	Monkey GTC (General Transcription Complex) ELISA Kit
E-EL-MK2125	Monkey IGFALS (Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) ELISA Kit
E-EL-MK2126	Monkey CKLF (Chemokine Like Factor) ELISA Kit
E-EL-MK2127	Monkey PAL (L-Phenylalanine Ammonia-Lyase) ELISA Kit
E-EL-MK2129	Monkey PF3 (Platelet Factor 3) ELISA Kit
E-EL-MK2130	Monkey CENPI (Centromere Protein I) ELISA Kit
E-EL-MK2131	Monkey VDR (Vitamin D Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-MK2132	Monkey CENPH (Centromere Protein H) ELISA Kit
E-EL-MK2133	Monkey THR α (Thyroid Hormone Receptor Alpha) ELISA Kit
E-EL-MK2134	Monkey LGALS3 (Lectin-Galactoside Binding Soluble 3) ELISA Kit
E-EL-MK2135	Monkey GPRC5A (G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit
E-EL-MK2136	Monkey CR1 (Complement Receptor type 1) ELISA Kit
E-EL-MK2137	Monkey MSLN (Mesothelin) ELISA Kit
E-EL-MK2138	Monkey TIEG1 (TGF Beta Inducible Early Response Gene 1) ELISA Kit
E-EL-MK2139	Monkey PTPRK (Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit
E-EL-MK2140	Monkey GKN1 (Gastrokine 1) ELISA Kit
E-EL-MK2141	Monkey GDF7 (Growth Differentiation Factor 7) ELISA Kit
E-EL-MK2142	Monkey CA125 (Carbohydrate Antigen 125) ELISA Kit
E-EL-MK2143	Monkey CK-19/KRT19 (Cytokeratin 19) ELISA Kit
E-EL-MK2144	Monkey ACTN2 (Actinin Alpha 2) ELISA Kit
E-EL-MK2145	Monkey EFNA1 (Ephrin A1) ELISA Kit
E-EL-MK2146	Monkey VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
E-EL-MK2147	Monkey PLA1 (Phospholipase A1) ELISA Kit
E-EL-MK2148	Monkey PGDS (Prostaglandin D Synthase) ELISA Kit
E-EL-MK2149	Monkey APP (Amyloid Precursor Protein) ELISA Kit
E-EL-MK2150	Monkey FAK (Focal Adhesion Kinase) ELISA Kit
E-EL-MK2151	Monkey BTA (Bladder Tumor Antigen) ELISA Kit

Cat.NO	Product Name
E-EL-MK2152	Monkey EAR3 (Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit
E-EL-MK2153	Monkey BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit
E-EL-MK2155	Monkey FOXP1 (Forkhead Box Protein P1) ELISA Kit
E-EL-MK2156	Monkey NXF1 (Nuclear RNA Export Factor 1) ELISA Kit
E-EL-MK2157	Monkey PEX2 (Peroxisomal Biogenesis Factor 2) ELISA Kit
E-EL-MK2159	Monkey ANG II R1-Ab (Angiotensin II Receptor 1 Antibody) ELISA Kit
E-EL-MK2160	Monkey β -EP (Beta-Endorphin) ELISA Kit
E-EL-MK2161	Monkey TSI (Thyroid Stimulating Immunoglobulin) ELISA Kit
E-EL-MK2162	Monkey PDP (Pyruvate Dehydrogenase Phosphatase) ELISA Kit
E-EL-MK2163	Monkey NT-ProANP (N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit
E-EL-MK2164	Monkey CLCF1 (Cardiotrophin Like Cytokine Factor 1) ELISA Kit
E-EL-MK2165	Monkey CDK9 (Cyclin Dependent Kinase 9) ELISA Kit
E-EL-MK2166	Monkey PHSA (Poly Serum Antigen) ELISA Kit
E-EL-MK2167	Monkey CCRL1 (Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
E-EL-MK2168	Monkey TAFI (Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit
E-EL-MK2169	Monkey PACAPR (Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit
E-EL-MK2171	Monkey CYP1A1 (Cytochrome P450 1A1) ELISA Kit
E-EL-MK2172	Monkey UCP-2 (Uncoupling Protein 2, Mitochondrial) ELISA Kit
E-EL-MK2174	Monkey TR (Thrombin Receptor) ELISA Kit
E-EL-MK2175	Monkey CHAc (Choline acetylase) ELISA Kit

Cat.NO	Product Name
E-EL-MK2176	Monkey CHD5 (Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
E-EL-MK2177	Monkey MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
E-EL-MK2179	Monkey PSAP (Prosaposin) ELISA Kit
E-EL-MK2180	Monkey D2R/DRD2 (Dopamine Receptor D2) ELISA Kit
E-EL-MK2181	Monkey PF4 (Platelet Factor 4) ELISA Kit
E-EL-MK2182	Monkey MUC4 (Mucin 4) ELISA Kit
E-EL-MK2183	Monkey F13B (Coagulation Factor XIII B Polypeptide) ELISA Kit
E-EL-MK2184	Monkey IGF2R (Insulin Like Growth Factor 2 Receptor) ELISA Kit
E-EL-MK2185	Monkey CCL14 (Chemokine C-C-Motif Ligand 14) ELISA Kit
E-EL-MK2188	Monkey PGD2S (Prostaglandin D2 Synthase, Brain) ELISA Kit
E-EL-MK2190	Monkey CSF (Colony Stimulating Factor) ELISA Kit
E-EL-MK2191	Monkey TSA (Thyroid Stimulating Antibody) ELISA Kit
E-EL-MK2192	Monkey RPL26 (Ribosomal Protein L26) ELISA Kit
E-EL-MK2193	Monkey LAB7-2/CD86 (B-Lymphocyte Activation Antigen B7-2) ELISA Kit
E-EL-MK2194	Monkey MOG (Myelin Oligodendrocyte Glycoprotein) ELISA Kit
E-EL-MK2196	Monkey PDIA2 (Protein Disulfide Isomerase A2) ELISA Kit
E-EL-MK2197	Monkey BAK1 (Bcl2 Antagonist/Killer 1) ELISA Kit
E-EL-MK2198	Monkey TSPAN30 (Tetraspanin 30) ELISA Kit
E-EL-MK2199	Monkey CTSL (Cathepsin L) ELISA Kit
E-EL-MK2200	Monkey MUC5B (Mucin 5 Subtype B) ELISA Kit

Cat.NO	Product Name
E-EL-MK2201	Monkey CYP26A1 (Cytochrome P450 26A1) ELISA Kit
E-EL-MK2202	Monkey TSP-1 (Thrombospondin-1) ELISA Kit
E-EL-MK2204	Monkey CGRP (Calcitonin Gene Related Peptide) ELISA Kit
E-EL-MK2205	Monkey AST (Aspartate Aminotransferase) ELISA Kit
E-EL-MK2206	Monkey TAA (Tumor Associated Antigen) ELISA Kit
E-EL-MK2207	Monkey XCR1 (XC-Chemokine Receptor 1) ELISA Kit
E-EL-MK2208	Monkey FN (Fibronectin) ELISA Kit
E-EL-MK2209	Monkey 6-K-PG (6-Keto-prostaglandin) ELISA Kit
E-EL-MK2210	Monkey LTA4H (Leukotriene A4 Hydrolase) ELISA Kit
E-EL-MK2211	Monkey ERK2 (Extracellular Signal Regulated Kinase 2) ELISA Kit
E-EL-MK2212	Monkey MMP-11 (Matrix Metalloproteinase 11) ELISA Kit
E-EL-MK2213	Monkey FGF19 (Fibroblast Growth Factor 19) ELISA Kit
E-EL-MK2214	Monkey CYPD (Cyclophilin D) ELISA Kit
E-EL-MK2215	Monkey PKC θ (Protein Kinase C Theta)ELISA Kit
E-EL-MK2216	Monkey Surv (Survivin) ELISA Kit
E-EL-MK2217	Monkey ApoD (Apolipoprotein D) ELISA Kit
E-EL-MK2218	Monkey P-GP (Permeability Glycoprotein) ELISA Kit
E-EL-MK2222	Monkey ADCY1 (Adenylate Cyclase 1, Brain) ELISA Kit
E-EL-MK2223	Monkey CA1 (Carbonic Anhydrase I) ELISA Kit
E-EL-MK2224	Monkey FLN α (Filamin A, Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-MK2225	Monkey RNASEH (Ribonuclease H)ELISA Kit
E-EL-MK2226	Monkey GUCA2A (Guanylate Cyclase Activator 2A) ELISA Kit
E-EL-MK2227	Monkey FcγBP (Fc Fragment of IgG Binding Protein) ELISA Kit
E-EL-MK2228	Monkey ANG2 (Angiopoietin 2) ELISA Kit
E-EL-MK2229	Monkey PDE4D (Phosphodiesterase 4D, cAMP Specific) ELISA Kit
E-EL-MK2230	Monkey PGF (Prostaglandin F) ELISA Kit
E-EL-MK2231	Monkey α Manase (α Mannosidase) ELISA Kit
E-EL-MK2232	Monkey Ep-CAM/CD326 (Epithelial Cell Adhesion Molecule) ELISA Kit
E-EL-MK2233	Monkey CYFRA21-1 (Cytokeratin Fragment Antigen 21-1) ELISA Kit
E-EL-MK2234	Monkey PRDX1 (Peroxiredoxin 1) ELISA Kit
E-EL-MK2236	Monkey GPX3 (Glutathione Peroxidase 3, Plasma) ELISA Kit
E-EL-MK2237	Monkey ATM (Ataxia Telangiectasia Mutated) ELISA Kit
E-EL-MK2238	Monkey NF-κB p65 (Nuclear Factor-κB p65) ELISA Kit
E-EL-MK2239	Monkey COL8α1 (Collagen Type VIII Alpha 1) ELISA Kit
E-EL-MK2241	Monkey SA (Sialic Acid) ELISA Kit
E-EL-MK2242	Monkey CEA (Carcinoembryonic Antigen) ELISA Kit
E-EL-MK2243	Monkey CSTA (Cystatin A) ELISA Kit
E-EL-MK2244	Monkey GKN2 (Gastrokine 2) ELISA Kit
E-EL-MK2245	Monkey PF4V1 (Platelet Factor 4 Variant 1) ELISA Kit
E-EL-MK2246	Monkey CTXI (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit

Cat.NO	Product Name
E-EL-MK2247	Monkey ANG (Angiostatin) ELISA Kit
E-EL-MK2248	Monkey CLA/CTAGE (Cutaneous Lymphocyte Associated Antigen) ELISA Kit
E-EL-MK2249	Monkey IL1R1 (Interleukin 1 Receptor Type I) ELISA Kit
E-EL-MK2250	Monkey MAF (Macrophage Activating Factor) ELISA Kit
E-EL-MK2251	Monkey REG1α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit
E-EL-MK2252	Monkey TAT (Thrombin/Antithrombin Complex) ELISA Kit
E-EL-MK2253	Monkey FASN (Fatty Acid Synthase) ELISA Kit
E-EL-MK2254	Monkey CTNNβ1 (Catenin, Beta 1) ELISA Kit
E-EL-MK2255	Monkey CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit
E-EL-MK2256	Monkey 6-OHDA (6-Hydroxydopamine) ELISA Kit
E-EL-MK2257	Monkey C/EBPα (CCAAT/Enhancer Binding Protein Alpha) ELISA Kit
E-EL-MK2258	Monkey ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-MK2259	Monkey Ta1 (Thymosin-α1) ELISA Kit
E-EL-MK2260	Monkey PS1 (Presenilin 1) ELISA Kit
E-EL-MK2261	Monkey PACAP-38 (Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
E-EL-MK2262	Monkey PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-MK2263	Monkey IgM (Immunoglobulin M) ELISA Kit
E-EL-MK2264	Monkey INS (Insulin) ELISA Kit
E-EL-MK2265	Monkey GLP-1, active form (Glucagon Like Peptide 1, active form) ELISA Kit
E-EL-MK2266	Monkey CIC(circulating immune complex)ELISA Kit

Cat.NO	Product Name
E-EL-RB0001	Rabbit ALB (Albumin) ELISA Kit
E-EL-RB0002	Rabbit ANP (Atrial Natriuretic Peptide) ELISA Kit
E-EL-RB0003	Rabbit HP (Haptoglobin) ELISA Kit
E-EL-RB0004	Rabbit PLG (Plasminogen) ELISA Kit
E-EL-RB0005	Rabbit CRP (C-reactive protein) ELISA Kit
E-EL-RB0006	Rabbit TRF (Transferrin) ELISA Kit
E-EL-RB0007	Rabbit IgG (Immunoglobulin G) ELISA Kit
E-EL-RB0008	Rabbit MYO (Myoglobin) ELISA Kit
E-EL-RB0009	Rabbit FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
E-EL-RB0010	Rabbit cTn-I/TNNI3 (cardiac Troponin I) ELISA Kit
E-EL-RB0011	Rabbit TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-RB0013	Rabbit IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-RB0014	Rabbit IL-6 (Interleukin 6) ELISA Kit
E-EL-RB0015	Rabbit RETN (Resistin) ELISA Kit
E-EL-RB0016	Rabbit D2D (D-Dimer) ELISA Kit
E-EL-RB0017	Rabbit ACE (Angiotensin I Converting Enzyme) ELISA Kit
E-EL-RB0018	Rabbit F1+2 (Prothrombin Fragment 1+2) ELISA Kit
E-EL-RB0019	Rabbit GSH-Px (Glutathione peroxidase) ELISA Kit
E-EL-RB0021	Rabbit GLP-1 (Glucagon Like Peptide 1) ELISA Kit
E-EL-RB0022	Rabbit ALP (Alkaline Phosphatase) ELISA Kit

Cat.NO	Product Name
E-EL-RB0025	Rabbit β Manase ELISA Kit
E-EL-RB0026	Rabbit PIGR/SC (Polymeric Immunoglobulin Receptor/Membrane Secretory Component) ELISA Kit
E-EL-RB0027	Rabbit plkB α (Phospho-Inhibitory Subunit Of NF Kappa B Alpha) ELISA Kit
E-EL-RB0028	Rabbit UGCG (UDP-glucose ceramide glucosyltransferase) ELISA Kit
E-EL-RB0029	Rabbit PARP4 (Poly ADP Ribose Polymerase 4) ELISA Kit
E-EL-RB0030	Rabbit FLOT2 (Flotillin 2) ELISA Kit
E-EL-RB0033	Rabbit UBE2C (Ubiquitin Conjugating Enzyme E2C) ELISA Kit
E-EL-RB0034	Rabbit β -Glu (Beta Glucosidase) ELISA Kit
E-EL-RB0035	Rabbit DSE (Dermatan Sulfate Epimerase) ELISA Kit
E-EL-RB0036	Rabbit EGF (Epidermal Growth Factor) ELISA Kit
E-EL-RB0037	Rabbit EMILIN 1 (Elastin Microfibril Interface Located Protein 1) ELISA Kit
E-EL-RB0039	Rabbit NUP107 (Nucleoporin 107kDa) ELISA Kit
E-EL-RB0040	Rabbit MAOA (Type A Monoamine Oxidase) ELISA Kit
E-EL-RB0041	Rabbit MMP-12 (Matrix Metalloproteinase 12) ELISA Kit
E-EL-RB0043	Rabbit S100B (S100 Calcium Binding Protein B) ELISA Kit
E-EL-RB0044	Rabbit EPCR (Endothelial Protein C Receptor) ELISA Kit
E-EL-RB0045	Rabbit F5 (Coagulation Factor V) ELISA Kit
E-EL-RB0046	Rabbit UFC (Urinary Free Cortisol) ELISA Kit
E-EL-RB0047	Rabbit C8b (Complement Component 8b) ELISA Kit
E-EL-RB0049	Rabbit CGRP (Calcitonin Gene Related Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-RB0050	Rabbit LUM (Lumican) ELISA Kit
E-EL-RB0051	Rabbit MRP (Multidrug Resistance-Associated Protein) ELISA Kit
E-EL-RB0052	Rabbit Gp49a (Glycoprotein 49A) ELISA Kit
E-EL-RB0053	Rabbit BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit
E-EL-RB0054	Rabbit TWEAK (Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit
E-EL-RB0055	Rabbit OSM (Oncostatin M) ELISA Kit
E-EL-RB0056	Rabbit CK (Creatine Kinase) ELISA Kit
E-EL-RB0057	Rabbit CKMB (Creatine Kinase MB Isoenzyme) ELISA Kit
E-EL-RB0059	Rabbit IL1R1 (Interleukin 1 Receptor Type I) ELISA Kit
E-EL-RB0060	Rabbit NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
E-EL-RB0061	Rabbit PP (Pancreatic Polypeptide) ELISA Kit
E-EL-RB0062	Rabbit MUC6 (Mucin 6) ELISA Kit
E-EL-RB0063	Rabbit MIgM (Membrane IgM) ELISA Kit
E-EL-RB0064	Rabbit COL8α1 (Collagen Type VIII Alpha 1) ELISA Kit
E-EL-RB0065	Rabbit GSTα4 (Glutathione S Transferase Alpha 4) ELISA Kit
E-EL-RB0066	Rabbit TGM3 (Transglutaminase 3, Epidermal) ELISA Kit
E-EL-RB0067	Rabbit MMP-10 (Matrix Metalloproteinase 10) ELISA Kit
E-EL-RB0068	Rabbit FSP (Fibroblast Surface Protein) ELISA Kit
E-EL-RB0069	Rabbit IL-12 (Interleukin 12) ELISA Kit
E-EL-RB0071	Rabbit PLCγ2 (Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit

Cat.NO	Product Name
E-EL-RB0073	Rabbit SCG2 (Secretogranin II) ELISA Kit
E-EL-RB0075	Rabbit GS (Gelsolin) ELISA Kit
E-EL-RB0076	Rabbit Collagenase II ELISA Kit
E-EL-RB0077	Rabbit BFP (Brain Finger Protein) ELISA Kit
E-EL-RB0078	Rabbit PRSS8 (Protease, Serine 8) ELISA Kit
E-EL-RB0079	Rabbit CAP1 (Adenylate Cyclase Associated Protein 1) ELISA Kit
E-EL-RB0081	Rabbit D β H (Dopamine- β -Hydroxylase) ELISA Kit
E-EL-RB0082	Rabbit Slit2 (Slit Homolog 2) ELISA Kit
E-EL-RB0083	Rabbit UGT2B4 (Uridine Diphosphate Glucuronosyltransferase 2 Family Polypeptide B4) ELISA Kit
E-EL-RB0084	Rabbit GP (Anti-Glycoprotein Antibody) ELISA Kit
E-EL-RB0085	Rabbit CNX (Calnexin) ELISA Kit
E-EL-RB0086	Rabbit CA153 (Carbohydrate Antigen 153) ELISA Kit
E-EL-RB0087	Rabbit MUSK (Muscle Skeletal Receptor Tyrosine Kinase) ELISA Kit
E-EL-RB0089	Rabbit BACE2 (β -site APP Cleaving Enzyme 2) ELISA Kit
E-EL-RB0090	Rabbit MTX (Methotrexate) ELISA Kit
E-EL-RB0091	Rabbit TE (Telomerase) ELISA Kit
E-EL-RB0092	Rabbit 5-NT (5-Nucleotidase) ELISA Kit
E-EL-RB0094	Rabbit TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit
E-EL-RB0095	Rabbit AChE (Acetylcholinesterase) ELISA Kit
E-EL-RB0096	Rabbit TF (Transfer Factor) ELISA Kit

Cat.NO	Product Name
E-EL-RB0097	Rabbit RNASET2 (Ribonuclease T2) ELISA Kit
E-EL-RB0098	Rabbit CDSN (Corneodesmosin) ELISA Kit
E-EL-RB0100	Rabbit CNTF (Ciliary Neurotrophic Factor) ELISA Kit
E-EL-RB0101	Rabbit GDNF (Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
E-EL-RB0102	Rabbit NUP35 (Nucleoporin 35kDa) ELISA Kit
E-EL-RB0105	Rabbit USP8 (Ubiquitin Specific Peptidase 8) ELISA Kit
E-EL-RB0106	Rabbit AIF1 (Allograft Inflammatory Factor 1) ELISA Kit
E-EL-RB0107	Rabbit GR (Glutathione Reductase) ELISA Kit
E-EL-RB0108	Rabbit SERT (Serotonin Transporter) ELISA Kit
E-EL-RB0109	Rabbit NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
E-EL-RB0110	Rabbit COL3α1 (Collagen Type III Alpha 1) ELISA Kit
E-EL-RB0111	Rabbit MFG-E8 (Milk Fat Globule EGF Factor 8) ELISA Kit
E-EL-RB0112	Rabbit C4 (Complement Component 4) ELISA Kit
E-EL-RB0113	Rabbit PDGFRα (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit
E-EL-RB0114	Rabbit PTX 3/TSG-14 (Pentraxin 3, Long) ELISA Kit
E-EL-RB0115	Rabbit 6-K-PG (6-Keto-prostaglandin) ELISA Kit
E-EL-RB0116	Rabbit PTHrP (Parathyroid Hormone Related Protein) ELISA Kit
E-EL-RB0117	Rabbit CT (Calcitonin) ELISA Kit
E-EL-RB0118	Rabbit FMOD (Fibromodulin) ELISA Kit
E-EL-RB0119	Rabbit TAP (Trypsinogen Activation Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-RB0120	Rabbit SFRP2 (Secreted Frizzled Related Protein 2) ELISA Kit
E-EL-RB0121	Rabbit LGALS3 (Lectin-Galactoside Binding Soluble 3) ELISA Kit
E-EL-RB0122	Rabbit F9 (Coagulation Factor IX) ELISA Kit
E-EL-RB0123	Rabbit β MSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-RB0124	Rabbit CABIN1 (Calcineurin Binding Protein 1) ELISA Kit
E-EL-RB0125	Rabbit MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit
E-EL-RB0126	Rabbit Pax (Paxillin) ELISA Kit
E-EL-RB0127	Rabbit EAR12 (Eosinophil Associated Ribonuclease A family Member 12) ELISA Kit
E-EL-RB0129	Rabbit GDF15 (Growth Differentiation Factor 15) ELISA Kit
E-EL-RB0130	Rabbit PTGIS (Prostaglandin I Synthase) ELISA Kit
E-EL-RB0131	Rabbit α 1-AGP (Alpha-1-Acid Glycoprotein) ELISA Kit
E-EL-RB0132	Rabbit Bcl-3 (B-cell Leukemia/Lymphoma 3) ELISA Kit
E-EL-RB0133	Rabbit PLVAP (Plasmalemma Vesicle Associated Protein) ELISA Kit
E-EL-RB0134	Rabbit PDGFRL (Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
E-EL-RB0135	Rabbit CDX2 (Caudal Type Homeobox Transcription Factor 2) ELISA Kit
E-EL-RB0137	Rabbit PSA/KLK3 (Prostate Specific Antigen) ELISA Kit
E-EL-RB0138	Rabbit PL (Pancreatic Lipase) ELISA Kit
E-EL-RB0139	Rabbit Bcl-10 (B-cell Leukemia/Lymphoma 10) ELISA Kit
E-EL-RB0140	Rabbit β 2-GP1 Ab IgG (β 2-Glycoprotein 1 Antibody IgG) ELISA Kit
E-EL-RB0141	Rabbit FABP6 (Fatty Acid Binding Protein 6, Ileal) ELISA Kit

Cat.NO	Product Name
E-EL-RB0143	Rabbit PPAR- δ (Peroxisome Proliferator Activated Receptor Delta) ELISA Kit
E-EL-RB0144	Rabbit NSE (Neuron-Specific Enolase) ELISA Kit
E-EL-RB0145	Rabbit ASK-1 (Apoptosis Signal Regulating Kinase 1) ELISA Kit
E-EL-RB0146	Rabbit M-AChR M1 (Muscarinic Acetylcholine Receptor M1) ELISA Kit
E-EL-RB0147	Rabbit MUC2 (Mucin 2) ELISA Kit
E-EL-RB0148	Rabbit CENPF (Centromere Protein F) ELISA Kit
E-EL-RB0149	Rabbit BCAN (Brevican) ELISA Kit
E-EL-RB0150	Rabbit FOXO3 (Forkhead Box Protein O3) ELISA Kit
E-EL-RB0152	Rabbit TSC22 (Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
E-EL-RB0153	Rabbit RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
E-EL-RB0155	Rabbit β 2-GP1 Ab IgA (β 2-Glycoprotein 1 antibody IgA) ELISA Kit
E-EL-RB0156	Rabbit MYD88 (Myeloid Differentiation Factor 88) ELISA Kit
E-EL-RB0157	Rabbit IgG2 (Immunoglobulin G2) ELISA Kit
E-EL-RB0158	Rabbit NRG-2 (Neuregulin 2) ELISA Kit
E-EL-RB0159	Rabbit NUP214 (Nucleoporin 214kDa) ELISA Kit
E-EL-RB0160	Rabbit FGF21(Fibroblast Growth Factor 21) ELISA Kit
E-EL-RB0161	Rabbit ATA (Angiotensinase A) ELISA Kit
E-EL-RB0162	Rabbit FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit
E-EL-RB0163	Rabbit MUC13 (Mucin 13) ELISA Kit
E-EL-RB0164	Rabbit α FP (Alpha-Fetoprotein) ELISA Kit

Cat.NO	Product Name
E-EL-RB0166	Rabbit SAA3 (Serum Amyloid A3) ELISA Kit
E-EL-RB0167	Rabbit TS/TYMS (Thymidylate Synthetase) ELISA Kit
E-EL-RB0169	Rabbit PIICP (Procollagen II C-Terminal ProPeptide) ELISA Kit
E-EL-RB0172	Rabbit MBP/MBL (Mannose Binding Protein/Mannose Binding Lectin) ELISA Kit
E-EL-RB0173	Rabbit GHAb (Anti-Growth Hormone Antibody) ELISA Kit
E-EL-RB0174	Rabbit γ -ECS (γ -Glutamyl Systeine Synthetase) ELISA Kit
E-EL-RB0175	Rabbit eIF3a (Eukaryotic Translation Initiation Factor 3a) ELISA Kit
E-EL-RB0176	Rabbit COL6 (Collagen Type VI) ELISA Kit
E-EL-RB0177	Rabbit PADI4 (Peptidyl Arginine Deiminase Type IV) ELISA Kit
E-EL-RB0178	Rabbit Neopterin ELISA Kit
E-EL-RB0179	Rabbit CD4 (Cluster of Differentiation 4) ELISA Kit
E-EL-RB0180	Rabbit ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit
E-EL-RB0182	Rabbit TCR (Tcell Receptor) ELISA Kit
E-EL-RB0184	Rabbit CLA/CTAGE (Cutaneous Lymphocyte Associated Antigen) ELISA Kit
E-EL-RB0186	Rabbit DPT (Dermatopontin) ELISA Kit
E-EL-RB0187	Rabbit SPOP (Speckle Type POZ Protein) ELISA Kit
E-EL-RB0188	Rabbit AAA (Anti-Actin Antibody) ELISA Kit
E-EL-RB0189	Rabbit MAP τ (Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
E-EL-RB0190	Rabbit TNNT1 (Troponin T Type 1, Slow Skeletal) ELISA Kit
E-EL-RB0192	Rabbit 16- α OHE-1 (16- α Hydroxyestrogen 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB0193	Rabbit IL-9 (Interleukin 9) ELISA Kit
E-EL-RB0194	Rabbit NRG-4 (Neuregulin 4) ELISA Kit
E-EL-RB0195	Rabbit EMA (Epithelial Membrane Antigen) ELISA Kit
E-EL-RB0196	Rabbit β 2-GP (Beta2 Glycoprotein) ELISA Kit
E-EL-RB0197	Rabbit HYOU1 (Hypoxia Up Regulated 1) ELISA Kit
E-EL-RB0198	Rabbit sTfR1 (Soluble Transferrin Receptor1) ELISA Kit
E-EL-RB0199	Rabbit CFI (Complement Factor I) ELISA Kit
E-EL-RB0201	Rabbit CA19-9 (Carbohydrate antigen19-9) ELISA Kit
E-EL-RB0202	Rabbit DCLK1 (Doublecortin Like Kinase 1) ELISA Kit
E-EL-RB0204	Rabbit EGR1 (Early Growth Response Protein 1) ELISA Kit
E-EL-RB0205	Rabbit CRY α B (Crystallin Alpha B) ELISA Kit
E-EL-RB0206	Rabbit Ntn1 (Netrin 1) ELISA Kit
E-EL-RB0207	Rabbit KAF (Keratinocyte Autocrine Factor) ELISA Kit
E-EL-RB0208	Rabbit CXCR1 (CXC-Chemokine Receptor 1) ELISA Kit
E-EL-RB0209	Rabbit PKC β 1 (Protein Kinase C Beta 1) ELISA Kit
E-EL-RB0210	Rabbit UCP-2 (Uncoupling Protein 2, Mitochondrial) ELISA Kit
E-EL-RB0211	Rabbit NRIP1 (Nuclear Receptor Interacting Protein 1) ELISA Kit
E-EL-RB0212	Rabbit TP5 (Thymopentin) ELISA Kit
E-EL-RB0213	Rabbit STAT4 (Signal Transducer and Activator of Transcription 4) ELISA Kit
E-EL-RB0214	Rabbit UBD (Ubiquitin D) ELISA Kit

Cat.NO	Product Name
E-EL-RB0215	Rabbit UCRP (Ubiquitin Cross Reactive Protein) ELISA Kit
E-EL-RB0216	Rabbit PGD2S (Prostaglandin D2 Synthase, Brain) ELISA Kit
E-EL-RB0217	Rabbit CTNN β 1 (Catenin, Beta 1) ELISA Kit
E-EL-RB0218	Rabbit MIP-5 (Macrophage Inflammatory Protein 5) ELISA Kit
E-EL-RB0219	Rabbit SSFA2 (Sperm Specific Antigen 2) ELISA Kit
E-EL-RB0221	Rabbit TSPAN30 (Tetraspanin 30) ELISA Kit
E-EL-RB0222	Rabbit DPD (Deoxypyridinoline) ELISA Kit
E-EL-RB0223	Rabbit ADAM9 (A Disintegrin And Metalloprotease 9) ELISA Kit
E-EL-RB0224	Rabbit Anti-CIIAb (Anti-Type II Collagen Antibody) ELISA Kit
E-EL-RB0225	Rabbit Flt3L(FMS Like Tyrosine Kinase 3 Ligand) ELISA Kit
E-EL-RB0226	Rabbit COL4 α 1 (Collagen Type IV Alpha 1) ELISA Kit
E-EL-RB0227	Rabbit TNNC2 (Troponin C Type 2) ELISA Kit
E-EL-RB0228	Rabbit ECFR (Eosinophil Chemotactic Factor Receptor) ELISA Kit
E-EL-RB0229	Rabbit MCP/CD46 (Membrane Cofactor Protein) ELISA Kit
E-EL-RB0230	Rabbit ADCY3 (Adenylate Cyclase 3) ELISA Kit
E-EL-RB0231	Rabbit DAO (Diamine Oxidase) ELISA Kit
E-EL-RB0232	Rabbit C3 (Complement Component 3) ELISA Kit
E-EL-RB0233	Rabbit SRP9 (Signal Recognition Particle 9kDa) ELISA Kit
E-EL-RB0234	Rabbit MAPKAPK3 (Mitogen Activated Protein Kinase Activated Protein Kinase 3) ELISA Kit
E-EL-RB0235	Rabbit GlyCAM1 (Glycosylation Dependent Cell Adhesion Molecule 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB0236	Rabbit α HSP (Alpha-Hemoglobin Stabilizing Protein) ELISA Kit
E-EL-RB0237	Rabbit PRSS33 (Protease, Serine 33) ELISA Kit
E-EL-RB0238	Rabbit MMP-25 (Matrix Metalloproteinase 25) ELISA Kit
E-EL-RB0240	Rabbit C3bR (Complement Fragment 3b Receptor) ELISA Kit
E-EL-RB0242	Rabbit ABCG2 (ATP Binding Cassette Transporter G2) ELISA Kit
E-EL-RB0243	Rabbit NPHN (Nephrin) ELISA Kit
E-EL-RB0245	Rabbit OXA (Orexin A) ELISA Kit
E-EL-RB0246	Rabbit AR (Androgen Receptor) ELISA Kit
E-EL-RB0247	Rabbit AQP-2 (Aquaporin 2) ELISA Kit
E-EL-RB0248	Rabbit MAP-2 (Microtubule Associated Protein 2) ELISA Kit
E-EL-RB0249	Rabbit NUP50 (Nucleoporin 50kDa) ELISA Kit
E-EL-RB0250	Rabbit Bcl2L11 (Bcl-2 Like Protein 11) ELISA Kit
E-EL-RB0251	Rabbit MAG Ab (Anti-Myelin Associated Glycoprotein Antibody) ELISA Kit
E-EL-RB0252	Rabbit OBTF1 (Octamer Binding Transcription Factor 1) ELISA Kit
E-EL-RB0253	Rabbit PYY (Peptide YY) ELISA Kit
E-EL-RB0254	Rabbit NT-4 (Neurotrophin 4) ELISA Kit
E-EL-RB0255	Rabbit OCLN (Occludin) ELISA Kit
E-EL-RB0256	Rabbit PCDHyA2 (Protocadherin Gamma A2) ELISA Kit
E-EL-RB0257	Rabbit ACA-IgM (Anti-Cardiolipin Antibody IgM) ELISA Kit
E-EL-RB0258	Rabbit BNP (Brain Natriuretic Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-RB0259	Rabbit PDAP1 (Pdgfa Associated Protein 1) ELISA Kit
E-EL-RB0260	Rabbit REN (Renin) ELISA Kit
E-EL-RB0263	Rabbit GSTpi (Glutathione S Transferases Pi) ELISA Kit
E-EL-RB0264	Rabbit CYTH1 (Cytohesin 1) ELISA Kit
E-EL-RB0265	Rabbit MRP1 (Multidrug Resistance Associated Protein 1) ELISA Kit
E-EL-RB0266	Rabbit AGC (Aggrecan) ELISA Kit
E-EL-RB0267	Rabbit FLN α (Filamin A, Alpha) ELISA Kit
E-EL-RB0268	Rabbit PK (Pyruvate Kinase, Liver And RBC) ELISA Kit
E-EL-RB0269	Rabbit LIF (Leukemia Inhibitory Factor) ELISA Kit
E-EL-RB0270	Rabbit OxLDL (Oxidized Low Density Lipoprotein) ELISA Kit
E-EL-RB0271	Rabbit CLU (Clusterin) ELISA Kit
E-EL-RB0273	Rabbit GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit
E-EL-RB0274	Rabbit ITIH4 (Inter Alpha-Globulin Inhibitor H4) ELISA Kit
E-EL-RB0277	Rabbit C-P (C-Peptide) ELISA Kit
E-EL-RB0278	Rabbit CFD (Complement Factor D) ELISA Kit
E-EL-RB0279	Rabbit PHLDA2 (Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit
E-EL-RB0280	Rabbit TUB α 3C (Tubulin Alpha 3c)ELISA Kit
E-EL-RB0282	Rabbit ACTg1 (Actin Gamma 1) ELISA Kit
E-EL-RB0283	Rabbit TSGF (Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
E-EL-RB0284	Rabbit PRDX5 (Peroxisredoxin 5) ELISA Kit

Cat.NO	Product Name
E-EL-RB0285	Rabbit CSTA (Cystatin A) ELISA Kit
E-EL-RB0286	Rabbit CYPD (Cyclophilin D) ELISA Kit
E-EL-RB0287	Rabbit APSA (Anti-Phosphatidyl Serine Antibody) ELISA Kit
E-EL-RB0288	Rabbit TREM-1 (Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
E-EL-RB0289	Rabbit SOX2 (Sex Determining Region Y Box Protein 2) ELISA Kit
E-EL-RB0290	Rabbit BLNK (B-Cell Linker Protein) ELISA Kit
E-EL-RB0291	Rabbit GUK1 (Guanylate Kinase 1) ELISA Kit
E-EL-RB0292	Rabbit ADIPOR1 (Adiponectin Receptor 1) ELISA Kit
E-EL-RB0293	Rabbit CRHBP (Corticotropin Releasing Hormone Binding Protein) ELISA Kit
E-EL-RB0294	Rabbit CCND2 (Cyclin-D2) ELISA Kit
E-EL-RB0295	Rabbit MAP2K1 (Mitogen Activated Protein Kinase Kinase 1) ELISA Kit
E-EL-RB0296	Rabbit F12 (Coagulation Factor XII) ELISA Kit
E-EL-RB0297	Rabbit GAD (Glutamic Acid Decarboxylase) ELISA Kit
E-EL-RB0298	Rabbit CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit
E-EL-RB0299	Rabbit NADPH (Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
E-EL-RB0301	Rabbit TR (Thrombin Receptor) ELISA Kit
E-EL-RB0302	Rabbit Ly6a (Lymphocyte Antigen 6 Complex Locus A) ELISA Kit
E-EL-RB0303	Rabbit SDF-4 (Stromal Cell Derived Factor 4) ELISA Kit
E-EL-RB0304	Rabbit PLGA (Plasminogen Activator) ELISA Kit
E-EL-RB0305	Rabbit PK2 (Prokineticin 2) ELISA Kit

Cat.NO	Product Name
E-EL-RB0306	Rabbit ACO1 (Aconitase 1) ELISA Kit
E-EL-RB0307	Rabbit TAGLN (Transgelin) ELISA Kit
E-EL-RB0308	Rabbit Smad3/MADH3 (Mothers Against Decapentaplegic Homolog 3) ELISA Kit
E-EL-RB0309	Rabbit C4a (Complement Component 4a) ELISA Kit
E-EL-RB0310	Rabbit MYL9 (Myosin Light Chain 9, Regulatory) ELISA Kit
E-EL-RB0311	Rabbit OCT (Ornithine Carbamoyl Transferase) ELISA Kit
E-EL-RB0313	Rabbit PROS (Protein S) ELISA Kit
E-EL-RB0315	Rabbit BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit
E-EL-RB0316	Rabbit VEGF-B (Vascular Endothelial Cell Growth Factor B) ELISA Kit
E-EL-RB0318	Rabbit THRSP (Thyroid Hormone Responsive) ELISA Kit
E-EL-RB0319	Rabbit NXF1 (Nuclear RNA Export Factor 1) ELISA Kit
E-EL-RB0320	Rabbit FE (Ferritin) ELISA Kit
E-EL-RB0321	Rabbit MEC/CCL28 (Mucosae Associated Epithelia Chemokine) ELISA Kit
E-EL-RB0322	Rabbit OTF2B (Octamer Transcription Factor 2B) ELISA Kit
E-EL-RB0323	Rabbit PKD2 (Protein Kinase D2)ELISA Kit
E-EL-RB0325	Rabbit HNRPA2B1 (Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit
E-EL-RB0326	Rabbit GP5 (Glycoprotein V, Platelet) ELISA Kit
E-EL-RB0327	Rabbit LATS (Long Acting Thyroid Stimulator) ELISA Kit
E-EL-RB0328	Rabbit GTEAGE (Glomerular Tissue Extract Advanced Glycosylation End Products) ELISA Kit
E-EL-RB0329	Rabbit UBE1L (Ubiquitin Activating Enzyme E1 Like Protein) ELISA Kit

Cat.NO	Product Name
E-EL-RB0330	Rabbit ANG II R2/AGTR2 (Angiotensin II Receptor 2) ELISA Kit
E-EL-RB0331	Rabbit FPB (Fibrinopeptide B) ELISA Kit
E-EL-RB0332	Rabbit β -LPH (β -Lipotropic Hormone) ELISA Kit
E-EL-RB0333	Rabbit cPLA2 (Cytosolic Phospholipase A2) ELISA Kit
E-EL-RB0334	Rabbit CRT (Calreticulin) ELISA Kit
E-EL-RB0336	Rabbit MANF (Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit
E-EL-RB0337	Rabbit IL-19 (Interleukin 19) ELISA Kit
E-EL-RB0339	Rabbit SAA2 (Serum Amyloid A2) ELISA Kit
E-EL-RB0340	Rabbit NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit
E-EL-RB0341	Rabbit PC (Plasma Anticoagulant Protein C) ELISA Kit
E-EL-RB0342	Rabbit SHH (Hedgehog Homolog, Sonic) ELISA Kit
E-EL-RB0343	Rabbit PRSS1 (Protease, Serine, 1) ELISA Kit
E-EL-RB0344	Rabbit OIT3 (Oncoprotein Induced Transcript 3) ELISA Kit
E-EL-RB0345	Rabbit CNTNAP2 (Contactin Associated Protein Like Protein 2) ELISA Kit
E-EL-RB0346	Rabbit FTH (Ferritin, Heavy Polypeptide) ELISA Kit
E-EL-RB0348	Rabbit LRP6 (Low Density Lipoprotein Receptor Related Protein) ELISA Kit
E-EL-RB0349	Rabbit MMP-1 (Matrix Metalloproteinase 1) ELISA Kit
E-EL-RB0350	Rabbit PP (Placental Protein) ELISA Kit
E-EL-RB0351	Rabbit SOX1 (Sex Determining Region Y Box Protein 1) ELISA Kit
E-EL-RB0352	Rabbit OVA sIgG (Ovalbumin Specific IgG) ELISA Kit

Cat.NO	Product Name
E-EL-RB0353	Rabbit LTD4 (Leukotriene D4) ELISA Kit
E-EL-RB0355	Rabbit rT3 (Reverse Triiodothyronine) ELISA Kit
E-EL-RB0356	Rabbit CLS (Cardiolipin Synthase) ELISA Kit
E-EL-RB0357	Rabbit PROZ (Protein Z) ELISA Kit
E-EL-RB0358	Rabbit PP13 (Placental Protein13) ELISA Kit
E-EL-RB0359	Rabbit SIGLEC2 (Sialic Acid Binding Ig Like Lectin 2) ELISA Kit
E-EL-RB0361	Rabbit EAR2 (Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit
E-EL-RB0362	Rabbit AQP-1 (Aquaporin 1) ELISA Kit
E-EL-RB0363	Rabbit PDXK (Pyridoxal Kinase) ELISA Kit
E-EL-RB0364	Rabbit CA4 (Carbonic Anhydrase IV) ELISA Kit
E-EL-RB0365	Rabbit GTC (General Transcription Complex) ELISA Kit
E-EL-RB0366	Rabbit OTF2A (Octamer Transcription Factor 2A) ELISA Kit
E-EL-RB0367	Rabbit IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) ELISA Kit
E-EL-RB0368	Rabbit CYFRA21-1 (Cytokeratin Fragment Antigen 21-1) ELISA Kit
E-EL-RB0369	Rabbit RAR α (Retinoic Acid Receptor Alpha) ELISA Kit
E-EL-RB0371	Rabbit PS (Plasma Anticoagulant Protein S) ELISA Kit
E-EL-RB0372	Rabbit α Manase (α Mannosidase) ELISA Kit
E-EL-RB0373	Rabbit TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit
E-EL-RB0374	Rabbit CDKN3 (Cyclin Dependent Kinase Inhibitor 3) ELISA Kit
E-EL-RB0375	Rabbit TF (Tissue Factor) ELISA Kit

Cat.NO	Product Name
E-EL-RB0376	Rabbit BMP-15 (Bone Morphogenetic Protein 15) ELISA Kit
E-EL-RB0377	Rabbit FGF23 (Fibroblast Growth Factor 23) ELISA Kit
E-EL-RB0378	Rabbit Sp100 (Sp100 Nuclear Antigen) ELISA Kit
E-EL-RB0379	Rabbit ANG2 (Angiopoietin 2) ELISA Kit
E-EL-RB0382	Rabbit UBR4 (Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit
E-EL-RB0383	Rabbit CNN1 (Calponin 1, Basic, Smooth Muscle) ELISA Kit
E-EL-RB0384	Rabbit F-TESTO (Free Testosterone) ELISA Kit
E-EL-RB0385	Rabbit TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit
E-EL-RB0386	Rabbit NPYR2 (Neuropeptide Y Receptor Y2) ELISA Kit
E-EL-RB0387	Rabbit A β 1-42 (Amyloid Beta Peptide 1-42) ELISA Kit
E-EL-RB0388	Rabbit CNTNAP4 (Contactin Associated Protein 4) ELISA Kit
E-EL-RB0390	Rabbit PGAM2 (Phosphoglycerate Mutase 2) ELISA Kit
E-EL-RB0391	Rabbit VCAM-1/CD106 (soluble Vascular Cell Adhesion Molecule 1) ELISA Kit
E-EL-RB0392	Rabbit CENPI (Centromere Protein I) ELISA Kit
E-EL-RB0393	Rabbit NT-ProANP (N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit
E-EL-RB0394	Rabbit β -actin (Actin Beta) ELISA Kit
E-EL-RB0395	Rabbit NNE (Non-Neuronal Enolase) ELISA Kit
E-EL-RB0396	Rabbit FUCA (Alpha-L-Fucosidase, Tissue) ELISA Kit
E-EL-RB0397	Rabbit CCL14 (Chemokine C-C-Motif Ligand 14) ELISA Kit
E-EL-RB0398	Rabbit SPT β N4 (Spectrin Beta, Non Erythrocytic 4) ELISA Kit

Cat.NO	Product Name
E-EL-RB0399	Rabbit SDC4 (Syndecan 4) ELISA Kit
E-EL-RB0400	Rabbit PCX (Podocalyxin) ELISA Kit
E-EL-RB0401	Rabbit HB μ (Hemoglobin, Mu) ELISA Kit
E-EL-RB0403	Rabbit HbH (Hemoglobin H) ELISA Kit
E-EL-RB0406	Rabbit PSPN (Persephin) ELISA Kit
E-EL-RB0407	Rabbit AP36 (Apelin 36) ELISA Kit
E-EL-RB0408	Rabbit CCSA-2 (Colon Cancer Specific Antigen-2) ELISA Kit
E-EL-RB0409	Rabbit P-GP (Permeability Glycoprotein) ELISA Kit
E-EL-RB0410	Rabbit ApoB (Apolipoprotein B) ELISA Kit
E-EL-RB0411	Rabbit CNTN1 (Contactin 1) ELISA Kit
E-EL-RB0412	Rabbit EK (Eepidermal Keratin) ELISA Kit
E-EL-RB0413	Rabbit NAGLU (N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit
E-EL-RB0414	Rabbit TACR1 (Tachykinin Receptor 1) ELISA Kit
E-EL-RB0415	Rabbit ANG (Angiostatin) ELISA Kit
E-EL-RB0416	Rabbit BLVRB (Biliverdin Reductase B) ELISA Kit
E-EL-RB0417	Rabbit CHD3 (Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
E-EL-RB0418	Rabbit PCP (Procollagen C Peptidase)ELISA Kit
E-EL-RB0419	Rabbit STMN1 (Stathmin 1) ELISA Kit
E-EL-RB0421	Rabbit NT (Neurotensin) ELISA Kit
E-EL-RB0423	Rabbit β -TG (β -Thromboglobulin) ELISA Kit

Cat.NO	Product Name
E-EL-RB0424	Rabbit β TG/PBP/CXCL7/NAP2 (Thromboglobulin, Beta) ELISA Kit
E-EL-RB0425	Rabbit DAT (Dopamine Transporter) ELISA Kit
E-EL-RB0427	Rabbit CASP7 (Caspase 7) ELISA Kit
E-EL-RB0428	Rabbit PAQR7 (Progesterin and AdipoQ Receptor Family Member VII) ELISA Kit
E-EL-RB0429	Rabbit ECM1 (Extracellular Matrix Protein 1) ELISA Kit
E-EL-RB0430	Rabbit FBN1 (Fibrillin 1) ELISA Kit
E-EL-RB0431	Rabbit PRB1 (Basic Salivary Proline Rich Protein 1) ELISA Kit
E-EL-RB0432	Rabbit SYP (Synaptophysin) ELISA Kit
E-EL-RB0433	Rabbit PDGFC (Platelet Derived Growth Factor C) ELISA Kit
E-EL-RB0434	Rabbit GasR (Gastrin Receptor) ELISA Kit
E-EL-RB0436	Rabbit D4R/DRD4 (Dopamine Receptor D4) ELISA Kit
E-EL-RB0437	Rabbit α 1BG (Alpha-1-B-Glycoprotein) ELISA Kit
E-EL-RB0438	Rabbit FOXO1 (Forkhead Box Protein O1) ELISA Kit
E-EL-RB0439	Rabbit FGL2 (Fibrinogen Like Protein 2) ELISA Kit
E-EL-RB0441	Rabbit MDC/CCL22 (Macrophage Derived Chemokine) ELISA Kit
E-EL-RB0442	Rabbit CNP (C-type Natriuretic Peptide) ELISA Kit
E-EL-RB0443	Rabbit LFA-1 (Lymphocyte Function Associated Antigen 1) ELISA Kit
E-EL-RB0444	Rabbit CLA (Collagen Autoantibody) ELISA Kit
E-EL-RB0445	Rabbit SOX9 (Sex Determining Region Y Box Protein 9) ELISA Kit
E-EL-RB0446	Rabbit CCL4L1 (Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB0447	Rabbit IL-12 p40 (Interleukin 12 p40) ELISA Kit
E-EL-RB0448	Rabbit LRG1 (Leucine Rich Alpha-2-Glycoprotein 1) ELISA Kit
E-EL-RB0449	Rabbit SNCy (Synuclein Gamma) ELISA Kit
E-EL-RB0450	Rabbit OPGL (Osteoprotegerin Ligand) ELISA Kit
E-EL-RB0451	Rabbit Cyr61 (Cysteine Rich Protein, Angiogenic Inducer 61) ELISA Kit
E-EL-RB0452	Rabbit F13 (Coagulation Factor X III) ELISA Kit
E-EL-RB0453	Rabbit MlgD (Membrane IgD) ELISA Kit
E-EL-RB0454	Rabbit ACTα2 (Actin Alpha 2, Smooth Muscle) ELISA Kit
E-EL-RB0456	Rabbit M-AChR M2 (Muscarinic Acetylcholine Receptor M2) ELISA Kit
E-EL-RB0457	Rabbit CNTN4 (Contactin 4) ELISA Kit
E-EL-RB0458	Rabbit PTF1α (Pancreas Specific Transcription Factor 1α) ELISA Kit
E-EL-RB0459	Rabbit MXD1 (MAX Dimerization Protein 1) ELISA Kit
E-EL-RB0460	Rabbit GLP-2 (Glucagon Like Peptide 2) ELISA Kit
E-EL-RB0461	Rabbit PGC (Pepsinogen C) ELISA Kit
E-EL-RB0462	Rabbit SIGLEC1 (Sialic Acid Binding Ig Like Lectin 1) ELISA Kit
E-EL-RB0463	Rabbit sFMC (Soluble Fibrin Monomer Complex) ELISA Kit
E-EL-RB0464	Rabbit BLC (B-Lymphocyte Chemoattractant) ELISA Kit
E-EL-RB0465	Rabbit FLNC (Filamin C, Gamma) ELISA Kit
E-EL-RB0466	Rabbit CKB (Creatine Kinase, Brain) ELISA Kit
E-EL-RB0467	Rabbit SIGLEC12 (Sialic Acid Binding Ig Like Lectin 12) ELISA Kit

Cat.NO	Product Name
E-EL-RB0468	Rabbit RPS6K β 1 (Ribosomal Protein S6 Kinase Beta 1) ELISA Kit
E-EL-RB0469	Rabbit GNLY (Granulysin) ELISA Kit
E-EL-RB0470	Rabbit PKC (Protein Kinase C) ELISA Kit
E-EL-RB0471	Rabbit CK-13/KRT13 (Cytokeratin 13) ELISA Kit
E-EL-RB0472	Rabbit PTD (Pentosidine) ELISA Kit
E-EL-RB0473	Rabbit CDK7 (Cyclin Dependent Kinase 7) ELISA Kit
E-EL-RB0474	Rabbit VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit
E-EL-RB0475	Rabbit Thyronine ELISA Kit
E-EL-RB0476	Rabbit CRLF1 (Cytokine Receptor Like Factor 1) ELISA Kit
E-EL-RB0477	Rabbit COL3 (Collagen Type III) ELISA Kit
E-EL-RB0478	Rabbit APP (Amyloid Precursor Protein) ELISA Kit
E-EL-RB0479	Rabbit FBN3 (Fibrillin 3) ELISA Kit
E-EL-RB0480	Rabbit ADH (Antidiuretic Hormone) ELISA Kit
E-EL-RB0481	Rabbit VIP (Vasoactive Intestinal Peptide) ELISA Kit
E-EL-RB0482	Rabbit IGFBP-7 (Insulin Like Growth Factor Binding Protein 7) ELISA Kit
E-EL-RB0483	Rabbit CHD5 (Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
E-EL-RB0484	Rabbit PCDH β 16 (Protocadherin Beta 16) ELISA Kit
E-EL-RB0485	Rabbit CD30 (Cluster of Differentiation 30) ELISA Kit
E-EL-RB0486	Rabbit RNF4 (Ring Finger Protein 4) ELISA Kit
E-EL-RB0487	Rabbit IL-10 (Interleukin-10) ELISA Kit

Cat.NO	Product Name
E-EL-RB0488	Rabbit PUMA (P53 Upregulated Modulator of Apoptosis) ELISA Kit
E-EL-RB0489	Rabbit FPR2 (Formyl Peptide Receptor 2) ELISA Kit
E-EL-RB0493	Rabbit BMP-3 (Bone Morphogenetic Protein 3) ELISA Kit
E-EL-RB0494	Rabbit GTF2A1 (General Transcription Factor IIA, Polypeptide 1) ELISA Kit
E-EL-RB0496	Rabbit CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
E-EL-RB0498	Rabbit AMY2 (Amylase Alpha 2, Pancreatic) ELISA Kit
E-EL-RB0499	Rabbit CTXI (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
E-EL-RB0500	Rabbit IGFALS (Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) ELISA Kit
E-EL-RB0501	Rabbit EP2 (Prostaglandin E Receptor 2) ELISA Kit
E-EL-RB0503	Rabbit FSHR (Follicle Stimulating Hormone Receptor) ELISA Kit
E-EL-RB0504	Rabbit SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit
E-EL-RB0505	Rabbit α -Glu (Alpha Glucosidase) ELISA Kit
E-EL-RB0506	Rabbit ELN (Elastin) ELISA Kit
E-EL-RB0507	Rabbit SPA (Pulmonary Surfactant Associated Protein A) ELISA Kit
E-EL-RB0508	Rabbit GKN2 (Gastrokine 2) ELISA Kit
E-EL-RB0509	Rabbit Surv (Survivin) ELISA Kit
E-EL-RB0510	Rabbit I κ B β (Inhibitory Subunit of NF Kappa B Beta) ELISA Kit
E-EL-RB0511	Rabbit PKR2 (Prokineticin Receptor 2) ELISA Kit
E-EL-RB0512	Rabbit CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) ELISA Kit
E-EL-RB0513	Rabbit sCKR (Soluble Cytokine Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-RB0514	Rabbit HSF2 (Heat Shock Transcription Factor 2) ELISA Kit
E-EL-RB0515	Rabbit β -CTx (Beta Crosslaps) ELISA Kit
E-EL-RB0516	Rabbit FGF19 (Fibroblast Growth Factor 19) ELISA Kit
E-EL-RB0517	Rabbit PF4V1 (Platelet Factor 4 Variant 1) ELISA Kit
E-EL-RB0518	Rabbit DKK1 (Dickkopf Related Protein 1) ELISA Kit
E-EL-RB0519	Rabbit SLC/CCL21 (Secondary Lymphoid Tissue Chemokine) ELISA Kit
E-EL-RB0520	Rabbit CTSC (Cathepsin C) ELISA Kit
E-EL-RB0523	Rabbit tPA (Plasminogen Activator, Tissue) ELISA Kit
E-EL-RB0524	Rabbit G6PD (Glucose 6 Phosphate Dehydrogenase) ELISA Kit
E-EL-RB0525	Rabbit NGB (Neuroglobin) ELISA Kit
E-EL-RB0526	Rabbit ATF-3 (Activating Transcription Factor 3) ELISA Kit
E-EL-RB0527	Rabbit ApoH (Apolipoprotein H) ELISA Kit
E-EL-RB0528	Rabbit ANG-R-Tie2 (Angiopoietin Receptor Tie2) ELISA Kit
E-EL-RB0529	Rabbit CXCR3 (Chemokine CXC Receptor 3) ELISA Kit
E-EL-RB0530	Rabbit BlysR (B Lymphocyte Stimulator Receptor) ELISA Kit
E-EL-RB0531	Rabbit BMP-1 (Bone Morphogenetic Protein 1) ELISA Kit
E-EL-RB0532	Rabbit CG α (Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit
E-EL-RB0534	Rabbit BACE1 (Beta-site APP Cleaving Enzyme 1) ELISA Kit
E-EL-RB0535	Rabbit Flt3 (FMS Like Tyrosine Kinase 3) ELISA Kit
E-EL-RB0536	Rabbit CG (Chorionic Gonadotropin) ELISA Kit

Cat.NO	Product Name
E-EL-RB0537	Rabbit CMLC-1 (Cardiac Myosin Light Chains 1) ELISA Kit
E-EL-RB0538	Rabbit STAT6 (Signal Transducer and Activator of Transcription 6) ELISA Kit
E-EL-RB0540	Rabbit PKN1 (Protein Kinase N1) ELISA Kit
E-EL-RB0541	Rabbit TLA (Thymus Leukemia Antigen) ELISA Kit
E-EL-RB0542	Rabbit TLN (Talin) ELISA Kit
E-EL-RB0543	Rabbit CA125 (Carbohydrate Antigen 125) ELISA Kit
E-EL-RB0544	Rabbit AMA IgA (Anti-Myelin antibody IgA) ELISA Kit
E-EL-RB0545	Rabbit GALR1 (Galanin Receptor 1) ELISA Kit
E-EL-RB0547	Rabbit AMBP (Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit
E-EL-RB0548	Rabbit PDH E1 (Pyruvate Dehydrogenase E1) ELISA Kit
E-EL-RB0549	Rabbit COLEC11 (Collectin Sub-Family Member 11) ELISA Kit
E-EL-RB0550	Rabbit NE/ELA2 (Elastase 2, Neutrophil) ELISA Kit
E-EL-RB0551	Rabbit TAP (Transporter Associated with Antigen Processing) ELISA Kit
E-EL-RB0552	Rabbit NP (Nandrolone Phenylpropionate) ELISA Kit
E-EL-RB0553	Rabbit Grb2 (Growth Factor Receptor Bound Protein 2) ELISA Kit
E-EL-RB0554	Rabbit CR (Calretinin) ELISA Kit
E-EL-RB0555	Rabbit DPP4 (Dipeptidyl Peptidase IV) ELISA Kit
E-EL-RB0556	Rabbit LDL (Low Density Lipoprotein) ELISA Kit
E-EL-RB0557	Rabbit ABCG1 (ATP Binding Cassette Transporter G1) ELISA Kit
E-EL-RB0558	Rabbit PGF (Placental Growth Factor) ELISA Kit

Cat.NO	Product Name
E-EL-RB0559	Rabbit c-myc (c-myc Oncogene Product) ELISA Kit
E-EL-RB0561	Rabbit MVP (Major Vault Protein) ELISA Kit
E-EL-RB0562	Rabbit GR β (Glucocorticoid Receptor Beta) ELISA Kit
E-EL-RB0565	Rabbit ApoM (Apolipoprotein M) ELISA Kit
E-EL-RB0566	Rabbit LFA-3 (Lymphocyte Function Associated Antigen 3) ELISA Kit
E-EL-RB0567	Rabbit C2 (Complement Component 2) ELISA Kit
E-EL-RB0568	Rabbit CA1 (Carbonic Anhydrase I) ELISA Kit
E-EL-RB0569	Rabbit PLA2G10 (Phospholipase A2, Group X) ELISA Kit
E-EL-RB0570	Rabbit ADM (Adrenomedullin) ELISA Kit
E-EL-RB0571	Rabbit MTHFR (Methylenetetrahydrofolate Reductase) ELISA Kit
E-EL-RB0572	Rabbit ARO (Aromatase) ELISA Kit
E-EL-RB0573	Rabbit GAS (N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
E-EL-RB0575	Rabbit GST α 1 (Glutathione S Transferase Alpha 1) ELISA Kit
E-EL-RB0576	Rabbit ME20M (Melanoma Associated ME20) ELISA Kit
E-EL-RB0577	Rabbit MUC1 (Mucin 1) ELISA Kit
E-EL-RB0578	Rabbit IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
E-EL-RB0579	Rabbit AGGF1 (Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit
E-EL-RB0580	Rabbit STAT1 (Signal Transducer and Activator of Transcription 1) ELISA Kit
E-EL-RB0582	Rabbit β -EPR (Beta-Endorphin Receptor) ELISA Kit
E-EL-RB0583	Rabbit ANG-R-Tie1 (Angiopoietin Receptor Tie1) ELISA Kit

Cat.NO	Product Name
E-EL-RB0584	Rabbit PTPRN (Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
E-EL-RB0586	Rabbit CA9 (Carbonic Anhydrase IX) ELISA Kit
E-EL-RB0587	Rabbit FGFR3 (Fibroblast Growth Factor Receptor 3) ELISA Kit
E-EL-RB0588	Rabbit TXLNa (Taxilin Alpha) ELISA Kit
E-EL-RB0589	Rabbit CCL16 (Chemokine C-C-Motif Ligand 16) ELISA Kit
E-EL-RB0590	Rabbit Extein ELISA Kit
E-EL-RB0592	Rabbit PDGFR β (Platelet Derived Growth Factor Receptor Beta) ELISA Kit
E-EL-RB0593	Rabbit TACR2 (Tachykinin Receptor 2) ELISA Kit
E-EL-RB0594	Rabbit CASP11 (Caspase 11) ELISA Kit
E-EL-RB0595	Rabbit AGR2 (Anterior Gradient Protein 2) ELISA Kit
E-EL-RB0597	Rabbit BTA (Bladder Tumor Antigen) ELISA Kit
E-EL-RB0598	Rabbit AMH (Anti-Mullerian Hormone) ELISA Kit
E-EL-RB0600	Rabbit TCGF- III (Tcell Growth Factor III) ELISA Kit
E-EL-RB0602	Rabbit Try2 (Trypsinogen II) ELISA Kit
E-EL-RB0603	Rabbit DAP6 (Death Associated Protein 6) ELISA Kit
E-EL-RB0605	Rabbit EGR3 (Early Growth Response Protein 3) ELISA Kit
E-EL-RB0606	Rabbit MMP-3 (Matrix Metalloproteinase 3) ELISA Kit
E-EL-RB0607	Rabbit CLDN4 (Claudin 4) ELISA Kit
E-EL-RB0608	Rabbit C3b (Complement Fragment 3b) ELISA Kit
E-EL-RB0609	Rabbit CG (Cryoglobulin) ELISA Kit

Cat.NO	Product Name
E-EL-RB0610	Rabbit SNRPD1 (Small Nuclear Ribonucleoprotein Polypeptide D1) ELISA Kit
E-EL-RB0611	Rabbit FAIM3 (Fas Apoptotic Inhibitory Molecule 3) ELISA Kit
E-EL-RB0612	Rabbit IGFBP-5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit
E-EL-RB0613	Rabbit FBG (Fibrinogen) ELISA Kit
E-EL-RB0614	Rabbit TPS (Tissue Polypeptide Specific Antigen) ELISA Kit
E-EL-RB0615	Rabbit TFPI2 (Tissue Factor Pathway Inhibitor 2) ELISA Kit
E-EL-RB0616	Rabbit PDGF-BB (Platelet Derived Growth Factor-BB) ELISA Kit
E-EL-RB0617	Rabbit FRS3 (Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
E-EL-RB0619	Rabbit CYP7A1 (Cytochrome P450 7A1) ELISA Kit
E-EL-RB0622	Rabbit MX1 (Myxovirus Resistance 1) ELISA Kit
E-EL-RB0623	Rabbit S100A8 (S100 Calcium Binding Protein A8/Calgranulin A) ELISA Kit
E-EL-RB0624	Rabbit MSTN (Myostatin) ELISA Kit
E-EL-RB0625	Rabbit Fcγ (Fc Fragment of IgG) ELISA Kit
E-EL-RB0626	Rabbit AVPR1A (Arginine Vasopressin Receptor 1A) ELISA Kit
E-EL-RB0627	Rabbit ADCY4 (Adenylate Cyclase 4) ELISA Kit
E-EL-RB0628	Rabbit PHB (Prohibitin) ELISA Kit
E-EL-RB0629	Rabbit PLAA (Phospholipase A2 Activating Protein) ELISA Kit
E-EL-RB0630	Rabbit CKLF (Chemokine Like Factor) ELISA Kit
E-EL-RB0632	Rabbit CX37 (Connexin 37) ELISA Kit
E-EL-RB0633	Rabbit PP14 (Placental Protein14) ELISA Kit

Cat.NO	Product Name
E-EL-RB0634	Rabbit TEP1 (Telomerase Associated Protein 1) ELISA Kit
E-EL-RB0635	Rabbit GMCSF/CSF2(Granulocyte Macrophage Colony Stimulating Factor 2) ELISA Kit
E-EL-RB0636	Rabbit PLP1 (Proteolipid Protein 1, Myelin) ELISA Kit
E-EL-RB0637	Rabbit GATA4 (GATA Binding Protein 4) ELISA Kit
E-EL-RB0638	Rabbit GDF5 (Growth Differentiation Factor 5) ELISA Kit
E-EL-RB0639	Rabbit ATF-1 (Activating Transcription Factor 1) ELISA Kit
E-EL-RB0640	Rabbit ANG III (Angiotensin III) ELISA Kit
E-EL-RB0641	Rabbit HB β (Hemoglobin Beta) ELISA Kit
E-EL-RB0642	Rabbit Alk-Smase (Alkaline Sphingomyelinase) ELISA Kit
E-EL-RB0643	Rabbit TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
E-EL-RB0644	Rabbit MIF (Macrophage Migration Inhibitory Factor) ELISA Kit
E-EL-RB0645	Rabbit α -SMA (α -Smooth muscle actin) ELISA Kit
E-EL-RB0646	Rabbit Dyn (Big Dynorphin) ELISA Kit
E-EL-RB0647	Rabbit EFNA5 (Ephrin A5) ELISA Kit
E-EL-RB0649	Rabbit DNM2 (Dynamin 2) ELISA Kit
E-EL-RB0651	Rabbit E-Cad (E-Cadherin) ELISA Kit
E-EL-RB0653	Rabbit IgG3 (Immunoglobulin G3) ELISA Kit
E-EL-RB0654	Rabbit TEF (Thyrotroph Embryonic Factor) ELISA Kit
E-EL-RB0655	Rabbit SELE (E-Selectin) ELISA Kit
E-EL-RB0656	Rabbit CASP3 (Caspase 3) ELISA Kit

Cat.NO	Product Name
E-EL-RB0658	Rabbit VEGF (Vascular Endothelial Cell Growth Factor) ELISA Kit
E-EL-RB0659	Rabbit IgE (Immunoglobulin E) ELISA Kit
E-EL-RB0660	Rabbit GAL (Galanin) ELISA Kit
E-EL-RB0661	Rabbit KLK7 (Kallikrein 7) ELISA Kit
E-EL-RB0662	Rabbit LIPC (Lipase, Hepatic) ELISA Kit
E-EL-RB0663	Rabbit α HBDH (α -Hydroxybutyrate Dehydrogenase) ELISA Kit
E-EL-RB0664	Rabbit Fc ϵ R I (Receptor I for the Fc Fragment of IgE) ELISA Kit
E-EL-RB0665	Rabbit SGK2 (Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit
E-EL-RB0666	Rabbit NKB (Neurokinins B) ELISA Kit
E-EL-RB0667	Rabbit RIPK2 (Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit
E-EL-RB0668	Rabbit PDP (Pyruvate Dehydrogenase Phosphatase) ELISA Kit
E-EL-RB0670	Rabbit C3a (Complement Component 3a) ELISA Kit
E-EL-RB0671	Rabbit CALU (Calumenin) ELISA Kit
E-EL-RB0672	Rabbit CTNN α 1 (Catenin, Alpha 1) ELISA Kit
E-EL-RB0673	Rabbit CYP26A1 (Cytochrome P450 26A1) ELISA Kit
E-EL-RB0674	Rabbit PAR4 (Protease Activated Receptor 4) ELISA Kit
E-EL-RB0675	Rabbit Smad1/MADH1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit
E-EL-RB0677	Rabbit CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
E-EL-RB0678	Rabbit CDK8 (Cyclin Dependent Kinase 8) ELISA Kit
E-EL-RB0679	Rabbit IFN- γ (Interferon Gamma) ELISA Kit

Cat.NO	Product Name
E-EL-RB0680	Rabbit CCRL1 (Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
E-EL-RB0681	Rabbit α 1ACT (Alpha1 Antichymotrypsin) ELISA Kit
E-EL-RB0682	Rabbit TSI (Thyroid Stimulating Immunoglobulin) ELISA Kit
E-EL-RB0683	Rabbit SEPCR/PROCR (Soluble Endothelial Protein C Receptor) ELISA Kit
E-EL-RB0684	Rabbit REG4 (Regenerating Islet Derived Protein 4) ELISA Kit
E-EL-RB0687	Rabbit IL-2R α /CD25 (Interleukin 2 Receptor Alpha) ELISA Kit
E-EL-RB0688	Rabbit BSP (Bone Sialoprotein) ELISA Kit
E-EL-RB0690	Rabbit RARS (Arginyl tRNA Synthetase) ELISA Kit
E-EL-RB0691	Rabbit PEX2 (Peroxisomal Biogenesis Factor 2) ELISA Kit
E-EL-RB0692	Rabbit TR (Targeted Ribonuclease) ELISA Kit
E-EL-RB0693	Rabbit OLAB (Anti-Oxidized Low Density Lipoprotein Antibody) ELISA Kit
E-EL-RB0694	Rabbit PON1 (Paraoxonase 1) ELISA Kit
E-EL-RB0695	Rabbit OSF (Osteoclast Stimulating Factor 1) ELISA Kit
E-EL-RB0696	Rabbit PKI γ (Protein Kinase Inhibitor Gamma) ELISA Kit
E-EL-RB0697	Rabbit ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit
E-EL-RB0699	Rabbit ACP3/PAP (Prostatic Acid Phosphatase 3) ELISA Kit
E-EL-RB0700	Rabbit FSTL1 (Follistatin Like Protein 1) ELISA Kit
E-EL-RB0701	Rabbit F11 (Coagulation Factor XI) ELISA Kit
E-EL-RB0702	Rabbit PKC ϵ (Protein Kinase C Epsilon) ELISA Kit
E-EL-RB0703	Rabbit GAD-Ab (Anti-Glutamic Acid Decarboxylase antibody) ELISA Kit

Cat.NO	Product Name
E-EL-RB0704	Rabbit DNM1 (Dynamin 1) ELISA Kit
E-EL-RB0705	Rabbit Cyt-C (Cytochrome C) ELISA Kit
E-EL-RB0706	Rabbit GPLD1 (Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Kit
E-EL-RB0707	Rabbit GDF3 (Growth Differentiation Factor 3) ELISA Kit
E-EL-RB0708	Rabbit BGN (Biglycan) ELISA Kit
E-EL-RB0709	Rabbit PGM (Phosphoglucomutase) ELISA Kit
E-EL-RB0710	Rabbit CBX3 (Chromobox Homolog 3) ELISA Kit
E-EL-RB0711	Rabbit GTF2B (General Transcription Factor IIB) ELISA Kit
E-EL-RB0712	Rabbit NOX1 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
E-EL-RB0713	Rabbit ApoA2 (Apolipoprotein A2) ELISA Kit
E-EL-RB0714	Rabbit SPTAN1 (Alpha-Fodrin) ELISA Kit
E-EL-RB0715	Rabbit DKK4 (Dickkopf Related Protein 4) ELISA Kit
E-EL-RB0716	Rabbit VEGFR1/Flt1 (Vascuoar Endothelial Growth Factor Receptor 1) ELISA Kit
E-EL-RB0717	Rabbit SDF-2 (Stromal Cell Derived Factor 2) ELISA Kit
E-EL-RB0718	Rabbit ANGPTL2 (Angiopoietin Like Protein 2) ELISA Kit
E-EL-RB0719	Rabbit PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit
E-EL-RB0720	Rabbit HCCS (Holocytochrome C Synthase) ELISA Kit
E-EL-RB0721	Rabbit DDT (D-Dopachrome Tautomerase) ELISA Kit
E-EL-RB0722	Rabbit IL-23 (Interleukin 23) ELISA Kit
E-EL-RB0723	Rabbit VDR (Vitamin D Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-RB0724	Rabbit VEGF-D (Vascular Endothelial Growth Factor D) ELISA Kit
E-EL-RB0725	Rabbit PLA1 (Phospholipase A1) ELISA Kit
E-EL-RB0726	Rabbit GEF (Guanylate Exchange Factor) ELISA Kit
E-EL-RB0728	Rabbit CDNF (Cerebral Dopamine Neurotrophic Factor) ELISA Kit
E-EL-RB0729	Rabbit NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) ELISA Kit
E-EL-RB0730	Rabbit Preptin ELISA Kit
E-EL-RB0732	Rabbit FTL (Ferritin, Light Polypeptide) ELISA Kit
E-EL-RB0733	Rabbit BANP (BTG3 Associated Nuclear Protein) ELISA Kit
E-EL-RB0734	Rabbit RNASEN (Ribonuclease III, Nuclear) ELISA Kit
E-EL-RB0735	Rabbit STK39 (Serine/Threonine Kinase 39) ELISA Kit
E-EL-RB0736	Rabbit GAL8 (Galectin 8) ELISA Kit
E-EL-RB0737	Rabbit CLDN14 (Claudin 14) ELISA Kit
E-EL-RB0738	Rabbit DBP (Vitamin D Binding Protein) ELISA Kit
E-EL-RB0739	Rabbit TP63 (Tumor Protein p63) ELISA Kit
E-EL-RB0740	Rabbit GPIIb/IIIa (Platelet Membrane Glycoprotein IIb/IIIa) ELISA Kit
E-EL-RB0741	Rabbit PDCN (Podocin) ELISA Kit
E-EL-RB0742	Rabbit ANGPTL3 (Angiopoietin Like Protein 3) ELISA Kit
E-EL-RB0743	Rabbit Sp17 (Sperm Protein 17) ELISA Kit
E-EL-RB0744	Rabbit SOCS3 (Suppressors of Cytokine Signaling 3) ELISA Kit
E-EL-RB0745	Rabbit HSP-90 (Heat Shock Protein 90) ELISA Kit

Cat.NO	Product Name
E-EL-RB0746	Rabbit ProGRP (Pro-Gastrin Releasing Peptide) ELISA Kit
E-EL-RB0747	Rabbit PA-IgM (Platelet Associated Antibody) ELISA Kit
E-EL-RB0748	Rabbit CFTR (Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit
E-EL-RB0750	Rabbit AIRA (Anti-Insulin Receptor Antibody) ELISA Kit
E-EL-RB0751	Rabbit PCI (Protein C Inhibitor) ELISA Kit
E-EL-RB0752	Rabbit GKRP (Glucokinase Regulatory Protein) ELISA Kit
E-EL-RB0753	Rabbit ATXN1 (Ataxin 1) ELISA Kit
E-EL-RB0754	Rabbit OSCAR (Osteoclast Associated Receptor) ELISA Kit
E-EL-RB0755	Rabbit C7 (Complement Component 7) ELISA Kit
E-EL-RB0756	Rabbit RHOA (Ras Homolog Gene Family, Member A) ELISA Kit
E-EL-RB0757	Rabbit APC (Activated Protein C) ELISA Kit
E-EL-RB0758	Rabbit MN (Metanephrine) ELISA Kit
E-EL-RB0759	Rabbit CASP 9 (Caspase 9) ELISA Kit
E-EL-RB0762	Rabbit PRKA α 1 (Protein Kinase, AMP Activated Alpha 1) ELISA Kit
E-EL-RB0765	Rabbit Cx (Connexin) ELISA Kit
E-EL-RB0767	Rabbit TNK2 (Tyrosine Kinase, Non Receptor 2) ELISA Kit
E-EL-RB0768	Rabbit PR3Ab (Proteinase 3 Antibody) ELISA Kit
E-EL-RB0769	Rabbit IL-11 (Interleukin 11) ELISA Kit
E-EL-RB0772	Rabbit PGA (Pepsinogen A) ELISA Kit
E-EL-RB0773	Rabbit α ZGP1 (Alpha-2-Glycoprotein 1, Zine Binding) ELISA Kit

Cat.NO	Product Name
E-EL-RB0774	Rabbit CLDN9 (Claudin 9) ELISA Kit
E-EL-RB0775	Rabbit γ GT1 (Gamma Glutamyltransferase 1) ELISA Kit
E-EL-RB0776	Rabbit ANG II R2-Ab (Angiotensin II Receptor 2 Antibody) ELISA Kit
E-EL-RB0777	Rabbit I κ B α (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
E-EL-RB0778	Rabbit PAP (Plasmin-Antiplasmin Complex) ELISA Kit
E-EL-RB0779	Rabbit CK-HMW (Cytokeratin High Molecular Weight) ELISA Kit
E-EL-RB0780	Rabbit HMGCS1 (3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit
E-EL-RB0781	Rabbit ECF/CCL11 (Eosinophil Chemotactic Factor) ELISA Kit
E-EL-RB0782	Rabbit GPR37 (G Protein Coupled Receptor 37) ELISA Kit
E-EL-RB0783	Rabbit MUC5AC (Mucin 5 Subtype AC) ELISA Kit
E-EL-RB0784	Rabbit N-MID-OT (N-MID Osteocalcin) ELISA Kit
E-EL-RB0785	Rabbit GUS β (β -glucuronidase) ELISA Kit
E-EL-RB0786	Rabbit ARR β 2 (Arrestin Beta 2) ELISA Kit
E-EL-RB0787	Rabbit DCN (Decorin) ELISA Kit
E-EL-RB0788	Rabbit SAP (Serum Amyloid P Component) ELISA Kit
E-EL-RB0789	Rabbit AGPA (Anti-Glucoprotein 210 Antibody) ELISA Kit
E-EL-RB0790	Rabbit PLC β 1/PI-PLC (Phospholipase C Beta 1, Phosphoinositide Specific) ELISA Kit
E-EL-RB0791	Rabbit GAL1 (Galectin 1) ELISA Kit
E-EL-RB0792	Rabbit PCT (Procalcitonin) ELISA Kit
E-EL-RB0793	Rabbit AATF (Apoptosis Antagonizing Transcription Factor) ELISA Kit

Cat.NO	Product Name
E-EL-RB0794	Rabbit DUB (Deubiquitinating Enzyme) ELISA Kit
E-EL-RB0795	Rabbit DCX (Doublecortin) ELISA Kit
E-EL-RB0796	Rabbit PLN (Phospholamban) ELISA Kit
E-EL-RB0798	Rabbit TTR (Transthyretin) ELISA Kit
E-EL-RB0799	Rabbit CCP (Complement Control Protein) ELISA Kit
E-EL-RB0800	Rabbit GUCA2B (Guanylate Cyclase Activator 2B) ELISA Kit
E-EL-RB0801	Rabbit CLC (Charcot Leyden Crystal Protein) ELISA Kit
E-EL-RB0803	Rabbit IGF2BP3 (Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit
E-EL-RB0804	Rabbit C3c (Complement Fragment 3c) ELISA Kit
E-EL-RB0805	Rabbit ANNA1(Anti-Neuronal Nuclear Antibody 1) ELISA Kit
E-EL-RB0806	Rabbit ApoC4 (Apolipoprotein C4) ELISA Kit
E-EL-RB0807	Rabbit CR2/CD21 (Complement Receptor 2) ELISA Kit
E-EL-RB0808	Rabbit FPA (Fibrinopeptide A) ELISA Kit
E-EL-RB0809	Rabbit RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
E-EL-RB0810	Rabbit α 1-AT (Alpha 1-Antitrypsin) ELISA Kit
E-EL-RB0811	Rabbit MART/Melan-A (Melanoma Marker A) ELISA Kit
E-EL-RB0812	Rabbit SCT (Secretin) ELISA Kit
E-EL-RB0813	Rabbit TJP1 (Tight Junction Protein 1) ELISA Kit
E-EL-RB0814	Rabbit VLA (Very Late Appearing Antigen) ELISA Kit
E-EL-RB0815	Rabbit IL-1ra (Interleukin 1 Receptor Antagonist) ELISA Kit

Cat.NO	Product Name
E-EL-RB0816	Rabbit UPP2 (Uridine Phosphorylase 2) ELISA Kit
E-EL-RB0817	Rabbit FOXC1 (Forkhead Box Protein C1) ELISA Kit
E-EL-RB0818	Rabbit CETP (Cholesteryl Ester Transfer Protein) ELISA Kit
E-EL-RB0820	Rabbit CCND1 (Cyclin-D1) ELISA Kit
E-EL-RB0821	Rabbit MAdCAM1 (Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit
E-EL-RB0822	Rabbit FASN (Fatty Acid Synthase) ELISA Kit
E-EL-RB0823	Rabbit GIP (Gastric Inhibitory Polypeptide) ELISA Kit
E-EL-RB0824	Rabbit PKBG (Protein Kinase B Gamma) ELISA Kit
E-EL-RB0825	Rabbit TPA (Tissue Polypeptide Antigen) ELISA Kit
E-EL-RB0826	Rabbit SRF/IF (Skin Reactive Factor/Inflammatory Factor) ELISA Kit
E-EL-RB0828	Rabbit L1CAM/CD171 (L1-Cell Adhesion Molecule) ELISA Kit
E-EL-RB0829	Rabbit HbC (Hemoglobin C) ELISA Kit
E-EL-RB0830	Rabbit P III CP (Procollagen III C-Terminal ProPeptide) ELISA Kit
E-EL-RB0831	Rabbit PRSS (Serine Protease) ELISA Kit
E-EL-RB0832	Rabbit LH (Luteinizing Hormone) ELISA Kit
E-EL-RB0833	Rabbit C/EBPδ (CCAAT/Enhancer Binding Protein, Delta) ELISA Kit
E-EL-RB0834	Rabbit TACE (TNF α Converting Enzyme) ELISA Kit
E-EL-RB0835	Rabbit CLDN1 (Claudin 1) ELISA Kit
E-EL-RB0836	Rabbit TAFI (Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit
E-EL-RB0837	Rabbit CAPF (Cuticular Active Peptide Factor) ELISA Kit

Cat.NO	Product Name
E-EL-RB0838	Rabbit IGF-1 (Insulin Like Growth Factor 1) ELISA Kit
E-EL-RB0839	Rabbit EL/LIPG (Endothelial Lipase) ELISA Kit
E-EL-RB0840	Rabbit Ext-N (N-Terminal Extein) ELISA Kit
E-EL-RB0841	Rabbit HSP-70/HSPA9 (Heat Shock Protein 70) ELISA Kit
E-EL-RB0842	Rabbit INSL5 (Insulin Like Protein 5) ELISA Kit
E-EL-RB0843	Rabbit AP 13 (Apelin 13) ELISA Kit
E-EL-RB0844	Rabbit TBP2 (Thioredoxin Binding Protein 2) ELISA Kit
E-EL-RB0845	Rabbit TGFβ1 (Transforming Growth Factor Beta Induced Protein) ELISA Kit
E-EL-RB0847	Rabbit RPA-70 (replication protein A) ELISA Kit
E-EL-RB0849	Rabbit DR-70TM (Tumor Marker DR-70 For Lung Cancer) ELISA Kit
E-EL-RB0850	Rabbit PDE4D (Phosphodiesterase 4D, cAMP Specific) ELISA Kit
E-EL-RB0851	Rabbit APAF-1 (Apoptosis Protease Activating Factor-1) ELISA Kit
E-EL-RB0853	Rabbit Gal (α-Galactoyl) ELISA Kit
E-EL-RB0854	Rabbit S100A9 (S100 Calcium Binding Protein A9) ELISA Kit
E-EL-RB0855	Rabbit S100A6 (S100 Calcium Binding Protein A6) ELISA Kit
E-EL-RB0856	Rabbit TSLP (Thymic Stromal Lymphopoietin) ELISA Kit
E-EL-RB0857	Rabbit SPB (Pulmonary Surfactant Associated Protein B) ELISA Kit
E-EL-RB0858	Rabbit bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit
E-EL-RB0859	Rabbit AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit
E-EL-RB0860	Rabbit CSF (Colony Stimulating Factor) ELISA Kit

Cat.NO	Product Name
E-EL-RB0861	Rabbit Gas6 (Growth Arrest Specific Protein 6) ELISA Kit
E-EL-RB0862	Rabbit TGF- β 3 (Transforming Growth Factor β 3) ELISA Kit
E-EL-RB0863	Rabbit MSLN (Mesothelin) ELISA Kit
E-EL-RB0864	Rabbit SDF-1/CXCL12 (Stromal Cell Derived Factor 1) ELISA Kit
E-EL-RB0865	Rabbit GM (Gastric Mucin) ELISA Kit
E-EL-RB0866	Rabbit ESF (Erythropoiesis Stimulating Factor) ELISA Kit
E-EL-RB0867	Rabbit dDAVP (1-Desamino-8D-Arginine Vasopressin) ELISA Kit
E-EL-RB0868	Rabbit GPI (Glucose 6 Phosphate Isomerase) ELISA Kit
E-EL-RB0869	Rabbit TSP-1 (Thrombospondin-1) ELISA Kit
E-EL-RB0870	Rabbit SIGLEC3 (Sialic Acid Binding Ig Like Lectin 3) ELISA Kit
E-EL-RB0871	Rabbit PIINP (Procollagen II N-Terminal ProPeptide) ELISA Kit
E-EL-RB0872	Rabbit CASP10 (Caspase 10) ELISA Kit
E-EL-RB0873	Rabbit GPC3 (Glypican 3) ELISA Kit
E-EL-RB0874	Rabbit GT-34 (Gastrin 34) ELISA Kit
E-EL-RB0875	Rabbit MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
E-EL-RB0876	Rabbit GTF2F1 (General Transcription Factor IIF, Polypeptide 1) ELISA Kit
E-EL-RB0878	Rabbit ANGPTL4 (Angiopoietin Like Protein 4) ELISA Kit
E-EL-RB0879	Rabbit CTAPIII (Connective Tissue Activating Peptide III) ELISA Kit
E-EL-RB0880	Rabbit MIP-3 β /ELC/CCL19 (Macrophage Inflammatory Protein 3 β) ELISA Kit
E-EL-RB0881	Rabbit PGDS (Prostaglandin D Synthase) ELISA Kit

Cat.NO	Product Name
E-EL-RB0882	Rabbit SNC α (Synuclein Alpha) ELISA Kit
E-EL-RB0883	Rabbit 6-OHDA (6-Hydroxydopamine) ELISA Kit
E-EL-RB0884	Rabbit CPN10 (Chaperonin 10) ELISA Kit
E-EL-RB0885	Rabbit PTHR2 (Parathyroid Hormone Receptor 2) ELISA Kit
E-EL-RB0886	Rabbit PKD1 (Protein Kinase D1) ELISA Kit
E-EL-RB0887	Rabbit CA3 (Carbonic Anhydrase III, Muscle Specific) ELISA Kit
E-EL-RB0888	Rabbit GAP43 (Growth Associated Protein 43) ELISA Kit
E-EL-RB0889	Rabbit sICAM-1 (Soluble Intercellular Adhesion Molecule-1) ELISA Kit
E-EL-RB0890	Rabbit NLR (NOD-Like Receptor9) ELISA Kit
E-EL-RB0892	Rabbit FAS/CD95/TNFRSF6 (Factor Related Apoptosis) ELISA Kit
E-EL-RB0893	Rabbit MAP2K6/MKK6 (Mitogen Activated Protein Kinase Kinase 6) ELISA Kit
E-EL-RB0894	Rabbit GPX1 (Glutathione Peroxidase 1) ELISA Kit
E-EL-RB0895	Rabbit EAR3 (Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit
E-EL-RB0896	Rabbit PHC (Primary Hepatic Carcinoma) ELISA Kit
E-EL-RB0897	Rabbit NUMA1 (Nuclear Mitotic Apparatus Protein 1) ELISA Kit
E-EL-RB0898	Rabbit PKM2 (Pyruvate Kinase, Muscle) ELISA Kit
E-EL-RB0899	Rabbit COL7 (Collagen Type VII) ELISA Kit
E-EL-RB0900	Rabbit ALCAM (Activated Leukocyte Cell Adhesion Molecule) ELISA Kit
E-EL-RB0901	Rabbit FBL (Fibrillarlin) ELISA Kit
E-EL-RB0902	Rabbit TDT (Terminal Deoxynucleotidyl Transferase) ELISA Kit

Cat.NO	Product Name
E-EL-RB0903	Rabbit A-CoA (Acetyl Coenzyme A) ELISA Kit
E-EL-RB0904	Rabbit TXLNb (Taxilin Beta) ELISA Kit
E-EL-RB0905	Rabbit NF (Neurofilament) ELISA Kit
E-EL-RB0906	Rabbit OCT4 (Octamer Binding Transcription Factor 4) ELISA Kit
E-EL-RB0908	Rabbit snRNP/Sm (Small Nuclear Ribonucleoprotein) ELISA Kit
E-EL-RB0909	Rabbit FABP4 (Fatty Acid-binding Protein 4, Adipocyte) ELISA Kit
E-EL-RB0910	Rabbit CLDN11 (Claudin 11) ELISA Kit
E-EL-RB0911	Rabbit M-AChR (Muscarinic Acetylcholine Receptor) ELISA Kit
E-EL-RB0912	Rabbit FGF-7/KGF (Fibroblast Growth Factor 7) ELISA Kit
E-EL-RB0913	Rabbit STAT2 (Signal Transducer and Activator of Transcription 2) ELISA Kit
E-EL-RB0914	Rabbit PGF2 α (Prostaglandin F2 α) ELISA Kit
E-EL-RB0915	Rabbit PRSS23 (Protease, Serine 23) ELISA Kit
E-EL-RB0916	Rabbit DNase-I (Deoxyribonuclease I) ELISA Kit
E-EL-RB0917	Rabbit NTRK3 (Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
E-EL-RB0918	Rabbit ASPN (Asporin) ELISA Kit
E-EL-RB0919	Rabbit ft3 (Free Triiodothyronine) ELISA Kit
E-EL-RB0920	Rabbit DHEA (Dehydroepiandrosterone) ELISA Kit
E-EL-RB0921	Rabbit GADL1 (Glutamate Decarboxylase Like Protein 1) ELISA Kit
E-EL-RB0922	Rabbit RAPGEF1 (Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit
E-EL-RB0923	Rabbit PHSA (Poly Serum Antigen) ELISA Kit

Cat.NO	Product Name
E-EL-RB0924	Rabbit GP6 (Glycoprotein VI, Platelet) ELISA Kit
E-EL-RB0925	Rabbit TFF3 (Trefoil Factor 3, Intestinal) ELISA Kit
E-EL-RB0926	Rabbit COL10 (Collagen Type X) ELISA Kit
E-EL-RB0928	Rabbit PAI2 (Plasminogen Activator Inhibitor 2)ELISA Kit
E-EL-RB0929	Rabbit PIAS1 (Protein Inhibitor of Activated STAT 1) ELISA Kit
E-EL-RB0931	Rabbit EGR4 (Early Growth Response Protein 4) ELISA Kit
E-EL-RB0932	Rabbit MYH1 (Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit
E-EL-RB0933	Rabbit ET-1 (Endothelin 1) ELISA Kit
E-EL-RB0934	Rabbit PTGS1 (Prostaglandin Endoperoxide Synthase 1) ELISA Kit
E-EL-RB0935	Rabbit PAEP (Progestagen-Associated Endometrial Protein) ELISA Kit
E-EL-RB0936	Rabbit GZMK (Granzyme K) ELISA Kit
E-EL-RB0937	Rabbit ERCC1 (Excision Repair Cross Complementing Rodent Repair Deficiency Complementation 1) ELISA Kit
E-EL-RB0939	Rabbit CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit
E-EL-RB0941	Rabbit PDGF-AB (Platelet Derived Growth Factor AB) ELISA Kit
E-EL-RB0942	Rabbit FCN1 (Ficolin 1) ELISA Kit
E-EL-RB0943	Rabbit EL (EuglobulinLysis) ELISA Kit
E-EL-RB0944	Rabbit EGR2 (Early Growth Response Protein 2) ELISA Kit
E-EL-RB0945	Rabbit Tn-I (Troponin I) ELISA Kit
E-EL-RB0946	Rabbit CDK2 (Cyclin Dependent Kinase 2) ELISA Kit
E-EL-RB0947	Rabbit SPARC (Secreted Protein Acidic And Rich In Cysteine) ELISA Kit

Cat.NO	Product Name
E-EL-RB0948	Rabbit AAP (Alanine Aminopeptidase) ELISA Kit
E-EL-RB0949	Rabbit NUP85 (Nucleoporin 85kDa) ELISA Kit
E-EL-RB0950	Rabbit CK-19/KRT19 (Cytokeratin 19) ELISA Kit
E-EL-RB0951	Rabbit COL11 α 1 (Collagen Type XI Alpha 1) ELISA Kit
E-EL-RB0952	Rabbit AT (AntiThrombin) ELISA Kit
E-EL-RB0953	Rabbit EIPS2 (Estrogen Induced Protein pS2) ELISA Kit
E-EL-RB0954	Rabbit AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) ELISA Kit
E-EL-RB0955	Rabbit MIA1 (Melanoma Inhibitory Activity Protein 1) ELISA Kit
E-EL-RB0956	Rabbit BMPR1A (Bone Morphogenetic Protein Receptor, type IA) ELISA Kit
E-EL-RB0957	Rabbit PGR (Progesterone Receptor) ELISA Kit
E-EL-RB0958	Rabbit PICK1 (Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
E-EL-RB0959	Rabbit EPHA3 (Ephrin Type A Receptor 3) ELISA Kit
E-EL-RB0960	Rabbit Anti-Surv (Anti-Survivin Antibody) ELISA Kit
E-EL-RB0961	Rabbit GPI (Glycophosphatidylinositol) ELISA Kit
E-EL-RB0962	Rabbit CD99 (Cluster of Differentiation 99) ELISA Kit
E-EL-RB0963	Rabbit PMN Elastase (Polymorphonuclear Elastase) ELISA Kit
E-EL-RB0965	Rabbit HO1 (Heme Oxygenase 1) ELISA Kit
E-EL-RB0966	Rabbit TAT (Thrombin/Antithrombin Complex) ELISA Kit
E-EL-RB0967	Rabbit C1r (Complement Component 1R) ELISA Kit
E-EL-RB0968	Rabbit FN (Fibronectin) ELISA Kit

Cat.NO	Product Name
E-EL-RB0969	Rabbit LRP (Lung Resistance-Related Protein) ELISA Kit
E-EL-RB0971	Rabbit MUC3 (Mucin 3) ELISA Kit
E-EL-RB0972	Rabbit NUP62 (Nucleoporin 62kDa) ELISA Kit
E-EL-RB0973	Rabbit FGFBP2 (Fibroblast Growth Factor Binding Protein 2) ELISA Kit
E-EL-RB0974	Rabbit β -EP (Beta-Endorphin) ELISA Kit
E-EL-RB0976	Rabbit EphA2 (Tyrosine protein kinase receptor A2) ELISA Kit
E-EL-RB0977	Rabbit CR1 (Complement Receptor type 1) ELISA Kit
E-EL-RB0978	Rabbit GST θ 1/GST τ 1 (Glutathione S Transferase Theta 1) ELISA Kit
E-EL-RB0979	Rabbit CENPA (Centromere Protein A) ELISA Kit
E-EL-RB0980	Rabbit CKLFSF5 (Chemokine Like Factor Superfamily 5) ELISA Kit
E-EL-RB0981	Rabbit ADCY2 (Adenylate Cyclase 2) ELISA Kit
E-EL-RB0982	Rabbit TI-Ag (Thymus Independent Antigen) ELISA Kit
E-EL-RB0983	Rabbit PKI β (Protein Kinase Inhibitor Beta) ELISA Kit
E-EL-RB0984	Rabbit GST ω 1 (Glutathione S Transferase Omega 1) ELISA Kit
E-EL-RB0985	Rabbit ANG (Angiogenin) ELISA Kit
E-EL-RB0986	Rabbit UBC (Urinary Bladder Cancer Antigen) ELISA Kit
E-EL-RB0987	Rabbit TFF1 (Trefoil Factor 1) ELISA Kit
E-EL-RB0988	Rabbit THAA (Thyroid Hormone Autoantibody) ELISA Kit
E-EL-RB0989	Rabbit MMP-17 (Matrix Metalloproteinase 17) ELISA Kit
E-EL-RB0990	Rabbit APT (Abnormal Prothrombin) ELISA Kit

Cat.NO	Product Name
E-EL-RB0991	Rabbit GDF2 (Growth Differentiation Factor 2) ELISA Kit
E-EL-RB0992	Rabbit PIIINP (N-Terminal Procollagen III Propeptide) ELISA Kit
E-EL-RB0993	Rabbit ADA (Adenosine Deaminase) ELISA Kit
E-EL-RB0994	Rabbit GRO α /CXCL1/NAP3 (Growth Regulated Oncogene Alpha) ELISA Kit
E-EL-RB0995	Rabbit HPF4(Heparin/Platelet Factor 4 Antibodies)ELISA Kit
E-EL-RB0996	Rabbit ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit
E-EL-RB0997	Rabbit α MSH (Alpha-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-RB0998	Rabbit NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) ELISA Kit
E-EL-RB0999	Rabbit α 2-AP (α 2-Antiplasmin) ELISA Kit
E-EL-RB1000	Rabbit MMP-7 (Matrix Metalloproteinase 7) ELISA Kit
E-EL-RB1001	Rabbit PCDH15 (Protocadherin 15) ELISA Kit
E-EL-RB1002	Rabbit LEP (Leptin) ELISA Kit
E-EL-RB1004	Rabbit ADAMTS4 (ADAM with Thrombospondin Type 1 Motif4) ELISA Kit
E-EL-RB1005	Rabbit DKK2 (Dickkopf Related Protein 2) ELISA Kit
E-EL-RB1006	Rabbit FST (Follistatin) ELISA Kit
E-EL-RB1007	Rabbit ASGR1 (Asialoglycoprotein Receptor 1) ELISA Kit
E-EL-RB1008	Rabbit IgM (Immunoglobulin M) ELISA Kit
E-EL-RB1010	Rabbit OJ/IleRS (Anti-OJ-Antibody) ELISA Kit
E-EL-RB1011	Rabbit PLA2G4D (Phospholipase A2, Group IV D) ELISA Kit
E-EL-RB1012	Rabbit MASP1 (Mannan Associated Serine Protease 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB1013	Rabbit MUC17 (Mucin 17) ELISA Kit
E-EL-RB1014	Rabbit SYK (Spleen Tyrosine Kinase) ELISA Kit
E-EL-RB1015	Rabbit TGM1 (Transglutaminase 1, Keratinocyte) ELISA Kit
E-EL-RB1016	Rabbit C5c (Complement C5 Convertase) ELISA Kit
E-EL-RB1017	Rabbit PSAP (Prosaposin) ELISA Kit
E-EL-RB1018	Rabbit MCP-4/CCL13 (Monocyte Chemotactic Protein 4) ELISA Kit
E-EL-RB1022	Rabbit ATM (Ataxia Telangiectasia Mutated) ELISA Kit
E-EL-RB1023	Rabbit AXIN2 (Axis Inhibition Protein 2) ELISA Kit
E-EL-RB1024	Rabbit ACA-IgA (Anti-Cardiolipin Antibody IgA) ELISA Kit
E-EL-RB1025	Rabbit ATF-4 (Activating Transcription Factor 4) ELISA Kit
E-EL-RB1026	Rabbit ASGR2 (Asialoglycoprotein Receptor 2) ELISA Kit
E-EL-RB1027	Rabbit PTN (Pleiotrophin) ELISA Kit
E-EL-RB1028	Rabbit EloA (Elongin A) ELISA Kit
E-EL-RB1029	Rabbit TGM2 (Transglutaminase 2, Tissue) ELISA Kit
E-EL-RB1032	Rabbit SPJ (Surfactant Associated Protein J) ELISA Kit
E-EL-RB1033	Rabbit Smad4/MADH4 (Mothers Against Decapentaplegic Homolog 4)ELISA Kit
E-EL-RB1034	Rabbit EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit
E-EL-RB1035	Rabbit TNNI2 (Troponin I Type 2, Fast Skeletal) ELISA Kit
E-EL-RB1036	Rabbit PPIL1 (Peptidylprolyl Isomerase/Cyclophilin Like 1) ELISA Kit
E-EL-RB1037	Rabbit BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit

Cat.NO	Product Name
E-EL-RB1038	Rabbit LT-E4 (Leukotriene E4) ELISA Kit
E-EL-RB1039	Rabbit ES-Ab (Anti-Endostatin Antibody) ELISA Kit
E-EL-RB1040	Rabbit BCDF (B-Cell Differentiation Factor) ELISA Kit
E-EL-RB1041	Rabbit SST (Somatostatin) ELISA Kit
E-EL-RB1043	Rabbit TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
E-EL-RB1044	Rabbit PLA2G12B (Phospholipase A2, Group XII B) ELISA Kit
E-EL-RB1045	Rabbit GRN (Granulin) ELISA Kit
E-EL-RB1046	Rabbit HLTF (Helicase Like Transcription Factor) ELISA Kit
E-EL-RB1049	Rabbit TM (Thrombomodulin) ELISA Kit
E-EL-RB1050	Rabbit OSMR (Oncostatin M Receptor) ELISA Kit
E-EL-RB1051	Rabbit CK-17/KRT17 (Cytokeratin 17) ELISA Kit
E-EL-RB1053	Rabbit LDHC (Lactate Dehydrogenase C) ELISA Kit
E-EL-RB1054	Rabbit ACAT2 (Acetyl Coenzyme A Acetyltransferase 2) ELISA Kit
E-EL-RB1055	Rabbit GHRH (Growth Hormone Releasing Hormone) ELISA Kit
E-EL-RB1056	Rabbit MT1E (Metallothionein 1E) ELISA Kit
E-EL-RB1058	Rabbit ADH1 (Alcohol Dehydrogenase 1) ELISA Kit
E-EL-RB1059	Rabbit CA2 (Carbonic Anhydrase II) ELISA Kit
E-EL-RB1060	Rabbit TPO (Thyroid Peroxidase) ELISA Kit
E-EL-RB1061	Rabbit Pro-MMP-10 (Pro-Matrix Metalloproteinase 10) ELISA Kit
E-EL-RB1062	Rabbit DPP6 (Dipeptidyl Peptidase VI) ELISA Kit

Cat.NO	Product Name
E-EL-RB1063	Rabbit Grb10 (Growth Factor Receptor Bound Protein 10) ELISA Kit
E-EL-RB1064	Rabbit HNRPU (Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit
E-EL-RB1065	Rabbit BMPR2 (Bone Morphogenetic Protein Receptor II) ELISA Kit
E-EL-RB1066	Rabbit PCDHβ2 (Protocadherin Beta 2) ELISA Kit
E-EL-RB1067	Rabbit ApoC2 (Apolipoprotein C2) ELISA Kit
E-EL-RB1068	Rabbit GFAP (Glial Fibrillary Acidic Protein) ELISA Kit
E-EL-RB1069	Rabbit MUC4 (Mucin 4) ELISA Kit
E-EL-RB1070	Rabbit Hpc (Hepcidin) ELISA Kit
E-EL-RB1071	Rabbit F13B (Coagulation Factor XIII B Polypeptide) ELISA Kit
E-EL-RB1072	Rabbit DAG1 (Dystrophin Associated Glycoprotein 1) ELISA Kit
E-EL-RB1073	Rabbit GARS (Glycyl tRNA Synthetase) ELISA Kit
E-EL-RB1074	Rabbit RBP4 (Retinol Binding Protein 4, Plasma) ELISA Kit
E-EL-RB1075	Rabbit FABP7 (Fatty Acid Binding Protein 7, Brain) ELISA Kit
E-EL-RB1076	Rabbit Cav-1 (Caveolin-1) ELISA Kit
E-EL-RB1078	Rabbit CEBPγ (CCAAT/Enhancer Binding Protein Gamma) ELISA Kit
E-EL-RB1079	Rabbit CDK4 (Cyclin Dependent Kinase 4) ELISA Kit
E-EL-RB1080	Rabbit HSPβ6/HSP-20 (Heat Shock Protein Beta 6) ELISA Kit
E-EL-RB1081	Rabbit ANG II R1-Ab (Angiotensin II Receptor 1 Antibody) ELISA Kit
E-EL-RB1082	Rabbit NES1 (Nesfatin 1) ELISA Kit
E-EL-RB1083	Rabbit RHOB (Ras Homolog Gene Family, Member B) ELISA Kit

Cat.NO	Product Name
E-EL-RB1084	Rabbit UCHL1 (Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
E-EL-RB1085	Rabbit CD40L/TNFSF5 (Cluster of Differentiation 40 Ligand) ELISA Kit
E-EL-RB1087	Rabbit ODF (Osteoclast Differentiation Factor) ELISA Kit
E-EL-RB1088	Rabbit MPIF-1/CCL23 (Myeloid Progenitor Inhibitory Factor 1) ELISA Kit
E-EL-RB1089	Rabbit DLK-1 (Delta Like 1 Homolog) ELISA Kit
E-EL-RB1091	Rabbit F13A1 (Coagulation Factor XIII A1 Polypeptide) ELISA Kit
E-EL-RB1093	Rabbit GROγ (Growth Regulated Oncogene Gamma) ELISA Kit
E-EL-RB1094	Rabbit NUP133 (Nucleoporin 133kDa) ELISA Kit
E-EL-RB1095	Rabbit f-Hb (Free Haemoglobin) ELISA Kit
E-EL-RB1096	Rabbit AP 17 (Apelin 17) ELISA Kit
E-EL-RB1097	Rabbit ANXA5 (Annexin A5) ELISA Kit
E-EL-RB1098	Rabbit ALDOC (Aldolase C, Fructose Bisphosphate) ELISA Kit
E-EL-RB1100	Rabbit TP/TXA2R (Thromboxane Receptor) ELISA Kit
E-EL-RB1101	Rabbit SHBG (Sex Hormone Binding Globulin) ELISA Kit
E-EL-RB1102	Rabbit C4BP (C4 Binding Protein) ELISA Kit
E-EL-RB1103	Rabbit M-AChR M3 (Muscarinic Acetylcholine Receptor M3) ELISA Kit
E-EL-RB1104	Rabbit PRELP (Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit
E-EL-RB1105	Rabbit AGA (Anti-Glycolipid Antibody) ELISA Kit
E-EL-RB1106	Rabbit ODC (Ornithine Decarboxylase) ELISA Kit
E-EL-RB1107	Rabbit PC/CPG (Choline Phosphoglyceride) ELISA Kit

Cat.NO	Product Name
E-EL-RB1108	Rabbit RND1 (Rho Family GTPase 1) ELISA Kit
E-EL-RB1109	Rabbit IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit
E-EL-RB1110	Rabbit TFR/CD71 (Transferrin Receptor) ELISA Kit
E-EL-RB1111	Rabbit EFNA4 (Ephrin A4) ELISA Kit
E-EL-RB1112	Rabbit TFF2 (Trefoil Factor 2) ELISA Kit
E-EL-RB1114	Rabbit MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit
E-EL-RB1115	Rabbit IGFBP-6 (Insulin Like Growth Factor Binding Protein 6) ELISA Kit
E-EL-RB1117	Rabbit ATGA/TGAB (Anti-Thyroid-Globulin Antibody) ELISA Kit
E-EL-RB1118	Rabbit C5b (Complement Fragment 5b) ELISA Kit
E-EL-RB1119	Rabbit DAF (Decay Accelerating Factor) ELISA Kit
E-EL-RB1120	Rabbit BMF (Bcl-2 Modifying Factor) ELISA Kit
E-EL-RB1121	Rabbit BAG3 (Bcl2 Associated Athanogene 3) ELISA Kit
E-EL-RB1122	Rabbit FOXP3 (Forkhead Box Protein P3) ELISA Kit
E-EL-RB1123	Rabbit CYPA (Cyclophilin A) ELISA Kit
E-EL-RB1125	Rabbit AZU1 (Azurocidin 1) ELISA Kit
E-EL-RB1126	Rabbit ARRβ1 (Arrestin Beta 1) ELISA Kit
E-EL-RB1127	Rabbit PS2 (Presenilin 2) ELISA Kit
E-EL-RB1128	Rabbit PPP1R1A (Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit
E-EL-RB1129	Rabbit BRAK/CXCL14 (Breast and Kidney Expressed Chemokine) ELISA Kit
E-EL-RB1130	Rabbit GJβ1 (Gap Junction Protein Beta 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB1131	Rabbit CENPE (Centromere Protein E) ELISA Kit
E-EL-RB1132	Rabbit F10 (Coagulation Factor X) ELISA Kit
E-EL-RB1133	Rabbit CCKAR (Cholecystokinin A Receptor) ELISA Kit
E-EL-RB1134	Rabbit PKC η (Protein Kinase C Eta) ELISA Kit
E-EL-RB1135	Rabbit CKMT1A (Creatine Kinase, Mitochondrial 1A) ELISA Kit
E-EL-RB1136	Rabbit CA12 (Carbonic Anhydrase XII) ELISA Kit
E-EL-RB1137	Rabbit AKA (Anti-Keratin Antibody) ELISA Kit
E-EL-RB1138	Rabbit ANG I (Angiotensin I) ELISA Kit
E-EL-RB1140	Rabbit CTSB (Cathepsin B) ELISA Kit
E-EL-RB1141	Rabbit CDK5 (Cyclin Dependent Kinase 5) ELISA Kit
E-EL-RB1142	Rabbit IL-8 (Interleukin 8) ELISA Kit
E-EL-RB1143	Rabbit LGALS3BP (Lectin Galactoside Binding, Soluble 3 Binding Protein) ELISA Kit
E-EL-RB1144	Rabbit CD97 (Cluster of Differentiation 97) ELISA Kit
E-EL-RB1145	Rabbit HSP-27/HSPB1 (Heat Shock Protein 27) ELISA Kit
E-EL-RB1147	Rabbit PDH β (Pyruvate Dehydrogenase Beta) ELISA Kit
E-EL-RB1148	Rabbit 15-LO (Arachidonate 15-Lipoxygenase) ELISA Kit
E-EL-RB1149	Rabbit MYL1 (Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit
E-EL-RB1150	Rabbit CX26 (Connexin 26) ELISA Kit
E-EL-RB1151	Rabbit AQP-3 (Aquaporin 3) ELISA Kit
E-EL-RB1152	Rabbit VWF (Von Willebrand Factor) ELISA Kit

Cat.NO	Product Name
E-EL-RB1153	Rabbit PRI (Placental Ribonuclease Inhibitor) ELISA Kit
E-EL-RB1154	Rabbit CD109 (Cluster of Differentiation 109) ELISA Kit
E-EL-RB1155	Rabbit MUC5B (Mucin 5 Subtype B) ELISA Kit
E-EL-RB1157	Rabbit C3c (Complement C3 Convertase) ELISA Kit
E-EL-RB1158	Rabbit PCDH1 (Protocadherin 1) ELISA Kit
E-EL-RB1160	Rabbit GATA2 (GATA Binding Protein 2) ELISA Kit
E-EL-RB1161	Rabbit IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit
E-EL-RB1162	Rabbit TNC (Tenascin C) ELISA Kit
E-EL-RB1165	Rabbit D5R/DRD5 (Dopamine Receptor D5) ELISA Kit
E-EL-RB1166	Rabbit TMPRSS2 (Transmembrane Protease, Serine 2) ELISA Kit
E-EL-RB1167	Rabbit NEFL (Neurofilament, Light Polypeptide) ELISA Kit
E-EL-RB1169	Rabbit GP39 (Glycoprotein 39, Cartilage) ELISA Kit
E-EL-RB1170	Rabbit BCAR1 (Breast Cancer Anti Estrogen Resistance 1) ELISA Kit
E-EL-RB1172	Rabbit TAF1 (TATA Box Binding Protein Associated Factor 1) ELISA Kit
E-EL-RB1173	Rabbit IL-15 (Interleukin 15) ELISA Kit
E-EL-RB1174	Rabbit IL-23 p40(Interleukin 23) ELISA Kit
E-EL-RB1175	Rabbit PROC (Protein C) ELISA Kit
E-EL-RB1177	Rabbit REG1 β (Regenerating Islet Derived Protein 1 Beta) ELISA Kit
E-EL-RB1178	Rabbit CLCF1 (Cardiotrophin Like Cytokine Factor 1) ELISA Kit
E-EL-RB1179	Rabbit MPF (Maturation Promoting Factor) ELISA Kit

Cat.NO	Product Name
E-EL-RB1180	Rabbit MYOG (Myogenin) ELISA Kit
E-EL-RB1181	Rabbit CDC23 (Cell Division Cycle Protein 23) ELISA Kit
E-EL-RB1182	Rabbit PAC3 (Platelet Associated Complement 3) ELISA Kit
E-EL-RB1183	Rabbit NUP155 (Nucleoporin 155kDa) ELISA Kit
E-EL-RB1184	Rabbit NEB (Noradrenaline Bitartas) ELISA Kit
E-EL-RB1185	Rabbit GA (Glycated Albumin) ELISA Kit
E-EL-RB1187	Rabbit CgA (Chromogranin A/Pancreastatin) ELISA Kit
E-EL-RB1188	Rabbit CTR (Calcitonin Receptor) ELISA Kit
E-EL-RB1190	Rabbit β 2-GP1 Ab IgM (β 2-Glycoprotein 1 Antibody IgM) ELISA Kit
E-EL-RB1191	Rabbit SEP1 (Selenoprotein 1) ELISA Kit
E-EL-RB1193	Rabbit PYGM (Glycogen Phosphorylase, Muscle) ELISA Kit
E-EL-RB1194	Rabbit AFP-L3 (Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit
E-EL-RB1195	Rabbit HJV/HFE2 (Hemojuvelin/Hemochromatosis Type 2) ELISA Kit
E-EL-RB1196	Rabbit IDU α (Alpha-L-Iduronidase) ELISA Kit
E-EL-RB1197	Rabbit Ras-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit
E-EL-RB1198	Rabbit SCF (Stem Cell Factor) ELISA Kit
E-EL-RB1199	Rabbit S100A11 (S100 Calcium Binding Protein A11) ELISA Kit
E-EL-RB1200	Rabbit BCL2L1/BCL-X (Bcl-2 Like Protein 1) ELISA Kit
E-EL-RB1201	Rabbit MRC1 (Mannose Receptor C Type 1) ELISA Kit
E-EL-RB1202	Rabbit TM (Thrombin)ELISA Kit

Cat.NO	Product Name
E-EL-RB1203	Rabbit ANXA2 (Annexin A2) ELISA Kit
E-EL-RB1204	Rabbit AKAP1 (A Kinase Anchor Protein 1) ELISA Kit
E-EL-RB1205	Rabbit TPP (Thrombus Precursor Protein) ELISA Kit
E-EL-RB1206	Rabbit α -EP (α -Endomorphin) ELISA Kit
E-EL-RB1207	Rabbit GLRX3 (Glutaredoxin 3) ELISA Kit
E-EL-RB1208	Rabbit IgA (Immunoglobulin A) ELISA Kit
E-EL-RB1209	Rabbit PIAS3 (Protein Inhibitor of Activated STAT3) ELISA Kit
E-EL-RB1210	Rabbit EPO (Erythropoietin) ELISA Kit
E-EL-RB1211	Rabbit LBP (Lipopolysaccharide Binding Protein) ELISA Kit
E-EL-RB1212	Rabbit THR α (Thyroid Hormone Receptor Alpha) ELISA Kit
E-EL-RB1213	Rabbit C/EBP β (CCAAT/Enhancer Binding Protein Beta) ELISA Kit
E-EL-RB1214	Rabbit 5-LO (Arachidonate 5-Lipoxygenase) ELISA Kit
E-EL-RB1215	Rabbit PG (Proteoglycan) ELISA Kit
E-EL-RB1216	Rabbit CVL/EZR (Cytovillin/Ezrin) ELISA Kit
E-EL-RB1217	Rabbit IL-2 (Interleukin 2) ELISA Kit
E-EL-RB1218	Rabbit MBLR (Mannose Binding Lectin Receptor) ELISA Kit
E-EL-RB1219	Rabbit CDK1 (Cyclin Dependent Kinase 1) ELISA Kit
E-EL-RB1220	Rabbit TU3A (TU3A Protein) ELISA Kit
E-EL-RB1221	Rabbit CYSLTR1 (Cysteinyl Leukotriene Receptor 1) ELISA Kit
E-EL-RB1222	Rabbit GnRH (Gonadotropin Releasing Hormone) ELISA Kit

Cat.NO	Product Name
E-EL-RB1223	Rabbit PRL (Prolactin) ELISA Kit
E-EL-RB1225	Rabbit C5a (Complement Component 5a) ELISA Kit
E-EL-RB1226	Rabbit Lptn/XCL1 (Lymphotactin) ELISA Kit
E-EL-RB1227	Rabbit HELIX- II (Type II Collagen Helical Peptide) ELISA Kit
E-EL-RB1228	Rabbit RLR (Retinoic Acid Inducible Gene/RIG-Like Receptor) ELISA Kit
E-EL-RB1229	Rabbit FZD8 (Frizzled Homolog 8) ELISA Kit
E-EL-RB1230	Rabbit CASP12 (Caspase 12) ELISA Kit
E-EL-RB1231	Rabbit TSA (Thyroid Stimulating Antibody) ELISA Kit
E-EL-RB1232	Rabbit ANGPTL1 (Angiopoietin Like Protein 1) ELISA Kit
E-EL-RB1233	Rabbit CX40 (Connexin 40) ELISA Kit
E-EL-RB1234	Rabbit TGF- β 2 (Transforming Growth Factor β 2) ELISA Kit
E-EL-RB1236	Rabbit ST2 (Syntenin 2) ELISA Kit
E-EL-RB1237	Rabbit Fc ϵ R II (Receptor II for the Fc Fragment of IgE) ELISA Kit
E-EL-RB1238	Rabbit AQP-4 (Aquaporin 4) ELISA Kit
E-EL-RB1239	Rabbit CA242(Pancreatic Carcinoma Markers-CA242)ELISA Kit
E-EL-RB1240	Rabbit PARC/CCL18 (Pulmonary Activation Regulated Chemokine) ELISA Kit
E-EL-RB1241	Rabbit HLE (Leukocyte Exastase)ELISA Kit
E-EL-RB1242	Rabbit CaSR (Calcium Sensing Receptor) ELISA Kit
E-EL-RB1243	Rabbit FRS2 (Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit
E-EL-RB1244	Rabbit COL6 α 3 (Collagen Type VI Alpha 3) ELISA Kit

Cat.NO	Product Name
E-EL-RB1245	Rabbit TLR2 (Toll-Like Receptor 2) ELISA Kit
E-EL-RB1246	Rabbit SM/SPH (Sphingomyelin) ELISA Kit
E-EL-RB1247	Rabbit ATF6 (Activating Transcription Factor 6) ELISA Kit
E-EL-RB1248	Rabbit 26S-PSM (26S Proteasome) ELISA Kit
E-EL-RB1249	Rabbit POMC (Pro-Opiomelanocortin) ELISA Kit
E-EL-RB1251	Rabbit Anti-AChR (Anti-Acetylcholine Receptor Antibody) ELISA Kit
E-EL-RB1252	Rabbit PFKP (Phosphofructokinase, Platelet) ELISA Kit
E-EL-RB1253	Rabbit PDGF-AA (Platelet Derived Growth Factor-AA) ELISA Kit
E-EL-RB1254	Rabbit BDNF(Brain Derived Neurotrophic Factor) ELISA Kit
E-EL-RB1255	Rabbit NFκB (Nuclear Factor Kappa B) ELISA Kit
E-EL-RB1256	Rabbit Collagenase I ELISA Kit
E-EL-RB1258	Rabbit COL14 (Collagen Type XIV) ELISA Kit
E-EL-RB1259	Rabbit GCLM (Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit
E-EL-RB1260	Rabbit Arg (Arginase) ELISA Kit
E-EL-RB1261	Rabbit GSK3β (Glycogen Synthase Kinase 3 Beta) ELISA Kit
E-EL-RB1263	Rabbit CCR1 (CC-Chemokine Receptor 1) ELISA Kit
E-EL-RB1264	Rabbit NPY (Neuropeptide Y) ELISA Kit
E-EL-RB1265	Rabbit PPAR-γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
E-EL-RB1266	Rabbit GROβ (Growth Regulated Oncogene Beta) ELISA Kit
E-EL-RB1267	Rabbit DDC (Dopamine Decarboxylase) ELISA Kit

Cat.NO	Product Name
E-EL-RB1268	Rabbit SCTR (Secretin Receptor) ELISA Kit
E-EL-RB1269	Rabbit EM2 (Endomorphin 2) ELISA Kit
E-EL-RB1271	Rabbit HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit
E-EL-RB1272	Rabbit MCRS1 (Microspherule Protein 1) ELISA Kit
E-EL-RB1273	Rabbit TERT (Telomerase Reverse Transcriptase) ELISA Kit
E-EL-RB1274	Rabbit SPD (Pulmonary Surfactant Associated Protein D) ELISA Kit
E-EL-RB1275	Rabbit PLA2G2D (Phospholipase A2, Group II D) ELISA Kit
E-EL-RB1276	Rabbit Rac-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit
E-EL-RB1277	Rabbit AFA (Anti-Filaggrin Antibody) ELISA Kit
E-EL-RB1278	Rabbit NFR κ B (Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit
E-EL-RB1279	Rabbit PACAPR (Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit
E-EL-RB1280	Rabbit VLDLR (Very Low Density Lipoprotein Receptor) ELISA Kit
E-EL-RB1283	Rabbit LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit
E-EL-RB1284	Rabbit PAFAHIB3 (Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit
E-EL-RB1285	Rabbit TR (Testosterone Receptor) ELISA Kit
E-EL-RB1286	Rabbit DMD (Dystrophin) ELISA Kit
E-EL-RB1287	Rabbit NTXI (Cross Linked N-telopeptide of Type I Collagen) ELISA Kit
E-EL-RB1288	Rabbit BAFF/CD257/TNFSF13B (B-Cell Activating Factor) ELISA Kit
E-EL-RB1289	Rabbit MFI2 (Melanotransferrin) ELISA Kit
E-EL-RB1290	Rabbit CYTIP (Cytohesin Interacting Protein) ELISA Kit

Cat.NO	Product Name
E-EL-RB1291	Rabbit GC (Glycocalicin) ELISA Kit
E-EL-RB1292	Rabbit GzmB (Granzyme B) ELISA Kit
E-EL-RB1293	Rabbit NUP153 (Nucleoporin 153kDa) ELISA Kit
E-EL-RB1294	Rabbit ACTN3 (Actinin Alpha 3) ELISA Kit
E-EL-RB1295	Rabbit ELA2B (Elastase 2B) ELISA Kit
E-EL-RB1296	Rabbit PCPE1 (Procollagen C Proteinase Enhancer 1) ELISA Kit
E-EL-RB1297	Rabbit CFB (Complement Factor B) ELISA Kit
E-EL-RB1298	Rabbit PAR2 (Protease Activated Receptor 2) ELISA Kit
E-EL-RB1299	Rabbit FcαR/CD89 (Receptor I for the Fc Fragment of IgA) ELISA Kit
E-EL-RB1300	Rabbit CRH (Corticotropin Releasing Hormone) ELISA Kit
E-EL-RB1301	Rabbit FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit
E-EL-RB1302	Rabbit Pro-MMP-13 (Pro-Matrix Metalloproteinase 13) ELISA Kit
E-EL-RB1303	Rabbit IGFBP-1 (Insulin Like Growth Factor Binding Protein 1) ELISA Kit
E-EL-RB1304	Rabbit NFKBIE (Nuclear Factor Of Kappa Light Polypeptide Gene Enhancer In B-cells Inhibitor Epsilon) ELISA Kit
E-EL-RB1305	Rabbit PDIA5 (Protein Disulfide Isomerase A5) ELISA Kit
E-EL-RB1306	Rabbit MAG (Myelin Associated Glycoprotein) ELISA Kit
E-EL-RB1307	Rabbit NGAL (Neutrophil Gelatinase Associated Lipocalin) ELISA Kit
E-EL-RB1308	Rabbit PAR3 (Protease Activated Receptor 3) ELISA Kit
E-EL-RB1309	Rabbit TSP-2 (Thrombospondin-2) ELISA Kit
E-EL-RB1310	Rabbit cADPRH/CD38 (Cyclic ADP Hydrolase) ELISA Kit

Cat.NO	Product Name
E-EL-RB1311	Rabbit CDT (Carbohydrate Deficient Transferrin) ELISA Kit
E-EL-RB1312	Rabbit TDP43 (Tar DNA Binding Protein 43KDa)ELISA Kit
E-EL-RB1313	Rabbit CA724 (Tumor Marker CA724) ELISA Kit
E-EL-RB1314	Rabbit TPP1 (Tripeptidyl Peptidase I) ELISA Kit
E-EL-RB1315	Rabbit MT1 (Metallothionein 1) ELISA Kit
E-EL-RB1317	Rabbit MGP (Matrix Gla Protein) ELISA Kit
E-EL-RB1318	Rabbit MSM β (Microseminoprotein Beta) ELISA Kit
E-EL-RB1320	Rabbit MMSAM (Melanoma Metastasis Surface Adhesion Molecule) ELISA Kit
E-EL-RB1321	Rabbit COL1 α 2 (Collagen Type I Alpha 2) ELISA Kit
E-EL-RB1322	Rabbit MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit
E-EL-RB1323	Rabbit TMS β 10 (Thymosin Beta 10) ELISA Kit
E-EL-RB1324	Rabbit SYCN (Syncollin) ELISA Kit
E-EL-RB1325	Rabbit APA (Anti-Phospholipid Antibody) ELISA Kit
E-EL-RB1328	Rabbit PRSS12 (Protease, Serine 12) ELISA Kit
E-EL-RB1329	Rabbit PRF1 (Perforin 1/Pore-Forming Protein) ELISA Kit
E-EL-RB1331	Rabbit SEMA3A (Semaphorin 3A) ELISA Kit
E-EL-RB1332	Rabbit CD59 (Protectin) ELISA Kit
E-EL-RB1333	Rabbit C5 (Complement Component 5) ELISA Kit
E-EL-RB1334	Rabbit TAF (TATA Box Binding Protein/TBP-Associated Factors) ELISA Kit
E-EL-RB1335	Rabbit MSR1 (Macrophage Scavenger Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB1336	Rabbit TCF20 (Transcription Factor 20) ELISA Kit
E-EL-RB1337	Rabbit MT (Melatonin) ELISA Kit
E-EL-RB1338	Rabbit CKLFSF4 (Chemokine Like Factor Superfamily 4) ELISA Kit
E-EL-RB1339	Rabbit PKI α (Protein Kinase Inhibitor Alpha) ELISA Kit
E-EL-RB1342	Rabbit SIGLEC10 (Sialic Acid Binding Ig Like Lectin 10) ELISA Kit
E-EL-RB1343	Rabbit VS (Versican/PG-M) ELISA Kit
E-EL-RB1345	Rabbit FSA (Fetal Sulfoslycoprotein Antigen) ELISA Kit
E-EL-RB1347	Rabbit MPA (Mycophenolic Acid) ELISA Kit
E-EL-RB1348	Rabbit UreC (Urease Related Protein C) ELISA Kit
E-EL-RB1349	Rabbit TAB (Thyroxine Antibody) ELISA Kit
E-EL-RB1350	Rabbit F XII a (Activated Coagulation Factor XII) ELISA Kit
E-EL-RB1351	Rabbit PDGFA (Platelet Derived Growth Factor Subunit A) ELISA Kit
E-EL-RB1352	Rabbit α -CTx (Alpha Crosslaps) ELISA Kit
E-EL-RB1355	Rabbit CRY1 (Cryptochrome 1) ELISA Kit
E-EL-RB1356	Rabbit OC/BGP (Osteocalcin/Bone Gla Protein) ELISA Kit
E-EL-RB1357	Rabbit GSTs (Glutathione S-Transferases) ELISA Kit
E-EL-RB1358	Rabbit c-jun ELISA Kit
E-EL-RB1359	Rabbit PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit
E-EL-RB1360	Rabbit LN (Laminin) ELISA Kit
E-EL-RB1361	Rabbit GLUT4 (Glucose Transporter 4) ELISA Kit

Cat.NO	Product Name
E-EL-RB1362	Rabbit CALD (Caldesmon) ELISA Kit
E-EL-RB1363	Rabbit CPP (Copeptin) ELISA Kit
E-EL-RB1364	Rabbit CVNPF (Cobra Venom Neuronal Protective Factor) ELISA Kit
E-EL-RB1365	Rabbit GDI1 (GDP Dissociation Inhibitor 1) ELISA Kit
E-EL-RB1366	Rabbit CSF3/GCSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit
E-EL-RB1367	Rabbit CKMT1B (Creatine Kinase, Mitochondrial 1B) ELISA Kit
E-EL-RB1369	Rabbit C1qR1 (Complement Component 1q Receptor) ELISA Kit
E-EL-RB1370	Rabbit COL1 α 1 (Collagen Type I Alpha 1) ELISA Kit
E-EL-RB1371	Rabbit PI (Proinsulin) ELISA Kit
E-EL-RB1372	Rabbit TP53 (Tumor Protein p53) ELISA Kit
E-EL-RB1373	Rabbit TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
E-EL-RB1374	Rabbit COL6 α 1 (Collagen Type VI Alpha 1) ELISA Kit
E-EL-RB1375	Rabbit CASP4 (Caspase 4) ELISA Kit
E-EL-RB1377	Rabbit CLDN3 (Claudin 3) ELISA Kit
E-EL-RB1378	Rabbit STAT6B (Signal Transducer and Activator of Transcription 6B) ELISA Kit
E-EL-RB1379	Rabbit sTNF-1 (Soluble Tumor Necrosis Factor 1) ELISA Kit
E-EL-RB1380	Rabbit TPO-Ab (Anti-Thyroid-Peroxidase Antibody) ELISA Kit
E-EL-RB1381	Rabbit A β 1-40 (Amyloid Beta Peptide 1-40) ELISA Kit
E-EL-RB1382	Rabbit CP (Ceruloplasmin) ELISA Kit
E-EL-RB1383	Rabbit MAGE (Melanoma-associated Antigen) ELISA Kit

Cat.NO	Product Name
E-EL-RB1384	Rabbit ELA3B (Elastase 3B) ELISA Kit
E-EL-RB1386	Rabbit GDF7 (Growth Differentiation Factor 7) ELISA Kit
E-EL-RB1387	Rabbit MASP2 (Mannan Associated Serine Protease 2) ELISA Kit
E-EL-RB1388	Rabbit MMP-8 (Matrix Metalloproteinase 8) ELISA Kit
E-EL-RB1389	Rabbit CTSL (Cathepsin L) ELISA Kit
E-EL-RB1390	Rabbit MMP-11 (Matrix Metalloproteinase 11) ELISA Kit
E-EL-RB1391	Rabbit ELMO2 (Engulfment And Cell Motility 2) ELISA Kit
E-EL-RB1392	Rabbit ALD (Aldolase) ELISA Kit
E-EL-RB1393	Rabbit CA6 (Carbonic Anhydrase VI) ELISA Kit
E-EL-RB1395	Rabbit Trx (Thioredoxin) ELISA Kit
E-EL-RB1397	Rabbit PLSCR1 (Phospholipid Scramblase 1) ELISA Kit
E-EL-RB1398	Rabbit RB1 (Retinoblastoma Protein 1) ELISA Kit
E-EL-RB1399	Rabbit F8 (Coagulation Factor VIII) ELISA Kit
E-EL-RB1400	Rabbit LIPF (Lipase, Gastric) ELISA Kit
E-EL-RB1401	Rabbit CX43 (Connexin 43) ELISA Kit
E-EL-RB1402	Rabbit Gal-6S (Galactose-6-Sulphatase) ELISA Kit
E-EL-RB1403	Rabbit ASB (Arylsulfatase B) ELISA Kit
E-EL-RB1404	Rabbit TL1A (Tumor Necrosis Factor Related Ligand 1A) ELISA Kit
E-EL-RB1405	Rabbit GSTκ1 (Glutathione S Transferase Kappa 1) ELISA Kit
E-EL-RB1406	Rabbit 11-DH-TXB2 (11-Dehydro-thromboxane B2) ELISA Kit

Cat.NO	Product Name
E-EL-RB1407	Rabbit PRDX1 (Peroxiredoxin 1) ELISA Kit
E-EL-RB1408	Rabbit Cys-C (Cystatin C) ELISA Kit
E-EL-RB1409	Rabbit CX31 (Connexin 31) ELISA Kit
E-EL-RB1410	Rabbit GCSFR (Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
E-EL-RB1411	Rabbit IL-5 (Interleukin 5) ELISA Kit
E-EL-RB1412	Rabbit PTPRO (Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
E-EL-RB1414	Rabbit SAP25 (Synaptosome Associated Protein 25) ELISA Kit
E-EL-RB1415	Rabbit Grb14 (Growth Factor Receptor Bound Protein 14) ELISA Kit
E-EL-RB1416	Rabbit CAM-Ab (Anti-Calmodulin Specific Antibody) ELISA Kit
E-EL-RB1417	Rabbit SOX18 (Sex Determining Region Y Box Protein 18) ELISA Kit
E-EL-RB1418	Rabbit LCA/CD45(Leukocyte Common Antigen) ELISA Kit
E-EL-RB1419	Rabbit AGRN (Agrin) ELISA Kit
E-EL-RB1420	Rabbit LR (Leukoregulin) ELISA Kit
E-EL-RB1421	Rabbit MCT (Mast Cell Tryptase) ELISA Kit
E-EL-RB1422	Rabbit FABP8 (Fatty Acid Binding Protein 8, Myelin) ELISA Kit
E-EL-RB1424	Rabbit GLUT1 (Glucose Transporter 1) ELISA Kit
E-EL-RB1425	Rabbit VCL (Vinculin) ELISA Kit
E-EL-RB1426	Rabbit CYB561 (Cytochrome b-561) ELISA Kit
E-EL-RB1427	Rabbit SEMA3B (Semaphorin 3B) ELISA Kit
E-EL-RB1428	Rabbit NAIP (Neuronal Apoptosis Inhibitory Protein) ELISA Kit

Cat.NO	Product Name
E-EL-RB1429	Rabbit Neo-T4 (Neonatal Thyroxine) ELISA Kit
E-EL-RB1430	Rabbit TA (Anti-Teichoic Acid Antibody) ELISA Kit
E-EL-RB1431	Rabbit EED (Embryonic Ectoderm Development) ELISA Kit
E-EL-RB1432	Rabbit STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit
E-EL-RB1433	Rabbit BAFFR (B Cell Activating Factor Receptor) ELISA Kit
E-EL-RB1434	Rabbit PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit
E-EL-RB1435	Rabbit LpPLA2 (Phospholipase A2, Lipoprotein Associated) ELISA Kit
E-EL-RB1436	Rabbit ICAM-4 (Intercellular Adhesion Molecule-4) ELISA Kit
E-EL-RB1437	Rabbit FGL1 (Fibrinogen Like Protein 1) ELISA Kit
E-EL-RB1438	Rabbit TNX (Tenascin X) ELISA Kit
E-EL-RB1439	Rabbit CPR (Cytochrome P450 Reductase) ELISA Kit
E-EL-RB1440	Rabbit SMM (Smooth Muscle Myosin) ELISA Kit
E-EL-RB1441	Rabbit HbA1c(Glycosylated Hemoglobin/Hemoglobin A1c)ELISA Kit
E-EL-RB1442	Rabbit LP-a (Lipoprotein a) ELISA Kit
E-EL-RB1443	Rabbit TRH (Thyrotropin Releasing Hormone) ELISA Kit
E-EL-RB1444	Rabbit CK-20/KRT20 (Cytokeratin 20) ELISA Kit
E-EL-RB1445	Rabbit RGS4 (Regulator of G Protein Signaling 4) ELISA Kit
E-EL-RB1446	Rabbit TBG (Thyroxine Binding Globulin) ELISA Kit
E-EL-RB1447	Rabbit CCSA-4 (Colon Cancer Specific Antigen-4) ELISA Kit
E-EL-RB1448	Rabbit ENG/CD105 (Endoglin) ELISA Kit

Cat.NO	Product Name
E-EL-RB1449	Rabbit CTX-III (Cross Linked C-telopeptide of Type III Collagen) ELISA Kit
E-EL-RB1450	Rabbit C4b (Complement Component 4b) ELISA Kit
E-EL-RB1451	Rabbit CYP21B (Cytochrome P450 21B) ELISA Kit
E-EL-RB1452	Rabbit TNPO1 (Transportin 1) ELISA Kit
E-EL-RB1453	Rabbit GAL9 (Galectin 9) ELISA Kit
E-EL-RB1454	Rabbit BMP-10 (Bone Morphogenetic Protein 10) ELISA Kit
E-EL-RB1455	Rabbit MAO (Monoamine Oxidase) ELISA Kit
E-EL-RB1456	Rabbit Exendin-4 ELISA Kit
E-EL-RB1457	Rabbit LFA-2/CD2 (Lymphocyte Function Associated Antigen 2) ELISA Kit
E-EL-RB1458	Rabbit Mly/CXCL9 (Monocyte Interferon Gamma Inducing Factor) ELISA Kit
E-EL-RB1459	Rabbit α 2-PI (α 2-Plasmin Inhibitor) ELISA Kit
E-EL-RB1460	Rabbit FG γ (Fibrinogen Gamma) ELISA Kit
E-EL-RB1461	Rabbit MUC20 (Mucin 20) ELISA Kit
E-EL-RB1462	Rabbit IAPP (Islet Amyloid Polypeptide) ELISA Kit
E-EL-RB1464	Rabbit DHEA-S7 (Dehydroepiandrosterone S7) ELISA Kit
E-EL-RB1465	Rabbit ATF7 (Activating Transcription Factor 7) ELISA Kit
E-EL-RB1467	Rabbit FABP1 (Fatty Acid Binding Protein 1, Liver) ELISA Kit
E-EL-RB1468	Rabbit HGF (Hepatocyte Growth Factor) ELISA Kit
E-EL-RB1469	Rabbit SOX3 (Sex Determining Region Y Box Protein 3) ELISA Kit
E-EL-RB1470	Rabbit TrxR1 (Thioredoxin Reductase 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB1471	Rabbit MIgA (Membrane IgA) ELISA Kit
E-EL-RB1472	Rabbit PA2G4 (Proliferation Associated Protein 2G4) ELISA Kit
E-EL-RB1473	Rabbit CYSLTR2 (Cysteinyl Leukotriene Receptor 2) ELISA Kit
E-EL-RB1474	Rabbit Ep-CAM/CD326 (Epithelial Cell Adhesion Molecule) ELISA Kit
E-EL-RB1475	Rabbit PTPRJ (Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit
E-EL-RB1476	Rabbit PKCθ (Protein Kinase C Theta)ELISA Kit
E-EL-RB1477	Rabbit IRAK3 (Interleukin 1 Receptor Associated Kinase 3) ELISA Kit
E-EL-RB1478	Rabbit AIG1 (Androgen Induced Protein 1) ELISA Kit
E-EL-RB1479	Rabbit CKLFSF7 (Chemokine Like Factor Superfamily 7) ELISA Kit
E-EL-RB1480	Rabbit PRKG1 (Protein Kinase, cGMP Dependent Type I) ELISA Kit
E-EL-RB1481	Rabbit HTRA1 (HtrA Serine, Peptidase 1) ELISA Kit
E-EL-RB1483	Rabbit GPER1 (G Protein Coupled Estrogen Receptor 1) ELISA Kit
E-EL-RB1486	Rabbit IL-4 (Interleukin 4) ELISA Kit
E-EL-RB1488	Rabbit FBLN4 (Fibulin 4) ELISA Kit
E-EL-RB1490	Rabbit FGFR1 (Fibroblast Growth Factor Receptor 1) ELISA Kit
E-EL-RB1491	Rabbit CASP14 (Caspase 14) ELISA Kit
E-EL-RB1492	Rabbit MT2 (Metallothionein 2) ELISA Kit
E-EL-RB1493	Rabbit SDH (Sorbitol Dehydrogenase) ELISA Kit
E-EL-RB1494	Rabbit ABCA1(ATP-binding cassette transporter A1) ELISA Kit
E-EL-RB1495	Rabbit α2-M (Alpha-2 Macroglobulin) ELISA Kit

Cat.NO	Product Name
E-EL-RB1496	Rabbit IL-1 α (Interleukin 1 Alpha) ELISA Kit
E-EL-RB1497	Rabbit PF3 (Platelet Factor 3) ELISA Kit
E-EL-RB1498	Rabbit TBP (TATA Binding Protein) ELISA Kit
E-EL-RB1499	Rabbit DLAT (Dihydrolipoyl Transacetylase) ELISA Kit
E-EL-RB1500	Rabbit LDHA (Lactate Dehydrogenase A) ELISA Kit
E-EL-RB1501	Rabbit ETRB (Endothelin Receptor B) ELISA Kit
E-EL-RB1502	Rabbit NT-3 (Neurotrophin-3) ELISA Kit
E-EL-RB1503	Rabbit BPI (Bactericidal/Permeability Increasing Protein) ELISA Kit
E-EL-RB1504	Rabbit ERK (Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-RB1505	Rabbit PDIA6 (Protein Disulfide Isomerase A6) ELISA Kit
E-EL-RB1506	Rabbit GP4 (Platelet Membrane Glycoprotein IV) ELISA Kit
E-EL-RB1507	Rabbit SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-RB1508	Rabbit GPR120 (G Protein Coupled Receptor 120) ELISA Kit
E-EL-RB1510	Rabbit ApoE(Apolipoprotein E)ELISA Kit
E-EL-RB1512	Rabbit TSP-4 (Thrombospondin-4) ELISA Kit
E-EL-RB1513	Rabbit ALK (Anaplastic Lymphoma Kinase) ELISA Kit
E-EL-RB1514	Rabbit LDHB (Lactate Dehydrogenase B) ELISA Kit
E-EL-RB1515	Rabbit AMY1 (Amylase Alpha 1, Salivary) ELISA Kit
E-EL-RB1516	Rabbit KERA (Keratocan) ELISA Kit
E-EL-RB1517	Rabbit β APP (β Amyloid Precursor Protein) ELISA Kit

Cat.NO	Product Name
E-EL-RB1518	Rabbit IGF2R (Insulin Like Growth Factor 2 Receptor) ELISA Kit
E-EL-RB1519	Rabbit CEP110 (Centrosomal Protein 110kDa) ELISA Kit
E-EL-RB1521	Rabbit DPYSL2 (Dihydropyrimidinase Like Protein 2) ELISA Kit
E-EL-RB1522	Rabbit FcγR II/CD32 (Receptor II for the Fc Fragment of IgG) ELISA Kit
E-EL-RB1523	Rabbit PR3 (Proteinase 3) ELISA Kit
E-EL-RB1524	Rabbit TSH (Thyroid Stimulating Hormone) ELISA Kit
E-EL-RB1525	Rabbit CaN (Calcineurin) ELISA Kit
E-EL-RB1526	Rabbit FRA (Fibronectin Related Antigen) ELISA Kit
E-EL-RB1527	Rabbit PGF (Prostaglandin F) ELISA Kit
E-EL-RB1528	Rabbit OLFM4 (Olfactomedin 4) ELISA Kit
E-EL-RB1529	Rabbit EPOR (Erythropoietin Receptor) ELISA Kit
E-EL-RB1530	Rabbit GAL12 (Galectin 12) ELISA Kit
E-EL-RB1531	Rabbit Rock-2 (Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit
E-EL-RB1532	Rabbit MYCBP (C-Myc Binding Protein) ELISA Kit
E-EL-RB1533	Rabbit LTA4H (Leukotriene A4 Hydrolase) ELISA Kit
E-EL-RB1534	Rabbit OPG (Osteoprotegerin) ELISA Kit
E-EL-RB1535	Rabbit TNKS1 (Tankyrase 1) ELISA Kit
E-EL-RB1536	Rabbit COMP (Cartilage Oligomeric Matrix Protein) ELISA Kit
E-EL-RB1537	Rabbit GSK3α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit
E-EL-RB1538	Rabbit PAF (Platelet Activating Factor) ELISA Kit

Cat.NO	Product Name
E-EL-RB1539	Rabbit IL-16 (Interleukin 16) ELISA Kit
E-EL-RB1540	Rabbit MMP-2 (Matrix Metalloproteinase 2) ELISA Kit
E-EL-RB1541	Rabbit RUNX2 (Runt-Related Transcription Factor 2) ELISA Kit
E-EL-RB1542	Rabbit ACTN2 (Actinin Alpha 2) ELISA Kit
E-EL-RB1543	Rabbit CDK9 (Cyclin Dependent Kinase 9) ELISA Kit
E-EL-RB1544	Rabbit NRP1 (Neuropilin 1) ELISA Kit
E-EL-RB1545	Rabbit PAPPA2 (Pappalysin 2) ELISA Kit
E-EL-RB1546	Rabbit PAL (L-Phenylalanine Ammonia-Lyase) ELISA Kit
E-EL-RB1547	Rabbit MK (Midkine) ELISA Kit
E-EL-RB1548	Rabbit CNTN2 (Contactin 2) ELISA Kit
E-EL-RB1549	Rabbit AAA (Anti-Albumin Antibody) ELISA Kit
E-EL-RB1550	Rabbit TXB2 (Thromboxane B2) ELISA Kit
E-EL-RB1551	Rabbit PDIA3 (Protein Disulfide Isomerase A3) ELISA Kit
E-EL-RB1554	Rabbit ITLN1 (Intelectin 1/Omentin) ELISA Kit
E-EL-RB1555	Rabbit ApoA1 (Apolipoprotein A1) ELISA Kit
E-EL-RB1556	Rabbit Grb7 (Growth Factor Receptor Bound Protein 7) ELISA Kit
E-EL-RB1558	Rabbit GM-CSF Ab (Anti-Granulocyte Macrophage Colony Stimulating Factor Antibody) ELISA Kit
E-EL-RB1559	Rabbit VIPR1 (Vasoactive Intestinal Peptide Receptor) ELISA Kit
E-EL-RB1560	Rabbit WARS (Tryptophanyl tRNA Synthetase) ELISA Kit
E-EL-RB1562	Rabbit CCNB (Cyclin-B) ELISA Kit

Cat.NO	Product Name
E-EL-RB1563	Rabbit TRY (Trypsin) ELISA Kit
E-EL-RB1564	Rabbit GKN1 (Gastroke 1) ELISA Kit
E-EL-RB1565	Rabbit CDK6 (Cyclin Dependent Kinase 6) ELISA Kit
E-EL-RB1566	Rabbit ApoF (Apolipoprotein F) ELISA Kit
E-EL-RB1567	Rabbit FG (Fibrinogen) ELISA Kit
E-EL-RB1569	Rabbit KIM-1 (Kidney Injury Molecule 1) ELISA Kit
E-EL-RB1570	Rabbit CK-18/KRT18 (Cytokeratin 18) ELISA Kit
E-EL-RB1571	Rabbit SSRP1 (Structure Specific Recognition Protein 1) ELISA Kit
E-EL-RB1572	Rabbit MYBPC3 (Myosin Binding Protein C, Cardiac) ELISA Kit
E-EL-RB1574	Rabbit IGFBP-2 (Insulin Like Growth Factor Binding Protein 2) ELISA Kit
E-EL-RB1575	Rabbit CCSA-3 (Colon Cancer Specific Antigen-3) ELISA Kit
E-EL-RB1576	Rabbit CMR (Costimulatory Molecules Receptor) ELISA Kit
E-EL-RB1577	Rabbit FETUA (Fetuin A) ELISA Kit
E-EL-RB1578	Rabbit PLA2R1 (Phospholipase A2 Receptor 1) ELISA Kit
E-EL-RB1580	Rabbit APLN (Apelin) ELISA Kit
E-EL-RB1581	Rabbit F7 (Coagulation Factor VII) ELISA Kit
E-EL-RB1582	Rabbit ELA4 (Elastase 4) ELISA Kit
E-EL-RB1584	Rabbit CRN (Corin) ELISA Kit
E-EL-RB1585	Rabbit PKCγ (Protein Kinase C Gamma) ELISA Kit
E-EL-RB1587	Rabbit OGP (Osteogenic Growth Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-RB1588	Rabbit CKLFSF8 (Chemokine Like Factor Superfamily 8) ELISA Kit
E-EL-RB1589	Rabbit GSTθ2/GSTt2 (Glutathione S Transferase Theta 2) ELISA Kit
E-EL-RB1590	Rabbit ES (Endostatin) ELISA Kit
E-EL-RB1591	Rabbit FABP5 (Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
E-EL-RB1592	Rabbit ADCY1 (Adenylate Cyclase 1, Brain) ELISA Kit
E-EL-RB1593	Rabbit CNTNAP1 (Contactin Associated Protein 1) ELISA Kit
E-EL-RB1594	Rabbit S100A7 (S100 Calcium Binding Protein A7) ELISA Kit
E-EL-RB1595	Rabbit Ntn4 (Netrin 4) ELISA Kit
E-EL-RB1597	Rabbit PLC (Phospholipase C) ELISA Kit
E-EL-RB1598	Rabbit MYH9 (Myosin Heavy Chain 9, Non Muscle) ELISA Kit
E-EL-RB1599	Rabbit Bax (Bcl-2 Associated X Protein) ELISA Kit
E-EL-RB1600	Rabbit D3R/DRD3 (Dopamine Receptor D3) ELISA Kit
E-EL-RB1601	Rabbit TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
E-EL-RB1602	Rabbit EFNA1 (Ephrin A1) ELISA Kit
E-EL-RB1603	Rabbit MT3 (Metallothionein 3) ELISA Kit
E-EL-RB1604	Rabbit ERβ (Estrogen Receptor Beta) ELISA Kit
E-EL-RB1605	Rabbit SA (Sialic Acid) ELISA Kit
E-EL-RB1606	Rabbit ACP5 (Tartrate Resistant Acid Phosphatase 5) ELISA Kit
E-EL-RB1607	Rabbit GUCA2A (Guanylate Cyclase Activator 2A) ELISA Kit
E-EL-RB1609	Rabbit MCP-3 (Monocyte Chemotactic Protein 3) ELISA Kit

Cat.NO	Product Name
E-EL-RB1610	Rabbit GPC4 (Glypican 4) ELISA Kit
E-EL-RB1611	Rabbit sgp130 (soluble Glucoprotein 130) ELISA Kit
E-EL-RB1612	Rabbit HSPG (Heparan Sulfate Proteoglycan) ELISA Kit
E-EL-RB1614	Rabbit AP62A (Anti-P62 Antibody) ELISA Kit
E-EL-RB1615	Rabbit MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
E-EL-RB1617	Rabbit PGRP-S (Peptidoglycan Recognition Protein-S) ELISA Kit
E-EL-RB1618	Rabbit PEDF (Pigment Epithelium Derived Factor) ELISA Kit
E-EL-RB1619	Rabbit ATXN2 (Ataxin 2) ELISA Kit
E-EL-RB1620	Rabbit SREBP (Sterol Regulatory Element Binding Proteins) ELISA Kit
E-EL-RB1621	Rabbit HPR (Haptoglobin Related Protein) ELISA Kit
E-EL-RB1622	Rabbit GDF11 (Growth Differentiation Factor 11) ELISA Kit
E-EL-RB1623	Rabbit Ponticulin ELISA Kit
E-EL-RB1625	Rabbit ANG II (Angiotensin II) ELISA Kit
E-EL-RB1626	Rabbit PDI (Protein Disulfide Isomerase) ELISA Kit
E-EL-RB1627	Rabbit ETF β (Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit
E-EL-RB1628	Rabbit PLAP/ALPP (Placental Alkaline Phosphatase) ELISA Kit
E-EL-RB1629	Rabbit PP (Pepsin) ELISA Kit
E-EL-RB1630	Rabbit GSP (Glycosylated Serum Protein) ELISA Kit
E-EL-RB1631	Rabbit PC4 (Positive Cofactor 4) ELISA Kit
E-EL-RB1632	Rabbit KLK6 (Kallikrein 6) ELISA Kit

Cat.NO	Product Name
E-EL-RB1633	Rabbit TGF- α (Transforming Growth Factor α) ELISA Kit
E-EL-RB1634	Rabbit Anti-AlaRS (Anti-Alanyl-tRNA Synthetase/Anti-PL12-Antibody) ELISA Kit
E-EL-RB1635	Rabbit GIF (Gastric Intrinsic Factor) ELISA Kit
E-EL-RB1636	Rabbit CgB (Chromogranin B) ELISA Kit
E-EL-RB1637	Rabbit ApoC1 (Apolipoprotein C1) ELISA Kit
E-EL-RB1638	Rabbit PICP (Procollagen I C-terminal Propeptide) ELISA Kit
E-EL-RB1639	Rabbit CHRM2 (Cholinergic Receptor, Muscarinic 2) ELISA Kit
E-EL-RB1640	Rabbit ET (Epitestosterone) ELISA Kit
E-EL-RB1641	Rabbit MUC21 (Mucin 21) ELISA Kit
E-EL-RB1642	Rabbit RELB (V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
E-EL-RB1643	Rabbit CYP1B1 (Cytochrome P450 1B1) ELISA Kit
E-EL-RB1644	Rabbit DLD (Dihydrolipoyl Dehydrogenase) ELISA Kit
E-EL-RB1645	Rabbit PS (Phosphatidylserine) ELISA Kit
E-EL-RB1646	Rabbit TTN (Titin) ELISA Kit
E-EL-RB1647	Rabbit AR (Amphiregulin) ELISA Kit
E-EL-RB1648	Rabbit HBF (Fetal Hemoglobin) ELISA Kit
E-EL-RB1649	Rabbit ANG1 (Angiopoietin 1) ELISA Kit
E-EL-RB1651	Rabbit SAG (S Antigen) ELISA Kit
E-EL-RB1652	Rabbit CX3CL1 (Chemokine C-X3-C-Motif Ligand) ELISA Kit
E-EL-RB1653	Rabbit GPX3 (Glutathione Peroxidase 3, Plasma) ELISA Kit

Cat.NO	Product Name
E-EL-RB1654	Rabbit ANG II R1/AGTR1 (Angiotensin II Receptor 1) ELISA Kit
E-EL-RB1655	Rabbit GOD (Glucose Oxidase) ELISA Kit
E-EL-RB1657	Rabbit TALLA-1 (T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit
E-EL-RB1658	Rabbit ANEA (Anti-Neutrophil Elastase Antibody) ELISA Kit
E-EL-RB1659	Rabbit PAH (Phenylalanine Hydroxylase) ELISA Kit
E-EL-RB1660	Rabbit SPARCL1 (SPARC Like Protein 1) ELISA Kit
E-EL-RB1661	Rabbit CD30L (Cluster of Differentiation 30 Lligand) ELISA Kit
E-EL-RB1662	Rabbit GNβ1 (G Protein Beta 1) ELISA Kit
E-EL-RB1664	Rabbit AP (Aprotinin) ELISA Kit
E-EL-RB1665	Rabbit GATA3 (GATA Binding Protein 3) ELISA Kit
E-EL-RB1666	Rabbit TACR3 (Tachykinin Receptor 3) ELISA Kit
E-EL-RB1669	Rabbit NCAM/CD56 (Neural Cell Adhesion Molecule) ELISA Kit
E-EL-RB1671	Rabbit GDF1 (Growth Differentiation Factor 1) ELISA Kit
E-EL-RB1672	Rabbit NT-ProBNP (N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
E-EL-RB1673	Rabbit SBEM (Small Breast Epithelial Mucin) ELISA Kit
E-EL-RB1674	Rabbit I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit
E-EL-RB1675	Rabbit CALLA (Common Acute Lymphocytic Leukaemia Antigen) ELISA Kit
E-EL-RB1676	Rabbit FGFBP1 (Fibroblast Growth Factor Binding Protein 1) ELISA Kit
E-EL-RB1677	Rabbit GPC2 (Glypican 2) ELISA Kit
E-EL-RB1678	Rabbit CLEC2C (C-Type Lectin Domain Family 2 Member C) ELISA Kit

Cat.NO	Product Name
E-EL-RB1679	Rabbit PLSCR2 (Phospholipid Scramblase 2) ELISA Kit
E-EL-RB1680	Rabbit Tgtp (T-cell Specific Gtpase) ELISA Kit
E-EL-RB1681	Rabbit CD28 (Cluster of Differentiation 28) ELISA Kit
E-EL-RB1682	Rabbit DEFα1 (Defensin Alpha 1, Neutrophil) ELISA Kit
E-EL-RB1683	Rabbit COL5α2 (Collagen Type V Alpha 2) ELISA Kit
E-EL-RB1685	Rabbit RPL10 (Ribosomal Protein L10) ELISA Kit
E-EL-RB1686	Rabbit TIEG1 (TGF Beta Inducible Early Response Gene 1) ELISA Kit
E-EL-RB1687	Rabbit IL-22 (Interleukin 22) ELISA Kit
E-EL-RB1688	Rabbit PLP (Parathyroid Hormone-Like Protein) ELISA Kit
E-EL-RB1690	Rabbit P-Cadherin (Cadherin, Placental) ELISA Kit
E-EL-RB1692	Rabbit 5-HIAA (5-Hydroxyindolacetic Acid) ELISA Kit
E-EL-RB1693	Rabbit PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit
E-EL-RB1694	Rabbit VIM (Vimentin) ELISA Kit
E-EL-RB1695	Rabbit SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit
E-EL-RB1696	Rabbit PRB2 (Basic Salivary Proline Rich Protein 2) ELISA Kit
E-EL-RB1697	Rabbit Ub (Ubiquitin) ELISA Kit
E-EL-RB1698	Rabbit PCI (Procollagen I) ELISA Kit
E-EL-RB1700	Rabbit IL-3 (Interleukin 3) ELISA Kit
E-EL-RB1701	Rabbit GLUT3 (Glucose Transporter 3) ELISA Kit
E-EL-RB1702	Rabbit ApoA5 (Apolipoprotein A5) ELISA Kit

Cat.NO	Product Name
E-EL-RB1703	Rabbit LSP (Liver Specific Lipoprotein) ELISA Kit
E-EL-RB1704	Rabbit sCD14 (Soluble Cluster of Differentiation14) ELISA Kit
E-EL-RB1705	Rabbit MCAM (Melanoma Cell Adhesion Molecule) ELISA Kit
E-EL-RB1706	Rabbit P-CK (Pan-Cytokeratin) ELISA Kit
E-EL-RB1707	Rabbit CFP (Complement Factor P) ELISA Kit
E-EL-RB1708	Rabbit VTN (Vitronectin) ELISA Kit
E-EL-RB1709	Rabbit VEGF-C (Vascular Endothelial Growth Factor C) ELISA Kit
E-EL-RB1710	Rabbit LT β (Lymphotoxin Beta) ELISA Kit
E-EL-RB1711	Rabbit F-Box Protein 2(FBXO2) ELISA Kit
E-EL-RB1713	Rabbit C1q-Ab (Anti-Complement 1q Antibody) ELISA Kit
E-EL-RB1714	Rabbit Smad7/MADH7 (Mothers Against Decapentaplegic Homolog 7) ELISA Kit
E-EL-RB1715	Rabbit PL (Placental Lactogen) ELISA Kit
E-EL-RB1716	Rabbit LT β R (Lymphotoxin Beta Receptor) ELISA Kit
E-EL-RB1717	Rabbit L-LDH (L-Lactate Dehydrogenase) ELISA Kit
E-EL-RB1718	Rabbit C/EBP α (CCAAT/Enhancer Binding Protein Alpha) ELISA Kit
E-EL-RB1719	Rabbit GR α (Glucocorticoid Receptor Alpha) ELISA Kit
E-EL-RB1720	Rabbit C9 (Complement Component 9) ELISA Kit
E-EL-RB1721	Rabbit PIR (Pirin) ELISA Kit
E-EL-RB1722	Rabbit PAR1 (Protease Activated Receptor 1) ELISA Kit
E-EL-RB1723	Rabbit NMP-22 (Nuclear Matrix Protein 22) ELISA Kit

Cat.NO	Product Name
E-EL-RB1724	Rabbit SEPP1 (Selenoprotein P) ELISA Kit
E-EL-RB1725	Rabbit EPHA1 (Ephrin Type A Receptor 1) ELISA Kit
E-EL-RB1726	Rabbit C4BP β (C4 Binding Protein Beta) ELISA Kit
E-EL-RB1727	Rabbit TLR9 (Toll-Like Receptor 9) ELISA Kit
E-EL-RB1728	Rabbit IgG1 (Immunoglobulin G1) ELISA Kit
E-EL-RB1729	Rabbit PRSS2 (Protease, Serine, 2) ELISA Kit
E-EL-RB1730	Rabbit OPN (Osteopontin) ELISA Kit
E-EL-RB1731	Rabbit ALDOA (Aldolase A, Fructose Bisphosphate) ELISA Kit
E-EL-RB1732	Rabbit GPRC5A (G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit
E-EL-RB1733	Rabbit PLA2G2A (Phospholipase A2, Group II A) ELISA Kit
E-EL-RB1734	Rabbit TMS β 4 (Thymosin Beta 4) ELISA Kit
E-EL-RB1736	Rabbit COL17 (Collagen Type XVII) ELISA Kit
E-EL-RB1737	Rabbit CAMK2 α (Calcium/Calmodulin Dependent Protein Kinase II Alpha) ELISA Kit
E-EL-RB1738	Rabbit FLN β (Filamin B, Beta) ELISA Kit
E-EL-RB1739	Rabbit VPI (Vasoactive Peptidase Inhibitors) ELISA Kit
E-EL-RB1740	Rabbit VEGFR3/Flt4 (Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit
E-EL-RB1741	Rabbit PRDX2 (Peroxiredoxin 2) ELISA Kit
E-EL-RB1742	Rabbit PSG1/SP1 (Pregnancy Specific Beta-1-Glycoprotein 1) ELISA Kit
E-EL-RB1743	Rabbit NRG-3 (Neuregulin 3) ELISA Kit
E-EL-RB1745	Rabbit PCK1 (Phosphoenolpyruvate Carboxykinase 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB1746	Rabbit INSR (Insulin Receptor) ELISA Kit
E-EL-RB1747	Rabbit CTSK (Cathepsin K) ELISA Kit
E-EL-RB1748	Rabbit Tyk-2 (Tyrosine Kinase 2) ELISA Kit
E-EL-RB1749	Rabbit AIM2 (Absent In Melanoma 2) ELISA Kit
E-EL-RB1750	Rabbit BECN1 (Beclin 1) ELISA Kit
E-EL-RB1752	Rabbit IFN- β (Interferon Beta) ELISA Kit
E-EL-RB1753	Rabbit EAR1 (Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit
E-EL-RB1754	Rabbit ACT α 1 (Actin Alpha 1, Skeletal Muscle) ELISA Kit
E-EL-RB1755	Rabbit TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit
E-EL-RB1756	Rabbit COL15 (Collagen Type XV) ELISA Kit
E-EL-RB1757	Rabbit ApoD (Apolipoprotein D) ELISA Kit
E-EL-RB1758	Rabbit MGST1 (Microsomal Glutathione S Transferase 1) ELISA Kit
E-EL-RB1759	Rabbit PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
E-EL-RB1760	Rabbit δ SIP- α (Delta Sleep Inducing Peptide α) ELISA Kit
E-EL-RB1762	Rabbit PTGES2 (Prostaglandin E Synthase 2) ELISA Kit
E-EL-RB1763	Rabbit GREM1 (Gremlin 1) ELISA Kit
E-EL-RB1764	Rabbit PP (Protein Phosphatase) ELISA Kit
E-EL-RB1766	Rabbit CENPB (Centromere Protein B) ELISA Kit
E-EL-RB1768	Rabbit SAA (Serum Amyloid A) ELISA Kit
E-EL-RB1769	Rabbit Salusin α (Salusin Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-RB1770	Rabbit GL β (Galactosidase Beta) ELISA Kit
E-EL-RB1771	Rabbit bTC (Betacellulin) ELISA Kit
E-EL-RB1773	Rabbit DT (Dopachrome Tautomerase) ELISA Kit
E-EL-RB1774	Rabbit FLRT2 (Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit
E-EL-RB1775	Rabbit GST α 5 (Glutathione S Transferase Alpha 5) ELISA Kit
E-EL-RB1776	Rabbit NAP (Neutrophil Alkaline Phosphatase) ELISA Kit
E-EL-RB1777	Rabbit MCP-2/CCL8 (Monocyte Chemotactic Protein 2) ELISA Kit
E-EL-RB1778	Rabbit c-fos ELISA Kit
E-EL-RB1779	Rabbit LMP7/PSMB9 (Low Molecular-weight Protein/proteasome Beta-type Subunit) ELISA Kit
E-EL-RB1780	Rabbit AGT (Angiotensinogen) ELISA Kit
E-EL-RB1781	Rabbit DBP (DNA Binding Protein) ELISA Kit
E-EL-RB1782	Rabbit CFH (Complement Factor H) ELISA Kit
E-EL-RB1783	Rabbit SEMA7A (Semaphorin 7A) ELISA Kit
E-EL-RB1784	Rabbit AIF (Apoptosis Inducing Factor) ELISA Kit
E-EL-RB1785	Rabbit PGD2-MOX (Prostaglandin D2-Methoxime) ELISA Kit
E-EL-RB1786	Rabbit TUB β (β -Tubulin) ELISA Kit
E-EL-RB1787	Rabbit PYGB (Glycogen Phosphorylase, Brain) ELISA Kit
E-EL-RB1788	Rabbit PLTP (Phospholipid Transfer Protein) ELISA Kit
E-EL-RB1789	Rabbit D2R/DRD2 (Dopamine Receptor D2) ELISA Kit
E-EL-RB1790	Rabbit LZM (Lysozyme) ELISA Kit

Cat.NO	Product Name
E-EL-RB1791	Rabbit ACD (Adrenocortical Dysplasia Homolog) ELISA Kit
E-EL-RB1792	Rabbit COL9α3 (Collagen Type IX Alpha 3) ELISA Kit
E-EL-RB1793	Rabbit VMA (Vanillylmandelic Acid) ELISA Kit
E-EL-RB1794	Rabbit SLPI (Secretory Leukocyte Protease Inhibitor) ELISA Kit
E-EL-RB1795	Rabbit IL-17 (Interleukin 17) ELISA Kit
E-EL-RB1798	Rabbit TDGF1 (Teratocarcinoma Derived Growth Factor 1) ELISA Kit
E-EL-RB1799	Rabbit HHIP (Hedgehog Interacting Protein) ELISA Kit
E-EL-RB1800	Rabbit CTAGE (Cutaneous T-Cell Lymphocyte Associated Antigen) ELISA Kit
E-EL-RB1801	Rabbit CENPH (Centromere Protein H) ELISA Kit
E-EL-RB1802	Rabbit TARC (Thymus Activation Regulated Chemokine) ELISA Kit
E-EL-RB1803	Rabbit AIRE (Autoimmune Regulator) ELISA Kit
E-EL-RB1804	Rabbit D1R/DRD1 (Dopamine Receptor D1) ELISA Kit
E-EL-RB1805	Rabbit CCA (Colon Cancer Antigen) ELISA Kit
E-EL-RB1806	Rabbit AChRab (Acetylcholine Receptor Antibody) ELISA Kit
E-EL-RB1808	Rabbit EMP (Erythrocyte Membrane Protein) ELISA Kit
E-EL-RB1809	Rabbit MCSFR (Macrophage Colony Stimulating Factor Receptor) ELISA Kit
E-EL-RB1811	Rabbit MYC (V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit
E-EL-RB1812	Rabbit STAT5A (Signal Transducer and Activator of Transcription 5A) ELISA Kit
E-EL-RB1813	Rabbit NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
E-EL-RB1814	Rabbit τPK1 (Tau-Protein Kinase) ELISA Kit

Cat.NO	Product Name
E-EL-RB1815	Rabbit NEF3 (Neurofilament 3) ELISA Kit
E-EL-RB1816	Rabbit MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
E-EL-RB1817	Rabbit FGα (Fibrinogen Alpha) ELISA Kit
E-EL-RB1818	Rabbit KNG1 (Kininogen 1) ELISA Kit
E-EL-RB1819	Rabbit CTGF (Connective Tissue Growth Factor) ELISA Kit
E-EL-RB1820	Rabbit GAL7 (Galectin 7) ELISA Kit
E-EL-RB1821	Rabbit EGFR (Epidermal Growth Factor Receptor) ELISA Kit
E-EL-RB1822	Rabbit ADAM10 (A Disintegrin And Metalloprotease 10) ELISA Kit
E-EL-RB1823	Rabbit CAPA (Cytotoxin Associated Protein A) ELISA Kit
E-EL-RB1824	Rabbit RCOR3 (Rest Corepressor 3) ELISA Kit
E-EL-RB1825	Rabbit LOXL4 (Lysyl Oxidase Like Protein 4) ELISA Kit
E-EL-RB1826	Rabbit MMRN1 (Multimerin 1) ELISA Kit
E-EL-RB1827	Rabbit Anti-MPO (Anti-Myeloperoxidase Antibody) ELISA Kit
E-EL-RB1828	Rabbit CAM (Calmodulin) ELISA Kit
E-EL-RB1829	Rabbit NUP160 (Nucleoporin 160kDa) ELISA Kit
E-EL-RB1830	Rabbit FOXO4 (Forkhead Box Protein O4) ELISA Kit
E-EL-RB1831	Rabbit PPOX (Protoporphyrinogen Oxidase) ELISA Kit
E-EL-RB1832	Rabbit TECK/CCL25 (Thymus Expressed Chemokine) ELISA Kit
E-EL-RB1833	Rabbit MAF (Macrophage Activating Factor) ELISA Kit
E-EL-RB1834	Rabbit F2 (Coagulation Factor II) ELISA Kit

Cat.NO	Product Name
E-EL-RB1835	Rabbit ApoA4 (Apolipoprotein A4) ELISA Kit
E-EL-RB1837	Rabbit Tn-T (Troponin T) ELISA Kit
E-EL-RB1839	Rabbit NOX5 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit
E-EL-RB1840	Rabbit EXT1 (Exostoses 1) ELISA Kit
E-EL-RB1841	Rabbit HSF1 (Heat Shock Transcription Factor 1) ELISA Kit
E-EL-RB1842	Rabbit RAG-1 (Recombination Activating Gene 1) ELISA Kit
E-EL-RB1843	Rabbit FcγBP (Fc Fragment of IgG Binding Protein) ELISA Kit
E-EL-RB1844	Rabbit NUP188 (Nucleoporin 188kDa) ELISA Kit
E-EL-RB1845	Rabbit LAP (Leucine Aminopeptidase) ELISA Kit
E-EL-RB1846	Rabbit IGF-2 (Insulin Like Growth Factor 2) ELISA Kit
E-EL-RB1847	Rabbit DPYD (Dihydropyrimidine Dehydrogenase) ELISA Kit
E-EL-RB1848	Rabbit PRG4 (Proteoglycan 4) ELISA Kit
E-EL-RB1849	Rabbit MUC7 (Mucin 7) ELISA Kit
E-EL-RB1850	Rabbit PLCγ1 (Phospholipase C Gamma 1) ELISA Kit
E-EL-RB1851	Rabbit ACP1 (Acid Phosphatase 1) ELISA Kit
E-EL-RB1852	Rabbit NPYR1 (Neuropeptide Y Receptor Y1) ELISA Kit
E-EL-RB1853	Rabbit ADAM8 (A Disintegrin And Metalloprotease 8) ELISA Kit
E-EL-RB1854	Rabbit Fbn (Fibrinolysin) ELISA Kit
E-EL-RB1855	Rabbit NPYR5 (Neuropeptide Y Receptor Y5) ELISA Kit
E-EL-RB1856	Rabbit HSP-60/HSPD1 (Heat Shock Protein 60) ELISA Kit

Cat.NO	Product Name
E-EL-RB1858	Rabbit NKA (Neurokinin A) ELISA Kit
E-EL-RB1859	Rabbit PDH α (Pyruvate Dehydrogenase Alpha) ELISA Kit
E-EL-RB1860	Rabbit CTACK (Cutaneous T-cell Attracting Chemokine) ELISA Kit
E-EL-RB1861	Rabbit PIK3C β (Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit
E-EL-RB1862	Rabbit CAD (Caspase Activated Deoxyribonuclease) ELISA Kit
E-EL-RB1863	Rabbit GAL2 (Galectin 2) ELISA Kit
E-EL-RB1864	Rabbit GPC1 (Glypican 1) ELISA Kit
E-EL-RB1865	Rabbit 12-HETE (12-Hydroxyeicosatetraenoic Acid) ELISA Kit
E-EL-RB1866	Rabbit PI3K2 α (Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit
E-EL-RB1867	Rabbit PKB (Protein Kinase B) ELISA Kit
E-EL-RB1870	Rabbit DSP (Desmoplakin) ELISA Kit
E-EL-RB1871	Rabbit CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit
E-EL-RB1873	Rabbit Pro-MMP-1 (pro-Matrix Metalloproteinase 1) ELISA Kit
E-EL-RB1875	Rabbit CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit
E-EL-RB1877	Rabbit MAT (Matrilysin) ELISA Kit
E-EL-RB1878	Rabbit PKC ζ (Protein Kinase C Zeta) ELISA Kit
E-EL-RB1879	Rabbit TAF2 (TATA Box Binding Protein Associated Factor 2) ELISA Kit
E-EL-RB1880	Rabbit PKC δ 1 (Protein Kinase C Delta) ELISA Kit
E-EL-RB1881	Rabbit CNTN3 (Contactin 3) ELISA Kit
E-EL-RB1882	Rabbit C8g (Complement Component 8g) ELISA Kit

Cat.NO	Product Name
E-EL-RB1883	Rabbit SUMO2 (Small Ubiquitin Related Modifier 2) ELISA Kit
E-EL-RB1884	Rabbit LAT (Linker For Activation of T-cell) ELISA Kit
E-EL-RB1885	Rabbit Bid (BH3 Interacting Domain Death Agonist) ELISA Kit
E-EL-RB1886	Rabbit GLUT2 (Glucose Transporter 2) ELISA Kit
E-EL-RB1887	Rabbit ANKRD1 (Ankyrin Repeat Domain 1) ELISA Kit
E-EL-RB1888	Rabbit PKP1 (Plakophilin 1) ELISA Kit
E-EL-RB1889	Rabbit AP (Allopregnanolone) ELISA Kit
E-EL-RB1890	Rabbit NF- κ B p65 (Nuclear Factor- κ B p65) ELISA Kit
E-EL-RB1891	Rabbit PDK4 (Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit
E-EL-RB1892	Rabbit TFR2 (Transferrin Receptor 2) ELISA Kit
E-EL-RB1893	Rabbit PDIA2 (Protein Disulfide Isomerase A2) ELISA Kit
E-EL-RB1894	Rabbit GAD2 (Glutamate Decarboxylase 2) ELISA Kit
E-EL-RB1895	Rabbit GST μ 2/GSTm2 (Glutathione S Transferase Mu 2, Muscle) ELISA Kit
E-EL-RB1896	Rabbit EPP (Erythrocyte Protoporphyrin) ELISA Kit
E-EL-RB1897	Rabbit TNPO2 (Transportin 2) ELISA Kit
E-EL-RB1898	Rabbit IL-13 (Interleukin 13) ELISA Kit
E-EL-RB1899	Rabbit ECE1 (Endothelin Converting Enzyme 1) ELISA Kit
E-EL-RB1900	Rabbit OBTF2 (Octamer Binding Transcription Factor 2) ELISA Kit
E-EL-RB1902	Rabbit PDGF (Platelet-Derived Growth Factor) ELISA Kit
E-EL-RB1903	Rabbit ConA (Concanavalin A) ELISA Kit

Cat.NO	Product Name
E-EL-RB1904	Rabbit BCR (B-Cell Receptor) ELISA Kit
E-EL-RB1905	Rabbit ACHAP (Acetylcholinesterase Associated Protein) ELISA Kit
E-EL-RB1906	Rabbit MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit
E-EL-RB1908	Rabbit PDGFB (Platelet Derived Growth Factor Subunit B) ELISA Kit
E-EL-RB1909	Rabbit PLB (Phospholipase B) ELISA Kit
E-EL-RB1910	Rabbit PLD2 (Phospholipase D2) ELISA Kit
E-EL-RB1911	Rabbit pERK (Phospho Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-RB1912	Rabbit Fc γ R3B (Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit
E-EL-RB1913	Rabbit MLCK/MYLK (Myosin Light Chain Kinase) ELISA Kit
E-EL-RB1914	Rabbit GDF10 (Growth Differentiation Factor 10) ELISA Kit
E-EL-RB1915	Rabbit PAPPA (Pregnancy Associated Plasma Protein A/ Pappalysin-1) ELISA Kit
E-EL-RB1917	Rabbit IgA2 (Immunoglobulin A2) ELISA Kit
E-EL-RB1918	Rabbit NUP98 (Nucleoporin 98kDa) ELISA Kit
E-EL-RB1919	Rabbit CREB (Cyclic AMP Response Element Binding Protein) ELISA Kit
E-EL-RB1920	Rabbit MS (Melatonin Sulfate) ELISA Kit
E-EL-RB1921	Rabbit ADIPOR2 (Adiponectin Receptor 2) ELISA Kit
E-EL-RB1922	Rabbit ApoB100 (Apolipoprotein B100) ELISA Kit
E-EL-RB1923	Rabbit MAPK (Mitogen Activated Protein Kinase) ELISA Kit
E-EL-RB1924	Rabbit GDF9 (Growth Differentiation Factor 9) ELISA Kit
E-EL-RB1926	Rabbit ANGPT4 (Angiopoietin 4) ELISA Kit

Cat.NO	Product Name
E-EL-RB1927	Rabbit Ank-B (Ankyrin B) ELISA Kit
E-EL-RB1928	Rabbit GRIK2 (Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit
E-EL-RB1929	Rabbit HB (Hemoglobin) ELISA Kit
E-EL-RB1930	Rabbit CT-1 (Cardiotrophin 1) ELISA Kit
E-EL-RB1931	Rabbit SHC/SLP-76 (SH2-Containing Protein) ELISA Kit
E-EL-RB1932	Rabbit DS (Dermatan Sulfate) ELISA Kit
E-EL-RB1933	Rabbit CLL1 (Collectin Liver 1) ELISA Kit
E-EL-RB1934	Rabbit PTPRB (Protein Tyrosine Phosphatase Receptor Type B) ELISA Kit
E-EL-RB1935	Rabbit SDF2L1 (Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
E-EL-RB1936	Rabbit CD23/FcεRII (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
E-EL-RB1937	Rabbit C1q (Complement 1q) ELISA Kit
E-EL-RB1938	Rabbit MFAP4 (Microfibrillar Associated Protein 4) ELISA Kit
E-EL-RB1940	Rabbit CEA (Carcinoembryonic Antigen) ELISA Kit
E-EL-RB1942	Rabbit IAA (Insulin Autoantibody) ELISA Kit
E-EL-RB1943	Rabbit IGFBP-3 (Insulin Like Growth Factor Binding Protein 3) ELISA Kit
E-EL-RB1944	Rabbit fFDFT1 (Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit
E-EL-RB1945	Rabbit AGA (Anti-Gliadin Antibody) ELISA Kit
E-EL-RB1946	Rabbit TNR (Tenascin R) ELISA Kit
E-EL-RB1947	Rabbit P27 (P27 Protein) ELISA Kit
E-EL-RB1949	Rabbit HSP-40 (Heat Shock Protein 40) ELISA Kit

Cat.NO	Product Name
E-EL-RB1950	Rabbit PXR (Pregnane X Receptor) ELISA Kit
E-EL-RB1951	Rabbit P2 (Myelin Protein 2) ELISA Kit
E-EL-RB1952	Rabbit NRG-1 (Neuregulin 1) ELISA Kit
E-EL-RB1953	Rabbit PI3K (Phosphoinositide-3-Kinase) ELISA Kit
E-EL-RB1955	Rabbit ON (Osteonectin) ELISA Kit
E-EL-RB1956	Rabbit ARHGD1 α (Rho GDP Dissociation Inhibitor Alpha) ELISA Kit
E-EL-RB1957	Rabbit DSG-E1 (Desmogleins 1) ELISA Kit
E-EL-RB1958	Rabbit IL-24 (Interleukin 24) ELISA Kit
E-EL-RB1959	Rabbit CSF2R α (Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit
E-EL-RB1960	Rabbit MRC1 (Mannose Receptor C Type 1 Like Protein 1) ELISA Kit
E-EL-RB1961	Rabbit Apo A1-Ab (Anti-Apolipoprotein A1 antibody) ELISA Kit
E-EL-RB1962	Rabbit EREG (Epiregulin) ELISA Kit
E-EL-RB1963	Rabbit CC16 (Clara Cell 16KD Protein) ELISA Kit
E-EL-RB1965	Rabbit LIFR (Leukemia Inhibitory Factor Receptor) ELISA Kit
E-EL-RB1966	Rabbit COL4 (Collagen Type IV) ELISA Kit
E-EL-RB1967	Rabbit PAIg (Platelet Associated Immunoglobulins) ELISA Kit
E-EL-RB1969	Rabbit VEGF165 (Vascular Endothelial Growth Factor165) ELISA Kit
E-EL-RB1970	Rabbit HLA-A (Leukocyte Antigen A) ELISA Kit
E-EL-RB1971	Rabbit TNNI1 (Troponin I Type 1, Slow Skeletal) ELISA Kit
E-EL-RB1972	Rabbit ERN1 (Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB1973	Rabbit FcγR3A (Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit
E-EL-RB1974	Rabbit IL-18 (Interleukin 18) ELISA Kit
E-EL-RB1975	Rabbit ASM (Acid Sphingomyelinase) ELISA Kit
E-EL-RB1976	Rabbit DEFα5 (Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
E-EL-RB1977	Rabbit KLK8 (Kallikrein 8) ELISA Kit
E-EL-RB1978	Rabbit MLZE (Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit
E-EL-RB1979	Rabbit CLEC4E (C-Type Lectin Domain Family 4 Member E) ELISA Kit
E-EL-RB1980	Rabbit LHRH/GnRH (Luteinizing Hormone Releasing Hormone) ELISA Kit
E-EL-RB1981	Rabbit NFYC (Nuclear Transcription Factor Y Gamma) ELISA Kit
E-EL-RB1982	Rabbit CD163 (Cluster of Differentiation 163) ELISA Kit
E-EL-RB1983	Rabbit MF/MPF (Maturation Promoting Factor) ELISA Kit
E-EL-RB1984	Rabbit cTnT/TNNT2 (Troponin T Type 2, Cardiac) ELISA Kit
E-EL-RB1985	Rabbit SPC (Pulmonary Surfactant Associated Protein C) ELISA Kit
E-EL-RB1986	Rabbit C4BPα (C4 Binding Protein Alpha) ELISA Kit
E-EL-RB1987	Rabbit D-LDH (D-Lactate Dehydrogenase) ELISA Kit
E-EL-RB1989	Rabbit DHEA-S (Dehydroepiandrosterone Sulfate) ELISA Kit
E-EL-RB1991	Rabbit SCFR/c-Kit (Stem Cell Factor Receptor) ELISA Kit
E-EL-RB1992	Rabbit SDC1 (Syndecan 1) ELISA Kit
E-EL-RB1993	Rabbit DPP10 (Dipeptidyl Peptidase X) ELISA Kit
E-EL-RB1994	Rabbit MT1M (Metallothionein 1M) ELISA Kit

Cat.NO	Product Name
E-EL-RB1995	Rabbit SF1 (Steroidogenic Factor 1) ELISA Kit
E-EL-RB1996	Rabbit PTH (Parathyroid Hormone) ELISA Kit
E-EL-RB1997	Rabbit MMP-9 (Matrix Metalloproteinase 9) ELISA Kit
E-EL-RB1998	Rabbit PON3 (Paraoxonase 3) ELISA Kit
E-EL-RB1999	Rabbit MR (Mannose Receptor) ELISA Kit
E-EL-RB2000	Rabbit IgA1 (Immunoglobulin A1) ELISA Kit
E-EL-RB2001	Rabbit sAC (Soluble Adenylate Cyclase) ELISA Kit
E-EL-RB2002	Rabbit ARSA (Arylsulfatase A) ELISA Kit
E-EL-RB2003	Rabbit ANGPTL6 (Angiopoietin Like Protein 6) ELISA Kit
E-EL-RB2004	Rabbit EPYC (Epiphycan) ELISA Kit
E-EL-RB2005	Rabbit GAL3 (Galectin 3) ELISA Kit
E-EL-RB2006	Rabbit LIPH (Lipase H) ELISA Kit
E-EL-RB2007	Rabbit α -SCA (α -sarcomeric Actin) ELISA Kit
E-EL-RB2008	Rabbit CCKBR (Cholecystokinin B Receptor) ELISA Kit
E-EL-RB2009	Rabbit Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit
E-EL-RB2011	Rabbit S100A13 (S100 Calcium Binding Protein A13) ELISA Kit
E-EL-RB2013	Rabbit PRB (Retinoblastoma Tumor Suppressor Protein) ELISA Kit
E-EL-RB2014	Rabbit GHITM (Growth Hormone Inducible Transmembrane Protein) ELISA Kit
E-EL-RB2015	Rabbit Pg (Progesterone) ELISA Kit
E-EL-RB2016	Rabbit TD-Ag (Thymus Dependent Antigen) ELISA Kit

Cat.NO	Product Name
E-EL-RB2017	Rabbit MCSP (Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit
E-EL-RB2018	Rabbit FADS2 (Fatty Acid Desaturase 2) ELISA Kit
E-EL-RB2019	Rabbit Anti-GM1 (Anti-Ganglioside M1 Antibody) ELISA Kit
E-EL-RB2020	Rabbit FBLN5 (Fibulin 5) ELISA Kit
E-EL-RB2022	Rabbit UBE3A (Ubiquitin Protein Ligase E3A) ELISA Kit
E-EL-RB2023	Rabbit NAGase (N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
E-EL-RB2024	Rabbit NUP37 (Nucleoporin 37kDa) ELISA Kit
E-EL-RB2025	Rabbit BAD (Bcl2 Associated Death Promoter) ELISA Kit
E-EL-RB2026	Rabbit COL13 (Collagen Type XIII) ELISA Kit
E-EL-RB2027	Rabbit TG (Thyroglobulin) ELISA Kit
E-EL-RB2028	Rabbit HRG (Histidine Rich Glycoprotein) ELISA Kit
E-EL-RB2029	Rabbit ALDH (Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit
E-EL-RB2030	Rabbit GH (Growth Hormone) ELISA Kit
E-EL-RB2032	Rabbit PDPN (Podoplanin) ELISA Kit
E-EL-RB2033	Rabbit Methylase ELISA Kit
E-EL-RB2034	Rabbit FDP (Fibrinogen Degradation Product) ELISA Kit
E-EL-RB2035	Rabbit TNF- β (Tumor Necrosis Factor Beta) ELISA Kit
E-EL-RB2036	Rabbit GATA1 (GATA Binding Protein 1) ELISA Kit
E-EL-RB2037	Rabbit PIK3C α (Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit
E-EL-RB2038	Rabbit RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit

Cat.NO	Product Name
E-EL-RB2039	Rabbit BAK1 (Bcl2 Antagonist/Killer 1) ELISA Kit
E-EL-RB2040	Rabbit PYGL (Glycogen Phosphorylase, Liver) ELISA Kit
E-EL-RB2041	Rabbit SIAE (Sialic Acid Acetylesterase) ELISA Kit
E-EL-RB2043	Rabbit LCN1 (Lipocalin 1) ELISA Kit
E-EL-RB2044	Rabbit FSH (Follicle-Stimulating Hormone) ELISA Kit
E-EL-RB2045	Rabbit MEPE (Matrix Extracellular Phosphoglycoprotein) ELISA Kit
E-EL-RB2047	Rabbit NGF (Nerve Growth Factor) ELISA Kit
E-EL-RB2048	Rabbit IFN- α (Interferon Alpha) ELISA Kit
E-EL-RB2049	Rabbit fT4 (Free Thyroxine) ELISA Kit
E-EL-RB2050	Rabbit MAN1A1 (Mannosidase Alpha Class 1A Member 1) ELISA Kit
E-EL-RB2051	Rabbit RNASEH (Ribonuclease H)ELISA Kit
E-EL-RB2052	Rabbit BALP (Bone Alkaline Phosphatase) ELISA Kit
E-EL-RB2053	Rabbit PTGES (Prostaglandin E Synthase, Microsomal)ELISA Kit
E-EL-RB2054	Rabbit RXFP1 (Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
E-EL-RB2056	Rabbit RPL26 (Ribosomal Protein L26) ELISA Kit
E-EL-RB2058	Rabbit LAB7-2/CD86 (B-Lymphocyte Activation Antigen B7-2) ELISA Kit
E-EL-RB2059	Rabbit VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
E-EL-RB2060	Rabbit ATF-5 (Activating Transcription Factor 5) ELISA Kit
E-EL-RB2061	Rabbit IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit
E-EL-RB2062	Rabbit PCIII (Procollagen III) ELISA Kit

Cat.NO	Product Name
E-EL-RB2063	Rabbit OFQ/N (Orphanin FQ/Nociceptin) ELISA Kit
E-EL-RB2064	Rabbit IL-6R (Interleukin 6 Receptor) ELISA Kit
E-EL-RB2066	Rabbit PRCP (Prolylcarboxypeptidase) ELISA Kit
E-EL-RB2067	Rabbit MTL (Motilin) ELISA Kit
E-EL-RB2069	Rabbit GSANA (Granulocyte Specific Anti-Nuclear Antibody) ELISA Kit
E-EL-RB2070	Rabbit TPPP (Tubulin Polymerization Promoting Protein) ELISA Kit
E-EL-RB2072	Rabbit FP (Prostaglandin F Receptor) ELISA Kit
E-EL-RB2073	Rabbit ERK3 (Extracellular Signal Regulated Kinase 3) ELISA Kit
E-EL-RB2074	Rabbit EG-VEGF/PK1 (Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit
E-EL-RB2075	Rabbit TSP1 (Thrombin Sensitive Protein 1) ELISA Kit
E-EL-RB2077	Rabbit GMFβ (Glia Maturation Factor, Beta) ELISA Kit
E-EL-RB2078	Rabbit EMMPRIN/CD147 (Extracellular Matrix Metalloproteinase Inducer) ELISA Kit
E-EL-RB2080	Rabbit FOXO2 (Forkhead Box Protein O2) ELISA Kit
E-EL-RB2081	Rabbit MBP (Myelin Basic Protein) ELISA Kit
E-EL-RB2082	Rabbit SELL (L-Selectin) ELISA Kit
E-EL-RB2083	Rabbit FHA (Filamentous Hemagglutinin) ELISA Kit
E-EL-RB2084	Rabbit SMAF (Specific Macrophage Arming Factor) ELISA Kit
E-EL-RB2085	Rabbit NCAD (Neural Cadherin) ELISA Kit
E-EL-RB2086	Rabbit VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
E-EL-RB2087	Rabbit LILRB4 (Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit

Cat.NO	Product Name
E-EL-RB2088	Rabbit FETUB (Fetuin B) ELISA Kit
E-EL-RB2089	Rabbit GBP4 (Guanylate Binding Protein 4) ELISA Kit
E-EL-RB2090	Rabbit CD3d (Cluster of Differentiation 3d) ELISA Kit
E-EL-RB2091	Rabbit CKLFSF6 (Chemokine Like Factor Superfamily 6) ELISA Kit
E-EL-RB2092	Rabbit GYPE (Glycophorin E) ELISA Kit
E-EL-RB2093	Rabbit ALT (Alanine Transaminase) ELISA Kit
E-EL-RB2094	Rabbit PFKM (Phosphofructokinase, Muscle) ELISA Kit
E-EL-RB2095	Rabbit AhR (Aryl Hydrocarbon Receptor) ELISA Kit
E-EL-RB2096	Rabbit PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit
E-EL-RB2097	Rabbit IL-20 (Interleukin 20) ELISA Kit
E-EL-RB2099	Rabbit NUP205 (Nucleoporin 205kDa) ELISA Kit
E-EL-RB2100	Rabbit CYP24A1 (Cytochrome P450 24A1) ELISA Kit
E-EL-RB2101	Rabbit GCK (Glucokinase) ELISA Kit
E-EL-RB2103	Rabbit GLTP (Glycolipid Transfer Protein) ELISA Kit
E-EL-RB2104	Rabbit SEMA4C (Semaphorin 4C) ELISA Kit
E-EL-RB2107	Rabbit BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit
E-EL-RB2108	Rabbit OXB (Orexin B) ELISA Kit
E-EL-RB2110	Rabbit MYO (Myosin) ELISA Kit
E-EL-RB2111	Rabbit LPL (Lipoprotein Lipase) ELISA Kit
E-EL-RB2113	Rabbit SNRPC (Small Nuclear Ribonucleoprotein Polypeptide C) ELISA Kit

Cat.NO	Product Name
E-EL-RB2115	Rabbit FGF6 (Fibroblast Growth Factor 6) ELISA Kit
E-EL-RB2117	Rabbit DPYSL3 (Dihydropyrimidinase Like Protein 3) ELISA Kit
E-EL-RB2118	Rabbit ADRP (Adipose Differentiation Related Protein) ELISA Kit
E-EL-RB2120	Rabbit PKP2 (Plakophilin 2) ELISA Kit
E-EL-RB2121	Rabbit α -La (Alpha-Lactalbumin) ELISA Kit
E-EL-RB2122	Rabbit GDH/GLDH (Glutamate dehydrogenase) ELISA Kit
E-EL-RB2123	Rabbit TYR (Tyrosinase) ELISA Kit
E-EL-RB2124	Rabbit C1INH (Complement 1 Inhibitor) ELISA Kit
E-EL-RB2125	Rabbit TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
E-EL-RB2126	Rabbit TAA (Tumor Associated Antigen) ELISA Kit
E-EL-RB2127	Rabbit BMG/ β 2-MG (Beta-2-Microglobulin) ELISA Kit
E-EL-RB2128	Rabbit PF4 (Platelet Factor 4) ELISA Kit
E-EL-RB2129	Rabbit α -GST (α -glutathione S-transferases) ELISA Kit
E-EL-RB2130	Rabbit CYP1A2 (Cytochrome P450 1A2) ELISA Kit
E-EL-RB2131	Rabbit CK-LMW (Cytokeratin Low Molecular Weight) ELISA Kit
E-EL-RB2132	Rabbit PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit
E-EL-RB2133	Rabbit CYPB (Cyclophilin B) ELISA Kit
E-EL-RB2134	Rabbit MCF (Macrophage Chemotatic Factor) ELISA Kit
E-EL-RB2135	Rabbit FGFR4 (Fibroblast Growth Factor Receptor 4) ELISA Kit
E-EL-RB2136	Rabbit GAL14 (Galectin 14) ELISA Kit

Cat.NO	Product Name
E-EL-RB2137	Rabbit ER α (Estrogen Receptor Alpha) ELISA Kit
E-EL-RB2138	Rabbit FSCN (Fascin) ELISA Kit
E-EL-RB2139	Rabbit EFNA3 (Ephrin A3) ELISA Kit
E-EL-RB2140	Rabbit SOD3 (Superoxide Dismutase 3, Extracellular) ELISA Kit
E-EL-RB2141	Rabbit CNR-1 (Cadherin Related Neuronal Receptor1) ELISA Kit
E-EL-RB2142	Rabbit UBAP2 (Ubiquitin Associated Protein 2) ELISA Kit
E-EL-RB2143	Rabbit EPHA10 (Ephrin Type A Receptor 10) ELISA Kit
E-EL-RB2144	Rabbit CACYBP (Calcyclin Binding Protein) ELISA Kit
E-EL-RB2146	Rabbit S100A12 (S100 Calcium Binding Protein A12) ELISA Kit
E-EL-RB2148	Rabbit CASP1 (Caspase 1) ELISA Kit
E-EL-RB2149	Rabbit CS (Chondroitin Sulfate) ELISA Kit
E-EL-RB2151	Rabbit MYH10 (Myosin Heavy Chain 10, Non Muscle) ELISA Kit
E-EL-RB2152	Rabbit RAG-2 (Recombination Activating Gene 2) ELISA Kit
E-EL-RB2153	Rabbit HA (Hemagglutinin) ELISA Kit
E-EL-RB2154	Rabbit PDGFD (Platelet Derived Growth Factor D) ELISA Kit
E-EL-RB2155	Rabbit SIGLEC8 (Sialic Acid Binding Ig Like Lectin 8) ELISA Kit
E-EL-RB2156	Rabbit PLP (Proteolipid Protein) ELISA Kit
E-EL-RB2157	Rabbit EJ/GlyRS (Anti-EJ-Antibody) ELISA Kit
E-EL-RB2158	Rabbit GHR (Growth Hormone Receptor) ELISA Kit
E-EL-RB2159	Rabbit RNASE (Ribonuclease) ELISA Kit

Cat.NO	Product Name
E-EL-RB2160	Rabbit TP (Thymidine Phosphorylase) ELISA Kit
E-EL-RB2161	Rabbit PTPRK (Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit
E-EL-RB2162	Rabbit IDE (Insulin Degrading Enzyme) ELISA Kit
E-EL-RB2163	Rabbit GATA5 (GATA Binding Protein 5) ELISA Kit
E-EL-RB2164	Rabbit PKR1 (Prokineticin Receptor 1) ELISA Kit
E-EL-RB2165	Rabbit GH2 (Growth Hormone 2) ELISA Kit
E-EL-RB2166	Rabbit GL α (Galactosidase Alpha) ELISA Kit
E-EL-RB2167	Rabbit ACA-IgG (Anti-Cardiolipin Antibody IgG) ELISA Kit
E-EL-RB2168	Rabbit VacA (Vacuolating Cytotoxin A) ELISA Kit
E-EL-RB2169	Rabbit TFPI (Tissue Factor Pathway Inhibitor) ELISA Kit
E-EL-RB2170	Rabbit CD42 (Cluster of Differentiation 42) ELISA Kit
E-EL-RB2171	Rabbit MMP-24 (Matrix Metalloproteinase 24) ELISA Kit
E-EL-RB2172	Rabbit CA5B (Carbonic Anhydrase VB, Mitochondrial) ELISA Kit
E-EL-RB2173	Rabbit SELP (P-Selectin) ELISA Kit
E-EL-RB2174	Rabbit TPO (Thrombopoietin) ELISA Kit
E-EL-RB2175	Rabbit TH (Tyrosine Hydroxylase) ELISA Kit
E-EL-RB2176	Rabbit PCSK9 (Proprotein Convertase 9) ELISA Kit
E-EL-RB2177	Rabbit SYNPO (Synaptopodin) ELISA Kit
E-EL-RB2178	Rabbit REG1 α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit
E-EL-RB2179	Rabbit CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit

Cat.NO	Product Name
E-EL-RB2180	Rabbit XCR1 (XC-Chemokine Receptor 1) ELISA Kit
E-EL-RB2181	Rabbit MTF (Microtransferrin) ELISA Kit
E-EL-RB2182	Rabbit α Actin (Alpha-Actin) ELISA Kit
E-EL-RB2183	Rabbit OTR (Oxytocin Receptor) ELISA Kit
E-EL-RB2184	Rabbit INSL3 (Insulin Like Protein 3) ELISA Kit
E-EL-RB2185	Rabbit REG3 γ (Regenerating Islet Derived Protein 3 Gamma) ELISA Kit
E-EL-RB2186	Rabbit GAP (GTPase Activating Protein) ELISA Kit
E-EL-RB2187	Rabbit iPLA2 (Phospholipidase A2, Calcium Independent) ELISA Kit
E-EL-RB2189	Rabbit FGF4 (Fibroblast Growth Factor 4) ELISA Kit
E-EL-RB2190	Rabbit CTSD (Cathepsin D) ELISA Kit
E-EL-RB2191	Rabbit AST (Aspartate Aminotransferase) ELISA Kit
E-EL-RB2192	Rabbit AGER (Advanced Glycosylation End Product Specific Receptor) ELISA Kit
E-EL-RB2194	Rabbit GDF6 (Growth Differentiation Factor 6) ELISA Kit
E-EL-RB2195	Rabbit FG β (Fibrinogen Beta) ELISA Kit
E-EL-RB2196	Rabbit CTX-II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit
E-EL-RB2197	Rabbit ACTH (Adrenocorticotrophic Hormone) ELISA Kit
E-EL-RB2198	Rabbit CHAc (Choline acetylase) ELISA Kit
E-EL-RB2199	Rabbit ADP/Acrp30 (Adiponectin) ELISA Kit
E-EL-RB2200	Rabbit RI (Ribonuclease Inhibitor) ELISA Kit
E-EL-RB2201	Rabbit FGF9 (Fibroblast Growth Factor 9) ELISA Kit

Cat.NO	Product Name
E-EL-RB2202	Rabbit GAD1 (Glutamate Decarboxylase 1) ELISA Kit
E-EL-RB2204	Rabbit HO2 (Heme Oxygenase 2) ELISA Kit
E-EL-RB2205	Rabbit BF (Factor B) ELISA Kit
E-EL-RB2206	Rabbit SDC3 (Syndecan 3) ELISA Kit
E-EL-RB2207	Rabbit CCND3 (Cyclin-D3) ELISA Kit
E-EL-RB2208	Rabbit EP3 (Prostaglandin E Receptor 3) ELISA Kit
E-EL-RB2209	Rabbit C/EBP ϵ (CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit
E-EL-RB2211	Rabbit TAF13 (TATA Box Binding Protein Associated Factor 13) ELISA Kit
E-EL-RB2212	Rabbit ApoC3 (Apolipoprotein C3) ELISA Kit
E-EL-RB2213	Rabbit FOXP1 (Forkhead Box Protein P1) ELISA Kit
E-EL-RB2214	Rabbit sIgA (Secretory Immunoglobulin A) ELISA Kit
E-EL-RB2216	Rabbit GP-73 (Golgi Protein 73) ELISA Kit
E-EL-RB2217	Rabbit PMP3 (Peroxisomal Membrane Protein 3) ELISA Kit
E-EL-RB2218	Rabbit BK (Bradykinin) ELISA Kit
E-EL-RB2219	Rabbit UPP1 (Uridine Phosphorylase 1) ELISA Kit
E-EL-RB2220	Rabbit GAL4 (Galectin 4) ELISA Kit
E-EL-RB2222	Rabbit REG3 α (Regenerating Islet Derived Protein 3 Alpha) ELISA Kit
E-EL-RB2223	Rabbit H-Cadherin (Cadherin, Heart) ELISA Kit
E-EL-RB2224	Rabbit CYP1A1 (Cytochrome P450 1A1) ELISA Kit
E-EL-RB2225	Rabbit MMP-13 (Matrix Metalloproteinase 13) ELISA Kit

Cat.NO	Product Name
E-EL-RB2226	Rabbit LDLR (Low Density Lipoprotein Receptor) ELISA Kit
E-EL-RB2227	Rabbit TNKS2 (Tankyrase 2) ELISA Kit
E-EL-RB2228	Rabbit MRC2 (Mannose Receptor C Type 2) ELISA Kit
E-EL-RB2229	Rabbit MASP (Mannan Binding Lectin Associated Serine Protease) ELISA Kit
E-EL-RB2230	Rabbit MOG (Myelin Oligodendrocyte Glycoprotein) ELISA Kit
E-EL-RB2231	Rabbit Des (Desmin) ELISA Kit
E-EL-RB2232	Rabbit ACVA (Activin A) ELISA Kit
E-EL-RB2233	Rabbit FAK (Focal Adhesion Kinase) ELISA Kit
E-EL-RB2234	Rabbit CPE (Carboxypeptidase E) ELISA Kit
E-EL-RB2235	Rabbit MSP (Macrophage Stimulating Protein) ELISA Kit
E-EL-RB2237	Rabbit MAU (Microalbuminuria) ELISA Kit
E-EL-RB2238	Rabbit TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) ELISA Kit
E-EL-RB2239	Rabbit NAP3 (Neutrophil Activating Protein 3) ELISA Kit
E-EL-RB2240	Rabbit ATF-2 (Activating Transcription Factor 2) ELISA Kit
E-EL-RB2241	Rabbit LTC4 (Leukotriene C4) ELISA Kit
E-EL-RB2242	Rabbit P-PKC (Phospho Protein Kinase C) ELISA Kit
E-EL-RB2243	Rabbit TARS (Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit
E-EL-RB2244	Rabbit SSTR1 (Somatostatin Receptor 1) ELISA Kit
E-EL-RB2245	Rabbit FBN2 (Fibrillin 2) ELISA Kit
E-EL-RB2246	Rabbit ETRA (Endothelin Receptor A) ELISA Kit

Cat.NO	Product Name
E-EL-RB2247	Rabbit TSTA (Tumor Specific Transplantation Antigen) ELISA Kit
E-EL-RB2248	Rabbit ORM2 (Orosomucoid 2) ELISA Kit
E-EL-RB2249	Rabbit ADRA1A (Adrenergic Receptor alpha-1A) ELISA Kit
E-EL-RB2250	Rabbit UBE2A (Ubiquitin Conjugating Enzyme E2A) ELISA Kit
E-EL-RB2251	Rabbit PT (Prothrombin) ELISA Kit
E-EL-RB2252	Rabbit VF (Visfatin) ELISA Kit
E-EL-RB2253	Rabbit ERK2 (Extracellular Signal Regulated Kinase 2) ELISA Kit
E-EL-RB2254	Rabbit GC (Glucagon) ELISA Kit
E-EL-RB2255	Rabbit FcαR (Fc Fragment of IgA Receptor) ELISA Kit
E-EL-RB2257	Rabbit LEPR (Leptin Receptor) ELISA Kit
E-EL-RB2258	Rabbit COX-2R (Cyclooxygenase-2 Receptor) ELISA Kit
E-EL-RB2259	Rabbit CNR1 (Cannabinoid Receptor 1, Brain) ELISA Kit
E-EL-RB2260	Rabbit Gzma (Granzyme A) ELISA Kit
E-EL-RB2261	Rabbit ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-RB2262	Rabbit Ta1 (Thymosin-α1) ELISA Kit
E-EL-RB2263	Rabbit PS1 (Presenilin 1) ELISA Kit
E-EL-RB2264	Rabbit PACAP-38 (Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
E-EL-RB2265	Rabbit PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-RB2266	Rabbit ApoB48 (Apolipoprotein B48) ELISA Kit
E-EL-RB2267	Rabbit hs-CRP (high-sensitivity C-reactive protein) ELISA Kit

Cat.NO	Product Name
E-EL-RB2268	Rabbit A1 α (Protein Kinase A1 Alpha)ELISA Kit
E-EL-RB2269	Rabbit A1 β (Protein Kinase A1 Beta)ELISA Kit
E-EL-RB2270	Rabbit A2 α (Protein Kinase A2 Alpha)ELISA Kit
E-EL-RB2271	Rabbit A2 β (Protein Kinase A2 Beta)ELISA Kit
E-EL-RB2272	Rabbit PTK6 (Protein Tyrosine Kinase 6) ELISA Kit
E-EL-RB2273	Rabbit cGMP(cyclic guanosine monophosphate)ELISA Kit
E-EL-RB2274	Rabbit INS (Insulin) ELISA Kit
E-EL-RB2275	Rabbit HMGB-1(High mobility group protein B1) ELISA Kit
E-EL-RB2276	Rabbit CAT(Catalase) ELISA kit
E-EL-GP0001	Guinea pig C3 (Complement Component 3) ELISA Kit
E-EL-GP0002	Guinea pig IgG (Immunoglobulin G) ELISA Kit
E-EL-GP0003	Guinea pig cTn-I/TNNI3 (cardiac Troponin I) ELISA Kit
E-EL-GP0004	Guinea pig MYO (Myoglobin) ELISA Kit
E-EL-GP0005	Guinea pig IFN- α (Interferon Alpha) ELISA Kit
E-EL-GP0006	Guinea pig NGF (Nerve Growth Factor) ELISA Kit
E-EL-GP0007	Guinea pig TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit
E-EL-GP0009	Guinea pig PEDF (Pigment Epithelium Derived Factor) ELISA Kit
E-EL-GP0010	Guinea pig CGRP (Calcitonin Gene Related Peptide) ELISA Kit
E-EL-GP0011	Guinea pig IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-GP0013	Guinea pig EGF (Epidermal Growth Factor) ELISA Kit

Cat.NO	Product Name
E-EL-GP0014	Guinea pig TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
E-EL-GP0015	Guinea pig NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit
E-EL-GP0016	Guinea pig IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) ELISA Kit
E-EL-GP0017	Guinea pig IFN- γ (Interferon Gamma) ELISA Kit
E-EL-GP0018	Guinea pig IgG1 (Immunoglobulin G1) ELISA Kit
E-EL-GP0019	Guinea pig IgG2 α (Immunoglobulin G2 α) ELISA Kit
E-EL-GP0021	Guinea pig GPLD1 (Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Kit
E-EL-GP0022	Guinea pig CTSD (Cathepsin D) ELISA Kit
E-EL-GP0023	Guinea pig VLA (Very Late Appearing Antigen) ELISA Kit
E-EL-GP0024	Guinea pig GKN2 (Gastroke 2) ELISA Kit
E-EL-GP0026	Guinea pig TE (Telomerase) ELISA Kit
E-EL-GP0027	Guinea pig GAL7 (Galectin 7) ELISA Kit
E-EL-GP0028	Guinea pig PDGFA (Platelet Derived Growth Factor Subunit A) ELISA Kit
E-EL-GP0029	Guinea pig MCSP (Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit
E-EL-GP0030	Guinea pig PP13 (Placental Protein13) ELISA Kit
E-EL-GP0031	Guinea pig TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-GP0032	Guinea pig IgE (Immunoglobulin E) ELISA Kit
E-EL-GP0034	Guinea pig SGK2 (Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit
E-EL-GP0035	Guinea pig NFR κ B (Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit
E-EL-GP0036	Guinea pig GKRP (Glucokinase Regulatory Protein) ELISA Kit

Cat.NO	Product Name
E-EL-GP0037	Guinea pig EG-VEGF/PK1 (Endocrine Gland Derived Vascular Endothelial Growth Factor)ELISA Kit
E-EL-GP0038	Guinea pig iPLA2 (Phospholipidase A2, Calcium Independent) ELISA Kit
E-EL-GP0040	Guinea pig IL-2 (Interleukin 2) ELISA Kit
E-EL-GP0041	Guinea pig GP (Anti-Glycoprotein Antibody) ELISA Kit
E-EL-GP0042	Guinea pig MlgD (Membrane IgD) ELISA Kit
E-EL-GP0043	Guinea pig IL-4 (Interleukin 4) ELISA Kit
E-EL-GP0044	Guinea pig PRSS8 (Protease, Serine 8) ELISA Kit
E-EL-GP0045	Guinea pig IL-6(Interleukin 6)ELISA Kit
E-EL-GP0047	Guinea pig CNTN2 (Contactin 2) ELISA Kit
E-EL-GP0048	Guinea pig IL-12 (Interleukin 12) ELISA Kit
E-EL-GP0050	Guinea pig CNR-1 (Cadherin Related Neuronal Receptor1) ELISA Kit
E-EL-GP0051	Guinea pig CTGF (Connective Tissue Growth Factor) ELISA Kit
E-EL-GP0052	Guinea pig HSF2 (Heat Shock Transcription Factor 2) ELISA Kit
E-EL-GP0053	Guinea pig IL-17 (Interleukin 17) ELISA Kit
E-EL-GP0054	Guinea pig Des (Desmin) ELISA Kit
E-EL-GP0055	Guinea pig F11 (Coagulation Factor XI) ELISA Kit
E-EL-GP0056	Guinea pig GAL3 (Galectin 3) ELISA Kit
E-EL-GP0058	Guinea pig OxLDL (Oxidized Low Density Lipoprotein) ELISA Kit
E-EL-GP0061	Guinea pig NSE (Neuron-Specific Enolase) ELISA Kit
E-EL-GP0062	Guinea pig β 2-GP1 (Beta2 Glycoprotein 1) ELISA Kit

Cat.NO	Product Name
E-EL-GP0063	Guinea pig CENPF (Centromere Protein F) ELISA Kit
E-EL-GP0065	Guinea pig SIGLEC3 (Sialic Acid Binding Ig Like Lectin 3) ELISA Kit
E-EL-GP0066	Guinea pig CFTR (Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit
E-EL-GP0067	Guinea pig MIP-5 (Macrophage Inflammatory Protein 5) ELISA Kit
E-EL-GP0068	Guinea pig CA1 (Carbonic Anhydrase I) ELISA Kit
E-EL-GP0070	Guinea pig Cx (Connexin) ELISA Kit
E-EL-GP0072	Guinea pig MIP-3 β /ELC/CCL19 (Macrophage Inflammatory Protein 3 β) ELISA Kit
E-EL-GP0073	Guinea pig C-P (C-Peptide) ELISA Kit
E-EL-GP0074	Guinea pig SOD3 (Superoxide Dismutase 3, Extracellular) ELISA Kit
E-EL-GP0076	Guinea pig GDF2 (Growth Differentiation Factor 2) ELISA Kit
E-EL-GP0077	Guinea pig PRELP (Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit
E-EL-GP0078	Guinea pig NOX5 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit
E-EL-GP0079	Guinea pig GAL4 (Galectin 4) ELISA Kit
E-EL-GP0080	Guinea pig CLA (Collagen Autoantibody) ELISA Kit
E-EL-GP0081	Guinea pig NA/NE (Noradrenaline/Norepinephrine) ELISA Kit
E-EL-GP0082	Guinea pig CDK2 (Cyclin Dependent Kinase 2) ELISA Kit
E-EL-GP0083	Guinea pig PAH (Phenylalanine Hydroxylase) ELISA Kit
E-EL-GP0084	Guinea pig MVP (Major Vault Protein) ELISA Kit
E-EL-GP0085	Guinea pig ARO (Aromatase) ELISA Kit
E-EL-GP0086	Guinea pig GLUT4 (Glucose Transporter 4) ELISA Kit

Cat.NO	Product Name
E-EL-GP0087	Guinea pig PI3K2 α (Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit
E-EL-GP0088	Guinea pig AGA (Anti-Glycolipid Antibody) ELISA Kit
E-EL-GP0089	Guinea pig GAL9 (Galectin 9) ELISA Kit
E-EL-GP0091	Guinea pig RHOA (Ras Homolog Gene Family, Member A) ELISA Kit
E-EL-GP0092	Guinea pig CD97 (Cluster of Differentiation 97) ELISA Kit
E-EL-GP0093	Guinea pig TAP (Transporter Associated with Antigen Processing) ELISA Kit
E-EL-GP0095	Guinea pig Fc α R/CD89 (Receptor I for the Fc Fragment of IgA) ELISA Kit
E-EL-GP0096	Guinea pig CTNN β 1 (Catenin, Beta 1) ELISA Kit
E-EL-GP0097	Guinea pig ANXA5 (Annexin A5) ELISA Kit
E-EL-GP0098	Guinea pig FGF23 (Fibroblast Growth Factor 23) ELISA Kit
E-EL-GP0099	Guinea pig CNTNAP2 (Contactin Associated Protein Like Protein 2) ELISA Kit
E-EL-GP0100	Guinea pig Tn-T (Troponin T) ELISA Kit
E-EL-GP0101	Guinea pig ABCG1 (ATP Binding Cassette Transporter G1) ELISA Kit
E-EL-GP0103	Guinea pig LFA-2/CD2 (Lymphocyte Function Associated Antigen 2) ELISA Kit
E-EL-GP0104	Guinea pig PROZ (Protein Z) ELISA Kit
E-EL-GP0105	Guinea pig Bax (Bcl-2 Associated X Protein) ELISA Kit
E-EL-GP0106	Guinea pig TLA (Thymus Leukemia Antigen) ELISA Kit
E-EL-GP0107	Guinea pig AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) ELISA Kit
E-EL-GP0108	Guinea pig FGF-7/KGF (Fibroblast Growth Factor 7) ELISA Kit
E-EL-GP0109	Guinea pig FMOD (Fibromodulin) ELISA Kit

Cat.NO	Product Name
E-EL-GP0110	Guinea pig CAM (Calmodulin) ELISA Kit
E-EL-GP0111	Guinea pig CLS (Cardiolipin Synthase) ELISA Kit
E-EL-GP0112	Guinea pig α ZGP1 (Alpha-2-Glycoprotein 1, Zine Binding) ELISA Kit
E-EL-GP0113	Guinea pig PLN (Phospholamban) ELISA Kit
E-EL-GP0114	Guinea pig Ponticulin ELISA Kit
E-EL-GP0115	Guinea pig PP (Pancreatic Polypeptide) ELISA Kit
E-EL-GP0116	Guinea pig LpPLA2 (Phospholipase A2, Lipoprotein Associated) ELISA Kit
E-EL-GP0117	Guinea pig TAP (Trypsinogen Activation Peptide) ELISA Kit
E-EL-GP0118	Guinea pig PDH E1 (Pyruvate Dehydrogenase E1) ELISA Kit
E-EL-GP0119	Guinea pig NRG-1 (Neuregulin 1) ELISA Kit
E-EL-GP0120	Guinea pig AP62A (Anti-P62 Antibody) ELISA Kit
E-EL-GP0121	Guinea pig HSP-40 (Heat Shock Protein 40) ELISA Kit
E-EL-GP0122	Guinea pig DLD (Dihydrolipoyl Dehydrogenase) ELISA Kit
E-EL-GP0123	Guinea pig CRP (C-reactive protein) ELISA Kit
E-EL-GP0124	Guinea pig GDH/GLDH (Glutamate dehydrogenase) ELISA Kit
E-EL-GP0125	Guinea pig ORM2 (Orosomucoid 2) ELISA Kit
E-EL-GP0126	Guinea pig DPYD (Dihydropyrimidine Dehydrogenase) ELISA Kit
E-EL-GP0127	Guinea pig ApoC4 (Apolipoprotein C4) ELISA Kit
E-EL-GP0128	Guinea pig UGCG (UDP-glucose ceramide glucosyltransferase) ELISA Kit
E-EL-GP0129	Guinea pig RBP4 (Retinol Binding Protein 4, Plasma) ELISA Kit

Cat.NO	Product Name
E-EL-GP0131	Guinea pig CFB (Complement Factor B) ELISA Kit
E-EL-GP0133	Guinea pig LDHC (Lactate Dehydrogenase C) ELISA Kit
E-EL-GP0134	Guinea pig OSCAR (Osteoclast Associated Receptor) ELISA Kit
E-EL-GP0135	Guinea pig BAK1 (Bcl2 Antagonist/Killer 1) ELISA Kit
E-EL-GP0136	Guinea pig GR α (Glucocorticoid Receptor Alpha) ELISA Kit
E-EL-GP0137	Guinea pig C1r (Complement Component 1R) ELISA Kit
E-EL-GP0138	Guinea pig MAT (Matrilysin) ELISA Kit
E-EL-GP0139	Guinea pig M-AChR (Muscarinic Acetylcholine Receptor) ELISA Kit
E-EL-GP0140	Guinea pig MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit
E-EL-GP0141	Guinea pig MMP-13 (Matrix Metalloproteinase 13) ELISA Kit
E-EL-GP0142	Guinea pig F-Box Protein 2(FBXO2) ELISA Kit
E-EL-GP0144	Guinea pig BK (Bradykinin) ELISA Kit
E-EL-GP0145	Guinea pig PLAA (Phospholipase A2 Activating Protein) ELISA Kit
E-EL-GP0146	Guinea pig AQP-2 (Aquaporin 2) ELISA Kit
E-EL-GP0147	Guinea pig β Manase ELISA Kit
E-EL-GP0148	Guinea pig MCAM (Melanoma Cell Adhesion Molecule) ELISA Kit
E-EL-GP0149	Guinea pig MART/Melan-A (Melanoma Marker A) ELISA Kit
E-EL-GP0150	Guinea pig C9 (Complement Component 9) ELISA Kit
E-EL-GP0152	Guinea pig ITLN1 (Intelectin 1/Omentin) ELISA Kit
E-EL-GP0153	Guinea pig UCP-2 (Uncoupling Protein 2, Mitochondrial) ELISA Kit

Cat.NO	Product Name
E-EL-GP0154	Guinea pig SOCS3 (Suppressors of Cytokine Signaling 3) ELISA Kit
E-EL-GP0155	Guinea pig SEMA3A (Semaphorin 3A) ELISA Kit
E-EL-GP0157	Guinea pig FSH (Follicle-Stimulating Hormone) ELISA Kit
E-EL-GP0158	Guinea pig TGF- α (Transforming Growth Factor α) ELISA Kit
E-EL-GP0159	Guinea pig ACTg1 (Actin Gamma 1) ELISA Kit
E-EL-GP0160	Guinea pig PMP3 (Peroxisomal Membrane Protein 3) ELISA Kit
E-EL-GP0161	Guinea pig RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
E-EL-GP0162	Guinea pig CX37 (Connexin 37) ELISA Kit
E-EL-GP0163	Guinea pig PTD (Pentosidine) ELISA Kit
E-EL-GP0165	Guinea pig IGFBP-2 (Insulin Like Growth Factor Binding Protein 2) ELISA Kit
E-EL-GP0166	Guinea pig WARS (Tryptophanyl tRNA Synthetase) ELISA Kit
E-EL-GP0167	Guinea pig COL7 (Collagen Type VII) ELISA Kit
E-EL-GP0168	Guinea pig HP (Haptoglobin) ELISA Kit
E-EL-GP0169	Guinea pig LH (Luteinizing Hormone) ELISA Kit
E-EL-GP0170	Guinea pig PDXK (Pyridoxal Kinase) ELISA Kit
E-EL-GP0171	Guinea pig OPN (Osteopontin) ELISA Kit
E-EL-GP0172	Guinea pig TDT (Terminal Deoxynucleotidyl Transferase) ELISA Kit
E-EL-GP0173	Guinea pig RUNX2 (Runt-Related Transcription Factor 2) ELISA Kit
E-EL-GP0174	Guinea pig PC/CPG (Choline Phosphoglyceride) ELISA Kit
E-EL-GP0175	Guinea pig HNRPU (Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit

Cat.NO	Product Name
E-EL-GP0176	Guinea pig PKCθ (Protein Kinase C Theta)ELISA Kit
E-EL-GP0177	Guinea pig LIF (Leukemia Inhibitory Factor) ELISA Kit
E-EL-GP0178	Guinea pig ET-1 (Endothelin 1) ELISA Kit
E-EL-GP0179	Guinea pig RAPGEF1 (Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit
E-EL-GP0180	Guinea pig EM2 (Endomorphin 2) ELISA Kit
E-EL-GP0181	Guinea pig CNR1 (Cannabinoid Receptor 1, Brain) ELISA Kit
E-EL-GP0182	Guinea pig C4 (Complement Component 4) ELISA Kit
E-EL-GP0183	Guinea pig NCAM/CD56 (Neural Cell Adhesion Molecule) ELISA Kit
E-EL-GP0185	Guinea pig D2R/DRD2 (Dopamine Receptor D2) ELISA Kit
E-EL-GP0186	Guinea pig PRKG1 (Protein Kinase, cGMP Dependent Type I) ELISA Kit
E-EL-GP0187	Guinea pig TFR/CD71 (Transferrin Receptor) ELISA Kit
E-EL-GP0188	Guinea pig PKC (Protein Kinase C) ELISA Kit
E-EL-GP0189	Guinea pig TAT (Thrombin/Antithrombin Complex) ELISA Kit
E-EL-GP0190	Guinea pig PMN Elastase (Polymorphonuclear Elastase) ELISA Kit
E-EL-GP0191	Guinea pig SREBP (Sterol Regulatory Element Binding Proteins) ELISA Kit
E-EL-GP0192	Guinea pig PAFAHIB3 (Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit
E-EL-GP0193	Guinea pig DCLK1 (Doublecortin Like Kinase 1) ELISA Kit
E-EL-GP0194	Guinea pig Trx (Thioredoxin) ELISA Kit
E-EL-GP0195	Guinea pig GAL8 (Galectin 8) ELISA Kit
E-EL-GP0196	Guinea pig sFMC (Soluble Fibrin Monomer Complex) ELISA Kit

Cat.NO	Product Name
E-EL-GP0197	Guinea pig CD163 (Cluster of Differentiation 163) ELISA Kit
E-EL-GP0198	Guinea pig MCRS1 (Microspherule Protein 1) ELISA Kit
E-EL-GP0199	Guinea pig BANP (BTG3 Associated Nuclear Protein) ELISA Kit
E-EL-GP0200	Guinea pig EPYC (Epiphycan) ELISA Kit
E-EL-GP0201	Guinea pig OSMR (Oncostatin M Receptor) ELISA Kit
E-EL-GP0202	Guinea pig CASP4 (Caspase 4) ELISA Kit
E-EL-GP0203	Guinea pig DBP (DNA Binding Protein) ELISA Kit
E-EL-GP0204	Guinea pig ATF-5 (Activating Transcription Factor 5) ELISA Kit
E-EL-GP0205	Guinea pig EPHA1 (Ephrin Type A Receptor 1) ELISA Kit
E-EL-GP0206	Guinea pig IL-18 (Interleukin 18) ELISA Kit
E-EL-GP0207	Guinea pig NKB (Neurokinin B) ELISA Kit
E-EL-GP0208	Guinea pig CTAPIII (Connective Tissue Activating Peptide III) ELISA Kit
E-EL-GP0209	Guinea pig ASPN (Asporin) ELISA Kit
E-EL-GP0210	Guinea pig TAF13 (TATA Box Binding Protein Associated Factor 13) ELISA Kit
E-EL-GP0211	Guinea pig VEGFR1/Flt1 (Vascuoar Endothelial Growth Factor Receptor 1) ELISA Kit
E-EL-GP0213	Guinea pig PR3 (Proteinase 3) ELISA Kit
E-EL-GP0214	Guinea pig PAR2 (Protease Activated Receptor 2) ELISA Kit
E-EL-GP0215	Guinea pig Grb14 (Growth Factor Receptor Bound Protein 14) ELISA Kit
E-EL-GP0216	Guinea pig GZMK (Granzyme K) ELISA Kit
E-EL-GP0217	Guinea pig CC16 (Clara Cell 16KD Protein) ELISA Kit

Cat.NO	Product Name
E-EL-GP0219	Guinea pig CFI (Complement Factor I) ELISA Kit
E-EL-GP0220	Guinea pig MUC3 (Mucin 3) ELISA Kit
E-EL-GP0221	Guinea pig FLNC (Filamin C, Gamma) ELISA Kit
E-EL-GP0222	Guinea pig COLEC11 (Collectin Sub-Family Member 11) ELISA Kit
E-EL-GP0223	Guinea pig CTNN α 1 (Catenin, Alpha 1) ELISA Kit
E-EL-GP0224	Guinea pig GH (Growth Hormone) ELISA Kit
E-EL-GP0225	Guinea pig α -CTx (Alpha Crosslaps) ELISA Kit
E-EL-GP0226	Guinea pig EAR3 (Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit
E-EL-GP0228	Guinea pig ASB (Arylsulfatase B) ELISA Kit
E-EL-GP0229	Guinea pig PTHrP (Parathyroid Hormone Related Protein) ELISA Kit
E-EL-GP0230	Guinea pig VF (Visfatin) ELISA Kit
E-EL-GP0231	Guinea pig CNTF (Ciliary Neurotrophic Factor) ELISA Kit
E-EL-GP0232	Guinea pig MRC1 (Mannose Receptor C Type 1 Like Protein 1) ELISA Kit
E-EL-GP0233	Guinea pig HHIP (Hedgehog Interacting Protein) ELISA Kit
E-EL-GP0234	Guinea pig ADCY3 (Adenylate Cyclase 3) ELISA Kit
E-EL-GP0235	Guinea pig ApoC1 (Apolipoprotein C1) ELISA Kit
E-EL-GP0236	Guinea pig CCKAR (Cholecystokinin A Receptor) ELISA Kit
E-EL-GP0237	Guinea pig Collagenase II ELISA Kit
E-EL-GP0238	Guinea pig ADRP (Adipose Differentiation Related Protein) ELISA Kit
E-EL-GP0239	Guinea pig CLDN11 (Claudin 11) ELISA Kit

Cat.NO	Product Name
E-EL-GP0240	Guinea pig COL1 α 1 (Collagen Type I Alpha 1) ELISA Kit
E-EL-GP0241	Guinea pig GTF2A1 (General Transcription Factor IIA, Polypeptide 1) ELISA Kit
E-EL-GP0242	Guinea pig ATF-1 (Activating Transcription Factor 1) ELISA Kit
E-EL-GP0243	Guinea pig TMS β 4 (Thymosin Beta 4) ELISA Kit
E-EL-GP0244	Guinea pig PGD2S (Prostaglandin D2 Synthase, Brain) ELISA Kit
E-EL-GP0245	Guinea pig β -actin (Actin Beta) ELISA Kit
E-EL-GP0246	Guinea pig REG3 α (Regenerating Islet Derived Protein 3 Alpha) ELISA Kit
E-EL-GP0247	Guinea pig Flt3 (FMS Like Tyrosine Kinase 3) ELISA Kit
E-EL-GP0248	Guinea pig BCAN (Brevican) ELISA Kit
E-EL-GP0249	Guinea pig MMP-24 (Matrix Metalloproteinase 24) ELISA Kit
E-EL-GP0250	Guinea pig MMP-25 (Matrix Metalloproteinase 25) ELISA Kit
E-EL-GP0251	Guinea pig CAPA (Cytotoxin Associated Protein A) ELISA Kit
E-EL-GP0252	Guinea pig VEGF-D (Vascular Endothelial Growth Factor D) ELISA Kit
E-EL-GP0253	Guinea pig CYP7A1 (Cytochrome P450 7A1) ELISA Kit
E-EL-GP0254	Guinea pig CFH (Complement Factor H) ELISA Kit
E-EL-GP0255	Guinea pig TPA (Tissue Polypeptide Antigen) ELISA Kit
E-EL-GP0256	Guinea pig HJV/HFE2 (Hemojuvelin/Hemochromatosis Type 2) ELISA Kit
E-EL-GP0257	Guinea pig SDF2L1 (Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
E-EL-GP0258	Guinea pig GAP (GTPase Activating Protein) ELISA Kit
E-EL-GP0259	Guinea pig Arg (Arginase) ELISA Kit

Cat.NO	Product Name
E-EL-GP0260	Guinea pig FABP5 (Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
E-EL-GP0261	Guinea pig ApoB100 (Apolipoprotein B100) ELISA Kit
E-EL-GP0262	Guinea pig bTC (Betacellulin) ELISA Kit
E-EL-GP0264	Guinea pig FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
E-EL-GP0266	Guinea pig AMA IgA (Anti-Myelin antibody IgA) ELISA Kit
E-EL-GP0267	Guinea pig MIP-2 (Macrophage Inflammatory Protein 2) ELISA Kit
E-EL-GP0268	Guinea pig TSC22 (Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
E-EL-GP0270	Guinea pig 6-K-PG (6-Keto-prostaglandin) ELISA Kit
E-EL-GP0271	Guinea pig SAA2 (Serum Amyloid A2) ELISA Kit
E-EL-GP0272	Guinea pig TF (Tissue Factor) ELISA Kit
E-EL-GP0273	Guinea pig GLTP (Glycolipid Transfer Protein) ELISA Kit
E-EL-GP0275	Guinea pig C3a (Complement Component 3a) ELISA Kit
E-EL-GP0276	Guinea pig NKA (Neurokinin A) ELISA Kit
E-EL-GP0277	Guinea pig N-MID-OT (N-MID Osteocalcin) ELISA Kit
E-EL-GP0278	Guinea pig CENPE (Centromere Protein E) ELISA Kit
E-EL-GP0279	Guinea pig USP8 (Ubiquitin Specific Peptidase 8) ELISA Kit
E-EL-GP0280	Guinea pig α 2-PI (α 2-Plasmin Inhibitor) ELISA Kit
E-EL-GP0281	Guinea pig DLAT (Dihydrolipoyl Transacetylase) ELISA Kit
E-EL-GP0282	Guinea pig G6PD (Glucose 6 Phosphate Dehydrogenase) ELISA Kit
E-EL-GP0283	Guinea pig FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit

Cat.NO	Product Name
E-EL-GP0284	Guinea pig CRT (Calreticulin) ELISA Kit
E-EL-GP0285	Guinea pig CTR (Calcitonin Receptor) ELISA Kit
E-EL-GP0286	Guinea pig Lptn/XCL1 (Lymphotactin) ELISA Kit
E-EL-GP0287	Guinea pig PON1 (Paraoxonase 1) ELISA Kit
E-EL-GP0289	Guinea pig CASP14 (Caspase 14) ELISA Kit
E-EL-GP0290	Guinea pig H-Cadherin (Cadherin, Heart) ELISA Kit
E-EL-GP0291	Guinea pig KAF (Keratinocyte Autocrine Factor) ELISA Kit
E-EL-GP0294	Guinea pig NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) ELISA Kit
E-EL-GP0295	Guinea pig IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit
E-EL-GP0296	Guinea pig ERK3 (Extracellular Signal Regulated Kinase 3) ELISA Kit
E-EL-GP0297	Guinea pig CFP (Complement Factor P) ELISA Kit
E-EL-GP0298	Guinea pig CENPI (Centromere Protein I) ELISA Kit
E-EL-GP0299	Guinea pig SPARCL1 (SPARC Like Protein 1) ELISA Kit
E-EL-GP0300	Guinea pig BAFFR (B Cell Activating Factor Receptor) ELISA Kit
E-EL-GP0301	Guinea pig CKLFSF6 (Chemokine Like Factor Superfamily 6) ELISA Kit
E-EL-GP0302	Guinea pig ABCG2 (ATP Binding Cassette Transporter G2) ELISA Kit
E-EL-GP0303	Guinea pig CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit
E-EL-GP0304	Guinea pig NAP3 (Neutrophil Activating Protein 3) ELISA Kit
E-EL-GP0305	Guinea pig DKK4 (Dickkopf Related Protein 4) ELISA Kit
E-EL-GP0306	Guinea pig PPAR- γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit

Cat.NO	Product Name
E-EL-GP0307	Guinea pig DEF α 5 (Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
E-EL-GP0308	Guinea pig GPX1 (Glutathione Peroxidase 1) ELISA Kit
E-EL-GP0309	Guinea pig PAIg (Platelet Associated Immunoglobulins) ELISA Kit
E-EL-GP0310	Guinea pig SDF-1/CXCL12 (Stromal Cell Derived Factor 1) ELISA Kit
E-EL-GP0311	Guinea pig THR α (Thyroid Hormone Receptor Alpha) ELISA Kit
E-EL-GP0312	Guinea pig PGF (Placental Growth Factor) ELISA Kit
E-EL-GP0313	Guinea pig COL1 α 2 (Collagen Type I Alpha 2) ELISA Kit
E-EL-GP0315	Guinea pig PDGF-BB (Platelet Derived Growth Factor-BB) ELISA Kit
E-EL-GP0316	Guinea pig Slit2 (Slit Homolog 2) ELISA Kit
E-EL-GP0318	Guinea pig FRS2 (Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit
E-EL-GP0319	Guinea pig BTA (Bladder Tumor Antigen) ELISA Kit
E-EL-GP0321	Guinea pig VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
E-EL-GP0322	Guinea pig FGFR4 (Fibroblast Growth Factor Receptor 4) ELISA Kit
E-EL-GP0324	Guinea pig TPS (Tryptase) ELISA Kit
E-EL-GP0325	Guinea pig CLDN14 (Claudin 14) ELISA Kit
E-EL-GP0326	Guinea pig CKLFSF5 (Chemokine Like Factor Superfamily 5) ELISA Kit
E-EL-GP0327	Guinea pig PLA1 (Phospholipase A1) ELISA Kit
E-EL-GP0329	Guinea pig SPJ (Surfactant Associated Protein J) ELISA Kit
E-EL-GP0330	Guinea pig PYGM (Glycogen Phosphorylase, Muscle) ELISA Kit
E-EL-GP0331	Guinea pig PRB2 (Basic Salivary Proline Rich Protein 2) ELISA Kit

Cat.NO	Product Name
E-EL-GP0332	Guinea pig CK-17/KRT17 (Cytokeratin 17) ELISA Kit
E-EL-GP0333	Guinea pig TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
E-EL-GP0334	Guinea pig DPP6 (Dipeptidyl Peptidase VI) ELISA Kit
E-EL-GP0335	Guinea pig tPA (Plasminogen Activator, Tissue) ELISA Kit
E-EL-GP0336	Guinea pig GSTθ1/GSTt1 (Glutathione S Transferase Theta 1) ELISA Kit
E-EL-GP0337	Guinea pig C5a (Complement Component 5a) ELISA Kit
E-EL-GP0338	Guinea pig Fcγ (Fc Fragment of IgG) ELISA Kit
E-EL-GP0340	Guinea pig AMY2 (Amylase Alpha 2, Pancreatic) ELISA Kit
E-EL-GP0341	Guinea pig ST3 (Stromelysin-3) ELISA Kit
E-EL-GP0342	Guinea pig ACTH (Adrenocorticotrophic Hormone) ELISA Kit
E-EL-GP0343	Guinea pig IGF-2 (Insulin Like Growth Factor 2) ELISA Kit
E-EL-GP0346	Guinea pig TNR (Tenascin R) ELISA Kit
E-EL-GP0347	Guinea pig LIPC (Lipase, Hepatic) ELISA Kit
E-EL-GP0348	Guinea pig PKD2 (Protein Kinase D2) ELISA Kit
E-EL-GP0349	Guinea pig C5 (Complement Component 5) ELISA Kit
E-EL-GP0350	Guinea pig CA153 (Carbohydrate Antigen 153) ELISA Kit
E-EL-GP0351	Guinea pig CKMT1B (Creatine Kinase, Mitochondrial 1B) ELISA Kit
E-EL-GP0352	Guinea pig HSP-60/HSPD1 (Heat Shock Protein 60) ELISA Kit
E-EL-GP0353	Guinea pig ES (Endostatin) ELISA Kit
E-EL-GP0354	Guinea pig Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit

Cat.NO	Product Name
E-EL-GP0356	Guinea pig TGF- β 3 (Transforming Growth Factor β 3) ELISA Kit
E-EL-GP0357	Guinea pig OCT (Ornithine Carbamoyl Transferase) ELISA Kit
E-EL-GP0358	Guinea pig PRL (Prolactin) ELISA Kit
E-EL-GP0359	Guinea pig LCN1 (Lipocalin 1) ELISA Kit
E-EL-GP0360	Guinea pig SCTR (Secretin Receptor) ELISA Kit
E-EL-GP0361	Guinea pig Grb7 (Growth Factor Receptor Bound Protein 7) ELISA Kit
E-EL-GP0362	Guinea pig CRN (Corin) ELISA Kit
E-EL-GP0363	Guinea pig LTA4H (Leukotriene A4 Hydrolase) ELISA Kit
E-EL-GP0364	Guinea pig MAGE (Melanoma-associated Antigen) ELISA Kit
E-EL-GP0365	Guinea pig NUP188 (Nucleoporin 188kDa) ELISA Kit
E-EL-GP0367	Guinea pig γ GT1 (Gamma Glutamyltransferase 1) ELISA Kit
E-EL-GP0368	Guinea pig DCX (Doublecortin) ELISA Kit
E-EL-GP0369	Guinea pig TSH (Thyroid Stimulating Hormone) ELISA Kit
E-EL-GP0370	Guinea pig IL-10 (Interleukin-10) ELISA Kit
E-EL-GP0371	Guinea pig ANG (Angiostatin) ELISA Kit
E-EL-GP0372	Guinea pig MMSAM (Melanoma Metastasis Surface Adhesion Molecule) ELISA Kit
E-EL-GP0373	Guinea pig ANG-R-Tie2 (Angiopoietin Receptor Tie2) ELISA Kit
E-EL-GP0374	Guinea pig CG α (Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit
E-EL-GP0375	Guinea pig CETP (Cholesteryl Ester Transfer Protein) ELISA Kit
E-EL-GP0376	Guinea pig CLU (Clusterin) ELISA Kit

Cat.NO	Product Name
E-EL-GP0377	Guinea pig Grb2 (Growth Factor Receptor Bound Protein 2) ELISA Kit
E-EL-GP0378	Guinea pig EIPS2 (Estrogen Induced Protein pS2) ELISA Kit
E-EL-GP0379	Guinea pig FASN (Fatty Acid Synthase) ELISA Kit
E-EL-GP0380	Guinea pig P III CP (Procollagen III C-Terminal ProPeptide) ELISA Kit
E-EL-GP0381	Guinea pig CYSLTR2 (Cysteinyl Leukotriene Receptor 2) ELISA Kit
E-EL-GP0382	Guinea pig cTnT/TNNT2 (Troponin T Type 2, Cardiac) ELISA Kit
E-EL-GP0383	Guinea pig GLP-1 (Glucagon Like Peptide 1) ELISA Kit
E-EL-GP0384	Guinea pig PLSCR1 (Phospholipid Scramblase 1) ELISA Kit
E-EL-GP0385	Guinea pig UBE3A (Ubiquitin Protein Ligase E3A) ELISA Kit
E-EL-GP0386	Guinea pig EloA (Elongin A) ELISA Kit
E-EL-GP0387	Guinea pig EAR2 (Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit
E-EL-GP0388	Guinea pig P-GP (Permeability Glycoprotein) ELISA Kit
E-EL-GP0389	Guinea pig NAIP (Neuronal Apoptosis Inhibitory Protein) ELISA Kit
E-EL-GP0390	Guinea pig PIK3C α (Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit
E-EL-GP0392	Guinea pig PCK1 (Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
E-EL-GP0393	Guinea pig ANG (Angiogenin) ELISA Kit
E-EL-GP0394	Guinea pig LDL (Low Density Lipoprotein) ELISA Kit
E-EL-GP0395	Guinea pig Fc γ R3A (Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit
E-EL-GP0396	Guinea pig SPT β N4 (Spectrin Beta, Non Erythrocytic 4) ELISA Kit
E-EL-GP0398	Guinea pig cADPRH/CD38 (Cyclic ADP Hydrolase) ELISA Kit

Cat.NO	Product Name
E-EL-GP0399	Guinea pig C5c (Complement C5 Convertase) ELISA Kit
E-EL-GP0400	Guinea pig TRH (Thyrotropin Releasing Hormone) ELISA Kit
E-EL-GP0401	Guinea pig CNTN1 (Contactin 1) ELISA Kit
E-EL-GP0402	Guinea pig HSPG (Heparan Sulfate Proteoglycan) ELISA Kit
E-EL-GP0403	Guinea pig BFP (Brain Finger Protein) ELISA Kit
E-EL-GP0404	Guinea pig MCP-4/CCL13 (Monocyte Chemotactic Protein 4) ELISA Kit
E-EL-GP0405	Guinea pig CDK4 (Cyclin Dependent Kinase 4) ELISA Kit
E-EL-GP0406	Guinea pig MAPK (Mitogen Activated Protein Kinase) ELISA Kit
E-EL-GP0407	Guinea pig CLEC4E (C-Type Lectin Domain Family 4 Member E) ELISA Kit
E-EL-GP0408	Guinea pig CDKN3 (Cyclin Dependent Kinase Inhibitor 3) ELISA Kit
E-EL-GP0409	Guinea pig AGT (Angiotensinogen) ELISA Kit
E-EL-GP0410	Guinea pig BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit
E-EL-GP0411	Guinea pig NPYR2 (Neuropeptide Y Receptor Y2) ELISA Kit
E-EL-GP0413	Guinea pig UCHL1 (Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
E-EL-GP0414	Guinea pig DAG1 (Dystrophin Associated Glycoprotein 1) ELISA Kit
E-EL-GP0416	Guinea pig LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit
E-EL-GP0417	Guinea pig PDP (Pyruvate Dehydrogenase Phosphatase) ELISA Kit
E-EL-GP0418	Guinea pig PKR1 (Prokineticin Receptor 1) ELISA Kit
E-EL-GP0419	Guinea pig PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit
E-EL-GP0420	Guinea pig CA2 (Carbonic Anhydrase II) ELISA Kit

Cat.NO	Product Name
E-EL-GP0421	Guinea pig β -TG (β -Thromboglobulin) ELISA Kit
E-EL-GP0422	Guinea pig SPOP (Speckle Type POZ Protein) ELISA Kit
E-EL-GP0423	Guinea pig IL-24 (Interleukin 24) ELISA Kit
E-EL-GP0424	Guinea pig TG (Thyroglobulin) ELISA Kit
E-EL-GP0426	Guinea pig Rac-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit
E-EL-GP0428	Guinea pig SLPI (Secretory Leukocyte Protease Inhibitor) ELISA Kit
E-EL-GP0429	Guinea pig AQP-3 (Aquaporin 3) ELISA Kit
E-EL-GP0430	Guinea pig CDK5 (Cyclin Dependent Kinase 5) ELISA Kit
E-EL-GP0431	Guinea pig MMP-3 (Matrix Metalloproteinase 3) ELISA Kit
E-EL-GP0433	Guinea pig TSLP (Thymic Stromal Lymphopoietin) ELISA Kit
E-EL-GP0434	Guinea pig GHAb (Anti-Growth Hormone Antibody) ELISA Kit
E-EL-GP0435	Guinea pig NUP153 (Nucleoporin 153kDa) ELISA Kit
E-EL-GP0436	Guinea pig GLUT2 (Glucose Transporter 2) ELISA Kit
E-EL-GP0437	Guinea pig BMG/ β 2-MG (Beta-2-Microglobulin) ELISA Kit
E-EL-GP0438	Guinea pig SEMA3B (Semaphorin 3B) ELISA Kit
E-EL-GP0439	Guinea pig HNRPA2B1 (Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit
E-EL-GP0440	Guinea pig GDF15 (Growth Differentiation Factor 15) ELISA Kit
E-EL-GP0442	Guinea pig EMMPRIN/CD147 (Extracellular Matrix Metalloproteinase Inducer) ELISA Kit
E-EL-GP0443	Guinea pig PACAPR (Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit
E-EL-GP0444	Guinea pig OGP (Osteogenic Growth Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-GP0445	Guinea pig LFA-3 (Lymphocyte Function Associated Antigen 3) ELISA Kit
E-EL-GP0446	Guinea pig PLC (Phospholipase C) ELISA Kit
E-EL-GP0447	Guinea pig PS (Plasma Anticoagulant Protein S) ELISA Kit
E-EL-GP0448	Guinea pig NF- κ B p65 (Nuclear Factor- κ B p65) ELISA Kit
E-EL-GP0449	Guinea pig ANGPTL4 (Angiopoietin Like Protein 4) ELISA Kit
E-EL-GP0450	Guinea pig ASGR1 (Asialoglycoprotein Receptor 1) ELISA Kit
E-EL-GP0451	Guinea pig PS (Phosphatidylserine) ELISA Kit
E-EL-GP0452	Guinea pig α Actin (Alpha-Actin) ELISA Kit
E-EL-GP0453	Guinea pig EphA2 (Tyrosine protein kinase receptor A2) ELISA Kit
E-EL-GP0454	Guinea pig CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
E-EL-GP0455	Guinea pig BMP-15 (Bone Morphogenetic Protein 15) ELISA Kit
E-EL-GP0457	Guinea pig SPC (Pulmonary Surfactant Associated Protein C) ELISA Kit
E-EL-GP0458	Guinea pig ACA-IgG (Anti-Cardiolipin Antibody IgG) ELISA Kit
E-EL-GP0459	Guinea pig D2D (D-Dimer) ELISA Kit
E-EL-GP0461	Guinea pig S100A12 (S100 Calcium Binding Protein A12) ELISA Kit
E-EL-GP0462	Guinea pig CYP1A2 (Cytochrome P450 1A2) ELISA Kit
E-EL-GP0463	Guinea pig DKK2 (Dickkopf Related Protein 2) ELISA Kit
E-EL-GP0464	Guinea pig LIFR (Leukemia Inhibitory Factor Receptor) ELISA Kit
E-EL-GP0465	Guinea pig CNTN4 (Contactin 4) ELISA Kit
E-EL-GP0466	Guinea pig FE (Ferritin) ELISA Kit

Cat.NO	Product Name
E-EL-GP0467	Guinea pig PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit
E-EL-GP0468	Guinea pig FBLN4 (Fibulin 4) ELISA Kit
E-EL-GP0470	Guinea pig MFG-E8 (Milk Fat Globule EGF Factor 8) ELISA Kit
E-EL-GP0471	Guinea pig SHBG (Sex Hormone Binding Globulin) ELISA Kit
E-EL-GP0472	Guinea pig TPP (Thrombus Precursor Protein) ELISA Kit
E-EL-GP0474	Guinea pig LDHB (Lactate Dehydrogenase B) ELISA Kit
E-EL-GP0475	Guinea pig PLP (Parathyroid Hormone-Like Protein) ELISA Kit
E-EL-GP0476	Guinea pig HSP-27/HSPB1 (Heat Shock Protein 27) ELISA Kit
E-EL-GP0477	Guinea pig NOX1 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
E-EL-GP0478	Guinea pig OBTF1 (Octamer Binding Transcription Factor 1) ELISA Kit
E-EL-GP0479	Guinea pig PTPRN (Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
E-EL-GP0480	Guinea pig MEPE (Matrix Extracellular Phosphoglycoprotein) ELISA Kit
E-EL-GP0481	Guinea pig F XII a (Activated Coagulation Factor XII) ELISA Kit
E-EL-GP0482	Guinea pig ALP (Alkaline Phosphatase) ELISA Kit
E-EL-GP0483	Guinea pig RNASEN (Ribonuclease III, Nuclear) ELISA Kit
E-EL-GP0484	Guinea pig MRP1 (Multidrug Resistance Associated Protein 1) ELISA Kit
E-EL-GP0485	Guinea pig PDGF-AB (Platelet Derived Growth Factor AB) ELISA Kit
E-EL-GP0486	Guinea pig ACA-IgM (Anti-Cardiolipin Antibody IgM) ELISA Kit
E-EL-GP0487	Guinea pig CG (Chorionic Gonadotropin) ELISA Kit
E-EL-GP0488	Guinea pig TGM1 (Transglutaminase 1, Keratinocyte) ELISA Kit

Cat.NO	Product Name
E-EL-GP0489	Guinea pig Try2 (Trypsinogen II) ELISA Kit
E-EL-GP0490	Guinea pig CABIN1 (Calcineurin Binding Protein 1) ELISA Kit
E-EL-GP0491	Guinea pig ADAM9 (A Disintegrin And Metalloprotease 9) ELISA Kit
E-EL-GP0492	Guinea pig ATXN2 (Ataxin 2) ELISA Kit
E-EL-GP0493	Guinea pig PRKA α 1 (Protein Kinase, AMP Activated Alpha 1) ELISA Kit
E-EL-GP0494	Guinea pig MUC13 (Mucin 13) ELISA Kit
E-EL-GP0495	Guinea pig NP (Nandrolone Phenylpropionate) ELISA Kit
E-EL-GP0496	Guinea pig CDNF (Cerebral Dopamine Neurotrophic Factor) ELISA Kit
E-EL-GP0497	Guinea pig CDT (Carbohydrate Deficient Transferrin) ELISA Kit
E-EL-GP0498	Guinea pig CA242 (Pancreatic Carcinoma Markers-CA242) ELISA Kit
E-EL-GP0499	Guinea pig SAP (Serum Amyloid P Component) ELISA Kit
E-EL-GP0500	Guinea pig CD109 (Cluster of Differentiation 109) ELISA Kit
E-EL-GP0501	Guinea pig PDGFR α (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit
E-EL-GP0502	Guinea pig CASP 9 (Caspase 9) ELISA Kit
E-EL-GP0503	Guinea pig PIINP (Procollagen II N-Terminal ProPeptide) ELISA Kit
E-EL-GP0504	Guinea pig PAR3 (Protease Activated Receptor 3) ELISA Kit
E-EL-GP0505	Guinea pig FGFBP1 (Fibroblast Growth Factor Binding Protein 1) ELISA Kit
E-EL-GP0506	Guinea pig NUP98 (Nucleoporin 98kDa) ELISA Kit
E-EL-GP0507	Guinea pig TFF2 (Trefoil Factor 2) ELISA Kit
E-EL-GP0508	Guinea pig GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit

Cat.NO	Product Name
E-EL-GP0509	Guinea pig LN (Laminin) ELISA Kit
E-EL-GP0510	Guinea pig TXLNb (Taxilin Beta) ELISA Kit
E-EL-GP0511	Guinea pig CLA/CTAGE (Cutaneous Lymphocyte Associated Antigen) ELISA Kit
E-EL-GP0512	Guinea pig gp130 (Glucoprotein 130) ELISA Kit
E-EL-GP0513	Guinea pig TPPP (Tubulin Polymerization Promoting Protein) ELISA Kit
E-EL-GP0515	Guinea pig DPP10 (Dipeptidyl Peptidase X) ELISA Kit
E-EL-GP0516	Guinea pig IGFBP-4 (Insulin Like Growth Factor Binding Protein 4) ELISA Kit
E-EL-GP0517	Guinea pig AP 13 (Apelin 13) ELISA Kit
E-EL-GP0518	Guinea pig S100B (S100 Calcium Binding Protein B) ELISA Kit
E-EL-GP0519	Guinea pig GATA3 (GATA Binding Protein 3) ELISA Kit
E-EL-GP0520	Guinea pig TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
E-EL-GP0521	Guinea pig PTH (Parathyroid Hormone) ELISA Kit
E-EL-GP0522	Guinea pig SEP1 (Selenoprotein 1) ELISA Kit
E-EL-GP0523	Guinea pig PAF (Platelet Activating Factor) ELISA Kit
E-EL-GP0524	Guinea pig CXCR3 (Chemokine CXC Receptor 3) ELISA Kit
E-EL-GP0525	Guinea pig P-PKC (Phospho Protein Kinase C) ELISA Kit
E-EL-GP0526	Guinea pig SRP9 (Signal Recognition Particle 9kDa) ELISA Kit
E-EL-GP0527	Guinea pig HA (Hemagglutinin) ELISA Kit
E-EL-GP0528	Guinea pig MYL9 (Myosin Light Chain 9, Regulatory) ELISA Kit
E-EL-GP0530	Guinea pig CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) ELISA Kit

Cat.NO	Product Name
E-EL-GP0531	Guinea pig TPP1 (Tripeptidyl Peptidase I) ELISA Kit
E-EL-GP0532	Guinea pig TAF1 (TATA Box Binding Protein Associated Factor 1) ELISA Kit
E-EL-GP0533	Guinea pig PAP (Plasmin-Antiplasmin Complex) ELISA Kit
E-EL-GP0534	Guinea pig NT-3 (Neurotrophin-3) ELISA Kit
E-EL-GP0535	Guinea pig ANGPTL2 (Angiopoietin Like Protein 2) ELISA Kit
E-EL-GP0536	Guinea pig TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
E-EL-GP0537	Guinea pig NAP (Neutrophil Alkaline Phosphatase) ELISA Kit
E-EL-GP0538	Guinea pig sICAM-1 (Soluble Intercellular Adhesion Molecule-1) ELISA Kit
E-EL-GP0540	Guinea pig GRN (Granulin) ELISA Kit
E-EL-GP0541	Guinea pig GPI (Glucose 6 Phosphate Isomerase) ELISA Kit
E-EL-GP0542	Guinea pig ADH1 (Alcohol Dehydrogenase 1) ELISA Kit
E-EL-GP0543	Guinea pig GUK1 (Guanylate Kinase 1) ELISA Kit
E-EL-GP0545	Guinea pig LEPR (Leptin Receptor) ELISA Kit
E-EL-GP0546	Guinea pig MYD88 (Myeloid Differentiation Factor 88) ELISA Kit
E-EL-GP0547	Guinea pig CPP (Copeptin) ELISA Kit
E-EL-GP0548	Guinea pig CASP10 (Caspase 10) ELISA Kit
E-EL-GP0550	Guinea pig FLOT2 (Flotillin 2) ELISA Kit
E-EL-GP0551	Guinea pig TJP1 (Tight Junction Protein 1) ELISA Kit
E-EL-GP0553	Guinea pig TP/TXA2R (Thromboxane Receptor) ELISA Kit
E-EL-GP0554	Guinea pig MSR1 (Macrophage Scavenger Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-GP0555	Guinea pig PHB (Prohibitin) ELISA Kit
E-EL-GP0556	Guinea pig NEF3 (Neurofilament 3) ELISA Kit
E-EL-GP0557	Guinea pig LP-a (Lipoprotein a) ELISA Kit
E-EL-GP0558	Guinea pig MAO (Monoamine Oxidase) ELISA Kit
E-EL-GP0559	Guinea pig SIGLEC12 (Sialic Acid Binding Ig Like Lectin 12) ELISA Kit
E-EL-GP0560	Guinea pig TNNI1 (Troponin I Type 1, Slow Skeletal) ELISA Kit
E-EL-GP0561	Guinea pig PYGB (Glycogen Phosphorylase, Brain) ELISA Kit
E-EL-GP0562	Guinea pig GT-34 (Gastrin 34) ELISA Kit
E-EL-GP0563	Guinea pig ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit
E-EL-GP0564	Guinea pig CX31 (Connexin 31) ELISA Kit
E-EL-GP0565	Guinea pig ATF-2 (Activating Transcription Factor 2) ELISA Kit
E-EL-GP0566	Guinea pig MYOG (Myogenin) ELISA Kit
E-EL-GP0567	Guinea pig CLL1 (Collectin Liver 1) ELISA Kit
E-EL-GP0568	Guinea pig MMP-8 (Matrix Metalloproteinase 8) ELISA Kit
E-EL-GP0569	Guinea pig GPRC5A (G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit
E-EL-GP0570	Guinea pig sIgA (Secretory Immunoglobulin A) ELISA Kit
E-EL-GP0571	Guinea pig SRF/IF (Skin Reactive Factor/Inflammatory Factor) ELISA Kit
E-EL-GP0572	Guinea pig SELE (E-Selectin) ELISA Kit
E-EL-GP0573	Guinea pig MCP/CD46 (Membrane Cofactor Protein) ELISA Kit
E-EL-GP0574	Guinea pig ANG III (Angiotensin III) ELISA Kit

Cat.NO	Product Name
E-EL-GP0575	Guinea pig Extein ELISA Kit
E-EL-GP0577	Guinea pig PKP2 (Plakophilin 2) ELISA Kit
E-EL-GP0578	Guinea pig PXR (Pregnane X Receptor) ELISA Kit
E-EL-GP0579	Guinea pig c-myc (c-myc Oncogene Product) ELISA Kit
E-EL-GP0580	Guinea pig GAL12 (Galectin 12) ELISA Kit
E-EL-GP0582	Guinea pig CKMT1A (Creatine Kinase, Mitochondrial 1A) ELISA Kit
E-EL-GP0583	Guinea pig PCDH1 (Protocadherin 1) ELISA Kit
E-EL-GP0584	Guinea pig TXB2 (Thromboxane B2) ELISA Kit
E-EL-GP0585	Guinea pig FGL2 (Fibrinogen Like Protein 2) ELISA Kit
E-EL-GP0586	Guinea pig AGGF1 (Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit
E-EL-GP0587	Guinea pig DSG-E1 (Desmogleins 1) ELISA Kit
E-EL-GP0588	Guinea pig CR2/CD21 (Complement Receptor 2) ELISA Kit
E-EL-GP0590	Guinea pig HB (Hemoglobin) ELISA Kit
E-EL-GP0591	Guinea pig SCT (Secretin) ELISA Kit
E-EL-GP0592	Guinea pig CD23/FcεRII (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
E-EL-GP0593	Guinea pig CYP1B1 (Cytochrome P450 1B1) ELISA Kit
E-EL-GP0594	Guinea pig LTC4 (Leukotriene C4) ELISA Kit
E-EL-GP0595	Guinea pig PDGFD (Platelet Derived Growth Factor D) ELISA Kit
E-EL-GP0596	Guinea pig CDX2 (Caudal Type Homeobox Transcription Factor 2) ELISA Kit
E-EL-GP0597	Guinea pig COL6 (Collagen Type VI) ELISA Kit

Cat.NO	Product Name
E-EL-GP0598	Guinea pig DCN (Decorin) ELISA Kit
E-EL-GP0599	Guinea pig TLR2 (Toll-Like Receptor 2) ELISA Kit
E-EL-GP0601	Guinea pig β 2-GP1 Ab IgG (β 2-Glycoprotein 1 Antibody IgG) ELISA Kit
E-EL-GP0602	Guinea pig TMS β 10 (Thymosin Beta 10) ELISA Kit
E-EL-GP0603	Guinea pig PLA2G4D (Phospholipase A2, Group IV D) ELISA Kit
E-EL-GP0604	Guinea pig EPO (Erythropoietin) ELISA Kit
E-EL-GP0605	Guinea pig NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
E-EL-GP0606	Guinea pig OLA β (Anti-Oxidized Low Density Lipoprotein Antibody) ELISA Kit
E-EL-GP0607	Guinea pig PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit
E-EL-GP0608	Guinea pig NGB (Neuroglobin) ELISA Kit
E-EL-GP0609	Guinea pig β -EP (Beta-Endorphin) ELISA Kit
E-EL-GP0610	Guinea pig Smad7/MADH7 (Mothers Against Decapentaplegic Homolog 7) ELISA Kit
E-EL-GP0611	Guinea pig TDGF1 (Teratocarcinoma Derived Growth Factor 1) ELISA Kit
E-EL-GP0612	Guinea pig COL9 α 3 (Collagen Type IX Alpha 3) ELISA Kit
E-EL-GP0613	Guinea pig SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-GP0614	Guinea pig TFPI (Tissue Factor Pathway Inhibitor) ELISA Kit
E-EL-GP0615	Guinea pig VCL (Vinculin) ELISA Kit
E-EL-GP0618	Guinea pig BMPR1A (Bone Morphogenetic Protein Receptor, type IA) ELISA Kit
E-EL-GP0619	Guinea pig ANXA2 (Annexin A2) ELISA Kit
E-EL-GP0620	Guinea pig C1q (Complement 1q) ELISA Kit

Cat.NO	Product Name
E-EL-GP0621	Guinea pig LOXL4 (Lysyl Oxidase Like Protein 4) ELISA Kit
E-EL-GP0622	Guinea pig DAT (Dopamine Transporter) ELISA Kit
E-EL-GP0623	Guinea pig LILRB4 (Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit
E-EL-GP0624	Guinea pig LZM (Lysozyme) ELISA Kit
E-EL-GP0625	Guinea pig CKLFSF7 (Chemokine Like Factor Superfamily 7) ELISA Kit
E-EL-GP0626	Guinea pig EPHA3 (Ephrin Type A Receptor 3) ELISA Kit
E-EL-GP0627	Guinea pig BACE1 (Beta-site APP Cleaving Enzyme 1) ELISA Kit
E-EL-GP0628	Guinea pig PLG (Plasminogen) ELISA Kit
E-EL-GP0629	Guinea pig ALT (Alanine Transaminase) ELISA Kit
E-EL-GP0630	Guinea pig GTC (General Transcription Complex) ELISA Kit
E-EL-GP0631	Guinea pig MYL1 (Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit
E-EL-GP0632	Guinea pig OIT3 (Oncoprotein Induced Transcript 3) ELISA Kit
E-EL-GP0633	Guinea pig CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit
E-EL-GP0634	Guinea pig GPX3 (Glutathione Peroxidase 3, Plasma) ELISA Kit
E-EL-GP0635	Guinea pig Collagenase I ELISA Kit
E-EL-GP0636	Guinea pig TACR3 (Tachykinin Receptor 3) ELISA Kit
E-EL-GP0637	Guinea pig MUC20 (Mucin 20) ELISA Kit
E-EL-GP0638	Guinea pig PTGES2 (Prostaglandin E Synthase 2) ELISA Kit
E-EL-GP0639	Guinea pig TNC (Tenascin C) ELISA Kit
E-EL-GP0640	Guinea pig GJβ1 (Gap Junction Protein Beta 1) ELISA Kit

Cat.NO	Product Name
E-EL-GP0641	Guinea pig TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
E-EL-GP0643	Guinea pig C/EBP α (CCAAT/Enhancer Binding Protein Alpha) ELISA Kit
E-EL-GP0645	Guinea pig PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
E-EL-GP0646	Guinea pig ADIPOR1 (Adiponectin Receptor 1) ELISA Kit
E-EL-GP0647	Guinea pig FOXP3 (Forkhead Box Protein P3) ELISA Kit
E-EL-GP0648	Guinea pig FBG (Fibrinogen) ELISA Kit
E-EL-GP0649	Guinea pig TBP2 (Thioredoxin Binding Protein 2) ELISA Kit
E-EL-GP0650	Guinea pig CNX (Calnexin) ELISA Kit
E-EL-GP0652	Guinea pig CS (Chondroitin Sulfate) ELISA Kit
E-EL-GP0653	Guinea pig GEF (Guanylate Exchange Factor) ELISA Kit
E-EL-GP0654	Guinea pig CCL4L1 (Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit
E-EL-GP0655	Guinea pig ACE (Angiotensin I Converting Enzyme) ELISA Kit
E-EL-GP0656	Guinea pig PTX 3/TSG-14 (Pentraxin 3, Long) ELISA Kit
E-EL-GP0657	Guinea pig TGM2 (Transglutaminase 2, Tissue) ELISA Kit
E-EL-GP0658	Guinea pig CK-LMW (Cytokeratin Low Molecular Weight) ELISA Kit
E-EL-GP0659	Guinea pig CD30L (Cluster of Differentiation 30 Lligand) ELISA Kit
E-EL-GP0660	Guinea pig PTN (Pleiotrophin) ELISA Kit
E-EL-GP0661	Guinea pig ACD (Adrenocortical Dysplasia Homolog) ELISA Kit
E-EL-GP0663	Guinea pig snRNP/Sm (Small Nuclear Ribonucleoprotein) ELISA Kit
E-EL-GP0664	Guinea pig MK (Midkine) ELISA Kit

Cat.NO	Product Name
E-EL-GP0667	Guinea pig ACA-IgA (Anti-Cardiolipin Antibody IgA) ELISA Kit
E-EL-GP0668	Guinea pig VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit
E-EL-GP0669	Guinea pig AQP-4 (Aquaporin 4) ELISA Kit
E-EL-GP0670	Guinea pig TSPAN30 (Tetraspanin 30) ELISA Kit
E-EL-GP0671	Guinea pig TNPO2 (Transportin 2) ELISA Kit
E-EL-GP0673	Guinea pig DPYSL2 (Dihydropyrimidinase Like Protein 2) ELISA Kit
E-EL-GP0674	Guinea pig RI (Ribonuclease Inhibitor) ELISA Kit
E-EL-GP0675	Guinea pig PI3K (Phosphoinositide-3-Kinase) ELISA Kit
E-EL-GP0676	Guinea pig PKM2 (Pyruvate Kinase, Muscle) ELISA Kit
E-EL-GP0677	Guinea pig CKB (Creatine Kinase, Brain) ELISA Kit
E-EL-GP0678	Guinea pig TACR2 (Tachykinin Receptor 2) ELISA Kit
E-EL-GP0679	Guinea pig GNβ1 (G Protein Beta 1) ELISA Kit
E-EL-GP0680	Guinea pig MGST1 (Microsomal Glutathione S Transferase 1) ELISA Kit
E-EL-GP0681	Guinea pig MMP-7 (Matrix Metalloproteinase 7) ELISA Kit
E-EL-GP0683	Guinea pig PROS (Protein S) ELISA Kit
E-EL-GP0684	Guinea pig ANP (Atrial Natriuretic Peptide) ELISA Kit
E-EL-GP0685	Guinea pig PUMA (P53 Upregulated Modulator of Apoptosis) ELISA Kit
E-EL-GP0686	Guinea pig TP53 (Tumor Protein p53) ELISA Kit
E-EL-GP0687	Guinea pig RANTES (Regulated On Activation, Normal T-Cell Expressed And Secreted) ELISA Kit
E-EL-GP0688	Guinea pig RAG-2 (Recombination Activating Gene 2) ELISA Kit

Cat.NO	Product Name
E-EL-GP0690	Guinea pig CgB (Chromogranin B) ELISA Kit
E-EL-GP0692	Guinea pig FTL (Ferritin, Light Polypeptide) ELISA Kit
E-EL-GP0693	Guinea pig MMP-2 (Matrix Metalloproteinase 2) ELISA Kit
E-EL-GP0695	Guinea pig POMC (Pro-Opiomelanocortin) ELISA Kit
E-EL-GP0696	Guinea pig ADH (Antidiuretic Hormone) ELISA Kit
E-EL-GP0697	Guinea pig INSR (Insulin Receptor) ELISA Kit
E-EL-GP0698	Guinea pig GSTpi (Glutathione S Transferases Pi) ELISA Kit
E-EL-GP0699	Guinea pig α -Glu (Alpha Glucosidase) ELISA Kit
E-EL-GP0701	Guinea pig fT4 (Free Thyroxine) ELISA Kit
E-EL-GP0702	Guinea pig SEMA7A (Semaphorin 7A) ELISA Kit
E-EL-GP0703	Guinea pig CaSR (Calcium Sensing Receptor) ELISA Kit
E-EL-GP0704	Guinea pig CA724 (Tumor Marker CA724) ELISA Kit
E-EL-GP0705	Guinea pig GnRH (Gonadotropin Releasing Hormone) ELISA Kit
E-EL-GP0706	Guinea pig C3b (Complement Fragment 3b) ELISA Kit
E-EL-GP0707	Guinea pig CVNPF (Cobra Venom Neuronal Protective Factor) ELISA Kit
E-EL-GP0709	Guinea pig CT-1 (Cardiotrophin 1) ELISA Kit
E-EL-GP0710	Guinea pig MT1 (Metallothionein 1) ELISA Kit
E-EL-GP0711	Guinea pig sCD14 (Soluble Cluster of Differentiation14) ELISA Kit
E-EL-GP0712	Guinea pig ANG2 (Angiopoietin 2) ELISA Kit
E-EL-GP0713	Guinea pig LT-E4 (Leukotriene E4) ELISA Kit

Cat.NO	Product Name
E-EL-GP0714	Guinea pig IL-12 p40 (Interleukin 12 p40) ELISA Kit
E-EL-GP0715	Guinea pig MUC17 (Mucin 17) ELISA Kit
E-EL-GP0716	Guinea pig PCP (Procollagen C Peptidase)ELISA Kit
E-EL-GP0717	Guinea pig NUP50 (Nucleoporin 50kDa) ELISA Kit
E-EL-GP0718	Guinea pig TCR (Tcell Receptor) ELISA Kit
E-EL-GP0719	Guinea pig AChE (Acetylcholinesterase) ELISA Kit
E-EL-GP0720	Guinea pig GDF6 (Growth Differentiation Factor 6) ELISA Kit
E-EL-GP0721	Guinea pig CKLFSF8 (Chemokine Like Factor Superfamily 8) ELISA Kit
E-EL-GP0722	Guinea pig MUC1 (Mucin 1) ELISA Kit
E-EL-GP0724	Guinea pig A β 1-42 (Amyloid Beta Peptide 1-42) ELISA Kit
E-EL-GP0725	Guinea pig fDFT1 (Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit
E-EL-GP0726	Guinea pig NT-ProANP (N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit
E-EL-GP0727	Guinea pig NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit
E-EL-GP0728	Guinea pig MAG Ab (Anti-Myelin Associated Glycoprotein Antibody) ELISA Kit
E-EL-GP0729	Guinea pig GTF2B (General Transcription Factor IIB) ELISA Kit
E-EL-GP0730	Guinea pig TSP-1 (Thrombospondin-1) ELISA Kit
E-EL-GP0732	Guinea pig ON (Osteonectin) ELISA Kit
E-EL-GP0733	Guinea pig SOX2 (Sex Determining Region Y Box Protein 2) ELISA Kit
E-EL-GP0734	Guinea pig α HSP (Alpha-Hemoglobin Stabilizing Protein) ELISA Kit
E-EL-GP0735	Guinea pig AGA (Anti-Gliadin Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-GP0736	Guinea pig TTR (Transthyretin) ELISA Kit
E-EL-GP0737	Guinea pig PK2 (Prokineticin 2) ELISA Kit
E-EL-GP0738	Guinea pig CHRM2 (Cholinergic Receptor, Muscarinic 2) ELISA Kit
E-EL-GP0742	Guinea pig CASP1 (Caspase 1) ELISA Kit
E-EL-GP0743	Guinea pig FGF9 (Fibroblast Growth Factor 9) ELISA Kit
E-EL-GP0744	Guinea pig PGF2 α (Prostaglandin F2 α) ELISA Kit
E-EL-GP0745	Guinea pig ALD (Aldolase) ELISA Kit
E-EL-GP0747	Guinea pig MIF (Macrophage Migration Inhibitory Factor) ELISA Kit
E-EL-GP0748	Guinea pig CCND2 (Cyclin-D2) ELISA Kit
E-EL-GP0749	Guinea pig AIRE (Autoimmune Regulator) ELISA Kit
E-EL-GP0750	Guinea pig FLN α (Filamin A, Alpha) ELISA Kit
E-EL-GP0751	Guinea pig FTH (Ferritin, Heavy Polypeptide) ELISA Kit
E-EL-GP0752	Guinea pig TFF3 (Trefoil Factor 3, Intestinal) ELISA Kit
E-EL-GP0753	Guinea pig HPR (Haptoglobin Related Protein) ELISA Kit
E-EL-GP0754	Guinea pig EP2 (Prostaglandin E Receptor 2) ELISA Kit
E-EL-GP0755	Guinea pig MUC5AC (Mucin 5 Subtype AC) ELISA Kit
E-EL-GP0757	Guinea pig PDH β (Pyruvate Dehydrogenase Beta) ELISA Kit
E-EL-GP0758	Guinea pig IL-13 (Interleukin 13) ELISA Kit
E-EL-GP0759	Guinea pig LBP (Lipopolysaccharide Binding Protein) ELISA Kit
E-EL-GP0760	Guinea pig TAB (Thyroxine Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-GP0761	Guinea pig CCND1 (Cyclin-D1) ELISA Kit
E-EL-GP0762	Guinea pig LEP (Leptin) ELISA Kit
E-EL-GP0763	Guinea pig CLC (Charcot Leyden Crystal Protein) ELISA Kit
E-EL-GP0764	Guinea pig CCRL1 (Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
E-EL-GP0767	Guinea pig SEMA4C (Semaphorin 4C) ELISA Kit
E-EL-GP0768	Guinea pig AXIN2 (Axis Inhibition Protein 2) ELISA Kit
E-EL-GP0769	Guinea pig GPR37 (G Protein Coupled Receptor 37) ELISA Kit
E-EL-GP0770	Guinea pig GzmA (Granzyme A) ELISA Kit
E-EL-GP0771	Guinea pig PTGES (Prostaglandin E Synthase, Microsomal) ELISA Kit
E-EL-GP0772	Guinea pig MN (Metanephrine) ELISA Kit
E-EL-GP0773	Guinea pig GATA2 (GATA Binding Protein 2) ELISA Kit
E-EL-GP0774	Guinea pig FABP6 (Fatty Acid Binding Protein 6, Ileal) ELISA Kit
E-EL-GP0775	Guinea pig Cav-1 (Caveolin-1) ELISA Kit
E-EL-GP0776	Guinea pig CA125 (Carbohydrate Antigen 125) ELISA Kit
E-EL-GP0777	Guinea pig NPYR5 (Neuropeptide Y Receptor Y5) ELISA Kit
E-EL-GP0778	Guinea pig ADCY1 (Adenylate Cyclase 1, Brain) ELISA Kit
E-EL-GP0779	Guinea pig CX40 (Connexin 40) ELISA Kit
E-EL-GP0781	Guinea pig CYPA (Cyclophilin A) ELISA Kit
E-EL-GP0782	Guinea pig EMILIN 1 (Elastin Microfibril Interface Located Protein 1) ELISA Kit
E-EL-GP0783	Guinea pig PLCβ1/PI-PLC (Phospholipase C Beta 1, Phosphoinositide Specific) ELISA Kit

Cat.NO	Product Name
E-EL-GP0784	Guinea pig NT-ProBNP (N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
E-EL-GP0785	Guinea pig FBLN5 (Fibulin 5) ELISA Kit
E-EL-GP0786	Guinea pig MMP-10 (Matrix Metalloproteinase 10) ELISA Kit
E-EL-GP0787	Guinea pig PF4 (Platelet Factor 4) ELISA Kit
E-EL-GP0788	Guinea pig GPIIb/IIIa (Platelet Membrane Glycoprotein IIb/IIIa) ELISA Kit
E-EL-GP0789	Guinea pig SSFA2 (Sperm Specific Antigen 2) ELISA Kit
E-EL-GP0790	Guinea pig UPP2 (Uridine Phosphorylase 2) ELISA Kit
E-EL-GP0792	Guinea pig GHRH (Growth Hormone Releasing Hormone) ELISA Kit
E-EL-GP0794	Guinea pig HSF1 (Heat Shock Transcription Factor 1) ELISA Kit
E-EL-GP0795	Guinea pig PIIINP (N-Terminal Procollagen III Propeptide) ELISA Kit
E-EL-GP0796	Guinea pig IGFBP-7 (Insulin Like Growth Factor Binding Protein 7) ELISA Kit
E-EL-GP0797	Guinea pig PSAP (Prosaposin) ELISA Kit
E-EL-GP0798	Guinea pig NUP160 (Nucleoporin 160kDa) ELISA Kit
E-EL-GP0799	Guinea pig PROC (Protein C) ELISA Kit
E-EL-GP0800	Guinea pig PYY (Peptide YY) ELISA Kit
E-EL-GP0802	Guinea pig ERα (Estrogen Receptor Alpha) ELISA Kit
E-EL-GP0803	Guinea pig ANG1 (Angiopoietin 1) ELISA Kit
E-EL-GP0804	Guinea pig FcαR (Fc Fragment of IgA Receptor) ELISA Kit
E-EL-GP0805	Guinea pig ALDOA (Aldolase A, Fructose Bisphosphate) ELISA Kit
E-EL-GP0806	Guinea pig GR (Glutathione Reductase) ELISA Kit

Cat.NO	Product Name
E-EL-GP0807	Guinea pig OBTF4 (Octamer Binding Transcription Factor 4) ELISA Kit
E-EL-GP0808	Guinea pig FcεR II (Receptor II for the Fc Fragment of IgE) ELISA Kit
E-EL-GP0809	Guinea pig DAP6 (Death Associated Protein 6) ELISA Kit
E-EL-GP0810	Guinea pig TFR2 (Transferrin Receptor 2) ELISA Kit
E-EL-GP0811	Guinea pig GH2 (Growth Hormone 2) ELISA Kit
E-EL-GP0812	Guinea pig CENPA (Centromere Protein A) ELISA Kit
E-EL-GP0813	Guinea pig ET (Epitestosterone) ELISA Kit
E-EL-GP0814	Guinea pig PAR1 (Protease Activated Receptor 1) ELISA Kit
E-EL-GP0815	Guinea pig IGF2BP3 (Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit
E-EL-GP0816	Guinea pig MIA1 (Melanoma Inhibitory Activity Protein 1) ELISA Kit
E-EL-GP0817	Guinea pig PLD2 (Phospholipase D2) ELISA Kit
E-EL-GP0819	Guinea pig PIGR/SC (Polymeric Immunoglobulin Receptor/ Membrane Secretory Component) ELISA Kit
E-EL-GP0820	Guinea pig MYC (V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit
E-EL-GP0822	Guinea pig EXT1 (Exostoses 1) ELISA Kit
E-EL-GP0823	Guinea pig MT (Melatonin) ELISA Kit
E-EL-GP0824	Guinea pig MLCK/MYLK (Myosin Light Chain Kinase) ELISA Kit
E-EL-GP0825	Guinea pig FBN3 (Fibrillin 3) ELISA Kit
E-EL-GP0827	Guinea pig cPLA2 (Cytosolic Phospholipase A2) ELISA Kit
E-EL-GP0828	Guinea pig PTGS1 (Prostaglandin Endoperoxide Synthase 1) ELISA Kit
E-EL-GP0830	Guinea pig 15-LO (Arachidonate 15-Lipoxygenase) ELISA Kit

Cat.NO	Product Name
E-EL-GP0832	Guinea pig RXFP1 (Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
E-EL-GP0833	Guinea pig MAP τ (Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
E-EL-GP0834	Guinea pig MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
E-EL-GP0835	Guinea pig NADPH (Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
E-EL-GP0836	Guinea pig PK (Pyruvate Kinase, Liver And RBC) ELISA Kit
E-EL-GP0837	Guinea pig CHD5 (Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
E-EL-GP0838	Guinea pig DPT (Dermatopontin) ELISA Kit
E-EL-GP0839	Guinea pig LIPH (Lipase H) ELISA Kit
E-EL-GP0840	Guinea pig Anti-CIIAb (Anti-Type II Collagen Antibody) ELISA Kit
E-EL-GP0841	Guinea pig P-CK (Pan-Cytokeratin) ELISA Kit
E-EL-GP0842	Guinea pig PDGF (Platelet-Derived Growth Factor) ELISA Kit
E-EL-GP0843	Guinea pig LAT (Linker For Activation of T-cell) ELISA Kit
E-EL-GP0844	Guinea pig OSM (Oncostatin M) ELISA Kit
E-EL-GP0846	Guinea pig VEGF-C (Vascular Endothelial Growth Factor C) ELISA Kit
E-EL-GP0847	Guinea pig MTX (Methotrexate) ELISA Kit
E-EL-GP0848	Guinea pig MANF (Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit
E-EL-GP0850	Guinea pig CR (Calretinin) ELISA Kit
E-EL-GP0851	Guinea pig FGF21(Fibroblast Growth Factor 21) ELISA Kit
E-EL-GP0852	Guinea pig CTSK (Cathepsin K) ELISA Kit
E-EL-GP0853	Guinea pig RGS4 (Regulator of G Protein Signaling 4) ELISA Kit

Cat.NO	Product Name
E-EL-GP0854	Guinea pig PCDHyA2 (Protocadherin Gamma A2) ELISA Kit
E-EL-GP0856	Guinea pig PTPRO (Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
E-EL-GP0857	Guinea pig VIP (Vasoactive Intestinal Peptide) ELISA Kit
E-EL-GP0858	Guinea pig ACP1 (Acid Phosphatase 1) ELISA Kit
E-EL-GP0859	Guinea pig FAIM3 (Fas Apoptotic Inhibitory Molecule 3) ELISA Kit
E-EL-GP0860	Guinea pig MAP2K6/MKK6 (Mitogen Activated Protein Kinase Kinase 6) ELISA Kit
E-EL-GP0861	Guinea pig PG (Proteoglycan) ELISA Kit
E-EL-GP0863	Guinea pig GP39 (Glycoprotein 39, Cartilage) ELISA Kit
E-EL-GP0864	Guinea pig MYH1 (Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit
E-EL-GP0865	Guinea pig CA5B (Carbonic Anhydrase VB, Mitochondrial) ELISA Kit
E-EL-GP0866	Guinea pig CDC23 (Cell Division Cycle Protein 23) ELISA Kit
E-EL-GP0867	Guinea pig FSHR (Follicle Stimulating Hormone Receptor) ELISA Kit
E-EL-GP0868	Guinea pig GST μ 2/GSTm2 (Glutathione S Transferase Mu 2, Muscle) ELISA Kit
E-EL-GP0869	Guinea pig TSGF (Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
E-EL-GP0870	Guinea pig IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit
E-EL-GP0871	Guinea pig PRDX5 (Peroxisome oxidoreductase 5) ELISA Kit
E-EL-GP0872	Guinea pig MRP (Multidrug Resistance-Associated Protein) ELISA Kit
E-EL-GP0873	Guinea pig Bid (BH3 Interacting Domain Death Agonist) ELISA Kit
E-EL-GP0874	Guinea pig TD-Ag (Thymus Dependent Antigen) ELISA Kit
E-EL-GP0875	Guinea pig AZU1 (Azurocidin 1) ELISA Kit

Cat.NO	Product Name
E-EL-GP0876	Guinea pig GAL2 (Galectin 2) ELISA Kit
E-EL-GP0877	Guinea pig CEP110 (Centrosomal Protein 110kDa) ELISA Kit
E-EL-GP0878	Guinea pig GREM1 (Gremlin 1) ELISA Kit
E-EL-GP0879	Guinea pig IL-6R (Interleukin 6 Receptor) ELISA Kit
E-EL-GP0880	Guinea pig BMP-1 (Bone Morphogenetic Protein 1) ELISA Kit
E-EL-GP0881	Guinea pig MAU (Microalbuminuria) ELISA Kit
E-EL-GP0882	Guinea pig PCDH15 (Protocadherin 15) ELISA Kit
E-EL-GP0883	Guinea pig ESF (Erythropoiesis Stimulating Factor) ELISA Kit
E-EL-GP0885	Guinea pig ERN1 (Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit
E-EL-GP0886	Guinea pig PTPRJ (Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit
E-EL-GP0888	Guinea pig PDGFR β (Platelet Derived Growth Factor Receptor Beta) ELISA Kit
E-EL-GP0889	Guinea pig GHITM (Growth Hormone Inducible Transmembrane Protein) ELISA Kit
E-EL-GP0890	Guinea pig ELA2B (Elastase 2B) ELISA Kit
E-EL-GP0891	Guinea pig α -GST (α -glutathione S-transferases) ELISA Kit
E-EL-GP0892	Guinea pig COL8 α 1 (Collagen Type VIII Alpha 1) ELISA Kit
E-EL-GP0893	Guinea pig CA4 (Carbonic Anhydrase IV) ELISA Kit
E-EL-GP0894	Guinea pig CHD3 (Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
E-EL-GP0895	Guinea pig LAP (Leucine Aminopeptidase) ELISA Kit
E-EL-GP0896	Guinea pig TU3A (TU3A Protein) ELISA Kit
E-EL-GP0897	Guinea pig MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-GP0898	Guinea pig GDNF (Glial Cell Line Derived Neurotrophic Factor) ELISA Kit
E-EL-GP0900	Guinea pig APT (Abnormal Prothrombin) ELISA Kit
E-EL-GP0901	Guinea pig α FP (Alpha-Fetoprotein) ELISA Kit
E-EL-GP0902	Guinea pig RNASET2 (Ribonuclease T2) ELISA Kit
E-EL-GP0903	Guinea pig BACE2 (β -site APP Cleaving Enzyme 2) ELISA Kit
E-EL-GP0905	Guinea pig E-Cad (E-Cadherin) ELISA Kit
E-EL-GP0906	Guinea pig GST ω 1 (Glutathione S Transferase Omega 1) ELISA Kit
E-EL-GP0907	Guinea pig PON3 (Paraoxonase 3) ELISA Kit
E-EL-GP0908	Guinea pig PKD1 (Protein Kinase D1) ELISA Kit
E-EL-GP0909	Guinea pig PGR (Progesterone Receptor) ELISA Kit
E-EL-GP0910	Guinea pig P2 (Myelin Protein 2) ELISA Kit
E-EL-GP0911	Guinea pig NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
E-EL-GP0912	Guinea pig D-LDH (D-Lactate Dehydrogenase) ELISA Kit
E-EL-GP0913	Guinea pig AR (Amphiregulin) ELISA Kit
E-EL-GP0914	Guinea pig HELIX- II (Type II Collagen Helical Peptide) ELISA Kit
E-EL-GP0915	Guinea pig MSM β (Microseminoprotein Beta) ELISA Kit
E-EL-GP0916	Guinea pig PFKP (Phosphofructokinase, Platelet) ELISA Kit
E-EL-GP0917	Guinea pig PAL (L-Phenylalanine Ammonia-Lyase) ELISA Kit
E-EL-GP0918	Guinea pig PCSK9 (Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit
E-EL-GP0919	Guinea pig GDF9 (Growth Differentiation Factor 9) ELISA Kit

Cat.NO	Product Name
E-EL-GP0922	Guinea pig F13 (Coagulation Factor X III) ELISA Kit
E-EL-GP0923	Guinea pig CYSLTR1 (Cysteinyl Leukotriene Receptor 1) ELISA Kit
E-EL-GP0924	Guinea pig MAN1A1 (Mannosidase Alpha Class 1A Member 1) ELISA Kit
E-EL-GP0925	Guinea pig APC (Activated Protein C) ELISA Kit
E-EL-GP0926	Guinea pig β TG/PBP/CXCL7/NAP2 (Thromboglobulin, Beta) ELISA Kit
E-EL-GP0927	Guinea pig STK39 (Serine/Threonine Kinase 39) ELISA Kit
E-EL-GP0928	Guinea pig CK (Creatine Kinase) ELISA Kit
E-EL-GP0929	Guinea pig PGF (Prostaglandin F) ELISA Kit
E-EL-GP0930	Guinea pig DAO (Diamine Oxidase) ELISA Kit
E-EL-GP0931	Guinea pig VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
E-EL-GP0932	Guinea pig F2 (Coagulation Factor II) ELISA Kit
E-EL-GP0934	Guinea pig PRG4 (Proteoglycan 4) ELISA Kit
E-EL-GP0935	Guinea pig APA (Anti-Phospholipid Antibody) ELISA Kit
E-EL-GP0937	Guinea pig RPL10 (Ribosomal Protein L10) ELISA Kit
E-EL-GP0938	Guinea pig MBLR (Mannose Binding Lectin Receptor) ELISA Kit
E-EL-GP0939	Guinea pig SNCy (Synuclein Gamma) ELISA Kit
E-EL-GP0940	Guinea pig TR (Thrombin Receptor) ELISA Kit
E-EL-GP0941	Guinea pig DPD (Deoxy pyridinoline) ELISA Kit
E-EL-GP0942	Guinea pig PKI γ (Protein Kinase Inhibitor Gamma) ELISA Kit
E-EL-GP0943	Guinea pig AP (Aprotinin) ELISA Kit

Cat.NO	Product Name
E-EL-GP0944	Guinea pig FOXO1 (Forkhead Box Protein O1) ELISA Kit
E-EL-GP0946	Guinea pig CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit
E-EL-GP0947	Guinea pig GLRX3 (Glutaredoxin 3) ELISA Kit
E-EL-GP0948	Guinea pig OPGL (Osteoprotegerin Ligand) ELISA Kit
E-EL-GP0949	Guinea pig IL-15 (Interleukin 15) ELISA Kit
E-EL-GP0950	Guinea pig TPO (Thrombopoietin) ELISA Kit
E-EL-GP0951	Guinea pig DHEA-S7 (Dehydroepiandrosterone S7) ELISA Kit
E-EL-GP0952	Guinea pig PKBG (Protein Kinase B Gamma) ELISA Kit
E-EL-GP0953	Guinea pig C4a (Complement Component 4a) ELISA Kit
E-EL-GP0954	Guinea pig PRF1 (Perforin 1/Pore-Forming Protein) ELISA Kit
E-EL-GP0955	Guinea pig Gas6 (Growth Arrest Specific Protein 6) ELISA Kit
E-EL-GP0956	Guinea pig IDE (Insulin Degrading Enzyme) ELISA Kit
E-EL-GP0957	Guinea pig AKAP1 (A Kinase Anchor Protein 1) ELISA Kit
E-EL-GP0958	Guinea pig SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit
E-EL-GP0959	Guinea pig PFKM (Phosphofructokinase, Muscle) ELISA Kit
E-EL-GP0960	Guinea pig GIP (Gastric Inhibitory Polypeptide) ELISA Kit
E-EL-GP0961	Guinea pig VIM (Vimentin) ELISA Kit
E-EL-GP0962	Guinea pig FGF4 (Fibroblast Growth Factor 4) ELISA Kit
E-EL-GP0963	Guinea pig TGF- β 2 (Transforming Growth Factor β 2) ELISA Kit
E-EL-GP0964	Guinea pig CD30 (Cluster of Differentiation 30) ELISA Kit

Cat.NO	Product Name
E-EL-GP0965	Guinea pig D5R/DRD5 (Dopamine Receptor D5) ELISA Kit
E-EL-GP0966	Guinea pig Tgtp (T-cell Specific Gtpase) ELISA Kit
E-EL-GP0967	Guinea pig GAL (Galanin) ELISA Kit
E-EL-GP0968	Guinea pig D3R/DRD3 (Dopamine Receptor D3) ELISA Kit
E-EL-GP0969	Guinea pig HbA1c (Glycosylated Hemoglobin/Hemoglobin A1c) ELISA Kit
E-EL-GP0970	Guinea pig TR (Testosterone Receptor) ELISA Kit
E-EL-GP0971	Guinea pig SPD (Pulmonary Surfactant Associated Protein D) ELISA Kit
E-EL-GP0972	Guinea pig PLGA (Plasminogen Activator) ELISA Kit
E-EL-GP0973	Guinea pig TM (Thrombomodulin) ELISA Kit
E-EL-GP0974	Guinea pig IAA (Insulin Autoantibody) ELISA Kit
E-EL-GP0975	Guinea pig ST2 (Syntenin 2) ELISA Kit
E-EL-GP0976	Guinea pig TPO-Ab (Anti-Thyroid-Peroxidase Antibody) ELISA Kit
E-EL-GP0977	Guinea pig MSLN (Mesothelin) ELISA Kit
E-EL-GP0978	Guinea pig MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
E-EL-GP0979	Guinea pig rT3 (Reverse Triiodothyronine) ELISA Kit
E-EL-GP0980	Guinea pig VEGF (Vascular Endothelial Cell Growth Factor) ELISA Kit
E-EL-GP0981	Guinea pig TARC (Thymus Activation Regulated Chemokine) ELISA Kit
E-EL-GP0984	Guinea pig PDE4D (Phosphodiesterase 4D, cAMP Specific) ELISA Kit
E-EL-GP0985	Guinea pig ALB (Albumin) ELISA Kit
E-EL-GP0986	Guinea pig MAG (Myelin Associated Glycoprotein) ELISA Kit

Cat.NO	Product Name
E-EL-GP0987	Guinea pig TRY (Trypsin) ELISA Kit
E-EL-GP0988	Guinea pig Ly6a (Lymphocyte Antigen 6 Complex Locus A) ELISA Kit
E-EL-GP0989	Guinea pig CPR (Cytochrome P450 Reductase) ELISA Kit
E-EL-GP0991	Guinea pig TAF2 (TATA Box Binding Protein Associated Factor 2) ELISA Kit
E-EL-GP0992	Guinea pig TPO (Thyroid Peroxidase) ELISA Kit
E-EL-GP0994	Guinea pig CT (Calcitonin) ELISA Kit
E-EL-GP0995	Guinea pig TECK/CCL25 (Thymus Expressed Chemokine) ELISA Kit
E-EL-GP0996	Guinea pig PLB (Phospholipase B) ELISA Kit
E-EL-GP0997	Guinea pig MTHFR (Methylenetetrahydrofolate Reductase) ELISA Kit
E-EL-GP0998	Guinea pig TTN (Titin) ELISA Kit
E-EL-GP0999	Guinea pig MCF (Macrophage Chemotatic Factor) ELISA Kit
E-EL-GP1000	Guinea pig FCN1 (Ficolin 1) ELISA Kit
E-EL-GP1001	Guinea pig GKN1 (Gastrokine 1) ELISA Kit
E-EL-GP1003	Guinea pig C/EBP δ (CCAAT/Enhancer Binding Protein, Delta) ELISA Kit
E-EL-GP1004	Guinea pig Gal (α -Galactoyl) ELISA Kit
E-EL-GP1005	Guinea pig MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit
E-EL-GP1006	Guinea pig SA (Sialic Acid) ELISA Kit
E-EL-GP1007	Guinea pig PTHR2 (Parathyroid Hormone Receptor 2) ELISA Kit
E-EL-GP1008	Guinea pig FOXP1 (Forkhead Box Protein P1) ELISA Kit
E-EL-GP1010	Guinea pig TM (Thrombin) ELISA Kit

Cat.NO	Product Name
E-EL-GP1011	Guinea pig ANEA (Anti-Neutrophil Elastase Antibody) ELISA Kit
E-EL-GP1012	Guinea pig TF (Transfer Factor) ELISA Kit
E-EL-GP1013	Guinea pig dDAVP (1-Desamino-8D-Arginine Vasopressin) ELISA Kit
E-EL-GP1014	Guinea pig EFNA3 (Ephrin A3) ELISA Kit
E-EL-GP1015	Guinea pig SDC4 (Syndecan 4) ELISA Kit
E-EL-GP1017	Guinea pig AFA (Anti-Filaggrin Antibody) ELISA Kit
E-EL-GP1018	Guinea pig IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit
E-EL-GP1019	Guinea pig PAC3 (Platelet Associated Complement 3) ELISA Kit
E-EL-GP1020	Guinea pig ODC (Ornithine Decarboxylase) ELISA Kit
E-EL-GP1021	Guinea pig CKLF (Chemokine Like Factor) ELISA Kit
E-EL-GP1022	Guinea pig GDF5 (Growth Differentiation Factor 5) ELISA Kit
E-EL-GP1023	Guinea pig ATF6 (Activating Transcription Factor 6) ELISA Kit
E-EL-GP1024	Guinea pig PICP (Procollagen I C-terminal Propeptide) ELISA Kit
E-EL-GP1025	Guinea pig GAL1 (Galectin 1) ELISA Kit
E-EL-GP1026	Guinea pig MMP-12 (Matrix Metalloproteinase 12) ELISA Kit
E-EL-GP1027	Guinea pig NTRK3 (Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
E-EL-GP1029	Guinea pig CCND3 (Cyclin-D3) ELISA Kit
E-EL-GP1030	Guinea pig NTXI (Cross Linked N-telopeptide of Type I Collagen) ELISA Kit
E-EL-GP1031	Guinea pig ATM (Ataxia Telangiectasia Mutated) ELISA Kit
E-EL-GP1032	Guinea pig PDGFRL (Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit

Cat.NO	Product Name
E-EL-GP1034	Guinea pig GC (Glycocalicin) ELISA Kit
E-EL-GP1036	Guinea pig GATA5 (GATA Binding Protein 5) ELISA Kit
E-EL-GP1037	Guinea pig FDP (Fibrinogen Degradation Product) ELISA Kit
E-EL-GP1039	Guinea pig TAFI (Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit
E-EL-GP1041	Guinea pig PIICP (Procollagen II C-Terminal ProPeptide) ELISA Kit
E-EL-GP1042	Guinea pig ACHAP (Acetylcholinesterase Associated Protein) ELISA Kit
E-EL-GP1043	Guinea pig ATXN1 (Ataxin 1) ELISA Kit
E-EL-GP1044	Guinea pig C3c (Complement C3 Convertase) ELISA Kit
E-EL-GP1045	Guinea pig LUM (Lumican) ELISA Kit
E-EL-GP1046	Guinea pig CaN (Calcineurin) ELISA Kit
E-EL-GP1047	Guinea pig TrxR1 (Thioredoxin Reductase 1) ELISA Kit
E-EL-GP1048	Guinea pig TSP-4 (Thrombospondin-4) ELISA Kit
E-EL-GP1049	Guinea pig ATF7 (Activating Transcription Factor 7) ELISA Kit
E-EL-GP1050	Guinea pig SPB (Pulmonary Surfactant Associated Protein B) ELISA Kit
E-EL-GP1051	Guinea pig PTF1 α (Pancreas Specific Transcription Factor 1 α) ELISA Kit
E-EL-GP1052	Guinea pig TNK2 (Tyrosine Kinase, Non Receptor 2) ELISA Kit
E-EL-GP1054	Guinea pig sCKR (Soluble Cytokine Receptor) ELISA Kit
E-EL-GP1055	Guinea pig IL-19 (Interleukin 19) ELISA Kit
E-EL-GP1056	Guinea pig MYH9 (Myosin Heavy Chain 9, Non Muscle) ELISA Kit
E-EL-GP1057	Guinea pig STAT1 (Signal Transducer and Activator of Transcription 1) ELISA Kit

Cat.NO	Product Name
E-EL-GP1059	Guinea pig PRDX2 (Peroxiredoxin 2) ELISA Kit
E-EL-GP1060	Guinea pig ACP3/PAP (Prostatic Acid Phosphatase 3) ELISA Kit
E-EL-GP1061	Guinea pig OPG (Osteoprotegerin) ELISA Kit
E-EL-GP1062	Guinea pig Cyt-C (Cytochrome C) ELISA Kit
E-EL-GP1063	Guinea pig TI-Ag (Thymus Independent Antigen) ELISA Kit
E-EL-GP1064	Guinea pig Dyn (Big Dynorphin) ELISA Kit
E-EL-GP1065	Guinea pig CAM-Ab (Anti-Calmodulin Specific Antibody) ELISA Kit
E-EL-GP1066	Guinea pig PCT (Procalcitonin) ELISA Kit
E-EL-GP1067	Guinea pig SOX3 (Sex Determining Region Y Box Protein 3) ELISA Kit
E-EL-GP1068	Guinea pig PKI β (Protein Kinase Inhibitor Beta) ELISA Kit
E-EL-GP1069	Guinea pig FGF19 (Fibroblast Growth Factor 19) ELISA Kit
E-EL-GP1070	Guinea pig CNTN3 (Contactin 3) ELISA Kit
E-EL-GP1071	Guinea pig GAD (Glutamic Acid Decarboxylase) ELISA Kit
E-EL-GP1072	Guinea pig DNM2 (Dynamin 2) ELISA Kit
E-EL-GP1073	Guinea pig α -SMA (α -Smooth muscle actin) ELISA Kit
E-EL-GP1074	Guinea pig TFPI2 (Tissue Factor Pathway Inhibitor 2) ELISA Kit
E-EL-GP1075	Guinea pig ERCC1 (Excision Repair Cross Complementing Rodent Repair Deficiency Complementation 1) ELISA Kit
E-EL-GP1076	Guinea pig MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit
E-EL-GP1077	Guinea pig β 2-GP (Beta2 Glycoprotein) ELISA Kit
E-EL-GP1078	Guinea pig IL-11 (Interleukin 11) ELISA Kit

Cat.NO	Product Name
E-EL-GP1080	Guinea pig BLNK (B-Cell Linker Protein) ELISA Kit
E-EL-GP1081	Guinea pig MASP2 (Mannan Associated Serine Protease 2) ELISA Kit
E-EL-GP1082	Guinea pig UreC (Urease Related Protein C) ELISA Kit
E-EL-GP1084	Guinea pig HB β (Hemoglobin Beta) ELISA Kit
E-EL-GP1085	Guinea pig SELP (P-Selectin) ELISA Kit
E-EL-GP1086	Guinea pig MMP-9 (Matrix Metalloproteinase 9) ELISA Kit
E-EL-GP1087	Guinea pig S100A6 (S100 Calcium Binding Protein A6) ELISA Kit
E-EL-GP1088	Guinea pig VWF (Von Willebrand Factor) ELISA Kit
E-EL-GP1089	Guinea pig PCPE1 (Procollagen C Proteinase Enhancer 1) ELISA Kit
E-EL-GP1090	Guinea pig Salusin α (Salusin Alpha) ELISA Kit
E-EL-GP1091	Guinea pig ApoE (Apolipoprotein E) ELISA Kit
E-EL-GP1092	Guinea pig STAT2 (Signal Transducer and Activator of Transcription 2) ELISA Kit
E-EL-GP1094	Guinea pig Methylase ELISA Kit
E-EL-GP1095	Guinea pig CD3d (Cluster of Differentiation 3d) ELISA Kit
E-EL-GP1096	Guinea pig SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit
E-EL-GP1097	Guinea pig α 1-AGP (Alpha-1-Acid Glycoprotein) ELISA Kit
E-EL-GP1098	Guinea pig CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
E-EL-GP1100	Guinea pig C/EBP β (CCAAT/Enhancer Binding Protein Beta) ELISA Kit
E-EL-GP1101	Guinea pig S100A11 (S100 Calcium Binding Protein A11) ELISA Kit
E-EL-GP1102	Guinea pig FG (Fibrinogen) ELISA Kit

Cat.NO	Product Name
E-EL-GP1104	Guinea pig MUC7 (Mucin 7) ELISA Kit
E-EL-GP1105	Guinea pig CAMK2 α (Calcium/Calmodulin Dependent Protein Kinase II Alpha) ELISA Kit
E-EL-GP1106	Guinea pig EL (EuglobulinLysis) ELISA Kit
E-EL-GP1107	Guinea pig TERT (Telomerase Reverse Transcriptase) ELISA Kit
E-EL-GP1108	Guinea pig SIGLEC1 (Sialic Acid Binding Ig Like Lectin 1) ELISA Kit
E-EL-GP1109	Guinea pig GFAP (Glial Fibrillary Acidic Protein) ELISA Kit
E-EL-GP1110	Guinea pig Fc γ R II/CD32 (Receptor II for the Fc Fragment of IgG) ELISA Kit
E-EL-GP1111	Guinea pig NRIP1 (Nuclear Receptor Interacting Protein 1) ELISA Kit
E-EL-GP1112	Guinea pig OTF2A (Octamer Transcription Factor 2A) ELISA Kit
E-EL-GP1113	Guinea pig NE/ELA2 (Elastase 2, Neutrophil) ELISA Kit
E-EL-GP1114	Guinea pig ALDOC (Aldolase C, Fructose Bisphosphate) ELISA Kit
E-EL-GP1115	Guinea pig FPR2 (Formyl Peptide Receptor 2) ELISA Kit
E-EL-GP1116	Guinea pig HLTF (Helicase Like Transcription Factor) ELISA Kit
E-EL-GP1117	Guinea pig OSF (Osteoclast Stimulating Factor 1) ELISA Kit
E-EL-GP1118	Guinea pig FAS/CD95/TNFRSF6 (Factor Related Apoptosis) ELISA Kit
E-EL-GP1119	Guinea pig pERK (Phospho Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-GP1120	Guinea pig MBP/MBL (Mannose Binding Protein/Mannose Binding Lectin) ELISA Kit
E-EL-GP1121	Guinea pig CDSN (Corneodesmosin) ELISA Kit
E-EL-GP1122	Guinea pig CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit
E-EL-GP1123	Guinea pig PKB (Protein Kinase B) ELISA Kit

Cat.NO	Product Name
E-EL-GP1124	Guinea pig SIGLEC10 (Sialic Acid Binding Ig Like Lectin 10) ELISA Kit
E-EL-GP1125	Guinea pig VEGF-B (Vascular Endothelial Cell Growth Factor B) ELISA Kit
E-EL-GP1126	Guinea pig PIAS1 (Protein Inhibitor of Activated STAT 1) ELISA Kit
E-EL-GP1129	Guinea pig SAA (Serum Amyloid A) ELISA Kit
E-EL-GP1130	Guinea pig CKLFSF4 (Chemokine Like Factor Superfamily 4) ELISA Kit
E-EL-GP1131	Guinea pig XCR1 (XC-Chemokine Receptor 1) ELISA Kit
E-EL-GP1132	Guinea pig NUP37 (Nucleoporin 37kDa) ELISA Kit
E-EL-GP1133	Guinea pig PGA (Pepsinogen A) ELISA Kit
E-EL-GP1135	Guinea pig GIP (Glucose Dependent Insulin Releasing Polypeptide) ELISA Kit
E-EL-GP1137	Guinea pig CRYαB (Crystallin Alpha B) ELISA Kit
E-EL-GP1138	Guinea pig TMPRSS2 (Transmembrane Protease, Serine 2) ELISA Kit
E-EL-GP1139	Guinea pig CALU (Calumenin) ELISA Kit
E-EL-GP1140	Guinea pig IFN-β (Interferon Beta) ELISA Kit
E-EL-GP1141	Guinea pig MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit
E-EL-GP1142	Guinea pig IgM (Immunoglobulin M) ELISA Kit
E-EL-GP1143	Guinea pig GATA4 (GATA Binding Protein 4) ELISA Kit
E-EL-GP1144	Guinea pig NFKBIE (Nuclear Factor Of Kappa Light Polypeptide Gene Enhancer In B-cells Inhibitor Epsilon) ELISA Kit
E-EL-GP1145	Guinea pig PRSS23 (Protease, Serine 23) ELISA Kit
E-EL-GP1147	Guinea pig IGFBP-1 (Insulin Like Growth Factor Binding Protein 1) ELISA Kit
E-EL-GP1149	Guinea pig MTF (Microtransferrin) ELISA Kit

Cat.NO	Product Name
E-EL-GP1150	Guinea pig CLDN4 (Claudin 4) ELISA Kit
E-EL-GP1152	Guinea pig NXF1 (Nuclear RNA Export Factor 1) ELISA Kit
E-EL-GP1153	Guinea pig SF1 (Steroidogenic Factor 1) ELISA Kit
E-EL-GP1154	Guinea pig FRS3 (Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
E-EL-GP1155	Guinea pig TH (Tyrosine Hydroxylase) ELISA Kit
E-EL-GP1156	Guinea pig ADCY4 (Adenylate Cyclase 4) ELISA Kit
E-EL-GP1157	Guinea pig PRI (Placental Ribonuclease Inhibitor) ELISA Kit
E-EL-GP1158	Guinea pig PT (Prothrombin) ELISA Kit
E-EL-GP1159	Guinea pig RAG-1 (Recombination Activating Gene 1) ELISA Kit
E-EL-GP1160	Guinea pig α Manase (α Mannosidase) ELISA Kit
E-EL-GP1161	Guinea pig GCSFR (Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
E-EL-GP1164	Guinea pig ApoA1 (Apolipoprotein A1) ELISA Kit
E-EL-GP1165	Guinea pig PARC/CCL18 (Pulmonary Activation Regulated Chemokine) ELISA Kit
E-EL-GP1166	Guinea pig HSP-90 (Heat Shock Protein 90) ELISA Kit
E-EL-GP1167	Guinea pig PDI (Protein Disulfide Isomerase) ELISA Kit
E-EL-GP1168	Guinea pig S100A8 (S100 Calcium Binding Protein A8/ Calgranulin A) ELISA Kit
E-EL-GP1169	Guinea pig Ext-N (N-Terminal Extein) ELISA Kit
E-EL-GP1170	Guinea pig BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit
E-EL-GP1171	Guinea pig CLEC2C (C-Type Lectin Domain Family 2 Member C) ELISA Kit
E-EL-GP1172	Guinea pig AGR2 (Anterior Gradient Protein 2) ELISA Kit

Cat.NO	Product Name
E-EL-GP1173	Guinea pig SEPP1 (Selenoprotein P) ELISA Kit
E-EL-GP1174	Guinea pig EREG (Epiregulin) ELISA Kit
E-EL-GP1175	Guinea pig PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit
E-EL-GP1176	Guinea pig ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit
E-EL-GP1177	Guinea pig MUC2 (Mucin 2) ELISA Kit
E-EL-GP1178	Guinea pig PRSS12 (Protease, Serine 12) ELISA Kit
E-EL-GP1179	Guinea pig PPAR- δ (Peroxisome Proliferator Activated Receptor Delta) ELISA Kit
E-EL-GP1180	Guinea pig NPY (Neuropeptide Y) ELISA Kit
E-EL-GP1181	Guinea pig INSL3 (Insulin Like Protein 3) ELISA Kit
E-EL-GP1182	Guinea pig CALD (Caldesmon) ELISA Kit
E-EL-GP1183	Guinea pig SYNPO (Synaptopodin) ELISA Kit
E-EL-GP1184	Guinea pig C/EBP ϵ (CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit
E-EL-GP1185	Guinea pig PIR (Pirin) ELISA Kit
E-EL-GP1186	Guinea pig CCKBR (Cholecystokinin B Receptor) ELISA Kit
E-EL-GP1187	Guinea pig PF3 (Platelet Factor 3) ELISA Kit
E-EL-GP1188	Guinea pig SDF-2 (Stromal Cell Derived Factor 2) ELISA Kit
E-EL-GP1189	Guinea pig M-AChR M2 (Muscarinic Acetylcholine Receptor M2) ELISA Kit
E-EL-GP1190	Guinea pig VDR (Vitamin D Receptor) ELISA Kit
E-EL-GP1191	Guinea pig IL-2R α (Interleukin 2 Receptor Alpha) ELISA Kit
E-EL-GP1192	Guinea pig DHEA-S (Dehydroepiandrosterone Sulfate) ELISA Kit

Cat.NO	Product Name
E-EL-GP1193	Guinea pig RIPK2 (Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit
E-EL-GP1194	Guinea pig PEX2 (Peroxisomal Biogenesis Factor 2) ELISA Kit
E-EL-GP1195	Guinea pig FOXC1 (Forkhead Box Protein C1) ELISA Kit
E-EL-GP1196	Guinea pig TAGLN (Transgelin) ELISA Kit
E-EL-GP1197	Guinea pig Hpc (Hepcidin) ELISA Kit
E-EL-GP1198	Guinea pig ApoF (Apolipoprotein F) ELISA Kit
E-EL-GP1199	Guinea pig PHLDA2 (Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit
E-EL-GP1200	Guinea pig TEF (Thyrotroph Embryonic Factor) ELISA Kit
E-EL-GP1201	Guinea pig GPR120 (G Protein Coupled Receptor 120) ELISA Kit
E-EL-GP1202	Guinea pig Grb10 (Growth Factor Receptor Bound Protein 10) ELISA Kit
E-EL-GP1203	Guinea pig OFQ/N (Orphanin FQ/Nociceptin) ELISA Kit
E-EL-GP1204	Guinea pig 11-DH-TXB2 (11-Dehydro-thromboxane B2) ELISA Kit
E-EL-GP1205	Guinea pig MUC5B (Mucin 5 Subtype B) ELISA Kit
E-EL-GP1206	Guinea pig TXLNa (Taxilin Alpha) ELISA Kit
E-EL-GP1207	Guinea pig MR (Mannose Receptor) ELISA Kit
E-EL-GP1208	Guinea pig CLDN3 (Claudin 3) ELISA Kit
E-EL-GP1209	Guinea pig PDGF-AA (Platelet Derived Growth Factor-AA) ELISA Kit
E-EL-GP1210	Guinea pig 5-NT (5-Nucleotidase) ELISA Kit
E-EL-GP1211	Guinea pig LR (Leukoregulin) ELISA Kit
E-EL-GP1212	Guinea pig S100A9 (S100 Calcium Binding Protein A9) ELISA Kit

Cat.NO	Product Name
E-EL-GP1213	Guinea pig α HBDH (α -Hydroxybutyrate Dehydrogenase) ELISA Kit
E-EL-GP1214	Guinea pig TREM-1 (Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
E-EL-GP1215	Guinea pig TSP-2 (Thrombospondin-2) ELISA Kit
E-EL-GP1216	Guinea pig EMA (Epithelial Membrane Antigen) ELISA Kit
E-EL-GP1217	Guinea pig EPCR (Endothelial Protein C Receptor) ELISA Kit
E-EL-GP1218	Guinea pig LATS (Long Acting Thyroid Stimulator) ELISA Kit
E-EL-GP1219	Guinea pig SEPCR/PROCR (Soluble Endothelial Protein C Receptor) ELISA Kit
E-EL-GP1220	Guinea pig EK (Eepidermal Keratin) ELISA Kit
E-EL-GP1221	Guinea pig CK-HMW (Cytokeratin High Molecular Weight) ELISA Kit
E-EL-GP1223	Guinea pig NRG-3 (Neuregulin 3) ELISA Kit
E-EL-GP1224	Guinea pig CASP3 (Caspase 3) ELISA Kit
E-EL-GP1225	Guinea pig GRIK2 (Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit
E-EL-GP1226	Guinea pig SNRPD1 (Small Nuclear Ribonucleoprotein Polypeptide D1) ELISA Kit
E-EL-GP1227	Guinea pig SDF-4 (Stromal Cell Derived Factor 4) ELISA Kit
E-EL-GP1228	Guinea pig AGPA (Anti-Glucoprotein 210 Antibody) ELISA Kit
E-EL-GP1229	Guinea pig ERK (Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-GP1231	Guinea pig CYB561 (Cytochrome b-561) ELISA Kit
E-EL-GP1233	Guinea pig PP (Protein Phosphatase) ELISA Kit
E-EL-GP1234	Guinea pig I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit
E-EL-GP1235	Guinea pig IgA1 (Immunoglobulin A1) ELISA Kit

Cat.NO	Product Name
E-EL-GP1237	Guinea pig TL1A (Tumor Necrosis Factor Related Ligand 1A) ELISA Kit
E-EL-GP1238	Guinea pig BSP (Bone Sialoprotein) ELISA Kit
E-EL-GP1239	Guinea pig PDAP1 (Pdgfa Associated Protein 1) ELISA Kit
E-EL-GP1240	Guinea pig CX43 (Connexin 43) ELISA Kit
E-EL-GP1241	Guinea pig VPI (Vasoactive Peptidase Inhibitors) ELISA Kit
E-EL-GP1242	Guinea pig GMFβ (Glia Maturation Factor, Beta) ELISA Kit
E-EL-GP1243	Guinea pig GLUT3 (Glucose Transporter 3) ELISA Kit
E-EL-GP1244	Guinea pig NUP85 (Nucleoporin 85kDa) ELISA Kit
E-EL-GP1245	Guinea pig NNE (Non-Neuronal Enolase) ELISA Kit
E-EL-GP1246	Guinea pig PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit
E-EL-GP1247	Guinea pig SELL (L-Selectin) ELISA Kit
E-EL-GP1248	Guinea pig CREB (Cyclic AMP Response Element Binding Protein) ELISA Kit
E-EL-GP1249	Guinea pig HYOU1 (Hypoxia Up Regulated 1) ELISA Kit
E-EL-GP1250	Guinea pig PDK4 (Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit
E-EL-GP1251	Guinea pig Bcl-10 (B-cell Leukemia/Lymphoma 10) ELISA Kit
E-EL-GP1252	Guinea pig AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit
E-EL-GP1253	Guinea pig PKP1 (Plakophilin 1) ELISA Kit
E-EL-GP1254	Guinea pig GP-73 (Golgi Protein 73) ELISA Kit
E-EL-GP1255	Guinea pig FcεR I (Receptor I for the Fc Fragment of IgE) ELISA Kit
E-EL-GP1256	Guinea pig CYTIP (Cytohesin Interacting Protein) ELISA Kit

Cat.NO	Product Name
E-EL-GP1257	Guinea pig IL-3 (Interleukin 3) ELISA Kit
E-EL-GP1258	Guinea pig PP (Placental Protein) ELISA Kit
E-EL-GP1259	Guinea pig UPP1 (Uridine Phosphorylase 1) ELISA Kit
E-EL-GP1260	Guinea pig Gal-6S (Galactose-6-Sulphatase) ELISA Kit
E-EL-GP1261	Guinea pig CG (Cryoglobulin) ELISA Kit
E-EL-GP1262	Guinea pig NPYR1 (Neuropeptide Y Receptor Y1) ELISA Kit
E-EL-GP1263	Guinea pig ELN (Elastin) ELISA Kit
E-EL-GP1264	Guinea pig GPC1 (Glypican 1) ELISA Kit
E-EL-GP1265	Guinea pig TNNI2 (Troponin I Type 2, Fast Skeletal) ELISA Kit
E-EL-GP1266	Guinea pig IL-23 p40(Interleukin 23) ELISA Kit
E-EL-GP1267	Guinea pig CSF (Colony Stimulating Factor) ELISA Kit
E-EL-GP1268	Guinea pig DAF (Decay Accelerating Factor) ELISA Kit
E-EL-GP1269	Guinea pig CENPB (Centromere Protein B) ELISA Kit
E-EL-GP1270	Guinea pig CSF2R α (Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit
E-EL-GP1271	Guinea pig GAS (N-Acetylgalactosamine 6-Sulfatase) ELISA Kit
E-EL-GP1272	Guinea pig GS (Gelsolin) ELISA Kit
E-EL-GP1273	Guinea pig REN (Renin) ELISA Kit
E-EL-GP1274	Guinea pig BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit
E-EL-GP1277	Guinea pig α MSH (Alpha-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-GP1278	Guinea pig DMD (Dystrophin) ELISA Kit

Cat.NO	Product Name
E-EL-GP1279	Guinea pig C4BP α (C4 Binding Protein Alpha) ELISA Kit
E-EL-GP1280	Guinea pig NUP205 (Nucleoporin 205kDa) ELISA Kit
E-EL-GP1281	Guinea pig BPI (Bactericidal/Permeability Increasing Protein) ELISA Kit
E-EL-GP1282	Guinea pig GM (Gastric Mucin) ELISA Kit
E-EL-GP1283	Guinea pig RNF4 (Ring Finger Protein 4) ELISA Kit
E-EL-GP1285	Guinea pig L-LDH (L-Lactate Dehydrogenase) ELISA Kit
E-EL-GP1286	Guinea pig CYTH1 (Cytohesin 1) ELISA Kit
E-EL-GP1287	Guinea pig CDK7 (Cyclin Dependent Kinase 7) ELISA Kit
E-EL-GP1288	Guinea pig BAFF/CD257/TNFSF13B (B-Cell Activating Factor) ELISA Kit
E-EL-GP1289	Guinea pig PC4 (Positive Cofactor 4) ELISA Kit
E-EL-GP1290	Guinea pig COMP (Cartilage Oligomeric Matrix Protein) ELISA Kit
E-EL-GP1291	Guinea pig IL-9 (Interleukin 9) ELISA Kit
E-EL-GP1292	Guinea pig COL6 α 3 (Collagen Type VI Alpha 3) ELISA Kit
E-EL-GP1294	Guinea pig AMY1 (Amylase Alpha 1, Salivary) ELISA Kit
E-EL-GP1295	Guinea pig NUP107 (Nucleoporin 107kDa) ELISA Kit
E-EL-GP1297	Guinea pig ATA (Angiotensinase A) ELISA Kit
E-EL-GP1298	Guinea pig EGR3 (Early Growth Response Protein 3) ELISA Kit
E-EL-GP1299	Guinea pig OVA sIgG (Ovalbumin Specific IgG) ELISA Kit
E-EL-GP1301	Guinea pig bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit
E-EL-GP1302	Guinea pig CD59 (Protectin) ELISA Kit

Cat.NO	Product Name
E-EL-GP1303	Guinea pig TRF (Transferrin) ELISA Kit
E-EL-GP1304	Guinea pig KLK6 (Kallikrein 6) ELISA Kit
E-EL-GP1306	Guinea pig CYFRA21-1 (Cytokeratin Fragment Antigen 21-1) ELISA Kit
E-EL-GP1307	Guinea pig CYPB (Cyclophilin B) ELISA Kit
E-EL-GP1308	Guinea pig PI (Proinsulin) ELISA Kit
E-EL-GP1309	Guinea pig SYP (Synaptophysin) ELISA Kit
E-EL-GP1310	Guinea pig ACVA (Activin A) ELISA Kit
E-EL-GP1311	Guinea pig TR (Targeted Ribonuclease) ELISA Kit
E-EL-GP1312	Guinea pig NT (Neurotensin) ELISA Kit
E-EL-GP1313	Guinea pig BDNF(Brain Derived Neurotrophic Factor) ELISA Kit
E-EL-GP1314	Guinea pig F12 (Coagulation Factor XII) ELISA Kit
E-EL-GP1315	Guinea pig PDIA3 (Protein Disulfide Isomerase A3) ELISA Kit
E-EL-GP1317	Guinea pig ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-GP1318	Guinea pig β MSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-GP1319	Guinea pig Ta1 (Thymosin- α 1) ELISA Kit
E-EL-GP1320	Guinea pig PS1 (Presenilin 1) ELISA Kit
E-EL-GP1321	Guinea pig PACAP-38 (Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
E-EL-GP1322	Guinea pig PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-GP1323	Guinea pig NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) ELISA Kit
E-EL-GP1324	Guinea pig IgA (Immunoglobulin A) ELISA Kit

Cat.NO	Product Name
E-EL-GP1325	Guinea pig OC/BGP (Osteocalcin) ELISA Kit
E-EL-GP1326	Guinea pig CTXI (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
E-EL-GP1351	Guinea pig CTX-II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit
E-EL-GP1355	Guinea pig CTX-III (Cross Linked C-telopeptide of Type III Collagen) ELISA Kit
E-EL-GP1356	Guinea pig INS (Insulin) ELISA Kit
E-EL-GP1357	Guinea pig IgG3(immunoglobulin G3)ELISA Kit
E-EL-GP1358	Guinea pig IgG2b(immunoglobulin G2b)ELISA Kit
E-EL-Ch0001	Chicken IgY (Immunoglobulin Y) ELISA Kit
E-EL-Ch0003	Chicken ORM1 (Orosomucoid 1) ELISA Kit
E-EL-Ch0005	Chicken α -Glu (Alphah Glucosidase) ELISA Kit
E-EL-Ch0006	Chicken AmAC-1 (Alternative Macrophage Activation-associated CC Chemokine 1) ELISA Kit
E-EL-Ch0007	Chicken AT- III (Anti-Thrombin III) ELISA Kit
E-EL-Ch0008	Chicken bFGF/FGF2 (Basic Fibroblast Growth Factor)
E-EL-Ch0010	Chicken CRP (C-reactive protein) ELISA Kit
E-EL-Ch0012	Chicken EGF (Epidermal growth factor) ELISA Kit
E-EL-Ch0013	Chicken CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit
E-EL-Ch0014	Chicken MCSP (Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit
E-EL-Ch0015	Chicken PLSCR2 (Phospholipid Scramblase 2) ELISA Kit
E-EL-Ch0016	Chicken MT (Melatonin) ELISA Kit
E-EL-Ch0017	Chicken FLN β (Filamin B, Beta) ELISA Kit

Cat.NO	Product Name
E-EL-Ch0018	Chicken IL-4 (Interleukin 4) ELISA Kit
E-EL-Ch0019	Chicken MYL9 (Myosin Light Chain 9, Regulatory) ELISA Kit
E-EL-Ch0020	Chicken C/EBP ϵ (CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit
E-EL-Ch0021	Chicken GUS β (β -glucuronidase) ELISA Kit
E-EL-Ch0022	Chicken HRG (Histidine Rich Glycoprotein) ELISA Kit
E-EL-Ch0023	Chicken GPC2 (Glypican 2) ELISA Kit
E-EL-Ch0024	Chicken GzmA (Granzyme A) ELISA Kit
E-EL-Ch0025	Chicken C1q-Ab (Anti-Complement 1q Antibody) ELISA Kit
E-EL-Ch0026	Chicken IFN- γ (Interferon Gamma) ELISA Kit
E-EL-Ch0027	Chicken ADAM10 (A Disintegrin And Metalloprotease 10) ELISA Kit
E-EL-Ch0028	Chicken TMPRSS2 (Transmembrane Protease, Serine 2) ELISA Kit
E-EL-Ch0029	Chicken EK (Eepidermal Keratin) ELISA Kit
E-EL-Ch0030	Chicken CLCF1 (Cardiotrophin Like Cytokine Factor 1) ELISA Kit
E-EL-Ch0031	Chicken PP (Pepsin) ELISA Kit
E-EL-Ch0032	Chicken DDT (D-Dopachrome Tautomerase) ELISA Kit
E-EL-Ch0033	Chicken C9 (Complement Component 9) ELISA Kit
E-EL-Ch0034	Chicken BAG3 (Bcl2 Associated Athanogene 3) ELISA Kit
E-EL-Ch0035	Chicken BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit
E-EL-Ch0036	Chicken PCX (Podocalyxin) ELISA Kit
E-EL-Ch0038	Chicken TXLNa (Taxilin Alpha) ELISA Kit

E-EL-Ch0039	Chicken SDC4 (Syndecan 4) ELISA Kit
E-EL-Ch0040	Chicken Grb14 (Growth Factor Receptor Bound Protein 14) ELISA Kit
E-EL-Ch0041	Chicken FcγR3B (Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit
E-EL-Ch0043	Chicken CKLFSF8 (Chemokine Like Factor Superfamily 8) ELISA Kit
E-EL-Ch0044	Chicken HbH (Hemoglobin H) ELISA Kit
E-EL-Ch0045	Chicken ASGR2 (Asialoglycoprotein Receptor 2) ELISA Kit
E-EL-Ch0046	Chicken GUCA2A (Guanylate Cyclase Activator 2A) ELISA Kit
E-EL-Ch0047	Chicken TSP-4 (Thrombospondin-4) ELISA Kit
E-EL-Ch0048	Chicken CYSLTR1 (Cysteinyl Leukotriene Receptor 1) ELISA Kit
E-EL-Ch0049	Chicken f-Hb (Free Haemoglobin) ELISA Kit
E-EL-Ch0050	Chicken PCDHβ16 (Protocadherin Beta 16) ELISA Kit
E-EL-Ch0051	Chicken MLCK/MYLK (Myosin Light Chain Kinase) ELISA Kit
E-EL-Ch0052	Chicken PLA2G4D (Phospholipase A2, Group IV D) ELISA Kit
E-EL-Ch0053	Chicken CCND2 (Cyclin-D2) ELISA Kit
E-EL-Ch0054	Chicken APT (Abnormal Prothrombin) ELISA Kit
E-EL-Ch0055	Chicken PDGF-BB (Platelet Derived Growth Factor-BB) ELISA Kit
E-EL-Ch0058	Chicken GPER1 (G Protein Coupled Estrogen Receptor 1) ELISA Kit
E-EL-Ch0059	Chicken PDGFD (Platelet Derived Growth Factor D) ELISA Kit
E-EL-Ch0060	Chicken ADP/Acrp30 (Adiponectin) ELISA Kit
E-EL-Ch0061	Chicken GAL (Galanin) ELISA Kit
E-EL-Ch0062	Chicken GAS (N-Acetylgalactosamine 6-Sulfatase) ELISA Kit

E-EL-Ch0063	Chicken CKMB (Creatine Kinase MB Isoenzyme) ELISA Kit
E-EL-Ch0064	Chicken MEC/CCL28 (Mucosae Associated Epithelia Chemokine) ELISA Kit
E-EL-Ch0065	Chicken AGA (Anti-Glycolipid Antibody) ELISA Kit
E-EL-Ch0066	Chicken STAT6B (Signal Transducer and Activator of Transcription 6B) ELISA Kit
E-EL-Ch0067	Chicken HPR (Haptoglobin Related Protein) ELISA Kit
E-EL-Ch0069	Chicken GREM1 (Gremlin 1) ELISA Kit
E-EL-Ch0070	Chicken CPN10 (Chaperonin 10) ELISA Kit
E-EL-Ch0071	Chicken AGA (Anti-Gliadin Antibody) ELISA Kit
E-EL-Ch0072	Chicken ALCAM (Activated Leukocyte Cell Adhesion Molecule) ELISA Kit
E-EL-Ch0073	Chicken MAT (Matrilysin) ELISA Kit
E-EL-Ch0074	Chicken SYK (Spleen Tyrosine Kinase) ELISA Kit
E-EL-Ch0076	Chicken TP (Thymidine Phosphorylase) ELISA Kit
E-EL-Ch0077	Chicken GATA3 (GATA Binding Protein 3) ELISA Kit
E-EL-Ch0079	Chicken GPLD1 (Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Kit
E-EL-Ch0080	Chicken COL9α3 (Collagen Type IX Alpha 3) ELISA Kit
E-EL-Ch0081	Chicken LCN2 (Lipocalin 2) ELISA Kit
E-EL-Ch0082	Chicken LTβR (Lymphotoxin Beta Receptor) ELISA Kit
E-EL-Ch0083	Chicken CA125 (Carbohydrate Antigen 125) ELISA Kit
E-EL-Ch0084	Chicken CKMT1B (Creatine Kinase, Mitochondrial 1B) ELISA Kit
E-EL-Ch0085	Chicken PDGFC (Platelet Derived Growth Factor C) ELISA Kit
E-EL-Ch0086	Chicken SDH (Sorbitol Dehydrogenase) ELISA Kit

E-EL-Ch0087	Chicken PLC γ 2 (Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit
E-EL-Ch0088	Chicken MCP-4/CCL13 (Monocyte Chemotactic Protein 4) ELISA Kit
E-EL-Ch0089	Chicken ERK (Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-Ch0090	Chicken F XII a (Activated Coagulation Factor XII) ELISA Kit
E-EL-Ch0092	Chicken NUP35 (Nucleoporin 35kDa) ELISA Kit
E-EL-Ch0093	Chicken CLDN3 (Claudin 3) ELISA Kit
E-EL-Ch0094	Chicken CD40L/TNFSF5 (Cluster of Differentiation 40 Ligand) ELISA Kit
E-EL-Ch0096	Chicken LT-E4 (Leukotriene E4) ELISA Kit
E-EL-Ch0097	Chicken Neopterin ELISA Kit
E-EL-Ch0098	Chicken PKC γ (Protein Kinase C Gamma) ELISA Kit
E-EL-Ch0099	Chicken CAM (Calmodulin) ELISA Kit
E-EL-Ch0100	Chicken GST α 1 (Glutathione S Transferase Alpha 1) ELISA Kit
E-EL-Ch0101	Chicken FSP (Fibroblast Surface Protein) ELISA Kit
E-EL-Ch0102	Chicken GRO α /CXCL1/NAP3 (Growth Regulated Oncogene Alpha) ELISA Kit
E-EL-Ch0103	Chicken AQP-1 (Aquaporin 1) ELISA Kit
E-EL-Ch0104	Chicken HbC (Hemoglobin C) ELISA Kit
E-EL-Ch0105	Chicken SREBP (Sterol Regulatory Element Binding Proteins) ELISA Kit
E-EL-Ch0106	Chicken PAH (Phenylalanine Hydroxylase) ELISA Kit
E-EL-Ch0108	Chicken snRNP/Sm (Small Nuclear Ribonucleoprotein) ELISA Kit
E-EL-Ch0110	Chicken P-Cadherin (Cadherin, Placental) ELISA Kit
E-EL-Ch0111	Chicken BMPR2 (Bone Morphogenetic Protein Receptor II) ELISA Kit

E-EL-Ch0112	Chicken NA/NE (Noradrenaline/Norepinephrine) ELISA Kit
E-EL-Ch0113	Chicken NUP50 (Nucleoporin 50kDa) ELISA Kit
E-EL-Ch0114	Chicken FGL2 (Fibrinogen Like Protein 2) ELISA Kit
E-EL-Ch0115	Chicken α -SMA (α -Smooth muscle actin) ELISA Kit
E-EL-Ch0116	Chicken IGF-1 (Insulin Like Growth Factor 1) ELISA Kit
E-EL-Ch0117	Chicken MMP-1 (Matrix Metalloproteinase 1) ELISA Kit
E-EL-Ch0118	Chicken ApoA1 (Apolipoprotein A1) ELISA Kit
E-EL-Ch0119	Chicken FGF19 (Fibroblast Growth Factor 19) ELISA Kit
E-EL-Ch0120	Chicken IL-2 (Interleukin 2) ELISA Kit
E-EL-Ch0121	Chicken ANP (Atrial Natriuretic Peptide) ELISA Kit
E-EL-Ch0122	Chicken PDH β (Pyruvate Dehydrogenase Beta) ELISA Kit
E-EL-Ch0123	Chicken REG3 γ (Regenerating Islet Derived Protein 3 Gamma) ELISA Kit
E-EL-Ch0124	Chicken MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
E-EL-Ch0126	Chicken GDI1 (GDP Dissociation Inhibitor 1) ELISA Kit
E-EL-Ch0127	Chicken CLDN1 (Claudin 1) ELISA Kit
E-EL-Ch0128	Chicken Ank-B (Ankyrin B) ELISA Kit
E-EL-Ch0129	Chicken CA242 (Pancreatic Carcinoma Markers-CA242) ELISA Kit
E-EL-Ch0130	Chicken IFN- α (Interferon Alpha) ELISA Kit
E-EL-Ch0131	Chicken CK-13/KRT13 (Cytokeratin 13) ELISA Kit
E-EL-Ch0132	Chicken PKR2 (Prokineticin Receptor 2) ELISA Kit
E-EL-Ch0133	Chicken FAK (Focal Adhesion Kinase) ELISA Kit

E-EL-Ch0134	Chicken BACE2 (β -site APP Cleaving Enzyme 2) ELISA Kit
E-EL-Ch0135	Chicken GCSFR (Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit
E-EL-Ch0136	Chicken CK-HMW (Cytokeratin High Molecular Weight) ELISA Kit
E-EL-Ch0137	Chicken ELMO2 (Engulfment And Cell Motility 2) ELISA Kit
E-EL-Ch0139	Chicken TFPI (Tissue Factor Pathway Inhibitor) ELISA Kit
E-EL-Ch0141	Chicken GPX3 (Glutathione Peroxidase 3, Plasma) ELISA Kit
E-EL-Ch0143	Chicken PEDF (Pigment Epithelium Derived Factor) ELISA Kit
E-EL-Ch0144	Chicken REG1 α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit
E-EL-Ch0145	Chicken PP13 (Placental Protein13) ELISA Kit
E-EL-Ch0146	Chicken EPCR (Endothelial Protein C Receptor) ELISA Kit
E-EL-Ch0147	Chicken PHLDA2 (Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit
E-EL-Ch0148	Chicken CAPF (Cuticular Active Peptide Factor) ELISA Kit
E-EL-Ch0149	Chicken SNCy (Synuclein Gamma) ELISA Kit
E-EL-Ch0150	Chicken Ub (Ubiquitin) ELISA Kit
E-EL-Ch0151	Chicken PMP3 (Peroxisomal Membrane Protein 3) ELISA Kit
E-EL-Ch0152	Chicken OIT3 (Oncoprotein Induced Transcript 3) ELISA Kit
E-EL-Ch0153	Chicken TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
E-EL-Ch0154	Chicken RAPGEF1 (Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit
E-EL-Ch0155	Chicken EPYC (Epiphycan) ELISA Kit
E-EL-Ch0156	Chicken PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit
E-EL-Ch0157	Chicken THR α (Thyroid Hormone Receptor Alpha) ELISA Kit

E-EL-Ch0159	Chicken ACTN3 (Actinin Alpha 3) ELISA Kit
E-EL-Ch0160	Chicken GLP-1 (Glucagon Like Peptide 1) ELISA Kit
E-EL-Ch0161	Chicken APSA (Anti-Phosphatidyl Serine Antibody) ELISA Kit
E-EL-Ch0162	Chicken CCSA-3 (Colon Cancer Specific Antigen-3) ELISA Kit
E-EL-Ch0163	Chicken EMMPRIN/CD147 (Extracellular Matrix Metalloproteinase Inducer) ELISA Kit
E-EL-Ch0165	Chicken GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit
E-EL-Ch0166	Chicken ACT α 2 (Actin Alpha 2,Smooth Muscle) ELISA Kit
E-EL-Ch0167	Chicken GasR (Gastrin Receptor) ELISA Kit
E-EL-Ch0168	Chicken DLAT (Dihydrolipoyl Transacetylase) ELISA Kit
E-EL-Ch0169	Chicken BlysR (B Lymphocyte Stimulator Receptor) ELISA Kit
E-EL-Ch0170	Chicken HSP-40 (Heat Shock Protein 40) ELISA Kit
E-EL-Ch0171	Chicken PICP (Procollagen I C-terminal Propeptide) ELISA Kit
E-EL-Ch0173	Chicken STAT2 (Signal Transducer and Activator of Transcription 2) ELISA Kit
E-EL-Ch0175	Chicken EMA (Epithelial Membrane Antigen) ELISA Kit
E-EL-Ch0178	Chicken CDK2 (Cyclin Dependent Kinase 2) ELISA Kit
E-EL-Ch0179	Chicken TNPO1 (Transportin 1) ELISA Kit
E-EL-Ch0180	Chicken CDK9 (Cyclin Dependent Kinase 9) ELISA Kit
E-EL-Ch0181	Chicken α HBDH (α -Hydroxybutyrate Dehydrogenase) ELISA Kit
E-EL-Ch0182	Chicken PRKG1 (Protein Kinase, cGMP Dependent Type I) ELISA Kit
E-EL-Ch0183	Chicken IGF-2 (Insulin Like Growth Factor 2) ELISA Kit
E-EL-Ch0184	Chicken F13 (Coagulation Factor X III) ELISA Kit

E-EL-Ch0185	Chicken δ SIP- α (Delta Sleep Inducing Peptide α) ELISA Kit
E-EL-Ch0186	Chicken NRP1 (Neuropilin 1) ELISA Kit
E-EL-Ch0187	Chicken STAT5A (Signal Transducer and Activator of Transcription 5A) ELISA Kit
E-EL-Ch0188	Chicken BCDF (B-Cell Differentiation Factor) ELISA Kit
E-EL-Ch0189	Chicken MASP (Mannan Binding Lectin Associated Serine Protease) ELISA Kit
E-EL-Ch0190	Chicken HYOU1 (Hypoxia Up Regulated 1) ELISA Kit
E-EL-Ch0191	Chicken MLZE (Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit
E-EL-Ch0192	Chicken PDGF (Platelet-Derived Growth Factor) ELISA Kit
E-EL-Ch0193	Chicken CaN (Calcineurin) ELISA Kit
E-EL-Ch0194	Chicken CPR (Cytochrome P450 Reductase) ELISA Kit
E-EL-Ch0196	Chicken CDKN3 (Cyclin Dependent Kinase Inhibitor 3) ELISA Kit
E-EL-Ch0197	Chicken CXCR1 (CXC-Chemokine Receptor 1) ELISA Kit
E-EL-Ch0198	Chicken Ntn1 (Netrin 1) ELISA Kit
E-EL-Ch0200	Chicken AIRA (Anti-Insulin Receptor Antibody) ELISA Kit
E-EL-Ch0201	Chicken γ GT1 (Gamma Glutamyltransferase 1) ELISA Kit
E-EL-Ch0202	Chicken GAL12 (Galectin 12) ELISA Kit
E-EL-Ch0203	Chicken BMF (Bcl-2 Modifying Factor) ELISA Kit
E-EL-Ch0204	Chicken GMF β (Glia Maturation Factor, Beta) ELISA Kit
E-EL-Ch0205	Chicken EPHA3 (Ephrin Type A Receptor 3) ELISA Kit
E-EL-Ch0206	Chicken α -CTx (Alpha Crosslaps) ELISA Kit
E-EL-Ch0207	Chicken CCND3 (Cyclin-D3) ELISA Kit

E-EL-Ch0208	Chicken DPT (Dermatopontin) ELISA Kit
E-EL-Ch0209	Chicken ADAM8 (A Disintegrin And Metalloprotease 8) ELISA Kit
E-EL-Ch0210	Chicken NEF3 (Neurofilament 3) ELISA Kit
E-EL-Ch0211	Chicken LUM (Lumican) ELISA Kit
E-EL-Ch0212	Chicken GP-73 (Golgi Protein 73) ELISA Kit
E-EL-Ch0213	Chicken COLEC11 (Collectin Sub-Family Member 11) ELISA Kit
E-EL-Ch0214	Chicken PDXK (Pyridoxal Kinase) ELISA Kit
E-EL-Ch0215	Chicken TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-Ch0216	Chicken GST ω 1 (Glutathione S Transferase Omega 1) ELISA Kit
E-EL-Ch0218	Chicken GDF11 (Growth Differentiation Factor 11) ELISA Kit
E-EL-Ch0219	Chicken CD30L (Cluster of Differentiation 30 Lligand) ELISA Kit
E-EL-Ch0220	Chicken TGM2 (Transglutaminase 2, Tissue) ELISA Kit
E-EL-Ch0221	Chicken FBL (Fibrillarin) ELISA Kit
E-EL-Ch0222	Chicken CA153 (Carbohydrate Antigen 153) ELISA Kit
E-EL-Ch0223	Chicken RGS4 (Regulator of G Protein Signaling 4) ELISA Kit
E-EL-Ch0224	Chicken C5b (Complement Fragment 5b) ELISA Kit
E-EL-Ch0225	Chicken GRO γ (Growth Regulated Oncogene Gamma) ELISA Kit
E-EL-Ch0227	Chicken RUNX2 (Runt-Related Transcription Factor 2) ELISA Kit
E-EL-Ch0228	Chicken IL-6(Interleukin 6) ELISA Kit
E-EL-Ch0229	Chicken FZD8 (Frizzled Homolog 8) ELISA Kit
E-EL-Ch0230	Chicken sgp130 (soluble Glucoprotein 130) ELISA Kit

E-EL-Ch0231	Chicken PKN1 (Protein Kinase N1) ELISA Kit
E-EL-Ch0232	Chicken β -LPH (β -Lipotropic Hormone) ELISA Kit
E-EL-Ch0233	Chicken Smad7/MADH7 (Mothers Against Decapentaplegic Homolog 7) ELISA Kit
E-EL-Ch0234	Chicken RHOA (Ras Homolog Gene Family, Member A) ELISA Kit
E-EL-Ch0235	Chicken TU3A (TU3A Protein) ELISA Kit
E-EL-Ch0236	Chicken MAN1A1 (Mannosidase Alpha Class 1A Member 1) ELISA Kit
E-EL-Ch0237	Chicken EPHA1 (Ephrin Type A Receptor 1) ELISA Kit
E-EL-Ch0238	Chicken CT-1 (Cardiotrophin 1) ELISA Kit
E-EL-Ch0239	Chicken GP (Anti-Glycoprotein Antibody) ELISA Kit
E-EL-Ch0240	Chicken NAP (Neutrophil Alkaline Phosphatase) ELISA Kit
E-EL-Ch0242	Chicken GP5 (Glycoprotein V, Platelet) ELISA Kit
E-EL-Ch0243	Chicken GP4 (Platelet Membrane Glycoprotein IV) ELISA Kit
E-EL-Ch0245	Chicken ON (Osteonectin) ELISA Kit
E-EL-Ch0246	Chicken PPAR- γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
E-EL-Ch0247	Chicken Gal (α -Galactoyl) ELISA Kit
E-EL-Ch0248	Chicken Anti-Surv (Anti-Survivin Antibody) ELISA Kit
E-EL-Ch0249	Chicken D2D (D-Dimer) ELISA Kit
E-EL-Ch0251	Chicken MPA (Mycophenolic Acid) ELISA Kit
E-EL-Ch0252	Chicken POMC (Pro-Opiomelanocortin) ELISA Kit
E-EL-Ch0253	Chicken ApoD (Apolipoprotein D) ELISA Kit
E-EL-Ch0254	Chicken IL-13 (Interleukin 13) ELISA Kit

E-EL-Ch0255	Chicken NUP153 (Nucleoporin 153kDa) ELISA Kit
E-EL-Ch0256	Chicken CKLFSF7 (Chemokine Like Factor Superfamily 7) ELISA Kit
E-EL-Ch0257	Chicken SRF/IF (Skin Reactive Factor/Inflammatory Factor) ELISA Kit
E-EL-Ch0258	Chicken LT β (Lymphotoxin Beta) ELISA Kit
E-EL-Ch0259	Chicken CDK7 (Cyclin Dependent Kinase 7) ELISA Kit
E-EL-Ch0260	Chicken TrxR1 (Thioredoxin Reductase 1) ELISA Kit
E-EL-Ch0261	Chicken MTX (Methotrexate) ELISA Kit
E-EL-Ch0263	Chicken GKRP (Glucokinase Regulatory Protein) ELISA Kit
E-EL-Ch0264	Chicken TACR1 (Tachykinin Receptor 1) ELISA Kit
E-EL-Ch0265	Chicken C-P (C-Peptide) ELISA Kit
E-EL-Ch0266	Chicken GTF2F1 (General Transcription Factor IIF, Polypeptide 1) ELISA Kit
E-EL-Ch0267	Chicken IDE (Insulin Degrading Enzyme) ELISA Kit
E-EL-Ch0268	Chicken CYPA (Cyclophilin A) ELISA Kit
E-EL-Ch0269	Chicken FABP1 (Fatty Acid Binding Protein 1, Liver) ELISA Kit
E-EL-Ch0270	Chicken ECFR (Eosinophil Chemotactic Factor Receptor) ELISA Kit
E-EL-Ch0271	Chicken SHBG (Sex Hormone Binding Globulin) ELISA Kit
E-EL-Ch0272	Chicken D1R/DRD1 (Dopamine Receptor D1) ELISA Kit
E-EL-Ch0273	Chicken MIF (Macrophage Migration Inhibitory Factor) ELISA Kit
E-EL-Ch0275	Chicken PLC γ 1 (Phospholipase C Gamma 1) ELISA Kit
E-EL-Ch0276	Chicken MYC (V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit
E-EL-Ch0277	Chicken CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit

E-EL-Ch0278	Chicken 6-OHDA (6-Hydroxydopamine) ELISA Kit
E-EL-Ch0279	Chicken ARHGDI α (Rho GDP Dissociation Inhibitor Alpha) ELISA Kit
E-EL-Ch0280	Chicken MSTN (Myostatin) ELISA Kit
E-EL-Ch0281	Chicken DBP (Vitamin D Binding Protein) ELISA Kit
E-EL-Ch0282	Chicken CGRP (Calcitonin Gene Related Peptide) ELISA Kit
E-EL-Ch0283	Chicken TG (Thyroglobulin) ELISA Kit
E-EL-Ch0284	Chicken TLN (Talin) ELISA Kit
E-EL-Ch0285	Chicken ANG I (Angiotensin I) ELISA Kit
E-EL-Ch0286	Chicken LN (Laminin) ELISA Kit
E-EL-Ch0288	Chicken SEMA4C (Semaphorin 4C) ELISA Kit
E-EL-Ch0289	Chicken LRP (Lung Resistance-Related Protein) ELISA Kit
E-EL-Ch0290	Chicken PDGFR α (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit
E-EL-Ch0291	Chicken TDT (Terminal Deoxynucleotidyl Transferase) ELISA Kit
E-EL-Ch0292	Chicken PKBG (Protein Kinase B Gamma) ELISA Kit
E-EL-Ch0293	Chicken CA5B (Carbonic Anhydrase VB, Mitochondrial) ELISA Kit
E-EL-Ch0294	Chicken Fc α R (Fc Fragment of IgA Receptor) ELISA Kit
E-EL-Ch0295	Chicken COL14 (Collagen Type XIV) ELISA Kit
E-EL-Ch0296	Chicken C1q (Complement 1q) ELISA Kit
E-EL-Ch0298	Chicken C1r (Complement Component 1R) ELISA Kit
E-EL-Ch0299	Chicken ACA-IgA (Anti-Cardiolipin Antibody IgA) ELISA Kit
E-EL-Ch0301	Chicken BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit

E-EL-Ch0302	Chicken GDH/GLDH (Glutamate dehydrogenase) ELISA Kit
E-EL-Ch0303	Chicken DAO (Diamine Oxidase) ELISA Kit
E-EL-Ch0304	Chicken MMP-25 (Matrix Metalloproteinase 25) ELISA Kit
E-EL-Ch0305	Chicken GDF5 (Growth Differentiation Factor 5) ELISA Kit
E-EL-Ch0306	Chicken PACAPR (Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit
E-EL-Ch0307	Chicken GEF (Guanylate Exchange Factor) ELISA Kit
E-EL-Ch0308	Chicken PAQR7 (Progestin and AdipoQ Receptor Family Member VII) ELISA Kit
E-EL-Ch0309	Chicken NUP214 (Nucleoporin 214kDa) ELISA Kit
E-EL-Ch0310	Chicken XCR1 (XC-Chemokine Receptor 1) ELISA Kit
E-EL-Ch0311	Chicken AATF (Apoptosis Antagonizing Transcription Factor) ELISA Kit
E-EL-Ch0312	Chicken DAT (Dopamine Transporter) ELISA Kit
E-EL-Ch0313	Chicken FGF6 (Fibroblast Growth Factor 6) ELISA Kit
E-EL-Ch0314	Chicken CK-19/KRT19 (Cytokeratin 19) ELISA Kit
E-EL-Ch0315	Chicken NRG-3 (Neuregulin 3) ELISA Kit
E-EL-Ch0316	Chicken ACD (Adrenocortical Dysplasia Homolog) ELISA Kit
E-EL-Ch0317	Chicken ARO (Aromatase) ELISA Kit
E-EL-Ch0318	Chicken Cyt-C (Cytochrome C) ELISA Kit
E-EL-Ch0319	Chicken CT (Calcitonin) ELISA Kit
E-EL-Ch0320	Chicken sICAM-1 (Soluble Intercellular Adhesion Molecule-1) ELISA Kit
E-EL-Ch0321	Chicken Pro-MMP-1 (pro-Matrix Metalloproteinase 1) ELISA Kit
E-EL-Ch0322	Chicken ANGPTL6 (Angiopoietin Like Protein 6) ELISA Kit

E-EL-Ch0323	Chicken GAP (GTPase Activating Protein) ELISA Kit
E-EL-Ch0324	Chicken TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
E-EL-Ch0325	Chicken PG (Proteoglycan) ELISA Kit
E-EL-Ch0326	Chicken LDHC (Lactate Dehydrogenase C) ELISA Kit
E-EL-Ch0327	Chicken D4R/DRD4 (Dopamine Receptor D4) ELISA Kit
E-EL-Ch0328	Chicken Cyr61 (Cysteine Rich Protein, Angiogenic Inducer 61) ELISA Kit
E-EL-Ch0329	Chicken GM (Gastric Mucin) ELISA Kit
E-EL-Ch0330	Chicken CD28 (Cluster of Differentiation 28) ELISA Kit
E-EL-Ch0331	Chicken PTHR2 (Parathyroid Hormone Receptor 2) ELISA Kit
E-EL-Ch0333	Chicken RAG-1 (Recombination Activating Gene 1) ELISA Kit
E-EL-Ch0334	Chicken DCX (Doublecortin) ELISA Kit
E-EL-Ch0335	Chicken CLC (Charcot Leyden Crystal Protein) ELISA Kit
E-EL-Ch0336	Chicken PSA/KLK3 (Prostate Specific Antigen) ELISA Kit
E-EL-Ch0337	Chicken CgA (Chromogranin A/Pancreastatin) ELISA Kit
E-EL-Ch0338	Chicken ADAM9 (A Disintegrin And Metalloprotease 9) ELISA Kit
E-EL-Ch0339	Chicken ECM1 (Extracellular Matrix Protein 1) ELISA Kit
E-EL-Ch0340	Chicken MIgD (Membrane IgD) ELISA Kit
E-EL-Ch0341	Chicken COL6 (Collagen Type VI) ELISA Kit
E-EL-Ch0342	Chicken TNNT1 (Troponin T Type 1, Slow Skeletal) ELISA Kit
E-EL-Ch0343	Chicken CAD (Caspase Activated Deoxyribonuclease) ELISA Kit
E-EL-Ch0344	Chicken S100A8 (S100 Calcium Binding Protein A8/Calgranulin A) ELISA Kit

E-EL-Ch0345	Chicken LRP6 (Low Density Lipoprotein Receptor Related Protein) ELISA Kit
E-EL-Ch0346	Chicken VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
E-EL-Ch0347	Chicken PK (Pyruvate Kinase, Liver And RBC) ELISA Kit
E-EL-Ch0348	Chicken MMP-11 (Matrix Metalloproteinase 11) ELISA Kit
E-EL-Ch0349	Chicken CG (Chorionic Gonadotropin) ELISA Kit
E-EL-Ch0350	Chicken IGF2R (Insulin Like Growth Factor 2 Receptor) ELISA Kit
E-EL-Ch0351	Chicken BMPR1A (Bone Morphogenetic Protein Receptor, type IA) ELISA Kit
E-EL-Ch0352	Chicken GL β (Galactosidase Beta) ELISA Kit
E-EL-Ch0353	Chicken PTPRO (Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
E-EL-Ch0354	Chicken IL1R1 (Interleukin 1 Receptor Type I) ELISA Kit
E-EL-Ch0355	Chicken E-Cad (E-Cadherin) ELISA Kit
E-EL-Ch0356	Chicken Tgtp (T-cell Specific Gtpase) ELISA Kit
E-EL-Ch0357	Chicken C/EBP β (CCAAT/Enhancer Binding Protein Beta) ELISA Kit
E-EL-Ch0358	Chicken LAP (Leucine Aminopeptidase) ELISA Kit
E-EL-Ch0359	Chicken GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit
E-EL-Ch0360	Chicken GRN (Granulin) ELISA Kit
E-EL-Ch0362	Chicken ADCY3 (Adenylate Cyclase 3) ELISA Kit
E-EL-Ch0363	Chicken NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
E-EL-Ch0364	Chicken GST α 5 (Glutathione S Transferase Alpha 5) ELISA Kit
E-EL-Ch0365	Chicken IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit
E-EL-Ch0366	Chicken IL-16 (Interleukin 16) ELISA Kit

E-EL-Ch0367	Chicken IgM (Immunoglobulin M) ELISA Kit
E-EL-Ch0368	Chicken LFA-2/CD2 (Lymphocyte Function Associated Antigen 2) ELISA Kit
E-EL-Ch0369	Chicken CgB (Chromogranin B) ELISA Kit
E-EL-Ch0370	Chicken PLAP/ALPP (Placental Alkaline Phosphatase) ELISA Kit
E-EL-Ch0371	Chicken AP36 (Apelin 36) ELISA Kit
E-EL-Ch0373	Chicken M-AChR M2 (Muscarinic Acetylcholine Receptor M2) ELISA Kit
E-EL-Ch0374	Chicken IgG3 (Immunoglobulin G3) ELISA Kit
E-EL-Ch0375	Chicken TAF (TATA Box Binding Protein/TBP-Associated Factors) ELISA Kit
E-EL-Ch0377	Chicken GATA1 (GATA Binding Protein 1) ELISA Kit
E-EL-Ch0378	Chicken TFR/CD71 (Transferrin Receptor) ELISA Kit
E-EL-Ch0379	Chicken GHAb (Anti-Growth Hormone Antibody) ELISA Kit
E-EL-Ch0380	Chicken KNG1 (Kininogen 1) ELISA Kit
E-EL-Ch0381	Chicken UBR4 (Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit
E-EL-Ch0382	Chicken COL6α1 (Collagen Type VI Alpha 1) ELISA Kit
E-EL-Ch0383	Chicken BANP (BTG3 Associated Nuclear Protein) ELISA Kit
E-EL-Ch0384	Chicken ALK (Anaplastic Lymphoma Kinase) ELISA Kit
E-EL-Ch0385	Chicken Aβ1-40 (Amyloid Beta Peptide 1-40) ELISA Kit
E-EL-Ch0386	Chicken GLTP (Glycolipid Transfer Protein) ELISA Kit
E-EL-Ch0388	Chicken GSTθ2/GSTt2 (Glutathione S Transferase Theta 2) ELISA Kit
E-EL-Ch0389	Chicken PK2 (Prokineticin 2) ELISA Kit
E-EL-Ch0390	Chicken NTRK2 (Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit

E-EL-Ch0391	Chicken FLNC (Filamin C, Gamma) ELISA Kit
E-EL-Ch0393	Chicken EXT1 (Exostoses 1) ELISA Kit
E-EL-Ch0396	Chicken VIPR1 (Vasoactive Intestinal Peptide Receptor) ELISA Kit
E-EL-Ch0398	Chicken CD99 (Cluster of Differentiation 99) ELISA Kit
E-EL-Ch0399	Chicken ARRβ2 (Arrestin Beta 2) ELISA Kit
E-EL-Ch0401	Chicken DPP10 (Dipeptidyl Peptidase X) ELISA Kit
E-EL-Ch0402	Chicken Smad4/MADH4 (Mothers Against Decapentaplegic Homolog 4)ELISA Kit
E-EL-Ch0403	Chicken CFH (Complement Factor H) ELISA Kit
E-EL-Ch0404	Chicken TWEAK (Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit
E-EL-Ch0405	Chicken PLP1 (Proteolipid Protein 1, Myelin) ELISA Kit
E-EL-Ch0406	Chicken CaSR (Calcium Sensing Receptor) ELISA Kit
E-EL-Ch0407	Chicken AChRab (Acetylcholine Receptor Antibody) ELISA Kit
E-EL-Ch0408	Chicken ERN1 (Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit
E-EL-Ch0409	Chicken TTR (Transthyretin) ELISA Kit
E-EL-Ch0410	Chicken CSF (Colony Stimulating Factor) ELISA Kit
E-EL-Ch0411	Chicken AP62A (Anti-P62 Antibody) ELISA Kit
E-EL-Ch0412	Chicken M-AChR M3 (Muscarinic Acetylcholine Receptor M3) ELISA Kit
E-EL-Ch0413	Chicken NKB (Neurokinins B) ELISA Kit
E-EL-Ch0414	Chicken PRF1 (Perforin 1/Pore-Forming Protein) ELISA Kit
E-EL-Ch0415	Chicken NUP155 (Nucleoporin 155kDa) ELISA Kit
E-EL-Ch0416	Chicken COL8α1 (Collagen Type VIII Alpha 1) ELISA Kit

E-EL-Ch0417	Chicken LGALS3 (Lectin-Galactoside Binding Soluble 3) ELISA Kit
E-EL-Ch0418	Chicken ATF-1 (Activating Transcription Factor 1) ELISA Kit
E-EL-Ch0419	Chicken ANNA1(Anti-Neuronal Nuclear Antibody 1) ELISA Kit
E-EL-Ch0420	Chicken PCPE1 (Procollagen C Proteinase Enhancer 1) ELISA Kit
E-EL-Ch0421	Chicken TP/TXA2R (Thromboxane Receptor) ELISA Kit
E-EL-Ch0422	Chicken MIgM (Membrane IgM) ELISA Kit
E-EL-Ch0423	Chicken CLA/CTAGE (Cutaneous Lymphocyte Associated Antigen) ELISA Kit
E-EL-Ch0424	Chicken TS/TYMS (Thymidylate Synthetase) ELISA Kit
E-EL-Ch0425	Chicken β TG/PBP/CXCL7/NAP2 (Thromboglobulin, Beta) ELISA Kit
E-EL-Ch0426	Chicken ATM (Ataxia Telangiectasia Mutated) ELISA Kit
E-EL-Ch0427	Chicken SAA2 (Serum Amyloid A2) ELISA Kit
E-EL-Ch0428	Chicken SDF-2 (Stromal Cell Derived Factor 2) ELISA Kit
E-EL-Ch0430	Chicken γ -ECS (γ -Glutamyl Systeine Synthetase) ELISA Kit
E-EL-Ch0431	Chicken PTD (Pentosidine) ELISA Kit
E-EL-Ch0432	Chicken ME20M (Melanoma Associated ME20) ELISA Kit
E-EL-Ch0433	Chicken PTH (Parathyroid Hormone) ELISA Kit
E-EL-Ch0434	Chicken BMP-15 (Bone Morphogenetic Protein 15) ELISA Kit
E-EL-Ch0435	Chicken ADA (Adenosine Deaminase) ELISA Kit
E-EL-Ch0436	Chicken AFA (Anti-Filaggrin Antibody) ELISA Kit
E-EL-Ch0437	Chicken α ZGP1 (Alpha-2-Glycoprotein 1, Zine Binding) ELISA Kit
E-EL-Ch0438	Chicken TI-Ag (Thymus Independent Antigen) ELISA Kit

E-EL-Ch0440	Chicken FGFR3 (Fibroblast Growth Factor Receptor 3) ELISA Kit
E-EL-Ch0441	Chicken CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit
E-EL-Ch0442	Chicken PGF2 α (Prostaglandin F2 α) ELISA Kit
E-EL-Ch0443	Chicken PRSS (Serine Protease) ELISA Kit
E-EL-Ch0444	Chicken CD97 (Cluster of Differentiation 97) ELISA Kit
E-EL-Ch0445	Chicken CKB (Creatine Kinase, Brain) ELISA Kit
E-EL-Ch0447	Chicken GPR37 (G Protein Coupled Receptor 37) ELISA Kit
E-EL-Ch0449	Chicken IL-9 (Interleukin 9) ELISA Kit
E-EL-Ch0450	Chicken AGC (Aggrecan) ELISA Kit
E-EL-Ch0451	Chicken ANEA (Anti-Neutrophil Elastase Antibody) ELISA Kit
E-EL-Ch0452	Chicken C2 (Complement Component 2) ELISA Kit
E-EL-Ch0453	Chicken GUK1 (Guanylate Kinase 1) ELISA Kit
E-EL-Ch0454	Chicken SIGLEC1 (Sialic Acid Binding Ig Like Lectin 1) ELISA Kit
E-EL-Ch0455	Chicken TXB2 (Thromboxane B2) ELISA Kit
E-EL-Ch0456	Chicken LTA4H (Leukotriene A4 Hydrolase) ELISA Kit
E-EL-Ch0457	Chicken PRCP (Prolylcarboxypeptidase) ELISA Kit
E-EL-Ch0458	Chicken Bcl2L11 (Bcl-2 Like Protein 11) ELISA Kit
E-EL-Ch0459	Chicken SAP25 (Synaptosome Associated Protein 25) ELISA Kit
E-EL-Ch0461	Chicken NOX5 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit
E-EL-Ch0462	Chicken ADIPOR1 (Adiponectin Receptor 1) ELISA Kit
E-EL-Ch0464	Chicken GHRH (Growth Hormone Releasing Hormone) ELISA Kit

E-EL-Ch0466	Chicken MRC2 (Mannose Receptor C Type 2) ELISA Kit
E-EL-Ch0467	Chicken BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit
E-EL-Ch0468	Chicken PRB (Retinoblastoma Tumor Suppressor Protein) ELISA Kit
E-EL-Ch0469	Chicken IL-23 p40(Interleukin 23) ELISA Kit
E-EL-Ch0470	Chicken SIGLEC3 (Sialic Acid Binding Ig Like Lectin 3) ELISA Kit
E-EL-Ch0471	Chicken GAL8 (Galectin 8) ELISA Kit
E-EL-Ch0472	Chicken Dyn (Big Dynorphin) ELISA Kit
E-EL-Ch0473	Chicken MT1E (Metallothionein 1E) ELISA Kit
E-EL-Ch0474	Chicken VLDLR (Very Low Density Lipoprotein Receptor) ELISA Kit
E-EL-Ch0476	Chicken TF (Transfer Factor) ELISA Kit
E-EL-Ch0477	Chicken MUC1 (Mucin 1) ELISA Kit
E-EL-Ch0478	Chicken CCA (Colon Cancer Antigen) ELISA Kit
E-EL-Ch0479	Chicken CTSC (Cathepsin C) ELISA Kit
E-EL-Ch0480	Chicken VF (Visfatin) ELISA Kit
E-EL-Ch0481	Chicken PRG4 (Proteoglycan 4) ELISA Kit
E-EL-Ch0483	Chicken GDF1 (Growth Differentiation Factor 1) ELISA Kit
E-EL-Ch0484	Chicken MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit
E-EL-Ch0485	Chicken TAA (Tumor Associated Antigen) ELISA Kit
E-EL-Ch0486	Chicken TPS (Tissue Polypeptide Specific Antigen) ELISA Kit
E-EL-Ch0487	Chicken TAF1 (TATA Box Binding Protein Associated Factor 1) ELISA Kit
E-EL-Ch0488	Chicken GDF15 (Growth Differentiation Factor 15) ELISA Kit

E-EL-Ch0489	Chicken RARS (Arginyl tRNA Synthetase) ELISA Kit
E-EL-Ch0490	Chicken TAF2 (TATA Box Binding Protein Associated Factor 2) ELISA Kit
E-EL-Ch0491	Chicken NSE (Neuron-Specific Enolase) ELISA Kit
E-EL-Ch0492	Chicken CCRL1 (Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit
E-EL-Ch0493	Chicken COL3 (Collagen Type III) ELISA Kit
E-EL-Ch0494	Chicken GHITM (Growth Hormone Inducible Transmembrane Protein) ELISA Kit
E-EL-Ch0495	Chicken GZMK (Granzyme K) ELISA Kit
E-EL-Ch0496	Chicken CCL14 (Chemokine C-C-Motif Ligand 14) ELISA Kit
E-EL-Ch0497	Chicken PRDX2 (Peroxiredoxin 2) ELISA Kit
E-EL-Ch0499	Chicken ETRB (Endothelin Receptor B) ELISA Kit
E-EL-Ch0500	Chicken GPC3 (Glypican 3) ELISA Kit
E-EL-Ch0501	Chicken NUP85 (Nucleoporin 85kDa) ELISA Kit
E-EL-Ch0503	Chicken PLD2 (Phospholipase D2) ELISA Kit
E-EL-Ch0504	Chicken CVNPF (Cobra Venom Neuronal Protective Factor) ELISA Kit
E-EL-Ch0505	Chicken MMP-17 (Matrix Metalloproteinase 17) ELISA Kit
E-EL-Ch0506	Chicken ProGRP (Pro-Gastrin Releasing Peptide) ELISA Kit
E-EL-Ch0507	Chicken EPOR (Erythropoietin Receptor) ELISA Kit
E-EL-Ch0508	Chicken ALDH (Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit
E-EL-Ch0509	Chicken NT-4 (Neurotrophin 4) ELISA Kit
E-EL-Ch0510	Chicken FGFR4 (Fibroblast Growth Factor Receptor 4) ELISA Kit
E-EL-Ch0511	Chicken ApoB (Apolipoprotein B) ELISA Kit

E-EL-Ch0512	Chicken 26S-PSM (26S Proteasome) ELISA Kit
E-EL-Ch0513	Chicken PS (Phosphatidylserine) ELISA Kit
E-EL-Ch0514	Chicken LDL (Low Density Lipoprotein) ELISA Kit
E-EL-Ch0515	Chicken TLA (Thymus Leukemia Antigen) ELISA Kit
E-EL-Ch0516	Chicken β MSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-Ch0517	Chicken THAA (Thyroid Hormone Autoantibody) ELISA Kit
E-EL-Ch0518	Chicken SEPP1 (Selenoprotein P) ELISA Kit
E-EL-Ch0520	Chicken ACT α 1 (Actin Alpha 1, Skeletal Muscle) ELISA Kit
E-EL-Ch0521	Chicken PHC (Primary Hepatic Carcinoma) ELISA Kit
E-EL-Ch0522	Chicken ANG (Angiogenin) ELISA Kit
E-EL-Ch0523	Chicken NCAM/CD56 (Neural Cell Adhesion Molecule) ELISA Kit
E-EL-Ch0524	Chicken GAL14 (Galectin 14) ELISA Kit
E-EL-Ch0525	Chicken ALT (Alanine Transaminase) ELISA Kit
E-EL-Ch0526	Chicken Collagenase I ELISA Kit
E-EL-Ch0528	Chicken MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
E-EL-Ch0529	Chicken VacA (Vacuolating Cytotoxin A) ELISA Kit
E-EL-Ch0530	Chicken UBE1L (Ubiquitin Activating Enzyme E1 Like Protein) ELISA Kit
E-EL-Ch0531	Chicken PI3K (Phosphoinositide-3-Kinase) ELISA Kit
E-EL-Ch0532	Chicken PRDX1 (Peroxiredoxin 1) ELISA Kit
E-EL-Ch0533	Chicken UBE2C (Ubiquitin Conjugating Enzyme E2C) ELISA Kit
E-EL-Ch0534	Chicken MMP-24 (Matrix Metalloproteinase 24) ELISA Kit

E-EL-Ch0535	Chicken MCSFR (Macrophage Colony Stimulating Factor Receptor) ELISA Kit
E-EL-Ch0536	Chicken ACVA (Activin A) ELISA Kit
E-EL-Ch0537	Chicken MYO (Myosin) ELISA Kit
E-EL-Ch0538	Chicken EP2 (Prostaglandin E Receptor 2) ELISA Kit
E-EL-Ch0539	Chicken IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-Ch0540	Chicken ANG II R2/AGTR2 (Angiotensin II Receptor 2) ELISA Kit
E-EL-Ch0541	Chicken TUB α 3C (Tubulin Alpha 3c)ELISA Kit
E-EL-Ch0542	Chicken MSLN (Mesothelin) ELISA Kit
E-EL-Ch0543	Chicken CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit
E-EL-Ch0545	Chicken NUP188 (Nucleoporin 188kDa) ELISA Kit
E-EL-Ch0546	Chicken CKLFSF4 (Chemokine Like Factor Superfamily 4) ELISA Kit
E-EL-Ch0547	Chicken ANG2 (Angiopoietin 2) ELISA Kit
E-EL-Ch0548	Chicken C4BP β (C4 Binding Protein Beta) ELISA Kit
E-EL-Ch0550	Chicken AQP-3 (Aquaporin 3) ELISA Kit
E-EL-Ch0551	Chicken COL11 α 1 (Collagen Type XI Alpha 1) ELISA Kit
E-EL-Ch0552	Chicken FPB (Fibrinopeptide B) ELISA Kit
E-EL-Ch0553	Chicken CTNN α 1 (Catenin, Alpha 1) ELISA Kit
E-EL-Ch0554	Chicken CALD (Caldesmon) ELISA Kit
E-EL-Ch0555	Chicken HBF (Fetal Hemoglobin) ELISA Kit
E-EL-Ch0556	Chicken PCIII (Procollagen III) ELISA Kit
E-EL-Ch0557	Chicken CDT (Carbohydrate Deficient Transferrin) ELISA Kit

E-EL-Ch0558	Chicken FOXO4 (Forkhead Box Protein O4) ELISA Kit
E-EL-Ch0560	Chicken COL6α3 (Collagen Type VI Alpha 3) ELISA Kit
E-EL-Ch0561	Chicken 16-α OHE-1 (16-α Hydroxyestrogen 1) ELISA Kit
E-EL-Ch0562	Chicken Pro-MMP-13 (Pro-Matrix Metalloproteinase 13) ELISA Kit
E-EL-Ch0563	Chicken CACYBP (Calcyclin Binding Protein) ELISA Kit
E-EL-Ch0564	Chicken TYR (Tyrosinase) ELISA Kit
E-EL-Ch0570	Chicken α1-AT (Alpha 1-Antitrypsin) ELISA Kit
E-EL-Ch0571	Chicken CENPB (Centromere Protein B) ELISA Kit
E-EL-Ch0572	Chicken HSP-60/HSPD1 (Heat Shock Protein 60) ELISA Kit
E-EL-Ch0573	Chicken DCN (Decorin) ELISA Kit
E-EL-Ch0574	Chicken MYH10 (Myosin Heavy Chain 10, Non Muscle) ELISA Kit
E-EL-Ch0575	Chicken MVP (Major Vault Protein) ELISA Kit
E-EL-Ch0576	Chicken MCP-2/CCL8 (Monocyte Chemotactic Protein 2) ELISA Kit
E-EL-Ch0577	Chicken BFP (Brain Finger Protein) ELISA Kit
E-EL-Ch0578	Chicken AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit
E-EL-Ch0579	Chicken MRP1 (Multidrug Resistance Associated Protein 1) ELISA Kit
E-EL-Ch0580	Chicken BMP-10 (Bone Morphogenetic Protein 10) ELISA Kit
E-EL-Ch0581	Chicken PAI2 (Plasminogen Activator Inhibitor 2) ELISA Kit
E-EL-Ch0583	Chicken GRIK2 (Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit
E-EL-Ch0584	Chicken Tyk-2 (Tyrosine Kinase 2) ELISA Kit
E-EL-Ch0585	Chicken ApoM (Apolipoprotein M) ELISA Kit

E-EL-Ch0586	Chicken UBE2A (Ubiquitin Conjugating Enzyme E2A) ELISA Kit
E-EL-Ch0587	Chicken EAR1 (Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit
E-EL-Ch0588	Chicken RPL10 (Ribosomal Protein L10) ELISA Kit
E-EL-Ch0589	Chicken Preptin ELISA Kit
E-EL-Ch0591	Chicken Bcl-10 (B-cell Leukemia/Lymphoma 10) ELISA Kit
E-EL-Ch0593	Chicken BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit
E-EL-Ch0594	Chicken RXFP1 (Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
E-EL-Ch0595	Chicken PSPN (Persephin) ELISA Kit
E-EL-Ch0596	Chicken CX40 (Connexin 40) ELISA Kit
E-EL-Ch0598	Chicken NAGase (N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
E-EL-Ch0599	Chicken CNTF (Ciliary Neurotrophic Factor) ELISA Kit
E-EL-Ch0600	Chicken APC (Activated Protein C) ELISA Kit
E-EL-Ch0602	Chicken EL (EuglobulinLysis) ELISA Kit
E-EL-Ch0603	Chicken AMBP (Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit
E-EL-Ch0604	Chicken CRYαB (Crystallin Alpha B) ELISA Kit
E-EL-Ch0605	Chicken PSAP (Prosaposin) ELISA Kit
E-EL-Ch0606	Chicken IL-5 (Interleukin 5) ELISA Kit
E-EL-Ch0607	Chicken SELE (E-Selectin) ELISA Kit
E-EL-Ch0608	Chicken PADI4 (Peptidyl Arginine Deiminase Type IV) ELISA Kit
E-EL-Ch0609	Chicken TCGF- III (Tcell Growth Factor III) ELISA Kit
E-EL-Ch0610	Chicken IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit

E-EL-Ch0611	Chicken FADS2 (Fatty Acid Desaturase 2) ELISA Kit
E-EL-Ch0613	Chicken IL-20 (Interleukin 20) ELISA Kit
E-EL-Ch0614	Chicken GLUT2 (Glucose Transporter 2) ELISA Kit
E-EL-Ch0615	Chicken PKI β (Protein Kinase Inhibitor Beta) ELISA Kit
E-EL-Ch0616	Chicken BF (Factor B) ELISA Kit
E-EL-Ch0617	Chicken SUMO2 (Small Ubiquitin Related Modifier 2) ELISA Kit
E-EL-Ch0618	Chicken DHEA-S7 (Dehydroepiandrosterone S7) ELISA Kit
E-EL-Ch0619	Chicken MS (Melatonin Sulfate) ELISA Kit
E-EL-Ch0620	Chicken IgA (Immunoglobulin A) ELISA Kit
E-EL-Ch0621	Chicken PF4 (Platelet Factor 4) ELISA Kit
E-EL-Ch0624	Chicken ARR β 1 (Arrestin Beta 1) ELISA Kit
E-EL-Ch0625	Chicken MOG (Myelin Oligodendrocyte Glycoprotein) ELISA Kit
E-EL-Ch0626	Chicken ALD (Aldolase) ELISA Kit
E-EL-Ch0627	Chicken FETUB (Fetuin B) ELISA Kit
E-EL-Ch0628	Chicken GTF2A1 (General Transcription Factor IIA, Polypeptide 1) ELISA Kit
E-EL-Ch0629	Chicken PGRP-S (Peptidoglycan Recognition Protein-S) ELISA Kit
E-EL-Ch0630	Chicken MBP/MBL (Mannose Binding Protein/Mannose Binding Lectin) ELISA Kit
E-EL-Ch0631	Chicken PDGFB (Platelet Derived Growth Factor Subunit B) ELISA Kit
E-EL-Ch0633	Chicken GnRH (Gonadotropin Releasing Hormone) ELISA Kit
E-EL-Ch0634	Chicken RIPK2 (Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit
E-EL-Ch0636	Chicken GAL2 (Galectin 2) ELISA Kit

E-EL-Ch0637	Chicken c-jun ELISA Kit
E-EL-Ch0638	Chicken MAP2K6/MKK6 (Mitogen Activated Protein Kinase Kinase 6) ELISA Kit
E-EL-Ch0639	Chicken P III CP (Procollagen III C-Terminal ProPeptide) ELISA Kit
E-EL-Ch0640	Chicken REG4 (Regenerating Islet Derived Protein 4) ELISA Kit
E-EL-Ch0641	Chicken β -EP (Beta-Endorphin) ELISA Kit
E-EL-Ch0642	Chicken C4BP α (C4 Binding Protein Alpha) ELISA Kit
E-EL-Ch0644	Chicken OPN (Osteopontin) ELISA Kit
E-EL-Ch0645	Chicken cTn-I/TNNI3 (cardiac Troponin I) ELISA Kit
E-EL-Ch0646	Chicken MPIF-1/CCL23 (Myeloid Progenitor Inhibitory Factor 1) ELISA Kit
E-EL-Ch0647	Chicken SRP9 (Signal Recognition Particle 9kDa) ELISA Kit
E-EL-Ch0648	Chicken SMAF (Specific Macrophage Arming Factor) ELISA Kit
E-EL-Ch0649	Chicken LFA-3 (Lymphocyte Function Associated Antigen 3) ELISA Kit
E-EL-Ch0650	Chicken IGFBP-1 (Insulin Like Growth Factor Binding Protein 1) ELISA Kit
E-EL-Ch0651	Chicken GSTs (Glutathione S-Transferases) ELISA Kit
E-EL-Ch0652	Chicken VEGF (Vascular Endothelial Cell Growth Factor) ELISA Kit
E-EL-Ch0653	Chicken PAF (Platelet Activating Factor) ELISA Kit
E-EL-Ch0654	Chicken Fc ϵ R I (Receptor I for the Fc Fragment of IgE) ELISA Kit
E-EL-Ch0657	Chicken GAL9 (Galectin 9) ELISA Kit
E-EL-Ch0658	Chicken ApoA4 (Apolipoprotein A4) ELISA Kit
E-EL-Ch0661	Chicken MAP2K1 (Mitogen Activated Protein Kinase Kinase 1) ELISA Kit
E-EL-Ch0662	Chicken HHIP (Hedgehog Interacting Protein) ELISA Kit

E-EL-Ch0663	Chicken LAB7-2/CD86 (B-Lymphocyte Activation Antigen B7-2) ELISA Kit
E-EL-Ch0664	Chicken PON3 (Paraoxonase 3) ELISA Kit
E-EL-Ch0665	Chicken GM-CSF Ab (Anti-Granulocyte Macrophage Colony Stimulating Factor Antibody) ELISA Kit
E-EL-Ch0666	Chicken NT-ProANP (N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit
E-EL-Ch0667	Chicken SNCα (Synuclein Alpha) ELISA Kit
E-EL-Ch0668	Chicken CD59 (Protectin) ELISA Kit
E-EL-Ch0670	Chicken MIP-1α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
E-EL-Ch0671	Chicken SAA3 (Serum Amyloid A3) ELISA Kit
E-EL-Ch0672	Chicken GTEAGE (Glomerular Tissue Extract Advanced Glycosylation End Products) ELISA Kit
E-EL-Ch0673	Chicken PRSS2 (Protease, Serine, 2) ELISA Kit
E-EL-Ch0674	Chicken FBN2 (Fibrillin 2) ELISA Kit
E-EL-Ch0675	Chicken FcγR3A (Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit
E-EL-Ch0676	Chicken UFC (Urinary Free Cortisol) ELISA Kit
E-EL-Ch0677	Chicken TNK2 (Tyrosine Kinase, Non Receptor 2) ELISA Kit
E-EL-Ch0678	Chicken Des (Desmin) ELISA Kit
E-EL-Ch0679	Chicken PCSK9 (Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit
E-EL-Ch0680	Chicken CPE (Carboxypeptidase E) ELISA Kit
E-EL-Ch0681	Chicken MCAM (Melanoma Cell Adhesion Molecule) ELISA Kit
E-EL-Ch0683	Chicken OPGL (Osteoprotegerin Ligand) ELISA Kit
E-EL-Ch0684	Chicken COL10 (Collagen Type X) ELISA Kit
E-EL-Ch0685	Chicken PGA (Pepsinogen A) ELISA Kit

E-EL-Ch0686	Chicken NF (Neurofilament) ELISA Kit
E-EL-Ch0687	Chicken F9 (Coagulation Factor IX) ELISA Kit
E-EL-Ch0688	Chicken PAEP (Progestagen-Associated Endometrial Protein) ELISA Kit
E-EL-Ch0689	Chicken SCG2 (Secretogranin II) ELISA Kit
E-EL-Ch0691	Chicken PCP (Procollagen C Peptidase)ELISA Kit
E-EL-Ch0692	Chicken BPI (Bactericidal/Permeability Increasing Protein) ELISA Kit
E-EL-Ch0694	Chicken RAR α (Retinoic Acid Receptor Alpha) ELISA Kit
E-EL-Ch0695	Chicken C3b (Complement Fragment 3b) ELISA Kit
E-EL-Ch0697	Chicken iPLA2 (Phospholipidase A2, Calcium Independent) ELISA Kit
E-EL-Ch0698	Chicken NMP-22 (Nuclear Matrix Protein 22) ELISA Kit
E-EL-Ch0700	Chicken CTGF (Connective Tissue Growth Factor) ELISA Kit
E-EL-Ch0701	Chicken THRSP (Thyroid Hormone Responsive) ELISA Kit
E-EL-Ch0702	Chicken FABP3 (Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit
E-EL-Ch0703	Chicken CNP (C-type Natriuretic Peptide) ELISA Kit
E-EL-Ch0704	Chicken HSP β 6/HSP-20 (Heat Shock Protein Beta 6) ELISA Kit
E-EL-Ch0705	Chicken TCF20 (Transcription Factor 20) ELISA Kit
E-EL-Ch0706	Chicken CMR (Costimulatory Molecules Receptor) ELISA Kit
E-EL-Ch0707	Chicken cADPRH/CD38 (Cyclic ADP Hydrolase) ELISA Kit
E-EL-Ch0708	Chicken CNX (Calnexin) ELISA Kit
E-EL-Ch0709	Chicken PRB2 (Basic Salivary Proline Rich Protein 2) ELISA Kit
E-EL-Ch0710	Chicken CCNB (Cyclin-B) ELISA Kit

E-EL-Ch0711	Chicken FBN1 (Fibrillin 1) ELISA Kit
E-EL-Ch0712	Chicken KLK6 (Kallikrein 6) ELISA Kit
E-EL-Ch0713	Chicken APAF-1 (Apoptosis Protease Activating Factor-1) ELISA Kit
E-EL-Ch0714	Chicken RB1 (Retinoblastoma Protein 1) ELISA Kit
E-EL-Ch0715	Chicken H-Cadherin (Cadherin, Heart) ELISA Kit
E-EL-Ch0716	Chicken CYP26A1 (Cytochrome P450 26A1) ELISA Kit
E-EL-Ch0717	Chicken FLOT2 (Flotillin 2) ELISA Kit
E-EL-Ch0718	Chicken ACA-IgG (Anti-Cardiolipin Antibody IgG) ELISA Kit
E-EL-Ch0719	Chicken CYP24A1 (Cytochrome P450 24A1) ELISA Kit
E-EL-Ch0720	Chicken TAB (Thyroxine Antibody) ELISA Kit
E-EL-Ch0721	Chicken F11 (Coagulation Factor XI) ELISA Kit
E-EL-Ch0724	Chicken MMP-2 (Matrix Metalloproteinase 2) ELISA Kit
E-EL-Ch0725	Chicken FBG (Fibrinogen) ELISA Kit
E-EL-Ch0727	Chicken HMGCS1 (3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit
E-EL-Ch0728	Chicken CASP7 (Caspase 7) ELISA Kit
E-EL-Ch0729	Chicken Sp17 (Sperm Protein 17) ELISA Kit
E-EL-Ch0730	Chicken CFB (Complement Factor B) ELISA Kit
E-EL-Ch0731	Chicken LDHA (Lactate Dehydrogenase A) ELISA Kit
E-EL-Ch0732	Chicken α HSP (Alpha-Hemoglobin Stabilizing Protein) ELISA Kit
E-EL-Ch0736	Chicken DAP6 (Death Associated Protein 6) ELISA Kit
E-EL-Ch0737	Chicken DSE (Dermatan Sulfate Epimerase) ELISA Kit

E-EL-Ch0738	Chicken TMSβ10 (Thymosin Beta 10) ELISA Kit
E-EL-Ch0739	Chicken DT (Dopachrome Tautomerase) ELISA Kit
E-EL-Ch0740	Chicken SOX1 (Sex Determining Region Y Box Protein 1) ELISA Kit
E-EL-Ch0741	Chicken MK (Midkine) ELISA Kit
E-EL-Ch0742	Chicken sTNF-1 (Soluble Tumor Necrosis Factor 1) ELISA Kit
E-EL-Ch0743	Chicken MAPKAPK3 (Mitogen Activated Protein Kinase Activated Protein Kinase 3) ELISA Kit
E-EL-Ch0744	Chicken IκBβ (Inhibitory Subunit of NF Kappa B Beta) ELISA Kit
E-EL-Ch0745	Chicken TFPI2 (Tissue Factor Pathway Inhibitor 2) ELISA Kit
E-EL-Ch0746	Chicken PAFAHIB3 (Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit
E-EL-Ch0748	Chicken CASP4 (Caspase 4) ELISA Kit
E-EL-Ch0749	Chicken SAG (S Antigen) ELISA Kit
E-EL-Ch0750	Chicken Lptn/XCL1 (Lymphotactin) ELISA Kit
E-EL-Ch0751	Chicken RI (Ribonuclease Inhibitor) ELISA Kit
E-EL-Ch0753	Chicken PKCθ (Protein Kinase C Theta)ELISA Kit
E-EL-Ch0754	Chicken HBμ (Hemoglobin, Mu) ELISA Kit
E-EL-Ch0755	Chicken IL-6R (Interleukin 6 Receptor) ELISA Kit
E-EL-Ch0757	Chicken IL-11 (Interleukin 11) ELISA Kit
E-EL-Ch0758	Chicken SPA (Pulmonary Surfactant Associated Protein A) ELISA Kit
E-EL-Ch0759	Chicken PLAA (Phospholipase A2 Activating Protein) ELISA Kit
E-EL-Ch0760	Chicken ASB (Arylsulfatase B) ELISA Kit
E-EL-Ch0761	Chicken Try2 (Trypsinogen II) ELISA Kit

E-EL-Ch0762	Chicken CETP (Cholesteryl Ester Transfer Protein) ELISA Kit
E-EL-Ch0763	Chicken RLR (Retinoic Acid Inducible Gene/RIG-Like Receptor) ELISA Kit
E-EL-Ch0764	Chicken RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
E-EL-Ch0765	Chicken AIG1 (Androgen Induced Protein 1) ELISA Kit
E-EL-Ch0766	Chicken LIPF (Lipase, Gastric) ELISA Kit
E-EL-Ch0767	Chicken GAL7 (Galectin 7) ELISA Kit
E-EL-Ch0768	Chicken TDGF1 (Teratocarcinoma Derived Growth Factor 1) ELISA Kit
E-EL-Ch0769	Chicken PKC ζ (Protein Kinase C Zeta) ELISA Kit
E-EL-Ch0770	Chicken ERK2 (Extracellular Signal Regulated Kinase 2) ELISA Kit
E-EL-Ch0771	Chicken NUP205 (Nucleoporin 205kDa) ELISA Kit
E-EL-Ch0772	Chicken PP (Protein Phosphatase) ELISA Kit
E-EL-Ch0773	Chicken IAA (Insulin Autoantibody) ELISA Kit
E-EL-Ch0774	Chicken TNPO2 (Transportin 2) ELISA Kit
E-EL-Ch0775	Chicken ATXN1 (Ataxin 1) ELISA Kit
E-EL-Ch0777	Chicken CA2 (Carbonic Anhydrase II) ELISA Kit
E-EL-Ch0778	Chicken PGD2S (Prostaglandin D2 Synthase, Brain) ELISA Kit
E-EL-Ch0779	Chicken α FP (Alpha-Fetoprotein) ELISA Kit
E-EL-Ch0780	Chicken MRP (Multidrug Resistance-Associated Protein) ELISA Kit
E-EL-Ch0781	Chicken AGRN (Agrin) ELISA Kit
E-EL-Ch0782	Chicken SCTR (Secretin Receptor) ELISA Kit
E-EL-Ch0783	Chicken TACE (TNF α Converting Enzyme) ELISA Kit

E-EL-Ch0785	Chicken CKMT1A (Creatine Kinase, Mitochondrial 1A) ELISA Kit
E-EL-Ch0787	Chicken CLU (Clusterin) ELISA Kit
E-EL-Ch0789	Chicken TREM-1 (Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
E-EL-Ch0790	Chicken ApoA5 (Apolipoprotein A5) ELISA Kit
E-EL-Ch0791	Chicken MYOG (Myogenin) ELISA Kit
E-EL-Ch0793	Chicken CKLFSF6 (Chemokine Like Factor Superfamily 6) ELISA Kit
E-EL-Ch0794	Chicken GIF (Gastric Intrinsic Factor) ELISA Kit
E-EL-Ch0795	Chicken pERK (Phospho Extracellular Signal Regulated Kinase) ELISA Kit
E-EL-Ch0796	Chicken NEB (Noradrenaline Bitartas) ELISA Kit
E-EL-Ch0797	Chicken BMG/ β 2-MG (Beta-2-Microglobulin) ELISA Kit
E-EL-Ch0798	Chicken Bcl-3 (B-cell Leukemia/Lymphoma 3) ELISA Kit
E-EL-Ch0799	Chicken MUC6 (Mucin 6) ELISA Kit
E-EL-Ch0800	Chicken Neo-T4 (Neonatal Thyroxine) ELISA Kit
E-EL-Ch0802	Chicken IFABP/FABP2 (Intestinal Fatty Acid Binding Protein) ELISA Kit
E-EL-Ch0804	Chicken IL-1ra (Interleukin 1 Receptor Antagonist) ELISA Kit
E-EL-Ch0805	Chicken PARC/CCL18 (Pulmonary Activation Regulated Chemokine) ELISA Kit
E-EL-Ch0806	Chicken ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit
E-EL-Ch0807	Chicken TSP1 (Thrombin Sensitive Protein 1) ELISA Kit
E-EL-Ch0808	Chicken L-LDH (L-Lactate Dehydrogenase) ELISA Kit
E-EL-Ch0809	Chicken 5-HIAA (5-Hydroxyindolacetic Acid) ELISA Kit
E-EL-Ch0811	Chicken PRELP (Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit

E-EL-Ch0812	Chicken NUMA1 (Nuclear Mitotic Apparatus Protein 1) ELISA Kit
E-EL-Ch0813	Chicken CD23/FcεRII (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit
E-EL-Ch0814	Chicken MXD1 (MAX Dimerization Protein 1) ELISA Kit
E-EL-Ch0815	Chicken UBAP2 (Ubiquitin Associated Protein 2) ELISA Kit
E-EL-Ch0816	Chicken PDHα (Pyruvate Dehydrogenase Alpha) ELISA Kit
E-EL-Ch0817	Chicken Flt3L(FMS Like Tyrosine Kinase 3 Ligand) ELISA Kit
E-EL-Ch0818	Chicken GDF6 (Growth Differentiation Factor 6) ELISA Kit
E-EL-Ch0819	Chicken TBG (Thyroxine Binding Globulin) ELISA Kit
E-EL-Ch0820	Chicken ALDOA (Aldolase A, Fructose Bisphosphate) ELISA Kit
E-EL-Ch0821	Chicken SGK2 (Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit
E-EL-Ch0822	Chicken APA (Anti-Phospholipid Antibody) ELISA Kit
E-EL-Ch0823	Chicken SOX3 (Sex Determining Region Y Box Protein 3) ELISA Kit
E-EL-Ch0824	Chicken AAA (Anti-Actin Antibody) ELISA Kit
E-EL-Ch0825	Chicken GTF2B (General Transcription Factor IIB) ELISA Kit
E-EL-Ch0827	Chicken cPLA2 (Cytosolic Phospholipase A2) ELISA Kit
E-EL-Ch0828	Chicken ALDOC (Aldolase C, Fructose Bisphosphate) ELISA Kit
E-EL-Ch0829	Chicken FGF9 (Fibroblast Growth Factor 9) ELISA Kit
E-EL-Ch0830	Chicken USP8 (Ubiquitin Specific Peptidase 8) ELISA Kit
E-EL-Ch0831	Chicken VCAM-1/CD106 (soluble Vascular Cell Adhesion Molecule 1) ELISA Kit
E-EL-Ch0832	Chicken MYH1 (Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit
E-EL-Ch0833	Chicken PC/CPG (Choline Phosphoglyceride) ELISA Kit

E-EL-Ch0834	Chicken Ntn4 (Netrin 4) ELISA Kit
E-EL-Ch0835	Chicken PT (Prothrombin) ELISA Kit
E-EL-Ch0836	Chicken FABP8 (Fatty Acid Binding Protein 8, Myelin) ELISA Kit
E-EL-Ch0837	Chicken FABP4 (Fatty Acid-binding Protein 4, Adipocyte) ELISA Kit
E-EL-Ch0838	Chicken OTF2A (Octamer Transcription Factor 2A) ELISA Kit
E-EL-Ch0839	Chicken plkB α (Phospho-Inhibitory Subunit Of NF Kappa B Alpha) ELISA Kit
E-EL-Ch0840	Chicken CX43 (Connexin 43) ELISA Kit
E-EL-Ch0841	Chicken IL-23 (Interleukin 23) ELISA Kit
E-EL-Ch0842	Chicken SYCN (Syncollin) ELISA Kit
E-EL-Ch0843	Chicken GzmB (Granzyme B) ELISA Kit
E-EL-Ch0844	Chicken LEP (Leptin) ELISA Kit
E-EL-Ch0845	Chicken C3c (Complement C3 Convertase) ELISA Kit
E-EL-Ch0846	Chicken LTD4 (Leukotriene D4) ELISA Kit
E-EL-Ch0847	Chicken FUCA (Alpha-L-Fucosidase, Tissue) ELISA Kit
E-EL-Ch0848	Chicken IL-18 (Interleukin 18) ELISA Kit
E-EL-Ch0849	Chicken SSRP1 (Structure Specific Recognition Protein 1) ELISA Kit
E-EL-Ch0850	Chicken DKK2 (Dickkopf Related Protein 2) ELISA Kit
E-EL-Ch0851	Chicken VIM (Vimentin) ELISA Kit
E-EL-Ch0852	Chicken COL1 α 1 (Collagen Type I Alpha 1) ELISA Kit
E-EL-Ch0853	Chicken 6-K-PG (6-Keto-prostaglandin) ELISA Kit
E-EL-Ch0854	Chicken IL-22 (Interleukin 22) ELISA Kit

E-EL-Ch0855	Chicken SOX9 (Sex Determining Region Y Box Protein 9) ELISA Kit
E-EL-Ch0856	Chicken LZM (Lysozyme) ELISA Kit
E-EL-Ch0857	Chicken DKK4 (Dickkopf Related Protein 4) ELISA Kit
E-EL-Ch0859	Chicken TLR9 (Toll-Like Receptor 9) ELISA Kit
E-EL-Ch0860	Chicken TPO (Thrombopoietin) ELISA Kit
E-EL-Ch0862	Chicken SAP (Serum Amyloid P Component) ELISA Kit
E-EL-Ch0863	Chicken CAP1 (Adenylate Cyclase Associated Protein 1) ELISA Kit
E-EL-Ch0864	Chicken IGFBP-3 (Insulin Like Growth Factor Binding Protein 3) ELISA Kit
E-EL-Ch0865	Chicken CX3CL1 (Chemokine C-X3-C-Motif Ligand) ELISA Kit
E-EL-Ch0866	Chicken sTfR1 (Soluble Transferrin Receptor1) ELISA Kit
E-EL-Ch0867	Chicken DBP (DNA Binding Protein) ELISA Kit
E-EL-Ch0868	Chicken CASP 9 (Caspase 9) ELISA Kit
E-EL-Ch0869	Chicken Ext-N (N-Terminal Extein) ELISA Kit
E-EL-Ch0870	Chicken KLK8 (Kallikrein 8) ELISA Kit
E-EL-Ch0871	Chicken PRSS33 (Protease, Serine 33) ELISA Kit
E-EL-Ch0873	Chicken ICAM-4 (Intercellular Adhesion Molecule-4) ELISA Kit
E-EL-Ch0874	Chicken HLE (Leukocyte Exastase)ELISA Kit
E-EL-Ch0875	Chicken PL (Pancreatic Lipase) ELISA Kit
E-EL-Ch0876	Chicken CYP21B (Cytochrome P450 21B) ELISA Kit
E-EL-Ch0877	Chicken dDAVP (1-Desamino-8D-Arginine Vasopressin) ELISA Kit
E-EL-Ch0878	Chicken FcαR/CD89 (Receptor I for the Fc Fragment of IgA) ELISA Kit

E-EL-Ch0879	Chicken GAL1 (Galectin 1) ELISA Kit
E-EL-Ch0880	Chicken NT-3 (Neurotrophin-3) ELISA Kit
E-EL-Ch0882	Chicken VEGFR1/Flt1 (Vascuoar Endothelial Growth Factor Receptor 1) ELISA Kit
E-EL-Ch0883	Chicken β -Glu (Beta Glucosidase) ELISA Kit
E-EL-Ch0884	Chicken AFP-L3 (Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit
E-EL-Ch0885	Chicken TE (Telomerase) ELISA Kit
E-EL-Ch0886	Chicken AMY2 (Amylase Alpha 2, Pancreatic) ELISA Kit
E-EL-Ch0887	Chicken MSR1 (Macrophage Scavenger Receptor 1) ELISA Kit
E-EL-Ch0888	Chicken PAR4 (Protease Activated Receptor 4) ELISA Kit
E-EL-Ch0889	Chicken PKD2 (Protein Kinase D2)NELISA Kit
E-EL-Ch0890	Chicken PHB (Prohibitin) ELISA Kit
E-EL-Ch0891	Chicken S100A7 (S100 Calcium Binding Protein A7) ELISA Kit
E-EL-Ch0893	Chicken ACO1 (Aconitase 1) ELISA Kit
E-EL-Ch0894	Chicken BALP (Bone Alkaline Phosphatase) ELISA Kit
E-EL-Ch0895	Chicken FLRT2 (Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit
E-EL-Ch0896	Chicken Anti-MPO (Anti-Myeloperoxidase Antibody) ELISA Kit
E-EL-Ch0897	Chicken ApoF (Apolipoprotein F) ELISA Kit
E-EL-Ch0898	Chicken TPO (Thyroid Peroxidase) ELISA Kit
E-EL-Ch0899	Chicken Fc γ R II/CD32 (Receptor II for the Fc Fragment of IgG) ELISA Kit
E-EL-Ch0900	Chicken TAP (Trypsinogen Activation Peptide) ELISA Kit
E-EL-Ch0901	Chicken sFMC (Soluble Fibrin Monomer Complex) ELISA Kit

E-EL-Ch0902	Chicken SOCS3 (Suppressors of Cytokine Signaling 3) ELISA Kit
E-EL-Ch0903	Chicken EphA2 (Tyrosine protein kinase receptor A2) ELISA Kit
E-EL-Ch0904	Chicken F7 (Coagulation Factor VII) ELISA Kit
E-EL-Ch0905	Chicken MX1 (Myxovirus Resistance 1) ELISA Kit
E-EL-Ch0906	Chicken ApoE (Apolipoprotein E) ELISA Kit
E-EL-Ch0908	Chicken ES-Ab (Anti-Endostatin Antibody) ELISA Kit
E-EL-Ch0909	Chicken CDSN (Corneodesmosin) ELISA Kit
E-EL-Ch0910	Chicken ITIH4 (Inter Alpha-Globulin Inhibitor H4) ELISA Kit
E-EL-Ch0911	Chicken GLUT1 (Glucose Transporter 1) ELISA Kit
E-EL-Ch0912	Chicken ATF-2 (Activating Transcription Factor 2) ELISA Kit
E-EL-Ch0913	Chicken DNase-I (Deoxyribonuclease I) ELISA Kit
E-EL-Ch0915	Chicken FRS2 (Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit
E-EL-Ch0916	Chicken MUC17 (Mucin 17) ELISA Kit
E-EL-Ch0917	Chicken EGR2 (Early Growth Response Protein 2) ELISA Kit
E-EL-Ch0918	Chicken Slit2 (Slit Homolog 2) ELISA Kit
E-EL-Ch0919	Chicken TMSβ4 (Thymosin Beta 4) ELISA Kit
E-EL-Ch0920	Chicken P27 (P27 Protein) ELISA Kit
E-EL-Ch0921	Chicken GPX1 (Glutathione Peroxidase 1) ELISA Kit
E-EL-Ch0922	Chicken OSMR (Oncostatin M Receptor) ELISA Kit
E-EL-Ch0923	Chicken DHEA (Dehydroepiandrosterone) ELISA Kit
E-EL-Ch0924	Chicken PDIA3 (Protein Disulfide Isomerase A3) ELISA Kit

E-EL-Ch0925	Chicken SCFR/c-Kit (Stem Cell Factor Receptor) ELISA Kit
E-EL-Ch0926	Chicken CSF3/GCSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit
E-EL-Ch0927	Chicken BNP (Brain Natriuretic Peptide) ELISA Kit
E-EL-Ch0928	Chicken UCRP (Ubiquitin Cross Reactive Protein) ELISA Kit
E-EL-Ch0929	Chicken GR α (Glucocorticoid Receptor Alpha) ELISA Kit
E-EL-Ch0931	Chicken C8b (Complement Component 8b) ELISA Kit
E-EL-Ch0932	Chicken PHSA (Poly Serum Antigen) ELISA Kit
E-EL-Ch0933	Chicken AR (Androgen Receptor) ELISA Kit
E-EL-Ch0934	Chicken RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit
E-EL-Ch0935	Chicken TSA (Thyroid Stimulating Antibody) ELISA Kit
E-EL-Ch0936	Chicken TPP (Thrombus Precursor Protein) ELISA Kit
E-EL-Ch0937	Chicken CASP12 (Caspase 12) ELISA Kit
E-EL-Ch0938	Chicken OBTF1 (Octamer Binding Transcription Factor 1) ELISA Kit
E-EL-Ch0939	Chicken CCSA-2 (Colon Cancer Specific Antigen-2) ELISA Kit
E-EL-Ch0940	Chicken TSPAN30 (Tetraspanin 30) ELISA Kit
E-EL-Ch0943	Chicken ANKRD1 (Ankyrin Repeat Domain 1) ELISA Kit
E-EL-Ch0944	Chicken CTR (Calcitonin Receptor) ELISA Kit
E-EL-Ch0945	Chicken ETRA (Endothelin Receptor A) ELISA Kit
E-EL-Ch0948	Chicken MlgA (Membrane IgA) ELISA Kit
E-EL-Ch0949	Chicken CYP1B1 (Cytochrome P450 1B1) ELISA Kit
E-EL-Ch0950	Chicken ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit

E-EL-Ch0951	Chicken RNASEH (Ribonuclease H)ELISA Kit
E-EL-Ch0952	Chicken DNM1 (Dynamin 1) ELISA Kit
E-EL-Ch0953	Chicken PP (Pancreatic Polypeptide) ELISA Kit
E-EL-Ch0954	Chicken VLA (Very Late Appearing Antigen) ELISA Kit
E-EL-Ch0956	Chicken PPAR- δ (Peroxisome Proliferator Activated Receptor Delta) ELISA Kit
E-EL-Ch0957	Chicken PTN (Pleiotrophin) ELISA Kit
E-EL-Ch0958	Chicken β -EPR (Beta-Endorphin Receptor) ELISA Kit
E-EL-Ch0959	Chicken RBP4 (Retinol Binding Protein 4, Plasma) ELISA Kit
E-EL-Ch0960	Chicken EFNA5 (Ephrin A5) ELISA Kit
E-EL-Ch0961	Chicken CTSL (Cathepsin L) ELISA Kit
E-EL-Ch0965	Chicken CRH (Corticotropin Releasing Hormone) ELISA Kit
E-EL-Ch0966	Chicken DNM2 (Dynamin 2) ELISA Kit
E-EL-Ch0968	Chicken BAD (Bcl2 Associated Death Promoter) ELISA Kit
E-EL-Ch0969	Chicken GLUT3 (Glucose Transporter 3) ELISA Kit
E-EL-Ch0970	Chicken MCP/CD46 (Membrane Cofactor Protein) ELISA Kit
E-EL-Ch0971	Chicken HNRPA2B1 (Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit
E-EL-Ch0973	Chicken Anti-CIIAb (Anti-Type II Collagen Antibody) ELISA Kit
E-EL-Ch0974	Chicken CD30 (Cluster of Differentiation 30) ELISA Kit
E-EL-Ch0975	Chicken SPARC (Secreted Protein Acidic And Rich In Cysteine) ELISA Kit
E-EL-Ch0976	Chicken ATGA/TGAB (Anti-Thyroid-Globulin Antibody) ELISA Kit
E-EL-Ch0977	Chicken P-GP (Permeability Glycoprotein) ELISA Kit

E-EL-Ch0978	Chicken NFRkB (Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit
E-EL-Ch0979	Chicken S100B (S100 Calcium Binding Protein B) ELISA Kit
E-EL-Ch0980	Chicken ANG II R2-Ab (Angiotensin II Receptor 2 Antibody) ELISA Kit
E-EL-Ch0981	Chicken HLTF (Helicase Like Transcription Factor) ELISA Kit
E-EL-Ch0982	Chicken ATA (Angiotensinase A) ELISA Kit
E-EL-Ch0983	Chicken PGM (Phosphoglucomutase) ELISA Kit
E-EL-Ch0984	Chicken EloA (Elongin A) ELISA Kit
E-EL-Ch0985	Chicken PA-IgM (Platelet Associated Antibody) ELISA Kit
E-EL-Ch0986	Chicken PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit
E-EL-Ch0987	Chicken CKLFSF5 (Chemokine Like Factor Superfamily 5) ELISA Kit
E-EL-Ch0988	Chicken CD109 (Cluster of Differentiation 109) ELISA Kit
E-EL-Ch0989	Chicken GP39 (Glycoprotein 39, Cartilage) ELISA Kit
E-EL-Ch0990	Chicken TSP-2 (Thrombospondin-2) ELISA Kit
E-EL-Ch0991	Chicken AZU1 (Azurocidin 1) ELISA Kit
E-EL-Ch0992	Chicken NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit
E-EL-Ch0993	Chicken PGC (Pepsinogen C) ELISA Kit
E-EL-Ch0994	Chicken SOX2 (Sex Determining Region Y Box Protein 2) ELISA Kit
E-EL-Ch0995	Chicken GSTθ1/GSTt1 (Glutathione S Transferase Theta 1) ELISA Kit
E-EL-Ch0996	Chicken CFI (Complement Factor I) ELISA Kit
E-EL-Ch0997	Chicken Gal-6S (Galactose-6-Sulphatase) ELISA Kit
E-EL-Ch0998	Chicken NPHN (Nephrin) ELISA Kit

E-EL-Ch1001	Chicken C4a (Complement Component 4a) ELISA Kit
E-EL-Ch1003	Chicken SSFA2 (Sperm Specific Antigen 2) ELISA Kit
E-EL-Ch1004	Chicken EGR1 (Early Growth Response Protein 1) ELISA Kit
E-EL-Ch1005	Chicken F5 (Coagulation Factor V) ELISA Kit
E-EL-Ch1006	Chicken PLC (Phospholipase C) ELISA Kit
E-EL-Ch1008	Chicken CLDN11 (Claudin 11) ELISA Kit
E-EL-Ch1010	Chicken TR (Targeted Ribonuclease) ELISA Kit
E-EL-Ch1012	Chicken IκBα (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit
E-EL-Ch1013	Chicken HTRA1 (HtrA Serine, Peptidase 1) ELISA Kit
E-EL-Ch1014	Chicken SMM (Smooth Muscle Myosin) ELISA Kit
E-EL-Ch1016	Chicken SOD3 (Superoxide Dismutase 3, Extracellular) ELISA Kit
E-EL-Ch1017	Chicken ApoC3 (Apolipoprotein C3) ELISA Kit
E-EL-Ch1018	Chicken TJP1 (Tight Junction Protein 1) ELISA Kit
E-EL-Ch1019	Chicken RAG-2 (Recombination Activating Gene 2) ELISA Kit
E-EL-Ch1020	Chicken LHRH/GnRH (Luteinizing Hormone Releasing Hormone) ELISA Kit
E-EL-Ch1021	Chicken RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit
E-EL-Ch1022	Chicken CSTA (Cystatin A) ELISA Kit
E-EL-Ch1023	Chicken Bax (Bcl-2 Associated X Protein) ELISA Kit
E-EL-Ch1024	Chicken EGFR (Epidermal Growth Factor Receptor) ELISA Kit
E-EL-Ch1025	Chicken PAR3 (Protease Activated Receptor 3) ELISA Kit
E-EL-Ch1026	Chicken LpPLA2 (Phospholipase A2, Lipoprotein Associated) ELISA Kit

E-EL-Ch1027	Chicken β 2-GP1 Ab IgG (β 2-Glycoprotein 1 Antibody IgG) ELISA Kit
E-EL-Ch1028	Chicken CHAc (Choline acetylase) ELISA Kit
E-EL-Ch1030	Chicken CTXI (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
E-EL-Ch1031	Chicken NUP62 (Nucleoporin 62kDa) ELISA Kit
E-EL-Ch1032	Chicken ADH1 (Alcohol Dehydrogenase 1) ELISA Kit
E-EL-Ch1034	Chicken Tn-I (Troponin I) ELISA Kit
E-EL-Ch1035	Chicken PC (Plasma Anticoagulant Protein C) ELISA Kit
E-EL-Ch1036	Chicken TSI (Thyroid Stimulating Immunoglobulin) ELISA Kit
E-EL-Ch1037	Chicken RPS6K β 1 (Ribosomal Protein S6 Kinase Beta 1) ELISA Kit
E-EL-Ch1038	Chicken CRN (Corin) ELISA Kit
E-EL-Ch1039	Chicken MMP-3 (Matrix Metalloproteinase 3) ELISA Kit
E-EL-Ch1040	Chicken EFNA3 (Ephrin A3) ELISA Kit
E-EL-Ch1041	Chicken BLNK (B-Cell Linker Protein) ELISA Kit
E-EL-Ch1042	Chicken PCI (Protein C Inhibitor) ELISA Kit
E-EL-Ch1043	Chicken NPY (Neuropeptide Y) ELISA Kit
E-EL-Ch1044	Chicken GPC4 (Glypican 4) ELISA Kit
E-EL-Ch1045	Chicken UBD (Ubiquitin D) ELISA Kit
E-EL-Ch1046	Chicken MMP-13 (Matrix Metalloproteinase 13) ELISA Kit
E-EL-Ch1047	Chicken GAL4 (Galectin 4) ELISA Kit
E-EL-Ch1048	Chicken PLA2R1 (Phospholipase A2 Receptor 1) ELISA Kit
E-EL-Ch1050	Chicken GRO β (Growth Regulated Oncogene Beta) ELISA Kit

E-EL-Ch1051	Chicken CNTN4 (Contactin 4) ELISA Kit
E-EL-Ch1052	Chicken ACTH (Adrenocorticotrophic Hormone) ELISA Kit
E-EL-Ch1053	Chicken PLA2G2A (Phospholipase A2, Group II A) ELISA Kit
E-EL-Ch1054	Chicken S100A12 (S100 Calcium Binding Protein A12) ELISA Kit
E-EL-Ch1056	Chicken GR β (Glucocorticoid Receptor Beta) ELISA Kit
E-EL-Ch1057	Chicken Smad3/MADH3 (Mothers Against Decapentaplegic Homolog 3) ELISA Kit
E-EL-Ch1058	Chicken NE/ELA2 (Elastase 2, Neutrophil) ELISA Kit
E-EL-Ch1060	Chicken 8-epi-PGF2 α (8-Epi-Prostaglandin F2 Alpha) ELISA Kit
E-EL-Ch1061	Chicken CYPD (Cyclophilin D) ELISA Kit
E-EL-Ch1062	Chicken FRA (Fibronectin Related Antigen) ELISA Kit
E-EL-Ch1063	Chicken PLA2G2D (Phospholipase A2, Group II D) ELISA Kit
E-EL-Ch1064	Chicken SPB (Pulmonary Surfactant Associated Protein B) ELISA Kit
E-EL-Ch1065	Chicken TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
E-EL-Ch1066	Chicken DPYSL2 (Dihydropyrimidinase Like Protein 2) ELISA Kit
E-EL-Ch1067	Chicken sCKR (Soluble Cytokine Receptor) ELISA Kit
E-EL-Ch1068	Chicken CASP1 (Caspase 1) ELISA Kit
E-EL-Ch1069	Chicken ACP5 (Tartrate Resistant Acid Phosphatase 5) ELISA Kit
E-EL-Ch1070	Chicken MTHFR (Methylenetetrahydrofolate Reductase) ELISA Kit
E-EL-Ch1071	Chicken CNTN3 (Contactin 3) ELISA Kit
E-EL-Ch1072	Chicken MUC2 (Mucin 2) ELISA Kit
E-EL-Ch1074	Chicken Ras-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit

E-EL-Ch1076	Chicken COL15 (Collagen Type XV) ELISA Kit
E-EL-Ch1077	Chicken C3a (Complement Component 3a) ELISA Kit
E-EL-Ch1078	Chicken TGF- α (Transforming Growth Factor α) ELISA Kit
E-EL-Ch1079	Chicken PIICP (Procollagen II C-Terminal ProPeptide) ELISA Kit
E-EL-Ch1080	Chicken DPP6 (Dipeptidyl Peptidase VI) ELISA Kit
E-EL-Ch1081	Chicken GATA5 (GATA Binding Protein 5) ELISA Kit
E-EL-Ch1082	Chicken MAG Ab (Anti-Myelin Associated Glycoprotein Antibody) ELISA Kit
E-EL-Ch1083	Chicken GYPE (Glycophorin E) ELISA Kit
E-EL-Ch1085	Chicken Gp49a (Glycoprotein 49A) ELISA Kit
E-EL-Ch1086	Chicken PSG1/SP1 (Pregnancy Specific Beta-1-Glycoprotein 1) ELISA Kit
E-EL-Ch1087	Chicken MYL1 (Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit
E-EL-Ch1089	Chicken CA1 (Carbonic Anhydrase I) ELISA Kit
E-EL-Ch1090	Chicken C3c (Complement Fragment 3c) ELISA Kit
E-EL-Ch1091	Chicken bTC (Betacellulin) ELISA Kit
E-EL-Ch1093	Chicken SLC/CCL21 (Secondary Lymphoid Tissue Chemokine) ELISA Kit
E-EL-Ch1095	Chicken MEPE (Matrix Extracellular Phosphoglycoprotein) ELISA Kit
E-EL-Ch1096	Chicken BDNF (Brain Derived Neurotrophic Factor) ELISA Kit
E-EL-Ch1097	Chicken Pro-MMP-10 (Pro-Matrix Metalloproteinase 10) ELISA Kit
E-EL-Ch1098	Chicken ELA2B (Elastase 2B) ELISA Kit
E-EL-Ch1099	Chicken PIK3C α (Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit
E-EL-Ch1100	Chicken FOXO2 (Forkhead Box Protein O2) ELISA Kit

E-EL-Ch1101	Chicken MANF (Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit
E-EL-Ch1103	Chicken PLP (Parathyroid Hormone-Like Protein) ELISA Kit
E-EL-Ch1104	Chicken PTGIS (Prostaglandin I Synthase) ELISA Kit
E-EL-Ch1105	Chicken S100A9 (S100 Calcium Binding Protein A9) ELISA Kit
E-EL-Ch1106	Chicken DAF (Decay Accelerating Factor) ELISA Kit
E-EL-Ch1107	Chicken NRG-2 (Neuregulin 2) ELISA Kit
E-EL-Ch1109	Chicken ARSA (Arylsulfatase A) ELISA Kit
E-EL-Ch1110	Chicken INSR (Insulin Receptor) ELISA Kit
E-EL-Ch1111	Chicken FMOD (Fibromodulin) ELISA Kit
E-EL-Ch1112	Chicken P2 (Myelin Protein 2) ELISA Kit
E-EL-Ch1113	Chicken Sp100 (Sp100 Nuclear Antigen) ELISA Kit
E-EL-Ch1114	Chicken Surv (Survivin) ELISA Kit
E-EL-Ch1115	Chicken OC/BGP (Osteocalcin/Bone Gla Protein) ELISA Kit
E-EL-Ch1116	Chicken A-CoA (Acetyl Coenzyme A) ELISA Kit
E-EL-Ch1117	Chicken COL13 (Collagen Type XIII) ELISA Kit
E-EL-Ch1118	Chicken GPRC5A (G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit
E-EL-Ch1119	Chicken FPR2 (Formyl Peptide Receptor 2) ELISA Kit
E-EL-Ch1120	Chicken DSG-E1 (Desmogleins 1) ELISA Kit
E-EL-Ch1121	Chicken NADPH (Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit
E-EL-Ch1122	Chicken CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) ELISA Kit
E-EL-Ch1124	Chicken CVL/EZR (Cytovillin/Ezrin) ELISA Kit

E-EL-Ch1125	Chicken REG3α (Regenerating Islet Derived Protein 3 Alpha) ELISA Kit
E-EL-Ch1126	Chicken NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) ELISA Kit
E-EL-Ch1127	Chicken GDF9 (Growth Differentiation Factor 9) ELISA Kit
E-EL-Ch1128	Chicken GAP43 (Growth Associated Protein 43) ELISA Kit
E-EL-Ch1131	Chicken GAD (Glutamic Acid Decarboxylase) ELISA Kit
E-EL-Ch1133	Chicken SDF-1/CXCL12 (Stromal Cell Derived Factor 1) ELISA Kit
E-EL-Ch1134	Chicken β-CTx (Beta Crosslaps) ELISA Kit
E-EL-Ch1136	Chicken LATS (Long Acting Thyroid Stimulator) ELISA Kit
E-EL-Ch1137	Chicken BLC (B-Lymphocyte Chemoattractant) ELISA Kit
E-EL-Ch1138	Chicken AP 13 (Apelin 13) ELISA Kit
E-EL-Ch1139	Chicken PLGA (Plasminogen Activator) ELISA Kit
E-EL-Ch1140	Chicken ASGR1 (Asialoglycoprotein Receptor 1) ELISA Kit
E-EL-Ch1141	Chicken SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit
E-EL-Ch1142	Chicken ENG/CD105 (Endoglin) ELISA Kit
E-EL-Ch1144	Chicken TNX (Tenascin X) ELISA Kit
E-EL-Ch1145	Chicken CXCR3 (Chemokine CXC Receptor 3) ELISA Kit
E-EL-Ch1146	Chicken F12 (Coagulation Factor XII) ELISA Kit
E-EL-Ch1147	Chicken GSANA (Granulocyte Specific Anti-Nuclear Antibody) ELISA Kit
E-EL-Ch1148	Chicken IL-12 (Interleukin 12) ELISA Kit
E-EL-Ch1149	Chicken FSTL1 (Follistatin Like Protein 1) ELISA Kit
E-EL-Ch1150	Chicken COL4α1 (Collagen Type IV Alpha 1) ELISA Kit

E-EL-Ch1151	Chicken AGR2 (Anterior Gradient Protein 2) ELISA Kit
E-EL-Ch1152	Chicken CKLF (Chemokine Like Factor) ELISA Kit
E-EL-Ch1153	Chicken D5R/DRD5 (Dopamine Receptor D5) ELISA Kit
E-EL-Ch1154	Chicken PDCN (Podocin) ELISA Kit
E-EL-Ch1155	Chicken AhR (Aryl Hydrocarbon Receptor) ELISA Kit
E-EL-Ch1156	Chicken MYD88 (Myeloid Differentiation Factor 88) ELISA Kit
E-EL-Ch1157	Chicken α 1ACT (Alpha1 Antichymotrypsin) ELISA Kit
E-EL-Ch1158	Chicken HSF2 (Heat Shock Transcription Factor 2) ELISA Kit
E-EL-Ch1159	Chicken c-myc (c-myc Oncogene Product) ELISA Kit
E-EL-Ch1160	Chicken PYGB (Glycogen Phosphorylase, Brain) ELISA Kit
E-EL-Ch1161	Chicken ATF-5 (Activating Transcription Factor 5) ELISA Kit
E-EL-Ch1162	Chicken M-AChR M1 (Muscarinic Acetylcholine Receptor M1) ELISA Kit
E-EL-Ch1163	Chicken L1CAM/CD171 (L1-Cell Adhesion Molecule) ELISA Kit
E-EL-Ch1164	Chicken RND1 (Rho Family GTPase 1) ELISA Kit
E-EL-Ch1165	Chicken MUC3 (Mucin 3) ELISA Kit
E-EL-Ch1166	Chicken PP14 (Placental Protein14) ELISA Kit
E-EL-Ch1168	Chicken KAF (Keratinocyte Autocrine Factor) ELISA Kit
E-EL-Ch1169	Chicken TRF (Transferrin) ELISA Kit
E-EL-Ch1170	Chicken CCND1 (Cyclin-D1) ELISA Kit
E-EL-Ch1171	Chicken LOXL4 (Lysyl Oxidase Like Protein 4) ELISA Kit
E-EL-Ch1172	Chicken D3R/DRD3 (Dopamine Receptor D3) ELISA Kit

E-EL-Ch1173	Chicken GLUT4 (Glucose Transporter 4) ELISA Kit
E-EL-Ch1174	Chicken ApoC2 (Apolipoprotein C2) ELISA Kit
E-EL-Ch1175	Chicken Mly/CXCL9 (Monocyte Interferon Gamma Inducing Factor) ELISA Kit
E-EL-Ch1177	Chicken Collagenase II ELISA Kit
E-EL-Ch1178	Chicken CLDN9 (Claudin 9) ELISA Kit
E-EL-Ch1179	Chicken VEGFR3/Flt4 (Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit
E-EL-Ch1180	Chicken SPARCL1 (SPARC Like Protein 1) ELISA Kit
E-EL-Ch1181	Chicken VCL (Vinculin) ELISA Kit
E-EL-Ch1182	Chicken SDF2L1 (Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit
E-EL-Ch1183	Chicken SDF-4 (Stromal Cell Derived Factor 4) ELISA Kit
E-EL-Ch1184	Chicken Thyronine ELISA Kit
E-EL-Ch1185	Chicken NAIP (Neuronal Apoptosis Inhibitory Protein) ELISA Kit
E-EL-Ch1186	Chicken PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit
E-EL-Ch1187	Chicken PTPRJ (Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit
E-EL-Ch1188	Chicken HO2 (Heme Oxygenase 2) ELISA Kit
E-EL-Ch1189	Chicken PMN Elastase (Polymorphonuclear Elastase) ELISA Kit
E-EL-Ch1190	Chicken GDF7 (Growth Differentiation Factor 7) ELISA Kit
E-EL-Ch1191	Chicken CRT (Calreticulin) ELISA Kit
E-EL-Ch1192	Chicken MFI2 (Melanotransferrin) ELISA Kit
E-EL-Ch1193	Chicken DUB (Deubiquitinating Enzyme) ELISA Kit
E-EL-Ch1194	Chicken CDNF (Cerebral Dopamine Neurotrophic Factor) ELISA Kit

E-EL-Ch1195	Chicken FcεR II (Receptor II for the Fc Fragment of IgE) ELISA Kit
E-EL-Ch1196	Chicken βAPP (β Amyloid Precursor Protein) ELISA Kit
E-EL-Ch1197	Chicken DEFα1 (Defensin Alpha 1, Neutrophil) ELISA Kit
E-EL-Ch1198	Chicken OGP (Osteogenic Growth Peptide) ELISA Kit
E-EL-Ch1200	Chicken CDX2 (Caudal Type Homeobox Transcription Factor 2) ELISA Kit
E-EL-Ch1201	Chicken PKP1 (Plakophilin 1) ELISA Kit
E-EL-Ch1202	Chicken SPOP (Speckle Type POZ Protein) ELISA Kit
E-EL-Ch1205	Chicken ACA-IgM (Anti-Cardiolipin Antibody IgM) ELISA Kit
E-EL-Ch1206	Chicken SAA (Serum Amyloid A) ELISA Kit
E-EL-Ch1207	Chicken PFKM (Phosphofructokinase, Muscle) ELISA Kit
E-EL-Ch1209	Chicken GDF3 (Growth Differentiation Factor 3) ELISA Kit
E-EL-Ch1210	Chicken ORM2 (Orosomucoid 2) ELISA Kit
E-EL-Ch1211	Chicken HSPG (Heparan Sulfate Proteoglycan) ELISA Kit
E-EL-Ch1213	Chicken ET (Epitestosterone) ELISA Kit
E-EL-Ch1214	Chicken AGER (Advanced Glycosylation End Product Specific Receptor) ELISA Kit
E-EL-Ch1215	Chicken P-CK (Pan-Cytokeratin) ELISA Kit
E-EL-Ch1216	Chicken LILRB4 (Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit
E-EL-Ch1217	Chicken CNR-1 (Cadherin Related Neuronal Receptor1) ELISA Kit
E-EL-Ch1218	Chicken FLNα (Filamin A, Alpha) ELISA Kit
E-EL-Ch1219	Chicken PIK3Cβ (Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit
E-EL-Ch1221	Chicken ANG II (Angiotensin II) ELISA Kit

E-EL-Ch1222	Chicken DKK1 (Dickkopf Related Protein 1) ELISA Kit
E-EL-Ch1223	Chicken CFD (Complement Factor D) ELISA Kit
E-EL-Ch1224	Chicken GPI (Glucose 6 Phosphate Isomerase) ELISA Kit
E-EL-Ch1226	Chicken TNF- β (Tumor Necrosis Factor Beta) ELISA Kit
E-EL-Ch1227	Chicken PCDH1 (Protocadherin 1) ELISA Kit
E-EL-Ch1228	Chicken NEFL (Neurofilament, Light Polypeptide) ELISA Kit
E-EL-Ch1229	Chicken MUC7 (Mucin 7) ELISA Kit
E-EL-Ch1230	Chicken ANG-R-Tie2 (Angiopoietin Receptor Tie2) ELISA Kit
E-EL-Ch1231	Chicken SLPI (Secretory Leukocyte Protease Inhibitor) ELISA Kit
E-EL-Ch1232	Chicken 12-HETE (12-Hydroxyeicosatetraenoic Acid) ELISA Kit
E-EL-Ch1233	Chicken CHD5 (Chromodomain Helicase DNA Binding Protein 5) ELISA Kit
E-EL-Ch1234	Chicken IL-8 (Interleukin 8) ELISA Kit
E-EL-Ch1235	Chicken C/EBP α (CCAAT/Enhancer Binding Protein Alpha) ELISA Kit
E-EL-Ch1236	Chicken C5 (Complement Component 5) ELISA Kit
E-EL-Ch1237	Chicken VWF (Von Willebrand Factor) ELISA Kit
E-EL-Ch1238	Chicken IRAK3 (Interleukin 1 Receptor Associated Kinase 3) ELISA Kit
E-EL-Ch1239	Chicken COX-2R (Cyclooxygenase-2 Receptor) ELISA Kit
E-EL-Ch1241	Chicken α -EP (α -Endomorphin) ELISA Kit
E-EL-Ch1242	Chicken CEBP γ (CCAAT/Enhancer Binding Protein Gamma) ELISA Kit
E-EL-Ch1243	Chicken VEGF-C (Vascular Endothelial Growth Factor C) ELISA Kit
E-EL-Ch1244	Chicken OTF2B (Octamer Transcription Factor 2B) ELISA Kit

E-EL-Ch1245	Chicken CR2/CD21 (Complement Receptor 2) ELISA Kit
E-EL-Ch1246	Chicken FTL (Ferritin, Light Polypeptide) ELISA Kit
E-EL-Ch1247	Chicken MASP1 (Mannan Associated Serine Protease 1) ELISA Kit
E-EL-Ch1248	Chicken CDK8 (Cyclin Dependent Kinase 8) ELISA Kit
E-EL-Ch1249	Chicken C3bR (Complement Fragment 3b Receptor) ELISA Kit
E-EL-Ch1250	Chicken UBC (Urinary Bladder Cancer Antigen) ELISA Kit
E-EL-Ch1252	Chicken MAGE (Melanoma-associated Antigen) ELISA Kit
E-EL-Ch1253	Chicken PFKP (Phosphofructokinase, Platelet) ELISA Kit
E-EL-Ch1255	Chicken AP (Allopregnanolone) ELISA Kit
E-EL-Ch1256	Chicken SST (Somatostatin) ELISA Kit
E-EL-Ch1257	Chicken FGγ (Fibrinogen Gamma) ELISA Kit
E-EL-Ch1258	Chicken TDP43 (Tar DNA Binding Protein 43KDa)ELISA Kit
E-EL-Ch1259	Chicken FSHR (Follicle Stimulating Hormone Receptor) ELISA Kit
E-EL-Ch1260	Chicken NOX1 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit
E-EL-Ch1261	Chicken HELIX- II (Type II Collagen Helical Peptide) ELISA Kit
E-EL-Ch1262	Chicken MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit
E-EL-Ch1263	Chicken EAR12 (Eosinophil Associated Ribonuclease A family Member 12) ELISA Kit
E-EL-Ch1264	Chicken HNRPU (Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit
E-EL-Ch1265	Chicken OTR (Oxytocin Receptor) ELISA Kit
E-EL-Ch1266	Chicken Apo A1-Ab (Anti-Apolipoprotein A1 antibody) ELISA Kit
E-EL-Ch1267	Chicken ADRA1A (Adrenergic Receptor alpha-1A) ELISA Kit

E-EL-Ch1268	Chicken ApoA2 (Apolipoprotein A2) ELISA Kit
E-EL-Ch1269	Chicken PRDX5 (Peroxisome oxidoreductase 5) ELISA Kit
E-EL-Ch1270	Chicken CTNNβ1 (Catenin, Beta 1) ELISA Kit
E-EL-Ch1271	Chicken PGD2-MOX (Prostaglandin D2-Methoxime) ELISA Kit
E-EL-Ch1272	Chicken GlyCAM1 (Glycosylation Dependent Cell Adhesion Molecule 1) ELISA Kit
E-EL-Ch1273	Chicken CENPA (Centromere Protein A) ELISA Kit
E-EL-Ch1274	Chicken GLP-2 (Glucagon Like Peptide 2) ELISA Kit
E-EL-Ch1275	Chicken TBP (TATA Binding Protein) ELISA Kit
E-EL-Ch1276	Chicken TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
E-EL-Ch1278	Chicken CTX-II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit
E-EL-Ch1279	Chicken TBP2 (Thioredoxin Binding Protein 2) ELISA Kit
E-EL-Ch1280	Chicken PAPPA (Pregnancy Associated Plasma Protein A/ Pappalysin-1) ELISA Kit
E-EL-Ch1281	Chicken CR1 (Complement Receptor type 1) ELISA Kit
E-EL-Ch1282	Chicken AVPR1A (Arginine Vasopressin Receptor 1A) ELISA Kit
E-EL-Ch1283	Chicken PDGFR (Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit
E-EL-Ch1284	Chicken PGDS (Prostaglandin D Synthase) ELISA Kit
E-EL-Ch1285	Chicken PAR1 (Protease Activated Receptor 1) ELISA Kit
E-EL-Ch1286	Chicken MUC13 (Mucin 13) ELISA Kit
E-EL-Ch1287	Chicken DAG1 (Dystrophin Associated Glycoprotein 1) ELISA Kit
E-EL-Ch1288	Chicken HCCS (Holocytochrome C Synthase) ELISA Kit
E-EL-Ch1290	Chicken CA3 (Carbonic Anhydrase III, Muscle Specific) ELISA Kit

E-EL-Ch1291	Chicken MCP-3 (Monocyte Chemotactic Protein 3) ELISA Kit
E-EL-Ch1292	Chicken CTAGE (Cutaneous T-Cell Lymphocyte Associated Antigen) ELISA Kit
E-EL-Ch1295	Chicken CFP (Complement Factor P) ELISA Kit
E-EL-Ch1297	Chicken CDK6 (Cyclin Dependent Kinase 6) ELISA Kit
E-EL-Ch1298	Chicken IgA1 (Immunoglobulin A1) ELISA Kit
E-EL-Ch1299	Chicken ANG (Angiostatin) ELISA Kit
E-EL-Ch1300	Chicken PLA2G12B (Phospholipase A2, Group XII B) ELISA Kit
E-EL-Ch1301	Chicken DPYD (Dihydropyrimidine Dehydrogenase) ELISA Kit
E-EL-Ch1302	Chicken EED (Embryonic Ectoderm Development) ELISA Kit
E-EL-Ch1303	Chicken REN (Renin) ELISA Kit
E-EL-Ch1304	Chicken FSCN (Fascin) ELISA Kit
E-EL-Ch1305	Chicken ApoB100 (Apolipoprotein B100) ELISA Kit
E-EL-Ch1306	Chicken MYH9 (Myosin Heavy Chain 9, Non Muscle) ELISA Kit
E-EL-Ch1307	Chicken HA (Hemagglutinin) ELISA Kit
E-EL-Ch1308	Chicken PCK1 (Phosphoenolpyruvate Carboxykinase 1) ELISA Kit
E-EL-Ch1309	Chicken CENPF (Centromere Protein F) ELISA Kit
E-EL-Ch1310	Chicken ANG-R-Tie1 (Angiopoietin Receptor Tie1) ELISA Kit
E-EL-Ch1312	Chicken FASN (Fatty Acid Synthase) ELISA Kit
E-EL-Ch1313	Chicken NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit
E-EL-Ch1314	Chicken α 2-PI (α 2-Plasmin Inhibitor) ELISA Kit
E-EL-Ch1315	Chicken CK-18/KRT18 (Cytokeratin 18) ELISA Kit

E-EL-Ch1316	Chicken PROS (Protein S) ELISA Kit
E-EL-Ch1317	Chicken PKM2 (Pyruvate Kinase, Muscle) ELISA Kit
E-EL-Ch1318	Chicken PTHrP (Parathyroid Hormone Related Protein) ELISA Kit
E-EL-Ch1319	Chicken PKCε (Protein Kinase C Epsilon) ELISA Kit
E-EL-Ch1320	Chicken SA (Sialic Acid) ELISA Kit
E-EL-Ch1321	Chicken PAL (L-Phenylalanine Ammonia-Lyase) ELISA Kit
E-EL-Ch1322	Chicken GC (Glycocalicin) ELISA Kit
E-EL-Ch1323	Chicken NTXI (Cross Linked N-telopeptide of Type I Collagen) ELISA Kit
E-EL-Ch1324	Chicken FN (Fibronectin) ELISA Kit
E-EL-Ch1325	Chicken HSP-27/HSPB1 (Heat Shock Protein 27) ELISA Kit
E-EL-Ch1326	Chicken PDK4 (Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit
E-EL-Ch1327	Chicken PROC (Protein C) ELISA Kit
E-EL-Ch1328	Chicken NRIP1 (Nuclear Receptor Interacting Protein 1) ELISA Kit
E-EL-Ch1329	Chicken VEGF-D (Vascular Endothelial Growth Factor D) ELISA Kit
E-EL-Ch1330	Chicken PTGES2 (Prostaglandin E Synthase 2) ELISA Kit
E-EL-Ch1331	Chicken UGT2B4 (Uridine Diphosphate Glucuronosyltransferase 2 Family Polypeptide B4) ELISA Kit
E-EL-Ch1333	Chicken MIP-3β/ELC/CCL19 (Macrophage Inflammatory Protein 3β) ELISA Kit
E-EL-Ch1334	Chicken EPP (Erythrocyte Protoporphyrin) ELISA Kit
E-EL-Ch1335	Chicken PRI (Placental Ribonuclease Inhibitor) ELISA Kit
E-EL-Ch1336	Chicken GH2 (Growth Hormone 2) ELISA Kit
E-EL-Ch1337	Chicken VTN (Vitronectin) ELISA Kit

E-EL-Ch1338	Chicken GATA4 (GATA Binding Protein 4) ELISA Kit
E-EL-Ch1339	Chicken FGFBP1 (Fibroblast Growth Factor Binding Protein 1) ELISA Kit
E-EL-Ch1340	Chicken IgG2 (Immunoglobulin G2) ELISA Kit
E-EL-Ch1341	Chicken CPP (Copeptin) ELISA Kit
E-EL-Ch1342	Chicken IGFBP-5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit
E-EL-Ch1343	Chicken COL3α1 (Collagen Type III Alpha 1) ELISA Kit
E-EL-Ch1344	Chicken CASP14 (Caspase 14) ELISA Kit
E-EL-Ch1345	Chicken SHC/SLP-76 (SH2-Containing Protein) ELISA Kit
E-EL-Ch1346	Chicken PLN (Phospholamban) ELISA Kit
E-EL-Ch1347	Chicken SEPCR/PROCR (Soluble Endothelial Protein C Receptor) ELISA Kit
E-EL-Ch1348	Chicken GCK (Glucokinase) ELISA Kit
E-EL-Ch1349	Chicken PGF (Prostaglandin F) ELISA Kit
E-EL-Ch1350	Chicken PDGFRβ (Platelet Derived Growth Factor Receptor Beta) ELISA Kit
E-EL-Ch1352	Chicken EREG (Epiregulin) ELISA Kit
E-EL-Ch1353	Chicken TALLA-1 (T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit
E-EL-Ch1354	Chicken FGF-7/KGF (Fibroblast Growth Factor 7) ELISA Kit
E-EL-Ch1355	Chicken ECF/CCL11 (Eosinophil Chemotactic Factor) ELISA Kit
E-EL-Ch1357	Chicken SIGLEC10 (Sialic Acid Binding Ig Like Lectin 10) ELISA Kit
E-EL-Ch1358	Chicken FABP7 (Fatty Acid Binding Protein 7, Brain) ELISA Kit
E-EL-Ch1359	Chicken AGPA (Anti-Glucoprotein 210 Antibody) ELISA Kit
E-EL-Ch1360	Chicken PGAM2 (Phosphoglycerate Mutase 2) ELISA Kit

E-EL-Ch1361	Chicken PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
E-EL-Ch1362	Chicken CSF2R α (Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit
E-EL-Ch1363	Chicken TNNI2 (Troponin I Type 2, Fast Skeletal) ELISA Kit
E-EL-Ch1364	Chicken TAGLN (Transgelin) ELISA Kit
E-EL-Ch1365	Chicken FSH (Follicle-Stimulating Hormone) ELISA Kit
E-EL-Ch1366	Chicken PRB1 (Basic Salivary Proline Rich Protein 1) ELISA Kit
E-EL-Ch1367	Chicken TACR2 (Tachykinin Receptor 2) ELISA Kit
E-EL-Ch1368	Chicken α Actin (Alpha-Actin) ELISA Kit
E-EL-Ch1369	Chicken PIAS1 (Protein Inhibitor of Activated STAT 1) ELISA Kit
E-EL-Ch1370	Chicken TGM3 (Transglutaminase 3, Epidermal) ELISA Kit
E-EL-Ch1371	Chicken CRY1 (Cryptochrome 1) ELISA Kit
E-EL-Ch1373	Chicken CNTN1 (Contactin 1) ELISA Kit
E-EL-Ch1374	Chicken TERT (Telomerase Reverse Transcriptase) ELISA Kit
E-EL-Ch1375	Chicken CD3d (Cluster of Differentiation 3d) ELISA Kit
E-EL-Ch1377	Chicken UreC (Urease Related Protein C) ELISA Kit
E-EL-Ch1379	Chicken CALLA (Common Acute Lymphocytic Leukaemia Antigen) ELISA Kit
E-EL-Ch1380	Chicken STAT4 (Signal Transducer and Activator of Transcription 4) ELISA Kit
E-EL-Ch1381	Chicken IL-2R α (Interleukin 2 Receptor Alpha) ELISA Kit
E-EL-Ch1382	Chicken SERT (Serotonin Transporter) ELISA Kit
E-EL-Ch1384	Chicken MIP-5 (Macrophage Inflammatory Protein 5) ELISA Kit
E-EL-Ch1385	Chicken CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit

E-EL-Ch1387	Chicken MTL (Motilin) ELISA Kit
E-EL-Ch1388	Chicken MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit
E-EL-Ch1389	Chicken SM/SPH (Sphingomyelin) ELISA Kit
E-EL-Ch1391	Chicken COL17 (Collagen Type XVII) ELISA Kit
E-EL-Ch1392	Chicken PRSS23 (Protease, Serine 23) ELISA Kit
E-EL-Ch1393	Chicken VDR (Vitamin D Receptor) ELISA Kit
E-EL-Ch1394	Chicken MMSAM (Melanoma Metastasis Surface Adhesion Molecule) ELISA Kit
E-EL-Ch1395	Chicken ERCC1 (Excision Repair Cross Complementing Rodent Repair Deficiency Complementation 1) ELISA Kit
E-EL-Ch1396	Chicken FBLN4 (Fibulin 4) ELISA Kit
E-EL-Ch1397	Chicken C/EBP δ (CCAAT/Enhancer Binding Protein, Delta) ELISA Kit
E-EL-Ch1398	Chicken ADCY2 (Adenylate Cyclase 2) ELISA Kit
E-EL-Ch1399	Chicken TNNC2 (Troponin C Type 2) ELISA Kit
E-EL-Ch1400	Chicken FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit
E-EL-Ch1401	Chicken NUP160 (Nucleoporin 160kDa) ELISA Kit
E-EL-Ch1402	Chicken ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit
E-EL-Ch1403	Chicken CYSLTR2 (Cysteinyl Leukotriene Receptor 2) ELISA Kit
E-EL-Ch1404	Chicken Ly6a (Lymphocyte Antigen 6 Complex Locus A) ELISA Kit
E-EL-Ch1405	Chicken CENPE (Centromere Protein E) ELISA Kit
E-EL-Ch1407	Chicken MYBPC3 (Myosin Binding Protein C, Cardiac) ELISA Kit
E-EL-Ch1408	Chicken β -actin (Actin Beta) ELISA Kit
E-EL-Ch1409	Chicken CEA (Carcinoembryonic Antigen) ELISA Kit

E-EL-Ch1410	Chicken CYP1A1 (Cytochrome P450 1A1) ELISA Kit
E-EL-Ch1411	Chicken tPA (Plasminogen Activator, Tissue) ELISA Kit
E-EL-Ch1412	Chicken FCN1 (Ficolin 1) ELISA Kit
E-EL-Ch1413	Chicken IL-15 (Interleukin 15) ELISA Kit
E-EL-Ch1414	Chicken LEPR (Leptin Receptor) ELISA Kit
E-EL-Ch1415	Chicken CAPA (Cytotoxin Associated Protein A) ELISA Kit
E-EL-Ch1417	Chicken PYGM (Glycogen Phosphorylase, Muscle) ELISA Kit
E-EL-Ch1418	Chicken AIF1 (Allograft Inflammatory Factor 1) ELISA Kit
E-EL-Ch1419	Chicken GDF10 (Growth Differentiation Factor 10) ELISA Kit
E-EL-Ch1420	Chicken PDAP1 (Pdgfa Associated Protein 1) ELISA Kit
E-EL-Ch1421	Chicken Salusin α (Salusin Alpha) ELISA Kit
E-EL-Ch1422	Chicken ABCG2 (ATP Binding Cassette Transporter G2) ELISA Kit
E-EL-Ch1423	Chicken F-Box Protein 2(FBXO2) ELISA Kit
E-EL-Ch1424	Chicken CCSA-4 (Colon Cancer Specific Antigen-4) ELISA Kit
E-EL-Ch1425	Chicken BLVRB (Biliverdin Reductase B) ELISA Kit
E-EL-Ch1426	Chicken LSP (Liver Specific Lipoprotein) ELISA Kit
E-EL-Ch1427	Chicken TUB β (β -Tubulin) ELISA Kit
E-EL-Ch1428	Chicken EAR3 (Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit
E-EL-Ch1429	Chicken PTF1 α (Pancreas Specific Transcription Factor 1 α) ELISA Kit
E-EL-Ch1430	Chicken SF1 (Steroidogenic Factor 1) ELISA Kit
E-EL-Ch1431	Chicken SDC3 (Syndecan 3) ELISA Kit

E-EL-Ch1432	Chicken DPD (Deoxypyridinoline) ELISA Kit
E-EL-Ch1434	Chicken AQP-4 (Aquaporin 4) ELISA Kit
E-EL-Ch1435	Chicken PAI-1 (Plasminogen Activator Inhibitor-1) ELISA Kit
E-EL-Ch1436	Chicken Tn-T (Troponin T) ELISA Kit
E-EL-Ch1437	Chicken PLG (Plasminogen) ELISA Kit
E-EL-Ch1438	Chicken DβH (Dopamine-β-Hydroxylase) ELISA Kit
E-EL-Ch1439	Chicken β2-GP (Beta2 Glycoprotein) ELISA Kit
E-EL-Ch1440	Chicken FE (Ferritin) ELISA Kit
E-EL-Ch1442	Chicken PLP (Proteolipid Protein) ELISA Kit
E-EL-Ch1443	Chicken OxLDL (Oxidized Low Density Lipoprotein) ELISA Kit
E-EL-Ch1444	Chicken CA724 (Tumor Marker CA724) ELISA Kit
E-EL-Ch1446	Chicken MYO (Myoglobin) ELISA Kit
E-EL-Ch1448	Chicken PF3 (Platelet Factor 3) ELISA Kit
E-EL-Ch1449	Chicken TP53 (Tumor Protein p53) ELISA Kit
E-EL-Ch1450	Chicken C1INH (Complement 1 Inhibitor) ELISA Kit
E-EL-Ch1451	Chicken CLEC4E (C-Type Lectin Domain Family 4 Member E) ELISA Kit
E-EL-Ch1452	Chicken CNTN2 (Contactin 2) ELISA Kit
E-EL-Ch1453	Chicken CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit
E-EL-Ch1454	Chicken ALB (Albumin) ELISA Kit
E-EL-Ch1455	Chicken OPG (Osteoprotegerin) ELISA Kit
E-EL-Ch1456	Chicken TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit

E-EL-Ch1457	Chicken GSK3 β (Glycogen Synthase Kinase 3 Beta) ELISA Kit
E-EL-Ch1459	Chicken Gas6 (Growth Arrest Specific Protein 6) ELISA Kit
E-EL-Ch1460	Chicken COL7 (Collagen Type VII) ELISA Kit
E-EL-Ch1461	Chicken IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit
E-EL-Ch1462	Chicken FG α (Fibrinogen Alpha) ELISA Kit
E-EL-Ch1463	Chicken TECK/CCL25 (Thymus Expressed Chemokine) ELISA Kit
E-EL-Ch1464	Chicken PKB (Protein Kinase B) ELISA Kit
E-EL-Ch1466	Chicken UCHL1 (Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit
E-EL-Ch1467	Chicken LIFR (Leukemia Inhibitory Factor Receptor) ELISA Kit
E-EL-Ch1469	Chicken ANGPTL4 (Angiopoietin Like Protein 4) ELISA Kit
E-EL-Ch1470	Chicken VMA (Vanillylmandelic Acid) ELISA Kit
E-EL-Ch1471	Chicken KIM-1 (Kidney Injury Molecule 1) ELISA Kit
E-EL-Ch1472	Chicken C4BP (C4 Binding Protein) ELISA Kit
E-EL-Ch1473	Chicken MUC4 (Mucin 4) ELISA Kit
E-EL-Ch1474	Chicken TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit
E-EL-Ch1475	Chicken CREB (Cyclic AMP Response Element Binding Protein) ELISA Kit
E-EL-Ch1476	Chicken D2R/DRD2 (Dopamine Receptor D2) ELISA Kit
E-EL-Ch1477	Chicken ANG II R1-Ab (Angiotensin II Receptor 1 Antibody) ELISA Kit
E-EL-Ch1478	Chicken ApoC1 (Apolipoprotein C1) ELISA Kit
E-EL-Ch1479	Chicken DR-70TM (Tumor Marker DR-70 For Lung Cancer) ELISA Kit
E-EL-Ch1480	Chicken AIM2 (Absent In Melanoma 2) ELISA Kit

E-EL-Ch1481	Chicken CNR1 (Cannabinoid Receptor 1, Brain) ELISA Kit
E-EL-Ch1482	Chicken PDP (Pyruvate Dehydrogenase Phosphatase) ELISA Kit
E-EL-Ch1483	Chicken PIIINP (N-Terminal Procollagen III Propeptide) ELISA Kit
E-EL-Ch1484	Chicken APP (Amyloid Precursor Protein) ELISA Kit
E-EL-Ch1486	Chicken F1+2 (Prothrombin Fragment 1+2) ELISA Kit
E-EL-Ch1487	Chicken ApoH (Apolipoprotein H) ELISA Kit
E-EL-Ch1488	Chicken MYCBP (C-Myc Binding Protein) ELISA Kit
E-EL-Ch1490	Chicken NRG-4 (Neuregulin 4) ELISA Kit
E-EL-Ch1491	Chicken PAC3 (Platelet Associated Complement 3) ELISA Kit
E-EL-Ch1492	Chicken UGCG (UDP-glucose ceramide glucosyltransferase) ELISA Kit
E-EL-Ch1493	Chicken EFNA4 (Ephrin A4) ELISA Kit
E-EL-Ch1494	Chicken fFDFT1 (Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit
E-EL-Ch1495	Chicken ERβ (Estrogen Receptor Beta) ELISA Kit
E-EL-Ch1496	Chicken SIAE (Sialic Acid Acetyltransferase) ELISA Kit
E-EL-Ch1497	Chicken CFTR (Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit
E-EL-Ch1498	Chicken MTF (Microtransferrin) ELISA Kit
E-EL-Ch1499	Chicken CLDN14 (Claudin 14) ELISA Kit
E-EL-Ch1500	Chicken BAFF/CD257/TNFSF13B (B-Cell Activating Factor) ELISA Kit
E-EL-Ch1501	Chicken SCT (Secretin) ELISA Kit
E-EL-Ch1502	Chicken PEX2 (Peroxisomal Biogenesis Factor 2) ELISA Kit
E-EL-Ch1504	Chicken TNNI1 (Troponin I Type 1, Slow Skeletal) ELISA Kit

E-EL-Ch1505	Chicken IL-3 (Interleukin 3) ELISA Kit
E-EL-Ch1506	Chicken PYGL (Glycogen Phosphorylase, Liver) ELISA Kit
E-EL-Ch1507	Chicken PLA2G10 (Phospholipase A2, Group X) ELISA Kit
E-EL-Ch1508	Chicken OCT (Ornithine Carbamoyl Transferase) ELISA Kit
E-EL-Ch1509	Chicken CX26 (Connexin 26) ELISA Kit
E-EL-Ch1510	Chicken MT2 (Metallothionein 2) ELISA Kit
E-EL-Ch1511	Chicken NKA (Neurokinin A) ELISA Kit
E-EL-Ch1512	Chicken CYTIP (Cytohesin Interacting Protein) ELISA Kit
E-EL-Ch1513	Chicken PDE4D (Phosphodiesterase 4D, cAMP Specific) ELISA Kit
E-EL-Ch1514	Chicken 5-NT (5-Nucleotidase) ELISA Kit
E-EL-Ch1515	Chicken ACTg1 (Actin Gamma 1) ELISA Kit
E-EL-Ch1516	Chicken CDC23 (Cell Division Cycle Protein 23) ELISA Kit
E-EL-Ch1517	Chicken GSTk1 (Glutathione S Transferase Kappa 1) ELISA Kit
E-EL-Ch1518	Chicken ANG II R1/AGTR1 (Angiotensin II Receptor 1) ELISA Kit
E-EL-Ch1520	Chicken CHD3 (Chromodomain Helicase DNA Binding Protein 3) ELISA Kit
E-EL-Ch1521	Chicken LBP (Lipopolysaccharide Binding Protein) ELISA Kit
E-EL-Ch1523	Chicken SIGLEC8 (Sialic Acid Binding Ig Like Lectin 8) ELISA Kit
E-EL-Ch1524	Chicken HGF (Hepatocyte Growth Factor) ELISA Kit
E-EL-Ch1525	Chicken LPL (Lipoprotein Lipase) ELISA Kit
E-EL-Ch1526	Chicken TP5 (Thymopentin) ELISA Kit
E-EL-Ch1527	Chicken sCD14 (Soluble Cluster of Differentiation14) ELISA Kit

E-EL-Ch1532	Chicken EPO (Erythropoietin) ELISA Kit
E-EL-Ch1533	Chicken MRC1 (Mannose Receptor C Type 1) ELISA Kit
E-EL-Ch1535	Chicken IGF2BP3 (Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit
E-EL-Ch1536	Chicken LAT (Linker For Activation of T-cell) ELISA Kit
E-EL-Ch1538	Chicken AKAP1 (A Kinase Anchor Protein 1) ELISA Kit
E-EL-Ch1539	Chicken TSGF (Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
E-EL-Ch1540	Chicken GKN2 (Gastrokine 2) ELISA Kit
E-EL-Ch1542	Chicken PDH E1 (Pyruvate Dehydrogenase E1) ELISA Kit
E-EL-Ch1543	Chicken Anti-AlaRS (Anti-Alanyl-tRNA Synthetase/Anti-PL12-Antibody) ELISA Kit
E-EL-Ch1544	Chicken DHEA-S (Dehydroepiandrosterone Sulfate) ELISA Kit
E-EL-Ch1545	Chicken FSA (Fetal Sulfoslycoprotein Antigen) ELISA Kit
E-EL-Ch1546	Chicken TH (Tyrosine Hydroxylase) ELISA Kit
E-EL-Ch1547	Chicken SPTAN1 (Alpha-Fodrin) ELISA Kit
E-EL-Ch1549	Chicken GIP (Glucose Dependent Insulin Releasing Polypeptide) ELISA Kit
E-EL-Ch1550	Chicken FDP (Fibrinogen Degradation Product) ELISA Kit
E-EL-Ch1552	Chicken CCKAR (Cholecystokinin A Receptor) ELISA Kit
E-EL-Ch1553	Chicken OSF (Osteoclast Stimulating Factor 1) ELISA Kit
E-EL-Ch1554	Chicken MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit
E-EL-Ch1557	Chicken INSL3 (Insulin Like Protein 3) ELISA Kit
E-EL-Ch1558	Chicken S100A6 (S100 Calcium Binding Protein A6) ELISA Kit
E-EL-Ch1559	Chicken GDNF (Glial Cell Line Derived Neurotrophic Factor) ELISA Kit

E-EL-Ch1560	Chicken FP (Prostaglandin F Receptor) ELISA Kit
E-EL-Ch1561	Chicken rT3 (Reverse Triiodothyronine) ELISA Kit
E-EL-Ch1562	Chicken IGFBP-6 (Insulin Like Growth Factor Binding Protein 6) ELISA Kit
E-EL-Ch1563	Chicken GKN1 (Gastrokeine 1) ELISA Kit
E-EL-Ch1564	Chicken TPO-Ab (Anti-Thyroid-Peroxidase Antibody) ELISA Kit
E-EL-Ch1565	Chicken CYB561 (Cytochrome b-561) ELISA Kit
E-EL-Ch1566	Chicken DPYSL3 (Dihydropyrimidinase Like Protein 3) ELISA Kit
E-EL-Ch1567	Chicken ETFβ (Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit
E-EL-Ch1568	Chicken BMP-1 (Bone Morphogenetic Protein 1) ELISA Kit
E-EL-Ch1569	Chicken LH (Luteinizing Hormone) ELISA Kit
E-EL-Ch1570	Chicken LIPH (Lipase H) ELISA Kit
E-EL-Ch1571	Chicken MMP-10 (Matrix Metalloproteinase 10) ELISA Kit
E-EL-Ch1572	Chicken PDGF-AB (Platelet Derived Growth Factor AB) ELISA Kit
E-EL-Ch1573	Chicken RNASET2 (Ribonuclease T2) ELISA Kit
E-EL-Ch1574	Chicken FOXO3 (Forkhead Box Protein O3) ELISA Kit
E-EL-Ch1576	Chicken TLR2 (Toll-Like Receptor 2) ELISA Kit
E-EL-Ch1577	Chicken OXB (Orexin B) ELISA Kit
E-EL-Ch1578	Chicken α Manase (α Mannosidase) ELISA Kit
E-EL-Ch1579	Chicken TTN (Titin) ELISA Kit
E-EL-Ch1580	Chicken CEP110 (Centrosomal Protein 110kDa) ELISA Kit
E-EL-Ch1581	Chicken GDF2 (Growth Differentiation Factor 2) ELISA Kit

E-EL-Ch1582	Chicken NCAD (Neural Cadherin) ELISA Kit
E-EL-Ch1583	Chicken DLK-1 (Delta Like 1 Homolog) ELISA Kit
E-EL-Ch1585	Chicken SCF (Stem Cell Factor) ELISA Kit
E-EL-Ch1586	Chicken TR (Thrombin Receptor) ELISA Kit
E-EL-Ch1587	Chicken CALU (Calumenin) ELISA Kit
E-EL-Ch1588	Chicken DPP4 (Dipeptidyl Peptidase IV) ELISA Kit
E-EL-Ch1589	Chicken SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit
E-EL-Ch1590	Chicken α 1BG (Alpha-1-B-Glycoprotein) ELISA Kit
E-EL-Ch1591	Chicken UPP2 (Uridine Phosphorylase 2) ELISA Kit
E-EL-Ch1592	Chicken FETUA (Fetuin A) ELISA Kit
E-EL-Ch1593	Chicken PLA1 (Phospholipase A1) ELISA Kit
E-EL-Ch1594	Chicken TAP (Transporter Associated with Antigen Processing) ELISA Kit
E-EL-Ch1595	Chicken NFYC (Nuclear Transcription Factor Y Gamma) ELISA Kit
E-EL-Ch1596	Chicken AT (AntiThrombin) ELISA Kit
E-EL-Ch1597	Chicken VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit
E-EL-Ch1598	Chicken DLD (Dihydrolipoyl Dehydrogenase) ELISA Kit
E-EL-Ch1599	Chicken SPD (Pulmonary Surfactant Associated Protein D) ELISA Kit
E-EL-Ch1600	Chicken NNE (Non-Neuronal Enolase) ELISA Kit
E-EL-Ch1601	Chicken EM2 (Endomorphin 2) ELISA Kit
E-EL-Ch1602	Chicken TXLNb (Taxilin Beta) ELISA Kit
E-EL-Ch1603	Chicken MAOA (Type A Monoamine Oxidase) ELISA Kit

E-EL-Ch1604	Chicken ABCG1 (ATP Binding Cassette Transporter G1) ELISA Kit
E-EL-Ch1605	Chicken GPI (Glycophosphatidylinositol) ELISA Kit
E-EL-Ch1606	Chicken TNKS1 (Tankyrase 1) ELISA Kit
E-EL-Ch1607	Chicken CTAPIII (Connective Tissue Activating Peptide III) ELISA Kit
E-EL-Ch1608	Chicken LTC4 (Leukotriene C4) ELISA Kit
E-EL-Ch1609	Chicken Cx (Connexin) ELISA Kit
E-EL-Ch1610	Chicken CASP11 (Caspase 11) ELISA Kit
E-EL-Ch1611	Chicken CD4 (Cluster of Differentiation 4) ELISA Kit
E-EL-Ch1612	Chicken FHA (Filamentous Hemagglutinin) ELISA Kit
E-EL-Ch1613	Chicken PGR (Progesterone Receptor) ELISA Kit
E-EL-Ch1614	Chicken CCL4L1 (Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit
E-EL-Ch1615	Chicken Methylase ELISA Kit
E-EL-Ch1616	Chicken SFRP2 (Secreted Frizzled Related Protein 2) ELISA Kit
E-EL-Ch1617	Chicken GNLY (Granulysin) ELISA Kit
E-EL-Ch1618	Chicken PLSCR1 (Phospholipid Scramblase 1) ELISA Kit
E-EL-Ch1620	Chicken Grb7 (Growth Factor Receptor Bound Protein 7) ELISA Kit
E-EL-Ch1621	Chicken IL-10 (Interleukin-10) ELISA Kit
E-EL-Ch1622	Chicken PCI (Procollagen I) ELISA Kit
E-EL-Ch1623	Chicken VEGF-B (Vascular Endothelial Cell Growth Factor B) ELISA Kit
E-EL-Ch1624	Chicken RPL26 (Ribosomal Protein L26) ELISA Kit
E-EL-Ch1626	Chicken PCDHyA2 (Protocadherin Gamma A2) ELISA Kit

E-EL-Ch1627	Chicken FOXO1 (Forkhead Box Protein O1) ELISA Kit
E-EL-Ch1628	Chicken CK-17/KRT17 (Cytokeratin 17) ELISA Kit
E-EL-Ch1630	Chicken CK-LMW (Cytokeratin Low Molecular Weight) ELISA Kit
E-EL-Ch1631	Chicken GARS (Glycyl tRNA Synthetase) ELISA Kit
E-EL-Ch1632	Chicken Grb2 (Growth Factor Receptor Bound Protein 2) ELISA Kit
E-EL-Ch1633	Chicken FAIM3 (Fas Apoptotic Inhibitory Molecule 3) ELISA Kit
E-EL-Ch1634	Chicken TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit
E-EL-Ch1635	Chicken BCAN (Brevican) ELISA Kit
E-EL-Ch1636	Chicken SPTβN4 (Spectrin Beta, Non Erythrocytic 4) ELISA Kit
E-EL-Ch1637	Chicken MGP (Matrix Gla Protein) ELISA Kit
E-EL-Ch1639	Chicken PLAU/uPA (Urokinase-Type Plasminogen Activator) ELISA Kit
E-EL-Ch1640	Chicken SPJ (Surfactant Associated Protein J) ELISA Kit
E-EL-Ch1641	Chicken AAA (Anti-Albumin Antibody) ELISA Kit
E-EL-Ch1642	Chicken CYPB (Cyclophilin B) ELISA Kit
E-EL-Ch1643	Chicken COL1α2 (Collagen Type I Alpha 2) ELISA Kit
E-EL-Ch1644	Chicken ET-1 (Endothelin 1) ELISA Kit
E-EL-Ch1645	Chicken COL4 (Collagen Type IV) ELISA Kit
E-EL-Ch1646	Chicken NAP3 (Neutrophil Activating Protein 3) ELISA Kit
E-EL-Ch1647	Chicken LIPC (Lipase, Hepatic) ELISA Kit
E-EL-Ch1648	Chicken Grb10 (Growth Factor Receptor Bound Protein 10) ELISA Kit
E-EL-Ch1649	Chicken MN (Metanephrine) ELISA Kit

E-EL-Ch1650	Chicken DCLK1 (Doublecortin Like Kinase 1) ELISA Kit
E-EL-Ch1651	Chicken IL-12 p40 (Interleukin 12 p40) ELISA Kit
E-EL-Ch1652	Chicken VPI (Vasoactive Peptidase Inhibitors) ELISA Kit
E-EL-Ch1653	Chicken Trx (Thioredoxin) ELISA Kit
E-EL-Ch1655	Chicken AGT (Angiotensinogen) ELISA Kit
E-EL-Ch1656	Chicken PROZ (Protein Z) ELISA Kit
E-EL-Ch1657	Chicken PPP1R1A (Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit
E-EL-Ch1658	Chicken PIAS3 (Protein Inhibitor of Activated STAT3) ELISA Kit
E-EL-Ch1659	Chicken PDIA6 (Protein Disulfide Isomerase A6) ELISA Kit
E-EL-Ch1660	Chicken ABCA1(ATP-binding cassette transporter A1) ELISA Kit
E-EL-Ch1661	Chicken PLCβ1/PI-PLC (Phospholipase C Beta 1, Phosphoinositide Specific) ELISA Kit
E-EL-Ch1662	Chicken CDK1 (Cyclin Dependent Kinase 1) ELISA Kit
E-EL-Ch1663	Chicken ApoC4 (Apolipoprotein C4) ELISA Kit
E-EL-Ch1664	Chicken FST (Follistatin) ELISA Kit
E-EL-Ch1667	Chicken TSH (Thyroid Stimulating Hormone) ELISA Kit
E-EL-Ch1669	Chicken HLA-A (Leukocyte Antigen A) ELISA Kit
E-EL-Ch1670	Chicken EJ/GlyRS (Anti-EJ-Antibody) ELISA Kit
E-EL-Ch1673	Chicken EPHA10 (Ephrin Type A Receptor 10) ELISA Kit
E-EL-Ch1674	Chicken PKCδ1 (Protein Kinase C Delta) ELISA Kit
E-EL-Ch1675	Chicken ADAMTS4 (ADAM with Thrombospondin Type 1 Motif4) ELISA Kit
E-EL-Ch1676	Chicken EGR3 (Early Growth Response Protein 3) ELISA Kit

E-EL-Ch1677	Chicken ACHAP (Acetylcholinesterase Associated Protein) ELISA Kit
E-EL-Ch1678	Chicken AChE (Acetylcholinesterase) ELISA Kit
E-EL-Ch1679	Chicken STAT1 (Signal Transducer and Activator of Transcription 1) ELISA Kit
E-EL-Ch1680	Chicken APLN (Apelin) ELISA Kit
E-EL-Ch1681	Chicken STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit
E-EL-Ch1683	Chicken OXA (Orexin A) ELISA Kit
E-EL-Ch1684	Chicken NUP133 (Nucleoporin 133kDa) ELISA Kit
E-EL-Ch1685	Chicken ALP (Alkaline Phosphatase) ELISA Kit
E-EL-Ch1686	Chicken PRL (Prolactin) ELISA Kit
E-EL-Ch1687	Chicken ANXA2 (Annexin A2) ELISA Kit
E-EL-Ch1688	Chicken TL1A (Tumor Necrosis Factor Related Ligand 1A) ELISA Kit
E-EL-Ch1689	Chicken CS (Chondroitin Sulfate) ELISA Kit
E-EL-Ch1690	Chicken GLRX3 (Glutaredoxin 3) ELISA Kit
E-EL-Ch1691	Chicken TIEG1 (TGF Beta Inducible Early Response Gene 1) ELISA Kit
E-EL-Ch1692	Chicken AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) ELISA Kit
E-EL-Ch1693	Chicken OCLN (Occludin) ELISA Kit
E-EL-Ch1694	Chicken ATF-4 (Activating Transcription Factor 4) ELISA Kit
E-EL-Ch1695	Chicken DMD (Dystrophin) ELISA Kit
E-EL-Ch1696	Chicken SIGLEC2 (Sialic Acid Binding Ig Like Lectin 2) ELISA Kit
E-EL-Ch1698	Chicken C1qR1 (Complement Component 1q Receptor) ELISA Kit
E-EL-Ch1699	Chicken S100A13 (S100 Calcium Binding Protein A13) ELISA Kit

E-EL-Ch1700	Chicken FOXP3 (Forkhead Box Protein P3) ELISA Kit
E-EL-Ch1701	Chicken AR (Amphiregulin) ELISA Kit
E-EL-Ch1703	Chicken LCA/CD45(Leukocyte Common Antigen) ELISA Kit
E-EL-Ch1704	Chicken INSL5 (Insulin Like Protein 5) ELISA Kit
E-EL-Ch1705	Chicken G6PD (Glucose 6 Phosphate Dehydrogenase) ELISA Kit
E-EL-Ch1706	Chicken HSP-70/HSPA9 (Heat Shock Protein 70) ELISA Kit
E-EL-Ch1708	Chicken NPYR5 (Neuropeptide Y Receptor Y5) ELISA Kit
E-EL-Ch1710	Chicken CNN1 (Calponin 1, Basic, Smooth Muscle) ELISA Kit
E-EL-Ch1711	Chicken ATF6 (Activating Transcription Factor 6) ELISA Kit
E-EL-Ch1712	Chicken FTH (Ferritin, Heavy Polypeptide) ELISA Kit
E-EL-Ch1714	Chicken ECE1 (Endothelin Converting Enzyme 1) ELISA Kit
E-EL-Ch1715	Chicken IFN- β (Interferon Beta)ELISA Kit
E-EL-Ch1716	Chicken LIF (Leukemia Inhibitory Factor) ELISA Kit
E-EL-Ch1718	Chicken PKI α (Protein Kinase Inhibitor Alpha) ELISA Kit
E-EL-Ch1719	Chicken C3 (Complement Component 3) ELISA Kit
E-EL-Ch1722	Chicken Ep-CAM/CD326 (Epithelial Cell Adhesion Molecule) ELISA Kit
E-EL-Ch1723	Chicken ANG1 (Angiopoietin 1) ELISA Kit
E-EL-Ch1725	Chicken F10 (Coagulation Factor X) ELISA Kit
E-EL-Ch1726	Chicken TPA (Tissue Polypeptide Antigen) ELISA Kit
E-EL-Ch1727	Chicken Alk-Smase (Alkaline Sphingomyelinase) ELISA Kit
E-EL-Ch1728	Chicken ANGPTL1 (Angiopoietin Like Protein 1) ELISA Kit

E-EL-Ch1729	Chicken DDC (Dopamine Decarboxylase) ELISA Kit
E-EL-Ch1730	Chicken MRC1 (Mannose Receptor C Type 1 Like Protein 1) ELISA Kit
E-EL-Ch1731	Chicken CMLC-1 (Cardiac Myosin Light Chains 1) ELISA Kit
E-EL-Ch1732	Chicken C4 (Complement Component 4) ELISA Kit
E-EL-Ch1733	Chicken LMP7/PSMB9 (Low Molecular-weight Protein/ proteasome Beta-type Subunit) ELISA Kit
E-EL-Ch1734	Chicken GSTpi (Glutathione S Transferases Pi) ELISA Kit
E-EL-Ch1735	Chicken fT4 (Free Thyroxine) ELISA Kit
E-EL-Ch1736	Chicken ELA3B (Elastase 3B) ELISA Kit
E-EL-Ch1737	Chicken NXF1 (Nuclear RNA Export Factor 1) ELISA Kit
E-EL-Ch1738	Chicken PDI (Protein Disulfide Isomerase) ELISA Kit
E-EL-Ch1739	Chicken NPYR2 (Neuropeptide Y Receptor Y2) ELISA Kit
E-EL-Ch1741	Chicken TRY (Trypsin) ELISA Kit
E-EL-Ch1742	Chicken PARP4 (Poly ADP Ribose Polymerase 4) ELISA Kit
E-EL-Ch1743	Chicken PS2 (Presenilin 2) ELISA Kit
E-EL-Ch1744	Chicken STMN1 (Stathmin 1) ELISA Kit
E-EL-Ch1745	Chicken ANGPT4 (Angiopoietin 4) ELISA Kit
E-EL-Ch1746	Chicken FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit
E-EL-Ch1747	Chicken CLEC2C (C-Type Lectin Domain Family 2 Member C) ELISA Kit
E-EL-Ch1748	Chicken MUSK (Muscle Skeletal Receptor Tyrosine Kinase) ELISA Kit
E-EL-Ch1749	Chicken TP63 (Tumor Protein p63) ELISA Kit
E-EL-Ch1751	Chicken PRSS1 (Protease, Serine, 1) ELISA Kit

E-EL-Ch1752	Chicken ASPN (Asporin) ELISA Kit
E-EL-Ch1753	Chicken AQP-2 (Aquaporin 2) ELISA Kit
E-EL-Ch1754	Chicken SELP (P-Selectin) ELISA Kit
E-EL-Ch1755	Chicken α -SCA (α -sarcomeric Actin) ELISA Kit
E-EL-Ch1756	Chicken PLB (Phospholipase B) ELISA Kit
E-EL-Ch1757	Chicken SIGLEC12 (Sialic Acid Binding Ig Like Lectin 12) ELISA Kit
E-EL-Ch1758	Chicken NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) ELISA Kit
E-EL-Ch1759	Chicken GC (Glucagon) ELISA Kit
E-EL-Ch1761	Chicken CTSK (Cathepsin K) ELISA Kit
E-EL-Ch1764	Chicken COL5 α 2 (Collagen Type V Alpha 2) ELISA Kit
E-EL-Ch1765	Chicken sAC (Soluble Adenylate Cyclase) ELISA Kit
E-EL-Ch1766	Chicken LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit
E-EL-Ch1767	Chicken IL-1 α (Interleukin 1 Alpha) ELISA Kit
E-EL-Ch1769	Chicken α -GST (α -glutathione S-transferases) ELISA Kit
E-EL-Ch1770	Chicken 11-DH-TXB2 (11-Dehydro-thromboxane B2) ELISA Kit
E-EL-Ch1771	Chicken CD163 (Cluster of Differentiation 163) ELISA Kit
E-EL-Ch1774	Chicken PAR2 (Protease Activated Receptor 2) ELISA Kit
E-EL-Ch1775	Chicken UBE3A (Ubiquitin Protein Ligase E3A) ELISA Kit
E-EL-Ch1776	Chicken PKI γ (Protein Kinase Inhibitor Gamma) ELISA Kit
E-EL-Ch1777	Chicken FBLN5 (Fibulin 5) ELISA Kit
E-EL-Ch1778	Chicken M-AChR (Muscarinic Acetylcholine Receptor) ELISA Kit

E-EL-Ch1779	Chicken BK (Bradykinin) ELISA Kit
E-EL-Ch1780	Chicken ELA4 (Elastase 4) ELISA Kit
E-EL-Ch1781	Chicken EAR2 (Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit
E-EL-Ch1782	Chicken AGGF1 (Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit
E-EL-Ch1783	Chicken NT-ProBNP (N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
E-EL-Ch1784	Chicken GP6 (Glycoprotein VI, Platelet) ELISA Kit
E-EL-Ch1785	Chicken PICK1 (Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
E-EL-Ch1786	Chicken PKP2 (Plakophilin 2) ELISA Kit
E-EL-Ch1787	Chicken CLDN4 (Claudin 4) ELISA Kit
E-EL-Ch1789	Chicken N-MID-OT (N-MID Osteocalcin) ELISA Kit
E-EL-Ch1790	Chicken MAO (Monoamine Oxidase) ELISA Kit
E-EL-Ch1792	Chicken GPIIb/IIIa (Platelet Membrane Glycoprotein IIb/IIIa) ELISA Kit
E-EL-Ch1793	Chicken ASK-1 (Apoptosis Signal Regulating Kinase 1) ELISA Kit
E-EL-Ch1794	Chicken OJ/IleRS (Anti-OJ-Antibody) ELISA Kit
E-EL-Ch1795	Chicken GATA2 (GATA Binding Protein 2) ELISA Kit
E-EL-Ch1796	Chicken Flt3 (FMS Like Tyrosine Kinase 3) ELISA Kit
E-EL-Ch1797	Chicken BID (BH3 Interacting Domain Death Agonist) ELISA Kit
E-EL-Ch1798	Chicken HSF1 (Heat Shock Transcription Factor 1) ELISA Kit
E-EL-Ch1799	Chicken Rock-2 (Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit
E-EL-Ch1800	Chicken C7 (Complement Component 7) ELISA Kit
E-EL-Ch1801	Chicken MMP-9 (Matrix Metalloproteinase 9) ELISA Kit

E-EL-Ch1802	Chicken GOD (Glucose Oxidase) ELISA Kit
E-EL-Ch1803	Chicken BAFFR (B Cell Activating Factor Receptor) ELISA Kit
E-EL-Ch1804	Chicken CNTNAP4 (Contactin Associated Protein 4) ELISA Kit
E-EL-Ch1805	Chicken MFAP4 (Microfibrillar Associated Protein 4) ELISA Kit
E-EL-Ch1806	Chicken PKC η (Protein Kinase C Eta) ELISA Kit
E-EL-Ch1807	Chicken SEMA3A (Semaphorin 3A) ELISA Kit
E-EL-Ch1808	Chicken LGALS3BP (Lectin Galactoside Binding, Soluble 3 Binding Protein) ELISA Kit
E-EL-Ch1809	Chicken CK (Creatine Kinase) ELISA Kit
E-EL-Ch1810	Chicken CG (Cryoglobulin) ELISA Kit
E-EL-Ch1811	Chicken CYTH1 (Cytohesin 1) ELISA Kit
E-EL-Ch1812	Chicken AMA IgA (Anti-Myelin antibody IgA) ELISA Kit
E-EL-Ch1813	Chicken EIPS2 (Estrogen Induced Protein pS2) ELISA Kit
E-EL-Ch1815	Chicken ACP3/PAP (Prostatic Acid Phosphatase 3) ELISA Kit
E-EL-Ch1816	Chicken Pg (Progesterone) ELISA Kit
E-EL-Ch1817	Chicken MBLR (Mannose Binding Lectin Receptor) ELISA Kit
E-EL-Ch1818	Chicken WARS (Tryptophanyl tRNA Synthetase) ELISA Kit
E-EL-Ch1819	Chicken RNF4 (Ring Finger Protein 4) ELISA Kit
E-EL-Ch1820	Chicken PR3 (Proteinase 3) ELISA Kit
E-EL-Ch1821	Chicken GTC (General Transcription Complex) ELISA Kit
E-EL-Ch1822	Chicken IL-19 (Interleukin 19) ELISA Kit
E-EL-Ch1823	Chicken TSP-1 (Thrombospondin-1) ELISA Kit

E-EL-Ch1824	Chicken PAPPA2 (Pappalysin 2) ELISA Kit
E-EL-Ch1825	Chicken FOXC1 (Forkhead Box Protein C1) ELISA Kit
E-EL-Ch1826	Chicken HSP-90 (Heat Shock Protein 90) ELISA Kit
E-EL-Ch1827	Chicken PIR (Pirin) ELISA Kit
E-EL-Ch1828	Chicken GSP (Glycosylated Serum Protein) ELISA Kit
E-EL-Ch1829	Chicken VIP (Vasoactive Intestinal Peptide) ELISA Kit
E-EL-Ch1830	Chicken OBTF2 (Octamer Binding Transcription Factor 2) ELISA Kit
E-EL-Ch1831	Chicken A β 1-42 (Amyloid Beta Peptide 1-42) ELISA Kit
E-EL-Ch1833	Chicken PTX 3/TSG-14 (Pentraxin 3, Long) ELISA Kit
E-EL-Ch1834	Chicken MUC5AC (Mucin 5 Subtype AC) ELISA Kit
E-EL-Ch1835	Chicken PR3Ab (Proteinase 3 Antibody) ELISA Kit
E-EL-Ch1836	Chicken ODC (Ornithine Decarboxylase) ELISA Kit
E-EL-Ch1837	Chicken TARC (Thymus Activation Regulated Chemokine) ELISA Kit
E-EL-Ch1838	Chicken BECN1 (Beclin 1) ELISA Kit
E-EL-Ch1839	Chicken C5c (Complement C5 Convertase) ELISA Kit
E-EL-Ch1840	Chicken Fbn (Fibrinolysin) ELISA Kit
E-EL-Ch1841	Chicken P-PKC (Phospho Protein Kinase C) ELISA Kit
E-EL-Ch1842	Chicken NES1 (Nesfatin 1) ELISA Kit
E-EL-Ch1843	Chicken UPP1 (Uridine Phosphorylase 1) ELISA Kit
E-EL-Ch1844	Chicken PKC β 1 (Protein Kinase C Beta 1) ELISA Kit
E-EL-Ch1846	Chicken TEP1 (Telomerase Associated Protein 1) ELISA Kit

E-EL-Ch1847	Chicken ADIPOR2 (Adiponectin Receptor 2) ELISA Kit
E-EL-Ch1848	Chicken GL α (Galactosidase Alpha) ELISA Kit
E-EL-Ch1849	Chicken SNRPD1 (Small Nuclear Ribonucleoprotein Polypeptide D1) ELISA Kit
E-EL-Ch1850	Chicken GT-34 (Gastrin 34) ELISA Kit
E-EL-Ch1851	Chicken ACTN2 (Actinin Alpha 2) ELISA Kit
E-EL-Ch1852	Chicken CA4 (Carbonic Anhydrase IV) ELISA Kit
E-EL-Ch1853	Chicken CTACK (Cutaneous T-cell Attracting Chemokine) ELISA Kit
E-EL-Ch1854	Chicken SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-Ch1857	Chicken PIGR/SC (Polymeric Immunoglobulin Receptor/ Membrane Secretory Component) ELISA Kit
E-EL-Ch1858	Chicken Ponticulin ELISA Kit
E-EL-Ch1859	Chicken OBTF4 (Octamer Binding Transcription Factor 4) ELISA Kit
E-EL-Ch1860	Chicken GA (Glycated Albumin) ELISA Kit
E-EL-Ch1861	Chicken PLVAP (Plasmalemma Vesicle Associated Protein) ELISA Kit
E-EL-Ch1862	Chicken AMH (Anti-Mullerian Hormone) ELISA Kit
E-EL-Ch1864	Chicken PAIg (Platelet Associated Immunoglobulins) ELISA Kit
E-EL-Ch1865	Chicken STK39 (Serine/Threonine Kinase 39) ELISA Kit
E-EL-Ch1866	Chicken GHR (Growth Hormone Receptor) ELISA Kit
E-EL-Ch1867	Chicken FG β (Fibrinogen Beta) ELISA Kit
E-EL-Ch1868	Chicken HB β (Hemoglobin Beta) ELISA Kit
E-EL-Ch1869	Chicken ADM (Adrenomedullin) ELISA Kit
E-EL-Ch1870	Chicken TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit

E-EL-Ch1871	Chicken SEMA3B (Semaphorin 3B) ELISA Kit
E-EL-Ch1872	Chicken PTPRN (Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit
E-EL-Ch1873	Chicken PUMA (P53 Upregulated Modulator of Apoptosis) ELISA Kit
E-EL-Ch1874	Chicken DSP (Desmoplakin) ELISA Kit
E-EL-Ch1875	Chicken LFA-1 (Lymphocyte Function Associated Antigen 1) ELISA Kit
E-EL-Ch1876	Chicken TSC22 (Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
E-EL-Ch1877	Chicken EL/LIPG (Endothelial Lipase) ELISA Kit
E-EL-Ch1878	Chicken α 2-AP (α 2-Antiplasmin) ELISA Kit
E-EL-Ch1879	Chicken BGN (Biglycan) ELISA Kit
E-EL-Ch1880	Chicken STAT6 (Signal Transducer and Activator of Transcription 6) ELISA Kit
E-EL-Ch1881	Chicken TPPP (Tubulin Polymerization Promoting Protein) ELISA Kit
E-EL-Ch1882	Chicken TCR (Tcell Receptor) ELISA Kit
E-EL-Ch1883	Chicken fT3 (Free Triiodothyronine) ELISA Kit
E-EL-Ch1884	Chicken MMP-7 (Matrix Metalloproteinase 7) ELISA Kit
E-EL-Ch1886	Chicken PCT (Procalcitonin) ELISA Kit
E-EL-Ch1889	Chicken MMP-8 (Matrix Metalloproteinase 8) ELISA Kit
E-EL-Ch1890	Chicken EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit
E-EL-Ch1891	Chicken PON1 (Paraoxonase 1) ELISA Kit
E-EL-Ch1892	Chicken IgA2 (Immunoglobulin A2) ELISA Kit
E-EL-Ch1894	Chicken TEF (Thyrotroph Embryonic Factor) ELISA Kit
E-EL-Ch1895	Chicken NAGLU (N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit

E-EL-Ch1896	Chicken HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit
E-EL-Ch1898	Chicken MCT (Mast Cell Tryptase) ELISA Kit
E-EL-Ch1899	Chicken ConA (Concanavalin A) ELISA Kit
E-EL-Ch1900	Chicken PTPRB (Protein Tyrosine Phosphatase Receptor Type B) ELISA Kit
E-EL-Ch1901	Chicken CYP7A1 (Cytochrome P450 7A1) ELISA Kit
E-EL-Ch1902	Chicken β 2-GP1 Ab IgM (β 2-Glycoprotein 1 Antibody IgM) ELISA Kit
E-EL-Ch1903	Chicken Cav-1 (Caveolin-1) ELISA Kit
E-EL-Ch1904	Chicken CLA (Collagen Autoantibody) ELISA Kit
E-EL-Ch1905	Chicken FGL1 (Fibrinogen Like Protein 1) ELISA Kit
E-EL-Ch1906	Chicken PA2G4 (Proliferation Associated Protein 2G4) ELISA Kit
E-EL-Ch1907	Chicken β 2-GP1 Ab IgA (β 2-Glycoprotein 1 antibody IgA) ELISA Kit
E-EL-Ch1908	Chicken ANGPTL2 (Angiopoietin Like Protein 2) ELISA Kit
E-EL-Ch1909	Chicken Extein ELISA Kit
E-EL-Ch1910	Chicken TACR3 (Tachykinin Receptor 3) ELISA Kit
E-EL-Ch1911	Chicken Anti-AChR (Anti-Acetylcholine Receptor Antibody) ELISA Kit
E-EL-Ch1912	Chicken CRLF1 (Cytokine Receptor Like Factor 1) ELISA Kit
E-EL-Ch1913	Chicken CCKBR (Cholecystokinin B Receptor) ELISA Kit
E-EL-Ch1914	Chicken MMRN1 (Multimerin 1) ELISA Kit
E-EL-Ch1915	Chicken ADCY1 (Adenylate Cyclase 1, Brain) ELISA Kit
E-EL-Ch1916	Chicken LRG1 (Leucine Rich Alpha-2-Glycoprotein 1) ELISA Kit
E-EL-Ch1917	Chicken Smad1/MADH1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit

E-EL-Ch1918	Chicken I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit
E-EL-Ch1919	Chicken TNKS2 (Tankyrase 2) ELISA Kit
E-EL-Ch1920	Chicken GPC1 (Glypican 1) ELISA Kit
E-EL-Ch1922	Chicken ANXA5 (Annexin A5) ELISA Kit
E-EL-Ch1923	Chicken VEGF165 (Vascular Endothelial Growth Factor165) ELISA Kit
E-EL-Ch1924	Chicken c-fos ELISA Kit
E-EL-Ch1925	Chicken MASP2 (Mannan Associated Serine Protease 2) ELISA Kit
E-EL-Ch1926	Chicken F2 (Coagulation Factor II) ELISA Kit
E-EL-Ch1927	Chicken CASP3 (Caspase 3) ELISA Kit
E-EL-Ch1928	Chicken BRAK/CXCL14 (Breast and Kidney Expressed Chemokine) ELISA Kit
E-EL-Ch1929	Chicken LR (Leukoregulin) ELISA Kit
E-EL-Ch1930	Chicken CAM-Ab (Anti-Calmodulin Specific Antibody) ELISA Kit
E-EL-Ch1931	Chicken TR (Testosterone Receptor) ELISA Kit
E-EL-Ch1932	Chicken PKR1 (Prokineticin Receptor 1) ELISA Kit
E-EL-Ch1935	Chicken ADRP (Adipose Differentiation Related Protein) ELISA Kit
E-EL-Ch1937	Chicken GAD1 (Glutamate Decarboxylase 1) ELISA Kit
E-EL-Ch1940	Chicken eIF3a (Eukaryotic Translation Initiation Factor 3a) ELISA Kit
E-EL-Ch1941	Chicken Pax (Paxillin) ELISA Kit
E-EL-Ch1942	Chicken FABP5 (Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
E-EL-Ch1943	Chicken MT1 (Metallothionein 1) ELISA Kit
E-EL-Ch1946	Chicken CHRM2 (Cholinergic Receptor, Muscarinic 2) ELISA Kit

E-EL-Ch1947	Chicken SBEM (Small Breast Epithelial Mucin) ELISA Kit
E-EL-Ch1948	Chicken CCR1 (CC-Chemokine Receptor 1) ELISA Kit
E-EL-Ch1949	Chicken SPC (Pulmonary Surfactant Associated Protein C) ELISA Kit
E-EL-Ch1950	Chicken MAU (Microalbuminuria) ELISA Kit
E-EL-Ch1951	Chicken CCP (Complement Control Protein) ELISA Kit
E-EL-Ch1953	Chicken IgG1 (Immunoglobulin G1) ELISA Kit
E-EL-Ch1954	Chicken MCF (Macrophage Chemotatic Factor) ELISA Kit
E-EL-Ch1955	Chicken NGB (Neuroglobin) ELISA Kit
E-EL-Ch1957	Chicken DS (Dermatan Sulfate) ELISA Kit
E-EL-Ch1959	Chicken TIMP-4 (Tissue Inhibitors of Metalloproteinase 4) ELISA Kit
E-EL-Ch1960	Chicken TNR (Tenascin R) ELISA Kit
E-EL-Ch1962	Chicken PCDHβ2 (Protocadherin Beta 2) ELISA Kit
E-EL-Ch1964	Chicken CK-20/KRT20 (Cytokeratin 20) ELISA Kit
E-EL-Ch1965	Chicken FOXP1 (Forkhead Box Protein P1) ELISA Kit
E-EL-Ch1966	Chicken OLAB (Anti-Oxidized Low Density Lipoprotein Antibody) ELISA Kit
E-EL-Ch1967	Chicken IL-24 (Interleukin 24) ELISA Kit
E-EL-Ch1968	Chicken GADL1 (Glutamate Decarboxylase Like Protein 1) ELISA Kit
E-EL-Ch1969	Chicken PKC (Protein Kinase C) ELISA Kit
E-EL-Ch1970	Chicken FPA (Fibrinopeptide A) ELISA Kit
E-EL-Ch1971	Chicken MUC20 (Mucin 20) ELISA Kit
E-EL-Ch1972	Chicken IGFALS (Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) ELISA Kit

E-EL-Ch1974	Chicken DEF α 5 (Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
E-EL-Ch1975	Chicken RHOB (Ras Homolog Gene Family, Member B) ELISA Kit
E-EL-Ch1976	Chicken REG1 β (Regenerating Islet Derived Protein 1 Beta) ELISA Kit
E-EL-Ch1977	Chicken CBX3 (Chromobox Homolog 3) ELISA Kit
E-EL-Ch1978	Chicken PDGFA (Platelet Derived Growth Factor Subunit A) ELISA Kit
E-EL-Ch1979	Chicken NUP107 (Nucleoporin 107kDa) ELISA Kit
E-EL-Ch1981	Chicken GCLM (Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit
E-EL-Ch1982	Chicken TFR2 (Transferrin Receptor 2) ELISA Kit
E-EL-Ch1983	Chicken CENPI (Centromere Protein I) ELISA Kit
E-EL-Ch1984	Chicken VS (Versican/Pg-M) ELISA Kit
E-EL-Ch1985	Chicken RELB (V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit
E-EL-Ch1987	Chicken PCDH15 (Protocadherin 15) ELISA Kit
E-EL-Ch1988	Chicken NF- κ B p65 (Nuclear Factor- κ B p65) ELISA Kit
E-EL-Ch1989	Chicken TF (Tissue Factor) ELISA Kit
E-EL-Ch1990	Chicken MAP-2 (Microtubule Associated Protein 2) ELISA Kit
E-EL-Ch1991	Chicken OSCAR (Osteoclast Associated Receptor) ELISA Kit
E-EL-Ch1992	Chicken GALR1 (Galanin Receptor 1) ELISA Kit
E-EL-Ch1993	Chicken BACE1 (Beta-site APP Cleaving Enzyme 1) ELISA Kit
E-EL-Ch1994	Chicken TAT (Thrombin/Antithrombin Complex) ELISA Kit
E-EL-Ch1995	Chicken NGF (Nerve Growth Factor) ELISA Kit
E-EL-Ch1996	Chicken PDPN (Podoplanin) ELISA Kit

E-EL-Ch1997	Chicken TNC (Tenascin C) ELISA Kit
E-EL-Ch1998	Chicken CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit
E-EL-Ch1999	Chicken NFκB (Nuclear Factor Kappa B) ELISA Kit
E-EL-Ch2000	Chicken GAL3 (Galectin 3) ELISA Kit
E-EL-Ch2001	Chicken EMP (Erythrocyte Membrane Protein) ELISA Kit
E-EL-Ch2002	Chicken RNASEN (Ribonuclease III, Nuclear) ELISA Kit
E-EL-Ch2003	Chicken IAPP (Islet Amyloid Polypeptide) ELISA Kit
E-EL-Ch2004	Chicken RCOR3 (Rest Corepressor 3) ELISA Kit
E-EL-Ch2006	Chicken MCRS1 (Microspherule Protein 1) ELISA Kit
E-EL-Ch2007	Chicken EGR4 (Early Growth Response Protein 4) ELISA Kit
E-EL-Ch2008	Chicken CRHBP (Corticotropin Releasing Hormone Binding Protein) ELISA Kit
E-EL-Ch2011	Chicken HbA1c (Glycosylated Hemoglobin/Hemoglobin A1c) ELISA Kit
E-EL-Ch2012	Chicken GSTα4 (Glutathione S Transferase Alpha 4) ELISA Kit
E-EL-Ch2013	Chicken EP3 (Prostaglandin E Receptor 3) ELISA Kit
E-EL-Ch2016	Chicken PXR (Pregnane X Receptor) ELISA Kit
E-EL-Ch2017	Chicken CDK4 (Cyclin Dependent Kinase 4) ELISA Kit
E-EL-Ch2018	Chicken FcγBP (Fc Fragment of IgG Binding Protein) ELISA Kit
E-EL-Ch2019	Chicken RNASE (Ribonuclease) ELISA Kit
E-EL-Ch2020	Chicken NTRK3 (Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
E-EL-Ch2022	Chicken GSTμ2/GSTm2 (Glutathione S Transferase Mu 2, Muscle) ELISA Kit
E-EL-Ch2023	Chicken PDGF-AA (Platelet Derived Growth Factor-AA) ELISA Kit

E-EL-Ch2024	Chicken SNRPC (Small Nuclear Ribonucleoprotein Polypeptide C) ELISA Kit
E-EL-Ch2025	Chicken AIRE (Autoimmune Regulator) ELISA Kit
E-EL-Ch2026	Chicken PTGS1 (Prostaglandin Endoperoxide Synthase 1) ELISA Kit
E-EL-Ch2027	Chicken PYY (Peptide YY) ELISA Kit
E-EL-Ch2028	Chicken CTSB (Cathepsin B) ELISA Kit
E-EL-Ch2029	Chicken PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit
E-EL-Ch2030	Chicken ELN (Elastin) ELISA Kit
E-EL-Ch2031	Chicken PDIA5 (Protein Disulfide Isomerase A5) ELISA Kit
E-EL-Ch2032	Chicken FGF4 (Fibroblast Growth Factor 4) ELISA Kit
E-EL-Ch2034	Chicken EFNA1 (Ephrin A1) ELISA Kit
E-EL-Ch2037	Chicken CP (Ceruloplasmin) ELISA Kit
E-EL-Ch2038	Chicken GNβ1 (G Protein Beta 1) ELISA Kit
E-EL-Ch2040	Chicken CABIN1 (Calcineurin Binding Protein 1) ELISA Kit
E-EL-Ch2041	Chicken AP 17 (Apelin 17) ELISA Kit
E-EL-Ch2042	Chicken PKD1 (Protein Kinase D1) ELISA Kit
E-EL-Ch2043	Chicken GAD2 (Glutamate Decarboxylase 2) ELISA Kit
E-EL-Ch2044	Chicken CYP1A2 (Cytochrome P450 1A2) ELISA Kit
E-EL-Ch2046	Chicken GFAP (Glial Fibrillary Acidic Protein) ELISA Kit
E-EL-Ch2047	Chicken ACAT2 (Acetyl Coenzyme A Acetyltransferase 2) ELISA Kit
E-EL-Ch2048	Chicken IGFBP-2 (Insulin Like Growth Factor Binding Protein 2) ELISA Kit
E-EL-Ch2049	Chicken SYNPO (Synaptopodin) ELISA Kit

E-EL-Ch2050	Chicken BCL2L1/BCL-X (Bcl-2 Like Protein 1) ELISA Kit
E-EL-Ch2051	Chicken GS (Gelsolin) ELISA Kit
E-EL-Ch2052	Chicken FG (Fibrinogen) ELISA Kit
E-EL-Ch2053	Chicken MPF (Maturation Promoting Factor) ELISA Kit
E-EL-Ch2054	Chicken PPOX (Protoporphyrinogen Oxidase) ELISA Kit
E-EL-Ch2055	Chicken CA19-9 (Carbohydrate antigen19-9) ELISA Kit
E-EL-Ch2056	Chicken AAP (Alanine Aminopeptidase) ELISA Kit
E-EL-Ch2058	Chicken CAMK2 α (Calcium/Calmodulin Dependent Protein Kinase II Alpha) ELISA Kit
E-EL-Ch2059	Chicken TGF- β 2 (Transforming Growth Factor β 2) ELISA Kit
E-EL-Ch2060	Chicken PTGES (Prostaglandin E Synthase, Microsomal) ELISA Kit
E-EL-Ch2061	Chicken CCL16 (Chemokine C-C-Motif Ligand 16) ELISA Kit
E-EL-Ch2063	Chicken GMCSF/CSF2(Granulocyte Macrophage Colony Stimulating Factor 2) ELISA Kit
E-EL-Ch2064	Chicken MDC/CCL22 (Macrophage Derived Chemokine) ELISA Kit
E-EL-Ch2066	Chicken GR (Glutathione Reductase) ELISA Kit
E-EL-Ch2067	Chicken TA (Anti-Teichoic Acid Antibody) ELISA Kit
E-EL-Ch2068	Chicken β -TG (β -Thromboglobulin) ELISA Kit
E-EL-Ch2069	Chicken FGFR1 (Fibroblast Growth Factor Receptor 1) ELISA Kit
E-EL-Ch2070	Chicken AMY1 (Amylase Alpha 1, Salivary) ELISA Kit
E-EL-Ch2072	Chicken MIA1 (Melanoma Inhibitory Activity Protein 1) ELISA Kit
E-EL-Ch2073	Chicken sIgA (Secretory Immunoglobulin A) ELISA Kit
E-EL-Ch2074	Chicken CR (Calretinin) ELISA Kit

E-EL-Ch2075	Chicken PF4V1 (Platelet Factor 4 Variant 1) ELISA Kit
E-EL-Ch2076	Chicken AST (Aspartate Aminotransferase) ELISA Kit
E-EL-Ch2078	Chicken F13A1 (Coagulation Factor XIII A1 Polypeptide) ELISA Kit
E-EL-Ch2079	Chicken MF/MPF (Maturation Promoting Factor) ELISA Kit
E-EL-Ch2080	Chicken MAPK (Mitogen Activated Protein Kinase) ELISA Kit
E-EL-Ch2082	Chicken NUP37 (Nucleoporin 37kDa) ELISA Kit
E-EL-Ch2083	Chicken MFG-E8 (Milk Fat Globule EGF Factor 8) ELISA Kit
E-EL-Ch2084	Chicken F-TESTO (Free Testosterone) ELISA Kit
E-EL-Ch2086	Chicken MSM β (Microseminoprotein Beta) ELISA Kit
E-EL-Ch2087	Chicken TPP1 (Tripeptidyl Peptidase I) ELISA Kit
E-EL-Ch2088	Chicken BMP-3 (Bone Morphogenetic Protein 3) ELISA Kit
E-EL-Ch2089	Chicken AKA (Anti-Keratin Antibody) ELISA Kit
E-EL-Ch2090	Chicken PIINP (Procollagen II N-Terminal ProPeptide) ELISA Kit
E-EL-Ch2091	Chicken MAP τ (Microtubule Associated Protein Tau/Tau Protein) ELISA Kit
E-EL-Ch2092	Chicken PRKA α 1 (Protein Kinase, AMP Activated Alpha 1) ELISA Kit
E-EL-Ch2093	Chicken NPYR1 (Neuropeptide Y Receptor Y1) ELISA Kit
E-EL-Ch2096	Chicken Rac-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit
E-EL-Ch2097	Chicken SEMA7A (Semaphorin 7A) ELISA Kit
E-EL-Ch2098	Chicken FGFBP2 (Fibroblast Growth Factor Binding Protein 2) ELISA Kit
E-EL-Ch2099	Chicken ESF (Erythropoiesis Stimulating Factor) ELISA Kit
E-EL-Ch2101	Chicken CLL1 (Collectin Liver 1) ELISA Kit

E-EL-Ch2102	Chicken TGM1 (Transglutaminase 1, Keratinocyte) ELISA Kit
E-EL-Ch2104	Chicken PRSS8 (Protease, Serine 8) ELISA Kit
E-EL-Ch2105	Chicken PDIA2 (Protein Disulfide Isomerase A2) ELISA Kit
E-EL-Ch2106	Chicken ERK3 (Extracellular Signal Regulated Kinase 3) ELISA Kit
E-EL-Ch2107	Chicken MT3 (Metallothionein 3) ELISA Kit
E-EL-Ch2108	Chicken C4b (Complement Component 4b) ELISA Kit
E-EL-Ch2109	Chicken GAD-Ab (Anti-Glutamic Acid Decarboxylase antibody) ELISA Kit
E-EL-Ch2110	Chicken LCN1 (Lipocalin 1) ELISA Kit
E-EL-Ch2111	Chicken ATF7 (Activating Transcription Factor 7) ELISA Kit
E-EL-Ch2113	Chicken PI3K2 α (Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit
E-EL-Ch2114	Chicken GJ β 1 (Gap Junction Protein Beta 1) ELISA Kit
E-EL-Ch2115	Chicken CYFRA21-1 (Cytokeratin Fragment Antigen 21-1) ELISA Kit
E-EL-Ch2116	Chicken GH (Growth Hormone) ELISA Kit
E-EL-Ch2117	Chicken RETN (Resistin) ELISA Kit
E-EL-Ch2118	Chicken HP (Haptoglobin) ELISA Kit
E-EL-Ch2120	Chicken SSTR1 (Somatostatin Receptor 1) ELISA Kit
E-EL-Ch2121	Chicken IL-17 (Interleukin 17) ELISA Kit
E-EL-Ch2122	Chicken FGF21(Fibroblast Growth Factor 21) ELISA Kit
E-EL-Ch2123	Chicken MSP (Macrophage Stimulating Protein) ELISA Kit
E-EL-Ch2124	Chicken CENPH (Centromere Protein H) ELISA Kit
E-EL-Ch2125	Chicken SOX18 (Sex Determining Region Y Box Protein 18) ELISA Kit

E-EL-Ch2126	Chicken ANG III (Angiotensin III) ELISA Kit
E-EL-Ch2127	Chicken PC4 (Positive Cofactor 4) ELISA Kit
E-EL-Ch2128	Chicken SYP (Synaptophysin) ELISA Kit
E-EL-Ch2129	Chicken MAF (Macrophage Activating Factor) ELISA Kit
E-EL-Ch2130	Chicken CA12 (Carbonic Anhydrase XII) ELISA Kit
E-EL-Ch2131	Chicken MMP-12 (Matrix Metalloproteinase 12) ELISA Kit
E-EL-Ch2132	Chicken TGFβ1 (Transforming Growth Factor Beta Induced Protein) ELISA Kit
E-EL-Ch2133	Chicken S100A11 (S100 Calcium Binding Protein A11) ELISA Kit
E-EL-Ch2134	Chicken EMILIN 1 (Elastin Microfibril Interface Located Protein 1) ELISA Kit
E-EL-Ch2135	Chicken PRSS12 (Protease, Serine 12) ELISA Kit
E-EL-Ch2136	Chicken αMSH (Alpha-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-Ch2138	Chicken Cys-C (Cystatin C) ELISA Kit
E-EL-Ch2139	Chicken NUP98 (Nucleoporin 98kDa) ELISA Kit
E-EL-Ch2140	Chicken CNTNAP2 (Contactin Associated Protein Like Protein 2) ELISA Kit
E-EL-Ch2141	Chicken ACP1 (Acid Phosphatase 1) ELISA Kit
E-EL-Ch2142	Chicken F13B (Coagulation Factor XIII B Polypeptide) ELISA Kit
E-EL-Ch2144	Chicken ASM (Acid Sphingomyelinase) ELISA Kit
E-EL-Ch2146	Chicken IGFBP-4 (Insulin Like Growth Factor Binding Protein 4) ELISA Kit
E-EL-Ch2147	Chicken PPIL1 (Peptidylprolyl Isomerase/Cyclophilin Like 1) ELISA Kit
E-EL-Ch2148	Chicken PLTP (Phospholipid Transfer Protein) ELISA Kit
E-EL-Ch2149	Chicken KLK7 (Kallikrein 7) ELISA Kit

E-EL-Ch2151	Chicken CA6 (Carbonic Anhydrase VI) ELISA Kit
E-EL-Ch2152	Chicken UCP-2 (Uncoupling Protein 2, Mitochondrial) ELISA Kit
E-EL-Ch2153	Chicken GMP-140 (Alpha-granular Membrane Protein) ELISA Kit
E-EL-Ch2154	Chicken MBP (Myelin Basic Protein) ELISA Kit
E-EL-Ch2155	Chicken TM (Thrombin)ELISA Kit
E-EL-Ch2156	Chicken TAF13 (TATA Box Binding Protein Associated Factor 13) ELISA Kit
E-EL-Ch2157	Chicken CASP10 (Caspase 10) ELISA Kit
E-EL-Ch2159	Chicken GIP (Gastric Inhibitory Polypeptide) ELISA Kit
E-EL-Ch2160	Chicken MART/Melan-A (Melanoma Marker A) ELISA Kit
E-EL-Ch2161	Chicken TGF-β3 (Transforming Growth Factor β3) ELISA Kit
E-EL-Ch2162	Chicken CA9 (Carbonic Anhydrase IX) ELISA Kit
E-EL-Ch2163	Chicken SEP1 (Selenoprotein 1) ELISA Kit
E-EL-Ch2164	Chicken ITLN1 (Intelectin 1/Omentin) ELISA Kit
E-EL-Ch2165	Chicken CX37 (Connexin 37) ELISA Kit
E-EL-Ch2166	Chicken GBP4 (Guanylate Binding Protein 4) ELISA Kit
E-EL-Ch2168	Chicken TM (Thrombomodulin) ELISA Kit
E-EL-Ch2169	Chicken MAG (Myelin Associated Glycoprotein) ELISA Kit
E-EL-Ch2170	Chicken CX31 (Connexin 31) ELISA Kit
E-EL-Ch2172	Chicken Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit
E-EL-Ch2173	Chicken GPR120 (G Protein Coupled Receptor 120) ELISA Kit
E-EL-Ch2174	Chicken TSTA (Tumor Specific Transplantation Antigen) ELISA Kit

E-EL-Ch2175	Chicken Anti-GM1 (Anti-Ganglioside M1 Antibody) ELISA Kit
E-EL-Ch2176	Chicken MT1M (Metallothionein 1M) ELISA Kit
E-EL-Ch2177	Chicken ES (Endostatin) ELISA Kit
E-EL-Ch2178	Chicken TAFI (Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit
E-EL-Ch2179	Chicken ODF (Osteoclast Differentiation Factor) ELISA Kit
E-EL-Ch2181	Chicken HO1 (Heme Oxygenase 1) ELISA Kit
E-EL-Ch2182	Chicken LDHB (Lactate Dehydrogenase B) ELISA Kit
E-EL-Ch2183	Chicken HJV/HFE2 (Hemojuvelin/Hemochromatosis Type 2) ELISA Kit
E-EL-Ch2185	Chicken IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit
E-EL-Ch2186	Chicken MGST1 (Microsomal Glutathione S Transferase 1) ELISA Kit
E-EL-Ch2187	Chicken EG-VEGF (Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit
E-EL-Ch2188	Chicken FABP6 (Fatty Acid Binding Protein 6, Ileal) ELISA Kit
E-EL-Ch2189	Chicken RPA-70 (replication protein A) ELISA Kit
E-EL-Ch2190	Chicken TRH (Thyrotropin Releasing Hormone) ELISA Kit
E-EL-Ch2191	Chicken CLS (Cardiolipin Synthase) ELISA Kit
E-EL-Ch2193	Chicken FBN3 (Fibrillin 3) ELISA Kit
E-EL-Ch2194	Chicken CTSD (Cathepsin D) ELISA Kit
E-EL-Ch2195	Chicken CTX-III (Cross Linked C-telopeptide of Type III Collagen) ELISA Kit
E-EL-Ch2196	Chicken NT (Neurotensin) ELISA Kit
E-EL-Ch2197	Chicken ATXN2 (Ataxin 2) ELISA Kit
E-EL-Ch2198	Chicken TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit

E-EL-Ch2199	Chicken CNTNAP1 (Contactin Associated Protein 1) ELISA Kit
E-EL-Ch2200	Chicken FGF23 (Fibroblast Growth Factor 23) ELISA Kit
E-EL-Ch2201	Chicken AP (Aprotinin) ELISA Kit
E-EL-Ch2202	Chicken PAP (Plasmin-Antiplasmin Complex) ELISA Kit
E-EL-Ch2203	Chicken FRS3 (Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
E-EL-Ch2204	Chicken SDC1 (Syndecan 1) ELISA Kit
E-EL-Ch2205	Chicken CD42 (Cluster of Differentiation 42) ELISA Kit
E-EL-Ch2206	Chicken SELL (L-Selectin) ELISA Kit
E-EL-Ch2207	Chicken MUC5B (Mucin 5 Subtype B) ELISA Kit
E-EL-Ch2208	Chicken LP-a (Lipoprotein a) ELISA Kit
E-EL-Ch2209	Chicken FAS/CD95/TNFRSF6 (Factor Related Apoptosis) ELISA Kit
E-EL-Ch2210	Chicken OSM (Oncostatin M) ELISA Kit
E-EL-Ch2211	Chicken PS (Plasma Anticoagulant Protein S) ELISA Kit
E-EL-Ch2212	Chicken TARS (Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit
E-EL-Ch2213	Chicken ATF-3 (Activating Transcription Factor 3) ELISA Kit
E-EL-Ch2215	Chicken NFKBIE (Nuclear Factor Of Kappa Light Polypeptide Gene Enhancer In B-cells Inhibitor Epsilon) ELISA Kit
E-EL-Ch2216	Chicken C8g (Complement Component 8g) ELISA Kit
E-EL-Ch2217	Chicken VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
E-EL-Ch2218	Chicken ANGPTL3 (Angiopoietin Like Protein 3) ELISA Kit
E-EL-Ch2221	Chicken Fcγ (Fc Fragment of IgG) ELISA Kit
E-EL-Ch2222	Chicken TD-Ag (Thymus Dependent Antigen) ELISA Kit

E-EL-Ch2223	Chicken CG α (Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit
E-EL-Ch2224	Chicken KERA (Keratocan) ELISA Kit
E-EL-Ch2226	Chicken AIF (Apoptosis Inducing Factor) ELISA Kit
E-EL-Ch2227	Chicken β Manase ELISA Kit
E-EL-Ch2228	Chicken TSLP (Thymic Stromal Lymphopoietin) ELISA Kit
E-EL-Ch2229	Chicken ADCY4 (Adenylate Cyclase 4) ELISA Kit
E-EL-Ch2231	Chicken α 1-AGP (Alpha-1-Acid Glycoprotein) ELISA Kit
E-EL-Ch2232	Chicken PTPRK (Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit
E-EL-Ch2233	Chicken NRG-1 (Neuregulin 1) ELISA Kit
E-EL-Ch2234	Chicken OFQ/N (Orphanin FQ/Nociceptin) ELISA Kit
E-EL-Ch2235	Chicken IDU α (Alpha-L-Iduronidase) ELISA Kit
E-EL-Ch2236	Chicken α 2-M (Alpha-2 Macroglobulin) ELISA Kit
E-EL-Ch2237	Chicken F8 (Coagulation Factor VIII) ELISA Kit
E-EL-Ch2238	Chicken ST2 (Syntenin 2) ELISA Kit
E-EL-Ch2239	Chicken cTnT/TNNT2 (Troponin T Type 2, Cardiac) ELISA Kit
E-EL-Ch2240	Chicken CDK5 (Cyclin Dependent Kinase 5) ELISA Kit
E-EL-Ch2241	Chicken MUC21 (Mucin 21) ELISA Kit
E-EL-Ch2242	Chicken MAdCAM1 (Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit
E-EL-Ch2244	Chicken AXIN2 (Axis Inhibition Protein 2) ELISA Kit
E-EL-Ch2245	Chicken BSP (Bone Sialoprotein) ELISA Kit
E-EL-Ch2246	Chicken PI (Proinsulin) ELISA Kit

E-EL-Ch2247	Chicken COMP (Cartilage Oligomeric Matrix Protein) ELISA Kit
E-EL-Ch2248	Chicken GUCA2B (Guanylate Cyclase Activator 2B) ELISA Kit
E-EL-Ch2249	Chicken MR (Mannose Receptor) ELISA Kit
E-EL-Ch2250	Chicken PP (Placental Protein) ELISA Kit
E-EL-Ch2251	Chicken D-LDH (D-Lactate Dehydrogenase) ELISA Kit
E-EL-Ch2252	Chicken IGFBP-7 (Insulin Like Growth Factor Binding Protein 7) ELISA Kit
E-EL-Ch2253	Chicken NP (Nandrolone Phenylpropionate) ELISA Kit
E-EL-Ch2254	Chicken OLFM4 (Olfactomedin 4) ELISA Kit
E-EL-Ch2256	Chicken C5a (Complement Component 5a) ELISA Kit
E-EL-Ch2257	Chicken HB (Hemoglobin) ELISA Kit
E-EL-Ch2258	Chicken ER α (Estrogen Receptor Alpha) ELISA Kit
E-EL-Ch2259	Chicken NLR (NOD-Like Receptor9) ELISA Kit
E-EL-Ch2260	Chicken τ PK1 (Tau-Protein Kinase) ELISA Kit
E-EL-Ch2261	Chicken BTA (Bladder Tumor Antigen) ELISA Kit
E-EL-Ch2262	Chicken ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-Ch2263	Chicken Ta1 (Thymosin- α 1) ELISA Kit
E-EL-Ch2264	Chicken PS1 (Presenilin 1) ELISA Kit
E-EL-Ch2265	Chicken PACAP-38 (Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
E-EL-Ch2266	Chicken PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-Ch2267	Chicken INS (Insulin) ELISA Kit
E-EL-S0001	Sheep RETN (Resistin) ELISA Kit

Cat.NO	Product Name
E-EL-S0002	Sheep S100B (S100 Calcium Binding Protein B) ELISA Kit
E-EL-S0003	Sheep SELP (P-Selectin) ELISA Kit
E-EL-S0004	Sheep β -EP (Beta-Endorphin) ELISA Kit
E-EL-S0005	Sheep IL-10 (Interleukin-10) ELISA Kit
E-EL-S0007	Sheep SAA (Serum Amyloid A) ELISA Kit
E-EL-S0008	Sheep AMPK (Phosphorylated Adenosine Monophosphate Activated Protein Kinase) ELISA Kit
E-EL-S0009	Sheep IFN- γ (Interferon Gamma) ELISA Kit
E-EL-S0010	Sheep IL-2 (Interleukin 2) ELISA Kit
E-EL-S0011	Sheep IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-S0012	Sheep CSTA (Cystatin A) ELISA Kit
E-EL-S0013	Sheep IL-12 (Interleukin 12) ELISA Kit
E-EL-S0014	Sheep BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit
E-EL-S0015	Sheep IL-2R (Interleukin 2 Receptor) ELISA Kit
E-EL-S0016	Sheep OPG (Osteoprotegerin) ELISA Kit
E-EL-S0017	Sheep IL-6 (Interleukin 6) ELISA Kit
E-EL-S0018	Sheep IL-18 (Interleukin 18) ELISA Kit
E-EL-S0019	Sheep CKLFSF5 (Chemokine Like Factor Superfamily 5) ELISA Kit
E-EL-S0020	Sheep TrxR1 (Thioredoxin Reductase 1) ELISA Kit
E-EL-S0021	Sheep AP 13 (Apelin 13) ELISA Kit
E-EL-S0022	Sheep DCX (Doublecortin) ELISA Kit

Cat.NO	Product Name
E-EL-S0023	Sheep CK-HMW (Cytokeratin High Molecular Weight) ELISA Kit
E-EL-S0024	Sheep TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-S0025	Sheep cTn-I/TNNI3 (cardiac Troponin I) ELISA Kit
E-EL-S0026	Sheep COL1 α 2 (Collagen Type I Alpha 2) ELISA Kit
E-EL-S0028	Sheep TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit
E-EL-S0029	Sheep CLA/CTAGE (Cutaneous Lymphocyte Associated Antigen) ELISA Kit
E-EL-S0030	Sheep Preptin ELISA Kit
E-EL-S0031	Sheep TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit
E-EL-S0033	Sheep FABP1 (Fatty Acid Binding Protein 1, Liver) ELISA Kit
E-EL-S0034	Sheep TMS β 4 (Thymosin Beta 4) ELISA Kit
E-EL-S0035	Sheep GATA4 (GATA Binding Protein 4) ELISA Kit
E-EL-S0036	Sheep IL-12 p40 (Interleukin 12 p40) ELISA Kit
E-EL-S0037	Sheep MANF (Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit
E-EL-S0038	Sheep PRELP (Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit
E-EL-S0039	Sheep CLDN14 (Claudin 14) ELISA Kit
E-EL-S0040	Sheep GHAb (Anti-Growth Hormone Antibody) ELISA Kit
E-EL-S0041	Sheep PHB (Prohibitin) ELISA Kit
E-EL-S0043	Sheep GAL7 (Galectin 7) ELISA Kit
E-EL-S0044	Sheep F13A1 (Coagulation Factor XIII A1 Polypeptide) ELISA Kit
E-EL-S0046	Sheep ApoE (Apolipoprotein E) ELISA Kit

Cat.NO	Product Name
E-EL-S0047	Sheep FAK (Focal Adhesion Kinase) ELISA Kit
E-EL-S0048	Sheep ATA (Angiotensinase A) ELISA Kit
E-EL-S0049	Sheep CK-19/KRT19 (Cytokeratin 19) ELISA Kit
E-EL-S0050	Sheep ANG II R2/AGTR2 (Angiotensin II Receptor 2) ELISA Kit
E-EL-S0051	Sheep IL-23 (Interleukin 23) ELISA Kit
E-EL-S0052	Sheep PARP4 (Poly ADP Ribose Polymerase 4) ELISA Kit
E-EL-S0053	Sheep PGAM2 (Phosphoglycerate Mutase 2) ELISA Kit
E-EL-S0054	Sheep SDC1 (Syndecan 1) ELISA Kit
E-EL-S0055	Sheep IL-27 (Interleukin 27) ELISA Kit
E-EL-S0056	Sheep SHBG (Sex Hormone Binding Globulin) ELISA Kit
E-EL-S0057	Sheep ESF (Erythropoiesis Stimulating Factor) ELISA Kit
E-EL-S0058	Sheep C4 (Complement Component 4) ELISA Kit
E-EL-S0059	Sheep ADAM9 (A Disintegrin And Metalloprotease 9) ELISA Kit
E-EL-S0060	Sheep ON (Osteonectin) ELISA Kit
E-EL-S0061	Sheep MBP (Myelin Basic Protein) ELISA Kit
E-EL-S0062	Sheep PRSS12 (Protease, Serine 12) ELISA Kit
E-EL-S0063	Sheep TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit
E-EL-S0064	Sheep TUB α 3C (Tubulin Alpha 3c)ELISA Kit
E-EL-S0065	Sheep CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit
E-EL-S0066	Sheep FTH (Ferritin, Heavy Polypeptide) ELISA Kit

Cat.NO	Product Name
E-EL-S0067	Sheep OTF2A (Octamer Transcription Factor 2A) ELISA Kit
E-EL-S0068	Sheep MUC1 (Mucin 1) ELISA Kit
E-EL-S0069	Sheep ALT (Alanine Transaminase) ELISA Kit
E-EL-S0070	Sheep Sp100 (Sp100 Nuclear Antigen) ELISA Kit
E-EL-S0071	Sheep PKI α (Protein Kinase Inhibitor Alpha) ELISA Kit
E-EL-S0072	Sheep IL-35 (Interleukin 35) ELISA Kit
E-EL-S0074	Sheep UBE2C (Ubiquitin Conjugating Enzyme E2C) ELISA Kit
E-EL-S0075	Sheep GC (Glucagon) ELISA Kit
E-EL-S0077	Sheep AMBP (Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit
E-EL-S0078	Sheep GCK (Glucokinase) ELISA Kit
E-EL-S0079	Sheep COL15 (Collagen Type XV) ELISA Kit
E-EL-S0080	Sheep PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit
E-EL-S0081	Sheep Fc ϵ R I (Receptor I for the Fc Fragment of IgE) ELISA Kit
E-EL-S0082	Sheep PCDH β 16 (Protocadherin Beta 16) ELISA Kit
E-EL-S0084	Sheep TXLNb (Taxilin Beta) ELISA Kit
E-EL-S0085	Sheep bTC (Betacellulin) ELISA Kit
E-EL-S0086	Sheep SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-S0087	Sheep BRAK/CXCL14 (Breast and Kidney Expressed Chemokine) ELISA Kit
E-EL-S0088	Sheep PLAP/ALPP (Placental Alkaline Phosphatase) ELISA Kit
E-EL-S0089	Sheep I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit

Cat.NO	Product Name
E-EL-S0090	Sheep GADL1 (Glutamate Decarboxylase Like Protein 1) ELISA Kit
E-EL-S0091	Sheep PR3Ab (Proteinase 3 Antibody) ELISA Kit
E-EL-S0092	Sheep GP (Anti-Glycoprotein Antibody) ELISA Kit
E-EL-S0093	Sheep STAT2 (Signal Transducer and Activator of Transcription 2) ELISA Kit
E-EL-S0094	Sheep FLN β (Filamin B, Beta) ELISA Kit
E-EL-S0095	Sheep Smad3/MADH3 (Mothers Against Decapentaplegic Homolog 3) ELISA Kit
E-EL-S0096	Sheep SIAE (Sialic Acid Acetyltransferase) ELISA Kit
E-EL-S0098	Sheep BCL2L1/BCL-X (Bcl-2 Like Protein 1) ELISA Kit
E-EL-S0099	Sheep CRN (Corin) ELISA Kit
E-EL-S0100	Sheep TXLNa (Taxilin Alpha) ELISA Kit
E-EL-S0102	Sheep BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit
E-EL-S0103	Sheep ASB (Arylsulfatase B) ELISA Kit
E-EL-S0104	Sheep CLC (Charcot Leyden Crystal Protein) ELISA Kit
E-EL-S0105	Sheep FRS3 (Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit
E-EL-S0106	Sheep TSP-1 (Thrombospondin-1) ELISA Kit
E-EL-S0107	Sheep BACE2 (β -site APP Cleaving Enzyme 2) ELISA Kit
E-EL-S0108	Sheep A β 1-40 (Amyloid Beta Peptide 1-40) ELISA Kit
E-EL-S0109	Sheep Pro-MMP-13 (Pro-Matrix Metalloproteinase 13) ELISA Kit
E-EL-S0110	Sheep VEGFR1/Flt1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit
E-EL-S0111	Sheep PDXK (Pyridoxal Kinase) ELISA Kit

Cat.NO	Product Name
E-EL-S0112	Sheep POMC (Pro-Opiomelanocortin) ELISA Kit
E-EL-S0113	Sheep PAR4 (Protease Activated Receptor 4) ELISA Kit
E-EL-S0114	Sheep GATA3 (GATA Binding Protein 3) ELISA Kit
E-EL-S0115	Sheep Cys-C (Cystatin C) ELISA Kit
E-EL-S0117	Sheep IGFBP-7 (Insulin Like Growth Factor Binding Protein 7) ELISA Kit
E-EL-S0119	Sheep GAL (Galanin) ELISA Kit
E-EL-S0120	Sheep TXB2 (Thromboxane B2) ELISA Kit
E-EL-S0121	Sheep β 2-GP (Beta2 Glycoprotein) ELISA Kit
E-EL-S0122	Sheep AGGF1 (Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit
E-EL-S0123	Sheep PC/CPG (Choline Phosphoglyceride) ELISA Kit
E-EL-S0124	Sheep FG γ (Fibrinogen Gamma) ELISA Kit
E-EL-S0125	Sheep DAO (Diamine Oxidase) ELISA Kit
E-EL-S0127	Sheep α Actin (Alpha-Actin) ELISA Kit
E-EL-S0128	Sheep C3c (Complement C3 Convertase) ELISA Kit
E-EL-S0129	Sheep ELA2B (Elastase 2B) ELISA Kit
E-EL-S0130	Sheep ApoB (Apolipoprotein B) ELISA Kit
E-EL-S0131	Sheep C8g (Complement Component 8g) ELISA Kit
E-EL-S0132	Sheep FGL2 (Fibrinogen Like Protein 2) ELISA Kit
E-EL-S0133	Sheep CGRP (Calcitonin Gene Related Peptide) ELISA Kit
E-EL-S0135	Sheep REG3 γ (Regenerating Islet Derived Protein 3 Gamma) ELISA Kit

Cat.NO	Product Name
E-EL-S0136	Sheep PL (Placental Lactogen) ELISA Kit
E-EL-S0138	Sheep Des (Desmin) ELISA Kit
E-EL-S0139	Sheep CK-LMW (Cytokeratin Low Molecular Weight) ELISA Kit
E-EL-S0140	Sheep CG α (Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit
E-EL-S0141	Sheep Ly6a (Lymphocyte Antigen 6 Complex Locus A) ELISA Kit
E-EL-S0142	Sheep PG (Proteoglycan) ELISA Kit
E-EL-S0143	Sheep FGF4 (Fibroblast Growth Factor 4) ELISA Kit
E-EL-S0145	Sheep AFA (Anti-Filaggrin Antibody) ELISA Kit
E-EL-S0146	Sheep MFI2 (Melanotransferrin) ELISA Kit
E-EL-S0148	Sheep TFF2 (Trefoil Factor 2) ELISA Kit
E-EL-S0149	Sheep α -SMA (α -Smooth muscle actin) ELISA Kit
E-EL-S0150	Sheep NUP155 (Nucleoporin 155kDa) ELISA Kit
E-EL-S0151	Sheep MAN1A1 (Mannosidase Alpha Class 1A Member 1) ELISA Kit
E-EL-S0152	Sheep PDGFR α (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit
E-EL-S0153	Sheep 5-NT (5-Nucleotidase) ELISA Kit
E-EL-S0156	Sheep PICK1 (Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit
E-EL-S0157	Sheep TP63 (Tumor Protein p63) ELISA Kit
E-EL-S0158	Sheep α 1BG (Alpha-1-B-Glycoprotein) ELISA Kit
E-EL-S0159	Sheep GAP43 (Growth Associated Protein 43) ELISA Kit
E-EL-S0160	Sheep SYP(Synaptophysin)ELISA Kit

Cat.NO	Product Name
E-EL-S0161	Sheep CKMB (Creatine Kinase MB Isoenzyme) ELISA Kit
E-EL-S0163	Sheep MCP-3 (Monocyte Chemotactic Protein 3) ELISA Kit
E-EL-S0164	Sheep TAT (Thrombin/Antithrombin Complex) ELISA Kit
E-EL-S0165	Sheep Salusin α (Salusin Alpha) ELISA Kit
E-EL-S0166	Sheep COX-2R (Cyclooxygenase-2 Receptor) ELISA Kit
E-EL-S0167	Sheep TUB β (β -Tubulin) ELISA Kit
E-EL-S0168	Sheep GMP-140 (Alpha-granular Membrane Protein) ELISA Kit
E-EL-S0169	Sheep PKM2 (Pyruvate Kinase, Muscle) ELISA Kit
E-EL-S0170	Sheep RXFP1 (Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit
E-EL-S0171	Sheep CYFRA21-1 (Cytokeratin Fragment Antigen 21-1) ELISA Kit
E-EL-S0172	Sheep GP5 (Glycoprotein V, Platelet) ELISA Kit
E-EL-S0173	Sheep UP (Proteinuria) ELISA Kit
E-EL-S0174	Sheep KIM-1 (Kidney Injury Molecule 1) ELISA Kit
E-EL-S0175	Sheep BMPR2 (Bone Morphogenetic Protein Receptor II) ELISA Kit
E-EL-S0177	Sheep CDC23 (Cell Division Cycle Protein 23) ELISA Kit
E-EL-S0178	Sheep 6-OHDA (6-Hydroxydopamine) ELISA Kit
E-EL-S0179	Sheep LRP6 (Low Density Lipoprotein Receptor Related Protein) ELISA Kit
E-EL-S0180	Sheep MMP-13 (Matrix Metalloproteinase 13) ELISA Kit
E-EL-S0181	Sheep TSTA (Tumor Specific Transplantation Antigen) ELISA Kit
E-EL-S0182	Sheep OSMR (Oncostatin M Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-S0183	Sheep CA1 (Carbonic Anhydrase I) ELISA Kit
E-EL-S0184	Sheep AGPA (Anti-Glucoprotein 210 Antibody) ELISA Kit
E-EL-S0186	Sheep PHSA (Poly Serum Antigen) ELISA Kit
E-EL-S0187	Sheep HSP-40 (Heat Shock Protein 40) ELISA Kit
E-EL-S0188	Sheep DEF β 2 (Defensin Beta 2) ELISA Kit
E-EL-S0189	Sheep STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit
E-EL-S0190	Sheep MMP-9 (Matrix Metalloproteinase 9) ELISA Kit
E-EL-S0191	Sheep LT-E4 (Leukotriene E4) ELISA Kit
E-EL-S0193	Sheep DMD (Dystrophin) ELISA Kit
E-EL-S0194	Sheep LIPF (Lipase, Gastric) ELISA Kit
E-EL-S0195	Sheep Smad1/MADH1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit
E-EL-S0196	Sheep STMN1 (Stathmin 1) ELISA Kit
E-EL-S0197	Sheep MIgD (Membrane IgD) ELISA Kit
E-EL-S0198	Sheep TD-Ag (Thymus Dependent Antigen) ELISA Kit
E-EL-S0199	Sheep IL-8 (Interleukin 8) ELISA Kit
E-EL-S0200	Sheep CD40L/TNFSF5 (Cluster of Differentiation 40 Ligand) ELISA Kit
E-EL-S0201	Sheep EAR1 (Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit
E-EL-S0202	Sheep ANG II (Angiotensin II) ELISA Kit
E-EL-S0203	Sheep ASGR2 (Asialoglycoprotein Receptor 2) ELISA Kit
E-EL-S0204	Sheep MUC21 (Mucin 21) ELISA Kit

Cat.NO	Product Name
E-EL-S0205	Sheep FGFR4 (Fibroblast Growth Factor Receptor 4) ELISA Kit
E-EL-S0206	Sheep MYH10 (Myosin Heavy Chain 10, Non Muscle) ELISA Kit
E-EL-S0207	Sheep FcαR/CD89 (Receptor I for the Fc Fragment of IgA) ELISA Kit
E-EL-S0208	Sheep MEC/CCL28 (Mucosae Associated Epithelia Chemokine) ELISA Kit
E-EL-S0209	Sheep GSTθ2/GSTt2 (Glutathione S Transferase Theta 2) ELISA Kit
E-EL-S0210	Sheep HPR (Haptoglobin Related Protein) ELISA Kit
E-EL-S0211	Sheep CASP11 (Caspase 11) ELISA Kit
E-EL-S0212	Sheep CYP7A1 (Cytochrome P450 7A1) ELISA Kit
E-EL-S0213	Sheep CXCR3 (Chemokine CXC Receptor 3) ELISA Kit
E-EL-S0214	Sheep VLDLR (Very Low Density Lipoprotein Receptor) ELISA Kit
E-EL-S0215	Sheep DUB (Deubiquitinating Enzyme) ELISA Kit
E-EL-S0217	Sheep TNKS1 (Tankyrase 1) ELISA Kit
E-EL-S0218	Sheep IgA2 (Immunoglobulin A2) ELISA Kit
E-EL-S0219	Sheep CHRM2 (Cholinergic Receptor, Muscarinic 2) ELISA Kit
E-EL-S0220	Sheep DAF (Decay Accelerating Factor) ELISA Kit
E-EL-S0221	Sheep RHOA (Ras Homolog Gene Family, Member A) ELISA Kit
E-EL-S0222	Sheep PTF1α (Pancreas Specific Transcription Factor 1α) ELISA Kit
E-EL-S0223	Sheep Cx (Connexin) ELISA Kit
E-EL-S0225	Sheep CP (Ceruloplasmin) ELISA Kit
E-EL-S0227	Sheep GROα/CXCL1/NAP3 (Growth Regulated Oncogene Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-S0228	Sheep CEBPγ (CCAAT/Enhancer Binding Protein Gamma) ELISA Kit
E-EL-S0229	Sheep DPP4 (Dipeptidyl Peptidase IV) ELISA Kit
E-EL-S0230	Sheep HHIP (Hedgehog Interacting Protein) ELISA Kit
E-EL-S0231	Sheep TG (Thyroglobulin) ELISA Kit
E-EL-S0232	Sheep ARHGDIα (Rho GDP Dissociation Inhibitor Alpha) ELISA Kit
E-EL-S0233	Sheep PGC (Pepsinogen C) ELISA Kit
E-EL-S0234	Sheep SDF-4 (Stromal Cell Derived Factor 4) ELISA Kit
E-EL-S0235	Sheep TFR2 (Transferrin Receptor 2) ELISA Kit
E-EL-S0236	Sheep G6PD (Glucose 6 Phosphate Dehydrogenase) ELISA Kit
E-EL-S0237	Sheep KLK7 (Kallikrein 7) ELISA Kit
E-EL-S0238	Sheep TR (Thrombin Receptor) ELISA Kit
E-EL-S0239	Sheep Ras-GAP1 (Rac-GTPase Activating Protein 1) ELISA Kit
E-EL-S0240	Sheep CKLFSF7 (Chemokine Like Factor Superfamily 7) ELISA Kit
E-EL-S0241	Sheep BF (Factor B) ELISA Kit
E-EL-S0243	Sheep RPS6Kβ1 (Ribosomal Protein S6 Kinase Beta 1) ELISA Kit
E-EL-S0244	Sheep PTPRO (Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit
E-EL-S0245	Sheep SPJ (Surfactant Associated Protein J) ELISA Kit
E-EL-S0246	Sheep PROZ (Protein Z) ELISA Kit
E-EL-S0247	Sheep DEFα5 (Defensin Alpha 5, Paneth Cell Specific) ELISA Kit
E-EL-S0248	Sheep ADH (Antidiuretic Hormone) ELISA Kit

Cat.NO	Product Name
E-EL-S0249	Sheep TP (Thymidine Phosphorylase) ELISA Kit
E-EL-S0250	Sheep S100A7 (S100 Calcium Binding Protein A7) ELISA Kit
E-EL-S0251	Sheep TPP1 (Tripeptidyl Peptidase I) ELISA Kit
E-EL-S0253	Sheep TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
E-EL-S0254	Sheep PDGF-AB (Platelet Derived Growth Factor AB) ELISA Kit
E-EL-S0255	Sheep HCCS (Holocytochrome C Synthase) ELISA Kit
E-EL-S0256	Sheep COL14 (Collagen Type XIV) ELISA Kit
E-EL-S0257	Sheep β TG/PBP/CXCL7/NAP2 (Thromboglobulin, Beta) ELISA Kit
E-EL-S0258	Sheep HO1 (Heme Oxygenase 1) ELISA Kit
E-EL-S0260	Sheep GAL14 (Galectin 14) ELISA Kit
E-EL-S0261	Sheep FGF-7/KGF (Fibroblast Growth Factor 7) ELISA Kit
E-EL-S0262	Sheep CCL14 (Chemokine C-C-Motif Ligand 14) ELISA Kit
E-EL-S0263	Sheep CA3 (Carbonic Anhydrase III, Muscle Specific) ELISA Kit
E-EL-S0265	Sheep CPE (Carboxypeptidase E) ELISA Kit
E-EL-S0266	Sheep NF- κ B p65 (Nuclear Factor- κ B p65) ELISA Kit
E-EL-S0267	Sheep NFYC (Nuclear Transcription Factor Y Gamma) ELISA Kit
E-EL-S0268	Sheep RAR α (Retinoic Acid Receptor Alpha) ELISA Kit
E-EL-S0269	Sheep C1qR1 (Complement Component 1q Receptor) ELISA Kit
E-EL-S0271	Sheep bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit
E-EL-S0272	Sheep MR (Mannose Receptor) ELISA Kit

Cat.NO	Product Name
E-EL-S0274	Sheep FRS2 (Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit
E-EL-S0276	Sheep XCR1 (XC-Chemokine Receptor 1) ELISA Kit
E-EL-S0277	Sheep CTX-III (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit
E-EL-S0278	Sheep IGFBP-1 (Insulin Like Growth Factor Binding Protein 1) ELISA Kit
E-EL-S0279	Sheep DEFα1 (Defensin Alpha 1, Neutrophil) ELISA Kit
E-EL-S0280	Sheep PDGFA (Platelet Derived Growth Factor Subunit A) ELISA Kit
E-EL-S0281	Sheep COL10 (Collagen Type X) ELISA Kit
E-EL-S0283	Sheep LTA (Lymphotoxin α)ELISA Kit
E-EL-S0284	Sheep CKMT1B (Creatine Kinase, Mitochondrial 1B) ELISA Kit
E-EL-S0285	Sheep MYD88 (Myeloid Differentiation Factor 88) ELISA Kit
E-EL-S0287	Sheep PKBG (Protein Kinase B Gamma) ELISA Kit
E-EL-S0288	Sheep AIF (Apoptosis Inducing Factor) ELISA Kit
E-EL-S0289	Sheep CYTH1 (Cytohesin 1) ELISA Kit
E-EL-S0290	Sheep BTA (Bladder Tumor Antigen) ELISA Kit
E-EL-S0291	Sheep C3a (Complement Component 3a) ELISA Kit
E-EL-S0292	Sheep AQP-3 (Aquaporin 3) ELISA Kit
E-EL-S0293	Sheep OSCAR (Osteoclast Associated Receptor) ELISA Kit
E-EL-S0294	Sheep COL4 (Collagen Type IV) ELISA Kit
E-EL-S0296	Sheep tPA (Plasminogen Activator, Tissue) ELISA Kit
E-EL-S0297	Sheep PDIA6 (Protein Disulfide Isomerase A6) ELISA Kit

Cat.NO	Product Name
E-EL-S0298	Sheep CTR (Calcitonin Receptor) ELISA Kit
E-EL-S0299	Sheep ACAT2 (Acetyl Coenzyme A Acetyltransferase 2) ELISA Kit
E-EL-S0300	Sheep F7 (Coagulation Factor VII) ELISA Kit
E-EL-S0301	Sheep AZU1 (Azurocidin 1) ELISA Kit
E-EL-S0302	Sheep OTR (Oxytocin Receptor) ELISA Kit
E-EL-S0303	Sheep CA12 (Carbonic Anhydrase XII) ELISA Kit
E-EL-S0304	Sheep NUMA1 (Nuclear Mitotic Apparatus Protein 1) ELISA Kit
E-EL-S0305	Sheep PMN Elastase (Polymorphonuclear Elastase) ELISA Kit
E-EL-S0306	Sheep MUC5AC (Mucin 5 Subtype AC) ELISA Kit
E-EL-S0307	Sheep CX31 (Connexin 31) ELISA Kit
E-EL-S0308	Sheep NEFL (Neurofilament, Light Polypeptide) ELISA Kit
E-EL-S0309	Sheep GH2 (Growth Hormone 2) ELISA Kit
E-EL-S0310	Sheep MYL1 (Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit
E-EL-S0311	Sheep DPYD (Dihydropyrimidine Dehydrogenase) ELISA Kit
E-EL-S0312	Sheep APLN (Apelin) ELISA Kit
E-EL-S0313	Sheep ApoF (Apolipoprotein F) ELISA Kit
E-EL-S0314	Sheep GP4 (Platelet Membrane Glycoprotein IV) ELISA Kit
E-EL-S0315	Sheep FOXP1 (Forkhead Box Protein P1) ELISA Kit
E-EL-S0316	Sheep ECM1 (Extracellular Matrix Protein 1) ELISA Kit
E-EL-S0317	Sheep PAR2 (Protease Activated Receptor 2) ELISA Kit

Cat.NO	Product Name
E-EL-S0318	Sheep HSP-70/HSPA9 (Heat Shock Protein 70) ELISA Kit
E-EL-S0319	Sheep UCP-2 (Uncoupling Protein 2, Mitochondrial) ELISA Kit
E-EL-S0320	Sheep CETP (Cholesteryl Ester Transfer Protein) ELISA Kit
E-EL-S0321	Sheep GzmB (Granzyme B) ELISA Kit
E-EL-S0322	Sheep CC16 (Clara Cell 16KD Protein) ELISA Kit
E-EL-S0324	Sheep PGD2-MOX (Prostaglandin D2-Methoxime) ELISA Kit
E-EL-S0325	Sheep AGT (Angiotensinogen) ELISA Kit
E-EL-S0326	Sheep TSA (Thyroid Stimulating Antibody) ELISA Kit
E-EL-S0327	Sheep Alk-Smase (Alkaline Sphingomyelinase) ELISA Kit
E-EL-S0328	Sheep FABP5 (Fatty Acid Binding Protein 5, Epidermal) ELISA Kit
E-EL-S0329	Sheep LTC4 (Leukotriene C4) ELISA Kit
E-EL-S0330	Sheep ENG/CD105 (Endoglin) ELISA Kit
E-EL-S0331	Sheep PKC ζ (Protein Kinase C Zeta) ELISA Kit
E-EL-S0332	Sheep CDX2 (Caudal Type Homeobox Transcription Factor 2) ELISA Kit
E-EL-S0334	Sheep FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit
E-EL-S0336	Sheep β MSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-S0338	Sheep LCN2 (Lipocalin 2) ELISA Kit
E-EL-S0339	Sheep SCFR/c-Kit (Stem Cell Factor Receptor) ELISA Kit
E-EL-S0340	Sheep ANGPTL3 (Angiopoietin Like Protein 3) ELISA Kit
E-EL-S0341	Sheep CYB561 (Cytochrome b-561) ELISA Kit

Cat.NO	Product Name
E-EL-S0342	Sheep PIK3C β (Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit
E-EL-S0343	Sheep CD3d (Cluster of Differentiation 3d) ELISA Kit
E-EL-S0344	Sheep TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit
E-EL-S0345	Sheep COL17 (Collagen Type XVII) ELISA Kit
E-EL-S0346	Sheep UBE3A (Ubiquitin Protein Ligase E3A) ELISA Kit
E-EL-S0347	Sheep FBN2 (Fibrillin 2) ELISA Kit
E-EL-S0348	Sheep cPLA2 (Cytosolic Phospholipase A2) ELISA Kit
E-EL-S0349	Sheep Gal-6S (Galactose-6-Sulphatase) ELISA Kit
E-EL-S0350	Sheep VIP (Vasoactive Intestinal Peptide) ELISA Kit
E-EL-S0351	Sheep PLA2R1 (Phospholipase A2 Receptor 1) ELISA Kit
E-EL-S0353	Sheep IgG1 (Immunoglobulin G1) ELISA Kit
E-EL-S0354	Sheep CYP1A2 (Cytochrome P450 1A2) ELISA Kit
E-EL-S0356	Sheep ApoC1 (Apolipoprotein C1) ELISA Kit
E-EL-S0357	Sheep Sp17 (Sperm Protein 17) ELISA Kit
E-EL-S0359	Sheep γ -ECS (γ -Glutamyl Systeine Synthetase) ELISA Kit
E-EL-S0360	Sheep Surv (Survivin) ELISA Kit
E-EL-S0361	Sheep NLR (NOD-Like Receptor9) ELISA Kit
E-EL-S0362	Sheep SA (Sialic Acid) ELISA Kit
E-EL-S0363	Sheep M-AChR M2 (Muscarinic Acetylcholine Receptor M2) ELISA Kit
E-EL-S0364	Sheep NT-4 (Neurotrophin 4) ELISA Kit

Cat.NO	Product Name
E-EL-S0365	Sheep IGF-1 (Insulin Like Growth Factor 1) ELISA Kit
E-EL-S0366	Sheep Thyronine ELISA Kit
E-EL-S0367	Sheep LTβR (Lymphotoxin Beta Receptor) ELISA Kit
E-EL-S0368	Sheep PRDX2 (Peroxisredoxin 2) ELISA Kit
E-EL-S0369	Sheep ATXN1 (Ataxin 1) ELISA Kit
E-EL-S0370	Sheep SDC4 (Syndecan 4) ELISA Kit
E-EL-S0371	Sheep MSTN (Myostatin) ELISA Kit
E-EL-S0372	Sheep COL11α1 (Collagen Type XI Alpha 1) ELISA Kit
E-EL-S0374	Sheep FRA (Fibronectin Related Antigen) ELISA Kit
E-EL-S0375	Sheep CALLA (Common Acute Lymphocytic Leukaemia Antigen) ELISA Kit
E-EL-S0376	Sheep GFAP (Glial Fibrillary Acidic Protein) ELISA Kit
E-EL-S0377	Sheep Grb2 (Growth Factor Receptor Bound Protein 2) ELISA Kit
E-EL-S0378	Sheep Pro-MMP-1 (pro-Matrix Metalloproteinase 1) ELISA Kit
E-EL-S0380	Sheep PTGIS (Prostaglandin I Synthase) ELISA Kit
E-EL-S0381	Sheep GPX3 (Glutathione Peroxidase 3, Plasma) ELISA Kit
E-EL-S0382	Sheep SNRPD1 (Small Nuclear Ribonucleoprotein Polypeptide D1) ELISA Kit
E-EL-S0383	Sheep ST2 (Syntenin 2) ELISA Kit
E-EL-S0384	Sheep SPC (Pulmonary Surfactant Associated Protein C) ELISA Kit
E-EL-S0385	Sheep GDF1 (Growth Differentiation Factor 1) ELISA Kit
E-EL-S0386	Sheep PLCγ1 (Phospholipase C Gamma 1) ELISA Kit

Cat.NO	Product Name
E-EL-S0387	Sheep CAPF (Cuticular Active Peptide Factor) ELISA Kit
E-EL-S0388	Sheep TSC22 (Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit
E-EL-S0389	Sheep CLDN4 (Claudin 4) ELISA Kit
E-EL-S0390	Sheep FGβ (Fibrinogen Beta) ELISA Kit
E-EL-S0391	Sheep MRP (Multidrug Resistance-Associated Protein) ELISA Kit
E-EL-S0392	Sheep TI-Ag (Thymus Independent Antigen) ELISA Kit
E-EL-S0393	Sheep RAG-1 (Recombination Activating Gene 1) ELISA Kit
E-EL-S0394	Sheep ANG III (Angiotensin III) ELISA Kit
E-EL-S0395	Sheep GSTθ1/GSTt1 (Glutathione S Transferase Theta 1) ELISA Kit
E-EL-S0396	Sheep EAR12 (Eosinophil Associated Ribonuclease A family Member 12) ELISA Kit
E-EL-S0397	Sheep ATF6 (Activating Transcription Factor 6) ELISA Kit
E-EL-S0398	Sheep CDKN3 (Cyclin Dependent Kinase Inhibitor 3) ELISA Kit
E-EL-S0399	Sheep C4a (Complement Component 4a) ELISA Kit
E-EL-S0400	Sheep HNRPA2B1 (Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit
E-EL-S0401	Sheep AATF (Apoptosis Antagonizing Transcription Factor) ELISA Kit
E-EL-S0402	Sheep MMP-25 (Matrix Metalloproteinase 25) ELISA Kit
E-EL-S0403	Sheep LCA/CD45(Leukocyte Common Antigen) ELISA Kit
E-EL-S0404	Sheep CDK1 (Cyclin Dependent Kinase 1) ELISA Kit
E-EL-S0405	Sheep PFKP (Phosphofructokinase, Platelet) ELISA Kit
E-EL-S0406	Sheep KERA (Keratocan) ELISA Kit

Cat.NO	Product Name
E-EL-S0408	Sheep MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
E-EL-S0409	Sheep AQP-1 (Aquaporin 1) ELISA Kit
E-EL-S0410	Sheep NP (Nandrolone Phenylpropionate) ELISA Kit
E-EL-S0411	Sheep AMA IgA (Anti-Myelin antibody IgA) ELISA Kit
E-EL-S0412	Sheep CNTN3 (Contactin 3) ELISA Kit
E-EL-S0413	Sheep CNTN4 (Contactin 4) ELISA Kit
E-EL-S0414	Sheep LZM (Lysozyme) ELISA Kit
E-EL-S0415	Sheep FBLN5 (Fibulin 5) ELISA Kit
E-EL-S0416	Sheep CD4 (Cluster of Differentiation 4) ELISA Kit
E-EL-S0417	Sheep AAA (Anti-Albumin Antibody) ELISA Kit
E-EL-S0418	Sheep ACT α 2 (Actin Alpha 2, Smooth Muscle) ELISA Kit
E-EL-S0419	Sheep MN (Metanephrine) ELISA Kit
E-EL-S0420	Sheep ET (Epitestosterone) ELISA Kit
E-EL-S0421	Sheep GDF11 (Growth Differentiation Factor 11) ELISA Kit
E-EL-S0422	Sheep GR α (Glucocorticoid Receptor Alpha) ELISA Kit
E-EL-S0423	Sheep PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit
E-EL-S0424	Sheep D4R/DRD4 (Dopamine Receptor D4) ELISA Kit
E-EL-S0425	Sheep CBX3 (Chromobox Homolog 3) ELISA Kit
E-EL-S0426	Sheep GL β (Galactosidase Beta) ELISA Kit
E-EL-S0427	Sheep TLR9 (Toll-Like Receptor 9) ELISA Kit

Cat.NO	Product Name
E-EL-S0428	Sheep Methylase ELISA Kit
E-EL-S0429	Sheep β -LPH (β -Lipotropic Hormone) ELISA Kit
E-EL-S0431	Sheep BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit
E-EL-S0432	Sheep NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit
E-EL-S0433	Sheep PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit
E-EL-S0434	Sheep SPARCL1 (SPARC Like Protein 1) ELISA Kit
E-EL-S0435	Sheep CKLFSF8 (Chemokine Like Factor Superfamily 8) ELISA Kit
E-EL-S0436	Sheep AGRN (Agrin) ELISA Kit
E-EL-S0437	Sheep TFPI (Tissue Factor Pathway Inhibitor) ELISA Kit
E-EL-S0438	Sheep ANEA (Anti-Neutrophil Elastase Antibody) ELISA Kit
E-EL-S0439	Sheep ApoH (Apolipoprotein H) ELISA Kit
E-EL-S0440	Sheep eIF3a (Eukaryotic Translation Initiation Factor 3a) ELISA Kit
E-EL-S0441	Sheep C7 (Complement Component 7) ELISA Kit
E-EL-S0442	Sheep MUC5B (Mucin 5 Subtype B) ELISA Kit
E-EL-S0443	Sheep OVA sIgG (Ovalbumin Specific IgG) ELISA Kit
E-EL-S0445	Sheep MAP2K1 (Mitogen Activated Protein Kinase Kinase 1) ELISA Kit
E-EL-S0446	Sheep TAP (Transporter Associated with Antigen Processing) ELISA Kit
E-EL-S0447	Sheep CASP7 (Caspase 7) ELISA Kit
E-EL-S0448	Sheep AQP-4 (Aquaporin 4) ELISA Kit
E-EL-S0449	Sheep GALR1 (Galanin Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-S0450	Sheep ApoA1 (Apolipoprotein A1) ELISA Kit
E-EL-S0452	Sheep CAPA (Cytotoxin Associated Protein A) ELISA Kit
E-EL-S0453	Sheep PRKA α 1 (Protein Kinase, AMP Activated Alpha 1) ELISA Kit
E-EL-S0455	Sheep ProGRP (Pro-Gastrin Releasing Peptide) ELISA Kit
E-EL-S0456	Sheep PMP3 (Peroxisomal Membrane Protein 3) ELISA Kit
E-EL-S0457	Sheep CYSLTR2 (Cysteinyl Leukotriene Receptor 2) ELISA Kit
E-EL-S0458	Sheep RB1 (Retinoblastoma Protein 1) ELISA Kit
E-EL-S0459	Sheep CTSL (Cathepsin L) ELISA Kit
E-EL-S0460	Sheep ER β (Estrogen Receptor Beta) ELISA Kit
E-EL-S0461	Sheep DHEA-S (Dehydroepiandrosterone Sulfate) ELISA Kit
E-EL-S0462	Sheep MUC4 (Mucin 4) ELISA Kit
E-EL-S0463	Sheep Pax (Paxillin) ELISA Kit
E-EL-S0464	Sheep TSI (Thyroid Stimulating Immunoglobulin) ELISA Kit
E-EL-S0465	Sheep Fc γ BP (Fc Fragment of IgG Binding Protein) ELISA Kit
E-EL-S0466	Sheep GPC4 (Glypican 4) ELISA Kit
E-EL-S0467	Sheep IGF2R (Insulin Like Growth Factor 2 Receptor) ELISA Kit
E-EL-S0468	Sheep GKRP (Glucokinase Regulatory Protein) ELISA Kit
E-EL-S0469	Sheep NT (Neurotensin) ELISA Kit
E-EL-S0470	Sheep PKC γ (Protein Kinase C Gamma) ELISA Kit
E-EL-S0471	Sheep β -actin (Actin Beta) ELISA Kit

Cat.NO	Product Name
E-EL-S0472	Sheep SBEM (Small Breast Epithelial Mucin) ELISA Kit
E-EL-S0473	Sheep VEGF-C (Vascular Endothelial Growth Factor C) ELISA Kit
E-EL-S0475	Sheep CS (Chondroitin Sulfate) ELISA Kit
E-EL-S0476	Sheep Flt3 (FMS Like Tyrosine Kinase 3) ELISA Kit
E-EL-S0477	Sheep PIINP (Procollagen II N-Terminal ProPeptide) ELISA Kit
E-EL-S0478	Sheep LEPR (Leptin Receptor) ELISA Kit
E-EL-S0479	Sheep PLA1 (Phospholipase A1) ELISA Kit
E-EL-S0480	Sheep CX43 (Connexin 43) ELISA Kit
E-EL-S0481	Sheep c-fos ELISA Kit
E-EL-S0483	Sheep UBAP2 (Ubiquitin Associated Protein 2) ELISA Kit
E-EL-S0485	Sheep PEX2 (Peroxisomal Biogenesis Factor 2) ELISA Kit
E-EL-S0486	Sheep GP6 (Glycoprotein VI, Platelet) ELISA Kit
E-EL-S0487	Sheep CTNNβ1 (Catenin, Beta 1) ELISA Kit
E-EL-S0488	Sheep CSF3/GCSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit
E-EL-S0489	Sheep SIGLEC8 (Sialic Acid Binding Ig Like Lectin 8) ELISA Kit
E-EL-S0490	Sheep IL-3 (Interleukin 3) ELISA Kit
E-EL-S0491	Sheep FADS2 (Fatty Acid Desaturase 2) ELISA Kit
E-EL-S0492	Sheep MSR1 (Macrophage Scavenger Receptor 1) ELISA Kit
E-EL-S0493	Sheep IFN-α (Interferon Alpha) ELISA Kit
E-EL-S0494	Sheep AMY1 (Amylase Alpha 1, Salivary) ELISA Kit

Cat.NO	Product Name
E-EL-S0495	Sheep IL-13 (Interleukin 13) ELISA Kit
E-EL-S0496	Sheep RNASEN (Ribonuclease III, Nuclear) ELISA Kit
E-EL-S0497	Sheep ANG II R1-Ab (Angiotensin II Receptor 1 Antibody) ELISA Kit
E-EL-S0498	Sheep KNG1 (Kininogen 1) ELISA Kit
E-EL-S0499	Sheep DAP6 (Death Associated Protein 6) ELISA Kit
E-EL-S0500	Sheep CT-1 (Cardiotrophin 1) ELISA Kit
E-EL-S0501	Sheep ABCG1 (ATP Binding Cassette Transporter G1) ELISA Kit
E-EL-S0502	Sheep BMF (Bcl-2 Modifying Factor) ELISA Kit
E-EL-S0505	Sheep ASM (Acid Sphingomyelinase) ELISA Kit
E-EL-S0506	Sheep Ext-N (N-Terminal Extein) ELISA Kit
E-EL-S0507	Sheep FcεR II (Receptor II for the Fc Fragment of IgE) ELISA Kit
E-EL-S0508	Sheep TREM-1 (Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit
E-EL-S0510	Sheep PXR (Pregnane X Receptor) ELISA Kit
E-EL-S0511	Sheep SREBP (Sterol Regulatory Element Binding Proteins) ELISA Kit
E-EL-S0514	Sheep SST (Somatostatin) ELISA Kit
E-EL-S0515	Sheep IL-2Rβ (Interleukin 2 Receptor Beta) ELISA Kit
E-EL-S0516	Sheep PTPRK (Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit
E-EL-S0517	Sheep CCKBR (Cholecystokinin B Receptor) ELISA Kit
E-EL-S0518	Sheep CKB (Creatine Kinase, Brain) ELISA Kit
E-EL-S0519	Sheep iPLA2 (Phospholipidase A2, Calcium Independent) ELISA Kit

Cat.NO	Product Name
E-EL-S0520	Sheep dDAVP (1-Desamino-8D-Arginine Vasopressin) ELISA Kit
E-EL-S0521	Sheep MEPE (Matrix Extracellular Phosphoglycoprotein) ELISA Kit
E-EL-S0522	Sheep CK-17/KRT17 (Cytokeratin 17) ELISA Kit
E-EL-S0523	Sheep CD97 (Cluster of Differentiation 97) ELISA Kit
E-EL-S0524	Sheep RAPGEF1 (Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit
E-EL-S0525	Sheep α ZGP1 (Alpha-2-Glycoprotein 1, Zine Binding) ELISA Kit
E-EL-S0526	Sheep TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit
E-EL-S0527	Sheep PUMA (P53 Upregulated Modulator of Apoptosis) ELISA Kit
E-EL-S0528	Sheep CYP24A1 (Cytochrome P450 24A1) ELISA Kit
E-EL-S0529	Sheep EIPS2 (Estrogen Induced Protein pS2) ELISA Kit
E-EL-S0530	Sheep Fcy (Fc Fragment of IgG) ELISA Kit
E-EL-S0531	Sheep ECFR (Eosinophil Chemotactic Factor Receptor) ELISA Kit
E-EL-S0532	Sheep CNR1 (Cannabinoid Receptor 1, Brain) ELISA Kit
E-EL-S0533	Sheep ANGPTL6 (Angiopoietin Like Protein 6) ELISA Kit
E-EL-S0535	Sheep OxLDL (Oxidized Low Density Lipoprotein) ELISA Kit
E-EL-S0536	Sheep Bid (BH3 Interacting Domain Death Agonist) ELISA Kit
E-EL-S0537	Sheep OC/BGP (Osteocalcin/Bone Gla Protein) ELISA Kit
E-EL-S0538	Sheep ITIH4 (Inter Alpha-Globulin Inhibitor H4) ELISA Kit
E-EL-S0539	Sheep MT1M (Metallothionein 1M) ELISA Kit
E-EL-S0540	Sheep FGF21(Fibroblast Growth Factor 21) ELISA Kit

Cat.NO	Product Name
E-EL-S0541	Sheep CLS (Cardiolipin Synthase) ELISA Kit
E-EL-S0543	Sheep HTRA1 (HtrA Serine, Peptidase 1) ELISA Kit
E-EL-S0544	Sheep ANGPT4 (Angiopoietin 4) ELISA Kit
E-EL-S0546	Sheep ADIPOR1 (Adiponectin Receptor 1) ELISA Kit
E-EL-S0547	Sheep FOXO4 (Forkhead Box Protein O4) ELISA Kit
E-EL-S0548	Sheep CNTNAP1 (Contactin Associated Protein 1) ELISA Kit
E-EL-S0549	Sheep CKMT1A (Creatine Kinase, Mitochondrial 1A) ELISA Kit
E-EL-S0550	Sheep MAF (Macrophage Activating Factor) ELISA Kit
E-EL-S0551	Sheep DNM2 (Dynamin 2) ELISA Kit
E-EL-S0552	Sheep CVL/EZR (Cytovillin/Ezrin) ELISA Kit
E-EL-S0554	Sheep C5 (Complement Component 5) ELISA Kit
E-EL-S0555	Sheep Tn-T (Troponin T) ELISA Kit
E-EL-S0556	Sheep AGR2 (Anterior Gradient Protein 2) ELISA Kit
E-EL-S0557	Sheep FOXO3 (Forkhead Box Protein O3) ELISA Kit
E-EL-S0558	Sheep NEF3 (Neurofilament 3) ELISA Kit
E-EL-S0559	Sheep gp130 (Glucoprotein 130) ELISA Kit
E-EL-S0561	Sheep GPRC5A (G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit
E-EL-S0562	Sheep AGC (Aggrecan) ELISA Kit
E-EL-S0563	Sheep TSP1 (Thrombin Sensitive Protein 1) ELISA Kit
E-EL-S0564	Sheep PCSK9 (Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit

Cat.NO	Product Name
E-EL-S0565	Sheep TP5 (Thymopentin) ELISA Kit
E-EL-S0566	Sheep MMP-10 (Matrix Metalloproteinase 10) ELISA Kit
E-EL-S0567	Sheep PYGL (Glycogen Phosphorylase, Liver) ELISA Kit
E-EL-S0568	Sheep β -Lg (β -Lactoglobulin) ELISA Kit
E-EL-S0569	Sheep CDK5 (Cyclin Dependent Kinase 5) ELISA Kit
E-EL-S0570	Sheep Mly/CXCL9 (Monocyte Interferon Gamma Inducing Factor) ELISA Kit
E-EL-S0571	Sheep GHITM (Growth Hormone Inducible Transmembrane Protein) ELISA Kit
E-EL-S0572	Sheep FAS/CD95/TNFRSF6 (Factor Related Apoptosis) ELISA Kit
E-EL-S0573	Sheep ANG I (Angiotensin I) ELISA Kit
E-EL-S0577	Sheep TF (Tissue Factor) ELISA Kit
E-EL-S0578	Sheep SFRP2 (Secreted Frizzled Related Protein 2) ELISA Kit
E-EL-S0579	Sheep GATA2 (GATA Binding Protein 2) ELISA Kit
E-EL-S0580	Sheep MGST1 (Microsomal Glutathione S Transferase 1) ELISA Kit
E-EL-S0581	Sheep ACTg1 (Actin Gamma 1) ELISA Kit
E-EL-S0582	Sheep CT (Calcitonin) ELISA Kit
E-EL-S0583	Sheep VDR (Vitamin D Receptor) ELISA Kit
E-EL-S0584	Sheep LAT (Linker For Activation of T-cell) ELISA Kit
E-EL-S0585	Sheep PCDH1 (Protocadherin 1) ELISA Kit
E-EL-S0586	Sheep NCAD (Neural Cadherin) ELISA Kit
E-EL-S0587	Sheep TPO-Ab (Anti-Thyroid-Peroxidase Antibody) ELISA Kit

Cat.NO	Product Name
E-EL-S0588	Sheep F-TESTO (Free Testosterone) ELISA Kit
E-EL-S0589	Sheep PP13 (Placental Protein13) ELISA Kit
E-EL-S0590	Sheep PAEP (Progesterone-Associated Endometrial Protein) ELISA Kit
E-EL-S0591	Sheep ACHAP (Acetylcholinesterase Associated Protein) ELISA Kit
E-EL-S0592	Sheep COL6α3 (Collagen Type VI Alpha 3) ELISA Kit
E-EL-S0594	Sheep MT (Melatonin) ELISA Kit
E-EL-S0595	Sheep SDF-2 (Stromal Cell Derived Factor 2) ELISA Kit
E-EL-S0596	Sheep ECE1 (Endothelin Converting Enzyme 1) ELISA Kit
E-EL-S0597	Sheep GHR (Growth Hormone Receptor) ELISA Kit
E-EL-S0598	Sheep Bcl-10 (B-cell Leukemia/Lymphoma 10) ELISA Kit
E-EL-S0599	Sheep FGF23 (Fibroblast Growth Factor 23) ELISA Kit
E-EL-S0601	Sheep GROβ (Growth Regulated Oncogene Beta) ELISA Kit
E-EL-S0602	Sheep NAGase (N-Acetyl Beta-D-Glucosaminidase) ELISA Kit
E-EL-S0603	Sheep MK (Midkine) ELISA Kit
E-EL-S0605	Sheep β-EPR (Beta-Endorphin Receptor) ELISA Kit
E-EL-S0606	Sheep PDE4D (Phosphodiesterase 4D, cAMP Specific) ELISA Kit
E-EL-S0607	Sheep CAM (Calmodulin) ELISA Kit
E-EL-S0608	Sheep EPCR (Endothelial Protein C Receptor) ELISA Kit
E-EL-S0609	Sheep ELN (Elastin) ELISA Kit
E-EL-S0610	Sheep COL7 (Collagen Type VII) ELISA Kit

Cat.NO	Product Name
E-EL-S0612	Sheep NT-ProANP (N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit
E-EL-S0613	Sheep NRIP1 (Nuclear Receptor Interacting Protein 1) ELISA Kit
E-EL-S0614	Sheep UBC (Urinary Bladder Cancer Antigen) ELISA Kit
E-EL-S0616	Sheep TFPI2 (Tissue Factor Pathway Inhibitor 2) ELISA Kit
E-EL-S0617	Sheep DPYSL2 (Dihydropyrimidinase Like Protein 2) ELISA Kit
E-EL-S0618	Sheep FcγR3B (Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit
E-EL-S0619	Sheep FcγR II/CD32 (Receptor II for the Fc Fragment of IgG) ELISA Kit
E-EL-S0620	Sheep PAPPA2 (Pappalysin 2) ELISA Kit
E-EL-S0621	Sheep SDH (Sorbitol Dehydrogenase) ELISA Kit
E-EL-S0622	Sheep SIGLEC10 (Sialic Acid Binding Ig Like Lectin 10) ELISA Kit
E-EL-S0624	Sheep CXCR1 (CXC-Chemokine Receptor 1) ELISA Kit
E-EL-S0625	Sheep VWF (Von Willebrand Factor) ELISA Kit
E-EL-S0626	Sheep BAFFR (B Cell Activating Factor Receptor) ELISA Kit
E-EL-S0627	Sheep HNRPU (Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit
E-EL-S0628	Sheep ALDH (Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit
E-EL-S0629	Sheep FCN1 (Ficolin 1) ELISA Kit
E-EL-S0630	Sheep CLDN1 (Claudin 1) ELISA Kit
E-EL-S0631	Sheep PAP (Plasmin-Antiplasmin Complex) ELISA Kit
E-EL-S0633	Sheep CPR (Cytochrome P450 Reductase) ELISA Kit
E-EL-S0634	Sheep Bcl-3 (B-cell Leukemia/Lymphoma 3) ELISA Kit

Cat.NO	Product Name
E-EL-S0635	Sheep BALP (Bone Alkaline Phosphatase) ELISA Kit
E-EL-S0636	Sheep GAD-Ab (Anti-Glutamic Acid Decarboxylase antibody) ELISA Kit
E-EL-S0637	Sheep PL (Pancreatic Lipase) ELISA Kit
E-EL-S0638	Sheep CTNN α 1 (Catenin, Alpha 1) ELISA Kit
E-EL-S0639	Sheep c-sis ELISA Kit
E-EL-S0640	Sheep DPT (Dermatopontin) ELISA Kit
E-EL-S0641	Sheep CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) ELISA Kit
E-EL-S0642	Sheep MIP-5 (Macrophage Inflammatory Protein 5) ELISA Kit
E-EL-S0643	Sheep STAT6 (Signal Transducer and Activator of Transcription 6) ELISA Kit
E-EL-S0644	Sheep ETF β (Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit
E-EL-S0645	Sheep REG3 α (Regenerating Islet Derived Protein 3 Alpha) ELISA Kit
E-EL-S0646	Sheep α FP (Alpha-Fetoprotein) ELISA Kit
E-EL-S0647	Sheep TTN (Titin) ELISA Kit
E-EL-S0649	Sheep GDF9 (Growth Differentiation Factor 9) ELISA Kit
E-EL-S0650	Sheep CA153 (Carbohydrate Antigen 153) ELISA Kit
E-EL-S0651	Sheep PA-IgM (Platelet Associated Antibody) ELISA Kit
E-EL-S0652	Sheep ADA (Adenosine Deaminase) ELISA Kit
E-EL-S0653	Sheep SNC α (Synuclein Alpha) ELISA Kit
E-EL-S0654	Sheep EPOR (Erythropoietin Receptor) ELISA Kit
E-EL-S0656	Sheep ANG II R1/AGTR1 (Angiotensin II Receptor 1) ELISA Kit

Cat.NO	Product Name
E-EL-S0657	Sheep CCL16 (Chemokine C-C-Motif Ligand 16) ELISA Kit
E-EL-S0659	Sheep THRSP (Thyroid Hormone Responsive) ELISA Kit
E-EL-S0660	Sheep CHAc (Choline acetylase) ELISA Kit
E-EL-S0661	Sheep GSP (Glycosylated Serum Protein) ELISA Kit
E-EL-S0662	Sheep F8 (Coagulation Factor VIII) ELISA Kit
E-EL-S0663	Sheep EK (Eepidermal Keratin) ELISA Kit
E-EL-S0664	Sheep PGD2S (Prostaglandin D2 Synthase, Brain) ELISA Kit
E-EL-S0665	Sheep AGA (Anti-Glycolipid Antibody) ELISA Kit
E-EL-S0668	Sheep GA (Glycated Albumin) ELISA Kit
E-EL-S0669	Sheep PKP2 (Plakophilin 2) ELISA Kit
E-EL-S0670	Sheep ATF5 (Activating Transcription Factor 5) ELISA Kit
E-EL-S0671	Sheep TPPP (Tubulin Polymerization Promoting Protein) ELISA Kit
E-EL-S0674	Sheep FSH (Follicle-Stimulating Hormone) ELISA Kit
E-EL-S0675	Sheep TGFβI (Transforming Growth Factor Beta Induced Protein) ELISA Kit
E-EL-S0677	Sheep CD30 (Cluster of Differentiation 30) ELISA Kit
E-EL-S0678	Sheep GPR120 (G Protein Coupled Receptor 120) ELISA Kit
E-EL-S0679	Sheep FcaR (Fc Fragment of IgA Receptor) ELISA Kit
E-EL-S0680	Sheep HRG (Histidine Rich Glycoprotein) ELISA Kit
E-EL-S0681	Sheep Ntn4 (Netrin 4) ELISA Kit
E-EL-S0682	Sheep SPOP (Speckle Type POZ Protein) ELISA Kit

Cat.NO	Product Name
E-EL-S0683	Sheep ES (Endostatin) ELISA Kit
E-EL-S0684	Sheep GPER1 (G Protein Coupled Estrogen Receptor 1) ELISA Kit
E-EL-S0685	Sheep β Manase ELISA Kit
E-EL-S0686	Sheep PLP (Parathyroid Hormone-Like Protein) ELISA Kit
E-EL-S0687	Sheep C3 (Complement Component 3) ELISA Kit
E-EL-S0688	Sheep GAD (Glutamic Acid Decarboxylase) ELISA Kit
E-EL-S0689	Sheep PDGFD (Platelet Derived Growth Factor D) ELISA Kit
E-EL-S0690	Sheep TSGF (Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit
E-EL-S0691	Sheep IGFBP-6 (Insulin Like Growth Factor Binding Protein 6) ELISA Kit
E-EL-S0692	Sheep SMM (Smooth Muscle Myosin) ELISA Kit
E-EL-S0693	Sheep LEP(Leptin) ELISA Kit
E-EL-S0694	Sheep SUMO2 (Small Ubiquitin Related Modifier 2) ELISA Kit
E-EL-S0696	Sheep L-LDH (L-Lactate Dehydrogenase) ELISA Kit
E-EL-S0697	Sheep FBG(Fibrinogen)ELISA Kit
E-EL-S0698	Sheep C1q (Complement 1q) ELISA Kit
E-EL-S0700	Sheep IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit
E-EL-S0702	Sheep F-Box Protein 2(FBXO2) ELISA Kit
E-EL-S0703	Sheep NXF1 (Nuclear RNA Export Factor 1) ELISA Kit
E-EL-S0704	Sheep NNE (Non-Neuronal Enolase) ELISA Kit
E-EL-S0705	Sheep IgG3 (Immunoglobulin G3) ELISA Kit

Cat.NO	Product Name
E-EL-S0706	Sheep NTRK3 (Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit
E-EL-S0707	Sheep WARS (Tryptophanyl tRNA Synthetase) ELISA Kit
E-EL-S0708	Sheep AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit
E-EL-S0710	Sheep MUC17 (Mucin 17) ELISA Kit
E-EL-S0711	Sheep ANG2 (Angiopoietin 2) ELISA Kit
E-EL-S0712	Sheep COL13 (Collagen Type XIII) ELISA Kit
E-EL-S0713	Sheep Cyt-C (Cytochrome C) ELISA Kit
E-EL-S0714	Sheep FBL (Fibrillarin) ELISA Kit
E-EL-S0716	Sheep SOX18 (Sex Determining Region Y Box Protein 18) ELISA Kit
E-EL-S0718	Sheep HSP-60/HSPD1 (Heat Shock Protein 60) ELISA Kit
E-EL-S0719	Sheep EFNA4 (Ephrin A4) ELISA Kit
E-EL-S0720	Sheep IκBβ (Inhibitory Subunit of NF Kappa B Beta) ELISA Kit
E-EL-S0721	Sheep TNX (Tenascin X) ELISA Kit
E-EL-S0722	Sheep PCP (Procollagen C Peptidase) ELISA Kit
E-EL-S0723	Sheep Anti-Surv (Anti-Survivin Antibody) ELISA Kit
E-EL-S0724	Sheep PROC (Protein C) ELISA Kit
E-EL-S0725	Sheep PKCδ1 (Protein Kinase C Delta) ELISA Kit
E-EL-S0726	Sheep TE (Telomerase) ELISA Kit
E-EL-S0727	Sheep β-TG (β-Thromboglobulin) ELISA Kit
E-EL-S0728	Sheep Bcl2L11 (Bcl-2 Like Protein 11) ELISA Kit

Cat.NO	Product Name
E-EL-S0729	Sheep DHEA-S7 (Dehydroepiandrosterone S7) ELISA Kit
E-EL-S0730	Sheep IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
E-EL-S0731	Sheep BCAR1 (Breast Cancer Anti Estrogen Resistance 1) ELISA Kit
E-EL-S0732	Sheep BNP (Brain Natriuretic Peptide) ELISA Kit
E-EL-S0733	Sheep HLE (Leukocyte Exostase) ELISA Kit
E-EL-S0734	Sheep PGM (Phosphoglucomutase) ELISA Kit
E-EL-S0735	Sheep SHC/SLP-76 (SH2-Containing Protein) ELISA Kit
E-EL-S0736	Sheep TPA (Tissue Polypeptide Antigen) ELISA Kit
E-EL-S0737	Sheep DLAT (Dihydrolipoyl Transacetylase) ELISA Kit
E-EL-S0738	Sheep CLEC2C (C-Type Lectin Domain Family 2 Member C) ELISA Kit
E-EL-S0739	Sheep NUP35 (Nucleoporin 35kDa) ELISA Kit
E-EL-S0740	Sheep NOX5 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit
E-EL-S0741	Sheep TALLA-1 (T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit
E-EL-S0743	Sheep Try2 (Trypsinogen II) ELISA Kit
E-EL-S0744	Sheep E-Cad (E-Cadherin) ELISA Kit
E-EL-S0745	Sheep GAD1 (Glutamate Decarboxylase 1) ELISA Kit
E-EL-S0746	Sheep δ SIP- α (Delta Sleep Inducing Peptide α) ELISA Kit
E-EL-S0747	Sheep β APP (β Amyloid Precursor Protein) ELISA Kit
E-EL-S0748	Sheep CDNF (Cerebral Dopamine Neurotrophic Factor) ELISA Kit
E-EL-S0749	Sheep REG1 α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit

Cat.NO	Product Name
E-EL-S0750	Sheep COL8 α 1 (Collagen Type VIII Alpha 1) ELISA Kit
E-EL-S0752	Sheep TLN (Talin) ELISA Kit
E-EL-S0753	Sheep EGR3 (Early Growth Response Protein 3) ELISA Kit
E-EL-S0754	Sheep APSA (Anti-Phosphatidyl Serine Antibody) ELISA Kit
E-EL-S0755	Sheep MYH1 (Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit
E-EL-S0756	Sheep NUP133 (Nucleoporin 133kDa) ELISA Kit
E-EL-S0757	Sheep PI3K (Phosphoinositide-3-Kinase) ELISA Kit
E-EL-S0758	Sheep A-CoA (Acetyl Coenzyme A) ELISA Kit
E-EL-S0759	Sheep AKA (Anti-Keratin Antibody) ELISA Kit
E-EL-S0760	Sheep DCLK1 (Doublecortin Like Kinase 1) ELISA Kit
E-EL-S0761	Sheep SPTAN1 (Alpha-Fodrin) ELISA Kit
E-EL-S0762	Sheep SOCS3 (Suppressors of Cytokine Signaling 3) ELISA Kit
E-EL-S0763	Sheep PK (Pyruvate Kinase, Liver And RBC) ELISA Kit
E-EL-S0764	Sheep CLL1 (Collectin Liver 1) ELISA Kit
E-EL-S0765	Sheep APA (Anti-Phospholipid Antibody) ELISA Kit
E-EL-S0766	Sheep MMP-3 (Matrix Metalloproteinase 3) ELISA Kit
E-EL-S0767	Sheep ER α (Estrogen Receptor Alpha) ELISA Kit
E-EL-S0768	Sheep GATA1 (GATA Binding Protein 1) ELISA Kit
E-EL-S0769	Sheep SAA2 (Serum Amyloid A2) ELISA Kit
E-EL-S0770	Sheep DS (Dermatan Sulfate) ELISA Kit

Cat.NO	Product Name
E-EL-S0771	Sheep GC (Glycocalicin) ELISA Kit
E-EL-S0772	Sheep C4BP β (C4 Binding Protein Beta) ELISA Kit
E-EL-S0773	Sheep UBR4 (Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit
E-EL-S0775	Sheep ABCG2 (ATP Binding Cassette Transporter G2) ELISA Kit
E-EL-S0776	Sheep LDLR (Low Density Lipoprotein Receptor) ELISA Kit
E-EL-S0777	Sheep D-LDH (D-Lactate Dehydrogenase) ELISA Kit
E-EL-S0778	Sheep NOS1/nNOS (Nitric Oxide Synthase 1, Neuronal) ELISA Kit
E-EL-S0779	Sheep PPP1R1A (Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit
E-EL-S0780	Sheep A β 1-42 (Amyloid Beta Peptide 1-42) ELISA Kit
E-EL-S0781	Sheep MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit
E-EL-S0782	Sheep AGER (Advanced Glycosylation End Product Specific Receptor) ELISA Kit
E-EL-S0783	Sheep LH (Luteinizing Hormone) ELISA Kit
E-EL-S0784	Sheep GARS (Glycyl tRNA Synthetase) ELISA Kit
E-EL-S0785	Sheep P-GP (Permeability Glycoprotein) ELISA Kit
E-EL-S0786	Sheep PGF2 α (Prostaglandin F2 α) ELISA Kit
E-EL-S0788	Sheep CD59 (Protectin) ELISA Kit
E-EL-S0789	Sheep ANP (Atrial Natriuretic Peptide) ELISA Kit
E-EL-S0790	Sheep CA9 (Carbonic Anhydrase IX) ELISA Kit
E-EL-S0791	Sheep MT2 (Metallothionein 2) ELISA Kit
E-EL-S0792	Sheep EFNA1 (Ephrin A1) ELISA Kit

Cat.NO	Product Name
E-EL-S0793	Sheep ALDOA (Aldolase A, Fructose Bisphosphate) ELISA Kit
E-EL-S0794	Sheep MLZE (Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit
E-EL-S0795	Sheep ADCY3 (Adenylate Cyclase 3) ELISA Kit
E-EL-S0796	Sheep ASPN (Asporin) ELISA Kit
E-EL-S0797	Sheep BPI (Bactericidal/Permeability Increasing Protein) ELISA Kit
E-EL-S0798	Sheep CKLFSF4 (Chemokine Like Factor Superfamily 4) ELISA Kit
E-EL-S0800	Sheep PRB (Retinoblastoma Tumor Suppressor Protein) ELISA Kit
E-EL-S0801	Sheep plkB α (Phospho-Inhibitory Subunit Of NF Kappa B Alpha) ELISA Kit
E-EL-S0802	Sheep CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit
E-EL-S0803	Sheep ADAM8 (A Disintegrin And Metalloprotease 8) ELISA Kit
E-EL-S0804	Sheep C/EBP β (CCAAT/Enhancer Binding Protein Beta) ELISA Kit
E-EL-S0805	Sheep GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit
E-EL-S0807	Sheep PKC η (Protein Kinase C Eta) ELISA Kit
E-EL-S0808	Sheep CX40 (Connexin 40) ELISA Kit
E-EL-S0811	Sheep HMGCS1 (3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit
E-EL-S0813	Sheep PDI (Protein Disulfide Isomerase) ELISA Kit
E-EL-S0814	Sheep DLK-1 (Delta Like 1 Homolog) ELISA Kit
E-EL-S0815	Sheep PEDF (Pigment Epithelium Derived Factor) ELISA Kit
E-EL-S0816	Sheep PIR (Pirin) ELISA Kit
E-EL-S0817	Sheep MMP-7 (Matrix Metalloproteinase 7) ELISA Kit

Cat.NO	Product Name
E-EL-S0819	Sheep C1INH (Complement 1 Inhibitor) ELISA Kit
E-EL-S0820	Sheep PC (Plasma Anticoagulant Protein C) ELISA Kit
E-EL-S0821	Sheep GPC3 (Glypican 3) ELISA Kit
E-EL-S0822	Sheep GSANA (Granulocyte Specific Anti-Nuclear Antibody) ELISA Kit
E-EL-S0823	Sheep ACA-IgG (Anti-Cardiolipin Antibody IgG) ELISA Kit
E-EL-S0824	Sheep LAP (Leucine Aminopeptidase) ELISA Kit
E-EL-S0825	Sheep CK-18/KRT18 (Cytokeratin 18) ELISA Kit
E-EL-S0826	Sheep TGM2 (Transglutaminase 2, Tissue) ELISA Kit
E-EL-S0827	Sheep PYY (Peptide YY) ELISA Kit
E-EL-S0829	Sheep TSLP (Thymic Stromal Lymphopoietin) ELISA Kit
E-EL-S0830	Sheep SGK2 (Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit
E-EL-S0831	Sheep PTH (Parathyroid Hormone) ELISA Kit
E-EL-S0832	Sheep APT (Abnormal Prothrombin) ELISA Kit
E-EL-S0833	Sheep GTF2B (General Transcription Factor IIB) ELISA Kit
E-EL-S0835	Sheep τ PK1 (Tau-Protein Kinase) ELISA Kit
E-EL-S0836	Sheep ANG (Angiogenin) ELISA Kit
E-EL-S0837	Sheep MTHFR (Methylenetetrahydrofolate Reductase) ELISA Kit
E-EL-S0838	Sheep CEA (Carcinoembryonic Antigen) ELISA Kit
E-EL-S0839	Sheep FE (Ferritin) ELISA Kit
E-EL-S0840	Sheep PKC (Protein Kinase C) ELISA Kit

Cat.NO	Product Name
E-EL-S0841	Sheep P-CK (Pan-Cytokeratin) ELISA Kit
E-EL-S0842	Sheep ALK (Anaplastic Lymphoma Kinase) ELISA Kit
E-EL-S0843	Sheep CCL4L1 (Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit
E-EL-S0844	Sheep ANGPTL1 (Angiopoietin Like Protein 1) ELISA Kit
E-EL-S0846	Sheep PHLDA2 (Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit
E-EL-S0847	Sheep EMP (Erythrocyte Membrane Protein) ELISA Kit
E-EL-S0848	Sheep INSR (Insulin Receptor) ELISA Kit
E-EL-S0849	Sheep IL-4 (Interleukin 4) ELISA Kit
E-EL-S0851	Sheep PAlg (Platelet Associated Immunoglobulins) ELISA Kit
E-EL-S0852	Sheep CASP14 (Caspase 14) ELISA Kit
E-EL-S0854	Sheep α 1-AGP (Alpha-1-Acid Glycoprotein) ELISA Kit
E-EL-S0855	Sheep RND1 (Rho Family GTPase 1) ELISA Kit
E-EL-S0856	Sheep NRG-1 (Neuregulin 1) ELISA Kit
E-EL-S0857	Sheep Grb7 (Growth Factor Receptor Bound Protein 7) ELISA Kit
E-EL-S0858	Sheep CRY1 (Cryptochrome 1) ELISA Kit
E-EL-S0860	Sheep CTSD (Cathepsin D) ELISA Kit
E-EL-S0861	Sheep RARS (Arginyl tRNA Synthetase) ELISA Kit
E-EL-S0862	Sheep ST3 (Stromelysin-3) ELISA Kit
E-EL-S0863	Sheep PTGES (Prostaglandin E Synthase, Microsomal) ELISA Kit
E-EL-S0864	Sheep FGF19 (Fibroblast Growth Factor 19) ELISA Kit

Cat.NO	Product Name
E-EL-S0865	Sheep CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit
E-EL-S0866	Sheep Hpc (Hepcidin) ELISA Kit
E-EL-S0867	Sheep AAP (Alanine Aminopeptidase) ELISA Kit
E-EL-S0868	Sheep ACE (Angiotensin I Converting Enzyme) ELISA Kit
E-EL-S0869	Sheep CRLF1 (Cytokine Receptor Like Factor 1) ELISA Kit
E-EL-S0870	Sheep NFRkB (Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit
E-EL-S0871	Sheep PAR1 (Protease Activated Receptor 1) ELISA Kit
E-EL-S0872	Sheep MSMβ (Microseminoprotein Beta) ELISA Kit
E-EL-S0873	Sheep M-AChR M1 (Muscarinic Acetylcholine Receptor M1) ELISA Kit
E-EL-S0875	Sheep CASP1 (Caspase 1) ELISA Kit
E-EL-S0876	Sheep fDFT1 (Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit
E-EL-S0877	Sheep ATGA/TGAB (Anti-Thyroid-Globulin Antibody) ELISA Kit
E-EL-S0878	Sheep 6-K-PG (6-Keto-prostaglandin) ELISA Kit
E-EL-S0879	Sheep BMG/β2-MG (Beta-2-Microglobulin) ELISA Kit
E-EL-S0880	Sheep DAT (Dopamine Transporter) ELISA Kit
E-EL-S0881	Sheep PKR1 (Prokineticin Receptor 1) ELISA Kit
E-EL-S0882	Sheep MAP-2 (Microtubule Associated Protein 2) ELISA Kit
E-EL-S0883	Sheep F12 (Coagulation Factor XII) ELISA Kit
E-EL-S0884	Sheep α1-AT (Alpha 1-Antitrypsin) ELISA Kit
E-EL-S0885	Sheep VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit

Cat.NO	Product Name
E-EL-S0886	Sheep CYP1B1 (Cytochrome P450 1B1) ELISA Kit
E-EL-S0887	Sheep PIGR/SC (Polymeric Immunoglobulin Receptor/ Membrane Secretory Component) ELISA Kit
E-EL-S0888	Sheep FPR2 (Formyl Peptide Receptor 2) ELISA Kit
E-EL-S0889	Sheep CDK7 (Cyclin Dependent Kinase 7) ELISA Kit
E-EL-S0890	Sheep HB μ (Hemoglobin, Mu) ELISA Kit
E-EL-S0891	Sheep ELA4 (Elastase 4) ELISA Kit
E-EL-S0893	Sheep UGCG (UDP-glucose ceramide glucosyltransferase) ELISA Kit
E-EL-S0894	Sheep EPHA1 (Ephrin Type A Receptor 1) ELISA Kit
E-EL-S0895	Sheep TERT (Telomerase Reverse Transcriptase) ELISA Kit
E-EL-S0896	Sheep sICAM-1 (Soluble Intercellular Adhesion Molecule-1) ELISA Kit
E-EL-S0897	Sheep CTAGE (Cutaneous T-Cell Lymphocyte Associated Antigen) ELISA Kit
E-EL-S0898	Sheep IL-22 (Interleukin 22) ELISA Kit
E-EL-S0899	Sheep β 2-GP1 Ab IgA (β 2-Glycoprotein 1 antibody IgA) ELISA Kit
E-EL-S0900	Sheep ODC (Ornithine Decarboxylase) ELISA Kit
E-EL-S0901	Sheep PLC γ 2 (Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit
E-EL-S0902	Sheep OSF (Osteoclast Stimulating Factor 1) ELISA Kit
E-EL-S0903	Sheep cTnT/TNNT2 (Troponin T Type 2, Cardiac) ELISA Kit
E-EL-S0904	Sheep CASP 9 (Caspase 9) ELISA Kit
E-EL-S0906	Sheep CR1 (Complement Receptor type 1) ELISA Kit
E-EL-S0907	Sheep HA (Hemagglutinin) ELISA Kit

Cat.NO	Product Name
E-EL-S0908	Sheep RNASEH (Ribonuclease H)ELISA Kit
E-EL-S0909	Sheep SOX9 (Sex Determining Region Y Box Protein 9) ELISA Kit
E-EL-S0910	Sheep TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14)ELISA Kit
E-EL-S0911	Sheep CKLFSF6 (Chemokine Like Factor Superfamily 6) ELISA Kit
E-EL-S0912	Sheep EGR1 (Early Growth Response Protein 1) ELISA Kit
E-EL-S0915	Sheep PRSS1 (Protease, Serine, 1) ELISA Kit
E-EL-S0916	Sheep MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit
E-EL-S0917	Sheep EMA (Epithelial Membrane Antigen) ELISA Kit
E-EL-S0918	Sheep Gal (α -Galactoyl) ELISA Kit
E-EL-S0920	Sheep MUC13 (Mucin 13) ELISA Kit
E-EL-S0921	Sheep HLA-A (Leukocyte Antigen A) ELISA Kit
E-EL-S0922	Sheep IgA (Immumoglobulin A) ELISA Kit
E-EL-S0923	Sheep LPL (Lipoprotein Lipase) ELISA Kit
E-EL-S0924	Sheep GasR (Gastrin Receptor) ELISA Kit
E-EL-S0925	Sheep Slit2 (Slit Homolog 2) ELISA Kit
E-EL-S0927	Sheep PPAR- γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
E-EL-S0928	Sheep PCI (Protein C Inhibitor) ELISA Kit
E-EL-S0929	Sheep GPIIb/IIIa (Platelet Membrane Glycoprotein IIb/IIIa) ELISA Kit
E-EL-S0930	Sheep FABP4 (Fatty Acid-binding Protein 4, Adipocyte) ELISA Kit
E-EL-S0931	Sheep PDPN (Podoplanin) ELISA Kit

Cat.NO	Product Name
E-EL-S0932	Sheep GAL3 (Galectin 3) ELISA Kit
E-EL-S0933	Sheep NA/NE (Noradrenaline/Norepinephrine) ELISA Kit
E-EL-S0934	Sheep FOXO1 (Forkhead Box Protein O1) ELISA Kit
E-EL-S0935	Sheep NAGLU (N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit
E-EL-S0936	Sheep TPP (Thrombus Precursor Protein) ELISA Kit
E-EL-S0937	Sheep PCIII (Procollagen III) ELISA Kit
E-EL-S0938	Sheep ALB (Albumin) ELISA Kit
E-EL-S0940	Sheep BFP (Brain Finger Protein) ELISA Kit
E-EL-S0941	Sheep UBD (Ubiquitin D) ELISA Kit
E-EL-S0942	Sheep C1q-Ab (Anti-Complement 1q Antibody) ELISA Kit
E-EL-S0943	Sheep NOX4 (Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit
E-EL-S0944	Sheep CPP (Copeptin) ELISA Kit
E-EL-S0945	Sheep IgG2 (Immunoglobulin G2) ELISA Kit
E-EL-S0946	Sheep BMP-1 (Bone Morphogenetic Protein 1) ELISA Kit
E-EL-S0947	Sheep FSA (Fetal Sulfoslycoprotein Antigen) ELISA Kit
E-EL-S0948	Sheep AKAP1 (A Kinase Anchor Protein 1) ELISA Kit
E-EL-S0949	Sheep Tyk-2 (Tyrosine Kinase 2) ELISA Kit
E-EL-S0950	Sheep SEMA3A (Semaphorin 3A) ELISA Kit
E-EL-S0951	Sheep PP (Placental Protein) ELISA Kit
E-EL-S0952	Sheep CNP (C-type Natriuretic Peptide) ELISA Kit

Cat.NO	Product Name
E-EL-S0953	Sheep GZMK (Granzyme K) ELISA Kit
E-EL-S0954	Sheep DHEA (Dehydroepiandrosterone) ELISA Kit
E-EL-S0955	Sheep NE/ELA2 (Elastase 2, Neutrophil) ELISA Kit
E-EL-S0956	Sheep NSE (Neuron-Specific Enolase) ELISA Kit
E-EL-S0957	Sheep MX1 (Myxovirus Resistance 1) ELISA Kit
E-EL-S0958	Sheep STAT4 (Signal Transducer and Activator of Transcription 4) ELISA Kit
E-EL-S0960	Sheep ERCC1 (Excision Repair Cross Complementing Rodent Repair Deficiency Complementation 1) ELISA Kit
E-EL-S0962	Sheep BAFF/CD257/TNFSF13B (B-Cell Activating Factor) ELISA Kit
E-EL-S0963	Sheep EP3 (Prostaglandin E Receptor 3) ELISA Kit
E-EL-S0964	Sheep FPB (Fibrinopeptide B) ELISA Kit
E-EL-S0965	Sheep GAP (GTPase Activating Protein) ELISA Kit
E-EL-S0966	Sheep Extein ELISA Kit
E-EL-S0967	Sheep MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit
E-EL-S0968	Sheep SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit
E-EL-S0969	Sheep TTR (Transthyretin) ELISA Kit
E-EL-S0970	Sheep ANXA2 (Annexin A2) ELISA Kit
E-EL-S0971	Sheep S100A8 (S100 Calcium Binding Protein A8/Calgranulin A) ELISA Kit
E-EL-S0972	Sheep TDP43 (Tar DNA Binding Protein 43KDa)ELISA Kit
E-EL-S0973	Sheep CX3CL1 (Chemokine C-X3-C-Motif Ligand) ELISA Kit
E-EL-S0974	Sheep FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit

Cat.NO	Product Name
E-EL-S0975	Sheep α -CTx (Alpha Crosslaps) ELISA Kit
E-EL-S0976	Sheep VASPIN (Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit
E-EL-S0978	Sheep AR (Amphiregulin) ELISA Kit
E-EL-S0979	Sheep MRP1 (Multidrug Resistance Associated Protein 1) ELISA Kit
E-EL-S0980	Sheep SELL (L-Selectin) ELISA Kit
E-EL-S0981	Sheep SIGLEC1 (Sialic Acid Binding Ig Like Lectin 1) ELISA Kit
E-EL-S0982	Sheep HO2 (Heme Oxygenase 2) ELISA Kit
E-EL-S0983	Sheep CaN (Calcineurin) ELISA Kit
E-EL-S0986	Sheep CPN10 (Chaperonin 10) ELISA Kit
E-EL-S0987	Sheep AIF1 (Allograft Inflammatory Factor 1) ELISA Kit
E-EL-S0988	Sheep CCND1 (Cyclin-D1) ELISA Kit
E-EL-S0989	Sheep PRF1 (Perforin 1/Pore-Forming Protein) ELISA Kit
E-EL-S0990	Sheep MASP (Mannan Binding Lectin Associated Serine Protease) ELISA Kit
E-EL-S0991	Sheep HB (Hemoglobin) ELISA Kit
E-EL-S0992	Sheep Neopterin ELISA Kit
E-EL-S0993	Sheep CA4 (Carbonic Anhydrase IV) ELISA Kit
E-EL-S0994	Sheep CASP4 (Caspase 4) ELISA Kit
E-EL-S0995	Sheep MYL9 (Myosin Light Chain 9, Regulatory) ELISA Kit
E-EL-S0996	Sheep CG (Chorionic Gonadotropin) ELISA Kit
E-EL-S0997	Sheep IL-24 (Interleukin 24) ELISA Kit

Cat.NO	Product Name
E-EL-S0998	Sheep AP (Allopregnanolone) ELISA Kit
E-EL-S0999	Sheep AQP-2 (Aquaporin 2) ELISA Kit
E-EL-S1000	Sheep COL3 (Collagen Type III) ELISA Kit
E-EL-S1001	Sheep ApoA5 (Apolipoprotein A5) ELISA Kit
E-EL-S1002	Sheep ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-S1003	Sheep Ta1 (Thymosin- α 1) ELISA Kit
E-EL-S1004	Sheep PS1 (Presenilin 1) ELISA Kit
E-EL-S1005	Sheep PACAP-38 (Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
E-EL-S1006	Sheep PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-S1007	Sheep NOS2/iNOS (Nitric Oxide Synthase 2, Inducible) ELISA Kit
E-EL-S1008	Sheep IL-5 (Interleukin 5) ELISA Kit
E-EL-S1010	Sheep GHRL(Ghrelin/Obestatin Preprohormone)ELISA Kit
E-EL-S1011	Sheep OB(Obestatin)ELISA Kit
E-EL-S1012	Sheep Pg (Progesterone) ELISA Kit
E-EL-S1227	Sheep PRL (Prolactin) ELISA Kit
E-EL-S1263	Sheep IgG (Immunoglobulin G) ELISA Kit
E-EL-S1264	Sheep INS (Insulin) ELISA Kit
E-EL-S1265	Sheep ADP(Adiponectin) ELISA Kit
E-EL-S1266	Sheep EDN3(Endothelin-3)ELISA kit
E-EL-S1267	Sheep IgM (Immunoglobulin M) ELISA Kit

Cat.NO	Product Name
E-EL-S1268	Sheep CAST (Calpastatin) ELISA Kit
E-EL-S1269	Sheep IL-7 (Interleukin 7) ELISA Kit
E-EL-S1270	Sheep GM-CSF (Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit
E-EL-S1271	Sheep IL-17(Interleukin-17)ELISA Kit
E-EL-C0001	Canine MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit
E-EL-C0002	Canine IL-8 (Interleukin 8) ELISA
E-EL-C0003	Canine IFN- γ (Interferon Gamma) ELISA
E-EL-C0004	Canine TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit
E-EL-C0005	Canine IL-6 (Interleukin 6) ELISA Kit
E-EL-C0006	Canine IL-10 (Interleukin 10) ELISA Kit
E-EL-C0007	Canine TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit
E-EL-C0008	Canine VEGF (Vascular Endothelial cell Growth Factor) ELISA Kit
E-EL-C0009	Canine Fbg (Fibrinogen) ELISA Kit
E-EL-C0012	Canine ADP/Acrp30 (Adiponectin) ELISA Kit
E-EL-C0013	Canine C3(Complement 3)ELISA Kit
E-EL-C0014	Canine TRF (Transferrin) ELISA Kit
E-EL-C0015	Canine HP (Haptoglobin) ELISA Kit
E-EL-C0016	Canine SERPINA1 (Serine/Cysteine Peptidase Inhibitor Clade A Member 1) ELISA Kit
E-EL-C0017	Canine IgE (Immunoglobulin E) ELISA Kit
E-EL-C0018	Canine IgG (Immunoglobulin G) ELISA Kit

Cat.NO	Product Name
E-EL-C0019	Canine PLG (Plasminogen) ELISA Kit
E-EL-C0020	Canine CRP (C-Reactive Protein) ELISA Kit
E-EL-C0021	Canine Microalbumin ELISA Kit
E-EL-C0022	Canine MYO (Myoglobin) ELISA Kit
E-EL-C0023	Canine ADORA1 (Adenosine A1 Receptor) ELISA Kit
E-EL-C0024	Canine NA/NE (Noradrenaline/Norepinephrine) ELISA Kit
E-EL-C0025	Canine ACCPA (anti-cyclic citrullinated peptide antibody) ELISA Kit
E-EL-C0026	Canine β MSH (beta-Melanocyte Stimulating Hormone) ELISA Kit
E-EL-C0027	Canine Ta1 (Thymosin- α 1) ELISA Kit
E-EL-C0028	Canine PS1 (Presenilin 1) ELISA Kit
E-EL-C0029	Canine PACAP-38 (Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit
E-EL-C0030	Canine PACAP-27 (Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit
E-EL-C0031	Canine KIM-1 (Kidney Injury Molecule 1) ELISA Kit
E-EL-C0043	Canine BNP (Brain Natriuretic Peptide) ELISA Kit
E-EL-C0044	Canine CKMB (Creatine Kinase MB Isoenzyme) ELISA Kit
E-EL-C0045	Canine α 1-AGP (Alpha-1-Acid Glycoprotein) ELISA Kit
E-EL-C0046	Canine BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit
E-EL-C0047	Canine EPO (Erythropoietin) ELISA Kit
E-EL-C0048	Canine IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit
E-EL-C0049	Canine IgA (Immunoglobulin A) ELISA Kit

Cat.NO	Product Name
E-EL-C0050	Canine IL-13 (Interleukin 13) ELISA Kit
E-EL-C0051	Canine IL-5 (Interleukin 5) ELISA Kit
E-EL-C0052	Canine IL-4 (Interleukin 4) ELISA Kit
E-EL-C0053	Canine MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit
E-EL-C0054	Canine MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit
E-EL-C0055	Canine RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit
E-EL-C0056	Canine TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit
E-EL-C0057	Canine α 1-AT (Alpha 1-Antitrypsin) ELISA Kit
E-EL-C0058	Canine CD40L/TNFSF5 (Cluster of Differentiation 40 Ligand) ELISA Kit
E-EL-C0059	Canine IL-18 (Interleukin 18) ELISA Kit
E-EL-C0060	Canine CLU (Clusterin) ELISA Kit
E-EL-C0061	Canine SCF (Stem Cell Factor) ELISA Kit
E-EL-C0062	Canine IL-1 β (Interleukin 1 Beta) ELISA Kit
E-EL-C0063	Canine IL-12 p40 (Interleukin 12 p40) ELISA Kit
E-EL-C0064	Canine Pg (Progesterone) ELISA Kit
E-EL-C0065	Canine IFN- α (Interferon Alpha) ELISA Kit
E-EL-C0092	Canine IFN- β (Interferon Beta) ELISA Kit
E-EL-C0116	Canine IL-12 (Interleukin 12) ELISA Kit
E-EL-C0135	Canine TSHA (Thyroid Stimulating Hormone, Alpha) ELISA Kit
E-EL-C0136	Canine PAI-1 (Plasminogen Activator Inhibitor-1) ELISA Kit

Cat.NO	Product Name
E-EL-C0205	Canine S-100 (Soluble Protein-100) ELISA Kit
E-EL-C0207	Canine DMD (Dystrophin) ELISA Kit
E-EL-C0208	Canine Neopterin ELISA Kit
E-EL-C0209	Canine D2D (D-Dimer) ELISA Kit
E-EL-C0210	Canine NT-ProBNP (N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
E-EL-C0211	Canine Cys-C (Cystatin C) ELISA Kit
E-EL-C0212	Canine cTn-I/TNNI3 (cardiac Troponin I) ELISA Kit
E-EL-C0213	Canine GLP-2 (glucagon-like peptide-2) ELISA Kit
E-EL-C0214	Canine SELE (E-Selectin) ELISA Kit
E-EL-C0215	Canine β -OHB (Beta-Hydroxybutyric Acid) ELISA Kit
E-EL-C0216	Canine INS (Insulin) ELISA Kit
E-EL-C0217	Canine GM-CSF (Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit
E-EL-C0218	Canine GAP43 (Growth Associated Protein 43) ELISA Kit
E-EL-C0219	Canine CKB (Creatine Kinase, Brain) ELISA Kit
E-EL-C0220	Canine CKM (Creatine Kinase, Muscle) ELISA Kit
E-EL-C0222	Canine AST (Aspartate Aminotransferase) ELISA Kit
E-EL-C0223	Canine Ang-I (Angiotensin I) ELISA Kit
E-EL-C0224	Canine Ang-II (Angiotensin II) ELISA Kit
E-EL-C0225	Canine Ang-III (Angiotensin III) ELISA Kit
E-EL-C0226	Canine SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit

Cat.NO	Product Name
E-EL-C0227	Canine SOD2 (Superoxide Dismutase 2, Mitochondrial) ELISA Kit
E-EL-C0228	Canine SOD3 (Superoxide Dismutase 3, Extracellular) ELISA Kit
E-EL-C0229	Canine IL-1 α (Interleukin 1 Alpha) ELISA Kit
E-EL-C0230	Canine TNF- β (Tumor Necrosis Factor Beta) ELISA Kit
E-EL-C0231	Canine TCC C5b-9 (Terminal Complement Complex C5b-9) ELISA Kit
E-EL-C0232	Canine IL-2 (Interleukin 2) ELISA
E-EL-C0233	Canine GT (Gastrin) ELISA Kit
E-EL-C0234	Canine ET-1 (Endothelin 1) ELISA Kit
E-EL-0001	CecA (Cecropin A) ELISA Kit
E-EL-0002	CecB (Cecropin B) ELISA Kit
E-EL-0003	CecD (Cecropin D) ELISA Kit
E-EL-0004	TSST-1 (Toxic Shock Syndrome Toxin) ELISA Kit
E-EL-0005	CsA (Cyclosporine A) ELISA Kit
E-EL-0006	VA (Vitamin A) ELISA Kit
E-EL-0007	VB1 (Vitamin B1) ELISA Kit
E-EL-0008	VB6 (Vitamin B6) ELISA Kit
E-EL-0009	FA/VB9 (Folic Acid/Vitamin B9) ELISA Kit
E-EL-0010	VB12 (Vitamin B12) ELISA Kit
E-EL-0011	VC (Vitamin C) ELISA Kit
E-EL-0012	VD (Vitamin D) ELISA Kit

Cat.NO	Product Name
E-EL-0013	VD2 (Vitamin D2) ELISA Kit
E-EL-0014	VD3 (Vitamin D3) ELISA Kit
E-EL-0015	25-HVD3 (25-Hydroxy Vitamin D3) ELISA Kit
E-EL-0016	DHVD3 (1,25-Dihydroxyvitamin D3) ELISA Kit
E-EL-0017	VD4 (Vitamin D4) ELISA Kit
E-EL-0018	VE (Vitamin E) ELISA Kit
E-EL-0019	VK1 (Vitamin K1) ELISA Kit
E-EL-0020	MEK (Methionine-Enkephalin) ELISA Kit
E-EL-0021	LEK (Leucine-Enkephalin) ELISA Kit
E-EL-0022	PGI2 (Prostacyclin) ELISA Kit
E-EL-0023	ACSDKP (N-Acetyl-Ser-Asp-Lys-Pro) ELISA Kit
E-EL-0024	HYD (Hydrocortisone) ELISA Kit
E-EL-0025	LPS (Lipopolysaccharides) ELISA Kit
E-EL-0026	GSH (Glutathione) ELISA Kit
E-EL-0028	8-OHdG (8-Hydroxydeoxyguanosine) ELISA Kit
E-EL-0029	OT (Oxytocin) ELISA Kit
E-EL-0030	Cortisol ELISA Kit
E-EL-0031	DHT (Dihydrotestosterone) ELISA Kit
E-EL-0032	HIS (Histamine) ELISA Kit
E-EL-0033	ST/5-HT (Serotonin/5-Hydroxytryptamine) ELISA Kit

Cat.NO	Product Name
E-EL-0034	PGE2 (Prostaglandin E2) ELISA Kit
E-EL-0035	β-Nph (Beta-Naphthol) ELISA Kit
E-EL-0036	HA (Hyaluronic Acid) ELISA Kit
E-EL-0037	E3 (Estriol) ELISA Kit
E-EL-0038	F-E3 (Free Estriol)ELISA Kit
E-EL-0039	17-OHP (17-Hydroxyprogesterone) ELISA Kit
E-EL-0040	3-NT (3-Nitrotyrosine) ELISA Kit
E-EL-0041	8-epi-PGF2α (8-Epi-Prostaglandin F2 Alpha) ELISA Kit
E-EL-0042	ADMA (Asymmetrical Dimethylarginine) ELISA Kit
E-EL-0045	EPI (Epinephrine/Adrenaline) ELISA Kit
E-EL-0046	DA (Dopamine) ELISA Kit
E-EL-0047	NA/NE (Noradrenaline/Norepinephrine) ELISA Kit
E-EL-0048	VB2 (Vitamin B2) ELISA Kit
E-EL-0049	FMN (Flavin mononucleotide) ELISA Kit
E-EL-0050	FAD (Flavin adenine dinucleotide) ELISA Kit
E-EL-0051	AA (Arachidonic Acid) ELISA Kit
E-EL-0052	PGE1 (Prostaglandin E1) ELISA Kit
E-EL-0053	LXA4 (Lipoxin A4) ELISA Kit
E-EL-0054	6-keto-PGF1a (6-keto-prostaglandin F1a) ELISA Kit
E-EL-0055	2,3-dinor-TXB2 (2,3-dinor Thromboxane B2) ELISA Kit

Cat.NO	Product Name
E-EL-0056	cAMP (Cyclic adenosine monophosphate) ELISA Kit
E-EL-0057	TXA2 (Thromboxane A2) ELISA Kit
E-EL-0058	Cr (creatinine) ELISA Kit
E-EL-0059	IP3 (Inositol Triphosphate) ELISA Kit
E-EL-0060	MDA (Malondialdehyde) ELISA Kit
E-EL-0061	LTB4 (Leukotriene B4) ELISA Kit
E-EL-0062	E1 (Estrone) ELISA Kit
E-EL-0063	fMet (Formylmethionine) ELISA Kit
E-EL-0064	Cotinine ELISA Kit
E-EL-0065	E2 (Estradiol) ELISA Kit
E-EL-0066	PGD2 (Prostaglandin D2) ELISA Kit
E-EL-0067	SP (Substance P) ELISA Kit
E-EL-0068	ASD (Androstenedione) ELISA Kit
E-EL-0069	Car (Carnosine) ELISA Kit
E-EL-0070	ALD (Aldosterone) ELISA Kit
E-EL-0071	GnRH (Gonadotropin-Releasing Hormone) ELISA Kit
E-EL-0072	T (Testosterone) ELISA Kit
E-EL-0073	TNI(Troponin I)ELISA Kit
E-EL-0074	ADT(Androsterone)ELISA Kit
E-EL-0075	5-HIAA(5-Hydroxyindoleacetic Acid)ELISA Kit

Cat.NO	Product Name
E-EL-0076	Bb(Bilirubin) ELISA Kit
E-EL-0078	HVA (homovanillic acid) ELISA kit
E-EL-0081	ACH(Acetylcholine)ELISA Kit
E-EL-0082	VH (Vitamin H) ELISA Kit
E-EL-0083	cGMP(Cyclic GMP)ELISA Kit