


Elabscience®



Monoclonal Antibodies

Focus on Lab Research

Service for Life Science

About Elabscience

Elabscience specializes in immunodiagnostic technology for life science community. We have complete platform for R&D and manufacture. At the same time, we have in house QC for every Antibody, endeavoring to keep your experiment results more consistent and precise. Through unremitting effort and development, our customers have spread in more than 100 countries all over the world.

Elabscience major products cover Proteins, Antibodies, ELISA Kits, CLIA Kits, Labeling Kits and related reagents.

Elabscience also offers custom services for our customers including Protein Expression and Purification Service, Custom Antibody Service, Custom ELISA Service, Molecular Biology Service, Cell Biology Service and Whole Research Project Solutions.





Monoclonal Antibody

A monoclonal antibody specifically recognizes a single epitope within an antigen, thus having following advantages compared with polyclonal antibody, low cross-reactivity, high specificity, excellent lot-to-lot stability.

Elabscience offers mouse, rat, rabbit and other species derived recombinant monoclonal antibodies, all of which are strictly verified, so the quality could be guaranteed.

Elabscience also provides custom monoclonal antibody service, designing and producing high quality antibodies according to your project requirements, to benefit your research at the greatest extent.

EGFR Monoclonal Antibody

Catalog number: E-AB-22007

Synonyms: EGFR, ERBB, ERBB1, HER1, Epidermal growth factor receptor, Proto-oncogene c-ErbB-1, Receptor tyrosine-protein kinase erbB-1


Host: Mouse

Reactivity: H


Applications: WB, IF

Recommended dilution: WB 1:500-1:2000, IF 1:50-1:100

Receptor tyrosine kinase binding ligands of the EGF family and activating several signaling cascades to convert extracellular cues into appropriate cellular responses. Known ligands include EGF, TGFA/TGF- α , amphiregulin, EPGN/epigen, BTC/betacellulin, EREG/epiregulin and HBEGF/heparin-binding EGF. Ligand binding triggers receptor homo- and/or heterodimerization and autophosphorylation on key cytoplasmic residues. The phosphorylated receptor recruits adapter proteins like GRB2 which in turn activates complex downstream signaling cascades. It activates at least 4 major downstream signaling cascades including the RAS-RAF-MEK-ERK, PI3 kinase-AKT, PL Cgamma-PKC and STATs modules. It may also activate the NF-kappa-B signaling cascade. It also directly phosphorylates other proteins like RGS16, activating its GTPase activity and probably coupling the EGF receptor signaling to the G protein-coupled receptor signaling. It also phosphorylates MUC1 and increases its interaction with SRC and CTNNB1/beta-catenin.



Western Blot analysis of Hela cells using E-AB-22007 EGFR Monoclonal Antibody at dilution of 1:1000



Immunofluorescence analysis of Hela cells using E-AB-22007 EGFR Monoclonal Antibody at dilution of 1:100



p53 Monoclonal Antibody

Catalog number: E-AB-22015

Synonyms: TP53, P53, Cellular tumor antigen p53, Antigen NY-CO-13, Phosphoprotein p53, Tumor suppressor p53


Host: Mouse

Reactivity: H


Applications: WB, IHC, IF

Recommended dilution: WB 1:500-1:2000, IHC 1:50-1:200, IF 1:50-1:100


P53 acts as a tumor suppressor in many tumor types; induces growth arrest or apoptosis depending on the physiological circumstances and cell type. Involved in cell cycle regulation as a trans-activator that acts to negatively regulate cell division by controlling a set of genes required for this process. One of the activated genes is an inhibitor of cyclin-dependent kinases. Apoptosis induction seems to be mediated either by stimulation of BAX and FAS antigen expression, or by repression of Bcl-2 expression, implicated in Notch signaling cross-over. Isoform 2 enhances the transactivation activity of isoform 1 from some but not all TP53-inducible promoters. Isoform 4 suppresses transactivation activity and impairs growth suppression mediated by isoform 1. Isoform 7 inhibits isoform 1-mediated apoptosis.



Immunohistochemistry of paraffin-embedded Human colon cancer tissue using E-AB-22015 P53 Monoclonal antibody at dilution of 1:200



Immunofluorescence analysis of HeLa cells using E-AB-22015 P53 Monoclonal Antibody at dilution of 1:100



Western Blot analysis of 293T cells using E-AB-22015 P53 Monoclonal Antibody at dilution of 1:1000 (1) and 1:2000 (2)



Monoclonal Antibodies

Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27413	14-3-3 Monoclonal Antibody	P27348	Mouse	H,R	ELISA,WB,IHC,IF
E-AB-27562	14-3-3 sigma Monoclonal Antibody	P31947	Mouse	H	ELISA,WB,IHC,IF
E-AB-27001	4EBP1 Monoclonal Antibody	Q13541	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-22080	α skeletal muscle actin Monoclonal Antibody	P68133	Mouse	H,M,R	WB,IP
E-AB-27633	α -Synuclein Monoclonal Antibody	P37840	Mouse	H,M,R	ELISA,WB
E-AB-27002	A1BG Monoclonal Antibody	P04217	Mouse	H	ELISA,WB,IHC,IF
E-AB-22032	ABCB5 Monoclonal Antibody	Q2M3G0	Mouse	H	WB
E-AB-27003	ABCD1 Monoclonal Antibody	P33897	Mouse	H	ELISA,WB,IF
E-AB-27005	ACPP Monoclonal Antibody	P15309	Mouse	H	ELISA,WB,IHC
E-AB-27018	ACTC1 Monoclonal Antibody	P68032	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27006	ADIPOQ Monoclonal Antibody	Q15848	Mouse	H	ELISA,WB
E-AB-27623	ADRM1 Monoclonal Antibody	Q16186	Mouse	H	ELISA,WB
E-AB-27414	ADRP Monoclonal Antibody	Q99541	Mouse	H	ELISA,WB,IHC
E-AB-27007	AEBP2 Monoclonal Antibody	Q6ZN18	Mouse	H	ELISA,WB
E-AB-27008	AFM Monoclonal Antibody	P43652	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-22033	AFP alpha 1 Fetoprotein Monoclonal Antibody	P02771	Mouse	H	WB
E-AB-27624	AGTR1 Monoclonal Antibody	P30556	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27009	AHCY Monoclonal Antibody	P23526	Mouse	H	ELISA,WB,IHC,IF
E-AB-27415	AIRE Monoclonal Antibody	O43918	Mouse	H	ELISA,WB
E-AB-27010	AK2 Monoclonal Antibody	P54819	Mouse	H	ELISA,WB,IHC,IF
E-AB-27416	AKAP3 Monoclonal Antibody	O75969	Mouse	H	ELISA,WB
E-AB-27012	AKT Monoclonal Antibody	P31749	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27011	AKT1/3 Monoclonal Antibody	P31749	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27013	Albumin Monoclonal Antibody	Q13740	Mouse	H,R	ELISA,WB,IHC,IF
E-AB-27014	ALDH1A1 Monoclonal Antibody	P00352	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27015	Aldolase C Monoclonal Antibody	P04075	Mouse	H	ELISA,WB,IHC,IF
E-AB-27016	ALK Monoclonal Antibody	Q9UM73	Mouse	H	ELISA,WB,IHC
E-AB-27625	ALOX5 Monoclonal Antibody	P09917	Mouse	H,Pig	ELISA,WB,IHC
E-AB-27017	Alpha 1 Antichymotrypsin Monoclonal Antibody	P01011	Mouse	H	ELISA,WB,IHC
E-AB-27019	Alpha galactosidase A Monoclonal Antibody	P06280	Mouse	H	ELISA,WB,IHC,IF



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-22042	Alpha Lactalbumin Monoclonal Antibody	P00709	Mouse	H	WB
E-AB-22042	Alpha Lactalbumin Monoclonal Antibody	P01009	Mouse	H,R	WB
E-AB-22042	Alpha Lactalbumin Monoclonal Antibody	P01023	Mouse	H	ELISA, WB, IHC, IF
E-AB-22010	AMACR Monoclonal Antibody	Q9UHK6	Mouse	H,M,R	WB, IHC
E-AB-27022	AMOT Monoclonal Antibody	Q4VCS5	Mouse	H	ELISA, WB, IHC
E-AB-27628	Amphiregulin Monoclonal Antibody	P15514	Mouse	H,R	ELISA, WB, IHC
E-AB-27023	Amylase alpha Monoclonal Antibody	P06278	Mouse	H,R	ELISA, WB, IHC, IF
E-AB-27025	Angiotensinogen Monoclonal Antibody	P01019	Mouse	H	ELISA, WB, IHC
E-AB-27629	Annexin A1 Monoclonal Antibody	P04083	Mouse	H,M	ELISA, WB, IHC, IF
E-AB-27026	Annexin A2 Monoclonal Antibody	P07355	Mouse	H,M	ELISA, WB, IHC, IF
E-AB-27563	Annexin A2 Monoclonal Antibody	P07355	Mouse	H	ELISA, WB, IHC, IF
E-AB-27630	Annexin A3 Monoclonal Antibody	P12429	Mouse	H	ELISA, WB, IHC
E-AB-27418	Annexin IV Monoclonal Antibody	P09525	Mouse	H	ELISA, WB, IHC, IF
E-AB-27419	Annexin V Monoclonal Antibody	P08758	Mouse	H	ELISA, WB, IHC, IF
E-AB-27027	ANP32B Monoclonal Antibody	Q92688	Mouse	H	ELISA, WB, IF
E-AB-27028	Antithrombin III Monoclonal Antibody	P01008	Mouse	H,M,R	ELISA, WB, IHC, IF
E-AB-27029	APOA1 Monoclonal Antibody	P02647	Mouse	H,M	ELISA, WB, IHC
E-AB-27420	APOA4 Monoclonal Antibody	P06727	Mouse	H	ELISA, WB
E-AB-27030	APOBEC3G Monoclonal Antibody	Q9HC16	Mouse	H,M,R	ELISA, WB, IHC, IF
E-AB-27031	APOD Monoclonal Antibody	P05090	Mouse	H	ELISA, WB
E-AB-27032	APOL1 Monoclonal Antibody	O14791	Mouse	H	ELISA, WB, IHC, IF
E-AB-27033	Apolipoprotein H Monoclonal Antibody	P02655	Mouse	H	ELISA, WB, IHC
E-AB-27024	APP Monoclonal Antibody	P05067	Mouse	H	ELISA, WB, IHC
E-AB-22025	Aquaporin 4 Monoclonal Antibody	P55087	Mouse	H,M,R	WB, IF
E-AB-27235	ARG1 Monoclonal Antibody	P07824	Mouse	H,M,R	ELISA, WB, IHC
E-AB-27631	ARHGDI A Monoclonal Antibody	P52565	Mouse	H,M,R	ELISA, WB
E-AB-27034	ARHGDI B Monoclonal Antibody	P52566	Mouse	H	ELISA, WB
E-AB-27632	ASNA1 Monoclonal Antibody	O43681	Mouse	H,M,R	ELISA, WB, IHC
E-AB-27035	ASS1 Monoclonal Antibody	P00966	Mouse	H	ELISA, WB, IHC, IF
E-AB-27036	ATF4 Monoclonal Antibody	P18848	Mouse	H,M	ELISA, WB, IHC, IF
E-AB-27037	ATG5 Monoclonal Antibody	Q9H1Y0	Mouse	H,M	ELISA, WB
E-AB-27038	ATP5A1 Monoclonal Antibody	P25705	Mouse	H,M,R	ELISA, WB, IHC, IF



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27039	ATP5C1 Monoclonal Antibody	P36542	Mouse	H	ELISA,WB,IHC
E-AB-27040	Aurora B Monoclonal Antibody	Q96GD4	Mouse	H,M	ELISA,WB,IHC
E-AB-27041	AZGP1 Monoclonal Antibody	P25311	Mouse	H	ELISA,WB,IHC
E-AB-22028	β -Catenin Monoclonal Antibody	P35222	Mouse	H,M,R	WB,IHC
E-AB-27640	β -Catenin Monoclonal Antibody	P35222	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27042	B23 Monoclonal Antibody	P06748	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27043	B2M Monoclonal Antibody	P61769	Mouse	H	ELISA,WB,IHC
E-AB-27634	B7-1 Monoclonal Antibody	P33681	Mouse	H,M	ELISA,WB,IHC
E-AB-27044	BAX Monoclonal Antibody	Q07812	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27421	BBS2 Monoclonal Antibody	Q9BXC9	Mouse	H	ELISA,WB,IF
E-AB-27045	BBS5 Monoclonal Antibody	Q8N3I7	Mouse	H,M	ELISA,WB,IHC
E-AB-27046	BCL2 Monoclonal Antibody	P10415	Mouse	H,R	ELISA,WB,IHC
E-AB-22004	Bcl-2 Monoclonal Antibody	P10415	Mouse	H,M,R	WB,IHC
E-AB-27635	BCL6 Monoclonal Antibody	P41182	Mouse	H	ELISA,WB,IHC
E-AB-27047	Bcl-XL Monoclonal Antibody	Q07817	Mouse	H	ELISA,WB,IHC,IF
E-AB-27048	BCS1L Monoclonal Antibody	Q9Y276	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27422	BDNF Monoclonal Antibody	P23560	Mouse	H,R	ELISA,WB,IF
E-AB-27049	BEST1 Monoclonal Antibody	O76090	Mouse	H	ELISA,WB
E-AB-27050	BID Monoclonal Antibody	P55957	Mouse	H	ELISA,WB,IHC,IF
E-AB-27423	BMI1 Monoclonal Antibody	P35226	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27642	BMP2 Monoclonal Antibody	P12643	Mouse	H,Pig	ELISA,WB,IHC
E-AB-27643	BMP5 Monoclonal Antibody	P22003	Mouse	H	ELISA,WB
E-AB-27424	BORIS Monoclonal Antibody	Q8NI51	Mouse	H,M	ELISA,WB,IF
E-AB-27425	BRCA1 Monoclonal Antibody	P38398	Mouse	H	ELISA,WB,IHC
E-AB-27051	BRD8 Monoclonal Antibody	Q9H0E9	Mouse	H	ELISA,WB,IHC,IF
E-AB-27052	BrdU Monoclonal Antibody		Mouse	H,M,R	ELISA,IHC,IF
E-AB-27644	BRE Monoclonal Antibody	Q9NXR7	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27053	BSA Monoclonal Antibody	P02769	Mouse	Bovine	ELISA,WB
E-AB-27427	BSAP Monoclonal Antibody	Q02548	Mouse	H	ELISA,WB,IHC
E-AB-27054	C1qC Monoclonal Antibody	P02747	Mouse	H	ELISA,WB
E-AB-27055	C3/C3b/C3c Monoclonal Antibody	P01024	Mouse	H	ELISA,WB,IHC,IF
E-AB-27056	C4 (gamma chain) Monoclonal Antibody	P0C0L4	Mouse	H	ELISA,WB



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27057	C9orf72 Monoclonal Antibody	Q96LT7	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-22034	CA IX Monoclonal Antibody	Q16790	Mouse	H	WB,IP
E-AB-27428	CA9 Monoclonal Antibody	Q16790	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27645	Calbindin-D28k Monoclonal Antibody	P05937	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27646	CALCOCO2 Monoclonal Antibody	Q13137	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27647	Calsequestrin 2 Monoclonal Antibody	O14958	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27648	CAMSAP2 Monoclonal Antibody	Q08AD1	Mouse	H,M,R	ELISA,WB
E-AB-27058	CAPG Monoclonal Antibody	P40121	Mouse	H	ELISA,WB,IHC,IF
E-AB-27649	CAPRIN1 Monoclonal Antibody	Q14444	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27059	CAPZA1 Monoclonal Antibody	P52907	Mouse	H,M,R	ELISA,WB,IF
E-AB-22024	Carcinoembryonic Antigen Monoclonal Antibody	P06731	Mouse	H	IHC
E-AB-27650	Cardiac Troponin I Monoclonal Antibody	P19429	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27651	Caspase 3 Monoclonal Antibody	P42574	Mouse	H,M,R	ELISA,WB
E-AB-27060	Caspase 8 Monoclonal Antibody	Q14790	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-22035	Caspase 9 Monoclonal Antibody	P55211	Mouse	H,M,R	WB,IP
E-AB-27061	Caspase 9/p35/p10 Monoclonal Antibody	P55211	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27062	Caveolin-1 Monoclonal Antibody	Q03135	Mouse	H	ELISA,WB,IHC
E-AB-27652	CBLB Monoclonal Antibody	Q13191	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27653	CBX3 Monoclonal Antibody	Q13185	Mouse	H,M	ELISA,WB,IF
E-AB-27063	CCBP2 Monoclonal Antibody	Q58D77	Mouse	H	ELISA,IHC
E-AB-27429	CCT3 Monoclonal Antibody	P49368	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27064	CD109 Monoclonal Antibody	Q6YHK3	Mouse	H	ELISA,WB
E-AB-27564	CD11a Monoclonal Antibody	P20701	Mouse	H	ELISA,IHC
E-AB-27065	CD11c Monoclonal Antibody	P20702	Mouse	H,M	ELISA,WB,IHC
E-AB-27066	CD13 Monoclonal Antibody	P15144	Mouse	H	ELISA,WB,IHC
E-AB-27067	CD138 Monoclonal Antibody	P18827	Mouse	H	ELISA,WB
E-AB-27068	CD138 Monoclonal Antibody	P18827	Mouse	H	ELISA,WB,IHC
E-AB-27069	CD14 Monoclonal Antibody	P08571	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27070	CD146 Monoclonal Antibody	P43121	Mouse	H	ELISA,WB,IHC
E-AB-27654	CD147 Monoclonal Antibody	P35613	Mouse	H	ELISA,WB,IHC
E-AB-22048	CD16 Monoclonal Antibody	P08637	Mouse	H	WB
E-AB-27430	CD19 Monoclonal Antibody	P15391	Mouse	H	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27071	CD1A Monoclonal Antibody	P06126	Mouse	H	ELISA,WB
E-AB-27431	CD1d Monoclonal Antibody	P15813	Mouse	H	ELISA,WB,IHC
E-AB-27655	CD2 Monoclonal Antibody	P06729	Mouse	H	ELISA,WB,IHC
E-AB-22011	CD20 Monoclonal Antibody	P11836	Mouse	H,M,R	IHC
E-AB-27072	CD20 Monoclonal Antibody	P11836	Mouse	H	ELISA,WB,IHC
E-AB-27432	CD200 Monoclonal Antibody	P41217	Mouse	H	ELISA,WB,IHC
E-AB-27073	CD206 Monoclonal Antibody	P22897	Mouse	H,M,R	ELISA,WB
E-AB-22023	CD21 Monoclonal Antibody	P20023	Mouse	H,M,R	IHC
E-AB-27074	CD22 Monoclonal Antibody	P20273	Mouse	H	ELISA,WB,IHC
E-AB-22012	CD23 Monoclonal Antibody	P06734	Mouse	H,M,R	IHC
E-AB-22040	CD23 Monoclonal Antibody	P06734	Mouse	H,M,R	WB
E-AB-27075	CD23 Monoclonal Antibody	P06734	Mouse	H,M	ELISA,WB,IHC
E-AB-27076	CD23 Monoclonal Antibody	P06734	Mouse	H,M	ELISA,WB,IHC
E-AB-27433	CD27 Monoclonal Antibody	P26842	Mouse	H	ELISA,WB,IHC
E-AB-27656	CD276 Monoclonal Antibody	Q5ZPR3	Mouse	H	ELISA,WB
E-AB-27077	CD3 Monoclonal Antibody	P07766	Mouse	H	ELISA,WB,IHC
E-AB-27565	CD3 delta Monoclonal Antibody	P04234	Mouse	H,M	ELISA,WB
E-AB-27078	CD3 gamma Monoclonal Antibody	P09693	Mouse	H	ELISA,WB,IHC
E-AB-27079	CD31 Monoclonal Antibody	P16284	Mouse	H	ELISA,WB,IHC,IF
E-AB-27080	CD34 Monoclonal Antibody	P28906	Mouse	H	ELISA,WB,IHC
E-AB-27657	CD36 Monoclonal Antibody	P16671	Mouse	H	ELISA,WB,IHC
E-AB-27081	CD38 Monoclonal Antibody	P28907	Mouse	H,M	ELISA,WB,IF
E-AB-27082	CD41 Monoclonal Antibody	P08514	Mouse	H	ELISA,WB
E-AB-27083	CD43 Monoclonal Antibody	P16150	Mouse	H	ELISA,WB,IHC
E-AB-27084	CD44 Monoclonal Antibody	P16070	Mouse	H	ELISA,WB,IHC,IF
E-AB-22009	CD45 Monoclonal Antibody	P08575	Mouse	H	WB
E-AB-27085	CD45 Monoclonal Antibody	P08575	Mouse	H	ELISA,WB,IHC
E-AB-27086	CD47 Monoclonal Antibody	Q08722	Mouse	H	ELISA,WB
E-AB-22013	CD68 Monoclonal Antibody	P34810	Mouse	H,M,R	IHC
E-AB-27087	CD68 Monoclonal Antibody	P34810	Mouse	H,M	ELISA,IHC
E-AB-27089	CD7 Monoclonal Antibody	P09564	Mouse	H,M	ELISA,IHC
E-AB-27090	CD7 Monoclonal Antibody	P09564	Mouse	H,M	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27091	CD71 Monoclonal Antibody	P02786	Mouse	H	ELISA,WB,IHC,IF
E-AB-27658	CD74 Monoclonal Antibody	P04233	Mouse	H	ELISA,WB
E-AB-27092	CD9 Monoclonal Antibody	P21926	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27659	CD99 Monoclonal Antibody	P14209	Mouse	H	ELISA,WB
E-AB-27660	CDC37 Monoclonal Antibody	Q16543	Mouse	H,M,R	ELISA,WB,IF
E-AB-27093	CDC6 Monoclonal Antibody	Q99741	Mouse	H,M,R	ELISA,WB,IF
E-AB-27435	CDH17 Monoclonal Antibody	Q12864	Mouse	H	ELISA,WB,IHC
E-AB-27094	CDK2 Monoclonal Antibody	P24941	Mouse	H	ELISA,WB,IHC
E-AB-27437	CDK6 Monoclonal Antibody	Q00534	Mouse	H,M,R	ELISA,WB,IF
E-AB-22020	CDX2 Monoclonal Antibody	Q99626	Mouse	H,M,R	WB,IHC
E-AB-27095	CDX2 Monoclonal Antibody	Q99626	Mouse	H,M,R	ELISA,WB
E-AB-27096	Ceruloplasmin Monoclonal Antibody	P00450	Mouse	H	ELISA,WB,IHC
E-AB-27438	CGB Monoclonal Antibody	P0DN86	Mouse	H	ELISA,WB,IHC
E-AB-27439	CHCHD2 Monoclonal Antibody	Q9Y6H1	Mouse	H,R	ELISA,WB,IHC,IF
E-AB-27440	CHD4 Monoclonal Antibody	Q14839	Mouse	H	ELISA,WB,IF
E-AB-27097	Chk1 Monoclonal Antibody	O14757	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27098	CHOP; GADD153 Monoclonal Antibody	P35638	Mouse	H	ELISA,WB
E-AB-27661	CHRNA5 Monoclonal Antibody	P30532	Mouse	H,M	ELISA,WB,IHC
E-AB-27099	CHRNA7 Monoclonal Antibody	P36544	Mouse	H	ELISA,WB,IHC
E-AB-27100	Chromogranin A Monoclonal Antibody	P10645	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27101	CISD2 Monoclonal Antibody	Q8N5K1	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-22019	CK16 Monoclonal Antibody	P08779	Mouse	H,M,R	IHC
E-AB-22039	CK17 Monoclonal Antibody	Q04695	Mouse	H	WB,IP
E-AB-22014	CK19 Monoclonal Antibody	P08727	Mouse	H,M,R	IHC
E-AB-22017	CK7 Monoclonal Antibody	P08729	Mouse	H,M,R	WB,IHC,IP
E-AB-27662	CK7 Monoclonal Antibody	P08729	Mouse	H	ELISA,WB,IHC
E-AB-22018	CK8 Monoclonal Antibody	P05787	Mouse	H,M,R	WB,IHC
E-AB-27116	CKM Monoclonal Antibody	P06732	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27102	Claudin 18 Monoclonal Antibody	P56856	Mouse	H,M	ELISA,IHC
E-AB-27103	CLEC11A Monoclonal Antibody	Q9Y240	Mouse	H	ELISA,WB,IHC
E-AB-27444	CLEC4G Monoclonal Antibody	Q6UXB4	Mouse	H	ELISA,WB
E-AB-27663	CLIC4 Monoclonal Antibody	Q9Y696	Mouse	H,R	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27445	CLPP Monoclonal Antibody	Q16740	Mouse	H	ELISA,WB,IHC
E-AB-27104	CLTB Monoclonal Antibody	P09497	Mouse	H,Pig	ELISA,WB,IHC,IF
E-AB-27664	CLTC Monoclonal Antibody	Q00610	Mouse	H,M,R	ELISA,WB
E-AB-27105	CLTCL1 Monoclonal Antibody	P53675	Mouse	H	ELISA,WB,IHC,IF
E-AB-27106	Clusterin Monoclonal Antibody	P10909	Mouse	H	ELISA,WB,IHC
E-AB-27107	CNN2 Monoclonal Antibody	Q99439	Mouse	H	ELISA,WB,IHC
E-AB-27108	CNPY2 Monoclonal Antibody	Q9Y2B0	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27109	Cofilin Monoclonal Antibody	P23528	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-22071	Collagen III Monoclonal Antibody	P02461	Mouse	H,M,R	WB
E-AB-27110	Complement factor B Monoclonal Antibody	P00751	Mouse	H	ELISA,WB,IHC
E-AB-27665	COPS2 Monoclonal Antibody	P61201	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27111	COTL1 Monoclonal Antibody	Q14019	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-22002	COX IV Monoclonal Antibody	P13073	Mouse	H,M,R	WB
E-AB-27666	COX2 Monoclonal Antibody	P35354	Mouse	H,M	ELISA,WB,IHC
E-AB-27112	COXIV Monoclonal Antibody	P13073	Mouse	H	ELISA,WB,IHC,IF
E-AB-27113	COXIV Monoclonal Antibody	P13073	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27667	CP110 Monoclonal Antibody	O43303	Mouse	H	ELISA,WB
E-AB-27114	CPT1A Monoclonal Antibody	P50416	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27115	CPT1C Monoclonal Antibody	Q8TCG5	Mouse	H	ELISA,WB,IHC
E-AB-27668	CRABP2 Monoclonal Antibody	P29373	Mouse	H,M,R	ELISA,WB
E-AB-27669	CRBN Monoclonal Antibody	Q96SW2	Mouse	H,M	ELISA,WB,IHC
E-AB-27117	CREST Monoclonal Antibody	O75177	Mouse	H	ELISA,WB
E-AB-27118	CRP Monoclonal Antibody	P02741	Mouse	H,M	ELISA,WB,IHC
E-AB-27447	CRYBB1 Monoclonal Antibody	P53674	Mouse	H,M,R	ELISA,WB
E-AB-27119	CSK Monoclonal Antibody	P41240	Mouse	H,M	ELISA,WB,IF
E-AB-27120	c-SRC Monoclonal Antibody	P12931	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27448	CTAGE1 Monoclonal Antibody	Q9HC47	Mouse	H	ELISA,WB,IHC,IF
E-AB-27449	CTNNA1 Monoclonal Antibody	P35221	Mouse	H,M	ELISA,WB,IHC
E-AB-27121	CUL4A Monoclonal Antibody	Q13619	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27597	CUL4B Monoclonal Antibody	Q13620	Mouse	H,R	ELISA,WB,IHC
E-AB-27450	CX3CL1 Monoclonal Antibody	P78423	Mouse	H	ELISA,WB,IHC
E-AB-27122	CXCL16 Monoclonal Antibody	Q9H2A7	Mouse	H,M	ELISA,WB



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27123	CXCL8 Monoclonal Antibody	P10145	Mouse	H	ELISA,WB,IHC
E-AB-27124	CXCR3B Monoclonal Antibody	P49682	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27125	CXCR4 Monoclonal Antibody	P61073	Mouse	H,M	ELISA,WB
E-AB-27126	CXCR7 Monoclonal Antibody	P25106	Mouse	H,M	ELISA,WB,IHC
E-AB-27127	CYB5R3 Monoclonal Antibody	P00387	Mouse	H	ELISA,WB,IHC,IF
E-AB-27670	Cyclin A2 Monoclonal Antibody	P20248	Mouse	H	ELISA,WB,IHC,IF
E-AB-27128	Cyclin D1 Monoclonal Antibody	P24385	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27671	Cyclin D3 Monoclonal Antibody	P30281	Mouse	H	ELISA,WB,IHC
E-AB-27566	Cyclophilin B Monoclonal Antibody	P23284	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27451	CYP2D6 Monoclonal Antibody	P10635	Mouse	H	ELISA,WB
E-AB-27672	Cystatin C Monoclonal Antibody	P01034	Mouse	H	ELISA,WB,IF
E-AB-27129	Cytochrome c Monoclonal Antibody	P99999	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27130	Cytoglobin Monoclonal Antibody	Q8WWM9	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27131	Cytokeratin 14 Monoclonal Antibody	Q6A162	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27132	Cytokeratin 15 Monoclonal Antibody	P19012	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-22053	Cytokeratin 18 Monoclonal Antibody	P05783	Mouse	H,M,R	WB,IF
E-AB-22054	Cytokeratin 18 Monoclonal Antibody	P05783	Mouse	H,M	WB
E-AB-22055	Cytokeratin 18 Monoclonal Antibody	P05783	Mouse	H,M,R	WB
E-AB-27133	Cytokeratin 18 Monoclonal Antibody	P05783	Mouse	H	ELISA,WB,IHC,IF
E-AB-27134	Cytokeratin 19 Monoclonal Antibody	P08727	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27135	Cytokeratin 20 Monoclonal Antibody	P35900	Mouse	H	ELISA,WB,IHC
E-AB-22038	Cytokeratin 6 Monoclonal Antibody	P02538	Mouse	H	WB
E-AB-27136	Cytokeratin 8 Monoclonal Antibody	P05787	Mouse	H	ELISA,WB,IHC
E-AB-27673	D2HGDH Monoclonal Antibody	Q8N465	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27598	DACH1 Monoclonal Antibody	Q9UI36	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27452	DAO Monoclonal Antibody	P14920	Mouse	H	ELISA,WB,IHC
E-AB-27453	DBC1 Monoclonal Antibody	O60477	Mouse	H,M	ELISA,WB,IF
E-AB-27137	DDB1 Monoclonal Antibody	Q16531	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27138	DEK Monoclonal Antibody	P35659	Mouse	H	ELISA,WB,IHC,IF
E-AB-27674	Dermatopontin Monoclonal Antibody	Q07507	Mouse	H,R	ELISA,WB
E-AB-27139	Desmin Monoclonal Antibody	P17661	Mouse	H,M	ELISA,WB,IHC
E-AB-27140	Desmocollin 2 Monoclonal Antibody	Q02487	Mouse	H,M	ELISA,WB



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27141	DGCR8 Monoclonal Antibody	Q8WYQ5	Mouse	H,M,R	ELISA,WB
E-AB-27142	DSE Monoclonal Antibody	Q9UL01	Mouse	H,M	ELISA,WB,IHC
E-AB-27143	Dynamin Monoclonal Antibody	Q13561	Mouse	H,M	ELISA,WB,IHC
E-AB-27675	EBI3 Monoclonal Antibody	Q14213	Mouse	H,R	ELISA,WB
E-AB-27676	E-cadherin Monoclonal Antibody	P12830	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27454	ECHS1 Monoclonal Antibody	P30084	Mouse	H	ELISA,WB,IHC,IF
E-AB-27144	ECM1 Monoclonal Antibody	Q16610	Mouse	H	ELISA,WB,IHC
E-AB-27455	EEA1 Monoclonal Antibody	Q15075	Mouse	H	ELISA,IF
E-AB-27456	EEF1B2 Monoclonal Antibody	P24534	Mouse	H	ELISA,WB,IF
E-AB-27145	EEF1D Monoclonal Antibody	P29692	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27146	EEF1D Monoclonal Antibody	P29692	Mouse	H,M,R	ELISA,WB,IHC
E-AB-22043	EFHD1 Monoclonal Antibody	Q9BUP0	Mouse	H,M,R	WB,IF
E-AB-22007	EGFR Monoclonal Antibody	P00533	Mouse	H	WB,IF
E-AB-22075	EGFR Monoclonal Antibody	P00533	Mouse	H	WB
E-AB-27677	EGFR Monoclonal Antibody	P00533	Mouse	H	ELISA,WB,IHC
E-AB-27678	EIF2A Monoclonal Antibody	Q9BY44	Mouse	H	ELISA,WB
E-AB-27147	EIF3D Monoclonal Antibody	O15371	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27148	EIF3M Monoclonal Antibody	Q7L2H7	Mouse	H	ELISA,WB,IF
E-AB-22076	EIF4A1 Monoclonal Antibody	P60842	Mouse	H,M,R	WB,IF
E-AB-27599	EMILIN1 Monoclonal Antibody	A0A0C4DFX3	Mouse	H	ELISA,WB,IHC,IF
E-AB-22088	eNOS Monoclonal Antibody	P29474	Mouse	H,M,R	WB
E-AB-27149	ENTPD1 Monoclonal Antibody	P49961	Mouse	H	ELISA,WB,IF
E-AB-27679	EPCAM Monoclonal Antibody	P16422	Mouse	H	ELISA,WB,IHC,IF
E-AB-27600	EPHA1 Monoclonal Antibody	P21709	Mouse	H	ELISA,WB,IHC
E-AB-27567	EPLIN Monoclonal Antibody	Q9UHB6	Mouse	H	ELISA,WB,IHC,IF
E-AB-27568	ERAB Monoclonal Antibody	Q99714	Mouse	H,M	ELISA,WB
E-AB-22036	ERCC1 Monoclonal Antibody	P07992	Mouse	H	WB
E-AB-27457	ERCC1 Monoclonal Antibody	P07992	Mouse	H	ELISA,WB,IF
E-AB-27150	ERK1/2 Monoclonal Antibody	P27361	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27680	Erp72 Monoclonal Antibody	P13667	Mouse	H	ELISA,WB,IHC,IF
E-AB-27151	ESR2 Monoclonal Antibody	Q92731	Mouse	H,M	ELISA,WB
E-AB-27152	EXOSC2 Monoclonal Antibody	Q13868	Mouse	H	ELISA,WB,IHC,IF



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27681	EZH2 Monoclonal Antibody	Q15910	Mouse	H,M,R	ELISA,WB
E-AB-27458	FABP3 Monoclonal Antibody	P05413	Mouse	H	ELISA,WB
E-AB-27459	FABP5 Monoclonal Antibody	Q01469	Mouse	H	ELISA,WB,IHC
E-AB-27153	Factor XII Monoclonal Antibody	P00748	Mouse	H	ELISA,WB,IHC,IF
E-AB-27682	Factor XIIIa Monoclonal Antibody	P00488	Mouse	H,M	ELISA,WB,IHC
E-AB-27460	FAF1 Monoclonal Antibody	Q9UNN5	Mouse	H,M,R	ELISA,WB
E-AB-27461	FAK Monoclonal Antibody	Q05397	Mouse	H,M,R	ELISA,WB,IF
E-AB-27154	FAM3C Monoclonal Antibody	Q92520	Mouse	H	ELISA,WB,IHC,IF
E-AB-27155	FAS Monoclonal Antibody	P25445	Mouse	H,M	ELISA,WB,IHC
E-AB-27683	Fascin Monoclonal Antibody	Q16658	Mouse	H	ELISA,WB,IHC,IF
E-AB-27156	Fetuin-A Monoclonal Antibody	P02765	Mouse	H	ELISA,WB,IHC,IF
E-AB-27157	FGF12 Monoclonal Antibody	P61328	Mouse	H	ELISA,WB
E-AB-27462	FGF18 Monoclonal Antibody	O76093	Mouse	H,M	ELISA,WB
E-AB-27158	FGFR1 Monoclonal Antibody	P11362	Mouse	H,M	ELISA,WB,IHC
E-AB-22031	FH Monoclonal Antibody	P07954	Mouse	H,M,R	WB,IF
E-AB-22083	FH Fumarase Monoclonal Antibody	P07954	Mouse	H,M,R	WB,IF
E-AB-27159	FHIT Monoclonal Antibody	P49789	Mouse	H	ELISA,WB,IHC,IF
E-AB-27160	Fibrinogen beta chain Monoclonal Antibody	P02675	Mouse	H	ELISA,WB
E-AB-27161	Fibrinogen gamma chain Monoclonal Antibody	P02679	Mouse	H	ELISA,WB,IHC,IF
E-AB-27162	Fibromodulin Monoclonal Antibody	Q06828	Mouse	H,M	ELISA,WB,IHC
E-AB-22077	Fibronectin Monoclonal Antibody	P02751	Mouse	H,M,R	WB
E-AB-27163	Fibronectin Monoclonal Antibody	P02751	Mouse	H,M	ELISA,WB,IF
E-AB-27164	Fibulin 5 Monoclonal Antibody	Q9UBX5	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27165	FKBP2 Monoclonal Antibody	P26885	Mouse	H,R	ELISA,WB,IHC
E-AB-27166	FKBP52 Monoclonal Antibody	Q02790	Mouse	H	ELISA,WB,IHC,IF
E-AB-27167	FKBP52 Monoclonal Antibody	Q02790	Mouse	H	ELISA,WB,IHC,IF
E-AB-27601	FKBPL Monoclonal Antibody	Q9UIM3	Mouse	H,M	ELISA,WB
E-AB-22084	Flotillin-1 Monoclonal Antibody	P02751	Mouse	M,R,(H)	WB
E-AB-27569	Follistatin Monoclonal Antibody	P19883	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27168	FOLR1 Monoclonal Antibody	P15328	Mouse	H	ELISA,WB,IHC
E-AB-27686	FOXN1 Monoclonal Antibody	O15353	Mouse	H	ELISA,WB
E-AB-27687	FOXO1 Monoclonal Antibody	Q12778	Mouse	H	ELISA,WB



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27688	FOXO3A Monoclonal Antibody	O43524	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27169	FRS2 Monoclonal Antibody	Q8WU20	Mouse	H	ELISA,WB,IHC,IF
E-AB-27170	FUS Monoclonal Antibody	P35637	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27464	FUT8 Monoclonal Antibody	Q9BYC5	Mouse	H	ELISA,WB,IHC
E-AB-27465	FUT9 Monoclonal Antibody	Q9Y231	Mouse	H	ELISA,WB,IHC
E-AB-27689	G6PD Monoclonal Antibody	P11413	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27690	GABARAPL1 Monoclonal Antibody	Q9H0R8	Mouse	H,M,R	ELISA,WB
E-AB-27171	Galectin-1 Monoclonal Antibody	P09382	Mouse	H,M	ELISA,WB,IHC
E-AB-22006	Galectin-3 Monoclonal Antibody	P17931	Mouse	H	WB,IHC
E-AB-27172	Galectin-3 Monoclonal Antibody	P17931	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27173	Gamma cystathionase Monoclonal Antibody	P32929	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27692	GAMT Monoclonal Antibody	Q14353	Mouse	H,M,R	ELISA,WB
E-AB-27174	Gastrin Monoclonal Antibody	P01350	Mouse	H	ELISA,IHC
E-AB-27602	GATA1 Monoclonal Antibody	P15976	Mouse	H,M	ELISA,WB
E-AB-27695	GATA3 Monoclonal Antibody	P23771	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27696	GDF10 Monoclonal Antibody	P55107	Mouse	H,R	ELISA,WB
E-AB-27697	GDI1 Monoclonal Antibody	P31150	Mouse	H,M,R	ELISA,WB
E-AB-27175	GDI2 Monoclonal Antibody	P50395	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27176	GEFT Monoclonal Antibody	Q86VW2	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27177	Gelsolin Monoclonal Antibody	P06396	Mouse	H	ELISA,WB,IHC
E-AB-22022	GFAP Monoclonal Antibody	P14136	Mouse	M,R	WB
E-AB-27178	GFAP Monoclonal Antibody	P14136	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27179	GKN1 Monoclonal Antibody	Q9NS71	Mouse	H,R	ELISA,WB,IHC
E-AB-27470	GLS Monoclonal Antibody	O94925	Mouse	H	ELISA,WB,IF
E-AB-27180	GLUT1 Monoclonal Antibody	P11166	Mouse	H,M	ELISA,WB,IHC
E-AB-27700	Glutamine synthetase Monoclonal Antibody	P15104	Mouse	H	ELISA,WB,IHC
E-AB-27181	GM2A Monoclonal Antibody	P17900	Mouse	H	ELISA,WB
E-AB-27603	GMF-beta Monoclonal Antibody	P60983	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27701	GNPAT Monoclonal Antibody	O15228	Mouse	H	ELISA,WB
E-AB-27182	GOSR2 Monoclonal Antibody	O14653	Mouse	H	ELISA,WB,IHC
E-AB-27471	GOT1 Monoclonal Antibody	P17174	Mouse	H	ELISA,WB,IHC,IF
E-AB-27702	GP73 Monoclonal Antibody	Q8NBJ4	Mouse	H	ELISA,WB,IHC,IF



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27183	GRP94 Monoclonal Antibody	P14625	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27703	GSDMD Monoclonal Antibody	P57764	Mouse	H	ELISA,WB,IF
E-AB-27184	GSNOR Monoclonal Antibody	P11766	Mouse	H	ELISA,WB,IHC,IF
E-AB-27473	GSTO1 Monoclonal Antibody	P78417	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27185	HADHA Monoclonal Antibody	P40939	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27188	HBE1 Monoclonal Antibody	P02100	Mouse	H	ELISA,WB
E-AB-27474	HBG1 Monoclonal Antibody	P69891	Mouse	H	ELISA,IHC
E-AB-27186	HDAC1 Monoclonal Antibody	Q13547	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27187	HDGF Monoclonal Antibody	P51858	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27707	Hemopexin Monoclonal Antibody	P02790	Mouse	H,R	ELISA,WB
E-AB-22008	HER2 Monoclonal Antibody	P04626	Mouse	H,M,R	WB,IHC
E-AB-22037	HER2 Monoclonal Antibody	P04626	Mouse	H	WB
E-AB-27570	HER2 Monoclonal Antibody	P04626	Mouse	H	ELISA,IHC
E-AB-27475	HEXA Monoclonal Antibody	P06865	Mouse	H	ELISA,WB
E-AB-27476	HEXIM1 Monoclonal Antibody	O94992	Mouse	H	ELISA,WB
E-AB-22086	Histone H2B Monoclonal Antibody	P33778	Mouse	H,M,R	WB
E-AB-22003	Histone H3 Monoclonal Antibody	P68431	Mouse	H,M,R,Y	WB,IF,IP
E-AB-22070	Histone H3 Monoclonal Antibody	P68431	Mouse	Z	WB
E-AB-22079	Histone H3 Monoclonal Antibody	P68431	Mouse	H,M,R	WB
E-AB-22081	Histone H3 Monoclonal Antibody	P68431	Mouse	H,M,R	WB
E-AB-27189	HLA class I (HLA-A) Monoclonal Antibody	P01892	Mouse	H	ELISA,WB,IHC
E-AB-27708	HLA-G Monoclonal Antibody	P17693	Mouse	H	ELISA,WB,IHC
E-AB-27478	HP Monoclonal Antibody	P00738	Mouse	H	ELISA,WB
E-AB-27479	HPSE Monoclonal Antibody	Q9Y251	Mouse	H	ELISA,WB,IHC,IF
E-AB-27709	HRPT2; CDC73 Monoclonal Antibody	Q6P1J9	Mouse	H,M,R	ELISA,WB
E-AB-27710	Hsc70 Monoclonal Antibody	P11142	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27190	HSP60 Monoclonal Antibody	P10809	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-22005	HSP70 Monoclonal Antibody	P34931	Mouse	H,M,R	WB,IF
E-AB-27481	HSP70 Monoclonal Antibody	P34931	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27191	HSP90 Monoclonal Antibody	P07900	Mouse	H,M,R	ELISA,WB,IHC
E-AB-22072	HSP90 β Monoclonal Antibody	P08238	Mouse	H,M,R	WB
E-AB-27192	HSPA2 Monoclonal Antibody	P54652	Mouse	H,M,R	ELISA,WB,IHC,IF



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27711	HVCN1 Monoclonal Antibody	Q96D96	Mouse	H	ELISA,WB,IHC
E-AB-27194	IA-2 Monoclonal Antibody	Q16849	Mouse	H	ELISA,WB,IHC,IF
E-AB-27195	ICAM-1 Monoclonal Antibody	P05362	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-22041	IDE Monoclonal Antibody	P14735	Mouse	H	WB
E-AB-27196	IDH1 Monoclonal Antibody	O75874	Mouse	H	ELISA,WB,IHC,IF
E-AB-27197	Iduronate 2 sulfatase Monoclonal Antibody	P22304	Mouse	H	ELISA,WB,IHC,IF
E-AB-27198	IFITM1 Monoclonal Antibody	P13164	Mouse	H	ELISA,WB,IHC,IF
E-AB-27200	IFITM2/3 Monoclonal Antibody	Q01628	Mouse	H	WB,IHC,IF
E-AB-27199	IFITM2 Monoclonal Antibody	Q01629	Mouse	H	ELISA,WB
E-AB-27484	IFT88 Monoclonal Antibody	Q13099	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27201	IGBP1 Monoclonal Antibody	P78318	Mouse	H,M	ELISA,WB,IHC
E-AB-27202	IGF1R beta chain Monoclonal Antibody	P08069	Mouse	H,R	ELISA,WB,IHC
E-AB-27712	IGHG4 Monoclonal Antibody	P01861	Mouse	H	ELISA,WB
E-AB-27713	IGHM Monoclonal Antibody	P01871	Mouse	H	ELISA,WB
E-AB-27714	IGSF2 Monoclonal Antibody	Q93033	Mouse	H,M	ELISA,WB,IHC
E-AB-27715	IkB alpha Monoclonal Antibody	P25963	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27717	IKBKG Monoclonal Antibody	Q9Y6K9	Mouse	H	ELISA,WB
E-AB-27203	IL10 Monoclonal Antibody	P22301	Mouse	H	ELISA,WB,IHC
E-AB-27204	IL12A Monoclonal Antibody	P29459	Mouse	H	ELISA,WB,IHC
E-AB-27487	IL12B Monoclonal Antibody	P29460	Mouse	H	ELISA,WB
E-AB-27571	IL16 Monoclonal Antibody	Q14005	Mouse	H	ELISA,WB
E-AB-27572	IL17 Monoclonal Antibody	Q16552	Mouse	H	ELISA,WB
E-AB-27718	IL17A Monoclonal Antibody	Q62386	Mouse	H	ELISA,WB
E-AB-27573	IL17D Monoclonal Antibody	Q8TAD2	Mouse	H	ELISA,WB
E-AB-27574	IL17F Monoclonal Antibody	Q96PD4	Mouse	H,M	ELISA,WB
E-AB-27205	IL19 Monoclonal Antibody	Q9UHD0	Mouse	H	ELISA,WB,IHC
E-AB-27489	IL1F10 Monoclonal Antibody	Q8WWZ1	Mouse	H	ELISA,WB
E-AB-27491	IL1F8 Monoclonal Antibody	Q9NZH7	Mouse	H	ELISA,WB,IHC
E-AB-27206	IL1R2 Monoclonal Antibody	P27930	Mouse	H,M	ELISA,WB,IF
E-AB-27207	IL2 Monoclonal Antibody	P60568	Mouse	H,M	ELISA,WB,IHC
E-AB-27208	IL20RB Monoclonal Antibody	Q6UXL0	Mouse	H	ELISA,WB
E-AB-27719	IL21R Monoclonal Antibody	Q9HBE5	Mouse	H	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27209	IL22RA2 Monoclonal Antibody	Q969J5	Mouse	H	ELISA,WB,IHC
E-AB-27492	IL23A Monoclonal Antibody	Q9NPF7	Mouse	H	ELISA,WB,IHC
E-AB-27575	IL27 Monoclonal Antibody	Q8NEV9	Mouse	H	ELISA,WB,IHC
E-AB-27210	IL28A Monoclonal Antibody	8IZJ0	Mouse	H	ELISA,WB,IHC
E-AB-27211	IL29 Monoclonal Antibody	Q8IU54	Mouse	H	ELISA,WB
E-AB-27495	IL33 Monoclonal Antibody	O95760	Mouse	H	ELISA,WB,IHC,IF
E-AB-27490	IL36RN Monoclonal Antibody	A0A024R518	Mouse	H	ELISA,WB
E-AB-27496	IL37 Monoclonal Antibody	Q9NZH6	Mouse	H	ELISA,WB,IHC
E-AB-27212	IL4 Monoclonal Antibody	P05112	Mouse	H	ELISA,WB,IHC
E-AB-27213	IL6 Monoclonal Antibody	P05231	Mouse	H	ELISA,WB,IHC
E-AB-27214	IL7 Monoclonal Antibody	P13232	Mouse	H	ELISA,WB,IHC,IF
E-AB-27215	IL9 Monoclonal Antibody	P15248	Mouse	H	ELISA,WB,IHC
E-AB-27498	IMP3 Monoclonal Antibody	Q9NV31	Mouse	H	ELISA,WB,IHC,IF
E-AB-27216	INHBA Monoclonal Antibody	P08476	Mouse	H	ELISA,WB,IHC
E-AB-27217	Inhibin beta A Monoclonal Antibody	P08476	Mouse	H	ELISA,WB
E-AB-27218	Insulin Monoclonal Antibody	P01308	Mouse	H,M	ELISA,IHC
E-AB-27576	Integrin alpha-3 Monoclonal Antibody	P26006	Mouse	H	ELISA,WB,IF
E-AB-27720	Integrin beta-1 Monoclonal Antibody	P05556	Mouse	H	ELISA,WB
E-AB-27219	Interferon alpha Monoclonal Antibody	P01562	Mouse	H	ELISA,WB
E-AB-27721	IRF4 Monoclonal Antibody	Q15306	Mouse	H,R	ELISA,WB,IHC
E-AB-27220	IRF5 Monoclonal Antibody	Q13568	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27499	ITM2C Monoclonal Antibody	Q9NQX7	Mouse	H	ELISA,IF
E-AB-27500	IVD Monoclonal Antibody	P26440	Mouse	H	ELISA,WB,IHC,IF
E-AB-27722	JAK1 Monoclonal Antibody	P23458	Mouse	H,R	ELISA,WB
E-AB-27221	JNK Monoclonal Antibody	P45983	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27723	JUN Monoclonal Antibody	P05412	Mouse	H,M,R	ELISA,WB
E-AB-27724	KChIP1 Monoclonal Antibody	Q9NZI2	Mouse	H	ELISA,WB,IHC
E-AB-27607	KCNIP4 Monoclonal Antibody	Q6PIL6	Mouse	H,R	ELISA,WB,IHC
E-AB-27504	KCNMB4 Monoclonal Antibody	Q86W47	Mouse	H,M	ELISA,WB,IHC
E-AB-27222	KCNN4 Monoclonal Antibody	O15554	Mouse	H,M,R	ELISA,WB
E-AB-27223	KEAP1 Monoclonal Antibody	Q14145	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-22027	Ki 67 Monoclonal Antibody	P46013	Mouse	H,M,R	IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-22026	Kif 7 Monoclonal Antibody	Q2M1P5	Mouse	H,M,R	IHC
E-AB-27224	KIFAP3 Monoclonal Antibody	Q92845	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27225	Kininogen 1 Monoclonal Antibody	P01042	Mouse	H	ELISA,WB,IHC,IF
E-AB-27226	KLF15 Monoclonal Antibody	Q9UIH9	Mouse	H	ELISA,WB
E-AB-27505	KLK3 Monoclonal Antibody	P07288	Mouse	H	ELISA,IHC,IF
E-AB-27228	KMO Monoclonal Antibody	O15229	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27506	L2HGDH Monoclonal Antibody	Q9H9P8	Mouse	H	ELISA,WB,IHC
E-AB-27507	LAMP2 Monoclonal Antibody	P13473	Mouse	H	ELISA,WB,IHC,IF
E-AB-27277	LAP3 Monoclonal Antibody	P28838	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27229	LBP Monoclonal Antibody	P18428	Mouse	H	ELISA,WB,IHC
E-AB-27230	LC3A/B Monoclonal Antibody	Q9H492	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27231	LCK Monoclonal Antibody	P06239	Mouse	H	ELISA,WB
E-AB-27509	LDHA Monoclonal Antibody	P00338	Mouse	H,M,Pig	ELISA,WB,IHC,IF
E-AB-27728	LDHB Monoclonal Antibody	P07195	Mouse	H,M,R	ELISA,WB
E-AB-27729	LDLR Monoclonal Antibody	P01130	Mouse	H	ELISA,WB,IHC
E-AB-27232	LGALS3BP Monoclonal Antibody	Q08380	Mouse	H	ELISA,WB,IHC,IF
E-AB-27233	LIN28 Monoclonal Antibody	Q9H9Z2	Mouse	H	ELISA,WB,IHC
E-AB-27234	LIPH Monoclonal Antibody	Q8WWY8	Mouse	H	ELISA,WB,IHC,IF
E-AB-27236	LONP1 Monoclonal Antibody	P36776	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27237	LPCAT1 Monoclonal Antibody	Q3TFD2	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27238	LYN Monoclonal Antibody	P07948	Mouse	H,M	ELISA,WB,IHC
E-AB-27730	Lysozyme Monoclonal Antibody	P61626	Mouse	H,M,R	ELISA,WB
E-AB-27239	MAD2L1 Monoclonal Antibody	Q13257	Mouse	H	ELISA,WB,IHC
E-AB-27510	MAGEA3 Monoclonal Antibody	P43357	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27731	MALT1 Monoclonal Antibody	Q9UDY8	Mouse	H	ELISA,WB,IHC,IF
E-AB-27240	Mammaglobin A Monoclonal Antibody	Q13296	Mouse	H	ELISA,IHC
E-AB-27511	MAOB Monoclonal Antibody	P27338	Mouse	H	ELISA,WB,IHC
E-AB-22030	MAP2 Monoclonal Antibody	P11137	Mouse	H,M,R	IHC
E-AB-27512	MAPK13 Monoclonal Antibody	O15264	Mouse	H,M	ELISA,WB
E-AB-27513	MBIP Monoclonal Antibody	Q9NS73	Mouse	H	ELISA,WB,IHC
E-AB-27241	MCL1 Monoclonal Antibody	Q07820	Mouse	H	ELISA,WB,IHC,IF
E-AB-27514	MCM2 Monoclonal Antibody	P49736	Mouse	H	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27515	Mcp1 Monoclonal Antibody	P13500	Mouse	M	ELISA,WB,IHC
E-AB-27242	MCRS1 Monoclonal Antibody	Q96EZ8	Mouse	H,M,R	ELISA,WB
E-AB-27733	MEA1 Monoclonal Antibody	Q16626	Mouse	H,R	ELISA,WB
E-AB-27243	Melanophilin Monoclonal Antibody	Q9BV36	Mouse	H	ELISA,WB,IHC,IF
E-AB-27244	MEST Monoclonal Antibody	Q5EB52	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27245	MGEA5 Monoclonal Antibody	O60502	Mouse	H	ELISA,WB,IHC,IF
E-AB-27734	MGP Monoclonal Antibody	P08493	Mouse	H,M	ELISA,WB,IHC
E-AB-27735	MICA Monoclonal Antibody	Q29983	Mouse	H	ELISA,WB,IHC
E-AB-27736	MKP-2 Monoclonal Antibody	Q13115	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27516	MLANA Monoclonal Antibody	Q16655	Mouse	H	ELISA,IHC
E-AB-27246	MME Monoclonal Antibody	P08473	Mouse	H	ELISA,WB,IHC
E-AB-27737	MMP2 Monoclonal Antibody	P08253	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27738	MMP3 Monoclonal Antibody	P08254	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27247	MMS19 Monoclonal Antibody	Q96T76	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27248	MPO Monoclonal Antibody	P05164	Mouse	H	ELISA,WB,IHC
E-AB-27249	MSH2 Monoclonal Antibody	P43246	Mouse	H	ELISA,WB,IHC
E-AB-27250	MSH6 Monoclonal Antibody	P52701	Mouse	H	ELISA,WB,IHC,IF
E-AB-27739	MSLN Monoclonal Antibody	Q13421	Mouse	H	ELISA,WB,IHC
E-AB-27251	MTA2 Monoclonal Antibody	O94776	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27252	MUC16 Monoclonal Antibody	Q8WX17	Mouse	H	ELISA,IHC
E-AB-27740	MUM1 Monoclonal Antibody	Q2TAK8	Mouse	H	ELISA,WB
E-AB-27253	MYH9 Monoclonal Antibody	P35579	Mouse	H,M	ELISA,WB
E-AB-27254	MYL3 Monoclonal Antibody	P08590	Mouse	H,R	ELISA,WB,IHC
E-AB-27743	MYOC Monoclonal Antibody	Q99972	Mouse	H,M,R	ELISA,WB
E-AB-27517	MYOD1 Monoclonal Antibody	P15172	Mouse	H	ELISA,WB
E-AB-27255	Myoglobin Monoclonal Antibody	P02144	Mouse	H,R	ELISA,WB,IHC
E-AB-27256	Myosin 2a Monoclonal Antibody	Q9UKX2	Mouse	H	ELISA,WB,IHC
E-AB-22021	Myosin Heavy Chain Monoclonal Antibody	Q7Z406	Mouse	H,M,R,Fruit Fly,Nematode	IF
E-AB-27257	Myosin Light Chain 2 Monoclonal Antibody	P10916	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27744	NAMPT Monoclonal Antibody	P43490	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27258	NAPRT1 Monoclonal Antibody	Q6XQN6	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27259	Napsin A Monoclonal Antibody	O96009	Mouse	H,M	ELISA,IHC,IF



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27260	N-cadherin Monoclonal Antibody	P19022	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27261	NCALD Monoclonal Antibody	P61601	Mouse	H	ELISA,WB,IHC
E-AB-27577	NCAM1 Monoclonal Antibody	P13591	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27262	NDUFA4L2 Monoclonal Antibody	Q9NRX3	Mouse	H	ELISA,WB,IHC,IF
E-AB-27263	NDUFS5 Monoclonal Antibody	O43920	Mouse	H	ELISA,WB,IHC
E-AB-27745	NEFM Monoclonal Antibody	P07197	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27519	NENF Monoclonal Antibody	Q9UMX5	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27520	NESP55 Monoclonal Antibody	O95467	Mouse	H	ELISA,WB,IHC
E-AB-27264	Nestin Monoclonal Antibody	P48681	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27265	Neuropilin 1 Monoclonal Antibody	O14786	Mouse	H	ELISA,WB,IHC,IF
E-AB-27266	NF200 Monoclonal Antibody	P12036	Mouse	H,Pig	ELISA,WB,IHC,IF
E-AB-27746	NFE2 Monoclonal Antibody	Q16621	Mouse	H	ELISA,WB
E-AB-27747	NFIL3 Monoclonal Antibody	Q16649	Mouse	H	ELISA,WB
E-AB-22016	NFkB p65 Monoclonal Antibody	Q04206	Mouse	H,M,R	WB,IHC,IP
E-AB-22066	NFkB p65 Monoclonal Antibody	Q04206	Mouse	H,M,R	WB,IP
E-AB-27267	NF-L Monoclonal Antibody	P07196	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27268	NINJ2 Monoclonal Antibody	Q9NZG7	Mouse	H	ELISA,WB,IHC
E-AB-27269	NOB1 Monoclonal Antibody	Q9ULX3	Mouse	H	ELISA,WB,IHC
E-AB-27748	NONO Monoclonal Antibody	Q15233	Mouse	H,M,R	ELISA,WB
E-AB-27521	NR1H2 Monoclonal Antibody	P55055	Mouse	H	ELISA,WB
E-AB-27270	NR1H3 Monoclonal Antibody	Q13133	Mouse	H,R	ELISA,WB,IF
E-AB-27749	NR2E3 Monoclonal Antibody	Q9Y5X4	Mouse	H	ELISA,WB
E-AB-27271	NR2F6 Monoclonal Antibody	P10588	Mouse	H	ELISA,WB,IHC
E-AB-27272	NR2F6 Monoclonal Antibody	P10588	Mouse	H,M	ELISA,WB,IHC
E-AB-27750	NRG1 Monoclonal Antibody	Q96IB3	Mouse	H,M,R	ELISA,WB
E-AB-22029	NSE Monoclonal Antibody	P09104	Mouse	H,M,R	WB,IHC
E-AB-27273	NSE Monoclonal Antibody	P09104	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27274	NTF2 Monoclonal Antibody	P61970	Mouse	H	ELISA,WB,IHC,IF
E-AB-27751	NUDT21 Monoclonal Antibody	O43809	Mouse	H,M	ELISA,WB,IF
E-AB-27275	NUMB Monoclonal Antibody	P49757	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27276	NUMBL Monoclonal Antibody	Q9Y6R0	Mouse	H	ELISA,WB,IHC
E-AB-27752	NUP155 Monoclonal Antibody	O75694	Mouse	H,R	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27753	Occludin Monoclonal Antibody	Q16625	Mouse	H,M,R	ELISA,WB,IHC
E-AB-22047	Oct1 Monoclonal Antibody	P14859	Mouse	H	WB
E-AB-27277	OCT4 Monoclonal Antibody	Q01860	Mouse	H,M	ELISA,WB,IHC
E-AB-27522	OGDH Monoclonal Antibody	Q02218	Mouse	H,R	ELISA,WB,IHC,IF
E-AB-27754	OGN Monoclonal Antibody	P20774	Mouse	H	ELISA,WB,IHC
E-AB-27755	OIP5 Monoclonal Antibody	O43482	Mouse	H	ELISA,WB
E-AB-27278	OPTN Monoclonal Antibody	Q96CV9	Mouse	H	ELISA,WB,IHC
E-AB-27279	ORAI1 Monoclonal Antibody	Q96D31	Mouse	H	ELISA,WB,IHC
E-AB-27280	ORM1/2 Monoclonal Antibody	P02763	Mouse	H	ELISA,WB,IHC,IF
E-AB-27281	ORM2 Monoclonal Antibody	P19652	Mouse	H	ELISA,WB,IF
E-AB-27756	OTUD7B Monoclonal Antibody	Q6GQQ9	Mouse	H	ELISA,WB
E-AB-27524	OXA1L Monoclonal Antibody	Q15070	Mouse	H	ELISA,WB,IHC
E-AB-27578	p120 Catenin Monoclonal Antibody	O60716	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27282	P21 Monoclonal Antibody	P38936	Mouse	H,M	ELISA,WB
E-AB-27757	P2RX4 Monoclonal Antibody	Q99571	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27758	P2RX5 Monoclonal Antibody	Q93086	Mouse	H,M,R	ELISA,WB
E-AB-27283	p38 MAPK Monoclonal Antibody	Q16539	Mouse	H	ELISA,WB,IHC
E-AB-27284	P4HA1 Monoclonal Antibody	P13674	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27285	p504S Monoclonal Antibody	Q9UHK6	Mouse	H,M	ELISA,WB
E-AB-22015	p53 Monoclonal Antibody	P04637	Mouse	H	WB,IHC,IF
E-AB-27286	p53 Monoclonal Antibody	P04637	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27287	p62 Monoclonal Antibody	Q13501	Mouse	H	ELISA,WB,IHC,IF
E-AB-27525	PA2G4 Monoclonal Antibody	Q9UQ80	Mouse	H	ELISA,WB,IHC,IF
E-AB-27759	PADI2 Monoclonal Antibody	Q9Y2J8	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27526	PAI-1 Monoclonal Antibody	P05121	Mouse	H	ELISA,WB,IHC,IF
E-AB-27288	PAIP1 Monoclonal Antibody	Q9H074	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27227	pan Ras Monoclonal Antibody	P01116	Mouse	H,M,R	ELISA,WB,IHC
E-AB-22074	Cleaved-PARP Monoclonal Antibody	P09874	Mouse	H	WB
E-AB-22078	Cleaved-PARP Monoclonal Antibody	P09874	Mouse	H	WB
E-AB-22073	PARP Monoclonal Antibody	P09874	Mouse	H	WB
E-AB-27290	PAX1 Monoclonal Antibody	P15863	Mouse	H,M	ELISA,WB
E-AB-27291	PAX8 Monoclonal Antibody	Q06710	Mouse	H	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27292	PAX8 Monoclonal Antibody	Q06710	Mouse	H	ELISA,WB,IHC
E-AB-27436	P-cadherin Monoclonal Antibody	P22223	Mouse	H	ELISA,WB,IHC
E-AB-27527	PCMT1 Monoclonal Antibody	P22061	Mouse	H	ELISA,WB,IF
E-AB-22001	PCNA Monoclonal Antibody	P12004	Mouse	H,M,R	WB,IHC
E-AB-27293	PCNA Monoclonal Antibody	P12004	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27294	PD-1 Monoclonal Antibody	Q15116	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27760	PDCD10 Monoclonal Antibody	Q9BUL8	Mouse	H	ELISA,WB,IHC
E-AB-27295	PDCD4 Monoclonal Antibody	Q53EL6	Mouse	H	ELISA,WB,IHC
E-AB-27604	PDGFRA Monoclonal Antibody	P16234	Mouse	H	ELISA,WB,IHC
E-AB-27296	PDH E1a Monoclonal Antibody	P08559	Mouse	H	ELISA,WB,IHC,IF
E-AB-27761	PDI Monoclonal Antibody	P07237	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27762	PDIA3 Monoclonal Antibody	P30101	Mouse	H,M,R	ELISA,WB
E-AB-27297	PD-L1 Monoclonal Antibody	Q9NZQ7	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27763	PEBP1 Monoclonal Antibody	P30086	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27764	Periostin Monoclonal Antibody	Q15063	Mouse	H,M,R	ELISA,WB
E-AB-27765	Peripherin Monoclonal Antibody	P41219	Mouse	H,R	ELISA,WB,IHC
E-AB-22067	Peroxiredoxin 1 Monoclonal Antibody	Q06830	Mouse	H,M,R	WB,IF
E-AB-27298	Peroxiredoxin 2 Monoclonal Antibody	P32119	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27299	PEX5 Monoclonal Antibody	P50542	Mouse	H,R	ELISA,WB,IHC
E-AB-27300	PFN2 Monoclonal Antibody	P35080	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27301	PGM1 Monoclonal Antibody	P36871	Mouse	H	ELISA,WB,IHC,IF
E-AB-27766	PGRMC1 Monoclonal Antibody	O00264	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27302	PGRMC2 Monoclonal Antibody	O15173	Mouse	H,M	ELISA,WB,IHC
E-AB-27767	PHF10 Monoclonal Antibody	Q8WUB8	Mouse	H,M	ELISA,WB,IHC
E-AB-27529	PHOX2B Monoclonal Antibody	Q99453	Mouse	H,M	ELISA,WB,IHC
E-AB-27304	PI3K p85(alpha) Monoclonal Antibody	P27986	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27305	Pikachurin Monoclonal Antibody	Q63HQ2	Mouse	H	ELISA,WB
E-AB-27769	PKC alpha Monoclonal Antibody	P17252	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27770	PKC gamma Monoclonal Antibody	P05129	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27771	PKC iota Monoclonal Antibody	P41743	Mouse	H,M,R	ELISA,WB
E-AB-27306	PKM2 Monoclonal Antibody	P14618	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27772	PLA2G1B Monoclonal Antibody	Q9BS22	Mouse	H	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27307	PLAA Monoclonal Antibody	Q9Y263	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27773	Plakoglobin Monoclonal Antibody	P14923	Mouse	H,R	ELISA,WB,IHC
E-AB-27308	PLAP Monoclonal Antibody	P05187	Mouse	H	ELISA,WB,IHC
E-AB-27774	PLEK Monoclonal Antibody	P08567	Mouse	H	ELISA,WB,IHC
E-AB-27775	PLG Monoclonal Antibody	P00747	Mouse	H	ELISA,WB,IHC
E-AB-27776	PLOD2 Monoclonal Antibody	O00469	Mouse	H	ELISA,WB,IHC
E-AB-27309	PLOD3 Monoclonal Antibody	O60568	Mouse	H	ELISA,WB,IHC,IF
E-AB-27531	PMS2 Monoclonal Antibody	P54278	Mouse	H	ELISA,WB,IHC
E-AB-27777	POMC Monoclonal Antibody	P01189	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27310	PPAR gamma Monoclonal Antibody	P37231	Mouse	H	ELISA,WB,IHC
E-AB-27311	PPARD Monoclonal Antibody	Q03181	Mouse	H,M	ELISA,WB
E-AB-27778	PPARGC1A Monoclonal Antibody	Q9UBK2	Mouse	H,M,R	ELISA,WB,IF
E-AB-27312	PR Monoclonal Antibody	P06401	Mouse	H	ELISA,IHC
E-AB-27313	PRDX4 Monoclonal Antibody	Q13162	Mouse	H	ELISA,WB,IHC,IF
E-AB-27314	Prealbumin Monoclonal Antibody	P02766	Mouse	H,M	ELISA,WB,IHC
E-AB-27315	Pre-IL18 Monoclonal Antibody	Q14116	Mouse	H	ELISA,WB,IF
E-AB-27316	Prohibitin Monoclonal Antibody	P35232	Mouse	H,R	ELISA,WB,IHC,IF
E-AB-27779	Prohibitin 2 Monoclonal Antibody	Q99623	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27579	Prostein Monoclonal Antibody	Q96JT2	Mouse	H	ELISA,WB,IHC
E-AB-27317	Protein C inhibitor Monoclonal Antibody	P05154	Mouse	H	ELISA,WB,IHC
E-AB-27780	PSAP Monoclonal Antibody	P07602	Mouse	H	ELISA,WB,IHC
E-AB-27318	P-selectin Monoclonal Antibody	P16109	Mouse	H	ELISA,WB,IHC
E-AB-27781	PSMB10 Monoclonal Antibody	P40306	Mouse	H	ELISA,WB
E-AB-27319	PSMD4 Monoclonal Antibody	P55036	Mouse	H	ELISA,WB,IHC,IF
E-AB-27320	PTEN Monoclonal Antibody	P60484	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27321	RAB18 Monoclonal Antibody	Q9NP72	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27322	RAB27A Monoclonal Antibody	P51159	Mouse	H	ELISA,WB,IHC,IF
E-AB-27782	RAB5A Monoclonal Antibody	P20339	Mouse	H	ELISA,WB,IHC,IF
E-AB-27323	RABEP2 Monoclonal Antibody	Q9H5N1	Mouse	H	ELISA,WB,IHC
E-AB-27324	RAC1 Monoclonal Antibody	P63000	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27325	RAC2 Monoclonal Antibody	P15153	Mouse	H,M	ELISA,WB,IHC
E-AB-27535	RACGAP1 Monoclonal Antibody	Q9H0H5	Mouse	H,M,R	ELISA,WB,IHC,IF



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27326	RBAP48 Monoclonal Antibody	Q09028	Mouse	H	ELISA,WB,IHC
E-AB-27327	RBBP9 Monoclonal Antibody	O75884	Mouse	H	ELISA,WB,IHC,IF
E-AB-27328	RBBP9 Monoclonal Antibody	O75884	Mouse	H	ELISA,WB,IHC
E-AB-27329	RBM15 Monoclonal Antibody	Q96T37	Mouse	H	ELISA,WB,IF
E-AB-27330	RBM4 Monoclonal Antibody	Q9BWF3	Mouse	H	ELISA,WB,IHC,IF
E-AB-27331	RBP4 Monoclonal Antibody	P02753	Mouse	H	ELISA,WB,IHC
E-AB-27538	RBPJ Monoclonal Antibody	Q06330	Mouse	H	ELISA,WB
E-AB-27539	RCAS1 Monoclonal Antibody	O00559	Mouse	H	ELISA,WB,IHC
E-AB-27332	RENALASE Monoclonal Antibody	Q5VYX0	Mouse	H	ELISA,WB,IHC
E-AB-27783	RGS4 Monoclonal Antibody	P49798	Mouse	H,M,R	ELISA,WB
E-AB-27333	RISC Monoclonal Antibody	Q9HB40	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27540	RND3 Monoclonal Antibody	P61587	Mouse	H	ELISA,WB
E-AB-27334	RNH1 Monoclonal Antibody	P13489	Mouse	H	ELISA,WB
E-AB-27335	RPL3 Monoclonal Antibody	P39023	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27541	RPS3 Monoclonal Antibody	P23396	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27336	RRM1 Monoclonal Antibody	P23921	Mouse	H	ELISA,WB,IHC,IF
E-AB-27542	RRM1 Monoclonal Antibody	P23921	Mouse	H,M	ELISA,WB,IHC
E-AB-27337	RSPO3 Monoclonal Antibody	Q9BXY4	Mouse	H	ELISA,WB,IHC,IF
E-AB-27338	RXRA Monoclonal Antibody	P19793	Mouse	H,M	ELISA,WB
E-AB-27784	RXRG Monoclonal Antibody	P48443	Mouse	H	ELISA,WB
E-AB-27339	S100A10 Monoclonal Antibody	P60903	Mouse	H	ELISA,WB,IHC
E-AB-27340	S100A11 Monoclonal Antibody	P31949	Mouse	H,M	ELISA,WB,IHC
E-AB-27786	S100A4 Monoclonal Antibody	P26447	Mouse	H,M	ELISA,WB,IHC
E-AB-27341	S100A6 Monoclonal Antibody	P06703	Mouse	H	ELISA,WB,IHC,IF
E-AB-27342	SALL4 Monoclonal Antibody	Q9UJQ4	Mouse	H	ELISA,WB
E-AB-27543	SDCCAG8 Monoclonal Antibody	Q86SQ7	Mouse	H	ELISA,WB,IHC
E-AB-27343	SEC5 Monoclonal Antibody	Q96KP1	Mouse	H,M	ELISA,WB,IF
E-AB-27344	SECTM1 Monoclonal Antibody	Q8WVN6	Mouse	H	ELISA,WB,IHC
E-AB-27345	SEPT2 Monoclonal Antibody	Q15019	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27346	SEPT8 Monoclonal Antibody	Q92599	Mouse	H	ELISA,WB
E-AB-27544	SERPINA10 Monoclonal Antibody	Q9UK55	Mouse	H,R	ELISA,WB,IHC,IF
E-AB-27787	SERPINA7 Monoclonal Antibody	P05543	Mouse	H	ELISA,WB



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27347	SERPINE2 Monoclonal Antibody	P07093	Mouse	H	ELISA,WB,IHC,IF
E-AB-27348	Serum amyloid P component Monoclonal Antibody	P02743	Mouse	H	ELISA,WB,IHC
E-AB-27545	Sestrin2 Monoclonal Antibody	P58004	Mouse	H,R	ELISA,WB,IHC,IF
E-AB-27349	SETDB1 Monoclonal Antibody	Q15047	Mouse	H	ELISA,WB
E-AB-27788	SFRP2 Monoclonal Antibody	Q96HF1	Mouse	H,R	ELISA,WB,IHC,IF
E-AB-27546	SGTA Monoclonal Antibody	O43765	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27350	SIRT1 Monoclonal Antibody	Q96EB6	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27789	SIRT2 Monoclonal Antibody	Q8IXJ6	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27790	SIX2 Monoclonal Antibody	Q9NPC8	Mouse	H,R	ELISA,WB
E-AB-27791	SKAP2 Monoclonal Antibody	O75563	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27351	SLAM Monoclonal Antibody	Q13291	mouse	H	ELISA,WB
E-AB-27792	SLC01B3 Monoclonal Antibody	Q0VGL9	Mouse	H	ELISA,WB,IHC
E-AB-27793	SLP76 Monoclonal Antibody	Q13094	Mouse	H	ELISA,WB
E-AB-27794	SMAD7 Monoclonal Antibody	O15105	Mouse	H,M,R	ELISA,WB
E-AB-27352	SMARCD1 Monoclonal Antibody	Q96GM5	Mouse	H	ELISA,WB
E-AB-27353	SMARCE1 Monoclonal Antibody	Q969G3	Mouse	H	ELISA,WB,IHC,IF
E-AB-27580	SMMHC Monoclonal Antibody	P35749	Mouse	H,M	ELISA,WB
E-AB-27581	SMMHC Monoclonal Antibody	P35749	Mouse	H	ELISA,WB
E-AB-27354	SMN Monoclonal Antibody	Q16637	Mouse	H,M,R	ELISA,WB
E-AB-27355	SMN Monoclonal Antibody	Q16637	Mouse	H	ELISA,WB,IHC,IF
E-AB-27356	SMN-Exon7 Monoclonal Antibody	Q16637	Mouse	H,M	ELISA,WB,IF
E-AB-27795	SMYD3 Monoclonal Antibody	Q9H7B4	Mouse	H	ELISA,WB,IHC
E-AB-27357	SNAP25 Monoclonal Antibody	P60880	Mouse	H,M,R	ELISA,WB,IF
E-AB-27358	SND1 Monoclonal Antibody	Q7KZF4	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27359	SNRPD2 Monoclonal Antibody	P62316	Mouse	H	ELISA,WB
E-AB-27796	SOD2 Monoclonal Antibody	P04179	Mouse	H,M,R	ELISA,WB
E-AB-27797	SOX2 Monoclonal Antibody	P48431	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27798	SPARC Monoclonal Antibody	P09486	Mouse	H	ELISA,WB
E-AB-27799	SRD5A1 Monoclonal Antibody	P18405	Mouse	H	ELISA,WB,IHC,IF
E-AB-27360	SRP9 Monoclonal Antibody	P49458	Mouse	H	ELISA,WB
E-AB-27361	SRPX2 Monoclonal Antibody	O60687	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27362	ST2 Monoclonal Antibody	Q01638	Mouse	H,M,R	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27547	ST6GALNAC6 Monoclonal Antibody	Q969X2	Mouse	H,M	ELISA,WB,IHC
E-AB-27582	Stanniocalcin 2 Monoclonal Antibody	O76061	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27800	STAT2 Monoclonal Antibody	P52630	Mouse	H	ELISA,WB
E-AB-27548	STAT3 Monoclonal Antibody	P40763	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27801	STAT5A Monoclonal Antibody	P42229	Mouse	H,M,R	ELISA,WB
E-AB-27802	STAT5B Monoclonal Antibody	P51692	Mouse	H	ELISA,WB
E-AB-27363	Stathmin Monoclonal Antibody	P16949	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27364	STIM1 Monoclonal Antibody	Q13586	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27365	STMN3 Monoclonal Antibody	Q9NZ72	Mouse	H	ELISA,WB
E-AB-27366	STOML2 Monoclonal Antibody	Q9UJZ1	Mouse	H,M	ELISA,WB,IHC
E-AB-27367	STXBP2 Monoclonal Antibody	Q15833	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27549	SULT1A1 Monoclonal Antibody	P50225	Mouse	H	ELISA,WB
E-AB-27550	SULT4A1 Monoclonal Antibody	Q9BR01	Mouse	H,M	ELISA,WB,IHC
E-AB-27803	SYCP3 Monoclonal Antibody	Q8IZU3	Mouse	H	ELISA,WB
E-AB-27368	Synaptophysin Monoclonal Antibody	P21579	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27804	Syntaxin 1A/1B Monoclonal Antibody	P61266	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27369	Syntaxin 6 Monoclonal Antibody	O43752	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27805	TBL2 Monoclonal Antibody	Q9Y4P3	Mouse	H	ELISA,WB,IHC
E-AB-27370	TBP Monoclonal Antibody	P20226	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27806	TCN2 Monoclonal Antibody	P20062	Mouse	H	ELISA,WB
E-AB-27371	TDP-43 Monoclonal Antibody	Q13148	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27372	TELO2 Monoclonal Antibody	Q9Y4R8	Mouse	H	ELISA,WB,IHC,IF
E-AB-27605	TEM1 Monoclonal Antibody	Q9HCU0	Mouse	H	ELISA,WB,IHC
E-AB-27807	TFF2 Monoclonal Antibody	Q03403	Mouse	H,R	ELISA,WB
E-AB-27373	TFIP11 Monoclonal Antibody	Q9UBB9	Mouse	H	ELISA,WB,IHC,IF
E-AB-27374	TGFB1 Monoclonal Antibody	Q15582	Mouse	H,M	ELISA,WB,IHC
E-AB-27808	TGOLN2 Monoclonal Antibody	O43493	Mouse	H,R	ELISA,WB
E-AB-27809	TH Monoclonal Antibody	P07101	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27375	Thyroglobulin Monoclonal Antibody	P01266	Mouse	H	ELISA,WB,IHC
E-AB-27810	TIM3 Monoclonal Antibody	Q8TDQ0	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27551	TIMM44 Monoclonal Antibody	O43615	Mouse	H,M	ELISA,WB,IHC
E-AB-27811	TINP1 Monoclonal Antibody	O95478	Mouse	H	ELISA,WB,IHC



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27376	TLE1 Monoclonal Antibody	Q04724	Mouse	H	ELISA,WB,IHC
E-AB-27377	TLE3 Monoclonal Antibody	Q04726	Mouse	H	ELISA,WB,IHC
E-AB-27812	TLR4 Monoclonal Antibody	O00206	Mouse	H,M	ELISA,WB,IHC
E-AB-27378	TMEM106B Monoclonal Antibody	Q9NUM4	Mouse	H	ELISA,WB,IHC
E-AB-27379	TMEM70 Monoclonal Antibody	Q9BUB7	Mouse	H,M	ELISA,WB,IHC
E-AB-27380	TNF alpha Monoclonal Antibody	P01375	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27553	TNFAIP1 Monoclonal Antibody	Q13829	Mouse	H	ELISA,WB
E-AB-27381	TNFR1 Monoclonal Antibody	P19438	Mouse	H,M	ELISA,WB,IHC
E-AB-27813	TNFSF9 Monoclonal Antibody	P41273	Mouse	H	ELISA,WB
E-AB-27382	TP63 Monoclonal Antibody	Q9H3D4	Mouse	H	ELISA,WB
E-AB-27383	TPI1 Monoclonal Antibody	P60174	Mouse	H	ELISA,WB
E-AB-27384	TPPII Monoclonal Antibody	P29144	Mouse	H	ELISA,WB,IHC,IF
E-AB-27814	TPRKB Monoclonal Antibody	Q9Y3C4	Mouse	H	ELISA,WB
E-AB-27385	TPSAB1 Monoclonal Antibody	Q15661	Mouse	H	ELISA,WB,IHC
E-AB-27386	TRAF3 Monoclonal Antibody	Q13114	Mouse	H,M	ELISA,WB,IF
E-AB-27387	Transcription Termination Factor Monoclonal Antibody	Q15361	Mouse	H,M	ELISA,WB
E-AB-22082	Transferrin Monoclonal Antibody	P02787	Mouse	H	WB
E-AB-27388	Transferrin Monoclonal Antibody	P02787	Mouse	H	ELISA,WB,IHC
E-AB-27389	Transgelin 2 Monoclonal Antibody	P37802	Mouse	H,M	ELISA,WB,IHC
E-AB-27390	Transgelin Monoclonal Antibody	Q01995	Mouse	H,R	ELISA,WB,IHC
E-AB-27391	Transketolase Monoclonal Antibody	P23254	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27554	TRAPPC9 Monoclonal Antibody	Q96Q05	Mouse	H,M	ELISA,WB,IHC
E-AB-27815	TRIM24 Monoclonal Antibody	O15164	Mouse	H,M	ELISA,WB,IHC
E-AB-27555	TRIM44 Monoclonal Antibody	96DX7	Mouse	H	ELISA,WB,IF
E-AB-27556	TRIP6 Monoclonal Antibody	Q15654	Mouse	H,M	ELISA,WB
E-AB-27816	TrkC Monoclonal Antibody	Q16288	Mouse	H,M,R	ELISA,WB
E-AB-27392	TRPC4AP Monoclonal Antibody	Q8TEL6	Mouse	H	ELISA,WB,IHC
E-AB-27393	TST Monoclonal Antibody	Q16762	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27394	TTF1 Monoclonal Antibody	P43699	Mouse	H	ELISA,WB,IHC
E-AB-27558	TTL Monoclonal Antibody	Q8NG68	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27817	TUBB3 Monoclonal Antibody	Q13509	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27819	TXN Monoclonal Antibody	P10599	Mouse	H	ELISA,WB



Cat. No.	Product Name	Swissprot	Host	Reactivity	Application
E-AB-27395	U2AF35 Monoclonal Antibody	Q01081	Mouse	H	ELISA,WB,IHC
E-AB-27396	UBC9 Monoclonal Antibody	P63279	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27397	UBE2C Monoclonal Antibody	O00762	Mouse	H	ELISA,WB,IHC,IF
E-AB-27606	UBE3A Monoclonal Antibody	Q05086	Mouse	H,M	ELISA,WB,IHC
E-AB-27398	UCHL1 Monoclonal Antibody	P09936	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27399	USP1 Monoclonal Antibody	O94782	Mouse	H	ELISA,WB,IHC
E-AB-27400	USP13 Monoclonal Antibody	Q92995	Mouse	H	ELISA,WB,IHC,IF
E-AB-27560	USP5 Monoclonal Antibody	P45974	Mouse	H	ELISA,WB
E-AB-27821	VAMP3 Monoclonal Antibody	Q15836	Mouse	H,M,R	ELISA,WB
E-AB-27401	VAPB Monoclonal Antibody	O95292	Mouse	H	ELISA,WB,IHC,IF
E-AB-27402	VCAM-1 Monoclonal Antibody	P19320	Mouse	H	ELISA,WB,IHC
E-AB-27403	VCP Monoclonal Antibody	P55072	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-27822	VDAC1 Monoclonal Antibody	P21796	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27823	VDAC2 Monoclonal Antibody	P45880	Mouse	H,M,R	ELISA,WB,IHC
E-AB-27404	Villin Monoclonal Antibody	P09327	Mouse	H,M	ELISA,WB,IHC,IF
E-AB-22085	Vimentin Monoclonal Antibody	P08670	Mouse	H,M,R	WB
E-AB-27405	Vimentin Monoclonal Antibody	P08670	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27406	Vinculin Monoclonal Antibody	P18206	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27561	Vitamin D binding protein Monoclonal Antibody	P02774	Mouse	H	ELISA,WB
E-AB-27824	VPS37A Monoclonal Antibody	Q8NEZ2	Mouse	H	ELISA,WB,IHC
E-AB-27825	VTN Monoclonal Antibody	P04004	Mouse	H	ELISA,WB,IHC
E-AB-27407	WTAP Monoclonal Antibody	Q15007	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-22044	XRCC4 Monoclonal Antibody	Q13426	Mouse	H	WB,IP
E-AB-27408	YY1 Monoclonal Antibody	P25490	Mouse	H,M,R	ELISA,WB,IHC,IF
E-AB-27409	ZAP70 Monoclonal Antibody	P43403	Mouse	H,M	ELISA,WB,IHC
E-AB-27410	ZAP70 Monoclonal Antibody	P43404	Mouse	H,M	ELISA,WB
E-AB-27826	ZC3HAV1 Monoclonal Antibody	Q7Z2W4	Mouse	H	ELISA,WB,IHC,IF
E-AB-27411	ZEB1 Monoclonal Antibody	P37275	Mouse	H	ELISA,WB
E-AB-27827	ZO-1 Monoclonal Antibody	Q07157	Mouse	H	ELISA,WB,IHC
E-AB-27412	Zyxin Monoclonal Antibody	Q15942	Mouse	H,M	ELISA,WB,IHC,IF

Elabscience Biotechnology Inc.

Web: www.elabscience.com

Tel: 1-240-252-7368

Fax: 1-240-252-7376

Order: orders@elabscience.com

Add: 14780 memorial drive, suite 216, houston, TX77079