

Elabscience®

Neuroscience Research

Focus on your research Service for life science




www.elabscience.com



Neuroscience

Neuroscience, or neurobiology, is a field that is devoted to the scientific study of the nervous system. Diseases of the brain, collectively called neurological diseases, can have widespread psychological and biological effects, causing harm to the diseased person. With increasing population of patients suffering from neurological diseases such as Alzheimer's and Parkinson's disease, effective prevention and treatment of neurological diseases have been the major medical challenges. Numerous research resource has been devoted to understanding the cellular and molecular mechanism within the nervous system.



Elabscience offers extensive collection of sensitive and specific ELISA kits and antibodies for detection of cytokines involved in cellular communication, inflammation, growth factors and more important biomarkers in neuroscience research. Most of Elabscience products have been used by researchers all around the world. An increasing number of citations are available in scientific publications. Each of our antibodies and ELISA kits is developed with elaborate design and strict validation.



Products & Services

Products

- ELISA Kits
- Primary/ Secondary Antibodies
- Flow Cytometry Antibodies
- Recombinant Proteins
- Biochemical Kits
- Labeling Kits
- Other Immunology Reagents

Services

- Custom Gene Synthesis Service
- Custom Peptide Synthesis Service
- Custom Protein Service
- Custom Antibody Service



Our Advantages





Publications in Neuroscience Research

Elabscience products have been cited by many top publications, e.g. Nature Medicine, Circulation, Cell Research, Gastroenterology, etc.

Olfactory Deprivation Hastens Alzheimer-Like Pathologies in a Human Tau-Overexpressed Mouse Model via Activation of cdk5.
Mol Neurobiol, 2016 Jan, 53(1): 391-401.

DOI: 10.1007/s12035-014-9007-z.


Interpretation:

Excessive secretion of corticosterone in htau mice after OBX led to calcium dyshomeostasis and caused higher activity of calpain, then activated cdk5, resulting in acceleration of taupathology.


Sample: Mouse plasma.

Product


Cat. No.: E-EL-M0349 (Mouse CORT (Corticosterone) ELISA Kit)



Typical Standard Curve of Mouse CORT ELISA Kit.



The recovery of Mouse CORT spiked at three different levels in samples throughout the range of the assay was evaluated in various matrices.



Mouse CORT was measured in biological fluids.

Neuroprotective Effect of Nanodiamond in Alzheimer's disease Rat Model: a Pivotal Role for Modulating NF- κ B and STAT3 Signaling.
Mol Neurobiol, 2017 Apr, 54(3): 1906-1918.

DOI: 10.1007/s12035-016-9762-0.


Interpretation:

This study aimed at investigating the neuroprotective effect of nanodiamond (ND). ND mitigated the increase of hippocampal beta-amyloid (A β 42) and beta-site amyloid precursor protein cleaving enzyme-1 (BACE1) together with down-regulation of phosphorylated tau protein.


Sample: Rat brain.

Product

Cat. No.: E-AB-20988 (Phospho-Tau (Ser235) Polyclonal Antibody)



Western Blot analysis of HuvEc cells using Phospho-Tau (Ser235) Polyclonal Antibody at dilution of 1:1000.



Western Blot analysis of Human Brain and HuvEc cells using Phospho-Tau (Ser235) Polyclonal Antibody at dilution of 1:1000.



Products for Popular Topics in Neuroscience Research

Elabscience provides 900+ ELISA Kits and 3500+ Antibodies regarding the Neuroscience Research on human, mouse and rat.

ELISA

Target	Cat.No.	Product Name	Test Method	Detection Range	Sensitivity
ADP/Acrp30	E-EL-H0004	Human ADP/Acrp30(Adiponectin) ELISA Kit	Sandwich	0.78--50 ng/mL	0.47 ng/mL
	E-EL-M0002	Mouse ADP/Acrp30(Adiponectin) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0329	Rat ADP/Acrp30(Adiponectin) ELISA Kit	Sandwich	46.88--3000 pg/mL	28.13 pg/mL
AGER	E-EL-H0295	Human AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit	Sandwich	78.13--5000 pg/mL	46.88 pg/mL
	E-EL-M0149	Mouse AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit	Sandwich	125--8000 pg/mL	75 pg/mL
	E-EL-R0643	Rat AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
BDNF	E-EL-H0010	Human BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-M0203	Mouse BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R1235	Rat BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
CNTF	E-EL-H0039	Human CNTF(Ciliary Neurotrophic Factor) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-M0289	Mouse CNTF(Ciliary Neurotrophic Factor) ELISA Kit	Sandwich	39.06--2500 pg/mL	23.44 pg/mL
	E-EL-R0207	Rat CNTF(Ciliary Neurotrophic Factor) ELISA Kit	Sandwich	39.06--2500 pg/mL	23.44 pg/mL
CRP	E-EL-H0043	Human CRP(C-Reactive Protein) ELISA Kit	Sandwich	0.39--25 ng/mL	0.23 ng/mL
	E-EL-M0053	Mouse CRP(C-Reactive Protein) ELISA Kit	Sandwich	0.03--25 ng/mL	0.02 ng/mL
	E-EL-R0022	Rat CRP(C-Reactive Protein) ELISA Kit	Sandwich	0.31--20 ng/mL	0.19 ng/mL
CX3CL1	E-EL-H0044	Human CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
	E-EL-M0267	Mouse CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	Sandwich	0.31--20 ng/mL	0.19 ng/mL
	E-EL-R0385	Rat CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
ECF/CCL11	E-EL-H0025	Human ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-M0464	Mouse ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R0368	Rat ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
FAS/CD95	E-EL-H0067	Human FAS/CD95(Factor Related Apoptosis) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-M0484	Mouse FAS/CD95(Factor Related Apoptosis) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R0373	Rat FAS/CD95(Factor Related Apoptosis) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
FASL/TNFSF6	E-EL-H0068	Human FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-M0028	Mouse FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R2373	Rat FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
G-CSF	E-EL-H0079	Human G-CSF(Granulocyte Colony-stimulating Factor) ELISA Kit	Sandwich	39.06--2500 pg/mL	23.44 pg/mL
	E-EL-M0031	Mouse G-CSF(Granulocyte Colony-stimulating Factor) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0454	Rat G-CSF(Granulocyte Colony Stimulating Factor) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL



Target	Cat.No.	Product Name	Test Method	Detection Range	Sensitivity
GDNF	E-EL-H1495	Human GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	Sandwich	0.31--20 ng/mL	0.19 ng/mL
	E-EL-M0553	Mouse GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	Sandwich	78.13--5000 pg/mL	46.88 pg/mL
	E-EL-R0420	Rat GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	Sandwich	78.13--5000 pg/mL	46.88 pg/mL
GM-CSF	E-EL-H0081	Human GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	Sandwich	7.81--500 pg/mL	4.69 pg/mL
	E-EL-M0032	Mouse GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0008	Rat GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
IFN- γ	E-EL-H0108	Human IFN- γ (Interferon Gamma) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-M0048	Mouse IFN- γ (Interferon Gamma) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0009	Rat IFN- γ (Interferon Gamma) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
IGF-1	E-EL-H0086	Human IGF-1(Insulin-like Growth Factor 1) ELISA Kit	Sandwich	1.56--100 ng/mL	0.94 ng/mL
	E-EL-M0034	Mouse IGF-1(Insulin-like Growth Factor 1) ELISA Kit	Sandwich	78.13--5000 pg/mL	46.88 pg/mL
	E-EL-R0010	Rat IGF-1(Insulin-like Growth Factor 1) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
IL-1 β	E-EL-H0149	Human IL-1 β (Interleukin 1 Beta) ELISA Kit	Sandwich	7.81--500 pg/mL	4.69 pg/mL
	E-EL-M0037	Mouse IL-1 β (Interleukin 1 Beta) ELISA Kit	Sandwich	7.81--500 pg/mL	4.69 pg/mL
	E-EL-R0012	Rat IL-1 β (Interleukin 1 Beta) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
IL-6	E-EL-H0102	Human IL-6(Interleukin 6) ELISA Kit	Sandwich	7.81--500 pg/mL	4.69 pg/mL
	E-EL-M0044	Mouse IL-6(Interleukin 6) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R0015	Rat IL-6(Interleukin 6) ELISA Kit	Sandwich	12.5--800 pg/mL	7.5 pg/mL
IL-8	E-EL-H0048	Human IL-8(Interleukin 8) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
IL-10	E-EL-H0103	Human IL-10(Interleukin 10) ELISA Kit	Sandwich	7.81--500 pg/mL	4.69 pg/mL
	E-EL-M0046	Mouse IL-10(Interleukin 10) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0016	Rat IL-10(Interleukin 10) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
LIF	E-EL-H0094	Human LIF(Leukemia Inhibitory Factor) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-M0040	Mouse LIF(Leukemia Inhibitory Factor) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0585	Rat LIF(Leukemia Inhibitory Factor) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
MCP-1	E-EL-H0020	Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-M0006	Mouse MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0633	Rat MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
MPIF2	E-EL-H0030	Human MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	Sandwich	78.13--5000 pg/mL	46.88 pg/mL
	E-EL-M0466	Mouse MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0356	Rat MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL



Target	Cat.No.	Product Name	Test Method	Detection Range	Sensitivity
MIP-1 α	E-EL-H0021	Human MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	Sandwich	23.44--1500 pg/mL	14.06 pg/mL
	E-EL-M0007	Mouse MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0602	Rat MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
MIP-1 β	E-EL-H0022	Human MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-M0008	Mouse MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-R0603	Rat MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
MMP-3	E-EL-H1446	Human MMP-3(Matrix Metalloproteinase 3) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
	E-EL-M0626	Mouse MMP-3(Matrix Metalloproteinase 3) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
	E-EL-R0619	Rat MMP-3(Matrix Metalloproteinase 3) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
NGF	E-EL-H1205	Human NGF(Nerve Growth Factor) ELISA Kit	Sandwich	15.63--1000 pg/mL	9.38 pg/mL
	E-EL-M0815	Mouse NGF(Nerve growth factor) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R0652	Rat NGF(Nerve Growth Factor) ELISA Kit	Sandwich	39.06--2500 pg/mL	23.44 pg/mL
NRP1	E-EL-H1247	Human NRP1(Neuropilin 1) ELISA Kit	Sandwich	78.13--5000 pg/mL	46.88 pg/mL
S100B	E-EL-H1297	Human S100B(S100 Calcium Binding Protein B) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-M1033	Mouse S100B(S100 Calcium Binding Protein B) ELISA Kit	Sandwich	78.13--5000 pg/mL	46.88 pg/mL
	E-EL-R0868	Rat S100B(S100 Calcium Binding Protein B) ELISA Kit	Sandwich	62.5--4000 pg/mL	37.5 pg/mL
SDF-1/CXCL12	E-EL-H0052	Human SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
	E-EL-M1094	Mouse SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
	E-EL-R0922	Rat SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
TARC	E-EL-H0026	Human TARC(Thymus Activation Regulated Chemokine) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-M0012	Mouse TARC(Thymus Activation Regulated Chemokine) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R0969	Rat TARC(Thymus Activation Regulated Chemokine) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
TGF- β 1	E-EL-H0110	Human TGF- β 1(Transforming Growth Factor Beta 1) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-M0051	Mouse TGF- β 1(Transforming Growth Factor Beta 1) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R0084	Rat TGF- β 1(Transforming Growth Factor Beta 1) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
TIMP-1	E-EL-H0184	Human TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
	E-EL-M0641	Mouse TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
	E-EL-R0540	Rat TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	Sandwich	0.16--10 ng/mL	0.09 ng/mL
TNF- α	E-EL-H0109	Human TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	Sandwich	7.81--500 pg/mL	4.69 pg/mL
	E-EL-M0049	Mouse TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R0019	Rat TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	Sandwich	78.13--5000 pg/mL	46.88 pg/mL
VEGF-A	E-EL-H0111	Human VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-M1292	Mouse VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL
	E-EL-R2603	Rat VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	Sandwich	31.25--2000 pg/mL	18.75 pg/mL



Antibody

Cat.No.	Product Name	Host	Reactivity	Applications
E-AB-30415	ADIPOQ Polyclonal Antibody	Rabbit	H,M,R	WB,ELISA
E-AB-10500	BDNF Polyclonal Antibody	Rabbit	H,M,R	WB,IHC-p,ELISA
E-AB-30189	CCL11 Polyclonal Antibody	Rabbit	H,M,R	WB,ELISA
E-AB-12309	CCL17 Polyclonal Antibody	Rabbit	H,M	WB,IHC-p,ELISA
E-AB-12307	CCL2 Polyclonal Antibody	Rabbit	H	IHC-p,ELISA
E-AB-10216	CCL24 Polyclonal Antibody	Rabbit	H	WB,ELISA
E-AB-33420	CCL26 Polyclonal Antibody	Rabbit	H	WB,IHC-p,ELISA
E-AB-33477	CCL4 Polyclonal Antibody	Rabbit	H	WB,ELISA
E-AB-30055	Cleaved-NOTCH2 (A1734) Polyclonal Antibody	Rabbit	H,M,R	WB,ELISA
E-AB-40235	CNTF Polyclonal Antibody	Rabbit	H	IHC-p
E-AB-33466	CRP Polyclonal Antibody	Rabbit	H,R	WB,IHC-p,ELISA
E-AB-36275	CSF2 Polyclonal Antibody	Rabbit	H	IHC-p,ELISA
E-AB-12464	CSF3 Polyclonal Antibody	Rabbit	H	IHC-p,ELISA
E-AB-30101	CX3CL1 Polyclonal Antibody	Rabbit	H	WB,IHC-p,ELISA
E-AB-32858	CXCL12 Polyclonal Antibody	Rabbit	H,M,R	WB,IHC-p,IF,ELISA
E-AB-30842	FAS Polyclonal Antibody	Rabbit	H	WB,IHC-p,ELISA
E-AB-31410	FASLG Polyclonal Antibody	Rabbit	H,M	WB,IHC-p,IF,ELISA
E-AB-70291	GDNF Polyclonal Antibody	Rabbit	M,R	IHC-p
E-AB-31763	IFNG Polyclonal Antibody	Rabbit	H	WB,IHC-p,ELISA
E-AB-64111	IGF1 Polyclonal Antibody	Rabbit	H,M,R	WB
E-AB-40248	IL1A Polyclonal Antibody	Rabbit	H	IHC-p
E-AB-33471	IL1B Polyclonal Antibody	Rabbit	H,M,R	WB,IHC-p,ELISA
E-AB-33415	IL4 Polyclonal Antibody	Rabbit	H	WB,IHC-p,ELISA
E-AB-30095	IL6 Polyclonal Antibody	Rabbit	H	WB,IHC-p,ELISA
E-AB-31801	IL8 Polyclonal Antibody	Rabbit	H	WB,IHC-p,ELISA
E-AB-33408	IL10 Polyclonal Antibody	Rabbit	H	WB,IHC-p,ELISA
E-AB-33377	LIF Polyclonal Antibody	Rabbit	H	WB,IHC-p,ELISA
E-AB-32054	MMP2 Polyclonal Antibody	Rabbit	H,M,R,MK	WB,IHC-p,IF,ELISA
E-AB-40220	MMP3 Polyclonal Antibody	Rabbit	H	IHC-p
E-AB-32239	NGF Polyclonal Antibody	Rabbit	H,M,R	WB,IHC-p,ELISA
E-AB-12815	NOTCH1 Polyclonal Antibody	Rabbit	H,R	WB,IHC-p,ELISA
E-AB-32211	NRP1 Polyclonal Antibody	Rabbit	H,M,R	WB,IHC-p,ELISA
E-AB-70136	RAGE Polyclonal Antibody	Rabbit	H,M,R	WB
E-AB-36018	S100B Polyclonal Antibody	Rabbit	H,M,R	WB,IHC-p
E-AB-22214	TGFB1 Monoclonal Antibody	Mouse	M,R	IHC-p
E-AB-33105	TIMP1 Polyclonal Antibody	Rabbit	H,M,R	WB,IHC-p,IF,ELISA
E-AB-22116	TNF alpha Monoclonal Antibody	Mouse	H,M,R	WB,IHC-p
E-AB-34220	VEGFA Polyclonal Antibody	Rabbit	H,M,R	IHC-p,ELISA



We are constantly developing new products to provide new tools for your Neuroscience Research. Please visit our website, www.elabscience.com for up-to-date details on all the products we have. For ordering, please email to orders@elabscience.com.

Elabscience Biotechnology Inc.

Web: www.elabscience.com

Tel: 1-240-252-7368

Fax: 1-240-252-7376

Email: techsupport@elabscience.com

76354126329727672839840583162797

1826472578

154273410927

00927513428

0483516109362

6382795

3162797

338912541083651730389

381384683092713038178

049438543508917815327

04958724170338410289161

76354126329727672839840

