

Name : SARS-CoV-2 ORF3a

Protein: P0DTC3

Length: 275 AA

## Sequence:

MDLFMRIFTI GTVTLKQGEI KDATPSDFVR ATATIPIQAS LPFGWLIVGV ALLAVFQSAS KIITLKKRWQ  
LALSKGVHVF CNLLLLFVTY YSHLLLVAAG LEAPFLYLYA LUYFLQSINF VRIIMRLWLC WKCRSKNPLL  
YDANYFLCWH TNCYDYCIPY NSVTSSIVIT SGDGTTPIS EHDYQIGGYT EKWESGVKDC VVLHSYFTSD  
YYQLYSTQLS TDTGVEHVTF FIYNKIVDEP EEHVQIHTID GSSGVVNPVM EPIYDEPTTT TSVPL

You can design the peptide sequences according to the Elabscience® peptide design principle or send the demand to [techsupport@elabscience.com](mailto:techsupport@elabscience.com). Elabscience® peptide expert team will provide you with professional technical support.