



Antibodies

Meeting Research Needs Across Multiple Fields

- ⌘ Virus Antibodies
- ⌘ Polyclonal Antibodies
- ⌘ Tags&Controls Antibodies
- ⌘ Acetyl/Phospho/Methyl Antibodies
- ⌘ Monoclonal Antibodies
- ⌘ Pathological Antibodies
- ⌘ Recombinant Antibodies



A Reliable Research Partner in Life Science and Medicine

Introduction to Elabscience® Antibody Products

Elabscience® has established advanced antibody development platforms that generate over 17,000 high-quality antibodies. Our antibody products support research in life science and medicine, including cell biology, cancer, neurobiology, epigenetics, infectious diseases, and immunology.

Product Features

Scientifically Sound

High specificity, high affinity

Strict Quality Control Systems

Standardized operating procedures ensure high batch-to-batch consistency

Cost-effective

An efficient supply chain allows for rapid turnaround time and competitive pricing

Technical Support


Professional, on-demand assistance for optimized procedural guidance




Elabscience® was awarded as "Antibody Supplier to Watch in 2021" by CiteAb

Antibodies by Applications


WB Antibodies




IHC Antibodies



IF Antibodies




IP Antibodies




ELISA Antibodies




FCM Antibodies



Neutralization Antibodies



Block Antibodies



Antibodies by Research Areas



Basic Biology and Cell Status

- Cell Biology
- Stem Cells
- Cell Adhesion
- Protein Trafficking and Clearance
- Developmental Biology
- Apoptosis/Tumor Suppressors/Cell Death



Immune System and Host Defense

- Immunology
- Adaptive Immunity
- Innate Immunity
- Costimulatory Molecules
- Inhibitory Molecules
- Complement



Cancer Research and Translational Medicine

- Cancer
- Angiogenesis
- Cancer Biomarkers
- Tags and Cell Markers
- Immuno-oncology



Neuroscience and Brain Diseases

- Neuroscience
- Neurodegeneration
- Neuroinflammation



Metabolism and Signal Networks

- Metabolism
- Signal Transduction
- Transcription Factors
- Epigenetics and Nuclear Signaling



Application and Systems Biology

- Cell Markers
- Cardiovascular
- Microbiology

Featured Product Lists

Cancer

Classification	Product Name	Cat. No.	Application	Reactivity	Host
Tumor Microenvironment	Integrin-β1/CD29 Polyclonal Antibody	E-AB-93359	WB, IF	Human, Mouse, Rat	Rabbit
	alpha Tubulin Monoclonal Antibody	E-AB-20036	WB, IHC-p, IF, IP	Human, Mouse, Rat	Mouse
	Cytokeratin-18 Polyclonal Antibody	E-AB-64048	IHC, IF	Human, Mouse, Rat	Rabbit
	JNK1/2/3 Polyclonal Antibody	E-AB-31853	WB, IHC-p, IF	Human, Mouse, Rat, Chicken	Rabbit
Tumor Proliferation	Phospho-Pan-Akt (Ser473) Polyclonal Antibody	E-AB-20802	WB, IHC-p, IF	Human, Mouse, Rat	Rabbit
	CEACAM5 Polyclonal Antibody	E-AB-60116	WB	Human, Mouse	Rabbit
	Phospho-ERK 1/2 (Tyr204) Polyclonal Antibody	E-AB-20869	WB, IHC-p, IF	Human, Mouse, Rat	Rabbit
	PD-L1/CD274 Polyclonal Antibody	E-AB-92367	WB, IHC	Human, Mouse	Rabbit
	RAPSN Polyclonal Antibody	E-AB-53167	IHC	Human, Mouse	Rabbit
Tumor Immune Checkpoint	CD27 Polyclonal Antibody	E-AB-11055	WB, IHC	Human	Rabbit
	CD40 Polyclonal Antibody	E-AB-16318	IHC	Human, Mouse	Rabbit
	TNF alpha Polyclonal Antibody	E-AB-33121	WB, IHC-p, IF	Human, Mouse, Rat	Rabbit
	IL4 Polyclonal Antibody	E-AB-62102	IHC	Human	Rabbit
	CD86 Polyclonal Antibody	E-AB-18593	WB, IHC	Human	Rabbit
Tumor Biomarkers	CgA Monoclonal Antibody	PA7034	IHC-P	Human	Rabbit
	CD34 Polyclonal Antibody	E-AB-60105	IHC, IF	Human, Rat	Rabbit
	CD4 Monoclonal Antibody	E-AB-22098	IHC-p, IF	Human, Mouse, Rat	Mouse
	PCNA Monoclonal Antibody	E-AB-22001	WB, IHC-p, IF	Human, Mouse, Rat	Mouse
	FGFR2 Polyclonal Antibody	E-AB-60590	IHC	Human, Mouse, Rat	Rabbit
	PD-L1 Polyclonal Antibody	E-AB-60391	IF	Human, Mouse, Rat	Rabbit

■ **Neurobiology**

Classification	Product Name	Cat. No.	Application	Reactivity	Host
Neural Lineage Markers	SOX2 Polyclonal Antibody	E-AB-18159	WB, IHC	Human, Mouse	Rabbit
	Occludin Polyclonal Antibody	E-AB-68357	WB	Human, Mouse, Rat	Rabbit
	GFAP Monoclonal Antibody	E-AB-70205	WB, IHC, IF	Mouse, Rat	Mouse
	HES1 Polyclonal Antibody	E-AB-70255	WB	Mouse, Rat	Rabbit
	PAX6 Polyclonal Antibody	E-AB-61653	IF	Human, Mouse, Rat	Rabbit
	TH Polyclonal Antibody	E-AB-70077	WB, IHC, IF	Human, Mouse, Rat	Rabbit
	CHAT Polyclonal Antibody	E-AB-70028	IHC	Mouse, Rat	Rabbit
Neurotransmission	CGRP Polyclonal Antibody	E-AB-93381	WB, IF	Human, Mouse, Rat	Rabbit
	VDAC1 Polyclonal Antibody	E-AB-52417	WB, IHC	Human, Mouse, Rat	Rabbit
	IP3R Polyclonal Antibody	E-AB-93302	WB, IF	Human, Mouse, Rat	Rabbit
	NMDAR1 Polyclonal Antibody	E-AB-93399	WB, IHC, IF	Human, Mouse, Rat	Rabbit
	GRIA2 Polyclonal Antibody	E-AB-65661	WB	Human, Mouse, Rat	Rabbit
	RYR2 Polyclonal Antibody	E-AB-32840	WB, IHC-p, IF	Human, Mouse, Rat	Rabbit
	Calmodulin 1/2/3 Polyclonal Antibody	E-AB-90310	WB	Human, Mouse, Rat	Rabbit
Neurogenesis	BMP2 Polyclonal Antibody	E-AB-70111	WB, IHC	Human, Mouse, Rat	Rabbit
	Ki67 Monoclonal Antibody	AN004340L	IHC	Human, Mouse, Rat	Rabbit
	Nestin Polyclonal Antibody	E-AB-63599	WB, IHC, IF	Human, Mouse, Rat	Rabbit
	NOTCH1 Polyclonal Antibody	E-AB-12815	WB, IHC	Human, Rat	Rabbit
	Vimentin Polyclonal Antibody	E-AB-70081	WB, IHC	Human, Mouse, Rat	Rabbit
	Recombinant SOX9 Monoclonal Antibody	E-AB-81453	WB, IHC-P	Human, Mouse	Rabbit
	Catenin beta Monoclonal Antibody	E-AB-22111	WB, IHC-p, IF	Human, Mouse, Rat	Mouse
Neurodegenerative Diseases	APP Polyclonal Antibody	E-AB-68269	IF	Human, Mouse, Rat	Rabbit
	AQP4 Monoclonal Antibody	E-AB-22025	WB, IHC-p, IF	Human, Mouse, Rat	Mouse
	HMOX1 Polyclonal Antibody	E-AB-10366	IHC	Human, Mouse, Rat	Rabbit
	HSD17B10 Polyclonal Antibody	E-AB-60967	IHC, IF	Human, Mouse, Rat	Rabbit
	NEFH Polyclonal Antibody	E-AB-16639	IHC	Human	Rabbit
	Tau Monoclonal Antibody	E-AB-22141	WB, IHC-p	Human, Mouse, Rat	Mouse
	TP53 Polyclonal Antibody	E-AB-52130	WB	Human	Rabbit
Neuroinflammation	CD11B Polyclonal Antibody	E-AB-60354	IF	Human, Mouse, Rat	Rabbit
	Caspase-1 Polyclonal Antibody	E-AB-40529	WB	Zebrafish	Rabbit
	CD163 Polyclonal Antibody	E-AB-70306	WB	Human, Mouse, Rat	Rabbit
	IL10 Polyclonal Antibody	E-AB-70235	IHC, IF	Mouse, Rat	Rabbit
	Phospho-P38 (Thr180/Tyr182) Polyclonal Antibody	E-AB-21027	WB, IHC-p, IF	Human, Mouse, Rat	Rabbit
	NOX2/gp91phox Polyclonal Antibody	E-AB-93080	WB	Human, Mouse, Rat	Rabbit
	P38 Polyclonal Antibody	D-AB-10194L	WB, IF	Human, Mouse	Rabbit
Neurotrophin	(KO Validated) FGF2 Polyclonal Antibody	E-AB-93292	WB	Human, Mouse, Rat	Rabbit
	BDNF Polyclonal Antibody	E-AB-18244	WB, IHC	Human, Mouse, Rat	Rabbit
	GTPBP1 Polyclonal Antibody	E-AB-91132	WB	Human, Mouse	Rabbit
	NGF Polyclonal Antibody	E-AB-10487	IHC	Human, Mouse, Rat	Rabbit
	NRG1 Polyclonal Antibody	E-AB-13305	WB	Human	Rabbit
	NTF4 Polyclonal Antibody	E-AB-60174	WB, IHC	Human, Rat	Rabbit

■ Programmed Cell Death

Classification	Product Name	Cat. No.	Application	Reactivity	Host
Ferroptosis	STAT3 Monoclonal Antibody	E-AB-22166	WB, IHC-p	Human, Mouse, Rat	Mouse
	YAP1 Polyclonal Antibody	E-AB-62344	IHC, IF	Human, Mouse, Rat	Rabbit
	Ceruloplasmin Polyclonal Antibody	E-AB-90643	IF	Human, Mouse, Rat	Rabbit
	RAF1 Polyclonal Antibody	E-AB-60020	WB, IHC, IF	Human, Mouse, Rat	Rabbit
	Fth1 Polyclonal Antibody	D-AB-10200L	WB	Mouse, Rat	Rabbit
	Heme Oxygenase 1 Polyclonal Antibody	E-AB-93300	WB, IF	Human, Mouse, Rat	Rabbit
Apoptosis	BAX Polyclonal Antibody	D-AB-10386L	WB	Human, Mouse, Rat	Rabbit
	BAD Polyclonal Antibody	E-AB-70084	IHC	Human, Mouse, Rat	Rabbit
	AIF1 Polyclonal Antibody	E-AB-10382	WB, IHC	Human, Mouse, Rat	Rabbit
	TRAF6 Polyclonal Antibody	E-AB-18251	WB, IHC	Human, Mouse, Rat	Rabbit
	AIF Polyclonal Antibody	E-AB-90029	WB, IHC, IF	Human, Mouse, Rat	Rabbit
	FAS Polyclonal Antibody	E-AB-13830	WB, IHC	Human	Rabbit
	FADD Polyclonal Antibody	E-AB-10318	WB, IHC	Human, Mouse	Rabbit
	Caspase-3 Polyclonal Antibody	E-AB-60017	WB, IF	Human, Mouse, Rat	Rabbit
Pyroptosis	AIM2 Polyclonal Antibody	E-AB-10974	IHC	Human	Rabbit
	CASP1 Polyclonal Antibody	E-AB-93322	WB, IF	Human, Mouse, Rat	Rabbit
	NLRP3 Polyclonal Antibody	E-AB-93112	WB, IF	Human, Mouse, Rat	Rabbit
	ASC / TMS1 Polyclonal Antibody	E-AB-60228	WB, IF	Human, Mouse, Rat	Rabbit
	IL18 Polyclonal Antibody	E-AB-14154	WB, IHC	Human	Rabbit
	IL1 beta Polyclonal Antibody	E-AB-40530	WB	Zebrafish	Rabbit
	CDK1 Polyclonal Antibody	E-AB-60019	WB, IHC, IF	Human, Mouse, Rat	Rabbit
Cuproptosis	PDHA1 Polyclonal Antibody	E-AB-52930	WB, IHC	Human, Mouse, Rat	Rabbit
	GLS Polyclonal Antibody	E-AB-19340	WB, IHC	Human, Mouse, Rat	Rabbit
	CDKN2A Polyclonal Antibody	E-AB-19334	IHC	Human	Rabbit
	DLD Polyclonal Antibody	E-AB-14947	IHC	Human, Mouse, Rat	Rabbit
	PDHB Polyclonal Antibody	E-AB-61420	WB	Human, Mouse, Rat	Rabbit
	DBT Polyclonal Antibody	E-AB-92231	WB, IHC, IF	Human, Mouse, Rat	Rabbit
	GCSH Polyclonal Antibody	E-AB-19143	IHC	Human, Mouse, Rat	Rabbit
Autophagy	ATG5 Polyclonal Antibody	E-AB-10814	WB, IHC	Human, Mouse, Rat	Rabbit
	Bcl-2 Polyclonal Antibody	E-AB-60012	WB, IHC, IF	Human, Mouse, Rat	Rabbit
	(KO Validated) LC3B Polyclonal Antibody	E-AB-63506	WB, IHC, IF	Human, Mouse, Rat	Rabbit
	BECN1 Monoclonal Antibody	E-AB-22139	WB, IHC-p	Human, Mouse, Rat	Mouse
	Recombinant LC3A Monoclonal Antibody	E-AB-81578	WB	Human, Mouse, Rat	Rabbit
	LC3A/LC3B Polyclonal Antibody	E-AB-61027	IF	Human, Mouse, Rat	Rabbit
	MAP1LC3A Monoclonal Antibody	E-AB-22136	WB, IHC-p, IF	Human, Mouse, Rat	Mouse
Necroptosis	TNF RI/TNFRSF1A Polyclonal Antibody	AN005930L	WB	Human, Rat	Rabbit
	TNFR2 Polyclonal Antibody	E-AB-60182	WB	Human, Mouse, Rat	Rabbit
	RIPK1 Polyclonal Antibody	E-AB-18284	WB, IHC	Human	Rabbit
	Phospho-MLKL (Ser358) Monoclonal Antibody	E-AB-21338	IHC-p	Human	Mouse
	MLKL Polyclonal Antibody	E-AB-18957	IHC	Human	Rabbit
	CASP8 Monoclonal Antibody	E-AB-22107	WB, IHC-p, IF	Human, Mouse, Rat	Mouse

■ Cell Metabolism

Cat. No.	Product Name	Cat. No.	Application	Reactivity	Host
Carbohydrate metabolism	ENO2 Polyclonal Antibody	E-AB-12823	WB	Human, Mouse, Rat	Rabbit
	LDHA Polyclonal Antibody	D-AB-10210L	WB, IHC	Human, Mouse, Rat	Rabbit
	Insulin Polyclonal Antibody	E-AB-70202	IHC	Mouse	Rabbit
	c-Myc Polyclonal Antibody	E-AB-30975	WB, IHC-p	Human, Mouse, Rat	Rabbit
Fat Metabolism	ACACA Polyclonal Antibody	E-AB-66498	WB, IHC	Human, Mouse, Rat	Rabbit
	PRKAB1 Polyclonal Antibody	E-AB-62178	IHC	Human, Mouse, Rat	Rabbit
Amino Acid Metabolism	NOS2 Polyclonal Antibody	E-AB-10495	IHC	Human	Rabbit
	NOS3 Monoclonal Antibody	E-AB-22088	WB	Human, Mouse, Rat	Mouse
Signaling Pathway	AKT1 Polyclonal Antibody	E-AB-30467	WB, IHC-p	Human, Mouse, Rat	Rabbit
	Pan-Akt Polyclonal Antibody	E-AB-30460	WB, IHC-p	Human, Mouse, Rat	Rabbit
Cellular Metabolic Stress Response	eNOS Polyclonal Antibody	E-AB-66416	IF	Human, Mouse, Rat	Rabbit
	NOX4 Polyclonal Antibody	E-AB-70215	WB, IHC	Human, Mouse, Rat	Rabbit
	GPX4 Polyclonal Antibody	E-AB-67390	WB, IF	Human, Mouse, Rat	Rabbit
	NQO1 Polyclonal Antibody	E-AB-16011	WB	Human, Mouse, Rat	Rabbit
	SOD1 Polyclonal Antibody	E-AB-40047	WB	Human	Rabbit

For more antibody products, please visit www.elabscience.com or contact local distributors.

Citations Featuring Our Product

Autonomous transposons tune their sequences to ensure somatic suppression, *Nature*, 2024.

Product: SRSF3 Polyclonal Antibody (E-AB-32966)

Cross-disorder and disease-specific pathways in dementia revealed by single-cell genomics, *Cell*, 2024.

Product: KCNH7 Polyclonal Antibody (E-AB-53444)

Maternal natural killer cells drive neuroimmune disorders in offspring through aberrant secretion of extracellular granzyme B, *Signal Transduction and Targeted Therapy*, 2025.

Product: DLAT Polyclonal Antibody (E-AB-66302)

JUN mediates the senescence associated secretory phenotype and immune cell recruitment to prevent prostate cancer progression, *Molecular Cancer*, 2024.

Product: EpCAM Antibody (E-AB-70132)

For more citation information, please visit www.elabscience.com.

Elabscience Bionovation Inc.

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: orders@elabscience.com; techsupport@elabscience.com

EF0059-2505v1