

Featured Recombinant Proteins

Empowering Biopharmaceutical Innovation

Elabscience® offers an extensive selection of over 6,000 proteins, including cytokines, growth factors, immune checkpoint proteins, CD antigens, Fc receptors, enzymes, and hormones, for R&D in pharmaceutical companies.



High Quality

Purity > 90%, endotoxin < 1.0 EU/μg, validated bioactivity



Strict Quality Control Systems

Standardized operating procedures ensure high inter-batch consistency



Versatile Tags

70+ types of tags are available, such as GST, His, Fc, SUMO, Trx, and more



Customization Services

One-stop solution from codon optimization to protein expression and purification



Quick Turnaround Time

An efficient supply chain ensures rapid product delivery



Technical Support

Professional, on-demand assistance for optimized experimental guidance

Functional Research Areas

Stem Cell culture

PD-1 GM-CSF 4-1BB
FGF-2 IL-6 Wnt3a
FGF-2 EGF
LIF Vitronectin

Immune Cell Culture

IL-21 TNF-α IL-15
CD47 IL-6 IFN-γ
M-CSF Fas Ligand
IL-18 NKp30 ICOS

Organoid Culture

Wnt-3a R-spondin 1
Noggin HGF MMP1
BMP-2 IGF-1 Activin A
Noggin SHH TGFB1

Diagnostic Reagents

HIV gp120 PCT AFP
Myoglobin Ferritin
CA125 PSA CA125
Roquin I Ferritin



Vaccine Development

SARS-CoV-2 Helicase
HCoV-229E Spike1
MERS-CoV Nucleoprotein VP1
HCV E2 HIV-1 gag Protein p24

Gene Therapy

VEGF165/VEGFA
BDNF IL-2 PD-1
FGF-2 IFN-γ HGF
BCMA Insulin

Experimental Validation



Featured Products

Category	Cat. No.	Product Name	Tag	Bio-activity	Endotoxin	Accession
Cytokines	PKSH033475	Recombinant Human VEGF165/VEGFA Protein	None	Active	< 0.01 EU/ μ g	P15692-4
	PDEH100712	Recombinant Human IL-6 Protein	N-His	Active	< 0.01 EU/ μ g	P05231
Growth Factors	PKSH032439	Recombinant Human FGF-2/FGFb Protein	None	Active	< 0.01 EU/ μ g	P09038-4
	PKSH031612	Recombinant Human IGFBP3/IBP3 Protein	N-His	Active	< 1.0 EU/ μ g	NP_000589.2
Immune Checkpoint Proteins	PKSH032861	Recombinant Human PD-1/PDCD1 Protein	C-His	Active	< 1.0 EU/ μ g	Q15116
	PKSH032026	Recombinant Human 4-1BB/TNFRSF9 Protein	C-Fc	Active	< 1.0 EU/ μ g	Q07011
CD Antigens	PKSH030305	Recombinant Human CD47 Protein	None	Active	< 1.0 EU/ μ g	NP_942088.1
Enzymes	PKSH031542	Recombinant Human MMP1 Protein	C-His	Active	< 1.0 EU/ μ g	NP_002412.1
Hormones	PKSH030373	Recombinant Human Insulin Receptor/INSR Protein	N-His-GST	Active	< 1.0 EU/ μ g	NP_000199.2
Viral Proteins	PKSR030474	Recombinant SARS-CoV-2 Spike Protein	C-mFc	Active	< 1.0 EU/ μ g	QHD43416.1
	PKSV030241	Recombinant MERS-CoV Spike/S1 Protein	C-His	Active	< 1.0 EU/ μ g	AFS88936.1

For more related products, please visit www.elabscience.com or contact local distributors.

Citations Featuring Our Products

- A targetable secreted neural protein drives pancreatic cancer metastatic colonization and HIF1a nuclear retention. *Cancer Discovery*, 2024; 14 (12): 2489-2508
Products: Recombinant Human AMIGO2 Protein (PKSH030569)
- Supersulphides provide airway protection in viral and chronic lung diseases. *Nature Communications*, 2023; 14 (1): 4476
Products: Recombinant Human Cathepsin B/CTSB Protein (PKSH031575)
- R-spondin 2 governs xenopus left-right body axis formation by establishing an FGF signaling gradient. *Nature Communications*, 2024; 15 (1): 1003
Products: Recombinant Human Wnt3a/Wnt-3a (PKSH033972)

For more citation information, please visit www.elabscience.com.

Elabscience Bionovation Inc.

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: orders@elabscience.com; techsupport@elabscience.com

EFL061-2504V1