

Research Protocol for Mouse T Cell Analysis

One-Stop Service, From Isolation to Detection



EasySort™ Mouse T Cell Isolation Kit

High-purity isolated T cells



Mouse CD3/CD28 T Cell Activation Beads

APC-mimetic magnetic beads enable efficient T cell activation



Recombinant Mouse Interleukin-2/ IL-2 Protein (active)


High-quality protein, stable expansion



CFSE Cell Division Tracker Kit

Accurate detection of cell proliferation





T Cell Function Research Experimental Procedure



T Cell Function Research Products

Product Name	Cat. No.	Size	Application
EasySort™ Mouse CD3+T Cell Isolation Kit	MIM001N	10/100/200 Assays	Cell Isolation
EasySort™ Mouse CD4+T Cell Isolation Kit	MIM002N	10/100/200 Assays	
EasySort™ Mouse CD8+T Cell Isolation Kit	MIM003N	10/100/200 Assays	
Mouse CD3/CD28 T Cell Activation Beads	MIM001A	0.2 mL/1 mL/1 mL*5	T Cell Activation
Recombinant Mouse Interleukin-2/IL-2 Protein (active)	PKSM041523	20/50 µg	T Cell Proliferation
CFSE Cell Division Tracker Kit	E-CK-A345	500/2000 Assays	Cell Proliferation Detection

Experimental Results



Elabscience Bionovation Inc.

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: orders@elabscience.com; techsupport@elabscience.com

EFL110-2505v1