
Elabscience®



Regenerative Medicine

Regenerative Medicine

Regenerative medicine is a branch of translational research in tissue engineering and molecular biology, which deals with replacing, engineering or regenerating human cells, tissues or organs to restore or establish normal functions. This field holds the promise of engineering damaged tissues and organs by stimulating the body's own repair mechanism to functionally heal previous irreparable tissues or organs.

Regenerative medicine also makes the possibility of culturing tissues and organs in the laboratory and implanting them when the body cannot heal itself. If cells of regenerated organs would be derived from the patients' own tissues or cells, this would potentially solve the problems of donor organs shortage and organ transplant rejection.

Some of the biomedical approaches within the field of regenerative medicine may involve the use of stem cells. Examples include the injection of stem cells or progenitor cells obtained through directed differentiation(cell therapy), the induction of regeneration by biologically active molecules administered alone or as a secretion by infused cells(immune modulation therapy), and transplantation of in vitro grown organs and tissues(tissue engineering).

The use of cord blood stem cells in treating conditions such as brain injury and type 1 diabetes has already been studied in humans, and earlier research is being conducted for treatment of stroke and hearing loss.



Contents

Cell therapy	 1
Tissue engineering	 5
Organ regeneration	 7
Therapeutic cloning	 8
Stem cells	 9

Cell therapy

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0910	Human TP53 (Tumor Protein p53) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H0109	Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0060	Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0114	Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0059	Human EGF (Epidermal Growth Factor) ELISA Kit	3.91250pg/mL	2.34pg/mL	Sandwich
E-EL-H0756	Human c-myc (c-myc Oncogene Product) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0108	Human IFN-γ (Interferon Gamma) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H2327	Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit	0.6340ng/mL	0.38ng/mL	Sandwich
E-EL-H0483	Human bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	12.5800pg/mL	7.5pg/mL	Sandwich
E-EL-H1629	Human EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H2092	Human CCND1 (Cyclin-D1) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0111	Human VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1299	Human SCF (Stem Cell Factor) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0666	Human CTNNβ1 (Catenin, Beta 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0084	Human HGF (Hepatocyte Growth Factor) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H0086	Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H2182	Human KRAS (V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0099	Human IL-2 (Interleukin 2) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H1133	Human ERα (Estrogen Receptor Alpha) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H0067	Human FAS/CD95 (Factor Related Apoptosis) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H2473	Human CDKN1A (Cyclin Dependent Kinase Inhibitor 1A)ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H1603	Human VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H0103	Human IL-10 (Interleukin 10) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H1030	Human IL-2Rα (Interleukin 2 Receptor Alpha) ELISA Kit	93.756000pg/mL	56.25pg/mL	Sandwich
E-EL-H0562	Human BAX (Bcl-2 Associated X Protein) ELISA Kit	0.7850ng/mL	0.47ng/mL	Sandwich

E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.13-5000pg/mt. 46.88pg/mt. Sandwich E-EL-H0930 Human SHH (Hedgehog Homolog, Sanic) ELISA Kit 78.13-5000pg/mt. 46.88pg/mt. Sandwich E-EL-H0742 Human CD4 (Cluster of Differentiation 4) ELISA Kit 0.31-20ng/mt. 0.19ng/mt. Sandwich E-EL-H0149 Human PCNA (Proifferating Cell Nuclear Antigen) ELISA Kit 0.31-20ng/mt. 0.19ng/mt. Sandwich E-EL-H0177 Human PCNA (Proifferating Cell Nuclear Antigen) ELISA Kit 0.31-20ng/mt. 0.19ng/mt. Sandwich E-EL-H0177 Human CSPS (Caspasa 3) ELISA Kit 0.31-20ng/mt. 0.19ng/mt. Sandwich E-EL-H0178 Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit 0.16-10ng/mt. 0.09ng/mt. Sandwich E-EL-H0035 Human SDA (Adenosine Dearminase) ELISA Kit 0.31-20ng/mt. 0.19ng/mt. Sandwich E-EL-H0035 Human SDA (Soluble Cluster of Differentiation 40 Ligand) 62.5-4000pg/mt. 37.5pg/mt. Sandwich E-EL-H0035 Human SDA (Soluble Cluster of Differentiation 44) ELISA Kit 0.16-10ng/mt. 0.09ng/mt. Sandwich E-EL-H02	Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0752 Human CD4 (Cluster of Differentiation 4) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H0149 Human IL-16 (Interleukin 1 Beta) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H02499 Human PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H0271 Human PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H0374 Human PCNAP (Claspase 3) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich ELISA Kit 0.31–20ng/mL 0.19ng/m	E-EL-H5490	Human CXCR4 (Chemokine C-X-C-Motif Receptor)ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H0149 Human IL-1β (Interleukin 1 Beta) ELISA Kit 7.81–500gg/mL 4.69gg/mL Sandwich E-EL-H2399 Human PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H0170 Human CASP3 (Caspase 3) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1846 Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H03052 Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H0313 Human SICAM-1/CD54 (Soluble Intercellular Adhesion Molecule 78.13–5000pg/mL 46.88pg/mL Sandwich E-EL-H0326 Human SCAM-1/CD54 (Soluble Intercellular Adhesion Molecule 78.13–5000pg/mL 46.88pg/mL Sandwich E-EL-H0305 Human SCD40 (Soluble Cluster of Differentiation 40 Ligand) 62.5–4000pg/mL 37.5pg/mL Sandwich E-EL-H2727 Human SCD44 (Soluble Cluster of Differentiation 44) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H3161 Human MUC1 (Mucin 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H3161 <t< td=""><td>E-EL-H0930</td><td>Human SHH (Hedgehog Homolog, Sonic) ELISA Kit</td><td>78.135000pg/mL</td><td>46.88pg/mL</td><td>Sandwich</td></t<>	E-EL-H0930	Human SHH (Hedgehog Homolog, Sonic) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H2399 Human PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H0017 Human CASP3 (Caspase 3) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1084 Human PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H0052 Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H0052 Human SICAM-1/CD54 (Soluble Intercellular Adhesion Molecule F-EL-H035 Human SICAM-1/CD54 (Soluble Intercellular Adhesion Molecule F-EL-H035 Human SCAM-1/CD54 (Soluble Intercellular Adhesion Molecule F-EL-H0362 Human ADA (Adenosine Deaminase) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H0362 Human SCD40 (Soluble Cluster of Differentiation 40 Ligand) E-EL-H0635 Human SCD40 (Soluble Cluster of Differentiation 40 Ligand) E-EL-H0616 Human MC1 (Mucin 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H1614 Human MC1 (Mucin 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H1614 Human BPAF (B-Raf Proto Oncogene Serine/Threonine Protein 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1988 EliSA Kit 1.50–10ng/mL 0.09ng/mL Sandwich E-EL-H1988 Human STAT3 (Signal Transducer and Activator of Transcription 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H079 Human FN (Fibronectin) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwich E-EL-H079 Human FN (Fibronectin) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwich E-EL-H0808 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit 0.16–10ng/mL 0.99ng/mL Sandwich E-EL-H0808 Human FR (Relinoblastoma Protein 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H099 Human FR (Relinoblastoma Protein 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H099 Human FR (Relinoblastoma Protein 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H099 Human FR (Relinoblastoma Protein 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H099 Human FR (Relinoblastoma Protein 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H099 Human FR (Relinoblastoma Protein 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H099 Human FR (E-EL-H0752	Human CD4 (Cluster of Differentiation 4) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0017 Human CASP3 (Caspase 3) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1848 Human PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H0052 Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H0135 Human SICAM-1/CD54 (Soluble Intercellular Adhesion Molecule 1) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H0262 Human ADA (Adenosine Deaminase) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H0281 Human SCD40 (Soluble Cluster of Differentiation 40 Ligand) 62.54000pg/mL 37.5pg/mL Sandwich E-EL-H0272 Human SCD44 (Soluble Cluster of Differentiation 44) ELISA Kit 62.54000pg/mL 37.5pg/mL Sandwich E-EL-H0161 Human MUC1 (Mucin 1) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich E-EL-H0141 Human SRAF (B-Raf Proto Oncogene Serine/Threonine Protein (Kinase) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H0198 Human STAT3 (Signal Transducer and Activator of Transcription (Sinase) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich	E-EL-H0149	Human IL-1β (Interleukin 1 Beta) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H1846 Human PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) 0.31–20ng/mL 0.19ng/mL Sandwich ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich 1) ELISA Kit 0.31–20ng/mL 10.19ng/mL Sandwich 1) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich 1) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich 1) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich 1 ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich 1 ELISA Kit 0.16–10ng/mL 0.19ng/mL Sandwich 1 ELISA Kit 0.16–10ng/mL 0.19ng/mL Sandwich 1 ELISA Kit 0.16–10ng/mL 0.94ng/mL Sandwich 1 ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich 1 ELISA Ki	E-EL-H2399	Human PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H10405 ELISA Kit 0.31-20ng/mL 0.19ng/mL Sandwich E-EL-H0052 Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich E-EL-H0135 Human SICAM-1/CD54 (Soluble Intercellular Adhesion Molecule To JELISA Kit 0.31-20ng/mL 0.19ng/mL Sandwich E-EL-H0262 Human ADA (Adenosine Deaminase) ELISA Kit 0.31-20ng/mL 0.19ng/mL Sandwich E-EL-H035 Human scD40 (Soluble Cluster of Differentiation 40 Ligand) 62.5-4000pg/mL 37.5pg/mL Sandwich E-EL-H0616 Human SCD44 (Soluble Cluster of Differentiation 44) ELISA Kit 62.5-4000pg/mL 37.5pg/mL Sandwich E-EL-H0616 Human MUC1 (Mucin 1) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich E-EL-H0616 Human MUC1 (Mucin 1) ELISA Kit 0.16-10ng/mL 0.19ng/mL Sandwich E-EL-H0616 Human BRAF (B-Raf Proto Oncogene Serine/Threonine Protein 0.31-20ng/mL 0.19ng/mL Sandwich E-EL-H0618 Human BRAF (B-Raf Proto Oncogene Serine/Threonine Protein 0.16-10ng/mL 0.09ng/mL Sandwich E-EL-H0818 Human BRAF (B-Raf Proto Oncogene Serine/Th	E-EL-H0017	Human CASP3 (Caspase 3) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0135 Human sICAM-1/CD54 (Soluble Intercellular Adhesion Molecule 1) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich 1) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 37.5pg/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 37.5pg/mL Sandwich 1 ELISA Kit 0.52–4000pg/mL 37.5pg/mL Sandwich 1 ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 0.31–20ng/mL 0.19ng/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 0.31–20ng/mL 0.30pg/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 0.31–20ng/mL 0.30pg/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 0.31–20ng/mL 0.34pg/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 0.34pg/mL 0.34pg/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 0.09ng/mL 0.34pg/mL Sandwich 1 ELISA Kit 0.31–20ng/mL 0.09ng/mL	E-EL-H1846		0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0262 Human ADA (Adenosine Deaminase) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H0035 Human sCD40 (Soluble Cluster of Differentiation 40 Ligand) 62.54000pg/mL 37.5pg/mL Sandwich E-EL-H0036 Human sCD40 (Soluble Cluster of Differentiation 44) ELISA Kit 62.54000pg/mL 37.5pg/mL Sandwich E-EL-H0616 Human MUC1 (Mucin 1) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich E-EL-H0616 Human BRAF (B-Raf Proto Oncogene Serine/Threonine Protein Rimsee) ELISA Kit 1.5e-L-H0616 Human BRAF (B-Raf Proto Oncogene Serine/Threonine Protein 0.3120ng/mL 0.19ng/mL Sandwich Kinase) ELISA Kit 1.5e-L-H0985 Beta) ELISA Kit 1.5e-L-H0985 Human POGFRB (Platelet Derived Growth Factor Receptor Beta) ELISA Kit 1.5e-L-H0995 Human STAT3 (Signal Transducer and Activator of Transcription 0.3120ng/mL 0.19ng/mL Sandwich 3) ELISA Kit 1.5e-L-H079 Human FN (Fibronectin) ELISA Kit 1.5e-100ng/mL 0.94ng/mL Sandwich 1.5e-L-H079 Kit 1.5e-100ng/mL 0.94ng/mL Sandwich 1.5e-L-H0048 Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Rit 1.5e-2000pg/mL 2.3.44pg/mL Sandwich 1.5e-L-H0048 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit 1.5e-2000pg/mL 9.38pg/mL Sandwich 1.5e-L-H0068 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich 1.5e-L-H0706 Human TERT (Telomerase Reverse Transcriptase) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich 1.5e-L-H0493 Human RB1 (Relinoblastoma Protein 1) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich 1.5e-L-H0493 Human AR (Androgen Receptor) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich 1.5e-L-H0493 Human AR (Androgen Receptor) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich 1.5e-L-H0493 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich 1.5e-L-H059 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich 1.5e-L-H059 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich 1.5e-L-H059 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwich 1.5e-L-H059 Human CASP8 (Caspase 8) ELISA	E-EL-H0052	Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0035	E-EL-H0135		78.135000pg/mL	46.88pg/mL	Sandwich
ELEL-H0035 ELISA Kit E-EL-H0272 Human sCD44 (Soluble Cluster of Differentiation 44) ELISA Kit E-EL-H0616 Human MUC1 (Mucin 1) ELISA Kit E-EL-H0616 Human MUC1 (Mucin 1) ELISA Kit E-EL-H1614 Human BRAF (B-Raf Proto Oncogene Serine/Threonine Protein Kinase) ELISA Kit E-EL-H0985 Human PDGFRB (Platelet Derived Growth Factor Receptor Beta) ELISA Kit E-EL-H02132 Human STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit E-EL-H0179 Human FN (Fibronectin) ELISA Kit E-EL-H0179 Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit E-EL-H0048 Human IL-8 (Interleukin 8) ELISA Kit E-EL-H0048 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit E-EL-H0706 Human TERT (Telomerase Reverse Transcriptase) ELISA Kit E-EL-H0707 Human RB1 (Retinoblastoma Protein 1) ELISA Kit E-EL-H0708 Human TERT (Telomerase Reverse Transcriptase) ELISA Kit E-EL-H0709 Human RB1 (Retinoblastoma Protein 1) ELISA Kit E-EL-H0709 Human RB1 (Retinoblastoma Protein 1) ELISA Kit E-EL-H1932 Human RB1 (Retinoblastoma Protein 1) ELISA Kit E-EL-H0943 Human AR (Androgen Receptor) ELISA Kit E-EL-H1933 Human AR (Androgen Receptor) ELISA Kit E-EL-H1939 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit E-EL-H1939 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit E-EL-H1939 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit E-EL-H1939 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit E-EL-H1938 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit E-EL-H1938 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit E-EL-H1938 Human NF-kB p105 (Nuclear Factor Kappa B p105) ELISA	E-EL-H0262	·	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0616 Human MUC1 (Mucin 1) ELISA Kit E-EL-H1614 Kinase) ELISA Kit E-EL-H1614 Human BRAF (B-Raf Proto Oncogene Serine/Threonine Protein Kinase) ELISA Kit E-EL-H0985 Human PDGFRβ (Platelet Derived Growth Factor Receptor Beta) ELISA Kit E-EL-H0985 Human STAT3 (Signal Transducer and Activator of Transcription 3.1–20ng/mL 0.19ng/mL Sandwich E-EL-H0179 Human STAT3 (Signal Transducer and Activator of Transcription 3.1–20ng/mL 0.19ng/mL Sandwich E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwich E-EL-H0079 Kit 2.3.44pg/mL Sandwich E-EL-H0079 Kit 3.1.25–2000pg/mL 18.75pg/mL Sandwich E-EL-H0048 Human IL-8 (Interleukin 8) ELISA Kit 31.25–2000pg/mL 18.75pg/mL Sandwich E-EL-H0048 Kit Sandwich Skit 31.25–2000pg/mL 3.38pg/mL Sandwich E-EL-H0068 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H0706 Human TERT (Telomerase Reverse Transcriptase) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H1932 Human RB1 (Retinoblastoma Protein 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H1939 Human AR (Androgen Receptor) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human PPARyC1α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit 0.31–2000pg/mL 18.75pg/mL Sandwich E-EL-H16659 Human CASP8 (Caspase 8) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H16659 Human CASP8 (Caspase 8) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich	E-EL-H0035		62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H1614 Kinase) ELISA Kit E-EL-H0985 Human PDGFRß (Platelet Derived Growth Factor Receptor Beta) ELISA Kit E-EL-H2132 Human STAT3 (Signal Transducer and Activator of Transcription 3.31–20ng/mL 0.19ng/mL Sandwich E-EL-H2132 Human STAT3 (Signal Transducer and Activator of Transcription 3.31–20ng/mL 0.19ng/mL Sandwich E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwich E-EL-H0079 Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA 39.06-2500pg/mL 23.44pg/mL Sandwich E-EL-H0048 Human IL-8 (Interleukin 8) ELISA Kit 31.25-2000pg/mL 18.75pg/mL Sandwich E-EL-H0068 Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H0068 Human FRSL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H0706 Human TERT (Telomerase Reverse Transcriptase) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H0493 Human RB1 (Retinoblastoma Protein 1) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H0493 Human AR (Androgen Receptor) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human Cocctivator 1 Alpha) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human COK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H1386 Human NF-KB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich	E-EL-H2272	Human sCD44 (Soluble Cluster of Differentiation 44) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H1614 Kinase) ELISA Kit E-EL-H0985 Human PDG-FRβ (Platelet Derived Growth Factor Receptor Beta) ELISA Kit E-EL-H2132 3) ELISA Kit E-EL-H2132 3) ELISA Kit E-EL-H0179 Human FN (Fibronectin) ELISA Kit E-EL-H0179 Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA 39.062500pg/mL 23.44pg/mL Sandwich E-EL-H0079 Kit E-EL-H0079 Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA 39.062500pg/mL 23.44pg/mL Sandwich E-EL-H0048 Human IL-8 (Interleukin 8) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwich E-EL-H0068 Kit E-EL-H0068 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H0706 Human TERT (Telomerase Reverse Transcriptase) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1932 Human RB1 (Retinoblastoma Protein 1) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H0493 Human AR (Androgen Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human PPARγC1α (Peroxisome Proliferator Activated Receptor 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1359 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1451 Human CASP8 (Caspase 8) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1386 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich	E-EL-H0616	Human MUC1 (Mucin 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0383 Beta) ELISA Kit 0.18-101g/mL 0.09ng/mL Sandwich E-EL-H2132 Human STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwich E-EL-H0079 Kit 39.062500pg/mL 23.44pg/mL Sandwich E-EL-H0048 Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit 31.252000pg/mL 23.44pg/mL Sandwich E-EL-H0048 Human IL-8 (Interleukin 8) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwich E-EL-H0068 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit 0.1610ng/mL 9.38pg/mL Sandwich E-EL-H0706 Human TERT (Telomerase Reverse Transcriptase) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1932 Human RB1 (Retinoblastoma Protein 1) ELISA Kit 0.1610ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human AR (Androgen Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.31	E-EL-H1614		0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2132 3) ELISA Kit 0.31-20ng/mL 0.19ng/mL Sandwich E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwich E-EL-H0079 Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit 39.06–2500pg/mL 23.44pg/mL Sandwich E-EL-H0048 Human IL-8 (Interleukin 8) ELISA Kit 31.25–2000pg/mL 18.75pg/mL Sandwich E-EL-H0068 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit 0.16–10ng/mL 9.38pg/mL Sandwich E-EL-H0706 Human TERT (Telomerase Reverse Transcriptase) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H1932 Human RB1 (Retinoblastoma Protein 1) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwich E-EL-H0493 Human AR (Androgen Receptor) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human PPARγC1α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit 0.31–20ng/mL 0.09ng/mL Sandwich E-EL-H1351 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.31–20ng/mL 18.75pg/mL Sandwich E-EL-H0659 Human CASP8 (Caspa	E-EL-H0985		0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0079 Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA (39.062500pg/mL 23.44pg/mL Sandwich Kit 31.252000pg/mL 18.75pg/mL Sandwich E-EL-H0068 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA (15.631000pg/mL 9.38pg/mL Sandwich Kit 0.1610ng/mL 0.09ng/mL Sandwich Human TERT (Telomerase Reverse Transcriptase) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1932 Human RB1 (Retinoblastoma Protein 1) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H0493 Human AR (Androgen Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human PPARγC1α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H2323 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwich E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1386 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich	E-EL-H2132		0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0079 KitKit39.062500pg/mL23.44pg/mLSandwichE-EL-H0048 E-EL-H0068 KitHuman IL-8 (Interleukin 8) ELISA Kit31.252000pg/mL18.75pg/mLSandwichE-EL-H0068 KitHuman FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit15.631000pg/mL9.38pg/mLSandwichE-EL-H0706 E-EL-H0706Human TERT (Telomerase Reverse Transcriptase) ELISA Kit0.1610ng/mL0.09ng/mLSandwichE-EL-H1932 E-EL-H0493Human RB1 (Retinoblastoma Protein 1) ELISA Kit0.1610ng/mL0.09ng/mLSandwichE-EL-H0493 Gamma Coactivator 1 Alpha) ELISA Kit0.3120ng/mL0.19ng/mLSandwichE-EL-H1359 Gamma Coactivator 1 Alpha) ELISA Kit0.1610ng/mL0.09ng/mLSandwichE-EL-H2323 E-EL-H2323 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit0.3120ng/mL0.19ng/mLSandwichE-EL-H1451 E-EL-H1451Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit31.252000pg/mL18.75pg/mLSandwichE-EL-H0659 E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit0.1610ng/mL0.09ng/mLSandwichE-EL-H1386 E-EL-H1386Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit0.1610ng/mL0.09ng/mLSandwich	E-EL-H0179	Human FN (Fibronectin) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H0068 Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit E-EL-H0706 Human TERT (Telomerase Reverse Transcriptase) ELISA Kit E-EL-H1932 Human RB1 (Retinoblastoma Protein 1) ELISA Kit E-EL-H0493 Human AR (Androgen Receptor) ELISA Kit E-EL-H0493 Human PPARγC1α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit E-EL-H1359 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit E-EL-H2323 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit E-EL-H0659 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL 0.09ng/mL Sandwich Sandwich Sandwich 0.1610ng/mL 0.09ng/mL Sandwich Sandwich	E-EL-H0079	, , , , , , , , , , , , , , , , , , , ,	39.062500pg/mL	23.44pg/mL	Sandwich
E-EL-H0068 KitKit15.631000pg/mL9.38pg/mLSandwichE-EL-H0706Human TERT (Telomerase Reverse Transcriptase) ELISA Kit0.1610ng/mL0.09ng/mLSandwichE-EL-H1932Human RB1 (Retinoblastoma Protein 1) ELISA Kit0.1610ng/mL0.09ng/mLSandwichE-EL-H0493Human AR (Androgen Receptor) ELISA Kit0.3120ng/mL0.19ng/mLSandwichE-EL-H1359Human PPARγC1α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit0.1610ng/mL0.09ng/mLSandwichE-EL-H2323Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit0.3120ng/mL0.19ng/mLSandwichE-EL-H1451Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit31.252000pg/mL18.75pg/mLSandwichE-EL-H0659Human CASP8 (Caspase 8) ELISA Kit0.1610ng/mL0.09ng/mLSandwichE-EL-H1386Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit0.1610ng/mL0.09ng/mLSandwich	E-EL-H0048	Human IL-8 (Interleukin 8) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1932 Human RB1 (Retinoblastoma Protein 1) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H0493 Human AR (Androgen Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human PPARγC1α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H2323 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwich E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1386 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich	E-EL-H0068	(15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H0493 Human AR (Androgen Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1359 Human PPARγC1α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H2323 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwich E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1386 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich	E-EL-H0706	Human TERT (Telomerase Reverse Transcriptase) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1359 Human PPARγC1α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit 0.3120ng/mL 0.09ng/mL Sandwich E-EL-H2323 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwich E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1386 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich	E-EL-H1932	Human RB1 (Retinoblastoma Protein 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1359 Gamma Coactivator 1 Alpha) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H2323 Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwich E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwich E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1386 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich	E-EL-H0493	Human AR (Androgen Receptor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H1451 Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwich E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1386 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich	E-EL-H1359		0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0659 Human CASP8 (Caspase 8) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich E-EL-H1386 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich	E-EL-H2323	Human CDK4 (Cyclin Dependent Kinase 4) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H1386 Human NF-κB p105 (Nuclear Factor Kappa B p105) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwich	E-EL-H1451	Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
	E-EL-H0659	Human CASP8 (Caspase 8) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0100 Human IL-3 (Interleukin 3) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwich	E-EL-H1386	Human NF-кВ p105 (Nuclear Factor Карра В p105) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
	E-EL-H0100	Human IL-3 (Interleukin 3) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0878	Human Flt3 (FMS Like Tyrosine Kinase 3) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2069	Human CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0010	Human BDNF (Brain Derived Neurotrophic Factor) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1495	Human GDNF (Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2214	Human PDGFRα (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0345	Human FGFR1 (Fibroblast Growth Factor Receptor 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H2239	Human JAK2 (Janus Kinase 2) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H0066	Human EPO (Erythropoietin) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0356	Human FGFR3 (Fibroblast Growth Factor Receptor 3) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1445	Human MMP-2 (Matrix Metalloproteinase 2) ELISA Kit	0.7850ng/mL	0.47ng/mL	Sandwich
E-EL-H1277	Human HIF-1α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H2407	Human ΙκΒα (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1607	Human CD34 (Cluster Of Differentiation 34) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0549	Human BCL2L1/BCL-X (Bcl-2 Like Protein 1) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2181	Human ALB (Albumin) ELISA Kit	3.13200µg/mL	1.88µg/mL	Sandwich
E-EL-H0101	Human IL-4 (Interleukin 4) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0012	Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0158	Human FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1655	Human mTOR (Mammalian Target of Rapamycin) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0081	Human GM-CSF (Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H1104	Human FOXP3 (Forkhead Box Protein P3) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0011	Human BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H0171	Human TF (Transferrin) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H0430	Human Anti-Surv (Anti-Survivin Antibody) ELISA Kit	1.56100mU/mL	0.94mU/mL	Sandwich
E-EL-H5610	Human CD19 (Cluster Of Differentiation 19) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0547	Human BCR (B-Cell Receptor) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0240	Human ARO (Aromatase) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1441	Human MMP-1 (Matrix Metalloproteinase 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0601	Human BRCA1 (Breast Cancer Susceptibility Protein 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H5634	Human TGFBR2 (TGF-beta receptor type 2) ELISA Kit	0.7850ng/mL	0.47ng/mL	Sandwich
E-EL-H1094	Human VIM (Vimentin) ELISA Kit	7.81500ng/mL	4.69ng/mL	Sandwich
E-EL-H1640	Human PECAM1/CD31 (Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1046	Human GP4/CD36 (Platelet Membrane Glycoprotein IV) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1592	Human IL7R (Interleukin 7 Receptor) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H0425	Human IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2086	Human ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H1205	Human NGF (Nerve Growth Factor) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich

Tissue engineering

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H2106	Human tPA (Plasminogen Activator, Tissue) ELISA Kit	0.6340ng/mL	0.38ng/mL	Sandwich
E-EL-H0828	Human CTGF (Connective Tissue Growth Factor) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H0184	Human TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0910	Human TP53 (Tumor Protein p53) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H2104	Human PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0111	Human VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1453	Human TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H1445	Human MMP-2 (Matrix Metalloproteinase 2) ELISA Kit	0.7850ng/mL	0.47ng/mL	Sandwich
E-EL-H0109	Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0060	Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H1163	Human ELN (Elastin) ELISA Kit	0.4730ng/mL	0.28ng/mL	Sandwich
E-EL-H1454	Human TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H2266	Human FBN1 (Fibrillin 1) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H1451	Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1046	Human GP4/CD36 (Platelet Membrane Glycoprotein IV) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H2339	Human TGF-β3 (Transforming Growth Factor β3) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H0103	Human IL-10 (Interleukin 10) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0483	Human bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	12.5800pg/mL	7.5pg/mL	Sandwich
E-EL-H0113	Human LEP (Leptin) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H1629	Human EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H2181	Human ALB (Albumin) ELISA Kit	3.13200µg/mL	1.88µg/mL	Sandwich
E-EL-H2182	Human KRAS (V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0706	Human TERT (Telomerase Reverse Transcriptase) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1846	Human PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0158	Human FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0002	Human ACE (Angiotensin I Converting Enzyme) ELISA Kit	0.7850ng/mL	0.47ng/mL	Sandwich
E-EL-H0108	Human IFN-γ (Interferon Gamma) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H0493	Human AR (Androgen Receptor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0084	Human HGF (Hepatocyte Growth Factor) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H0011	Human BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H2473	Human CDKN1A(Cyclin Dependent Kinase Inhibitor 1A)ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich

Organ regeneration

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H2092	Human CCND1 (Cyclin-D1) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2242	Human REG1 α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0666	Human CTNNβ1 (Catenin, Beta 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0910	Human TP53 (Tumor Protein p53) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H1766	Human CFTR (Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H2243	Human REG4 (Regenerating Islet Derived Protein 4) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0111	Human VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich

Therapeutic cloning

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0060	Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2182	Human KRAS (V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0109	Human TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0111	Human VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0059	Human EGF (Epidermal Growth Factor) ELISA Kit	3.91250pg/mL	2.34pg/mL	Sandwich
E-EL-H0064	Human ET-1 (Endothelin 1) ELISA Kit	1.25-80pg/mL	0.75pg/mL	Sandwich
E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H1629	Human EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H0002	Human ACE (Angiotensin I Converting Enzyme) ELISA Kit	0.7850ng/mL	0.47ng/mL	Sandwich
E-EL-H0910	Human TP53 (Tumor Protein p53) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H0099	Human IL-2 (Interleukin 2) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich

Stem cells

E-EL-H1299 Human SCF (Stem Cell Factor) ELISA Kit 31.25-2000pg/mL 18.75pg/mL Sandwi E-EL-H0910 Human TP53 (Tumor Protein p53) ELISA Kit 78.13-5000pg/mL 46.88pg/mL Sandwi E-EL-H0102 Human IL-6 (Interleukin 6) ELISA Kit 7.81-500pg/mL 4.69pg/mL Sandwi E-EL-H0059 Human EGF (Epidermal Growth Factor) ELISA Kit 3.91-250pg/mL 2.34pg/mL Sandwi E-EL-H0930 Human SHH (Hedgehog Homolog, Sonic) ELISA Kit 78.13-5000pg/mL 46.88pg/mL Sandwi E-EL-H0930 Human SHH (Hedgehog Homolog, Sonic) ELISA Kit 62.5-4000pg/mL 46.88pg/mL Sandwi E-EL-H0930 Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit 62.5-4000pg/mL 46.88pg/mL Sandwi E-EL-H0109 Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit 7.81-500pg/mL 4.69pg/mL Sandwi E-EL-H01012 Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit 31.25-2000pg/mL 18.75pg/mL Sandwi E-EL-H0102 Human GCT4 (Octamer-binding protein 4) ELISA Kit 0.16-10ng/mL 0.09ng/mL Sandwi E-EL-H0084 Human IGF-1 (Insulin-Like Growth Factor) ELISA Kit 1.56-100ng/mL 0.94ng/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor) ELISA Kit 0.16-10ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.31-20ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human cmyc (c-myc Oncogene Product) ELISA Kit 0.31-20ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.31-20ng/mL 0.19ng/mL Sandwi E-EL-H0079 Human CMX2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit 0.31-20ng/mL 0.19ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56-100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 1.56-100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human IFN-y (Iriterferon Gamma) ELISA Kit 1.56-100ng/mL 0.93ng/mL Sandwi E-EL-H0114 Human IFN-y (Iriterferon Gamma) ELISA Kit 1.56-2000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.25-2000pg/mL 18.75pg/mL Sandwi	Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0102 Human TP53 (Tumor Protein p53) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi E-EL-H0102 Human IL-6 (Interleukin 6) ELISA Kit 7.81500pg/mL 4.69pg/mL Sandwi E-EL-H0059 Human EGF (Epidermal Growth Factor) ELISA Kit 3.91250pg/mL 2.34pg/mL Sandwi E-EL-H0300 Human SHH (Hedgehog Homolog, Sonic) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi E-EL-H0272 Human sCD44 (Soluble Cluster of Differentiation 44) ELISA Kit 62.54000pg/mL 37.5pg/mL Sandwi E-EL-H0109 Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit 7.81500pg/mL 4.69pg/mL Sandwi E-EL-H0109 Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit 7.81500pg/mL 4.69pg/mL Sandwi E-EL-H0112 Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H0012 Human OCT4 (Octamer-binding protein 4) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0084 Human HGF (Hepatocyte Growth Factor) ELISA Kit 1258000pg/mL 75pg/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H1607 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0179 Human CDXV2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit 0.3120ng/mL 0.38ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.1610ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human FN (Fibronectin) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.25200pg/mL 46.88pg/mL Sandwi	E-EL-H0483	Human bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	12.5800pg/mL	7.5pg/mL	Sandwich
E-EL-H0102 Human IL-6 (Interleukin 6) ELISA Kit 7.81–500pg/mL 4.69pg/mL Sandwi E-EL-H0059 Human EGF (Epidermal Growth Factor) ELISA Kit 3.91–250pg/mL 2.34pg/mL Sandwi E-EL-H0930 Human SHH (Hedgehog Homolog, Sonic) ELISA Kit 78.13–5000pg/mL 46.88pg/mL Sandwi E-EL-H0272 Human sCD44 (Soluble Cluster of Differentiation 44) ELISA Kit 62.5–4000pg/mL 37.5pg/mL Sandwi E-EL-H0109 Human TNF-q (Tumor Necrosis Factor Alpha) ELISA Kit 7.81–500pg/mL 4.69pg/mL Sandwi E-EL-H0101 Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit 31.25–2000pg/mL 18.75pg/mL Sandwi E-EL-H0012 Human OCT4 (Octamer-binding protein 4) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwi E-EL-H0084 Human HGF (Hepatocyte Growth Factor) ELISA Kit 125–8000pg/mL 75pg/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwi E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0600 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0759 Human c-Myc (C-glin Dependent Kinase Inhibitor 2A) ELISA Kit 0.63–40ng/mL 0.38ng/mL Sandwi E-EL-H0179 Human BC-KN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit 0.63–40ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human BC-(B-cell Leukemia/Lymphoma 2) ELISA Kit 0.16–10ng/mL 0.99ng/mL Sandwi E-EL-H0114 Human FN (Fibronectin) ELISA Kit 15.63–1000ng/mL 0.99ng/mL Sandwi E-EL-H0108 Human IFN-y (Interferon Gamma) ELISA Kit 15.63–1000ng/mL 0.99ng/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.25–2000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.13–5000pg/mL 46.88pg/mL Sandwi	E-EL-H1299	Human SCF (Stem Cell Factor) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0059 Human EGF (Epidermal Growth Factor) ELISA Kit 3.91~250pg/mL 2.34pg/mL Sandwi E-EL-H0930 Human SHH (Hedgehog Homolog, Sonic) ELISA Kit 78.13~5000pg/mL 46.88pg/mL Sandwi E-EL-H0272 Human sCD44 (Soluble Cluster of Differentiation 44) ELISA Kit 62.5~4000pg/mL 37.5pg/mL Sandwi E-EL-H0109 Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit 7.81~500pg/mL 4.69pg/mL Sandwi E-EL-H0012 Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit 31.25~2000pg/mL 18.75pg/mL Sandwi E-EL-H0012 Human CCT4 (Octamer-binding protein 4) ELISA Kit 0.16~10ng/mL 0.09ng/mL Sandwi E-EL-H0084 Human HGF (Hepatocyte Growth Factor) ELISA Kit 125~8000pg/mL 75pg/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit 1.56~100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.16~10ng/mL 0.99ng/mL Sandwi E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.31~20ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.31~20ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.31~20ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.31~20ng/mL 0.19ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56~100ng/mL 0.94ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56~100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 1.56~100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human FN (Fibronectin) ELISA Kit 1.56~100ng/mL 0.99ng/mL Sandwi E-EL-H0114 Human FN-γ (Interferon Gamma) ELISA Kit 1.56~3000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.25~2000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.25~2000pg/mL 46.88pg/mL Sandwi E-EL-H0111 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.13~5000pg/mL 46.88pg/mL Sandwi	E-EL-H0910	Human TP53 (Tumor Protein p53) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H0930 Human SHH (Hedgehog Homolog, Sonic) ELISA Kit 78.13–5000pg/mL 46.88pg/mL Sandwi E-EL-H2272 Human sCD44 (Soluble Cluster of Differentiation 44) ELISA Kit 62.5–4000pg/mL 37.5pg/mL Sandwi E-EL-H0109 Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit 7.81–500pg/mL 4.69pg/mL Sandwi E-EL-H0012 Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit 31.25–2000pg/mL 18.75pg/mL Sandwi E-EL-H0012 Human OCT4 (Octamer-binding protein 4) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwi E-EL-H0084 Human HGF (Hepatocyte Growth Factor) ELISA Kit 125–8000pg/mL 75pg/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwi E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0179 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit 0.31–20ng/mL 0.94ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwi E-EL-H0110 Human BC-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 1.56–100ng/mL 0.99ng/mL Sandwi E-EL-H01010 Human IFN-y (Interferon Gamma) ELISA Kit 15.63–1000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 11.25–2000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 11.25–2000pg/mL 18.75pg/mL Sandwi	E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H2272 Human sCD44 (Soluble Cluster of Differentiation 44) ELISA Kit 62.54000pg/mL 37.5pg/mL Sandwi E-EL-H0109 Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit 7.81500pg/mL 4.69pg/mL Sandwi E-EL-H0012 Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H0022 Human OCT4 (Octamer-binding protein 4) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0084 Human HGF (Hepatocyte Growth Factor) ELISA Kit 1258000pg/mL 75pg/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.3120ng/mL 0.38ng/mL Sandwi E-EL-H0179 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA (1.56100ng/mL 0.94ng/mL Sandwi E-EL-H01179 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 1.56100ng/mL 0.09ng/mL Sandwi E-EL-H0108 Human IFN-y (Interferon Gamma) ELISA Kit 1.5631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 18.75pg/mL Sandwi	E-EL-H0059	Human EGF (Epidermal Growth Factor) ELISA Kit	3.91250pg/mL	2.34pg/mL	Sandwich
E-EL-H0109 Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit 7.81500pg/mL 4.69pg/mL Sandwi E-EL-H0012 Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H1022 Human OCT4 (Octamer-binding protein 4) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0084 Human HGF (Hepatocyte Growth Factor) ELISA Kit 1258000pg/mL 75pg/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H1495 ELISA Kit 0.1610ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0600 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0104 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.38ng/mL Sandwi E-EL-H0114 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H0108 Human IFN-γ (Interferon Gamma) ELISA Kit 1.56100ng/mL 0.39ng/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A 1.562000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A 1.562000pg/mL 18.75pg/mL Sandwi E-EL-H0490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H0930	Human SHH (Hedgehog Homolog, Sonic) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H0012 Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H1022 Human OCT4 (Octamer-binding protein 4) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0084 Human HGF (Hepatocyte Growth Factor) ELISA Kit 1258000pg/mL 75pg/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0327 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit 0.3120ng/mL 0.38ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.1610ng/mL 0.99ng/mL Sandwi E-EL-H0108 Human IFN-y (Interferon Gamma) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 9.38pg/mL Sandwi E-EL-H0400 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H2272	Human sCD44 (Soluble Cluster of Differentiation 44) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H1022 Human OCT4 (Octamer-binding protein 4) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwi E-EL-H0084 Human HGF (Hepatocyte Growth Factor) ELISA Kit 125–8000pg/mL 75pg/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwi E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H2327 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit 0.63–40ng/mL 0.38ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwi E-EL-H0108 Human IFN-y (Interferon Gamma) ELISA Kit 15.63–1000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 1.25–2000pg/mL 18.75pg/mL Sandwi E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.13–5000pg/mL 46.88pg/mL Sandwi	E-EL-H0109	Human TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0084 Human HGF (Hepatocyte Growth Factor) ELISA Kit 125–8000pg/mL 75pg/mL Sandwi E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H2327 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit 0.6340ng/mL 0.38ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0108 Human IFN-y (Interferon Gamma) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H0012	Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0086 Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwi E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwi E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.31–20ng/mL 0.19ng/mL Sandwi E-EL-H2327 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit 0.63–40ng/mL 0.38ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56–100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.16–10ng/mL 0.09ng/mL Sandwi E-EL-H0108 Human IFN-y (Interferon Gamma) ELISA Kit 15.63–1000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 1.25–2000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.13–5000pg/mL 46.88pg/mL Sandwi	E-EL-H1022	Human OCT4 (Octamer-binding protein 4) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1607 Human CD34 (Cluster Of Differentiation 34) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H2327 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA (it 0.6340ng/mL 0.38ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0108 Human IFN-γ (Interferon Gamma) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 1.252000pg/mL 18.75pg/mL Sandwi E-EL-H0111 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H0084	Human HGF (Hepatocyte Growth Factor) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H1495 Human GDNF (Glial Cell Line Derived Neurotrophic Factor) E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit E-EL-H2327 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA (kit E-EL-H0179 Human FN (Fibronectin) ELISA Kit E-EL-H0179 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit E-EL-H0108 Human IFN-y (Interferon Gamma) ELISA Kit E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit E-EL-H0111 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit T8.135000pg/mL O.19ng/mL O.19ng/mL Sandwi D.19ng/mL O.19ng/mL O.38ng/mL Sandwi	E-EL-H0086	Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H1495 ELISA Kit E-EL-H0756 Human c-myc (c-myc Oncogene Product) ELISA Kit E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit E-EL-H2327 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA (Kit E-EL-H0179 Human FN (Fibronectin) ELISA Kit E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit E-EL-H0108 Human IFN-γ (Interferon Gamma) ELISA Kit E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) E-EL-H0111 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit D.19ng/mL D.19	E-EL-H1607	Human CD34 (Cluster Of Differentiation 34) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0060 Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi E-EL-H2327 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA 0.6340ng/mL 0.38ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0108 Human IFN-y (Interferon Gamma) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H1495	,	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2327 Human CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA (0.6340ng/mL 0.38ng/mL Sandwi E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0108 Human IFN-γ (Interferon Gamma) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H0756	Human c-myc (c-myc Oncogene Product) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0179 Human FN (Fibronectin) ELISA Kit 1.56100ng/mL 0.94ng/mL Sandwi E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0108 Human IFN-γ (Interferon Gamma) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) E-EL-H0111 ELISA Kit 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H0060	Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0114 Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi E-EL-H0108 Human IFN-γ (Interferon Gamma) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H2327		0.6340ng/mL	0.38ng/mL	Sandwich
E-EL-H0108 Human IFN-γ (Interferon Gamma) ELISA Kit 15.631000pg/mL 9.38pg/mL Sandwi E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H0179	Human FN (Fibronectin) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H0111 Human VEGF-A (Vascular Endothelial Cell Growth Factor A) 31.252000pg/mL 18.75pg/mL Sandwi E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H0114	Human Bcl-2 (B-cell Leukemia/Lymphoma 2) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0111 ELISA Kit E-EL-H5490 Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit 78.135000pg/mL 46.88pg/mL Sandwi	E-EL-H0108	Human IFN-γ (Interferon Gamma) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
	E-EL-H0111	,	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H5496 Human PROM1 (Prominin-1) ELISA Kit 0.3120ng/mL 0.19ng/mL Sandwi	E-EL-H5490	Human CXCR4 (Chemokine C-X-C-Motif Receptor) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
	E-EL-H5496	Human PROM1 (Prominin-1) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0666 Human CTNNβ1 (Catenin, Beta 1) ELISA Kit 0.1610ng/mL 0.09ng/mL Sandwi	E-EL-H0666	Human CTNNβ1 (Catenin, Beta 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H2139 Human SOX2 (Sex Determining Region Y Box Protein 2) ELISA 0.1610ng/mL 0.09ng/mL Sandwi	E-EL-H2139	,	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0011 Human BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit 62.54000pg/mL 37.5pg/mL Sandwi	E-EL-H0011	Human BMP-2 (Bone Morphogenetic Protein 2) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity 1	Test Method
E-EL-H0239	Human INHA (Inhibin A) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H0062	Human ENG (Endoglin) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0103	Human IL-10 (Interleukin 10) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0706	Human TERT (Telomerase Reverse Transcriptase) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0052	Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0010	Human BDNF (Brain Derived Neurotrophic Factor) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0100	Human IL-3 (Interleukin 3) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H0305	Human ALCAM (Activated Leukocyte Cell Adhesion Molecule) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H5461	Human FGF8 (Fibroblast Growth Factor 8) ELISA Kit	12.5800pg/mL	7.5pg/mL	Sandwich
E-EL-H0356	Human FGFR3 (Fibroblast Growth Factor Receptor 3) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H2092	Human CCND1 (Cyclin-D1) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0878	Human Flt3 (FMS Like Tyrosine Kinase 3) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0135	Human sICAM-1/CD54 (Soluble Intercellular Adhesion Molecule 1) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H1603	Human VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H2473	Human CDKN1A (Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H0094	Human LIF (Leukemia Inhibitory Factor) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H2182	Human KRAS (V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0067	Human FAS/CD95 (Factor Related Apoptosis) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0079	Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit	39.062500pg/mL	23.44pg/mL	Sandwich
E-EL-H1945	Human NOG (Noggin) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H2304	Human TP63 (Tumor Protein p63) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0985	Human PDGFRβ (Platelet Derived Growth Factor Receptor Beta) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H2132	Human STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H1361	Human PPAR-γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0552	Human FGF4 (Fibroblast Growth Factor 4) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H0039	Human CNTF (Ciliary Neurotrophic Factor) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H2339	Human TGF-β3 (Transforming Growth Factor β3) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H0149	Human IL-1β (Interleukin 1 Beta) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0345	Human FGFR1 (Fibroblast Growth Factor Receptor 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H2239	Human JAK2 (Janus Kinase 2) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H1349	Human SOX9 (Sex Determining Region Y Box Protein 9) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2214	Human PDGFRα (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0099	Human IL-2 (Interleukin 2) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H0158	Human FGFR2 (Fibroblast Growth Factor Receptor 2) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1659	Human HLA-A (Leukocyte Antigen A) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H1629	Human EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit	1258000pg/mL	75pg/mL	Sandwich
E-EL-H0224	Human 5-NT (5-Nucleotidase) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0752	Human CD4 (Cluster of Differentiation 4) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0171	Human TF (Transferrin) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H0262	Human ADA (Adenosine Deaminase) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H1030	Human IL-2Rα (Interleukin 2 Receptor Alpha) ELISA Kit	93.756000pg/mL	56.25pg/mL	Sandwich
E-EL-H1896	Human NT-3 (Neurotrophin-3) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H1614	Human BRAF (B-Raf Proto Oncogene Serine/Threonine Protein Kinase) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0013	Human BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1205	Human NGF (Nerve Growth Factor) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H1932	Human RB1 (Retinoblastoma Protein 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0092	Human FGF7 (Fibroblast Growth Factor 7) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1588	Human TPO (Thrombopoietin) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0985	Human PDGFRβ (Platelet Derived Growth Factor Receptor Beta) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0017	Human CASP3 (Caspase 3) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0048	Human IL-8 (Interleukin 8) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1888	Human GFAP (Glial Fibrillary Acidic Protein) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0562	Human BAX (Bcl-2 Associated X Protein) ELISA Kit	0.7850ng/mL	0.47ng/mL	Sandwich
E-EL-H5462	Human FGF10 (Fibroblast Growth Factor 10) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0072	Human AAP (Alanine Aminopeptidase) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H1640	Human PECAM1/CD31 (Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0880	Human GATA2 (GATA Binding Protein 2) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2246	Human PTPN11 (Protein Tyrosine Phosphatase, Non Receptor Type 11) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0075	Human Flt3L (FMS Like Tyrosine Kinase 3 Ligand) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H0195	Human NCAD (Neural Cadherin) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H5610	Human CD19 (Cluster Of Differentiation 19) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2157	Human VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0101	Human IL-4 (Interleukin 4) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H2069	Human CTLA4 (Cytotoxic T Lymphocyte Associated Antigen 4) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1094	Human VIM (Vimentin) ELISA Kit	7.81500ng/mL	4.69ng/mL	Sandwich
E-EL-H1587	Human TGF-β2 (Transforming Growth Factor β2) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0879	Human GATA1 (GATA Binding Protein 1) ELISA Kit	0.6340ng/mL	0.38ng/mL	Sandwich
E-EL-H1386	Human NF-кВ p105 (Nuclear Factor Карра В p105) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1347	Human OPN (Osteopontin) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0057	Human DKK1 (Dickkopf Related Protein 1) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0081	Human GM-CSF (Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit	7.81500pg/mL	4.69pg/mL	Sandwich
E-EL-H0068	Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H1277	Human HIF-1α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0925	Human RUNX2 (Runt-Related Transcription Factor 2) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0035	Human sCD40L (Soluble Cluster of Differentiation 40 Ligand) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H0259	Human ACTα2 (Actin Alpha 2, Smooth Muscle) ELISA Kit	3.13200ng/mL	1.88ng/mL	Sandwich
E-EL-H0659	Human CASP8 (Caspase 8) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0035	Human sCD40L (Soluble Cluster of Differentiation 40 Ligand) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H2086	Human ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0530	Human ABCG2 (ATP Binding Cassette Transporter G2) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1117	Human FGF9 (Fibroblast Growth Factor 9) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H1592	Human IL7R (Interleukin 7 Receptor) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H0097	Human MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H0425	Human IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2075	Human CK-19/KRT19 (Cytokeratin 19) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H0066	Human EPO (Erythropoietin) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H1846	Human PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H5636	Human CD3e (Cluster of Differentiation 3e) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H0435	Human IGF-2 (Insulin Like Growth Factor 2) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2130	Human STAT1 (Signal Transducer and Activator of Transcription 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0513	Human MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H0505	Human MAP2K1 (Mitogen Activated Protein Kinase Kinase 1) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H0547	Human BCR (B-Cell Receptor) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H1539	Human TLR4 (Toll-Like Receptor 4) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H5559	Human JAK3 (Janus Kinase 3) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H2247	Human OSM (Oncostatin M) ELISA Kit	15.631000pg/mL	9.38pg/mL	Sandwich
E-EL-H2399	Human PCNA (Proliferating Cell Nuclear Antigen) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H2000	Human GSK3β (Glycogen Synthase Kinase 3 Beta) ELISA Kit	0.6340ng/mL	0.38ng/mL	Sandwich
E-EL-H2403	Human MCAM (Melanoma Cell Adhesion Molecule) ELISA Kit	0.6340ng/mL	0.38ng/mL	Sandwich
E-EL-H0493	Human AR (Androgen Receptor) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H5022	Human IL-11 (Interleukin 11) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H5634	Human TGFBR2 (TGF-beta receptor type 2) ELISA Kit	0.7850ng/mL	0.47ng/mL	Sandwich
E-EL-H1451	Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit	31.252000pg/mL	18.75pg/mL	Sandwich
E-EL-H2168	Human VWF (Von Willebrand Factor) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H1104	Human FOXP3 (Forkhead Box Protein P3) ELISA Kit	0.3120ng/mL	0.19ng/mL	Sandwich
E-EL-H1133	Human ERα (Estrogen Receptor Alpha) ELISA Kit	62.54000pg/mL	37.5pg/mL	Sandwich
E-EL-H2171	Human HBβ (Hemoglobin Beta) ELISA Kit	1.56100ng/mL	0.94ng/mL	Sandwich
E-EL-H2203	Human ITGB3 (Integrin Beta 3) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H5483	Human HES1 (Transcription factor HES-1) ELISA Kit	78.135000pg/mL	46.88pg/mL	Sandwich
E-EL-H1046	Human GP4/CD36 (Platelet Membrane Glycoprotein IV) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich
E-EL-H2407	Human IκBα (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit	0.1610ng/mL	0.09ng/mL	Sandwich





International Quality and Management Certification ISO9001:2008



Elabscience Biotechnology Inc.

Web: www.elabscience.com

Tel: 1-240-252-7368 Fax: 1-240-252-7376

Email: orders@elabscience.com