

Catalog Number	Product Name
E-BC-D001	Angiotensin Converting Enzyme 2 (ACE2) Inhibitor Screening Kit
E-BC-D002	Cyclooxygenase-2 (COX-2) Inhibitor Screening Kit
E-BC-D003	Matrix Metalloproteinase 3 (MMP-3) Inhibitor Screening Kit
E-BC-D004	Neuraminidases (NA) Inhibitor Screening Assay Kit
E-BC-D007	Dipeptidyl Peptidase IV (DPP4) Inhibitor Screening Assay Kit
E-BC-D008	Ferroptosis Suppressor Protein-1 (FSP-1) Inhibitor Screening Kit
E-BC-D011	3-Hydroxy-3-methylglutaryl-CoA Reductase (HMGR) Inhibitor Screening Kit
E-BC-D012	Monoamine Oxidase A (MAO-A) Inhibitor Screening Kit
E-BC-D016	Angiotensin Converting Enzyme 1 (ACE1) Inhibitor Screening Kit
E-BC-D017	Alpha-Glucosidase Inhibitor Screening Kit
E-BC-D018	Alpha-Amylase Inhibitor Screening Kit
E-BC-D019	Xanthine Oxidase (XOD) Inhibitor Screening Kit
E-BC-D022	Acetylcholinesterase (AChE) Inhibitor Screening Kit
E-BC-E001	Animal Tissue Mitochondrial Extraction Assay Kit
E-BC-E002	Total Protein Extraction Assay Kit
E-BC-E005	Plant Mitochondrial Extraction Assay Kit(Enzyme Method)
E-BC-E006	Cell Mitochondrial Extraction Assay Kit
E-BC-F001	Hydrogen Peroxide (H ₂ O ₂) Fluorometric Assay Kit
E-BC-F003	Lipid Peroxide (LPO) Fluorometric Assay Kit
E-BC-F004	ATP/ADP Ratio Chemiluminescence Assay kit
E-BC-F005	Reactive Oxygen Species (ROS) Fluorometric Assay Kit (Red)
E-BC-F006	Catalase (CAT) Activity Fluorometric Assay Kit
E-BC-F007	Malondialdehyde (MDA) Fluorometric Assay Kit
E-BC-F008	Mitochondrial Superoxide Fluorometric Assay Kit
E-BC-F009	Adenosine Diphosphate (ADP) Fluorometric Assay Kit
E-BC-F012	Peroxidase (POD) Activity Fluorometric Assay Kit
E-BC-F013	Myeloperoxidase (MPO) Peroxidation Activity Fluorometric Assay Kit
E-BC-F018	Uric Acid (UA) Fluorometric Assay Kit
E-BC-F019	Xanthine Oxidase (XOD) Activity Fluorometric Assay Kit
E-BC-F020	Nitric Oxide Synthase (NOS) Activity Fluorometric Assay Kit
E-BC-F021	Nitric Oxide (NO) Fluorometric Assay Kit
E-BC-F022	Ferroptosis Suppressor Protein-1 (FSP-1) Activity Fluorometric Assay Kit
E-BC-F023	Cytochrome P450 2D6 (CYP2D6) Activity Fluorometric Assay Kit
E-BC-F024	Xanthine/Hypoxanthine Fluorometric Assay Kit
E-BC-F025	N-Acetylcysteine (NAC) Fluorometric Assay Kit
E-BC-F032	Total Cholesterol And Cholesteryl Ester Fluorometric Assay Kit

E-BC-F033	Triglyceride (TG) Fluorometric Assay Kit
E-BC-F035	High-Density Lipoprotein Cholesterol (HDL-C) Fluorometric Assay Kit
E-BC-F037	Glucose (GLU) Fluorometric Assay Kit
E-BC-F038	Alanine Aminotransferase (ALT/GPT) Activity Fluorometric Assay Kit
E-BC-F039	Free Fatty Acids (NEFA/FFA) Fluorometric Assay Kit
E-BC-F041	Glucose Uptake Fluorometric Assay Kit
E-BC-F042	Sucrose Fluorometric Assay Kit
E-BC-F043	Aspartate Aminotransferase (AST/GOT) Activity Fluorometric Assay Kit
E-BC-F044	Aldehyde Dehydrogenase (ALDH) Fluorometric Activity Assay Kit
E-BC-F045	Total Glutathione (T-GSH) And Reduced Glutathione (GSH) Fluorometric Assay Kit
E-BC-F046	Acetyl Coenzyme A (Acetyl-CoA) Fluorometric Assay Kit
E-BC-F047	α -Ketoglutarate (α -KG) Fluorometric Assay Kit
E-BC-F048	γ -Aminobutyric Acid (GABA) Fluorometric Assay Kit
E-BC-F049	Lactulose Fluorometric Assay Kit
E-BC-F050	Cyclooxygenase (COX) Activity Fluorometric Assay Kit
E-BC-F051	Histone Deacetylase (HDAC) Activity Fluorometric Assay Kit
E-BC-F052	Cathepsin D (CatD) Activity Fluorometric Assay Kit
E-BC-F053	Dipeptidyl Peptidase IV (DPP4) Activity Fluorometric Assay Kit
E-BC-F054	Angiotensin I Converting Enzyme 2 (ACE2) Activity Fluorescence Assay Kit
E-BC-F055	Choline/Acetylcholine Fluorometric Assay Kit
E-BC-F056	Sirtuin 1 (SIRT-1) Activity Assay Kit
E-BC-F058	Pyruvate Fluorometric Assay Kit
E-BC-F059	Cathepsin B Activity Fluorometric Assay Kit
E-BC-F061	Matrix Metalloproteinase 3 (MMP-3) Activity Fluorometric Assay Kit
E-BC-F062	Singlet Oxygen(1O_2) Fluorometric Assay Kit
E-BC-F063	Cellular Intracellular Oxygen Fluorometric Assay Kit
E-BC-F064	Mitochondrial Permeability Transition Pore (mPTP) Fluorometric Assay Kit
E-BC-F065	Adipogenesis Fluorometric Assay Kit
E-BC-F066	Cystine Uptake Fluorometric Assay Kit
E-BC-F067	Fatty Acid Uptake Fluorometric Assay Kit
E-BC-F069	Extracellular Acidification Rate (ECAR) Fluorometric Assay Kit
E-BC-F070	Enhanced Oxygen Consumption Rate (OCR) Fluorometric Assay Kit
E-BC-F071	Dual Luciferase Reporter Gene Luminescence Assay Kit (Flash Type)
E-BC-F072	Firefly Luciferase Reporter Gene Luminescence Assay Kit (Flash Type)
E-BC-F073	Intracellular Acidification Rate Fluorometric Assay Kit
E-BC-F074	Dual Luciferase Reporter Gene Luminescence Assay Kit (Glow Type)
E-BC-F075	Firefly Luciferase Reporter Gene Luminescence Assay Kit (Glow Type)
E-BC-F076	β -Secretase 1 (BACE 1) Activity Fluorometric Assay Kit
E-BC-F077	Lipoxygenase (LOX) Activity Fluorometric Assay Kit

E-BC-F078	Mitochondrial Stress Fluorometric Assay Kit
E-BC-F079	Sphingomyelin (SM) Fluorometric Assay Kit
E-BC-F080	Fructose Fluorometric Assay Kit
E-BC-F081	Glycerol-3-phosphate (G3P) Fluorometric Assay Kit
E-BC-F082	Fatty Acid Amide Hydrolase (FAAH) Activity Fluorometric Assay Kit
E-BC-F083	Neuraminidase (NA) Activity Fluorometric Assay Kit
E-BC-F084	Glycolysis Stress Fluorometric Assay Kit
E-BC-F086	Cellular Senescence Fluorometric Assay Kit
E-BC-F087	Fructose-6-Phosphate (F-6-P) Fluorometric Assay Kit
E-BC-F089	Cholesterol Uptake Fluorometric Assay Kit
E-BC-F091	Tryptophan (Trp) Fluorometric Assay Kit
E-BC-F093	Monoamine Oxidase (MAO) Activity Fluorometric Assay Kit
E-BC-F094	Glutamate Oxidase (GLOD) Activity Fluorometric Assay Kit
E-BC-F100	Fluo-4 Calcium Fluorometric Assay Kit
E-BC-F101	Cell Ferrous Iron (Fe ²⁺) Fluorometric Assay Kit
E-BC-F102	Cellular Cuprous Fluorometric Assay Kit
E-BC-F103	Cell Hydrogen Sulfide (H ₂ S) Fluorometric Assay Kit
E-BC-F106	Lead Ion (Pb ²⁺) Fluorometric Assay Kit
E-BC-F107	Intracellular Cuprous Ion (Cu ⁺) Fluorometric Assay Kit
E-BC-F108	Cell Ferrous (Fe ²⁺) Fluorometric Assay Kit
E-BC-F140	Dihydroorotate Dehydrogenase (DHODH) Activity Fluorometric Assay Kit
E-BC-F151	Dihydroxyacetone Phosphate (DHAP) Fluorometric Assay Kit
E-BC-F152	6-Phosphogluconic Acid (6-PGA) Fluorometric Assay Kit
E-BC-F153	Citrate Fluorometric Assay Kit
E-BC-F154	Phosphocreatine (PCr) Fluorometric Assay Kit
E-BC-F170	NADPH Fluorometric Assay Kit
E-BC-F171	Glycogen Fluorometric Assay Kit (Glycosidase Method)
E-BC-F200	Cell Viability Chemiluminescence Assay Kit
E-BC-F201	Enhanced ATP Chemiluminescence Assay Kit
E-BC-F202	Lysosomal Activity Fluorometric Assay Kit
E-BC-F300	ATP Chemiluminescence Assay Kit (Double Reagent)
E-BC-K001-M	Inhibition And Production Of Superoxide Anionic Colorimetric Assay Kit (WST-1 Method)
E-BC-K002-M	D-Lactic Acid/Lactate(LA)Colorimetric Assay Kit
E-BC-K003-M	Angiotensin Converting Enzyme (ACE1) Activity Assay Kit
E-BC-K003-S	Angiotensin Converting Enzyme (ACE1) Activity Assay Kit
E-BC-K004-M	Free Cholesterol (FC) Colorimetric Assay Kit
E-BC-K005-M	β-Amylase Activity Assay Kit
E-BC-K006-M	α-Amylase And β-Amylase Activity Assay Kit
E-BC-K007-M	α-Amylase Activity Assay Kit

E-BC-K008-M	Monoamine Oxidase (MAO) Activity Assay Kit
E-BC-K008-S	Monoamine Oxidase (MAO) Activity Assay Kit
E-BC-K009-M	Alkaline Phosphatase (ALP) Activity Assay Kit (PNPP Method)
E-BC-K010-M	Acid Phosphatase (ACP) Activity Assay Kit (PNPP Method)
E-BC-K011-M	Glucose-6-phosphate (G6P) Colorimetric Assay Kit
E-BC-K013-M	Non-esterified Free Fatty Acids (NEFA/FFA) Colorimetric Assay Kit
E-BC-K013-S	Non-esterified Free Fatty Acids (NEFA/FFA) Colorimetric Assay Kit
E-BC-K014	Non-esterified Free Fatty Acids (NEFA) Colorimetric Assay Kit
E-BC-K016-M	Uric Acid (UA) Colorimetric Assay Kit
E-BC-K016-S	Uric Acid (UA) Colorimetric Assay Kit
E-BC-K018-S	D-Xylose Colorimetric Assay Kit
E-BC-K019-M	Total Superoxide Dismutase (T-SOD) Activity Assay Kit (Hydroxylamine Method)
E-BC-K019-S	Total Superoxide Dismutase (T-SOD) Activity Assay Kit (Hydroxylamine Method)
E-BC-K020-M	Total Superoxide Dismutase (T-SOD) Activity Assay Kit (WST-1 Method)
E-BC-K022-M	CuZn/Mn Superoxide Dismutase (CuZn-SOD/Mn-SOD) Activity Assay Kit (Hydroxylamine Method)
E-BC-K022-S	CuZn/Mn Superoxide Dismutase (CuZn-SOD/Mn-SOD) Activity Assay Kit (Hydroxylamine Method)
E-BC-K025-M	Malondialdehyde (MDA) Colorimetric Assay Kit (TBA Method)
E-BC-K025-S	Malondialdehyde (MDA) Colorimetric Assay Kit (TBA Method)
E-BC-K027-M	Malondialdehyde (MDA) Colorimetric Assay Kit (Plant Samples)
E-BC-K027-S	Malondialdehyde (MDA) Colorimetric Assay Kit (Plant Samples)
E-BC-K030-M	Reduced Glutathione (GSH) Colorimetric Assay Kit
E-BC-K030-S	Reduced Glutathione (GSH) Colorimetric Assay Kit
E-BC-K031-M	Catalase (CAT) Activity Assay Kit
E-BC-K031-S	Catalase (CAT) Activity Assay Kit
E-BC-K033-M	Vitamin E (VE) Colorimetric Assay Kit
E-BC-K033-S	Vitamin E (VE) Colorimetric Assay Kit
E-BC-K034-M	Vitamin C (VC) Colorimetric Assay Kit
E-BC-K034-S	Vitamin C (VC) Colorimetric Assay Kit
E-BC-K035-M	Nitric Oxide (NO) Colorimetric Assay Kit
E-BC-K035-S	Nitric Oxide (NO) Colorimetric Assay Kit
E-BC-K041-M	Maltase Activity Assay Kit
E-BC-K042-M	Hydroxyl Free Radical Scavenging Capacity Colorimetric Assay Kit (Fenton Method)
E-BC-K043-S	L-Lactic Acid (LA) Colorimetric Assay Kit (Whole Blood Samples)
E-BC-K044-M	L-Lactic Acid (LA) Colorimetric Assay Kit
E-BC-K044-S	L-Lactic Acid (LA) Colorimetric Assay Kit
E-BC-K046-M	Lactate Dehydrogenase (LDH) Activity Assay Kit
E-BC-K046-S	Lactate Dehydrogenase (LDH) Activity Assay Kit
E-BC-K052-M	Cholinesterase (CHE) Activity Assay Kit
E-BC-K052-S	Cholinesterase (CHE) Activity Assay Kit

E-BC-K055-M	Total Amino Acids (T-AA) Colorimetric Assay Kit
E-BC-K055-S	Total Amino Acids (T-AA) Colorimetric Assay Kit
E-BC-K056-M	Glucose-6-Phosphate Dehydrogenase (G-6-PD) Activity Assay Kit
E-BC-K057-M	Albumin (ALB) Colorimetric Assay Kit (Bromocresol Green Method)
E-BC-K057-S	Albumin (ALB) Colorimetric Assay Kit (Bromocresol Green Method)
E-BC-K061-M	Hydroxyproline (HYP) Colorimetric Assay Kit (Alkali Hydrolysis Method)
E-BC-K061-S	Hydroxyproline (HYP) Colorimetric Assay Kit (Alkali Hydrolysis Method)
E-BC-K062-M	Hydroxyproline (HYP) Colorimetric Assay Kit (Acid Hydrolysis Method)
E-BC-K062-S	Hydroxyproline (HYP) Colorimetric Assay Kit
E-BC-K064-M	β -N-Acetylglucosaminidase (NAG) Activity Assay Kit
E-BC-K068-M	Sialic Acid (SA) Colorimetric Assay Kit
E-BC-K068-S	Sialic Acid (SA) Colorimetric Assay Kit
E-BC-K069-M	Glycosylated Serum Protein (GSP) Colorimetric Assay Kit
E-BC-K070-S	Nitrite Colorimetric Assay Kit
E-BC-K071-M	Total Iron Binding Capacity (TIBC) Colorimetric Assay Kit
E-BC-K071-S	Total Iron Binding Capacity (TIBC) Colorimetric Assay Kit
E-BC-K073-S	Glycogen Colorimetric Assay Kit (Liver/Muscle Samples)
E-BC-K074-M	Myeloperoxidase (MPO) Activity Assay Kit
E-BC-K074-S	Myeloperoxidase (MPO) Activity Assay Kit
E-BC-K083-M	α -Ketoglutarate Dehydrogenase (α -KGDH) Activity Assay Kit
E-BC-K091-M	Alkaline Phosphatase (ALP) Activity Assay Kit
E-BC-K091-S	Alkaline Phosphatase (ALP) Activity Assay Kit
E-BC-K094-M	Acid Phosphatase (ACP) Activity Assay Kit
E-BC-K094-S	Acid Phosphatase (ACP) Activity Assay Kit
E-BC-K096-M	Glutathione Peroxidase (GSH-Px) Activity Assay Kit
E-BC-K096-S	Glutathione Peroxidase (GSH-Px) Activity Assay Kit
E-BC-K097-M	Total Glutathione (T-GSH)/Oxidized Glutathione (GSSG) Colorimetric Assay Kit
E-BC-K097-S	Total Glutathione (T-GSH)/Oxidized Glutathione (GSSG) Colorimetric Assay Kit
E-BC-K099-M	Glutathione Reductase (GR) Activity Assay Kit
E-BC-K099-S	Glutathione Reductase (GR) Activity Assay Kit
E-BC-K102-M	Hydrogen Peroxide (H ₂ O ₂) Colorimetric Assay Kit
E-BC-K102-S	Hydrogen Peroxide (H ₂ O ₂) Colorimetric Assay Kit
E-BC-K103-M	Calcium (Ca) Colorimetric Assay Kit
E-BC-K109-M	Total Cholesterol (TC) Colorimetric Assay Kit (Single Reagent, COD-PAP Method)
E-BC-K109-S	Total Cholesterol (TC) Colorimetric Assay Kit (Single Reagent, COD-PAP Method)
E-BC-K1100-M	Formate Colorimetric Assay Kit
E-BC-K1102-M	α -Hydroxybutyrate Dehydrogenase (α -HBDH) Activity Colorimetric Assay Kit
E-BC-K1103-M	Tissue Total Cholesterol (TC) Colorimetric Assay Kit
E-BC-K1104-M	Tissue Free Cholesterol (FC) Colorimetric Assay Kit

E-BC-K1105-M	Biotin Colorimetric Assay Kit
E-BC-K1106-M	Oxalate (Oxalic Acid) Colorimetric Assay Kit (Extraction method)
E-BC-K1107-M	Carnitine Palmitoyl Transferase-I (CPT-I) Activity Colorimetric Assay Kit
E-BC-K1109-M	Fatty Acid Synthase (FAS) Activity Colorimetric Assay Kit
E-BC-K1150-M	6-Phosphogluconate Dehydrogenase (6-PGDH) Activity Colorimetric Assay Kit
E-BC-K1151-M	Transketolase(TKT) Activity Colorimetric Assay Kit
E-BC-K1154-M	Glycogen Synthase(GCS) Activity Colorimetric Assay Kit
E-BC-K1155-M	Sorbitol Dehydrogenase (SDH) Activity Colorimetric Assay Kit
E-BC-K1156-M	Sorbitol Colorimetric Assay Kit
E-BC-K1157-M	Fructokinase (FRK) Activity Colorimetric Assay Kit
E-BC-K117-M	Protein Carbonyl Colorimetric Assay Kit (Tissue And Serum Samples)
E-BC-K117-S	Protein Carbonyl Colorimetric Assay Kit (Tissue And Serum Samples)
E-BC-K118-S	Glutamic Acid Colorimetric Assay Kit
E-BC-K1200-M	Amylose Colorimetric Assay Kit
E-BC-K1201-M	Amylopectin Colorimetric Assay Kit
E-BC-K1202-M	Soluble Starch Synthase (SSS) Activity Colorimetric Assay Kit
E-BC-K1203-M	Granule-Bound Starch Synthase (GBSS) Activity Colorimetric Assay Kit
E-BC-K1204-M	Starch Branching Enzyme (SBE) Activity Colorimetric Assay Kit
E-BC-K1205-M	Starch Phosphorylase (SP) Activity Colorimetric Assay Kit
E-BC-K1206-M	Starch Debranching Enzyme (DBE) Activity Colorimetric Assay Kit
E-BC-K1207-M	Sucrose Synthase- II (SS- II) Activity Colorimetric Assay Kit
E-BC-K1208-M	Sucrose Synthase- I (SS- I) Activity Colorimetric Assay Kit
E-BC-K1209-M	Starch Colorimetric Assay Kit
E-BC-K1210-M	Phytic Acid Colorimetric Assay Kit
E-BC-K1211-M	Phytase Activity Colorimetric Assay Kit
E-BC-K1217-M	Cellulase (CL) Activity Colorimetric Assay Kit
E-BC-K1223-M	β -1,3-Glucanase Activity Colorimetric Assay Kit
E-BC-K1224-M	Chitinase Activity Colorimetric Assay Kit
E-BC-K1225-M	Dehydroascorbic Acid (DHA) Colorimetric Assay Kit
E-BC-K122-S	H ⁺ K ⁺ -ATPase Activity Assay Kit
E-BC-K125-M	Choline Acetyltransferase (ChAT) Activity Assay Kit (Tissue Samples)
E-BC-K125-S	Choline Acetyltransferase (ChAT) Activity Assay Kit (Tissue Samples)
E-BC-K126-M	γ -Glutamyl Transferase (GGT/ γ -GT) Activity Assay Kit
E-BC-K1300-M	Alanine Aminotransferase (ALT/GPT) Activity Colorimetric Assay Kit
E-BC-K1301-M	Aspartate Aminotransferase (AST/GOT) Activity Colorimetric Assay Kit (Aspartate Substrate Method)
E-BC-K1302-M	Lactate Dehydrogenase (LDH) Activity Colorimetric Assay Kit (Lactate Substrate Method)
E-BC-K1304-M	Alkaline Phosphatase (ALP) Activity Colorimetric Assay Kit(Method of NPP-AMP)
E-BC-K130-M	Pyruvic Acid Colorimetric Assay Kit
E-BC-K130-S	Pyruvic Acid Colorimetric Assay Kit

E-BC-K131-M	Lactase Activity Assay Kit
E-BC-K134-M	Fructose Colorimetric Assay Kit
E-BC-K135-M	Nitric Oxide (NO) Colorimetric Assay Kit (Nitrate Reductase Method)
E-BC-K136-M	Total Antioxidant Capacity (T-AOC) Colorimetric Assay Kit
E-BC-K136-S	Total Antioxidant Capacity (T-AOC) Colorimetric Assay Kit
E-BC-K137-M	Zinc (Zn) Colorimetric Assay Kit
E-BC-K138-F	Reactive Oxygen Species (ROS) Fluorometric Assay Kit (Green)
E-BC-K139-M	Iron Colorimetric Assay Kit
E-BC-K139-S	Iron Colorimetric Assay Kit
E-BC-K1401-M	Tyrosine (Tyr) Colorimetric Assay Kit
E-BC-K1402-M	Tyrosinase Activity Colorimetric Assay Kit
E-BC-K1403-M	Protein Thiol Colorimetric Assay Kit
E-BC-K145-M	Blood Ammonia Colorimetric Assay Kit
E-BC-K145-S	Blood Ammonia Colorimetric Assay Kit
E-BC-K149-M	Mitochondrial Complex I (NADH-CoQ Reductase) Activity Assay Kit
E-BC-K150-M	Mitochondrial Complex II Activity Assay Kit
E-BC-K151-M	Mitochondrial Complex III (Coenzyme Q-Cytochrome C Reductase) Activity Assay Kit
E-BC-K152-M	Mitochondrial Complex IV (Cytochrome C Oxidase) Activity Assay Kit
E-BC-K153-M	Mitochondrial Complex V (F ₀ F ₁ -ATPase/ATP Synthase) Activity Assay Kit
E-BC-K157-M	ATP Colorimetric Assay Kit
E-BC-K157-S	ATP Colorimetric Assay Kit
E-BC-K161-S	Sucrose Colorimetric Assay Kit
E-BC-K162-M	Magnesium (Mg) Colorimetric Assay Kit
E-BC-K162-S	Magnesium (Mg) Colorimetric Assay Kit
E-BC-K165-M	Biuret Protein Colorimetric Assay Kit
E-BC-K165-S	Biuret Protein Colorimetric Assay Kit
E-BC-K168-M	Bradford Protein Colorimetric Assay Kit
E-BC-K168-S	Bradford Protein Colorimetric Assay Kit
E-BC-K171-M	Total Carbonyl Colorimetric Assay Kit
E-BC-K171-S	Total Carbonyl Colorimetric Assay Kit
E-BC-K174-M	Acetylcholinesterase (AChE) Activity Assay Kit
E-BC-K176-M	Lipid Peroxide (LPO) Colorimetric Assay Kit
E-BC-K177-M	Proline (Pro) Colorimetric Assay Kit
E-BC-K177-S	Proline (Pro) Colorimetric Assay Kit
E-BC-K178-M	Citrate Synthase (CS) Activity Assay Kit
E-BC-K181-M	Total Bile Acid (TBA) Colorimetric Assay Kit
E-BC-K183-M	Urea (BUN) Colorimetric Assay Kit (Urease Method)
E-BC-K183-S	Urea (BUN) Colorimetric Assay Kit (Urease Method)
E-BC-K188-M	Creatinine (Cr) Colorimetric Assay Kit (Sarcosine Oxidase Method)

E-BC-K189-M	Chlorine (Cl) Colorimetric Assay Kit
E-BC-K196-M	5'-Nucleotidase (5'-NT) Activity Assay Kit
E-BC-K197-M	Adenosine Deaminase (ADA) Activity Assay Kit
E-BC-K205-M	Low-Density Lipoprotein Cholesterol (LDL-C) Colorimetric Assay Kit (Double Reagents)
E-BC-K205-S	Low-Density Lipoprotein Cholesterol (LDL-C) Colorimetric Assay Kit (Double Reagents)
E-BC-K207-M	Sodium (Na) Colorimetric Assay Kit
E-BC-K207-S	Sodium (Na ⁺) Colorimetric Assay Kit
E-BC-K212-M	Ca ²⁺ -ATPase Activity Assay Kit
E-BC-K212-S	Ca ²⁺ -ATPase Activity Assay Kit
E-BC-K219-M	Total Antioxidant Capacity (T-AOC) Colorimetric Assay Kit (ABTS, Enzyme Method)
E-BC-K221-M	High-Density Lipoprotein Cholesterol (HDL-C) Colorimetric Assay Kit (Double Reagents)
E-BC-K225-M	Total Antioxidant Capacity (T-AOC) Colorimetric Assay Kit (FRAP Method)
E-BC-K227-M	Peroxidase (POD) Activity Assay Kit (Plant samples)
E-BC-K227-S	Peroxidase (POD) Activity Assay Kit (Plant Samples)
E-BC-K234-M	Glucose (Glu) Colorimetric Assay Kit (GOD-POD Method)
E-BC-K234-S	Glucose (Glu) Colorimetric Assay Kit (GOD-POD Method)
E-BC-K235-M	Alanine Aminotransferase (ALT/GPT) Activity Assay Kit
E-BC-K235-S	Alanine Aminotransferase (ALT/GPT) Activity Assay Kit
E-BC-K236-M	Aspartate Aminotransferase (AST/GOT) Activity Assay Kit
E-BC-K236-S	Aspartate Aminotransferase (AST/GOT) Activity Assay Kit
E-BC-K245-M	Phosphorus (Pi) Colorimetric Assay Kit (Phospho Molybdate Method)
E-BC-K245-S	Phosphorus (Pi) Colorimetric Assay Kit (Phospho Molybdate Method)
E-BC-K252-M	Urine Protein Colorimetric Assay Kit
E-BC-K259-M	Polyphenol Oxidase (PPO) Activity Assay Kit
E-BC-K259-S	Polyphenol Oxidase (PPO) Activity Assay Kit
E-BC-K261-M	Triglyceride (TG) Colorimetric Assay Kit (Single Reagent, GPO-PAP Method)
E-BC-K261-S	Triglyceride (TG) Colorimetric Assay Kit (Single Reagent, GPO-PAP Method)
E-BC-K265-M	Total Sulfhydryl Group/Total Thiol (-SH) Colorimetric Assay Kit
E-BC-K271-M	Total Antioxidant Capacity (T-AOC) Colorimetric Assay Kit (ABTS, Chemical Method)
E-BC-K278-S	Glutathione-S-Transferase (GST) Activity Assay Kit
E-BC-K279-M	Potassium (K) Turbidimetric Assay Kit
E-BC-K284-M	Plant Flavonoids Colorimetric Assay Kit
E-BC-K284-S	Plant Flavonoids Colorimetric Assay Kit
E-BC-K298-F	Thiobarbituric Acid Reactants (TBARS) Fluorometric Assay Kit
E-BC-K298-M	Thiobarbituric Acid Reactants (TBARS) Colorimetric Assay Kit
E-BC-K300-M	Copper (Cu ²⁺) Colorimetric Assay Kit (Complexing Method)
E-BC-K318-M	BCA Protein Colorimetric Assay Kit
E-BC-K329-S	Urea (BUN) Colorimetric Assay Kit (Diacetyl Oxime Method)
E-BC-K340-M	Glycerol Colorimetric Assay Kit

E-BC-K351-M	Citric Acid (CA) Colorimetric Assay Kit
E-BC-K351-S	Citric Acid (CA) Colorimetric Assay Kit
E-BC-K352-M	Cysteine (Cys) Colorimetric Assay Kit
E-BC-K353-M	Ascorbate Peroxidase (APX) Activity Assay Kit
E-BC-K353-S	Ascorbate Peroxidase (APX) Activity Assay Kit
E-BC-K354-M	Total Phenols Colorimetric Assay Kit (Plant samples)
E-BC-K354-S	Total Phenols Colorimetric Assay Kit (Plant Samples)
E-BC-K355-M	Hydrogen Sulfide (H ₂ S) Colorimetric Assay Kit (Indirect Method)
E-BC-K500-M	Pyruvate Decarboxylase (PDC) Activity Assay Kit
E-BC-K508-M	Acetyl-CoA Carboxylase (ACC) Activity Assay Kit
E-BC-K520-M	Glucose Oxidase (GOD) Activity Assay Kit
E-BC-K522-S	Phenylalanine Ammonia Lyase (PAL) Activity Assay Kit
E-BC-K524-M	Diamine Oxidase (DAO) Activity Assay Kit
E-BC-K527-M	Hydroxyl Free Radical Scavenging Capacity Assay Kit
E-BC-K539-M	Na ⁺ K ⁺ -ATPase Activity Assay Kit
E-BC-K548-M	Thioredoxin Reductase (TrxR) Activity Assay Kit
E-BC-K554-M	NADP-Isocitrate Dehydrogenase (NADP-IDH) Activity Assay Kit
E-BC-K555-M	NADP-Malic Enzyme (NADP-ME) Activity Colorimetric Assay Kit
E-BC-K556-M	NAD-Malic Enzyme (NAD-ME) Activity Colorimetric Assay Kit
E-BC-K558-M	Creatine Kinase (CK) Activity Assay Kit
E-BC-K558-S	Creatine kinase (CK) Activity Assay Kit
E-BC-K561-M	NAD-Malate Dehydrogenase (NAD-MDH) Activity Assay Kit
E-BC-K562-M	NADP-Malate Dehydrogenase (NADP-MDH) Activity Assay Kit
E-BC-K565-M	Aldehyde Dehydrogenase (ALDH) Activity Assay Kit
E-BC-K568-M	Leucine Aminopeptidase (LAP) Activity Assay Kit
E-BC-K579-M	Carboxylesterase (CES) Activity Colorimetric Assay Kit
E-BC-K605-M	Phosphoenolpyruvate Carboxykinase (PEPCK) Activity Colorimetric Assay Kit
E-BC-K608-M	Pyruvate Carboxylase (PC) Activity Assay Kit
E-BC-K610-M	Hexokinase (HK) Activity Assay Kit
E-BC-K611-M	Pyruvate Kinase (PK) Activity Assay Kit
E-BC-K612-M	Phosphofructokinase (PFK) Activity Assay Kit
E-BC-K613-M	Phosphoenolpyruvate Carboxylase (PEPC) Activity Assay Kit
E-BC-K630-M	α-Galactosidase (α-Gal) Activity Assay Kit
E-BC-K631-M	β-Galactosidase (β-GAL) Activity Assay Kit
E-BC-K649-M	Succinate Dehydrogenase (SDH) Activity Assay Kit
E-BC-K650-M	Pyruvate Dehydrogenase (PDH) Activity Assay Kit
E-BC-K651-M	NAD-Isocitrate Dehydrogenase (NAD-IDH) Activity Assay Kit
E-BC-K656-M	Aconitase (ACO) Activity Assay Kit
E-BC-K659-M	Glutamate Synthase (GOGAT) Activity Assay Kit

E-BC-K660-M	Glutaminase (GLS) Activity Assay Kit
E-BC-K664-M	Glutamine Synthetase (GS) Activity Colorimetric Assay Kit
E-BC-K692-S	Glycolate Oxidase Activity Assay Kit
E-BC-K751-M	Sucrase Activity Assay Kit
E-BC-K753-M	Lactose Colorimetric Assay Kit
E-BC-K754-M	Trehalose Colorimetric Assay Kit (Enzyme Method)
E-BC-K755-M	Total Carbohydrate Colorimetric Assay Kit
E-BC-K756-M	Maltose Colorimetric Assay Kit
E-BC-K757-M	Revertose Colorimetric Assay Kit (DNS Method)
E-BC-K759-M	Glutamate Dehydrogenase (GDH) Activity Assay Kit
E-BC-K760-M	Total Bilirubin (TBIL) Colorimetric Assay Kit
E-BC-K761-M	Direct Bilirubin (DBIL) Colorimetric Assay Kit
E-BC-K762-M	α -L-Fucosidase (AFU) Activity Assay Kit
E-BC-K763-M	Alcohol Dehydrogenase (ADH) Activity Assay Kit
E-BC-K766-M	Lactate Dehydrogenase (LDH) Activity Assay Kit (WST-8 Method)
E-BC-K768-M	Glucokinase (GCK) Activity Colorimetric Assay Kit
E-BC-K769-M	Acetaldehyde Colorimetric Assay Kit
E-BC-K770-M	Phosphoglucose Isomerase (PGI) Activity Colorimetric Assay Kit
E-BC-K771-M	Lactate Dehydrogenase (LDH) Cytotoxicity Colorimetric Assay Kit
E-BC-K772-M	Total Iron Colorimetric Assay Kit
E-BC-K773-M	Ferrous Iron Colorimetric Assay Kit
E-BC-K774-M	ATP Colorimetric Assay Kit (Enzyme Method)
E-BC-K775-M	Cell Copper (Cu ²⁺) Colorimetric Assay Kit (Complexing Method)
E-BC-K776-M	Selenium Colorimetric Assay Kit
E-BC-K777-M	Ammonia Colorimetric Assay Kit
E-BC-K778-M	Phosphoglycerate Kinase (PGK) Activity Colorimetric Assay Kit
E-BC-K779-M	Acyl-CoA Synthetase (ACS) Activity Colorimetric Assay Kit
E-BC-K781-M	Phospholipase C (PLC) Activity Colorimetric Assay Kit
E-BC-K782-M	Total Phosphodiesterase (PDEs) Activity Colorimetric Assay Kit
E-BC-K783-M	Phospholipase D (PLD) Activity Colorimetric Assay Kit
E-BC-K784-M	Fatty Acid Oxidation (FAO) Colorimetric Assay Kit
E-BC-K785-M	β -Hydroxybutyrate (Ketone Body) Colorimetric Assay Kit
E-BC-K786-M	Lipase (LPS) Activity Assay Kit
E-BC-K790-M	Glyceraldehyde-3-Phosphate Dehydrogenase (GAPDH) Activity Assay Kit
E-BC-K791-M	Glycerol Kinase (GK) Activity Colorimetric Assay Kit
E-BC-K792-M	Free Fatty Acids (NEFA/FFA) Colorimetric Assay Kit (Enzyme Method)
E-BC-K793-M	Monoacylglycerol Lipase (MAGL) Activity Colorimetric Assay Kit
E-BC-K794-M	Squalene Synthase (SQS) Activity Colorimetric Assay Kit
E-BC-K795-M	3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase (HMGR) Activity Colorimetric Assay Kit

E-BC-K796-M	Phosphatidylcholine (PC) Colorimetric Assay Kit
E-BC-K798-M	Lipolysis Colorimetric Assay Kit
E-BC-K800-M	Glutathione-S-Transferase (GST) Activity Assay Kit (DTNB Method)
E-BC-K801-M	Total Antioxidant Status (TAS) Colorimetric Assay Kit
E-BC-K802-M	Total Oxidant Status (TOS) Colorimetric Assay Kit
E-BC-K803-M	NADP ⁺ /NADPH Colorimetric Assay Kit (WST-8)
E-BC-K804-M	NAD ⁺ /NADH Colorimetric Assay Kit (WST-8)
E-BC-K805-M	Xanthine Oxidase (XOD) Activity Assay Kit
E-BC-K806-M	NADH Oxidase (NOX) Activity Assay Kit
E-BC-K807-M	DPPH Free Radical Scavenging Capacity Colorimetric Assay Kit
E-BC-K808-M	Cytochrome P450 Reductase (CPR) Activity Assay Kit
E-BC-K809-M	Cell Glutathione Peroxidase (GPX) Activity Assay Kit
E-BC-K810-M	Aldo-Keto Reductase (ARKs) Activity Assay Kit
E-BC-K811-M	Cystathionine- β -synthase (CBS) Activity Colorimetric Assay Kit
E-BC-K812-M	Human Immunoglobulin G (IgG) Turbidimetric Assay Kit
E-BC-K813-M	Sialic Acid (SA) Colorimetric Assay Kit (Enzyme Method)
E-BC-K814-M	Enhanced Cell Malondialdehyde (MDA) Colorimetric Assay Kit
E-BC-K815-M	NADPH Oxidase (NAO) Activity Colorimetric Assay Kit
E-BC-K816-M	Dihydrofolate Reductase (DHFR) Activity Colorimetric Assay Kit
E-BC-K817-M	Nitric Oxide (NO) Colorimetric Assay Kit
E-BC-K820-M	Tannins Colorimetric Assay Kit
E-BC-K821-M	α -Glucosidase (α -GC) Activity Assay Kit
E-BC-K822-M	β -Glucosidase (β -GC) Activity Assay Kit
E-BC-K823-M	Protein Tyrosine Phosphatase-1B (PTP-1B) Activity Colorimetric Assay Kit
E-BC-K826-M	ATPase/GTPase Activity Colorimetric Assay kit
E-BC-K831-M	ATPase Activity Assay Kit
E-BC-K832-M	Phosphate Colorimetric Assay Kit(Malachite Green Method)
E-BC-K833-M	Adenosine Colorimetric Assay Kit
E-BC-K834-M	Cell Mitochondrial Complex I (NADH-CoQ Reductase) Activity Assay Kit
E-BC-K835-M	Cell Mitochondrial Complex II Activity Assay Kit
E-BC-K836-M	Cell Mitochondrial Complex III (Coenzyme Q-Cytochrome C Reductase) Activity Assay Kit
E-BC-K837-M	Cell Mitochondrial Complex IV (Cytochrome C Oxidase) Activity Assay Kit
E-BC-K838-M	Cell Mitochondrial Complex V (F ₀ F ₁ -ATPase/ATP Synthase) Activity Assay Kit
E-BC-K839-M	Alanine (Ala) Colorimetric Assay Kit
E-BC-K840-M	Transglutaminases (TGs) Activity Colorimetric Assay Kit
E-BC-K841-M	S-Adenosyl-Homocysteine Hydrolase (SAHH) Activity Colorimetric Assay Kit
E-BC-K843-M	Trypsin Activity Assay Kit
E-BC-K844-M	Pepsin Activity Assay Kit
E-BC-K845-M	Chymotrypsin Activity Colorimetric Assay Kit

E-BC-K846-M	Phenylalanine (Phe) Colorimetric Assay Kit
E-BC-K847-M	Branched Chain Amino Acids (BCAA) Colorimetric Assay Kit
E-BC-K848-M	Arginase Activity Colorimetric Assay Kit
E-BC-K849-M	Aspartate (Asp) Colorimetric Assay Kit
E-BC-K850-M	L-Arginine (L-Arg) Colorimetric Assay Kit
E-BC-K851-M	Tyrosine Ammonia-lyase (TAL) Activity Assay Kit
E-BC-K852-M	γ -Aminobutyric Acid (GABA) Colorimetric Assay Kit
E-BC-K853-M	Glutamine (Gln) Colorimetric Assay Kit
E-BC-K854-M	Homocysteine (Hcy) Colorimetric Assay Kit
E-BC-K855-S	Plant Root Activity Assay Kit
E-BC-K856-M	Adenosine Diphosphate Glucose Prophosphorylase(AGP) Activity Assay Kit
E-BC-K857-M	Uridine Diphosphate Glucose Pyrophosphorylase (UGP) Activity Assay Kit
E-BC-K858-M	Plant Nitrate Nitrogen Colorimetric Assay Kit
E-BC-K859-M	Plant Ammonium Nitrogen Colorimetric Assay Kit
E-BC-K860-M	Proline Dehydrogenase (ProDH) Activity Assay Kit
E-BC-K861-M	Lipoxygenase (LOX) Activity Assay Kit
E-BC-K863-M	Nitrite Reductase(NiR) Activity Colorimetric Assay Kit
E-BC-K864-M	Dehydroascorbic Acid Colorimetric Assay Kit
E-BC-K865-S	Lactase Activity Colorimetric Assay Kit
E-BC-K866-M	Plant Soluble Sugar Colorimetric Assay Kit
E-BC-K867-M	Pectinase Activity Colorimetric Assay Kit
E-BC-K868-M	Chitinase Activity Colorimetric Assay Kit
E-BC-K871-M	Tartrate Resistant Acid Phosphatase (TRAP) Activity Assay Kit
E-BC-K872-M	Fluoride Resistant Acid Phosphatase (FRAP) Activity Assay Kit
E-BC-K880-M	Cell Total Iron Colorimetric Assay Kit
E-BC-K881-M	Cell Ferrous Iron Colorimetric Assay Kit
E-BC-K883-M	Glutathione Peroxidase 4 (GPX4) Activity Assay Kit
E-BC-K884-M	NETosis Colorimetric Assay Kit
E-BC-K885-M	Coenzyme Q10 (CoQ10) Colorimetric Assay Kit
E-BC-K886-M	Protease Activity Colorimetric Assay Kit
E-BC-K891-M	Ethanol Colorimetric Assay Kit
E-BC-K892-S	Oxalate Colorimetric Assay Kit
E-BC-K893-M	Histamine Colorimetric Assay Kit
E-BC-K894-M	Hydrogen sulfide (H ₂ S) Colorimetric Assay Kit (Direct Method)
E-BC-K895-M	Nitrate Colorimetric Assay Kit
E-BC-K896-M	Uricase Activity Colorimetric Assay Kit
E-BC-K897-M	Uric Acid (UA) Colorimetric Assay Kit (Enzyme Method)
E-BC-K899-M	Neutrophil Elastase (NE) Activity Colorimetric Assay Kit
E-BC-K900-M	Creatine Colorimetric Assay Kit(Enzyme Method)

E-BC-K901-M	Oxaloacetate Colorimetric Assay Kit
E-BC-K902-M	Succinic Acid Colorimetric Assay Kit
E-BC-K903-M	Glutamic Acid (Glu) Colorimetric Assay Kit
E-BC-K905-M	Malic Acid Colorimetric Assay Kit
E-BC-K906-M	Succinyl-CoA Synthetase (SCS) Activity Assay Kit
E-BC-K907-M	Coenzyme A (CoA) Colorimetric Assay Kit
E-BC-K908-M	Fumarase (FUM) Activity Colorimetric Assay Kit
E-BC-K909-M	Asparaginase (ASNase) Activity Colorimetric Assay Kit
E-BC-K910-M	Acetate Kinase (Ack) Activity Colorimetric Assay Kit
E-BC-K912-M	α -Ketoglutarate Dehydrogenase (α -KGDH) Activity Colorimetric Assay Kit
E-BC-K940-M	Hexavalent Chromium Ion (Cr^{6+}) Colorimetric Assay Kit
E-BC-K945-M	Total Iron Binding Capacity (TIBC) And Serum Iron Assay Kit (Direct Method)
E-BC-L001	ALTD lysis buffer
E-BC-L002	CAMT lysis buffer