



# Cardiovascular Antibodies

**Elabscience<sup>®</sup>**

Focus On Lab Research Service For Life Science

# About Elabscience

---

Elabscience specializes in immunodiagnostic technology for the life science community. We have complete R&D and manufacture platform and have good quality control for every kit. We endeavor to reduce lot to lot variation to keep your experiment results more consistent and precise. Through unremitting efforts and development, our customers have spread in more than 100 countries and districts all over the world.

Elabscience's major products are Proteins, Antibodies, ELISA kits, CLIA Kits, labeling kits and related reagents. Elabscience also provides custom services for our customers including protein expression and purification service, antibody customized service, ELISA customized service, cell biology service, molecular biology service, and whole research project solutions.



# Product introduction

---






# The cardiovascular system

Cardiovascular system is an organ system that permits blood to circulate and transport nutrients (such as amino acids and electrolytes), oxygen, carbon dioxide, hormones, and blood cells to and from the cells in the body to provide nourishment and help in fighting diseases, stabilize temperature and pH, and maintain homeostasis. The study of the blood flow is called hemodynamics. The study of the properties of the blood flow is called hemorheology. Much research aims to understand the cardiovascular system during disease development with angiogenesis, vasculature, atherosclerosis and heart research.

## Atherosclerosis








# Product details

## Angiogenesis



Cat. No.	Product Name	Swissprot	Reactivity	Application
EPM12352	LAMR1 RPSA Monoclonal Antibody	P08865	H M	ELISA WB IHC
ESAP14485	RPSA Polyclonal Antibody	P08865	H M R	ELISA WB IHC
ENT0112	ADAMTS-1 Polyclonal Antibody	Q9UHI8	H	WB ELISA
ESAP16247	ADAMTS2 Polyclonal Antibody	O95450	H	ELISA IHC
EAP0618	CTSK Polyclonal Antibody	P43235	H M R	WB
ESAP11295	CYR61 Polyclonal Antibody	O00622	H M R	ELISA WB IHC
ESAP11500	ITGA2B Polyclonal Antibody	P08514	H M	ELISA WB IHC
ESAP16712	IFT74 Polyclonal Antibody	Q96LB3	H	ELISA IHC
ENT1025	COL4A2 Polyclonal Antibody	P08572	H M	WB IHC IF ELISA
EAP0851	PTK2 Polyclonal Antibody	Q05397	H M R	WB IHC
ENP0108	FAK (Phospho-Y861) Polyclonal Antibody	Q05397	H M R Mk	WB ELISA



# Angiogenesis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENP0629	FAK (Phospho-Y576) Polyclonal Antibody	Q05397	H M R	WB IHC ELISA
ENP0649	FAK (Phospho-B17Y407) Polyclonal Antibody	Q05397	H M R Mk	WB IHC ELISA
ENP0739	FAK (Phospho-Y397) Polyclonal Antibody	Q05397	H M R	WB IHC ELISA
ENT1660	FAK Polyclonal Antibody	Q05397	H M R	WB IHC IF ELISA
ESAP12224	PTK2 (phospho-Tyr576/Tyr577) Polyclonal Antibody	Q05397	H	WB IF
ENM0137	Fibronectin Monoclonal Antibody	P02751	H M R	WB
ENT1733	FN1 Polyclonal Antibody	P02751	H M R	WB IHC ELISA
ESAP13758	HSPG2 Polyclonal Antibody	P98160	H M	ELISA WB IHC
ESAP13535	HAS1 Polyclonal Antibody	Q92839	H M	ELISA WB IHC
ESAP11528	LAMB3 Polyclonal Antibody	Q13751	H M	ELISA WB IHC
EAP0828	ITGAV Polyclonal Antibody	P06756	H M R	WB
ENC092	Cleaved-Integrin $\alpha$ V HC (K889) Polyclonal Antibody	P06756	H	WB ELISA
ENT2365	Integrin $\alpha$ V Polyclonal Antibody	P06756	H M	WB IHC ELISA
ENP0524	Integrin $\beta$ 3 (Phospho-Y785) Polyclonal Antibody	P05106	H M R	WB ELISA
ENT2370	Integrin $\beta$ 3 Polyclonal Antibody	P05106	H M R	WB IHC ELISA
ENC091	Cleaved-Integrin $\alpha$ 7 LC (E959) Polyclonal Antibody	Q13683	H Mk	WB ELISA
ESAP16050	ITGA7 Polyclonal Antibody	Q13683	H M R	ELISA IHC
ENT2531	Laminin $\gamma$ -1 Polyclonal Antibody	P11047	H M R Mk	WB IHC IF ELISA
ENT2528	Laminin $\beta$ -1 Polyclonal Antibody	P07942	H M	WB IHC IF ELISA
ENT2527	Laminin $\alpha$ -5 Polyclonal Antibody	O15230	H M	WB IHC ELISA
ENK5379	Laminin $\gamma$ -2 Polyclonal Antibody	Q13753	H R	WB ELISA
ESAP12751	MADCAM1 Polyclonal Antibody	Q13477	H M	ELISA WB
ENC099	Cleaved-MMP-1 22k (F100) Polyclonal Antibody	P03956	H	WB ELISA
ESAP14082	MMP1 Polyclonal Antibody	P03956	H M	ELISA WB IHC
ENT2797	MMP-19 Polyclonal Antibody	Q99542	H M	WB IHC ELISA



# Angiogenesis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENC063	Cleaved-MMP-17 (Q129) Polyclonal Antibody	Q9ULZ9	H	WB ELISA
ESAP12775	MMP17 Polyclonal Antibody	Q9ULZ9	H	ELISA WB
ENC059	Cleaved-MMP-10 (F99) Polyclonal Antibody	P09238	H	WB ELISA
ESAP16113	MMP10 Polyclonal Antibody	P09238	H M	ELISA WB IHC
ESAP16114	MMP11 Polyclonal Antibody	P24347	H M	ELISA WB IHC
ENC060	Cleaved-MMP-12 (G106) Polyclonal Antibody	P39900	H	WB ELISA
ESAP14487	MMP12 Polyclonal Antibody	P39900	H	ELISA IHC
ESAP10179	MMP13 Polyclonal Antibody	P45452	H M R	ELISA WB IHC
EAP1022	MMP14 Polyclonal Antibody	P50281	H M R	WB IHC
ENC061	Cleaved-MMP-14 (Y112) Polyclonal Antibody	P50281	H M R	WB ELISA
ESAP14489	MMP15 Polyclonal Antibody	P51511	H M	ELISA WB
ENT2912	MMP-16 Polyclonal Antibody	P51512	H M R	WB IHC ELISA
ENT2798	MMP-2 Polyclonal Antibody	P08253	H M R Mk	WB IHC IF ELISA
ENC066	Cleaved-MMP-3 (F100) Polyclonal Antibody	P08254	H	WB ELISA
ENT4465	MMP-3 Polyclonal Antibody	P08254	H M R	WB IHC IF ELISA
ENT2663	MMP-7 Polyclonal Antibody	P09237	H M R Mk	WB IHC ELISA
ESAP16110	MMP8 Polyclonal Antibody	P22894	H	ELISA WB IHC
ENT1892	MMP-9 Polyclonal Antibody	P14780	H	WB IHC IF ELISA
ENT3467	OPN Polyclonal Antibody	P10451	H M R	WB IHC ELISA
EAP1448	PKP2 Polyclonal Antibody	Q99959	H M R	WB IHC IF
ESAP15664	ROBO4 Polyclonal Antibody	Q8WZ75	H	ELISA IHC
ESAP14655	SPARC Polyclonal Antibody	P09486	H M R	ELISA WB IHC
ENT4658	TIMP-1 Polyclonal Antibody	P01033	H M R	WB ELISA
ENT4659	TIMP-2 Polyclonal Antibody	P16035	H M R	WB IHC IF ELISA
EAP0433	TIMP3 Polyclonal Antibody	P35625	H M R	WB IHC IF





# Angiogenesis

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP2700	TIMP4 Polyclonal Antibody	Q99727	H	WB IHC IF
ESAP12404	ADM Polyclonal Antibody	P35318	H M	ELISA WB IHC
ESAP14056	ANG Polyclonal Antibody	P03950	H	ELISA WB
EAP0764	ANGPTL4 Polyclonal Antibody	Q9BY76	H M	WB IHC IF
ESAP13396	CHRD2 Polyclonal Antibody	Q6WN34	H	ELISA WB IHC
ESAP16723	CRIM1 Polyclonal Antibody	Q9NZV1	H M R	ELISA IHC
ESAP16552	RCAN2 Polyclonal Antibody	Q14206	H M R	ELISA IHC
ESAP11344	EDIL3 Polyclonal Antibody	O43854	H M	ELISA WB IHC
ESAP16823	GLMN Polyclonal Antibody	Q92990	H	ELISA IHC
ESAP12492	CCL15 Polyclonal Antibody	Q16663	H	ELISA IHC
ENT2821	MOX-2 Polyclonal Antibody	P50222	H M R	WB IHC ELISA
EPM12392	Nestin Monoclonal Antibody	P48681	H M R;	ELISA WB IHC IF
ESAP14526	NPY Polyclonal Antibody	P01303	H M R	ELISA WB
ENT0335	ARP-1 Polyclonal Antibody	P24468	H M R	WB IHC IF ELISA
EAP3588	PROK2 Polyclonal Antibody	Q9HC23	H	IHC
EAP0582	PGF Polyclonal Antibody	P49763	H M R	WB
ESAP16404	S100A4 Polyclonal Antibody	P26447	H	ELISA WB IHC
ENT5299	PDGF-D Polyclonal Antibody	Q9GZP0	H	WB ELISA
ESAP17226	SMOC2 Polyclonal Antibody	Q9H3U7	H M	ELISA WB IHC
EAP1802	ANTXR1 Polyclonal Antibody	Q9H6X2	H M R	WB IHC IF
EGM0307	VEGFA Polyclonal Antibody	P15692	R	WB IHC
ENT1731	Flt3-L Polyclonal Antibody	P49771	H	WB IHC ELISA
EAP0472	IFNB1 Polyclonal Antibody	P01574	H	WB IHC IF
ENK5470	IL-1RII Polyclonal Antibody	P27930	H	WB ELISA
EPM12328	IL1R2 Monoclonal Antibody	P27930	H M;	ELISA WB IF



# Angiogenesis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT2343	IL-8 Polyclonal Antibody	P10145	H	WB IHC ELISA
EAP2753	IL17F Polyclonal Antibody	Q96PD4	H	WB IHC IF
ENK5487	IL-33 Polyclonal Antibody	O95760	H	WB ELISA
EAP1176	CSF1R Polyclonal Antibody	P07333	H M R	WB
ENP0734	c-Fms (Phospho-Y561) Polyclonal Antibody	P07333	H M R	WB IHC ELISA
ENT0881	c-Fms Polyclonal Antibody	P07333	H	WB IHC IF ELISA
ESAP14742	IL1RL1 Polyclonal Antibody	Q01638	H	ELISA WB IHC
EAP2023	ELF1 Polyclonal Antibody	P32519	H M	WB IHC IF
ENT1589	Ephrin-A1 Polyclonal Antibody	P20827	H M R	WB IHC IF ELISA
ENC052	Cleaved-Ephrin-A2 (N188) Polyclonal Antibody	O43921	H M	WB ELISA
ENT1590	Ephrin-A2 Polyclonal Antibody	O43921	H M	WB ELISA
ENT1591	Ephrin-A3 Polyclonal Antibody	P52797	H M R	WB IHC IF ELISA
EAP2322	EPHA1 Polyclonal Antibody	P21709	H	WB IHC IF
EAP2227	EPHA2 Polyclonal Antibody	P29317	H M R	WB
ENP0098	EphA2 (Phospho-Y588) Polyclonal Antibody	P29317	H M	WB ELISA
ENT1575	EphA2 Polyclonal Antibody	P29317	H M	WB ELISA
ENP0550	EphA2/3/4 (Phospho-Y588/596) Polyclonal Antibody	P29317/P29320/ P54764	H M R	WB IF ELISA
ENT1576	EphA2/3/4 Polyclonal Antibody	P29317/P29320/ P54764	H R	WB IF ELISA
ENP0506	EphA2/5 (Phospho-Y594) Polyclonal Antibody	P29317/P54756	H M	WB ELISA
ENT1578	EphA3 Polyclonal Antibody	P29320	H M R	WB ELISA
EAP3827	EPHA4 Polyclonal Antibody	P54764	M R	WB
ENP0507	EphA4 (Phospho-Y596) Polyclonal Antibody	P54764	H M	WB ELISA
ENT1579	EphA5 Polyclonal Antibody	P54756	H M R	WB ELISA
ENP0508	EphA7 (Phospho-Y791) Polyclonal Antibody	Q15375	H M R	WB ELISA
ESAP13265	EPHA7 Polyclonal Antibody	Q15375	H M R	ELISA IHC



# Angiogenesis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT1582	EphB1 Polyclonal Antibody	P54762	H M R Mk	WB ELISA
ENP0551	EphB1/2 (Phospho-Y594/604) Polyclonal Antibody	P54762/P29323	H M R	WB IF ELISA
ENT1583	EphB1/2 Polyclonal Antibody	P54762/P29323	H M R	WB IHC IF ELISA
ENT1586	EphB2 Polyclonal Antibody	P29323	H	WB IHC ELISA
ENP0356	Ephrin-B1 (Phospho-Y317) Polyclonal Antibody	P98172	H M R	WB ELISA
ENT1594	Ephrin-B1 Polyclonal Antibody	P98172	H M R Mk	WB ELISA
ENP0284	Ephrin-B1/2 (Phospho-Y330) Polyclonal Antibody	P98172/P52799	H M R	WB ELISA
ENT1595	Ephrin-B1/2 Polyclonal Antibody	P98172/P52799	H M R	WB IHC ELISA
ENP0812	Ephrin-B1/2/3 (Phospho-Y324) Polyclonal Antibody	P98172/P52799/ Q15768	H M R	WB IHC ELISA
ESAP14326	EFNB2 Polyclonal Antibody	P52799	H M	ELISA IHC
ENT1596	Ephrin-B3 Polyclonal Antibody	Q15768	H M R	WB IHC ELISA
ESAP13610	KIDINS220 Polyclonal Antibody	Q9ULH0	H R	ELISA WB
ENK5409	Inhibin $\beta$ -A Polyclonal Antibody	P08476	H M R	WB ELISA
EPM12339	INHBA-Specific Monoclonal Antibody	P08476	H Pig;	ELISA WB IHC
ESAP10533	INHBA Polyclonal Antibody	P08476	H M R	ELISA WB IHC
ENK5408	Inhibin $\beta$ -B Polyclonal Antibody	P09529	H M R	WB ELISA
ENT0198	ALK-1 Polyclonal Antibody	P37023	H M R	WB IHC IF ELISA
ESAP14153	ANGPT1 Polyclonal Antibody	Q15389	H M R	ELISA IHC
ESAP10253	ANGPT2 Polyclonal Antibody	O15123	H M R	ELISA WB IHC
ESAP16267	ANGPTL2 Polyclonal Antibody	Q9UKU9	H M	ELISA WB IHC
EAP0128	ANGPTL3 Polyclonal Antibody	Q9Y5C1	H M R	WB
ENT0498	BMP-2 Polyclonal Antibody	P12643	H M R	WB IHC ELISA
ESAP14207	BMP4 Polyclonal Antibody	P12644	H M R	ELISA WB IHC
ENT0501	BMP-5 Polyclonal Antibody	P22003	H M	WB IHC ELISA
ESAP15098	BMP15 Polyclonal Antibody	O95972	H	ELISA WB IHC





# Angiogenesis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP15425	NOV Polyclonal Antibody	P48745	H M R	ELISA WB IHC
ESAP14252	CDK7 Polyclonal Antibody	P50613	H M R	ELISA WB IHC
ESAP12525	CTGF Polyclonal Antibody	P29279	H M R	ELISA WB IHC
ESAP13399	TSPYL2 Polyclonal Antibody	Q9H2G4	H	ELISA IHC
ESAP16705	TSPYL2 Polyclonal Antibody	Q9H2G4	H	ELISA WB
ESAP12476	FGF1 Polyclonal Antibody	P05230	H M R	ELISA IHC
ESAP15833	FGF2 Polyclonal Antibody	P09038	H M R	ELISA WB IHC
ESAP12477	FGF4 Polyclonal Antibody	P08620	H M	ELISA IHC
ENT5190	FGF-5 Polyclonal Antibody	P12034	H	WB ELISA
ENT5169	FGF-6 Polyclonal Antibody	P10767	H M R	WB ELISA
ESAP15836	FGF8 Polyclonal Antibody	P55075	H M	ELISA WB
ESAP10302	FGF9 Polyclonal Antibody	P31371	H M R	ELISA WB IHC
ENT5547	FGF-11 Polyclonal Antibody	Q92914	H M R	WB ELISA
ESAP10303	FGF12 Polyclonal Antibody	P61328	H M R	ELISA WB
ESAP16795	FGF16 Polyclonal Antibody	O43320	H M R	ELISA IHC
ESAP15328	FGF17 Polyclonal Antibody	O60258	H M R	ELISA WB
ESAP15330	FGF21 Polyclonal Antibody	Q9NSA1	H M	ELISA IHC
EAP0814	FGFR1 Polyclonal Antibody	P11362	H M R	WB IHC
ENP0518	Flg (Phospho-Y154) Polyclonal Antibody	P11362	H M R	WB IF ELISA
ENP0528	Flg (Phospho-Y654) Polyclonal Antibody	P11362	H M R	WB IF ELISA
ENT1717	Flg Polyclonal Antibody	P11362	H M R	WB IHC ELISA
EPM12279	FGFR1 Monoclonal Antibody	P11362	H M;	ELISA WB IHC
ENP0452	Flg/Bek (Phospho-Y463/466) Polyclonal Antibody	P11362/P21802	H M R	WB ELISA
ENT1719	Flg/Bek Polyclonal Antibody	P11362/P21802	H M R	WB ELISA
EAP0815	FGFR2 Polyclonal Antibody	P21802	H M R	WB IHC IF



# Angiogenesis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT1701	FGFR-3 Polyclonal Antibody	P22607	H M R	WB IHC ELISA
ENP0509	FGFR-4 (Phospho-Y642) Polyclonal Antibody	P22455	H M R	WB ELISA
ESAP12656	GDF1 Polyclonal Antibody	P27539	H M	ELISA WB
EAP3333	GDF5 Polyclonal Antibody	P43026	H M R	WB IHC
ESAP14367	GDF9 Polyclonal Antibody	O60383	H M	ELISA WB IHC
EAP3895	IGFBP4 Polyclonal Antibody	P22692	H M R	WB IHC
ENT5226	SHIP-2 Polyclonal Antibody	O15357	H M R	WB ELISA
ESAP13587	INPL1 Polyclonal Antibody	O15357	H M R	ELISA IHC
ESAP12478	FGF7 Polyclonal Antibody	P21781	H M R	ELISA IHC
ESAP15474	LIMS1 Polyclonal Antibody	P48059	H M	ELISA WB IHC
ENT3073	Neuropilin-1 Polyclonal Antibody	O14786	H M R	WB IHC ELISA
ESAP17061	PDGFC Polyclonal Antibody	Q9NRA1	H M R	ELISA IHC
EAP0635	SDC2 Polyclonal Antibody	P34741	H M R	WB
ENT4626	TGF $\alpha$ Polyclonal Antibody	P01135	H M R Mk	WB IHC ELISA
ESAP10795	TGFB1 Polyclonal Antibody	P01137	H M R	ELISA WB IHC
ENP1191	TGF $\beta$ RI (Phospho-S165) Polyclonal Antibody	P36897	H M R	IF ELISA
ENT4627	TGF $\beta$ RI Polyclonal Antibody	P36897	H M R	WB IF ELISA
ESAP13151	TGFBR1 Polyclonal Antibody	P36897	H M R	ELISA IHC
ENT4629	TGF $\beta$ RII Polyclonal Antibody	P37173	H M R	WB IHC ELISA
ESAP10410	TNFRSF12A Polyclonal Antibody	Q9NP84	H M	ELISA WB IHC
ENT4871	VEGF-B Polyclonal Antibody	P49765	H M R	WB IHC ELISA
EAP1028	VEGFC Polyclonal Antibody	P49767	H M R	WB IHC
ENT4872	VEGF-D Polyclonal Antibody	O43915	H M R Mk	WB IHC ELISA
ENP0678	Flt-1 (Phospho-Y1333) Polyclonal Antibody	P17948	H M R	WB IHC ELISA
ENP1063	Flt-1 (Phospho-Y1048) Polyclonal Antibody	P17948	H M R	IHC ELISA




# Angiogenesis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENP1199	FIt-1 (Phospho-Y1213) Polyclonal Antibody	P17948	H M R	WB ELISA
ENT1729	FIt-1 Polyclonal Antibody	P17948	H M R	IHC ELISA
EAP3846	FLT4 Polyclonal Antibody	P35916	H M	WB IHC IF
EPM12142	AMOT Monoclonal Antibody	Q4VCS5	H;	ELISA IHC WB
ESAP13332	BAI1 Polyclonal Antibody	O14514	H M	ELISA IHC
ESAP16637	BAI2 Polyclonal Antibody	O60241	H M	ELISA IHC
EAP4091	CST7 Polyclonal Antibody	O76096	H	WB
ESAP16767	AIMP1 Polyclonal Antibody	Q12904	H M	ELISA WB IHC
ENT2311	IL-12A Polyclonal Antibody	P29459	H M R	WB IHC ELISA
EPM12622	IL12B Monoclonal Antibody	P29460	H Pig;	ELISA WB
ESAP16895	JAG1 Polyclonal Antibody	P78504	H M R	ELISA IHC
EAP0990	EEF1D Polyclonal Antibody	P29692	H M R	WB IHC IF
EAP1882	LECT1 Polyclonal Antibody	O75829	H M	WB IHC IF
ESAP10677	SERPINF1 Polyclonal Antibody	P36955	H M	ELISA WB IHC
ESAP11708	RNH1 Polyclonal Antibody	P13489	H	ELISA WB IHC
ENT4235	SEMA4A Polyclonal Antibody	Q9H3S1	H M R Mk	WB IHC IF ELISA
ESAP17244	SPON1 Polyclonal Antibody	Q9HCB6	H M R	ELISA WB IHC
ESAP17296	TNFSF15 Polyclonal Antibody	O95150	H M R	ELISA WB IHC
ENT4853	Vasohibin Polyclonal Antibody	Q7L8A9	H M R	WB ELISA





# Atherosclerosis



Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP10159	ANXA5 Polyclonal Antibody	P08758	H M R	ELISA WB IHC
ESAP14112	BAG1 Polyclonal Antibody	Q99933	H	ELISA WB IHC
EAP2654	BTG1 Polyclonal Antibody	P62324	H M R	IHC IF
ENP0912	Calmodulin (Phospho-T80/S82) Polyclonal Antibody	P62158	H M R	WB IHC IF ELISA
ENT0612	Calmodulin Polyclonal Antibody	P62158	H M R	WB IHC IF ELISA
ENT0721	CD130 Polyclonal Antibody	P40189	H M	WB IHC IF ELISA
ESAP12717	IRF1 Polyclonal Antibody	P10914	H M R	ELISA WB IHC
ENT2286	IGFBP3 Polyclonal Antibody	P17936	H M	WB IHC ELISA
ESAP13937	SORT1 Polyclonal Antibody	Q99523	H M R	ELISA WB IHC
EAP3281	STK4 Polyclonal Antibody	Q13043	H M R	WB IHC IF
ESAP16665	STK3/STK4 Polyclonal Antibody	Q13043	H M	ELISA WB IHC
EAP3523	SIRT6 Polyclonal Antibody	Q8N6T7	H M	WB
EAP0161	SIRT7 Polyclonal Antibody	Q9NRC8	H M R	WB IHC IF



## Atherosclerosis

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP0875	TNFSF10 Polyclonal Antibody	P50591	H M	WB IHC
ENT4721	TRAIL Polyclonal Antibody	P50591	H	WB IHC IF ELISA
EAP1238	pan-AKT Polyclonal Antibody	P31749	H M R	WB
ESAP10871	ACOT1 Polyclonal Antibody	Q86TX2	H	ELISA IHC
ENC016	Adiponectin Polyclonal Antibody	Q15848	H	WB ELISA
ESAP16251	ADIPOR2 Polyclonal Antibody	Q86V24	H M	ELISA IHC
EAP0098	AKT2 Polyclonal Antibody	P31751	H M R	WB
ENP0009	Akt2 (Phospho-S474) Polyclonal Antibody	P31751	H M R	WB IHC ELISA
ENT0187	Akt2 Polyclonal Antibody	P31751	H M R	WB IHC IF ELISA
ESAP12208	AKT1/AKT2/AKT3 (phospho-Tyr315/316/312) Polyclonal Antibody	P31749	H M R	WB IHC IF
ENT0225	Amylin Polyclonal Antibody	P10997	H	WB IHC IF ELISA
ESAP17007	ADRB3 Polyclonal Antibody	P13945	H	ELISA IHC
ESAP10637	BDNF Polyclonal Antibody	P23560	H M R	ELISA WB IHC
ENT0638	Cardiotrophin-1 Polyclonal Antibody	Q16619	H M R	WB IHC ELISA
ENP0040	C/EBP $\beta$ (Phospho-T235) Polyclonal Antibody	P17676	H M R	WB IHC IF ELISA
ENT0553	C/EBP $\beta$ Polyclonal Antibody	P17676	H M R	WB IHC ELISA
EAP3512	MLXIPL Polyclonal Antibody	Q9NP71	H M R	WB
EAP3772	CKM Polyclonal Antibody	P06732	H M R	WB
ESAP15799	CKM Polyclonal Antibody	P06732	H M R	ELISA WB IHC
ESAP13457	ELOVL6 Polyclonal Antibody	Q9H5J4	H M	ELISA WB IHC
ESAP15349	FTO Polyclonal Antibody	Q9C0B1	H	ELISA WB IHC
ESAP11429	GCKR Polyclonal Antibody	Q14397	H M R	ELISA WB IHC
ENT1928	Glut1 Polyclonal Antibody	P11166	H M R	WB IHC ELISA
ESAP12665	SLC2A1 Polyclonal Antibody	P11166	H M R	ELISA WB IHC
ESAP12670	SLC2A11 Polyclonal Antibody	Q9BYW1	H	ELISA IHC



# Atherosclerosis

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP0444	UCP3 Polyclonal Antibody	P55916	H M R	WB IHC IF
EAP0595	USF1 Polyclonal Antibody	P22415	H M R	WB IHC IF
ESAP16474	S1PR3 Polyclonal Antibody	Q99500	H M	ELISA WB IHC
EAP4165	FBLN5 Polyclonal Antibody	Q9UBX5	H M	WB
ESAP12594	FBLN5 Polyclonal Antibody	Q9UBX5	H M R	ELISA IHC
ESAP13560	HSD11B2 Polyclonal Antibody	P80365	H	ELISA IHC
EPM12465	RISC Monoclonal Antibody	Q9HB40	H M;	ELISA WB IF
ESAP11711	SCPEP1 Polyclonal Antibody	Q9HB40	H M R	ELISA WB IHC
ENT4253	SERCA1 Polyclonal Antibody	O14983	H M R	WB IHC ELISA
ESAP13665	ATP2A1 Polyclonal Antibody	O14983	H M	ELISA WB IHC
ENT0132	Adenosine A3-R Polyclonal Antibody	P33765	H	WB IHC IF ELISA
ENT0360	AR- $\beta$ 2 Polyclonal Antibody	P07550	H	WB IHC IF ELISA
EAP0359	ALDH2 Polyclonal Antibody	P05091	H M R	WB IHC IF
EAP1133	EPOR Polyclonal Antibody	P19235	H M	WB
ENP0552	EpoR (Phospho-Y368) Polyclonal Antibody	P19235	H M R	WB IF ELISA
ENP1205	EpoR (Phospho-Y426) Polyclonal Antibody	P19235	H M R	WB ELISA
ESAP12582	EPOR Polyclonal Antibody	P19235	H M R	ELISA WB IHC
ESAP12676	UTS2R Polyclonal Antibody	Q9UKP6	H	ELISA WB
ESAP16863	HSPB6 Polyclonal Antibody	O14558	H M R	ELISA IHC
ENT2253	HSP27 Polyclonal Antibody	P04792	H	WB IHC ELISA
ESAP12123	HSPB1 (Phospho-Ser15) Polyclonal Antibody	P04792	H	WB IHC IF
ESAP12170	HSPB1 (Phospho-Ser78) Polyclonal Antibody	P04792	H	WB IHC IF
ESAP12171	HSPB1 (Phospho-Ser82) Polyclonal Antibody	P04792	H	WB IHC
ESAP11619	NT5C1A Polyclonal Antibody	Q9BX13	H M	ELISA IHC
EAP4036	SLC6A2 Polyclonal Antibody	P23975	H M R	WB



# Atherosclerosis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP14591	PRDX3 Polyclonal Antibody	P30048	H M R	ELISA WB IHC
ESAP14593	PRDX6 Polyclonal Antibody	P30041	H M R	ELISA WB IHC
ESAP13751	PDE5A Polyclonal Antibody	O76074	H	ELISA WB IHC
EAP0239	TNFRSF1B Polyclonal Antibody	P20333	H	WB IHC IF
ENT4688	TNF-R2 Polyclonal Antibody	P20333	H M R	WB IHC IF ELISA
ESAP10805	TNFRSF1B Polyclonal Antibody	P20333	H	ELISA WB IHC
EAP0238	TYMP Polyclonal Antibody	P19971	H M R	WB IHC IF
ENK5414	PD-ECGF Polyclonal Antibody	P19971	H	WB ELISA
ESAP13182	ABCA2 Polyclonal Antibody	Q9BZC7	H M R	ELISA IHC
ESAP12917	ABCG1 Polyclonal Antibody	P45844	H	ELISA WB IHC
ESAP12918	ABCG5 Polyclonal Antibody	Q9H222	H	ELISA IHC
ESAP15806	APOB Polyclonal Antibody	P04114	H M R	ELISA IHC
ENT0273	ApoE Polyclonal Antibody	P02649	H M R	WB IHC ELISA
ESAP14887	APOM Polyclonal Antibody	O95445	H M R	ELISA WB
ENT0276	ApoF Polyclonal Antibody	Q13790	H	WB ELISA
ESAP14885	APOA4 Polyclonal Antibody	P06727	H	ELISA IHC
ESAP13322	PLIN2 Polyclonal Antibody	Q99541	H M	ELISA WB IHC
EAP0395	CETP Polyclonal Antibody	P11597	H M	WB
ENP0366	HBP1 (Phospho-S402) Polyclonal Antibody	O60381	H	WB IHC ELISA
ENT2103	HBP1 Polyclonal Antibody	O60381	H M R Mk	WB IHC ELISA
EAP3711	ALOX5(Phospho-S271) Polyclonal Antibody	P09917	H	WB
ESAP10155	ALOX5 Polyclonal Antibody	P09917	H M R	ELISA IHC
ESAP13303	ALOX15 Polyclonal Antibody	P16050	H M R	ELISA WB IHC
EAP1995	ALOX15B Polyclonal Antibody	O15296	H M	WB IHC IF
EAP3692	ACAT2 Polyclonal Antibody	Q9BWD1	H M	WB



# Atherosclerosis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP14063	APOA1 Polyclonal Antibody	P02647	H	ELISA IHC
EPM12149	APOA1 Monoclonal Antibody	P02647	H;	ELISA WB IHC
ESAP14886	APOC1 Polyclonal Antibody	P02654	H	ELISA WB
ESAP12958	APOC4 Polyclonal Antibody	P55056	H	ELISA IHC
ESAP16277	APOL1 Polyclonal Antibody	O14791	H	ELISA WB IHC
EAP2170	LDLRAP1 Polyclonal Antibody	Q5SW96	H M	WB IHC IF
EAP3777	CNR1 Polyclonal Antibody	P21554	H M R	WB
ENT0687	CB1 Polyclonal Antibody	P21554	H M R Mk	WB IHC IF ELISA
ESAP12514	PTGS1 Polyclonal Antibody	P23219	H	ELISA IHC
ENT1067	Cox-1 Polyclonal Antibody	P23219	H	WB IHC ELISA
EAP2331	PTGS1 Polyclonal Antibody	P23219	H M R	WB IHC
ESAP16839	HDLBP Polyclonal Antibody	Q00341	H M R	ELISA WB IHC
EAP0520	HMGCR Polyclonal Antibody	P04035	H M R	WB IHC IF
ESAP16687	B3GAT1 Polyclonal Antibody	Q9P2W7	H M R	ELISA WB IHC
EAP4124	LIPG Polyclonal Antibody	Q9Y5X9	H M R	WB IHC
EAP0668	CES1 Polyclonal Antibody	P23141	H M R	WB IHC
ESAP15516	MSR1 Polyclonal Antibody	P21757	H M	ELISA IHC
ESAP17020	NPC1 Polyclonal Antibody	O15118	H	ELISA IHC
ESAP17065	PLIN1 Polyclonal Antibody	O60240	H M R	ELISA IHC
ENT5176	cPLA2 Polyclonal Antibody	P47712	H M R	WB ELISA
ESAP11264	PLA2G4A Polyclonal Antibody	P47712	H M R	ELISA WB IHC
EAP2624	SOAT1 Polyclonal Antibody	P35610	H	WB
ESAP14863	AUP1 Polyclonal Antibody	Q9Y679	H M R	ELISA WB IHC
EAP0097	CD40LG Polyclonal Antibody	P29965	H	WB IHC IF
ENT1251	Cystatin C Polyclonal Antibody	P01034	H M R Mk	WB ELISA





# Atherosclerosis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP12448	CST3 Polyclonal Antibody	P01034	H	ELISA WB IHC
ENK5359	Autotaxin Polyclonal Antibody	Q13822	H M R	WB ELISA
ENT1853	Gas6 Polyclonal Antibody	Q14393	H	WB IHC ELISA
ENT5045	NPAS2 Polyclonal Antibody	Q99743	H R	WB ELISA
ESAP13993	PTX3 Polyclonal Antibody	P26022	H M	ELISA WB IHC
EAP1023	TNFSF11 Polyclonal Antibody	O14788	H M R	WB IHC IF
ENT4001	RANKL Polyclonal Antibody	O14788	H M R	WB IF ELISA
EAP0751	TBXAS1 Polyclonal Antibody	P24557	H M R	WB IHC IF
ESAP12920	TNIP3 Polyclonal Antibody	Q96KP6	H M	ELISA WB
ENM0065	$\beta$ -Catenin Monoclonal Antibody	P35222	H M R	WB IHC
ENP0047	Catenin- $\beta$ (Phospho-S37) Polyclonal Antibody	P35222	H M R	WB IHC IF ELISA
ENP0048	Catenin- $\beta$ (Phospho-T41/S45) Polyclonal Antibody	P35222	H M R	WB ELISA
ENP0049	Catenin- $\beta$ (Phospho-Y654) Polyclonal Antibody	P35222	H M R	WB IHC ELISA
ENP1204	Catenin- $\beta$ (Phospho-S552) Polyclonal Antibody	P35222	H M R	WB ELISA
ENT0672	Catenin- $\beta$ Polyclonal Antibody	P35222	H M R	WB IHC IF ELISA
ESAP12150	CTNNB1 (Phospho-Ser33) Polyclonal Antibody	P35222	H M R	WB IHC
ESAP12151	CTNNB1 (Phospho-Ser37) Polyclonal Antibody	P35222	H M R	WB IHC
ENT0934	CKR-2 Polyclonal Antibody	P41597	M R	WB IF ELISA
ENT5635	CD192 Polyclonal Antibody	P41597	H	WB ELISA
ESAP15850	CCR2 Polyclonal Antibody	P41597	H M	ELISA WB IHC
ENT5290	CKR-4 Polyclonal Antibody	P51679	H M R	WB ELISA
ENP0289	CKR-5 (Phospho-S336) Polyclonal Antibody	P51681	H M R	WB ELISA
ENP0290	CKR-5 (Phospho-S349) Polyclonal Antibody	P51681	H	WB ELISA
ESAP16669	CCR5 Polyclonal Antibody	P51681	H M	ELISA IHC
ESAP16162	PECAM1 Polyclonal Antibody	P16284	H M	ELISA WB IHC



# Atherosclerosis

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP0022	CD40 Polyclonal Antibody	P25942	H M R	WB IHC
ENP0379	CD40 (Phospho-T254) Polyclonal Antibody	P25942	H M	WB ELISA
ESAP13381	CD59 Polyclonal Antibody	P13987	H	ELISA WB IHC
ESAP11108	NT5E Polyclonal Antibody	P21589	H M R	ELISA WB IHC
ESAP14598	SELPLG Polyclonal Antibody	Q14242	H	ELISA IHC
ESAP16320	CD163 Polyclonal Antibody	Q86VB7	H	ELISA IHC
EPM12566	CD200 Monoclonal Antibody	P41217	H;	ELISA WB IHC
ESAP15141	CD200R1 Polyclonal Antibody	Q8TD46	H	ELISA IHC
ESAP12503	CEBPD Polyclonal Antibody	P49716	H M R	ELISA WB IHC
ENT0556	C/EBP $\delta/\epsilon$ Polyclonal Antibody	P49716/Q15744	H M R	WB IHC IF ELISA
EAP1006	GJA4 Polyclonal Antibody	P35212	H M R	WB
EAP0900	GJA1 Polyclonal Antibody	P17302	H M R	WB
ENP0071	Connexin 43 (Phospho-S368) Polyclonal Antibody	P17302	H M R	WB IHC ELISA
ENT1046	Connexin 43 Polyclonal Antibody	P17302	H M R	WB IHC IF ELISA
ENK5354	Fractalkine Polyclonal Antibody	P78423	H	WB ELISA
ESAP12535	CX3CL1 Polyclonal Antibody	P78423	H M R	ELISA IHC
EAP1138	CXCR3 Polyclonal Antibody	P49682	H M R	WB IHC
ENT1461	EDG-2 Polyclonal Antibody	Q92633	H R	WB IHC IF ELISA
ENT1471	EDG-7 Polyclonal Antibody	Q9UBY5	H	WB IHC IF ELISA
ESAP16477	LPAR3 Polyclonal Antibody	Q9UBY5	H M R	ELISA IHC
ESAP15946	FPR2 Polyclonal Antibody	P25090	H	ELISA IHC
ESAP10500	HMOX1 Polyclonal Antibody	P09601	H M R	ELISA IHC
ESAP16685	ICAM1 Polyclonal Antibody	P05362	H	ELISA WB IHC
ENP0638	IKK $\gamma$ (Phospho-S31) Polyclonal Antibody	Q9Y6K9	H	WB IHC ELISA
ENT2305	IKK $\gamma$ Polyclonal Antibody	Q9Y6K9	H Mk	WB IHC IF ELISA



# Atherosclerosis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT2303	IKK $\beta$ Polyclonal Antibody	O14920	H M R	WB IHC ELISA
ENP0141	IKK $\alpha/\beta$ (Phospho-S176/177) Polyclonal Antibody	O15111/O14920	H M R	WB IHC ELISA
ENP0142	IKK $\alpha/\beta$ (Phospho-S180/181) Polyclonal Antibody	O15111/O14920	H M R	WB IHC ELISA
ENT2301	IKK $\alpha/\beta$ Polyclonal Antibody	O15111/O14920	H M R	WB IHC ELISA
ENP0140	IKK $\alpha$ (Phospho-T23) Polyclonal Antibody	O15111	H M R	WB IHC ELISA
ENT2300	IKK $\alpha$ Polyclonal Antibody	O15111	H M R	WB IHC ELISA
ENT2322	IL-1 $\beta$ Polyclonal Antibody	P01584	H M R	WB IHC ELISA
ENT5263	IL-1RI Polyclonal Antibody	P14778	H	WB ELISA
ESAP14131	IL1RAP Polyclonal Antibody	Q9NPH3	H M R	ELISA IHC
ENP0143	Integrin $\beta$ 1 (Phospho-T788) Polyclonal Antibody	P05556	H M R	WB IHC ELISA
ESAP14427	ITGB1 Polyclonal Antibody	P05556	H M R	ELISA WB IHC
ENT2389	IRAK-1 Polyclonal Antibody	P51617	H	WB IHC ELISA
ESAP10543	IRAK4 Polyclonal Antibody	Q9NWZ3	H M	ELISA IHC
ESAP12437	LOX Polyclonal Antibody	P28300	H M R	ELISA WB IHC
ENT5553	M-CSF Polyclonal Antibody	P09603	H	WB ELISA
ESAP15848	CCL2 Polyclonal Antibody	P13500	H	ELISA IHC
ENT2928	MyD88 Polyclonal Antibody	Q99836	H M R	WB IHC ELISA
ENH006	Mono-Methyl-NF $\kappa$ B-p65 (K314/K315) Polyclonal Antibody	Q04206	H M R	WB ELISA
ENK018	NF $\kappa$ B-p65 (Acetyl-K310) Polyclonal Antibody	Q04206	H M	WB IHC IF ELISA
ENK034	NF $\kappa$ B-p65 (Acetyl-K218) Polyclonal Antibody	Q04206	H M R	WB ELISA
ENK064	NF $\kappa$ B-p65 (Acetyl-K122) Polyclonal Antibody	Q04206	H M R	WB ELISA
ENK065	NF $\kappa$ B-p65 (Acetyl-K221) Polyclonal Antibody	Q04206	H M R	WB ELISA
ENK5428	NF $\kappa$ B-p65 Polyclonal Antibody	Q04206	H M R	WB ELISA
ENM0053	NF $\kappa$ B p65 Monoclonal Antibody	Q04206	H M R	WB IHC IP
ENP0187	NF $\kappa$ B-p65 (Phospho-S276) Polyclonal Antibody	Q04206	H M R	WB IHC IP ELISA



# Atherosclerosis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENP0188	NFκB-p65 (Phospho-S468) Polyclonal Antibody	Q04206	H M R	WB IHC IP ELISA
ENP0190	NFκB-p65 (Phospho-S529) Polyclonal Antibody	Q04206	H M	WB IHC ELISA
ENP0191	NFκB-p65 (Phospho-S536) Polyclonal Antibody	Q04206	H M R Mk	WB IHC IP ELISA
ENP0192	NFκB-p65 (Phospho-T254) Polyclonal Antibody	Q04206	H M R	WB IHC IP ELISA
ENP0847	NFκB-p65 (Phospho-S311) Polyclonal Antibody	Q04206	H M	WB IHC IF ELISA
ENT3102	NFκB-p65 Polyclonal Antibody	Q04206	H M R	WB ELISA
ESAP12056	RELA (Phospho-Thr254) Polyclonal Antibody	Q04206	H M R	WB IHC IF
ENK050	NFκB-p65 (Acetyl-K314/K315) Polyclonal Antibody	Q04206	H M R	WB ELISA
ESAP12124	RELA (Phospho-Thr505) Polyclonal Antibody	Q04206	H	WB IHC
ENK5382	Cryopyrin Polyclonal Antibody	Q96P20	H	WB ELISA
ESAP16727	NLRP3 Polyclonal Antibody	Q96P20	H M R	ELISA IHC
ESAP17033	OSM Polyclonal Antibody	P13725	H	ELISA IHC
EAP0759	POSTN Polyclonal Antibody	Q15063	H M R	WB
ENM0112	Peroxiredoxin 1 Monoclonal Antibody	Q06830	H M R	WB IF
EAP0840	PDGFRA Polyclonal Antibody	P16234	H M R	WB
EAP3966	PDGFRA(Phospho-Y754) Polyclonal Antibody	P16234	H	WB
EAP3968	PDGFRA(Phospho-Y1018) Polyclonal Antibody	P16234	H	WB
ENP0441	PDGFR-α (Phospho-Y849) Polyclonal Antibody	P16234	H M R	WB IHC ELISA
ENT3636	PDGFR-α Polyclonal Antibody	P16234	H M R	WB IHC IF ELISA
EPM12783	PDGFRA Monoclonal Antibody	P16234	H M;	ELISA
ENP0994	PDGFR-β (Phospho-Y1021) Polyclonal Antibody	P09619	H M R	IHC ELISA
ENT3639	PDGFR-β Polyclonal Antibody	P09619	H M R	WB IHC IF ELISA
ESAP16718	PLA2G4B Polyclonal Antibody	P0C869	H M	ELISA WB
ESAP13805	TNFRSF11A Polyclonal Antibody	Q9Y6Q6	H M	ELISA WB
EAP1495	CCL5 Polyclonal Antibody	P13501	H M	WB




# Atherosclerosis

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT4301	SIRP- $\alpha$ 1 Polyclonal Antibody	P78324	H M R	WB IHC ELISA
ENT4225	SDF-1 Polyclonal Antibody	P48061	H M R	WB IHC ELISA
ESAP13949	STAB1 Polyclonal Antibody	Q9NY15	H M	ELISA IHC
ESAP10754	SOD1 Polyclonal Antibody	P00441	H M R	ELISA WB
EAP1579	MAP3K7IP1 Polyclonal Antibody	Q15750	H M R	WB
ESAP11547	TAB1 Polyclonal Antibody	Q15750	H M	ELISA WB IHC
EAP3288	MAP3K7(Phospho-S439) Polyclonal Antibody	O43318	H M R	WB
ENP0424	Tak1 (Phospho-T187) Polyclonal Antibody	O43318	H M R	WB IHC ELISA
ENT4536	Tak1 Polyclonal Antibody	O43318	H M R	WB IHC ELISA
EAP0718	TIRAP Polyclonal Antibody	P58753	H M R	WB IHC
ESAP14944	TLR2 Polyclonal Antibody	O60603	H M	ELISA WB IHC
ESAP14088	TLR4 Polyclonal Antibody	O00206	H M R	ELISA WB IHC
ESAP16610	TLR6 Polyclonal Antibody	Q9Y2C9	H	ELISA IHC
EGM0321	TNF- $\alpha$ Polyclonal Antibody	P01375	H R M	WB IHC
EAP0448	TNFRSF1A Polyclonal Antibody	P19438	H M	WB IF
ESAP13998	TNIP2 Polyclonal Antibody	Q8NFB5	H M	ELISA WB IHC
ESAP17303	TPST2 Polyclonal Antibody	O60704	H	ELISA IHC
ENT4716	TRADD Polyclonal Antibody	Q15628	H M Mk	WB IHC ELISA
ENT4717	TRAF2 Polyclonal Antibody	Q12933	H	WB IHC ELISA
EPM12536	VCAM-1 Monoclonal Antibody	P19320	H;	ELISA WB IHC
ESAP13706	NPR2 Polyclonal Antibody	P20594	H M R	ELISA WB IHC
EAP1459	TFCP2 Polyclonal Antibody	Q12800	H M R	WB
ENT5587	CD42d Polyclonal Antibody	P40197	H M R	WB ELISA





## Blood



Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP14841	ORM1 Polyclonal Antibody	P02763	H	ELISA WB IHC
EGM0344	beta Amyloid Polyclonal Antibody	P05067	H R	WB IHC
ENP0855	Amyloid- $\beta$ (Phospho-T743) Polyclonal Antibody	P05067	H M R	WB IHC IF ELISA
ESAP12134	APP (Phospho-668) Polyclonal Antibody	P05067	H M R	WB
EAP1715	APP(Phospho-T743) Polyclonal Antibody	P05067	H M R	WB
ESAP15826	BDKRB2 Polyclonal Antibody	P30411	H M R	ELISA WB IHC
ESAP14279	CRP Polyclonal Antibody	P02741	H	ELISA WB IHC
ESAP15060	ORM2 Polyclonal Antibody	P19652	H M	ELISA WB IHC
ESAP16597	A4GALT Polyclonal Antibody	Q9NPC4	H R	ELISA IHC
ESAP16543	ART4 Polyclonal Antibody	Q93070	H	ELISA WB IHC
ESAP14183	B2M Polyclonal Antibody	P61769	H	ELISA WB IHC
ESAP13334	SLC4A1 Polyclonal Antibody	P02730	H	ELISA IHC
EAP0482	ABO Polyclonal Antibody	P16442	H M R	WB IHC IF



## Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP16120	MUC5AC Polyclonal Antibody	P98088	H	ELISA IHC
ESAP15096	KEL Polyclonal Antibody	P23276	H	ELISA IHC
EAP4070	WRB Polyclonal Antibody	O00258	H M R	WB IHC
ESAP11218	BCAM Polyclonal Antibody	P50895	H M	ELISA WB IHC
ESAP16807	FUT1 Polyclonal Antibody	P19526	H	ELISA WB IHC
EAP3859	GLA Polyclonal Antibody	P06280	H M R	WB IHC IF
ENT5638	CD235a Polyclonal Antibody	P02724	H	WB ELISA
ESAP11444	GYPC Polyclonal Antibody	P04921	H	ELISA IHC
ESAP14398	HP Polyclonal Antibody	P00738	H M R	ELISA WB IHC
EAP3883	HBA2 Polyclonal Antibody	P69905	M R	WB
EAP3885	HLA-DRB1 Polyclonal Antibody	P04229	H	WB IHC IF
ESAP15421	ICAM4 Polyclonal Antibody	Q14773	H	ELISA WB IHC
EAP2156	LILRB4 Polyclonal Antibody	Q8NHJ6	H M	WB
ESAP15527	MB Polyclonal Antibody	P02144	H M R	ELISA WB IHC
ENT5284	CD241 Polyclonal Antibody	Q02094	H	WB ELISA
ENT5088	RhCG Polyclonal Antibody	Q9UBD6	H	WB ELISA
ENT5626	CD240d Polyclonal Antibody	Q02161	H	WB ELISA
ESAP13859	SEMA7A Polyclonal Antibody	O75326	H M	ELISA WB IHC
ESAP16514	ADRA1A Polyclonal Antibody	P35348	H M R	ELISA WB IHC
ESAP16516	ADRA1B Polyclonal Antibody	P35368	H M R	ELISA WB IHC
ENT0372	AT1 Polyclonal Antibody	P30556	H M R	WB IHC IF ELISA
ESAP10241	AGTR1 Polyclonal Antibody	P30556	H R	ELISA WB
ESAP12949	AGTR2 Polyclonal Antibody	P50052	H M R	ELISA IHC
ESAP10960	ACE2 Polyclonal Antibody	Q9BYF1	H	ELISA IHC
ESAP12948	AGT Polyclonal Antibody	P01019	H M R	ELISA WB IHC



## Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT5163	Apelin Polyclonal Antibody	Q9ULZ1	H M R	WB ELISA
ENT1606	ERAP1 Polyclonal Antibody	Q9NZ08	H M R	WB IHC IF ELISA
EAP0498	NPPA Polyclonal Antibody	P01160	H M R	WB IHC IF
ENT0359	AR-β1 Polyclonal Antibody	P08588	H M R	WB IHC IF ELISA
ENK041	BMAL1 (Acetyl-K538) Polyclonal Antibody	O00327	H M R	WB ELISA
ENK5423	BMAL1 Polyclonal Antibody	O00327	H M R	WB ELISA
ESAP12481	NPPB Polyclonal Antibody	P16860	H	ELISA IHC
EAP2013	CRY2 Polyclonal Antibody	Q49AN0	H M	WB IHC IF
EPM12739	EMILIN1 Monoclonal Antibody	Q9Y6C2	H;	ELISA
ENT2070	GRK 5 Polyclonal Antibody	P34947	H M R	WB IHC ELISA
EAP3867	GNB3 Polyclonal Antibody	P16520	H M R	WB
EAP2012	CRY1 Polyclonal Antibody	Q16526	H M	WB IHC
ENT0978	Clock Polyclonal Antibody	O15516	H M R	WB IHC ELISA
EAP2176	LAP3 Polyclonal Antibody	P28838	H M R	WB IHC
ENC096	Cleaved-Kininogen-1 HC (K380) Polyclonal Antibody	P01042	H	WB ELISA
EPM12348	Kininogen 1 Monoclonal Antibody	P01042	H;	ELISA WB IHC
EAP1177	ABCC1 Polyclonal Antibody	P33527	H M R	WB IHC
ESAP12454	NPR1 Polyclonal Antibody	P16066	H M R	ELISA IHC
EAP3951	NPR3 Polyclonal Antibody	P17342	H M R	WB
ESAP16361	NPY1R Polyclonal Antibody	P25929	H R	ELISA WB IHC
EAP1909	P2RX4 Polyclonal Antibody	Q99571	H	WB IHC IF
ENP0764	Per2 (Phospho-S662) Polyclonal Antibody	O15055	H M	WB IHC ELISA
ESAP17064	PER2 Polyclonal Antibody	O15055	H	ELISA IHC
EAP1310	KLKB1 Polyclonal Antibody	P03952	H M R	WB IHC IF
ENC102	Cleaved-Plasma Kallikrein HC (R390) Polyclonal Antibody	P03952	H	WB IHC ELISA



# Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP10697	PRCP Polyclonal Antibody	P42785	H M	ELISA WB IHC
ESAP16360	NPPC Polyclonal Antibody	P23582	H M R	ELISA IHC
ESAP14967	REN Polyclonal Antibody	P00797	H	ELISA IHC
ESAP10726	RGS2 Polyclonal Antibody	P41220	H M R	ELISA IHC
ESAP15643	RGS5 Polyclonal Antibody	O15539	H M R	ELISA IHC
EPM12464	RENALASE Monoclonal Antibody	Q5VYX0	H;	ELISA WB IHC
ENT4826	Urotensin II Polyclonal Antibody	O95399	H Mk	WB ELISA
ESAP16517	SERPINF2 Polyclonal Antibody	P08697	H	ELISA IHC
EAP0483	ADORA2A Polyclonal Antibody	P29274	H	WB
EAP0471	SERPINC1 Polyclonal Antibody	P01008	H M R	WB IHC IF
ESAP14157	ANXA4 Polyclonal Antibody	P09525	H M R	ELISA WB IHC
ESAP14159	ANXA8 Polyclonal Antibody	P13928	H M R	ELISA WB IHC
ESAP14191	APOH Polyclonal Antibody	P02749	H	ELISA WB IHC
ENC020	Cleaved-C1s HC (R437) Polyclonal Antibody	P09871	H	WB ELISA
ESAP13869	SERPING1 Polyclonal Antibody	P05155	H	ELISA WB IHC
EAP3786	CPB2 Polyclonal Antibody	Q961Y4	H M R	IHC
EAP1590	PROCR Polyclonal Antibody	Q9UNN8	H M R	WB IHC
EAP0560	F12 Polyclonal Antibody	P00748	H M R	WB
ENC085	Cleaved-Factor XII HC (R372) Polyclonal Antibody	P00748	H	WB ELISA
EPM12274	Factor XII Monoclonal Antibody	P00748	H M;	ELISA WB IF IHC
ENC081	Cleaved-Factor VII LC (R212) Polyclonal Antibody	P08709	H	WB IHC ELISA
ESAP14686	F7 Polyclonal Antibody	P08709	H M	ELISA WB IHC
ESAP14688	F9 Polyclonal Antibody	P00740	H M R	ELISA WB
EAP0417	F10 Polyclonal Antibody	P00742	H M R	WB IHC IF
ENC082	Cleaved-Factor X/ Factor X LC (A41) Polyclonal Antibody	P00742	H	WB ELISA



## Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENC083	Cleaved-Factor Xa activated HC (I235) Polyclonal Antibody	P00742	H M R	WB ELISA
ENT1654	Factor XIII B Polyclonal Antibody	P05160	H	IHC IF ELISA
ESAP14343	FGL1 Polyclonal Antibody	Q08830	H M R	ELISA WB IHC
ESAP15339	FGB Polyclonal Antibody	P02675	H M R	ELISA WB IHC
EAP1502	FGG Polyclonal Antibody	P02679	H M R	WB
ESAP12593	FBLN1 Polyclonal Antibody	P23142	H	ELISA IHC
ENC077	Cleaved-Thrombin APII (R327) Polyclonal Antibody	P00734	H	WB ELISA
ENT3862	Prothrombin Polyclonal Antibody	P00734	H M R	WB IHC IF ELISA
ESAP11431	GGCX Polyclonal Antibody	P38435	H M R	ELISA WB IHC
EAP3103	GNAQ Polyclonal Antibody	P50148	H M	WB
ENT2189	HNF-4 $\alpha$ Polyclonal Antibody	P41235	H M R	WB IHC ELISA
ENP0951	HNF4- $\alpha$ (Phospho-S313) Polyclonal Antibody	P41235	H M R	WB IHC IF ELISA
ENT2190	HNF-4 $\alpha$ / $\gamma$ Polyclonal Antibody	P41235/Q14541	H M R	WB IHC ELISA
EAP3888	HRG Polyclonal Antibody	P04196	M	WB IHC
ESAP14445	KLK1 Polyclonal Antibody	P06870	H	ELISA IHC
EAP1892	MMRN1 Polyclonal Antibody	Q13201	H	IHC IF
ESAP17045	F2RL2 Polyclonal Antibody	O00254	H	ELISA IHC
ENT5221	Protein C Polyclonal Antibody	P04070	H	WB ELISA
ESAP14741	SERPINA5 Polyclonal Antibody	P05154	H	ELISA WB IHC
ESAP10839	PROS1 Polyclonal Antibody	P07225	H	ELISA IHC
ESAP13276	THBS1 Polyclonal Antibody	P07996	H M	ELISA WB IHC
ENC104	Cleaved-Thrombin R (S42) Polyclonal Antibody	P25116	H	WB ELISA
ENT4645	Thrombin R Polyclonal Antibody	P25116	H	WB ELISA
ENT5194	TF Polyclonal Antibody	P13726	H	WB ELISA
ESAP10788	F3 Polyclonal Antibody	P13726	H	ELISA WB IHC





# Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP0512	TFPI Polyclonal Antibody	P10646	H M R	WB IHC IF
ESAP14690	TFPI2 Polyclonal Antibody	P48307	H	ELISA IHC
ESAP14738	VWF Polyclonal Antibody	P04275	H	ELISA IHC
EAP1799	ALAS2 Polyclonal Antibody	P22557	H	WB IHC IF
EAP2738	AHSP Polyclonal Antibody	Q9NZD4	H	WB
EAP3737	BPGM Polyclonal Antibody	P07738	H M R	WB IHC
EAP0973	EPB41 Polyclonal Antibody	P11171	H M R	WB IHC IF
ENP0468	4.1R (Phospho-Y660) Polyclonal Antibody	P11171	H M	WB IHC ELISA
ESAP16773	EPO Polyclonal Antibody	P01588	H M R	ELISA WB IHC
EAP3839	FECH Polyclonal Antibody	P22830	H M R	WB IHC
EAP0200	HK1 Polyclonal Antibody	P19367	H M R	WB IHC IF
EPM12310	Hemoglobin epsilon-Specific Monoclonal Antibody	P02100	H;	ELISA WB
EGM0317	Hepcidin Polyclonal Antibody	P81172	H M	IHC
EAP2613	MPP1 Polyclonal Antibody	Q00013	H M	WB
ESAP17089	PLG Polyclonal Antibody	P00747	H R	ELISA WB IHC
ENT3569	PAI-1 Polyclonal Antibody	P05121	H M R	WB IHC ELISA
ESAP11743	SERPINB1 Polyclonal Antibody	P30740	H M	ELISA IHC
ESAP14551	SERPINE2 Polyclonal Antibody	P07093	H M R	ELISA WB IHC
ENK5413	TPO Polyclonal Antibody	P40225	H	WB ELISA
ESAP14580	PLAUR Polyclonal Antibody	Q03405	H	ELISA WB IHC
EAP0918	PLAU Polyclonal Antibody	P00749	M R	WB
EAP3694	ACO1 Polyclonal Antibody	P21399	H M R	WB
ENP0754	IRP-1 (Phospho-S711) Polyclonal Antibody	P21399	H M R Mk	WB IHC ELISA
EAP3707	ALAD Polyclonal Antibody	P13716	H	WB IHC
EAP1798	ALAS1 Polyclonal Antibody	P13196	H	WB IHC IF



## Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP10940	AMBP Polyclonal Antibody	P02760	H	ELISA WB IHC
ENT5612	BST-2 Polyclonal Antibody	Q10589	H	WB ELISA
ESAP14222	CA1 Polyclonal Antibody	P00915	H M R	ELISA WB
ESAP10327	CA2 Polyclonal Antibody	P00918	H M R	ELISA WB IHC
ENT0782	CD9 Polyclonal Antibody	P21926	H M R	WB IHC IF ELISA
ENK5370	L-Selectin Polyclonal Antibody	P14151	H M R	WB ELISA
EAP2070	RARRES2 Polyclonal Antibody	Q99969	H M R	WB IHC IF
ENT3721	pICln Polyclonal Antibody	P54105	H	WB ELISA
ESAP12534	CXCR5 Polyclonal Antibody	P32302	H	ELISA IHC
ENT1412	DREAM Polyclonal Antibody	Q9Y2W7	H M	WB IHC IF ELISA
ENT5562	Dkk-1 Polyclonal Antibody	O94907	H M R	WB ELISA
EPM12609	HBG1-Specific Monoclonal Antibody	P69891	H;	ELISA IHC
ESAP14845	AHSG Polyclonal Antibody	P02765	H	ELISA WB IHC
ESAP11395	FETUB Polyclonal Antibody	Q9UGM5	H	ELISA WB
ESAP15348	FSTL1 Polyclonal Antibody	Q12841	H M R	ELISA IHC
EPM12295	Gastrin Monoclonal Antibody	P01350	H;	ELISA IHC
ENT1900	Ghrelin Polyclonal Antibody	Q9UBU3	H M R	WB IHC ELISA
ENT1913	GLCNE Polyclonal Antibody	Q9Y223	H M R	WB IHC ELISA
ESAP15388	GSTP1 Polyclonal Antibody	P09211	H M R	ELISA WB IHC
EAP2264	ITLN1 Polyclonal Antibody	Q8WWA0	H	WB IHC IF
ENT2449	Kanadaplin Polyclonal Antibody	Q9BWU0	H M	WB IHC IF ELISA
EAP2235	KLF3 Polyclonal Antibody	P57682	H M R	WB IHC IF
ESAP14464	LYZ Polyclonal Antibody	P61626	H	ELISA WB IHC
ESAP16942	LTB4R2 Polyclonal Antibody	Q9NPC1	H	ELISA IHC
ESAP13650	MAPKAPK2 Polyclonal Antibody	P49137	H M	ELISA WB IHC



## Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT5123	Chymase Polyclonal Antibody	P23946	H	WB ELISA
ENC100	Cleaved-MPO 89k (A49) Polyclonal Antibody	P05164	H	WB ELISA
EPM12376	MPO Monoclonal Antibody	P05164	H;	ELISA WB IHC
EAP3936	MLF1 Polyclonal Antibody	P58340	H	WB IF
EAP0938	NOX4 Polyclonal Antibody	Q9NPH5	H M R	WB IHC IF
ESAP15099	BPI Polyclonal Antibody	P17213	H	ELISA IHC
EAP0390	GHRL Polyclonal Antibody	Q9UBU3	H M R	WB
ENT3466	OPG Polyclonal Antibody	O00300	H M	WB IHC IF ELISA
ESAP16373	PTAFR Polyclonal Antibody	P25105	H	ELISA IHC
EGM0324	PAFAH Polyclonal Antibody	Q13093	H R M	WB IHC
ESAP15033	PGK1 Polyclonal Antibody	P00558	H M R	ELISA WB IHC
EAP2058	PPID Polyclonal Antibody	Q08752	H M R	WB IHC IF
ENT4380	SP-B Polyclonal Antibody	P07988	H	WB ELISA
EAP3419	PDPN Polyclonal Antibody	Q86YL7	H	IHC
EAP0278	SH2D1A Polyclonal Antibody	O60880	H M R	WB IHC IF
ESAP13368	SIGLEC7 Polyclonal Antibody	Q9Y286	H	ELISA WB IHC
ESAP12450	TF Polyclonal Antibody	P02787	H	ELISA WB IHC
EAP1436	UROD Polyclonal Antibody	P06132	H M R	WB IHC IF
ESAP16227	UCN Polyclonal Antibody	P55089	H M R	ELISA IHC
ESAP13372	CD36 Polyclonal Antibody	P16671	H	ELISA WB IHC
ENK5393	CD42b Polyclonal Antibody	P07359	H	WB ELISA
ESAP11448	GP1BA Polyclonal Antibody	P07359	H M	ELISA WB IHC
ENT5593	P-Selectin Polyclonal Antibody	P16109	H M R	WB ELISA
ENK5525	CD63 Polyclonal Antibody	P08962	H	WB ELISA
EAP1036	PRKG1 Polyclonal Antibody	Q13976	H M R	WB IHC



## Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT1073	Cox-2 Polyclonal Antibody	P35354	H	WB IHC IF ELISA
ENP0084	Dok-1 (Phospho-Y362) Polyclonal Antibody	Q99704	H M R	WB IF ELISA
ENT1393	Dok-1 Polyclonal Antibody	Q99704	H M R	WB IHC IF ELISA
ESAP12183	DOK1 (phospho-Tyr398) Polyclonal Antibody	Q99704	H M R	WB IHC
ESAP14311	DOK1 Polyclonal Antibody	Q99704	H M R	ELISA WB
ENM0164	eNOS Monoclonal Antibody	P29474	H M R	WB
ENP0196	NOS3 (Phospho-T494) Polyclonal Antibody	P29474	H M R	WB ELISA
ENP0514	NOS3 (Phospho-S1177) Polyclonal Antibody	P29474	H M R	WB IF ELISA
ENT3171	NOS3 Polyclonal Antibody	P29474	H M R	WB ELISA
EAP0374	F11R Polyclonal Antibody	Q9Y624	H	WB IHC IF
ENK5479	JAM-A Polyclonal Antibody	Q9Y624	H R	WB ELISA
ESAP15446	F11R Polyclonal Antibody	Q9Y624	H	ELISA WB IHC
ESAP10409	LAMP2 Polyclonal Antibody	P13473	H	ELISA WB IHC
ESAP12389	MT-CO2 Polyclonal Antibody	P00403	H M R	ELISA WB IHC
ENT3630	PDGF-A Polyclonal Antibody	P04085	H M R	WB IHC ELISA
ESAP13526	GPR171 Polyclonal Antibody	O14626	H	ELISA WB IHC
EAP3984	PLEK Polyclonal Antibody	P08567	H M R	WB IHC
ESAP15688	SERPINB8 Polyclonal Antibody	P50452	H	ELISA IHC
ENP0614	Syk (Phospho-Y348) Polyclonal Antibody	P43405	H M R Mk	WB IHC ELISA
ENP0615	Syk (Phospho-Y525) Polyclonal Antibody	P43405	H M R Mk	WB IHC ELISA
ESAP10780	SYK Polyclonal Antibody	P43405	H M R	ELISA WB IHC
ESAP12225	SYK (phospho-Tyr323) Polyclonal Antibody	P43405	H	WB IHC IF
EPM12127	AFM Monoclonal Antibody	P43652	H M;	ELISA WB IHC
ESAP15752	SERPINA1 Polyclonal Antibody	P01009	H	ELISA WB IHC
ESAP12393	SERPINA3 Polyclonal Antibody	P01011	H	ELISA WB IHC





## Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP14019	AFP Polyclonal Antibody	P02771	H	ELISA WB IHC
EPM12141	Alpha-2-macroglobulin Monoclonal Antibody	P01023	H;	ELISA IHC WB
ESAP10116	ALB Polyclonal Antibody	P02768	H	ELISA WB IHC
ENT0201	ALS Polyclonal Antibody	P35858	H	WB ELISA
ESAP16520	ANGPT4 Polyclonal Antibody	Q9Y264	H	ELISA WB IHC
ESAP16521	ACE Polyclonal Antibody	P12821	H M R	ELISA WB IHC
EPM12217	Ceruloplasmin Monoclonal Antibody	P00450	H;	ELISA WB IHC
EAP3023	C8B Polyclonal Antibody	P07358	H M R	WB
EAP1860	GPLD1 Polyclonal Antibody	P80108	H M R	WB IHC IF
ESAP10494	CFHR1 Polyclonal Antibody	Q03591	H	ELISA WB IHC
EPM12251	Cytoglobin Monoclonal Antibody	Q8WWM9	H M;	ELISA WB IHC
ESAP13448	SLC11A2 Polyclonal Antibody	P49281	H M R	ELISA IHC
ESAP16758	SLC11A2 Polyclonal Antibody	P49281	H M R	ELISA IHC
ENT1653	Factor VIII Polyclonal Antibody	P00451	H M	WB IHC ELISA
EAP0607	FTL Polyclonal Antibody	P02792	H M R	WB IHC IF
EAP0279	FTH1 Polyclonal Antibody	P02794	H M R	WB IHC IF
EAP2035	HABP2 Polyclonal Antibody	Q14520	H M R	WB IHC IF
ENT5216	Kallistatin Polyclonal Antibody	P29622	H	WB ELISA
ESAP14987	LTF Polyclonal Antibody	P02788	H	ELISA WB IHC
ESAP12744	LRP1 Polyclonal Antibody	Q07954	H	ELISA WB IHC
ENT2656	MASP-2 Polyclonal Antibody	O00187	H	WB ELISA
ESAP13740	PAPPA Polyclonal Antibody	Q13219	H M	ELISA IHC
ESAP16384	PDGFB Polyclonal Antibody	P01127	H M R	ELISA IHC
ESAP15604	PLTP Polyclonal Antibody	P55058	H M	ELISA WB
ESAP13780	PON1 Polyclonal Antibody	P27169	H M R	ELISA WB IHC




## Blood

Cat. No.	Product Name	Swissprot	Reactivity	Application
EPM12446	Prealbumin/transthyretin Monoclonal Antibody	P02766	H;	ELISA WB IHC
EPM12463	RBP4 Monoclonal Antibody	P02753	H;	ELISA WB IHC
EPM12480	Serum amyloid P component Monoclonal Antibody	P02743	H;	ELISA WB IHC
EAP2180	SERPINA10 Polyclonal Antibody	Q9UK55	H	WB IHC IF
EAP1584	TCN2 Polyclonal Antibody	P20062	H M	WB
EAP2701	CLEC3B Polyclonal Antibody	P05452	H	WB IHC IF
ESAP16429	PLAT Polyclonal Antibody	P00750	H	ELISA WB IHC
EAP1671	TFRC Polyclonal Antibody	P02786	H M R	WB
ENP0729	CD71 (Phospho-S24) Polyclonal Antibody	P02786	H M	WB IHC ELISA
ENT0775	CD71 Polyclonal Antibody	P02786	H M	WB IHC ELISA
ESAP17109	TFRC Polyclonal Antibody	P02786	H M R	ELISA WB IHC
EAP0059	NAMPT Polyclonal Antibody	P43490	H M R	WB IHC IF
EGM0319	Visfatin Polyclonal Antibody	P43490	H R	WB IHC
ESAP14731	NAMPT Polyclonal Antibody	P43490	H M R	ELISA WB IHC
ESAP11068	VTN Polyclonal Antibody	P04004	H	ELISA IHC



# Cardiovascular markers

## Heart failure



Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP0928	SELE Polyclonal Antibody	P16581	H M R	WB IHC IF
EAP2258	TEK Polyclonal Antibody	Q02763	H M R	WB IHC IF
ENP0437	Tie-2 (Phospho-Y992) Polyclonal Antibody	Q02763	H M	WB ELISA
ENP0110	Flk-1 (Phospho-Y1175) Polyclonal Antibody	P35968	H M R	WB IHC ELISA
ENP0111	Flk-1 (Phospho-Y1214) Polyclonal Antibody	P35968	H M R	WB IHC IF ELISA
ENT1723	Flk-1 Polyclonal Antibody	P35968	H M	WB IHC IF ELISA
ESAP12217	KDR (phospho-Tyr1059) Polyclonal Antibody	P35968	H M	WB
ENT4869	VE-Cadherin Polyclonal Antibody	P33151	H M	WB IHC ELISA
ESAP13003	CDH5 Polyclonal Antibody	P33151	H	ELISA IHC
ESAP10325	TNNI3 Polyclonal Antibody	P19429	H M R	ELISA WB IHC
ENK5362	Troponin T-C Polyclonal Antibody	P45379	H M R	WB ELISA
ENT1861	GATA-4 Polyclonal Antibody	P43694	H M R	WB IHC IF ELISA
EAP1510	NKX2-5 Polyclonal Antibody	P52952	H	WB



## Cardiovascular markers

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT5053	α-SMA Polyclonal Antibody	P62736	H M R	WB IHC ELISA
EAP1968	TAGLN Polyclonal Antibody	Q01995	H	WB IHC IF
EAP4034	SMTN Polyclonal Antibody	P53814	H M	WB
ESAP16199	MYH11 Polyclonal Antibody	P35749	H M R	ELISA IHC

## Hypoxia



Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP0345	CA3 Polyclonal Antibody	P07451	H M R	WB IHC IF
ENM0076	CA IX Monoclonal Antibody	Q16790	H	WB IP
EPM12562	CA9 Monoclonal Antibody	Q16790	H M Pig	ELISA WB IHC IF
ENT2479	KIR6.2 Polyclonal Antibody	Q14654	H M R	WB IHC IF ELISA
ESAP15575	OS9 Polyclonal Antibody	Q13438	H M	ELISA WB I



# Hypoxia

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP1487	PGRMC1 Polyclonal Antibody	O00264	H M R	WB
ESAP11549	WDR83 Polyclonal Antibody	Q9BRX9	H M R	ELISA WB IHC
ESAP13807	HIF1AN Polyclonal Antibody	Q9NWT6	H M	ELISA WB IHC
ENT2133	HIF-1 $\alpha$ Polyclonal Antibody	Q16665	H M R	WB IHC ELISA
ENT2134	HIF-1 $\beta$ Polyclonal Antibody	P27540	H M R	WB IHC ELISA
ENT5325	EPAS-1 Polyclonal Antibody	Q99814	H M R	WB ELISA
ESAP10852	ABCB6 Polyclonal Antibody	Q9NP58	H M R	ELISA WB IHC
ESAP11314	DDIT4 Polyclonal Antibody	Q9NX09	H M R	ELISA IHC
ESAP16974	ERRF1 Polyclonal Antibody	Q9UJM3	H	ELISA WB IHC
EAP3904	IREB2 Polyclonal Antibody	P48200	H M	IF
ENT2507	KV1.5 Polyclonal Antibody	P22460	H M R	WB ELISA
EAP2054	PFKFB3 Polyclonal Antibody	Q16875	H M R	WB IHC IF
ESAP17174	SESN2 Polyclonal Antibody	P58004	H M	ELISA WB IHC
EAP0286	EGLN1 Polyclonal Antibody	Q9GZT9	H M R	WB IHC IF
ESAP10704	EGLN3 Polyclonal Antibody	Q9H6Z9	H M R	ELISA IHC





# Heart



Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP3153	MAP3K5(Phospho-S1033) Polyclonal Antibody	Q99683	H M R	WB
EAP3154	MAP3K5(Phospho-T918) Polyclonal Antibody	Q99683	H M R	WB
ENP0021	ASK 1 (Phospho-S83) Polyclonal Antibody	Q99683	H	WB IHC ELISA
ENP0022	ASK 1 (Phospho-S966) Polyclonal Antibody	Q99683	H M R	WB IHC ELISA
ENT0371	ASK 1 Polyclonal Antibody	Q99683	H	WB IHC ELISA
ESAP15022	BNIP3 Polyclonal Antibody	Q12983	H M	ELISA WB IHC
ESAP16619	BNIP3L Polyclonal Antibody	O60238	H M	ELISA IHC
ENP0781	CaMKII $\beta/\gamma/\delta$ (Phospho-T287) Polyclonal Antibody	Q13554/ Q13555/Q13557	H M R	WB IHC ELISA
ENT5063	CaMKII $\delta$ Polyclonal Antibody	Q13557	H M R	WB ELISA
ENP0279	CaMKII $\alpha/\beta/\delta$ (Phospho-T305) Polyclonal Antibody	Q9UQM7/ Q13554/Q13557	H M R	WB IHC ELISA
ENP0042	CaMKII $\alpha/\delta$ (Phospho-T286) Polyclonal Antibody	Q9UQM7/ Q13557	H M R	WB IHC IF ELISA
ENT0623	CaMKII $\alpha/\delta$ Polyclonal Antibody	Q9UQM7/ Q13557	H M R	WB IHC IF ELISA
EAP3771	CIDEA Polyclonal Antibody	O60543	H	IF



# Heart

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP12723	KLF5 Polyclonal Antibody	Q13887	H M	ELISA WB IHC
ESAP11580	MTFP1 Polyclonal Antibody	Q9UDX5	H M	ELISA WB IHC
ESAP10193	SMAD4 Polyclonal Antibody	Q13485	H M R	ELISA WB IHC
ESAP14854	AMBRA1 Polyclonal Antibody	Q9C0C7	H M	ELISA IHC
ESAP10976	ATG5 Polyclonal Antibody	Q9H1Y0	H M R	ELISA WB IHC
ESAP14881	ATG3 Polyclonal Antibody	Q9NT62	H M R	ELISA WB IHC
EAP0130	ATG7 Polyclonal Antibody	O95352	H M R	WB IHC
ESAP14880	ATG10 Polyclonal Antibody	Q9H0Y0	H M	ELISA WB IHC
ESAP10977	ATG12 Polyclonal Antibody	O94817	H	ELISA IHC
ESAP13247	ATG4A Polyclonal Antibody	Q8WYN0	H	ELISA WB IHC
ENT0394	Atg4b Polyclonal Antibody	Q9Y4P1	H M	WB IHC ELISA
ESAP14935	ATG4D Polyclonal Antibody	Q86TL0	H M	ELISA IHC
ESAP15023	ATG9A Polyclonal Antibody	Q7Z3C6	H R	ELISA WB IHC
EAP3679	ATG9B Polyclonal Antibody	Q674R7	H	WB
EAP0677	ATG16L1 Polyclonal Antibody	Q676U5	H M R	WB IHC IF IP
ESAP14934	ATG4C Polyclonal Antibody	Q96DT6	H M	ELISA WB IHC
ESAP15442	IRGM Polyclonal Antibody	A1A4Y4	H	ELISA WB IHC
EPM12354	LC3A/B Monoclonal Antibody	Q9H492/Q9GZQ8	H Pig;	ELISA IHC WB
ESAP13627	MAP1LC3C Polyclonal Antibody	Q9BXW4	H	ELISA IHC
EAP1486	MAP1LC3A Polyclonal Antibody	Q9H492	H M R	WB IHC IF
EAP2238	MAP1LC3B Polyclonal Antibody	Q9GZQ8	H M R	WB IHC IF
EAP3529	RAB24 Polyclonal Antibody	Q969Q5	H M R	WB IHC
EAP3350	SQSTM1 Polyclonal Antibody	Q13501	H M R	WB IHC
ESAP15972	SUMO1 Polyclonal Antibody	P63165	H M R	ELISA WB IHC
ESAP17291	VMP1 Polyclonal Antibody	Q96GC9	H M R	ELISA IHC



# Heart

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP3318	UVRAG Polyclonal Antibody	Q9P2Y5	H M R	IHC
ESAP15114	CASQ2 Polyclonal Antibody	O14958	H M R	ELISA IHC
ESAP13348	CACNA1C Polyclonal Antibody	Q13936	H	ELISA IHC
EAP3858	GJA5 Polyclonal Antibody	P36382	H	WB IHC IF
ESAP12969	GJA9 Polyclonal Antibody	P57773	H	ELISA WB
ENT1052	Connexin-45 Polyclonal Antibody	P36383	H M R	WB IHC ELISA
ENT1326	Desmin Polyclonal Antibody	P17661	H M R	WB IHC IF ELISA
ESAP15248	DSG2 Polyclonal Antibody	Q14126	H	ELISA IHC
ENT2475	KIR3.1 Polyclonal Antibody	P48549	H M R	WB IF ELISA
ESAP12693	HCN1 Polyclonal Antibody	O60741	H M R	ELISA WB IHC
ESAP16028	HCN2 Polyclonal Antibody	Q9UL51	H M R	ELISA IHC
EAP2562	KCND2 Polyclonal Antibody	Q9NZV8	H M R	WB
EAP2041	KCND3 Polyclonal Antibody	Q9UK17	H M R	WB IHC
EAP2175	KCNIP2 Polyclonal Antibody	Q9NS61	H M	WB
EAP0311	KCNE1 Polyclonal Antibody	P15382	H M	WB
EAP3391	KCNE2 Polyclonal Antibody	Q9Y6J6	H M	WB
EAP1150	KCNH2 Polyclonal Antibody	Q12809	H M R	WB IHC
ESAP13604	KCNG2 Polyclonal Antibody	Q9UJ96	H M R	ELISA WB IHC
ESAP15455	KCNQ1 Polyclonal Antibody	P51787	H M R	ELISA WB IHC
ESAP13102	SCN5A Polyclonal Antibody	Q14524	H M R	ELISA IHC
ENT2970	Na <sup>+</sup> CP-pan Polyclonal Antibody	P35498/Q99250/9NY46/ P35499/Q14524/9UQD0/ Q15858/Q9Y5Y9/Q9UI33	H M R	WB IHC IF ELISA
ENT2967	Na <sup>+</sup> CP type IX $\alpha$ Polyclonal Antibody	Q15858	H M R	WB IHC ELISA
ESAP16409	SCN9A Polyclonal Antibody	Q15858	H	ELISA IHC
ESAP16655	CACNA1E Polyclonal Antibody	Q15878	H M R	ELISA IHC



# Heart

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP11701	RGS4 Polyclonal Antibody	P49798	H M R	ELISA IHC
ENT3956	Rad GTPase Polyclonal Antibody	P55042	H	WB IHC IF ELISA
ENT0221	AMPK $\gamma$ 1/2/3 Polyclonal Antibody	P54619/Q9UGJ0/ Q9UGI9	H M R	WB IF ELISA
ENT4196	RyR-2 Polyclonal Antibody	Q92736	H M R	IHC IF ELISA
ESAP17149	ACTN2 Polyclonal Antibody	P35609	H M	ELISA WB IHC
ENT0102	Actinin- $\alpha$ 2/3 Polyclonal Antibody	P35609/Q08043	H M R	WB IHC IF ELISA
ENT0101	Actinin- $\alpha$ 1/2/3/4 Polyclonal Antibody	P12814/P35609/ Q08043/O43707	H M R	WB IHC ELISA
ENT4882	Vinculin Polyclonal Antibody	P18206	H M R	WB IHC IF ELISA
ESAP10243	AK1 Polyclonal Antibody	P00568	H M R	ELISA WB IHC
ESAP15126	CAT Polyclonal Antibody	P04040	H M R	ELISA WB IHC
ESAP16592	CAV3 Polyclonal Antibody	P56539	H	ELISA WB IHC
EPM12729	CKB-Specific Monoclonal Antibody	P12277	H;	ELISA
ESAP15372	GLO1 Polyclonal Antibody	Q04760	H M R	ELISA WB IHC
ENT1930	Glut4 Polyclonal Antibody	P14672	H M R	WB IHC ELISA
ENT1091	CPTI-M Polyclonal Antibody	Q92523	H	WB IHC ELISA
EAP2687	PYGB Polyclonal Antibody	P11216	H M R	WB IHC IF
ESAP15005	HK2 Polyclonal Antibody	P52789	H M R	ELISA WB IHC
ESAP17307	MED12 Polyclonal Antibody	Q93074	H M	ELISA IHC
ESAP14963	SIGMAR1 Polyclonal Antibody	Q99720	H M R	ELISA WB IHC
ENP0711	PFK-2 car (Phospho-S483) Polyclonal Antibody	O60825	H M R	WB IHC ELISA
ENT3681	PFK-2 car Polyclonal Antibody	O60825	H M R	WB IHC ELISA
EAP2199	PINK1 Polyclonal Antibody	Q9BXM7	H	WB IHC IF
EAP1515	PPARD Polyclonal Antibody	Q03181	H M R	WB IHC IF
ENT3832	PP2C $\kappa$ Polyclonal Antibody	Q8N3J5	H	WB IHC ELISA
ESAP15623	PYGM Polyclonal Antibody	P11217	H M R	ELISA WB IHC



# Heart

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP14022	PKM2 Polyclonal Antibody	P14618	H M R	ELISA WB IHC
EAP3982	PKLR Polyclonal Antibody	P30613	H M R	WB IHC
ESAP10529	IMMT Polyclonal Antibody	Q16891	H M	ELISA WB IHC
ENT3189	Nrf2 Polyclonal Antibody	Q16236	H M R	WB IHC IF ELISA
ESAP16418	SLC2A4RG Polyclonal Antibody	Q9NR83	H	ELISA WB
ESAP16189	SIRT5 Polyclonal Antibody	Q9NXA8	H M R	ELISA WB IHC
EAP4103	SUCLA2 Polyclonal Antibody	Q9P2R7	H M R	WB
ENC067	Cleaved-Notch 1 (V1754) Polyclonal Antibody	P46531	H M R	WB IHC ELISA
ENT0258	APC Polyclonal Antibody	P25054	H M R	WB IHC IF ELISA
EAP0254	APEX1 Polyclonal Antibody	P27695	H M R	WB IHC IF
ENK022	Ref-1 (Acetyl-K6) Polyclonal Antibody	P27695	H	WB ELISA
ENK023	Ref-1 (Acetyl-K7) Polyclonal Antibody	P27695	H	WB ELISA
ENT4039	Ref-1 Polyclonal Antibody	P27695	H M R	WB IHC ELISA
ESAP10963	ANKRD1 Polyclonal Antibody	Q15327	H M R	ELISA WB IHC
ENK5528	CD292 Polyclonal Antibody	P36894	H	WB ELISA
EPM12201	CD34 Monoclonal Antibody	P28906	H;	ELISA WB IHC
ESAP13414	CRELD1 Polyclonal Antibody	Q96HD1	H M R	ELISA WB IHC
EAP1832	CSRP3 Polyclonal Antibody	P50461	H M	WB IHC IF
ESAP16322	DDR2 Polyclonal Antibody	Q16832	H M	ELISA IHC
EAP0489	ETS1 Polyclonal Antibody	P14921	H M R	WB IHC IF
ENP0878	Ets-1 (Phospho-T38) Polyclonal Antibody	P14921	H M R	WB IHC IF ELISA
EAP3090	ETS2 Polyclonal Antibody	P15036	H	WB IHC IF
EAP3463	ZFPM2 Polyclonal Antibody	Q8WW38	H M	WB
EAP0117	GATA2 Polyclonal Antibody	P23769	H M R	WB
ESAP12634	GATA6 Polyclonal Antibody	Q92908	H M R	ELISA IHC





# Heart

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP2132	HAND2 Polyclonal Antibody	P61296	H M R	WB IHC IF
EAP1602	HEXIM1 Polyclonal Antibody	O94992	H M R	WB IHC IF
ENP0133	HSF1 (Phospho-S303) Polyclonal Antibody	Q00613	H M	WB IHC IF ELISA
ENT2236	HSF1 Polyclonal Antibody	Q00613	H M	WB IHC IF ELISA
ENT5330	KLF13 Polyclonal Antibody	Q9Y2Y9	H M	WB ELISA
ENP0806	DDR1 (Phospho-Y513) Polyclonal Antibody	Q08345	H M R	WB IHC ELISA
EAP1780	MTOR(Phospho-S2448) Polyclonal Antibody	P42345	H	WB
ENP0176	mTOR (Phospho-S2448) Polyclonal Antibody	P42345	H M R	WB IHC IF ELISA
ENP0384	mTOR (Phospho-T2446) Polyclonal Antibody	P42345	H M R	WB ELISA
ENT2913	mTOR Polyclonal Antibody	P42345	H M R	WB IHC ELISA
EAP0705	TFAM Polyclonal Antibody	Q00059	H M	WB IHC IF
ENP0164	MEF-2 (Phospho-T319) Polyclonal Antibody	Q02078	H M R	WB IHC ELISA
ENP0860	MEF-2 (Phospho-S408) Polyclonal Antibody	Q02078	H M R	WB IHC IF ELISA
ENT2698	MEF-2 Polyclonal Antibody	Q02078	H M R	WB IHC ELISA
ENT2701	MEF-2B Polyclonal Antibody	Q02080	H M	WB IHC ELISA
EAP1763	MEF2C(Phospho-S396) Polyclonal Antibody	Q06413	H	WB IF
ENT2702	MEF-2C Polyclonal Antibody	Q06413	H M	WB IHC ELISA
EAP2062	MAP2K5 Polyclonal Antibody	Q13163	H M R	WB IHC IF
ENT2895	MRTF-A Polyclonal Antibody	Q969V6	H M	WB IHC ELISA
ENT2905	Msx-2 Polyclonal Antibody	P35548	H M	WB IHC ELISA
EAP0447	NFATC1 Polyclonal Antibody	O95644	H M R	WB IHC IF
ESAP13034	NFAT5 Polyclonal Antibody	O94916	H	ELISA IHC
ESAP13050	NOTCH1 Polyclonal Antibody	P46531	H R	ELISA WB IHC
ENC071	Cleaved-Notch 4 (V1432) Polyclonal Antibody	Q99466	H Mk	WB ELISA
ESAP16359	NOTCH4 Polyclonal Antibody	Q99466	H	ELISA IHC



# Heart

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP13052	NOTCH3 Polyclonal Antibody	Q9UM47	H M R	ELISA WB IHC
ESAP14971	PITX2 Polyclonal Antibody	Q99697	H R	ELISA WB
ESAP12883	SRF Polyclonal Antibody	P11831	H M	ELISA WB IHC
EAP2514	SFRP2 Polyclonal Antibody	Q96HF1	H M R	WB IHC IF
ENT4311	Skp2 p45 Polyclonal Antibody	Q13309	H	WB IHC IF ELISA
EAP3149	SMAD6 Polyclonal Antibody	O43541	H M R	WB
ESAP14644	SNAI1 Polyclonal Antibody	O95863	H M	ELISA WB IHC
ENP0568	SNAI 1 (Phospho-S246) Polyclonal Antibody	O95863	H M Mk	WB IHC IF ELISA
ENT4351	SNAI 1 Polyclonal Antibody	O95863	H M Mk	WB IHC IF ELISA
ESAP17228	SIK1 Polyclonal Antibody	P57059	H M R	ELISA WB IHC
ENP0250	Stat3 (Phospho-S727) Polyclonal Antibody	P40763	H M R Mk	WB IHC IP ELISA
ENP0251	Stat3 (Phospho-Y705) Polyclonal Antibody	P40763	H M R	WB IHC IP ELISA
ENT4445	Stat3 Polyclonal Antibody	P40763	H M R	WB IHC IF ELISA
EAP0352	STIP1 Polyclonal Antibody	P31948	H M R	WB IHC IF
ESAP13661	TBX5 Polyclonal Antibody	Q99593	H M R	ELISA WB IHC
ESAP11798	TBX18 Polyclonal Antibody	O95935	H M	ELISA IHC
ENT4566	TBX15/18 Polyclonal Antibody	O95935/Q96SF7	H M	WB IHC IF ELISA
EAP1974	TEAD1 Polyclonal Antibody	P28347	H M R	WB IHC
EAP3308	TWIST1 Polyclonal Antibody	Q15672	H M R	WB IHC IF
EAP4283	VANGL2 Polyclonal Antibody	Q9ULK5	H M R	WB IHC
EPM12138	alpha cardiac Actin-specific Monoclonal Antibody	P68032	H;	ELISA WB IHC
ENM0149	α skeletal muscle actin Monoclonal Antibody	P68133	H M R	WB IP
ENP0874	Caldesmon (Phospho-S789) Polyclonal Antibody	Q05682	H M R	WB IHC IF ELISA
ESAP14213	CALD1 Polyclonal Antibody	Q05682	H R	ELISA WB IHC
EAP0706	TNNC1 Polyclonal Antibody	P63316	H M R	WB IHC IF



# Heart

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT0619	Calponin 2 Polyclonal Antibody	Q99439	H M	WB IF ELISA
ENP0070	Cofilin (Phospho-S3) Polyclonal Antibody	P23528	H M R	WB IHC ELISA
ENT1006	Cofilin Polyclonal Antibody	P23528	H M R	WB IHC IF ELISA
EAP3742	CALCB Polyclonal Antibody	P10092	H M R	WB IHC
ESAP15039	FXD1 Polyclonal Antibody	O00168	H M R	ELISA WB IHC
ENK5427	MYH-pan Polyclonal Antibody	P12882/9UKX2/11055/ Q9Y623/13533/12883/ P13535	H M R	WB ELISA
ESAP14504	MYH7 Polyclonal Antibody	P12883	H M R	ELISA IHC
ESAP13024	TRIM63 Polyclonal Antibody	Q969Q1	H M	ELISA WB IHC
ENT2926	MYBPC1 Polyclonal Antibody	Q00872	H M R	WB ELISA
EPM12384	Myosin 2a Monoclonal Antibody	Q9UKX2	H;	ELISA WB IHC
ENT2939	MYLIP Polyclonal Antibody	Q8WY64	H M R	WB IHC ELISA
ESAP10594	MYLK2 Polyclonal Antibody	Q9H1R3	H M R	ELISA WB
ESAP14996	MYLK3 Polyclonal Antibody	Q32MK0	H M R	ELISA IHC
EAP2718	MYOT Polyclonal Antibody	Q9UBF9	H	WB IHC IF
ESAP16125	MYL9 Polyclonal Antibody	P24844	H M R	ELISA WB IHC
ESAP11583	MYL3 Polyclonal Antibody	P08590	H M R	ELISA WB IHC
ESAP14501	MYOZ1 Polyclonal Antibody	Q9NP98	H M	ELISA WB IHC
EAP2741	MYOZ2 Polyclonal Antibody	Q9NPC6	H	WB IHC IF
ENT5103	NCX1 Polyclonal Antibody	P32418	H M R	WB ELISA
EAP3377	PDE4DIP Polyclonal Antibody	Q5VU43	H M R	WB IHC
EAP0241	ATP2A2 Polyclonal Antibody	P16615	H M R	WB IHC IF
EAP4041	SRI Polyclonal Antibody	P30626	H M R	IHC
ESAP17265	TMSB4X Polyclonal Antibody	P62328	H M R	ELISA IHC
EAP3085	ENPEP Polyclonal Antibody	Q07075	H M R	WB IHC



# Heart

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP2113	CAMK1 Polyclonal Antibody	Q14012	H M R	IHC IF
EAP0989	CAMK2B Polyclonal Antibody	Q13554	H M R	WB
EAP3743	CAMK2G Polyclonal Antibody	Q13555	H M R	WB IHC IF
ENT0625	CaMKII $\beta/\gamma$ Polyclonal Antibody	Q13555/Q13554	H M R	WB IF ELISA
ESAP16566	DGKZ Polyclonal Antibody	Q13574	H M R	ELISA IHC
EAP4337	FBXO32 Polyclonal Antibody	Q969P5	H M R	WB
ESAP15336	FHL1 Polyclonal Antibody	Q13642	H M R	ELISA WB IHC
EAP0006	GDF15 Polyclonal Antibody	Q99988	H R	WB
EAP2670	GRK6 Polyclonal Antibody	P43250	H M R	WB IHC IF
ENT2145	Histone deacetylase 1 Polyclonal Antibody	Q13547	H M R	WB IHC IF ELISA
EAP0821	HDAC2 Polyclonal Antibody	Q92769	H M R	WB IHC IF
ESAP12135	HDAC2 (Phospho-Ser394) Polyclonal Antibody	Q92769	H M R	WB IHC IF
EAP0876	HDAC3 Polyclonal Antibody	O15379	H M R	WB IHC IF
ENT2114	HDAC3 Polyclonal Antibody	O15379	H M R	WB IHC IF ELISA
EAP0042	HDAC4 Polyclonal Antibody	P56524	H M R	WB IHC IF
ENP0125	HDAC4 (Phospho-S632) Polyclonal Antibody	P56524	H M R	WB ELISA
ESAP12214	HDAC4/HDAC5/HDAC9 (phospho-Ser246/259/220) Polyclonal Antibody	P56524/Q9UQL6/Q9UKV0	H	WB IHC
EAP2231	HDAC5 Polyclonal Antibody	Q9UQL6	H M	WB IHC IF
ENP0126	HDAC5 (Phospho-S498) Polyclonal Antibody	Q9UQL6	H M	WB IHC IF ELISA
ENT3217	NY-CO-9 Polyclonal Antibody	Q9UQL6	H M R	WB IHC IF ELISA
ENP0299	HDAC5/9 (Phospho-S259/220) Polyclonal Antibody	Q9UQL6/Q9UKV0	H M	WB IHC ELISA
ENT2117	HDAC5/9 Polyclonal Antibody	Q9UQL6/Q9UKV0	H M	WB IHC ELISA
EAP2288	HDAC7 Polyclonal Antibody	Q8WUI4	H M R	WB IHC IF
ENP0127	HDAC8 (Phospho-S39) Polyclonal Antibody	Q9BY41	H M R	WB IHC ELISA
ESAP14392	HDAC8 Polyclonal Antibody	Q9BY41	H M R	ELISA WB IHC




# Heart

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP0436	HDAC9 Polyclonal Antibody	Q9UKV0	H M R	WB
ENT2149	Histone deacetylase 9 Polyclonal Antibody	Q9UKV0-5 6 7	H	WB IHC IF ELISA
ENP0014	AP-1/Jun D (Phospho-S73/100) Polyclonal Antibody	P05412/P17535	H M R	WB IHC ELISA
ENT0246	AP-1/Jun D Polyclonal Antibody	P05412/P17535	H M R	WB IHC ELISA
ESAP12505	JUND Polyclonal Antibody	P17535	H M R	ELISA WB IHC
EAP3366	KLF4 Polyclonal Antibody	O43474	H M R	WB IHC
EAP2699	KLF10 Polyclonal Antibody	Q13118	H	IHC IF
EAP3500	KLF15 Polyclonal Antibody	Q9UIH9	H	WB IHC IF
ESAP16233	NFATC2 Polyclonal Antibody	Q13469	H M	ELISA IHC
ENP0339	PKC $\alpha$ (Phospho-Y657) Polyclonal Antibody	P17252	H M R	WB ELISA
ENP0703	PKC $\alpha$ (Phospho-T638) Polyclonal Antibody	P17252	H M R	WB IHC ELISA
EAP3216	PKC $\alpha$ /beta II(Phospho-T638/641) Polyclonal Antibody	P17252/P05771	H	WB
EAP3217	PKC $\alpha$ /beta II(Phospho-T497) Polyclonal Antibody	P17252/P05771	H	WB
ENP0227	PKC (Phospho-T497) Polyclonal Antibody	P17252/P05771/ Q05655/Q02156/ P05129/P24723/ Q04759/Q05513	H M R	WB IHC IF ELISA
ENT3752	PKC Polyclonal Antibody	P17252/P05771/ Q05655/Q02156/ P05129/P24723/ Q04759/Q05513	H M R	WB IHC IF ELISA
ENP0695	Krs-1/2 (Phospho-T183) Polyclonal Antibody	Q13188/Q13043	H M	WB IHC ELISA
EAP1972	TCAP Polyclonal Antibody	O15273	H	WB IHC IF
ESAP14718	TPM1 Polyclonal Antibody	P09493	H M R	ELISA WB IHC
ESAP16466	TRPC3 Polyclonal Antibody	Q13507	H M R	ELISA WB IHC
EAP2406	WISP2 Polyclonal Antibody	O76076	H M R	WB
ESAP14737	BMPR2 Polyclonal Antibody	Q13873	H M	ELISA IHC





## Lipids and lipoproteins



Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP14762	ABHD5 Polyclonal Antibody	Q8WTS1	H M R	ELISA WB IHC
EAP1016	ADIPOQ Polyclonal Antibody	Q15848	H M R	WB
EGM0310	Acrp 30 Polyclonal Antibody	Q15848	R	WB IHC
ESAP12403	ADIPOQ Polyclonal Antibody	Q15848	H M	ELISA WB
ESAP14142	ADIPOR1 Polyclonal Antibody	Q96A54	H M	ELISA WB IHC
ESAP15173	CIDEC Polyclonal Antibody	Q96AQ7	H M R	ELISA WB IHC
ESAP12915	ABCD2 Polyclonal Antibody	Q9UBJ2	H M R	ELISA IHC
EAP2373	ACAA1 Polyclonal Antibody	P09110	H M R	WB IHC
ESAP14765	ACAA2 Polyclonal Antibody	P42765	H M R	ELISA WB IHC
EAP3690	ACADL Polyclonal Antibody	P28330	H M	WB
EAP0678	ACADM Polyclonal Antibody	P11310	H M R	WB IHC
ESAP14770	ACADS Polyclonal Antibody	P16219	H M R	ELISA WB IHC
ESAP12921	ACADVL Polyclonal Antibody	P49748	H	ELISA WB IHC



## Lipids and lipoproteins

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENP0595	ACC $\alpha$ (Phospho-S80) Polyclonal Antibody	Q13085	H M R	WB IHC ELISA
ENT0074	ACC $\alpha$ Polyclonal Antibody	Q13085	H M R	WB IHC ELISA
ESAP14774	ACOT11 Polyclonal Antibody	Q8WXI4	H M	ELISA WB IHC
ESAP10873	ACOX1 Polyclonal Antibody	Q15067	H M R	ELISA WB IHC
ESAP10875	ACOX2 Polyclonal Antibody	Q99424	H M	ELISA WB IHC
ESAP14776	ACOX3 Polyclonal Antibody	O15254	H	ELISA IHC
ESAP10877	ACSM5 Polyclonal Antibody	Q6NUN0	H M	ELISA IHC
ENT0092	ACSL6 Polyclonal Antibody	Q9UKU0	H M R	WB IHC ELISA
ENT0087	ACOT2 Polyclonal Antibody	P49753	H	WB IHC IF ELISA
ENT0088	ACOT4 Polyclonal Antibody	Q8N9L9	H	WB IHC IF ELISA
ENT0089	ACOT8 Polyclonal Antibody	O14734	H	WB IHC IF ELISA
EAP1796	AGPAT2 Polyclonal Antibody	O15120	H	WB IHC IF
ESAP12935	AGRP Polyclonal Antibody	O00253	H	ELISA IHC
EAP4273	ALG1 Polyclonal Antibody	Q9BT22	H	IHC IF
EAP0362	PRKAA1 Polyclonal Antibody	Q13131	H M R	WB IHC IF
EAP1785	AMPK $\alpha$ (Phospho-T172) Polyclonal Antibody	Q13131	H M R	WB
ENP0010	AMPK $\alpha$ 1 (Phospho-S496) Polyclonal Antibody	Q13131	H M R	WB IHC ELISA
ENT0215	AMPK $\alpha$ 1 Polyclonal Antibody	Q13131	H M R	WB IHC IF ELISA
ESAP16519	PRKAA1 Polyclonal Antibody	Q13131	H M R	ELISA WB IHC
ENP0575	AMPK $\alpha$ 1/2 (Phospho-T183/172) Polyclonal Antibody	Q13131/P54646	H M R Mk	WB IHC ELISA
ESAP14860	PRKAA2 Polyclonal Antibody	P54646	H M R	ELISA IHC
ENP0011	AMPK $\beta$ 1 (Phospho-S182) Polyclonal Antibody	Q9Y478	H M R	WB IHC ELISA
ENT0218	AMPK $\beta$ 1 Polyclonal Antibody	Q9Y478	H M R Mk	WB IHC ELISA
EAP2061	PRKAB2 Polyclonal Antibody	O43741	H M R	WB IHC IF
ENT0220	AMPK $\gamma$ 1 Polyclonal Antibody	P54619	H M R	WB ELISA



## Lipids and lipoproteins

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP15777	PRKAG1 Polyclonal Antibody	P54619	H	ELISA WB
EAP1306	FABP3 Polyclonal Antibody	P05413	H M R	WB IHC IF
ENP0463	C/EBP $\alpha$ (Phospho-T230) Polyclonal Antibody	P49715	H M R	WB ELISA
ENP0597	C/EBP $\alpha$ (Phospho-S21) Polyclonal Antibody	P49715	H M R	WB IHC ELISA
ENT0551	C/EBP $\alpha$ Polyclonal Antibody	P49715	H M R	WB IHC ELISA
EAP1302	CPT1A Polyclonal Antibody	P50416	H M R	WB IHC IF
EPM12236	CPT1C Monoclonal Antibody	Q8TCG5	H;	ELISA WB IHC
EAP1110	CPT2 Polyclonal Antibody	P23786	H M R	WB IHC IF
EAP2659	CR Polyclonal Antibody	P43155	H M	WB IHC IF
ESAP11271	CROT Polyclonal Antibody	Q9UKG9	H M R	ELISA WB IHC
EAP2381	CYP2C18 Polyclonal Antibody	P33260	H M R	WB IHC
ENT1224	CYP2U1 Polyclonal Antibody	Q7Z449	H M R	WB IHC IF ELISA
ENT1221	CYP2J2 Polyclonal Antibody	P51589	H Mk	WB IHC IF ELISA
ENT1232	CYP4A11/22 Polyclonal Antibody	Q02928/Q5TCH4	H	WB ELISA
EAP0344	ECI1 Polyclonal Antibody	P42126	H M R	WB IHC IF
ENT5229	DCI Polyclonal Antibody	P42126	H	WB ELISA
ESAP11319	DEGS1 Polyclonal Antibody	O15121	H M R	ELISA WB IHC
ESAP13443	DGAT1 Polyclonal Antibody	O75907	H M R	ELISA WB IHC
ESAP11340	ECHS1 Polyclonal Antibody	P30084	H M R	ELISA WB IHC
EAP2526	EHHADH Polyclonal Antibody	Q08426	H M R	WB IHC IF
ENT3610	PBFE Polyclonal Antibody	Q08426	H R	WB IHC ELISA
ESAP11346	EHHADH Polyclonal Antibody	Q08426	H M	ELISA WB IHC
ENT1538	ELOVL4 Polyclonal Antibody	Q9GZR5	H M	WB ELISA
EAP0684	EPHX2 Polyclonal Antibody	P34913	H M R	WB IHC
EAP0035	FABP4 Polyclonal Antibody	P15090	H M R	WB IHC IF



## Lipids and lipoproteins

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP11375	FADS1 Polyclonal Antibody	O60427	H M R	ELISA WB IHC
ESAP15317	FAR2 Polyclonal Antibody	Q96K12	H M	ELISA WB IHC
ESAP13117	SLC27A1 Polyclonal Antibody	Q6PCB7	H	ELISA IHC
ESAP15694	SLC27A2 Polyclonal Antibody	O14975	H M R	ELISA IHC
EAP2666	FABP5 Polyclonal Antibody	Q01469	H M	WB IHC IF
EAP0106	FASN Polyclonal Antibody	P49327	H M R	WB IHC IF
ENT1683	Fatty Acid Synthase Polyclonal Antibody	P49327	H	WB ELISA
EAP1555	GPD1 Polyclonal Antibody	P21695	H M R	WB IHC
ESAP13523	GNPAT Polyclonal Antibody	O15228	H M R	ELISA IHC
EAP1858	GPAM Polyclonal Antibody	Q9HCL2	H M	WB IHC
EAP1332	HADHA Polyclonal Antibody	P40939	H M R	WB IHC IF
EAP2743	HAO1 Polyclonal Antibody	Q9UJM8	H M	WB IHC
ENK5386	17 $\beta$ -HSD4 Polyclonal Antibody	P51659	H M R	WB ELISA
ESAP15408	HSD17B4 Polyclonal Antibody	P51659	H M R	ELISA WB IHC
ESAP16890	FABP2 Polyclonal Antibody	P12104	H M R	ELISA WB IHC
ENP1200	Ob-R (Phospho-Y1141) Polyclonal Antibody	P48357	H M	WB ELISA
ENT3226	Ob-R Polyclonal Antibody	P48357	H M Mk	WB IHC ELISA
ESAP16350	LCN2 Polyclonal Antibody	P80188	H	ELISA WB IHC
ESAP12673	LRP2 Polyclonal Antibody	P98164	H M R	ELISA IHC
ESAP13694	NDUFAB1 Polyclonal Antibody	O14561	H M	ELISA IHC
ENT1415	DRS-1 Polyclonal Antibody	O75521	H	WB IHC ELISA
ENT2383	iPLA2 $\gamma$ Polyclonal Antibody	Q9NP80	H M	WB ELISA
ESAP13084	PPARA Polyclonal Antibody	Q07869	H M R	ELISA WB IHC
ESAP13846	SCD Polyclonal Antibody	O00767	H	ELISA WB IHC
ESAP17176	SGPL1 Polyclonal Antibody	O95470	H M R	ELISA IHC



## Lipids and lipoproteins

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT0095	ACSVL6 Polyclonal Antibody	Q9Y2P5	H	WB IHC IF ELISA
ESAP13894	SLC27A5 Polyclonal Antibody	Q9Y2P5	H M	ELISA WB IHC
ESAP13181	AADA4L4 Polyclonal Antibody	Q5VUY2	H	ELISA WB
ENT0043	ABCA6 Polyclonal Antibody	Q8N139	H	WB ELISA
ESAP14758	ABCG4 Polyclonal Antibody	Q9H172	H M	ELISA WB
EAP2585	PNPLA2 Polyclonal Antibody	Q96AD5	H M R	WB
ESAP10922	AKR1D1 Polyclonal Antibody	P51857	H M R	ELISA WB IHC
EAP2303	PRKAG1 Polyclonal Antibody	P54619	H M R	WB IHC IF
ESAP16776	MSMO1 Polyclonal Antibody	Q15800	H	ELISA IHC
ENT1188	CYP11B1/2 Polyclonal Antibody	P15538/P19099	H	WB ELISA
ESAP13296	CAV1 Polyclonal Antibody	Q03135	H M R	ELISA WB IHC
EAP2658	CLPS Polyclonal Antibody	P04118	H M	WB IHC IF
ENT1166	CYB5R3 Polyclonal Antibody	P00387	H	WB IHC ELISA
ESAP14292	CYP1B1 Polyclonal Antibody	Q16678	H M R	ELISA WB
ENT1207	CYP2A7 Polyclonal Antibody	P20853	H	WB IHC ELISA
ENT1203	CYP2A13 Polyclonal Antibody	Q16696	H	WB IHC IF ELISA
ENT1212	CYP2C8 Polyclonal Antibody	P10632	H	WB ELISA
ENT1213	CYP2C8/9/18/19 Polyclonal Antibody	P10632/P11712/ P33260/P33261	H	WB IHC IF ELISA
ENT1222	CYP2R1 Polyclonal Antibody	Q6VVX0	H M Mk	WB IHC ELISA
ENT1223	CYP2S1 Polyclonal Antibody	Q96SQ9	H M	WB IHC IF ELISA
ENT1225	CYP2W1 Polyclonal Antibody	Q8TAV3	H	WB IHC IF ELISA
ENT1231	CYP3A7 Polyclonal Antibody	P24462	H	WB ELISA
ENT1233	CYP4B1 Polyclonal Antibody	P13584	H	WB IHC ELISA
ENT1235	CYP4F2 Polyclonal Antibody	P78329	H	WB IHC ELISA
EAP4305	CYP4F12 Polyclonal Antibody	Q9HCS2	H M R	WB





## Lipids and lipoproteins

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT0095	ACSVL6 Polyclonal Antibody	Q9Y2P5	H	WB IHC IF ELISA
ESAP13894	SLC27A5 Polyclonal Antibody	Q9Y2P5	H M	ELISA WB IHC
ESAP13181	AADACL4 Polyclonal Antibody	Q5VUY2	H	ELISA WB
ENT0043	ABCA6 Polyclonal Antibody	Q8N139	H	WB ELISA
ESAP14758	ABCG4 Polyclonal Antibody	Q9H172	H M	ELISA WB
EAP2585	PNPLA2 Polyclonal Antibody	Q96AD5	H M R	WB
ESAP10922	AKR1D1 Polyclonal Antibody	P51857	H M R	ELISA WB IHC
EAP2303	PRKAG1 Polyclonal Antibody	P54619	H M R	WB IHC IF
ESAP16776	MSMO1 Polyclonal Antibody	Q15800	H	ELISA IHC
ENT1188	CYP11B1/2 Polyclonal Antibody	P15538/P19099	H	WB ELISA
ESAP13296	CAV1 Polyclonal Antibody	Q03135	H M R	ELISA WB IHC
EAP2658	CLPS Polyclonal Antibody	P04118	H M	WB IHC IF
ENT1166	CYB5R3 Polyclonal Antibody	P00387	H	WB IHC ELISA
ESAP14292	CYP1B1 Polyclonal Antibody	Q16678	H M R	ELISA WB
ENT1207	CYP2A7 Polyclonal Antibody	P20853	H	WB IHC ELISA
ENT1203	CYP2A13 Polyclonal Antibody	Q16696	H	WB IHC IF ELISA
ENT1212	CYP2C8 Polyclonal Antibody	P10632	H	WB ELISA
ENT1213	CYP2C8/9/18/19 Polyclonal Antibody	P10632/P11712/ P33260/P33261	H	WB IHC IF ELISA
ENT1222	CYP2R1 Polyclonal Antibody	Q6VVX0	H M Mk	WB IHC ELISA
ENT1223	CYP2S1 Polyclonal Antibody	Q96SQ9	H M	WB IHC IF ELISA
ENT1225	CYP2W1 Polyclonal Antibody	Q8TAV3	H	WB IHC IF ELISA
ENT1231	CYP3A7 Polyclonal Antibody	P24462	H	WB ELISA
ENT1233	CYP4B1 Polyclonal Antibody	P13584	H	WB IHC ELISA
ENT1235	CYP4F2 Polyclonal Antibody	P78329	H	WB IHC ELISA
EAP4305	CYP4F12 Polyclonal Antibody	Q9HCS2	H M R	WB



## Lipids and lipoproteins

Cat. No.	Product Name	Swissprot	Reactivity	Application
ENT1228	CYP3A4 Polyclonal Antibody	P08684	H	WB IHC ELISA
ENT1229	CYP3A4/5 Polyclonal Antibody	P08684/P20815	H	WB IHC IF ELISA
EAP3803	CYP3A5 Polyclonal Antibody	P20815	H	WB
ESAP11294	CYP4A11 Polyclonal Antibody	Q02928	H M	ELISA WB IHC
ENT1189	CYP17A1 Polyclonal Antibody	P05093	H	WB IHC IF ELISA
EAP3987	POR Polyclonal Antibody	P16435	H M R	WB IHC
ESAP11323	DHCR7 Polyclonal Antibody	Q9UBM7	H M	ELISA IHC
EAP2020	DLST Polyclonal Antibody	P36957	H M R	WB IHC
EAP2575	FDFT1 Polyclonal Antibody	P37268	H M R	WB IHC IF
EAP1575	FDPS Polyclonal Antibody	P14324	H M R	WB IHC
ESAP13507	GALC Polyclonal Antibody	P54803	H M	ELISA WB
ENT2107	HCCS Polyclonal Antibody	P53701	H M Mk	WB IHC IF ELISA
ESAP13544	HMGCS1 Polyclonal Antibody	Q01581	H	ELISA WB IHC
ESAP16847	HMGCS2 Polyclonal Antibody	P54868	H M	ELISA WB IHC
ENP0636	HSL (Phospho-S552) Polyclonal Antibody	Q05469	H M R	WB IHC ELISA
ENT2240	HSL Polyclonal Antibody	Q05469	H M R	WB IHC ELISA
ESAP13581	IMPAD1 Polyclonal Antibody	Q9NX62	H M R	ELISA IHC
ESAP11499	INSIG2 Polyclonal Antibody	Q9Y5U4	H M R	ELISA IHC
EAP0412	LIPC Polyclonal Antibody	P11150	H	WB IHC IF
ENT2126	Hepatic Lipase Polyclonal Antibody	P11150	H	WB IHC ELISA
EAP2111	LIPF Polyclonal Antibody	P07098	H M R	WB
EPM12358	LIPH Monoclonal Antibody	Q8WWY8	H;	ELISA WB
ENT2586	LPD lipase Polyclonal Antibody	Q6XZB0	H	WB IHC ELISA
ESAP15012	OLR1 Polyclonal Antibody	P78380	H	ELISA WB IHC
EAP3145	LRP5 Polyclonal Antibody	O75197	H M R	WB IHC



## Lipids and lipoproteins

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP16084	LRPAP1 Polyclonal Antibody	P30533	H M R	ELISA WB IHC
EAP2043	LSS Polyclonal Antibody	P48449	H M	WB
EAP0878	NR1H3 Polyclonal Antibody	Q13133	H M R	WB
ENT5143	LXRα Polyclonal Antibody	Q13133	H M R	WB ELISA
EPM12398	NR1H3 Monoclonal Antibody	Q13133	H R;	ELISA WB IF
EPM12656	NR1H2 Monoclonal Antibody	P55055	H Pig;	ELISA WB
ENT2674	MC2-R Polyclonal Antibody	Q01718	H M	WB IF ELISA
ESAP16097	MC4R Polyclonal Antibody	P32245	H M	ELISA IHC
EAP4179	STARD3 Polyclonal Antibody	Q14849	H M R	WB
ESAP16981	MGLL Polyclonal Antibody	Q99685	H	ELISA IHC
ESAP16988	MVD Polyclonal Antibody	P53602	H M	ELISA IHC
ESAP15524	MVK Polyclonal Antibody	Q03426	H M R	ELISA WB IHC
EAP1386	NPC2 Polyclonal Antibody	P61916	H M R	WB IHC IF
EAP2677	OGDH Polyclonal Antibody	Q02218	H M R	WB IHC IF
ESAP15585	PNLIP Polyclonal Antibody	P16233	H M	ELISA WB IHC
ESAP15507	PCK2 Polyclonal Antibody	Q16822	H M	ELISA WB IHC
ESAP17057	PCSK9 Polyclonal Antibody	Q8NBP7	H	ELISA WB
EAP1917	PLA2G2D Polyclonal Antibody	Q9UNK4	H	WB IHC IF
EAP3197	PLCG1 Polyclonal Antibody	P19174	H	WB
ENP0234	PLC γ1 (Phospho-Y783) Polyclonal Antibody	P19174	H M R Mk	WB IHC ELISA
ENP0607	PLC γ1 (Phospho-Y771) Polyclonal Antibody	P19174	H M R Mk	WB IHC ELISA
ENT3792	PLC γ1 Polyclonal Antibody	P19174	H M R Mk	WB IHC ELISA
ESAP12098	PLCG1 (Phospho-Tyr783) Polyclonal Antibody	P19174	H M R	WB
ESAP12215	PLCG1 (phospho-Tyr771) Polyclonal Antibody	P19174	H M R	WB IF
ENT3783	PLA1A Polyclonal Antibody	Q53H76	H	WB IHC ELISA




## Lipids and lipoproteins

Cat. No.	Product Name	Swissprot	Reactivity	Application
EAP4195	PLA2R1 Polyclonal Antibody	Q13018	H M	WB
EAP0919	PLCG2 Polyclonal Antibody	P16885	H M R	WB IHC IF
ESAP12127	PLCG2 (Phospho-Tyr753) Polyclonal Antibody	P16885	H M R	WB IF
ENT3622	PC-PLD4 Polyclonal Antibody	Q96BZ4	H M R	WB ELISA
ESAP16831	HCAR2 Polyclonal Antibody	Q8TDS4	H M R	ELISA IHC
ESAP15681	SCAP Polyclonal Antibody	Q12770	H M R	ELISA IHC
ESAP14386	SCARB1 Polyclonal Antibody	Q8WTV0	H	ELISA IHC
ESAP17160	SEC14L2 Polyclonal Antibody	O76054	H M	ELISA WB IHC
ESAP15682	DHCR24 Polyclonal Antibody	Q15392	H M R	ELISA IHC
ESAP11047	SPTLC1 Polyclonal Antibody	O15269	H	ELISA WB IHC
ENT4411	SREBP-1 Polyclonal Antibody	P36956	H R	WB IHC ELISA
ESAP16209	SREBF2 Polyclonal Antibody	Q12772	H M R	ELISA IHC
ESAP17250	STARD5 Polyclonal Antibody	Q9NSY2	H M	ELISA WB IHC
ESAP11784	SULT2B1 Polyclonal Antibody	O00204	H	ELISA WB IHC
ESAP17288	TM7SF2 Polyclonal Antibody	O76062	H M	ELISA WB IHC
ENT0147	AI-BP Polyclonal Antibody	Q8NCW5	H M R	WB IF ELISA
ESAP16624	APOL6 Polyclonal Antibody	Q9BWW8	H	ELISA IHC
ESAP16275	APOA5 Polyclonal Antibody	Q6Q788	H	ELISA IHC
ENT5550	ApoC-III Polyclonal Antibody	P02656	H	WB ELISA
EPM12151	APOD Monoclonal Antibody	P05090	H;	ELISA WB
ESAP16527	APOL2 Polyclonal Antibody	Q9BQE5	H	ELISA WB IHC
ENT0981	Clusterin Polyclonal Antibody	P10909	H M R	WB IHC IF ELISA
ENT1167	Cyclin A Polyclonal Antibody	P78396/P20248	H M R Mk	WB IHC IF ELISA
EAP3035	CCNA2 Polyclonal Antibody	P20248	H M R	WB IHC
ESAP10569	LRP12 Polyclonal Antibody	Q9Y561	H M	ELISA WB IHC
EAP0529	PON2 Polyclonal Antibody	Q15165	H M R	WB
EAP1918	PLSCR1 Polyclonal Antibody	O15162	H	WB IHC IF



## Vasculature



Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP13262	DDAH1 Polyclonal Antibody	O94760	H M R	ELISA WB IHC
ESAP14315	ECE1 Polyclonal Antibody	P42892	H M R	ELISA IHC
ESAP14317	EGFL7 Polyclonal Antibody	Q9UHF1	H M	ELISA IHC
EAP0125	EDN1 Polyclonal Antibody	P05305	H M R	WB IHC
EAP4288	JAM2 Polyclonal Antibody	P57087	H M R	WB
ENT5115	JAM-B Polyclonal Antibody	P57087	H M R	WB ELISA
ESAP10573	LYVE1 Polyclonal Antibody	Q9Y5Y7	H	ELISA WB IHC
ENT2648	MARCKS Polyclonal Antibody	P29966	H M R	WB IHC ELISA
EAP0391	SOD2 Polyclonal Antibody	P04179	H M R	WB IHC IF
ENT4649	Tie-1 Polyclonal Antibody	P35590	H M R	WB ELISA
EAP3355	TNFSF18 Polyclonal Antibody	Q9UNG2	H M R	WB
ESAP10433	ELN Polyclonal Antibody	P15502	H	ELISA WB IHC
ESAP14314	ECE2 Polyclonal Antibody	O60344	H	ELISA IHC



## Vasculature

Cat. No.	Product Name	Swissprot	Reactivity	Application
ESAP11364	EDN2 Polyclonal Antibody	P20800	H	ELISA WB IHC
ENT2666	MaxiK $\alpha$ Polyclonal Antibody	Q12791	H M R	WB IHC ELISA
ENT3067	Neuromedin-S Polyclonal Antibody	Q5H8A3	H	WB IHC ELISA
ESAP13323	GPR182 Polyclonal Antibody	O15218	H M	ELISA WB IHC
ESAP14269	CNGA2 Polyclonal Antibody	Q16280	H M R	ELISA WB IHC
ENT1885	GCS- $\alpha$ -1 Polyclonal Antibody	Q02108	H M R	WB IHC ELISA
EAP1457	PNMT Polyclonal Antibody	P11086	H M R	WB IHC
EAP0665	PTGIR Polyclonal Antibody	P43119	H M	WB IHC
ENT5574	VPAC1 Polyclonal Antibody	P32241	H	WB ELISA



The background of the image is a blurred, microscopic view of biological structures. It features a prominent, thick, yellowish, fibrous structure that curves across the frame. Surrounding and interspersed with this are numerous smaller, bright red, oval-shaped structures, which appear to be cells or organelles. The overall color palette is dominated by deep reds and bright yellows against a dark, almost black background.

**Elabscience<sup>®</sup>**

---

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Website: [www.elabscience.com](http://www.elabscience.com)