

Elabscience®

# Popular Antibodies

Focus on your research Service for life science



[www.elabscience.com](http://www.elabscience.com)



## About Elabscience

Elabscience specializes in immunodiagnostic technology for life science community. We have complete platform for R&D and manufacture. At the same time, we have in house QC for each antibody, endeavoring to keep your experiment results more consistent and precise. Most of Elabscience products have been used by researchers all around the world. An increasing number of citations are available in scientific publications.

Elabscience major products cover Proteins, Antibodies, ELISA Kits, CLIA Kits, Labeling Kits, Biochemical Kits and Other Reagents.

Elabscience also offers custom services including Protein Service, Antibody Service, Peptide Service, Gene Synthesis Service and Whole Research Project Solutions.




### >> Our Advantages

- Rich range of products
- Superior quality with high stability
- Detailed product information
- Fast delivery and reliable supply chain
- Competitive price
- Strict quality control standards
- Professional and efficient technical support

### >> Antibody Category

- Polyclonal Antibodies
- Monoclonal Antibodies
- Tags & Controls Antibodies
- Phosphorylation Antibodies
- Secondary Antibodies
- Flow Cytometry Antibodies



## Why choose Elabscience Antibodies?

1. Specialization: Specialize for Breast cancer, Liver cancer and Kidney injury.
2. Variable categories: Monoclonal / Polyclonal / Tag & Control / Flow Cytometry / Acetyl / Phospho / Methyl Antibodies.
3. Wide range of reactivity: Human, Mouse, Rat, Monkey, etc.
4. Superior quality: High specificity and affinity, good consistency and stability.
5. Detailed product parameters: All antibodies are with detailed information and authentic data to bring you comfortable experience during the whole process of ordering and experiments.
6. Rich applications: WB/IHC/IF/IP/FC/ELISA, etc.
7. Strict quality control: Standardized operation procedures and multiple validation to ensure product quality.

## Top Selling Antibodies

Product Name	Cat.No.	Reactivity	Application	Host
Pan-Akt Polyclonal Antibody	E-AB-30465	H,M,R	WB,IHC,ELISA	Rabbit
ALPL Polyclonal Antibody	E-AB-40142	H,M,R	WB,IHC	Rabbit
AMPK alpha2 Polyclonal Antibody	E-AB-40256	H,R	WB,IHC	Rabbit
Phospho-AMPK alpha1/2 (Thr183/172) Polyclonal Antibody	E-AB-21121	H,M,R,Mk	WB,IHC,ELISA	Rabbit
AMPK alpha1/2 Polyclonal Antibody	E-AB-30491	H,M,R,Mk	WB,IHC,ELISA	Rabbit
ANXA5 Polyclonal Antibody	E-AB-40114	H,M,R	WB,IHC	Rabbit
ANXA6 Polyclonal Antibody	E-AB-40260	H,M,R	WB,IHC	Rabbit



Product Name	Cat.No.	Reactivity	Application	Host
ATXN3 Polyclonal Antibody	E-AB-40246	H,M,R	WB,IHC	Rabbit
BCL-x Polyclonal Antibody	E-AB-40057	H,M,R	WB,IHC	Rabbit
BMP2 Polyclonal Antibody	E-AB-70111	H,M,R	WB,IHC	Rabbit
BRC1 Polyclonal Antibody	E-AB-40282	H	IHC	Rabbit
CAMP Polyclonal Antibody	E-AB-40068	H,M,R	WB,IHC	Rabbit
Cleaved-CASP1 (D210) Polyclonal Antibody	E-AB-30002	H	WB,IHC,ELISA	Rabbit
CASP1 Polyclonal Antibody	E-AB-30752	H,M,R	WB,IHC,ELISA	Rabbit
Cleaved-CASP3 p17 (D175) Polyclonal Antibody	E-AB-30004	H,M,R	WB,IHC,ELISA	Rabbit
CASP3 Polyclonal Antibody	E-AB-30756	H,M,R	WB,IHC,IF,ELISA	Rabbit
Catenin beta Polyclonal Antibody	E-AB-70005	H,M,R	WB,IHC,IF	Rabbit
CBX5 Polyclonal Antibody	E-AB-40261	H	IHC	Rabbit
CCL2 Polyclonal Antibody	E-AB-40274	H	IHC	Rabbit
CD14 Polyclonal Antibody	E-AB-70113	H,M,R	WB,IHC	Rabbit
CD31 Polyclonal Antibody	E-AB-30820	H,M	WB,IHC,IF,ELISA	Rabbit
CD54 Polyclonal Antibody	E-AB-70046	H,M,R	WB,IHC,IF	Rabbit
CDC25A Polyclonal Antibody	E-AB-70140	H,M,R	WB,IHC	Rabbit
CDK2 Polyclonal Antibody	E-AB-40095	H,M	WB,IHC	Rabbit
CDK4 Polyclonal Antibody	E-AB-40286	H,M,R	WB	Rabbit
c-Fos Monoclonal Antibody	E-AB-70204	H,M,R	WB,IHC,IF	Mouse
CK-18 Polyclonal Antibody	E-AB-70095	H,M,R	WB,IHC	Rabbit
CK-19 Polyclonal Antibody	E-AB-70083	H,M,R	WB,IHC	Rabbit
CK-20 Polyclonal Antibody	E-AB-40067	H,M,R	WB,IHC	Rabbit
Connexin 43 Polyclonal Antibody	E-AB-70097	H,M,R	WB,IHC	Rabbit
COX4I1 Polyclonal Antibody	E-AB-70109	H,M,R	WB,IHC,IF	Rabbit
COX5B Polyclonal Antibody	E-AB-40083	H,M,R	WB,IHC	Rabbit
Phospho-CREB1 (Ser133) Polyclonal Antibody	E-AB-20849	H,M,R	WB,IHC,IP,ELISA	Rabbit
CRP Polyclonal Antibody	E-AB-40277	R,M	IHC	Rabbit
CST5 Polyclonal Antibody	E-AB-40105	H,M	WB,IHC	Rabbit
Cyclin A Polyclonal Antibody	E-AB-31069	H,M,R,Mk	WB,IHC,IF,ELISA	Rabbit
Cyclin B1 Polyclonal Antibody	E-AB-70114	H,M,R	WB,IHC,IF	Rabbit
Cyclin E1 Polyclonal Antibody	E-AB-31074	H,M,R	WB,IHC,IF,ELISA	Rabbit
CYCS Polyclonal Antibody	E-AB-70033	H,M,R	WB,IHC,IF	Rabbit
Cytochrome c1 Polyclonal Antibody	E-AB-40271	H,M,R	WB,IHC	Rabbit
E-Cadherin Polyclonal Antibody	E-AB-40285	H	IHC	Rabbit
EEF1E1 Polyclonal Antibody	E-AB-40118	H,M,R	WB,IHC	Rabbit
EGFR Polyclonal Antibody	E-AB-40283	H	IHC	Rabbit
Phospho-EIF2 alpha (Ser51) Polyclonal Antibody	E-AB-20864	H,M,R	WB,IHC,ELISA	Rabbit



Product Name	Cat.No.	Reactivity	Application	Host
EIF2 alpha Polyclonal Antibody	E-AB-31295	H,M,R,Mk	WB,IHC,ELISA	Rabbit
ENO1 Polyclonal Antibody	E-AB-40064	H,M,R	WB,IHC	Rabbit
EpCAM Polyclonal Antibody	E-AB-70132	H,M,R	WB,IHC,IF	Rabbit
ER alpha Polyclonal Antibody	E-AB-70088	H,M,R	WB,IHC,IF	Rabbit
ERK 1/2 Polyclonal Antibody	E-AB-31374	H,M,R	WB,IHC,ELISA	Rabbit
ER beta Polyclonal Antibody	E-AB-31382	H,M,R	WB,IF,ELISA	Rabbit
FH Polyclonal Antibody	E-AB-40106	H,M,R	WB,IHC	Rabbit
FOXO3 Polyclonal Antibody	E-AB-31469	H,M,R	WB,IHC,IF,ELISA	Rabbit
GAP43 Polyclonal Antibody	E-AB-70039	H,M,R	WB,IHC	Rabbit
GFAP Polyclonal Antibody	E-AB-70040	H,M,R	WB,IHC,IF	Rabbit
Glucocorticoid Receptor Polyclonal Antibody	E-AB-70157	H,M,R	WB,IHC	Rabbit
GLUT-1 Polyclonal Antibody	E-AB-31556	H,M,R	WB,IHC,ELISA	Rabbit
GLUT-3 Polyclonal Antibody	E-AB-70154	H,M,R	WB,IHC	Rabbit
GLUT-4 Polyclonal Antibody	E-AB-31558	H,M,R	WB,IHC,ELISA	Rabbit
GSK3 beta Polyclonal Antibody	E-AB-70043	H,M,R	WB,IHC,IF	Rabbit
HIF1 alpha Polyclonal Antibody	E-AB-31662	H,M,R	WB,IHC,ELISA	Rabbit
Histone H2A.X Polyclonal Antibody	E-AB-31679	H,M,R	WB,IHC,IF,ELISA	Rabbit
Histone H3 Monoclonal Antibody	E-AB-22003	H,M,R,Yeast	WB,IF,IP	Mouse
HMG20B Polyclonal Antibody	E-AB-40258	H	WB	Rabbit
HMGB1 Polyclonal Antibody	E-AB-70044	H,M,R	WB,IHC,IF	Rabbit
HMGB2 Polyclonal Antibody	E-AB-40199	H	WB,IHC	Rabbit
HSPA5 Polyclonal Antibody	E-AB-31742	H,M,R,Mk	WB,IHC,IF,ELISA	Rabbit
HSP60 Polyclonal Antibody	E-AB-70103	H,M,R	WB,IHC	Rabbit
HSP70 Polyclonal Antibody	E-AB-40007	H,M,R	WB,IHC	Rabbit
HSPA1A Polyclonal Antibody	E-AB-70102	H,M,R	WB,IHC	Rabbit
IFITM3 Polyclonal Antibody	E-AB-40263	H,R	WB,IHC	Rabbit
IHH Polyclonal Antibody	E-AB-40281	H	IHC	Rabbit
IL2 Polyclonal Antibody	E-AB-40275	H,M	IHC	Rabbit
IL6 Polyclonal Antibody	E-AB-70050	H,M,R	WB,IHC,IF	Rabbit
ILF2 Polyclonal Antibody	E-AB-40139	H,R	WB,IHC	Rabbit
IMP3 Polyclonal Antibody	E-AB-40171	H	WB,IHC	Rabbit
IRF3 Polyclonal Antibody	E-AB-40145	H,R	WB,IHC	Rabbit
ITGAX Polyclonal Antibody	E-AB-70017	H,M,R	WB,IHC,IF	Rabbit
Phospho-IkB alpha (Ser32/S36) Polyclonal Antibody	E-AB-20911	H,M,R,Mk	WB,IHC,IF,ELISA	Rabbit
IkB alpha Polyclonal Antibody	E-AB-31839	H,M,R	WB,ELISA	Rabbit
JAG1 Polyclonal Antibody	E-AB-30148	H,M,R	WB,ELISA	Rabbit
JAK1 Polyclonal Antibody	E-AB-31843	H,M,R	WB,IHC,ELISA	Rabbit



Product Name	Cat.No.	Reactivity	Application	Host
KI67 Polyclonal Antibody	E-AB-31869	H	IHC,IF,ELISA	Rabbit
Lamin B1 Polyclonal Antibody	E-AB-40257	H,M,R	WB,IHC	Rabbit
LC3A/B Polyclonal Antibody	E-AB-70053	H,M,R	WB,IHC	Rabbit
LMNA Polyclonal Antibody	E-AB-40252	H	WB	Rabbit
MAPK13 Polyclonal Antibody	E-AB-40183	H	WB,IHC	Rabbit
MCM3 Polyclonal Antibody	E-AB-40253	H	WB	Rabbit
MMP2 Polyclonal Antibody	E-AB-32054	H,M,R,Mk	WB,IHC,IF,ELISA	Rabbit
MMP9 Polyclonal Antibody	E-AB-31531	H	WB,IHC,IF,ELISA	Rabbit
Phospho-MTOR (Ser2448) Polyclonal Antibody	E-AB-20929	H,M,R	WB,IHC,IF,ELISA	Rabbit
MUC5AC Polyclonal Antibody	E-AB-40276	H	IHC	Rabbit
MYD88 Polyclonal Antibody	E-AB-32136	H,M,R	WB,IHC,ELISA	Rabbit
N-cadherin Polyclonal Antibody	E-AB-70061	H,M,R	WB,IHC	Rabbit
NF $\kappa$ B-p65 Polyclonal Antibody	E-AB-32233	H,M,R	WB,IHC,ELISA	Rabbit
Phospho-NF $\kappa$ B-p65 (Ser276) Polyclonal Antibody	E-AB-20937	H,M,R	WB,IHC,IP,ELISA	Rabbit
Phospho-NF $\kappa$ B-p65 (Ser536) Polyclonal Antibody	E-AB-20941	H,M,R,Mk	WB,IHC,IP,ELISA	Rabbit
NGF Polyclonal Antibody	E-AB-32239	H,M,R	WB,IHC,ELISA	Rabbit
NOS1 Polyclonal Antibody	E-AB-70065	H,M,R	WB,IHC	Rabbit
NOS2 Polyclonal Antibody	E-AB-32264	H,M,R	WB,ELISA	Rabbit
SPP1 Polyclonal Antibody	E-AB-32430	H,M,R	WB,IHC,ELISA	Rabbit
p16 Polyclonal Antibody	E-AB-40264	H,R	IHC	Rabbit
Phospho-P38 (Thr180/Tyr182) Polyclonal Antibody	E-AB-21027	H,M,R	WB,IHC,ELISA	Rabbit
Di-Methyl-P53 (Lys370) Polyclonal Antibody	E-AB-30091	H	WB,ELISA	Rabbit
P53 Monoclonal Antibody	E-AB-22015	H	WB,IHC,IF	Mouse
p70 S6 kinase alpha Polyclonal Antibody	E-AB-70159	H,M,R	WB,IHC	Rabbit
PAFAH Polyclonal Antibody	E-AB-40018	H,M,R	WB,IHC	Rabbit
Cleaved-PARP1 (D214) Polyclonal Antibody	E-AB-30080	H,M	WB,ELISA	Rabbit
PDIA5 Polyclonal Antibody	E-AB-40121	H,M	WB,IHC	Rabbit
POSTN Polyclonal Antibody	E-AB-40279	H	WB,IHC	Rabbit
PPARG Polyclonal Antibody	E-AB-32647	H,M,R	WB,IHC,ELISA	Rabbit
PPM1A Polyclonal Antibody	E-AB-40254	M,R	WB	Rabbit
PCNA Polyclonal Antibody	E-AB-40056	H,M,R	WB,IHC	Rabbit
PTEN Polyclonal Antibody	E-AB-40284	H,M	WB,IHC	Rabbit
RAN Polyclonal Antibody	E-AB-40251	H,M,R	WB	Rabbit
SDHA Polyclonal Antibody	E-AB-40089	H,R	WB,IHC	Rabbit
SMAD3 Polyclonal Antibody	E-AB-40050	H,R	WB,IHC	Rabbit
SNAI1 Polyclonal Antibody	E-AB-32931	H,M,Mk	WB,IHC,IF,ELISA	Rabbit
SPTAN1 Polyclonal Antibody	E-AB-40268	H,M,R	WB,IHC	Rabbit



Product Name	Cat.No.	Reactivity	Application	Host
STAT1 Polyclonal Antibody	E-AB-40052	H,R	WB,IHC	Rabbit
STAT3 Polyclonal Antibody	E-AB-32982	H,M,R	WB,IHC,IF,ELISA	Rabbit
Tau Polyclonal Antibody	E-AB-70075	H,M,R	WB,IHC	Rabbit
TIMP1 Polyclonal Antibody	E-AB-33105	H,M,R	WB,ELISA	Rabbit
TRIP4 Polyclonal Antibody	E-AB-40162	H,M,R	WB,IHC	Rabbit
UQCRCF51 Polyclonal Antibody	E-AB-40082	H,M	WB,IHC	Rabbit
Vimentin Polyclonal Antibody	E-AB-70081	H,M,R	WB,IHC	Rabbit

## Tag & Control Antibodies

\* The antibody marked with "\*" in the column of Reactivity is suitable for all kinds of samples that contain the Tag proteins.

Product Name	Cat.No.	Reactivity	Application	Host
Flag-Tag Monoclonal Antibody	E-AB-20006	*	WB, IF, IP	Mouse
Flag-Tag Polyclonal Antibody	E-AB-20003	*	WB, ELISA	Rabbit
GFP Monoclonal Antibody	E-AB-20086	*	WB, IP	Mouse
HA-Tag Monoclonal Antibody	E-AB-20008	*	WB, IP, IF	Mouse
HA-Tag Polyclonal Antibody	E-AB-20002	*	WB, ELISA	Rabbit
His-Tag Monoclonal Antibody	E-AB-20009	*	WB, IP, IF	Mouse
His-Tag Polyclonal Antibody	E-AB-20001	*	WB, ELISA	Rabbit
Myc-Tag Monoclonal Antibody	E-AB-20007	*	WB, IP, IF	Mouse
GAPDH Monoclonal Antibody	E-AB-20032	H, M, R, Mk	WB, IHC	Mouse
GAPDH Polyclonal Antibody	E-AB-20059	H, M, R, Mk	WB, IHC	Rabbit
alpha Tubulin Polyclonal Antibody	E-AB-20069	H, M, R	WB, IHC, IF, ELISA	Rabbit
alpha Tubulin Monoclonal Antibody	E-AB-20036	H, M, R	WB, IP	Mouse
beta actin Monoclonal Antibody	E-AB-20031	H, M, R, Mk	WB, IHC	Mouse
beta actin Polyclonal Antibody	E-AB-20058	H, M, R, Mk	WB, IHC	Rabbit
beta Tubulin Monoclonal Antibody	E-AB-20033	H, M, R, Mk	WB, IHC	Mouse
Lamin B1 Monoclonal Antibody	E-AB-20076	H, M, R	WB, IP	Mouse

We are constantly adding new products to provide a complete product catalogue for your research. Please visit our website, [www.elabscience.com](http://www.elabscience.com) for the most up-to-date details on all the products we have.

For ordering, please email us : [orders@elabscience.com](mailto:orders@elabscience.com).

**Elabscience Biotechnology Inc.**

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-240-252-7368

Fax: 1-240-252-7376

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

