

## Purified Anti-Human CD10 Antibody[B-E3]

Catalog No. E-AB-F1220A  
Storage  
Applications FCM

Reactivity  
Host  
Isotype

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

### Images

### Immunogen Information

### Product Information

#### For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit [www.elabscience.com](http://www.elabscience.com).

#### Focus on your research Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.