

## Purified Anti-Human CD15 Antibody[BRA-4F1]

**Catalog No.** E-AB-F1222A  
**Storage**  
**Applications** FCM

**Reactivity**  
**Host**  
**Isotype** Mouse IgM, κ

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

### Images

### Immunogen Information

### Product Information

#### For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit [www.elabscience.com](http://www.elabscience.com).

**Focus on your research  
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.