


IFNAR1 Polyclonal Antibody


Catalog No.	E-AB-40306	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human Duodenum using IFNAR1 Polyclonal Antibody at dilution of 1:200.



Immunohistochemistry of paraffin-embedded Human liver using IFNAR1 Polyclonal Antibody at dilution of 1:200.

Immunogen Information

Immunogen	Recombinant Human Interferon alpha/beta receptor 1 protein
GeneID	3454
Swissprot	P17181
Synonyms	Alpha type antiviral protein, Antiviral protein, interferon receptor 1, Interferon-alpha receptor

Product Information

Buffer	PBS with 0.05% Proclin300 and 50% glycerol, pH7.4.
Purify	Antigen Affinity Purification
Dilution	IHC 1:100-1:300

Background

IFNAR1 (Interferon Alpha And Beta Receptor Subunit 1) is a Protein Coding gene. Diseases associated with IFNAR1 include Septicemic Plague and Hepatitis C. Among its related pathways are Type I Interferon Signaling Pathways and Immune response IFN alpha/beta signaling pathway. GO annotations related to this gene include interferon receptor activity and type I interferon receptor activity. An important paralog of this gene is ENSG00000249624. The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. The encoded protein also functions as an antiviral factor.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.