


DNMT1 Polyclonal Antibody


Catalog No.	E-AB-40321	Reactivity	H,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human ovary using DNMT1 Polyclonal Antibody at dilution of 1:200.



Immunohistochemistry of paraffin-embedded Rat brain using DNMT1 Polyclonal Antibody at dilution of 1:200.

Immunogen Information

Immunogen	Recombinant Rat DNA (cytosine-5)-methyltransferase 1 protein
GeneID	84350
Swissprot	Q9Z330
Synonyms	ADCADN,AIM,CXXC finger protein 9,DNA MTase HsaI, DNMT 1, DNMT, Dnmt1, DNMT1, Dnmt1 o

Product Information

Buffer	PBS with 0.05% Proclin300 and 50% glycerol, pH7.4.
Purify	Antigen Affinity Purification
Dilution	IHC 1:100-1:300

Background

DNA (cytosine-5-)-methyltransferase 1 has a role in the establishment and regulation of tissue-specific patterns of methylated cytosine residues. Aberrant methylation patterns are associated with certain human tumors and developmental abnormalities. Two transcript variants encoding different isoforms have been found for this gene.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.