

## Rabbit Anti-Human IgG(H+L)(Cyanine3 conjugated)

Catalog No. E-AB-1013

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

### Product Information

**Immunogen:** Human IgG

**Purification Method:** Affinity purification

**Conjugation:** CY3

**Concentration:** 0.2 mg/mL(specific antibody purification≥99%),CY3:Ab=3-9:1

**Applications:** IF

**Dilution:** Immunofluorescence: 1:50-100

Optimal working dilutions must be determined by end user. Dilute before when it's needed, and the residual should be discarded.

**Amax:** 554nm

**E<sub>max</sub>:** 568-574nm

**Buffer:** PBS with 50% glycerol, PH7.4

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.

### For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit [www.elabscience.com](http://www.elabscience.com).

Focus on your research  
Service for life science

Applications: Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA- Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFCM: Intracellular Staining for Flow Cytometry; WB- Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation