

Goat Anti-Mouse IgG

Catalog No. E-AB-1035

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Immunogen: Mouse IgG

Purification Method: Affinity purification

Conjugation: Unconjugated

Concentration: 1 mg/mL (specific antibody purification ≥99%)

Buffer: PBS, PH7.4

Storage: Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research
Service for life science

Applications: Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA- Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFCM: Intracellular Staining for Flow Cytometry; WB- Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation