

Goat Anti-Rabbit IgG(H+L)(Elab Fluor® 488 conjugated)

Catalog No. E-AB-1055

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Immunogen: Rabbit IgG

Purification Method: Affinity purification

Conjugation: AF488

Concentration: 0.5 mg/mL

Applications: IF

Dilution: Immunofluorescence: 1:100-300

Optimal working dilutions must be determined by end user. Dilute before when it's needed, and the residual should be discarded.

Amax: 493nm

Emax: 519nm

Buffer: PBS with 0.05% sodium azide and 50% glycerol, PH7.4

Storage: Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications: Activ- Activation; Block- Blocking; Separation- Cell Separation ; Cell Sep-Neg- Cell Separation by Negative Selection; FA- Functional Assay; Neut- Neutralization; Stim- Stimulation; FCM- Flow Cytometry; ICFM: Intracellular Staining for Flow Cytometry; WB- Western Blotting; IHC- Immunohistochemistry; IF- Immunofluorescence; IP- Immunoprecipitation