


PHLDA3 Polyclonal Antibody


Catalog No.	E-AB-19111	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC,ELISA	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using PHLDA3 Polyclonal Antibody at dilution of 1:110(×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using PHLDA3 Polyclonal Antibody at dilution of 1:110(×200)

Immunogen Information

Immunogen	Fusion protein of human PHLDA3
Gene Accession	BC014390
Swissprot	Q9Y5J5
Synonyms	PHLA3,PHLDA 3,phlda3,TDAG51/Ipl homologue 1,TIH 1,TIH1

Product Information

Buffer	PBS with 0.05% NaN3 and 40% Glycerol,pH7.4
Purify	Antigen affinity purification
Dilution	IHC 1:100-1:300, ELISA 1:5000-1:10000

Background

PHLDA3 (Pleckstrin Homology Like Domain Family A Member 3) is a Protein Coding gene. Among its related pathways are Apoptosis and Autophagy and PI3K / Akt Signaling. GO annotations related to this gene include phosphatidylinositol-4,5-bisphosphate binding and phosphatidylinositol-3-phosphate binding. An important paralog of this gene is PHLDA1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.