


ENPEP Polyclonal Antibody

Catalog No.	E-AB-19620	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC,ELISA	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ENPEP Polyclonal Antibody at dilution of 1:65(×200)

Immunogen Information

Immunogen	Synthetic peptide of human ENPEP
Gene Accession	NP001968
Swissprot	Q07075
Synonyms	Aminopeptidase A,AMPE,AP-A,APA,bp-1,bp1,CD249,EAP,Enpep,Gp160,Ly51

Product Information

Buffer	PBS with 0.05% NaN3 and 40% Glycerol,pH7.4
Purify	Antigen affinity purification
Dilution	IHC 1:50-1:300, ELISA 1:5000-1:10000

Background

ENPEP (Glutamyl Aminopeptidase) is a Protein Coding gene. Diseases associated with ENPEP include Neurogenic Hypertension and Gestational Choriocarcinoma. Among its related pathways are Metabolism of proteins and Peptide hormone metabolism. GO annotations related to this gene include peptidase activity and aminopeptidase activity. An important paralog of this gene is ANPEP.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.