


RFP-Tag Monoclonal Antibody

Catalog No.	E-AB-20016	Reactivity	
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Mouse
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western Blot analysis of 1µg RFP fusion protein using RFP-Tag Monoclonal Antibody at dilution of 1) 1:3000 2) 1:5000 3) 1:10000.

Immunogen Information

Immunogen	Recombinant Protein
Synonyms	DsRed,GFP like chromoprotein,red fluorescent protein,RFP

Product Information

Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Purify	Protein A purification
Clone No.	Clone:1L6
Dilution	WB 1:5000-1:10000

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.