

T7-Tag Monoclonal Antibody

Catalog No.	E-AB-20017	Reactivity	
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Mouse
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



www.elabscience.com

Western Blot analysis of 0.5ug T7 fusion protein using T7-Tag Monoclonal Antibody at dilution of 1) 1:3000 2) 1:5000 3) 1:10000.

Immunogen Information

Immunogen	Synthetic Peptide
Synonyms	10, Major capsid protein 10A, Major capsid protein, Minor capsid protein 10B, T7 epitope tag

Product Information

Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Purify	Protein A purification
Clone No.	Clone:1M4
Dilution	WB 1:5000-1:10000

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.