

StrepTag II Monoclonal Antibody

Catalog No. E-AB-20042

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG
Clone	Clone:1L2
Purification	Protein A purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:2000-5000
-----------	-------------

Data



www.elabscience.com

Western Blot analysis of recombinant protein using
StrepTag II Monoclonal Antibody at dilution of
1:5000.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only