


BOLA1 Polyclonal Antibody


Catalog No.	E-AB-52476	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC,ELISA	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using BOLA1 Polyclonal Antibody at dilution of 1:50(×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using BOLA1 Polyclonal Antibody at dilution of 1:50(×200)

Immunogen Information

Immunogen	Full length fusion protein
Gene Accession	BC063405
Swissprot	Q9Y3E2
Synonyms	BolA Family Member 1,BolA-Like Protein 1,BolA Homolog 1 (E. Coli),BolA-Like 1 (E. Coli),CGI-143,HBolA

Product Information

Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol,pH7.4
Purify	Antigen affinity purification
Dilution	IHC 1:40-1:200, ELISA 1:5000-1:10000

Background

BOLA1 (BolA Family Member 1) is a Protein Coding gene. An important paralog of this gene is BOLA2.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.