


GATD3B Polyclonal Antibody

Catalog No.	E-AB-52500	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB,ELISA	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of 293T cell lysates using GATD3B Polyclonal Antibody at dilution of 1:250

Immunogen Information

Immunogen	Fusion protein of human GATD3B
Gene Accession	BC003587
Swissprot	A0A0B4J2D5
Synonyms	ES1 Protein Homolog, Mitochondrial, Uncharacterized Protein GATD3B

Product Information

Calculated MW	28 kDa
Observed MW	Refer to figures
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol,pH7.4
Purify	Antigen affinity purification
Dilution	WB 1:500-1:2000, ELISA 1:5000-1:10000

Background

GATD3B (Glutamine Amidotransferase Like Class 1 Domain Containing 3B) is a Protein Coding gene. Among its related pathways are Gastric Cancer Network 1. An important paralog of this gene is GATD3A.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.