


## NADK2 Polyclonal Antibody

<b>Catalog No.</b>	E-AB-52506	<b>Reactivity</b>	H,R
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.	<b>Host</b>	Rabbit
<b>Applications</b>	WB,ELISA	<b>Isotype</b>	IgG

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

### Images



Western blot analysis of HEPG2 HeLa and A549 cell lysates using NADK2 Polyclonal Antibody at dilution of 1:600

### Immunogen Information

<b>Immunogen</b>	Fusion protein of human NADK2
<b>Gene Accession</b>	BC062567
<b>Swissprot</b>	Q4G0N4
<b>Synonyms</b>	C5orf33,CE033,Chromosome 5 open reading frame 33 ,FLJ30596,Hypothetical protein LOC133686,MGC 43298,NADKD1,OTTHUMP00000161732,UPF0465 protein C5orf33

### Product Information

<b>Calculated MW</b>	49 kDa
<b>Observed MW</b>	Refer to figures
<b>Buffer</b>	PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol,pH7.4
<b>Purify</b>	Antigen affinity purification
<b>Dilution</b>	WB 1:500-1:2000, ELISA 1:5000-1:10000

### Background

This gene encodes a mitochondrial kinase that catalyzes the phosphorylation of NAD to yield NADP. Mutations in this gene result in 2,4-dienoyl-CoA reductase deficiency. Alternative splicing results in multiple transcript variants.

#### For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit [www.elabscience.com](http://www.elabscience.com).

**Focus on your research  
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.