


DTYMK Polyclonal Antibody


Catalog No.	E-AB-52710	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB,IHC,ELISA	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of 293T cell and Human fetal liver tissue lysates using DTYMK Polyclonal Antibody at dilution of 1:550



Immunohistochemistry of paraffin-embedded Human tonsil tissue using DTYMK Polyclonal Antibody at dilution of 1:70(×200)

Immunogen Information

Immunogen	Fusion protein of human DTYMK
Gene Accession	BC001827
Swissprot	P23919
Synonyms	CDC8,dTMP kinase,Dtymk,FLJ44192,KTHY,MGC198617,Thymidylate kinase,TMPK,TYMK

Product Information

Calculated MW	24 kDa
Observed MW	Refer to figures
Buffer	PBS with 0.05% NaN3 and 40% Glycerol,pH7.4
Purify	Antigen affinity purification
Dilution	WB 1:500-1:2000, IHC 1:50-1:300, ELISA 1:5000-1:10000

Background

DTYMK (Deoxythymidylate Kinase) is a Protein Coding gene. Among its related pathways are superpathway of pyrimidine deoxyribonucleotides de novo biosynthesis and Purine metabolism (KEGG). GO annotations related to this gene include kinase activity and thymidylate kinase activity.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.