


UBD Polyclonal Antibody


Catalog No.	E-AB-52823	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC,ELISA	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human tonsil tissue using UBD Polyclonal Antibody at dilution of 1:50(×200)



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using UBD Polyclonal Antibody at dilution of 1:50(×200)

Immunogen Information

Immunogen	Fusion protein of human UBD
Gene Accession	BC012472
Swissprot	O15205
Synonyms	Diubiquitin,FAT10,GABBR1,UBD 3,Ubd,UBD,Ubiquitin D,Ubiquitin like protein FAT10,Ubiquitin-like protein FAT10

Product Information

Buffer	PBS with 0.05% NaN3 and 40% Glycerol,pH7.4
Purify	Antigen affinity purification
Dilution	IHC 1:50-1:300, ELISA 1:5000-1:10000

Background

UBD (Ubiquitin D) is a Protein Coding gene. Diseases associated with UBD include Hepatocellular Carcinoma. Among its related pathways are TRAF Pathway and fMLP Pathway. GO annotations related to this gene include proteasome binding. An important paralog of this gene is UBB.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-
 Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-
 Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.