


GLG1 Polyclonal Antibody

Catalog No.	E-AB-65865	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using GLG1 Polyclonal Antibody at dilution of 1:3000.

Immunogen Information

Immunogen	Recombinant fusion protein of human GLG1 (NP_036333.2).
GeneID	2734
Swissprot	Q92896
Synonyms	GLG1,CFR-1,ESL-1,MG-160,MG160

Product Information

Calculated MW	134kDa/135kDa/137kDa
Observed MW	135kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:1000

Background

GLG1 (Golgi Glycoprotein 1) is a Protein Coding gene. Diseases associated with GLG1 include Ecthyma and Trichinosis. Among its related pathways are Cell surface interactions at the vascular wall and Apoptotic Pathways in Synovial Fibroblasts. Gene Ontology (GO) annotations related to this gene include signaling receptor binding.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.